

VERSION 1

JAN 18, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.6qpvr6xjpvmk/v1

Protocol Citation: Yan Chen, Jennifer Gin, Christopher J Petzold 2023. Alkaline-SDS cell lysis of microbes with acetone protein precipitation for proteomic sample preparation in 96-well plate format. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr6xjpvmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Dec 09, 2021

Last Modified: Jan 18, 2023

PROTOCOL integer ID:
 55810

🌐 Alkaline-SDS cell lysis of microbes with acetone protein precipitation for proteomic sample preparation in 96-well plate format V.1

🔗 Forked from [Acetone-Salt Protein Extraction for bacteria and fungi - 1.7mL tube](#)


Yan Chen¹, Jennifer Gin¹, Christopher J Petzold¹

¹Lawrence Berkeley National Laboratory

LBNL omics

Agile BioFoundry

[1 more workspace](#) ↓



Jennifer Gin

Lawrence Berkeley National Laboratory

ABSTRACT

This high-throughput protocol details the steps to extract protein from Gram-negative bacteria, Gram-positive bacteria, or non-filamentous fungi in 96-well plate format for quantitative proteomic workflows. This protocol uses a bench-top automated liquid dispenser but the volumes and times also apply to manual and multi-channel pipetter use. This protocol is designed for lab-based, culture conditions and synthetic community experiments where complex sample matrices are minimized. Additional sample preservation and/or protein extraction methods may be required for environmental samples (e.g., feces, soil) to minimize protein degradation and maintain sample integrity.

This protocol works best as part of a high-throughput proteomic sample preparation workflow with:

[Automated Protein Quantitation with the Biomek-FX liquid handler system](#)

and

[Automated Protein Normalization and Tryptic Digestion on a Biomek-NX Liquid Handler System](#)

Keywords: Proteomics,
Sample preparation, Bacteria,
Fungi, Protein extraction

GUIDELINES

- All centrifuge steps use an Eppendorf 5810R centrifuge.
- You can use Qiagen Lysis Buffer P2 (Qiagen, Cat. #19052) as the alkaline-SDS cell lysis buffer.

Note

Increase reagent volumes accordingly when processing larger amounts of biomass. This protocol works in larger formats (e.g., 1.7 mL Eppendorf tubes) too.

MATERIALS

PCR Plate 96-well non-skirted Thermo Fisher Scientific Catalog #AB0600

Hydrochloric acid Sigma Catalog #320331

Sodium Hydroxide (200 mM)

1% Sodium Dodecyl Sulfate (SDS)

Ammonium bicarbonate LC-MS grade [VWR Scientific](#) Catalog #[BJ40867-50G](#)

Benzonase nuclease Millipore Catalog #70746

Acetone Sigma Catalog #179124

LC-MS grade Methanol VWR Scientific, Catalog #BJLC230-2.5

Optional:

- Qiagen Lysis Buffer P2 (Qiagen, Cat. #19052) in place of NaOH and SDS

SAFETY WARNINGS

! Acetone is used in this protocol so please follow the appropriate safety guidelines for handling and disposing of non-halogenated solvents at your institution.

Sodium Hydroxide is a HIGHLY CORROSIVE CHEMICAL and contact can severely irritate and burn the skin and eyes with possible eye damage. Inhaling Sodium Hydroxide can irritate the lungs.

Wear gloves and appropriate PPE for safety and to minimize contamination of samples.

BEFORE START INSTRUCTIONS

For this protocol you will need:

- a bench-top automated liquid dispenser (*e.g.*, Formulatrix Mantis) or manual/multi-channel pipetters
- an Eppendorf 5810R centrifuge with S-4-104 rotor or similar centrifuge


Mix at least 3 mL of NaOH/SDS buffer for final concentrations of:

- 200 mM NaOH
- 1% SDS

or use Qiagen Lysis Buffer P2 (Qiagen, Cat.#19052)

Cell lysis

6m

- 1 Start with  10 μL of cells per well a non-skirted PCR plate (Thermo Scientific, Cat.#AB0600).





Plate of cell pellets

2 Add  25 μL of alkaline-SDS cell lysis buffer (200 mM NaOH, 1% SDS) to each well.

1m

3 Resuspend the cell pellet in lysis buffer on a plate mixer.

5m



Note

Ensure proper mixing of the cell lysis buffer and the cell pellet.







Plate mixer



Cells resuspended in lysis buffer



Neutralization and Benzonase treatment

8m

- 4 Add  2.75 μL  1 Molarity (M) Hydrochloric acid (Sigma Cat.#320331) to each well.

2m



- 5 Add  25 μL 100 mM Ammonium bicarbonate (VWR Scientific Cat.#BJ40867-50G) and  0.5 μL Benzonase nuclease (Millipore Cat.#70746) to each well.

3m



- 6 Mix thoroughly on the plate mixer.

3m






Cell lysate after benzonase treatment

Salt-Acetone protein precipitation

8m

- 7 Add  200 μ L 100% Acetone (Sigma Cat.#179124) to each well and let sit at room temperature for 5 minutes. 5m





Protein precipitation after mixing in 100% acetone

8 Centrifuge at  4000 rpm, 00:02:00 .


2m



9 Remove supernatant.

1m






Protein pellet after centrifugation and removal of supernatant

Wash and resuspend protein


8m

- 10 Wash protein pellet twice using  150 μL 80% Acetone (Sigma Cat.#179124) .

6m


Note

Increase the number of washes depending on the sample matrix background.

- 11 Add  60 μL of **1M** 100 millimolar (mM) Ammonium bicarbonate in 10% Methanol to each well to resuspend protein pellet.

2m





Protein pellet resuspended in 60 μ L of Ammonium bicarbonate in 10% Methanol

- 12 Store at $-20\text{ }^{\circ}\text{C}$ until ready for [Automated Protein Quantitation with the Biomek-FX liquid handler system](#).