

# Determination of Ash in Biomass

## Laboratory Analytical Procedure (LAP)

**Issue Date: 7/17/2005**

A. Sluiter, B. Hames, R. Ruiz, C. Scarlata,  
J. Sluiter, and D. Templeton

**Technical Report**  
**NREL/TP-510-42622**  
**January 2008**

NREL is operated by Midwest Research Institute • Battelle Contract No. DE-AC36-99-GO10337




# Determination of Ash in Biomass

## Laboratory Analytical Procedure (LAP)

Issue Date: 7/17/2005

A. Sluiter, B. Hames, R. Ruiz, C. Scarlata,  
J. Sluiter, and D. Templeton

*Technical Report*  
NREL/TP-510-42622  
January 2008



**National Renewable Energy Laboratory**  
1617 Cole Boulevard, Golden, Colorado 80401-3393  
303-275-3000 • [www.nrel.gov](http://www.nrel.gov)

Operated for the U.S. Department of Energy  
Office of Energy Efficiency and Renewable Energy  
by Midwest Research Institute • Battelle

Contract No. DE-AC36-99-GO10337

## DISCLAIMER

These Standard Biomass Analytical Methods (“Methods”) are provided by the National Renewable Energy Laboratory (“NREL”), which is operated by the Midwest Research Institute (“MRI”) for the Department Of Energy.

Access to and use of these Methods shall impose the following obligations on the user. The user is granted the right, without any fee or cost, to use, copy, modify, alter, enhance and distribute these Methods for any purpose whatsoever, except commercial sales, provided that this entire notice appears in all copies of the Methods. Further, the user agrees to credit NREL/MRI in any publications that result from the use of these Methods. The names NREL/MRI, however, may not be used in any advertising or publicity to endorse or promote any products or commercial entity unless specific written permission is obtained from NREL/MRI. The user also understands that NREL/MRI is not obligated to provide the user with any support, consulting, training or assistance of any kind with regard to the use of these Methods or to provide the user with any updates, revisions or new versions.

THESE METHODS ARE PROVIDED BY NREL/MRI "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL NREL/MRI BE LIABLE FOR ANY SPECIAL, INDIRECT OR CONSEQUENTIAL DAMAGES OR ANY DAMAGES WHATSOEVER, INCLUDING BUT NOT LIMITED TO CLAIMS ASSOCIATED WITH THE LOSS OF DATA OR PROFITS, WHICH MAY RESULT FROM AN ACTION IN CONTRACT, NEGLIGENCE OR OTHER TORTIOUS CLAIM THAT ARISES OUT OF OR IN CONNECTION WITH THE ACCESS, USE OR PERFORMANCE OF THESE METHODS.

# Procedure Title: Determination of Ash in Biomass

## Laboratory Analytical Procedure

### 1. Introduction

- 1.1 The amount of inorganic material in biomass, either structural or extractable, should be measured as part of the total composition. Structural ash is inorganic material that is bound in the physical structure of the biomass, while extractable ash is inorganic material that can be removed by washing or extracting the material. Extractable ash can be the result of soil remaining in the biomass. Refer to LAP “Total Mass Closure” for information on where to utilize this procedure.
- 1.2 This procedure is substantially similar to ASTM Standard Method Number E1755-01 “Standard Method for the Determination of Ash in Biomass”
- 1.3 This test method covers the determination of ash, expressed as the percentage of residue remaining after dry oxidation at 550 to 600°C. All results are reported relative to the 105°C oven dry weight of the sample.

### 2. Scope

- 2.1 This procedure is applicable to hard and soft woods, herbaceous materials, agricultural residues, wastepaper, and solid fraction process samples.
- 2.2 All analyses should be performed in accordance with an appropriate laboratory specific Quality Assurance Plan (QAP).

### 3. Terminology

- 3.1 *Oven dry weight (ODW)*- the weight of biomass mathematically corrected for the amount of moisture present in the sample at the time of weighing
- 3.2 *Prepared biomass*- biomass prepared according to LAP “Preparation of Samples for Biomass Compositional Analysis”.
- 3.3 *Ash*- The inorganic residue left after dry oxidation at 575°C.

### 4. Significance and Use

- 4.1 The ash content is a measure of the mineral content and other inorganic matter in biomass and is used in conjunction with other procedures to determine the total composition of biomass samples.
- 4.2 This procedure is used, in conjunction with other procedures to determine the chemical composition of biomass samples, see LAP “Summative Mass Closure for Biomass Samples”.

### 5. Interferences

- 5.1 None

### 6. Apparatus

- 6.1 Muffle furnace, equipped with a thermostat, set to  $575 \pm 25$  °C or equipped with optional ramping program

- 6.2 Analytical balance, accurate to 0.1 mg.
- 6.3 Desiccator containing desiccant
- 6.4 Ashing crucibles, 50 mL, porcelain, silica, or platinum
- 6.5 Porcelain markers, high temperature, or equivalent crucible marking method
- 6.6 Ashing burner, ignition source, tongs, and clay triangle with stand
- 6.7 Convection drying oven, with temperature control of  $105 \pm 3^{\circ}\text{C}$ , optional

## **7. Reagents and materials**

- 7.1 Reagents
  - 7.1.1 None
- 7.2 Materials
  - 7.2.1 None

## **8. ES&H Considerations and Hazards**

- 8.1 Use appropriate safety measures when handling open flame.
- 8.2 When placing crucibles in a furnace or removing them, use appropriate personal protective equipment, including heat resistant gloves.
- 8.3 Follow all applicable NREL chemical handling procedures

## **9. Sampling, Test Specimens and Test Units**

- 9.1 Care must be taken to ensure a representative sample is taken for analysis.
- 9.2 LAP “Preparation of Samples for Biomass Compositional Analysis” should be performed prior to this analysis
- 9.3 LAP “Determination of Solids in Biomass” should be performed at the same time that samples for ash analysis are weighed out. Alternately,  $105^{\circ}\text{C}$  oven dry biomass may be used, eliminating the need for total solids corrections.

## **10. Procedure**

- 10.1 Using a porcelain marker, mark an appropriate number of crucibles with identifiers, and place them in the muffle furnace at  $575 \pm 25^{\circ}\text{C}$  for a minimum of four hours. (Marking crucibles with a porcelain marker will permanently mark them, so a generic identifier is recommended.) Remove the crucibles from the furnace directly into a desiccator. If using a furnace set to  $575 \pm 25^{\circ}\text{C}$ , cool for a specific period of time, one hour is recommended. Record the cool time. Weigh the crucibles to the nearest 0.1 mg and record this weight.
- 10.2 Place the sample back into the muffle furnace at  $575 \pm 25^{\circ}\text{C}$  and dry to constant weight. Constant weight is defined as less than  $\pm 0.3$  mg change in the weight upon one hour of re-heating the crucible.
- 10.3 Weigh 0.5 to 2.0 g, to the nearest 0.1 mg, of a test specimen into the tared crucible. Record the sample weight. If the sample being analyzed is a  $105^{\circ}\text{C}$  dried test specimen, the sample should be stored in a desiccator until use. If air dry samples are used, LAP “Determination of Total Solids in Biomass” should be performed at the same time, to accurately measure the percent solids for correction. Each sample should be analyzed in duplicate, at minimum.

- 10.4 Ash the samples. Use step 10.4.1 if using a muffle furnace set to  $575 \pm 25$  °C. Use step 10.4.2 if using a muffle furnace with a ramping program.
- 10.4.1 Ash the samples using a muffle furnace set to  $575 \pm 25$  °C.
- 10.4.1.1 Using an ashing burner and clay triangle with stand, place the crucible over the flame until smoke appears. Immediately ignite the smoke and allow the sample to burn until no more smoke or flame appears. Allow the crucible to cool before placing it in the muffle furnace. Alternately, a furnace with a temperature ramping function may be used to avoid pre-ignition.
- 10.4.1.2 Place the crucibles in the muffle furnace at  $575 \pm 25$  °C for  $24 \pm 6$  hours. When handling the crucible, protect the sample from drafts to avoid mechanical loss of sample.
- 10.4.1.3 Carefully remove the crucible from the furnace directly into a desiccator and cool for a specific amount of time, equal to the initial cool time of the crucibles. Weigh the crucibles and ash to the nearest 0.1 mg and record the weight.
- 10.4.1.4 Place the sample back into the muffle furnace at  $575 \pm 25$  °C and ash to constant weight. Constant weight is defined as less than  $\pm 0.3$  mg change in the weight upon one hour of re-heating the crucible. When allowing samples to cool in a desiccator, it is necessary to maintain the initial cool time.
- 10.4.2 Ash the samples using a muffle furnace equipped with a ramping program.
- Furnace Temperature Ramp Program: Ramp from room temperature to 105 °C  
 Hold at 105°C for 12 minutes  
 Ramp to 250 °C at 10°C / minute  
 Hold at 250 °C for 30 minutes  
 Ramp to 575 °C at 20 °C / minute  
 Hold at 575 °C for 180 minutes  
 Allow temperature to drop to 105 °C  
 Hold at 105 °C until samples are removed
- 10.4.2.1 Place the crucibles in the muffle furnace and begin the ramping program. When handling the crucible, protect the sample from drafts to avoid mechanical loss of sample.
- 10.4.2.2 Carefully remove the crucible from the furnace directly into a desiccator and cool. Weigh the crucibles and ash to the nearest 0.1 mg and record the weight.
- 10.4.2.3 Place the sample back into the muffle furnace at  $575 \pm 25$  °C and ash to constant weight. Constant weight is defined as less than  $\pm 0.3$  mg change in the weight upon one hour of re-heating the crucible

## 11. Calculations

- 11.1 If an air dry sample was used, calculate the oven dry weight (ODW) of the sample, using the average total solids content as determined by the LAP “Standard Method for the Determination of Total Solids in Biomass”.

$$ODW = \frac{\text{Weight}_{\text{air dry sample}} \times \% \text{ Total solids}}{100}$$

11.2 Calculate and record the percentage ash on an ODW basis.

$$\% \text{ Ash} = \frac{\text{Weight}_{\text{crucible plus ash}} - \text{Weight}_{\text{crucible}}}{\text{ODW}_{\text{sample}}} \times 100$$

11.3 To report or calculate the relative percent difference (RPD) between two samples, use the following calculation

$$RPD = \left( \frac{X_1 - X_2}{X_{\text{mean}}} \right) \times 100$$

Where:

$X_1$  and  $X_2$  = measured values

$X_{\text{mean}}$  = the mean of  $X_1$  and  $X_2$

11.4 To report or calculate the root mean square deviation (RMS deviation) or the standard deviation (st dev) of the samples, use the following calculations.

First find the root mean square (RMS), of the sample using

$$RMS = x_m = \text{mean} = \sqrt{\frac{\sum_1^n x^2}{n}}$$

Then find the root mean square deviation, or standard deviation, using

$$RMS\text{deviation} = \sigma = \text{stdev} = \sqrt{\frac{\sum_1^n (x_i - x_m)^2}{n}}$$

Where:

$x_m$ =the root mean square of all x values in the set

$n$ =number of samples in set

$x_i$ =a measured value from the set

## 12. Report Format

12.1 Report ash as a percentage of the ODW of the sample. Report the percent solids, if applicable. Standard deviation and relative percent difference may also be reported.

## 13. Precision and Bias

13.1 *Round robin testing* – For a report documenting an international round robin test of biomass analysis methods, including this procedure, see Milne et al., 1992.

## 14. Quality Control

Reported Significant Figures or decimal places: Determined by data quality objectives and laboratory specific Quality Assurance Plan.

- 14.1 Replicates: Run all samples and method verification standards, if applicable, in duplicate, at minimum.
- 14.2 Blank: An empty aluminum dish or crucible should be run through the analysis. The dish should be weighed empty, ashed, and reweighed. The difference in weight must be less than the equivalent of a 0.5% error.
- 14.3 Relative percent difference criteria: Each sample must reproduce ash content  $\pm 0.5$  wt %.
- 14.4 Method verification standard (MVS): A MVS, such as a QA/QC sample, should be run in duplicate with every batch.
- 14.5 Sample size: 1-4 grams. If there is insufficient sample, the result should be flagged and the lack of precision noted.
- 14.6 Sample storage: If applicable, oven dried samples should be stored in a desiccator until ready to use.
- 14.7 Standard storage: Not applicable
- 14.8 Standard preparation: Not applicable
- 14.9 Definition of a batch: Any number of samples that are analyzed and recorded together. The maximum size of a batch will be limited by equipment constraints
- 14.10 Control charts: MVS or a QA/QC material should be control charted to verify reproducibility
- 14.11 Others: Biomass can rapidly gain or lose moisture when in contact with air. During the weighing steps, minimize the amount of time the sample and crucibles are exposed to the air.

## 15. Appendices

- 15.1 None

## 16. References

- 16.1 Moore, W., and D. Johnson. 1967. *Procedures for the Chemical Analysis of Wood and Wood Products*. Madison, WI: U.S. Forest Products Laboratory, U.S. Department of Agriculture.
- 16.2 NREL CAT Task Laboratory Analytical Procedure #005 "Standard Method for Ash in Biomass".
- 16.3 TAPPI Test Method T211, "Ash in Wood and Pulp." *In Tappi Test Methods*. Atlanta, GA: Technical Association of the Pulp and Paper Industry.
- 16.4 Milne, T. A.; Chum, H. L.; Agblevor, F. A.; Johnson, D. K. (1992). "Standardized Analytical Methods" Biomass & Bioenergy. Proceedings of International Energy Agency Bioenergy Agreement Seminar", 2-3 April 1992, Edinburgh, U.K. Vol. 2(1-6), 1992; pp. 341-366
- 16.5 ASTM E1755-01 "Standard Method for the Determination of Ash in Biomass" In *2003 Annual Book of ASTM Standards, Volume 11.05*. Philadelphia, PA: American Society for Testing and Materials, International.