

ezRAD Modifications to Illumina TruSeq DNA Kit Protocol

Modified from Ethan Ford's TruSeq RNA Kit protocol

http://ethanomics.files.wordpress.com/2012/02/truseq_rna_kit_protocol_third2.pdf

Perform End Repair

1. Use 13.3 μ l of End Repair Mix instead of 40 μ l.
2. Add 53.3 μ l of Ampure XP Beads instead of 160 μ l.
3. Place on magnetic rack for 5 minutes instead of 15 minutes.
4. Add 10 μ l Resuspension Buffer instead 17.5 μ l.
5. Incubate at room temperature for 5 minutes instead of 2 minutes.
6. Transfer 9 μ l of the supernatant to a new 0.2 ml PCR tube instead of 15 μ l.

Adenylate 3' Ends

1. Add 4.17 μ l A-Tailing Mix to sample instead of 12.5 μ l.

Ligate Adapters

1. Change recipe to : 0.83 μ l DNA Ligase Mix, 0.83 μ l Resuspension Buffer, 0.83 μ l Adapter Index .
2. Add 1.67 μ l Stop Ligase Mix instead of 5 μ l.
3. Add 20 μ l of AMPure XP beads instead of 42.5 μ l.
4. Place on magnetic rack for 5 minutes instead of 2 minutes.
5. Add 27 μ l Resuspension Buffer instead of 52.5 μ l.
6. Place in magnetic rack for 5 minutes instead of 2 minutes.
7. Transfer 26 μ l of the supernatant to a new 0.2 ml PCR tube instead of 50 μ l.
8. Add 26 μ l of well-mixed AMPure XP beads instead of 50 μ l.
9. Place on magnetic rack for 5 minutes instead of 2 minutes.

10. Add 25 μ l Resuspension Buffer instead of 22.5 μ l.
11. Place in magnetic rack for 5 minutes instead of 2 minutes.
12. Transfer 23 μ l of the supernatant to a new 0.2 ml PCR tube.

Proceed directly to Enrich DNA Fragments

1. Modify PCR mix as follows: 10 μ l of ligated DNA, 1.67 μ l PCR Primer Cocktail, 8.33 μ l PCR Master Mix.
2. Use 15-18 cycles of PCR instead of 10.

Return to Purify Ligation Products

1. Add 5 μ l of 4X loading buffer to each sample and to ladder instead of 7 μ l.