SIEMENS

Dimension Vista[®] System

Flex[®] reagent cartridge

PN 782423.001 2017-10-17 E

REF K6423A

PBNP

Dimension Vista® System

Dimension Vista® 500, Dimension Vista® 1000T, Dimension Vista® 1500, Dimension Vista® 3000T Systems

Flex® reagent cartridge

PBNP

See shaded sections: Updated information from 2017-04 version.

Issue Date 2017-10-17

N-terminal Pro-Brain Natriuretic Peptide

Intended Use: The PBNP method is an *in vitro* diagnostic assay for the quantitative measurement of N-terminal pro-brain natriuretic peptide (NT-proBNP) in human serum and plasma on the Dimension Vista® System. In individuals suspected of having congestive heart failure (CHF), measurements of NT-proBNP are used as an aid in the diagnosis and assessment of severity. The test is further indicated for the risk stratification of patients with acute coronary syndrome and heart failure.

Summary: Left ventricular dysfunction can occur as part of coronary heart disease, arterial hypertension, valvular disease, and primary myocardial disease. If the left ventricular dysfunction remains untreated and is progressive, the potential for mortality is high, e.g., due to sudden cardiac death.

Chronic cardiac insufficiency is a clinical syndrome caused by impairment of the cardiac pumping function. Based on the symptoms, the severity of cardiac insufficiency is classified in stages I-IV by the New York Heart Association (NYHA). Clinical tests and imaging procedures are used to diagnose left ventricular dysfunction.¹

The significance of natriuretic peptides in the control of cardiovascular system function has been established. Initial studies reveal that natriuretic peptides can be used for diagnosing clinical problems associated with left ventricular dysfunction.² The following natriuretic peptides have been reported in the literature: atrial natriuretic peptide (ANP), brain natriuretic peptide (BNP), and C-type natriuretic peptide (CNP).^{3,4}

ANP and BNP have natriuretic and diuretic properties. As antagonists of the renin-angiotensin-aldosterone system, they influence the electrolyte and fluid balance in an organism.^{5,6} In subjects with left ventricular dysfunction, serum and plasma concentrations of BNP increase, as do the concentrations of the biologically inactive prohormone, proBNP. ProBNP consists of 108 amino acids. It is secreted mainly by the left ventricle of the heart and, in this process, is cleaved into physiologically active BNP (77 – 108), and the N-terminal fragment NT-proBNP (1 – 76).⁴

Studies indicate that NT-proBNP can be used in diagnostic and prognostic applications.^{7 – 13} The concentration of NT-proBNP in plasma indicates the prognosis for the left ventricular dysfunction. It is also useful in assigning symptoms to cardiac or non-cardiac causes.

Determination of NT-proBNP helps to identify subjects with left ventricular dysfunction. Changes in NT-proBNP concentration can be used to evaluate the success of treatment in patients with left ventricular dysfunction. There are indications that NT-proBNP, due to its functions, is suitable for use in assessing vascular remodeling, and therefore contributes to the establishment of individualized rehabilitation procedures.^{10,11}

Fisher et al. found that congestive heart failure patients with elevated NT-proBNP values had a one year mortality rate of 53% compared to 11% in patients with lower values.¹² In the Gusto IV study which involved more than 6800 patients it was shown that NT-proBNP was the strongest independent predictor of one year mortality in patients with acute coronary syndrome.¹³

The European Society of Cardiology Task Force for the Diagnosis and Treatment of Chronic Heart Failure recommends in their guidelines that natriuretic peptides including NT-proBNP "may be most useful clinically as a rule out test due to consistent and very high negative predictive values."

NT-proBNP levels are increased in patients with unstable angina and following myocardial infarction.^{14,15} Although not diagnostic for these conditions, several studies indicate NT-proBNP measurements provide prognostic information for the short- and long-term risk stratification of patients with unstable angina or myocardial infarction.^{15,16,17}

Principles of Procedure: The PBNP method is a one-step sandwich chemiluminescent immunoassay based on LOCI® technology. LOCI® reagents include two synthetic bead reagents and a biotinylated monoclonal antibody fragment which recognize an epitope located in the N-terminal part of proBNP. The first bead reagent (Sensibeads) is coated with streptavidin and contains photosensitive dye. The second bead reagent (Chemibeads) is coated with a second antibody specific for a second independent epitope on NT-proBNP and contains chemiluminescent dye. Sample is incubated with Chemibeads and biotinylated antibody to form a particle/NT-proBNP/biotinylated antibody sandwich. Sensibeads then are added and bind to the biotin to form a bead-aggregated immunocomplex. Illumination of the complex by light at 680 nm generates singlet oxygen from Sensibeads, which diffuses to the Chemibeads and triggers a chemiluminescent reaction. The resulting chemiluminescent signal is measured at 612 nm and is directly proportional to the concentration of NT-proBNP in the sample.^{18,19,20}

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 2	Liquid	Biotinylated antibody	8 µg/mL	Sheep, monoclonal
3 – 4	Liquid	NT-proBNP chemibeads	150 µg/mL	Sheep, monoclonal
7 – 8	Liquid	Streptavidin Sensibeads	1400 µg/mL	Recombinant <i>E. coli</i>
9 – 12	Liquid	Assay buffer		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

Risk and Safety:

H317
P280, P272, P302 + P352, P333 + P313, P501

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. Contaminated work clothing should not be allowed out of the workplace. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Dispose of contents and container in accordance with all local, regional, and national regulations.

Contains: 5-chloro-2-methyl-3(2h)-isothiazolone mixture with 2-methyl-3(2h)-isothiazolone.

Safety data sheets (MSDS/SDS) available on siemens.com/healthcare

Precautions: Used LOCI® reaction vessels contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open well stability: 5 days for wells 1 – 12

Specimen Collection and Handling: Recommended specimen types: Serum or plasma samples (lithium heparin, sodium heparin, and EDTA).

Samples and controls stabilized with azide cannot be used.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.^{21,22,23} Each laboratory must determine its appropriate sample timing/ collection protocol.

Follow the instructions provided with your specimen collection device for use and processing.²³

Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.²⁴ Specimens should be free of particulate matter.

If clotting time is increased due to thrombolytic or anticoagulant therapy, the use of plasma specimens will allow for faster sample processing and reduce the risk of interference by particulate matter.^{21,24}

Samples are stable for three days when stored at 20 – 25 °C or at 2 – 8 °C or for 12 months when frozen at -20 °C or below.²⁵ Repetitive freezing and thawing of specimens should be avoided.

Ensure that patient samples, calibrators, and controls are equilibrated at ambient temperature (22 – 28 °C) before testing. Samples containing precipitates must be centrifuged before performing the assay.

Procedure

Materials Provided

PBNP Flex® reagent cartridge, Cat. No. K6423A

Materials Required But Not Provided

PBNP CAL, Cat. No. KC676A

Quality Control products such as:

- MAS® CardiImmune® proBNP Control / MAS® CardiImmune® TL/XL MoniTrol Cardiac TL
- CardiImmune® TL Liquid Assayed Cardiac Marker Control
- MAS® and CardiImmune® are registered trademarks of Medical Analysis Systems, Inc. Camarillo, CA.
- Bio-Rad Liquichek™ Cardiac Markers Controls
- Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, California 92618.
- Audit™ MicroFD™ BNP Control
- Audit™ and MicroFD™ are trademarks of Audit™ MicroControls, Inc. Las Vegas, NV 89102.

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista® System. For details of this processing, refer to your Dimension Vista® Operator's Guide.

Test Conditions

Sample Volume (delivered to the vessel)	8 µL
Chemibead Reagent Volume	20 µL
Biotinylated Antibody Reagent Volume	20 µL
Sensibead Volume	15 µL
Assay Buffer Volume	100 µL
Temperature	37.0 °C
Reaction time	10 minutes
Wavelength	680 and 612 nm
Type of Measurement	Chemiluminescence

Calibration

Calibration Material	PBNP Calibrator, Cat. No. KC676A
Calibration Scheme	5 levels, n = 3
Units	pg/mL
Typical Calibration Levels	Level 1 (Calibrator A): 0 pg/mL Level 2 (Calibrator B): 250 pg/mL Level 3 (Calibrator C): 1500 pg/mL Level 4 (Calibrator D): 12000 pg/mL Level 5 (Calibrator E): 36750 pg/mL
Calibration Frequency	Every 30 days for any one lot Calibration interval may be extended based on acceptable verification of calibration.
A new calibration is required:	<ul style="list-style-type: none"> • For each new lot of Flex® reagent cartridges • After major maintenance or service, if indicated by quality control results • As indicated in laboratory quality control procedures • When required by government regulations

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known NT-proBNP concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of PBNP in pg/mL using the calculation scheme described in your Dimension Vista® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 5 – 35000 pg/mL

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 35000 pg/mL should be repeated on dilution.
 - Manual Dilution: Dilute 5 fold with Reagent grade water to obtain results within reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.
 - Autodilution (AD): Instrument does a 1:5 dilution in the aliquot tray and the normal 8 µL sample is still used for over range samples. Refer to your Dimension Vista® Operator's Guide.
- Samples with results less than 5 pg/mL will be reported as "less than 5 pg/mL" by the instrument.

Limitations of Procedure

Concentrations of natriuretic peptides may be elevated in patients with acute myocardial infarction, patients that are candidates for renal dialysis, and patients that have undergone renal dialysis. Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This assay has been designed to minimize interference from heterophilic antibodies. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.^{26,27}

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in PBNP results. Refer to your Dimension Vista® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

Specimens that contain biotin at a concentration of 100 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to falsely depressed results for patient samples.

Results from patients taking biotin supplements or receiving high-dose biotin therapy should be interpreted with caution due to possible interference with this test.

Maximum Observed Repeatability

The expected maximum observed standard deviations for repeatability (within-run precision) using n = 5 replicates at the following analyte concentrations are:

PBNP concentration	Acceptable SD Maximum
125 pg/mL	6.5 pg/mL
450 pg/mL	21.4 pg/mL

A system malfunction may exist if the acceptable SD maximum is exceeded.

Expected Values

NT-proBNP concentrations in the Reference Group are shown in the following tables. Each laboratory should establish its own reference interval for the PBNP method as determined on the Dimension Vista® System. The recommended medical decision thresholds, by age group, are:

Patients < 75 years: 125 pg/mL
 Patients ≥ 75 years: 450 pg/mL

Reference Study Group

NT-proBNP concentrations were determined in 318 individuals without congestive heart failure (163 women and 155 men). This population included apparently healthy individuals and individuals with diabetes, hypertension, and pulmonary disease. The statistics for NT-proBNP concentrations in the reference study group are shown in the following table.

Reference Study Group

All	<55 yrs	55 – 64 yrs	65 – 74 yrs	≥ 75 yrs
Mean	45.9	54.6	175.5	451.3
SD	53.6	62.3	197.6	987.1
Median	27.2	28.8	106.2	173.5
95th Percentile	131.2	206.6	880.8	1615.7
% < 125 pg/mL	94%	80%	61%	-
% < 450 pg/mL	-	-	-	79%
N	163	15	18	122

Males	<55 yrs	55 – 64 yrs	65 – 74 yrs	≥ 75 yrs
Mean	39.0	56.7	219.4	520.8
SD	57.8	58.2	242.4	1080.9
Median	16.6	25.4	123.6	173.5
95th Percentile	168.8	154.7	880.8	1969.7
% < 125 pg/mL	93%	83%	55%	-
% < 450 pg/mL	-	-	-	73%
N	76	6	11	62

Females	<55 yrs	55 – 64 yrs	65 – 74 yrs	≥ 75 yrs
Mean	51.9	53.2	106.6	379.5
SD	49.2	68.3	60.8	883.1
Median	39.4	28.8	87.7	167.3
95th Percentile	124.0	206.6	215.2	1453.2
% < 125 pg/mL	95%	78%	71%	-
% < 450 pg/mL	-	-	-	85%
N	87	9	7	60

Disease Study Group

Blood samples were obtained from 269 patients diagnosed with congestive heart failure (CHF). The population included 78 women and 191 men. The descriptive statistics and New York Heart Association (NYHA) functional classes are provided below. Each laboratory should establish a reference range that represents the patient population that is to be evaluated. In addition, laboratories should be aware of their institution's current practice for the evaluation of CHF.

CHF Population – All

	<55 yrs	55 – 64 yrs	65 – 74 yrs	≥75 yrs
Mean	6131.8	4455.5	6168.1	8142.7
SD	13470.3	9296.0	10950.5	12699.1
Median	1550.9	1639.5	2541.3	3780.5
95th Percentile	28099.6	19481.5	17768.6	26923.1
% > 125 pg/mL	89%	92%	95%	-
% > 450 pg/mL	-	-	-	94%
N	72	72	59	66

CHF Population – Males

	<55 yrs	55 – 64 yrs	65 – 74 yrs	≥75 yrs
Mean	7771.4	4954	6537.2	8827.3
SD	15797.2	10536.7	12383.9	12856.6
Median	2526.1	1699.8	2541.3	5529.2
95th Percentile	31237.1	19481.5	17768.6	18634.8
% > 125 pg/mL	94%	94%	93%	-
% > 450 pg/mL	-	-	-	96%
N	48	51	43	49

CHF Population – Females

	<55 yrs	55 – 64 yrs	65 – 74 yrs	≥75 yrs
Mean	2852.6	3244.9	5176.1	6169.7
SD	5774.0	5200.7	5727.8	12397.4
Median	452.3	1516.9	2270.8	1125.3
95th Percentile	15569.3	7903.1	16889.9	45824.4
% > 125 pg/mL	79%	86%	100%	-
% > 450 pg/mL	-	-	-	88%
N	24	21	16	17

CHF Population – All

NYHA Functional Class	All CHF				
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	6184.5	3941.9	3648.8	6643.1	11089.2
SD	11737.8	8461.3	4897.0	13567.0	14428.2
Median	2282.2	1161.0	1091.7	2541.3	4537.7
5th Percentile	101.1	100.4	46.6	143.1	74.6
95th Percentile	22331.9	14184.1	14238.5	17497.1	28099.6
% > Cutoff	92%	93%	90%	94%	88%
Minimum	5.1	52.8	24.5	5.1	39.6
Maximum	91446.4	52327.0	18634.8	91446.4	78045.2
N	269	45	60	124	40

CHF Population – Males

NYHA Functional Class	All CHF				
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	7012.1	4567.5	4402.2	7298.4	14351.2
SD	12991.9	9204.3	5167.4	14789.5	17258.9
Median	3162.6	1639.5	2230.9	3442.2	6015.9
5th Percentile	126.1	136.3	46.5	126.1	101.1
95th Percentile	24004.7	19481.5	14238.5	17497.1	31237.1
% > Cutoff	94%	97%	93%	95%	91%
Minimum	5.1	72.5	44.7	5.1	74.6
Maximum	91446.4	52327.0	18634.8	91446.4	78045.2
N	191	37	40	91	23


CHF Population – Females

NYHA Functional Class	All CHF				
	All CHF	NYHA I	NYHA II	NYHA III	NYHA IV
Mean	4157.8	1048.3	2142.1	4836.0	6676.0
SD	7561.7	1636.1	4008.9	9366.5	7864.9
Median	950.2	225.4	705.3	1286.7	3780.5
5th Percentile	84.3	52.8	24.5	143.1	39.6
95th Percentile	22331.9	4960.2	11020.8	26923.1	22997.0
% > Cutoff	87%	75%	85%	94%	82%
Minimum	24.5	52.8	24.5	126.8	39.6
Maximum	45824.4	4960.2	15893.8	45824.4	22997.0
N	78	8	20	33	17


These results show that there is a relationship between the severity of the clinical signs and symptoms of CHF and the median NT-proBNP concentration, demonstrating that the Dimension Vista® PBNP method can be used as an aid in the diagnosis of all degrees of CHF severity including asymptomatic patients.

Interpretation of Results

The Receiver Operator Characteristics (ROC) Curve presents the clinical sensitivity and specificity at various cutoffs for the 269 patients diagnosed with CHF and 318 subjects without CHF. The ROC curve for the Dimension Vista® PBNP assay is shown below. The area under ROC curve (AUC) for the Dimension Vista® PBNP assay is 0.928 with a 95% confidence interval of 0.907 to 0.949.

A box and whiskers plot of the clinical study population is presented below. Recommended clinical thresholds are 125 pg/mL for patients younger than 75 years and 450 pg/mL for patients 75 years and older. Six disease group samples with values above the assay range are not displayed in the plot.

Sensitivity and Specificity vs. Age and Gender

The tables below show the clinical sensitivity and specificity of the Dimension Vista® PBNP assay using a cutoff of 125 pg/mL for patients younger than 75 years and 450 pg/mL for patients 75 years or older.

	Males	
	<75 yrs	≥75 yrs
Sensitivity	94% (133/142)	96% (47/49)
95% Confidence Interval	90 – 98	90 – 100
Specificity	87% (81/93)	73% (45/62)
95% Confidence Interval	80 – 94	61 – 84

	Females	
	<75 yrs	≥75 yrs
Sensitivity	89% (54/61)	88% (15/17)
95% Confidence Interval	81 – 97	73 – 100
Specificity	92% (95/103)	85% (51/60)
95% Confidence Interval	87 – 97	76 – 94

Age-matched ROC analysis

An age-matched ROC analysis of the clinical data was performed via the weighted method described in Kondratovich (2002).²⁹ The resulting AUC is 0.929 with a 95% confidence interval of 0.909 to 0.950.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista® System.

Precision^{29,c}

Material	Mean pg/mL	Standard Deviation (%CV)	
		Repeatability	Within-Lab
Serum Pool 1	120	1.72 (1.4)	2.06 (1.7)
Serum Pool 2	438	4.75 (1.1)	4.75 (1.1)
Serum Pool 3	880	9.05 (1.0)	12.24 (1.4)
Serum Pool 4	5075	69.06 (1.4)	91.97 (1.8)
Internal QC	32633	763.39 (2.3)	1140.27 (3.5)

c. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Method Comparison^{30,d}

Serum and Plasma Equivalency

Comparison testing results on the Dimension Vista® System for lithium heparin, sodium heparin and EDTA plasma samples versus serum samples is provided as follows:

Comparative Specimen	Slope	Intercept pg/mL	Correlation Coefficient	n
Li Heparin	1.01	-2.0	0.9984	53
Na Heparin	1.00	4.8	0.9987	57
EDTA	1.00	0.5	0.9997	59

d. Model equation by Passing-Bablok regression is: [results for Dimension Vista® System] = slope x [serum results] + intercept.

Regression Statistics^e

Comparative Method	Slope	Intercept pg/mL	Correlation Coefficient	n
Dimension Vista® PBNP (polyclonal)	1.02	-3.2	0.996	160 ^f

e. CLSI/NCCLS EP9-A2 was used. The method used to fit the regression line was Passing-Bablok Statistics.

f. The range of 160 values in the correlation study was 12 – 29519 pg/mL.

Specificity

HIL interference

The PBNP method was evaluated for interference to hemoglobin, icterus and lipemia (HIL) according to CLSI/NCCLS EP7-A2.³¹ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	NT-proBNP pg/mL	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer)	134	
	1000 mg/dL [mmol/L] ^a	2044	< 10
Bilirubin (unconjugated)	60 mg/dL [1026 μmol/L]	132	
		2019	< 10
Bilirubin (conjugated)	60 mg/dL [1026 μmol/L]	134	
		2066	< 10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	131	
		2028	< 10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* Analyte results should not be corrected based on this bias.

g. Système International d'Unités [SI units] are in brackets.

Non Interfering Substances

The following substances do not interfere with the PBNP method when present in serum at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at NT-proBNP concentration of 135 and 2000 pg/mL.

Compound	Concentration	SI Units
Abciximab	21 μg/mL	441 μmol/L
Acetaminophen	0.025 mg/dL	1.66 μmol/L
Allopurinol	2.5 mg/dL	184 μmol/L
Amikacin	15 mg/dL	256 μmol/L
Amiodarone	20 μg/mL	31 μmol/L
Amiodipine Besylate	4 μg/mL	7.1 μmol/L
Ampicillin	5.3 mg/dL	152 μmol/L
Ascorbic Acid	5 mg/dL	227 μmol/L
Atenolol	1.0 mg/dL	37.6 μmol/L
Atorvastatin	32 mg/dL	26.5 μmol/L
Caffeine	6 mg/dL	308 μmol/L
Captopril	5.0 mg/dL	230 μmol/L
Carbamazepine	3 mg/dL	127 μmol/L
Chloramphenicol	5 mg/dL	155 μmol/L
Chlordiazepoxide	1 mg/dL	33.3 μmol/L
Chlorpromazine	0.2 mg/dL	6.27 μmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 μmol/L
Cinnarizine	3.0 mg/dL	81.4 μmol/L
Clopidogrel bisulfate	30 μg/mL	71.4 μmol/L
Creatinine	30 mg/dL	2652 μmol/L
Cyclosporine A	4000 ng/mL	3.3 μmol/L
Dextran 40	6000 mg/dL	1500 μmol/L
Diazepam	0.5 mg/dL	17.6 μmol/L
Diclofenac	60 μg/mL	188.6 μmol/L
Digitoxin	60 ng/mL	78.6 nmol/L
Digoxin	5 ng/mL	6.15 nmol/L
Diltiazem	120 μg/mL	266.1 μmol/L
Dipyridamole	30 μg/mL	59.4 μmol/L
Disopyramide	40 μg/mL	117 μmol/L
Dopamine	65 mg/dL	3.4 mmol/L
Enalapril maleate	16 μg/mL	32.5 μmol/L
Erythromycin	6 mg/dL	81.6 μmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 μmol/L
Furosemide	6 mg/dL	181 μmol/L
Gentamicin	12 mg/dL	251 μmol/L

Compound	Concentration	SI Units
Heparin	3 U/mL	3000 U/L
Hydralazine	20 µg/mL	101.7 µmol/L
Hydrochlorothiazide	20 µg/mL	67.2 µmol/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Indomethacin	16 µg/mL	36.9 µmol/L
Isosorbide dinitrate	6.0 mg/dL	254 µmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lisinopryl	16 µg/mL	36.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Lovastatin	16 µg/mL	39.6 µmol/L
L-thyroxine	60 µg/dL	0.77 µmol/L
Methyldopa	2.5 mg/dL	118 µmol/L
Milrinone lactate	2.4 µg/mL	11.4 µmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nifedipine	6.0 mg/dL	173 µmol/L
Nitrofurantoin	40 µg/mL	168 µmol/L
Nitroglycerin	0.16 µg/mL	0.577 µmol/L
Oxazepam	12 µg/mL	41.8 µmol/L
Oxytetracycline	100 µg/mL	217.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Probenecid	200 µg/mL	700.9 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Propranolol	0.15 mg/dL	5.1 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	11.1 g/dL	111 g/L
Quinidine	20 µg/mL	61.6 µmol/L
Rheumatoid Factor	500 IU/mL	500 IU/mL
Salicylic Acid	60 mg/dL	4.34 mmol/L
Simvastatin	32 µg/mL	76.5 µmol/L
Sulfamethoxazole	320 µg/mL	1263.3 µmol/L
Theophylline	4 mg/dL	222 µmol/L
Trimethoprim	64 µg/mL	220.7 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Verapamil	16 mg/dL	0.33 µmol/L
Warfarin	4 µg/mL	130 µmol/L

Hook effect

One step sandwich immunoassays are susceptible to a high dose "hook effect", where an excess of antigen prevents simultaneous binding of the capture and detection antibodies to a single analyte molecule.³² Such samples must be diluted and reassayed prior to reporting the results. The PBNP method shows no hook effect up to 400,000 pg/mL.

Cross-reactivity

The following substances have no significant cross-reactivity (less than 1%) at the concentrations indicated when added to samples containing 0 and approximately 125 pg/mL NT-proBNP:

Substance	Concentration
ANP ₂₈	3.1 µg/mL
NT-proANP ₁₋₃₀ (preproANP ₂₆₋₅₅)	3.5 µg/mL
NT-proANP ₃₁₋₆₇ (preproANP ₅₆₋₉₂)	1.0 ng/mL
NT-proANP ₇₉₋₉₈ (preproANP ₁₀₄₋₁₂₃)	1.0 ng/mL
BNP ₃₂ (Natreacor®)	3.5 µg/mL
CNP ₃₂	2.2 µg/mL
DNP	1.0 ng/mL
VNP	1.0 ng/mL
Adrenomedullin	1.0 ng/mL
Aldosterone	0.6 ng/mL
Angiotensin I	0.6 ng/mL
Angiotensin II	0.6 ng/mL
Angiotensin III	1.0 ng/mL
Endothelin	20 pg/mL
Renin	50 ng/mL
Urodilatin	3.5 µg/mL
Arg-Vasopressin	1.0 ng/mL

Analytical Sensitivity: ≤ 5 pg/mL

The analytical sensitivity of the PBNP method is ≤ 5 pg/mL. Analytical sensitivity is defined as the concentration at two standard deviations (n = 20) above a sample devoid of NT-proBNP such as the PBNP Level 1 Calibrator (0 pg/mL).

Limit of Detection: 0.8 pg/mL

The Limit of Detection (LoD) for ANALYTE is 0.8 pg/mL, determined consistent with CLSI guideline EP17-A and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 60 determinations, with 5 blank and 5 low level samples. LOD is the lowest concentration that can be detected reliably.

Functional Sensitivity: Evaluation of the 20% CV limit

Functional sensitivity was evaluated by determining the total imprecision of natural NT-proBNP serum samples. Total imprecision was measured by a two replicate per day, twenty day ANOVA study. The functional sensitivity/ imprecision profile for the Dimension Vista® PBNP method corresponds to a coefficient of variation (CV) of 20% at an NT-proBNP level of ≤ 30 pg/mL.

Dilutional Linearity

A patient plasma sample with a high NT-proBNP concentration was diluted across the assay range. Additional plasma samples were diluted linearly to confirm dilutional accuracy at the low end of the assay range. All samples were diluted with low level patient plasma pools. The following table summarizes the results of these studies:

	Sample A	Sample B
Starting Concentration (pg/mL)	41812	1266
Regression Statistics		
Slope	1.004	1.003
Intercept	-55.8	-3.5
Correlation Coefficient	0.999	0.999
N	5	5
% Recovery		
Average	99.9	99.5
Range	98 – 101	98 – 100

Manufactured under license from Roche Diagnostics.

US 5786163 and foreign equivalents

US Patent Application: 09/890442 and foreign equivalents

Symbols Key: See Adjacent Panel.

Bibliography: See Adjacent Panel.

Dimension Vista®, LOCI®, Dimension®, and Flex® are trademarks of Siemens Healthcare Diagnostics.

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©2008 Siemens Healthcare Diagnostics

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Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía:

1. Remme WJ, Swedberg K. Guidelines for the diagnosis and treatment of chronic heart failure, Task Force for the Diagnosis and Treatment of Chronic Heart Failure, European Society of Cardiology. Eur Heart J 2001; 22: 1527- 1560.
2. Richards AM, Nicholls GM, Yandle TG, Frampton C, Espiner EA, Turner JG, et al. Plasma N-terminal pro-brain natriuretic peptide and adrenomedullin: New neurohormonal predictors of left ventricular function and prognosis after myocardial infarction. Circulation 1998; 97: 1921 – 1929.
3. de Bold AJ. Atrial natriuretic factor: A hormone produced by the heart. Science 1985; 230:767 – 770.
4. Valli N, Gobinet A, Bordernave L. Review of 10 years of the clinical use of brain natriuretic peptide in cardiology. J Lab Clin Med 1999; 134: 437-444.
5. de Bold AJ, Boerenstein HB, Varess AT, Sonnenberg H. A rapid and potent natriuretic response to intravenous injection of atrial extracts in rats. Life Sci 1981; 28: 89-94.
6. Epstein M, Loutzenhiser R, Friedland E, Aceto RM, Camargo MJF, Atlas SA. Relationship of increased plasma natriuretic factor and renal sodium handling during immersion-induced central hypervolemia in normal humans. J Clin Invest 1987; 79: 738-745.
7. Struthers AD. How to use natriuretic peptide levels for diagnosis and prognosis. The European Society of Cardiology. Eur Heart J 1999;20:1374-1375.
8. Hunt PJ, Richards AM, Nicholls MG, Yandle TG, Doughty RN, Espiner EA. Immunoreactive aminoterminal pro-brain natriuretic peptide (NT-proBNP): A new marker for cardiac impairment. Clin Endocrinol 1997; 47:287-296.
9. Talwar S, Squire IB, Davies JE, Barnett DB, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and the ECG in the assessment of left ventricular systolic dysfunction in a high risk population. Eur Heart J 1999; 20: 1736-1744.
10. Darbar D, Davidson NC, Gillespie N, Choy A-MJ, Lang CC, Shyr Y, et al. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. A J Cardiol 1996; 78: 284-287.
11. McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, et al. Biochemical detection of left-ventricular systolic dysfunction. Lancet 1998; 351: 9-13.
12. Fisher C et al. NT-proBNP predicts prognosis in patients with chronic heart failure. Heart 2003; 89:879-881.
13. James SK et al. NT-proBNP and other risk markers for the separate prediction of mortality and subsequent myocardial infarction in patients with unstable coronary artery disease. GUSTO IV Substudy. Circulation 2003; 108:275-281.
14. Talwar S, Squire I, Downie P, Davies J, Ng LL. Plasma N-terminal pro-brain natriuretic peptide and cardiotrophin 1 are raised in unstable angina. Heart 2000; 84: 421-424.
15. Jernberg T, Stridsberg M, Venge P, Lindahl B. N-Terminal pro-brain natriuretic peptide on admission for early risk stratification of patients with chest pain and no ST-segment elevation. J Am Coll Cardiol 2002; 40: 437-445.
16. Zeller M, Cottin Y, Laurent Y, Danchin N, L'Huillier I, et al. N-Terminal pro-brain natriuretic peptide levels in patients with non-ST-elevation myocardial infarction. Cardiology 2004; 102: 37-40.
17. Richards AM, Nicholls MG, Espiner E, Lainchbury J, et al. B-type natriuretic peptides and ejection fraction for prognosis after myocardial infarction. Circulation 2003; 107: 2786-2792.
18. Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et al. Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogeneous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
19. Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et al. Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl.Acad.Sci. USA, Vol 91, pp.5426-5430, June 1994 Biochemistry.
20. Ullman EF. Homogeneous Immunoassays, The Immunoassay Handbook, 2nd ed., David Wild editor, 2001 192-194.
21. Burtis CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Biology, Fourth Edition, Elsevier Saunders, St. Louis, MO 2006; pp. 41-58 (specimen collection).
22. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
23. Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
24. Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
25. Elecsys® proBNP Assay Package Insert Sheet, Roche Diagnostics GmbH, December 2003.
26. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. Clin Chem 1992; 38: 1737-1742.
27. Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7; 942-956.
28. Kondratovich, M (2002). Matched Receiver Operating Characteristic (ROC) analysis and propensity scores. Proceedings of the 2002 Joint Statistical Meeting, Biopharmaceutical Section, New York, NY.
29. Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline- Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA, 19087-1898 USA, 2004.
30. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline- Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
31. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
32. Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the "hook effect" in two-site immunoradiometric assays. Anal Biochem 982; 127: 308 – 315.

**Symbols Key
Symbolschlüssel
Explication des Symboles
Interpretazione simboli
Clave de los Símbolos**

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleiddokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabbicante / Fabricante
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
	Level / Konzentration / Niveau / Livello / Nivel

2014-10_EFRS

Made in: USA
Siemens Healthcare Diagnostics Inc.
500 GBC Drive
Newark, DE 19714 U.S.A.
siemens.com/healthcare

Siemens Healthcare Diagnostics Ltd.
Sir William Siemens Sq.
Frimley, Camberley, UK GU16 8QD

**Global Siemens
Headquarters**
Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

**Global Siemens
Healthcare Headquarters**
Siemens AG
Healthcare
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens.com/healthcare

Global Division
Siemens Healthcare
Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5097
USA
siemens.com/healthcare

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All rights reserved.