DADE BEHRING
Dimension Vista® System
Flex® reagent cartridge

REF K7046

hsCRP

K7046G00E0841 2008-02

DADE BEHRING

Dimension Vista[®] System

Flex[®] reagent cartridge

hsCRP

Issue Date 2008-02

CardioPhase[®] high sensitivity CRP

Intended Use: The CardioPhase[®] high sensitivity CRP (hsCRP) method is an *in vitro* diagnostic test for the quantitative measurement of C-reactive protein (CRP) in human serum and heparinized plasma on the Dimension Vista[®] System.

High sensitivity CRP measurements may be used for evaluation of conditions thought to be associated with inflammation, in otherwise healthy individuals and as an independent risk marker for the identification and stratification of individuals at risk for future cardiovascular disease*. Measurements of hsCRP, when used in conjunction with traditional clinical laboratory evaluation of acute coronary syndromes, may be useful as an independent marker of prognosis for recurrent events, in patients with stable coronary disease or acute coronary syndromes.

Summary: CRP is one of the 'acute-phase' proteins, the serum or plasma levels of which rise during general, nonspecific response to infectious and non-infectious inflammatory processes. CRP is synthesized in the liver and is normally present as a trace constituent of serum or plasma.

In various disease states resulting in tissue injury, infection or acute inflammation, CRP values may rise above normal to 20 to 500 mg/L¹. As elevated CRP values are always associated with pathological changes, the CRP method provides useful information for the diagnosis, therapy and monitoring of inflammatory processes and associated diseases^{3,4}. Increases in CRP values are non-specific and should not be interpreted without a complete clinical history. Studies have shown that measurement of CRP by high sensitivity methods is a strong independent predictor of risk for future cardiovascular and peripheral vascular disease⁵. High sensitivity CRP measurements have also been shown to add to the predictive value of other markers used to assess the risk of cardiovascular and peripheral vascular disease^{6,7,8}. Elevated CRP values determined by high sensitivity CRP methods may be indicative of the prognosis of individuals with acute coronary syndromes, and may be useful in the management of such individuals⁵.

Principles of Procedure: Polystyrene particles coated with monoclonal antibodies specific to human CRP are aggregated when mixed with samples containing CRP. These aggregates scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the respective protein in the sample. The result is evaluated by comparison with a standard of known concentration.

Reagents

Wells ^{a,b}	Form	Ingredient	Concentration ^c	Source
1 - 8	Liquid	hsCRP Supplement Reagent: Phosphate buffer; Polidocanol	1.9 g/L	
9 - 12	Liquid	High sensitivity CRP Reagent: Polystyrene particles; Monoclonal antibodies	1 g/L 13 mg/L	Mouse

a) Wells are numbered consecutively from the wide end of the cartridge.

b) Contain preservatives.

c) Nominal value per well in a cartridge.

Precautions:

Contains sodium azide (< 0.1 %) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Contains human source material.

Each donor or donor unit was tested and found to be negative for human immunodeficiency virus (HIV) 1 and 2, hepatitis B virus (HBV) and hepatitis C virus (HCV) using either tests found to be in conformance with the In Vitro Diagnostic Directive in the EU or FDA approved tests. Because no known test can offer complete assurance of the absence of infectious agents, all human derived products should be handled with appropriate caution.

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

Reagent Preparation:

All reagents are liquid and ready to use.

Store at: 2 - 8 °C.

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 90 days.

Open well stability: 21 days for wells 1 - 12

Specimen Collection and Handling:

Recommended specimen types: serum or heparinized plasma.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture⁹.

Follow the instructions provided with your specimen collection device for use and processing¹⁰.

For serum, complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection¹¹.

Samples should be as fresh as possible (stored for no more than seven days at 2 - 8 °C) or stored frozen. Samples can be stored at below -20 °C for up to eight months¹², if they are frozen within 24 hours after collection and if repeated freeze-thaw cycles are avoided. Lipemic or frozen samples, which become turbid after thawing, must be clarified by centrifugation (10 minutes at approximately 15,000 x g) prior to testing. Specimens should be free of particulate matter.

Procedure

Materials Provided

CardioPhase[®] hsCRP Flex[®] reagent cartridge, Cat. No. K7046

Materials Required But Not Provided

PROT2 CAL, Cat. No. KC780
System Diluent, Cat. No. KS804
Quality Control material, such as:
hsCRP CON L, Cat. No. KC755
hsCRP CON H, Cat. No. KC757

Test Steps

Sampling, reagent delivery, mixing, and processing are automatically performed by the Dimension Vista[®] System. For details of this processing, refer to your Dimension Vista[®] Operator's Guide.

	Test Conditions
Initial Sample Dilution	1:1 (neat)
Sample Volume (delivered to the cuvette)	Cuvette 1.37 µL
Diluent Volume	118 µL
High Sensitivity CRP Reagent	27.3 µL
Temperature	37.0 °C
Reaction time	5 minutes 50 seconds
Wavelength	840 nm
Type of Measurement	Nephelometric

	Calibration
Calibration Material	PROT2 CAL, Cat. No. KC780
Calibration Scheme	7 levels, n = 3
Units	mg/dL [mg/L] ^a (mg/dL x 0.1) = [mg/L]
Typical Calibration Levels	0.012, 0.025, 0.056, 0.111, 0.222, 0.50, 1.0 mg/dL [0.12, 0.25, 0.56, 1.11, 2.22, 5.0, 10.0 mg/L]

To obtain calibrator levels that span the measuring range, PROT2 CAL is diluted automatically with System Diluent by the instrument to the following dilutions:

Level 1: 1:162 dilution	
Level 2: 1:81 dilution	
Level 3: 1:36 dilution	
Level 4: 1:18 dilution	
Level 5: 1:9 dilution	
Level 6: 1:4 dilution	
Level 7: 1:2 dilution	
Calibration Frequency	Every 45 days for any one lot
Calibration interval may be extended based on acceptable verification of calibration	

A new calibration is required:

- For each new lot of Flex[®] reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

d) Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze a minimum of two levels of a Quality Control (QC) material with known C-reactive protein concentrations, e. g. hsCRP CON L or H.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of C-reactive protein in mg/dL [mg/L] using the calculation scheme described in your Dimension Vista[®] Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.016 - 0.95 mg/dL [0.16 - 9.50 mg/L]

This is the measuring range for undiluted samples that are automatically processed by the instrument. If the readings obtained are outside the initial measuring range, the method can be repeated using a higher dilution of the sample.

Refer to your Dimension Vista[®] Operator's Guide for information on repeat measurements using other dilutions.

- Samples with results in excess of 0.95 mg/dL [9.50 mg/L] can be repeated on a higher dilution.
- Samples with results less than 0.016 mg/dL [0.16 mg/L] will be reported as "less than 0.016 mg/dL" by the instrument.

Limitations of Procedure

Turbidity and particles in the samples may interfere with the determination. Therefore, samples containing particles must be centrifuged prior to testing. Lipemic or turbid samples, which cannot be clarified by centrifugation (10 minutes at approximately 15,000 x g), must not be used.

Due to matrix effects, inter-laboratory survey samples and control samples may yield results that differ from those obtained with other methods. It may therefore be necessary to assess these results in relation to method-specific target values.

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in CardioPhase[®] hsCRP results. Refer to your Dimension Vista[®] Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported. Patient samples may contain heterophilic antibodies that could react in immunoassays to give falsely elevated or depressed results. This method has been designed to minimize interference from heterophilic antibodies¹³. Nevertheless, complete elimination of this interference from all patient specimens cannot be guaranteed. A test result that is inconsistent with the clinical picture and patient history should be interpreted with caution.

AHA/CDC Expert Panel Recommendations¹⁴:

The AHA/CDC Expert Panel has the following recommendations for limitations on the use of hsCRP assays:

- hsCRP levels should not be substituted for assessment of traditional cardiovascular risk factors.
- Application of management guidelines for acute coronary syndromes should not be dependent on hsCRP levels.
- When using the assay for risk assessment, patients with persistently unexplained, marked elevation of hsCRP (> 10 mg/L) after repeated testing should be evaluated for non-cardiovascular etiologies.
- The expert panel recommends against screening of the entire adult population for hsCRP as a public health measure.
- Patients with evidence of active infection, systemic inflammatory processes or trauma should not be tested for cardiovascular disease risk assessment until these conditions have abated.
- Application of secondary prevention measures should not depend on hsCRP determination, but rather an array of risk factors (global risk assessment).
- Serial measurements of hsCRP should not be used to monitor effects of treatment.
- Two separate hsCRP measurements (optimally two weeks apart) should be obtained before performing risk assessment, due to within-subject hsCRP variability.


Interpretation of Results

Recent medical events resulting in tissue injury, infection or inflammation, which may cause elevated CRP levels should also be considered when interpreting results.

Cardiovascular and Peripheral Vascular Disease

Risk Assessment:

Several studies have examined the utility of hsCRP measurements for cardiovascular disease risk prediction. In a subset consisting of 280 individuals participating in the Physician's Health Study, baseline samples were tested for a number of biomarkers, including hsCRP. The population was followed for nine years to assess the development of peripheral arterial disease (PAD). The following graph shows the relative risk (95 % confidence interval) of developing future PAD as a function of quartile of hsCRP¹³.


The AHA/CDC Scientific Statement⁶ provides the following risk assessment guidelines:

Risk	hsCRP (mg/L)
Low	< 1.0
Average	1.0 - 3.0
High	> 3.0

Risk Stratification:

For patients with acute coronary syndromes, measurement of hsCRP may provide prognostic information. A value of CRP >10 mg/L in the early period (6 - 24 hours after onset of symptoms), has been shown to be indicative of an increased risk for short term (30 days - 1 year) recurrent cardiac events⁵.

A sub-study of 447 patients in the CAPTURE trial examined the clinical implications of elevated levels of CRP for risk stratification in patients with unstable angina. As shown in the following graph, patients with a CRP >10 mg/L experienced a higher event rate (mortality or MI) than patients with a CRP <10 mg/L¹⁴.

Expected Values: less than 0.30 mg/dL [3.00 mg/L]

This reference interval applies to serum samples from healthy adults^{15,16}.

The reference intervals are affected by many factors that may differ for each population studied. Each laboratory should establish its own expected values for *CardioPhase*[®] hsCRP as performed on the Dimension Vista[®] System.

Maximum Observed Repeatability

The expected maximum observed standard deviations (SD) for repeatability (within-run precision) using n = 5 replicates at the following nominal CRP concentrations are:

CRP Concentration	Acceptable SD Maximum
0.10 mg/dL [1 mg/L]	0.03 mg/dL [0.26 mg/L]
0.30 mg/dL [3 mg/L]	0.06 mg/dL [0.64 mg/L]

A system malfunction may exist if the acceptable SD maximum is exceeded.

Specific Performance Characteristics

The following data represent typical performance for the Dimension Vista[®] System.

Precision^{17,a}

Material	Mean		Standard Deviation mg/dL [mg/L] (% CV)					
	mg/dL	[mg/L]	Repeatability		Within-Lab			
hsCRP CON L	0.115	[1.15]	0.006	[0.06]	(4.8)	0.006	[0.06]	(5.4)
hsCRP CON H	0.346	[3.46]	0.014	[0.14]	(4.0)	0.015	[0.15]	(4.4)
Serum pool	0.239	[2.39]	0.013	[0.13]	(5.2)	0.013	[0.13]	(5.2)
Serum pool	0.650	[6.50]	0.033	[0.33]	(5.0)	0.037	[0.37]	(5.7)
Serum pool	0.772	[7.72]	0.036	[0.36]	(4.6)	0.042	[0.42]	(5.4)

e) CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Method Comparison¹⁸

Comparative Method	Regression Statistics ^d		Correlation Coefficient	n
	Slope	Intercept mg/dL [mg/L]		
<i>CardioPhase</i> [®] /hsCRP on the BN ProSpec [®] System	1.044	0.0003 [0.003]	0.995	133 ^f

f) CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was Passing Bablok.
g) The range of CRP values in the correlation study was 0.0169 mg/dL to 0.8922 mg/dL [0.169 mg/L to 8.922 mg/L].

Specificity

HIL Interference

The *CardioPhase*[®] hsCRP method was evaluated for interference according to CLSI/NCCLS EP7-A2¹⁹. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10 % is considered interference.

Substance Tested	Substance Concentration	CRP mg/dL [mg/L]	Bias* %
Hemoglobin (hemolysate)	1000 mg/dL [0.155 mmol/L]	0.093 [0.93]	+1
Bilirubin (unconjugated)	60 mg/dL [1026 µmol/L]	0.097 [0.97]	-3
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	0.097 [0.97]	+1
Lipemia	Refer to "Specimen Collection and Handling" section		

* Analyte results should not be corrected based on this bias.

Non Interfering Substances

The following substances do not interfere with the *CardioPhase*[®] hsCRP method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10 % at CRP concentration of 0.023 mg/dL to 0.591 mg/dL [0.23 mg/L to 5.91 mg/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ammonium heparin	3 U/mL	3000 U/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	5 ng/mL	6.15 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium chloride	2.3 mg/dL	3.2 mmol/L
Lithium heparin	3 U/mL	3000 U/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein, Albumin	6 g/dL	60 g/L
Rheumatoid Factors	500 IU/mL	500 IU/mL
Salicylic acid	60 mg/dL	4.34 mmol/L
Sodium heparin	3 U/mL	3000 U/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric acid	20 mg/dL	1190 µmol/L
Valproic acid	50 mg/dL	3467 µmol/L

The *CardioPhase*[®] hsCRP method shows no hook effect up to 130.22 mg/dL [1302.2 mg/L].

Recovery

Recovery of protein reference material ERM[®]-DA470 (CRM 470) ranged from 96.6 - 104.1 % with a mean recovery of 99.2 %.

Limit of Quantitation: 0.016 mg/dL [0.16 mg/L]

The limit of quantitation represents the lower limit of the reportable range for CRP.

Bibliography: See adjacent panel.

Symbols Key: See adjacent panel.

Dimension Vista[®], Dimension[®], Flex[®] and BN ProSpec[®] are trademarks of Siemens Healthcare Diagnostics.

CardioPhase[®] is a registered trademark in the USA, in Germany and other countries.

ERM[®] is the registered trademark for European Reference Materials.

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Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía

1. Kushner I, Rzewnicki DL. The acute phase response: general aspects. *Baillieres Clin Rheumatol* 1994;8:513-30.
2. Dixon JS, Bird HA, Sitton NG, et al. C-reactive protein in the serial assessment of disease activity in rheumatoid arthritis. *Scand J Rheum* 1984;13:39-44.
3. Hind CRK, Pepys MB. The role of serum C-reactive (CRP) measurement in clinical practice. *Int Med* 1984;5:112-51.
4. van Leeuwen M, van Rijswijk MH. Acute phase proteins in monitoring of inflammatory disorders. *Baillieres Clin Rheumatol* 1994;8:531-52.
5. Pearson TA, Mensah GA, Alexander RW, et al. Markers of Inflammation and Cardiovascular Disease: application to clinical and public health practice: A statement for healthcare professionals from the Centers for Disease Control and Prevention and the American Heart Association. *Circulation* 2003;107:499-511.
6. Ridker PM, Cushman M, Stampfer MJ, et al. Plasma concentration of C-Reactive Protein and risk of developing peripheral vascular Disease. *Circulation* 1998;97:425-8.
7. Ridker PM, Glynn RJ, Hennekens CH. C-reactive protein adds to the predictive value of total and HDL cholesterol in determining risk of first myocardial infarction. *Circulation* 1998;97:2007-11.
8. Ridker PM, Rifai N, Rose L, et al. Comparison of C-reactive protein and low-density lipoprotein cholesterol levels in the prediction of first cardiovascular events. *N Engl J Med* 2002;347:1557-65.
9. Clinical and Laboratory Standards Institute (CLSI). Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition. CLSI document H3-A5 (ISBN 1-56238-515-1). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
10. Clinical and Laboratory Standards Institute (CLSI). Tubes and Additives for Venous Blood Specimen Collection; Approved Standard-Fifth Edition. CLSI document H1-A5 (ISBN 1-56238-519-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
11. Clinical and Laboratory Standards Institute (CLSI). Procedures for the Handling and Processing of Blood Specimens; Approved Guideline-Third Edition. CLSI document H18-A3 (ISBN 1-56238-555-0). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
12. Kricka LJ. Human anti-animal antibody interference in immunological assays. *ClinChem* 1999;45:942-956.
13. Graph adapted from Ridker PM, Stampfer MJ, Rifai N. Novel risk factors for systemic atherosclerosis: a comparison of C-reactive protein, fibrinogen, homocysteine, lipoprotein(a), and standard cholesterol screening as predictors of peripheral arterial disease. *JAMA* 2001;285:2481-5.
14. Graph adapted from Heeschen C, Hamm C, Bruemmer J, et al. Predictive value of C-reactive protein and troponin T in patients with unstable angina: A comparative analysis. *J Am Coll Cardiol* 2000;35:1535-42.
15. Tracy RP, Lemaitre RN, Psaty BM, et al. Relationship of C-reactive protein to risk of cardiovascular disease in the elderly. Results from the Cardiovascular Health Study and the Rural Health Promotion Project. *Arterioscler Thromb Vasc Biol* 1997;17:1121-7.
16. Macy EM, Hayes TE, Tracy RP. Variability in the measurement of C-reactive protein in healthy subjects: implications for reference intervals and epidemiological applications. *Clin Chem* 1997;43:52-8.
17. Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 (ISBN 1-56238-542-9). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
18. Clinical and Laboratory Standards Institute (CLSI). Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition. CLSI document EP9-A2 (ISBN 1-56238-472-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
19. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition. CLSI document EP7-A2 (ISBN 1-56238-584-4). CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
20. Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. *Clin Chim Acta* 2005;355:S407.

**Symbols Key
 Symbolschlüssel
 Explication des Symboles
 Interpretazione simboli
 Clave de los Símbolos**

	Do not reuse / Nicht zur Wiederverwendung / Ne pas réutiliser / Non riutilizzare / No reutilizar
EXP CCYY-MM-DD 	Use By / Verwendbar bis / Utiliser jusque / Utilizzare entro / Fecha de caducidad
LOT	Batch Code / Chargenbezeichnung / Code du lot / Codice del lotto / Código de lote
REF	Catalogue Number / Bestellnummer / Référence du catalogue / Numero di catalogo / Número de catálogo
	Caution, consult accompanying documents / Achtung, Begleitdokumente beachten / Attention voir notice d'instructions / Attenzione, vedere le istruzioni per l'uso / Atención, ver instrucciones de uso
	Manufacturer / Hersteller / Fabricant / Fabricante / Fabricante
EC REP	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandataire dans la Communauté européenne / Mandatario nella Comunità Europea / Representante autorizado en la Comunidad Europea
	Contains sufficient for <n> tests / Inhalt ausreichend für <n> Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para <n> ensayos
IVD	In Vitro Diagnostic Medical Device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Producto sanitario para diagnóstico in vitro
	Temperature Limitation / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Limite de temperatura
	Consult Instructions for Use / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso
	Non-sterile / Nicht steril / Non stérile / Non sterile / No estéril
CE	CE Mark / CE Zeichen / Marquage CE / Marchio CE / Marca CE
CONTENTS	Contents / Inhalt / Contenu / Contenuto / Contenido
	Reconstitution Volume / Rekonstitutionsvolumen / Volume de reconstitution / Volume di ricostituzione / Volumen de reconstitución
LEVEL	Level / Konzentration / Niveau / Livello / Nivel

2006-06_EFIGS