EAA™

CE

ENDOTOXIN ACTIVITY ASSAY

For In Vitro Diagnostic Use Only
Pour diagnostic in vitro uniquement
Nur Für Die In-vitro-Diagnostik
Pour diagnostic in vitro uniquement
Только для диагностики In Vitro
Yalnızca In Vitro Teşhis Kullanımı İçindir

EN

FR

DE

IT

RU

TR

SPECTRAL | **UN**
Medical

Endotoxin Activity Assay (EAA™)

For In Vitro Diagnostic Use Only

EN

1. Intended Use

The EAA™ is a rapid in vitro diagnostic test that utilizes a specific monoclonal antibody to measure the endotoxin activity in EDTA whole blood specimens.

The information, when used in conjunction with microbial cultures and other relevant diagnostic tests (i.e., ultrasound, bronchoscopy or CAT scan) and clinical information (clinical signs, and/or radiographic test results) can:

- Aid in the risk assessment of patients for progression to severe sepsis.
- Direct the use of specific anti-endotoxin therapy

2. Endotoxin Activity Assay – Introduction

The EAA test assesses the level of endotoxin activity in the blood specimen of a specific patient in a semi-quantitative way, distinguishing 3 groups – low, intermediate and high endotoxin activity.

Endotoxin or lipopolysaccharide (LPS) is a major cell wall constituent of Gram-negative bacteria and the primary Gram-negative bacterial product responsible for severe sepsis¹⁻². Elevated endotoxin levels cause changes in the expression of more than 300 genes and activated macrophages, neutrophils, endothelial cells and the coagulation cascade producing the sepsis cascade³.

This assay uses the biological response of the neutrophils in a patient's blood to an immunological complex of endotoxin and exogenous antibody as a measure of the endotoxin activity in the patient⁴.

The assay reacts specifically with the Lipid A moiety of LPS of Gram-negative bacteria and does not cross-react with cell wall constituents of Gram-positive bacteria and other microorganisms.

3. Contents

3.1 EAA™ Kit

Name	Quantity
EAAST-20	1
EAA™ Trays	20
EAA™ Quality Control Trays	1
EAA™ Solution Vials	21
Instructions for Use	1

3.2 EAA™-Trays

Name	Quantity	Description
Tube 1	2	Stabilizers and luminol-zymosan
Tube 2	2	Murine monoclonal anti-endotoxin antibody, stabilizers and luminol-zymosan
Tube 3	2	Murine monoclonal anti-endotoxin antibody, stabilizers and luminol-zymosan
LPS Max Tube	1	Endotoxin (E. coli O55:B5) with stabilizers
Aliquot Tube	1	Empty tube

3.3 EAA™ Quality Control Trays

Name	Quantity	Description
Tube 1	2	Stabilizers and luminol-zymosan
Low Control Tube	1	Murine monoclonal antibody, stabilizers and luminol-zymosan
High Control Tube	1	Murine monoclonal anti-endotoxin antibody, stabilizers and luminol-zymosan
Tube 3	1	Murine monoclonal anti-endotoxin antibody, stabilizers and luminol-zymosan
QC LPS Max Tube	1	Endotoxin (E. coli "O55:B5") with stabilizers
QC Aliquot Tube	1	Empty tube

EN

3.4 EAA™ Solution

Name	Quantity	Description
EAA™ Solution	8 mL/Vial	HBSS and Heparin Solution

3.5 Accessories

Not supplied with the kit.

- Combipipettor
- 500 µL micropipette
- Sterile combipipette tips capable of delivering 40uL and 1000uL
- 1000uL micropipette tips
- EDTA blood collection tubes
- Timer
- Luminometer (Berthold Smartline TL Tube or equivalent)
- Incubating shaker capable of incubating samples at 37°C and agitating at 100 rpm
- Vortexer

4. Precautions

- Use general laboratory precautions when handling and storing the unopened assay materials.
- The EAA™ contains lyophilized beaded reagents in borosilicate glass test tubes. Under certain conditions (damaged packaging, improper sealing of the pouches or improper storage) the beads may shrink in the presence of excess humidity. Do not use tubes that contain shrunken beads as test results may be inaccurate.
- For manual pipetting of samples and controls, use individual pipette tips to eliminate carryover of residual.
- Do not use after the expiry date specified on the label.
- Use recommended universal precautions⁵ for handling used solutions and human specimens.

- Handle or dispose of used EAA™ solution and all human blood products using universal precautions.
- Do not pipette solutions by mouth.
- Clean up spills immediately with a 0.5% sodium hypochlorite solution.

5. Principle

A basal activity measurement (Tube 1) in the absence of the specific anti-endotoxin antibody measures the non-specific oxidative burst of the patient's neutrophils. An additional control measurement including the specific anti-endotoxin antibody and an excess of exogenous endotoxin (Tube 3) measures the maximum oxidative burst of the patient's neutrophils. The test measurement (Tube 2) includes the specific antibody to measure the actual level of endotoxin activity. The EAA result is calculated by normalizing the chemiluminescence in the test sample (Tube 2) against the maximum chemiluminescence (Tube 3), correcting both measurements for the basal activity chemiluminescence (Tube 1).

5.1 Specimen Collection and Preparation

Only use EDTA anti-coagulated blood collection tubes to obtain whole blood samples for this procedure. Guidelines recommended by the NCCLS⁶ should be followed when collecting, transporting and processing patient samples. Samples must be collected in blood collection tubes containing EDTA anti-coagulant by venipuncture or through an indwelling catheter. If sampling through an indwelling arterial or venous cannula or catheter, please flush the lines per your facility's policy prior to collecting blood in the EDTA anti-coagulated tube. Blood samples may be stored for up to 180 minutes at an ambient temperature range of 18°C to 25°C prior to analysis. Blood must be thoroughly mixed for at least 20-30 seconds by gentle inversion immediately prior to analysis.

A minimum of 1.5 mL of patient whole blood is suggested to run the patient assay (1.0 mL required). A minimum of 1.5 mL of whole blood is suggested to run the QC assay (1.0 mL required). Blood for QC test may be an independent whole blood sample.

5.2 Instruction for use with Berthold Smartline TL Tube Luminometer

NOTE: Prior to starting the EAA™ or EAA™ QC Test, ensure that Berthold SmartLine TL Tube Luminometer is “ON”, and the heating incubator shaker is “ON”, and warmed to the setting 37°C ± 1.0°C.

5.2.1 EAA™ Test Procedure

Good laboratory practice suggests that the EAA be performed in duplicate for each patient. The average of the duplicates is used to report the EAA™ result. Each laboratory should develop and apply acceptable CV's for the acceptance of replicates. You can perform up to 3 patient samples in one run.

1. Set up the EAA™ tubes for each patient sample in the supplied tube racks and remove the caps.
2. Combipipette 1.0 mL of the EAA™ Solution into Tubes 1, 2 and 3 (in duplicate).
3. After mixing the patient blood sample by gentle inversion, pipette 0.5 mL aliquots directly into the Aliquot Tube and the LPS MAX Tube. Vortex the LPS MAX tube.
4. Place the rack(s) in the incubating shaker, close the lid and incubate for 10 minutes.
5. Following the 10-minute incubation, open the lid, vortex the Aliquot tube then combipipette 40 µL of blood into Tubes 1 and 2 (in duplicate) using a sterile combipipette tip.
6. Vortex the LPS Max Tube, then combipipette 40 µL of blood into Tubes 3 (in duplicate) using a sterile combipipette tip.
7. Vortex Tubes 1, 2 and 3 and place back into the rack(s) in the incubating shaker, close the lid and start the motion (set at 100 rpm) and the timer (set at 14 minutes).
8. Insert the EAA labeled chip card into the luminometer and press Start.
9. Following the 14-minute incubation, follow the instructions on the luminometer display for reading the EAA tubes in the proper order.

10. Vortex each tube gently, open the sample drawer and place the tube in the sample holder, close the sample drawer and wait for an RLU reading. Repeat for the next tube.
11. After all tubes are read, the EAA results are printed automatically on the attached printout.
12. Repeat for each patient sample if applicable.

5.2.1.1 Interpretation of Results

Range (EAA™ result)	Interpretation of Results
0.0 to 0.39	Low endotoxin activity level represents a low risk for progression to severe sepsis.
0.40 to 0.59	Intermediate endotoxin activity level represents an elevated risk for severe sepsis.
≥ 0.60	High endotoxin activity level represents a high risk for developing severe sepsis

5.2.2 EAA™ QC Test Procedure

It is recommended that the QC test be run once per EAA™ kit box. However, each laboratory may produce their own QC testing regimen according to local standards.

1. Place the QC tubes into a supplied rack and remove the caps.
2. Combipipette 1.0 mL of the EAA™ Solution into each of the tubes.
3. Pipette 0.5 mL whole blood aliquots directly into the QC Aliquot Tube and the QC LPS MAX Tube. Vortex the QC LPS MAX tube.
4. Place the rack in the incubating shaker, close the lid and incubate for 10 minutes.
5. Following the 10-minute incubation, open the lid, vortex the QC Aliquot tube then combipipette 40 µL of blood into both Tube 1's and the Low Control Tube using a sterile combipipette tip.
6. Vortex the QC LPS MAX Tube, then combipipette 40 µL of blood

into the High Control Tube first and then both Tube 3's using a sterile combipipette tip.

7. Vortex all tubes and place back into the tube rack in the incubating shaker and start the motion (set at 100 rpm) and the timer (set at 14 minutes).
8. Insert the labeled EAA QC chip card into the luminometer and press Start.
9. Following the 14-minute incubation, follow the instructions on the luminometer display for reading the EAA QC tubes in the proper order.
10. Vortex each tube gently, open the sample drawer and place the tube in the sample holder, close the sample drawer and wait for an RLU reading. Repeat for the next tube.
11. After all tubes are read, the EAA™ QC results are printed automatically on the attached printout.

5.2.2.1 Interpretation of Results

EAA™ Solution QC is confirmed with the following results:

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If the QC test results for low or high control do not provide the desired results, please repeat the test. If the repeated QC test results remain outside acceptable limits, subsequent EAA results may be invalid. Please notify your EAA™ representative as soon as possible.

5.3 Limitations

5.3.1 Analytical

- As with all immunoassays or immunochemical reactions, careful technique is required.

- The user should be alert to the possible effects on results of potential interferences from medications or unknown endogenous substances. See Interfering Substances below.

5.3.2 Clinical

- The results of the EAA™ should be evaluated in context of all laboratory findings and the total clinical status of the patient. In cases where the laboratory results do not agree with the clinical picture or history, additional tests should be performed.
- The findings of positive endotoxin results in patients without severe sepsis have been reported previously in published studies.
- Results obtained on patients outside of the critical care setting should be interpreted with caution as they may not be of a significant clinical value.

5.4 Clinical Performance

The Multi-Center Endotoxin Detection in Critical Illness (MEDIC) trial was a multi-center, prospective observational study, performed in 10 Intensive Care Units (ICUs) of academic hospitals in North America and Europe. The presence of endotoxemia was evaluated on the first day of the patients' ICU stay to determine the odds of developing severe sepsis within 24 hours of ICU admission. The target population for the risk assessment study included all eligible patients enrolled in the MEDIC trial on first day of ICU admission who had evaluable samples, N=857. Patient demographics are shown in the following table.

Demographics and Baseline Features of Analyzed Patients

Characteristic	
Age: mean (SD) median [IQR]* years	60 (±17) 62 [49, 74]
Gender (% male)	58.9%
Race (% Non-Caucasian)	15.8%
APACHEII: mean (SD) median [IQR]*	15.2 (±9,5) 14 [8, 21]

Characteristic	
Days in ICU: mean (SD) median [IQR]*	5.4 (\pm 10,9) 2.0 [1, 5]
Days in Hospital: mean (SD) median [IQR]*	23.3 (\pm 27,0) 14.0 [7, 29]
ICU Mortality	13.3%
Hospital Mortality	19.9%

The critical care patient population is the most appropriate clinical environment for EAA™ testing. The utility of the EAA™ test has only been established for the critical care population. The relationship between endotoxin levels and the risk for developing severe sepsis in critically ill patients is intended as a supplement to currently available clinical information to establish risk for severe sepsis on the first day of ICU stay.

5.5 Performance Characteristics

5.5.1 Interfering Substances

- Triglyceride levels up to 1000 mg/dL do not interfere with EAA measurements. At triglyceride concentrations greater than 1500 mg/dL, a decrease in EAA results has been demonstrated. Grossly lipemic samples should be avoided for use in the EAA™.
- Hemoglobin at 100 mg/dL does not interfere with EAA measurements, while levels of 500 mg/dL attenuated low EAA™ measurements by less than EAA result of 0.03 and attenuated elevated EAA readings by an average of 15%. Grossly hemolyzed samples should be avoided for use in the EAA™.
- Exogenous infusion of endotoxin-free albumin into patients at levels which increase plasma albumin by more than 10 g/L above the upper limit of normal may lower EAA results. Samples from patients with hyper-protein anomalies or receiving albumin therapies should have protein determinations performed to ensure that elevated albumin levels are not a source of interference with the EAA™.

5.5.2 Precision

Whole blood samples, neat and with added exogenous endotoxin, were analyzed in eight replicates on each of three instruments. Both the within-run and total precision was calculated according to NCCLS procedures in EP5-A.

Sample	N	Mean EAA™	Within Run		Total Precision	
			1SD	%CV	1SD	%CV
1	24	0.11	0.015	14	0.023	22
2	24	0.20	0.030	15	0.029	14
3	24	0.30	0.043	15	0.042	14
4	24	0.50	0.059	12	0.064	13
5	24	0.52	0.036	7	0.034	6
6	24	0.59	0.050	8	0.046	8

5.5.3 Analytical Sensitivity

A low EAA patient sample was used to estimate the sensitivity of the EAA method. Twenty-four EAA determinations were made on a single patient sample with a low EAA level; eight measurements on each of three instruments. The mean EAA level of the patient was 0.11 with a SD of 0.023. Sensitivity, by the two standard deviation method, is thus estimated at EAA result of 0.046.

5.5.4 Analytical Specificity

Endotoxin from the Gram-negative bacteria, *E. coli* O55:B5, *P. aeruginosa*, *K. pneumoniae*, *S. enteritidis*, *E. coli* O127:B8 and *S. marcescens*, and *S. flexneri* induced an increase in EAA™ levels of 0.10 to 0.52 when 100 pg/mL of the endotoxin was added to a blood sample with EAA result of 0.33. In addition, preparations of *V. cholerae* LPS produced statistically significant increases in EAA results. No reactivity is seen with up to 2000 pg/mL lipoteichoic acid extracts from the following strains of Gram-positive bacteria: *S. mutans*, *S. pyogenes*,

S. sanguis, S. faecalis, S. aureus and B. subtilis. Similarly, there was no reactivity of yeast mannan up to 2000 pg/mL and cell wall extracts of C. albicans and Aspergillus fumigatus.

5.5.5 Linearity of Results

The assay has been optimized to be highly sensitive to low levels of endotoxin and to resist a high dose hook effect induced by grossly high endotoxin levels. Therefore, the curve is curvilinear in shape. When patient samples contain low levels of endotoxin, small increases in endotoxin will result in significant increases in EAA results. Dilutions are not required for samples that have generated high EAA results.

6. Storage and Stability

Store the EAA™ kit and its components (in the original sealed package) at 2-25°C. Do not use the EAA™ kit beyond the expiration date.

7. References:

1. Van Deventer SJH, et al. Endotoxemia: an early predictor of septicemia in febrile patients. The Lancet (March 19) 1988; 605-609.
2. Parrillo JE, Parker MM, Natanson C, et al. Septic shock in humans: advances in the understanding of pathogenesis, cardiovascular dysfunction, and therapy. Ann Int Med 1990; 113: 227-242.
3. Zhao B, Bowden RAS, Stavchansky SA, et al. Human endothelial cell response to gram-negative lipopolysaccharide assessed with cDNA microarrays. Am J Cell Physiol 2001; 281: C1587-C1595.
4. Romaschin AD, et al. A rapid assay of endotoxin in whole blood using autologous neutrophil dependent chemiluminescence. J. of Immunology. Methods 1998; 212: 169-185.
5. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # M29-A Protection of Laboratory Worker from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissues; Approved Guideline, 1997.

6. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # H3-#a Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture – Third Edition; Approved Standard, 1991.

FOR IN VITRO DIAGNOSTIC USE ONLY

Manufactured by:

Spectral Medical Inc.
135 The West Mall
Toronto, ON, Canada
M9C 1C2
+ 1.416.626.3233
FAX: + 1.416.626.7383
www.spectraldx.com

Authorized Representative for Europe:

Dr. Rolf Keck
Ritterstraße 25
D-77977 Rust/Baden
Germany
+ 49.782.286.5111
FAX: + 49.782.286.5113

Endotoxin Activity Assay (EAA™) is a registered trademark of Spectral Medical Inc.

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7. References:

1. Van Deventer SJH, et al. Endotoxemia: an early predictor of septicemia in febrile patients. *The Lancet* (March 19) 1988; 605-609.
2. Parrillo JE, Parker MM, Natanson C, et al. Septic shock in humans: advances in the understanding of pathogenesis, cardiovascular dysfunction, and therapy. *Ann Int Med* 1990; 113: 227-242.
3. Zhao B, Bowden RAS, Stavchansky SA, et al. Human endothelial cell response to gram-negative lipopolysaccharide assessed with cDNA microarrays. *Am J Cell Physiol* 2001; 281: C1587-C1595.
4. Romaschin AD, et al. A rapid assay of endotoxin in whole blood using autologous neutrophil dependent chemiluminescence. *J. of Immunology. Methods* 1998; 212: 169-185.
5. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # M29-A Protection of Laboratory Worker from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissues; Approved Guideline, 1997.
6. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # H3-#a Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture – Third Edition; Approved Standard, 1991.

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7. Literatur:

1. Van Deventer SJH, et al. Endotoxemia: an early predictor of septicemia in febrile patients. *The Lancet* (March 19) 1988; 605-609.
2. Parrillo JE, Parker MM, Natanson C, et al. Septic shock in humans: advances in the understanding of pathogenesis, cardiovascular dysfunction, and therapy. *Ann Int Med* 1990; 113: 227-242.
3. Zhao B, Bowden RAS, Stavchansky SA, et al. Human endothelial cell response to gram-negative lipopolysaccharide assessed with cDNA microarrays. *Am J Cell Physiol* 2001; 281: C1587-C1595.
4. Romaschin AD, et al. A rapid assay of endotoxin in whole blood using autologous neutrophil dependent chemiluminescence. *J. of Immunology. Methods* 1998; 212: 169-185.
5. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # M29-A Protection of Laboratory Worker from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissues; Approved Guideline, 1997.

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7. Bibliografia:

1. Van Deventer SJH, et al. Endotoxemia: an early predictor of septicemia in febrile patients. The Lancet (March 19) 1988; 605-609.
2. Parrillo JE, Parker MM, Natanson C, et al. Septic shock in humans: advances in the understanding of pathogenesis, cardiovascular dysfunction, and therapy. Ann Int Med 1990; 113: 227-242.
3. Zhao B, Bowden RAS, Stavchansky SA, et al. Human endothelial cell response to gram-negative lipopolysaccharide assessed with cDNA microarrays. Am J Cell Physiol 2001; 281: C1587-C1595.
4. Romaschin AD, et al. A rapid assay of endotoxin in whole blood using autologous neutrophil dependent chemiluminescence. J. of Immunology. Methods 1998; 212: 169-185.
5. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # M29-A Protection of Laboratory Worker from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissues; Approved Guideline, 1997.
6. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # H3-#a Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture – Third Edition; Approved Standard, 1991.

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Spectral Medical Inc.
135 The West Mall
Toronto, ON, Canada M9C 1C2
+ 1.416.626.3233
FAX: + 1.416.626.7383
www.spectraldx.com

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7. Referanslar:

1. Van Deventer SJH, et al. Endotoxemia: an early predictor of septicemia in febrile patients. The Lancet (March 19) 1988; 605-609.
2. Parrillo JE, Parker MM, Natanson C, et al. Septic shock in humans: advances in the understanding of pathogenesis, cardiovascular dysfunction, and therapy. Ann Int Med 1990; 113: 227-242.
3. Zhao B, Bowden RAS, Stavchansky SA, et al. Human endothelial cell response to gram-negative lipopolysaccharide assessed with cDNA microarrays. Am J Cell Physiol 2001; 281: C1587-C1595.

4. Romaschin AD, et al. A rapid assay of endotoxin in whole blood using autologous neutrophil dependent chemiluminescence. J. of Immunology. Methods 1998; 212: 169-185.
5. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # M29-A Protection of Laboratory Worker from Instrument Biohazards and Infectious Disease Transmitted by Blood, Body Fluids and Tissues; Approved Guideline, 1997.
6. NCCLS (National Committee for Clinical Laboratory Standards, 711 East Lancaster Avenue, Villanova, PA 19085) Document # H3-#a Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture – Third Edition; Approved Standard, 1991.

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Endotoksin Aktivite Tahlili (EAA™) iSpectral Medical Inc.'in tescilli ticari markasıdır.