

Appendix 4: Laboratory experiment result capture form*

EPHI: put PODTEC barcode here

1. Stool sample number	<input type="radio"/> 1 (first sample) <input type="radio"/> 2 <input type="radio"/> 3	<i>see middle part of code on barcode sticker on aliquot pot</i>
2. Aliquot number	--	<i>see last part of code on barcode sticker on aliquot pot</i>
3. Initials of laboratory technician conducting test: _____		
4. Remarks on stool appearance:	<input type="radio"/> looks normal <input type="radio"/> overgrown with fungi <input type="radio"/> gas formation has caused lid to drop off <input type="radio"/> other, specify: _____	<i>tick only one</i>
5. Experiment number:	<input type="radio"/> 1a (sample strategy: 3 different stools) → Q10 <input type="radio"/> 1b (sample strategy: 3 aliquots of 1 stool) → Q10 <input type="radio"/> 2 (storage conditions) → Q6 <input type="radio"/> 3a (processing conditions) → Q7 <input type="radio"/> 3b (maximum mixture storage time) → Q8 <input type="radio"/> 3c (optimum amount of stool processed) → Q9	<i>tick only one</i>
6. a. Sample storage temperature	<input type="radio"/> Room (20-22 ⁰ C) <input type="radio"/> Incubator (35-37 ⁰ C) <input type="radio"/> Refrigerator (2-6 ⁰ C) <input type="radio"/> Freezer (-20 ⁰ C)	<i>tick only one option</i>
b. Sample storage time	<input type="radio"/> 4-6 hours <input type="radio"/> 24 hours <input type="radio"/> 72 hours (3 days) <input type="radio"/> 120 hours (5 days)	<i>tick only one option</i> <i>After filling this question: → Q10</i>

7. a. Sample + SR buffer mixture processing time	<input type="radio"/> 10' + 10' <input type="radio"/> 10' + 5' <input type="radio"/> 15' <input type="radio"/> 10' + 10' vortex	<i>tick only one option</i>
8. a. What was the contact time between sample and SR after initial mixing?	<input type="radio"/> 2 hours <input type="radio"/> 5 hours <input type="radio"/> 24 hours	<i>tick only one option</i>
b. At what temperature was the sample/SR mixture stored?	<input type="radio"/> Room (20-22 ⁰ C) <input type="radio"/> Refrigerator (2-6 ⁰ C)	<i>tick only one option</i> <i>After filling this question: → Q10</i>
9. Amount of stool tested	<input type="radio"/> 0.3 g <input type="radio"/> 0.5 g <input type="radio"/> 0.8 g <input type="radio"/> 1.0 g <input type="radio"/> 1.2 g	<i>tick only one option</i>
10. Date and time of stool processing for Xpert-Ultra started	___ / ___ / ____ ___ : ___ DD / MM / YYYY HH : MM	
11. Date and time of Xpert-Ultra started	___ / ___ / ____ ___ : ___ DD / MM / YYYY HH : MM <input type="checkbox"/> test not started, reason: _____	
12. Xpert-Ultra result	<input type="radio"/> MTB trace detected <input type="radio"/> RIF resistant <input type="radio"/> MTB detected very low <input type="radio"/> RIF susceptible <input type="radio"/> MTB detected low <input type="radio"/> RIF indeterminate <input type="radio"/> MTB detected medium <input type="radio"/> MTB detected high <input type="radio"/> MTB not detected <input type="radio"/> Error code: ___ <input type="radio"/> Invalid <input type="radio"/> No result	<i>tick only one option</i> <i>for both the MTB</i> <i>and the RIF result</i>

13. Cycle threshold (Ct) values	Probe A: --	Probe E: --
	Probe B: --	Probe IS6110/IS1081: --
	Probe C: --	SPC: --
	Probe D: --	