


Protocol

Protocol for production and purification of SARS-CoV-2 3CL^{pro}



3CL^{pro} protease from SARS-CoV-2 is a primary target for COVID-19 antiviral drug development. Here, we present a protocol for 3CL^{pro} production in *Escherichia coli*. We describe steps to purify 3CL^{pro}, expressed as a fusion with the *Saccharomyces cerevisiae* SUMO protein, with yields up to 120 mg L⁻¹ following cleavage. The protocol also provides isotope-enriched samples suitable for nuclear magnetic resonance (NMR) studies. We also present methods to characterize 3CL^{pro} by mass spectrometry, X-ray crystallography, heteronuclear NMR, and a Förster-resonance-energy-transfer-based enzyme assay.

Publisher's note: Undertaking any experimental protocol requires adherence to local institutional guidelines for laboratory safety and ethics.

Luca Mazzei,
Rebecca Greene-
Cramer, Khushboo
Bafna, ..., Catherine
Ann Royer, Stefano
Ciurli, Gaetano T.
Montelione

stefano.ciurli@unibo.it
(S.C.)
monteg3@rpi.edu
(G.T.M.)

Highlights

Production and purification of SARS-CoV-2 3CL^{pro} with native N- and C-termini

3CL^{pro} inhibition assay and description for data processing/analysis

Structural characterization by X-ray crystallography and NMR spectroscopy

Detailed troubleshooting and alternative approaches

Mazzei et al., STAR Protocols
4, 102326
June 16, 2023 © 2023 The
Author(s).
<https://doi.org/10.1016/j.xpro.2023.102326>



Protocol

Protocol for production and purification of SARS-CoV-2 3CL^{pro}

Luca Mazzei,^{1,5,6} Rebecca Greene-Cramer,^{2,3,5,6} Khushboo Bafna,^{2,3} Aleksandar Jovanovic,¹ Anna De Falco,^{2,3} Thomas B. Acton,^{2,3} Catherine Ann Royer,^{2,4} Stefano Ciurli,^{1,*} and Gaetano T. Montelione^{2,3,7,*}

¹Department of Pharmacy and Biotechnology, University of Bologna, 40127 Bologna, Italy

²Center for Biotechnology and Interdisciplinary Sciences, Rensselaer Polytechnic Institute, Troy, NY 12180, USA

³Department of Chemistry and Chemical Biology, Rensselaer Polytechnic Institute, Troy, NY 12180, USA

⁴Department of Biology, Rensselaer Polytechnic Institute, Troy, NY 12180, USA

⁵These authors contributed equally

⁶Technical contacts: luca.mazzei2@unibo.it, greenr9@rpi.edu

⁷Lead contact

*Correspondence: stefano.ciurli@unibo.it (S.C.), monteg3@rpi.edu (G.T.M.)
<https://doi.org/10.1016/j.xpro.2023.102326>

SUMMARY

3CL^{pro} protease from SARS-CoV-2 is a primary target for COVID-19 antiviral drug development. Here, we present a protocol for 3CL^{pro} production in *Escherichia coli*. We describe steps to purify 3CL^{pro}, expressed as a fusion with the *Saccharomyces cerevisiae* SUMO protein, with yields up to 120 mg L⁻¹ following cleavage. The protocol also provides isotope-enriched samples suitable for nuclear magnetic resonance (NMR) studies. We also present methods to characterize 3CL^{pro} by mass spectrometry, X-ray crystallography, heteronuclear NMR, and a Förster-resonance-energy-transfer-based enzyme assay. For complete details on the use and execution of this protocol, please refer to Bafna et al.¹

BEFORE YOU BEGIN

Severe Acute Respiratory Syndrome CoronaVirus-2 3-Chymotrypsin-Like protease (SARS-CoV-2 3CL^{pro}), also referred to as SARS-CoV-2 Main protease (M^{pro}) or Non-structural Protein 5 (NSP5), is an essential viral cysteine protease required to hydrolyze the SARS CoV-2 1ab polyprotein at several sites during viral replication.^{2–4} It is recognized as a target for antiviral drug development. Structural and biochemical characterization of 3CL^{pro} inhibition to identify potential candidates for SARS-CoV-2 treatment and drug development rely on the production of large quantities of native and active 3CL^{pro}. Protein – ligand binding at the molecular level, as well as the potency of enzyme inhibition, can be investigated using various techniques, such as heteronuclear NMR, biomolecular crystallography, and biochemical and kinetic assays.

In this work, we present a detailed and integrated protocol for 3CL^{pro} production, purification, and structural characterization, as well as for biochemical investigation of enzyme inhibition. High yields of native and active 3CL^{pro} (30–120 mg per liter of broth, depending on the type of broth used for fermentation) are routinely obtained using *Escherichia coli* cells. The protein is initially expressed as a fusion with an N-terminal His-SUMO tag (His₆-SUMO-3CL^{pro}), then purified to produce the native target upon tag hydrolysis by SUMO protease (for which a production and purification protocol is also described). Protocols to produce isotopically enriched (e.g., ¹⁵N, ¹³C, ²H-enriched) 3CL^{pro} samples, to structurally characterize the protein using heteronuclear NMR, as well as to crystalize 3CL^{pro}






Figure 1. pGTM_CoV2_NSP5_004_SUMO construct design

(A) Multiple cloning site of the original vector pET15_SUMO2_NESG showing BsaI restriction sites used for cloning the CDS of 3CL^{pro}. Maps are generated using SnapGene® software (from Insightful Science; available at snapgene.com). The synthetic gene for 3CL^{pro} was obtained from Genscript, Inc. with standard Genscript codon optimization for *E. coli* expression.

(B) Final expression vector pGTM_CoV2_NSP5_004_SUMO.

(C) Schematic representation of the expressed His₆-SUMO-3CL^{pro}.

and determine the X-ray crystal structure at a resolution better than 2 Å, are also presented. These protocols provide a useful starting point for structural investigations of protein – inhibitor interactions. A detailed FRET-based enzyme assay for estimating IC₅₀ values of 3CL^{pro} inhibitors is also reported.

SARS-CoV-2 3CL^{pro} is a high priority target for drug development and biophysical studies, and its production has been described previously as a component of other STAR Protocols.^{5,6} This work aims to provide, in a single and integrated set of protocols, a complete toolbox to produce large amounts of native and active SARS-CoV-2 3CL^{pro} and characterize the binding of drug candidates to this crucial viral target.

Expression plasmids

His₆-SUMO-3CL^{pro} expression plasmid

Plasmid pGTM_CoV2_NSP5_004_SUMO (Figure 1) was obtained by cloning the synthetic coding sequence (CDS) for SARS-CoV-2 3CL^{pro} (NSP5) into the pET15_SUMO2_NESG expression vector,⁷ which contains the CDS for an N-terminal hexa-His tag upstream to a SUMO tag and SUMO protease cleavage site. The plasmid has a pMB1 origin of replication and carries the ampicillin resistance gene. Plasmid pGTM_CoV2_NSP5_004_SUMO is available from AddGene (AddGene ID: 190062); the expressed protein construct will be called hereafter His₆-SUMO-3CL^{pro}. The related plasmid pGTM_CoV2_NSP5_C145A_001_SUMO, with active-site mutation Cys145Ala, is also available from AddGene (AddGene ID: 192349).




Figure 2. His₈-MBP-SUMO^{Pro} construct design

(A) Multiple cloning site of the original vector pET15_8His_MBP_TEV_NESG. Maps are generated using SnapGene® software (from Insightful Science; available at snapgene.com). The synthetic gene for yeast SUMO^{Pro} was obtained from Genscript, Inc. with standard Genscript codon optimization for *E. coli* expression.
 (B) Final expression vector pGTM_YR375_SUMO_Protease_001.
 (C) Schematic representation of the expressed His₈-MBP-SUMO^{Pro} construct.

His₈-MBP-SUMO^{Pro} expression plasmid

Plasmid pGTM_YR375_SUMO_Protease_001 (Figure 2) was obtained by cloning the CDS for residues 347–621 of the Ubiquitin-like-specific protease 1 from *Saccharomyces cerevisiae* [Uniprot ID: Q02724 (ULP1_YEAST)] into the pET15_8His_MbpTEV_NESG expression vector,⁷ which contains the CDS for an N-terminal octa-His tag upstream to Maltose Binding Protein Solubility Enhancing Tag (MBP tag) and TEV protease cleavage site. This plasmid has a pMB1 origin of replication and carries the ampicillin resistance gene. Plasmid pGTM_YR375_SUMO_Protease_001 is available from AddGene (AddGene ID: 190063) and the expressed protein construct will be called hereafter His₈-MBP-SUMO^{Pro}.

Preparation of reagent stock solutions

⌚ Timing: 1 h

1. Prepare 50 mg mL⁻¹ ampicillin stock prepared in H₂O (or ²H₂O).
 - a. Weigh 2.5 g (or 0.050 g) of ampicillin sodium salt (Sigma Aldrich).
 - b. Dissolve it in a final volume of 50 mL (or 1.0 mL) using 50 % v/v ethanol solution (or ²H₂O).
 - c. Sterilize using 0.22 μm syringe filters and store at –20 °C.
2. Prepare 1 M isopropyl-β-d-thiogalactopyranoside (IPTG) in H₂O (or ²H₂O).
 - a. Weigh 2.38 g (or 0.238 g) of IPTG (molar mass 238.30 g mol⁻¹) (Sigma Aldrich).
 - b. Dissolve it in a final volume of 10 mL (or 1.0 mL) using MilliQ water (or ²H₂O).
 - c. Sterilize using 0.22 μm syringe filters and store at –20 °C as 1 mL aliquots.

3. Prepare 1 M magnesium chloride (MgCl_2) in H_2O (or $^2\text{H}_2\text{O}$).
 - a. Weigh 95.2 g (or 0.30 g) of MgCl_2 (molar mass 95.21 g mol^{-1}) (Sigma Aldrich).
 - b. Dissolve it in a final volume of 1 L (or 3.0 mL) using MilliQ water (or $^2\text{H}_2\text{O}$).
 - c. Filter the solution using $0.45 \mu\text{m}$ syringe filters and store at room temperature ($20\text{--}25 \text{ }^\circ\text{C}$).
4. Prepare 40 mg mL^{-1} deoxyribonuclease I from bovine pancreas (DNase I).
 - a. Weigh 40.0 mg of DNase I (Sigma Aldrich).
 - b. Dissolve it in a final volume of 1 mL using MilliQ water.
 - c. Store at $-20 \text{ }^\circ\text{C}$ as $30 \mu\text{L}$ aliquots.
5. Prepare 1 M dithiothreitol (DTT).
 - a. Weigh 1.54 g of DTT (molar mass $154.25 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - b. Dissolve it in a final volume of 10 mL using MilliQ water.
 - c. Store at $-20 \text{ }^\circ\text{C}$ as 1 mL aliquots.
6. Prepare 0.1 M Tris(2-carboxyethyl)phosphine (TCEP).
 - a. Weigh 286.7 mg of TCEP-HCl (molar mass $286.65 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - b. Dissolve it in a final volume of 10 mL using MilliQ water.
 - c. Store at $-20 \text{ }^\circ\text{C}$ as 1 mL aliquots.
7. Prepare 0.1 M NiSO_4 .
 - a. Weigh 1.55 g of anhydrous NiSO_4 (molar mass $154.76 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - b. Dissolve it in a final volume of 100 mL using MilliQ water.
 - c. Filter the solution using $0.45 \mu\text{m}$ syringe filters and store at room temperature.

For production using MJ9 medium, also prepare Trace Element Stock in H_2O or $^2\text{H}_2\text{O}$, as appropriate:

8. Prepare 60 mL of trace elements stock.⁸
 - a. Weigh 1.626 g of $\text{FeCl}_3 \times 6\text{H}_2\text{O}$ (molar mass $270.30 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - b. Weigh 0.24 g of $\text{ZnSO}_4 \times 7\text{H}_2\text{O}$ (molar mass $287.56 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - c. Weigh 0.42 g of $\text{CoCl}_2 \times 6\text{H}_2\text{O}$ (molar mass $237.93 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - d. Weigh 0.42 g of $\text{Na}_2\text{MoO}_4 \times 2\text{H}_2\text{O}$ (molar mass $241.95 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - e. Weigh 0.48 g of $\text{CuSO}_4 \times 5\text{H}_2\text{O}$ (molar mass $249.68 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - f. Weigh 0.12 g of H_3BO_3 (molar mass 61.83 g mol^{-1}) (Sigma Aldrich).
 - g. Weigh 0.3 g of MnSO_4 (molar mass $169.02 \text{ g mol}^{-1}$) (Sigma Aldrich).
 - h. Dissolve in a final volume of 54 mL of MilliQ water.
 - i. Add 6 mL of HCl (6 N) (Sigma Aldrich).

Note: Trace element stock for MJ9 medium can be stored for several months in frozen aliquots at $-20 \text{ }^\circ\text{C}$.

Note: Trace Element Stock is not required when using the LB or M9 medium protocols

Preparation of LB medium for bacterial growth

⌚ Timing: 4 h

9. Prepare Lysogeny Broth (LB) medium (calculations refer to 1 L final volume).
 - a. Weigh 25 g of LB broth powder (Miller) (Sigma Aldrich).
 - b. Dissolve the powder in a final volume of 1 L adding MilliQ water.
 - c. Sterilize in autoclave ($121 \text{ }^\circ\text{C}$, 20 min).
10. Prepare LB-agar plates for bacterial cell culture (calculations refer to 250 mL final volume).
 - a. Weigh 6.25 g of LB broth powder (Miller) (Sigma Aldrich) and 3.75 g of agar (Sigma Aldrich).
 - b. Add 250 mL of MilliQ water and dissolve the powder mixture (agar will completely dissolve only after heating during autoclave sterilization – see next point).
 - c. Sterilize in autoclave ($121 \text{ }^\circ\text{C}$, 20 min).

- d. Allow the sterilized LB-agar mixture to cool down to $\sim 50^{\circ}\text{C}$.
- e. Add $500\ \mu\text{L}$ of $50\ \text{mg mL}^{-1}$ ampicillin stock solution.

△ CRITICAL: Antibiotics must be added after cooling down to $\sim 50^{\circ}\text{C}$ to prevent thermal degradation.

- f. Pour the LB-agar into the Petri dishes under aseptic conditions inside the laminar flow cabinet.
- g. Wait until the agar solidifies and seal the plate with parafilm.
- h. Store the LB-agar plates at 4°C upside down, to prevent aqueous vapor condensation on the surface of the LB-Agar.

Alternatives: to obtain higher yields (up to $120\ \text{mg L}^{-1}$ fermentation), other rich media (*i.e.* Super Broth, Research Products International) can be used in place of LB.

Preparation of medium for production of ^{15}N -, ^{13}C - or ^{15}N , ^{13}C -enriched 3CL^{pro}

⌚ Timing: 1 h

11. MJ9 medium in H_2O . For ^{15}N - (or ^{15}N , ^{13}C -) enriched protein samples production in MJ9 medium, prepare MJ9 minimal medium as follows (calculations refer to 1 L final volume):
 - a. Weigh 6 g of K_2HPO_4 (J.T.Baker).
 - b. Weigh 9 g of KH_2PO_4 (J.T.Baker).
 - c. Weigh 1.5 g of $(\text{NH}_4)_2\text{SO}_4$ or $(^{15}\text{NH}_4)_2\text{SO}_4$ (Sigma Aldrich).
 - d. Weigh 0.5 g of sodium citrate dihydrate (J.T.Baker).
 - e. Weigh between 3 and 5 g of D-Glucose or D-Glucose ($\text{U-}^{13}\text{C}_6$ 99 %) (Sigma Aldrich, Cambridge Isotope).
 - f. Dissolve the powder in a final volume of 1 L adding MilliQ water.
 - g. Adjust the pH to 6.8.
 - h. Sterile filter.
 - i. Add 1 mL of vitamin stock (MEM Vitamin Solution 100 \times – Sigma Aldrich M6895).
 - j. Add 1 g of $\text{MgSO}_4 \times 7\text{H}_2\text{O}$ (J.T. Baker 2500-01).
 - k. Add 1 mL of trace elements stock.
 - l. Add 1 mL of stock ampicillin.

Alternatives: M9 medium. For ^{15}N - (or ^{15}N , ^{13}C -) enriched protein sample production in M9 medium, prepare M9 minimal medium as follows (calculations refer to 375 mL final volume):

- m. Weigh 2.25 g of Na_2HPO_4 (Sigma Aldrich).
- n. Weigh 1.12 g of KH_2PO_4 (J.T. Baker).
- o. Weigh 0.47 g of $(\text{NH}_4)_2\text{SO}_4$ or $(^{15}\text{NH}_4)_2\text{SO}_4$ (Sigma Aldrich).
- p. Weigh 0.09 g of MgSO_4 (Sigma Aldrich).
- q. Weigh 0.19 g of NaCl (Sigma Aldrich).
- r. Weigh between 1.2 and 1.9 g of D-Glucose or D-Glucose ($\text{U-}^{13}\text{C}_6$ 99 %) (Sigma Aldrich, Cambridge Isotope).
- s. Dissolve the powder in a final volume of 1 L by adding MilliQ water.
- t. Sterile filter.

Preparation of media for production of perdeuterated ^{15}N , ^{13}C , ^2H -enriched 3CL^{pro}

⌚ Timing: 2 h

12. For production of perdeuterated ^{15}N , ^{13}C , ^2H -enriched protein samples, prepare solutions as follows:

- a. For MJ9 minimal medium A (calculations refer to 50 mL final volume in H₂O).
 - i. Weigh 0.3 g of K₂HPO₄ (J.T.Baker).
 - ii. Weigh 0.45 g of KH₂PO₄ (J.T.Baker).
 - iii. Weigh 0.075 g of (¹⁵NH₄)₂SO₄ (Sigma Aldrich).
 - iv. Weigh 0.025 g of sodium citrate dihydrate (J.T.Baker).
 - v. Weigh between 0.15 and 0.25 g of D-Glucose (U-¹³C₆ 99 %) (Cambridge Isotope).
 - vi. Dissolve the powder in a final volume of 50 mL adding MilliQ water.
 - vii. Adjust the pH to 6.8.
 - viii. Sterile filter.
 - ix. Add 50 μL of vitamin stock (MEM Vitamin Solution 100× – Sigma Aldrich M6895).
 - x. Add 0.05 g of MgSO₄ × 7H₂O (J.T. Baker 2500-01).
 - xi. Add 50 μL of trace elements stock (1000×).
- b. For MJ9 minimal medium B (calculations refer to 100 mL final volume in ²H₂O).
 - i. Weigh 0.6 g of K₂HPO₄ (J.T.Baker).
 - ii. Weigh 0.9 g of KH₂PO₄ (J.T.Baker).
 - iii. Weigh 0.15 g of (¹⁵NH₄)₂SO₄ (Sigma Aldrich).
 - iv. Weigh 0.05 g of sodium citrate dihydrate (J.T.Baker).
 - v. Weigh between 0.3 and 0.5 g of D-Glucose (¹³C₆ 99 %) (Cambridge Isotope).
 - vi. Dissolve the powder in a final volume of 100 mL adding ²H₂O (CortecNet).
 - vii. Adjust the pH to 6.8.
 - viii. Sterile filter.
 - ix. Add 100 μL of vitamin stock (MEM Vitamin Solution 100× – Sigma Aldrich M6895).
 - x. Add 0.1 g of MgSO₄ × 7H₂O (J.T. Baker 2500-01).
 - xi. Add 100 μL of trace elements stock (1000×).
- c. For MJ9 minimal medium C (calculations refer to 500 mL final volume in ²H₂O).
 - i. Weigh 3.0 g of K₂HPO₄ (J.T.Baker).
 - ii. Weigh 4.5 g of KH₂PO₄ (J.T.Baker).
 - iii. Weigh 0.75 g of (¹⁵NH₄)₂SO₄ (Sigma Aldrich).
 - iv. Weigh 0.25 g of Sodium citrate dihydrate (J.T.Baker).
 - v. Weigh between 1.5 and 2.5 g of D-Glucose (U-¹³C₆ 99 %; 1,2,3,4,5,6,6-D7, 97–98 %) (Sigma Aldrich, Cambridge Isotope).
 - vi. Dissolve the powder in a final volume of 500 mL adding ²H₂O (CortecNet).
 - vii. Adjust the pH to 6.8.
 - viii. Sterile filter.
 - ix. Add 500 μL of vitamin stock (MEM Vitamin Solution 100× – Sigma Aldrich M6895).
 - x. Add 0.5 g of MgSO₄ × 7H₂O (J.T. Baker 2500-01).
 - xi. Add 500 μL of trace elements stock (1000×).

SUMO protease expression and purification

⌚ Timing: 6 days

Cell transformation and DNA extraction are performed following the protocol described in the section [transformation of *Escherichia coli* cells](#) using plasmid pGTM_YR375_SUMO_Protease_001 instead of plasmid pGTM_CoV2_NSP5_004_SUMO. Expression of His₈-MBP-SUMO^{PRO}, as well as cell lysis, soluble extract recovery, and a first step of immobilized metal affinity chromatography (IMAC) are performed following steps 1–3 as described in the section [expression and purification of SARS-CoV-2 3CL^{PRO}](#), using lysis, washing, and elution buffers at pH 7.5 instead of pH 8.0. The His₈-MBP-SUMO^{PRO} eluted from the IMAC is buffer-exchanged in a storage buffer containing 50 mM Tris-HCl, 150 mM NaCl, 10 % glycerol, 0.5 mM DTT, at adjusted pH 7.5; the expressed protease retains its proteolytic activity with no need for cleavage of the His₈-MBP tag by TEV protease. It can be concentrated up to 1 mg mL⁻¹ and stored in 1 mL aliquots at –80 °C. It will be used at a ratios (w/w) of 1:1000 to 1:100 to cleave the His₆-SUMO tag from His₆-SUMO-3CL^{PRO}, to obtain native 3CL^{PRO}.

Note: His₈-MBP-SUMO^{Pro} has a pI of ~ 10, and its purification protocol uses a buffer at pH 7.5. However, for His₆-SUMO-3CL^{Pro} (pI ~ 5.8) and 3CL^{Pro} (pI ~ 5.9), the purification protocol (described below) used a buffer at pH 8.0, which is farther from the isoelectric points of these constructs.

Note: Specific activity of His₈-MBP-SUMO^{Pro} may vary depending on the protein batch preparation, and the optimal ratio and cleavage temperature should be empirically determined by performing cleavage assays on small scale samples. In our experience, undergoing multiple freeze/thaw cycles can reduce protease activity. Therefore, we recommend freezing the purified SUMO^{Pro} in small (e.g. 50 µL) aliquots, which can be stored for up to 3 months. Individual aliquots should be thawed immediately before needed in the 12–16 h cleavage/dialysis.

Transformation of *Escherichia coli* cells

⌚ Timing: 18 h

13. Transformation of XL10 Gold cells by heat shock for DNA propagation purposes⁹ (Day 1).

Note: Additions of DNA and liquid medium to the cell suspension, as well as transferring and plating procedures, must be carried out working under aseptic conditions inside the laminar flow cabinet or using a Bunsen burner to prevent cell culture contamination.

- a. Pre-warm 1 mL of Super Optimal broth with Catabolite repression (SOC) medium at 37 °C in a 1.5 mL sterile microcentrifuge tube.
 - b. Thaw one aliquot (50–100 µL) of XL10 Gold *Escherichia coli* cells on ice.
 - c. Pipette ca. 50 ng of the pGTM_CoV2_NSP5_004_SUMO plasmid directly into the bacterial suspension.
 - d. Mix by gently flicking the tube and incubate on ice for 30 min.
 - e. Perform the heat shock by transferring the suspension at 42 °C (use either a dry or water bath) for 60 s and then on ice for 2 min.
 - f. Add the pre-warmed SOC medium to the suspension and incubate for 45 min at 37 °C under agitation (at 225 rpm).
 - g. Plate 100 µL of the suspension on a previously prepared LB-agar Petri dish with 100 µg mL⁻¹ ampicillin.
 - h. Incubate at 37 °C for 16 h under agitation (at 225 rpm).
- #### 14. Plasmid propagation and DNA extraction from XL10 Gold cells (Days 2 and 3).
- a. Pick a single colony from the plates prepared in step 13 and containing transformed XL10 Gold cells, inoculate it in 5 mL of SOC or LB medium with 100 µg mL⁻¹ ampicillin and incubate at 37 °C for 16 h under agitation (at 225 rpm).
 - b. Carry out the DNA extraction by using a NucleoSpin Plasmid (NoLid) kit (Macherey-Nagel) and following manufacturer's instructions (<https://www.mn-net.com/media/pdf/45/51/02/Instruction-NucleoSpin-Plasmid.pdf>, page 15).
 - c. Plasmid DNA can be stored in aliquots at –20 °C.

Alternatives: Other commercially available plasmid extraction kits may also be used.

- #### 15. Transformation of BL21(DE3) *Escherichia coli* cells by heat shock for protein expression purposes is carried out using the pGTM_CoV2_NSP5_004_SUMO plasmid extracted from XL10 Gold cells in the previous step, and following the protocol outlined in step 13 of this section using instead BL21(DE3) *Escherichia coli* cells. This will provide transformed BL21(DE3) cells in a plate that can be used for the next part of the protocol.

KEY RESOURCES TABLE

REAGENT or RESOURCE	SOURCE	IDENTIFIER
Bacterial strains		
XL10-Gold Ultracompetent Cells	Agilent Technologies	200317
BL21(DE3) Competent Cells	Agilent Technologies	200131
Chemicals, peptides, and recombinant proteins		
Ulp1 Sumo Protease	This work	Described above
Ampicillin sodium salt	Sigma-Aldrich	A9518
Dithiothreitol (DTT)	Sigma-Aldrich	D9779
Isopropyl-β-D-thiogalactopyranoside (IPTG)	Sigma-Aldrich	I6758
Tris(2-carboxyethyl)phosphine hydrochloride (TCEP)	Sigma-Aldrich	C4706
Super Optimal Broth with catabolite repression (SOC) medium	Sigma-Aldrich	S1797
LB broth powder (Miller)	Sigma-Aldrich	L3522
Super broth powder	Research Products International	S33050
Agar	Sigma-Aldrich	A1296
Dnase I	Sigma-Aldrich	DN25
Magnesium chloride (MgCl ₂)	Sigma-Aldrich	M8266
Sodium chloride (NaCl)	Sigma-Aldrich	S9888
Trizma base (Tris)	Sigma-Aldrich	T6066
Imidazole	Sigma-Aldrich	I2399
Ethylenediaminetetraacetic acid (EDTA)	Sigma-Aldrich	E9884
Iron(III) chloride hexahydrate (FeCl ₃ × 6H ₂ O)	Sigma-Aldrich	236489
Zinc sulfate heptahydrate (ZnSO ₄ × 7H ₂ O)	Sigma-Aldrich	221376
Cobalt(II) chloride hexahydrate (CoCl ₂ × 6H ₂ O)	Sigma-Aldrich	255599
Sodium molybdate dihydrate (Na ₂ MoO ₄ × 2H ₂ O)	Sigma-Aldrich	M1003
Copper(II) sulfate pentahydrate (CuSO ₄ × 5H ₂ O)	Sigma-Aldrich	209198
Boric acid (H ₃ BO ₃)	Sigma-Aldrich	B0394
Manganese(II) sulfate (MnSO ₄)	Sigma-Aldrich	M7506
Hydrochloric acid (HCl)	Sigma-Aldrich	XX0628
Sodium phosphate monobasic (NaH ₂ PO ₄)	Sigma-Aldrich	74092
Potassium phosphate monobasic (KH ₂ PO ₄)	J.T.Baker	3246-05
Potassium phosphate dibasic (K ₂ HPO ₄)	J.T.Baker	3252-01
Sodium citrate dihydrate (NaC ₆ H ₅ O ₇ × 2H ₂ O)	J.T.Baker	3646-01
Magnesium sulfate heptahydrate (MgSO ₄ × 7H ₂ O)	J.T.Baker	2500-01
MEM Vitamin Solution (100×)	Sigma-Aldrich	M6895
¹⁵ N-ammonium sulfate (¹⁵ NH ₄) ₂ SO ₄	Sigma-Aldrich	299286
¹² C D-Glucose	Sigma-Aldrich	G7520
D-Glucose (U- ¹³ C ₆ , 99 %; 1,2,3,4,5,6,6-D ₇ , 97–98 %)	Cambridge Isotope	CLM-1396-5
Deuterium oxide (² H ₂ O)	CortecNet	CD5251P1000
Nickel(II) sulfate (NiSO ₄)	Sigma-Aldrich	656895
Critical commercial assays		
NucleoSpin Plasmid NoLid	Macherey-Nagel	740499.50
Dabsyl-KTSAVLQ/SGFRKME-EDANS	GenScript Inc	Custom Synthesis
5(2Aminoethylamino)1Naphthalene sulfonic acid (EDANS)	MedChemExpress	ICN15474001
Recombinant DNA		
pGTM_CoV2_NSP5_004_SUMO	This work – available from Addgene	AddGene ID: 190062
pGTM_CoV2_NSP5_C145A_001_SUMO	This work – available from Addgene	AddGene ID: 192349
pGTM_YR375_SUMO_Protease_001	This work – available from Addgene	AddGene ID: 190063

(Continued on next page)

<i>Continued</i>		
REAGENT or RESOURCE	SOURCE	IDENTIFIER
<i>Other</i>		
Sonicator 3000	Misonix	S-3000/R2225
Amicon Ultra Centrifugal Filter Unit (MWCO 10000 Da)	Merck Millipore	UFC90110
5-mL HiTrap Chelating HP	Cytiva	17040901
HiPrep 26/10 Desalting	Cytiva	17508701
HiLoad® 16/60 Superdex® 75	Cytiva	28988946/17104401
Infinite M1000 PRO microplate reader	TECAN	30063849/140100256
Jupiter C4 column	Phenomenex	00F-4167-B0
ESI-Q-TOF spectrometer	Xevo, Waters	N/A
Shigemi 5 mm Symmetrical NMR microtube	Shigemi, INC	BMS-005B
Bruker INOVA 700 MHz NMR spectrometer	Bruker	N/A
EasyXtal 15-well crystallization plate	NeXtal	132006 132006
Protein Crystallization Cabinet	RUMED	PKS-E100
Magnetic Cryo Vials	Molecular Dimensions	MD7-402
SPINE Sample Changer Basket	Molecular Dimensions	MD7-510
Magnetic Mounted LithoLoops (0.20 mm)	Molecular Dimensions	MD7-136
Magnetic Cryo wand	Molecular Dimensions	MD7-411
<i>Software and algorithms</i>		
NMRPipe	Delaglio, F. et al. (1995) ¹⁰	https://www.ibbr.umd.edu/nmrpipe/
POKY	Lee, W. et al. (2021) ¹¹	https://sites.google.com/view/pokynmr
NMRFAM-SPARKY	Lee, W. et al. (2015) ¹²	https://nmrfam.wisc.edu/nmrfam-sparky-distribution/
MXCuBE	Gabalinho, J. et al. (2010) ¹³ Oscarsson, M. et al. (2019) ¹⁴	https://mxcube.github.io/mxcube/
XDS	Kabsch, W. (2010) ¹⁵	https://strucbio.biologie.uni-konstanz.de/xdswiki/index.php/Xds
CCP4	Winn, M. D. et al. (2011) ¹⁶	https://www.ccp4.ac.uk/
AIMLESS	Evans, P. R. and Murshudov, G. N. (2013) ¹⁷	https://www.ccp4.ac.uk/html/aimless.html
REFMAC	Murshudov, G. N. et al. (2011) ¹⁸	https://www.ccp4.ac.uk/html/refmac5/description.html
COOT	Emsley P. and Cowtan K. (2004) ¹⁹	https://www2.mrc-lmb.cam.ac.uk/personal/pemsley/coot/
Magellan™ v 7.2		https://lifesciences.tecan.com/software-magellan
GraphPad Prism		https://www.graphpad.com
Microsoft Excel		https://www.microsoft.com/en-us/microsoft-365/excel

MATERIALS AND EQUIPMENT

Lysis buffer at pH 8.0

Prepare fresh before every use by adding DNase and MgCl₂ to 30 mL of loading/washing buffer (see below)

Reagent	Final concentration	Amount for 30 mL
Dnase 40 mg mL ⁻¹	40 µg mL ⁻¹	30 µL
MgCl ₂ 1 M	10 mM	300 µL

Loading / Washing Buffer at pH 8.0 (for immobilized metal affinity chromatography, IMAC)

Reagent	Final concentration	Amount for 1 L
Tris-HCl	20 mM	2.42 g
NaCl	300 mM	17.53 g
Imidazole	10 mM	0.68 g
Dithiothreitol (DTT) 1 M	1 mM	1 mL

Prepare 1 L, adjust to pH 8.0, and store at room temperature (20–25 °C) for up to three days. Pass through 0.22 µm filter and degas before use

Elution buffer at pH 8.0 (for IMAC)

Reagent	Final concentration	Amount for 0.5 L
Tris-HCl	20 mM	1.21 g
NaCl	300 mM	8.76 g
Imidazole	300 mM	10.21 g
DTT 1 M	1 mM	0.5 mL

Prepare 0.5 L, adjust to pH 8.0, and store at room temperature (20–25 °C) for up to three days. Pass through 0.2 µm filter and degas before use

Size exclusion chromatography (SEC) buffer at pH 7.5

Reagent	Final concentration	Amount for 0.5 L
Tris-HCl	20 mM	1.21 g
NaCl	50 mM	1.46 g
EDTA	1 mM	0.146 g
TCEP 0.1 M	1 mM	5 mL

Prepare 0.5 L, adjust to pH 7.5, and store at 4 °C for up to one week. Pass through 0.2 µm filter and degas before use

Förster resonance energy transfer (FRET) buffer at pH 7.3

Reagent	Final concentration	Amount for 0.1 L
Tris-HCl	20 mM	0.242 g
NaCl	100 mM	0.292 g
EDTA	1 mM	0.0292 g
DTT 1 M	1 mM	1 mL
BSA*	0.1 mg mL ⁻¹	

Prepare 0.1 L, adjust to pH 7.3, and store at 4 °C

*To be added just before setting up the assay

STEP-BY-STEP METHOD DETAILS

Expression and purification of SARS-CoV-2 3CL^{Pro}

⌚ Timing: 4 days

This section details all the experimental steps needed to express His₆-SUMO-3CL^{Pro} from BL21(DE3) *Escherichia coli* cells harboring the pGTM_CoV2_NSP5_004_SUMO (or pGTM_CoV2_NSP5_C145A_SUMO) plasmid, and to obtain 3CL^{Pro} (or C145A-3CL^{Pro}) in a pure and native form.

1. Expression of His₆-SUMO-3CL^{Pro} by induction with IPTG of transformed BL21(DE3) cells (days 1 and 2):

Note: The transferring of cell cultures, as well as the addition of reagents and the collection of cell samples for the measurement of cellular growth must be carried out working under aseptic conditions inside the laminar flow cabinet or using a Bunsen burner to prevent cell culture contamination.

- a. Prepare a starter culture of transformed BL21(DE3) cells by picking and inoculating a single colony (from plates prepared at step 13 in the [before you begin](#) section, which contain transformed BL21(DE3) cells) in 50 mL of LB broth with $100 \mu\text{g mL}^{-1}$ ampicillin added. Incubate the resulting suspension for 16 h at 37°C under agitation (at 225 rpm).
- b. Transfer the starter culture to 1 L of LB broth with $100 \mu\text{g mL}^{-1}$ ampicillin added and incubate the suspension at 37°C under agitation (225 rpm).
- c. Measure the optical density at 600 nm (OD_{600}) every 30 min until it reaches ca. 0.6 units (approximately 1.5–2 h).
- d. Add IPTG to a final concentration of 1 mM to induce protein expression.
- e. Incubate for 16 h at 17°C under agitation (225 rpm).

Alternatives: To prepare isotopically-enriched samples for NMR experiments using MJ9 minimal medium optimized for isotope-enrichment,⁸ at the end of point a above, transfer 25 mL of the obtained starter culture to 1 L of MJ9 minimal medium added with $50 \mu\text{g mL}^{-1}$ ampicillin. The resulting suspension should be incubated at 37°C under stirring (225 rpm) and treated as described from point c onwards.

Alternatives: While the best results are usually obtained using MJ9 medium, M9 minimal medium can also be used. In this case, once OD_{600} reaches ca. 0.6 units (point c above), harvest the cells by centrifugation for 20 min at $7000 \times g$ at room temperature and resuspend the pellet in 375 mL of M9 minimal medium containing ^{15}N ammonium sulphate (and/or ^{13}C -enriched glucose) added with $50 \mu\text{g mL}^{-1}$ ampicillin. Incubate the resulting suspension for 20 min at 17°C under agitation (225 rpm), then induce and further incubate as described in points d and e.

Alternatives: To prepare ^{15}N , ^{13}C , ^2H -enriched samples for NMR experiments using MJ9 minimal medium optimized for isotope-enrichment,⁸ fully replace step 1 a-e with the following procedure:

- f. Prepare a starter culture of transformed BL21(DE3) cells by picking and inoculating a single colony (from plates prepared at step 13 in the [before you begin](#) section) in 1.0 mL of LB medium added with $100 \mu\text{g mL}^{-1}$ ampicillin. Incubate the resulting suspension at 37°C until OD_{600} reaches 0.6 units (approximately 1.5–2 h).
- g. Harvest the cells by centrifugation ($7000 \times g$) at room temperature.
- h. Resuspend the cell pellet in 3.0 mL of MJ9 medium A added with $50 \mu\text{g mL}^{-1}$ ampicillin.
- i. Incubate the cell suspension in a 15 mL Falcon Tube at 37°C under agitation (225 rpm) until OD_{600} reaches 0.5 units (approximately 1.5–2 h).

△ CRITICAL: OD_{600} must not exceed 0.65 units to ensure that the culture is still in the log phase.

- j. Prepare four separate dilutions (1:40, 1:100, 1:200, 1:400) of the cells grown in MJ9 medium A in a final volume of 20 mL using MJ9 medium B (with $50 \mu\text{g mL}^{-1}$ ampicillin prepared in $^2\text{H}_2\text{O}$ added).
- k. Incubate the cell suspensions for 16 h at 37°C under agitation (225 rpm).
- l. Measure the OD_{600} of the cultures and select the one with OD_{600} closest to 0.5 units.
- m. Harvest the cells of the selected culture by centrifugation ($7000 \times g$) at room temperature ($20\text{--}25^\circ\text{C}$) for 5 min.

- n. Resuspend the cell pellet in 500 mL of MJ9 medium C with 50 $\mu\text{g mL}^{-1}$ ampicillin prepared in $^2\text{H}_2\text{O}$ added.
 - o. Incubate the cell suspension at 37°C under agitation (225 rpm) until OD_{600} reaches 0.35 units, then further incubate for 45 min at 17 °C under agitation (225 rpm).
 - p. Induce protein expression with 1 mM IPTG (using IPTG stock prepared in $^2\text{H}_2\text{O}$).
 - q. Incubate the culture for 12–16 h at 17 °C under agitation (225 rpm).
2. Cell lysis and soluble extract recovery (Day 3).
- a. Harvest the cells by centrifuging the culture for 30 min at 7000 $\times g$ at 4 °C.
 - b. Gently resuspend the harvested cells in 20 mL lysis buffer.
 - c. Lyse cells by sonication (50 % power, 5 s/5 s on/off, for a total on time of 10 min) on ice.
 - d. Collect a 10 μL aliquot of the lysate (total lysate, TL) in a 1.5 mL microcentrifuge tube for SDS-PAGE analysis; add the appropriate amount of loading dye (usually 5 μL or what is recommended by manufacturer) and freeze at -20 °C.
 - e. Centrifuge the lysate for 30 min at 17000 $\times g$ at 4 °C to separate the soluble extract from the insoluble proteins and cell debris: collect the soluble extract. Collect 10 μL aliquots each of the soluble extract (S) and the insoluble pellet (I) in 1.5 mL microcentrifuge tubes for SDS-PAGE analysis, add the loading dye to each tube, and freeze at -20 °C. Discard the remaining insoluble portion according to standard biosafety protocols.
 - f. Filter the soluble extract using 0.45 μm syringe filters and store at 4 °C.

Alternatives: Lysis can also be carried out by using a French press (SLM Aminco). The harvested cells can be resuspended in ca. 25 mL of lysis buffer and homogenized by using a Potter-Elvehjem tissue grinder to obtain a homogeneous cell suspension and to disrupt any solid salt precipitate possibly formed during cellular growth. The resuspended cells are lysed performing three lysis cycles using the French press at 20,000 psi (1 psi = 6.9 kPa).

△ CRITICAL: While lysing, both procedures cause an abrupt temperature increase in the cell suspension. It is highly recommended to keep the lysate on ice and to use a container with good thermal conductivity during the procedure to prevent thermal protein denaturation that would negatively affect final protein yield.

3. Immobilized metal affinity chromatography (IMAC) (Day 3).

Note: The expressed SARS-CoV-2-3CL^{Pro} contains a His₆-SUMO tag at its N-terminus that allows for protein purification using immobilized nickel affinity chromatography carried out using an ÄKTAprime system with a 5 mL HiTrap Chelating HP column (Cytiva) working at a flow rate of 4 mL min⁻¹.

Note: Due to the high amount of His₆-SUMO-3CL^{Pro} expressed from 1 L of cell culture, it is recommended to split the soluble fraction and perform the IMAC in two steps to prevent column overloading and possible loss of protein. Two or more columns placed in series can also be used to achieve high amounts of protein recovery and speed up the purification procedure.

- a. Wash the 5 mL HiTrap Chelating HP column with five column volumes (CVs) of MilliQ water and regenerate with 0.1 M NiSO₄ following the manufacturer's protocol <https://cdn.cytivalifesciences.com/api/public/content/digi-11225-pdf>.
- b. Equilibrate the column with 5 CVs of washing buffer.
- c. Load the soluble extract onto the column either through a superloop (Cytiva) (using washing buffer) or directly with pump A.
- d. Wash the column with washing buffer and collect the unbound material (hereafter called Flow-Through, FT). Monitor the absorbance at 280 nm (Abs_{280}) and wait until the baseline returns to zero (typically ca. 10–20 CVs) to ensure complete removal of the FT.




Figure 3. Immobilized metal (nickel) affinity chromatography (IMAC)

(A) Representative chromatogram showing the entire process of protein loading, flow through collection, and protein elution.

(B) SDS-PAGE showing total lysate (TL), soluble extract (S), flowthrough (FT) and fractions eluted in peak 1 (P1) of the IMAC. His₆-SUMO-3CL^{P^{ro}} (molar mass of ca. 46 kDa) is detected in lanes corresponding to peak P1.

- e. Collect a 10 μ L aliquot of FT in a 1.5 mL microcentrifuge tube for SDS-PAGE analysis, add the loading dye, and freeze at -20°C . Store the rest of FT at 4°C .
 - f. Elute the bound proteins by increasing the imidazole concentration by setting a linear gradient from 0 to 100 % of elution buffer in ten CVs and collecting eluate in 2 mL fractions. A representative chromatogram is shown in Figure 3A.
 - g. Collect a 10 μ L aliquot of each fraction in 1.5 mL microcentrifuge tubes for SDS-PAGE analysis, add the loading dye and freeze at -20°C .
 - h. Once the elution is complete, wash the column with five CVs of elution buffer, MilliQ water, and 20 % v/v ethanol each to return the column in storing conditions. Store the column at 4°C .
 - i. Analyze the total lysate (TL), the soluble extract (S), as well as the FT and the fractions eluted from the IMAC, for protein presence and purity by running an SDS-PAGE loading 5 μ L of the samples previously collected and boiled for 5 min at 90°C (a representative protein gel is shown in Figure 3B).
 - j. Pool the fractions containing His₆-SUMO-3CL^{P^{ro}} and store at 4°C .
4. Cleavage of the His₆-SUMO tag by SUMO protease and buffer exchange (Day 3).

Note: The proteolytic cleavage of the N-terminal His₆-SUMO tag is performed by the addition of His₈-MBP-SUMO^{P^{ro}} to the expressed His₆-SUMO-3CL^{P^{ro}}. General protocols for SUMO protease cleavage recommend carrying out cleavage in the presence of imidazole concentrations lower than 150 mM to prevent adverse effects of imidazole on the activity of SUMO protease. His₆-SUMO-3CL^{P^{ro}} elutes at ca. 180 mM imidazole concentration from the IMAC, therefore a buffer exchange may be useful to decrease the concentration of imidazole in the protein solution prior to treating with SUMO protease. However, uncontrolled cleavage of His₆-SUMO tag from His₆-SUMO-3CL^{P^{ro}}, probably by trace amounts of native *E. coli* proteases, has been observed prior to the addition of SUMO protease. This can be minimized by carrying out cleavage in imidazole immediately after purifying the His₆-SUMO-3CL^{P^{ro}} fusion, with minimal to no negative effects on the production of pure and native 3CL^{P^{ro}} (as described in this protocol). Nevertheless, if SUMO protease activity issues are observed, a buffer exchange into the Loading / Washing Buffer can be carried out before treatment with SUMO protease by performing points c-e before points a and b of the following procedure. Buffer exchange can be carried out using an ÄKTAprime system and a HiPrep 26/10 Desalting column (Cytiva) working at a flow rate of 10 mL min^{-1} , or by dialysis. In any case, buffer exchange is mandatory

to decrease the concentration of imidazole in the protein prior to the following protocol for removal of the cleaved tag and the SUMO protease from the native 3CL^{Pro} by reverse IMAC.

Note: Incubation times for complete cleavage of N-terminal His₆-SUMO tag by SUMO protease may vary depending on protease activity. As described in the [before you begin](#) section, the optimal ratio and cleavage temperature (4 or 26 °C) should be empirically determined by performing cleavage assays on small-scale samples. In this protocol a 1:100 (w/w) ratio of SUMO protease generally provides > 99 % cleavage after 3 h when carried out in a pH 8.0 buffer containing reducing agent (e.g. the Elution Buffer) at 26 °C.

- a. Collect a 10 µL aliquot of the pre-cleaved His₆-SUMO-3CL^{Pro} in a 1.5 mL microcentrifuge tube for SDS-PAGE analysis, add the loading dye, and freeze at -20 °C.
- b. Add 1:100 (w/w) ratio of His₈-MBP-SUMO^{Pro} (ULP1) to the protein solution (estimated by absorbance at 280 nm and considering $\epsilon_{280} = 100,855 \text{ M}^{-1} \text{ cm}^{-1}$ and $32,890 \text{ M}^{-1} \text{ cm}^{-1}$ for His₈-MBP-SUMO^{Pro} and His₆-SUMO-3CL^{Pro}, respectively) in Elution Buffer (or after buffer exchange in Loading / Washing Buffer to reduce imidazole concentrations, if preferred) and incubate for 3 h at 26 °C to ensure complete cleavage of the N-terminal tag.
- c. During this hydrolysis incubation, wash the HiPrep 26/10 Desalting column (Cytiva) with five CVs of MilliQ water and equilibrate it with five CVs of Loading / Washing Buffer.
- d. At the end of the hydrolysis incubation, load up to 10 mL of protein solution onto the HiPrep 26/10 Desalting column using a sample loop of the appropriate volume, and carry out the buffer exchange and protein elution using Loading / Washing Buffer, monitoring Abs₂₈₀ and stopping sample collection when Abs₂₈₀ returns to baseline.

Note: Imidazole elution strongly affects the conductance of the eluate. Therefore, it is recommended to monitor the change of this parameter value during elution to ensure the correct separation from the protein solution (to read more about the basis of desalting procedures and technical details using HiPrep 26/10 Desalting visit manufacturer's web site at <https://www.cytivalifesciences.com/>).

- e. Complete washing of the Desalting column with five CVs of washing buffer, MilliQ water, and 20 % v/v ethanol each to return the column in storing conditions. Store the column at 4 °C.

Alternatives: Cleavage of His₆-SUMO tag from 3CL^{Pro} and buffer exchange can also be carried simultaneously by dialyzing the protein eluted from the IMAC in the presence of a 1:100 (w/w) ratio SUMO protease, for 12–16 h at 4 °C, against Loading / Washing Buffer prepared without imidazole (i.e., 20 mM Tris pH 8.0, 300 mM NaCl, 1 mM DTT).

5. Removal of the cleaved tag and SUMO protease from cleaved 3CL^{Pro} solution (Day 4).

Note: The cleaved His₆-SUMO tag and His₈-MBP-SUMO^{Pro} both present N-terminal His-tags allowing for separation from the native 3CL^{Pro} through an additional step of IMAC. Cleaved 3CL^{Pro} will not interact with the column resin and will therefore be collected in the flow-through, while His₆-SUMO tag and His₈-MBP-SUMO^{Pro} will bind to the column and will be eluted at a higher imidazole concentration. Separation is carried out using an ÄKTAprime system and a 5 mL HiTrap Chelating HP column (Cytiva) loaded with Ni²⁺, working at a flow rate of 4 mL min⁻¹.

- a. Wash the 5 mL HiTrap Chelating HP column (Cytiva) used the day before with five CVs of MilliQ water and equilibrate it with five CVs of Loading / Washing Buffer.
- b. Load the protein solution onto the column using a superloop (Cytiva) or pump A (while using Loading / Washing Buffer).

- c. Wash the column with Loading / Washing Buffer while monitoring the absorbance at 280 nm (Abs_{280}) to recover the cleaved, native 3CL^{PRO} in the flowthrough. Store the protein solution at 4 °C.
- d. Collect a 10 μ L aliquot of protein solution in a 1.5 mL microcentrifuge tube for SDS-PAGE analysis, add the loading dye, and freeze at -20 °C.
- e. Elute the bound proteins remaining on the column using five CVs of 100 % Elution Buffer until the Abs_{280} baseline returns to zero.
- f. Collect a 10 μ L aliquot of the eluted proteins in a 1.5 mL microcentrifuge tube for SDS-PAGE analysis, add the loading dye and freeze at -20 °C.
- g. Wash the column with five CVs of MilliQ water and five CVs of 20 % v/v ethanol to return the column in storing conditions. Store the column at 4 °C.

Alternatives: To minimize cleavage of His₆-SUMO tag from 3CL^{PRO} by native *E. coli* proteases, an on-column cleavage of His₆-SUMO tag can be performed directly on the 5 mL Hi-Trap Chelating HP column prior to eluting with imidazole. This useful shortcut allows to perform IMAC purification, cleavage of His₆-SUMO tag from 3CL^{PRO}, and removal of cleaved His₆-SUMO tag and His₈-MBP-SUMO^{PRO} from native 3CL^{PRO} in one step. In the on-column cleavage protocol, follow points a-e of step 3 (*Immobilized metal affinity chromatography*). Then add 1–2 mL of 1 mg mL⁻¹ His₈-MBP-SUMO^{PRO} (1 mL is sufficient for lower yields, e. g., batches fermented in minimal medium, but 2 mL may be needed for higher yield batches). Run one CV of Loading / Washing Buffer to fully load the protease onto the column and allow cleavage by storing the column 12–16 h at 4 °C. The next morning, elute the cleaved 3CL^{PRO} with Loading / Washing Buffer and use the Elution Buffer to elute any remaining uncleaved His₆-SUMO-3CL^{PRO}, the His₈-MBP-SUMO^{PRO}, and the cleaved His₆-SUMO tag. In this protocol, the target product is also not exposed to high concentrations of imidazole.

6. Size exclusion chromatography (SEC) (Day 4).

Note: A final size exclusion chromatographic (SEC) step provides a fully monodisperse and pure 3CL^{PRO} sample in a buffer suitable for downstream structural and biochemical studies. This step is carried out using an ÄKTAprime system and a HiLoad® 16/60 Superdex® 75 column (Cytiva), working at a flow rate of 1 mL min⁻¹.

- a. Concentrate the protein solution to a final volume \leq 5 mL using an Amicon Ultra Centrifugal Filter Unit (MWCO 10,000 Da) (Merck Millipore).

△ CRITICAL: Although not strictly necessary, this concentration step is highly recommended to expedite the following SEC. During this step it is important to monitor the concentration of the protein solution (by measuring the absorbance at 280 nm and considering $\epsilon_{280} = 32,890 \text{ M}^{-1} \text{ cm}^{-1}$), which should not exceed 0.8 mM to avoid possible precipitation. If the desired volume of 5 mL is not achievable, split the protein solution into two or more aliquots to be independently treated in the next purification step.

- b. Wash the column with one CV of MilliQ water and equilibrate it with one CV of SEC buffer.
- c. Load up to 5 mL of protein solution using a sample loop of the appropriate volume and carry out the protein elution using SEC buffer (a representative chromatogram is shown in [Figure 4A](#)), collecting the eluate in 1 mL fractions and store at 4 °C.
- d. Collect a 10 μ L aliquot of each fraction in 1.5 mL microcentrifuge tubes for SDS-PAGE analysis, add the loading dye, and freeze at -20 °C.
- e. Repeat points c - d for each 5 mL of the remaining protein solution.
- f. Upon complete elution of the last protein sample wash the column with one CV of MilliQ water and one CV 20 % v/v ethanol. Store the column at 4 °C.




Figure 4. Size exclusion chromatography of native 3CL^{Pro}

(A) Representative SEC chromatogram and (B) corresponding SDS-PAGE from elution of the SEC column. Lanes C, FT, P1, and P2 correspond respectively to the pre-loaded cleaved protein solution, molecular weight standards, native 3CL^{Pro} eluted in peak P1 (at a molar mass of ca. 33 kDa for the monomer), and the His₆-SUMO tag eluted in peak P2 (at a molar mass of ca. 12 kDa).

- g. Analyze all collected fractions for protein presence and purity by running an SDS-PAGE gel loading 5 μL of the samples previously collected and boiled for 5 min at 90°C. Pool the fractions containing native 3CL^{Pro}. A representative SDS-PAGE is shown in Figure 4B.
- h. Native 3CL^{Pro} can be concentrated up to 25 mg mL⁻¹ (corresponding to 0.8 mM monomer). At these concentrations, native 3CL^{Pro} can be flash-frozen in liquid nitrogen and stored for up to one month at -80 °C.

Note: SEC buffer here reported is one among the several buffers in which native 3CL^{Pro} can be stored. Buffers most compatible with other downstream studies can be chosen as well.

Alternatives: In our experience, the cleaved His₆-SUMO tag and the His₈-MBP-SUMO^{Pro} can also be effectively separated from the native 3CL^{Pro} directly through the SEC procedure, rather than using the HiTrap column. To pursue this alternative, follow the purification procedure until point 4b (*Cleavage of the His₆-SUMO tag by SUMO protease*), then skip step 4 points c - e (*Buffer exchange*) and step 5, and go directly to step 6.

Characterization of SARS-CoV-2 3CL^{Pro} by mass spectrometry

⌚ Timing: 3 h

7. Sample separation by HPLC.
 - a. Dilute a 1 mg mL⁻¹ solution of 3CL^{Pro} 1:10 (v/v) in MilliQ water.
 - b. Load the diluted protein sample onto a Jupiter C4 column (Phenomenex) working at a flow rate of 0.4 mL min⁻¹ using mobile phases A (H₂O/formic acid; 100/0.1; v/v) and B (acetonitrile/formic acid; 100/0.1; v/v) to perform a gradient to enhance the separation of the eluting species.
8. Analyze the eluted species using an ESI-Q-TOF spectrometer (Xevo, Waters) using the software provided with the instrument.

Note: The resulting deconvoluted spectrum shows a single signal corresponding to an experimental molar mass of 33,796 Da (Figure 5). The experimental value perfectly matches that estimated by the amino acid protein sequence (UniProt entry P0DTD1, residues 3264–3569) and establishes that the current protocol for the expression and purification of SARS-CoV-2 3CL^{Pro} is a reliable procedure to obtain large quantities of native 3CL^{Pro}.




Figure 5. Characterization of purified 3CL^{Pro}

(A) SDS-PAGE gel showing high quality of SARS-CoV-2 3CL^{Pro} sample. Lane a – molecular weight standards; Lane b – SARS-CoV-2 3CL^{Pro} post SEC purification (sample prepared as described above in "Size exclusion chromatography (SEC)" section) (B) ESI-MS profile of SARS-CoV-2 3CL^{Pro} showing the signal corresponding to an experimental molecular mass of 33,796 Da.

FRET (Förster resonance energy transfer) assay for the characterization of 3CL^{Pro} inhibition by IC₅₀ determination

⌚ Timing: 2 h

Note: 3CL^{Pro} inhibition properties are characterized by performing FRET assays using peptide substrate KTSAVLQ/SGFRKME (designed with the 3CL^{Pro} cleavage site Q/S) labeled with a DABCYL and EDANS FRET pair (fluorophore/quencher) on N and C termini, respectively (DABCYL-KTSAVLQ/SGFRKME-EDANS) (GenScript). Excitation at 360 nm excites the donor EDANS which transfers its energy non-radiatively due to the close proximity of the acceptor (quencher), DABCYL. Upon enzymatic cleavage of the labeled FRET peptide by 3CL^{Pro}, energy transfer no longer occurs between the fluorophore and quencher, resulting in an increase in fluorescence that can be detected at the emission maximum of the fluorophore EDANS (460 nm). This rate of increase in fluorescence is directly proportional to the enzyme activity. FRET assays are carried out at increasing concentrations of potential 3CL^{Pro} inhibitors to screen candidate molecules and evaluate their IC₅₀ values (*i.e.*, the inhibitor concentration providing a 50 % decrease of enzyme activity).

FRET assays are carried out using black opaque 96-well microplates (Corning) containing a final volume of 100 μ L per well. [Figures 6](#) and [7](#) provide the detailed amounts of reagents placed in each well. Experiments should be run in triplicate and results are monitored and analyzed using a TECAN Infinite M1000 microplate reader and Magellan™ software (or comparable multi-well fluorimeter).

9. Instrument protocol setup:

Use the Magellan™ software of TECAN Infinite M1000 microplate reader with the following parameters: (i) Excitation wavelength: 360 nm; (ii) Emission wavelength: 460 nm. Fluorescence intensity is recorded every 3 s for a total of 15–60 min at room temperature (20–25 °C).

10. Calibration curve for fluorescence-to-activity conversion.

Note: Raw data collected during the FRET experiments indicate fluorescence intensity due to substrate proteolysis and are reported in terms of Relative Fluorescence Units (RFU). Preparation of a calibration curve for the conversion of the fluorescence data to cleaved peptide concentration (μ M) using 0.10 to 20 μ M EDANS is described below:




Figure 6. Materials preparation for calibration curve

(A) Schematic of 96 microwell plate set-up for FRET calibration curve.

(B) Table of reagents for each well.

- Prepare a stock solution of 500 μM EDANS (MedChem Express) in DMSO.
- Prepare 100 μM and 10 μM dilutions of the 500 μM stock solution using DMSO.
- Add FRET buffer to each well according to [Figure 6](#).
- Add the appropriate volume of EDANS to each well so that the final concentrations range from 0.10 to 20 μM.
- Measure emission at 460 nm using the Magellan™ software of TECAN Infinite M1000 microplate reader.
- Plot the measured fluorescence as a function of EDANS concentration in Microsoft Excel or GraphPad Prism.
- Perform a linear fit of the data.

Note: The resulting best-fit equation will be used to interpolate experimental fluorescence data obtained in the following experiment and convert from RFU per second to substrate hydrolysis rates ($\mu\text{M s}^{-1}$).

Alternatives: A calibration curve can also be set up by incubating 200 nM 3CL^{Pro} with increasing concentrations of FRET substrate (0.5–50 μM) and monitoring the reaction until the fluorescence signal reaches plateaus for all concentrations, thus indicating that the total amount of the FRET substrate has been cleaved by 3CL^{Pro}. Final fluorescence values can be plotted as a function of the FRET substrate concentration and fit linearly to generate a calibration curve as described in step 10.

Note: For certain FRET studies it is necessary to calibrate the fluorimeter for inner-filter effects arising from UV absorbance by the substrate (or inhibitor) at high concentrations.²⁰ However, at the concentrations of substrate described in this protocol, no inner filter effect corrections are required.

- 3CL^{Pro}, FRET substrate and inhibitor preparation:
 - Prepare a 10 mM peptide substrate stock solution in DMSO, and store it as 60 μL aliquots at –20 °C.
 - Prepare 5 mM peptide substrate samples from the 10 mM stock by diluting 35 μL with an equal volume of FRET buffer.
 - Prepare inhibitor stocks dissolving them in DMSO and diluting to 10.0, 1.0, and 0.1 mM concentrations in DMSO.
 - Gently thaw 30 μL of purified 3CL^{Pro} on ice. Estimate the protein concentration using Lambert-Beer equation, with $\epsilon_{280} = 32,890 \text{ A.U. M}^{-1} \text{ cm}^{-1}$, and dilute the 3CL^{Pro} sample to a final concentration of 10 μM using FRET Buffer.
- Setup of the reaction mixtures in the 96-well microplates:




Figure 7. Materials preparation and data processing for IC_{50} determination

(A) Schematic of 96-well plate set-up for dose response assay to determine IC_{50} .

(B) Table of reagents for each well (final DMSO concentration = 4 %).

(C) Representative experimental raw data showing the increase of fluorescence intensity as a function of reaction time.

(D) Representative dose response curve for Narlaprevir, error bars represent SD's of triplicate measurements.

- a. In each well of the 96-wells microplate add:
 - i. FRET buffer as described in Figures 7A and 7B (90 μL of FRET Buffer will be added to every well except for column 2, where a negative control experiment in the absence of enzyme will be set up using 94 μL of FRET Buffer).
 - ii. 4 μL of 10 μM 3CL^{Pro} stock solution (resulting in 400 nM 3CL^{Pro} final concentration) except for the negative control experiment in column 2.
 - iii. 0–2.0 μL inhibitor, corresponding to 0–200 μM concentration (refer to Figures 7A and 7B for details).
 - iv. Add the appropriate amount of DMSO to maintain the DMSO concentration constant in all wells.
 - v. Let the plate equilibrate for 10–15 min at room temperature (20–25 °C).
 - vi. Initiate the reaction by adding 4 μL of 5 mM peptide substrate stock (resulting in 200 μM substrate final concentration), and mix by pipetting up and down 2–3 times.

Note: The enzymatic reaction begins as soon as the substrate is added to each well. Carefully pipette the substrate into each well, moving as quickly as possible to initiate the FRET data collection (below) with minimal dead time.

- b. Place the 96-well microplate in the Infinite M1000 plate reader (TECAN) and press START on the Magellan software immediately to begin data collection.
- c. Let the reactions run for 10–20 min to measure initial enzyme velocities.

Note: Measurements are carried out at saturating substrate concentration (*i.e.* $[\text{S}] \gg K_m$ and reaction rate close to V_{max}). Final concentration of DMSO = 4 %.

Note: Inhibitors vary in their solubility. Therefore, the assay may be adjusted to use final concentrations of DMSO even lower than 1 % for appropriately soluble substrates and inhibitors.

13. Data export and analysis.
 - a. Export the data from Magellan in Microsoft Excel .xls or in .csv format.
 - b. Analyze the data using GraphPad or other kinetic analysis software:
 - i. Plot the fluorescence intensity vs. time data (Figure 7C) for each inhibitor concentration.
 - ii. Measure the slope of the linear portion of each curve (usually within the first 5 min of reaction) by performing a linear fit on that region.

Note: Each obtained slope value corresponds to the initial velocity in the presence of that inhibitor concentration (V_i). The slope determined for the positive control measurements in column 1 (Figure 7B), carried out in the absence of inhibitor, corresponds to the initial velocity of the non-inhibited 3CL^{Pro} (V_0).

- iii. Calculate the average \pm standard deviation values for the triplicate measurements carried out at each inhibitor concentration to obtain, for each concentration, an averaged V_i (or V_0) value.
- iv. Using the averaged values, calculate residual activity (%) at each inhibitor concentration using the formula

$$\text{Residual activity (\%)} = 100 \times \frac{V_i}{V_0} \quad (\text{Equation 1})$$

- v. Plot residual activity (%) data as a function of increasing inhibitor concentration and fit them by using the following equation:

$$\text{Residual activity (\%)} = \frac{100}{1+x/IC_{50}} \quad (\text{Equation 2})$$

- vi. Transform the x-axis visualization on a Log₁₀ scale to show the plotted experimental data and resulting fit following the typical sigmoidal behavior (Figure 7D).

Note: The IC₅₀ value estimated by Equation 2 corresponds to the inhibitor concentration at which residual activity is 50 %.

Alternatives: This protocol can also be adopted as a method for preliminary screening of a panel of molecules as potential inhibitors by using a single concentration (e.g., 20 μ M) of each.

NMR spectroscopy of 3CL^{Pro}

⌚ Timing: Approximately 4 h

14. Sample preparation for the ¹⁵N NMR experiment.
 - a. Gently thaw a 250 μ L aliquot of purified 0.8 mM ¹⁵N-enriched 3CL^{Pro} on ice.
 - b. Add 25 μ L deuterium oxide (CortecNet or Cambridge Isotopes) to the protein solution.
 - c. transfer the clear, transparent sample to the Shigemi 5 mm symmetrical NMR microtube (SHIGEMI, INC).
15. Data recording and analysis.
 - a. Record ¹H, ¹⁵N BEST-TROSY spectra of the protein at 298 K on a Bruker 700 MHz NMR spectrometer using 256 increments, 24 scans, and 1.2 s relaxation delay.




Figure 8. 2D ^1H , ^{15}N BEST-TROSY NMR spectrum of SARS-CoV-2 3CL^{Pro} at 298 K

(A) Wild type and (B) [Cys145Ala]- SARS-CoV-2 3CL^{Pro} mutant. No significant sample degradation was observed over a 1-week period at room temperature.

- b. Sweep widths of 13 ppm and 40 ppm and offsets of 4.7 ppm and 118 ppm are used for the ^1H and ^{15}N dimension, respectively. The total acquisition time for each experiment is 135 min.
- c. Process the spectra using the NMRPipe software.¹⁰
- d. Analyze the obtained spectra (Figure 8) using the POKY¹¹ or NMRFAM SPARKY.¹²

Alternatives: Alternative NMR software for processing include Bruker TopSpin 4.2.0, and earlier releases, the CcpNmr software suite, or software available in NMRBox software collection.^{21,22} An extensive review of NMR data processing is presented in Stern and Hoch (1996).²³

Sequence-specific resonance assignments for SARS-CoV-2 3CL^{Pro} are available in the BMRB entry 50780²⁴ and can be used to validate the spectrum and to label the observed peaks following the nearest neighbor criterion.

Crystallization and structural determination of SARS-CoV-2 3CL^{Pro}

⌚ Timing: Approximately 2 weeks

Note: SARS-CoV-2 3CL^{Pro} can be crystallized at 3–8 mg mL⁻¹ through vapor diffusion technique (hanging drop method) using a 20–60 mM ammonium acetate buffer at pH 7.0 and 20–30 % PEG4000 as a precipitant. This protocol reports a representative crystallization procedure using a 5 mg mL⁻¹ sample.

16. Crystallization of SARS-CoV-2 3CL^{Pro}.

- a. Prepare the crystallization cocktail by filling each well (0.5 mL) of the EasyXtal 15-well crystallization plate (Quiagen) with the correct volumes of 1 M ammonium acetate at pH 7.0 (AMA 1 M), 50 % PEG4000 (PEG 50 %), and MilliQ water according to the scheme in Table 1, as shown in Figure 9A.

Note: Due to different viscosities of the three solutions added, it is recommended to mix the content of each well using a micropipette.

- b. Gently thaw 50 μL of purified 5 mg mL⁻¹ 3CL^{Pro} on ice.
- c. Place a screw-in crystallization support upside down and pipette 2 μL of 3CL^{Pro} solution at its center.

Table 1. Preparation of the crystallization plate

[Ammonium acetate at pH 7.0] (mM)

P E G 4 0 0 0 (%)	20	30	40	50	60
	AMA 1 M = 10 μ L	AMA 1 M = 15 μ L	AMA 1 M = 20 μ L	AMA 1 M = 25 μ L	AMA 1 M = 30 μ L
	PEG 50 % = 200 μ L	PEG 50 % = 200 μ L	PEG 50 % = 200 μ L	PEG 50 % = 200 μ L	PEG 50 % = 200 μ L
	H ₂ O = 290 μ L	H ₂ O = 285 μ L	H ₂ O = 280 μ L	H ₂ O = 275 μ L	H ₂ O = 270 μ L
	AMA 1 M = 10 μ L	AMA 1 M = 15 μ L	AMA 1 M = 20 μ L	AMA 1 M = 25 μ L	AMA 1 M = 30 μ L
	PEG 50 % = 250 μ L	PEG 50 % = 250 μ L	PEG 50 % = 250 μ L	PEG 50 % = 250 μ L	PEG 50 % = 250 μ L
	H ₂ O = 240 μ L	H ₂ O = 235 μ L	H ₂ O = 230 μ L	H ₂ O = 225 μ L	H ₂ O = 220 μ L
	AMA 1 M = 10 μ L	AMA 1 M = 15 μ L	AMA 1 M = 20 μ L	AMA 1 M = 25 μ L	AMA 1 M = 30 μ L
	PEG 50 % = 300 μ L	PEG 50 % = 300 μ L	PEG 50 % = 300 μ L	PEG 50 % = 300 μ L	PEG 50 % = 300 μ L
	H ₂ O = 190 μ L	H ₂ O = 185 μ L	H ₂ O = 180 μ L	H ₂ O = 175 μ L	H ₂ O = 170 μ L

- d. Recover 2 μ L of the precipitation cocktail present in well A1 and pipette it onto the already deposited protein solution.
- e. Reverse the crystallization support in the original position and screw it onto the corresponding A1 well (Figure 9B).
- f. Repeat point c for all the wells present in the crystallization plate.
- g. Store the crystallization plate at 20 °C in a PKS-E100 Protein Crystallization Cabinet (RUMED) to avoid temperature oscillations as well as vibrations during crystallization.
- h. Monitor the evolution of each crystallization drop under a light microscope after regular time intervals (1 day, 2 days, 4 days, 1 week and then every week), noting any phase transition occurring (precipitation, liquid-liquid separation, crystal formation, or clear drop).

Alternatives: The same experimental setup used to crystallize native 3CL^{Pro} can be also successfully followed, after minor changes, to crystallize the protein in the presence of potential inhibitors by using either co-crystallization or crystal soaking. In the first case, appropriate amounts of test molecule should be added in the crystallization cocktail (point 16a) prior to drop deposition, in order to obtain a co-crystallization drop containing both the protein and the molecule. To perform the crystal soaking, procedure described at point 17c (see below), the protocol should be modified by exposing a crystal of native 3CL^{Pro} in a cryo-protectant droplet that also contains the appropriate amount of ligand. During soaking the ligand may diffuse into the preformed crystal through solvent channels and bind to the protein. Exposure time of the protein to the ligand needs to be optimized to achieve complete binding (usually from a few minutes to 12–16 h).

17. Crystal storage and shipping.
 - a. Fill a SPINE Sample Changer Basket (Molecular Dimensions) with ten magnetic Cryo-Vials (Molecular Dimensions) and place it in a foam Dewar topped up with liquid nitrogen.
 - b. Select the desired 3CL^{Pro} crystal candidate(s) (Figure 9C) to be shipped to the synchrotron facility among those grown in one or more wells by carefully observing the crystallization drops under the light microscope.
 - c. Open the screw-in crystallization support, fish the crystal from the crystallization drop using a cryo-loop (Molecular Dimensions) of the correct size mounted on a magnetic Cryo wand (Molecular Dimensions) (Figure 9D), and expose it in a new cryo-protectant drop composed of a slightly higher concentration of the crystallization cocktail.
 - d. Fish the crystal from the new drop and quickly flash freeze it by depositing the cryo-loop in the magnetic cryo-vial stored in the basket immersed in liquid nitrogen.
 - e. Repeat the same procedure for each crystal to be shipped, load the filled baskets in a dry shipper properly equilibrated with liquid nitrogen and ship the material to the synchrotron.

△ CRITICAL: The liquid contained in the drops tends to evaporate quickly. This phenomenon also causes the disruption of the crystal order, which in turn negatively affects data quality.

Table 2. Data collection, processing and refinement statistics (PDB code 8CDC)

	8CDC
Data collection and processing	
Wavelength (Å)	0.9762
Detector	EIGER 16M
Crystal-to-Detector distance (mm)	180.413
Oscillation angle (degrees)	0.1
Number of images	3600
Space group	C121
Unit cell (a, b, c, Å)	114.1 53.5 44.6
Resolution range (Å) ^a	55.57 - 1.54 (1.57 - 1.54)
Total number of reflections ^a	254428 (12272)
Unique reflections ^a	38825 (1885)
Multiplicity ^a	6.6 (6.5)
Completeness ^a (%)	100.0 (100.0)
R _{sym} ^{a,b} (%)	4.7 (96.1)
R _{pim} ^{a,c,2θ} (%)	3.0 (62.1)
Mean I half-set correlation CC(1/2) ^a	0.999 (0.758)
Mean I/σ(I) ^a	15.3 (1.6)
Refinement statistics	
Number of monomers in the asymmetric unit	1
R _{factor} ^d (%)	18.5
R _{free} ^d (%)	23.9
Cruickshank's DPI for coordinate error ^e based on R _{factor} (Å)	0.097
Wilson plot B-factor (Å ²)	22.7
Average all atom B-factor (Å ²)	28.1
RMS (bonds) ^d	0.013
RMS (angles) ^d	1.81
Total number of atoms	2574
Total number of water molecules	294
Solvent content (%)	36.40
Matthews Coefficient (Å ³ /Da)	1.93
Ramachandran plot	
Most favored regions (%)	97.0
Additionally allowed regions (%)	2.6
Disallowed regions (%)	0.4

^aHighest resolution bin in parentheses;

^b $R_{sym} = \sum_{hkl} \sum_j |I_j - \langle I \rangle| / \sum_{hkl} \sum_j I_j$, where I is the intensity of a reflection, and $\langle I \rangle$ is the mean intensity of all symmetry related reflections j;

^c $R_{p.i.m.} = \sum_{hkl} \{ [1/(N-1)]^{1/2} \sum_j |I_j - \langle I \rangle| \} / \sum_{hkl} \sum_j I_j$, where I is the intensity of a reflection, and $\langle I \rangle$ is the mean intensity of all symmetry related reflections j, and N is the multiplicity²⁸²⁷²⁸²⁷²⁸²⁸²⁸²⁸²⁷²⁶²⁵²⁶²⁶²⁷²⁷²⁷²⁷²⁷²⁷²⁷²⁶²⁵²⁴²³²²²²²¹¹⁸;

^dTaken from REFMAC; R_{free} is calculated using 5% of the total reflections that were randomly selected and excluded from refinement;

^e $DPI = R_{factor} \cdot D_{max} \cdot compl^{-1/3} \sqrt{\frac{N_{atoms}}{N_{refl} - N_{params}}}$, where N_{atoms} is the number of the atoms included in the refinement, N_{refl}

is the number of the reflections included in the refinement, D_{max} is the maximum resolution of reflections included in the refinement, $compl$ is the completeness of the observed data, and for isotropic refinement, $N_{params} \approx 4N_{atoms}$

Crystal fishing and flash-cooling procedures must be therefore carried out as quickly as possible. If additional time is needed between two consecutive crystal fishing steps, place the drop support back onto the corresponding well to minimize drop drying.

18. Crystallography data collection, processing, and structural determination.

Note: Diffraction data used for the structural determination of the X-ray crystal structure of native SARS-CoV-2 3CL^{Pro} were collected on one crystal at 100 K using synchrotron X-ray

Table 3. Yields of pure native 3CL^{PRO} obtained in different media

Medium	Yield (mg L ⁻¹)
Super broth	122
LB	82
MJ9 (minimal medium)	28

Protein was spectrophotometrically quantified assuming MW = 33,796 Da and $\epsilon_{280} = 32,890 \text{ A.U. M}^{-1} \text{ cm}^{-1}$

radiation at the EMBL P13 beamline of the Petra III storage ring, c/o DESY, Hamburg (Germany).²⁵

- Record diffraction data using the collection strategy suggested by the MXCuBE software.^{13,14}
- Process X-ray diffraction images by using XDS.¹⁵
- Scale and reduce processed data by using AIMLESS,¹⁷ included in the suite CCP4.¹⁶
- Carry out structural refinement and model building using the REFMAC5 software¹⁸ and COOT,²⁶ respectively, both included in the suite CCP4.¹⁶ Use the structure factor amplitudes obtained by AIMLESS as experimental data and the PDB entry 7ALH as a starting model until refinement convergence. A representative picture of the final X-ray crystal structure of 3CL^{PRO} obtained in this work (and deposited in the Protein Data Bank under PDB code 8CDC) is shown in Figure 10. Data collection, processing and refinement statistics are reported in Table 2.

Alternatives: Other data collection facilities may be used. Other software can be used for analyzing the crystallographic data, such as PHENIX.²⁷

Note: Crystals grown in experimental conditions described belong to space group C121, as does the starting model selected for the refinement (PDB code 7ALH). If the space groups of




Figure 9. Experimental set up for the crystallization of SARS-CoV-2 3CL^{PRO} by hanging drop and crystal recovery and storage

(A and B) (A) Preparation of the crystallization plate, and (B) crystallization drops deposition.

(C and D) (C) 3CL^{PRO} typically grows as bundles of sword-shaped crystals which (D) are fished using cryo-loops, flash-cooled in liquid nitrogen for storage, and shipped to the synchrotron facility for data collection.




Figure 10. Active site region the X-ray crystal structure of SARS-CoV-2 3CL^{Pro}

The atomic model is shown superimposed on the final 2Fo - Fc electron density map contoured at 1 σ and colored gray. The carbon, nitrogen, oxygen, and sulfur atoms are gray, blue, red, and yellow, respectively. The figure is generated using PyMol (The PyMOL Molecular Graphics System, v. 1.8 Schrödinger, LLC.).

experimental data and the starting model do not match, a molecular replacement procedure must be carried out to obtain initial crystallographic phases to be used in the following refinement and model building steps.

Sample quality control

In addition to crystallography and NMR data, if available, sample quality control includes SDS-PAGE, MALDI-TOF mass spectrometry, and specific enzyme activity measurements. The expected specific activity is 55 Units / μ mole of enzyme, where a Unit is 1.0 μ mol substrate cleaved / min, assayed using the FRET buffer conditions, as described in "materials and equipment" section, using 0.4 μ M 3CL^{Pro}, 200 μ M DABCYL-KTSAVLQ/SGFRKME-EDANS substrate at 25°C.

EXPECTED OUTCOMES

This paper provides a detailed protocol for the expression, purification, and production of tens- of-milligram quantities of highly pure (> 99 %) SARS-CoV-2 3CL^{Pro} with native N and C termini and in a biochemically-active form, as demonstrated by the bio-structural characterization and activity assays reported. Yields of pure SARS CoV-2 3CL^{Pro} obtained, as high at \sim 120 mg / mL, depend on the medium used for fermentation (Table 3).

3CL^{Pro} samples produced by this protocol are stable in NMR measurements for up to 1 week at room temperature and provide excellent HSQC-TROSY spectra suitable for structure, dynamic, ligand binding studies. Moreover, the set-up crystallization protocol provides high quality 3CL^{Pro} X-ray crystal structures at resolution better than 2 Å, thus suitable for biophysical and drug discovery research using X-ray crystallography. In the paper a standardized assay for estimating IC₅₀ for new potential 3CL^{Pro} inhibitors is also presented, thus complementing structural information with biochemical insights of catalytic and inhibition mechanism of this target.

LIMITATIONS

A percentage of the obtained 3CL^{Pro} crystals (\sim 20 %) do not diffract at high resolution. The protocol section that describes the characterization of enzyme inhibition by FRET assays assumes solubility of

the screening molecules in 4 % DMSO. For potent inhibitors which can be studied at lower concentrations, lower DMSO concentrations may be used if solubility allows. Each potential inhibitor should be evaluated for solubility prior to setting up the assay. Moreover, IC_{50} obtained from the biochemical assays is a preliminary measurement of inhibitor potency. A more accurate investigation of the inhibition should be provided by determining thermodynamic and kinetic parameters for the inhibition mechanism.

TROUBLESHOOTING

Problem 1

Following the 12–16 h cell growth on LB Petri dishes, no colonies are visible.

Potential solution

After checking each step of the transformation procedure was carried out correctly (e.g., correct antibiotic, plasmid concentration, duration of heat shock), the concentration of the cell suspension can be increased via centrifugation. At the end of section 13f in the “before you begin” section, centrifuge the cell suspension for 1 min at $2500 \times g$ and decant half of the supernatant. Resuspend the pellet in the remaining supernatant for a more concentrated cell culture.

Problem 2

During cell lysis (step 2. Cell lysis and soluble extract recovery), heat generated by sonication and French Press, can cause aggregation or misfolding of 3CL^{Pro}.

Potential solution

Sonication should be done in short blasts, with the sample cooled in an ice bath in a container with good thermal conductivity. Lysis through French press should be also performed constantly maintaining the sample in an ice bath. Microfluidizer or homogenizer can be used in place of the aforementioned techniques, though lysis by these methods may be less complete resulting in lower final yields.

Problem 3

The expressed His₆-SUMO-3CL^{Pro} undergoes a first step of purification through an immobilized nickel affinity chromatography (IMAC) (step 3). Due to the high amount of protein expressed from 1 L of cell culture and a column capacity of usually 10 mg mL^{-1} of resin (indicating that a 5 mL column will efficiently bind up to 50 mg of protein), column overloading issues can be encountered. Such an event is clearly detectable by running an SDS-PAGE of the flow-through, which would show large amounts of unbound His₆-SUMO-3CL^{Pro}.

Potential solution

It is strongly recommended to split the soluble lysate in two or more aliquots and perform the IMAC in successive steps. As an alternative, two or more columns placed in series can also be used to achieve high amounts of protein recovery and speed up the purification procedure.

Problem 4

Uncontrolled cleavage of His₆-SUMO tag from His₆-SUMO-3CL^{Pro} probably by trace amounts of native *E. coli* proteases has been observed at the end of step 3 (IMAC), just before the addition of His₈-MBP-SUMO^{Pro}.

Potential solution

Uncontrolled cleavage can be minimized by carrying out cleavage in imidazole immediately after purifying the His₆-SUMO-3CL^{Pro}, with low to no negative effects on the production of pure and native 3CL^{Pro}, or by on-column cleavage, as suggested by this protocol.

Problem 5

The on-column cleavage protocol set up to perform IMAC purification, SUMO protease cleavage of His₆-SUMO tag from 3CL^{Pro}, and removal of cleaved His₆-SUMO tag and His₈-MBP-SUMO^{Pro} from native 3CL^{Pro} in one step is more efficient and quicker. However, in some cases subsequent elution of the proteins reveals incomplete cleavage of the His₆-SUMO-3CL^{Pro}, resulting in reduced final yields.

Potential solution

His₆-SUMO-3CL^{Pro} fusion that elutes from the column can be dialyzed to remove imidazole and undergo a second on-column cleavage procedure.

Problem 6

Incomplete cleavage of the His₆-SUMO-3CL^{Pro} fusion is observed (*i.e.*, full target fusion is visible in the final elution from the IMAC).

Potential solution

Collect the elution fractions containing the His₆-SUMO-3CL^{Pro} fusion, perform an 12–16 h dialysis/buffer exchange and simultaneous cleavage in the presence of His₈-MBP-SUMO^{Pro} and repeat an IMAC step by collecting the flow-through.

Problem 7

No increasing fluorescence signal is observed after combining enzyme and substrate during the FRET Assay control (described in step 12).

Potential solution

3CL^{Pro} is very sensitive to oxidation and can lose enzymatic activity due to oxidation of the active site cysteine even though reducing agents are present. During purification it is advantageous to work quickly and at 4°C as much as possible to minimize this effect. Samples can also be protected from O₂ by storing them under nitrogen or argon gas.

Problem 8

Quality of ¹⁵N HSQC spectrum is not optimal.

Potential solution

Prepare a ²H, ¹⁵N -enriched sample to obtain a higher quality ¹⁵N HSQC spectrum.

Problem 9

Reservoir solutions for crystallization trials includes mixing of ammonium acetate, 50 % PEG4000, and water, which have different viscosity and may prevent the final solution to be homogeneous.

Potential solution

Carefully mix the reservoir solution present in each well using a micropipette.

RESOURCE AVAILABILITY

Lead contact

Further information and requests for resources and reagents should be directed to and will be fulfilled by the lead contact, G.T. Montelione, monteg3@rpi.edu.

Materials availability

Plasmids generated in this study have been deposited to AddGene, and are available under the Uniform Biological Material Transfer Agreement (“UBMTA”):

Plasmid pGTM_CoV2_NSP5_004_SUMO. AddGene ID: 190062.

Plasmid pGTM_CoV2_NSP5_C145A_001_SUMO. AddGene ID: 192349.

Plasmid pGTM_YR375_SUMO_Protease_001. AddGene ID: 190063.

Plasmid sequencing data is available at: <https://github.rpi.edu/RPIBioinformatics/SARS-CoV-2-3CLpro>.

Data and code availability

Datasets generated during this study, including the NMR free-induction decay (FID), and crystallographic processed and reduced reflections, are available at: Github: RPI Bioinformatics SARS-CoV-2-3CL^{PRO} (<https://github.rpi.edu/RPIBioinformatics/SARS-CoV-2-3CLpro>) The Protein Data Bank accession number for the crystal structure of CoV-2-3CL^{PRO} described here is 8CDC

ACKNOWLEDGMENTS

We thank N. Dube, B. Harish, L. Klang, T. Ramelot, B. Shurina, and L. Spaman for helpful discussions, suggestions, and comments on the manuscript. This research was supported by National Institutes of Health grants 1R35-GM141818 to G.T.M. and 1U19 AI171443-01 (to S. Chanda, A. Garcia-Sastre, and G.T.M. as coPI).

AUTHOR CONTRIBUTIONS

These authors contributed equally, L.M., R.G.-C.; these senior authors contributed equally, S.C., G.T.M.; conceptualization, L.M., R.G.-C., K.B., A.J., T.B.A., C.A.R., S.C., A.D.F., G.T.M.; biochemical investigations, L.M., R.G.-C., K.B., A.J., A.D.F.; NMR investigations, S.C.; crystallography investigations, L.M.; visualization, L.M., R.G.-C., S.C., T.B.A., G.T.M.; writing, reviewing, and editing, L.M., R.G.-C., K.B., A.J., A.D.F., T.B.A., C.A.R., S.C., G.T.M.; funding acquisition, S.C., G.T.M.

DECLARATION OF INTERESTS

G.T.M. is a founder and advisor to Nexomics Biosciences, Inc. This does not represent a conflict of interest with respect to this study.

REFERENCES

- Bafna, K., White, K., Harish, B., Rosales, R., Ramelot, T.A., Acton, T.B., Moreno, E., Kehrer, T., Miorin, L., Royer, C.A., et al. (2021). Hepatitis C virus drugs that inhibit SARS-CoV-2 papain-like protease synergize with remdesivir to suppress viral replication in cell culture. *Cell Rep.* 35, 109133. <https://doi.org/10.1016/j.celrep.2021.109133>.
- Kim, J.C., Spence, R.A., Currier, P.F., Lu, X., and Denison, M.R. (1995). Coronavirus protein processing and RNA synthesis is inhibited by the cysteine proteinase inhibitor E64d. *Virology* 208, 1–8. <https://doi.org/10.1006/viro.1995.1123>.
- Ziebuhr, J., Snijder, E.J., and Gorbalenya, A.E. (2000). Virus-encoded proteinases and proteolytic processing in the Nidovirales. *J. Gen. Virol.* 81, 853–879. <https://doi.org/10.1099/0022-1317-81-4-853>.
- Anand, K., Palm, G.J., Mesters, J.R., Siddell, S.G., Ziebuhr, J., and Hilgenfeld, R. (2002). Structure of coronavirus main proteinase reveals combination of a chymotrypsin fold with an extra alpha-helical domain. *EMBO J.* 21, 3213–3224. <https://doi.org/10.1093/emboj/cdf327>.
- Ertem, F.B., Guven, O., Buyukdag, C., Gocencler, O., Ayan, E., Yuksel, B., Gul, M., Usta, G., Cakilkaya, B., Johnson, J.A., et al. (2022). Protocol for structure determination of SARS-CoV-2 main protease at near-physiological-temperature by serial femtosecond crystallography. *STAR Protoc.* 3, 101158. <https://doi.org/10.1016/j.xpro.2022.101158>.
- Ihssen, J., Faccio, G., Yao, C., Sirec, T., and Spitz, U. (2021). Fluorogenic in vitro activity assay for the main protease M(pro) from SARS-CoV-2 and its adaptation to the identification of inhibitors. *STAR Protoc.* 2, 100793. <https://doi.org/10.1016/j.xpro.2021.100793>.
- Acton, T.B., Xiao, R., Anderson, S., Aramini, J., Buchwald, W.A., Ciccocanti, C., Conover, K., Everett, J., Hamilton, K., Huang, Y.J., et al. (2011). Preparation of protein samples for NMR structure, function, and small-molecule screening studies. *Methods Enzymol.* 493, 21–60. <https://doi.org/10.1016/B978-0-12-381274-2.00002-9>.
- Jansson, M., Li, Y.C., Jendeborg, L., Anderson, S., Montelione, G.T., and Nilsson, B. (1996). High-level production of uniformly (15)N- and (13)C-enriched fusion proteins in *Escherichia coli*. *J. Biomol. NMR* 7, 131–141. <https://doi.org/10.1007/BF00203823>.
- Bergmans, H.E., Van Die, I.M., and Hoekstra, W.P. (1981). Transformation in *Escherichia coli*: stages in the process. *J. Bacteriol.* 146, 564–570. <https://doi.org/10.1128/jb.146.2.564-570.1981>.
- Delaglio, F., Grzesiek, S., Vuister, G.W., Zhu, G., Pfeifer, J., and Bax, A. (1995). NMRPipe: a multidimensional spectral processing system based on UNIX pipes. *J. Biomol. NMR* 6, 277–293. <https://doi.org/10.1007/BF00197809>.
- Lee, W., Rahimi, M., Lee, Y., and Chiu, A. (2021). POKY: a software suite for multidimensional NMR and 3D structure calculation of biomolecules. *Bioinformatics* 37, 3041–3042. <https://doi.org/10.1093/bioinformatics/btab180>.
- Lee, W., Tonelli, M., and Markley, J.L. (2015). NMRFAM-SPARKY: enhanced software for biomolecular NMR spectroscopy. *Bioinformatics* 31, 1325–1327. <https://doi.org/10.1093/bioinformatics/btu830>.
- Gabadinho, J., Beteva, A., Guizarro, M., Rey-Bakaikoa, V., Spruce, D., Bowler, M.W., Brockhauser, S., Flot, D., Gordon, E.J., Hall, D.R., et al. (2010). MxCuBE: a synchrotron beamline control environment customized for macromolecular crystallography experiments.

- J. Synchrotron Radiat. 17, 700–707. <https://doi.org/10.1107/S0909049510020005>.
14. Oscarsson, M., Beteva, A., Flot, D., Gordon, E., Guijarro, M., Leonard, G., McSweeney, S., Monaco, S., Mueller-Dieckmann, C., Nanao, M., et al. (2019). MXCuBE2: the dawn of MXCuBE Collaboration. *J. Synchrotron Radiat.* 26, 393–405. <https://doi.org/10.1107/S1600577519001267>.
 15. Kabsch, W. (2010). Xds. *Acta Crystallogr. D* 66, 125–132. <https://doi.org/10.1107/S0907444909047337>.
 16. Winn, M.D., Ballard, C.C., Cowtan, K.D., Dodson, E.J., Emsley, P., Evans, P.R., Keegan, R.M., Krissinel, E.B., Leslie, A.G.W., McCoy, A., et al. (2011). Overview of the CCP4 suite and current developments. *Acta Crystallogr. D Biol. Crystallogr.* 67, 235–242. <https://doi.org/10.1107/s0907444910045749>.
 17. Evans, P.R., and Murshudov, G.N. (2013). How good are my data and what is the resolution? *Acta Crystallogr. D Biol. Crystallogr.* 69, 1204–1214. <https://doi.org/10.1107/S0907444913000061>.
 18. Murshudov, G.N., Vagin, A.A., and Dodson, E.J. (1997). Refinement of macromolecular structures by the maximum-likelihood method. *Acta Crystallogr. D* 53, 240–255. <https://doi.org/10.1107/S0907444996012255>.
 19. Emsley, P., and Cowtan, K. (2004). Coot: model-building tools for molecular graphics. *Acta Crystallogr. D Biol. Crystallogr.* 60, 2126–2132. <https://doi.org/10.1107/S0907444904019158>.
 20. Liu, Y., Kati, W., Chen, C.-M., Tripathi, R., Molla, A., and Kohlbrenner, W. (1999). Use of a fluorescence plate reader for measuring kinetic parameters with inner filter effect correction. *Anal. Biochem.* 267, 331–335.
 21. Skinner, S.P., Fogh, R.H., Boucher, W., Ragan, T.J., Mureddu, L.G., and Vuister, G.W. (2016). CcpNmr AnalysisAssign: a flexible platform for integrated NMR analysis. *J. Biomol. NMR* 66, 111–124. <https://doi.org/10.1007/s10858-016-0060-y>.
 22. Maciejewski, M.W., Schuyler, A.D., Gryk, M.R., Moraru, I.I., Romero, P.R., Ulrich, E.L., Eghbalnia, H.R., Livny, M., Delaglio, F., and Hoch, J.C. (2017). NMRbox: a resource for biomolecular NMR computation. *Biophys. J.* 112, 1529–1534. <https://doi.org/10.1016/j.bpj.2017.03.011>.
 23. Stern, A.S., and Hoch, J.C. (1996). *NMR Data Processing* (Wiley).
 24. Cantrelle, F.X., Boll, E., Brier, L., Moschidi, D., Belouzard, S., Landry, V., Leroux, F., Dewitte, F., Landrieu, I., Dubuisson, J., et al. (2021). NMR spectroscopy of the main protease of SARS-CoV-2 and fragment-based screening identify three protein hotspots and an antiviral fragment. *Angew. Chem. Int. Ed. Engl.* 60, 25428–25435. <https://doi.org/10.1002/anie.202109965>.
 25. Cianci, M., Bourenkov, G., Pompidor, G., Karpics, I., Kallio, J., Bento, I., Roessle, M., Cipriani, F., Fiedler, S., and Schneider, T.R. (2017). P13, the EMBL macromolecular crystallography beamline at the low-emittance PETRA III ring for high- and low-energy phasing with variable beam focusing. *J. Synchrotron Radiat.* 24, 323–332. <https://doi.org/10.1107/S1600577516016465>.
 26. Emsley, P., and Cowtan, K. (2004). Coot: model-building tools for molecular graphics. *Acta Crystallogr. D* 60, 2126–2132. <https://doi.org/10.1107/S0907444904019158>.
 27. Liebschner, D., Afonine, P.V., Baker, M.L., Bunkóczi, G., Chen, V.B., Croll, T.I., Hintze, B., Hung, L.-W., Jain, S., McCoy, A.J., et al. (2019). Macromolecular structure determination using X-rays, neutrons and electrons: recent developments in Phenix. *Acta Crystallogr. D Struct. Biol.* 75, 861–877. <https://doi.org/10.1107/s2059798319011471>.
 28. Weiss, M.S. (2001). Global indicators of X-ray data quality. *J. Appl. Crystallogr.* 34, 130–135. <https://doi.org/10.1107/S0021889800018227>.