cambridge.org/par

## **Research Article**

**Cite this article:** Souza DC *et al* (2020). Use of polyclonal IgY antibodies to detect serum immune complexes in patients with active hookworm infection. *Parasitology* **147**, 715–720. https://doi.org/10.1017/ S003118202000220

Received: 16 October 2019 Revised: 14 January 2020 Accepted: 5 February 2020 First published online: 13 February 2020

#### Key words:

Active infection; hookworm; immune complexes; IgY

### Author for correspondence:

Julia Maria Costa-Cruz, E-mail: costacruz@ufu.br

## Use of polyclonal IgY antibodies to detect serum immune complexes in patients with active hookworm infection

Dayane Costa de Souza<sup>1</sup>, Lucas Silva de Faria<sup>1</sup>, José Eduardo Neto de Sousa<sup>1</sup>, Camila de Almeida Lopes<sup>1</sup>, Vanessa da Silva Ribeiro<sup>1</sup>, Vivian Jordania da Silva<sup>2</sup>, Raphaella Paula Ribeiro<sup>3</sup>, Élida Mara Leite Rabelo<sup>2</sup>, Veridiana Melo Rodrigues Ávila<sup>4</sup>, Álvaro Ferreira Júnior<sup>3</sup> and Julia Maria Costa-Cruz<sup>1</sup>

<sup>1</sup>Laboratório de Diagnóstico de Parasitoses, Instituto de Ciências Biomédicas, Universidade Federal de Uberlândia – UFU, Av. Amazonas s/n, bloco 4C, Uberlândia, Minas Gerais, 38400-902, Brasil; <sup>2</sup>Laboratório de Parasitologia Molecular, Instituto de Ciências Biológicas, Universidade Federal de Minas Gerais – UFMG, Belo Horizonte, Minas Gerais, 31270-901, Brasil; <sup>3</sup>Programa de Pós-Graduação em Sanidade e Produção Animal nos Trópicos, Universidade de Uberaba – Uniube, Campus Aeroporto. Av. Nenê Sabino, sala 2D05, Uberaba, Minas Gerais, 38055-500, Brasil and <sup>4</sup>Laboratório de Bioquímica e Toxinas Animais, Instituto de Genética e Bioquímica, Universidade Federal de Uberlândia – UFU, Av. Pará 1720, Uberlândia, Minas Gerais, 38400-902, Brasil

#### Abstract

Definitive diagnosis of hookworm infection is usually based on the microscopic detection of eggs in a stool sample; however, several cases display a low or irregular egg output. Serodiagnosis can be a useful tool to identify these cases, but conventional tests do not differentiate past from active infections. The aim of this study was to obtain and apply egg yolk polyclonal immunoglobulin (IgY) antibodies to detect immune complexes (ICs) in serum samples from patients infected with hookworm. Hens were immunized with Ancylostoma ceylanicum saline extract, their eggs were collected and then IgY antibodies were extracted and purified. Antibody purity was tested by 12% sodium dodecyl sulphate polyacrylamide gel electrophoresis and specificity was assessed by immunoblotting and immunofluorescence. IgY production was evaluated by kinetics enzyme-linked immunosorbent assay (ELISA). Sandwich ELISA tested the ability of IgY to detect ICs in serum samples, from which diagnostic parameters were calculated. Antibody responses increased steadily from day 7 to 42. In the immunoblotting assay, IgY recognized two protein complexes. The immunofluorescence assay showed no staining in control samples. The sandwich ELISA presented a very high diagnostic value, with a sensitivity of 90% and a specificity of 86.7%. Our pioneer strategy highlights the potential use of egg yolk IgY as a diagnostic test to detect active hookworm infection.

## Introduction

Soil-transmitted helminthiasis (STH) is the most prevalent neglected tropical disease caused by several species of hookworm (*Necator americanus, Ancylostoma duodenale* and *Ancylostoma ceylanicum*), roundworm (*Ascaris lumbricoides*) and whipworm (*Trichuris trichiura*). Approximately 1.4 billion people worldwide are estimated to be infected with at least one of the main STHs, and more than 740 million people are infected with hookworm (Pullan *et al.*, 2014; Prieto-Pérez *et al.*, 2016). Infections with *N. americanus* and *A. duodenale* remain a major public health problem in several low- and middle-income countries (George *et al.*, 2016).

Diagnosis is mainly based on microscopic examination of faeces, but this method presents sensitivity problems in cases with low parasitic burden, intermittent shedding of eggs, single-sex infections or sexually immature forms (Elsemore *et al.*, 2017). In these cases, other diagnostic tests such as enzyme-linked immunosorbent assay (ELISA) must be applied.

The IgG assay can present false-positive reactions with other parasitic helminths and may not distinguish between past and active infections. To overcome this problem, diagnostic tests based on chicken egg yolk immunoglobulin (IgY) are being developed (Nie *et al.*, 2014). The use of egg yolk IgY was first proposed by Leslie and Clem (1969). IgY reacts with more epitopes on a mammalian antigen, amplifying the signal and improving the sensitivity of ELISA by 10–100 times (Ohnishi *et al.*, 2000). IgY antibodies have many advantages over IgG antibodies, including cost-effectiveness, strong avidity, scalable productivity, low assay background and high yield (Li *et al.*, 2015).

Due to the advantages mentioned above, polyclonal IgY antibodies have been widely used in applications related to parasitic diseases such as toxoplasmosis (Ferreira Júnior *et al.*, 2012), schistosomiasis (Lei *et al.*, 2012), trichinellosis (Liu *et al.*, 2013), strongyloidiasis (Faria *et al.*, 2019) and ascariasis (Lopes *et al.*, 2019). In view of this, we developed the first sensitive, specific and noninvasive IgY sandwich-based ELISA to detect immune complexes (ICs) in serum, indicating active hookworm infection.

© The Author(s), 2020. Published by Cambridge University Press



## Material and methods

## Serum samples

Serum samples from 90 individuals were tested and divided into three groups: group 1 consisted of 30 patients infected with hookworm; group 2 included 30 patients who had other parasitic diseases, including infections with *Giardia duodenalis* (n = 3), *Entamoeba histolytica/díspar* (n = 3), *Ascaris* sp. (n = 7), *Enterobius vermicularis* (n = 4), *Hymenolepis nana* (n = 3), *Schistosoma mansoni* (n = 4) and *Taenia* sp. (n = 6); and group 3 consisted of 30 apparently healthy volunteers. Individuals from all groups provided three faecal samples that were analysed according to Baermann (1917), Lutz (1919) and Moraes (1948).

### Preparation of Ancylostoma ceylanicum saline extract

Total saline extract (SE) from *A. ceylanicum* was produced using approximately 100 000 infective larvae (L3), as previously described by Gonzaga *et al.* (2013). The protein concentration of SE was determined according to Lowry *et al.* (1951), then proteins were subjected to 12% sodium dodecyl sulphate polyacrylamide gel electrophoresis (SDS-PAGE) (Laemmli, 1970) for the characterization of the antigenic profiles after silver staining (Friedman, 1982).

## Production of polyclonal IgY antibodies

Four 25-week-old laying hens (*Gallus gallus domesticus*, Hisex lineage) were immunized in the pectoral muscle with SE using Freund's adjuvant, as described previously by Schwarzkopf *et al.* (2001). Each hen was firstly immunized with  $100 \mu g$  s.E. diluted in  $250 \mu L$  of phosphate-buffered saline (PBS; 0.01 M, pH 7.2) and  $250 \mu L$  of complete Freund's adjuvant (Sigma Aldrich, USA). Hens then received three booster immunizations at 14-day intervals. These were performed with SE, as in the first immunization, but used incomplete Freund's adjuvant (Sigma Aldrich, USA). The control group was immunized with PBS. All testing groups contained two hens each.

The hens were maintained under standard specific pathogenfree conditions and were monitored daily for adverse effects and inflammatory processes. Chicken eggs were collected daily from before the first immunization until after the last boost.

## Characterization of polyclonal IgY antibodies

Polyclonal IgY antibodies were extracted from egg yolk using the water-soluble method described by Akita and Nakai (1993). IgY was purified by affinity chromatography (HiTrapIgY Purification HP 5 mL; GE Healthcare, USA) at a flow rate of 3 mL min<sup>-1</sup>. IgY antibodies were eluted using a linear gradient of 0.02 M sodium phosphate buffer (pH 7.2) in an ÄKTA prime liquid chromatography system (GE Healthcare, USA), then the eluted fraction was dialysed with ultrapure water (Amicon YM 30; Sigma-Aldrich Co., USA) using a 30 kDa cut-off membrane, and the protein concentration was determined at 280 nm (BioDrop, UK).

IgY was analysed by 12% SDS-PAGE under reducing conditions (5% 2-mercaptoethanol; Sigma-Aldrich Co., USA) and the gel was stained with Coomassie Blue (Vetec, Brazil). Purified IgY from eggs produced on day 42 was tested for target specificity by immunoblotting, where  $250 \,\mu\text{g}$  of SE diluted in PBS (final volume of  $200 \,\mu\text{L}$ ) was transferred to a nitrocellulose membrane after 12% SDS-PAGE according to Towbin *et al.* (1979). Strips were cut, and non-specific binding sites were blocked with PBS containing 0.05% Tween 20 and 5% skimmed milk (PBS-T-M) for 2 h at room temperature. Membranes were individually incubated

with anti-hookworm IgY (15µg) diluted in 1% PBS-T-M (final volume of  $500 \,\mu$ L) overnight at 4°C. Then,  $500 \,\mu$ L of HRP-conjugated anti-chicken IgY (whole molecule: Sigma-Aldrich Co., USA) diluted (1: 15000) in 1% PBS-T-M was used as the second antibody. Specific IgY binding was visualized by staining the membranes with 3,3-diaminobenzidine tetrahydrochloride (DAB SigmaFast tablets; Sigma-Aldrich Co., USA) and 30% hydrogen peroxide (Merck, Brazil) diluted in trisbuffered saline (0.02 M, pH 7.4). Between each step, the strips were washed six times with PBS-T. The reaction was stopped by adding distilled water, and positive reactions were identified by the appearance of defined bands. The relative molecular weights of the recognized protein fractions were determined by comparison with molecular markers (REC001; RBC, Taiwan).

#### Kinetic evaluation of IgY production

The reactivity and kinetics of anti-hookworm IgY production were monitored by indirect ELISA. In this procedure, samples were used at a final volume of  $50 \,\mu\text{L}$  per well. Briefly, 96-well lowaffinity polystyrene microplates (Greiner Bio-one, Germany) were coated with SE (10  $\mu$ g mL<sup>-1</sup>) diluted in carbonate bicarbonate buffer (0.06 M, pH 9.6) and incubated overnight at 4°C. After incubation, microplates were washed three times for 5 min each with PBS-T then blocked with 1% PBS-T-M for 1 h at 37°C. IgY  $(2 \mu \text{g mL}^{-1} \text{ in PBS-T-M})$  was added in duplicate and incubated for 1 h at 37°C. Washing was followed by incubation with anti-IgY antibodies conjugated to horseradish peroxidase (Sigma-Aldrich Co., USA) diluted 1:10 000 in 1% PBS-T-M for 45 min at 37°C. After washing, the reaction was revealed by o-phenylenediamine (Sigma-Aldrich Co., USA) with 0.03% hydrogen peroxide (Merck, Brazil) diluted in 0.1 M citrate phosphate buffer (pH 5.5) for 15 min. The reaction was stopped by adding  $25 \mu L$  of 2N H<sub>2</sub>SO<sub>4</sub> (Vetec, Brazil). The optical density (OD) was determined at 492 nm in an ELISA reader (Thermo plate, China). Data were expressed as the ELISA index (EI) as follows: EI = OD/cut-off of the negative controls (control group group immunized with PBS) + three standard deviations. To establish the cut-off values for the control group, the OD values from the first 6 weeks were used (Faria et al., 2019; Lopes et al., 2019).

#### Immunofluorescence assay

IgY was tested for target specificity by immunofluorescence. Prior to staining, parasites were added to tissue freezing medium (Sakura, Netherlands) and frozen at -20°C. Ancylostoma ceylanicum larvae were sectioned at  $2\mu$ m thickness with a cryomicrotome (CM 1850 UV; Leica, Germany), placed on glass slides and air dried for 30 min. Sections were incubated with antihookworm IgY  $(35 \mu g)$  diluted in PBS (final volume of  $50 \mu L$ ) for 30 min at 37°C. After six washes of 5 min each with PBS, sections were incubated at 37°C for 30 min with 50  $\mu$ L of the secondary antibody, rabbit anti-chicken IgY conjugated with fluorescein isothiocyanate (FITC; Sigma-Aldrich Co., USA) diluted 1:300 in 3% Evans blue (Vetec, Rio de Janeiro, Brazil) for counterstaining. Possible cross-reactivity was tested using anti-hookworm IgY antibodies on sections of Strongyloides venezuelensis filariform larvae under the same conditions. For the negative control, sections were incubated with IgY from control hens. After washing, the slides were mounted with glycerol/PBS (pH 9.0) and coverslips. FITC fluorescence emission was detected using an LSM 510 confocal microscope (Meta; Carl Zeiss, Germany).

# IgY-based sandwich ELISA to detect immune complexes in serum samples

Detection of circulating ICs by sandwich ELISA using anti-hookworm IgY was performed as previously described (Gonçalves et al., 2012) with some modifications. Preliminary tests with all the reagents were carried out to determine the high signal/noise ratio conditions and samples were diluted to a final volume of  $50 \,\mu$ L. Briefly, high-binding 96-well microplates (Corning-Costar, USA) were coated with purified IgY ( $10 \mu g$ mL<sup>-1</sup>; collected during the sixth week) in 0.06 M carbonate bicarbonate buffer (pH 9.6) and incubated overnight at 4°C. Serum samples were diluted 1:80 in PBS-T and incubated for 45 min at 37°C. Subsequently, goat anti-human IgG (Fc specific) HRP-conjugated antibody (Sigma-Aldrich Co., USA) diluted 1:2000 in PBS-T containing 1% bovine serum albumin was added and incubated for 45 min at 37°C. Between each incubation step, plates were washed three times for 5 min each with PBS-T. The following steps of the reaction were as previously described. Serum samples were tested in duplicate. The cut-off point was established from the receiver operating characteristic (ROC) curve using OD values from the negative and other parasite groups as negative controls.

## Statistical analysis

Analyses were performed using GraphPad software package 6.0 (GraphPad Software Inc., San Diego, USA). Statistical analysis was performed by calculating the lower limits of positivity (cut-off value) established for optimal sensitivity (Se) and specificity (Sp) using the ROC curve. The diagnostic parameters calculated were area under the curve (AUC), Se, Sp and positive likelihood ratio (LR+) according Ribeiro et al. (2013). Confidence intervals (CI) of 95% were provided for Se, Sp and AUC statistical calculations.

### Results

## Antigen preparation and IgY collection

The IgY yield at day 42 was  $4750 \,\mu g \text{ mL}^{-1}$  before purification. Samples were purified by affinity chromatography and pooled to give an adequate volume and concentration for the following experiments. The pooled samples reached a protein concentration of  $700 \,\mu g \text{ mL}^{-1}$ . In experiments with immunized hens, the final IgY concentration was found to vary. To overcome this variation, a standard concentration was established for each test.


After the second booster immunization, the levels of antihookworm IgY reached a peak (Fig. 1A). The antibody response increased steadily from day 7 to 42. The electrophoretic profiles of the s.E. showed bands with apparent molecular weights varying from 7 to 240 kDa (Fig 1B-2).

SDS-PAGE analysis of the recovered polyclonal IgY after HPLC fractionation confirmed the integrity of the antibodies, meaning they were not degraded (Fig 1B-3). This finding indicates that hens, as a host for the production of anti-hookworm IgY, show a good ability to rapidly and efficiently generate specific polyclonal IgY in a non-invasive way.

Immunoblot analysis indicated that IgY recognized two protein complexes of high (200 kDa) and medium (40–50 kDa) molecular weights (Fig. 1B-4). These bands were not recognized by control group IgY (Fig. 1B-5).

## Immunofluorescence

To confirm reactivity and determine recognition patterns, IgY antibodies were separately incubated with *A. ceylanicum* and *S. venezuelensis*, and the bound antibodies were visualized by



**Fig. 1.** Anti-hookworm IgY characterization. A: Kinetics of IgY antibodies production. \* indicate immunizations. Electrophoretical profiles in SDS-PAGE 12%. B: 1 – molecular weight, 2 – *A. ceylanicum* saline extract, 3 – IgY antibodies after HPLC fractionation. Immunoblotting to evaluate antigenic recognition of IgY antibodies from 4 – hens immunized with *A. ceylanicum* s.E. and 5 – hens immunized with PBS.

indirect fluorescence. The staining pattern observed in *A. ceylanicum* sections was not uniform, but rather distributed in small, interspersed textural patches with no apparent regularity. No fluorescence was detected for *S. venezuelensis* sections or the negative control (Fig. 2).


# IgY-based sandwich ELISA for the detection of immune complexes

Positivity in group 1 reached 90%. The cross-reactivity observed in group 2 (26.6%, n = 8) was due to samples harbouring other parasites including *G. duodenalis* (1/8), *E. histolytica/ dispar* (1/8), *E. vermicularis* (2/8), *H. nana* (1/8), *S. mansoni* (2/8) and *Taenia* sp. (1/8), whereas no reactivity was observed in group 3 (Fig. 3A).

Test performance, indicated by AUC (0.9539), almost reached the maximum value (1.00) of efficiency. The LR conferred a very high diagnostic value for this test (6.75), representing a high probability of a true positive case of active hookworm infection, and ROC curve analysis (Fig. 3B) confirmed high-efficiency values for the ELISA assays.

### Discussion

Diagnostic tests with high sensitivity for detecting neglected helminth tropical diseases are necessary; however, the analysis of faecal specimens by microscopy remains the most commonly used technique (Truscott *et al.*, 2016). This test is limited by egg release and the distribution of eggs in stool samples, which might provide



**Fig. 2.** Immunofluorescence antibody test using fractionated polyclonal IgY antibodies in *A. ceylanicum* sections and other parasitic sample (*S. venezuelensis* section) to evaluate the possibility of cross reaction. *Ancylostoma ceylanicum* and *S. venezuelensis* sections were incubated with anti-*A. ceylanicum* IgY or with IgY from the control group (immunized with PBS). Column A (green): anti-chicken IgY FITC conjugated; B (red): Evans blue counterstained; C (merged): from column 1 to 2.

false-negative results, especially in low-intensity infections and following treatment (Nikolay *et al.*, 2014). Therefore, improved diagnostic tools are needed to facilitate the detection of active hookworm infections, considering that a lack of specificity, cross-reactivity and the inability to distinguish between past and current infections are limitations of many immunodiagnostic tests.

Given the feasibility of using congeneric species to diagnose patients with parasites that share epitopes, such as *Taenia solium* and *Taenia saginata* (Oliveira *et al.*, 2007) and *Strongyloides stercoralis* and *S. venezuelensis* (Feliciano et al., 2010; Faria *et al.*, 2019), we evaluated the use of anti-*A. ceylanicum* IgY to detect hookworm infections caused by *A. duodenale* or *N. americanus* in Brazil, as *A. ceylanicum* is closely related to these species and therefore shares various antigenic determinants. In addition, unlike the major hookworm species, *A. ceylanicum* to parasitize laboratory animals has made it a useful model for the study of human hookworm infection (Hu *et al.*, 2013).


Tests to detect antibodies often present cross-reactivity with other parasitic nematodes, especially in tropical countries where infection with STH is common. Furthermore, detection of IgG antibodies may be attributed to current infection, previous infection or exposure to the parasite and not always to an active infection (Werkman *et al.*, 2018). Therefore, the presence of antigen and antibody immunocomplexes may serve as markers of disease activity, while antibody detection may be indicative of a resolved disease with residual antibodies (Gonçalves *et al.*, 2016).

Polyclonal antibodies produced in chickens (IgY) offer several important advantages over antibodies produced in mammals. For example, production in chickens avoids interference in immunological assays caused by the complement system (Larsson *et al.*, 1992), chicken antibodies do not react with the rheumatoid factor, they have poor cross-reactivity to mammalian IgGs due to immunological differences (Amro *et al.*, 2018) and they react with more epitopes on a mammalian antigen due to evolutionary differences between mammals and birds, thereby amplifying the signal (Ohnishi *et al.*, 2000).

IgY has been successfully used in different immunoassay formats, mainly ELISA and lateral flow immunoassays, for the detection of veterinary drugs, herbicides, mycotoxins, drugs of abuse, proteins and bacteria (Zhang *et al.*, 2015) for both diagnostic (Ferreira Júnior *et al.*, 2012; Faria *et al.*, 2019) and immunotherapeutic (Thirumalai *et al.*, 2019) purposes.

In the present study, IgY production reached the cut-off point after the first booster immunization (day 21), as chickens transfer maternal antibodies to the egg yolk to protect the developing embryos through passive immunity. Thus, the lag time observed in the first 2 weeks may be attributed to the time taken for specific antibodies to be produced and transferred from the serum to the egg yolk (Esmailnejad *et al.*, 2019).

In this study, the detection of circulating ICs in serum was significant (90%) in group 1. ICs are important because they neutralize disease pathogenesis and are rapidly eliminated from the bloodstream by the innate immune system; therefore, they act as a biomarker of active infection. In addition, ICs are produced at increased frequency in parasitic diseases (Ohyama *et al.*, 2016). The absence of detection in some individuals from group 1 (10%) may be due to the IC size, as smaller ICs generally remain in the circulation or are rapidly cleared or dissociated, intermediate-sized ICs can deposit in tissue while larger ones are generally rapidly cleared by the phagocytic system (Rojko *et al.*, 2014).



The positive detection of ICs in group 2 (26.6%, n = 8) can be explained by the lack of sensitivity of faecal tests, as the tested individuals all live in areas where hookworms are prevalent, so they could possibly be coinfected with this parasite. Parasitological investigation of a single stool specimen detects about 30% of infections, and the sensitivity is increased to 50% if three faecal samples are tested (Sudré *et al.*, 2006).

To our best of our knowledge, this is the first report of an IC detection assay based on polyclonal IgY to detect active hookworm infection in human serum samples. The IgY obtained here was successfully tested and can be used as a direct reagent in serological assays. However, it must be continuously tested to demonstrate its potential use, as novel diagnostic methods are expected to improve epidemiological studies and control efforts for the prevention and treatment of hookworm infection.

**Financial support.** The study was supported by the Coordenação de Aperfeiçoamento de Pessoal de Nível Superior (CAPES – NO 3200601204P-8), Fundação de Amparo à Pesquisa do Estado de Minas Gerais (FAPEMIG – CBB-PPM-00396-1) and Conselho Nacional de Desenvolvimento Científico e Tecnológico (CNPq – 303843/2015-2 and CNPq–404816/2016-9).

#### Conflicts of interest. None.

**Ethical standards.** The evaluation of the serological diagnosis was performed according to the ethical guidelines of the Brazilian Health Ministry and was approved by the Comitê de Ética em Pesquisas com Seres Humanos (CEP) of the Universidade Federal de Uberlândia (CAAE: 48492315.8.00005152 and 04/2008). All procedures using experimental animals (hens, *Gallus gallus domesticus*) were conducted after being approved by the *Comitê de Ética na Utilização de Animais* from the *Universidade Federal de Uberlândia*, State of Minas Gerais, Brazil (CEUA/UFU No 097/2014).

#### References

Akita EM and Nakai S (1993) Comparison of four purification methods for the production of immunoglobulins from eggs laid by hens immunized **Fig. 3.** ELISA to detect ICs in serum samples from patients with hookworm active infection (group 1, n = 43), other parasitic infections (group 2, n = 30) and apparently healthy individuals (group 3, n = 30) using anti-hookworm IgY. A: scatter dot plots. Dotted line indicates cut off, samples were considered positive when OD >0.4, horizontal bars indicate median. B: ROC curve indicating AUC, sensitivity (Se), specificity (Sp) and positive LR (95% CI indicated).

with an enterotoxigenic *E. coli* strain. *Journal of Immunological Methods* **160**, 207–214.

- Amro WA, Al-Qaisi W and Al-Razem F (2018) Production and purification of IgY antibodies from chicken egg yolk. *Journal of Genetic Engineering and Biotechnology* 16, 99–103.
- Baermann G (1917) Eine Einfache Methode zur Auffindung von Ankylostomum (Nematoden) Larven in Erdproben. In Mededeel. Mith H. Geneesk. Batavia. Lab Weltevreden Feestbundel 57, 41–47.
- Elsemore DA, Geng J, Cote J, Hanna R, Lucio-Forster A and Bowman DD (2017) Enzyme-linked immunosorbent assays for coproantigen detection of *Ancylostoma caninum* and *Toxocara canis* in dogs and *Toxocara cati* in cats. *Journal of Veterinary Diagnostic Investigation* **29**, 645–653.
- Esmailnejad A, Abdi-Hachesoo B, Hosseini Nasab E and Shakoori M (2019) Production, purification, and evaluation of quail immunoglobulin Y against Salmonella typhimurium and Salmonella enteritidis. Molecular Immunology 107, 79–83.
- Faria LS, de Souza DLN, Ribeiro RP, de Sousa JEN, Borges IP, Ávila VMR, Ferreira-Júnior Á, Goulart LR and Costa-Cruz JM (2019) Highly specific and sensitive anti-Strongyloides venezuelensis Igy antibodies applied to the human strongyloidiasis immunodiagnosis. Parasitology International 72, 101933.
- Feliciano ND, Gonzaga HT, Gonçalves-Pires MdR, Gonçalves AL, Rodrigues RM, Ueta MT and Costa-Cruz JM (2010) Hydrophobic fractions from *Strongyloides venezuelensis* for use in the human immunodiagnosis of strongyloidiasis. *Diagnostic Microbiology and Infectious Diseases* 67, 153–161.
- Ferreira Júnior Á, Santiago FM, Silva MV, Ferreira FB, Macêdo Júnior AG, Mota CM, Faria MS, Silva Filho HH, Silva DA, Cunha-Júnior JP, Mineo JR and Mineo TW (2012) Production, characterization and applications for *Toxoplasma gondii*-specific polyclonal chicken egg yolk immunoglobulins. *PLoS ONE* 7, e40391.
- Friedman RD (1982) Comparison of four different silver-staining techniques for salivary protein detection in alkaline polyacrylamide gels. *Analytical Biochemistry* 126, 346–349.
- George S, Levecke B, Kattula D, Velusamy V, Roy S, Geldhof P, Sarkar R and Kang G (2016) Molecular identification of hookworm isolates in humans, dogs and soil in a tribal area in Tamil Nadu, India. *PLoS Neglected Tropical Diseases* **10**, e0004891.

- Gonçalves AL, Nunes DS, Gonçalves-Pires MR, Ueta MT and Costa-Cruz JM (2012) Use of larval, parasitic female and egg antigens from *Strongyloides venezuelensis* to detect parasite-specific IgG and immune complexes in immunodiagnosis of human strongyloidiasis. *Parasitology* 139, 956–961.
- Gonçalves AL, de Araújo KC, Carvalho EF, Ueta MT and Costa-Cruz JM (2016) Specific IgG and immune complex responses to parthenogenetic females and eggs of nematode *Strongyloides venezuelensis* for the diagnosis of immunosuppression in infected rats. *Journal of Helminthology* **90**, 342–346.
- Gonzaga HT, Vila-Verde C, Nunes DS, Ribeiro VS, Cunha-Júnior JP and Costa-Cruz JM (2013) Ion-exchange protocol to obtain antigenic fractions with potential for serodiagnosis of strongyloidiasis. *Parasitology* 140, 69–75.
- Hu Y, Ellis BL, Yiu YY, Miller MM, Urban JF, Shi LZ and Aroian RV (2013) An extensive comparison of the effect of anthelmintic classes on diverse nematodes. *PLoS ONE* 8, e70702.
- Laemmli UK (1970) Cleavage of structural proteins during the assembly of the head of bacteriophage T4. *Nature* 227, 680–685.
- Larsson A, Wejåker PE, Forsberg PO and Lindahl T (1992) Chicken antibodies: a tool to avoid interference by complement activation in ELISA. *Journal of Immunology Methods* 156, 79–83.
- Lei JH, Guan F, Xu H, Chen L, Su BT, Zhou Y, Wang T, Li YL and Liu WQ (2012) Application of an immunomagnetic bead ELISA based on IgY for detection of circulating antigen in urine of mice infected with *Schistosoma japonicum*. Veterinary Parasitology 187, 196–202.
- Leslie GA and Clem LW (1969) Phylogeny of immunoglobulin structure and function 3. Immunoglobulins of the chicken. *Journal of Experimental Medicine* 130, 1337–1352.
- Li X, Wang L, Zhen Y, Li S and Xu Y (2015) Chicken egg yolk antibodies (IgY) as non-antibiotic production enhancers for use in swine production: a review. *Journal of Animal Science and Biotechnology* **26**, 40.
- Liu LN, Jing FJ, Cui J, Fu GY and Wang ZQ (2013) Detection of circulating antigen in serum of mice infected with *Trichinella spiralis* by an IgY-IgM mAb sandwich ELISA. *Experimental Parasitology* 133, 150–155.
- Lopes CA, de Faria LS, de Sousa JEN, Borges IP, Ribeiro RP, Bueno LL, Rodrigues Ávila VM, Ferreira-Júnior Á and Costa-Cruz JM (2019) Anti-Ascaris suum immunoglobulin Y as a novel biotechnological tool for the diagnosis of human ascariasis. Journal of Helminthology 94, e71.
- Lowry OH, Rosebrough NJ, Farr AL and Randall RJ (1951) Protein measurement with the Folin phenol reagent. *The Journal of Biological Chemistry* 193, 265–275.
- Lutz A (1919) O Schistosomum mansoni e a schistosomatose segundo observaçõões feitas no Brasil. Memórias do Instituto Oswaldo Cruz 11, 121–125.
- Moraes RG (1948) Contribuição para o estudo do *Strongyloides Stercoralis* e da estrongiloidíase no Brasil. *Revista do Serviço Especial de Saúde Pública* 1, 507–624.
- Nie G, Wang T, Lu S, Liu W, Li Y and Lei J (2014) Detection of *Clonorchis* sinensis circulating antigen in sera from Chinese patients by immunomagnetic bead ELISA based on IgY. *PLoS ONE* 9, e113208.
- Nikolay B, Brooker SJ and Pullan RL (2014) Sensitivity of diagnostic tests for human soil-transmitted helminth infections: a meta-analysis in the absence of a true gold standard. *International Journal for Parasitology* 44, 765–774.
- Ohnishi T, Kakimoto K, Hashida S, Fujii M, Hirono S, Nishiyama K, Amita Y, Ishikawa E, Tsubouchi H and Daikuhara Y (2000) Development of highly sensitive enzyme-linked immunosorbent assays for hepatocyte

growth factor/scatter factor (HGF/SF): determination of HGF/SF in serum and urine from normal human subjects. *Journal of Immunology Methods* **244**, 163–173.

- Ohyama K, Huy NT, Yoshimi H, Kishikawa N, Nishizawa JE, Roca Y, Revolvo Guzmán RJ, Velarde FU, Kuroda N and Hirayama K (2016) Proteomic profile of circulating immune complexes in chronic Chagas disease. *Parasite Immunology* **38**, 609–617.
- Oliveira HB, Machado GA, Cabral DD and Costa-Cruz JM (2007) Application of *Taenia saginata* Metacestodes as an alternative antigen for the serological diagnosis of human neurocysticercosis. *Parasitology Research* **101**, 1007–1013.
- Prieto-Pérez L, Pérez-Tanoira R, Cabello-Úbeda A, Petkova-Saiz E and Górgolas-Hernandéz-Mora M (2016) Geohelminths. Enfermedades Infecciosas y Microbiología Clínica 34, 384–389.
- Pullan RL, Smith JL, Jasrasaria R and Brooker SJ (2014) Global numbers of infection and disease burden of soil transmitted helminth infections in 2010. Parasites and Vectors 7, 37.
- Ribeiro VdS, Araújo TG, Gonzaga HT, Nascimento R, Goulart LR and Costa-Cruz JM (2013) Development of specific scFv antibodies to detect neurocysticercosis antigens and potential applications in immunodiagnosis. *Immunology Letters* 156, 59–67.
- Rojko JL, Evans MG, Price SA, Han B, Waine G, DeWitte M, Haynes J, Freimark B, Martin P, Raymond JT, Evering W, Rebelatto MC, Schenck E and Horvath C (2014) Formation, clearance, deposition, pathogenicity, and identification of biopharmaceutical-related immune complexes: review and case studies. *Toxicologic Pathology* 42, 725–764.
- Schwarzkopf, C, Staak, C, Behn, I and Erhard M (2001) Immunisation. In Schade R, Behn I, Erhard M, Hlinak A and Staak C (eds), *Chickens Egg* Yolk Antibodies, Production and Applications-IgY-Technology. Berlin, Heidelberg: Springer-Verlag, pp. 25–64.
- Sudré AP, Macedo HW, Peralta RHS and Peralta JM (2006) Diagnóstico da estrongiloidíase humana: importância e técnicas. *Revista de Patologia Tropical* 35, 173–184.
- Thirumalai D, Visaga Ambi S, Vieira-Pires RS, Xiaoying Z, Sekaran S and Krishnan U (2019) Chicken egg yolk antibody (IgY) as diagnostics and therapeutics in parasitic infections – A review. *International Journal of Biological Macromolecules* 136, 755–763.
- Towbin H, Staehelin T and Gordon J (1979) Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. Proceedings of the National Academy of Sciences of the United States of America 76, 4350–4354.
- **Truscott JE, Turner HC, Farrell SH and Anderson RM** (2016) Soil-transmitted helminths: mathematical models of transmission, the impact of mass drug administration and transmission elimination criteria. *Advances in Parasitology* **94**, 133–198.
- Werkman M, Wright JE, Truscott JE, Easton AV, Oliveira RG, Toor J, Ower A, Ásbjörnsdóttir KH, Means AR, Farrell SH, Walson JL and Anderson RM (2018) Testing for soil-transmitted helminth transmission elimination: analysing the impact of the sensitivity of different diagnostic tools. PLoS Neglected Tropical Diseases 12, e0006114.
- Zhang H, Mi T, Khan OY, Sheng Y, Eremin SA, Beier RC, Zhang S, Shen J and Wang Z (2015) Fluorescence polarization immunoassay using IgY antibodies for detection of valnemulin in swine tissue. *Analytical and Bioanalytical Chemistry* 407, 7843–7848.