



Identification of Potential Genomic Alterations Using Pan-Cancer Cell-Free DNA Next-Generation Sequencing in Patients With Gastric Cancer

Boyeon Kim , M.D., Ph.D.^{1,2,*}, Yoonjung Kim , M.D., Ph.D.^{3,*}, Jae Yong Cho , M.D., Ph.D.⁴, and Kyung-A Lee , M.D., Ph.D.³

¹Division of Biotechnology, Invites BioCore Co. Ltd., Yongin, Korea; ²Genome Service Development, Invites Genomics Co. Ltd., Jeju, Korea; ³Department of Laboratory Medicine, Gangnam Severance Hospital, Yonsei University College of Medicine, Seoul, Korea; ⁴Division of Medical Oncology, Department of Internal Medicine, Gangnam Severance Hospital, Yonsei University College of Medicine, Seoul, Korea

Background: Molecular cancer profiling may lead to appropriate trials for molecularly targeted therapies. Cell-free DNA (cfDNA) is a promising diagnostic and/or prognostic biomarker in gastric cancer (GC). We characterized somatic genomic alterations in cfDNA of patients with GC.

Methods: Medical records and cfDNA data of 81 patients diagnosed as having GC were reviewed. Forty-nine and 32 patients were tested using the OncoPrint Pan-Cancer Cell-Free Assay on the Ion Torrent platform and AlphaLiquid 100 kit on the Illumina platform, respectively.

Results: Tier I or II alterations were detected in 64.2% (52/81) of patients. Biomarkers for potential targeted therapy were detected in 55.6% of patients (45/81), and clinical trials are underway. *ERBB2* amplification is actionable and was detected in 4.9% of patients (4/81). Among biomarkers showing potential for possible targeted therapy, *TP53* mutation (38.3%, 35 variants in 31 patients, 31/81) and *FGFR2* amplification (6.2%, 5/81) were detected the most.

Conclusions: Next-generation sequencing of cfDNA is a promising technique for the molecular profiling of GC. Evidence suggests that cfDNA analysis can provide accurate and reliable information on somatic genomic alterations in patients with GC, potentially replacing tissue biopsy as a diagnostic and prognostic tool. Through cfDNA analysis for molecular profiling, it may be possible to translate the molecular classification into therapeutic targets and predictive biomarkers, leading to personalized treatment options for patients with GC in the future.

Key Words: Biomarkers, Cell-free nucleic acids, Diagnostics, Genomics, High-throughput nucleotide sequencing, Mutation, Stomach neoplasm

Received: May 3, 2023

Revision received: June 27, 2023

Accepted: September 12, 2023

Published online: October 30, 2023

Corresponding author:

Kyung-A Lee, M.D., Ph.D.
Department of Laboratory Medicine,
Gangnam Severance Hospital, Yonsei
University College of Medicine,
211 Eonju-Ro, Gangnam-Gu, Seoul 06273,
Korea
E-mail: kal1119@yuhs.ac

Co-corresponding author:

Jae Yong Cho, M.D., Ph.D.
Division of Medical Oncology, Department
of Internal Medicine, Gangnam Severance
Hospital, Yonsei University College of
Medicine, 211 Eonju-Ro, Gangnam-Gu,
Seoul 06273, Korea
E-mail: chojy@yuhs.ac

*These authors contributed equally to this study as co-first authors.



© Korean Society for Laboratory Medicine

This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<https://creativecommons.org/licenses/by-nc/4.0>) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

INTRODUCTION

Gastric cancer (GC) is the fifth most common cancer and the

fourth most common cause of cancer-related deaths worldwide [1], highlighting the need for a better understanding of its molecular characteristics. Advances in molecular analysis have

enabled the identification of patient subsets with differing genomic alterations despite the same histologic diagnosis in GC [2]. In 2014, The Cancer Genome Atlas (TCGA) suggested four subtypes of GC and gastroesophageal junction cancers: Epstein–Barr virus, microsatellite instability (MSI), chromosomal instability, and genomically stable [3]. The molecular and clinical characteristics of these subtypes differ, emphasizing the importance of molecular profiling in patients with GC. Kim, *et al.* [4] classified the molecular subtypes of GC in Korean patients according to TCGA system and suggested that comprehensive analysis of genomic alterations may aid in the design of more appropriate clinical trials.

According to Surveillance, Epidemiology, and End Results Program (SEER) data from the National Cancer Institute spanning 2011–2020, 33.9% of patients with GC are diagnosed in stage IV [5]. Data from the National Cancer Information Center spanning 2016–2020 indicated that 11.0% of patients are diagnosed in a late stage [6]. The National Comprehensive Cancer Network recommends the use of a validated comprehensive next-generation sequencing (NGS)-based cell-free DNA (cfDNA) genomic profiling assay for patients with advanced or metastatic GC who may be unable to undergo a tissue biopsy for disease progression monitoring [7]. cfDNA analysis has gained attention for its capability of sufficiently identifying the molecular profile of cancer, facilitating the development of appropriate trials for molecularly targeted therapies [8]. The Food and Drug Administration (FDA) recently approved *PIK3CA* and *EGFR* mutation testing of cfDNA in breast cancer and non-small cell lung cancer, respectively, highlighting the clinical utility of cfDNA analysis [9, 10].

Genetic analysis in patients with GC has revealed several biomarkers for targeted therapy in GC. For example, the identification of *ERBB2* amplification led to the first targeted treatment approach for GC [11]. The VIKTORY (targeted agent eValuation In gastric cancer basket KOREa) trial in 2019 was the first and largest platform study in GC, showing the feasibility of tumor profiling and its clinical utility [12]. Recently, the fibroblast growth factor receptor (FGFR) pathway has been reported as a promising treatment target, and clinical trials on FGFR pathway inhibition are being implemented [13, 14]. As suggested by several studies [4, 12, 15], GC molecular profiling aids not only in patient disease monitoring and treatment selection but also in the identification of candidate biomarkers for targeted therapies.

We aimed to identify well-defined targets or biological pathways that could potentially be used for molecular targeted ther-

apy in patients with GC by analyzing cfDNA data obtained using the Ion Torrent S5 XL (Thermo Fisher Scientific, Waltham, MA, USA) and Illumina NextSeq-550 (Illumina, San Diego, CA, USA) platforms. In addition, we demonstrated the advantages of cfDNA analysis in patients with GC.

MATERIALS AND METHODS

Study population

We conducted cfDNA analysis of 149 patients with solid cancer between September 2019 and February 2022. Among these patients, 81 were diagnosed as having advanced or metastatic GC by oncologists. In this retrospective study, cfDNA data and the medical records of these patients were reviewed. Forty of the patients were treatment-naïve, whereas 41 had previously undergone chemotherapy but had experienced recurrence or progression of GC. For the cfDNA analysis, all patients provided written informed consent for sample collection and genetic analysis. The study was approved by the Institutional Review Board of Gangnam Severance Hospital, Seoul, Korea (IRB No. 3-2022-0136).

Nucleic acid extraction from plasma

Peripheral blood samples were collected into EDTA blood collection tubes (Vacutainer K2EDTA; BD Biosciences, Franklin Lakes, NJ, USA). Plasma was isolated by centrifugation at 2,000 × *g* for 10 mins followed by high-spin centrifugation at 16,000 × *g*, 4°C for 10 mins. Plasma aliquots were stored at –80°C [16]. cfDNA was extracted using the MagMAX Cell-Free Total Nucleic Acid Kit (Thermo Fisher Scientific) according to the manufacturer's instructions. The cfDNA concentration and size distribution were assessed using a Qubit 3.0 Fluorometer (Thermo Fisher Scientific) for samples tested using Ion Torrent and a 2200 TapeStation Instrument (Agilent Technologies, Santa Clara, CA, USA) with the Agilent High Sensitivity D1000 ScreenTape System or Genomic DNA ScreenTape System for samples tested using Illumina. For assays A and B, approximately 30 µg of fragmented DNA from each plasma sample (4 mL) was used.

Cell-free DNA NGS assays

Genomic alterations that may serve as cancer biomarkers were analyzed using two NGS panels and sequencing platforms. The OncoPrint Pan-Cancer Cell-Free Assay (Thermo Fisher Scientific) using Ion Torrent semi-conductor sequencing technology and the AlphaLiquid 100 kit (IMBDx, Seoul, Korea) using sequencing-by-synthesis technology were used. The respective target genes

are listed in Table 1. In this paper, the OncoPrint Pan-Cancer Cell-Free Assay is referred to as “assay A” and the AlphaLiquid 100 assay as “assay B”. For assay A, *BRCA1*, *BRCA2*, and *MYCN* were additionally tested using the Ion AmpliSeq HD Custom Panel (Thermo Fisher Scientific).

Assay A

A DNA library of 49 samples was prepared. Templating and sequencing were performed using the Ion 540 Kit on the Ion Chef and on Ion Torrent. Alignment to the hg19 human reference genome and variant calling were performed using the Torrent Suite software v.5.12 and Ion Reporter v.5.14 (Thermo Fisher Scientific), respectively. The Torrent Suite software was used to assess molecular coverage depth and read coverage depth at the target base.

Assay B

A DNA library of 32 samples was prepared. Hybrid capture-selected libraries were sequenced on an Illumina in the 150-bp

paired-end mode. All sequencing reads were demultiplexed and converted into unmapped bam files, which store extracted unique molecular identifier (UMI) sequences. The unmapped bam files were aligned to the reference genome (GRCh38) using the Burrows–Wheeler Aligner MEM algorithm (0.7.17-r1188) [17]. GeneFuse [18] and SVICT [19] were used for gene fusion detection. Sequencing data were analyzed using AlphaLiquid Analysis Browser IMBER v.1.2 (IMBDx). Assay B can detect MSI in cfDNA. MSI was analyzed according to the manufacturer’s protocol.

Data interpretation and reporting

Both assays allow analyzing single-nucleotide variants (SNVs), insertions/deletions (indels), copy number variations (CNVs), and fusions. According to the manufacturers’ specifications, the limits of detection (LODs) for SNVs and indels were set at an allele frequency (AF) of 0.1%. For LOD validation, we conducted tests using Seracare circulating tumor DNA (ctDNA) complete samples with AFs of 0.5% and 2.5%, and Seracare ctDNA Muta-

Table 1. Target genes evaluated by the two assays used for cfDNA analysis in this study

Mutation type	Genes	
	OncoPrint Pan-Cancer Cell-Free Assay	AlphaLiquid 100 kit
SNVs, small indels	<i>AKT1, ALK, APC, AR, ARAF, BRAF, BRCA1*, BRCA2*, CHEK2, CTNNB1, DDR2, EGFR, ERBB2, ERBB3, ESR1, FBXW7, FGFR1, FGFR2, FGFR3, FGFR4, FLT3, GNA11, GNAQ, GNAS, HRAS, IDH1, IDH2, KIT, KRAS, MAP2K1, MAP2K2, MET, MTOR, NRAS, NTRK1, NTRK3, PDGFRA, PIK3CA, PTEN, RAF1, RET, ROS1, SF3B1, SMAD4, SMO, TP53</i>	<i>ABL1, AKT1, AKT2, ALK, APC, AR, ARAF, ARID1A, ATM, BRAF, BRCA1, BRCA2, BTK, CBL, CCND1, CCND2, CCNE1, CD274, CDH1, CDK4, CDK6, CDKN2A, CEBPA, CSF1R, CTNNB1, DDR2, DPYD, EGFR, ERBB2, ERBB3, ESR1, FBXW7, FGFR1, FGFR2, FGFR3, FLT3, GATA3, GNA11, GNAQ, GNAS, HRAS, IDH1, IDH2, IGF1R, JAK2, JAK3, KDM6A, KDR, KEAP1, KIT, KRAS, MAP2K1, MAP2K2, MAPK1, MAPK3, MDM2, MET, MLH1, MPL, MSH2, MSH6, MTOR, MYC, MYCN, NF1, NF2, NFE2L2, NOTCH1, NPM1, NRAS, NTRK1, NTRK2, NTRK3, PDCD1LG2, PDGFRA, PDGFRB, PIK3CA, PIK3R1, PMS2, PPP2R1A, PTEN, PTPN11, RAF1, RB1, RET, RHEB, RHOA, RIT1, RNF43, ROS1, RUNX1, SETD2, SMAD4, SMO, STAG2, STK11, TCF7L2, TERT, TOP2A, TP53, TSC1, TSC2, U2AF1, UGT1A1, VHL</i>
CNVs	<i>CCND1, CCND2, CCND3, CDK4, CDK6, EGFR, ERBB2, FGFR1, FGFR2, FGFR3, MET, MYC, MYCN*</i>	<i>ABL1, AKT1, AKT2, ALK, APC, AR, ARAF, ARID1A, ATM, BRAF, BRCA1, BRCA2, BTK, CBL, CCND1, CCND2, CCNE1, CD274, CDH1, CDK4, CDK6, CDKN2A, CEBPA, CSF1R, CTNNB1, DDR2, DPYD, EGFR, ERBB2, ERBB3, ESR1, FBXW7, FGFR1, FGFR2, FGFR3, FLT3, GATA3, GNA11, GNAQ, GNAS, HRAS, IDH1, IDH2, IGF1R, JAK2, JAK3, KDM6A, KDR, KEAP1, KIT, KRAS, MAP2K1, MAP2K2, MAPK1, MAPK3, MDM2, MET, MLH1, MPL, MSH2, MSH6, MTOR, MYC, MYCN, NF1, NF2, NFE2L2, NOTCH1, NPM1, NRAS, NTRK1, NTRK2, NTRK3, PDCD1LG2, PDGFRA, PDGFRB, PIK3CA, PIK3R1, PMS2, PPP2R1A, PTEN, PTPN11, RAF1, RB1, RET, RHEB, RHOA, RIT1, RNF43, ROS1, RUNX1, SETD2, SMAD4, SMO, STAG2, STK11, TCF7L2, TERT, TOP2A, TP53, TSC1, TSC2, U2AF1, UGT1A1, VHL</i>
Fusions	<i>ALK, BRAF, ERG, ETV1, FGFR1, FGFR2, FGFR3, MET, NTRK1, NTRK3, RET, ROS1</i>	<i>ALK, BCR, BRAF, EGFR, FGFR2, FGFR3, NTRK1, NTRK2, RET, ROS1</i>
Exon skipping	<i>MET exon 14 skipping</i>	<i>MET exon 14 skipping</i>

*Genes additionally tested using Ion AmpliSeq HD Custom Panels.

Abbreviations: cfDNA, cell-free DNA; SNV, single-nucleotide variant; indel, insertion/deletion; CNV, copy number variation.

tion V2 samples with an AF of 0.5%. We were unable to confirm variants with an AF < 0.5%. Deveson, *et al.* [20] suggested that variants with an AF > 0.5% can be detected with high sensitivity, precision, and reproducibility. Considering these findings, candidate variants were considered valid only if they met two criteria: (1) a minimum coverage of >5 molecular counts (assay A) or >5 fragment counts (assay B), and (2) an AF of $\geq 0.5\%$ of total reads covering a given base in a single sample. All genomic variants were manually checked using Integrative Genomics Viewer.

For CNVs and fusions, we adapted the manufacturers' specifications. For CNV reporting, assay A instructions recommend the sample-level QC criteria of the median absolute pairwise difference to be < 0.4. The CNV ratio call thresholds had been empirically derived by the manufacturer using plasma samples from healthy donors with a normal CNV status [21]. To pass the bioinformatic QC threshold for CNVs, the CNV ratio for a copy number gain must be > 1.15, and called copy number gains were visually reviewed from a copy number plot generated by Ion Reporter [22]. Assay A does not report copy number loss. For CNV reporting in assay B, at least ≥ 2.6 copies were considered. For fusion reporting, ≥ 2 reads with mapping quality ≥ 10 should be satisfied.

For interpretation and reporting, tiers I–III of both fusions and CNVs were reported according to the guidelines of the Association for Molecular Pathology, American Society of Clinical Oncology, and College of American Pathologists [23]. To exclude germline variants, we conducted Sanger sequencing on patients' peripheral blood mononuclear cells (PBMCs) for variants classified as tier I or tier II with AFs of 40%–60%.

Statistical analysis

Statistical analyses were conducted and graphs generated using GraphPad Prism software v.9.3 for Windows (GraphPad Software, San Diego, CA, USA) and Microsoft Excel 2013 (Microsoft, Redmond, WA, USA). Data are presented as medians and interquartile ranges.

RESULTS

Quality metrics

Sample-level QCs of assays A and B are summarized in Supplemental Data Tables S1 and S2, respectively. Sample DNA concentrations are presented in Supplemental Data Fig. S1. The median DNA concentrations measured using the Qubit 3.0 Fluorometer and 2200 TapeStation were 3.2 ng/ μ L and 2.8 ng/ μ L, respectively. The median cfDNA inputs in assays A and B were

32.0 ng and 30.9 ng, respectively.

For assay A, the median percentage on-target reads was 97%, median uniformity of base coverage was 99%, median overall mapped reads was 19,995,625, median read coverage was 60,380, and median molecular coverage was 6,201. For assay B, the median mapped reads in the high- and mid-output flow cells were 84,960,287 and 128,087,414, respectively. The median on-target read ratios were 64.5 for the high-output flow cell and 73.5 for the mid-output flow cell. The median on-target mean depths in the high- and mid-output flow cells were 13,021 and 20,875, respectively. Fragment uniformity was 97% for both flow cells.

Genomic alterations in GC

Among the 81 patients with GC, 48 were male and 33 were female, and the median age was 58 yrs (range, 25–82 yrs). At the time of diagnosis, three patients were in TNM stage IIb, one in IIIa, three in IIIb, two in IIIc, and 72 in stage IV. All patients were diagnosed as having advanced or metastatic GC by an oncologist. Patient characteristics and cfDNA analysis results are provided in Supplemental Data Table S3. Among the 81 patients, 40 patients were treatment-naïve. Tier I or II alterations were detected in 80% (32/40) of treatment-naïve patients and in 48.8% (20/41) of patients who were not treatment-naïve.

Among 32 samples tested with assay B, 30 samples showed microsatellite stability (MSS), one MSI-high, and one MSI-low. Six of the 32 patients were also tested for MSI in tissue; the results were all MSS, which agreed with the results for cfDNA. Tier I or II alterations were detected in 52 of 81 patients (64.2%). Assays A and B detected tier I or II variants in 28/49 patients (57.1%) and 24/32 patients (75.0%), respectively. Four CNVs were reported as tier I. Seventy-four SNVs, 21 indels, 23 CNVs, and one fusion were reported as tier II; 163 SNVs, 12 indels, and seven CNVs were reported as tier III. Patient No. 29 reportedly had an *EGFR-SEPT14* fusion, which is considered a noteworthy fusion in GC (Supplemental Data Table S3) [24]. The tier I and II mutations detected in the 81 patients are shown in Fig. 1.

The most commonly detected mutations were *TP53* mutations, which accounted for 38.3% (31/81 patients). Other frequently detected mutations included *APC* (13.6%, 11/81 patients), *KRAS* (9.9%, 8/81 patients), *ARID1A* (9.4%, 3/32 patients), *FBXW7* (4.9%, 4/81 patients), and *BRAF* (3.7%, 3/81 patients) mutations. CNVs were most frequently detected in *FGFR2*, accounting for 6.2% (5/81 patients), followed by *ERBB2* (4.9%, 4/81 patients), *CCNE1* (9.4%, 3/32 patients), *KRAS* (9.4%, 3/32 patients), *EGFR* (3.7%, 3/81 patients), and *MET*

(3.7%, 3/81 patients) (Table 2).

Molecular characterization for targeted therapy in clinical trials

The patients were categorized based on specific biomarkers for targeted therapy and markers undergoing clinical trials. Eleven biomarkers were selected from Nakamura, *et al.* [15] and Lee, *et al.* [12] (Table 3). Among these, *ERBB2* amplification is actionable and can be treated with a combination of trastuzumab and chemotherapy, as approved by the FDA for patients with GC [7, 15]. In our cfDNA data, *ERBB2* amplification was detected in four of 81 patients (4.9%). One of 49 patients (2.0%) showed *ERBB2* amplification in assay A and three of 32 patients (9.4%) in assay B.

The remaining 10 biomarkers are not currently actionable but show potential for possible targeted therapy, with clinical trials underway. Details on these 10 biomarkers in phase II or III clinical trials that are “active, not recruiting” or “recruiting” are listed in Table 3. Among the 52 patients with significant tier I or II variants, 45 patients (45/52, 86.5%; 45/81, 55.6%) had noteworthy biomarkers that could be detected using liquid biopsy.

DISCUSSION

We report on our experience with using two different NGS platforms for cfDNA analysis and demonstrated the effectiveness of liquid biopsy for molecular profiling of patients with GC. Using assays A and B, tier I or II variants were identified in 28/49 (57.1%) and 24/32 (75.0%) patients, respectively. Overall, tier I or II variants were detected in 64.2% (52/81) of patients with GC.

To increase the detection of true-positive mutations from highly fragmented and low quantities of circulating tumor DNA, UMIs were applied in both assays A and B. In assay B, the UniqSeq protocol and a series of IMBdx in-house filtering steps were used to eliminate false calls, and a machine learning model was used to distinguish true variants from false calls [25]. Both assays were able to detect significant variants in more than half of the patients, with assay B being developed later and containing more target genes than assay A.

Based on our routine laboratory practice, both NGS assays have their pros and cons. Assay A is developed for Ion Torrent, a PCR-based enrichment method. It has a two-day sample-to-report workflow [26], and cfDNA concentrations as low as 5 ng may be sufficient for evaluation [27]. However, because of the characteristics of the amplicon-based method, it may be difficult

Table 2. Commonly detected mutations in our study

Gene	Frequency	N by mutation type
<i>TP53</i>	38.3% (31/81)	31 SNVs and 4 indels
<i>APC</i>	13.6% (11/81)	5 SNVs and 6 indels
<i>KRAS</i>	9.9% (8/81)	8 SNVs
<i>ARID1A</i>	9.4% (3/32)	1 SNV and 3 indels
<i>FBXW7</i>	4.9% (4/81)	4 SNVs
<i>BRAF</i>	3.7% (3/81)	4 SNVs
<i>FGFR2*</i>	6.2% (5/81)	5 amplifications
<i>ERBB2*</i>	4.9% (4/81)	4 amplifications
<i>CCNE1*</i>	9.4% (3/32)	3 amplifications
<i>KRAS*</i>	9.4% (3/32)	3 amplifications
<i>EGFR*</i>	3.7% (3/81)	3 amplifications
<i>MET*</i>	3.7% (3/81)	3 amplifications

*CNVs were detected in these genes.

Abbreviations: SNV, single-nucleotide variant; CNV, copy number variation.

to detect unexpected fusions, and the fact that MSI information is not provided can be a disadvantage. Assay B is developed for Illumina, a hybrid capture-based method. This method takes longer than the PCR-based enrichment method but is effective for targeting large genomic regions [28]. Furthermore, assay B provides MSI information and can detect unexpected fusions such as *EGFR-SEPT14* [24].

We selected actionable potential biomarkers from Nakamura, *et al.* [15] and Lee, *et al.* [12], who reported biomarkers currently in clinical trials for solid cancers, including GC. These biomarkers are not only associated with major signaling pathways in various cancer types but also hold therapeutic relevance. They are associated with the phosphoinositide 3-kinase/Akt [29] and epidermal growth factor receptor (EGFR) [30] signaling pathways as therapeutic targets in human cancers. As targets for oncogenic mutations, biomarkers associated with the RAS/RAF/MEK/ERK (MAPK) signaling pathway were included [31]. Preclinical studies have highlighted the potential therapeutic value of targeting cyclin E1 (*CCNE1*), AT-rich interactive domain-containing protein 1A (*ARID1A*), and *TP53* [32-35].

Both our data and GENIE v11.0-public data in the cBioPortal database (<http://cbioportal.org>) indicate that *TP53* mutations are the most frequent mutations in patients with GC. We detected 11 tier II and five tier III variants in *APC*, which is more frequently mutated in intestinal-type GC [36]. Of the 14 patients with *APC* mutations, six patients had no information regarding the histologic type according to the Lauren classification, one patient had a diffuse type, and the remaining seven patients

Table 3. Patient subgroup according to established and potential biomarkers of targeted therapies in gastric cancer

Biomarker*	Established/in clinical trials [†]	N	Specifications
<i>ERBB2</i> amplification	First line: trastuzumab+chemotherapy Third line: trastuzumab, deruxtecan NCT05190445, NCT05152147, NCT02465060, NCT04143711	4 (4/81, 4.9%)	Assay A: 1/49, 2.0% Assay B: 3/32, 9.4%
<i>FGFR2</i> amplification	NCT05019794, NCT04189445	5 (5/81, 6.2%)	Assay A: 1/49, 2.0% Assay B: 4/32, 12.5%
<i>FGFR1</i> amplification	NCT05019794, NCT04189445	1 (1/81, 1.2%)	Assay A: 0/49, 0.0% Assay B: 1/32, 3.1%
<i>EGFR</i> amplification	NCT04077255, NCT04739202	3 (3/81, 3.7%)	Assay A: 2/49, 4.1% Assay B: 1/32, 3.1%
<i>CCNE1</i> amplification [†]	NCT05252416	3 (3/32, 9.4%)	-
<i>RAS</i> mutation or amplification [†]	NCT02465060	12 (9/81, 11.1%; 3/32, 9.4%)	Assay A: 5/49, 10.2% Assay B: 4/32, 12.5% 9 mutations (8 KRAS, 1 NRAS) 3 KRAS amplifications
<i>TP53</i> mutation	NCT03641313	31 (31/81, 38.3%)	Assay A: 17/49, 34.7% Assay B: 14/32, 43.8% 35 mutations in 31 patients
<i>PIK3CA</i> mutation or amplification [†]	NCT04526470, NCT04739202, NCT02465060	2 (2/81, 2.5%)	Assay A: 1/49, 2.0% Assay B: 1/32, 3.1% 2 mutations No amplification
<i>ARID1A</i> mutation [†]	NCT05379972	3 (3/32, 9.4%)	4 mutations from 3 patients
<i>MET</i> amplification	NCT04923932, NCT05620628, NCT05439993, NCT03993873, NCT02465060	3 (3/81, 3.7%)	Assay A: 3/49, 6.1% Assay B: 0/32, 0.0%
<i>TSC2</i> null [†]	NCT02465060	1 (1/32, 3.1%)	p.Phe471LeufsTer14 from #69 patient

*Biomarkers first described in Nakamura, *et al.* [15] and Lee, *et al.* [12].

[†]*CCNE1* amplification, *KRAS* amplification, *PIK3CA* amplification, *ARID1A* mutation, and *TSC2* null were only evaluated in the 32 patients tested using assay B (AlphaLiquid 100 kit).

[‡]Phase II or III clinical trials with “active, not recruiting” or “recruiting” status are listed.

Abbreviation: NCT, national clinical trial.

had an intestinal type of GC. This finding may be attributed to the high seroprevalence of *Helicobacter pylori* in Korea, which is closely related to the incidence of GC in Asia [37]. According to Rahman, *et al.* [38], the incidence of *H. pylori* infection is higher in intestinal GC than in diffuse GC.

FGFR2 amplification was the most frequently detected CNV in our study, followed by *ERBB2* amplification. *FGFR2* amplification can be attributed to specific characteristics of our patient cohort. All of our patients were diagnosed as having advanced or metastatic GC. According to Kim, *et al.* [39], patients harboring *FGFR2* amplification have significantly worse survival. All five patients with *FGFR2* amplification in our study were diagnosed as having metastatic GC, emphasizing the association between *FGFR2* amplification and advanced disease. Among them, pa-

tient No. 9 exhibited the highest *FGFR2* copy number (57.2) and had a history of recurrent GC. These findings suggest a potential association between *FGFR2* amplification and poor survival outcomes in patients with GC.

This study had several limitations. First, we only used cfDNA data from patients with advanced or metastatic GC, and data from patients with early-stage GC were not included. Moreover, not all cfDNA data could be directly compared with tissue biopsy data. As obtaining tissue biopsies from patients with advanced or metastatic GC is challenging, not all patients could undergo tissue biopsy. Nevertheless, although limited by the retrospective nature of this study, we included data from patients who underwent tissue biopsy in Supplemental Data Table S3 to provide additional insights. Second, we detected significant mutations in

55.6% of patients; however, the data were obtained using two different NGS assays with different gene panels. Direct comparison of the same patients using both assays was not feasible because our laboratory initially used assay A and later transitioned to assay B. Nonetheless, it is worthwhile to share our experience with using two different commercial assays for cfDNA testing considering that each assay uses a different NGS platform. Furthermore, MSI detection is limited to assay B; therefore, MSI status could only be evaluated in the 32 samples tested with assay B. Finally, we could not distinguish between somatic variants and clonal hematopoiesis of indeterminate potential based on our data; however, tier I or II variants with an AF of 40%–60% were subjected to Sanger sequencing to assess whether they were germline or somatic in origin. As shown in Supplemental Data Table S4 and Supplemental Data Fig. S2, we conducted Sanger sequencing on four variants using patients' PBMCs, which confirmed that they were somatic.

Despite these limitations, our study highlighted the advantages of liquid biopsy. Obtaining tissue samples, especially from patients with advanced or metastatic GC, can be challenging. Even when tissue is obtained, there is a risk of missing certain mutations because of tissue heterogeneity. In our study, notable biomarkers were found in 45 of 81 patients (55.6%), demonstrating the effectiveness of cfDNA analysis. Compared with other genome profiling studies in Korean patients [12], our cfDNA analysis demonstrated the advantage of liquid biopsy. While the VIKTORY trial relied on results from 715 tissue samples, we solely used cfDNA extracted from plasma (Supplemental Data Table S5). The prevalence of *TP53* mutation, *RAS* or *PIK3CA* mutation/amplification, and *MET*, *FGFR2*, *EGFR*, *FGFR1*, or *CCNE1* amplifications in our study was comparable to that reported by Lee, *et al.* [12]. Moreover, among the 42 patients who underwent tissue biopsy in our study, 16.7% (7/42) exhibited HER2-positivity (Supplemental Data Table S3). As ctDNA exists at a very low level in cfDNA, there may be some discrepancies between tissue and cfDNA results; however, patient Nos. 8 and 73 showed HER2-positivity in both their tissues and cfDNA samples (Supplemental Data Table S3). Patient Nos. 68 and 76 were unable to undergo tissue biopsy, but their cfDNA analysis revealed *ERBB2* amplification (Supplemental Data Table S3), demonstrating the benefits of cfDNA analysis in such cases.

Our study suggests that the genetic profiling of GC can be useful for finding not only novel molecular pathways but also biomarkers that may be targeted by existing therapies effective in other cancer types. We detected an *EGFR-SEPT14* fusion that has been reported in glioblastoma, non-small cell lung cancer,

and colorectal cancer [24, 40]. This fusion can be a novel target for patients with GC that may be treatable with EGFR tyrosine kinase inhibitors, which is effective in other cancer types [24].

Our study provided insights into the strengths and weaknesses of two cfDNA assays using different NGS platforms for use in patients with GC. Our research can help clinicians and researchers make better-informed decisions about which method may be best suited for their particular needs and goals. Furthermore, our research demonstrated the effectiveness of cfDNA analysis in providing oncologists with tumor genetic information to better understand the disease progression and select suitable treatment options for patients with GC.

In conclusion, our findings suggest that cfDNA analysis holds great promise as an alternative for tissue biopsy as a diagnostic tool. Molecular profiles obtained through cfDNA analysis can be translated into therapeutic targets and predictive biomarkers, which may allow for personalized treatment in the future. Our study highlights the immense potential of cfDNA analysis in improving the management of patients with GC.

SUPPLEMENTARY MATERIALS

Supplementary materials can be found via <https://doi.org/10.3343/alm.2023.0187>

ACKNOWLEDGEMENTS

None.

AUTHOR CONTRIBUTIONS

Kim B and Kim Y analyzed the data and wrote the manuscript. Cho JY acquired funding and administrated the project. Lee KA managed the study and supervised and discussed the results. All authors reviewed and approved the final version of the manuscript.

CONFLICTS OF INTEREST

None declared.

RESEARCH FUNDING

This study was supported by a faculty research grant from Yonsei University College of Medicine (6-2019-0077) and the Ministry of Health and Welfare, Korea (HI14C1277).

REFERENCES

1. Sung H, Ferlay J, Siegel RL, Laversanne M, Soerjomataram I, Jemal A, et al. Global Cancer Statistics 2020: GLOBOCAN estimates of incidence and mortality worldwide for 36 cancers in 185 countries. *CA Cancer J Clin* 2021;71:209-49.
2. Cristescu R, Lee J, Nebozhyn M, Kim KM, Ting JC, Wong SS, et al. Molecular analysis of gastric cancer identifies subtypes associated with distinct clinical outcomes. *Nat Med* 2015;21:449-56.
3. Cancer Genome Atlas Research Network. Comprehensive molecular characterization of gastric adenocarcinoma. *Nature* 2014;513:202-9.
4. Kim Y, Cho MY, Kim J, Kim SN, Oh SC, Lee KA. Profiling cancer-associated genetic alterations and molecular classification of cancer in Korean gastric cancer patients. *Oncotarget* 2017;8:69888-905.
5. National Cancer Institute. SEER*Explorer: an interactive website for SEER cancer statistics [Internet]. Surveillance Research Program. <https://seer.cancer.gov/statistics-network/explorer/> (Updated on Jun 2023).
6. National Cancer Information Center. 5-year cancer relative survival rates by SEER summary stage. <https://www.cancer.go.kr/lay1/S1T648C652/contents.do> (Updated on Jan 2023).
7. National Comprehensive Cancer Network Guidelines for Gastric Cancer (Version 1.2023). https://www.nccn.org/professionals/physician_gls/pdf/gastric.pdf. (Updated on Mar 2023).
8. De Mattos-Arruda L and Siravegna G. How to use liquid biopsies to treat patients with cancer. *ESMO Open* 2021;6:100060.
9. Martínez-Saéz O, Chic N, Pascual T, Adamo B, Vidal M, González-Farré B, et al. Frequency and spectrum of *PIK3CA* somatic mutations in breast cancer. *Breast Cancer Res* 2020;22:45.
10. Kwapisz D. The first liquid biopsy test approved. Is it a new era of mutation testing for non-small cell lung cancer? *Ann Transl Med* 2017;5:46.
11. Grávalos C, Gómez-Martín C, Rivera F, Alés I, Queralt B, Márquez A, et al. Phase II study of trastuzumab and cisplatin as first-line therapy in patients with HER2-positive advanced gastric or gastroesophageal junction cancer. *Clin Transl Oncol* 2011;13:179-84.
12. Lee J, Kim ST, Kim K, Lee H, Kozarewa I, Mortimer PGS, et al. Tumor genomic profiling guides patients with metastatic gastric cancer to targeted treatment: the VIKTORY umbrella trial. *Cancer Discov* 2019;9:1388-405.
13. Hur JY, Chao J, Kim K, Kim ST, Kim KM, Klempner SJ, et al. High-level FGFR2 amplification is associated with poor prognosis and Lower response to chemotherapy in gastric cancers. *Pathol Res Pract* 2020;216:152878.
14. Merz V, Zecchetto C, Simionato F, Cavaliere A, Casalino S, Pavarana M, et al. A phase II trial of the FGFR inhibitor pemigatinib in patients with metastatic esophageal-gastric junction/gastric cancer trastuzumab resistant: the FiGhTeR trial. *Ther Adv Med Oncol* 2020;12:1758835920937889.
15. Nakamura Y, Kawazoe A, Lordick F, Janjigian YY, Shitara K. Biomarker-targeted therapies for advanced-stage gastric and gastro-oesophageal junction cancers: an emerging paradigm. *Nat Rev Clin Oncol* 2021;18:473-87.
16. Shin S, Woo HI, Kim JW, Kim Y, Lee KA. Clinical practice guidelines for pre-analytical procedures of plasma epidermal growth factor receptor variant testing. *Ann Lab Med* 2022;42:141-9.
17. Li H and Durbin R. Fast and accurate short read alignment with Burrows-Wheeler transform. *Bioinformatics* 2009;25:1754-60.
18. Chen S, Liu M, Huang T, Liao W, Xu M, Gu J. GeneFuse: detection and visualization of target gene fusions from DNA sequencing data. *Int J Biol Sci* 2018;14:843-8.
19. Gawroński AR, Lin YY, McConeghy B, LeBihan S, Asghari H, Koçkan C, et al. Structural variation and fusion detection using targeted sequencing data from circulating cell free DNA. *Nucleic Acids Res* 2019;47:e38.
20. Deveson IW, Gong B, Lai K, LoCoco JS, Richmond TA, Schageman J, et al. Evaluating the analytical validity of circulating tumor DNA sequencing assays for precision oncology. *Nat Biotechnol* 2021;39:1115-28.
21. Analytical validation of the OncoPrint Pan-Cancer Cell-Free Assay. White Paper. ThermoFisher Scientific. COL33100 0619, 2019. <https://assets.thermofisher.com/TFS-Assets/CSD/Reference-Materials/oncomine-pan-cancer-assay-white-paper.pdf>.
22. So MK, Park JH, Kim JW, Jang JH. Analytical validation of a pan-cancer panel for cell-free assay for the detection of *EGFR* mutations. *Diagnostics (Basel)* 2021;11:1022.
23. Li MM, Datto M, Duncavage EJ, Kulkarni S, Lindeman NI, Roy S, et al. Standards and guidelines for the interpretation and reporting of sequence variants in cancer: a joint consensus recommendation of the Association for Molecular Pathology, American Society of Clinical Oncology, and College of American Pathologists. *J Mol Diagn* 2017;19:4-23.
24. Kim B, Kim Y, Park I, Cho JY, Lee KA. Detection of *EGFR-SEPT14* fusion in cell-free DNA of a patient with advanced gastric cancer: a case report. *World J Clin Cases* 2021;9:2884-9.
25. Kim S, Lim Y, Kang JK, Kim HP, Roh H, Kim SY, et al. Dynamic changes in longitudinal circulating tumour DNA profile during metastatic colorectal cancer treatment. *Br J Cancer* 2022;127:898-907.
26. Li Y, Lea K, Kshatriya P, Cao R, Gu J, Schageman J, et al. PO-086 an efficient Ion Torrent™ next generation sequencing workflow for liquid biopsy research to assess cell-free total nucleic acid. *ESMO Open* 2018;3:A260.
27. Shah M, Takayasu T, Zorofchian Moghadamtousi S, Arevalo O, Chen M, Lan C, et al. Evaluation of the OncoPrint Pan-Cancer Cell-Free Assay for analyzing circulating tumor DNA in the cerebrospinal fluid in patients with central nervous system malignancies. *J Mol Diagn* 2021;23:171-80.
28. Samorodnitsky E, Jewell BM, Hagopian R, Miya J, Wing MR, Lyon E, et al. Evaluation of hybridization capture versus amplicon-based methods for whole-exome sequencing. *Hum Mutat* 2015;36:903-14.
29. He Y, Sun MM, Zhang GG, Yang J, Chen KS, Xu WW, et al. Targeting PI3K/Akt signal transduction for cancer therapy. *Signal Transduct Target Ther* 2021;6:425.
30. Levantini E, Maroni G, Del Re M, Tenen DG. EGFR signaling pathway as therapeutic target in human cancers. *Semin Cancer Biol* 2022;85:253-75.
31. Degirmenci U, Wang M, Hu J. Targeting aberrant RAS/RAF/MEK/ERK signaling for cancer therapy. *Cells* 2020;9.
32. Min A, Kim JE, Kim YJ, Lim JM, Kim S, Kim JW, et al. Cyclin E overexpression confers resistance to the CDK4/6 specific inhibitor palbociclib in gastric cancer cells. *Cancer Lett* 2018;430:123-32.
33. Pang W, Li Y, Guo W, Shen H. Cyclin E: a potential treatment target to reverse cancer chemoresistance by regulating the cell cycle. *Am J Transl Res* 2020;12:5170-87.
34. Gu Y, Zhang P, Wang J, Lin C, Liu H, Li H, et al. Somatic ARID1A mutation stratifies patients with gastric cancer to PD-1 blockade and adjuvant chemotherapy. *Cancer Immunol Immunother* 2023;72:1199-208.
35. Qadir J, Majid S, Khan MS, Rashid F, Wani MD, Bhat SA. Implication of ARID1A undercurrents and PDL1, TP53 overexpression in advanced gastric cancer. *Pathol Oncol Res* 2021;27:1609826.
36. Fang DC, Luo YH, Yang SM, Li XA, Ling XL, Fang L. Mutation analysis of APC gene in gastric cancer with microsatellite instability. *World J Gastroenterol* 2002;8:787-91.

37. Lee JH, Choi KD, Jung HY, Baik GH, Park JK, Kim SS, et al. Seroprevalence of *Helicobacter pylori* in Korea: a multicenter, nationwide study conducted in 2015 and 2016. *Helicobacter* 2018;23:e12463.
38. Rahman R, Asombang AW, Ibdah JA. Characteristics of gastric cancer in Asia. *World J Gastroenterol* 2014;20:4483-90.
39. Kim HS, Kim JH, Jang HJ. Pathologic and prognostic impacts of FGFR2 amplification in gastric cancer: a meta-analysis and systemic review. *J Cancer* 2019;10:2560-7.
40. Zheng S, Li H, Feng J, Jiang C, Lin Y, Xie Y, et al. Complete remission in leptomeningeal metastasis of NSCLC with rare *EGFR-SEPT14* fusion treated with osimertinib combined with intrathecal chemotherapy with pemetrexed. *Anticancer Drugs* 2022;33:e795-8.