


Diagnostic performance of the Cobas MTB and Logix Smart MTB for diagnosing pulmonary and extrapulmonary tuberculosis: a cross-sectional study of diagnostic tests

Miguel Hueda-Zavaleta , Juan Carlos Gomez de la Torre , Diana Minchón-Vizconde ,
Claudia Barletta-Carrillo , Cesar Copaja-Corzo , Gustavo Tapia-Sequeiros ,
Cinthya Flores , Cristian Piscoche , Cecilia Miranda , Ada Mendoza  and
Vicente A. Benites-Zapata 

Abstract

Background: Molecular tests have contributed to reducing the mortality rate through early and accurate diagnosis of tuberculosis (TB). This is due to their low processing complexity and diagnostic accuracy superior to conventional methods.

Objective: To evaluate the diagnostic performance of Cobas MTB and Logix Smart MTB compared to Xpert MTB/RIF Ultra for pulmonary tuberculosis (PTB) and extrapulmonary tuberculosis (EPTB).

Design: A cross-sectional study of diagnostic tests was carried out in a clinical laboratory in Lima, Peru.

Methods: All pulmonary and extrapulmonary samples from patients with presumptive TB who had been subjected to smear microscopy, Xpert MTB/RIF Ultra, Cobas MTB, Logix Smart MTB, and mycobacteria culture were included to determine their diagnostic performance.

Results: A total of 175 samples were included, 102 (58.3%) of pulmonary origin and 73 (41.7%) of extrapulmonary origin. Among the total samples, 19 (10.9%) had positive cultures (all were pulmonary samples), 48 (27.4%) had positive Xpert MTB/RIF Ultra results, 45 (25.7%) had positive Cobas MTB results, and 36 (20.6%) had positive Logix Smart MTB results. The agreement between Cobas MTB and Logix Smart with the Xpert MTB/RIF Ultra was 97.1% and 93.8%, respectively. Compared to Xpert MTB/RIF Ultra, the area under the curve/receiver operating characteristic and sensitivity of the Cobas MTB and Logix Smart methods were 0.95 and 91.7%, and 0.90 and 81.0%, respectively.

Conclusion: Cobas MTB and Logix Smart presented adequate performance for diagnosing pulmonary and extrapulmonary tuberculosis comparable to the Xpert MTB/RIF Ultra.

Keywords: extrapulmonary, molecular diagnostic techniques, *Mycobacterium tuberculosis*, pulmonary, tuberculosis

Received: 9 May 2024; revised manuscript accepted: 18 November 2024.

Ther Adv Infect Dis

2024, Vol. 11: 1–13

DOI: 10.1177/
20499361241304516

© The Author(s), 2024.
Article reuse guidelines:
sagepub.com/journals-
permissions

Correspondence to:
Vicente A. Benites-Zapata
Unidad de Investigación
para la Generación y
Síntesis de Evidencias en
Salud, Universidad San
Ignacio de Loyola, Lima,
Campus 2, avenida La
Fontana 750, La Molina,
Lima, Peru
vbenites@usil.edu.pe

Miguel Hueda-Zavaleta
Diana Minchón-Vizconde
Gustavo Tapia-Sequeiros
Diagnóstico, tratamiento
e investigación de
enfermedades infecciosas
y tropicales, Universidad
Privada de Tacna, Tacna,
Peru

**Juan Carlos Gomez de la
Torre**
Laboratorio clínico Roe,
Lima, Peru

Facultad de Medicina
Humana, Universidad de
Piura, Lima, Peru

Claudia Barletta-Carrillo
Cinthya Flores
Cristian Piscoche
Cecilia Miranda
Ada Mendoza
Sequence Reference Lab,
Lima, Peru

Cesar Copaja-Corzo
Unidad de Investigación
para la Generación y
Síntesis de Evidencias en
Salud, Universidad San
Ignacio de Loyola, Lima,
Peru

Servicio de Infectología,
Hospital Nacional Edgardo
Rebagliati Martins,
EsSalud. Lima, Peru

Introduction

At the level of the Americas, Peru is the second country with the highest number of estimated tuberculosis (TB) cases.¹ The use of molecular tests to diagnose TB has contributed to a 3.4% reduction in the global mortality rate per year and a 1.9% reduction in the incidence of new cases per year.² TB control requires early and timely diagnosis. Currently, the World Health Organization (WHO) recommends molecular tests such as Xpert MTB/RIF Ultra and Truenat MTB.³ This is due to its low processing complexity and sensitivity and specificity superior to smear microscopy and culture for the rapid diagnosis of TB.⁴

Molecular tests have good diagnostic performance for the detection of TB. Xpert MTB/RIF Ultra has the highest sensitivity (91%) for pulmonary tuberculosis (PTB).⁵ However, its diagnostic performance is variable in extrapulmonary tuberculosis (EPTB), with a sensitivity of 50%–95.8% for lymph node samples,^{6,7} 71.4%–90.9% for cerebrospinal fluid,^{8,9} and 47.6%–84.2% for pleural fluid.^{10–12}

It is necessary to identify new diagnostic methods that allow infections to be identified earlier than conventional methods, mainly in groups that are difficult to diagnose (such as patients with EPTB and HIV coinfections), enabling early diagnosis of resistance to anti-TB drugs and with lower costs, for universal implementation.

It is unknown whether there are differences in the molecular diagnosis platforms for TB, specifically the three molecular methods marketed in our country: Xpert MTB/Rif Ultra, Cobas MTB, and Logix Smart MTB. Therefore, the objective of the present study was to determine the diagnostic performance of Cobas MTB and Logix Smart MTB compared to Xpert MTB/RIF Ultra for the rapid diagnosis of PTB and EPTB.

Materials and methods

Research setting and design

A cross-sectional study of diagnostic tests was carried out. To write the manuscript, we followed the guidelines of the Standards for Reporting of Diagnostic Accuracy Studies (STARD) for diagnostic tests (Supplemental File).¹³ The study was carried out in the private clinical laboratory ROE,

at its locations in the metropolitan city of Lima, Peru.

Population and eligibility criteria

All pulmonary samples (sputum, bronchial aspirate, and bronchoalveolar lavage) and extrapulmonary samples (pleural fluid, urine, cerebrospinal fluid, breast biopsy, soft tissue abscess, ascitic fluid, pericardial fluid, serum, bone biopsy, and bone marrow) were collected from patients with presumptive TB who underwent smear microscopy, Xpert MTB/RIF Ultra and mycobacterial culture in Löwenstein–Jensen media, which were stored in a biobank at -70°C and subsequently and simultaneously subjected to Cobas MTB and Logix Smart MTB to determine their diagnostic performance.

Reference standard

For this study, we used the Xpert MTB/RIF Ultra (Cepheid, Sunnyvale, CA, USA) as the reference standard. This automated, nested polymerase chain reaction (PCR) test runs on the GeneXpert platform, detecting *Mycobacterium tuberculosis* (MTB) and rifampicin resistance. The process begins with sample collection (e.g., sputum), followed by mixing with a reagent that liquefies and inactivates the bacteria. The sample is loaded into the Xpert cartridge, which is inserted into the GeneXpert system for automated DNA extraction, amplification, and detection. The test targets the IS6110 and IS1081 multicopy sequences and mutations in the *rpoB* gene for rifampicin resistance. With a detection limit of 15.6 colony-forming units (CFU)/mL, results available in approximately 90 minutes, and a sensitivity and specificity in patients with positive smear microscopy of 90% and 96%, respectively.¹⁴

Index test

The first index test was the Cobas MTB test (Roche Diagnostics, Mannheim, Germany) is performed on the Cobas 5800/6800/8800 systems and is designed to detect *M. tuberculosis* complex DNA in respiratory samples. Initially, samples are liquefied and inactivated using Cobas Microbial Inactivation Solution (MIS). After incubation at room temperature, the samples undergo sonication to ensure proper liquefaction. The system then automates nucleic acid

extraction using magnetic glass particles, where bacterial DNA is bound and purified. This is followed by real-time PCR amplification using dual-target primers (16S rRNA and *esx* genes: *esxJ*, *esxK*, *esxM*, *esxP*, *esxW*), ensuring accurate detection. The results are automatically interpreted by the system's software, with a detection limit of 7.6–8.8 CFU/mL. The test is highly specific and requires the use of additional modules, such as the Cobas MTB-RIF/INH, for detecting rifampicin and isoniazid resistance by identifying mutations in the RRDRs and *inhA* promoter regions, respectively.^{3,15}

The second index test, Logix Smart MTB (Co-Diagnostics, Inc., South Salt Lake, UT, USA), is a qualitative real-time PCR test that detects *M. tuberculosis* by targeting the IS6110 and MPB64 genes from respiratory samples. With a detection limit of 2–8 copies/mL and a sensitivity and specificity of 98% and 99%, respectively, it has been validated for use with the QIAamp DNA Mini Kit and can be run on the CoDx Box thermal cycler or other compatible platforms. The test processes respiratory fluid samples, undergoing nucleic acid extraction followed by real-time PCR amplification to identify the *M. tuberculosis* complex. Fluorescent dyes monitor the amplification of these genetic markers, ensuring precise detection. The Logix Smart MTB kit includes ready-to-use reagents like a master mix, positive, and negative controls. The entire process, from extraction to result interpretation, takes 40–60 min depending on the equipment, with results typically available in 2 h. Designed for streamlined processing, the test ensures reliable detection of both pulmonary and extrapulmonary tuberculosis with minimal manual intervention.¹⁶

The details of the analytical process used for the detection of Logix Smart MTB, Cobas MTB, and Xpert MTB/RIF Ultra are described in Table 1.

Statistical analysis

The data were entered into a Microsoft Excel v.2016 spreadsheet and subsequently exported for analysis in Stata v17 software (StataCorp, College Station, TX, USA) and R studio v2023.03.1 (Posi Software, PBC, MA, USA).

The performance of the reference tests for diagnosing pulmonary and extrapulmonary TB was calculated using the area under the curve (AUC) of the receiver operating characteristic (ROC), sensitivity, specificity, positive and negative likelihood ratios, and positive and negative predictive values. Finally, we measured the proportion of absolute agreement using Cohen's Kappa coefficient.

Ethics

The protocol was approved by the research ethics committee of the Faculty of Health Sciences at the Private University of Tacna (code: 94/FAC). The ethics committee waived the requirement for informed consent due to the nature of the study, which prevented contacting patients from thawed samples stored for laboratory quality control purposes. Similarly, confidentiality of the database was ensured, which was used exclusively for research purposes.

Results

One hundred seventy-five samples from patients with presumptive TB were included, of which 58.3% ($n=102$) were pulmonary samples, and 41.7% ($n=73$) were extrapulmonary samples. The median age was 57.8 years (IQR: 39.3–71.0). MTB was identified by Xpert MTB/RIF Ultra in up to 27.4% ($n=48$) of the samples studied, the majority of which were pulmonary samples (85.4%). Regarding the extrapulmonary samples, the majority were pleural fluid (32.9%), urine (28.8%), cerebrospinal fluid (9.6%), and breast biopsy (8.2%) (Table 2). We observed that 27.4% were positive for Xpert MTB/RIF Ultra, with a median cycle threshold (Ct) of 17.4 (IQR: 16.2–20.4), 25.7% were positive for Cobas MTB, with a median Ct of 29.8 (IQR: 25.6–31.9); and 20.6% ($n=36$) were positive for Logix Smart MTB, with a median Ct of 29.4 (IQR: 26.8–36.8). In contrast, only 7.9% and 10.9% of samples were positive for Ziehl-Neelsen staining and mycobacterial culture, respectively. Only 1/48 patients with a positive result of the Xpert MTB/RIF Ultra presented resistance to rifampin, which was also detected by Cobas MTB, who presented resistance to rifampicin and isoniazid simultaneously. Additionally, one sample was resistant to isoniazid but not to rifampicin according to Cobas

Table 1. Characteristics of three molecular methods (Logix Smart MTB, Cobas MTB, and Xpert MTB/RIF Ultra) and their processing for the diagnosis of tuberculosis.

Characteristics	Logix Smart MTB	Cobas MTB ^a	Xpert MTB/RIF Ultra
Number of samples per processing	32 concentrated samples	20 concentrated samples (maximum capacity of 90 samples)	16 concentrated samples (by GeneXpert® XVI platform)
Type of conventional samples	Sputum, airway fluids	Sputum, sputum sediments, or BAL with NALC-NaOH	Sputum or Sputum sediments
Detection target	Presence of IS6110 and MPB64 genes	Presence of <i>esx</i> genes and 16S rARN gene/Mutation in <i>rpoB</i> , <i>katG</i> and <i>inhA</i> genes	Presence of IS6110 and IS1081 genes/ <i>rpoB</i> gene mutation
Detection limit	8 copies/μL	7.6 UFC/mL	15.6 UFC/mL
Minimum sample quantity	200 μL	1000 μL	500 μL
Sample pretreatment	Liquefaction with reagents inside the plate + proteinase K	Inactivation + liquefaction with MIS and sonication	Liquefaction with sample reagent (composed of sodium hydroxide and isopropanolol)
Extraction type	Magnetic beads in HollySys Extraction Kit ^b	Magnetic beads in Cobas 6800	Sonication
Treatment and extraction time	30 min	150 min	20 min
Conservation temperature	HollySys between 4–30°C. For Logix Smart, below –15°C	2–8°C	2–28°C
Preparation time for amplification	40 min	20 min	1 min
Type of PCR	qPCR	qPCR	semi-nested qPCR
Amplification time	80 min	180 min	90 min
Results information	Qualitative presence of mycobacteria	Qualitative presence of mycobacteria. Determines resistance to R and H	Qualitative presence of mycobacteria. Determines resistance to R

^aIf positive samples exist another detection kit for rifampicin and isoniazid that meets the same parameters would be run.

^bComponents on the HollySys plate are Guanidine Hydrochloride, Sodium Chloride, SDS, Tris-HCl, Isopropanolol, and absolute ethanol.

BAL, bronchoalveolar lavage; H, isoniazid; MIS, Cobas® Microbial Inactivating Solution; MTB, *Mycobacterium tuberculosis*; NALC-NaOH, N-acetyl-L-cysteine-sodium hydroxide; qPCR, Real-time polymerase chain reaction; R, rifampin; rRNA, ribosomal ribonucleic acid.

MTB, which was not detected using Xpert MTB/RIF Ultra (Table 3).

Accuracy of Cobas MTB and Logix Smart for TB diagnosis compared to positive molecular result of Xpert MTB/RIF Ultra

In pulmonary samples, smear microscopy presented an AUC/ROC of 0.70, with a sensitivity of 40.0% and specificity of 100.0% for the diagnosis of TB. Mycobacterial culture presented an AUC/ROC of 0.72, with a sensitivity of 45.2%

and specificity of 100%. In comparison, the Cobas MTB molecular test showed better performance with an AUC/ROC of 0.96, presenting a sensitivity and specificity of 92.9% and 100.0%, respectively. In contrast, the Logix Smart MTB showed an AUC/ROC of 0.89, with a sensitivity of 78.4% and a specificity of 100.0%. Regarding the performance of molecular tests in extrapulmonary samples, Logix Smart MTB showed the best performance with a sensitivity and specificity of 100.0%, while Cobas MTB presented an AUC/ROC of 0.91,

Table 2. Main characteristics of the sample origin according to type and result of Xpert MTB/RIF Ultra.

Variable	Total (n = 175)	Positive (n = 48)	Negative (n = 127)
Age, years (n = 99) ^a	57.8 (39.3–71.0)	52.7 (31.8–65.7)	59.7 (40.2–72.0)
Sample origin			
Pulmonary (%)	102 (58.3)	41 (85.4)	61 (48.0)
Bronchoalveolar lavage	33 (32.4)	13 (31.7)	20 (32.8)
Bronchial aspirate	33 (32.4)	16 (39.0)	17 (27.9)
Sputum	36 (35.3)	12 (29.3)	24 (39.3)
Extrapulmonary (%)	73 (41.7)	7 (14.6)	66 (52.0)
Cerebrospinal fluid	7 (9.6)	0 (0.0)	7 (10.6)
Urine	21 (28.8)	3 (42.8)	18 (27.3)
Pleural fluid	24 (32.9)	1 (14.3)	23 (34.8)
Ascitic fluid	5 (6.8)	0 (0.0)	5 (7.6)
Serum	1 (1.4)	0 (0.0)	1 (1.5)
Pericardial fluid	2 (2.7)	0 (0.0)	2 (3.0)
Bone biopsy	1 (1.4)	1 (14.3)	0 (0.0)
Breast biopsy	6 (8.2)	0 (0.0)	6 (9.1)
Bone marrow biopsy	1 (1.4)	0 (0.0)	1 (1.5)
Soft tissue abscess	4 (5.4)	1 (14.3)	3 (4.6)
Intrathoracic abscess	1 (1.4)	1 (14.3)	0 (0.0)
^a Median and interquartile range.			

with a sensitivity and specificity of 83.3% and 97.7%, respectively (Table 4; Figure 1). Of the seven extrapulmonary samples that tested positive by Xpert MTB/RIF Ultra (three urine, one ascitic fluid, one bone biopsy, one soft tissue abscess, and one intrathoracic abscess), Logix Smart MTB was positive in 100% of them, while Cobas MTB was positive in 83.3% (failing to detect one positive urine sample), while none of them were positive in culture or smear (Table 4; Figure 1).

Concordance of molecular tests for the diagnosis of pulmonary and extrapulmonary TB

The percentage of agreement between the Cobas MTB and the Xpert MTB/RIF Ultra was 97.6% for pulmonary samples and 95.9% for extrapulmonary samples, with Cohen's Kappa

index of 0.95 and 0.81, respectively. On the other hand, the Logix Smart MTB presented agreement of 93.1% for pulmonary samples and 95.5% for extrapulmonary samples, with Cohen's Kappa index of 0.83 and 0.81, respectively (Table 5).

We also analyzed the positive molecular results, in which the mycobacterial culture was negative (Table 6). Thirty-one samples were positive according to molecular methods but negative according to mycobacterial culture, of which 23/31 were pulmonary, and 8/31 were extrapulmonary. Among these 23 pulmonary samples, 23/23 were detected by Xpert MTB/RIF Ultra (Ct range: 16–26), while Cobas MTB detected 20/23 samples (Ct range: 19–35) and Logix Smart 19/31 samples (Ct range: 24–37) with five samples with indeterminate results. Among the

Table 3. Laboratory results of conventional and molecular tests for the diagnosis of pulmonary and extrapulmonary tuberculosis.

Diagnostic method	Total	Positive (%)	Negative (%)
Smear (%)	126	10 (7.9)	116 (92.1)
Mycobacteria culture (%)	175	19 (10.9)	156 (89.1)
Xpert MTB/RIF Ultra (%)	175	48 (27.4)	127 (72.6)
Threshold cycle (Ct) ^a	–	17.4 (16.2–20.4)	–
Rifampicin resistance (%)	48	1 (2.0)	–
Cobas MTB (%)	175	45 (25.7)	130 (74.3)
Threshold cycle (Ct) ^a	–	29.8 (25.6–31.9)	–
Isoniazid resistance (%)	45	2 (4.4)	–
Rifampicin resistance (%)	45	1 (2.2)	–
Logix Smart MTB (%) ^b	175	36 (20.6)	124 (70.9)
Threshold cycle (Ct) ^a	–	29.4 (26.8–36.8)	–

^aMedian and interquartile range.
^bLogix Smart presented 15 samples with indeterminate results, which were excluded for this analysis.
Ct, cycle threshold; MTB, *Mycobacterium tuberculosis*.

eight extrapulmonary samples, 7/8 were detected by Logix Smart MTB (Ct range: 25–37), 6/8 were positive by Cobas MTB (Ct range: 24–38), and 6/8 were positive by Xpert MTB/ RIF Ultra (Ct range: 16–26) (Table 6).

Finally, 100% of samples positive by smear microscopy were also positive by Xpert MTB/ RIF Ultra, Cobas MTB, and Logix Smart MTB methods. While 100%, 100%, and 94.7% of samples were positive by mycobacterial culture, they were also positive by Xpert MTB/RIF Ultra, Cobas MTB, and Logix Smart MTB, respectively.

Discussion

The study evaluated the diagnostic performance of Cobas MTB and Logix Smart MTB for the diagnosis of PTB and EPTB. We observed that both molecular tests had adequate diagnostic performance on total samples when compared to the Xpert MTB/RIF Ultra. We also observed that the performance of Logix Smart was superior (AUC/ ROC: 0.97) to Cobas MTB for the diagnosis of EPTB, although it was inferior for the diagnosis

of PTB (AUC/ROC: 0.91). These findings demonstrate that both tests meet the minimum diagnostic performance requirements for PTB and EPTB proposed by the WHO.¹⁷

The use of molecular tests such as Xpert MTB/ RIF for the initial diagnosis of PTB has been shown to reduce mortality compared with smear microscopy, although not statistically significant.⁴ Currently, this test is recommended by the WHO because it has the best diagnostic performance, so we decided to compare its performance with the Cobas MTB and the Logix Smart MTB. It should also be noted that according to the WHO, a defined case of TB is through the identification of MTB either by culture or by some molecular method.¹⁸

Xpert MTB/RIF Ultra has the advantage over Xpert MTB/RIF of slightly higher sensitivity (88% vs 83%), especially when smear microscopy is negative (63% vs 43%).^{18,19} It also has a lower detection limit of 16 CFU per milliliter (compared to 114CFU/mL with Xpert MTB/RIF), and it offers better precision in detecting rifampicin resistance due to the incorporation of a

Table 4. The area under the ROC curve and cut-off points of conventional and molecular tests compared to Xpert MTB/RIF Ultra.

Variable	AUC/ROC (CI 95%)	Sensitivity	Specificity	LHR +	LHR-	VPP	VPN	Prevalence
All samples								
Smear	0.685 [0.592–0.777]	37.04	100.0	-	0.629	100.0	85.30	21.0
Mycobacteria culture	0.697 [0.628–0.767]	39.58	100.0	-	0.604	100.0	81.40	27.0
Cobas MTB	0.954 [0.914–0.994]	91.67	99.21	116.4	0.084	97.80	96.90	27.0
Logix Smart MTB	0.896 [0.835–0.957]	80.95	98.31	47.76	0.193	94.40	93.50	26.0
Pulmonary samples								
Smear	0.700 [0.602–0.798]	40.00	100.0	-	0.600	100.0	81.30	28.0
Mycobacteria culture	0.726 [0.650–0.802]	45.24	100.0	-	0.547	100.0	78.50	33.0
Cobas MTB	0.964 [0.924–0.999]	92.86	100.0	-	0.071	100.0	96.60	33.0
Logix Smart MTB	0.891 [0.824–0.959]	78.38	100.0	-	0.216	100.0	90.80	32.0
Extrapulmonary samples								
Smear	-	-	-	-	-	-	-	-
Mycobacteria culture	-	-	-	-	-	-	-	-
Cobas MTB	0.905 [0.740–0.999]	83.33	97.67	35.83	0.170	83.30	97.70	12.0
Logix Smart MTB	0.974 [0.939–0.999]	100.0	94.87	19.5	0.000	100.0	100.0	16.0

AUC, area under the curve; MTB, *Mycobacterium tuberculosis*; ROC, receiver operating characteristic; LHR, Likelihood ratio; 95% CI, 95% confidence interval.




Figure 1. Radar graph summarizing the diagnostic performance of conventional and molecular tests in pulmonary and extrapulmonary samples compared to Xpert MTB/RIF Ultra.

Table 5. Percentage of agreement of the Cobas MTB and Logix Smart MTB test with the Xpert MTB/RIF Ultra.

Comparator test	Sample origin	Agreement (95% CI)	Cohen's Kappa index (95% CI)
Cobas MTB vs Xpert MTB/RIF Ultra	Total	97.14 (94.65–99.64)	0.926 (0.86–0.99)
	Pulmonary	97.62 (94.9–100)	0.945 (88.3–100)
	Extrapulmonary	95.92 (90.1–100)	0.810 (0.54–100)
Logix Smart MTB vs Xpert MTB/RIF Ultra	Total	93.75 (88.02–99.48)	0.830 (0.72–0.94)
	Pulmonary	93.10 (86.3–99.8)	0.831 (0.71–0.95)
	Extrapulmonary	95.45 (84.1–100)	0.807 (0.54–100)
Cobas MTB vs Logix Smart MTB	Total	96.25 (90.93–100)	0.897 (0.80–0.98)
	Pulmonary	95.69 (89.4–100)	0.893 (0.79–0.99)
	Extrapulmonary	97.73 (87.1–100)	0.906 (0.69–100)

95% CI, 95% confidence interval; MTB, *Mycobacterium tuberculosis*.

Table 6. List of patients in whom MTB was identified by some molecular method, but not by culture.

Patient	Sample origin	Xpert MTB/RIF Ultra (Ct)	Cobas MTB (Ct)	Logix Smart MTB (Ct)
1	Pulmonary	17	31	Indeterminate
2	Pulmonary	16	23	Indeterminate
3	Extrapulmonary	26	Negative	Indeterminate
4	Pulmonary	16	19	Indeterminate
5	Pulmonary	17	31	Indeterminate
6	Pulmonary	20	30	Indeterminate
7	Extrapulmonary	17	32	27
8	Extrapulmonary	18	34	29
9	Pulmonary	21	34	Negative
10	Pulmonary	24	32	36
11	Pulmonary	26	35	Negative
12	Pulmonary	16	22	28
13	Pulmonary	23	36	Negative
14	Extrapulmonary	16	25	28
15	Pulmonary	16	20	27
16	Extrapulmonary	20	30	34
17	Pulmonary	16	27	27
18	Pulmonary	16	28	29
19	Pulmonary	25	35	Negative
20	Pulmonary	27	Negative	Negative
21	Pulmonary	16	25	24
22	Pulmonary	24	Negative	Negative
23	Pulmonary	21	32	35
24	Extrapulmonary	16	24	25
25	Pulmonary	19	32	32
26	Pulmonary	17	27	37
27	Pulmonary	19	32	34
28	Pulmonary	20	Negative	Negative
29	Pulmonary	20	30	34
30	Extrapulmonary	Negative	38	34
31	Extrapulmonary	Negative	Negative	37

Ct, cycle threshold; MTB, *Mycobacterium tuberculosis*.

melting temperature analysis using four probes that identify mutations in the *rpoB* gene.¹⁹

Cobas MTB showed a sensitivity of 89.2% (95% CI: 81.7%–93.9%) and specificity of 98.6%, compared to respiratory samples with positive cultures.²⁰ At the same time, the reported sensitivity of Cobas MTB in smear-negative samples was 63.0%–81.8%.

Unlike Xpert MTB-RIF Ultra, Cobas MTB-RIF/INH has the advantage of being able to confirm true MDR-TB since Xpert MTB-RIF Ultra only estimates MDR based on rifampicin resistance. A study evaluated its performance for this objective, showing a sensitivity of 77.8% and 90% for detecting resistance to isoniazid and rifampicin, respectively.²¹ The cases of false susceptibility to isoniazid in the Cobas MTB-RIF/INH were due to nontarget cobas mutations. Nevertheless, in terms of performance, Cobas MTB has shown greater analytical sensitivity than Xpert MTB/RIF and is similar to Genotype MTBDR plus for resistance to rifampicin and isoniazid.²²

According to the manufacturer's instructions, Cobas MTB-RIF/INH is only validated for respiratory samples; however, in our study, we observed adequate performance for extrapulmonary samples, similar to the Xpert, although inferior to the Logix Smart MTB.

Compared to other medium-complexity molecular methods supported by the WHO, the Cobas MTB-RIF/INH assay requires 3.7 h to process 94 samples and has the advantage of being able to simultaneously analyze 384 samples on the Cobas 6800 and 960 samples on the Cobas 8000. In addition, PCR can be used for MTB, hepatitis B, hepatitis C, and HIV simultaneously. One of its main limitations in sample preparation with Cobas MTB was sonication, which was the reason why it received a lower score regarding sample preparation by users (86% approval).^{23,24} Another disadvantage is that DNA extraction must be repeated if resistance testing is required, which can lead to increased costs and time, mainly in settings with high incidences of DR-TB. In addition, this technique requires more space and has special electrical requirements than the Xpert MTB/RIF Ultra.

In our study, we observed that Logix Smart MTB presented the best diagnostic performance for

extrapulmonary samples, although according to the supplier, this method has only validated for pulmonary samples.¹⁶ To our knowledge, there are no studies that have evaluated the diagnostic performance of Logix Smart MTB in this context; however, in pulmonary samples, it presented a lower performance, with a significant proportion of indeterminate or false negative results, most of which were detected using Xpert MTB/RIF Ultra or Cobas MTB, possibly due to a lower mycobacterial load, since some of these discordant results presented high Ct in the Cobas MTB. As a limitation, Logix Smart has only been validated with the QIAGEN and with the CoDx Box thermal cycler. Although it can provide results in less than 2 h, it is a method of moderate complexity and is not completely automated as is the case for the Xpert MTB/RIF Ultra. Furthermore, this test does not detect resistance to anti-TB drugs, which may be a major limitation.

This research is part of the third pillar of "The End TB strategy"²⁵ since although Xpert MTB/RIF Ultra is a test with high sensitivity and specificity for the diagnosis of TB, its implementation has not yet been universal, probably due to associated costs; Cobas MTB and Logix Smart MTB could be approximately 30–50% less expensive, and have the advantage of being able to be mass-processed, and with comparable performance. Likewise, although conventional tests like smear microscopy are more economical, they have low sensitivity in identifying TB cases.²⁶

This study has some strengths. It is the only study that compares the performance of the Xpert MTB/RIF Ultra, Cobas MTB, and Logix Smart MTB for diagnosing PTB and EPTB, even though some of these methods are out of license. In addition to being developed in a high-performance reference laboratory in a country with high TB endemicity, such as Peru, the generalizability of the results is further supported by the inclusion of a diverse range of samples (both pulmonary and extrapulmonary) and the use of different molecular methods. These factors provide a broader context that enhances the applicability of the results to different clinical and geographical settings.

Some of the limitations of this study are the following. Sample processing was carried out after freezing and thawing a limited volume of samples, and it is unknown whether this process could have affected the performance of the tests. It was also

not possible to know relevant clinical information such as sex, the presence of symptoms, clinical evolution, history of TB, and HIV infection or treatment with antiretrovirals because we did not have access to medical records and were unable to interview the patients to collect this information. Furthermore, the limited number of samples, due to the limited budget available, could lead to inaccuracy in the reported results. More studies with a greater number of samples analyzed are needed, especially in specific scenarios and with clinical information at the time of diagnosis and its outcomes. Additionally, a sample size calculation was not performed, which could also affect the precision of the results. Future studies should ensure proper sample size calculations to strengthen statistical validity and the detection of significant diagnostic differences.

Conclusion

We found that the Cobas MTB and Logix Smart MTB tests demonstrated adequate performance for tuberculosis diagnosis when compared to the Xpert MTB/RIF Ultra. For pulmonary samples, Cobas MTB showed a higher diagnostic yield, while for extrapulmonary samples, the Logix Smart MTB test was the most suitable. Based on our results, we conclude that molecular tests complement each other, and their selection should be guided by the type of sample being analyzed rather than a one-size-fits-all approach.

Declarations

Ethics approval and consent to participate

The protocol was approved by the research ethics committee of the Faculty of Health Sciences at the Private University of Tacna (code: 94/FAC). The ethics committee waived the requirement for informed consent due to the nature of the study, which prevented contacting patients from thawed samples stored for laboratory quality control purposes. Similarly, confidentiality of the database was ensured, which was used exclusively for research purposes.

Consent for publication

Not applicable.

Author contributions

Miguel Hueda-Zavaleta: Conceptualization; Formal analysis; Methodology; Resources;

Supervision; Writing – original draft; Writing – review & editing.

Juan Carlos Gomez de la Torre: Conceptualization; Investigation; Resources; Supervision; Writing – review & editing.

Diana Minchón-Vizconde: Data curation; Writing – original draft.

Claudia Barletta-Carrillo: Investigation; Writing – original draft.

Cesar Copaja-Corzo: Data curation; Methodology; Writing – original draft.

Gustavo Tapia-Sequeiros: Data curation; Formal analysis; Writing – original draft.

Cinthya Flores: Investigation; Writing – original draft.

Cristian Piscoche: Investigation; Writing – original draft.

Cecilia Miranda: Investigation; Writing – original draft.

Ada Mendoza: Investigation; Writing – original draft.

Vicente A. Benites-Zapata: Methodology; Writing – review & editing.

Acknowledgements

None.

Funding

The authors disclosed receipt of the following financial support for the research, authorship, and/or publication of this article: This study received funding through the Private University of Tacna, Tacna, Peru (323-2022-UPT-CU).

Competing interests

JCGdlT is a speaker at Cepheid, manufacturer of the GeneXpert platforms. The other authors declare that they have no competing interests.

Availability of data and materials

The datasets used and analyzed during the current study are available from the corresponding author on reasonable request.

ORCID iDs

Miguel Angel Hueda Zavaleta  <https://orcid.org/0000-0002-8049-7787>

Juan Carlos Gomez de la Torre  <https://orcid.org/0000-0003-4566-2027>
 Diana Julissa Minchon Vizconde  <https://orcid.org/0000-0003-2256-7609>
 Claudia Barletta-Carrillo  <https://orcid.org/0000-0001-9982-5368>
 Cesar Copaja-Corzo  <https://orcid.org/0000-0002-3497-0158>
 Gustavo Tapia Sequeiros  <https://orcid.org/0000-0002-9277-2731>
 Cinthya Gisele Flores Flores  <https://orcid.org/0009-0009-3452-3276>
 Nilver Cristian Piscoche Botello  <https://orcid.org/0000-0002-9907-8013>
 Cecilia Miranda  <https://orcid.org/0009-0000-7904-5735>
 Ada Patricia Mendoza Farro  <https://orcid.org/0009-0008-7202-856X>
 Vicente A. Benites-Zapata  <https://orcid.org/0000-0002-9158-1108>

Supplemental material

Supplemental material for this article is available online.

References

1. Pan American Health Organization. *Tuberculosis in the Americas, 2018*. Washington, DC: PAHO, <https://iris.paho.org/handle/10665.2/49510> (2018, accessed 24 September 2023).
2. Floyd K, Glaziou P, Zumla A, et al. The global tuberculosis epidemic and progress in care, prevention, and research: an overview in year 3 of the End TB era. *Lancet Respir Med* 2018; 6(4): 299–314.
3. World Health Organization. *Module 3: Diagnosis—Rapid Diagnostics for Tuberculosis Diagnosis, 2021 Update*. Geneva, Switzerland: WHO, <https://www.who.int/publications/item/9789240029415> (2021, accessed 31 March 2023).
4. Di Tanna GL, Miotto P, Migliori GB, et al. Effect of Xpert MTB/RIF on clinical outcomes in routine care settings: individual patient data meta-analysis. *Lancet Glob Health* 2019; 7(2): e191–e199.
5. Zifodya JS, Kreniske JS, House PS, et al. Xpert Ultra versus Xpert MTB/RIF for pulmonary tuberculosis and rifampicin resistance in adults with presumptive pulmonary tuberculosis. *Cochrane Database Syst Rev* 2021; 2: CD009593.
6. Bisognin F, Lombardi G, Lombardo D, et al. Improvement of *Mycobacterium tuberculosis* detection by Xpert MTB/RIF Ultra: a head-to-head comparison on Xpert-negative samples. *PLoS One* 2018; 13(8): e0201934.
7. Mekkaoui L, Vandersmissen G, Courcol R, et al. Performance of Xpert MTB/RIF Ultra for diagnosis of pulmonary and extra-pulmonary tuberculosis, one year of use in a multi-centric hospital laboratory in Brussels, Belgium. *PLoS One* 2021; 16(4): e0249734.
8. Donovan J, Phu Huong Lan N, Dinh Ple T, et al. Xpert MTB/RIF Ultra versus Xpert MTB/RIF for the diagnosis of tuberculous meningitis: a prospective, randomised, diagnostic accuracy study. *Lancet Infect Dis* 2020; 20(3): 299–310.
9. Huang M, Gao Y, Wang X, et al. Diagnostic accuracy of Xpert MTB/RIF Ultra for tuberculous meningitis in a clinical practice setting of China. *Diagn Microbiol Infect Dis* 2021; 100(1): 115306.
10. Wang G, Chen L, Wu X, et al. Xpert MTB/RIF Ultra improved the diagnosis of paucibacillary tuberculosis: a prospective cohort study. *J Infect* 2019; 78(4): 311–316.
11. Perez-Risco D, Montuenga LM, Saiz-Montero M, et al. Evaluation of the Xpert MTB/RIF Ultra assay for direct detection of *Mycobacterium tuberculosis* complex in smear-negative extrapulmonary samples. *J Clin Microbiol* 2018; 56(9): e00659-18.
12. Rindi L. Rapid molecular diagnosis of extra-pulmonary tuberculosis by Xpert/RIF Ultra. *Front Microbiol* 2022; 13: 817661.
13. Cohen JF, Korevaar DA, Altman DG, et al. STARD 2015 guidelines for reporting diagnostic accuracy studies: explanation and elaboration. *BMJ Open* 2016; 6(11): e012799.
14. Cepheid. Xpert MTB RIF Ultra Brochure. Sunnyvale, United States, <https://cepheid.widen.net/view/pdf/s6lwefaviz/Cepheid-Xpert-MTB-RIF-Ultra-Brochure-CE-IVD-3098-Spanish.pdf?t.download=true&u=escyfu> (2020, accessed 25 September 2023).
15. Roche Diagnostics. Cobas MTB. Basel, Switzerland, <https://diagnostics.roche.com/global/en/products/params/cobas-mtb.html#productSpecs> (2022, accessed 31 March 2022).
16. Co-Diagnostics. Logix Smart TM MTB. United States: Co-Dx, www.codiagnostics.com (2022, accessed 31 March 2022).

17. World Health Organization. High priority target product profiles for new tuberculosis diagnostics: report of a consensus meeting. Geneva, Switzerland: WHO, <https://www.who.int/publications/i/item/WHO-HTM-TB-2014.18> (2014, accessed 23 September 2023).
18. World Health Organization. *Treatment of tuberculosis: guidelines*. 4th ed. Geneva: WHO, <https://www.ncbi.nlm.nih.gov/books/NBK138748/> (2010, accessed 31 March 2023).
19. Zhang M, Xue M and He J. Diagnostic accuracy of the new Xpert MTB/RIF Ultra for tuberculosis disease: a preliminary systematic review and meta-analysis. *Int J Infect Dis* 2020; 90: 35–45.
20. Horita N, Miyazawa N, Yoshiyama T, et al. Sensitivity and specificity of Cobas TaqMan MTB real-time polymerase chain reaction for culture-proven *Mycobacterium tuberculosis*: meta-analysis of 26999 specimens from 17 studies. *Sci Rep* 2015; 5: 18113.
21. Nadarajan D, Schlager D, Clark AE, et al. Evaluation of the Roche Cobas MTB and MTB-RIF/INH assays in samples from Germany and Sierra Leone. *J Clin Microbiol* 2021; 59(5): e02983-20.
22. Aono A, Yano S, Yoshimura A, et al. Clinical evaluation of the Cobas® MTB-RIF/INH reagent and the cobas® 6800 for the detection of isoniazid and rifampicin resistance. *Tuberculosis* 2022; 134: 102199.
23. De Vos M, Collins S, Frahm N, et al. Comparative analytical evaluation of four centralized platforms for the detection of *Mycobacterium tuberculosis* complex and resistance to rifampicin and isoniazid. *J Clin Microbiol* 2021; 59(3): e02168-20.
24. David A, Shetty S, Anthony RM, et al. Feasibility, ease-of-use, and operational characteristics of World Health Organization–recommended moderate-complexity automated nucleic acid amplification tests for the detection of tuberculosis and resistance to rifampicin and isoniazid. *J Mol Diagn* 2023; 25(1): 46.
25. World Health Organization. The End TB Strategy. Geneva, Switzerland: WHO, <https://www.who.int/publications/i/item/WHO-HTM-TB-2015.19> (2015, accessed 24 September 2023).
26. Figueredo LJA, Miranda SS, Santos LBD, et al. Cost analysis of smear microscopy and the Xpert assay for tuberculosis diagnosis: average turnaround time. *Rev Soc Bras Med Trop* 2020; 53: e20200314.

Visit Sage journals online
[journals.sagepub.com/
 home/tai](https://journals.sagepub.com/home/tai)

 Sage journals