

Published in final edited form as:

Nat Protoc. 2009 ; 4(12): 1771–1783. doi:10.1038/nprot.2009.184.

High-throughput haplotype determination over long distances by Haplotype Fusion PCR and Ligation Haplotyping

Daniel J. Turner¹ and Matthew E. Hurlles¹

Daniel J. Turner: djt@sanger.ac.uk; Matthew E. Hurlles: meh@sanger.ac.uk

¹Wellcome Trust Sanger Institute, Wellcome Trust Genome Campus, Cambridge, UK.

Abstract

When combined with Haplotype Fusion PCR (HF-PCR), Ligation Haplotyping is a robust, high-throughput method for empirical determination of haplotypes, which can be applied to assaying both sequence and structural variation over long distances. Unlike alternative approaches to haplotype determination, such as allele-specific PCR and long PCR, HF-PCR and Ligation Haplotyping do not suffer from mispriming or template switching errors. In this method, HF-PCR is used to juxtapose DNA sequences from single molecule templates, that contain single nucleotide polymorphisms (SNPs) or paralogous sequence variants (PSVs) separated by several kilobases. HF-PCR employs an emulsion-based fusion PCR reaction, which can be performed rapidly, and in a 96-well format. Subsequently, a ligation-based assay is performed on the HF-PCR products to determine haplotypes. Products are resolved by capillary electrophoresis. Once optimized, the method is rapid to perform, taking a day and a half to generate phased haplotypes from genomic DNA.

Introduction

Knowledge of the combination of SNPs on the same parental chromosome is extremely valuable to several areas of research: for example, haplotypes can be used i) in evolutionary and population-genetic studies, to detect the influence of selection 1 and population migration 2; ii) for linkage disequilibrium mapping 3; and iii) in association studies, to investigate whether particular combinations of SNPs confer greater susceptibility to diseases (e.g. systemic lupus erythematosus 4). Unlike genotypes, for which several well-established assays are available, there are few reliable empirical methods for the determination of the phase of haplotype for diploid organisms, particularly over long distances, due to the need of a haplotyping assay to function on individual strands of DNA, without giving rise to artefacts.

Commonly used experimental approaches to haplotype determination involve allele-specific PCR 10, long PCR 11 or a combination of both 12. However, optimizing a PCR reaction so that it is truly allele-specific is laborious, as it is difficult to avoid mispriming. Additionally, the risk of 'jumping PCR' increases with increasing amplicon length, creating artificial recombinants by template switching 13. The efficiency and robustness of long PCR also decreases with increasing amplicon length 14. Consequently, approaches for haplotype determination that are based on allele-specific PCR or long PCR are unsuitable over

Corresponding author: Daniel J. Turner Wellcome Trust Sanger Institute, Wellcome Trust Genome Campus, Hinxton, Cambridge, UK, CB10 1SA Tel: +44 (0)1223 494763 Fax: +44 (0)1223 494919 djt@sanger.ac.uk.

Author contributions MEH and DJT designed this study. All experiments were performed by DJT. DJT and MEH performed the analysis and wrote the manuscript.

Competing Financial Interests The authors declare that they have no competing financial interests

distances of several kilobases, and so statistical haplotyping approaches, which reconstruct haplotypes from genotype information, have gained prominence (for a review of statistical haplotyping methods see ⁶).

Over short distances, statistical haplotyping methods are generally reliable^{7, 8}, but over longer distances, or in parts of the genome with low linkage disequilibrium, statistical approaches suffer from switch errors, rendering them inaccurate. Thus there is a need for a reliable empirical method for the determination of haplotype status, whose accuracy does not decrease with increasing distance ⁹. Such a method is Ligation Haplotyping ¹⁵.

The first stage of Ligation Haplotyping uses a high-throughput version of Haplotype Fusion PCR (HF-PCR) to generate condensed haplotypes, by juxtaposing DNA sequences from the same single molecule, thus retaining phase information ¹⁴ (Fig. 1a). The reaction is carried out in a water-in-oil emulsion – each aqueous compartment is independent from the others, allowing millions of individual reactions to take place in a single PCR tube ^{17,18}. If we know the diameter of these aqueous compartments, we can calculate the total number, and can then perform the PCR reaction with a low compartment : template ratio. In this way, the majority of occupied compartments can be made to contain a single target locus, which ensures that PSVs or SNPs from single molecules become fused together. The HF-PCR reaction generates a condensed haplotype, in which PSVs or SNPs that lie several kb apart in the genomic DNA are brought to within 100-150 bases of one another.

Many emulsification methods are not suitable for a high throughput assay because they use large reaction volumes, which makes them expensive, chiefly because of the large quantity of polymerase used ^{14,17,19}, and also because the techniques used to create the emulsions are only suitable for small numbers of simultaneous reactions. This led us to develop a method of emulsification that is effective on low reaction volumes, allowing us to use 2 units of polymerase per reaction instead of more than 35 units for other published methods ^{17,19}. Additionally, 96 emulsions can be generated simultaneously, in 150 seconds, and so this method is both relatively inexpensive and rapid. .

Both long PCR and HF-PCR require the nucleotides of interest to be present on the same template strand, and so for both techniques it is essential to have a good quality, high molecular weight DNA template. If long PCR is attempted on low molecular weight templates, the frequency of template switching will increase ¹³, leading to false positives. Because single template molecules are amplified independently during Ligation Haplotyping, template switching is prevented: low molecular weight DNA will result in no amplicon being produced. Unlike long PCR, the amplification efficiency of HF-PCR does not decrease over long distances – during the initial stages of the reaction, the two target regions are amplified independently within the aqueous compartment. The fusion reaction then begins once the concentration of these amplicons has increased to a sufficient level ¹⁴. The length of haplotype that can be condensed by HF-PCR should be limited only by the size of the template DNA used in the emulsion. We investigated the shearing effect of emulsion preparation on genomic DNA 100-150kb in length by pulsed-field gel electrophoresis, and were unable to detect any reduction in fragment size¹⁴, which indicates that haplotype determination in this way should be possible over a longer range than the 25kb offered by long PCR.

The high-throughput Ligation Haplotyping assay itself is performed on the condensed haplotypes generated by HF-PCR (Fig. 1b). There are several assays that exploit the high specificity of ligases at the 3' side of a ligation junction, such as the ligase detection reaction ²⁰, GoldenGate assay ²¹ and SNPlex assay ²². In Ligation Haplotyping, we make use of the ligase's ability to be highly discerning at both the 5' and 3' sides of the ligation junction:

allele-specific oligos for each SNP anneal to the condensed haplotype at either side of a central bridging oligo, and are ligated (Fig. 1). The assay described here was designed so that this bridging oligo is the same oligo that was used as the fusion primer in the prior reaction. However, in cases where it would be preferable to use PCR primers that do not about the PSVs or SNPs directly, it is possible to use a bridging oligo that differs from this PCR primer.

The product is linearly amplified by cycles of denaturation and annealing / ligation. Unlike exponential amplification, during this reaction, the ligation product itself does not serve as a template in subsequent cycles. Thus errors caused by misligation are not propagated – each is a de novo event, and this helps to maintain a high level of specificity during the reaction. Each allele-specific oligos used in the ligation reaction has a characteristic length so that, once ligated to the central bridging oligo, the size of the products reveals the identity of the haplotypes present. In addition, these oligos have universal sequences at their non-ligating termini, to facilitate simultaneous PCR amplification of all products with fluorescently labelled primers. The fluorescent labels allow the PCR-amplified ligation products to be identified by capillary electrophoresis (Fig. 1c).

The HF-PCR reaction should be amenable to multiplexing. For a given primer pair, there is only one target locus in every ~1000 aqueous compartments, and so additional sets of primers will not compete for resources, because amplification will proceed in separate compartments. Additionally, as with other well-established ligation based assays, the Ligation Haplotyping reaction should also be compatible with high levels of multiplexing. In the example presented in this manuscript, oligos were designed so that haplotypes generate peaks differing by 5 bases. This means that a window of 20 nucleotides must be reserved for each set of PSVs / SNPs. To increase the multiplexing capability of peak visualization, it is possible to reduce the difference in size between different haplotypes, to 2 bases, reducing the overall window size to 8 nucleotides. Additionally, it would be possible to use PCR primers with different fluorescent labels. Up to 4 different fluorophores can be used within each dye set, and in this way, overlapping size windows can be used. The assay could easily be adapted for use on addressable arrays, which will increase throughput and decrease cost further.

Experimental design

HF-PCR

Short regions are amplified more efficiently in the HF-PCR, and so it is necessary to select primer pairs that target regions of < 150bp. Primers must be designed for the two SNP- or PSV-containing loci separately, choosing similar annealing / melting temperatures for all primers using an oligonucleotide calculator (such as <http://www.cnr.berkeley.edu/~zimmer/oligoTMcalc.html>). A melting temperature of approximately 65°C is appropriate. The specificity of primers should be checked by in silico PCR (<http://genome.ucsc.edu/cgi-bin/hgPcr?command=start>).

Once all four primers have been designed, the ‘fusion primer’ is designed by appending the reverse complementary sequence of one primer from one locus onto the 5’ end of one primer from the second locus (Fig. 1a). For locus 1, the forward and reverse PCR primers are F1 and R1 respectively. For locus 2, they are F2 and R2 respectively. A fusion primer that would join the two loci together is F2’-R1 (i.e. the R1 sequence with the reverse complement of F2 attached to its 5’ end). In the PCR reaction, only primers F1, R2 and F2’-R1 are used. In the initial rounds, only locus 1 can amplify, but the strand that is generated by extension of primer F1 acquires the primer F2 sequence at its 3’ end, and so the whole amplicon is able to prime on locus 2. The F2’-R1 fusion primer is present at a lower

concentration than primers F1 and R2, and becomes depleted more rapidly. After a few rounds of PCR, amplification of the fused loci proceeds using primers F1 and R2.

There is no requirement for the 3' termini of primers at either locus to be adjacent to the SNPs of interest, or for the bridging oligo to be the same as the F2'-R1 fusion primer. The main consideration here is that the further the 3' terminus of F2' is from the locus 2 SNP, and the further the R1 is from the locus 1 SNP, the further apart the SNPs will be in the condensed haplotype, and hence the longer the bridging oligo will need to be in the ligation reaction, as it is the bridging oligo that must abut both SNPs.

Emulsion formation

When prepared singly by dropwise addition of aqueous phase to oil phase while stirring at 1000rpm, the emulsions produced contain aqueous compartments which, when viewed by light microscopy and compared to fluorescent marker beads, have an average diameter of approximately 15 μ m²³. If we assume that these compartments are spherical, their volume is 1.77×10^{-15} m³. If we also assume that aqueous compartments are densely packed in the emulsion, then there should be around 5.7×10^5 aqueous compartments per microlitre. If the speed of stirring is increased, the average diameter of aqueous compartments decreases. Increasing the speed of vortexing when preparing emulsions in 96-well format produces the same effect.

For the same mass of DNA in the emulsion, increasing the number of compartments by reducing their volume reduces the frequency of compartments that are occupied by more than one template, and so the background signal is minimised. But because only aqueous compartments that contain both polymerase and template can produce amplification products, increasing the number of compartments reduces the overall efficiency of the reaction. Thus it is necessary to find a compromise between efficient amplification and low background signal. We found that this was achieved by making emulsions with an average aqueous compartment diameter of approximately 5-10 μ M.

Emulsions made using 200ng genomic DNA contain approximately 30,000 copies of each autosomal target locus. Consequently only 1 in around 900 aqueous compartments contain a suitable template molecule for fusion PCR in the *IRF5* assay, and so in most compartments, amplicons are derived from a single template molecule. When made in 96-well format, the volume of aqueous and oil phases are reduced, but the final concentration of reagents is unchanged. Thus, only 50ng of template DNA is required per well, and the ratio of 1 template molecule in 900 aqueous compartments is maintained.

Ligations

The allele-specific oligos used in the ligation reaction are designed to be complementary to the HF-PCR product, and to have an annealing temperature of approximately 65°C. This is achieved simply by obtaining the sequence of the strand of the fused template that is complementary to the bridging oligo (see Fig 1b.), measuring the annealing temperature of the strand between the SNP and the end of the strand in an oligonucleotide calculator (e.g. <http://www.cnr.berkeley.edu/~zimmer/oligoTmcalc.html>) and reducing the length from the non-ligating end until the desired T_m is achieved. This is performed for both SNPs at both sides of the bridging oligo. Allele-specific oligos designed in this way may differ in length by one or two nucleotides, because the SNPs themselves will influence the T_m of the oligo. This is not problematic, since it is necessary to make the allele specific oligos from each locus differ from each other in length anyway. In the example given here, we added non-complementary 'spacer nucleotides' to one oligo from each pair, so that the oligos at the 5' end of the bridging oligo differed in length by 5 nucleotides and those at the 3' end of the

bridging oligo differed by 10 nucleotides. The universal annealing sequence is then added to the non-ligating end of the allele-specific oligos.

Thus the overall length of the allele-specific oligos are governed by the universal sequence, the spacer region and their T_m , so will differ from assay to assay. The length of the bridging oligo is governed by the distance between the SNPs. When allele-specific oligos ligate to the central bridging oligo, four different combinations are possible, and the combinations differ in length by 5 nucleotides.

Prior to PCR amplification, it is necessary to digest unligated oligonucleotides using exonucleases, and so it is necessary to add blocking groups to the non-ligating end of each allele-specific oligo. In this way, after successful ligation, products are protected from digestion at both ends, whereas unligated oligos are not. The exact nature of the blocking groups is not critical, but we have found suitable modifications to be four 2' O-methyl uracil bases added to the 3' end 5'-specific oligos, and a single amino-C6 modification added to the 5' end of 3'-specific oligos.

Template DNA

For a successful HF-PCR reaction, it is essential that the strands of template DNA are longer than the distance between the SNPs that are being assayed. Commercial DNA preparations typically have a fragment size in excess of 100kb, whereas biological specimens can be shorter, depending upon the extraction method used. If in doubt, fragment size and DNA integrity can be assayed by pulsed-field gel electrophoresis 14. Buffers that contain EDTA should be avoided, because this can inhibit PCR reactions.

Controls

To monitor reagent contamination, it is necessary to prepare reaction blanks for all enzymatic reactions. Here, an equivalent volume of ultrapure water is used instead of the DNA sample. The formation of reaction products indicates contamination.

Materials

Reagents

Oil phase for preparation of emulsions singly—

- Span 80 (Sigma cat. no. S6760)
- Tween 80 (Sigma cat. no. P8074)
- Triton X-100 (Sigma cat. no. T9284)
- Light mineral oil (Sigma cat. no. M5904)

Oil phase for preparation of emulsions in 96-well format—

- Silicone polyether / cyclopentasiloxane (Dow Corning cat. no. DC5225C)
- Cyclopentasiloxane / trimethylsiloxysilicate (Dow Corning cat. no. DC749)
- AR20 silicone oil (Aldrich cat. no. 10838)
- 10 μ M fluorescent beads (Beckman Coulter cat. no. 6601329)

Aqueous phase for emulsions—

Phusion DNA polymerase (Finnzymes; NEB cat. no. F-530L) **CRITICAL** – avoid the use of polymerases that do not generate blunt ended products.

5X Phusion HF buffer (Finnzymes; NEB cat. no. F-518S)
 50mM MgCl₂ (Finnzymes; supplied with NEB cat. no. F-518S)
 dNTPs (NEB cat. no. N0447L)
 Genomic DNA

HF-PCR primers for *IRF5* SNPs (Sigma)—

IRF5_F1 5' CGGGATGAAGACTGGAGTAGG
 IRF5_R2 5' GACAAGGAGGAGTAAGCAAGGAAC
 IRF5_F2'-R1 5' GGGTGCCTACAGCAGGGTTCTGACCCTGGCAGGTCC

Disruption of emulsions—

diethyl ether (Sigma cat. no. 309958) or hexane (Sigma cat. no. 34859) – CAUTION:
 highly flammable
 Proteinase K (Sigma cat. no. P4850)

Ligation Haplotyping reaction—

10X T4 ligase buffer (supplied with NEB cat. no. M0202S)
 T4 polynucleotide kinase (NEB cat. no. M0201L)
Thermus thermophilus (*Tth*) ligase (ABGene cat. no. AB-0325)
 10x *Tth* ligase buffer (supplied with ABGene cat. no AB-0325)

Ligation oligos for *IRF5* SNPs—

IRF5_F2'-R1 (Sigma, as above)
 IRF5_5T (IDT;
 TGAGCATCACCAATGTGACCTATCCTATCCGTGCCAGCAAGATCCAATCTAG
 A-2'UOMe-2'UOMe-2'UOMe-2'UOMe
 IRF5_5C (IDT;
 CGAGCATCACCAATGTGACCGTGCCAGCAAGATCCAATCTAGA-2'UOMe-2'
 UOMe-2'UOMe-2'UOMe)
 IRF5_3G (IDT; AmC6-GGGTTCCTAAGGGTTGGAGAAAGCGAGCTCGGGG)
 IRF5_3T (IDT; AmC6-
 GGGTTCCTAAGGGTTGGATTCGGAAAGCGAGCTCGGGT)

Ligation cleanup—

Lambda exonuclease (NEB cat. no. M0262L)
 10x Lambda exonuclease buffer (supplied with NEB cat. no. M0262L)
E. coli exonuclease I (USB cat. no. 70073Z)

Reamplification PCR—

Phusion DNA polymerase (Finnzymes; NEB cat. no. F-530L) CRITICAL – avoid the
 use of polymerases that do not generate blunt ended products.
 5X Phusion HF buffer, containing 7.5mM MgCl₂ (Finnzymes; NEB cat. no. F-518S)

50mM MgCl₂ (Finnzymes; supplied with NEB cat. no. F-518S)

dNTPs (NEB cat. no. N0447L)

MLPAF (Sigma; 5' 6FAM GGGTTCCCTAAGGGTTGGA 3')

MLPAR (Sigma; 5' GTGCCAGCAAGATCCAATCTAGA 3')

Visualizing reaction products—

Agarose (Invitrogen cat. no. 16500100)

10x TBE solution (Invitrogen cat. no. 15581028)

10mg / ml ethidium bromide (Sigma cat. no. E1510) – CAUTION: toxic – wear nitrile gloves when handling, and avoid formation of aerosols.

Loading dye (Qiagen cat. no. 239901)

Low molecular weight DNA ladder (NEB cat. no. N3233L)

8µl formamide (Sigma cat. no. F5786) – CAUTION: toxic – handle with care and avoid breathing vapour.

1µl GeneScan ROX 350 size standard (Applied Biosystems cat. no. 401735).

Equipment

Parafilm M (Sigma cat. no. P7793)

50ml Falcon Tubes (VWR cat. no. 734-0448)

Thermal cycler with 96-well heating blocks (e.g. MJ Research product no. PTC-225)

100ml and 1l measuring cylinders (e.g. Sigma cat. nos. Z324183 and Z650641 respectively)

Thin-walled 200µl PCR tubes (VWR cat. no. 732-0548)

3mm tungsten carbide beads (Qiagen cat. no. 69997).

96-well poly propylene PCR plates (Greiner cat. no. 652290)

Microseal 'A' plate seal (Bio-Rad cat. no. MSA-5001)

Microseal 'F' plate seal (Bio-Rad cat. no. MSA-1001)

Vortex mixer, fitted with a microtitre plate adapter (e.g. VWR cat. no. 444-0486 and VWR cat. no. 444-5919).

Plate centrifuge (e.g. Eppendorf product no. 5804, fitted with A-2-DWP rotor)

2ml polypropylene round-bottomed Cryogenic Vial (Corning cat. no. 430289)
CRITICAL: the shape of vial influences the stability of the resulting emulsion.

8mm × 3mm magnetic bar with pivot ring (VWR International cat. no. 442-4500)
CRITICAL: the size of magnetic stirrer and the presence / absence of a pivot ring influences the stability of the resulting emulsion.

Stirring plate, capable of 1,000 rpm (e.g. IKAMAG Midi MR1; Sigma cat. no. Z403733)

250ml conical flask (e.g. Sigma cat. no. 308811)

Horizontal electrophoresis gel tank (e.g. Sigma cat. no. Z338818)

Electrophoresis power supply (e.g. Sigma cat. no. Z654361)
 UV transilluminator (e.g. Sigma cat. no. Z363839)
 Fluorescence microscope (e.g. Olympus IX71)
 Glass microscope slides (e.g. Cole-Parmer cat. no. WZ-48500-00)
 ABI 3100 capillary electrophoresis sequencer with 36cm capillary array

Procedure

Optimization of PCRs (total time: 7 hours)

1. Optimize PCRs in solution before attempting them in emulsion. For locus 1, this is achieved using primer F1 and F2'-R1 and for locus 2 using F2 and R2 as forward and reverse primers respectively. For each locus, set up 12 identical PCR reactions as follows, adding the polymerase last. Each reaction will be run with a different annealing temperature, in order to find the most suitable one for the fusion reaction. Also include a reaction in which water is used instead of the DNA template, to monitor reagent contamination. Keep reagents on ice, apart from the polymerase, which should be kept at -20°C until it is needed:

5× Phusion HF buffer	10 μl (\rightarrow 1× final conc.)
50mM MgCl ₂	1 μl (\rightarrow 2.5mM TOTAL final conc.)
10mM dNTP mix	1.25 μl (\rightarrow 250 μM each final conc.)
10 μM forward primer	1.5 μl (\rightarrow 300nM final conc.)
10 μM reverse primer	1.5 μl (\rightarrow 300nM final conc.)
50ng / μl genomic DNA (or water)	2 μl (\rightarrow 2ng / μl final conc.)
2 U / μl Phusion polymerase	0.5 μl (\rightarrow 0.02 U / μl final conc.)
water	32.25 μl
Total volume	50 μl per reaction

2. Mix well by vortexing and spin down. Transfer 50 μl of each reaction to a 200 μl thin walled PCR tube and spin down.

3. Transfer tubes to a thermal cycler and perform standard PCR conditions with a gradient for the annealing temperature. The upper and lower range should be set at 5 $^{\circ}\text{C}$ higher and lower than the mean T_m for both primers. For primers designed as described in the Experimental Design section, suitable cycling conditions are:

98 $^{\circ}\text{C}$ for 30 seconds
 33 cycles of 98 $^{\circ}\text{C}$ for 10 seconds
 60-70 $^{\circ}\text{C}$ for 30 seconds
 72 $^{\circ}\text{C}$ for 15 seconds
 then 72 $^{\circ}\text{C}$ for 5 minutes
 4 $^{\circ}\text{C}$ indefinitely

4. Prepare 1l of 1x TBE buffer by adding 900ml deionized water to 100ml 10x TBE in a 1l measuring cylinder, and mixing well.

5. Add 150ml of 1x TBE to 3g of agarose in a 250ml conical flask, and heat in a microwave oven on full power until the liquid just begins to boil.
6. Remove the conical flask from the microwave and allow to cool, mixing frequently, until the flask can be held comfortably in a gloved hand.
7. Add 6 μ l of 10mg / ml ethidium bromide to the gel and mix well. Pour into the gel cassette, with a suitable comb, and allow to solidify. Submerge the gel in 1x TBE, also containing 0.4 μ g/ml ethidium bromide.
8. After PCR, take 5 μ l of each reaction from Step 3, mix with 2 μ l of loading dye in a 200 μ l PCR tube and load each reaction into a separate well of the gel from Step 7. Also include a well containing low molecular weight DNA ladder, for size estimation.
9. Run the gel at 6V / cm, until the yellow dye is close to the bottom of the gel (takes approximately 45 minutes for a 14cm gel).
10. Visualise the gel by UV light on a transilluminator.
11. Identify the highest annealing temperature at which the amplification of both loci produces a strong band, representing a single product of the intended size. Use this temperature as the optimal annealing temperature in subsequent reactions.

TROUBLESHOOTING

12. Take the PCR reactions for each locus that were performed at the optimal annealing temperature, and prepare a 10x dilution with water.
13. Perform a PCR using 1 μ l of both diluted templates from Step 12 to confirm that the fusion product is formed successfully. Prepare the following, in a 200 μ l thin-walled PCR tube:

5x Phusion HF buffer	10 μ l (\rightarrow 1x final conc.)
50mM MgCl ₂	1 μ l (\rightarrow 2.5mM TOTAL final conc.)
10mM dNTP mix	1.25 μ l (\rightarrow 250 μ M each final conc.)
10 μ M F1 primer	1.5 μ l (\rightarrow 300nM final conc.)
10 μ M R2 primer	1.5 μ l (\rightarrow 300nM final conc.)
diluted template 1	1 μ l
diluted template 2	1 μ l
2 U / μ l Phusion polymerase	0.5 μ l (\rightarrow 0.02 U / μ l final conc.)
water	32.25 μ l
Total volume	50 μ l per reaction

Mix well and spin down

14. Transfer to a thermal cycler and begin the PCR cycle described in Step 3, but using the optimal annealing temperature selected in Step 11. Increase the extension time, if necessary, so that it is suitable for amplification of the full length fusion product. As an approximate guide, 20 seconds extension time is sufficient for a 300bp amplicon.
15. Run the PCR product on a 2% agarose gel, as described in Steps 4-10, alongside the locus 1 and 2 amplicons from Step 3, to confirm that a single fusion product of the expected size is produced. TROUBLESHOOTING

HF-PCR

16. Emulsions can be prepared and HF-PCR carried out either singly or in 96-well format. Once the vortexing speed has been optimized, the 96-well format is simpler and costs less per reaction, but this format uses a different oil phase, which makes the post-emulsion cleanup slightly lengthier.

A. Preparation of emulsions and HF-PCR singly (Timing: 1 hour to make emulsions, 1.5 hours to perform PCR)

- i. In a 100ml measuring cylinder, add:
 - 4.5ml Span 80
 - 400 μ l Tween 80
 - 50 μ l Triton X-100
- ii. Make the volume up to 100ml with light mineral oil.
- iii. Cover the top of the measuring cylinder with Parafilm and mix thoroughly by vortexing / shaking.
- iv. Transfer to 2x 50ml Falcon Tubes

PAUSE POINT—This oil phase can be stored indefinitely at room temperature (20°C).

CRITICAL STEP—the components of this oil phase are very viscous, so are difficult to pipette accurately. For the Triton X-100, it is helpful to use a larger volume pipette (e.g. P1000) than would normally be used, as smaller volume pipettes can struggle to draw the desired volume of liquid into the pipette tip. When drawing liquids up, it is also essential to leave the pipette tip in the liquids for several seconds, to ensure that the intended volume has been taken, and to dispense slowly, so that no liquid is left behind. It is also possible to rinse tips out with light mineral oil, collecting all liquid in the 100ml measuring cylinder.

- v. Prepare 100 μ M aqueous phase for each reaction, adding the polymerase last. Keep reagents on ice, apart from the polymerase, which should be kept at -20°C until it is needed:

5 \times Phusion HF buffer	20 μ l (\rightarrow 1 \times final conc.)
50mM MgCl ₂	2 μ l (\rightarrow 2.5mM TOTAL final conc.)
10mM dNTP mix	2.5 μ l (\rightarrow 250 μ M each final conc.)
10 μ M F1 primer	10 μ l (\rightarrow 1 μ M final conc.)
10 μ M R2 primer	10 μ l (\rightarrow 1 μ M final conc.)
1 μ M F2'-R1 primer	1 μ l (\rightarrow 10nM final conc.)
50ng / μ l genomic DNA	4 μ l (\rightarrow 2ng / μ l final conc.)
2 U / μ l Phusion polymerase	8 μ l (\rightarrow 0.16 U / μ l final conc.)
water	42.5 μ l
Total volume	100 μ l
Mix well by vortexing and spin down.	

- vi. Add 200 μ l oil phase to a round bottomed Cryogenic Vial containing an 8mm \times 3mm magnetic bar, and position in the centre of a stirring plate set to 1,000 rpm. This is most easily done by holding the tube in a foam rack, making sure that the bottom of the tube

is in contact with the stirrer. Move the tube until the sound of the spinning magnetic bar is quietest; this indicates that the magnetic bar is in the centre of the tube.

vii. Add the 100 μ l of aqueous phase from step v. dropwise, one drop every five seconds, to the centre of the cryogenic vial, starting a timer as the first drop is added.

viii. Continue stirring for a total of five minutes since the first drop of aqueous phase was added and stop stirring.

ix. This emulsification method creates 300 μ l of emulsion, which is more than can fit into a single 200 μ l PCR tube, so either it can be divided between two tubes, by pipetting 125 μ l into each, creating duplicate reactions, or the surplus can be discarded. Overlay reactions with one drop / 30 μ l of light mineral oil.

x. Transfer to a thermal cycler and begin the following PCR program, using the optimal annealing temperature established in Step 11 (if using an MJ cycler, enter 100 μ l as the reaction volume).

98°C for 30 seconds

33 cycles of 98°C for 10 seconds

e.g. 65°C for 30 seconds

72°C for 30 seconds

then 72°C for 5 minutes

4°C indefinitely

B. Preparation of emulsions and HF-PCR in 96-well format (Timing: 4 hours for initial optimization, 1 hour subsequently, and 1.5 hours to perform PCR)

i. In a 100ml measuring cylinder, add:

40ml Silicone polyether / cyclopentasiloxane

30ml Cyclopentasiloxane / trimethylsiloxysilicate

30ml AR20 silicone oil

ii. Cover the top of the measuring cylinder with Parafilm and mix thoroughly by vortexing / shaking.

iii. Transfer to 2x 50ml Falcon Tubes

CRITICAL STEP

This oil phase may have a cloudy appearance, which will clear over time. Mix thoroughly before use.

PAUSE POINT

This oil phase can be stored for six months at room temperature.

iv. To allow optimisation of the emulsion preparation, prepare mock aqueous phase:

5x HF buffer	5 μ l
15 μ m fluorescent beads	1 μ l

water	19 μ l
Total volume	25 μ l
Mix thoroughly and spin down	

v. Transfer the mock aqueous phase from Step 16Biv to a well in the centre of a 96-well polypropylene PCR plate, along with a 3mm tungsten carbide bead, and add 50 μ l of oil phase from Step 16Biii.

vi. Apply a Microseal 'A' plate seal, pressing firmly, and centrifuge briefly in a plate centrifuge with the plate seal downmost. Holding the 'short spin' button for 10 seconds is sufficient. The purpose is to get the bead and both oil and aqueous phases to the correct part of the well.

vii. Vortex the plate for 150 seconds, inverted, at speed 5 on a Vortex Genie 2, fitted with a microtitre plate adapter.

At the same setting, different vortexers will operate at slightly different speeds, and so it is necessary to optimize the vortexing step to determine the settings required to obtain aqueous compartments that are 5-10 μ M in diameter.

viii. Revert the plate and centrifuge extremely briefly. Again, holding the 'short spin' button for 10 seconds is sufficient. The purpose is to get the bead and emulsions to the bottom of the wells before removing the plate seal. **CRITICAL STEP:** Excessive centrifugation may cause separation of the emulsions.

ix. Using a pipette tip, spread 10 μ l of the emulsion thinly on a clean microscope slide and inspect the emulsion by fluorescence microscopy. Estimate the average diameter of aqueous compartments by comparison with the fluorescent beads. This should be 5-10 μ m (Fig. 2 a-d); if it is not, repeat Steps 16Bv-ix using different vortexing settings: a lower vortexing speed will increase the diameter of aqueous compartments, whereas a faster speed will reduce the diameter. Note the setting on the vortexer at Step 16Bvii necessary to obtain 5-10 μ m average diameter of aqueous compartments, as this setting will be used for all subsequent emulsification. **TROUBLESHOOTING**

x. For HF-PCR, prepare aqueous phase as follows, adding the polymerase last. Keep reagents on ice, apart from the polymerase, which should be kept at -20°C until it is needed:

5 \times Phusion HF buffer	5 μ l (\rightarrow 1 \times final conc.)
50mM MgCl ₂	0.5 μ l (\rightarrow 2.5mM TOTAL final conc.)
2.5mM dNTP mix	2.5 μ l (\rightarrow 250 μ M each final conc.)
10 μ M F1 primer	2.5 μ l (\rightarrow 1 μ M final conc.)
10 μ M R2 primer	2.5 μ l (\rightarrow 1 μ M final conc.)
0.5 μ M F2'-R1 primer	0.5 μ l (\rightarrow 10nM final conc.)
50ng / μ l genomic DNA	1 μ l (\rightarrow 2ng / μ l final conc.)
2 U / μ l Phusion polymerase	2 μ l (\rightarrow 0.16 U / μ l final conc.)
water	9.5 μ l
Total volume	25 μ l
Mix well by vortexing and spin down.	

- x. Add one 3mm tungsten carbide bead to each well of a 96-well polypropylene PCR plate, and then add 50 μ l of oil phase from Step 16Biii and 25 μ l of aqueous phase from Step 16Bx to each well.
- xii. Seal the plate firmly with a Microseal 'A' plate seal, and centrifuge briefly, as described in Step 16Bvi.
- xiii. Vortex the plate for 150 seconds, inverted, at the setting determined by the optimization (Steps 16Bv-ix).
- xiv. Revert the plate and centrifuge extremely briefly as described in Step 16Bviii.
- xv. Transfer to a thermal cycler and begin the following PCR program, using the optimal annealing temperature established in Step 11 (if using an MJ cycler, enter 100 μ l as the reaction volume).
 - 98°C for 30 seconds
 - 33 cycles of 98°C for 10 seconds
 - e.g. 65°C for 30 seconds
 - 72°C for 30 seconds
 - then 72°C for 5 minutes
 - 4°C indefinitely
- xvi. Proceed to Step 17 as soon as possible, ideally within an hour of the PCR reaching 4°C, as emulsions may begin to separate if left for longer.

Disruption of emulsions (timing: 1 hour)

Transfer emulsions to a clean plate / tube using a multichannel pipette, taking care to leave any separated aqueous phase behind.

CRITICAL STEP

For emulsions in plate format, it is helpful to add 50 μ l of 1x Phusion HF PCR buffer to each well before transfer, to increase the volume.

CRITICAL STEP

Some separation of emulsions can occur during PCR. This will be evident as clear liquid at the bottom of wells / tubes. Fusion PCR that occurred in this aqueous phase will not necessarily result in the joining of loci from the same template molecule, and so it is essential to transfer intact emulsion to a clean plate / tube and to discard the separated aqueous phase. It is not necessary to recover the entire emulsion, and it is preferable to discard some than to retain this aqueous phase.

17. In a fume cupboard, add an equal volume of diethyl ether or hexane to emulsions, seal plates firmly with a foil seal and vortex, until a homogenous mixture is achieved.
18. For tubes, centrifuge at 13,000g, and for plates at 3000g, both for 3 minutes at 20°C.
19. In a fume cupboard, remove and discard the upper (solvent) layer using a multichannel pipette.
20. Repeat steps 18. to 20 until the solvent and aqueous phases separate cleanly.
21. Leave the plate / tubes uncovered for 15 minutes in a fume cupboard, to allow remaining ether / hexane to evaporate.

22. Add 0.8 units of Proteinase K to the recovered aqueous phase from Step 22 and incubate at 56°C for 1h to digest the polymerase.

23. Incubate at 95°C for 10 minutes to denature the Proteinase K.

24. Make the volume up to 200µl with water.

PAUSE POINT

reactions can be stored at –20°C for several weeks.

Preparation of ligation oligonucleotides (timing: 1 hour)

26. Phosphorylate unblocked 5' termini in a 200µl thin-walled PCR tube. Add:

10µM IRF5_5T	5µl (→ 1µM final conc.)
10µM IRF5_5C	5µl (→ 1µM final conc.)
10µM IRF5_FusF4	10µl (→ 2µM final conc.)
10× T4 ligase buffer	5µl (→ 1× final conc.)
10U / µl T4 polynucleotide kinase	1µl (→ 0.2U / µl final conc.)
water	24µl
Total volume	50µl

Mix thoroughly and spin down

27. Incubate for 30 minutes at 37°C

28. Denature the kinase by heating to 65°C for 20 minutes

PAUSE POINT: phosphorylated oligos can be stored in 10µl aliquots at –20°C for several weeks.

29. Add 5µl of 10µM primers IRF5_3G and IRF5_3T, along with water to a total volume of 100µl.

Ligation Haplotyping Reaction (total time 5 hours: ligation haplotyping reaction - 2 hours; digestion of surplus oligos - 1.5 hours; amplification of ligated products - 1.5 hours)

30. Set ligation reactions up in 200µl thin-walled PCR tubes as follows:

2µl primer mix	0.5µl
50 U / µl <i>Tth</i> ligase	0.4µl (→ 0.02 U / µl final conc.)
water	15.1µl
Total volume	20µl

Mix by vortexing and spin down

31. Transfer tubes to a thermal cycler and begin following ligation program: 95°C for 120 seconds

20 cycles of 95°C for 30 seconds

64°C for 240 seconds

then 4°C indefinitely

32. Add 2µl of 0.8 U / µl Proteinase K to each reaction, mix and spin down, and incubate as described in steps 23 and 24.

33. To digest surplus oligos, add:

10× Lambda exonuclease buffer	5µl (→ 1× final conc.)
5 U / µl Lambda exonuclease	0.2µl (→ 0.02 U / µl final conc.)
10 U / µl <i>E. coli</i> exo I	0.05µl (→ 0.001 U / µl final conc.)
water	22.75µl
Total volume	50µl
Mix thoroughly and spin down	

34. Incubate tubes at 37°C for an hour, followed by 65°C for 20 minutes to denature the exonucleases.

35. To amplify ligated oligos and to incorporate a fluorescent label, perform a final PCR reaction, in a 200µl thin-walled PCR tube. Use a proofreading DNA polymerase to avoid the addition of adenosine to the 3' end of PCR products, to give single peaks after capillary electrophoresis:

5× Phusion HF	10µl (→ 1× final conc.)
2.5mM dNTP mix	4µl (→ 200µM each final conc.)
10µM primer MLP AF	1.5µl (→ 300µM final conc.)
10µM primer MLP AR	1.5µl (→ 300µM final conc.)
Ligation product	2.5µl
2 U / µl Phusion polymerase	0.5µl (→ 0.02 U / µl final conc.)
water	20µl
Total volume	50µl
Mix thoroughly and spin down	

36. Transfer tubes to a thermal cycler and begin the PCR program:

98°C for 30 seconds

33 cycles of 98°C for 10 seconds

63°C for 30 seconds

72°C for 30 seconds

then 72°C for 5 minutes

4°C indefinitely

PAUSE POINT

Reactions can be stored at -20°C for several weeks.

Visualizing reaction products (timing: 1 hour)

37. Mix 1 μ l of PCR product from Step 36 with 8 μ l formamide and 1 μ l GeneScan ROX 350 size standard in a 96-well polypropylene PCR plate and spin down.

38. Heat to 96°C for 3 minutes and transfer immediately to ice

39. Spin plate down again, load into capillary sequencer and run, following the manufacturer's instructions.

Timing

1 Optimization of PCRs	7 hours
16 HF-PCR	
A. Preparation of emulsions singly	1 hour to make emulsions, 1.5 hours to perform PCR
B. Preparation of emulsions in 96-well format	4 hours for initial optimization 1 hour subsequently, and 1.5 hours to perform PCR
17 Disruption of emulsions	1 hour
26 Preparation of ligation oligonucleotides	1 hour
30 Ligation Haplotyping reaction	2 hours
33 digestion of surplus oligos	1.5 hours
35 amplification of ligated products	1.5 hours
37 Visualizing reaction products	1 hour

Step number	Problem	Possible reason	Solution
11	More than one amplicon per PCR	Primers insufficiently specific	Redesign primers and confirm specificity as described in steps 1-10
11	No amplicons produced by PCR	PCR master mix prepared incorrectly, or reagents have deteriorated	Repeat PCR with fresh reagents step 1
11		Temperature gradient range is too high	Repeat gradient PCRs using a lower range of annealing temperatures – step 3
11	No annealing temperature that works for both loci	Annealing temperature of primer pairs does not match	Adjust length of primers to give more similar melting temperatures step 1
15	No fusion PCR product	Incorrectly designed fusion primer	Check sequence of fusion primer
16 B ix	Aqueous compartments too large	Vortexing speed too low	Repeat optimization with higher vortexing speed
16 B ix	Aqueous compartments too small	Vortexing speed too high	Repeat optimization with lower vortexing speed
16 B x	Emulsion is not stable	Incorrect volumes of reagents used in oil	The components of the oil phase are very viscous. Make fresh oil phase,

Step number	Problem	Possible reason	Solution
	before PCR	phase	taking care to pipette the components slowly and accurately
16 B xii	Leakage of emulsion from plate	Plate sealed with insufficient pressure	Make sure that plate seal is Microfilm 'A' and that it is pressed into wells firmly
16 B xvi	Emulsion separates during PCR	Incorrect volumes of aqueous and oil phases used	The oil phase is very viscous. Repeat plate, taking care to pipette the oil phase slowly and accurately
16 B xvi		Incorrect volumes of reagents used in oil phase	See above
21	Some traces of emulsion remain after solvent steps	Inadequate mixing of solvent and emulsion	Repeat solvent washes, taking care to mix very thoroughly. However, it is not essential to obtain a completely clean aqueous phase
39	More than 2 peaks observed per individual	Aqueous phase recovered along with emulsion	Repeat HF-PCR, taking care only to recover intact emulsion
39	No peaks visible	Sequencing run failed	Repeat sequencing run
39	Only size standard peaks visible	No ligation in haplotyping reaction due to failure of phosphorylation	Repeat oligonucleotide preparation using fresh T4 ligase buffer. Make sure that all solids have dissolved before use.
39		No ligation in haplotyping reaction due to failure of ligase	Repeat oligonucleotide preparation using fresh T4 ligase buffer. Make sure that all solids have dissolved before use. Warm buffer to 37°C and vortex thoroughly. Inspect visually
39	Peaks too large (i.e. go off the Y axis scale on visualisation software)	Too much DNA used in sequencing reaction	Repeat sequencing with less DNA
39	Peaks differ in size by > 10x	Suboptimal concentration of oligos used in ligation reaction	Double concentration of oligos that generate small peaks, halve concentration of oligos that generate large peaks

Anticipated results

The allele-specific oligos for each SNP at each position differ in length, as described in the 'Experimental Design - Ligations' section, differ in length. In the example given here, the two oligonucleotides at SNP 1 were designed to differ by 5 bases, and the two oligonucleotides at SNP 2 were designed to differ by 10 bases. After the fusion ligation reaction, products have a different size, depending upon the two SNPs present in the condensed haplotype template.

Ligation products are amplified by PCR, using one fluorescently labelled primer, and are resolved by capillary electrophoresis. This generates a peak, the size of which indicates the SNPs present. For autosomal haplotypes, an individual who is homozygous at the

interrogated positions generates a single peak, whereas heterozygous individuals generate two peaks, representing the haplotypes on each parental chromosome. (Fig. 1).

Different allele-specific oligos will ligate with different efficiencies during the ligation haplotyping reaction, and this will give rise to different (Fig. 3). The concentration of allele-specific oligos given in steps 28 and 29 can be optimised to produce more uniform peak heights: the concentration of bridging oligo is kept constant, and for low peak heights the concentration of the relevant allele-specific oligo is doubled, whereas for large peak heights, the concentration is halved. This is usually sufficient, though the concentration can be halved or doubled once more, if necessary.

For the example given (Fig. 3), the 114bp GC peak is much larger than the 124bp GT peak, indicating that the concentration of the 5T oligo should be increased. In practice, it is not necessary to have identical peak heights to determine haplotypes unambiguously. It should be noted that the fragment size estimated by the sequencing software may not correspond exactly to the true size. It is helpful to perform the Ligation Haplotyping reaction on fusion PCR products that were generated in solution rather than emulsion, and to run the Ligation Haplotyping reaction products on the capillary sequencer. This will generate peaks of all four sizes, and will reveal any inaccuracy in the size calls.

***IRF5* SNPs**

The expression levels of several unique *IRF5* isoforms is controlled by two SNPs: presence of the T allele of SNP rs2004640, in exon 1 of *IRF5*, permits expression of these isoforms by creating a splice site that is absent from the G allele. The T allele of the second SNP, rs2280714 is associated with higher expression levels. The effect is acts in *cis*, increasing the likelihood that individuals with the T SNP at both positions on the same parental chromosome will develop systemic lupus erythematosus (SLE) (Fig. 4a) 4,24.

To test the method, we obtained statistically derived haplotypes from the 30 Centre d'Etude du Polymorphisme Humain (CEPH) HapMap trios 4,24, and compared these to the results obtained from Ligation Haplotyping. Both sets of results were in agreement in every case 15. Additionally, our results showed haplotype inheritance to be Mendelian, and the frequencies of diplotypes obeyed the Hardy-Weinberg equilibrium, as expected (Fig. 4b). Not all of the possible SNP combinations were observed in these samples: none of the individuals had the TC haplotype, and because of the low frequency (5.6%) of the GT haplotype, there were no GTGT homozygotes (Figs. 4b and 4c).

Acknowledgments

This work was funded by the Wellcome Trust [grant number 077014/Z/05/Z]. The authors would like to thank Chris Tyler-Smith for his work on the primary Ligation Haplotyping paper, and for his invaluable comments on this manuscript. We would also like to thank Robert Graham for statistically derived *IRF5* haplotype data and Oxford Journals for permission to reuse figures from the primary research paper.

References

1. Sabeti PC, et al. Detecting recent positive selection in the human genome from haplotype structure. *Nature*. 2002; 419:832–837. [PubMed: 12397357]
2. Beerli P, Felsenstein J. Maximum likelihood estimation of a migration matrix and effective population sizes in n subpopulations by using a coalescent approach. *Proc Natl Acad Sci U S A*. 2001; 98:4563–4568. [PubMed: 11287657]
3. Myers S, Bottolo L, Freeman C, McVean G, Donnelly P. A fine-scale map of recombination rates and hotspots across the human genome. *Science*. 2005; 310:321–324. [PubMed: 16224025]

4. Sigurdsson S, et al. Polymorphisms in the tyrosine kinase 2 and interferon regulatory factor 5 genes are associated with systemic lupus erythematosus. *Am J Hum Genet.* 2005; 76:528–537. [PubMed: 15657875]
5. Kleinjan DA, van Heyningen V. Long-range control of gene expression: emerging mechanisms and disruption in disease. *Am J Hum Genet.* 2005; 76:8–32. [PubMed: 15549674]
6. Salem RM, Wessel J, Schork NJ. A comprehensive literature review of haplotyping software and methods for use with unrelated individuals. *Hum Genomics.* 2005; 2:39–66. [PubMed: 15814067]
7. Marchini J, et al. A comparison of phasing algorithms for trios and unrelated individuals. *Am J Hum Genet.* 2006; 78:437–450. [PubMed: 16465620]
8. Consortium TIH. A haplotype map of the human genome. *Nature.* 2005; 437:1299–1320. [PubMed: 16255080]
9. Lin S, Cutler DJ, Zwick ME, Chakravarti A. Haplotype inference in random population samples. *Am J Hum Genet.* 2002; 71:1129–1137. [PubMed: 12386835]
10. Lo YM, et al. Direct haplotype determination by double ARMS: specificity, sensitivity and genetic applications. *Nucleic Acids Res.* 1991; 19:3561–3567. [PubMed: 1852606]
11. McDonald OG, Krynetski EY, Evans WE. Molecular haplotyping of genomic DNA for multiple single-nucleotide polymorphisms located kilobases apart using long-range polymerase chain reaction and intramolecular ligation. *Pharmacogenetics.* 2002; 12:93–99. [PubMed: 11875363]
12. Michalatos-Beloin S, Tishkoff SA, Bentley KL, Kidd KK, Ruano G. Molecular haplotyping of genetic markers 10 kb apart by allele-specific long-range PCR. *Nucleic Acids Res.* 1996; 24:4841–4843. [PubMed: 8972876]
13. Paabo S, Irwin DM, Wilson AC. DNA damage promotes jumping between templates during enzymatic amplification. *J Biol Chem.* 1990; 265:4718–4721. [PubMed: 2307682]
14. Turner DJ, et al. Assaying chromosomal inversions by single-molecule haplotyping. *Nat Methods.* 2006; 3:439–445. [PubMed: 16721377]
15. Turner DJ, Tyler-Smith C, Hurles ME. Long-range, high-throughput haplotype determination via haplotype-fusion PCR and ligation haplotyping. *Nucleic Acids Res.* 2008; 36:e82. [PubMed: 18562465]
16. Yon J, Fried M. Precise gene fusion by PCR. *Nucleic Acids Res.* 1989; 17:4895. [PubMed: 2748349]
17. Dressman D, Yan H, Traverso G, Kinzler KW, Vogelstein B. Transforming single DNA molecules into fluorescent magnetic particles for detection and enumeration of genetic variations. *Proc Natl Acad Sci U S A.* 2003; 100:8817–8822. [PubMed: 12857956]
18. Tawfik DS, Griffiths AD. Man-made cell-like compartments for molecular evolution. *Nat Biotechnol.* 1998; 16:652–656. [PubMed: 9661199]
19. Margulies M, et al. Genome sequencing in microfabricated high-density picolitre reactors. *Nature.* 2005; 437:376–380. [PubMed: 16056220]
20. Barany F. Genetic disease detection and DNA amplification using cloned thermostable ligase. *Proc Natl Acad Sci U S A.* 1991; 88:189–193. [PubMed: 1986365]
21. Fan JB, et al. Highly parallel SNP genotyping. *Cold Spring Harb Symp Quant Biol.* 2003; 68:69–78. [PubMed: 15338605]
22. Tobler AR, et al. The SNPlex genotyping system: a flexible and scalable platform for SNP genotyping. *J Biomol Tech.* 2005; 16:398–406. [PubMed: 16522862]
23. Ghadessy FJ, Ong JL, Holliger P. Directed evolution of polymerase function by compartmentalized self-replication. *Proc Natl Acad Sci U S A.* 2001; 98:4552–4557. [PubMed: 11274352]
24. Graham RR, et al. A common haplotype of interferon regulatory factor 5 (IRF5) regulates splicing and expression and is associated with increased risk of systemic lupus erythematosus. *Nat Genet.* 2006; 38:550–555. [PubMed: 16642019]





Figure 1. Ligation haplotyping

a) Haplotype fusion PCR. The two SNPs to be interrogated can lie several kilobases apart. Primer pairs F1 and R1, and F2 and R2 are designed for the separate loci, as for a typical PCR. The reverse complementary sequence of F2 is appended to the 5' terminus of primer R1, creating the 'fusion primer' F2'-R1. Three of these primers are used in the emulsion PCR reaction: F1, R2, and the fusion primer, F2'-R1, shown as black and black / indigo arrows. The fusion primer is used at a lower concentration than the other primers, so that in the initial rounds of amplification, where locus 1 is amplified exponentially, the fusion primer quickly becomes depleted. At the same time, locus 2 is amplified linearly, by extension of the R2 primer. Because the tail on the fusion primer matches the reverse complement of F2, the amplicon generated by this primer can prime on locus 2. Primers F1 and R2 can then amplify the fused template exponentially for the remainder of the reaction.

b) Ligation Haplotyping. During the ligation reaction, the bridging oligo (in this case, the same oligo as the fusion primer) anneals to the fused template, spanning the junction. Allele-specific oligos, with characteristic sizes, anneal to the positions at either end of the bridging oligo, and when the oligo matches the template, a ligation reaction occurs. For each haplotype, a differently sized product is generated. The allele-specific oligos are tailed with universal primer sequences (Tag 1 and Tag 2), so that all ligation products can be amplified together, with as little bias as possible.

c) Ligated products are amplified by PCR, using a fluorescently labelled primer, which anneals to the Tag 1 sequence, and these fluorescent amplicons are analyzed by capillary electrophoresis.

Figure adapted with permission from Oxford Journals, NAR 2008

**Figure 2. Preparation of emulsions**

We prepared emulsions by vortexing and stirring as described, but we also included 10µm fluorescent microspheres in the aqueous phase, to allow us to estimate the diameter of aqueous compartments (green). Aqueous compartments prepared by vortexing were typically 5-10µm (panel a), whereas those present in stirred emulsions were larger, approximately 15µm in diameter (panel b). Regardless of the method of emulsion preparation, the diameter of aqueous compartments tends to decrease slightly during the PCR reaction (see panel a compared to c, and b compared to d).

Figure adapted with permission from Oxford Journals, NAR 2008




Figure 3. Optimization of Ligation Haplotyping reactions

Peak heights may differ considerably after the Ligation Haplotyping reaction. To generate peaks with a more even size, the oligo concentration in the ligation reaction should be optimized. The 114bp GC peak is much larger than the 124bp GT peak – the 3G oligo is common to both peaks, indicating that the concentration of the 5T oligo should be increased. Figure adapted with permission from Oxford Journals, NAR 2008




Figure 4. Ligation haplotyping on genomic DNA

a) The G allele of rs2004640 creates *IRF5* splice variants that are initiated at exon 1A. The T allele causes expression of isoforms from exon 1B. At the second locus, rs2280714, the T allele is associated with elevated gene expression compared to the C allele, and the effect acts in *cis*. Consequently, the T-T haplotype is associated with a heightened risk of developing systemic lupus erythematosus.

b) Haplotype inheritance was observed to be Mendelian, and the frequencies of diplotypes obeyed the Hardy-Weinberg equilibrium. None of the individuals had the TC haplotype. The GT haplotype was present at a low frequency, and no GTGT homozygotes were observed.

c) Five diplotypes were observed in the CEPH individuals used in this study, and examples of each are shown. *IRF5* is on Chromosome 7, so individuals have two haplotypes, and heterozygotes will generate two peaks.

Figure adapted with permission from Oxford Journals, NAR 2008