

NIH Public Access

Author Manuscript

Vet Pathol. Author manuscript; available in PMC 2010 November 1.

Published in final edited form as:

Vet Pathol. 2009 November ; 46(6): 1230-1240. doi:10.1354/vp.09-VP-0002-S-FL.

A Data Capture Tool for Mouse Pathology Phenotyping

Beth A. Sundberg¹, Paul N. Schofield², Michael Gruenberger², and John P. Sundberg¹

¹ The Jackson Laboratory, Bar Harbor, Maine, USA

² Department of Physiology, Development, and Neuroscience, University of Cambridge, Cambridge, UK

Abstract

The Mouse Disease Information System (MoDIS) is a free Microsoft Access database (http://research.jax.org/faculty/sundberg/index.html) designed by veterinary pathologists to aid veterinary pathologists in data acquisition, analysis, and coordination of tissue sample archives. Linking the system to the Mouse Anatomy and Mouse Pathology Ontologies provides controlled vocabulary (and spelling) for organ, tissue, and diagnosis. Severity scores provide a quantitative assessment of all lesions to enable quantitative trait locus analysis for large scale studies. Individual diagnoses can be verified for their definition by online linkage to Pathbase.net. Histological images can be accessed from Pathbase using the Mouse Pathology Ontology (MPATH) directly for comparison to slides being viewed at the time of data entry providing the user with a reference and a "virtual second opinion."

Keywords

Access database; medical records; pathology training; Pathbase; Mouse Tumor Biology Database; Mouse Anatomy Ontology; Mouse Pathology Ontology

Introduction

Maintenance of records on mouse breeding colonies and related research is highly specialized and can easily become overwhelming if a systematic and proven method is not utilized. A number of very sophisticated systems are available today for this purpose (Jackson Laboratory colony management system, http://colonymanagement.jax.org/index.html; MausDB, http://www.helmholtz-muenchen.de/en/ieg/downloads/index.html; and MusDB¹¹). Likewise, sophisticated database systems are available for physiologic phenotyping data (blood chemistry, body weights, body composition, etc.). The same is true for medical records on mouse necropsy and diagnosis, both to organize all data generated and to access various specimens collected and stored. Generation of data and case materials is valueless over the long term unless they can be retrieved and analyzed on demand. Computerized databases and spreadsheets provide invaluable tools for maintaining large amounts of data that can be processed easily to generate tables and graphics for summary, interpretation, and presentation. ^{14,15} More importantly, these data can be downloaded into spreadsheets for quantitative trait locus (QTL) genetic mapping with the appropriate software. ^{1,2,4,6,8,9,13,16}

Today's researcher must keep extensive records on mouse colony, experimental test results, medical records, specimens, digital images, and other information. Many different methods

Corresponding author: John P. Sundberg, D.V.M., Ph.D., Dipl. A.C.V.P., The Jackson Laboratory, 600 Main Street, Bar Harbor, ME 04609-1500 USA, john.sundberg@jax.org, Phone: 207-288-6410, FAX: 207-288-6073/6078.

are often used to obtain and maintain this information, including paper records on forms in notebooks, computerized spreadsheets, computerized databases, automated laboratory equipment printouts, and specially designed computerized laboratory information management systems (LIMS). Additional information may then be created in the form of reports, statistics, graphs, posters, formal presentations, drafts of papers, and finally as publications. This information is frequently scattered through many locations and in the possession of various people involved in the project. Sometimes this proliferation of information can create great difficulties in retrieving specific data if a key person is no longer present. Hence, it is highly desirable that the researcher establish a record-keeping system at the beginning of the project that provides a consistent framework for use by all participants. Ancillary materials, such as glass microscope slides, digital images, paraffin blocks, etc. can then be filed under a common identifier (case/accession number) to facilitate retrieval.

The project record-keeping system will depend on the resources available to the individual researcher. Ideally, a single record-keeping system should be used by everyone in the same research laboratory so that materials from different individuals' projects can be easily identified, sorted, and retrieved, even after individuals leave the research group. Assessment of the available facilities and establishment of protocols for information recording at the beginning of the project (or even when establishing a research laboratory) will minimize problems later. Large automated (computerized) LIMS may already be in place and need to be used. For many researchers, however, these systems are not generally available, may only be available at a prohibitive cost, or cannot be adapted to a particular investigator's needs.

We describe here the use of a freely-available data capture tool built on a commonly used database software program (Microsoft Access) that focuses on collection of data by and for a pathologist. Referenced coding for organs and disease processes as well as summative diagnosis makes the system rapid and generates a report that can be printed and signed to yield a regulatory report or it can be emailed directly to the client/collaborator. Data can be downloaded in a structured form and sent to other databases for integration with genetic and phenotype data.

The overall strategy and potential use-cases of the system that we present here, called the Mouse Disease Information System (MoDIS), were previously described, ¹⁸ but in this complementary paper we provide detailed documentation on the structure and use of MoDIS for the target user group, predominantly veterinary pathologists.

Materials and Methods

Software

The database described here evolved over 20 years and was originally built on dBASE III Plus (Ashton-Tate, Torrance, CA)¹⁴ then migrated onto FoxPro 2.6 (Microsoft Corp., Redmond, WA)¹⁵ to its current form in Microsoft Access.¹⁸ Coding is built on the Anatomical Dictionary for the Adult Mouse (http://www.informatics.jax.org/searches/AMA_form.shtml; The Jackson Laboratory, Bar Harbor, ME) and mouse pathology ontology, MPATH, (http://www.pathbase.net; Pathbase European Consortium, University of Cambridge, Cambridge, UK).

Hardware

The Microsoft Access software is designed for computers using the Microsoft Windows operating system (Windows 2000 or XP, or VISTA and Microsoft Access XP 2003 or 2007). Other operating systems and hardware may be used with third party middleware solutions such as virtualization software.

The Mouse Disease Information System (MoDIS)

Standardization

Unique Identifiers: Large-scale data collection operations, such as diagnostic laboratories, clinics in hospitals, or drug safety studies in industry, assign a specific identifier for each individual animal. This is often referred to as a *case* or *accession* number and is utilized to organize all the materials generated by the case. In a relational database, this unique identifying number is used to tie all the database information together for one animal. Within the database, fields can be set up to record other identifiers for the same animal from other sources that use different numbers and bind them together. In this way all data can be collated unambiguously.

MoDIS uses a two-part unique case identifier. The beginning portion indicates the year of the necropsy. Accessions are then numbered numerically, or groups of numbers can be assigned to particular projects, programs, etc. For example, we use five digit codes. Final assignment of the number has the form: 99-1, 00-7123, 08-12555, etc.

Assignment of numbers to animals is usually done in one of two manners. Since this laboratory runs multiple types of research and diagnostic programs, most of which focus on necropsies of mice from a variety of sources, we assign numbers to the mice when they come to necropsy. A lab that imports mice or uses only mice from within their own colony may find that it is easier to assign the number to the animal when it is born or purchased. In fact, some of the colony management programs assign numbers to mice when they are weaned, which are then used throughout the life of the mouse until it is necropsied. This is particularly true when using some of the computerized software for colony breeding management. If a number is assigned when the mouse enters the colony/research laboratory, that number should be maintained for all data collected.

Complications in record keeping arise when specimens are shipped to collaborators or other services within the institution that have their own accession system. For example, specimens collected for ultrastructural studies are assigned a unique identifier and then sent to the Electron Microscopy (EM) Service for processing. The EM Service assigns its own accession number. When the investigator arrives to review the material with the electron microscopist, confusion is inevitable unless the numbers and data are cross referenced. Solutions include adding the service's accession number to a field in the database for maintaining a cross reference or submitting the investigator's unique identifier and signalment for the case (animal type [mouse, rat], age, and sex), as well as information on the study to the service prior to examination. This is the type of data flow problem that should be noted and solved before the project begins.

Computer programs can be used to generate labels for transparencies, folders, images, negative sleeves, etc. Using the case accession number and project identifier on the labels helps with filing and permits rapid identification if more data need to be looked up.

Project identifiers: A project name or acronym is defined and entered for each mouse record within that project. This project name allows searching through the database or output spreadsheet to find all mice or other data collected for that particular project. If the project involves collaborators or another investigator, separate fields can be set up under *submitter/collaborator* to permit rapid searches for all work that involves this individual or group. Often multiple studies are ongoing with a particular investigator and summaries of all work are needed. This approach enables the program to set up the summaries quickly and, most importantly, reproducibly.

<u>Collection of information from the colony:</u> It may be possible to bring computers into the animal rooms to collect information directly. The use of bar codes and readers, balances tied

to institutional networks to automatically record weights, etc. have become commonplace in mouse rooms. Specialized programs are becoming available for maintaining breeding records this way. If these data acquisition systems are not available, some type of paper form must be used for data collection. These data can then be entered later.

In this laboratory individual mouse information is normally collected into the database at the point the mouse is necropsied. An example of the data collection format is shown in Figure 1. The following signalment data are used:

<u>Case/accession number (necropsy number)</u>: A unique identifier is linked to an individual animal.

<u>Project identifier (clinical number)</u>: This defines the specific project so all cases related to that project can be downloaded for analysis.

<u>Submitter/collaborator (from)</u>: This identifies the primary research contact associated with this mouse (principal investigator, collaborator, other investigator, or student).

<u>Location (building and room)</u>: Where the mouse was housed (animal room number) can be critical for finding mice with similar problems or if an untoward response was identified that might have an environmental basis.

<u>Strain</u>: Strain identification abbreviations/symbols are very specific and will determine correct controls and what spontaneous diseases are common in that strain. Users should strictly adhere to the International Committee on Standardized Genetic Nomenclature for Mice standards of assigned strain and mutation symbols (http://www.informatics.jax.org/mgihome/nomen/index.shtml).¹⁷

<u>Stock No.</u>: The Jackson Laboratory assigns a production stock number associated with the strain or mutant stock. This serves as an internal control for the correct strain and mutant gene symbols.

<u>Genotype</u>: This would identify the specific gene and allelic mutation under investigation as well as differentiate mutant from control mice.

Line: This refers to the laboratory code used to record a designator for a breeding line.

<u>Pedigree numbers</u>: This may be assigned by the colony management program and serve as their accession number or it may be a unique number assigned by the person maintaining the animals. Often pedigrees are used to track the previous generation so there are fields for the sire and dam pedigree numbers.

<u>Other identifiers (other ID)</u>: This refers to the chip number (These are tracking tools implanted into individual mice⁵), cage number, special coding, ear tags, etc. used to collate a particular individual with other identifiers that may have been used previously.

Sex: F (female), M (male), FS (neutered [spayed] female), MC (castrated male)

<u>Birth date</u>: This information is used to calculate the age when subtracted from date of necropsy.

<u>Death/necropsy date (submission date)</u>: As stated above this is used to calculate age when the mouse died.

<u>Mating</u>: This is the type of mating scheme used: paired (male and female), trio (male and two females).

<u>Code, Source, and Category</u>: Fields that can be customized to further describe the case. Included are categories such as "breeder" or "pup" and sources such as "importation" or

"research". Codes may be created to indicate mutant versus control, study groups, or any other desired information.

<u>Construct name</u>: This would refer to the specific information on the genetic construct used to make a genetically engineered mouse.

Identification of other material associated with each mouse (cross-reference

information): Fields are set up to record in a yes/no fashion if the following material has been collected for each mouse: (gross) photographs, photomicrographs, electron micrographs, histology slides, frozen tissue, radiographs, *in situ* hybridization, immunohistochemistry, gene mapping data, etc. These items are physically stored or filed by project and case number in order to be easily accessible.

Import signalment from a spreadsheet: A simple import system exists for creating a new case and collecting some signalment data from a fixed format spreadsheet (Figure 2). This allows importation of mice at necropsy that have data in a different database. Figure 2 shows a sample of the fields imported.

<u>Collection of Necropsy Information</u>: Necropsy information is added in three phases, 1) directly after initial necropsy, 2) as test results are received, and 3) when the final diagnosis is made and the report finalized. The initial history, clinical signs, and gross description are entered as a searchable text field of variable length. Later on, comments about the final diagnosis are also entered as a searchable text field. The pathologist may enter as many different diagnoses (D_x) as needed to describe his/her findings (Figure. 3). Each diagnosis consists of several fields:

 $\underline{D_x Order}$: This field is used to control the order the diagnosis are listed on the printed report.

<u>Organ or Tissue</u>: This specifies the exact anatomical organ or structure within an organ that is examined.

 \underline{D}_x : One diagnosis associated with the organ or tissue examined.

<u>Disease name</u>: The disease or syndrome associated with the findings is listed here. This may be a well-defined disease or a morphologic diagnosis.

<u>Severity score</u>: Values range from 0 (no lesions) to 4 (extreme severity). This is reported as an adjective in the printed report (mild, moderate, severe, or extreme) but the numerical scores can be downloaded for quantitative trait analyses for haplotype association gene mapping or other specific uses.

The organ, diagnosis, and disease fields may be entered in a free format. However, it is recommended that the ontologies associated with each field are used to create consistency (accurate spelling, fixed vocabulary, and coding) and later numerical correlation of results.

The Anatomical Dictionary for the Adult Mouse⁷ is used to provide term choices for the organ field. Each organ and tissue has a unique code associated with it (for example: Molar MA: 00000350). A link is provided on the data-entry screen to the full hierarchical anatomy listing to help new users locate the correct terms. Once a term has become familiar, it is easily entered by typing the beginning of the term and pressing tab or enter.

The MPATH mouse pathology ontology¹² is used to provide term choices for the diagnosis (D_x) field. Each term has a unique code associated with it (for example: MPATH 445: Papilloma). A link is provided on the data entry screen to the full hierarchical pathology term listing to help new users locate the correct terms. This also provides a consensus definition for

the term by a panel of pathologists (Joint European Experimental Pathology Expert Group, JEEPEG). Once a term has become familiar, it is easily entered by typing the beginning of the term until a unique letter combination is entered at which point the full term and code appear automatically. For example, "bi" entry autofills to bile deposition MPATH: 35.

The *disease name* terms are user-generated and reflect the area or tradition of pathology terms familiar to the user. These can be summative diagnoses. Each new term is given a unique number (Figure. 4). The continuous expert curation of MPATH by a committee of veterinary and physician pathologists (JEEPEG) means that requests for the addition of new terms, synonyms, and definitions can be considered and, if approved, will be integrated into MPATH in an appropriate location in the ontology.

Images: Pathologists often obtain multiple images (photomicrographs and gross photographs) that may be stored in digital formats. It is inefficient to store these as part of a database such as MoDIS due to the large file sizes. Instead, MoDIS provides a method of recording hyperlinks to the location of digital files that are associated with specific diagnoses and cases.

Using terms for teaching and reference: Links are provided on the edit screen to the online Pathbase website (http://www.pathbase.net; Figure. 5). By following the links, examples of the MPATH terms, definitions, and representative annotated photomicrographs may be obtained (Figure. 6). This provides the pathologist or investigator verification of their interpretation and examples of similar if not identical lesions in other mutant mice, effectively providing a "virtual second opinion".

<u>Collection of Special Test results:</u> Any tests that may have been run, such as microbiological tests on selected tissues, are also recorded for each mouse along with the results. Results may be recorded as the number of negative or positive cultures identified or by entering a numeric value (Figure. 7). A separate field indicates the units for the value. A detected organism (parasite, bacteria, etc.) may be recorded. Space has been allocated to maintain a short description, serology number, and comment on each test if necessary. Standardized tables are set up with simple codes for the commonly used tests.

Report Generation: This database is used to output several standardized reports on the individual mice and projects as needed. The database is easily accessed by searching on a range of criteria. More importantly, printed final reports (Figure. 8) can easily be generated. If these reports are signed and dated by the pathologist of record, they become legal documents. It is always a concern that medical records can be subpoenaed or that research material may become important in patent processing. A signed, dated, and filed hard copy attached to the original hand written worksheet provides an unalterable reference. This can become important in legal situations but also provides a record that can be checked for accuracy if a particular data point stands out as apparently being entered incorrectly. For example, when doing retrospective case series studies of frequency and epidemiology of particular types of cancer affecting inbred strains of mice, exceptionally young mice may show up in the tabulation of case information. Checking the original case sheet may indicate that the wrong year was entered for the birth or necropsy date. For mice with an average life expectancy in a research or breeding colony of about a year, rarely two or more, this alters the data significantly.¹⁹ Maintenance of the original record allows for verification of accuracy of data entry, assuming the original worksheet was carefully prepared or the information can be followed back from there to obtain the correct result.

<u>Study Analysis:</u> Many scientists and technicians are familiar with the use of spreadsheets as a laboratory analysis tool. Our database system is used to output a Microsoft Excel spreadsheet of selected data for analysis. All database fields can be searched for and selected for output

(Figure 9). This spreadsheet is then used for additional data collection of project specific data that do not need to be kept in the larger historical database. Analysis tools in the spreadsheet package include summarization, statistics, and graphs. Final copies of tables and graphs can be created that are of publication quality. Graphs and tables can be copied to word processing software for inclusion in papers or to presentation software such as Microsoft PowerPoint, Adobe Illustrator, or Indesign. Spreadsheets can be used as input to statistical programs (JMP[®] Statistical Discovery Software http://www.jmp.com/) and gene mapping programs for further analysis.

Documentation: An extensive user guide is included, which contains instructions, examples with illustrations, and suggestions. The technical documentation provides step-by-step installation instructions, details about multiple users, how to use various operating systems, and problem-solving tips.

Data output: Apart from output to the documentary pathology report (see above), structured data can also be output as character (or comma) separated values (CSV) standard file format for upload into other databases. The formal constraints resulting from the use of standard ontologies mean that coded diagnoses are unambiguous and computable within any environment using the same standard description frameworks (MA and MPATH). Future developments of MoDIS will enable XML output into a defined schema for pathologic diagnosis allowing data transmission and upload into larger databases and integration with other phenotype and genotype data. Extensive documentation is provided online which details how this can be done.

Multiple Users, Security, and Data Tracking: The Access software can be used individually or with a server for multiple simultaneous users, which is how we use it. It is also possible, as long as the server is on the institutional network, to do data entry and manipulation off site. Security and data manipulation abilities are controlled by user passwords. User permissions are controlled by the system's administrator. There are three levels of access: administrator (full), user (add and edit), and guest (read-only, search, export, and print). Anyone may be given add/edit permission and data can then be entered by designated technicians, research administrative assistants (secretaries), and/or the pathologists. Guest users may not make any additions or changes. They have full abilities to view, search, and print. There is no tracking of updates to the database fields by user or time stamp. Several date fields are provided to track the completion status of a case: "submission date", "preliminary date", "telephone date", and "final date". We use these to indicate when a case was necropsied, when the signalment was entered, when the pathologist telephoned the investigator or otherwise reported results, and when the case's data entry was complete. The database contains fields for identification of the pathologist allowing records input by multiple users to be discriminated. In addition, there are the fields, "clinical number" and "code", in which codes, grant numbers, project numbers, etc. can be entered to allow filtering by project ID and, for example, Principal Investigator.

Discussion

Inability to maintain organized records limits the speed at which a study can be finalized and reports/manuscripts generated. Data collection based on a single accession number incorporated into a relational database enables centralized storage of all information on an individual, pooling of information in the form of summary reports, or copying to compatible programs. Spreadsheets, statistical packages, colony management, and gene mapping software that are compatible with the database make rapid analysis possible. Many marker analysis programs today accept data from standard spreadsheet programs as input. Many different programs of these types exist, and it would be impossible to identify all of them here.^{4,13} Use of a single accession number for all materials collected from an individual provides ease of

storage and retrieval of specimens. Retrospective studies can be done using both the stored data and finding the original material for verification, photography, etc. Linking to an online reference database (Pathbase) provides a "virtual second opinion" for those pathologists not familiar with the complexity of mouse pathology.

The population of major mouse phenotype databases such as Europhenome¹⁰ and Mouse Phenome Database³ will be facilitated by the development of data capture tools such as MoDIS which can capture pathology data and upload into central databases. Standardization of semantics and data structure is an essential prerequisite for the seamless integration of this data and we believe that the principles followed in the development of MoDIS provide a useful model for expert-friendly tools which gather and structure the data needed for phenotype analysis. It is important that the tool be familiar and intuitive to the expert collecting information in this case the pathologist and as a consequence the data is uploaded in a curated form which requires no further restructuring before integration into the central database.

The software and documentation for MoDIS can be downloaded from (http://research.jax.org/faculty/sundberg/index.html).¹⁸

Acknowledgments

This work was supported by grants from the Ellison Medical Foundation, the National Institutes of Health (CA34196, CA89713, and AG025707), and the Commission of the European Community, Framework Programme 6 contract nos. LSHG-CT-2006-037811; CASIMIR, and EUMODIC; LSHG-CT-2006-037188.

Abbreviations

JEEPEG	Joint European Experimental Pathology Expert Group
LIMS	laboratory information management systems
MGI	Mouse Genome Informatics
MoDIS	Mouse Disease Information System
MPD	Mouse Phenome Database
MTB	Mouse Tumor Biology Database
QTL	quantitative trait locus

References

- Beamer WG, Shultz KL, Churchill GA, Frankel WN, Baylink DJ, Rosen CJ, Donahue LR. Quantitative trait loci for bone density in C57BL/6J and CAST/EiJ inbred mice. Mammalian Genome 1999;10:1043–1049. [PubMed: 10556421]
- Beamer WG, Shultz KL, Donahue LR, Churchill GA, Sen S, Wergedal JR, Baylink DJ, Rosen CJ. Quantitative trait loci for femoral and lumbar vertebral bone mineral density in C57BL6J and C3H/ HeJ inbred strains of mice. J Bone Mineral Res 2001;16:1195–1206.
- Bogue MA, Grubb SC, Maddatu TP, Bult CJ. Mouse Phenome Database (MPD). Nucleic Acids Res 2007;35:D643–649. [PubMed: 17151079]
- Churchill GA, Doerge RW. Empirical threshold values for quantitative trait mapping. Genetics 1994;138:963–971. [PubMed: 7851788]
- Dell'omo G, Shore RF, Lipp HP. An automated system, based on microchips, for monitoring individual activity in wild small mammals. J Exp Zool 1998;280:97–99. [PubMed: 9437856]
- Farmer MA, Sundberg JP, Bristol IJ, Churchill GA, Li R, Elson CO, Leiter EH. A major quantitative trait locus on chromosome 3 controls colitis severity in IL-10-deficient mice. Proc Natl Acad Sci USA 2001;98:13820–13825. [PubMed: 11707574]

- 7. Hayamizu TF, Mangan M, Corradi JP, Kadin JA, Ringwald M. The adult mouse anatomical dictionary: a tool for annotating and integrating data. Genome Biol 2005;6:R29. [PubMed: 15774030]
- Ishimori N, Li R, Walsh KA, Korstanje R, Rollins JA, Petkov P, Pletcher MT, Wiltshire T, Donahue LR, Rosen CJ, Beamer WG, Churchill GA, Paigen B. Quantitative trait loci that determine BMD in C57BL/6J and 129S1/SvImJ inbred mice. J Bone Miner Res 2006;1:105–112. [PubMed: 16355279]
- Mahler M, Most C, Schmidtke S, Sundberg JP, Li R, Hedrich HJ, Churchill GA. Genetics of colitis susceptibility in IL-10-deficient mice: backcross versus F2 results contrasted by principal component analysis. Genomics 2002;80:274–282. [PubMed: 12213197]
- Mallon AM, Blake A, Hancock JM. EuroPhenome and EMPReSS: online mouse phenotyping resource. Nucleic Acids Res 2008;36:D715–718. [PubMed: 17905814]
- 11. Masuya H, Nakai Y, Motegi H, Niinaya N, Kida Y, Kaneko Y, Aritake H, Suzuki N, Ishii J, Koorikawa K, Suzuki T, Inoue M, Kobayashi K, Toki H, Wada Y, Kaneda H, Ishijima J, Takahashi KR, Minowa O, Noda T, Wakana S, Gondo Y, Shiroishi T. Development and implementation of a database system to manage a large-scale mouse ENU-mutagenesis program. Mamm Genome 2004;15:404–411. [PubMed: 15170230]
- Schofield, PN.; Bard, JBL.; Rozell, B.; Sundberg, JP. Computational pathology: challenges in the informatics of phenotype description in mutant mice. In: Sundberg, JP.; Ichiki, T., editors. Handbook on genetically engineered mice. CRC Press; Boca Raton: 2005. p. 61-81.
- Sen S, Churchill GA. A statistical framework for quantitative trait mapping. Genetics 2001;159:371– 387. [PubMed: 11560912]
- Sundberg BA, Sundberg JP. A database system for small diagnostic laboratories. Lab Anim 1990;19:55–58.
- Sundberg, BA.; Sundberg, JP. Medical record keeping for project analysis. In: Sundberg, JP.; Boggess, D., editors. Systematic approach to evaluation of mouse mutations. CRC Press; Boca Raton: 2000. p. 47-55.
- Sundberg JP, Boggess D, Silva KA, McElwee KJ, King LE, Li R, Churchill G, Cox GA. Major locus on mouse chromosome 17 and minor locus on chromosome 9 are linked with alopecia areata in C3H/ HeJ mice. J Invest Dermatol 2003;120:771–775. [PubMed: 12713579]
- 17. Sundberg JP, Schofield PN. A mouse by any other name. J Invest Dermatol. in press.
- Sundberg JP, Sundberg BA, Schofield PN. Integrating mouse anatomy and pathology ontologies into a diagnostic/phenotyping database: tools for record keeping and teaching. Mammalian Genome 2008;19:413–419. [PubMed: 18797968]
- Yuan R, Tsaih S-T, Petkova SB, deEvsikova CM, Xing S, Marion MA, Bogue MA, Mills KD, Peters LL, Bult CJ, Rosen CJ, Sundberg JP, Harrison DE, Churchill GA, Paigen B. Aging in inbred strains of mice: study design and interim report on median lifespans and circulating IGF1 levels. Aging Cell. submitted.

-

		Maintenance Funct	In the second						
Iouse Disease Information Sy		1	tions:						
dd or Edit Mouse Phint RFC Sheets Searches	Add Spreadsheet Data	Update Information Tables	Quit						
/ider Add/Edit Mouse	Path Base Terms	Mouse Anatomy Ter	ms						
arge Font Add/Edit Mouse	1	Disease Description Te	ms						
Data Entry for cases								-	•
				er out of g	ight click in field, in the funnel icon varacters. See o	type "text" in the	except Dx or S e filter for box and iter. This will not x help with search ests	press tab. Cli ind the * ? an	
Contraction of the second s	itrain by Stock #	659 Stock # that m			Name • Org BW 💌	anis - N - F 0 0	 Resul - Uni 0 19.85 g 0 	Chari	
GENOTYPE:	1000 C	Animals: 1 Sex F	Line #:	-		0	0		-
Pedigree Number	Contraction of the second s	Pedigree.							
Dam Pedigree:		2 Pedigree:							
	Code			•					
Date Born: 3/7/2005 Age in days	and the second se	Clinical No.: 01€		•					
Date Bom 3/7/2005 Age in days Construct Name:	795 Maling	Clinical No: 01€	LLS						
Date Born 3/7/2005 Age in days Construct Name: IT Color Photos IT BAY Photo IT Mic	795 Maxing 1	Elinical Na: 01€ Frozen Tissue I⊄ Wax	LLS						
Date Bonx 377/2005 Age in days Construct Name: □ Color Photos □ BAV Photo □ Mic 17 Histology □ X-tays □ SEM □ in	795 Maxing moscopic Photos I situ Clin Path	Clinical No: 01€ Frozen Tissue I Wax RNA	LLS Blocks IT Wet Ti	ioue 4					
Date Bonx 3/7/2005 Age in days Construct Name: Color Photos IT B/W Photo IT Mic PHistology IT X-tays IT SEM IT in Dix Order Organ or Tissue Search	795 Maxing moscopic Photos I situ Clin Path	Elinical Na: 01€ Frozen Tissue I⊄ Wax	Blocks F Wet Ti	toue		Disease Terms		evenity Score	_
Oate Born: 2/7/2005 Age in days Construct Name: If Construct Name: If Color Photos BAV Photo Meto If Hotogy K-stays SEM Fin Dx Order Organ or Ticsue Seatch 0.0 (eye MA.0000261) Fin	795 Mating moscopic Photos in Cin Path PathBase Anatomy	Clinical No: 01€ Frozen Tissue I Wax RNA	Blocks F Wet Ti	ioue 4	14		- M	PATH 3	1
Date Bonx 3/7/2005 Age in days Construct Name: Color Photos IF BAW Photo IF Mic IF Histology IF X-tays IF SEM IF in Dx Order Organ or Tissue Search	795 Mating moscopic Photos in Cin Path PathBase Anatomy	Clinical No: 01€ Frozen Tissue IF Wax RNA Search Disease Tem	Blocks F Wet Ti	toue			or view disease d	PATH 3 of Images	1
Oate Born: 2/7/2005 Age in days Construct Name: If Construct Name: If Color Photos BAV Photo Meto If Hotogy K-stays SEM Fin Dx Order Organ or Ticsue Seatch 0.0 (eye MA.0000261) Fin	795 Mating moscopic Photos in Cin Path PathBase Anatomy	Clinical No: 01€ Frozen Tissue IF Wax RNA Search Disease Tem	ELLS Blocks IT Wet Ti ns Dx 2000 MA degenerativ	toue	Add new o	isease Edit	CI View disease d	PATH 3 M. Images PATH 3	1
Date Bork 3/7/2005 Age in days Construct Name	795 Mating moscopic Photos 1 altu Clin Path PathBase Anatomy 090	Clinical No.: [01 € Frozen Tissue IF Wax RNA Search Disease Tem	ELLS Blocks IT Wet Ti ns Dx 2000 MA degenerativ	issue 4 (i in PathBase DX 0 e change MPATH: mation MPATH:21	Add new o	isease Edit	or view disease d	PATH 3 M. Images PATH 3	
Date Bork 3/7/2005 Age in days Construct Name	795 Mating moscopic Photos 1 altu Clin Path PathBase Anatomy 090	Clinical No.: [01 € Frozen Tissue IF Wax RNA Search Disease Tem	ELLS Blocks IT Wet Ti ns Dx 2000 MA degenerativ	in PathBase DX: 0 in PathBase DX: 0 in echange MPATH: mation MPATH:21	14 Add new o	isease Edit	or view disease d	PATH 3 M. Images PATH 3	
Date Bork 3/7/2005 Age in days Construct Name	795 Mating: noscopic Photos F I situ F Clin Path F PathBase Anatomy 090 with JPS:00006	Clinical No: [01 € Frozen Tissue IF Wax RNA Search Disease Tem	LLS Blocks IT Wet Tr ms Dx Seatt MA degenerativ MA acute inflam	in PathBase DX: 0 in PathBase DX: 0 in echange MPATH: mation MPATH:21	14 Add new o	isease Edit	or view disease d	PATH 3 ef. Images PATH 3 ef. Images PATH 1	
Oate Bork 3/7/2005 Age in days Construct Name: F Color Photos BAV Photo Main Color Photos BAV Photo Main BAV Main Valoage Xayar SEM In In Div Order Orgen or Tissue Seatch Seatch Direase name: [phrink:bubit JPS:00] Oil (molar MA:000056) Direase name: [primps body periodor Direase name: [primps body periodor Oil (save MA:0001459) Direase name: [primps body periodor	TPS Making Toscopic Photos F I PathBase Anatom UB0 Whis JPS:00006 PS:00019	Clinical No. (01 E Frozen Tissue Wax RNA Search Disease Tem	LLS Blocks IT Wet Tr ms Dx Seatt MA degenerativ MA acute inflam	ih PahBase DX:0 e charge MPATH: mation MPATH:21 MPATH:134	Add new of	isease Edit isease Edit isease Edit	M M Or view disease d M Or view disease d M Or view disease d Or view disease d	PATH 3 ef. Images PATH 3 ef. Images PATH 1	
Date Bork 3/7/2005 Age in days Construct Name:	795 Maxing moscopic Photos I stul ClnPath EathBase Acatemu 090 Whit JPS:00006 PS:00019 pilor BOTH: Comments: do alegoesia min Hud, love in pale ani	Ckricel No. [61 E Fracen Tissue Wax RNA Search Disease Tem	Blocks IT Wet Tr Blocks IT Wet Tr MA degenerativ MA scute inflam MA hyperplasia	ih PahBase DX:0 e charge MPATH: mation MPATH:21 MPATH:134	Add new o	isease Edit isease Edit isease Edit	M M Or view disease d M Or view disease d M Or view disease d Or view disease d	PATH 3 ef. Images PATH 3 ef. Images PATH 1	
Construct Name: Colo Photos T B/W Photo T Mic P Histology T X-lays T SEM T in Dx Oxder Organ or Tissue <u>Seasch</u> 0.0 Gree MA:0000261 Disease name: [phthrsis bubsi JFS:001 0.0 [mota: MA:0000350 Disease name: [phthrsis bubsi JFS:001 0.0 [mota: MA:0000350 Disease name: [ptersiseous] teison JI Symptome / Chrisal Signe / Gross Descr From LG: Cage ID: C3H-T18: E + Tag Symptome / Chrisal Signe / Gross Descr From LG: Cage ID: C3H-T18: E + Tag renthal areas conclaim Mid to moderate a (1.5 cm, 0.6 cm) are entraged and comb	795 Maxing rescoppic Photos F atu Chn Pah Esthill ace Anatoms 980 whit JPS:00006 PS:00019 point BDTH: Comments: du alspesito jale ar content la lage emounts	Clinical Na. [01 E Triazen Tissue Weak RNA Search Disease Tem stal arsol and betrine hanss of discoktreat of dark	LLS Blocks Wet Ti ms Dx Seal MA degenerativ MA acute inflam MA pyperplasia agnosis Memo / Co	ih PahBase DX:0 e charge MPATH: mation MPATH:21 MPATH:134	Add new o	isease <u>Edk</u> isease <u>Edk</u> isease <u>Edk</u> isease <u>Edk</u>	W M or view disease d W M or view disease d W M or view disease d larger,	PATH 3. PATH 3. PATH 3. M. Images PATH 1. M. Images ets the pathol	

Figure 1. Example showing the user interface and a data entry screen

The MoDIS entry form requests basic information on the mouse (signalment) to be entered in a systematic manner to enable data retrieval. Strain and stock number can be verified with the dropdown menus.

A	B	C	D	E	F	G	H	1	J	K	L	M	N
Expt. Week {ignored}	Expt. Date {ignored}		Case Number	Ear notch		Cage ID	Sex	birth date	PIXIMIS (Y/N)	Micronuclei (Y/N)	HSC or Bone	EKG (Y/N)	Submission date
1/3/2008	Mon.	00069A24E7 00069A24E8		L R	A/J A/J	CS-AJ-M1-01 CS-AJ-M1-02		01/23/05 01/24/05	Y Y	Y Y	HSC HSC	Y Y	1/5/2008 1/5/2008
▶ N Sheet1	0/		1995				10016	UKL.		3			Jean N.
	1000	New Cases Before		o be a	dded	Check and	ledit	this inf	ormation	- =	×		
	cas		any that	t shoul	Id not I	be importe	d. Cl	ick ADE		CORDS to			
		e number 0 8 - 7000 0			mber if	it has one		th Date /23/2008		bmission_dat 1/5/2008	0		
	Stra	ain			1					Stock 64			
		en Denement			nto el fre		alda			1	_		
	Fro	ss Descriptions m: CROSS S ronuclei:'Y: H	ECTION.	AL: Ca	ge ID: 0	CS-AJ-M1-0		Contraction of the					
	Sut JP:	omitter Sundberg/El	# Anii lison	mals P 1 J	atholoc PS	list[Histolog	y Wa	x block][r	Frozen Tis	Sue Clinical 01-ELL			
	DC ME vio rul	NOT COM SSAGE S lation is ca e violations	PLETE TATING used by are cau	THE A IT CAI one of ised by	DD IF NNOT r more y: the s	YOU REC APPEND A case num	LL T bers	HE RE already F; a dat	CORDS. existing e is inval	A key validation			

Figure 2.

An illustration of the procedure designed to download data from mice in a study in which each animal has a chip embedded subcutaneously for identification. MoDIS includes a sample Microsoft Excel spreadsheet for importing data (top). Information in specific fixed columns can be transferred by loading the spreadsheet into this form (bottom), which allows user editing of the data prior to finishing the download. This type of transfer eliminates entering this data again with the potential for errors.

	Disease Terms Dx Search PathBase D	X Ontology List of Disease Terms	Severity Score
0.0 eye MA:0000261	MA_degenerative change MPA	TH:14	MPATH 3
Disease name: phthisis bulbi JPS:00090		Add new disease Edit or v	iew disease def. Images
0.0 molar MA:0000350	MA acute inflammation MPATH	1:213	MPATH 3
Disease name: foreign body periodontitis JPS:00006		Add new disease Edit or v	iew disease def. Images
0.0 bone MA:0001459	MA hyperplasia MPATH:134		MPATH 1
Disease name: fibroosseous lesion JPS:00019		Add new disease Edit or v	iew disease def. 🛛 Images 🖕
Symptoms / Clinical Signs / Gross Description:	Diagnosis Memo / Comments:	SHIFT - F2 to Zoom fields larg	ger.
From: LG; Cage ID: C3H-F-18; Ear Tag: BOTH; Comments: dorsal and ventral areas contain mild to moderate amounts of alopecia. Uterine horn			
(1.5 cm, 0.6 cm) are enlarged and contain fluid, liver is pale and discolore gingival sulci is moderately swollen and contain large amounts of dark	ed,		
(1.5 cm, 0.6 cm) are enlarged and contain fluid, liver is pale and discolore gingival sulci is moderately swollen and contain large amounts of dark material			inalize case sets the pathologist
[1,5 cm, 0,6 cm] are enlarged and contain fluid, liver is pale and discolore gingival sulci is moderately swollen and contain large amounts of dark material "athologist JPS [] Preliminary Date: 11/19/2007 Telephone Date:		1 1 1	inalize case sets the pathologist o the default, final date to today, saves, and prints this case.

Figure 3.

Entering the diagnoses, reviewing findings at necropsy, and notes on the diagnoses. The first field lists the organ (eye) with an anatomy code that is automatically added as the organ or tissue name is written out. The second field on the right enters the disease process (degenerative changes) with an MPATH code that again is automatically coded. The next field is a disease severity code (3). The line below lists the disease name (phthisis bulbi) with a code provided by the pathologist of record. The lower left field provides gross features at the time of necropsy and the blank field to the right provides space for a description or comments on the specific lesions.

l erms to	r Disease Descriptions	
ID number JPS	ADD a new term Clicking the ADD button will create a new ID number for you and allow you to enter the disease name and definition. Then click the SAVE new term button. Use DON'T SAVE to change your mind.	
Disease name Ley	dig cell hyperplasia	
The name has a lim	nit of 100 characters. A change to an existing disease name will NOT change diagnoses where it	
was used previously		
was used previous Definition of the		
was used previous Definition		0
was used previous Definition of the term:		(

Figure 4.

Adding a disease name is followed by clicking on *Add a new term* to enter a local definition. The term is automatically assigned a code based on the initials of the pathologist. This definition may later be reviewed by the MPATH pathology panel and a consensus arrived at before it is formally accepted into the public ontology.




Figure 5.

Links to the Pathbase.net website (http://www.pathbase.net; buttons with blue underlined text) will open a searchable web page (as another window) showing the hierarchy of ontology terms, using the system's default browser.

e Edit View Favo	eulep.pdn.cam.ac.uk/Search_Path	ente unite hub	• 4 X Google	م
ogle 8 -		earch 🛉 🚇 • 🎕 🖉 🐑 🚇 • 🛛 🖏 • 🖂 Find • 🏷 Check • 🕤 AutoFill •	🖏 + 🍚 Sign In + 🛛 🏶	Convert 💌 🔂 Sele
🕼 🔝 Pathbase			©•⊜•⊗	Eage 🔹 🍈 Tools
PB 13834	Gene: N/A Pathology: MPATH 458 - normal	Normal optic nerve		
PB 13835	Gene: N/A Pathology: MPATH 458 - normal	Normal optic nerve		
PB 13836	Gene: N/A Pathology: MPATH 458 - normal	Normal retina near optic nerve		

Figure 6.

Sample Pathbase page illustrating the normal optic nerve and retina.

	Name -	. (Organis 🗸	1 -	F.	Resul -	Uni 🗸	Charg -	1000	Description	Se -
	BW			0	0	19.85	9				
*		-		0	0						Contraction of the second
4	MSKI MSPL MTHO MTRA MUTE MUTEC MVAG MVAGF MVEN		Skin Spleen Thoracic Trachea Uterus Uterine o Vagina Vaginal Ventral i	conter fluid	nts						,

Figure 7.

This field illustrates addition of test results to the record. BW is the code for body weight. The animal's weight, 19.85 g is entered. The drop down menu provides lists of other tests that can be automatically entered using abbreviations.

INSTITUTION NAME Institution address line

J.P.S., DVM Dinlomate ACVP Jane Doe

	D VAL, D I DIMALE I C VI	Acepear cu rassiscant	
From JP Sundberg/Ellison		Bldg Rm. 5-G3B	Date 5/11/2007
Male 718 Days of	d Born 5/23/2005 Pair	Line N/A	
Strain SWR/J			Necropsy Number
Stock#689 Genotype	+/+		07 - 10099
Pedigree #'s SH-537	Sire: SC-002	9	Clinical Number 01-ELLS
Dam 1: SC-1548	Dam 2: SC-1	572	01-ELLS
Other ID 00069A24E8	Code LG	Source R Category B	
Construct name N/A			

Items recorded include: Color photographs Photomicrographs Electronmicrographs Radiographs In Situ Prep. Histology Frozen Tissue Wax blocks Wet Tissue Clinical Pathology RNA

History and Clinical Signs -- Gross Description

From: LG; Cage ID: SWR-M-06; Ear Tag: RIGHT; Comments: Gingival sulci are severly swollen and contain small amounts of dark material between molars. Lungs contain small red foci throughout lobes.

Special Tests		#Neg #	Pos	Serology	Organism
Results	Charge				Description
Bodyweight		0	0		
28.02 g					
Microbiology.					
Lung		1	0	Pasteure	lla pneumotropica

Diagnosis - Comments

Molars. Broken hair fibers are packed between the mucous membranes and the molars. The area around the molars is ulcerated and acutely inflammated. This is hair impaction around the tooth, a common grooming abnormality seen in many strains of mice.

Lungs. Two descrete bronchoalveaolar adenomas are present. There are no other associated abnormalities in the lungs.

	Morphologic Diagnoses
	molar MA:0000350
	severe acute inflammation MPATH: 213
	foreign body periodontitis JPS:00006
	lung MA:0000415
0	moderate bronchioloalveolar adenoma MPATH:272
8	From Pathologist: JPS Date Telephoned 11/28/2007 Preliminary Report 11/19/2007 Final Report 12/6/2007

Figure 8.

Final printed report.

Search by Mu	Help	Click the Filter By Form search or to remove the	this form is first opened it shows a button and fill in the search criteria. Click e search. Note - a case will be repeated test. Click the help button for more instruct	the Funnel	button to activate	e the	nly.
Submitter	💌 Bld	lg. Rm	Submission date	Dx Order:	Score		
Strain			Stock # Assigned	Organ or		,	
Genotype		▼ # Animals	Sex "F' Line #	Tissue			•
Pedigree Number		Sire Pedigree			1	-	_
Dam Pedigree		Dam 2 Pedigree		Diagnos	is "degenerative change MPA1		-
Other ID		Code Sou	urce Category		and go in ear		
Date Born	Age in days	Mating Clinica	al No. ["01-ELLS"	Disease			
Construct Name		and a second second second		name:			-
層 Histology 層 X-rays 同	🖩 SEM 🕅 in situ 🗐 (ee 層 Wax Blocks 層 Wet Tissue ents:	Test name Organism: Number n	egative:	•	
層 Histology 層 Xrays Symptoms / Clinical Signs	副 SEM III in situ III (/ Gross Description: 1	Clin Path 🔲 RNA Diagnosis Memo / Comm	ents:	Organism Number n Number p Serology: Results: Units:	egative:		
m Histology m Xrays li Symptoms / Clinical Signs	🖩 SEM 🕅 in situ 🗐 (Clin Path 📧 RNA		Organism Number n Number p Serology: Results: Units: Charge:	: iegative: iositive: 		
爾 Histology 南 X-rays 斯 Symptoms / Clinical Signs	IF SEM IFF in shu IFF (C /Gross Description: 1 http://www.com/opension.com/opension/ http://www.com/opension/ http://wwwwww.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://www.com/opension/ http://wwww.com/opension/ http://www.com/opension/ http://wwww.com/ope	Clin Path 🔲 RNA Diagnosis Memo / Comm	ents:	Organism Number n Number p Serology: Results: Units:	i iegative: iositive: 		
F Histology F Xrays Symptoms / Clinical Signs Pathologist: F F Show in spreadsheet (c	IF SEM IFF in shu IFF (C / Gross Description: 1 http://www.sec.org/files/file	Clin Path III RNA Diagnosis Memo / Comm	ents:	Organism Number n Number p Serology Results: Units: Charge: Descriptio	i iegative: iositive: 		
E Histology E X-rays II Symptoms / Clinical Signs Pathologist: F F Show in spreadsheet fo pook for / Or /	IF SEM IFF in shu IFF (C / Gross Description: 1 http://www.sec.org/files/file	Clin Path III RNA Diagnosis Memo / Comm	ents:	Organism Number n Number p Serology Results: Units: Charge: Descriptio	i iegative: iositive: 		
Histology Krays Rymptoms / Clinical Signs Pathologist: Pathologist: F Show in spreadsheet fo pok for Or Search Form	IF SEM IF in shu IF (/ Gross Description: 1 heliminary Date: mmat Print this case	Clin Path III RNA Diagnosis Memo / Comm	Final Date:	Organism Number m Number p Serology: Results: Units: Charge: Descriptic Comment	i iegative: iositive: 		niss
Histology (Xrays Symptoms / Clinical Signs Pathologist: F Show in spreadsheet to ook for Or Search Form	T SEM IF in shu IF 0 / Gross Description: 4 heliminary Date: Print this case -1 Sex -37 Age - F 620 s	Clin Path III RNA Diagnosis Memo / Comm Telephone Date:	ents: Final Date:	Organism: Number n Number p Serology: Results: Units: Charge: Descriptic Comment: Score - 3	nr.		niss

Figure 9.

Top: Example of the search form for entering multiple criteria. A search on the "degenerative change MPATH:14" diagnosis, "01-ELLS" clinical no., and "F" sex resulted in the datasheet (bottom) listing 76 cases, which may be printed or exported in Microsoft Excel and several other formats.