

Performance and Clinical Evaluation of a Sensitive Multiplex Assay for the Rapid Detection of Common *NPM1* Mutations

Michael Hafez,* Fei Ye,[†] Keith Jackson,* Zhe Yang,[†]
Judith E. Karp,[‡] Emmanuel Labourier,[†]
and Christopher D. Gocke*[‡]

From the Departments of Pathology,* and Oncology,[†] Johns Hopkins University School of Medicine, Baltimore, Maryland; and Asuragen Inc.,[‡] Austin, Texas

Determination of *NPM1* mutation status has become essential for the molecular classification of acute myeloid leukemias (AML). Methods with high clinical sensitivity and specificity adapted to the molecular laboratory workflow are required for the diagnosis, prognosis, and monitoring of AML with normal karyotype. We report here the development and evaluation of a novel, streamlined, RNA-based assay for the rapid multiplex detection of common *NPM1* mutations in a 96-well assay format. Using synthetic transcripts and total RNA from leukemic cell lines, we show that the assay can specifically detect *NPM1* wild-type and mutants A, B, D, or J transcripts in the same reaction. Dilution experiments indicate an assay dynamic range >4 log units with an analytical sensitivity of approximately 0.01%. Evaluation of 69 clinical specimens at initial diagnosis resulted in 100% agreement with reference methods. Of patients with AML with normal karyotype, 53% carried one of four different mutations. The assay was also combined with other laboratory-developed tests to simultaneously detect *NPM1* mutant transcripts and fusion transcripts resulting from t(8;21) or inv(16) in a single reaction well. Overall, these results show that the assay is a versatile and specific tool for the screening of *NPM1* mutations in patients with AML. Its high analytical sensitivity further suggests potential utility for the monitoring of residual disease in AML with normal karyotype. (J Mol Diagn 2010, 12:629–635; DOI: 10.2353/jmoldx.2010.090219)

The nucleophosmin gene (*NPM1*) is thought to be the most frequently mutated gene in *de novo* acute myeloid leukemias (AMLs), particularly those with normal cytogenetics. In several series, approximately 25 to 35% of adult AMLs demonstrated mutation in one allele of *NPM1* (located on chromosome 5q35) at the time of diagnosis.^{1–3} The nucleophosmin protein is a nuclear shuttle protein

important in centrosome duplication, ribosome biogenesis, and regulation of transcription factors and also functions as a molecular chaperone.^{4–7} Mutations in *NPM1* in AML occur characteristically in the last exon of the gene, exon 12, and, although heterogeneous, are most commonly 4-bp insertions that result in frame shifts. These mutations typically do not truncate the protein but create a new nuclear export signal and remove nucleolar-localizing tryptophan residues.⁸ These changes result in the aberrant cytoplasmic localization of the mutant and remaining wild-type *NPM1* protein in a dominant-negative fashion.⁸

NPM1 mutation has important prognostic implications for patients with AML. Patients with *NPM1* mutations typically do not carry cytogenetic abnormalities and generally respond better to induction chemotherapy than those without mutations.¹ Several reports identified a prognostic interrelationship between *NPM1* mutation, internal tandem duplication mutation in the *FLT3* tyrosine kinase, and, perhaps, mutation in the CCAAT/enhancer binding protein α gene.^{9–12} These findings suggest that molecular stratification of AML cases may be important for prognosis and for therapeutic decisions, underlining the need for timely diagnostic testing for these mutations.

A number of methods of identifying *NPM1* mutations have been published. Because almost all cases of *NPM1* mutation result in length changes to exon 12, PCR amplification of the exon followed by sequencing or sizing by capillary electrophoresis or by high-performance liquid chromatography have all been used.^{13,14} Both genomic DNA and mRNA have served as substrates for these assays. Locked nucleic acid-clamping of the wild-type allele with chip-based detection of mutants has also

Supported in part by the National Institutes of Health/National Cancer Institute (Small Business Innovation Research grant R43CA130501).

M.H. and F.Y. contributed equally to this work.

Accepted for publication March 5, 2010.

F.Y., Z.Y., and E.L. were employees of Asuragen, Inc. (the manufacturer of the Signature *NPM1* Mutations assay) at the time of this work. None of the other authors declare any relevant financial relationships.

Supplemental material for this article can be found on <http://jmd.amjpathol.org>.

Address reprint requests to Christopher D. Gocke, M.D., Department of Pathology, Johns Hopkins Medical Institutions, 600 N. Wolfe St., Park SB202, Baltimore, MD. E-mail: cgocke1@jhmi.edu.

been used.³ The mutation spectrum in AML is fairly wide, with almost 40 different variants identified. However, the common mutations are relatively few, with the A mutation (a duplication of a TCTG tetranucleotide at position 956 to 959 of the reference sequence) accounting for 75 to 80% of mutant-bearing AML, the B mutation for approximately 10%, and the D mutation for approximately 5%; other mutations are very rare.¹⁵ Because the heterozygous mutations result in 4-nucleotide repeats or insertions in a background of wild-type sequences, sensitive and specific detection of *NPM1* mutations using standard PCR methods can be challenging.

We report here the rapid detection of common *NPM1* mutations using total RNA purified from cultured cells, bone marrow, or peripheral blood using a laboratory-developed test based on the Signature *NPM1* Mutations Research Use Only kit. The assay uses multiplex RT-PCR (RT-PCR) in combination with fluorescent bead-based detection to simultaneously identify transcripts for *NPM1* mutations and wild-type targets. The aim of this study was to evaluate the assay workflow and performance in a clinical setting.

Materials and Methods

RNA Specimens

Peripheral blood or bone marrow from patients with AML, acute lymphoid leukemia, chronic myeloid leukemia, acute bilineage leukemia, or myelodysplastic syndrome was collected as part of standard clinical care and following a protocol approved by The Johns Hopkins Medical Institutions Institutional Review Board. Total RNA was isolated from peripheral blood or bone marrow using the RNeasy Mini Kit according to the manufacturer's instructions for blood (Qiagen, Valencia, CA).

OCI-AML3 cells (DSMZ GmbH, Braunschweig, Germany) were cultured according to the supplier's recommendations, and total RNA was isolated with the RNeasy Mini Kit according to the manufacturer's instructions for cultured cells (Qiagen). When indicated, total RNA was diluted in purified HL-60 total RNA (Applied Biosystems, Foster City, CA), keeping the concentration of total RNA constant. Synthetic RNAs corresponding to mutations A, B, D, or J were prepared using standard *in vitro* transcription methods and diluted in HL-60 total RNA, keeping the concentration of total RNA constant. Purified total RNA was quantified using a NanoDrop ND1000 (NanoDrop Technologies, Wilmington, DE).

Signature Assays

The Signature *NPM1* Mutations Research Use Only kit (Asuragen Inc., Austin, TX) was evaluated in combination with Signature General Purpose Reagents (Asuragen Inc.). In brief, total RNA (5 μ l) was reverse-transcribed into cDNA in a 20- μ l reaction for 45 minutes at 42°C. The resulting cDNA (5 μ l) was amplified by multiplex PCR using 1 μ l of Signature *NPM1* Mutations Primer Set and 2.5 units of AmpliTaq Gold (Applied Biosystems) in a

20- μ l reaction (45 cycles consisting of 94°C for 30 seconds, 55°C for 30 seconds, and 72°C for 30 seconds). The PCR products (5 μ l) were then hybridized to 5 μ l of Signature *NPM1* Mutations Probe Set in a 50- μ l reaction for 30 minutes at 52°C. After addition of 25 μ l of Reporter Solution, the probe-bound PCR products were detected on a Luminex 200 system (Luminex Corporation, Austin, TX) with subsequent analysis of the median fluorescent intensity (MFI) signal outputs for bead ID #28 (*NPM1* wild-type) and #45 (*NPM1* mutant) in Excel (Microsoft, Richmond, WA). Positive and negative controls provided with the kit were included in each run.

For co-amplification and codetection of specific AML fusion transcripts, the same protocol was followed except that 1 μ l of Signature AML/ETO, *inv(16)*-ADE Primers (analyte-specific reagent; Asuragen Inc.) was added to the multiplex PCR reaction, and 5 μ l of Signature AML/ETO, *inv(16)*-ADE Probes (analyte-specific reagent) was added to the hybridization reaction. MFI signals corresponding to AML/ETO, *inv16*-A, *inv16*-B, and *inv16*-E fusion transcripts were detected on bead ID #14, #55, #16, and #41, respectively. Positive and negative controls were included in each run; the positive control consisted of pooled, *in vitro* transcribed RNAs for the pertinent targets in a background of total RNA and thus also controlled for the reverse transcription step. A GeneAmp PCR System 9700 (Applied Biosystems) was used for the RT, PCR, and hybridization steps.

Other Assays

Standard PCR was performed using 5 μ l of cDNA generated with the Signature reagents (Asuragen Inc.), pairs of primers flanking the *NPM1* exon 12 mutation site (Integrated DNA Technologies, Inc., Coralville, IA), and 35 PCR cycles (94°C for 30 seconds, 55°C for 30 seconds, and 72°C for 60 seconds) followed by a 10-minute elongation step at 72°C. Bidirectional sequencing of the PCR products using the BigDye terminator method was performed by ACGT Inc. (Wheeling, IL). For fragment size analysis, one of the PCR primers was 5'-modified with a FAM dye. The resulting labeled PCR products were diluted 1:50 in water, and 1 μ l was mixed with 0.5 μ l of GeneScan 500 ROX Size Standard and 13.5 μ l of Hi-Di Formamide (Applied Biosystems). After heat denaturation at 95°C for 2 minutes, analysis was performed on a 3130xl Genetic Analyzer (Applied Biosystems) using a 36-cm capillary with POP-7 Polymer and GeneMapper Software v4.0 (Applied Biosystems). Specific run parameters were prerun at 15 kV for 180 seconds, with injection at 1.2 kV for 23 seconds, and run at 15 kV for 30 minutes at a temperature of 60°C.

Results

Assay Design and Workflow

The Signature *NPM1* Mutations Research Use Only assay consists of four consecutive steps. Total RNA is first




Figure 1. Analytical performance evaluation. **A:** Wild-type and mutant sequences detected by the assay. **Bold** indicates insertion mutations. The relative representation of each mutation in the total mutation-positive population is indicated.^{16,17} **B:** Representative examples of signal (MFI) obtained with a serial dilution of a mutant A-positive cell line RNA (OCI-AML3) in wild-type cell line RNA (HL-60). **C:** Representative examples of signal (MFI) obtained with a serial dilution of mutant A *in vitro* transcript (mut-A) in wild-type cell line RNA (HL-60). **D:** Representative examples of data obtained with a serial dilution of four different *in vitro* transcripts (mut-A, mut-B, mut-D, or mut-J) in wild-type cell line RNA (HL-60). The graph shows the ratio of mutant to wild-type signal according to the initial copy number input of each *in vitro* transcript. WT, wild-type.

reverse-transcribed (RT step); then cDNA corresponding to *NPM1* mutation A, B, D, or J (Figure 1A)^{15,16} is amplified by PCR using gene-specific primers (PCR step). Labeled amplified products are directly hybridized to target-specific capture probes covalently bound to fluorescent microspheres (hybridization step). After incubation with a reporter molecule, the reaction mixture is directly analyzed by flow cytometry to detect the mutant transcripts on a single bead type (detection step). *NPM1* wild-type transcripts are co-amplified in each sample and concurrently detected on a second bead type to serve as endogenous internal controls. The MFI detected on at least 50 beads of each type is reported to determine the presence or absence of the respective *NPM1* wild-type or mutant targets in the RNA specimen. For this study, the cutoffs recommended in the package insert were used: a specimen was reported positive if the mutant bead generated a signal ≥ 400 MFI; a specimen was reported negative if the mutant bead generated a signal < 400 MFI with a wild-type bead signal ≥ 1000 MFI.

Positive and negative control samples included in the kit are processed with each batch of specimens to validate the RT, PCR, hybridization, and detection steps. The entire procedure is performed on 96-well plates and requires seven pipetting steps per reaction including two transfer steps, one from the RT plate to the PCR plate and one from the PCR plate to the hybridization plate. With one kit, up to 24 reactions can be performed in less than 6 hours with about 1.5 hours of hands-on time. The assay requires a Luminex 100 or 200 system and RNA isolation kit reagents; both already commonly found and validated in clinical laboratories. Overall, the Signature assay design is compatible with the molecular laboratory workflow, and we undertook evaluation of its analytical and

clinical performance for the testing of newly diagnosed leukemia clinical specimens.

Evaluation of Analytical Performance

Assay performance was first evaluated in triplicate using total RNA purified from well characterized cell lines isolated from patients with leukemia. An input mass titration experiment showed that *NPM1* mutant A transcripts were reproducibly detected by the assay in 0.1 to 1000 ng of total RNA from the cell line OCI-AML3 heterozygous for the *NPM1* mutation A (Supplemental Figure 1A, see <http://jmd.amjpathol.org>). To assess the analytical sensitivity of the assay, the OCI-AML3 total RNA was serially diluted in a background of total RNA purified from the *NPM1* wild-type cell line HL-60 (Figure 1B and Supplemental Figure 1B, see <http://jmd.amjpathol.org> for all data). *NPM1* mutation A transcripts were reproducibly detected at least twofold above the positive cutoff (400 MFI) in as low as 0.1 ng of OCI-AML3 RNA in 1000 ng of HL-60 RNA (0.01% mass equivalent). Above 10% of OCI-AML3 RNA, the mutant *NPM1* signal reached a plateau with maximum MFI signals between 10,000 and 12,000 MFI. In contrast, no *NPM1* mutant signal ≥ 400 MFI was observed using up to 1000 ng of HL-60 total RNA. Analytical specificity was also confirmed using synthetic *NPM1* mutant transcripts and total RNA samples purified from translocation positive cell lines and whole blood or bone marrow clinical specimens (see below).

To further assess analytical sensitivity, a synthetic *NPM1* mutation A transcript (mut-A) was prepared by *in vitro* transcription, quantified by absorbance, and serially diluted at a specific copy number concentration in a

constant background of excess *NPM1* wild-type transcripts provided by HL-60 total RNA (10 to 50 million copies/ μ g). The mut-A transcript was reproducibly detected in as low as 2000 copies per RT reaction (Figure 1C and Supplemental Figure 1C, see <http://jmd.amjpathol.org> for all data). The MFI response curve was similar to the OCI-AML3 dilution experiment with signal saturation greater than 200,000 mut-A copies per RT reaction. Experiments with mut-B, mut-D, and mut-J synthetic transcripts confirmed that the Signature assay detects all four mutations with similar efficiency and analytical sensitivity (Figure 1D). Together, these data suggest an assay limit of detection at or just below 0.01% dilution or 2000 copies of mutant transcripts.

The assay performance was next compared with that of standard methods based on PCR amplification and size fractionation analysis. The cDNA products from the same RT reactions were either analyzed with the Signature procedure or amplified by PCR using a pair of primers flanking the *NPM1* exon 12 mutation site followed by amplicon detection by capillary electrophoresis (Figure 2A). Although *NPM1* mutant transcripts were detected in all specimens with the Signature assay, the presence of a robust peak specific for *NPM1* mutant PCR products could only be detected at more than 10 ng of OCI-AML3 total RNA in a background of 1000 ng of HL-60 RNA or 1% equivalent (Figure 2B, top right). Likewise, at least 200,000 copies of synthetic mut-A transcript had to be spiked in HL-60 total RNA for positive *NPM1* mutant detection using the sizing assay (Figure 2B, bottom right). Similar results were obtained for the other synthetic *NPM1* mutant transcripts (Supplemental Figure 2, see <http://jmd.amjpathol.org>).

Evaluation of Clinical Specimens

The Signature *NPM1* Mutations assay was evaluated using 69 archived total RNA samples isolated from peripheral blood or bone marrow from leukemia patients at initial diagnosis (Supplemental Table 1 at <http://jmd.amjpathol.org>). The clinical set consisted of 48 AML specimens: 17 with normal karyotype (NK), 2 positive for *inv(16)*, 4 positive for *t(8;21)*, and 25 with various chromosomal abnormalities other than *inv(16)* or *t(8;21)*, in-

Table 1. Summary of Clinical Results

	Signature <i>NPM1</i> result	
	Positive (Mutant)	Negative (WT)
AML normal karyotype ($n = 17$)	9	8
AML <i>t(8;21)</i> or <i>inv(16)</i> ($n = 6$)	0	6
Other AML ($n = 25$)	1	24
Non-AML ($n = 21$)	0	21
Total ($n = 69$)	10	59
Mean <i>NPM1</i> mutant signal (MFI)	9849	42
Minimum <i>NPM1</i> mutant signal (MFI)	8014	0
Maximum <i>NPM1</i> mutant signal (MFI)	11,416	87
Mean RNA concentration (ng/ μ l)	239	257
Minimum RNA concentration (ng/ μ l)	54	8
Maximum RNA concentration (ng/ μ l)	732	825

WT, wild-type.

cluding 2 with acute bilineage leukemia. In addition, 21 non-AML leukemia samples consisting of 9 acute lymphoid leukemia, 9 chronic myeloid leukemia and 3 myelodysplastic syndrome were included. Approximately 53% (9 of 17) of the AML NK specimens were positive for *NPM1* mutation with the Signature assay (Table 1). Of the 52 other specimens, only 1 AML specimen with a deletion of chromosome 7q22 was positive for *NPM1* mutation. All positive signals were at least 20-fold above the cutoff (400 MFI) for the various RNA inputs tested (Table 1). In contrast, the *NPM1* mutant probe signal was well below the cutoff for the 59 negative specimens evaluated, with a mean value of 42 MFI and a SD of 26 MFI. Mixing experiments with representative clinical specimens confirmed that an analytical sensitivity of at least 0.01% could be reached by dilution of *NPM1*-positive RNA in a background of *NPM1*-negative RNA (Supplemental Figure 3, see <http://jmd.amjpathol.org>).

The same specimen set was also independently analyzed with the capillary electrophoresis assay described above (Table 2). There was 100% agreement between the two methods (95% confidence intervals, 94.7 to 100%). The 10 positive specimens and approximately 15% of the negative specimen set were further confirmed by bidirectional sequencing (Table 2). The four *NPM1* mutant transcripts A, B, D, and J were detected among




Figure 2. Comparison with a size-based assay. **A:** Summary of the positive (Pos) and negative (Neg) calls obtained with serial dilutions of a mutant A-positive cell line RNA (OCI-AML3) or *in vitro* transcript (mut-A) in wild-type cell line RNA (HL-60). A Pos/Neg call indicates that the results were not reproducible (below or at the limit of detection of the method). **B:** Representative examples of capillary electrophoresis traces. The wild-type (WT)- and mutant (M)-specific peaks are indicated. The small peaks at the beginning and end of each trace correspond to molecular weight markers.

Table 2. Method Comparison

	Signature NPM1 result	
	Positive (Mutant)	Negative (WT)
Positive by CE	10	0
Negative by CE	0	59
No. cases sequenced	10	9
Sequencing results	6A, 1B, 1D, 2J	9 wild type

WT, wild-type.

the *NPM1*-positive specimens, with 60% (6 of 10) of the specimens containing mutation A.

A protocol based on the Signature NPM1 Mutations assay and using additional Signature reagents was also developed to combine detection of *NPM1* mutations with detection of AML/ETO and *inv(16)* type A, D, or E fusion transcripts in a single multiplex assay (see *Materials and Methods*). Serial dilution of total RNA purified from OCI-AML3 or translocation-positive cell lines in a background of HL-60 total RNA showed that the individual *NPM1* assay and the combined assay are specific for *NPM1* mutations, AML/ETO and *inv(16)* targets with an analytical sensitivity of at least 0.1% for each target (Supplemental Figure 4, see <http://jmd.amjpathol.org>). Although the *NPM1* mutant signal was reduced by up to 30%, the combined assay still detected 0.01% or 2000 copies of *NPM1* mutant transcripts (Table 3). A subset of the clinical set, namely six AML NK, six AML positive for *inv(16)* or t(8:21), seven other AML, and seven non-AML, was also evaluated with this laboratory-developed test. All of the expected targets were specifically detected in the 26 specimens with 100% agreement between Signature and CE/sequencing for *NPM1* status and 100% agreement between Signature and fluorescence in situ hybridization/sequencing for fusion transcript detection (Table 4 and Supplemental Table 1, see <http://jmd.amjpathol.org>). Mixing experiments with clinical specimens also demonstrated that both *NPM1* mutations and AML fusion transcripts could be simultaneously detected in the same RNA sample (data not shown).

Discussion

AML with mutated *NPM1* is a leukemia entity with distinct pathological, molecular and prognostic features that was recently included in the 2008 World Health Organization classification of myeloid neoplasms. With a high frequency in translocation-negative AML (50 to 60%), *NPM1* mutations are also useful markers for the monitoring of minimal residual disease. However, because of the broad spectrum and nature of each mutation (small repeated or inserted sequences with almost 40 variants identified to date), sensitive molecular detection of multiple *NPM1* mutations in a single assay can be challenging. There is currently a need for improved methods with appropriate mutation coverage but without compromising analytical performance or laboratory efficiency.

The Signature technology is a 96-well assay format that enables rapid analysis of several RNA targets in a single reaction with automatic detection by flow cytometry. This step is performed using the Luminex 100 IS or 200 systems, a platform already validated in more than 4700 laboratories worldwide and compatible with a broad assay menu including 42 US Food and Drug Administration cleared tests (<http://www.luminexcorp.com/company/ceo.html> last accessed March 8, 2010). The streamlined and multiplexed workflow of the *NPM1* assay provides detection of an endogenous control, *NPM1* wild-type, and four *NPM1* mutations: A, B, D, and J (Figure 1). Depending on the study report, these four mutations represent at least 90 to 95% of all *NPM1*-positive cases.^{15,16} Although other synthetic mutant sequences were not tested in our study, detection of additional mutations with the Signature assay had not been reported at the time this manuscript was submitted (F. Ye and E. Labourier, unpublished results). Another advantage of the Signature technology platform is that the multiplexing can be increased using beads with a different identity, up to 100, for detection of the various analytes. For example, co-detection of *NPM1* mutations and *NPM1* wild-type as an internal control together with common AML-specific fusion transcripts in a single assay (Table 3) may be beneficial for the rapid

Table 3. Comparison of Analytical Sensitivity for the Signature NPM1 Mutations Assay Alone or Combined with the Signature AML/ETO, *inv(16)*-ADE Primers and Probes

	Signature NPM1 (MFI)*		Signature NPM1 + Signature AML/ETO, <i>inv(16)</i> ADE (MFI)*					
	NPM1-WT	NPM1-Mut	NPM1-WT	NPM1-Mut	Inv16-A	Inv16-D	Inv16-E	AML1/ETO
No RNA	87	87	42	81	86	64	43	57
100% OCI-AML3	13,019	12,604	9745	10,719	110	67	83	64
HL-60 + 10% OCI-AML3	11,523	11,839	9745	10,263	65	44	77	92
HL-60 + 1% OCI-AML3	14,023	9097	10,579	7623	134	78	52	37
HL-60 + 0.1% OCI-AML3	14,025	4455	9881	3408	25	73	30	36
HL-60 + 0.01% OCI-AML3	13,745	1088	10,166	785	159	81	81	74
100% HL-60	12,465	159	9236	77	2	66	24	17
HL-60 + 2E+07 copies	10,980	12,416	10,388	11,638	97	87	102	144
HL-60 + 2E+06 copies	10,245	10,718	9108	10,432	115	134	127	97
HL-60 + 2E+05 copies	11,593	8566	10,564	8151	83	87	100	57
HL-60 + 2E+04 copies	11,511	4158	10,044	3033	112	79	140	82
HL-60 + 2E+03 copies	10,588	866	10,129	664	107	122	77	97
HL-60 + 2E+02 copies	11,287	322	10,352	190	88	104	78	91

*Values above the cutoff are in bold.
 Mut, mutant; WT, wild-type.

Table 4. Codetection of *NPM1* Mutations and AML-Specific Fusion Transcripts in Seven Representative Clinical Specimens

	MFI signal*					
	NPM1-WT	NPM1-Mut	Inv16-A	Inv16-D	Inv16-E	AML1/ETO
NPM1 positive [†]	12,567	11,284	39	53	18	67
NPM1 negative [†]	10,415	107	35	79	59	33
No RNA [†]	36	50	28	29	63	50
AML-NK, 645 ng [‡]	13,094	11,095	29	48	41	57
AML-NK, 2000 ng	9800	64	48	85	51	16
AML t(8;21), 400 ng [‡]	3269	60	14	99	58	8210
AML t(7;11), 2330 ng	8367	69	89	68	61	36
AML del(7), 268 ng [‡]	11,856	11,240	50	65	86	29
CML t(9;22), 1815 ng	8560	84	30	61	84	60
ALL t(9;22), 1210 ng	9035	59	84	38	73	41

Mut, mutant; WT, wild-type.

*Obtained by combination of the Signature NPM1 Mutations assay with the Signature AML/ETO, inv(16)-ADE Primers and Probes (see *Materials and Methods*). Values above the cutoff are in bold.

[†]Controls provided with the Signature NPM1 Mutations kit.

[‡]These specimens were confirmed to be positive by bidirectional sequencing.

molecular characterization of patients with newly diagnosed AML.

An input range study showed that the Signature NPM1 assay is compatible with a broad range of total RNA purified from the OCI-AML3 cell line (Supplemental Figure 1A, see <http://jmd.amjpathol.org>). The saturation of mutant signal and lack of wild-type signal interference observed above a 10-ng input of OCI-AML3 RNA, as well as above 10% of OCI-AML3 RNA dilution or above 200,000 copies of mut-A synthetic transcript (Figure 1 and Supplemental Figure 1, see <http://jmd.amjpathol.org>), suggest that the assay is extremely robust. Indeed, evaluation of 69 clinical specimens using 5 μ l of total RNA corresponded to a broad range of input RNA (40 to 4000 ng per RT reaction) and resulted in 100% agreement with reference methods. Although a single commercially available RNA isolation method was used, independent studies showed that the Signature NPM1 assay is also compatible with other affinity column- or phenol-based extraction methods (F. Ye and E. Labourier, unpublished results). In addition, a recent multisite validation study to determine *FLT3* and *NPM1* mutation status in total nucleic acid samples (DNA and RNA) extracted with the Qiagen BioRobot EZ1 workstation further confirmed the robustness of the Signature assay and its compatibility with a wide range of clinical specimens (G. J. Tsongalis, Ph.D., personal communication).

In the present study, the specific mutant signals obtained with 10 positive clinical specimens and RNA inputs between 270 and 1200 ng were all robust and on average 25-fold above the cutoff set at 400 MFI (Table 1). Although the number of specimens tested is too low to fully assess clinical sensitivity (95% confidence intervals, 72.2 to 100%), we obtained 100% overall agreement with two independent reference methods (Table 2). The proportion of *NPM1* mutant-positive cases within the AML NK population was as expected (9 of 17 or 53%). A recent large series of patients with *NPM1* mutations showed that approximately 15% carried a karyotypic abnormality¹⁷; we observed a single patient with mutated *NPM1* and del(7q22). The overlapping biological, pathological, immunophenotypic, and clinical features of such cases,

independent of the presence of additional cytogenetic abnormalities, support the provisional diagnostic entity “AML with *NPM1* mutation” in the World Health Organization Classification.¹⁸ Importantly, the MFI signals obtained with the 59 negative specimens were all very low with a mean value at 42 MFI, almost 10-fold below the 400 MFI cutoff (Table 1). The tight distribution of these negative signals (SD = 26 MFI) is also indicative of a robust clinical specificity (95% confidence intervals, 93.9 to 100%) at that cutoff. Based on the normal-like distribution of negative signals in this sample set, a recalculated cutoff at approximately 200 MFI (mean of negative signals + 6.11 SD) would still correspond to a probability of approximately 1 in 1 billion of obtaining a false positive.

The assay analytical sensitivity was thoroughly investigated using a mutant A-positive cell line and synthetic mutant transcripts corresponding to each of the four targets detected by the assay. Together, these data show that the Signature NPM1 assay specifically detects *NPM1* mutations A, B, D, and J with excellent analytical sensitivity, 2 to 3 orders of magnitude better than a standard capillary electrophoresis (CE)-based PCR method (Figures 1 and 2 and Supplemental Figures 1 and 2, see <http://jmd.amjpathol.org>). In the presence of excess *NPM1* wild-type transcripts the assay could reproducibly detect the equivalent of one positive cell in a background of 10,000 negative cells (0.01%). Recent reports showed that the relative or absolute quantitative monitoring of *NPM1* mutation A expression levels in AML NK is an accurate method to measure treatment response and detect potential relapse.^{19,20} Although the Signature NPM1 Mutations assay is not a truly quantitative assay, analysis of the mutant/wild-type MFI ratio showed a semi-quantitative linear response for all four *NPM1* mutations as the relative copy number of mutant transcripts was decreased by up to 4 log units (Figure 1D and Supplemental Figure 1, see <http://jmd.amjpathol.org>). The excellent analytical sensitivity (~0.01%) and dynamic range (~4 log units) of this multiplex assay warrant further evaluation to determine its potential clinical utility for monitoring of residual disease during treatment of patients with AML.

References

1. Falini B, Mecucci C, Tiacci E, Alcalay M, Rosati R, Pasqualucci L, La Starza R, Diverio D, Colombo E, Santucci A, Bigerna B, Pacini R, Pucciarini A, Liso A, Vignetti M, Fazi P, Meani N, Pettrossi V, Saglio G, Mandelli F, Lo-Coco F, Pelicci PG, Martelli MF, GIMEMA Acute Leukemia Working Party: Cytoplasmic nucleophosmin in acute myelogenous leukemia with a normal karyotype. *N Engl J Med* 2005, 352:254–266
2. Suzuki T, Kiyoi H, Ozeki K, Tomita A, Yamaji S, Suzuki R, Koderu Y, Miyawaki S, Asou N, Kuriyama K, Yagasaki F, Shimazaki C, Akiyama H, Nishimura M, Motoji T, Shinagawa K, Takeshita A, Ueda R, Kinoshita T, Emi N, Naoe T: Clinical characteristics and prognostic implications of NPM1 mutations in acute myeloid leukemia. *Blood* 2005, 106:2854–2861
3. Thiede C, Creutzig E, Illmer T, Schaich M, Heise V, Ehninger G, Landt O: Rapid and sensitive typing of NPM1 mutations using LNA-mediated PCR clamping. *Leukemia* 2006, 20:1897–1899
4. Szebeni A, Olson MO: Nucleolar protein B23 has molecular chaperone activities. *Protein Sci* 1999, 8:905–912
5. Major nucleolar proteins shuttle between nucleus and cytoplasm. *Cell* 1989, 56:379
6. Okuda M, Horn HF, Tarapore P, Tokuyama Y, Smulian AG, Chan PK, Knudsen ES, Hofmann IA, Snyder JD, Bove KE, Fukasawa K: Nucleophosmin/B23 is a target of CDK2/cyclin E in centrosome duplication. *Cell* 2000, 103:127–140
7. Colombo E, Marine JC, Danovi D, Falini B, Pelicci PG: Nucleophosmin regulates the stability and transcriptional activity of p53. *Nat Cell Biol* 2002, 4:529–533
8. Falini B, Bolli N, Shan J, Martelli MP, Liso A, Pucciarini A, Bigerna B, Pasqualucci L, Mannucci R, Rosati R, Gorello P, Diverio D, Roti G, Tiacci E, Cazzaniga G, Biondi A, Schnittger S, Haferlach T, Hiddemann W, Martelli MF, Gu W, Mecucci C, Nicoletti I: Both the carboxyl terminus NES motif and mutated tryptophan(s) are crucial for aberrant nuclear export of nucleophosmin leukemic mutants in NPMc⁺ AML. *Blood* 2006, 107:4514–4523
9. Schnittger S, Schoch C, Kern W, Mecucci C, Tschulik C, Martelli MF, Haferlach T, Hiddemann W, Falini B: Nucleophosmin gene mutations are predictors of favorable prognosis in acute myelogenous leukemia with a normal karyotype. *Blood* 2005, 106:3733–3739
10. Verhaak RG, Goudswaard CS, van Putten W, Bijl MA, Sanders MA, Hagens W, Uitterlinden AG, Erpelinck CA, Delwel R, Lowenberg B, Valk PJ: Mutations in nucleophosmin (NPM1) in acute myeloid leukemia (AML): association with other gene abnormalities and previously established gene expression signatures and their favorable prognostic significance. *Blood* 2005, 106:3747–3754
11. Dohner K, Schlenk RF, Habdank M, Scholl C, Rucker FG, Corbacioglu A, Bullinger L, Frohling S, Dohner H: Mutant nucleophosmin (NPM1) predicts favorable prognosis in younger adults with acute myeloid leukemia and normal cytogenetics: interaction with other gene mutations. *Blood* 2005, 106:3740–3746
12. Schlenk RF, Dohner K, Krauter J, Frohling S, Corbacioglu A, Bullinger L, Habdank M, Spath D, Morgan M, Benner A, Schlegelberger B, Heil G, Ganser A, Dohner H, German-Austrian Acute Myeloid Leukemia Study Group: Mutations and treatment outcome in cytogenetically normal acute myeloid leukemia. *N Engl J Med* 2008, 358:1909–1918
13. Boissel N, Renneville A, Biggio V, Philippe N, Thomas X, Cayuela JM, Terre C, Tigaud I, Castaigne S, Raffoux E, De Botton S, Fenaux P, Dombret H, Preudhomme C: Prevalence, clinical profile, and prognosis of NPM mutations in AML with normal karyotype. *Blood* 2005, 106:3618–3620
14. Ammatuna E, Noguera NI, Zangrilli D, Curzi P, Panetta P, Bencivenga P, Amadori S, Federici G, Lo-Coco F: Rapid detection of nucleophosmin (NPM1) mutations in acute myeloid leukemia by denaturing HPLC. *Clin Chem* 2005, 51:2165–2167
15. Falini B, Nicoletti I, Martelli MF, Mecucci C: Acute myeloid leukemia carrying cytoplasmic/mutated nucleophosmin (NPMc⁺ AML): biologic and clinical features. *Blood* 2007, 109:874–885
16. Falini B, Sportoletti P, Martelli MP: Acute myeloid leukemia with mutated NPM1: diagnosis, prognosis and therapeutic perspectives. *Curr Opin Oncol* 2009, 21:573–581
17. Haferlach C, Mecucci C, Schnittger S, Kohlmann A, Mancini M, Cuneo A, Testoni N, Rege-Cambrin G, Santucci A, Vignetti M, Fazi P, Martelli MP, Haferlach T, Falini B: AML with mutated NPM1 carrying a normal or aberrant karyotype show overlapping biologic, pathologic, immunophenotypic, and prognostic features. *Blood* 2009, 114:3024–3032
18. Swerdlow SH, Campo E, Harris NL, Jaffe ES, Pileri SA, Stein H, Thiele J, Vardiman JW (Editors): WHO Classification of Tumours of Haematopoietic and Lymphoid Tissues. Lyon, IARC, 2008, pp. 120–121
19. Papadaki C, Dufour A, Seibl M, Schneider S, Bohlander SK, Zellmeier E, Mellert G, Hiddemann W, Spiekermann K: Monitoring minimal residual disease in acute myeloid leukaemia with NPM1 mutations by quantitative PCR: clonal evolution is a limiting factor. *Br J Haematol* 2009, 144:517–523
20. Bacher U, Badbaran A, Fehse B, Zabelina T, Zander AR, Kroger N: Quantitative monitoring of NPM1 mutations provides a valid minimal residual disease parameter following allogeneic stem cell transplantation. *Exp Hematol* 2009, 37:135–142