

Molecular Probes for Diagnosis of Clinically Relevant Bacterial Infections in Blood Cultures[∇]

Wendy L. J. Hansen, Judith Beuving, Cathrien A. Bruggeman, and Petra F. G. Wolffs*

Department of Medical Microbiology, Maastricht University Medical Center, and the Care and Public Health Research Institute (CAPHRI), Maastricht, Netherlands

Received 16 March 2010/Returned for modification 9 June 2010/Accepted 11 October 2010

Broad-range real-time PCR and sequencing of the 16S rRNA gene region is a widely known method for the detection and identification of bacteria in clinical samples. However, because of the need for sequencing, such identification of bacteria is time-consuming. The aim of our study was to develop a more rapid 16S real-time PCR-based identification assay using species- or genus-specific probes. The Gram-negative bacteria were divided into *Pseudomonas* species, *Pseudomonas aeruginosa*, *Escherichia coli*, and other Gram-negative species. Within the Gram-positive species, probes were designed for *Staphylococcus* species, *Staphylococcus aureus*, *Enterococcus* species, *Streptococcus* species, and *Streptococcus pneumoniae*. The assay also included a universal probe within the 16S rRNA gene region for the detection of all bacterial DNA. The assay was evaluated with a collection of 248 blood cultures. In this study, the universal probe and the probes targeting *Pseudomonas* spp., *P. aeruginosa*, *E. coli*, *Streptococcus* spp., *S. pneumoniae*, *Enterococcus* spp., and *Staphylococcus* spp. all had a sensitivity and specificity of 100%. The probe specific for *S. aureus* showed eight discrepancies, resulting in a sensitivity of 100% and a specificity of 93%. These data showed high agreement between conventional testing and our novel real-time PCR assay. Furthermore, this assay significantly reduced the time needed for identification. In conclusion, using pathogen-specific probes offers a faster alternative for pathogen detection and could improve the diagnosis of bloodstream infections.

Bloodstream infections are a major cause of death in the world and need a thorough and adequate therapeutic strategy. Inadequate antibiotic therapy is associated with higher mortality rates, the appearance of antibiotic resistance, and longer hospitalization lengths (14). Conventional identification and susceptibility testing have several limitations, such as lack of rapidity and sensitivity. The current gold standard, i.e., blood culture, usually requires 6 to 12 h of incubation before growth is detected and a further 24 to 48 h for the definitive identification of the infectious agent and its susceptibility to antibiotics (1, 9). Routine diagnostics already use molecular techniques for the direct detection of viral and bacterial pathogens. However, most in-house assays are targeted against one specific bacterium and/or virus and do not offer broad-range pathogen detection. Recently, several PCR assays have been developed targeting a panel of the most relevant bacterial and fungal bloodstream pathogens, which can be performed directly with blood, such as SeptiFast (Roche Diagnostics GmbH, Mannheim, Germany), SepsiTTest (Molzym GmbH & Co. KG, Bremen, Germany), and VYOO (SIRS-Lab GmbH, Jena, Germany), or using positive blood cultures, such as the microarray-based system Prove-it Sepsis (Mobidiag, Helsinki, Finland).

As discussed in our previous work, direct detection in whole blood is hampered by several factors, such as the presence of PCR inhibitors and background DNA, low bacterial load, in-

sufficient sensitivity, and difficulty of establishing an assay capable of detecting a wide range of pathogens (5). In contrast, molecular testing of growth-positive blood cultures does not require highly sensitive assays because of the presence of a high bacterial load. Furthermore, until now, the use of culturing remains essential to determine the microorganism's antimicrobial profile. Therefore, the role of blood cultures remains important for the detection and identification of causative bacterial agents. Molecular testing of blood cultures, possibly in combination with conventional testing, could enable more rapid identification and, consequently, more rapid diagnosis and start of correct therapy. Molecular approaches such as broad-range real-time PCR and sequencing of the 16S rRNA gene region are widely known methods for the detection and identification of bacteria in clinical samples (4, 7, 10, 15). However, because of the need for sequencing, the identification of bacteria is time-consuming.

The aim of our study was to develop a more rapid 16S real-time PCR-based identification assay using species- or genus-specific probes. The assay is particularly intended for the identification of positive blood cultures, for which Gram staining results are known. In this proof-of-concept study, priority was given to the genera or species most frequently found in blood cultures and/or those that could direct the choice of a suitable antibiotic therapy. Therefore, we selected a panel of eight species- or genus-specific probes. The Gram-negative bacteria were divided into *Pseudomonas* spp., *Pseudomonas aeruginosa*, *Escherichia coli*, and other Gram-negative species. Within the Gram-positive species, probes were designed for *Staphylococcus* spp., *Staphylococcus aureus*, *Enterococcus* spp., *Streptococcus* spp., and *Streptococcus pneumoniae*.

* Corresponding author. Mailing address: Department of Medical Microbiology, Maastricht University Medical Center, P. Debyelaan 25, Maastricht 6202 AZ, Netherlands. Phone: 31-43-3876669. Fax: 31-43-3876643. E-mail: p.wolffs@mumc.nl.

[∇] Published ahead of print on 20 October 2010.

Hence, a first indication about the causative microorganism is given after 2 h, while confirmation and precise identification can be achieved with sequencing. Consequently, multiple species can be detected in samples with polymicrobial infections. The present paper reports a retrospective study performed on blood cultures obtained from patients with suspected bloodstream infections. Results of this new multiprobe assay were compared with conventional blood culture findings.

MATERIALS AND METHODS

Clinical samples. A total of 248 blood cultures were collected at the Maastricht University Medical Center (MUMC; Maastricht, Netherlands). All samples were analyzed with standard conventional testing. Blood specimens drawn from patients suspected of having bloodstream infections were incubated in blood culture bottles (Plus Aerobic [product no. 442192] and Plus Anaerobic [product no. 442193]; BD Diagnostic Systems, Sparks, MD) and monitored for microbial growth in the Bactec 9240 automated blood culture device (BD). When growth was detected, Gram staining was performed. A small aliquot of each blood culture (1 ml) was requested for the novel molecular assay. Two separate assays were developed for Gram-negative and Gram-positive bacteria. Hence, further analysis was based on the results of Gram staining.

Conventional bacterial identification. Regarding the Gram-positive cocci (GPCs), to discern *Staphylococcus* spp. from other GPCs, a catalase test was performed by adding 1 colony to a drop of 3% H₂O₂. For the identification of *Staphylococcus* spp., catalase-positive strains were tested for coagulase and DNase production. If both tests were negative, the strain was identified as coagulase-negative *Staphylococcus* (CoNS). To discern *Enterococcus* spp. from other catalase-negative GPCs, bile esculin, Tellur diagnostic tablets (Rosco Diagnostica, Taastrup, Denmark), and an API 20 Strep test (bioMérieux SA, Marcy l'Etoile, France) were used, according to the manufacturer's guidelines. Optochin susceptibility (OXOID) was used to differentiate *Streptococcus pneumoniae* from the other viridans group streptococci, which were further identified by API 20 Strep. In the case of beta-hemolytic streptococci, latex agglutination was performed using the Prolex streptococcal grouping latex kit (product code PL.030; Bio Trading).

Bacterial strains. Reference strains were used to validate the specificity of the assay, including *S. aureus* ATCC 25923 and ATCC 29213, *Staphylococcus epidermidis* ATCC 12228 and ATCC 14990, *S. pneumoniae* ATCC 49619, *P. aeruginosa* ATCC 27853, *Enterococcus faecalis* ATCC 29212, *Enterococcus faecium* BM 4147, and *Escherichia coli* ATCC 35218.

Multiprobe assay. An aliquot (0.1 ml) of blood culture was 1:100 diluted in 0.9% NaCl. Dilutions were centrifuged for 5 min at 12,000 rpm, and the bacterial pellet was resuspended in 100 μ l nuclease-free water. The primers and the universal bacterial TaqMan probe have been described previously (8). The probes for *P. aeruginosa*, *Pseudomonas* spp., *E. coli*, *Staphylococcus* spp., *S. aureus*, *Enterococcus* spp., *Streptococcus* spp., and *S. pneumoniae* were designed by using the BLAST tool and ClustalW software. Multiple sequence alignments were made and are partly shown in Fig. 1. The designed probe sequences are given in Table 1. All primers and probes were tested for specificity and cross-reactions both manually and with use of the NCBI-BLAST tool (<http://www.ncbi.nlm.nih.gov/BLAST>). Primer and probe matrices were performed to determine optimal concentrations. Each test contained a 5- μ l purified sample and a 20- μ l reaction mixture. The reaction mixture contained 12.5 μ l of TaqMan environmental master mix 2.0 (Applied Biosystems, Foster City, CA), 0.9 μ M forward primer, 0.6 μ M reverse primer, and 0.2 μ M each probe. There were three separate reactions. The first reaction included the universal probe and the *P. aeruginosa* probe. The second reaction included the probes targeting *E. coli* and *Pseudomonas* spp. The third reaction included the *Staphylococcus* probe, the *S. aureus* probe, and the *Enterococcus* probe, the fourth and final reaction included the *Streptococcus* probe and the *S. pneumoniae* probe. Reactions were performed using the ABI Prism 7000 real-time PCR system (Applied Biosystems, Foster City, CA), and optimal thermal cycling conditions used were as follows: 2 min at 50°C, initial denaturation at 95°C for 15 min, 42 cycles of denaturation for 15 s at 95°C, and annealing at 60°C for 1 min. The cycle threshold (C_T) value, the cycle number at which amplicon fluorescence exceeded the preset detection threshold, was recorded for all samples. The cutoff value to consider a PCR result as positive was set to a C_T value of 35.

DNA sequencing. Samples with discrepant results were further analyzed with sequencing. The DNA was purified using the MSB PCRapace PCR cleanup kit (Invitex, Berlin, Germany) and dissolved in 30 μ l NASBA H₂O. Cyclic sequencing was performed using BigDye 3.0 (Applied Biosystems, Foster City, CA). The resulting sequence product was purified and separated with the 3730 DNA Analyzer (Applied Biosystems, Foster City, CA).

RESULTS

Design of the species- or genus-specific probes. The specificity of the probes was evaluated using BLAST. *In silico* analysis revealed that nearly only cross-reactivity was found with not clinically relevant microorganisms such as *Geobacillus* spp. and *Lactobacillus plantarum*. *In silico* cross-reactivity with clinically relevant microorganisms was tested *in vitro* with reference strains. For example, the *in silico* *E. coli* probe cross-reacted with *Acinetobacter* spp. and *Enterobacter* spp. However, *in vitro* tests showed no positive signals for these microorganisms. Furthermore, cross-reactions were found with *Escherichia albertii* and *Shigella* spp. Based on the literature, these microorganisms occur only rarely in blood cultures, and therefore, the risk of false identification was minimal. In a few cases such as with the *Streptococcus* and the *Staphylococcus* probes, cross-reactivity with *Bacillus cereus* was found. The species conferring cross-reactivity in PCR were then sequenced to check for mismatches between the probe-target hybrids. Sequencing results showed that this occurred because of at least two mismatches at the binding site of both probes (data not shown). However, in all these cases, the cross-reacting fluorescent signals generated by the mismatched probes were weak compared to that of the positive template control and uncharacteristically appeared over 5 cycles later than signals from 100% matched probes (data not shown).

Evaluation of the multiprobe assay by testing clinical blood culture samples. From a total of 248 blood cultures, the presence or absence of bacterial DNA was determined in 232 growth-positive and 16 growth-negative samples using the universal probe, yielding a sensitivity and specificity of 100% (Table 2). Our multiprobe assay was performed after Gram staining of the positive blood cultures. Considering the growth-positive blood cultures, results from the multiprobe assay were in accordance with conventional identification in 222 cases (96%). The specific probes targeting *Pseudomonas* spp., *P. aeruginosa*, *E. coli*, *Streptococcus* spp., *S. pneumoniae*, and *Staphylococcus* spp. all had a sensitivity and specificity of 100%. The majority of Gram-positive blood cultures, 85 of 145 (59%), contained a coagulase-negative *Staphylococcus* species. In 31 blood cultures, the causative agent was identified as *S. aureus*. Regarding the Gram-positive staphylococcal blood cultures, the multiprobe assay was in conflict with culture results in nine cases. The nine discrepancies were further analyzed by coagulase testing, specific *S. aureus* real-time PCR, and sequencing (Table 3; Fig. 2). One out of nine cases, clinical isolate (CI) 7, tested coagulase positive and was positive for two targets specific for *S. aureus*, *i.e.*, *femA* and *sa442*. Sequencing of the PCR product also confirmed the results of the multiprobe assay. The sequences of the eight other clinical isolates showed three mismatches with the reference *S. aureus* sequence, as shown in Fig. 2A. Hence, the remaining eight cases were confirmed as discordant. Consequently, sensitivity and specificity of 100% and 93%, respectively, and positive and




FIG. 1. Multiple sequence alignments of the amplified part of the bacterial 16S rRNA gene. Sequences of the bacteria included in the multiprobe assay are shown. The arrows represent the forward and reverse primer targeting the 16S rRNA gene. The sequences of the probes are in boldface and underlined. a, *Staphylococcus* probe; b, *Enterococcus* probe; c, *Pseudomonas* probe; d, *S. aureus* probe; e, *E. coli* probe; f, *Streptococcus* probe; g, universal probe; h, *S. pneumoniae* probe; i, *P. aeruginosa* probe.

negative predictive values of 79% and 100%, respectively, were achieved when sequencing was considered the gold standard. The probe specific for *Enterococcus* spp. showed one conflicting result. One blood culture was determined to contain *Aero-*

coccus viridans, while our multiprobe assay identified the infectious agent as an *Enterococcus* sp. (Table 3). Further analysis revealed that the partial sequence of CI 10 was completely similar to the sequence of a reference *Enterococcus*

Streptococcus species	AGGTCCCAGCGTTGTCCGGATTTATTGGGCGTAAAGCGAGCGCAGGCGG	250
Enterococcus species	AGGTGGCAAGCGTTGTCCGGATTTATTGGGCGTAAAGCGAGCGCAGGCGG	250
Staphylococcus aureus	AGGTGGCAAGCGTTATCCGGAATTATTGGGCGTAAAGCGCGCTAGGCGG	249
Staphylococcus species	AGGTGGCAAGCGTTATCCGGAATTATTGGGCGTAAAGCGCGCTAGGCGG	249
Pseudomonas aeruginosa	AGGGTGAAGCGTTAATCGGAATTACTGGGCGTAAAGCGCGCTAGGTGG	249
Pseudomonas species	AGGGTGAAGCGTTAATCGGAATTACTGGGCGTAAAGCGCGCTAGGTGG	249
Escherichia coli	AGGGTGAAGCGTTAATCGGAATTACTGGGCGTAAAGCGCACGCAGGCGG *** * ***** ** * ***** ** * ** *	249
Streptococcus pneumoniae	TTAGATAAGTCTGAAGTAAAGGCTG <u>TGGCTTAACCATAGTAGGCT-TTG</u>	299 h
Streptococcus species	TTAGAAAAGTCTGAAGTAAAGGCGAGTGGCTCAACCATTGTAGGCT-TTG	299
Enterococcus species	TTTCTTAAGTCTGATGTGAAAGCCCCGGCTCAACCGGGGAGGGTCATTG	300
Staphylococcus aureus	TTTTTTAAGTCTGATGTGAAAGCCACGGCTCAACCGTGGAGGGTCATTG	299
Staphylococcus species	TTTTTTAAGTCTGATGTGAAAGCCACGGCTCAACCGTGGAGGGTCATTG	299
Pseudomonas aeruginosa	TTCAGCAAGTTGGATGTGAAATCGCCGGGCTCAACCTGGGAAGTGCAT <u>CC</u>	299 i
Pseudomonas species	TTTGTTAAGTTGGATGTGAAATCCCCGGGCTCAACCTGGGAAGTGCATTC	299
Escherichia coli	TTTGTTAAGTCAGATGTGAAATCCCCGGGCTCAACCTGGGAAGTGCATCT ** ***** ** * ** * ***** ** * *	299
Streptococcus pneumoniae	<u>GAAACTGTTTAACTTGAGTGCAGAAGAGGGAGAGTGAAT-CCATGTGTAG</u>	348 h
Streptococcus species	GAAACTGTTTAACTTGAGTGCAGAAGAGGGAGAGTGAATTCATGTGTAG	349
Enterococcus species	GAAACTGGGAGACTTGAGTGCAGAAGAGGAGAGTGAATTCATGTGTAG	350
Staphylococcus aureus	GAAACTGGAAAAGTCTGAGTGCAGAAGAGGAAAGTGAATTCATGTGTAG	349
Staphylococcus species	GAAACTGGAAAAGTCTGAGTGCAGAAGAGGAAAGTGAATTCATGTGTAG	349
Pseudomonas aeruginosa	<u>AAAAGTACTGAGCTAGAGTACG</u> GTAGAGGGTGGTGAATT-CCTGTGTAG	348 i
Pseudomonas species	AAAAGTACTGACTAGAGTATGGTAGAGGGTGGTGAATTCCTGTGTAG	349
Escherichia coli	GATACTGGCAAGCTTGAGTCTCGTAGAGGGGGTGAATTCAGGTGTAG * ** * ** * ** * ** * ** * ** * ** * ** * ** *	349
Streptococcus pneumoniae	CGGTGAAATGCGTAGATATATGGAGGAACACCGGTGGCGAAAGCGGCTCT	398
Streptococcus species	CGGTGAAATGCGTAGATATATGGAGGAACACCGGTGGCGAAAGCGGCTCT	399
Enterococcus species	CGGTGAAATGCGTAGATATATGGAGGAACACCGTGGCGAAGGCGGCTCT	400
Staphylococcus aureus	CGGTGAAATGCGCAGAGATATGGAGGAACACCGTGGCGAAGGCGACTTT	399
Staphylococcus species	CGGTGAAATGCGCAGAGATATGGAGGAACACCGTGGCGAAGGCGGCTTT	399
Pseudomonas aeruginosa	CGGTGAAATGCGTAGATATAGGAAGGAACACCGTGGCGAAGGCGACCAC	398
Pseudomonas species	CGGTGAAATGCGTAGATATAGGAAGGAACACCGTGGCGAAGGCGACCAC	399
Escherichia coli	CGGTGAAATGCGTAGAGATCTGGAGGAATACCGTGGCGAAGGCGGCC ***** ** * ** * ***** ** * ** *	399
Streptococcus pneumoniae	CTGGCTTGTAAGTACGCTGAGGCTCGAAAGCGTGGGAGCA <u>AAACAGGAT</u>	448
Streptococcus species	CTGGCTGTCACTGACGCTGAGGCTCGAAAGCGTGGGAGCG <u>AAACAGGAT</u>	449

FIG. 1—Continued.

strain, as shown in Fig. 2B. Comparison of the reference *A. viridans* and *Enterococcus* strains showed five mismatches. Consequently, the probe specific for *Enterococcus* spp. yielded a sensitivity and specificity of 100% when sequencing was considered the gold standard.

DISCUSSION

New strategies for the detection and identification of bacterial pathogens in blood are continuously under investigation worldwide. Current diagnostic tools are hampered by several

TABLE 1. Probes designed for use in multiplex PCR

Probe	Sequence (5'-3') ^a
<i>Pseudomonas</i> species.....	NED-CCTTCCTCCCAACTTAAAGTGCTT-MGB
<i>P. aeruginosa</i>	JOE-CCAAAACACTACTGAGCTAGAGTACG-BHQ1
<i>E. coli</i>	JOE-GGAGTAAAGTTAATACCTTTGCTCATT-BHQ1
<i>Staphylococcus</i> species	NED-AATCTTCCGCAATGGGCGAAAAGC-MGB
<i>S. aureus</i>	FAM-AGATGTGCACAGTTACTTTACACATAT-BHQ1
<i>Enterococcus</i> species	JOE-TCCTTGTCTTCTCTAACAACAGAG-BHQ1
<i>Streptococcus</i> species.....	NED-CCAGAAAGGGACSGTAACT-MGB
<i>S. pneumoniae</i>	JOE-CCAAAAGCCTACTATGGTTAAGCCA-BHQ1

^a NED, fluorescent label (Applied Biosystems); MGB, minor groove binder; JOE, 6-carboxy-4,5-dichloro-2,7-dimethoxyfluorescein; BHQ1, black hole quencher 1; FAM, 6-carboxyfluorescein.

factors, such as a lack of sensitivity, a long time to results, the presence of PCR inhibitory compounds in certain sample materials, and the occurrence of fastidious and nonculturable microorganisms. Rapid universal 16S rRNA gene-based PCR

is often combined with sequence analysis, and a considerable number of commercial assays have been generated and tested recently, i.e., SeptiFast (Roche Diagnostics GmbH, Mannheim, Germany), Prove-it Sepsis (Mobidiag, Helsinki,

TABLE 2. Bacterial isolates in 232 positive blood cultures and results from universal and specific real-time PCR

Pathogen	Result of blood culture (no. [%]) ^a								
	UNI	PSEU	PSEUAE	ECOLI	STAPH	STAU	ENTE	STREPT	STREPN
Gram-negative pathogens									
<i>Acinetobacter lwoffii/haemolyticus</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Bacteroides fragilis</i>	2 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Citrobacter koseri</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Enterobacter cloacae</i>	3 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Escherichia coli</i>	43 (100)	0	0	43	NA	NA	NA	NA	NA
Gram-negative rod	2 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Klebsiella oxytoca</i>	5 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Klebsiella pneumoniae</i>	7 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Moraxella catarrhalis</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Morganella morganii</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Neisseria meningitidis</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Prevotella buccae</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Proteus mirabilis</i>	2 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Pseudomonas aeruginosa</i>	10 (100)	10	10	0	NA	NA	NA	NA	NA
<i>Pseudomonas oryzihabitans</i>	1 (100)	1	0	0	NA	NA	NA	NA	NA
<i>Serratia marcescens</i>	4 (100)	0	0	0	NA	NA	NA	NA	NA
Gram-positive pathogens									
<i>Aerococcus viridans</i>	1 (100)	NA	NA	NA	0	0	1	0	0
Coagulase-negative <i>Staphylococcus</i> spp.	85 (100)	NA	NA	NA	85	9	0	0	0
<i>Corynebacterium</i> spp.	2 (100)	NA	NA	NA	0	0	0	0	0
<i>Enterococcus avium</i>	1 (100)	NA	NA	NA	0	0	1	0	0
<i>Enterococcus faecalis</i>	3 (100)	NA	NA	NA	0	0	3	0	0
<i>Enterococcus faecium</i>	4 (100)	NA	NA	NA	0	0	4	0	0
<i>Lactobacillus</i> spp.	1 (100)	NA	NA	NA	0	0	0	0	0
<i>Propionibacterium acnes</i>	1 (100)	NA	NA	NA	0	0	0	0	0
<i>Staphylococcus aureus</i>	31 (100)	NA	NA	NA	31	31	0	0	0
<i>Streptococcus agalactiae</i>	1 (100)	NA	NA	NA	0	0	0	1	0
<i>Streptococcus</i> group A and/or D	2 (100)	NA	NA	NA	0	0	0	2	0
<i>Streptococcus milleri</i>	3 (100)	NA	NA	NA	0	0	0	3	0
<i>Streptococcus oralis</i>	1 (100)	NA	NA	NA	0	0	0	1	0
<i>Streptococcus pneumoniae</i>	6 (100)	NA	NA	NA	0	0	0	6	6
<i>Streptococcus pyogenes</i>	1 (100)	NA	NA	NA	0	0	0	1	0
<i>Streptococcus sanguis</i>	1 (100)	NA	NA	NA	0	0	0	1	0
Mixed									
<i>Alcaligenes faecalis</i> + <i>Propionibacterium</i> spp.	1 (100)	0	0	0	0	0	0	0	0
<i>Klebsiella oxytoca</i> + <i>Serratia marcescens</i>	1 (100)	0	0	0	NA	NA	NA	NA	NA
<i>Streptococcus salivarius</i> + <i>Streptococcus viridans</i>	1 (100)	NA	NA	NA	0	0	0	1	0
Totals									
Concordant		11	10	43	116	31	8	16	6
Discordant		0	0	0	0	9	1	0	0

^a UNI, universal probe; PSEU, *Pseudomonas* spp.; PSEUAE, *P. aeruginosa*; ECOLI, *E. coli*; STAPH, *Staphylococcus* spp.; STAU, *S. aureus*; ENTE, *Enterococcus* spp.; STREPT, *Streptococcus* spp.; STREPN, *S. pneumoniae*; NA, not applicable.

TABLE 3. Clinical isolates with discordant blood culture and PCR results

CI ^a	Specific probe(s)	Result of blood culture ^a		
		<i>femA/sa442</i>	Coagulase activity	Sequence
1	STAPH, STAU	N/N	N	Unidentified
2	STAPH, STAU	N/N	N	<i>Staphylococcus hominis</i>
3	STAPH, STAU	N/N	N	<i>S. hominis</i>
4	STAPH, STAU	N/N	N	<i>Staphylococcus pasteurii</i>
5	STAPH, STAU	N/N	N	<i>S. hominis</i>
6	STAPH, STAU	N/N	N	<i>Staphylococcus schleiferi</i>
7	STAPH, STAU	P/P	P	<i>S. aureus</i>
8	STAPH, STAU	N/N	N	<i>S. pasteurii</i>
9	STAPH, STAU	N/N	N	<i>S. epidermidis</i>
10	ENTE	NA	NA	<i>Enterococcus</i> spp.

^a Clinical isolates (CI) 1 to 9 are CoNS species, while CI 10 is *A. viridans*. STAPH, *Staphylococcus* spp.; STAU, *S. aureus*; ENTE, *Enterococcus* spp; N, negative; P, positive; NA, not applicable.

Finland), and SepsiT_{est} (Molz_{ym} GmbH & Co. KG, Bremen, Germany). The latter includes a universal 16S rRNA gene-based PCR assay combined with sequence analysis, while the two other commercial assays use a panel of probes for the detection of a range of bacterial and mycotic pathogens. Virtually all studies conducted with the SeptiFast and SepsiT_{est} assays found similar results, i.e., an overall agreement of 77 to 83% (2, 3, 6, 12, 13, 15, 16). Both assays are validated for whole-blood samples and can be completed in 4 to 6 h. Tissari et al. performed an observational study comparing conventional culture with Prove-it Sepsis. The assay, based on a microarray platform and performed with blood cultures, had a clinical sensitivity of 94.7% and a specificity of 98.8% (11).

In our study, the 16S rRNA gene was used for the design of eight species- or genus-specific probes. The sequence similarities found in *in silico* analysis (BLAST) were derived mainly from not clinically relevant microorganisms such as *Geobacil-*




FIG. 2. Partial sequence alignments of the 16S rRNA gene, derived from the clinical isolates with discordant results. (A) Region where the probe specific for *S. aureus* binds. The sequences of the eight discrepancies are shown in the rectangle. The nucleotides underlined in the rectangle represent the nucleotides which distinguish *S. aureus* from *Staphylococcus* (*Staph*) spp. The sequence alignment of clinical isolate (CI) 7 is the exact reverse complement of the *S. aureus* probe sequence, which implicates that this isolate, which was identified as CoNS, is in fact *S. aureus*. CI 1 was not shown because the sequence could not be determined. GenBank accession numbers of reference strains are as follows: *S. aureus*, FN433596; *Staphylococcus* spp., GQ222246; and *S. epidermidis*, GQ911565. (B) Region where the probe specific for *Enterococcus* spp. binds. The nucleotides underlined in the rectangle represent the nucleotides which distinguish *Enterococcus* spp. from *A. viridans*. Conventional testing identified CI 10 as *Aerococcus viridans*, whereas our assay and sequencing classified it as an *Enterococcus* sp. GenBank accession numbers of reference strains are as follows: *Enterococcus* spp., AB489105; *A. viridans*, M58797. *, reference sequences.

lus spp., *Paenibacillus* spp., *L. plantarum*, and *Lactobacillus fermentum*. Remaining similarities were, if available, tested *in vitro* with reference strains. A limited number of microorganisms showed cross-reactivity. For example, streptococcal and staphylococcal probes weakly cross-reacted with *Bacillus cereus*. In these cases, sequencing determined the correct identification. Additionally, as the assay is always preceded by Gram staining, in most cases, morphology data obtained by Gram staining was in contradiction with these false-positive signals and targeted the samples for confirmatory sequencing.

Our multiprobe assay was designed for use with blood culture material and can be completed within 2 h. Use of whole-blood samples could significantly improve the turnaround time of our assay. However, since culturing still remains essential to determine the microorganism's antimicrobial profile, we chose to use blood cultures instead of whole blood. In this study, we found an overall agreement of 97% between conventional testing and our multiprobe assay. One of the strengths of our approach is that the assay can be extended by adding more probes for other bacterial pathogens to the identification panel. The grouping of Gram-negative and Gram-positive bacteria was based on the identifications needed for antimicrobial susceptibility testing. The assay was designed to use a universal probe for detection instead of Sybr green, because it is known that Sybr green can generate false-positive signals because of the presence of background human DNA or the formation of primer-dimers.

Analysis of the conflicting results showed that, in some cases, our assay was in accordance with sequencing results. The reference sequence derived from a *Staphylococcus* sp. and a *S. epidermidis* strain showed seven and three mismatches, respectively, with the reference sequence derived from an *S. aureus* strain. The sequences derived from the discrepant clinical isolates only showed three or fewer mismatches with the *S. aureus* sequence. This indicated that within the group of staphylococci, more homologous sequences can cause false positives because of the less efficient binding capacity of the *S. aureus*-specific probe.

Overall, these results showed strong agreement between conventional testing and our novel, real-time PCR assay. Furthermore, this assay significantly reduced the time needed for identification in comparison to that of routine diagnostics. In conclusion, using pathogen-specific probes offers a faster alternative for pathogen detection and could improve the diagnosis of bloodstream infections. The assay will be implemented in a clinical trial investigating the impact of earlier pathogen identification in combination with susceptibility testing with the choice of therapy.

REFERENCES

1. Beekmann, S. E., D. J. Diekema, K. C. Chapin, and G. V. Doern. 2003. Effects of rapid detection of bloodstream infections on length of hospitalization and hospital charges. *J. Clin. Microbiol.* **41**:3119–3125.
2. Casalta, J. P., F. Gouriet, V. Roux, F. Thuny, G. Habib, and D. Raoult. 2009. Evaluation of the LightCycler SeptiFast test in the rapid etiologic diagnosis of infectious endocarditis. *Eur. J. Clin. Microbiol. Infect. Dis.* **28**:569–573.
3. Dierkes, C., B. Ehrenstein, S. Siebig, H. J. Linde, U. Reischl, and B. Salzberger. 2009. Clinical impact of a commercially available multiplex PCR system for rapid detection of pathogens in patients with presumed sepsis. *BMC Infect. Dis.* **9**:126.
4. Duan, Q. J., S. Q. Shang, and Y. D. Wu. 2009. Rapid diagnosis of bacterial meningitis in children with fluorescence quantitative polymerase chain reaction amplification in the bacterial 16S rRNA gene. *Eur. J. Pediatr.* **168**:211–216.
5. Hansen, W. L., C. A. Bruggeman, and P. F. Wolffs. 2009. Evaluation of new preanalysis sample treatment tools and DNA isolation protocols to improve bacterial pathogen detection in whole blood. *J. Clin. Microbiol.* **47**:2629–2631.
6. Mancini, N., D. Clerici, R. Diotti, M. Perotti, N. Ghidoli, D. De Marco, B. Pizzorno, T. Emrich, R. Burioni, F. Ciceri, and M. Clementi. 2008. Molecular diagnosis of sepsis in neutropenic patients with haematological malignancies. *J. Med. Microbiol.* **57**:601–604.
7. Mignard, S., and J. P. Flandrois. 2006. 16S rRNA sequencing in routine bacterial identification: a 30-month experiment. *J. Microbiol. Methods* **67**:574–581.
8. Nadkarni, M. A., F. E. Martin, N. A. Jacques, and N. Hunter. 2002. Determination of bacterial load by real-time PCR using a broad-range (universal) probe and primers set. *Microbiology* **148**:257–266.
9. Peters, R. P., M. A. van Agtmael, S. A. Danner, P. H. Savelkoul, and C. M. Vandembroucke-Grauls. 2004. New developments in the diagnosis of bloodstream infections. *Lancet Infect. Dis.* **4**:751–760.
10. Qian, Q., Y. W. Tang, C. P. Kolbert, C. A. Torgerson, J. G. Hughes, E. A. Vetter, W. S. Harmsen, S. O. Montgomery, F. R. Cockerill III, and D. H. Persing. 2001. Direct identification of bacteria from positive blood cultures by amplification and sequencing of the 16S rRNA gene: evaluation of BACTEC 9240 instrument true-positive and false-positive results. *J. Clin. Microbiol.* **39**:3578–3582.
11. Tissari, P., A. Zumla, E. Tarkka, S. Mero, L. Savolainen, M. Vaara, A. Aittakorpi, S. Laakso, M. Lindfors, H. Piiparinen, M. Maki, C. Carder, J. Huggett, and V. Gant. 2010. Accurate and rapid identification of bacterial species from positive blood cultures with a DNA-based microarray platform: an observational study. *Lancet* **375**:224–230.
12. Varani, S., M. Stanzani, M. Paolucci, F. Melchionda, G. Castellani, L. Nardi, M. P. Landini, M. Baccarani, A. Pession, and V. Sambri. 2009. Diagnosis of bloodstream infections in immunocompromised patients by real-time PCR. *J. Infect.* **58**:346–351.
13. von Lilienfeld-Toal, M., L. E. Lehmann, A. D. Raadts, C. Hahn-Ast, K. S. Orloff, G. Marklein, I. Purr, G. Cook, A. Hoeft, A. Glasmacher, and F. Stuber. 2009. Utility of a commercially available multiplex real-time PCR assay to detect bacterial and fungal pathogens in febrile neutropenia. *J. Clin. Microbiol.* **47**:2405–2410.
14. Wallet, F., S. Nseir, L. Baumann, S. Herwegh, B. Sendid, M. Boulo, M. Roussel-Delvallez, A. V. Durocher, and R. J. Courcol. 18 August 2009. Preliminary clinical study using a multiplex real-time PCR test for the detection of bacterial and fungal DNA directly in blood. *Clin. Microbiol. Infect.* doi:10.1111/j.1469-0691.2009.02940.x.
15. Wellinghausen, N., A. J. Kochem, C. Disque, H. Muhl, S. Gebert, J. Winter, J. Matten, and S. G. Sakka. 2009. Diagnosis of bacteremia in whole-blood samples by use of a commercial universal 16S rRNA gene-based PCR and sequence analysis. *J. Clin. Microbiol.* **47**:2759–2765.
16. Westh, H., G. Lisby, F. Breyse, B. Boddington, M. Chomarat, V. Gant, A. Goglio, A. Raglio, H. Schuster, F. Stuber, H. Wissing, and A. Hoeft. 2009. Multiplex real-time PCR and blood culture for identification of bloodstream pathogens in patients with suspected sepsis. *Clin. Microbiol. Infect.* **15**:544–551.