



Published in final edited form as:

Methods Enzymol. 2009 ; 460: 347–355. doi:10.1016/S0076-6879(09)05217-3.

MONITORING NF- κ B MEDIATED CHEMOKINE TRANSCRIPTION IN TUMORIGENESIS

Jinming Yang and Ann J. Richmond

Department of Cancer Biology, and Veterans Affairs Medical Center, Vanderbilt University School of Medicine, Nashville, Tennessee, USA

Abstract

Chemokine-receptor signaling plays an important role in the inflammatory response often associated with tumor growth and metastasis. The NF- κ B pathway, essential for transcription of chemokines/chemokine receptors and other key inflammatory modulators, has emerged as a potential target for tumor therapy. Here we describe an efficient approach to monitor drugs that target the NF- κ B signaling as related to tumor growth and metastasis *in vivo*. For bioluminescence imaging, the firefly luciferase (Fluc) reporter has the advantage of stable signaling, while Gaussia luciferase (Gluc) provides very sensitive signaling based on secretion of Gluc. We introduce the use of monitoring intratumoral Gluc, which rapidly diffuses into the blood circulation and urine. The peripheral Gluc assay may complement bioluminescence imaging and provide a kinetic, noninvasive, real-time read-out of NF- κ B activity by directly determining Gluc reporter activity in blood or urine samples from tumor-bearing mice.

1. INTRODUCTION

Nuclear factor- κ B (NF- κ B), a key signal transduction pathway in chemokine–chemokine receptor expression, inflammation, and cancer, is important target for drug discovery and development. It has been increasingly important to develop noninvasive, high-resolution, *in vivo* imaging to elucidate mechanisms that identify and validate drugs that target the NF- κ B pathway. Recent advances in technology allow visualization of signal transduction as related to biological processes *in vivo* (Phair and Misteli, 2001). Use of fluorescent proteins revolutionizes static microscopy images by providing the ability to make dynamic recordings of protein–protein interaction in living animals. More recently, a Gaussia luciferase that possesses a natural secretory signal, allowing secretion into the cell microenvironment, offers great promise for real-time *ex vivo* monitoring of NF- κ B signaling in tumor development and progression in conscious animals.

2. DEVELOPMENT OF NF- κ B REPORTER MODEL FOR TUMORS

To monitor NF- κ B signaling in tumorigenesis, we constructed a luciferase reporter vector that expresses a transcription factor–mediated reporter protein. To generate the NF- κ B–mediated transcription, firefly-luciferase reporter vector, we inserted a commercial vector with a transcription blocker that is composed of adjacent polyadenylation and transcription pause sites for reducing background transcription. An NF- κ B consensus sequence was fused to a TATA-like promoter, followed by a firefly- or Gaussia-luciferase reporter gene for monitoring NF- κ B signaling (Yang *et al.*, 2007). These vectors are designated as NF- κ B-Fluc or NF- κ B-Gluc, respectively. This system enables signal amplification of a transcription factor–mediated signal.

To create NF- κ B–promoter reporter cells, human melanoma cells (Hs294T) were transfected with the linearized NF- κ B-Gluc vector and stably transfected cells were selected with 2 mg/

ml G418. Gaussia luciferase in the cultured medium or cell lysate was characterized based on the catalytic reaction with its substrate, coelenterazine. We subcutaneously inoculated nude mice with these melanoma NF- κ B-Gluc reporter cells and allowed tumor xenografts to grow over 14 days. The following approaches are used to monitor NF- κ B signaling changes in real time during tumor progression and in response to molecular targeting therapy.

3. BIOLUMINESCENT IMAGING OF INTRATUMOR SIGNALING OF ANESTHETIZED MICE

3.1. Firefly luciferase reporter

Bioluminescence images can be taken using the IVIS 200 Imaging System (Xenogen Imaging Technologies). The system is composed of an imaging chamber, gas anesthesia system, and a highly sensitive charge-coupled device (CCD) camera. The CCD camera is cryogenically cooled for highly efficiency photon detection and displays a wide-range signal image. It takes approximately 3 min to anesthetize small mammals such as mice using the gas anesthesia system that is connected to an oxygen cylinder and isoflurane tank with flow rates set at scale. Fresh luciferin solution is prepared by dissolving luciferin powder (15 mg/ml) in phosphate-buffered saline (PBS) (137 mM NaCl, 2.7 mM KCl, 4.3 mM Na₂HPO₄, 1.47 mM KH₂PO₄, and adjusted to a final pH of 7.4). Mice are injected intraperitoneally with 150 μ g of luciferin per gram body weight and transferred into the image chamber. The image field is set according to the number of mice to be imaged. To acquire the live images, the focus is adjusted to 1.5cm for a subcutaneous (SC) tumor (Fig. 17.1), or to 1 cm for deep organ tumors such as a tumor metastatic to the liver (Fig. 17.2). The exposure time is usually 1 min for SC tumors and 3 min for tumors metastatic to the liver. The peak time of image intensity is around 10 min postinjection of substrate (luciferin), although this time is variable depending on the anatomy of the tumor location. The image file is saved as TIF format and the image intensity is quantitated using the Living Image software 3.0 from Xenogen Imaging Technologies (Fig. 17.1).

3.2. Gaussian luciferase reporter

To produce a bioluminescence image reporter for Gaussia luciferase activity, native coelenterazine (50 to 100 μ g/mouse) injected intravenously gives an accurate read-out of NF- κ B activity (unpublished). Coelenterazine is easily dissolved into either ethanol or methanol. To make up a stock of native coelenterazine, 4 mg coelenterazine is dissolved in 1 ml methanol. This stock solution can be stored at -70 °C for 2 weeks. However, for the most accurate reproducible comparative data, freshly mixed coelenterazine is recommended. The coelenterazine stock should be diluted 1:20 (v/v) with a buffer (150 mmol/l NaCl, 2 mmol/l KCl, 1 mmol/l MgCl₂, 10 mmol/l Na₂HPO₄, 2 mmol/l KH₂PO₄, pH 7.4) prior to intravenous injection. The injected mouse will be immediately subjected to bioluminescence imaging for the reason that the Gaussia luciferase signal rapidly fades relative to firefly luciferase. Thus, Gaussia luciferase imaging allows analysis of only one mouse at a time, while firefly luciferase imaging can be performed on a maximum 5 mice for one time point (Fig. 17.1).

4. CELL-BASED ASSAYS FOR KINASE AND TRANSCRIPTIONAL ACTIVITY *IN VITRO*

In a nonstimulated cell, I κ B α is tightly complexed with NF- κ B to hold NF- κ B in the cytoplasm and keep it in a biologically inactive form. In this manner, I κ B α serves as the brake for the NF- κ B signal transduction cascade. When ligand binding activates receptors at the cell surface, the signal cascade is triggered through activation of I κ B kinase (IKK). The activated IKK phosphorylates I κ B α protein, which enables ubiquitination followed by

degradation of I κ B α and the eventual nuclear translocation of NF- κ B (RelA/p50). The classical method for determining IKK activity has relied on an *in vitro* kinase assay where the IKK complex is immunoprecipitated from cell lysate and the activated IKK catalyzes the transfer of radiolabeled phosphate to a purified I κ B α protein. In contrast, the newly developed I κ B-Gluc reporter plasmid or NF- κ B-Gluc reporter plasmid reflects cellular I κ B α levels or transcriptional NF- κ B activity in intact cells without interruption. When cells are transfected to stably express the reporter plasmid, Gaussia luciferase is naturally secreted into the culture medium. Addition of the IKK inhibitor, BMS-345541 inhibits NF- κ B-Gluc and increase I κ B-Gluc activity. Gaussia luciferase as a reporter accurately reflects NF- κ B activity based on luminescence imaging (Fig. 17.3), or by measuring 10-s photon counts using a luminometer following a reaction of 10 μ l of medium with the substrate (100 micromoles/l coelenterazine) in buffer containing 500 mmol/l NaCl, 2 mmol/l KCl, 10 mmol/l MgCl₂, 10 mmol/l Na₂HPO₄, 2 mmol/l KH₂PO₄, 1 mmol/l EDTA and adjusted to a pH of 7.8.

5. PERIPHERAL SPYING OF INTRATUMORAL SIGNALING OF CONSCIOUS MICE

Bioluminescence imaging studies have shown that Gaussia luciferase as a reporter gives a 200-fold brighter signal than that of firefly luciferase (Tannous *et al.*, 2005). However, the Gaussia luciferase signal is less stable during bioluminescent process than the firefly luciferase luminescence imaging. Thus, Gaussia luciferase is not an ideal reporter for luminescence imaging. Rather, Gluc is a small monomer enzyme that is rapidly secreted from cells (Verhaegent and Christopoulos, 2002) and much more sensitive than the other established secretory marker, alkaline phosphatase (Wurdinger *et al.*, 2008). We demonstrate that Gaussia luciferase expressed in melanoma cells is secreted into culture medium *in vitro*; *in vivo* Gluc diffuses into blood circulation and subsequently is excreted into the urine in melanoma-bearing mice. This has been confirmed by reconstitution of Gluc diffusion through the vascular system (Fig. 17.4). Thus, Gluc may serve as an ideal reporter of intratumoral NF- κ B activity by following the release of Gluc into the blood and/or urine of tumor-bearing animals. During the course of exposure to drug treatment, the NF- κ B activity in cultured tumor cells or in mouse tumor xenografts can be continuously monitored by measuring the Gaussia luciferase activity in cultured medium or in blood and/or urine samples. Consequently, Gluc can play an important role toward predicting cancer drug efficacy *in vitro* and *in vivo*. This system has been validated by well-established IKK inhibitors and IKK stimuli, suggesting that this approach is feasible for screening the effect of cancer drugs on the NF- κ B pathway.

To prepare blood samples, 5 μ l blood were withdrawn using pipette tip from a small incision at the tail tip of conscious mice. The 5 μ l of blood was mixed into 10 μ l buffer (500 mmol/l NaCl, 2 mmol/l KCl, 10 mmol/l MgCl₂, 10 mmol/l Na₂HPO₄, 2 mmol/l KH₂PO₄, 1 mmol/l EDTA, pH 7.8) and stored at 4 °C for measuring Gluc activity within 3 weeks. Sampling of urine is easily done since the mouse readily offers approximately 100 μ l of urine triggered by the fear response to being “caught” by the gentle hand of the investigator. Of notice, the daily average mouse voiding frequency is 16 times and the urine volume per void is approximately 160 μ l; thus, urine samples were also collected from a clean glass jar where mouse was temporally kept. The urine Gluc activity was determined between 0 and 24 h without loss of Gluc activity. Gluc activity was measured using the Monolight TM 3010 luminometer (BD Biosciences Pharmingen, San Diego, CA). Ten microliters of sample were quickly mixed well with 20 μ l of 100 μ M coelenterazine in the buffer above and 10-s photon counts were acquired.

In conclusion, bioluminescence imaging of intratumoral NF- κ B signaling in living animals can provide useful information about the molecular basis of tumorigenesis and molecular response to therapeutic agents. Firefly luciferase is superior to Gaussia luciferase as a reporter for quantitative molecular imaging. Advantages of Gaussia luciferase are that it is naturally secreted and extremely sensitive, suggesting it is extremely useful as a peripheral marker for real-time monitoring of intratumoral molecular signaling in conscious mice.

Acknowledgments

We thank colleagues at the Richmond laboratory for insightful discussions and Richard Baeza of Vanderbilt University Institute of Imaging Science for excellent technical assistance. Work from the authors' program was supported through funding by the Department of Veterans Affairs through a VA Merit Award (AR) and a VA Senior Research Career Scientist Award (AR), National Institutes of Health grants (CA 098807) (A.R.), the Vanderbilt Ingram Cancer Center support grant (CA 68485), and the Skin Disease Research Center grant (SP30 AR 41943).

REFERENCES

- Phair RD, Misteli T. Kinetic modelling approaches to *in vivo* imaging. *Nat. Rev. Mol. Cell. Biol.* 2001; 2:898–907. [PubMed: 11733769]
- Tannous BA, Kim DE, Fernandez JL, Weissleder R, Breakefield XO. Codon-optimized Gaussia luciferase cDNA for mammalian gene expression in culture and *in vivo*. *Mol. Ther.* 2005; 11:435–443. [PubMed: 15727940]
- Verhaegent M, Christopoulos TK. Recombinant Gaussia luciferase. Overexpression, purification, and analytical application of a bioluminescent reporter for DNA hybridization. *Anal. Chem.* 2002; 74:4378–4385. [PubMed: 12236345]
- Wurdinger T, Badr C, Pike L, de Kleine R, Weissleder R, Breakefield XO, Tannous BA. A secreted luciferase for *ex vivo* monitoring of *in vivo* processes. *Nat. Methods.* 2008; 5:171–173. [PubMed: 18204457]
- Yang J, Pan WH, Clawson GA, Richmond A. Systemic targeting inhibitor of kappaB kinase inhibits melanoma tumor growth. *Cancer Res.* 2007; 67:3127–3134. [PubMed: 17409419]




Figure 17.1. Intratumoral NF- κ B reporter models. (A) NF- κ B-Fluc reporter model. Mouse melanocytes null for INK4A/ARF were genetically engineered with the NF- κ B-Fluc reporter (lanes 1 to 5), CMV promoter driven oncogenic H-RAS^{V12} expression (lanes 2 to 5), and/or a Tet-On inducible I κ B α (S32A/S36A) superrepressor expression vector (lanes 3 and 5). These cells were inoculated into nude mice to develop xenografts. Tumors were photographed (upper panels) and intratumoral NF- κ B activities were determined by quantitative luminescent imaging (lower panels), illustrating that H-RAS^{V12} induced NF- κ B activation *in vivo* was inhibited with the I κ B α superrepressor. (B) NF- κ B-Gluc reporter model. NF- κ B-Gluc was stably expressed in human melanoma Hs294T cells (1×10^6) and these cells were

subcutaneously inoculated into nude mice. The individual bioluminescent image was taken immediately following intravenously injection of 100 μg native coelenterazine per mouse.




Figure 17.2.

Melanoma metastasis to the liver. (A) NF- κ B-Fluc reporter model of metastasis. 2×10^5 of mouse melanocytes (NF- κ B-Fluc⁺, H-RAS^{V12+}, INK4A/ARF^{-/-}) were intravenously injected into each nude mouse. Sixty days after injection, the bioluminescent images were obtained (upper panel), mice were euthanized and (B) histological analysis (H&E staining) of organs were performed, indicating melanoma metastasis to the liver. (C) 2×10^5 of GFP-tagged melanocytes (H-RAS^{V12+}, INK4A/ARF^{-/-}) were intravenously injected into nude mice. Thirty days after cell injection, frozen sections were examined under the fluorescent microscope (10 \times magnification). H, hepatocytes; M, melanoma cells; RBC, red blood cells.




Figure 17.3.

Luminescent imaging of cultured dish. Hs294T $\text{I}\kappa\text{B-Gluc}$ reporter cells or Hs294T $\text{NF-}\kappa\text{B-Gluc}$ reporter cells were treated with the IKK inhibitor, BMS-345541 (0 to 10 mM) for 36 h. The plate containing the cultured cells and medium was subjected to bioluminescent imaging immediately following addition of 100 mM of coelenterazine. Two independent experiments were performed and similar results were achieved. Data show an increase in the level of $\text{I}\kappa\text{B}\alpha$ (upper panels) and reduced $\text{NF-}\kappa\text{B}$ transcriptional activity (lower panels) upon treatment with BMS-345541.




Figure 17.4. Reconstitution of Gluc diffusion via vascular system. Conditioned medium containing the secreted Gluc (2×10^8 RLU) was intraperitoneally injected into adult mice. The Gluc activity in the blood or the urine was detected over the course of 4 h. The average value of each time point shown was from six animals. Experiments were repeated twice with similar results. The peak time for Gluc diffusion from the intraperitoneal compartment into blood is 1h and into urine is 1.5 h.