

DNA amplification of a further exon of Duchenne muscular dystrophy locus increase possibilities for deletion screening

A.Speer, A.Rosenthal, H.Billwitz, R.Hanke, S.M.Forrest, D.Love, K.E.Davies and Ch.Coutelle

Molecular Genetics Group, Institute of Molecular Medicine, John Radcliffe Hospital, Headington, Oxford OX3 9DU, UK

Submitted May 16, 1989

The data of Chamberlain et al (1) allow the detection of 76% of deletions in the region Cf56a/Cf23a identified by hybridization in our patients. We have generated two oligonucleotides allowing the amplification of a further exon which is included in the 3.4kb hybridization fragment of Cf56a. This exon is deleted in about 10% of our patients.

5' cccaaaatatttttag CTCCTACTCAGACTGTTACTCTGGTGACACAACCTGTGGTTACT
AAGGAACTGCCATCTCAAACTAGAAATGCCATCTTCCTTGATGTTGGAGGTACCTGCTCTGGC
AGATTTCACCGGGCTTGACAGAACTTACCGACTGGCTTTCTCTGCCTTGATCAAGTTATAAAAT
CACAGAGGGTGATGGTGGGTGACCTTGAGGATATCAACGAGATGATCATCAAGCAGAAGgtatga
gaaaaaatgataaaagttggcagaagtttttcttttaaataagattttccaccaatcactttac
tctcctagaccattgccaccagttcttaggcaactgtttctctctcagcaaacacattactctc
actattcagcctaagtataatcaag 3'

Fig. 1: Nucleotide sequence of the exon included in the 3.4kb hybridization fragment of Cf56a and a part of the flanking introns. Exon sequences are in upper case, and intron sequences are in lower case. Primer sequences are underlined.




Fig. 2: Exon-containing PstI fragments detected by Cf56a and their deletion in patients D22, D101 and D132.




Fig. 3: Detection of DNA amplification products of 409 bp (PstI-3.4 kb/Cf56a) and 506 bp (HindIII-1.2+3.8 kb/probe 8, equivalent of Cf56a) (1) in the deletion cases and one normal individual.

REFERENCES

1. Chamberlain, J.S. et al. (1989) Nucl Acids Research, 16 (23): 11141-11156.