## Shotgun Isotope Array for Rapid, Substrate-Specific Detection of Microorganisms in a Microbial Community<sup>⊽</sup>†

Tomohiro Tobino,<sup>1</sup>\* Futoshi Kurisu,<sup>2</sup> Ikuro Kasuga,<sup>1</sup> and Hiroaki Furumai<sup>2</sup>

Department of Urban Engineering, Graduate School of Engineering, The University of Tokyo, 7-3-1 Hongo, Bunkyo-ku, Tokyo 113-8656, Japan,<sup>1</sup> and Research Center for Water Environment Technology, Graduate School of Engineering, The University of Tokyo, 7-3-1 Hongo, Bunkyo-ku, Tokyo 113-8656, Japan<sup>2</sup>

Received 19 January 2011/Accepted 15 August 2011

The shotgun isotope array method has been proposed to be an effective new tool for use in substrate-specific microbe exploration without any prior knowledge of the community composition. Proof of concept was demonstrated by detection of acetate-degrading microorganisms in activated sludge and further verified by independent stable isotope probing (SIP).

Isotope tracer techniques can be combined with molecular biological analyses to provide a means for elucidating the functions of specific microbes identified within a complex microbial community (2, 12). Nucleic acid-based stable isotope probing (SIP) (7, 10) and fluorescence *in situ* hybridization followed by microautoradiography (5, 9) have been widely used in functional studies of naturally occurring microbial communities (3, 12). The isotope array approach has also been reported as an alternative technique in which <sup>14</sup>C-labeled 16S rRNAs are detected by direct hybridization to oligonucleotide microarrays (1, 4). By the use of nanoscale secondary-ion mass spectrometry, incorporation and quantification of isotopes in microbial cells, together with their phylogenetical identities, can be visualized at the single-cell or subcellular level (6).

Here, we propose a novel method, namely, a shotgun isotope array approach (Fig. 1), which has potential advantages compared to the isotope array technique. In this approach, a <sup>14</sup>C-labeled compound is used as a tracer substrate, and extracted DNA is hybridized to a shotgun array (also known as metagenomic array [11]) that consists of genomic DNA fragment probes obtained by shotgun cloning of the sample to be analyzed (14). Sequences of probes with positive radio signals are then read to obtain information on the microorganisms involved in the assimilation of the tracer substrate. This shotgun array offers several advantages over oligonucleotide arrays, such as (i) freedom from the requirement for probe design and selection, (ii) applicability to any given sample, and (iii) the ability of the probe set to reflect the community composition of the sample, enabling unknown microorganisms to be detected. Proof of concept was demonstrated by hybridization of genomic DNA extracted from activated sludge grown in the presence of [<sup>14</sup>C]acetate with a membrane array prepared from the sludge DNA. The hybridization results were further verified by independent SIP.

An activated sludge sample (2,200 mg of suspended solids per liter) was collected from a bench-top conventional activated sludge process reactor that treats municipal wastewater in Japan. In a glass vial, 27 ml of the sludge sample was incubated under anoxic conditions (100 mg of N liter<sup>-1</sup> nitrate) with 660 mg liter<sup>-1</sup> sodium acetate containing 1.7 mCi [1-<sup>14</sup>C]sodium acetate (Moravek Biochemicals, Brea, CA) at room temperature on a shaker. During incubation, small subsamples were taken to monitor <sup>14</sup>C-labeled substrate uptake by the use of a liquid scintillation counter. Target <sup>14</sup>C-labeled DNA was extracted after 18 h and sonicated to obtain fragments averaging 400 bp.

Random genomic DNA fragment probes were prepared by shotgun cloning of the sludge DNA followed by PCR amplification and were manually spotted onto a nylon membrane. The membrane array consisted of 96 fragment probes (~2,000 bp in length) and both positive- and negative-control probes. Target <sup>14</sup>C-labeled DNA was hybridized to the membrane array in a plastic bag with 1.5 ml of hybridization buffer (digoxigenin [DIG] Easyhyb; Roche) and mixed gently at 55°C for 16 h. After washing was performed, radio signals on the membrane were detected using an imaging plate (MS-2010; Fujifilm, Tokyo, Japan) and an image reader (FLA-9000; Fujifilm). Spots that showed a signal-to-noise ratio (SNR) of >3 were considered to represent positive signals. Partial sequences (approximately 700 bp from one end) were determined for all the positive probes and four negative probes and searched in the DDBJ/EMBL/GenBank database.

SIP of the sludge sample was conducted using [1,2-<sup>13</sup>C]sodium acetate (99 atom%; Icon Isotopes, Summit, NJ) and unlabeled sodium acetate under the same conditions as and in parallel with the [<sup>14</sup>C]acetate incubation described above. Subsamples (2 ml) were taken every 6 h, and DNA was extracted from each subsample. Gradient density centrifugation was carried out essentially as described previously (8), and 16 to 18 density fractions were collected per tube. The copy numbers of the five positive and four negative probe sequences in each density fraction at different sample times were quantified by real-time PCR, using a primer pair that was designed to specifically amplify the probe sequence in question (see Table S2 and Fig. S2 in the supplemental material for further details of the materials and methods used).

<sup>\*</sup> Corresponding author. Mailing address: Department of Urban Engineering, Graduate School of Engineering, The University of Tokyo, 7-3-1 Hongo, Bunkyo-ku, Tokyo 113-8656, Japan. Phone: 81-3-5841-0669. Fax: 81-3-5841-6244. E-mail: tobino@env.t.u-tokyo.ac.jp.

<sup>†</sup> Supplemental material for this article may be found at http://aem .asm.org/.

<sup>&</sup>lt;sup>v</sup> Published ahead of print on 26 August 2011.



FIG. 1. Schematic diagram of shotgun isotope array approach.

**Detection of acetate-degrading microorganisms by shotgun isotope array.** [<sup>14</sup>C]acetate in the supernatant was consumed constantly, and 19% of the initial amount was incorporated into biomass after 18 h (see Fig. S1 in the supplemental material). The extracted DNA was hybridized to a shotgun array, and 20 of 96 probes gave positive signals (Fig. 2A and B). Some of the probe sequences that gave positive signals had a high level of similarity with sequences of nitrate-reducing bacteria, such as *Thauera*, *Azoarcus*, *Acidovorax*, and *Poralomonas* (Table 1; for a complete result set, see Table S1 in the supplemental material). However, the similarity values (66.9 to 89.8% for those with e-values of <1E-5) were not high enough for reliable phylogenetic identification. Positive identification should soon become more accurate as a result of the rapidly occurring growth of genome databases. Moreover, even without positive identification, the probe sequences can be used to monitor microorganisms of interest.

Table 1 summarizes the SIP results for each probe, together with the closest matching microorganisms found in the public databases (see Fig. S3 in the supplemental material for the SIP profiles). Four of the five positive probes gave consistent results with the shotgun isotope array approach, such that assimilation of the labeled compound was seen as the rise in the levels of the specific probe sequences in higher-density fractions over time. However, probe 12 gave a positive result in the hybridization assay, but substrate assimilation was not clearly observed in the SIP. That signal, as seen in the shotgun isotope array, might have arisen from <sup>14</sup>C-labeled DNA cross-hybridizing with the probe. Wu et al. reported that the threshold for discrimination by the use of PCR-amplified gene fragments was 80 to 85% sequence similarity for ~500-bp DNA-DNA hybridization on glass slides (13). In our preliminary study, higher specificity for discriminating instances of 85 to 90% sequence similarity could be achieved by the use of 2,000-bp DNA fragment probes under highly stringent conditions (hybridization performed at 75°C; data not shown). Further study is necessary to elucidate the specificity of shotgun fragment probes to ensure reliable microbial detection. As for the four negative probe sequences tested, three of them showed no clear substrate assimilation in the SIP. One probe, probe-59, was judged to be negative in the hybridization assay but showed a weak peak shift in the SIP (see Fig. S3 in the supplemental material). Insufficient labeling might result in negative signals in shotgun isotope array analysis, although these signals could remain detectable in the SIP due to the high sensitivity of PCR as a means of DNA sequence detection.

Further experiments are required to assess the application of the shotgun isotope array method. One promising option is to use tools of greater sensitivity, such as a  $\beta$ -Imager, in autoradiographic detection (4, 10). Higher sensitivity would enable not only the detection of less-abundant species but also the use of hybridization conditions of greater stringency to achieve



FIG. 2. Hybridization results. (A) Autoradiograph of shotgun array after hybridization with <sup>14</sup>C-labeled DNA in activated sludge. Arrows indicate positive controls. (B) Signal-to-noise ratios (SNR) of the probes. Asterisks above the bars indicate SNRs of >3. The probes used for verification in SIP are shown with closed bars, and arrows indicate each probe number.

Probe name	Hybridization result <sup>a</sup>	Closest match in databases	Gene feature(s) or encoded protein	E-value	% similarity (no. of identical nucleotides/total no. of nucleotides)	Assimilation of tracer compound in SIP <sup>c</sup>
probe-04	р	<i>Thauera</i> sp. MZ1T, complete genome	2-Isopropylmalate synthase	0	82.5 (611/741)	+
probe-07	р	Acidovorax sp. JS42, complete genome	Transposase, IS4 family protein	1.0E-137	78.0 (499/640)	+
probe-12	р	Methylococcus capsulatus strain Bath, complete genome	Glutathione S-transferase domain protein	1.0E-09	66.9 (192/287)	_
probe-16	р	Polaromonas naphthalenivorans CJ2, complete genome	RNA polymerase, sigma-24 subunit, ECF subfamily	2.0E-109	80.1 (370/462)	+
probe-18	р	Polaromonas naphthalenivorans CJ2 plasmid pPNAP05, complete sequence	NA <sup>b</sup>	0	81.5 (567/696)	+
probe-36	n	Plasmodium knowlesi strain H	DNA repair protein RAD50	6.2E-01	96.3 (26/27)	_
probe-49	n	Azoarcus sp. BH72, complete genome	Penicillin-binding protein	2.2E+00	96.2 (25/26)	-
probe-59	n	Paracoccus denitrificans PD1222 chromosome 2, complete sequence	hydrogenase-1 expression HyaE	1.0E-59	69.0 (460/667)	+
probe-86	n	Pedobacter heparinus DSM 2366, complete genome	Na <sup>+</sup> /H <sup>+</sup> antiporter	2.2E+00	90.3 (28/31)	-

TABLE 1. Verification of shotgun isotope array detection by SIP

<sup>a</sup> In [<sup>14</sup>C]acetate hybridization. p, positive; n, negative.

<sup>b</sup> NA, not available.

<sup>c</sup> Data indicate whether an increase in copy numbers in high-density fractions was observed (+, observed; -, not observed).

higher specificity. In addition, miniaturization with a glass slide array format using a larger number of probes would enable broader coverage for detection of microbe species.

**Nucleotide sequence accession numbers.** The nucleotide sequences used in this study can be found in the DDBJ/EMBL/ GenBank under accession numbers AB609074 to AB609122.

We express our gratitude to the staff of the Radioisotope Center for facilitating radioisotope analyses.

The study was supported, in part, by the CREST (Core Research for Evolutional Science and Technology) of the Japan Science and Technology Agency (JST).

## REFERENCES

- Adamczyk, J., et al. 2003. The isotope array, a new tool that employs substrate-mediated labeling of rRNA for determination of microbial community structure and function. Appl. Environ. Microbiol. 69:6875–6887.
- Dumont, M. G., and J. C. Murrell. 2005. Stable isotope probing—linking microbial identity to function. Nat. Rev. Microbiol. 3:499–504.
- Friedrich, M. W. 2006. Stable-isotope probing of DNA: insights into the function of uncultivated microorganisms from isotopically labeled metagenomes. Curr. Opin. Biotechnol. 17:59–66.
- Hesselsoe, M., et al. 2009. Isotope array analysis of Rhodocyclales uncovers functional redundancy and versatility in an activated sludge. ISME J. 3:1349– 1364.

- Lee, N., et al. 1999. Combination of fluorescent in situ hybridization and microautoradiography—a new tool for structure-function analyses in microbial ecology. Appl. Environ. Microbiol. 65:1289–1297.
- Li, T., et al. 2008. Simultaneous analysis of microbial identity and function using NanoSIMS. Environ. Microbiol. 10:580–588.
- Manefield, M., A. S. Whiteley, R. I. Griffiths, and M. J. Bailey. 2002. RNA stable isotope probing, a novel means of linking microbial community function to phylogeny. Appl. Environ. Microbiol. 68:5367–5373.
- Neufeld, J. D., et al. 2007. DNA stable-isotope probing. Nat. Protoc. 2:860– 866.
- Ouverney, C. C., and J. A. Fuhrman. 1999. Combined microautoradiography-16S rRNA probe technique for determination of radioisotope uptake by specific microbial cell types in situ. Appl. Environ. Microbiol. 65:1746–1752.
- Radajewski, S., P. Ineson, N. R. Parekh, and J. C. Murrell. 2000. Stableisotope probing as a tool in microbial ecology. Nature 403:646–649.
  Sebat, J. L., F. S. Colwell, and R. L. Crawford. 2003. Metagenomic profiling:
- Sebat, J. L., F. S. Colwell, and R. L. Crawford. 2003. Metagenomic profiling: microarray analysis of an environmental genomic library. Appl. Environ. Microbiol. 69:4927–4934.
- Wagner, M., P. H. Nielsen, A. Loy, J. L. Nielsen, and H. Daims. 2006. Linking microbial community structure with function: fluorescence in situ hybridization-microautoradiography and isotope arrays. Curr. Opin. Biotechnol. 17:83–91.
- Wu, L. Y., et al. 2001. Development and evaluation of functional gene arrays for detection of selected genes in the environment. Appl. Environ. Microbiol. 67:5780–5790.
- Yokoi, T., et al. 2007. 'FloraArray' for screening of specific DNA probes representing the characteristics of a certain microbial community. FEMS Microbiol. Lett. 273:166–171.