

Published in final edited form as:

*Bioanalysis*. 2012 January ; 4(2): 157–166. doi:10.4155/bio.11.303.

## Determination of lumefantrine in small-volume human plasma by LC–MS/MS: using a deuterated lumefantrine to overcome matrix effect and ionization saturation

Liusheng Huang<sup>1,\*</sup>, Xiaohua Li<sup>2</sup>, Florence Marzan<sup>1</sup>, Patricia S Lizak<sup>1</sup>, and Francesca T Aweeka<sup>1</sup>

<sup>1</sup>Drug Research Unit, Department of Clinical Pharmacy, School of Pharmacy, University of California, San Francisco, CA 94110, USA

<sup>2</sup>Drug Studies Unit, Department of Bioengineering & Therapeutic Sciences, School of Pharmacy, University of California, San Francisco, CA, USA

### Abstract

**Background**—To support pediatric study, a method to determine lumefantrine (LF) with small sample volume is needed. Matrix effect (ME) is a daunting issue in LF quantification in human plasma with LC MS/MS.

**Results**—Here we report an LC–MS/MS method with a deuterated LF as the internal standard (IS). Plasma sample (25 µl) was acidified with 5% formic acid prior to extraction with ethyl acetate. The recovery was over 80%. The absolute ME was within the range of 100 ± 8% for both LF and the IS, but cumulative ME was observed via large variation of IS as the IS. The linear signal. The cumulative ME and ionization saturation were overcome with the co-eluting LF-D<sub>9</sub> range of calibration curve was 50–20,000 ng/ml.

**Conclusion**—ME and ionization saturation was overcome with a deuterated IS. The method utilized a small sample volume, suitable for pediatric study with capillary tube blood collection method.

Lumefantrine (LF), previously named benflumetol, is a highly hydrophobic aromatic compound (Figure 1). It has a logP of 8.67 (25°C) and logD of 5.6 (pH 1–4, 25°C) and 7.0 (pH 7, 25°C), calculated using Advanced Chemistry Development Software V11.02. LF is used in combination with artemether under the brand name of Co-Artem<sup>®</sup> for antimalaria therapy [1]. The peak concentration of LF in patients is usually over 10,000 ng/ml reached

© 2012 Future Science Ltd.

\*Author for correspondence: Tel.: +1 415 502 2594, Fax: +1 415 476 0307, huangl@sfgsom.ucsf.edu.

For reprint orders, please contact reprints@future-science.com

Supplementary data

Supplementary data accompanies this paper and can be found at [www.future-science.com/doi/full/10.4155/BIO.11.303](http://www.future-science.com/doi/full/10.4155/BIO.11.303)

### Ethical conduct of research

The authors state that they have obtained appropriate institutional review board approval or have followed the principles outlined in the Declaration of Helsinki for all human or animal experimental investigations. In addition, for investigations involving human subjects, informed consent has been obtained from the participants involved.

### Financial & competing interests disclosure

This work was supported by NIH/NIAID International Maternal Pediatric Adolescent AIDS Clinical Trials (IMPAACT) Group (grant number 1U01AI068632), ACTG (1U01AI068636–02), UCSF-GIVI Center for AIDS Research (CFAR) and Novartis Pharm. Co. The authors have no other relevant affiliations or financial involvement with any organization or entity with a financial interest in or financial conflict with the subject matter or materials discussed in the manuscript apart from those disclosed.

No writing assistance was utilized in the production of this manuscript.

6–8 h after dosing [2]. Numerous methods have been reported for determination of LF in human plasma, most of these methods are based on HPLC-UV where a large plasma volume ( $\geq 200$   $\mu\text{l}$ ) was used [3–8]. Methods for quantifying LF on dry blood spots were also reported with relative small sample volume (50–100  $\mu\text{l}$ ) but the LLOQ was high (130–150 ng/ml) [9,10]. For pediatric studies, a sensitive method with small sample volume and acceptable LLOQ is required. Early attempts on application of LC–MS method for determination of LF in human plasma failed due to matrix effect (ME) [11]. Recently, several LC–MS/MS methods were reported, but the issue of ME has not been addressed sufficiently [12–15]. Using protein precipitation for sample preparation, Hodel *et al.* reported the first LC–MS/MS method for LF, but the ME was significant: 22.7–29.1% (100% = no ME) [12], while César *et al.* claimed the ME in their method was within the 85 to 115% range, but the calibration curve required a weighted quadratic regression [15]. None of these methods utilized a co-eluting internal standard (IS), which is the most effective for overcoming ME. A stable isotope-labeled analyte is generally considered to be the best IS. Because the analyte and IS are generally co-eluted and ionized identically, ME will affect both analyte and IS to the same extent and ionization saturation occurring at high analyte concentrations will also affect the IS, which in turn, compensates for those effects on the analyte. Here we report an LC–MS/MS method using a deuterated LF (LF-D<sub>9</sub>) as the IS, with a focus on ME and ionization saturation. The method was validated according to AIDS Clinical Trials Group (ACTG) guidelines and crossvalidated with an HPLC-UV method.

## Experimental

### Chemicals

LF (Figure 1) was purchased from AK Scientific Inc. (CA, USA). LF-D<sub>9</sub> was obtained as a gift from Novartis Pharma. Co. (NJ, USA). Acetonitrile (MeCN), water (HPLC grade) and ammonium formate were obtained from Fisher Scientific (NJ, USA). Ethyl acetate and formic acid (FA) was purchased from Sigma-Aldrich (MO, USA). All solvents were HPLC grade and chemicals were ACS reagents. Human plasma (EDTA as the anticoagulant) was purchased from Biological Specialty Co. (PA, USA).

### LC–MS/MS conditions

The LC–MS/MS system consisted of twin PE 200 micro-LC pumps, PE 200 autosampler (Perkin-Elmer; CT, USA), and the API 2000 triple quadrupole MS system (AB Sciex, ON, Canada). The instrument was placed in a room controlled with an air conditioner. Chromatographic separation was achieved on a Zorbax C<sub>8</sub> column (50  $\times$  2.1 mm, 5  $\mu\text{m}$ ; Agilent Technologies Inc., CA, USA). The LC setting was as follows: solvent A was aqueous ammonium formate 10 mM at pH 4.0. Solvent B was MeCN with FA 0.1%. The gradient program consisted of linear segments with 50% B (0–1 min), 50–100% B (1–4 min), 100% B (4–6 min), 100–50% B (6–6.1 min) and 50% B (6.1–8 min). The flow rate was 0.4 ml/min. Injection volume was 10  $\mu\text{l}$ . The retention times for LF and IS were both 3.5 min. ESI<sup>+</sup> was used for ion source and multiple reaction monitoring mode was chosen for quantification. The precursor–product ion pairs were  $m/z$  528 $\rightarrow$ 510 for LF and  $m/z$  537 $\rightarrow$ 519 for the IS. The optimized acquisition parameters were as follows: Turbo (Heater) set at 400°C; curtain gas: 25 psi; nebulizer gas (gas 1): 40 psi; auxiliary (turbo) gas (gas 2): 70 psi; collision-activated dissociation gas: 4; IonSpray Voltage: 4000 V. All gas lines were supplied from a liquid nitrogen tank (industrial grade). The optimized parameters for LF and IS are as follows: declustering potential: 56; entrance potential: 9.5; focusing potential: 370; collision cell entrance potential: 22; collision energy: 29; collision cell exit potential: 24. The scan time was set at 250 ms for each transition and pause between mass ranges was set at 5 ms. Data were processed with Analyst 1.4.2. (Danaher Co., Washington, DC, USA).

### Preparation of standard & validation (quality control) samples

Primary stock solutions (1 mg/ml) and working solutions of LF and IS were both prepared in MeCN:water (1:1) containing 0.5% FA. Calibration standards were prepared by spiking the LF working solutions into plasma to obtain concentrations of 50, 100, 500, 1000, 5000, 10,000, 15,000 and 20,000 ng/ml. Validation samples (also called quality control [QC]) were prepared at four different concentrations: 50 ng/ml (LLOQ), 120 ng/ml (low), 1500 ng/ml (medium) and 17,000 ng/ml (high). Calibration standards and validation samples were prepared from two different stock solutions made with separately weighted LF. The stock solutions, standards, QC samples and the IS solution were stored at  $-70^{\circ}\text{C}$  before uses.

### Sample preparation

A 25  $\mu\text{l}$  aliquot of plasma sample was mixed with 25  $\mu\text{l}$  100 ng/ml LF-D<sub>9</sub> and 100  $\mu\text{l}$  5% aqueous FA. The mixture was extracted with 900  $\mu\text{l}$  ethyl acetate by vortex mixing for 20 s followed by rotating on a tube rotor for 30 min. After centrifugation at  $15,000 \times g$  for 2 min, the organic phase was transferred into a glass tube (12  $\times$  75 mm) with a fine-tip disposable pipet, dried in  $40^{\circ}\text{C}$  water bath with purging of  $\text{N}_2$  and re-dissolved with 200  $\mu\text{l}$  MeCN:water (1:1, 0.5% FA) by vortex mixing for 10 s. The reconstituted sample was transferred to an autosampler vial. Injection volume was 10  $\mu\text{l}$ .

### Method validation procedure

The method validation was conducted according to the ACTG guidelines [16], which were developed based on US FDA guidelines. The calibration curve should contain at least six nonzero points with back-calculated values within the range of 85 to 115% of nominal values, except for LLOQ in which 80–120% is allowed. Intra- and inter-day accuracy and precision for low-, medium- and high-validation samples should be within  $\pm 15\%$ , but for LLOQ validation samples this can be  $\pm 20\%$ . Intra-day accuracy and precision were determined with at least five replicates of validation sample at each concentration (LLOQ, low, medium and high) along with a set of calibrators. The same experiment was performed on at least five different days to determine inter-day accuracy and precision. Stability tests were performed at ambient temperature,  $-70^{\circ}\text{C}$ , four freeze–thaw cycles, and in autosampler vials with low and high QC samples. The determined concentrations were compared with nominal values. Plasma samples at 60,000 ng/ml were diluted by four-, eight- and 12-fold to test sample integrity upon dilution. Stock solution in MeCN 50% containing FA 0.5% was diluted to 250 ng/ml and injected onto LC–MS/MS. The peak area was compared with that from freshly made stock solution. The following 12 potential concomitant drugs were tested for interference: nevirapine, lopinavir, ritonavir, zidovudine, lamivudine, efavirenz, chloroquine, sulfamethoxazole, trimethoprim, artemether, dihydroartemisinin and tenofovir, each drug spiked at 5000 ng/ml in medium validation sample (LF 1500 ng/ml). The ME was evaluated with three different approaches.

**Infusion experiment**—A 10  $\mu\text{g}/\text{ml}$  solution of LF was infused constantly at a flow rate of 10  $\mu\text{l}/\text{min}$  into the LC elute via a ‘T’ connector and the LC elute was directed into the MS source. A 10  $\mu\text{l}$  aliquot of blank plasma extract was injected onto LC column and signal for LF was monitored for 60 min when LC ran in isocratic mode at 0.4 ml/min (65% solvent B). A 10  $\mu\text{l}$  aliquot of 10  $\mu\text{g}/\text{ml}$  LF in mobile phase solvent was injected as a reference. The LC was also operated manually according to the final gradient method and LF signal was monitored for 29 min. A 10  $\mu\text{l}$  aliquot of mobile-phase solvent was injected as reference.

**The approach established in ACTG guidelines**—Six different lots of blank plasma were processed and injected into the LC–MS/MS. Medium QCs were prepared in triplicate

with each of six plasma lots, processed and injected into LC–MS/MS, their concentrations were compared with nominal concentration (1500 ng/ml).

**The approach recommended in ACTG guidelines**—This approach is based on the work of Matuszewski and co-workers [17]. For recovery and ME, four sets of samples were prepared. Set 1a, aliquots (25 µl each) of LF solutions at low (120 ng/ml), medium (1500 ng/ml) and high (17,000 ng/ml) levels were diluted with 175 µl MeCN:water (1:1) containing 0.5% FA and analyzed in six replicates (n = 6). A 25 µl aliquot of IS (100 ng/ml) was also diluted with 175 µl mobile-phase solvent and analyzed in six replicates.

Set 1b, aliquots (25 µl each) of LF solutions at low, medium and high levels were mixed with 25 µl IS (100 ng/ml), diluted with 150 µl MeCN:water (1:1) containing 0.5% FA, and analyzed in six replicates (n = 6).

Set 2, 25 µl aliquots of blank plasma from six different lots were processed with liquid–liquid extraction (LLE) and reconstituted with 200 µl mobile phase solvent containing LF and IS to make the final concentration the same level as those of the set 1b samples.

Set 3, a 25 µl aliquot of each plasma sample at the low (120 ng/ml), medium (1500 ng/ml) and high (17000 ng/ml) levels was processed following the normal procedure. These samples were prepared in the six different lots of plasma used in set 2.

The recovery (RE), ME and process efficiency (PE) of LF and IS were assessed by comparing the peak areas of LF and IS with the following formulae:

$$\begin{aligned} \text{Recovery (RE)} &= \frac{100 \times \text{Peak area of pre-extraction-spiked sample (set 3)}}{\text{Peak area of postextraction-spiked sample (set 2)}} \\ \text{Matrix effect (ME)} &= \frac{100 \times \text{Peak area of postextraction-spiked sample (set 2)}}{\text{Peak area of clean sample (set 1b)}} \\ \text{Process efficiency (PE)} &= \frac{100 \times \text{Peak area of pre-extraction-spiked sample (set 3)}}{\text{Peak area of clean sample (set 1b)}} \end{aligned}$$

### Crossvalidation

Ninety-six clinical samples from a previous study were reanalyzed with this assay and the results were compared with those from an HPLC–UV method, based on the criteria as described by Fast *et al.* [18]. At least 67% of reanalysis should be within 20% compared with the reference value, which is the mean of the initial and reanalyzed value:

$$\text{Dev}\% = \frac{100 \times (C_{\text{This assay}} - C_{\text{HPLC-UV}})}{[(C_{\text{This assay}} + C_{\text{HPLC-UV}})/2]}$$

All clinical samples were from studies approved by the institutional review board at University of California, San Francisco, USA.

### Clinical application to pediatric patients

The developed method was used to determine LF concentration from pediatric patients. Blood samples were collected with capillary tubes (200 µl capacity) and ~100 µl or less plasma was obtained. A total of 260 samples were analyzed with this method.

## Results & discussion

### Method development

Co-eluting stable isotope-labeled analyte as the IS is considered as the best IS to compensate for signal variation due to ME. Here we used a LF-D<sub>9</sub> as the IS. When a Waters symmetry<sup>®</sup> C<sub>18</sub> column (150 × 4.6 mm, 5 μm) was used, LF-D<sub>9</sub> was only partially co-eluted with LF. The ME on two compounds can be different if they are not co-eluted [19]. After switching to a Zorbax Eclipse C<sub>8</sub> column (50 × 2.1 mm, 5 μm), the IS was co-eluted with LF (Supplementary Figure 1).

In reverse-phase LC–MS/MS methods, ME was commonly seen in the early elution phase due to polar substances such as salts, and late phase due to hydrophobic substances such as phospholipids. In some cases, the highly hydrophobic substances do not affect the immediate sample but affect the subsequent samples analysis, especially in isocratic elution method. Here LC conditions were optimized to elute LF with gradient elution and clean the column with 100% organic mobile phase so as to achieve optimal separation and reduce hydrophobic matrix carried over to subsequent analysis.

### Sample preparation

Sample preparation is an important step in method development to overcome ME. A recently published LC–MS method utilized simple protein precipitation, ME reached 22.7% suggesting ion suppression >77% [12]. Combination of protein precipitation and SPE was also used for preparation of plasma sample of LF [6–7], but it is expensive and labor intensive. LLE is a superior method for sample preparation. It gives cleaner samples compared with protein precipitation and, in some cases, better samples than SPE.

To optimize sample preparation methods, SPE and LLE with different conditions were tested (Supplementary Table 1). SPE with a C<sub>8</sub> column yielded higher RE than a HLB column, probably due to hydrophilic interaction with the charged amine group on LF during elution phase. LLE with ethyl acetate was better than methyl t-butyl ether. Acidification of the sample prior to extraction was important in order to achieve high RE. This is in consistence with the previous finding [3]: extraction with 1% acetic acid in ethyl acetate afforded much higher RE than ethyl acetate alone. It is also justified via the observation that LF is poorly soluble in water and most organic solvents, but soluble in acidified organic solvents. On the other hand, LF is highly bound to plasma protein (>99%) [20], disruption of protein binding during extraction is critical for high extraction efficiency. Previous methods required protein precipitation before SPE to release LF from plasma protein [5–7]. However, from this study, protein precipitation prior to SPE was not necessary if the plasma sample was acidified before SPE, suggesting that acidified LF was able to release easily from plasma protein. To acidify the plasma sample, FA was better than TFA and at least 0.5% FA in final solution was required to achieve high extraction RE. Munjal *et al.* used a simple protein precipitation method with 0.1% acetic acid in MeCN [14]; the RE ranged from 50–65%, probably due to insufficient acidification. The final conditions were as follows: plasma sample (25 μl) and the IS (25 μl) were acidified with 5% FA 100 μl and extracted with ethyl acetate 900 μl.

### Calibration curve

Without LF-D<sub>9</sub> as the IS, the calibration curve was not linear in the range of 50 to 20,000 ng/ml. Quadratic fitting with 1/× weighting factor was required (Figure 2a). To check whether the signal saturation was caused by multiplier or ionization saturation, a lower abundant ion pair 528/291 was selected. Signal saturation was still observed with the weaker ion pair, suggesting ionization saturation occurred in the ion source. When using LF-D<sub>9</sub> as

the IS, the signal of the IS was decreased with the increase of LF concentration (Figure 2C). Since LF and LF-D<sub>9</sub> were co-eluted and had the same affinity for protons, both signals were reduced at the saturated ion source and the ratio remained unaffected. This consequently extended the linear range for quantification of LF (Figure 2b).

### Method validation

Using a 1/ $\times$  weighted linear regression, the calibration range of this method was 50 to 20,000 ng/ml. The calibration equation over 6 days (mean  $\pm$  SD) was  $y = (0.018 \pm 0.008)\times + (0.098 \pm 0.109)$ ,  $r = 0.9984 \pm 0.0014$ , where  $y$  represents peak area ratio (LF/LF-D<sub>9</sub>),  $\times$  is LF concentration and  $r$  is the correlation coefficient. Based on the minimal (>50 ng/ml) and maximal (>10,000 ng/ml) concentrations of LF from previous pharmacokinetic studies [1,2], an LLOQ at 50 ng/ml would be sufficient. The S/N at 50 ng/ml (LLOQ) was 27.

Representative multiple reaction monitoring ion chromatograms of blank sample, LLOQ and blank sample after ULOQ were shown in Figure 3. Notably, LF signal increased in the blank sample injected after ULOQ (dash line) when compared with the same blank injected in the beginning of assay (grey line), indicating carryover. But the carryover was generally less than 20% of LLOQ, which is not significant. Intra-day precisions ( $n = 5$ ) over 6 days ranged from 1.3 to 10% and inter-day precision ( $n = 30$ ) ranged from 5.3 to 6.1%. The intra- and inter-day accuracies were also within the acceptance limit (Table 1).

LF is a very stable compound. There was no significant change of concentration for stock solution in 0.5% FA stored at  $-70^{\circ}\text{C}$  for 18 months, and no significant degradation for plasma sample stored at  $-70^{\circ}\text{C}$  for 1 year and at ambient temperature for 3 days (Table 2). In addition, sample dilution did not affect LF quantification, and the commonly concomitant drugs were not interfered with LF quantification either (Supplementary Tables 2 & 3).

### ME & RE

Evaluation of ME in LC-MS/MS methods is complicated, since in most cases the ME is not displayed through an 'observable' peak, but indirectly affects analyte peak through changing its ionization. ACTG guidelines recommend evaluating ME with a complicated approach, but the criteria have not yet been finalized. Here, ME was evaluated with four different approaches.

Firstly, an infusion experiment was performed by infusion of LF 10  $\mu\text{g}/\text{ml}$  constantly while the LC ran at 0.4 ml/min at isocratic mode (65%B) (Figure 4). LF eluted at 2.1 min (green line). When 10  $\mu\text{l}$  of blank plasma extract from LLE was injected, no significant signal change was observed at the retention time of LF (2.1 min) and, thereafter, there is only one ion enhancement peak for 1 h (blue line). However, injection of blank plasma extract from protein precipitation with MeCN resulted in ion suppression at the retention time of LF and significant signal changes observed thereafter (red line), suggesting sample preparation with simple protein precipitation presented ME in the injected sample analysis and ME may also carryover to the following sample analysis. This also justified the need for column wash after each injection. Following this, an infusion experiment was also done with the final gradient elution method implemented with a wash phase. No significant signal change occurred at the retention time of LF when the blank plasma extract from LLE was injected (Supplementary Figure 2).

Second, six different sources of blank plasma were spiked with LF at medium QC level (1500 ng/ml), and then processed and analyzed in triplicate. The deviation from nominal value (dev%) and variation from the mean (CV%) are all within 15% (Supplementary Table 4).

Third, comprehensive evaluation of ME was performed according to ACTG guidelines based on the strategies proposed by Matuszewski *et al.* [16,17]. Six different lots of human plasma were used and the raw data were provided in supplemental materials (Supplementary Table 5). The most important issue during method validation is estimation of relative ME, which can be performed through comparison of %CV in repetitive injections of neat standard (set 1b) and post-extraction spiked samples from different lots of the matrix (set 2). When comparing the %CV from sets 1b, 2 and 3, the difference between sets 1b and 2 is 0, -0.1 and 0.4% for LF at low, medium and high concentration levels, respectively. These values are 1.0, 0.1 and -0.4% for IS, respectively. All are less than 5%, which indicates no significant relative ME presents (Tables 3 & 4). Alternatively, the variability of slopes of standard line fitted through low, medium and high validation samples from at least five different lots of the matrix can be used as a good indicator of relative ME. Slopes of lines connecting low, medium and high LF samples from each lot of plasma were calculated and these are presented in Table 5. The CV% from set 3 is 4.8% (<5%), confirming that no significant ME on quantification was presented. Notably, ACTG guidelines recommend CV % of slopes should be <5%. Based on our experience with other methods, this value over 5% does not necessarily make the method invalid. Absolute ME was evaluated with mean peak area values from sets 1b and 2. A value of 100% means no ME. At low, medium and high LF concentration levels, the MEs were 97.9–100.4% for LF and 92.1–97.1% for IS, respectively (Table 6). These results indicate that ME in the method is negligible; however, effect of co-eluting substance was significant depending on concentration of LF and co-eluting matrix. High-QC concentration of LF significantly decreased the signal of LF-D<sub>9</sub> (-30.4%), indicating ionization saturation. Surprisingly, the presence of LF-D<sub>9</sub> reduced LF signal significantly at the low QC concentration while ionization saturation was not expected. Probably due to low water solubility of LF, equilibrium existed on the droplet surface during ionization even at low LF concentration [21]. Presence of LF-D<sub>9</sub> shifted the equilibration point. It would be interesting to test non-aqueous mobile phase elution. We also noticed detection signal varied significantly from day to day, as demonstrated by IS levels, suggesting that matrix accumulated in the column affected subsequent analysis. Use of co-eluting LF-D<sub>9</sub> as the IS was critical to compensate effectively for the cumulative ME.

### Crossvalidation

An HPLC-UV method for LF has been developed in our laboratory for a previous drug–drug interaction study on healthy subjects [7]. Here we reanalyzed 96 of the clinical samples. When the results were compared with those from the HPLC-UV method, only one of 96 samples has >20% deviation from the reference value (Table 7), suggesting the two methods are equivalent. Notably, the IS levels were varied significantly during the last batch of analysis (Supplementary Figure 3), but the results were still in good agreement with the reference values.

### Clinical application

This method was applied to clinical samples from pediatric patients from whom ~200  $\mu$ l blood was collected by capillary tubes. So far, 259 samples were analyzed; 27 samples, including a few predose samples were below LLOQ and eight samples were above 10,000 ng/ml with the highest concentration at 28,000 ng/ml. The precision of QC samples during the period of sample analysis was 9.35–10.8%. The detailed analysis of data will be reported elsewhere.

### Conclusion

A robust LC–MS/MS method was developed for determination of LF in human plasma. With LF-D<sub>9</sub> as the IS, this method was significantly improved on the issues of ME and

signal saturation when compared with previously published methods. Through LLE of sample and optimization of LC conditions, absolute ME was minimized to 92–100%. Cumulative matrix was still observed, causing significant variation of signals over days of use. The co-eluting LF-D<sub>9</sub> as the IS was able to compensate effectively for the cumulative ME. Furthermore, ionization saturation was also compensated by the deuterated IS, and this extended the linear range of calibration curve to 50–20,000 ng/ml. Only 25 µl plasma volume is required for analysis and is suitable for pharmacokinetic study of LF in pediatric patients.

## Future perspective

In respect to the difficulty of blood sample collection and limited sample volume from pediatric patients, we believe the less the sample volume, the better the method. A quantification method using less than 50 µl plasma will facilitate blood collection with capillary tubes. In the future, simple finger pricking and blood collection with capillary tube combined with a small plasma volume quantification method or dried blood spot quantification method will likely be the standard for pediatric pharmacokinetic studies.

## Supplementary Material

Refer to Web version on PubMed Central for supplementary material.

## Key Terms

<b>Matrix effect</b>	Concomitant presence of other substances in the ion source affects ionization of the analyte leading to enhanced or suppressed signal of the analyte
<b>Stable isotope-labeled analyte</b>	Having the same chemical structure with the analyte except for certain atoms (usually H and/or C) substituted with corresponding stable isotopes ( <sup>2</sup> H or/and <sup>13</sup> C)
<b>Ionization saturation</b>	Occurred when high concentration of compounds presented in ion source competing for ionizing agent (H <sup>+</sup> , NH <sub>4</sub> <sup>+</sup> )

## References

Papers of special note have been highlighted as:

- of interest
  - of considerable interest
1. White NJ, van Vugt M, Ezzet F. Clinical pharmacokinetics and pharmacodynamics of artemether-lumefantrine. *Clin Pharmacokinet.* 1999; 37:105–125. [PubMed: 10496300]
  2. German P, Parikh P, Lawrence J, et al. Lopinavir/ritonavir affects pharmacokinetic exposure of artemether/lumefantrine in HIV-uninfected healthy volunteers. *J Acquir Immune Defic Syndr.* 2009; 51:424–429. [PubMed: 19506482]
  3. Zeng M, Lu Z, Yang S, et al. Determination of benflumetol in human plasma by reversed-phase high-performance liquid chromatography with ultraviolet detection. *J Chromatogr B.* 1996; 681:299–306.
  4. Khalil IF, Abildrup U, Alifrangis LH, et al. Measurement of lumefantrine and its metabolite in plasma by high performance liquid chromatography with ultraviolet detection. *J Pharm Biomed Anal.* 2011; 54:168–172. [PubMed: 20832961]



5. Lindegårdh N, Annerberg A, Blessborn D, Bergqvist Y, Day N, White NJ. Development and validation of a bioanalytical method using automated solid-phase extraction and LC–UV for the simultaneous determination of lumefantrine and its desbutyl metabolite in plasma. *J Pharm Biomed Anal.* 2005; 37:1081–1088. [PubMed: 15862688]
6. Annerberg A, Singtoroj T, Tipmanee P, White NJ, Day NP, Lindegårdh N. High-throughput assay for the determination of lumefantrine in plasma. *J Chromatogr B.* 2005; 822:330–333.
7. Huang L, Lizak PS, Jayewardene AL, Marzan F, Lee MT, Aweeka FT. A modified method for determination of lumefantrine in human plasma by HPLC–UV and combination of protein precipitation and solid-phase extraction: application to a pharmacokinetic study. *Anal Chem Insights.* 2010; 5:15–23. Describes the reference method for crossvalidation. [PubMed: 20448843]
8. Wong RPM, Salman S, Ilett KF, Siba PM, Mueller I, Davis TME. Desbutyl-lumefantrine is a metabolite of lumefantrine with potent in vitro antimalarial activity that may influence artemether–lumefantrine treatment outcome. *Antimicrob Agents Chemother.* 2011; 55:1194–1198. [PubMed: 21199927]
9. Blessborn D, Romsing S, Annerberg A, et al. Development and validation of an automated solid-phase extraction and liquid chromatographic method for determination of lumefantrine in capillary blood on sampling paper. *J Pharm Biomed Anal.* 2007; 45:282–287. [PubMed: 17719735]
10. Ntale M, Ogwal-Okeng JW, Mahindi M, Gustafsson LL, Beck O. A field-adapted sampling and HPLC quantification method for lumefantrine and its desbutyl metabolite in whole blood spotted on filter paper. *J Chromatogr B.* 2008; 876:261–265.
11. Ashley EA, Stepniewska K, Lindegårdh N, et al. Pharmacokinetic study of artemether–lumefantrine giving once daily for the treatment of uncomplicated multidrug-resistant falciparum malaria. *Trop Med Int Health.* 2007; 12:201–208. [PubMed: 17300626]
12. Hodel EM, Zanolari B, Mercier T, et al. A single LC–tandem mass spectrometry method for the simultaneous determination of 14 antimalarial drugs and their metabolites in human plasma. *J Chromatogr B.* 2009; 877:867–886. First reported LC–MS/MS method for lumefantrine: protein precipitation was used and matrix effect was significant.
13. Wahajuddin W, Singh SP, Jian GK. Determination of lumefantrine in rat plasma by liquid–liquid extraction using LC–MS/MS with electrospray ionization: assay development, validation and application to a pharmacokinetic study. *J Chromatogr B.* 2009; 877:1133–1139.
14. Munjal V, Paliwal N, Chaurasia BK, Varshney B, Ahmed T, Paliwal J. LC-tandem mass spectrometry method for quantification of lumefantrine in human plasma and its application to bioequivalence study. *Chromatographia.* 2010; 71:505–510. Low recovery from protein precipitation.
15. César IC, de Aquino Ribeiro AJ, Teixeira LS, et al. Liquid chromatography-tandem mass spectrometry for the simultaneous quantitation of artemether and lumefantrine in human plasma: application for a pharmacokinetic study. *J Pharm Biomed Anal.* 2011; 54:114–120. Quadratic fitting of calibration curve and protein precipitation. [PubMed: 20719459]
16. ACTG guidelines for method development and validation based on (and including) FDA guidelines dated by 2001. Version 2 (2005).
17. Matuszewski BK, Constanzer ML, Chavez-Eng CM. Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC–MS/MS. *Anal Chem.* 2003; 75:3019–3030. [PubMed: 12964746]
18. Fast DM, Kelley M, Viswanathan CT, et al. 2010 White paper on recent issues in regulated bioanalysis and global harmonization of bioanalytical guidance. *AAPS J.* 2009; 11(2):238. [PubMed: 19381839]
19. Lindegårdh N, Annerberg A, White NJ, Day NP. Development and validation of a liquid chromatographic-tandem mass spectrometric method for determination of piperazine in plasma: stable-isotope-labeled internal standard does not always compensate for matrix effects. *J Chromatogr B.* 2008; 862:227–236.
20. Colussi D, Parisot C, Legay F, Lefevre G. Binding of artemether and lumefantrine to plasma proteins and erythrocytes. *Eur J Pharm Sci.* 1999; 9:9–16. [PubMed: 10493991]

21. Trufelli H, Palma P, Famiglini G, Cappiello A. An overview of matrix effects in liquid chromatography-mass spectrometry. *Mass Spectrom Rev.* 2011; 30:491–509. [PubMed: 21500246]

## Executive summary

### Method development


- Acidifying plasma samples prior to sample preparation is important to achieve a high recovery.
- Matrix effect (ME) was reduced to minimal through sample preparation (liquid-liquid extraction) and LC separation (gradient elution followed with 100% organic washing phase).
- Stable isotope-labeled lumefantrine used as the internal standard compensated for ME and ionization saturation.

### Method validation


- Cumulative ME observed from day to day was compensated by a co-eluting stable isotope-labeled internal standard.

### Application


- A small plasma sample volume (25  $\mu$ l) used in the method was suitable for pediatric studies where blood sample collection with capillary tubes may be used.




**Figure 1.**  
Lumefantrine and the deuterated lumefantrine.



**Figure 2. Calibration standard curves with and without deuterated internal standard**  
**(A)** No IS, calibration curve fitted with  $1/x$  weighted quadratic regression. **(B)** Deuterated LF (LF-D<sub>9</sub>) as the IS, calibration curve fitted with  $1/x$  weighted linear regression. **(C)** IS signal intensity at the different calibrator concentrations. The LF-D<sub>9</sub> signal decreased at higher LF concentration, indicating signal saturation at the high end of calibrators. The ratio of LF/LF-D<sub>9</sub> still increased in a linear mode.  
 LF: Lumefantrine; IS: Internal standard.



**Figure 3.** Chromatograms of (A) lumefantrine and (B) the internal standard (deuterated lumefantrine) for blank plasma, LLOQ, and the blank plasma after ULOQ.



**Figure 4. Matrix effect of plasma extracts from different sample preparations evaluated via an infusion experiment**

The peak at 2.1 min was lumefantrine. A 10  $\mu\text{g}/\text{ml}$  of lumefantrine infused constantly into the LC elute via a 'T' joint and directed into MS source while the LC ran at 0.4 ml/min at isocratic mode (65%B). Signal was recorded after injected 10  $\mu\text{l}$  of blank plasma extract from L-L extraction (blue line) and PPT (red line). A 10  $\mu\text{l}$  of 10  $\mu\text{g}/\text{ml}$  lumefantrine in mobile phase solvent injected as reference (green line). L-L: Liquid-liquid extraction; PPT: Protein precipitation.

Table 1

Intra- and inter-day precision and accuracy.

Nominal (ng/ml)	Mean (ng/ml)	SD	RSD (%)	%dev	n
<i>Intra-day</i>					
50	51.0-58.5	0.79-2.80	1.3-5.2	2.1-17	5
120	119-131	4.2-8.1	3.5-6.3	-0.5-9.3	5
1500	1560-1646	35.6-160	2.2-10	4.0-9.7	5
17,000	16,660-18,220	709-1260	4.3-7.1	-2.0-7.2	5
<i>Inter-day</i>					
50	54.16	3.21	5.9	8.3	30
120	127.2	6.7	5.3	6	30
1500	1607	88.7	5.5	7.1	30
17,000	17,593	1067	6.1	3.5	30

%dev: Deviation from nominal value.

RSD: Relative standard deviation.



**Table 2**

Stability of lumefantrine (n = 3).

Stability test conditions	Conc. (ng/ml)	% change	SD
Four freeze–thaw cycles, plasma	120	1.4	7
	170,000	1.2	1249
RT (22–28°C), 3 days, plasma	120	11.9	4
	170,000	–1.6	404
RT (22°C), 72 h, autosampler	120	4.7	3.0
	170,000	–3.3	814
–70°C, 1 year, plasma	120	10	5
	170,000	7.0	796
–70°C, 1.5 year stock in 0.5% FA, 50% MeCN		–0.53	0.1

FA: Formic acid; MeCN: Acetonitrile; RT: Room temperature.

**Table 3**

Precision (%CV) of peak areas and peak area ratios in set 1–3.

	Low QC (conc.: 120 ng/ml)	Med. QC (conc.: 1500 ng/ml)	High QC (conc.: 17,000 ng/ml)	Column
<i>Lumefantrine peak area</i>				
Set 1b	2.8	0.8	1	A
Set 2	2.8	0.7	1.4	B
Set 3	4.3	3.3	2.5	C
<i>IS peak area</i>				
Set 1b	1.5	1.2	2.6	D
Set 2	2.5	1.3	2.2	E
Set 3	4.3	2	5.1	F
<i>Peak area ratio</i>				
Set 1b	2.2	1.7	2.4	G
Set 2	4.1	1.7	2.2	H
Set 3	5.4	2.1	4.8	I

n = 6.

IS: Internal standard; QC: Quality control.

**Table 4**

Relative matrix effect.

	Low QC (conc.: 120 ng/ml)	Med. QC (conc.: 1500 ng/ml)	High QC (conc.: 17,000 ng/ml)	%CV dif.
B-A	0	-0.1	0.4	<5
E-D	1	0.1	-0.4	<5
H-G	1.9	0	-0.2	<5
I-H	1.3	0.4	2.6	<5

n = 6.

IS: Internal standard; QC: Quality control.

**Table 5**

Precision of slopes of lines fitted through low, medium and high quality control samples (n = 6).

Plasma lot #	Slope		
	Set 1b	Set 2	Set 3
1	0.0079	0.0075	0.0067
2	0.0075	0.0078	0.0075
3	0.0074	0.0076	0.0073
4	0.0076	0.0078	0.0067
5	0.0073	0.0080	0.0071
6	0.0076	0.0076	0.0068
Mean	0.00755	0.00772	0.00702
SD	0.00021	0.00018	0.00034
%CV	2.75	2.38	4.80

**Table 6**

Matrix effect, recovery and process efficiency (n = 6).

	Low conc. (120 ng/ml)	Med conc. (1500 ng/ml)	High conc. (17,000 ng/ml)
<b>LF mean peak area × 10e4</b>			
Set 1b	4.72	57	492
Set 2	4.62	57.2	491
Set 3	3.89	46.4	409
<b>IS mean peak area × 10e4</b>			
Set 1	5.42	4.96	3.84
Set 2	4.99	4.81	3.72
Set 3	4.41	4.28	3.44
<b>ME</b>			
LF	97.9	100.4	99.7
IS	92.1	96.9	97.1
<b>RE</b>			
LF	84.3	81	83.4
IS	88.3	88.9	92.3
<b>PE</b>			
LF	82.5	81.3	83.1
IS	81.4	86.2	89.6

IS: Internal standard; LF: Lumefantrine; ME: Matrix effect; PE: Process efficiency; RE: Recovery.

**Table 7**

Comparison of lumefantrine concentrations in clinical samples using the HPLC-UV method (conc. 1) versus LC-MS/MS method (conc. 2).

Sample #	Conc. 1	Conc. 2	%dev
S1	9303	10,900	15.8
S2	9192	10,700	15.2
S3	10,623	11,500	7.93
S4	10,539	11,600	9.59
S5	11,820	12,800	7.96
S6	13,151	15,400	15.8
S7	11,217	11,500	2.49
S8	8236	8400	1.97
S9	3826	4440	14.9
S10	2289	1980	-14.5
S11	1356	1380	1.79
S12	940	1000	6.17
S13	699	735	5.09
S14	437	451	3.15
S15	352	351	-0.28
S16	310	308	-0.65
S17	5339	5800	8.28
S18	4762	4860	2.04
S19	5055	5410	6.79
S20	5374	5820	7.97
S21	4952	5340	7.54
S22	5777	6280	8.34
S23	4641	4750	2.33
S24	3047	3210	5.22
S25	1356	1230	-9.74
S26	632	613	-2.99
S27	311	272	-13.3
S28	247	254	2.98
S29	184	203	9.78
S30	119	111	-6.74
S31	100	89.4	-11.1
S32	103	99.4	-3.31
S33	13,893	14,900	7.00
S34	13,306	15,700	16.5
S35	13,359	14,700	9.56
S36	12,457	13,800	10.2

Sample #	Conc. 1	Conc. 2	%dev
S37	14,193	15,800	10.7
S38	18,966	21,300	11.6
S39	16,825	18,300	8.40
S40	11,167	12,900	14.4
S41	8082	8660	6.91
S42	4734	4910	3.64
S43	2728	2820	3.30
S44	1766	1900	7.34
S45	1272	1420	11.0
S46	808	841	3.98
S47	623	665	6.52
S48	524	540	3.00
S49	17,221	18,900	9.29
S50	16,701	16,400	-1.82
S51	16,471	17,000	3.16
S52	15,334	15,900	3.63
S53	14,836	15,200	2.42
S54	14,949	15,200	1.66
S55	12,217	15,100	21.1
S56	9168	10,900	17.3
S57	4512	4720	4.50
S58	1810	1920	5.92
S59	984	1060	7.43
S60	700	762	8.52
S61	498	539	7.91
S62	328	336	2.41
S63	272	288	5.86
S64	223	219	-1.91
S65	9508	9490	-0.19
S66	9680	8430	-13.8
S67	9061	9200	1.53
S68	8903	8800	-1.16
S69	8652	8550	-1.19
S70	15,601	14,900	-4.60
S71	14,278	14,400	0.85
S72	10,467	10,600	1.26
S73	4348	4380	0.74
S74	2612	2550	-2.40
S75	1508	1640	8.39

Sample #	Conc. 1	Conc. 2	%dev
S76	1167	1230	5.22
S77	910	951	4.44
S78	516	535	3.65
S79	491	476	-3.10
S80	270	268	-0.91
S81	10,022	10,600	5.61
S82	9792	9760	-0.33
S83	9293	9850	5.82
S84	8641	9200	6.27
S85	8002	9120	13.1
S86	19,647	16,900	-15.0
S87	15101	13,800	-9.00
S88	9589	8960	-6.78
S89	2484	2460	-0.98
S90	1250	1230	-1.64
S91	728	719	-1.18
S92	430	434	0.93
S93	332	315	-5.15
S94	279	265	-5.05
S95	196	185	-5.77
S96	191	184	-3.73

%dev: % deviation.