

Fractional Carbon Dioxide Laser Resurfacing

William M. Ramsdell, M.D.¹

¹ Private practice, Austin, Texas

Semin Plast Surg 2012;26:125–130.

Address for correspondence and reprint requests William M. Ramsdell, M.D., 102 Westlake Drive, Suite 100, Austin Texas, 78746 (e-mail: wmr@centexderm.com).

Abstract

Keywords

- ▶ fractional
- ▶ ablative
- ▶ collagen
- ▶ CO₂

Currently available ablative fractional CO₂ lasers provide excellent results and diminish down time with fewer complications than previous generation CO₂ lasers. Mechanisms of action, treatment parameters, as well as pre- and postoperative care will be discussed.

The sophisticated lasers used today might seem unimaginable to Dr. Leon Goldman who first studied laser-cutaneous interactions in the 1960s.¹ In 1983, Anderson and Parrish fundamentally changed our usage of lasers with the groundbreaking concept of photothermolysis via pulsed radiation.² During the 1990s, the ability to apply high energy with short-pulse durations made what is now termed “traditional” CO₂ very popular. Never before had a device been able to produce such changes in the appearance and function of skin. Many physicians consider the fully ablative traditional CO₂ laser to be the gold standard. Certainly, no other device can achieve the same results. Those results, however, come with side effects such as prolonged erythema and delayed hypopigmentation for some patients. Dr. Jeffrey Dover believes that today there are few instances where traditionally ablative CO₂ laser resurfacing is indicated (Personal Communication, January 16, 2012).

In 2004, Manstein et al developed a technique of creating thermal injury to skin in a fractional pattern in an attempt to achieve results similar to fully ablative resurfacing while maintaining an acceptable side-effect profile.³ Following initial development of nonablative fractional lasers, ablative fractional CO₂ lasers entered the market in hopes of delivering even better results.

Mechanisms of Actions

The benefits of CO₂ laser resurfacing are undeniable, whether fully ablative or fractional. The mechanisms by which such improvement is achieved have been the subject of much debate.⁴ Although fully ablative traditional and fractional CO₂ lasers have distinctive characteristics, much of the re-

search on fully ablative CO₂ resurfacing is applicable to fractional resurfacing.⁵ All CO₂ lasers have a wavelength of 10,600 nm, heavily absorbed by water. To achieve ablation without excessive thermal damage, a fluence of 5 J/cm² must be delivered within a pulse duration of < 1 ms, the generally agreed upon thermal relaxation time of skin.

The denaturation temperature of collagen is conventionally stated to be 66.8°C, although this varies.⁴ Once denatured by laser-generated heat, collagen rapidly contracts as fibers shrink to one-third of their length.⁶ The shrinkage of collagen is the primary mechanism of skin tightening, although vaporization of intracellular water and ablation contribute as well. Next, a wound-healing phase is initiated characterized by extremely high levels of collagenases (matrix metalloproteinases), which degrade the fragmented collagenous matrix.⁷ Rapid reconstitution of the epidermis from adjacent epidermal cells contrasts with healing after traditional resurfacing in which new epidermis is derived from cells that migrate from adnexal structures. A prolonged period of dermal neocollagenesis of up to at least 6 months follows.⁸

Although tissue ablation certainly plays a role in effectiveness, the depth of deep rhytides exceeds the depth of ablation.^{9,10} Residual thermal damage represents the true depth of injury and the primary factor leading to effectiveness.⁴ Residual thermal damage is increased with higher fluence,^{11,12} whereas increased density results in greater depth of ablation.

Patient Selection

As with any cosmetic procedure, patient selection is important. Realistic expectations and a thorough understanding of

the procedure and especially postoperative care are crucial. Although most patients require only one treatment, patients should be prepared for the possibility that more than one treatment may be necessary for best results. Contraindications to treatment include active cutaneous infection, diseases that inhibit healing, such as scleroderma, and a history of any diseases that Koebnerize such as vitiligo and psoriasis. Recent isotretinoin therapy or a history of radiation at the treatment site demands caution and nonaggressive parameters.

Preoperative Care

Regardless of history, all patients receive preoperative antiviral prophylaxis with valacyclovir 500 mg twice daily beginning one day preoperatively and continued for 10 days. Most laser surgeons also prescribe an antibiotic such as cefadroxil, dicloxacillin, doxycycline, or ciprofloxacin during the same period.¹³ Others believe that antibiotics for Gram+ organisms only promote Gram infection. Fluconazole 200 mg is given on the day of surgery. Prophylactic facial cleansing with chlorhexidine twice daily for 3 days preoperatively may further decrease the incidence of postoperative infection. Some surgeons prescribe intranasal mupirocin preoperatively. Preconditioning the skin with tretinoin in an attempt to speed healing and reduce the incidence of postoperative milia and acne is of questionable benefit and promotes increased postoperative erythema.¹⁴⁻¹⁸

Treatment

Fractional resurfacing begins with appropriate laser precautions and sedation of choice. Although general or tumescent anesthesia may be used, the liberal application of topical anesthesia 1 to 2 hours before surgery is usually adequate. Kilmer believes that a hydrophilic topical anesthetic such as EMLA (eutectic mixture of lidocaine and prilocaine; AstraZeneca Pharmaceuticals, LP, London, UK) may alter the response to the laser ensuring fewer side effects.¹⁹ Because the primary target of CO₂ lasers is water, this makes intuitive sense. However, Naouri et al believe the water content of topical anesthetics irrelevant.²⁰ Certainly, topical anesthesia should have very little effect on a narrow, deep ablative beam. Whether more superficial resurfacing is affected by topical anesthesia is open to debate. In addition to topical anesthesia, forced cold air provided by a Zimmer Cooler (LaserMed, Monroe, CT) or similar device is helpful.

Currently, there are several excellent fractional CO₂ devices on the market.²¹ This discussion will be limited to the specifics of using the Lumenis Ultrapulse Encore fractional CO₂ laser for photoaging and scarring of the face and neck. Two delivery devices offer dual depth ablation, Active FX and Deep FX, via separate handpieces. The Active FX handpiece has a collimated 1.3-mm spot size delivered in a nonsequential array (Cool Scan mode) to minimize thermal damage. Available energy varies from 2 mJ (150 mJ/cm²) to 225 mJ (169 J/cm²). Power settings range from 1 W (4.4 Hz) to 60 W (600 Hz). A computer pattern generator delivers seven patterns in nine different size options

and nine density options. It is the ability to alter the density that distinguishes fractional lasers from traditionally fully ablative ones. For density settings, 1 to 3 with the square pattern, density varies between 55% and 82%. Densities 4 to 9 are fully ablative (100%) although the depth is much less than with traditional CO₂ lasers.

Active FX creates a shallow, broad ablative crater that can extend into superficial papillary dermis. One or two passes are commonly used. Typical energy settings for facial skin range from 100 to 125 mJ. Density settings from 2 to 4 are typical, but may be increased for difficult cases such as acne scarring. Periorbital skin treatment parameters range from 60 to 90 mJ, density 2 to 3. Treatment of neck, chest, and extremities is possible at reduced settings such as 70 to 80 mJ, density 1 to 2. Opinion varies regarding whether to debride between laser passes for traditional resurfacing.²²⁻²⁵ For fractional resurfacing, it is not necessary.

For deep rhytides, scarring and for maximum collagen regeneration, Deep FX is used initially or exclusively. The Deep FX handpiece produces a narrower, noncollimated beam at 0.12 mm with four pattern and six size options. One pulse or 2-stacked pulses can be used. Power ranges up to 600 Hz with a single pulse and to 300 Hz with double pulsing. Energy settings between 2.5 mJ and 50 mJ are available. At 2.5 mJ, the ablation depth is quite shallow, ~76 μm, which extends through the epidermis and possibly to superficial papillary dermis depending on site. Although the depth can be increased to 1600 μ or more, these settings are much too powerful to be clinically appropriate.

A Deep FX energy setting between 15 to 22.5 mJ is typical for the face. At 22.5 mJ, ablation extends to ~675 μm, which correlates with deep dermis (depending on skin thickness). This depth will encompass most cases of solar elastosis. Density settings of 5 to 25% are available. This translates for a square scan area of 10 × 10 mm, from 196 spots at density 5% to 841 spots at 25%. Deep FX at aggressive settings can punch all of the way through eyelid skin. Safe settings for skin within the bony orbit include settings from 8 to 10 mJ and density 5 to 15%. Settings of 8 to 10 mJ, density 5 to 15% are appropriate for off of the face sites such as neck, chest, and extremities. The author uses one pulse only for all applications except scarring, for which double pulsing is sometimes used.

Oni et al documented that, on facial skin, a single pulse at 15 mJ extended to an average 416 μm, which correlates with mid-dermis.²⁶ Interestingly, double pulsing at 15 mJ was histologically equivalent to a single pulse at 30 mJ, extending 881 μm and 854 μm, respectively (deep reticular dermis). The relevance of this in terms of clinical outcome has yet to be determined. Little data exists in the peer-reviewed literature regarding how best to combine Active FX with Deep FX.

Because the precise boundaries of the scanned imprints on the skin can be difficult to see, skip areas, or worse, pulse stacking can occur with Deep FX resurfacing. Because treated skin will turn red or urticate several minutes after treatment, untreated areas can be recognized and treated. If Active FX treatment is planned following Deep FX, areas missed with Deep FX will usually blend in.



Figure 1 A 66-year-old woman (A) Before and (B) after Active FX resurfacing and transconjunctival lower blepharoplasty. Active FX settings: two passes, 100 mj, 60 W, pattern 1, size 9, density 6; periorbital areas: two passes, 80 mj, 45 W, pattern 5, size 6, density 4 and density 3; feathering with handpiece held tangentially: 80 mj, 45 W, pattern 1, size 9, density 3.

Postoperative Care

Meticulous postoperative care is absolutely crucial to avoid infection and scarring. Cold soaks of white vinegar 1 tsp to 1 cup water or aluminum acetate (Domeboro, Bayer Health-Care, Morristown, NJ) soaks reduce erythema, edema, and provide antisepsis. The skin should be kept moist with white petrolatum or Aquaphor® (Beiersdorf, Wilton, CT). Other topicals are strongly discouraged. The concept of promoting rapid healing via use of autologous platelet-rich plasma has recently been introduced.²⁷

Occlusive dressings such as Silon II (BioMed Sciences, Inc, Allentown, PA) have been shown to speed healing, reduce inflammation, and possibly improve results following traditional resurfacing.^{4,28,29} The dressings are cumbersome and not necessary for most fractional resurfacing. The author still uses Silon II with aggressive treatments for severe acne scarring. Reepithelialization usually takes 3 to 6 days. At that point, a light moisturizer and sunscreen may be used. Mild to moderate erythema persists for several weeks. Topical vitamin C can be used to encourage neocollagenesis and reduce erythema.^{30,31}



Figure 2 Combination ice-pick and box-car-type scarring. (A) Before and (B) after punch excisions followed by two passes of Active FX at 125 mj, density 6. Periorbital areas: one pass Active FX at 80 mj, density 3.



Figure 3 Box-car-type acne scarring. (A) Before and (B) after fractional resurfacing. One pass of Deep FX at 22.5 mj, 10% density, one pulse, followed by Active FX at 125 mj, density 5. Periorbital areas: one pass Active FX at 80 mj, density 2.

Efficacy

Efficacy can be difficult to assess because histologic changes do not always correlate with changes in appearance. Moreover, changes in appearance are difficult to quantify and subjective for both patient and physician. Virtually all patients show at least modest improvement (►Fig. 1). In general, thin skin responds better than thick. Young skin is capable of producing more collagen than old skin. Rahman

found moderate or better improvement in 83% of patients.³² Ortiz documented 50 to 75% improvement at 3 months with 25 to 50% improvement long term.³³ Improvement of eyelid rhytides by 53% has been documented as well.³⁴

Acne scarring is more difficult to treat and realistic expectations should be stressed during the consultation. After 2 to 3 treatments for acne scarring, 26 to 50% improvement in texture and atrophy has been noted.³⁵ The same study documented 68% improvement in the depth of scars as measured by



Figure 4 Rolling scars. (A) Before and (B) after treatment, showing only modest improvement despite aggressive settings. Punch excisions were followed by both Deep FX and Active FX. Deep FX: 125 mj, density 20%, two pulses to most severely affected areas. Active FX: 1 pass at 125 mj, density 7 to entire face except periorbital areas: 80 mj, density 4. Silon II was placed as a wound dressing for 3 days to speed healing and reduce inflammation.



Figure 5 Facelift scar. (A) Before and (B) after fractional resurfacing. Two passes of Active FX at 90 mj, density 4.

the Primos three-dimensional optical system. Ice pick-type scars should undergo punch excisions prior to laser therapy (►Fig. 2). Boxcar-type scarring responds well (►Fig. 3). Most challenging to treat are rolling scars (►Fig. 4). These patients may benefit from a facelift superior to the superficial muscular aponeurotic system (SMAS) first. This stretches the skin, thereby effacing some of the rolling depressions. The increased tension placed on fibroblasts may enhance subsequent collagen formation.⁷ Laser settings for rolling scars should be aggressive and many patients will require more than one treatment. Regardless of the etiology of scarring, thin scars respond better than thick ones (►Fig. 5).

Conclusion

Fractional CO₂ resurfacing represents one of the most useful devices produced by 50 years of research into laser-cutaneous interactions. Patients can expect excellent results with few complications. Additional research will further improve our ability to provide the best of care for our patients.

References

- Goldman L, Wilson RG, Hornby PM, Meyer RG. Radiation from a Q-switched ruby laser. Effect of repeated impacts of power output of 10 megawatts on a tattoo of a man. *J Invest Dermatol* 1965;44:69-71
- Anderson RR, Parrish JA. Selective photothermolysis: precise microsurgery by selective absorption of pulsed radiation. *Science* 1983;220(4596):524-527
- Manstein D, Herron GS, Sink RK, Tanner H, Anderson RR. Fractional photothermolysis: a new concept for cutaneous remodeling using microscopic patterns of thermal injury. *Lasers Surg Med* 2004;34(5):426-438
- Ross EV, McKinlay JR, Anderson RR. Why does carbon dioxide resurfacing work? A review. *Arch Dermatol* 1999;135(4):444-454
- Reilly MJ, Cohen M, Hokugo A, Keller GS. Molecular effects of fractional carbon dioxide laser resurfacing on photodamaged human skin. *Arch Facial Plast Surg* 2010;12(5):321-325
- Gardner ES, Reinisch L, Stricklin GP, Ellis DL. In vitro changes in non-facial human skin following CO₂ laser resurfacing: a comparison study. *Lasers Surg Med* 1996;19(4):379-387
- Fisher GJ, Varani J, Voorhees JJ. Looking older: fibroblast collapse and therapeutic implications. *Arch Dermatol* 2008;144(5):666-672
- Orringer JS, Kang S, Johnson TM, et al. Connective tissue remodeling induced by carbon dioxide laser resurfacing of photodamaged human skin. *Arch Dermatol* 2004;140(11):1326-1332
- Fitzpatrick R, Smith S, Sriprachya-anunt S. Depth of vaporization with the UltraPulse CO₂ laser. *Proc SPIE* 1997;2970:392-394
- Ross EV, Grossman MC, Duke D, Grevelink JM. Long-term results after CO₂ laser skin resurfacing: a comparison of scanned and pulsed systems. *J Am Acad Dermatol* 1997;37(5 Pt 1):709-718
- Bernstein EF, Brown DB, Kenkel J, Burns AJ. Residual thermal damage resulting from pulsed and scanned resurfacing lasers. *Dermatol Surg* 1999;25(10):739-744
- Saluja R, Khoury J, Detwiler SP, Goldman MP. Histologic and clinical response to varying density settings with a fractionally scanned carbon dioxide laser. *J Drugs Dermatol* 2009;8(1):17-20
- Ross EV, Amesbury EC, Barile A, Proctor-Shipman L, Feldman BD. Incidence of postoperative infection or positive culture after facial laser resurfacing: a pilot study, a case report, and a proposal for a rational approach to antibiotic prophylaxis. *J Am Acad Dermatol* 1998;39(6):975-981
- Lowe NE, Lask G, Griffin ME. Laser skin resurfacing. Pre and post treatment guidelines. *Dermatol Surg* 1995;21:1017-1019
- Goldman MP. Pre- and postoperative care of the laser patient. *Int J Aesthetic Restorative Surg* 1997;5:46-49
- Gilchrist BA. Treatment of photodamage with topical tretinoin: an overview. *J Am Acad Dermatol* 1997;36(3 Pt 2):S27-S36
- Hung VC, Lee JY, Zitelli JA, Hebda PA. Topical tretinoin and epithelial wound healing. *Arch Dermatol* 1989;125(1):65-69
- Nyirady J, Grossman RM. Use of tretinoin in precosmetic and postcosmetic procedures. A review. *Cos Dermatol* 2003;16(5):39-47
- Kilmer SL, Chotzen VA, Silva SK, McClaren ML. Safe and effective carbon dioxide laser skin resurfacing of the neck. *Lasers Surg Med* 2006;38(7):653-657

- 20 Naouri M, Atlan M, Perrodeau E, et al. High-resolution ultrasound imaging to demonstrate and predict efficacy of carbon dioxide fractional resurfacing laser treatment. *Dermatol Surg* 2011;37(5):596–603
- 21 Waibel J, Beer K, Narurkar V, Alster T. Preliminary observations on fractional ablative resurfacing devices: clinical impressions. *J Drugs Dermatol* 2009;8(5):481–485
- 22 Ross EV, Mowlavi A, Barnette D, Glatter RD, Grevelink JM. The effect of wiping on skin resurfacing in a pig model using a high energy pulsed CO₂ laser system. *Dermatol Surg* 1999;25(2):81–88
- 23 David L, Ruiz-Esparza J. Fast healing after laser skin resurfacing. The minimal mechanical trauma technique. *Dermatol Surg* 1997;23(5):359–361
- 24 Ruiz-Esparza J, Barba Gomez JM, Gomez de la Torre OL, et al. Erythema after laser skin resurfacing. *Dermatol Surg* 1998;24:31–34
- 25 Niamtu J III. To debride or not to debride? That is the question: rethinking char removal in ablative CO₂ laser skin resurfacing. *Dermatol Surg* 2008;34(9):1200–1211
- 26 Oni G, Robbins D, Bailey S, Brown SA, Kenkel JM. An in vivo histopathological comparison of single and double pulsed modes of a fractionated CO₂ laser. *Lasers Surg Med* 2012;44(1):4–10
- 27 Na JI, Choi JW, Choi HR, et al. Rapid healing and reduced erythema after ablative fractional carbon dioxide laser resurfacing combined with the application of autologous platelet-rich plasma. *Dermatol Surg* 2011;37(4):463–468
- 28 Batra RS, Ort RJ, Jacob C, Hobbs L, Arndt KA, Dover JS. Evaluation of a silicone occlusive dressing after laser skin resurfacing. *Arch Dermatol* 2001;137(10):1317–1321
- 29 Goldman MP, Roberts TL III, Skover G, Lettieri JT, Fitzpatrick RE. Optimizing wound healing in the face after laser abrasion. *J Am Acad Dermatol* 2002;46(3):399–407
- 30 Alster TS, West TB. Effect of topical vitamin C on postoperative carbon dioxide laser resurfacing erythema. *Dermatol Surg* 1998;24(3):331–334
- 31 Pinnel SR, Murad S, Darr D. Induction of collagen synthesis by ascorbic acid. A possible mechanism. *Arch Dermatol* 1987;123(12):1684–1686
- 32 Rahman Z, MacFalls H, Jiang K, et al. Fractional deep dermal ablation induces tissue tightening. *Lasers Surg Med* 2009;41(2):78–86
- 33 Ortiz AE, Tremaine AM, Zachary CB. Long-term efficacy of a fractional resurfacing device. *Lasers Surg Med* 2010;42(2):168–170
- 34 Kotlus BS. Dual-depth fractional carbon dioxide laser resurfacing for periocular rhytidosis. *Dermatol Surg* 2010;36(5):623–628
- 35 Chapas AM, Brightman L, Sukal S, et al. Successful treatment of acneiform scarring with CO₂ ablative fractional resurfacing. *Lasers Surg Med* 2008;40(6):381–386