

# NIH Public Access

Author Manuscript

AJR Am J Roentgenol. Author manuscript; available in PMC 2013 September 01.

#### Published in final edited form as:

AJR Am J Roentgenol. 2012 September ; 199(3): W392–W401. doi:10.2214/AJR.11.7255.

# Assessing the Standalone Sensitivity of Computer-aided Detection (CADe) with Cancer Cases from the Digital Mammographic Imaging Screening Trial (DMIST)

# Elodia B. Cole, MS,

Medical University of South Carolina Department of Radiology and Radiological Sciences 96 Jonathan Lucas St., MSC 323 Charleston, SC 29425

# Zheng Zhang, PhD,

Brown University Center for Statistical Science Box G-S121-7 121 S. Main St. Providence, RI 02912

# Helga S. Marques, MS,

Brown University Center for Statistical Science Box G-S121-7 121 S. Main St. Providence, RI 02912

# Robert M. Nishikawa, PhD,

Department of Radiology and the Committee on Medical Physics The University of Chicago 5841 South Maryland Ave., MC2026 Chicago, IL 60637

# R. Edward Hendrick, PhD,

University of Colorado – Denver, School of Medicine Department of Radiology, C278 12700 E. 19<sup>th</sup> Ave. Aurora, CO 80045

# Martin J. Yaffe, PhD,

Sunnybrook Health Sciences Centre 2075 Bayview Ave., S-Wing, Room S657 Toronto, ON M4N 3M5 CANADA

# Wittaya Padungchaichote, MD,

Lopburi Cancer Center Lopburi 15000, Thailand

# Cherie Kuzmiak, DO,

University of North Carolina Department of Radiology CB #5120 170 Manning Dr. Chapel Hill, NC 27599

# Jatuporn Chayakulkheeree, MD,

Memorial Hospital Chulalongkorn University Bangkok 10330, Thailand

# Emily F. Conant, MD,

University of Pennsylvania Department of Radiology/1 Silverstein 3400 Spruce St. Philadelphia, PA 19104

# Laurie Fajardo, MD,

Department of Radiology Room 3970 JPP University of Iowa Hospitals and Clinics 200 Hawkins Dr. Iowa City, IA 52242

# Janet Baum, MD,

Cambridge Health Alliance Dept of Radiology 1493 Cambridge St Cambridge, MA 02139

<sup>&</sup>lt;sup>\*</sup>Corresponding Author: Elodia B. Cole M.S., Department of Radiology and Radiological Science, Medical University of South Carolina, 96 Jonathan Lucas St., Suite 210, MSC 323, Charleston, SC 29425-3230. Phone: (843) 792-6965, Fax: (843) 792-1889, coleeb@musc.edu.

#### Constantine Gatsonis, PhD, and

Brown University Center for Statistical Science Box G-S121-7 121 S. Main St. Providence, RI 02912

#### Etta Pisano, MD

Medical University of South Carolina Department of Radiology and Radiological Sciences 96 Jonathan Lucas St., MSC 323 Charleston, SC 29425

# Abstract

**Objective**—To assess the sensitivities and false detection rates of two CADe systems when applied to digital or screen-film mammograms in detecting the known breast cancer cases from the DMIST breast cancer screening population.

**Materials and Methods**—Available screen-film and digital mammograms of 161 breast cancer cases from DMIST were analyzed by two CADe systems, iCAD SecondLook (iCAD) and R2 ImageChecker (R2). Three experienced breast imaging radiologists reviewed the CADe marks generated for each available cancer case, recording the number and locations of CADe marks and whether each CADe mark location corresponded with the known location of the cancer.

**Results**—For the 161 cancer cases included in this study, the sensitivities of the DMIST reader without CAD were 0.43 (69/161, 95% CI 0.35 to 0.51) for digital and 0.41 (66/161, 95% CI 0.33 to 0.49) for film-screen mammography. The sensitivities of iCAD were 0.74 (119/161, 95% CI 0.66 to 0.81) for digital and 0.69 (111/161, 95% CI 0.61 to 0.76) for screen-film mammogram, both significantly higher than the DMIST study sensitivities (p < 0.0001 for both). The average number of false CADe marks per case of iCAD was 2.57 (SD 1.92) for digital and 3.06(SD 1.72) for screen-film mammography. The sensitivity of R2 was 0.74 (119/161, 95% CI 0.66 to 0.81) for digital, and 0.60 (97/161, 95% CI 0.52 to 0.68) for screen-film mammography, both significantly higher than the DMIST study sensitivities (p < 0.0001 for both). The average number of false CADe marks per case of 1.52 to 0.68) for screen-film mammography, both significantly higher than the DMIST study sensitivities (p < 0.0001 for both). The average number of false case of R2 was 2.07 (SD 1.57) for digital and 1.52(SD 1.45) for screen-film mammogram.

**Conclusions**—Our results suggest the use of CADe in interpretation of digital and screen-film mammograms could lead to improvements in cancer detection.

# INTRODUCTION

In the Digital Mammographic Imaging Screening Trial (DMIST), the reported sensitivities for both digital and screen-film mammography were 0.41 based on 455 days of follow-up.<sup>1</sup> While low, similar sensitivities for digital and screen-film mammography have been previously reported, ranging from 0.40 to 0.77 for digital mammography and from 0.31 to 0.71 for screen-film mammography in recently published results from six other controlled laboratory studies,<sup>1-12,13</sup> while the National Breast Cancer Surveillance Consortium (NBCSC) reported an average screening sensitivity of 81.2% in their analysis of screening mammography records for over 3.8 million women between 1996 and 2007<sup>14</sup> based on 365 days of follow-up.<sup>3</sup> There is wide variation in reported radiologist sensitivities in these studies that is likely dependent not only on the skills of the participating radiologists, but also on the image case sets, the number of rounds of screening, and definitions of sensitivity<sup>15-18</sup>. Even so, there is a need for development of tools that improve the detection of breast cancer; thus the impetus for computer-aided detection software development. We use the abbreviation CADe to refer to computer-aided detection software, as opposed to the more ubiquitous CAD, which has been used to reference both computer-aided detection algorithms and computer-aided diagnostic algorithms. While the former is designed to prompt the radiologist to regions of a mammogram that may have been overlooked, the

latter, now frequently referred to as CADx in the literature, is designed to aid radiologists in distinguishing cancerous lesions from non-cancerous lesions.

It has been noted previously that the cases used to train and ultimately test a CADe system are the primary drivers on how well that technology will perform in cancer detection.<sup>19</sup>

While CADe publications have gone into great detail about the types and sizes of lesions included, the complexity of the dataset is hard to determine simply from these metrics. CADe systems were designed specifically to reduce the number of cancers missed by radiologists by marking suspicious regions of interest, prompting radiologists to inspect the annotated region. While the gold standard methodology for evaluation of a CADe system should involve radiologists reviewing clinical cases without and then with CADe, it is also useful to assess the cancer detection capability of the CADe system alone, without a radiologist reader. The study objectives detailed in this paper were: i) to estimate the sensitivity of each CADe system to known cancers when applied to digital and screen-film mammography; ii) to compare the standalone sensitivities of each CADe system to the original DMIST study sensitivities for digital and screen-film interpretation; iii) to report the average number of false positive marks per case for each CADe system; and, iv) to estimate the effect of tumor and patient characteristics on sensitivity.

# MATERIALS AND METHODS

Appropriate institutional review board approval was obtained prior to conducting this study. Informed consent and HIPAA consent were obtained from all evaluable subjects prior to image acquisition for DMIST. All data were handled in a HIPAA compliant manner. Digital and screen-film mammograms of the available DMIST pathologically-proven cancer cases, i.e., those detected within 455 days of acquisition of the paired digital and screen-film mammograms were included in this study. The screen-film and digital mammograms were obtained from the American College of Radiology Imaging Network (ACRIN) image archive. A total of 329 case pairs, including both digital and screen-film images, were available for inclusion in this study out of a total of 335 DMIST cancer cases.

#### CADe Systems Tested

The CADe systems for screen-film and digital mammography of two manufacturers were tested in this study: iCAD Secondlook v1.4 (digital and screen-film), R2 ImageChecker Cenova v1.0 (digital), and R2 ImageChecker v 8.0 (screen-film). The results from the default CAD sensitivity settings for each system were reported as the main results. In addition to the default CAD sensitivity setting for R2 screen-film, two additional CAD sensitivity settings were applied to a single digitized film set, and results are included in the appendix. While any film mammogram can be processed using any film CADe system, not all digital mammograms are suitable for processing with all digital CADe systems. At the time of this study, the iCAD Secondlook v1.4 used in the study was compatible with images from the Fuji CR, Hologic Selenia, and GE Senographe 2000D systems. The R2 ImageChecker Cenova CADe system tested was capable of analyzing images from the Fischer SenoScan, Fuji CR, Hologic Selenia, and GE Senographe 2000D digital mammography systems. The maximum available number of cases for each of our two digital CADe systems was limited by these compatibility constraints. In total, 245 cases were available for the iCAD digital system and all 329 cases were available for the R2 digital system.

#### **Film Preparation**

De-identified screen-film mammograms were digitized using manufacturer specified digitizer. CADe reports were printed to paper. Twenty-seven of the 329 available screen-film mammograms (8.2%) were rejected by the iCAD screen-film digitizer and therefore could not be digitized, resulting in 302 screen-film cases analyzed by iCAD. Twenty-six cases of 329 (7.9%) available screen-film mammograms were rejected by the R2 film digitizer and therefore could not be digitized, resulting in 303 screen-film cases analyzed by R2.

#### **Digital Image Preparation**

Digital mammograms were processed with CADe and displayed on manufacturer specific mammography workstations. The application of CADe algorithms was not successful for all digital images. In total, there were 179 evaluable digital cases for iCAD and 227 digital cases for R2. The CADe-generated electronic reports were displayed on a Hologic Secureview 6.0 workstation for R2 (Figure 1) and a Sectra IDS5.MX review workstation for iCAD (Figure 2).

#### **Common Cases by Both Vendors**

There are 161 common cases that were fully evaluated by both vendors, on both film screen and digital mammograms. We focused our data analysis on those common cases, and presented the patient and lesion characteristics of those cases on Table 1.

#### **CADe Scoring Procedure**

A radiologist with 26 years of mammography experience annotated acetate overlays with the screen-film mammograms for each case (Figure 3), denoting the locations of known cancers as seen on DMIST images and/or validated by DMIST pathology reports. The annotating radiologist specified whether cancer was visible for each modality.

Three different radiologists reviewed the generated CADe marks and compared these to the known cancer locations depicted on the overlays. All three readers completed breast imaging fellowships and had an average of 4 years (range 2-8 years) experience in mammography. An average of 48.3% (range 25%-70%) of their time was spent interpreting screening mammograms in their respective clinical practices. One radiologist scored the iCAD digital cases, a second radiologist scored the R2 digital cases, and a third radiologist scored all digitized screen-film CADe for both R2 and iCAD. The radiologists recorded the number of CADe marks generated, the location of each CADe mark, and determine whether the location of each CADe mark corresponded to a known cancer as recorded on pathology reports and the case overlay. Marks that did not correspond to a known cancer location were recorded as false positives by the radiologists.

#### **Statistical Methods**

All statistics were performed at the case level to allow for comparison to DMIST results. The study objectives were: i) to estimate the sensitivity of each CADe system to known cancers when applied to digital and screen-film mammography; ii) to compare the standalone sensitivities of each CADe system to the original DMIST study sensitivities for digital and screen-film interpretation; iii) to report the average number of false positive marks per case for each CADe system; and, iv) to estimate the effect of tumor and patient characteristics on sensitivity. Sensitivity and the associated 95% exact confidence interval for each CADe system were calculated and McNemar's test was used to test the difference between digital and screen-film within each CAD system. McNemar's test and 95% exact confidence intervals were used for the second objective. Mean, standard deviation, median,

and number of false positive marks per case were calculated for the third objective and paired t-test were used to test the difference between digital and screen-film within each CAD system. To address the fourth objective, sensitivities for subsets defined by age, breast density, menopausal status, cancer histology, tumor size, and lesion type were calculated by simple counts as well as univariate logistic regression-based estimates. 95% confidence intervals and p-values for comparing covariate effects on sensitivity were drawn from the same logistic regression model.

# RESULTS

We present the results for the default CADe sensitivity settings for screen-film and digital for iCAD and R2 for the 161 common cancer cases.

#### **Standalone Performance**

There were 161 common cases across both screen-film and digital modalities for each CAD vendor, allowing for modality comparison within each vendor. The overall sensitivity of iCAD digital was 119/161 (0.74) [95% CI: (0.66, 0.81)] and 111/161 (0.69) [95% CI: (0.61, 0.76)] for screen-film. This difference in sensitivity between iCAD digital and screen-film 0.05 [95% CI: (-0.03, 0.13) was not statistically significant (p=0.26). The average number of false CADe marks per case for iCAD digital ( $2.57\pm1.92$ ) was significantly lower than that for iCAD screen-film ( $3.06\pm1.72$ ) (p=0.0008). The overall sensitivity of R2 digital was 119/161 (0.74) [95% CI: (0.66, 0.81)] and 97/161 (0.60) [95% CI: (0.52, 0.68)] for R2 screen-film. The difference in sensitivity between R2 digital and R2 screen-film 0.14 [95% CI: (0.05, 0.23) was statistically significant (p=0.003). The average number of false CADe marks per case for R2 digital ( $2.07\pm1.57$ ) was significantly higher than for R2 screen-film ( $1.52\pm1.45$ ) (p <0.0001) (Table 2).

#### CADe Detection of Cancers Not Detected in the DMIST Study

In standalone mode, each CADe system for both screen-film and digital mammography was able to detect significantly more cancer cases than were detected by the original clinical radiologist readers in DMIST for the same cases. Of the 161 common cases, the original DMIST study has sensitivity of 0.43 (69/161, 95% CI 0.35 to 0.51) for digital and 0.41 (66/161, 95% CI 0.33 to 0.49) for screen-film. For digital mammogram, both CADe systems' sensitivity for digital (0.74) was significantly higher than the DMIST's sensitivity (p < 0.0001). For screen-film, sensitivities for iCAD(0.69) and R2(0.60) were both significantly higher than the DMIST's sensitivity (p < 0.0001 for both) (Table 3).

In addition, CADe system has shown ability to detect the cases that were not detected by any method in the DMIST study. Of the 46/161 cases that were detected by neither screen-film nor digital readers in the primary DMIST study, 25 (54.3%) were marked by the iCAD system when applied to digital, and 18(39.1%) when applied to screen-film. The comparable numbers for R2 system is 27(58.7%) for digital and 12 (26.1%) for screen-film.

#### Lesion and Subject Characteristics Impact on Sensitivity

CADe's sensitivity on lesion's histology type (DCIS or invasive), size and type (mass, calcification, asymmetric density and arch distortion), and the subject's age, breast density, and menopausal status was performed (Table 4) and compared using univariate logistic regression models (Table 5). CADe sensitivity did not appear to depend on subject's characteristics, although the sensitivity of R2 screen-film was lower on pre-menopausal women (49%) than on post-menopausal women (64%), but the difference was not significant (p=0.09). Lesion histology type and size didn't influence the sensitivity significantly, but lesion type seemed to have some effect. In general, CADe had the highest

sensitivity for calcification (83% for both iCAD and R2 digital, 73% for both iCAD and R2 screen-film), but R2 screen-film has a low sensitivity of 55% for mass, significantly lower than the 73% for calcification on the same system (p = 0.038). In addition, asymmetric density has low sensitivity of 50% on both iCAD and R2 digital, which are significantly lower than the 83% for calcification (p=0.017).

### DISCUSSION

In this study, we evaluated the standalone sensitivity of commercial CADe systems to breast cancers that were found at screening or follow-up in the DMIST study for both digital and screen-film mammography. We found that the CADe systems tested were able to detect significantly more cancers than were found on initial screening by the original radiologists in DMIST for both digital and screen-film mammography; for each manufacturer, standalone sensitivity was higher for digital than for screen-film; and the standalone sensitivities of CADe systems was not influenced by lesion characteristics (histology, lesion size, lesion type) or subject characteristics (age, breast density, menopausal status).

The overall average standalone CADe sensitivities reported here, ranging from 0.60 to 0.74, are low compared to other retrospective standalone CADe studies. Bolivar et al reported a digital R2 CADe sensitivity of 93% for a study that included cancers seen on both MLO and CC views.<sup>20</sup> The et al reported iCAD CADe sensitivity of 94% for cancers seen in at least one screening view.<sup>21</sup> Yang et al reported a sensitivity of 96.1% for cases that included only a single cancer per subject detected by radiologists on digital mammograms.<sup>22</sup> Our study differed in that we included mammographically occult lesions. Besides mammographic visibility of lesions, the characteristics of lesions and subjects have been shown to impact radiologist sensitivity. In our study, there was no statistical difference in CADe sensitivity for their study based on histopathology or tumor size.<sup>21</sup> Bolivar et al reported no difference in digital CADe sensitivity based on breast density, lesion type, or histopathology, but they did find a significant difference based on tumor size for masses.<sup>20</sup>

One unique benefit of using the DMIST dataset for CADe evaluation was the ability to assess the performance of CADe algorithms for both digital and screen-film mammography. In our study, we used a single CADe system for both digital and screen-film with iCAD, finding no difference in standalone CADe sensitivities between iCAD digital and screen-film (p=0.26). For R2 digital and screen-film cases, two different versions of R2 CADe algorithms were used. There was a significantly higher standalone CADe sensitivity for R2 digital than for R2 screen-film (p=0.003).

There were a few of limitations to this study. First, the methodology for scoring CADe mark localization of known cancer locations could have introduced bias. Three different radiologists conducted the scoring for screen-film (both iCAD and R2), iCAD digital, and R2 digital. While we attempted to minimize this potential bias by providing the readers with annotated overlays that directly matched the films generated by a fourth radiologist, the radiologists that scored the digital CAD did have to use their judgment as to whether the CADe mark shown on the computer screen overlapped the annotation on the overlay. Second, because our study included only cancer cases, we could not assess specificity. While higher sensitivity was realized under certain conditions, there is an unknown impact on specificity. We are in the process of completing analysis of a retrospective radiologist reader study that will allow us to assess the impact of CADe on radiologist sensitivity and specificity with digital mammography. Third, the results provide the performance of CADe at a given point in time. The algorithms tested were current at the time this study was

performed in 2008-2009 timeframe. Certainly newer versions of these algorithms would be available today resulting in perhaps different results. For reference we do provide results from R2 CADe at different sensitivity settings in the appendix. Standalone CADe sensitivity assessment is a critical first step in determining if a CADe algorithm is ready for clinical assessment. Our results show that CADe marked between 26.1%-58.7% more cancers than were detected by human readers alone, suggesting that the use of CADe with both screen-film and digital mammography could lead to improved breast cancer detection.

# Appendix

#### Appendix:

# Results

The following tables represent the full R2 screen-film dataset including all three sensitivity settings tested. R2-1 is the default sensitivity setting, with R2-2 having higher sensitivity than the default setting, and R2-3 having the highest sensitivity setting allowed by the system.

#### Table A-1

Patient and Lesion Characteristics for the R2 CADe Screen-Film Full Dataset

		SCREEN-FILM	
	R2-1(N=303) n (%)	R2-2(N=300) n (%)	R2-3(N=291) n (%)
AGE (years)			
Age<50	64 (21.1)	65 (21.7)	61 (21.0)
50<=Age<=64	150 (49.5)	147 (49.0)	145 (49.8)
Age>=65	89 (29.4)	88 (29.3)	85 (29.1)
BREAST DENSITY			
Dense Breasts	150 (49.5)	148 (49.3)	146 (50.2)
Fatty Breasts	153 (50.5)	152 (50.7)	145 (49.8)
MENOPAUSAL STATUS			
Pre or Peri-Menopausal	90 (30.1)	91 (30.7)	87 (30.3)
Post-Menopausal	209 (69.9)	205 (69.3)	200 (69.7)
HISTOLOGY			
DCIS	88 (29.1)	86 (28.8)	83 (28.6)
Invasive	214 (70.9)	213 (71.2)	207 (71.4)
TUMOR SIZE (mm)			
Size<=5	43 (19.7)	42 (19.4)	39 (18.6)
5 <size<=10< th=""><th>56 (25.7)</th><th>56 (25.9)</th><th>55 (26.2)</th></size<=10<>	56 (25.7)	56 (25.9)	55 (26.2)
Size>10	119 (54.6)	118 (54.6)	116 (55.2)
LESION TYPE			
Mass	129 (48.9)	127 (48.7)	123 (48.2)
Asym Density	21 (8.0)	21 (8.0)	21 (8.2)
Calcification	99 (37.5)	98 (37.5)	96 (37.6)

		SCREEN-FILM	
	R2-1(N=303) n (%)	R2-2(N=300) n (%)	R2-3(N=291) n (%)
Arch Distortion	15 (5.7)	15 (5.7)	15 (5.9)

#### Table A-2

Sensitivity Estimates and False Positive Marks per Case for the R2 CADe Screen-Film Full Dataset

			Screen-Film	
		R2-1	R2-2	R2-3
	n/N	176/303	190/300	184/291
Sensitivity	Estimate	0.58	0.63	0.63
	95% CI	(0.52,0.64)	(0.58,0.69)	(0.57,0.69)
	Ν	303	300	291
	Mean	1.5	2.0	2.0
False Positives	SD	1.42	1.65	1.64
	Range	(0, 7)	(0,8)	(0,8)
	Median	1	2	2

 $n = number \ of \ correctly \ detected \ cancers$ 

N = total number of cancer cases analyzable for each CADe system

95% CI = 95% exact confidence interval

SD =standard deviation

**NIH-PA Author Manuscript** 

**NIH-PA** Author Manuscript

ကို
∢
Ð
q
a'

	-
	<u> </u>
	5
	Υ.
	2
	$\overline{\sigma}$
	Ξ.
	$\overline{\sigma}$
	_
- 1	
	7
	<u> </u>
. [	L
	H
- 5	
۰.	<u> </u>
	ЪĽ.
	÷
	Ξ.
	e
	(1)
	<u> </u>
	5
	<u> </u>
<u>ر</u>	$\boldsymbol{n}$
	· •
	(1)
	ñ
	1
- 7	)
	-
	<b>~</b> 1
	. N
	~
	-
	e
- 5	
	Ξ.
	ч.
	$\circ$
	<u> </u>
	_
	cn.
	S
	S
	ICS
•	stics
•	Istics
•	ristics
•	eristics
•	eristics
•	steristics
•	cteristics
•	acteristics
•	racteristics
•	aracteristics
•	<b>aracteristics</b>
•	haracteristics
-	Characteristics
	Characteristics
	<b>I</b> Characteristics
	n Characteristics
	on Characteristics
	10n Characteristics
	sion Characteristics
	esion Characteristics
	esion Characteristics
	Lesion Characteristics
	Lesion Characteristics
	d Lesion Characteristics
	id Lesion Characteristics
	nd Lesion Characteristics
	and Lesion Characteristics
	and Lesion Characteristics
	t and Lesion Characteristics
	nt and Lesion Characteristics
	ent and Lesion Characteristics
	lent and Lesion Characteristics
	tient and Lesion Characteristics
5	atient and Lesion Characteristics
	atient and Lesion Characteristics
5	Patient and Lesion Characteristics
5	Patient and Lesion Characteristics
	y Patient and Lesion Characteristics
5	by Patient and Lesion Characteristics
	by Patient and Lesion Characteristics
	s by Patient and Lesion Characteristics
	es by Patient and Lesion Characteristics
	les by Patient and Lesion Characteristics
3	ties by Patient and Lesion Characteristics
	ities by Patient and Lesion Characteristics
5	vities by Patient and Lesion Characteristics
	<b>Exities by Patient and Lesion Characteristics</b>
- - - -	tivities by Patient and Lesion Characteristics
5 	itivities by Patient and Lesion Characteristics
5 - - - -	sitivities by Patient and Lesion Characteristics
	nsitivities by Patient and Lesion Characteristics
	insitivities by Patient and Lesion Characteristics
	ensitivities by Patient and Lesion Characteristics
5	Sensitivities by Patient and Lesion Characteristics

						SCREEN-F	ILM			
			R2-1			R2-2			R2-3	
		Se	nsitivity	95% CI	Sei	ısitivity	IJ %56	Se	nsitivity	13 %Se
	Age<50	0.53	(34/ 64)	(0.41, 0.65)	0.58	(38/ 65)	(0.46, 0.70)	0.59	(36/61)	(0.46, 0.71)
AGE(years)	50 Age 64	0.59	(88/ 150)	(0.51, 0.66)	0.63	(93/ 147)	(0.55, 0.71)	0.64	(93/ 145)	(0.56, 0.72)
	Age 65	0.61	(54/ 89)	(0.50, 0.70)	0.67	(29/ 88)	(0.57, 0.76)	0.65	(55/85)	(0.54, 0.74)
BREAST	Dense	0.60	(90/ 150)	(0.52, 0.68)	0.64	(95/ 148)	(0.56, 0.71)	0.65	(95/ 146)	(0.57, 0.72)
DENSITY	Fatty	0.56	(86/ 153)	(0.48, 0.64)	0.63	(95/ 152)	(0.55, 0.70)	0.61	(89/ 145)	(0.53, 0.69)
MENOPAUSAL	Pre/Peri	0.51	(46/ 90)	(0.41, 0.61)	0.58	(53/91)	(0.48, 0.68)	0.57	(50/87)	(0.47, 0.67)
STATUS	Post	0.61	(127/ 209)	(0.54, 0.67)	0.65	(134/205)	(0.59, 0.72)	0.66	(131/200)	(0.59, 0.72)
	DCIS	0.57	(50/ 88)	(0.46, 0.67)	0.63	(54/86)	(0.52, 0.72)	0.61	(51/83)	(0.51, 0.71)
IBOTOISH	Invasive	0.59	(126/214)	(0.52, 0.65)	0.64	(136/213)	(0.57, 0.70)	0.64	(133/207)	(0.57, 0.70)
	Size 5	0.53	(23/ 43)	(0.39, 0.68)	0.64	(27/42)	(0.49, 0.77)	0.67	(26/39)	(0.51, 0.80)
TUMOR SIZE(mm)	5 <size 10<="" th=""><th>0.48</th><th>(27/ 56)</th><th>(0.36, 0.61)</th><th>0.55</th><th>(31/56)</th><th>(0.42, 0.68)</th><th>0.56</th><th>(31/55)</th><th>(0.43, 0.69)</th></size>	0.48	(27/ 56)	(0.36, 0.61)	0.55	(31/56)	(0.42, 0.68)	0.56	(31/55)	(0.43, 0.69)
	Size>10	0.66	(78/ 119)	(0.57, 0.74)	0.69	(81/118)	(0.60, 0.76)	0.68	(79/ 116)	(0.59, 0.76)
	Mass	0.60	(78/ 129)	(0.52, 0.69)	0.66	(84/127)	(0.57, 0.74)	0.67	(83/123)	(0.59, 0.75)
LESION	Asym Density	0.52	(11/21)	(0.32, 0.72)	0.57	(12/21)	(0.36, 0.76)	0.57	(12/21)	(0.36, 0.76)
TYPE	Calcification	0.73	(72/ 99)	(0.63, 0.81)	0.77	(75/98)	(0.67, 0.84)	0.76	(73/96)	(0.67, 0.84)
	Arch Distortion	0.53	(8/ 15)	(0.29, 0.76)	0.67	(10/ 15)	(0.41, 0.85)	0.67	(10/ 15)	(0.41, 0.85)

Cole et al.

# Table A-4

Univariate Logistic Regression Models by Patient and Lesion Characteristics for the R2 CADe Screen-film Full Datasets

						SCREEN-FILI	M			
			R2-1			R2-2			R2-3	
	-	OR	95% CI(OR)	p-value	OR	95% CI(OR)	p-value	OR	95% CI(OR)	p-value
AGE										
02 <sup></sup> •	50 Age 64	1.25	(0.70, 2.26)	0.4538	1.22	(0.67, 2.22)	0.5071	1.24	(0.67, 2.29)	0.4882
Age<50	Age 65	1.36	(0.71, 2.61)	0.3520	1.45	(0.74, 2.81)	0.2768	1.27	(0.65, 2.51)	0.4844
50 Age 64	Age 65	1.09	(0.64, 1.86)	0.7600	1.18	(0.68, 2.06)	0.5575	1.03	(0.59, 1.79)	0.9308
BREAST DENSI	IY									
Dense	Fatty	0.86	(0.54, 1.35)	0.5038	0.93	(0.58, 1.49)	0.7615	0.85	(0.53, 1.37)	0.5142
MENOPAUSAL 5	STATUS									
Pre/Peri	Post	1.48	(0.90, 2.44)	0.1218	1.35	(0.82, 2.25)	0.2417	1.40	(0.84, 2.35)	0.1961
HISTOLOGY										
DCIS	Invasive	1.09	(0.66, 1.80)	0.7408	1.05	(0.62, 1.76)	0.8629	1.13	(0.67, 1.91)	0.6540
TUMOR SIZE (m	(un									
L L	5 <size 10<="" th=""><th>0.81</th><th>(0.37, 1.79)</th><th>0.6031</th><th>0.69</th><th>(0.30, 1.57)</th><th>0.3743</th><th>0.65</th><th>(0.28, 1.52)</th><th>0.3149</th></size>	0.81	(0.37, 1.79)	0.6031	0.69	(0.30, 1.57)	0.3743	0.65	(0.28, 1.52)	0.3149
c azic	Size>10	1.65	(0.81, 3.36)	0.1638	1.22	(0.58, 2.55)	0.6048	1.07	(0.49, 2.31)	0.8681
5 <size 10<="" th=""><th>Size&gt;10</th><th>2.04</th><th>(1.07, 3.90)</th><th>0.0302</th><th>1.77</th><th>(0.92, 3.40)</th><th>0.0889</th><th>1.65</th><th>(0.85, 3.20)</th><th>0.1359</th></size>	Size>10	2.04	(1.07, 3.90)	0.0302	1.77	(0.92, 3.40)	0.0889	1.65	(0.85, 3.20)	0.1359
LESION TYPE										
	Asym Density	0.72	(0.28, 1.82)	0.4855	0.68	(0.27, 1.75)	0.4252	0.64	(0.25, 1.65)	0.3578
Mass	Calcification	1.74	(0.99, 3.07)	0.0542	1.67	(0.92, 3.02)	0.0912	1.53	(0.84, 2.79)	0.1663
	Arch Distortion	0.75	(0.26, 2.19)	0.5949	1.02	(0.33, 3.18)	0.9676	0.96	(0.31, 3.01)	0.9494
Acres Doneiter	Calcification	2.42	(0.92, 6.35)	0.0718	2.45	(0.92, 6.53)	0.0744	2.38	(0.89, 6.36)	0.0838
And Density	Arch Distortion	1.04	(0.28, 3.92)	0.9550	1.50	(0.38, 5.95)	0.5641	1.50	(0.38, 5.95)	0.5640
Calcification	Arch Distortion	0.43	(0.14, 1.30)	0.1335	0.61	(0.19, 1.98)	0.4133	0.63	(0.20, 2.03)	0.4398

AJR Am J Roentgenol. Author manuscript; available in PMC 2013 September 01.

OR = odds ratio

# REFERENCES

- Pisano ED, Gatsonis C, Hendrick E, et al. Diagnostic performance of digital versus film mammography for breast-cancer screening. N Engl J Med. 2005; 353:1773–1783. [PubMed: 16169887]
- Lewin JM, D'Orsi CJ, Hendrick RE, et al. Clinical comparison of full-field digital mammography and screen-film mammography for detection of breast cancer. AJR Am J Roentgenol. 2002; 179:671–677. [PubMed: 12185042]
- Lewin JM, Hendrick RE, D'Orsi CJ, et al. Comparison of full-field digital mammography with screen-film mammography for cancer detection: results of 4,945 paired examinations. Radiology. 2001; 218:873–880. [PubMed: 11230669]
- 4. Skaane P, Balleyguier C, Diekmann F, et al. Breast lesion detection and classification: comparison of screen-film mammography and full-field digital mammography with soft-copy reading--observer performance study. Radiology. 2005; 237:37–44. [PubMed: 16100086]
- Skaane P, Kshirsagar A, Stapleton S, Young K, Castellino RA. Effect of computer-aided detection on independent double reading of paired screen-film and full-field digital screening mammograms. AJR Am J Roentgenol. 2007; 188:377–384. [PubMed: 17242245]
- 6. Skaane P, Skjennald A. Screen-film mammography versus full-field digital mammography with soft-copy reading: randomized trial in a population-based screening program--the Oslo II Study. Radiology. 2004; 232:197–204. [PubMed: 15155893]
- Skaane P, Young K, Skjennald A. Population-based mammography screening: comparison of screen-film and full-field digital mammography with soft-copy reading--Oslo I study. Radiology. 2003; 229:877–884. [PubMed: 14576447]
- Pisano ED, Gatsonis CA, Yaffe MJ, et al. American College of Radiology Imaging Network digital mammographic imaging screening trial: objectives and methodology. Radiology. 2005; 236:404– 412. [PubMed: 15961755]
- Pisano ED, Hendrick RE, Yaffe MJ, et al. Diagnostic accuracy of digital versus film mammography: exploratory analysis of selected population subgroups in DMIST. Radiology. 2008; 246:376–383. [PubMed: 18227537]
- Heddson B, Ronnow K, Olsson M, Miller D. Digital versus screen-film mammography: a retrospective comparison in a population-based screening program. Eur J Radiol. 2007; 64:419– 425. [PubMed: 17383841]
- Del Turco MR, Mantellini P, Ciatto S, et al. Full-field digital versus screen-film mammography: comparative accuracy in concurrent screening cohorts. AJR Am J Roentgenol. 2007; 189:860–866. [PubMed: 17885057]
- Vigeland E, Klaasen H, Klingen TA, Hofvind S, Skaane P. Full-field digital mammography compared to screen film mammography in the prevalent round of a population-based screening programme: the Vestfold County Study. Eur Radiol. 2008; 18:183–191. [PubMed: 17680246]
- 13. Skaane P. Studies comparing screen-film mammography and full-field digital mammography in breast cancer screening: updated review. Acta Radiol. 2009; 50:3–14. [PubMed: 19037825]
- Breast Cancer Surveillance Consortium. [Accessed 08/17, 2010] Performance Measures for 3,884,059 Screening Mammography Examinations from 1996 to 2997 by Age. Available at:http:// breastscreening.cancer.gov/data/performance/screening/perf\_age.html
- Balleyguier C, Kinkel K, Fermanian J, et al. Computer-aided detection (CAD) in mammography: does it help the junior or the senior radiologist? Eur J Radiol. 2005; 54:90–96. [PubMed: 15797297]
- Hukkinen K, Vehmas T, Pamilo M, Kivisaari L. Effect of computer-aided detection on mammographic performance: experimental study on readers with different levels of experience. Acta Radiol. 2006; 47:257–263. [PubMed: 16613306]
- Paquerault S, Samuelson FW, Petrick N, Myers KJ, Smith RC. Investigation of reading mode and relative sensitivity as factors that influence reader performance when using computer-aided detection software. Acad Radiol. 2009; 16:1095–1107. [PubMed: 19523855]

- Karssemeijer N, Bluekens AM, Beijerinck D, et al. Breast cancer screening results 5 years after introduction of digital mammography in a population-based screening program. Radiology. 2009; 253:353–358. [PubMed: 19703851]
- 19. Nishikawa RM, Giger ML, Doi K, et al. Effect of case selection on the performance of computeraided detection schemes. Med Phys. 1994; 21:265–269. [PubMed: 8177159]
- 20. Bolivar AV, Gomez SS, Merino P, et al. Computer-aided detection system applied to full-field digital mammograms. Acta Radiol. 2010; 51:1086–1092. [PubMed: 20883182]
- The JS, Schilling KJ, Hoffmeister JW, Friedmann E, McGinnis R, Holcomb RG. Detection of breast cancer with full-field digital mammography and computer-aided detection. AJR Am J Roentgenol. 2009; 192:337–340. [PubMed: 19155392]
- Yang SK, Moon WK, Cho N, et al. Screening mammography-detected cancers: sensitivity of a computer-aided detection system applied to full-field digital mammograms. Radiology. 2007; 244:104–111. [PubMed: 17507722]

Cole et al.



#### Figure 1.

Representative image of CAD marks and digital mammographic image R2.  $\Delta$  symbol indicates calcifications, while masses are circumscribed by a freeform shape.

**NIH-PA** Author Manuscript

Cole et al.



#### Figure 2.

Screen shot of CC views of digital mammogram with iCAD report overlay. Square represents calcifications and circle represents mass finding.



#### Figure 3.

Example of screen-film mammogram with annotated acetate overlay denoting location of cancer.

Patient and Lesion Characteristics for the Common Cases (N=161)

	n (%)
AGE (years)	
Age<50	30 (18.6)
50<=Age<=64	81 (50.3)
Age>=65	50 (31.1)
BREAST DENSITY	
Dense Breasts	75 (46.6)
Fatty Breasts	86 (53.4)
MENOPAUSAL STATUS	
Pre or Peri-Menopausal	43 (26.7)
Post-Menopausal	114 (70.8)
Missing	4 (2.5)
HISTOLOGY	
DCIS	46 (28.6)
Invasive	115 (71.4)
TUMOR SIZE (mm)	
Size<=5	21 (13.0)
5 <size<=10< th=""><th>31 (19.3)</th></size<=10<>	31 (19.3)
Size>10	69 (42.9)
Missing	40 (24.8)
LESION TYPE	
Mass	65 (40.4)
Asym Density	12 (7.5)
Calcification	60 (37.3)
Arch Distortion	8 (5.0)
Missing	16 (9.9)

Sensitivity and False Positive Marks for the Common Cases

		iCAD	R2
Sensitivity	DG	0.74(119/161) 95% CI: (0.66,0.81)	0.74(119/161) 95% CI: (0.66,0.81)
	SF	0.69(111/161) 95% CI: (0.61,0.76)	0.60(97/161) 95% CI: (0.52,0.68)
	DGvSF	0.05 95% CI: (-0.03,0.13) p=0.2559	0.14 95% CI: (0.05,0.23) p=0.0026
False Positive	DG	$2.57 \pm 1.92$	$2.07 \pm 1.57$
	SF	$3.06 \pm 1.72$	$1.52\pm1.45$
	DGvSF	-0.49 95% CI: (-0.77,-0.21)p=0.0008	0.55 95% CI: ( 0.33, 0.78) p<0.0001

Comparison of Sensitivities of CADe system and DMIST Study for the Common Cases

	iC	AD	R	2
	DG	SF	DG	SF
CAD	0.74(119/161)	0.69(111/161)	0.74(119/161)	0.60(97/161)
	95% CI:	95% CI:	95% CI:	95% CI:
	(0.66,0.81)	(0.61,0.76)	(0.66,0.81)	(0.52,0.68)
DMIST	0.43(69/161)	0.41(66/161)	0.43(69/161)	0.41(66/161)
	95% CI:	95% CI:	95% CI:	95% CI:
	(0.35,0.51)	(0.33,0.49)	(0.35,0.51)	(0.33,0.49)
CAD v DMIST	0.31	0.28	0.31	0.19
	95% CI:	95% CI:	95% CI:	95% CI:
	(0.22,0.40)	(0.19,0.37)	(0.22,0.40)	(0.10,0.28)
	p<.0001	p<.0001	p<.0001	p<.0001

**NIH-PA** Author Manuscript

Cases
Common
the
for
Characteristics
Lesion
and
Patient
by
Sensitivities

			iC	AD					R	5		
		DG			SF			DG			SF	
	ser	nsitivity	95% CI	sen	sitivity	95% CI	sen	sitivity	95% CI	sen	sitivity	95% CI
AGE(years)												
Age<50	0.73	(22/ 30)	(0.55, 0.86)	0.63	(19/30)	(0.45, 0.78)	0.80	(24/30)	(0.62, 0.91)	0.53	(16/30)	(0.36, 0.70)
50 Age 64	0.72	(58/ 81)	(0.61, 0.80)	0.73	(59/81)	(0.62, 0.81)	0.74	(60/ 81)	(0.63, 0.82)	0.62	(50/81)	(0.51, 0.72)
Age 65	0.78	(39/ 50)	(0.64, 0.87)	0.66	(33/50)	(0.52, 0.78)	0.70	(35/50)	(0.56, 0.81)	0.62	(31/50)	(0.48, 0.74)
<b>BREAST DENSI</b>	Υ.											
Dense	0.76	(57/75)	(0.65, 0.84)	0.67	(50/75)	(0.55, 0.76)	0.77	(58/75)	(0.67, 0.85)	0.60	(45/75)	(0.49, 0.70)
Fatty	0.72	(62/ 86)	(0.62, 0.81)	0.71	(61/ 86)	(0.61, 0.80)	0.71	(61/86)	(0.61, 0.80)	0.60	(52/ 86)	(0.50, 0.70)
MENOPAUSAL S	TATUS											
Pre or Peri- Menopausal	0.72	(31/ 43)	(0.57, 0.83)	0.63	(27/43)	(0.48, 0.76)	0.74	(32/43)	(0.59, 0.85)	0.49	(21/43)	(0.34, 0.63)
Post- Menopausal	0.75	(85/ 114)	(0.66, 0.82)	0.71	(81/114)	(0.62, 0.79)	0.74	(84/114)	(0.65, 0.81)	0.64	(73/114)	(0.55, 0.72)
HISTOLOGY												
DCIS	0.74	(34/ 46)	(0.59, 0.85)	0.59	(27/46)	(0.44, 0.72)	0.74	(34/46)	(0.59, 0.85)	0.63	(29/ 46)	(0.48, 0.76)
Invasive	0.74	(85/ 115)	(0.65, 0.81)	0.73	(84/115)	(0.64, 0.80)	0.74	(85/115)	(0.65, 0.81)	0.59	(68/ 115)	(0.50, 0.68)
<b>TUMOR SIZE (m</b>	<b>m</b> )											
Size 5	0.71	(15/21)	(0.49, 0.87)	0.76	(16/21)	(0.54, 0.90)	0.71	(15/21)	(0.49, 0.87)	0.48	(10/21)	(0.28, 0.68)
5 <size 10<="" th=""><th>0.74</th><th>(23/ 31)</th><th>(0.56, 0.87)</th><th>0.68</th><th>(21/31)</th><th>(0.50, 0.82)</th><th>0.81</th><th>(25/31)</th><th>(0.63, 0.91)</th><th>0.55</th><th>(17/31)</th><th>(0.37, 0.71)</th></size>	0.74	(23/ 31)	(0.56, 0.87)	0.68	(21/31)	(0.50, 0.82)	0.81	(25/31)	(0.63, 0.91)	0.55	(17/31)	(0.37, 0.71)
Size>10	0.75	(52/ 69)	(0.64, 0.84)	0.71	(49/ 69)	(0.59, 0.80)	0.75	(52/ 69)	(0.64, 0.84)	0.67	(46/ 69)	(0.55, 0.77)
LESION TYPE												
Mass	0.69	(45/ 65)	(0.57, 0.79)	0.69	(45/ 65)	(0.57, 0.79)	0.72	(47/ 65)	(0.60, 0.82)	0.55	(36/ 65)	(0.43, 0.67)
Asym Density	0.50	(6/ 12)	(0.24, 0.76)	0.75	(9/12)	(0.45, 0.92)	0.50	(6/12)	(0.24, 0.76)	0.75	(9/ 12)	(0.45, 0.92)
Calcification	0.83	(50/ 60)	(0.72, 0.91)	0.73	(44/ 60)	(0.61, 0.83)	0.83	(50/ 60)	(0.72, 0.91)	0.73	(44/ 60)	(0.61, 0.83)
Arch Distortion	0.88	(2/ 8)	(0.46, 0.98)	0.63	(5/ 8)	(0.28, 0.87)	0.75	(6/8)	(0.38, 0.94)	0.63	(5/ 8)	(0.28, 0.87)

Univariate Logistic Regression Models by Patient and Lesion Characteristics for the Common Cases

ICAD          ICAD         ICAD         ICAD         ICAD         ICAD         ICAD         ICAD         ICAD         ICAD   <														
$ \                                   $					iC	4D					R	5		
				DG			SF			DG			SF	
AE(cons)         0.35 (0.64, 3.78)         0.353 (0.51 (0.51, 19)         0.5196           Age 65         129         0.45, 3.68)         0.633 (0.41, 2.8)         0.64, 3.78)         0.3323         0.71         (0.26, 1.9)         0.5196           Age 65         129         0.45, 3.68)         0.633 (1.12)         0.44, 2.89)         0.8087         0.83 (0.20, 1.72)         0.3281           S0 Age 64         Age 65         141         (0.62, 3.21)         0.4184         0.72         (0.34, 1.55)         0.4064         0.82         0.3204           BREAST DENCIT         2.2.2         0.40, 1.66         0.5736         1.42         (0.66, 3.05)         0.3204         0.92         (0.37, 1.79)         0.510           BREAST DENCIT         2.2.2         0.40, 1.66         0.5736         1.42         (0.66, 3.05)         0.3204         0.96         0.43, 2.15)         0.3579           MENOPAUSAL         2.1.13         (0.52, 2.50)         0.7336         1.43         0.72         0.3504         0.43, 2.15)         0.9259           MENOPAUSAL         Post         1.13         (0.52, 2.50)         0.736         1.45         (0.66, 3.05)         0.3504         0.46, 2.18)         0.3504           M			OR	95% CI(OR)	p-value	OR	95% CI(OR)	p-value	OR	95% CI(OR)	p-value	OR	95% CI(OR)	p-value
	AGE(years)													
TREAD TotalRef1.20.431.12(0.44, 28)0.80870.808(0.20, 172)0.328150 Age 651.41(0.62, 3.21)0.41840.72(0.34, 1.55)0.40640.82(0.37, 1.79)0.6121BREAST DENSITYAge 651.41(0.62, 3.21)0.41840.72(0.34, 1.55)0.40640.82(0.37, 1.79)0.6121BREAST DENSITYAge 651.14(0.62, 3.28)0.55000.72(0.35, 1.46)0.35700.3570BREAST DENSITYAntity0.82(0.40, 1.66)0.77361.22(0.63, 2.38)0.56000.72(0.35, 1.46)0.3570BREAST DENSITY0.82(0.40, 1.66)0.77361.22(0.69, 3.05)0.56000.72(0.35, 1.46)0.3570MENOPAUSAL STATUSAntisty0.82(0.40, 1.56)0.77361.45(0.43, 2.15)0.93510.9351MENOPAUSAL STATUSAntisty0.920.73661.91(0.93, 3.91)0.07761.000.920.9351MENOPAUSAL STATUSAntisty0.920.9330.93010.07761.000.920.93310.97610.9251MENOPAUSAL STATUSAntisty0.920.9330.93010.97170.97310.97310.97310.97310.97310.9731MENOPAUSAL STATUSAntisty0.920.93330.8250.660.920.92310.97410.97310.944530.9351MENOPAUSALFree1.120.9333<		50 Age 64	0.92	(0.36, 2.35)	0.8570	1.55	(0.64, 3.78)	0.3323	0.71	(0.26, 1.99)	0.5196	1.41	(0.61, 3.29)	0.4246
50 age 64         Age 65         141 $(0.62, 3.21)$ $0.4184$ $0.72$ $(0.34, 1.55)$ $0.4064$ $0.82$ $(0.37, 1.79)$ $0.6121$ BIEAST DENSITY           Dense         Faity $0.82$ $(0.40, 1.66)$ $0.5736$ $1.22$ $(0.63, 2.38)$ $0.5600$ $0.72$ $(0.37, 1.79)$ $0.6121$ MENDATISAL         Automation of the state	Age<50	Age 65	1.29	(0.45, 3.68)	0.6354	1.12	(0.44, 2.89)	0.8087	0.58	(0.20, 1.72)	0.3281	1.43	(0.57, 3.57)	0.4466
BIREAST DENSITY           BIREAST DENSITY           Dense         Fatty         0.82         (0.40, 1.66)         0.5736         1.22         (0.63, 2.38)         0.753, 1.46)         0.35716           MENOPAUSAL STATUS           PrePeri         Post         1.113         (0.52, 2.50)         0.7536         1.45         (0.63, 2.38)         0.956.00         0.72         (0.43, 2.15)         0.9259           PrePeri         1.00         (0.73, 2.30)         0.7336         0.956.00         0.72, 1.46         0.99259           PrePeri         1.00         (0.73, 2.13)         (0.95, 2.13)         0.95600         0.735, 1.46         0.95269           PrePeri         1.13         (0.52, 2.30)         0.7336         0.7176           Size 10         1.12         (0.41, 3.65)         0.7176           Size 10          (0.41, 3.65)	50 Age 64	Age 65	1.41	(0.62, 3.21)	0.4184	0.72	(0.34, 1.55)	0.4064	0.82	(0.37, 1.79)	0.6121	1.01	(0.49, 2.09)	0.9752
Dense         Fatty         0.82         (0.40, 1.66)         0.5736         1.22         (0.63, 2.38)         0.5600         0.72         (0.35, 1.46)         0.3570           MENOPAUSAL         Fatt         1.13         (0.57, 2.50)         0.5736         1.45         (0.69, 3.05)         0.3204         0.96         (0.35, 1.45)         0.9259           Pre/Peri         Post         1.13         (0.52, 2.50)         0.7536         1.45         (0.69, 3.05)         0.3204         0.96         (0.43, 2.15)         0.9259           HISTOLOGY         Invasive         1.00         (0.46, 2.18)         1.0000         1.91         (0.92, 3.01)         0.776         1.00         (0.45, 2.15)         0.9259           TOMOR SIZE         Invasive         1.03         0.823, 3.83         0.8255         0.66         0.19, 2.30         0.7176         1.00         0.44, 3.55         0.9259           Size-10         1.12         0.33, 3.83         0.8255         0.66         0.19, 2.30         0.7176         0.44, 3.55         0.441, 3.65         0.7176           Size-10         1.22         0.41, 3.65         0.7176         0.733         0.733         0.73         0.26, 2.09         0.5656           Size-10         1.22 </th <th>BREAST DEN</th> <th>SITY</th> <th></th>	BREAST DEN	SITY												
$\begin{tabular}{ c c c c c c c c c c c c c c c c c c c$	Dense	Fatty	0.82	(0.40, 1.66)	0.5736	1.22	(0.63, 2.38)	0.5600	0.72	(0.35, 1.46)	0.3570	1.02	(0.54, 1.92)	0.9520
Pre/Peri         Post         1:13 $(0.52, 2.50)$ $0.7536$ $1.45$ $(0.69, 3.05)$ $0.3204$ $0.96$ $(0.43, 2.15)$ $0.9259$ HISTOLOCY         Invasive         1.01 $(0.73, 2.15)$ $0.735$ $1.47$ $(0.93, 3.91)$ $0.0776$ $1.00$ $0.46, 2.18$ $1.0000$ DCIS         Invasive         1.00 $1.91$ $(0.93, 3.91)$ $0.0776$ $1.00$ $0.46, 2.18$ $1.0000$ TUMOR SIZE         Invasive         1.05 $0.33, 3.98$ $0.8255$ $0.66$ $(0.19, 2.30)$ $0.776$ $1.00$ $0.4413$ Size         10         1.15 $(0.33, 3.98)$ $0.8255$ $0.66$ $(0.19, 2.30)$ $0.717$ $0.73$ $0.73$ $0.4413$ Size         10         1.12 $(0.33, 3.98)$ $0.8255$ $0.66$ $0.1943$ $0.73$ $0.73$ $0.73$ $0.73$ $0.4413$ Size         10         1.16 $0.7375$ $0.732$ $0.73$ $0.743$ $0.73$ $0.735$ $0.7165$ $0.7176$	MENOPAUSA	T STATUS												
HISTOLOGY           DCIS         Invasive         1.00         (0.46, 2.18)         1.0000         1.91         (0.93, 3.91)         0.0776         1.00         (0.46, 2.18)         1.0000           TUMOR SIZE (mm           Size 10         1.15         (0.33, 3.98)         0.8255         0.66         (0.19, 2.30)         0.5108         1.67         (0.46, 2.18)         1.0000           Size 10         1.15         (0.33, 3.98)         0.8255         0.66         (0.19, 2.30)         0.5108         1.67         (0.44, 3.65)         0.7176           Size 10         1.15         (0.33, 3.98)         0.8255         0.66         (0.19, 2.30)         0.5108         1.67         (0.41, 3.65)         0.7176           Size 10         1.12         (0.44, 2.83)         0.7177         (0.77         (0.25, 2.37)         0.6435         1.22         (0.41, 3.65)         0.7176           Size>10         1.02         (0.44, 2.83)         0.7177         (0.77, 2.91)         0.7413         0.73         0.264, 20)         0.562, 20)         0.562, 20)         0.562, 20)         0.562, 20)         0.562, 20)         0.562, 20)         0.562, 20)         0.562, 20)         0.562, 20)	Pre/Peri	Post	1.13	(0.52, 2.50)	0.7536	1.45	(0.69, 3.05)	0.3204	0.96	(0.43, 2.15)	0.9259	1.87	(0.92, 3.79)	0.0852
	HISTOLOGY													
	DCIS	Invasive	1.00	(0.46, 2.18)	1.0000	1.91	(0.93, 3.91)	0.0776	1.00	(0.46, 2.18)	1.0000	0.85	(0.42, 1.72)	0.6469
Figs         5-Size 10         1.15         (0.33, 3.98)         0.8255         0.66         0.19, 2.30)         0.5108         1.67         (0.45, 6.12)         0.413           Size>10         1.22         (0.41, 3.65)         0.7177         0.77         (0.25, 2.37)         0.6435         1.22         (0.41, 3.65)         0.7176           5-Size 10         1.22         (0.41, 3.65)         0.7177         (0.77         (0.25, 2.37)         0.6435         1.22         (0.41, 3.65)         0.7176           5-Size 10         Size>10         1.06         (0.40, 2.82)         0.9007         1.17         (0.47, 2.91)         0.7413         0.73         0.7176           LESION TYPE         Image         0.44         0.13, 1.55)         0.2029         1.33         (0.33, 5.45)         0.7413         0.1340         0.1340           Mass         Calcification         2.22         (0.94, 5.25)         0.0686         1.22         (0.56, 2.66)         0.6131         1.91         (0.14, 1.34)         0.1340           Mass         Calcification         2.22         (0.94, 5.25)         0.0686         1.22         (0.56, 2.66)         0.6131         1.91         (0.14, 1.34)         0.1341           Mass         Calcification	<b>TUMOR SIZE</b>	(uuu)												
Matrix 5         Size>10         1.22 $(0.41, 3.65)$ $0.7177$ $(0.77)$ $(0.72, 2.37)$ $(0.6435)$ $(1.3, 65)$ $0.7176$ S         Size>10         1.05 $(0.40, 2.82)$ $0.9007$ $1.17$ $(0.74, 2.91)$ $0.73$ $(0.26, 2.09)$ $0.7176$ LESION TYPE         Asym Density $0.44$ $(0.13, 1.55)$ $0.9007$ $1.17$ $(0.47, 2.91)$ $0.73$ $(0.26, 2.09)$ $0.5626$ Mass         Asym Density $0.44$ $(0.13, 1.55)$ $0.2029$ $1.33$ $(0.33, 5.45)$ $0.6890$ $0.38$ $0.11, 1.34)$ $0.1340$ Mass         Calcification $3.11$ $(0.36, 2.55)$ $0.0697$ $1.91$ $(0.30, 4.57)$ $0.1431$ Mass         Calcification $3.11$ $(0.36, 2.55)$ $0.0697$ $1.15$ $0.21, 6.23$ $0.9767$ Mass         Calcification $3.01$ $0.016, 3.401$ $0.0697$ $1.15$ $0.01, 6.3401$ $0.016, 8.757$ $0.0168^{*}$ Mass         Calcification $3.00$ $0.332$	Cta <sub>n</sub> E	5 <size 10<="" th=""><th>1.15</th><th>(0.33, 3.98)</th><th>0.8255</th><th>0.66</th><th>(0.19, 2.30)</th><th>0.5108</th><th>1.67</th><th>(0.45, 6.12)</th><th>0.4413</th><th>1.34</th><th>(0.44, 4.06)</th><th>0.6095</th></size>	1.15	(0.33, 3.98)	0.8255	0.66	(0.19, 2.30)	0.5108	1.67	(0.45, 6.12)	0.4413	1.34	(0.44, 4.06)	0.6095
\$	c azic	Size>10	1.22	(0.41, 3.65)	0.7177	0.77	(0.25, 2.37)	0.6435	1.22	(0.41, 3.65)	0.7176	2.20	(0.82, 5.93)	0.1193
	5 <size 10<="" th=""><th>Size&gt;10</th><th>1.06</th><th>(0.40, 2.82)</th><th>0.9007</th><th>1.17</th><th>(0.47, 2.91)</th><th>0.7413</th><th>0.73</th><th>(0.26, 2.09)</th><th>0.5626</th><th>1.65</th><th>(0.69, 3.92)</th><th>0.2591</th></size>	Size>10	1.06	(0.40, 2.82)	0.9007	1.17	(0.47, 2.91)	0.7413	0.73	(0.26, 2.09)	0.5626	1.65	(0.69, 3.92)	0.2591
	LESION TYPI	БJ												
Mass         Calcification         2.22 $(0.94, 5.25)$ $0.0686$ $1.22$ $(0.56, 2.66)$ $0.6131$ $1.91$ $(0.80, 4.57)$ $0.1431$ Arch Distortion $3.11$ $(0.36, 26.99)$ $0.3032$ $0.74$ $(0.16, 3.40)$ $0.6997$ $1.15$ $(0.21, 6.23)$ $0.8721$ Arch Distortion $3.11$ $(0.36, 26.99)$ $0.3032$ $0.74$ $(0.16, 3.40)$ $0.6997$ $1.15$ $(0.21, 6.23)$ $0.8721$ Asym Density         Calcification $5.00$ $(1.34, 18.71)$ $0.0168*$ $0.92$ $(0.22, 3.82)$ $0.9049$ $5.00$ $(1.34, 18.71)$ $0.0168*$ Asym Density         Arch Distortion $7.00$ $(0.65, 75.73)$ $0.1093$ $0.56$ $(0.08, 3.86)$ $0.5522$ $3.00$ $(0.42, 21.30)$ $0.2719$ Asym Density         Arch Distortion $1.40$ $(0.15, 12.67)$ $0.7647$ $0.61$ $(0.13, 2.83)$ $0.5722$ $3.00$ $(0.42, 21.30)$ $0.2719$		Asym Density	0.44	(0.13, 1.55)	0.2029	1.33	(0.33, 5.45)	0.6890	0.38	(0.11, 1.34)	0.1340	2.42	(0.60, 9.75)	0.2153
	Mass	Calcification	2.22	(0.94, 5.25)	0.0686	1.22	(0.56, 2.66)	0.6131	1.91	(0.80, 4.57)	0.1431	2.21	(1.04, 4.70)	$0.0384^{*}$
Asym Density         Calcification         5.00         (1.34, 18.71)         0.0168*         0.92         (0.22, 3.82)         0.9049         5.00         (1.34, 18.71)         0.0168*           Asym Density         Arch Distortion         7.00         (0.65, 75.73)         0.1093         0.56         (0.08, 3.86)         0.5522         3.00         (0.42, 21.30)         0.2719           Calcification         Arch Distortion         1.40         (0.15, 12.67)         0.7647         0.61         (0.13, 2.83)         0.5242         0.60         (0.11, 3.41)         0.5647		Arch Distortion	3.11	(0.36, 26.99)	0.3032	0.74	(0.16, 3.40)	0.6997	1.15	(0.21, 6.23)	0.8721	1.34	(0.30, 6.09)	0.7026
Asymptotism         Arch Distortion         7.00         (0.65, 75.73)         0.1093         0.56         (0.08, 3.86)         0.5522         3.00         (0.42, 21.30)         0.2719           Calcification         Arch Distortion         1.40         (0.15, 12.67)         0.7647         0.61         (0.13, 2.83)         0.5522         0.60         (0.11, 3.41)         0.5647	A cum Doncity	Calcification	5.00	(1.34, 18.71)	$0.0168^{*}$	0.92	(0.22, 3.82)	0.9049	5.00	(1.34, 18.71)	$0.0168^{*}$	0.92	(0.22, 3.82)	0.9051
Calcification Arch Distortion 1.40 (0.15, 12.67) 0.7647 0.61 (0.13, 2.83) 0.5242 0.60 (0.11, 3.41) 0.5647		Arch Distortion	7.00	(0.65, 75.73)	0.1093	0.56	(0.08, 3.86)	0.5522	3.00	(0.42, 21.30)	0.2719	0.56	(0.08, 3.86)	0.5525
	Calcification	Arch Distortion	1.40	(0.15, 12.67)	0.7647	0.61	(0.13, 2.83)	0.5242	0.60	(0.11, 3.41)	0.5647	0.61	(0.13, 2.83)	0.5244

AJR Am J Roentgenol. Author manuscript; available in PMC 2013 September 01.

OR: odds ratio