EVALUATION OF *MYCOBACTERIUM TUBERCULOSIS* **DRUG SUSCEPTIBILITY IN CLINICAL SPECIMENS FROM NIGERIA USING GENOTYPE MTBDRPLUS AND MTBDRSL ASSAYS**

Michael Felkel1,#, Robert Exner1,#, Regina Schleucher1,#, Helga Lay2 , Ingo B. Autenrieth² , Volkhard A. J. Kempf³ and Julia-Stefanie Frick2,#,*

¹Institute of Tropical Medicine, University Hospital of Tübingen, Tübingen, Germany ²Institute of Medical Microbiology and Hygiene, University Hospital of Tübingen, Tübingen, Germany ³Institute of Medical Microbiology and Hygiene, Johann-Wolfgang Goethe University, Frankfurt, Germany

Received: September 2, 2013; Accepted: September 15, 2013

The incidence of tuberculosis (TB) and especially multidrug-resistant TB (MDR) continues to increase alarmingly worldwide, and reliable and fast diagnosis of MDR is essential for the adequate treatment of patients. In contrast to the standard culture methods, nucleid acid amplification tests (NAATs) provide information about presence of *Mycobacterium tuberculosis* complex (MTBC) DNA and a potential resistance pattern within hours. We analyzed specimens of 110 patients from Nigeria comparing culture-based drug susceptibility testing (DST) to NAAT assays detecting isoniazid (INH), rifampicin (RMP) (GenoType MTBDRplus), and ethambutol (EMB) (GenoType MTBDRsl) resistance. Compared to DST, the GenoType MTBDRplus and MTBDRsl showed a specificity of 100% (86.3–100) and a sensitivity of 86% (42.1–99.6%) for detection of INH and a specificity of 100% (86.3–100) and a sensitivity of 83% (35.9–99.6%) for detection of RMP, and a sensitivity 100% (47.8–100%) for EMB resistance. However, in two strains, the NAAT assays provided false susceptible results as the mutations causing resistance were in genomic regions not covered by the probes of the GenoType MTBDRplus assay. We show that, in combination to DST, application of the GenoType MTBDRplus and GenoType MTBDRsl assays might be a useful additional tool to allow a rapid and safe diagnosis of MDR and extensively drugresistant (XDR) MTBC.

Keywords: *Mycobacterium tuberculosis* complex, drug susceptibility testing, GenoType, Nigeria

Introduction

The incidence of tuberculosis (TB) and especially multidrug-resistant TB (MDR) continues to increase alarmingly worldwide. Nigeria has a population of about 150 million people and ranks among the world's high-burden countries. The World Health Organization (WHO) report 2009 stated an incidence of TB at 311/100,000, a prevalence of 772,000 cases, and the rate of new sputum smear positive cases was approximately 131/100,000 per year. Studies on drug susceptibility profiles of *Mycobacterium tuberculosis* (MTB) in Nigeria are scarcely reported, and the prevalence of MDR strains is numeralized between 2% and 53% in new cases of TB [1, 2]. Treatment on MDR TB involves second line drugs which are costly, scarce, and toxic and require longer duration of therapy [2]. Therefore, a rapid diagnosis of MDR TB is essential to medicate patients adequate to the drug resistance pattern.

In 2006, the outbreak of rapidly fatal TB among hospitalized HIV patients in South Africa [3] led to the definition of extensively drug-resistant tuberculosis (XDR TB) as TB resistant to fluoroquinolones (FQs) and an injectable second line drug (amikacin, kanamycin, or capreomycin) [3, 4].

The GenoType MTBDRplus assay (Hain) allows identification of isoniazid (INH) and rifampicin (RMP) resistant TB stains directly from clinical isolates. This assay detects a broad variety of *rpoB* gene mutations and is described to identify more than 95% of all RMP resistant strains [5, 6]. The mutations causing isoniazid (INH) resistance are located in several genes and regions like in 50%–95% in the codon 315 of the *katG* gene [5–8], in 20%–35% in the *inhA* regulatory region [5, 6, 8, 9] and in 10%–15% in the *ahpC*–*oxyR* intergenic region [5, 6, 9, 10], often in conjunction with *katG* mutations outside of the codon 315 [11].

^{*}Corresponding author: Dr. Julia-Stefanie Frick; Institute of Medical Microbiology and Hygiene, University of Tübingen, Elfriede-Aulhorn-Str.6, D-72076 Tübingen, Germany; Tel.: 0049 7071 29 82352; Fax: 0049 7071 29 5440; E-mail: julia-stefanie.frick@med.uni-tuebingen.de

[#] These authors contributed equally to this work and share first authorship.

In order to rapidly detect second line drug resistance, the GenoType MTBDRsl test has been developed. This assay detects mutations in the *gyrA*, *rrs*, and *embB* genes and, therefore, resistance to fluoroquinolone (FQ) , amikacin (AM)/capreomycin (CM), and ethambutol (EMB). The sensitivity for detection of FQ resistance is described as about 75%–90% with specificity of 100%, for aminoglycoside resistance sensitivity, and specificity was 100% [4]; other studies report on sensitivity of 83.3% for AM and 86.8% for CM with a specificity of 100% for detection of AM and 99.1% of CM resistance [12]. For detection of EMB resistance, the sensitivity is indicated with 59%–64% and specificity 100% [4, 12]. Within the present study, we analyzed 110 specimens from Nigerian patients using the MTBDRplus and MTBDRsl tests. The major scientific objective was to determine the accuracy of these tests for direct detection of MTB and for determination of resistance pattern of MTB directly from the patients' samples.

Materials and methods

Clinical specimens

A total of 110 sputum specimens sent to our laboratory from a tertiary hospital in Hospital Nigeria were processed by conventional N-acetyl-L cystein-NaOH method (final NaOH concentration 1%). After decontamination, the concentrated sediment was suspended in 1.0 to 1.5 ml sterile phosphate. After inoculation of solid and liquid media for the growth detection and the GenoType MTBDRplus assay (Hain), the leftover sediment was stored at −20 °C.

DST

DST (drug susceptibility testing) with INH, RMP, EMB, STR, and PZA was performed by the BACTEC MGIT 960 method (Becton Dickinson, Heidelberg, Germany). Tests were performed with the standard critical concentrations of INH (0.1 μg/ml), RMP (1.0 μg/ml), EMB (5.0 μg/ml), STR $(1.0 \mu g/ml)$, and PZA $(100 \mu g/ml)$.

Identification of the MTBC strains from clinical specimens

For identification and differentiation of the MTBC strains isolated from the grown cultures, we used the genotype MTBC assay as recommended by the manufacturer (Hain Lifescience, Nehren, Germany).

GenoType MTBDRplus assay

From all specimens, the strip assay was performed as recommended by the manufacturer (Hain Lifescience, Nehren, Germany).

GenoType MTBDRsl assay

From the leftover sediment of the specimens that turned out in the DST to harbor RMP and INH resistant strains, the strip assay was performed as recommended by the manufacturer (Hain Lifescience, Nehren, Germany).

Sequencing

Sequencing of the key regions, involved in the development of resistance, was performed by Hain Lifescience (Nehren, Germany) as published previously [5].

Statistical analysis

Sensitivity and specificity were calculated from cross tabulation, comparing the strip assays (MTBDRplus and MT-BDRsl) (Hain Lifescience, Nehren, Germany) with gold standard culture techniques.

Results

Host characteristics

Of the 110 specimens, in 40, we were able to either cultivate MTBC or to detect MTBC complex DNA. Of these patients, 18 were female and 22 male, the mean age of the female patients was 37.27 years (14–82 years), and the mean age of the male patients was 33.28 years (18–82 years). Six females and four males were tested HIV positive, and seven females and ten males had a history of TB treatment or were undergoing TB treatment *(Table 1)*.

Table 1. Patient characteristics

Frequency of MTBC detected by culture or GenoType MTBDRplus assay

We analyzed the frequency of MTBC in the specimens using two different methods. The specimens were inoculated in liquid and on three solid media and incubated at 37 °C for at least 7 weeks. The solid media were checked weekly for growth of mycobacteria; growth in the liquid media was monitored by the BACTEC MGIT 960 method.

In cultures from patients without history of MTBC treatment, we detected, in 19 (20.3%) samples, growth

Table 2A. Frequency of MTBC detected by culture or Genotype MTBDRplus strip assay in patients without history of MTBC treatment

		Strip assay MTBDR plus (HAIN)		
Culture				
	_			

Table 2B. Frequency of MTBC detected by culture or Genotype MTBDRplus strip assay in patients undergoing MTBC treatment

of *M. tuberculosis* complex; 74 samples (79.56%) were negative for growth of *M. tuberculosis* up to 10 weeks of incubation. Using the GenoType MTBDRplus assay, 22 samples (23.6%) were tested positive for presence of *M. tuberculosis* complex DNA and 71 (76.3%) were tested negative. Four specimens showed positive detection of *M. tuberculosis* complex DNA by the GenoType MTBDRplus assay but negative results of culture; one specimen was culture positive but GenoType MTBDRplus negative. Compared to the results from the culture growth, the GenoType MTBDRplus assay showed specificity for detection of *M. tuberculosis* complex DNA of 95% (74.0– 99.9%) and a sensitivity of 95% (86.7–98.5%) *(Table 2A)*.

In patients undergoing MTBC treatment $(n = 16)$, we detected, in 12 (75%) samples, growth of *M. tuberculosis* complex; four samples (25%) were negative for growth of *M. tuberculosis* up to 10 weeks of incubation. Using the GenoType MTBDRplus assay, 14 samples (87.5%) were tested positive for presence of *M. tuberculosis* complex DNA and 2 (12.5%) were tested negative. From the 16 specimens, four showed positive detection of *M. tuberculosis* complex DNA by the GenoType MTBDRplus

> assay but negative results of culture, and two specimens were culture positive but GenoType MTBDRplus negative. Therefore, compared to the results from culture as gold standard technique, the GenoType MTBDRplus assay showed a sensitivity for detection of *M. tuberculosis* complex of 83% (51.6–97.9%) in patients undergoing treatment of MTBC *(Table 2B)*.

> For differentiation and identification of the grown MTBC strains, we used the GenoType MTBC. We identified 21 (67%) *M. tuberculosis* and 11 (34%) *M. africanum* strains *(Fig. 1A)*.

Frequency of drug resistance and MDR detected by culture methods

Using the BACTEC MGIT 960 method, we determined the drug resistance and susceptibility pattern of the cultivated MTB strains. Drug susceptibility testing (DST) revealed 21 (65.6%) of INH, RMP, PZA, STR, and EMB sensitive strains. Two strains (6.25%) showed a mono-resistance for either INH or STR. One strain (3.125%) was RMP and PZA resistant. We isolated five MDR MTB

Fig. 1. (A) Differentiation and identification of the MTBC strains obtained from culture. (B) Drug resistance pattern detected by culture based drug susceptibility testing (DST)

stains, of which one strain (3.125%) was resistant to INH, RMP, and STR and four strains (12.5%) were resistant to INH, RMP, STR, EMB, and PZA *(Fig. 1B)*. The MDR strains and the RMP/PZA resistant strain were sent to the National Reference Laboratory in Borstel for verification of the resistance pattern and for additional testing of second line drugs. All six strains were susceptible to protionamide, ofloxacin, PAS, cycloserine/terizidone, amikacin, and capreomycin**.**

Frequency of drug resistance and MDR detected by usage of the GenoType MTBDRplus and MTBDRsl assays

In order to analyze INH and RMP directly from the sediments of decontaminated sputum samples, we used the GenoType MTBDRplus assay and compared the results with the results from DST, which was used as gold standard technique. Twenty-six specimens (81.25%) were tested INH susceptible, six (18.75%) were tested INH resistant; this corresponds to a specificity of 100% (86.3–100%) and a sensitivity of 86% (42.1–99.6%). Twenty-seven specimens (84.4%) were tested RMP susceptible, five (15.6%) were tested RMP resistant, corresponding to a specificity of 100% (86.3–100%) and a sensitivity of 83% (35.9–99.6%) *(Table 3)*. In case of INH and/or RMP resistance, we additionally performed the GenoType MTBDRsl assay, which detects for fluoroquinolone (FQ), amikacin/capreomycin (AM/CM) , and EMB resistance. In five from five tested samples, EMB resistance was detected (sensitivity 100%) (47.8–100%), whereas no FQ or AM/CM resistance was detectable (specificity 100%) (54–100%) *(Table 3)*. This

is in line with the DST pattern from the liquid cultures, performed by the German National Reference Centre for Mycobacteria (Borstel, Germany).

Resistance to RMP and INH in *M. tuberculosis* is most often attributed to mutations in the *rpoB*, *katG*, and *inhA* genes. The GenoType MTBDRplus assay detects INH and RMP resistance-associated mutations in *rpoB*, *katG*, and *inhA* genes. The GenoType MTBDRsl assay detects mutations *gyrA*, *rrs*, and *embB* genes leading to resistance to FQ, AM/CM, and EMB. We found different RMP, INH, and EMB pattern; in case of INH resistance, most of the assays showed *katG* mutation specific bands with probe MUT 1, whereas in RMP and EMB resistance, pattern was diverse *(Table 4)*.

In sample no. 9, the GenoType MTBDRplus resulted in RMP susceptibility, whereas the DTS showed an RMP resistance of the strain. In sample no. 16, the GenoType MTBDRplus resulted in INH susceptibility, but the DTS showed an INH resistance of the strain. In both strains, sequencing of the *rpoB* gene and, respectively, the *inhA* promoter region and *katG* gene did not show any resistance-associated mutations. This might lead to the conclusion that mutations in regions not covered by the probes of the GenoType MTBDRplus assay might cause the resistant phenotype.

Discussion

Within the present study, we analyzed 110 respiratory specimens from Nigerian patients for infection with MTB and analyzed the MTB drug resistance pattern. Only few publications that deal with the susceptibility pattern of bacteria

Table 3. Sensitivity and specificity of the strip assay GenoType MTBDR plus assay and the GenoType MTBDR sl assay in comparison to gold standard method culture

			MTBDR plus				MTBDR sl					
			INH		RIF		ETB		Chinolone		Aminoglycoside/ capreomycin	
			sen	res	sen	res	sen	res	sen	res	sen	res
Culture	Resistance	sen	25(32)	0(32)	26(32)	0(32)		θ	6(6)	0	6(6)	$\bf{0}$
		res	(32)	6(32)	(32)	5(32)	5(5)	θ	$_{0}$	0	0	

Table 4. Detection of INH and RMP resistance-associated mutations in *rpoB, katG,* and *inhA* genes by the strip assay GenoType MTBDR plus assay and mutations in the *gyrA, rrs,* and *embB* genes leading to resistance to FQ, AM/CM, and EMB by the Geno-Type MTBDR sl assay

European Journal of Microbiology and Immunology 3 (2013) 4

of the *M. tuberculosis* complex in Nigeria are available. In our study, 68% of the isolated stains were pan-susceptible to all five drugs (INH, RMP, STR, PZA, and EMB) tested in the DST. This is in line with other studies from Nigeria, reporting on 59–62% pan-susceptible isolates [1, 13]. However, data collected in a tertiary hospital in Nigeria report on pan-susceptibility in only 46% of the isolated TB strains [2]. We found 6.25% of INH or STR mono-resistant strains. Other publications report on a streptomycin resistance rate of 11% but did not find INH mono-resistant strains [2]. The low number of isolates ranging from 32 (this study) to 50 [1, 2] might account for that effect. However, we identified 16% of the isolates as MDR TB stains; these strains were at least resistant to INH, RMP, STR, and EMB, and one was additionally resistant to PZA. The MTB prevalence described in the literature ranges between 4% and 31% [1, 2]. It is suggested that the incidence of MDR tuberculosis is associated to the type of hospital. Tertiary hospitals show a higher rate of MDR strains as compared to other hospitals, as patients treated in these hospitals might represent a selection of referred difficult cases that may have failed or defaulted from previous treatment [1, 2]. In our study, the specimens were obtained from patients treated in a tertiary hospital. However, for a better understanding of the epidemiological situation of MDR and XDR tuberculosis in Nigeria further studies are necessary, including large patient cohorts from different parts of the country.

Interestingly, our data revealed a sensitivity and specificity of 95% of the GenoType MTBDRplus assay on detection of *M. tuberculosis* complex DNA in the specimens of patients without history of TB treatment. Other commercial nucleid acid amplification tests (NAATs) for diagnosis of tuberculosis directly from respiratory specimens report on sensitivity between 85% and 92% and specificity between 96% and 97% [14, 15]. Therefore, the GenoType MTBDRplus assay seems to reliably detect *M. tuberculosis* complex DNA and INH and RMP resistance. Certainly, NAATs should be applied only in combination with gold standard culture methods.

As a rapid diagnosis of MDR or even XDR tuberculosis is of a high importance for the patient outcome and of a high epidemiological importance, we evaluated the GenoType MTBDRplus and GenoType MTBDRsl assays on the patient's specimens. These tests allow information about the MTB resistance pattern within 1 day. The time to positivity of MTB cultures can take up to 7 weeks, and the conventional DST testing takes additional 2 weeks after cultivation of the bacterium. We found a sensitivity of 86% and a specificity of 100% for INH resistance and a sensitivity of 83% and specificity of 100% for RMP resistance. Other studies report on sensitivity for INH resistance of 90% and 98% for RMP resistance [5]. We had two strains that were tested resistant to INH or, respectively, RMP in the DST and susceptible in the GenoType MTBDRplus. In both strains, sequencing of the *rpoB* gene and respectively the *inhA* promoter region and *katG* gene did not show any resistance associated mutations. This might lead to the

conclusion that mutations in regions not covered by the probes of the GenoType MTBDRplus assay might cause the resistant phenotype. Analysis of EMB, AM/CM, and FQ resistance using the GenoType MTBDRsl corresponded to the results obtained by DST; however, the number of isolates was low, and we did not detect any resistance to FQ or AM/CM. As mutations in loci other than those covered by the GenoType MTBDRplus and GenoType MTBDRsl assays might lead to false sensitive results, the NAAT-based methods need to be confirmed by conventional DST, based on culture techniques. Therefore, the stripe assays GenoType MTBDRplus and GenoType MT-BDRsl cannot replace gold standard culture methods but represent an additional tool for the first identification of drug resistant *M. tuberculosis* complex strains, resulting in a fast, adapted and more sufficient therapy.

Taken together, our data represent an important addition to the rare epidemiological data concerning resistance pattern of MTB in Nigeria and showed that the application of the GenoType MTBDRplus and GenoType MTBDRsl assays might be a useful additional tool to allow a rapid and safe diagnosis of MDR and XDR tuberculosis.

References

- 1. Ani AE, Idoko J, Dalyop YB, Pitmang SL: Drug resistance profile of *Mycobacterium tuberculosis* isolates from pulmonary tuberculosis patients in Jos, Nigeria. Trans R Soc Trop Med Hyg 103, 67–71 (2009)
- 2. Kehinde AO, Obaseki FA, Ishola OC, Ibrahim KD: Multidrug resistance to *Mycobacterium tuberculosis* in a tertiary hospital. J Natl Med Assoc 99, 1185–1189 (2007)
- 3. Gandhi NR, Moll A, Sturm AW, Pawinski R, Govender T, Lalloo U, Zeller K, Andrews J, Friedland G: Extensively drug-resistant tuberculosis as a cause of death in patients coinfected with tuberculosis and HIV in a rural area of South Africa. Lancet 368, 1575–1580 (2006)
- 4. Kiet VS, Lan NT, An DD, Dung NH, Hoa DV, van Vinh CN, Chinh NT, Farrar J, Caws M: Evaluation of the MTBDRsl test for detection of second-line-drug resistance in *Mycobacterium tuberculosis*. J Clin Microbiol 48, 2934–2939 (2010)
- 5. Hillemann D, Rusch-Gerdes S, Richter E: Evaluation of the GenoType MTBDRplus assay for rifampin and isoniazid susceptibility testing of *Mycobacterium tuberculosis* strains and clinical specimens. J Clin Microbiol 45, 2635–2640 (2007)
- 6. Telenti A, Honore N, Bernasconi C, March J, Ortega A, Heym B, Takiff HE, Cole ST: Genotypic assessment of isoniazid and rifampin resistance in *Mycobacterium tuberculosis*: a blind study at reference laboratory level. J Clin Microbiol 35, 719–723 (1997)
- 7. Mokrousov I, Narvskaya O, Otten T, Limeschenko E, Steklova L, Vyshnevskiy B: High prevalence of KatG Ser315Thr substitution among isoniazid-resistant *Mycobacterium tuberculosis* clinical isolates from northwestern Russia, 1996 to 2001. Antimicrob Agents Chemother 46, 1417–1424 (2002)
- 8. Musser JM, Kapur V, Williams DL, Kreiswirth BN, van SD, van Embden JD: Characterization of the catalase-per-

oxidase gene (katG) and inhA locus in isoniazid-resistant and -susceptible strains of *Mycobacterium tuberculosis* by automated DNA sequencing: restricted array of mutations associated with drug resistance. J Infect Dis 173, 196–202 (1996)

- 9. Piatek AS, Telenti A, Murray MR, El-Hajj H, Jacobs WR, Jr., Kramer FR, Alland D: Genotypic analysis of *Mycobacterium tuberculosis* in two distinct populations using molecular beacons: implications for rapid susceptibility testing. Antimicrob Agents Chemother 44, 103–110 (2000)
- 10. Kelley CL, Rouse DA, Morris SL: Analysis of ahpC gene mutations in isoniazid-resistant clinical isolates of *Mycobacterium tuberculosis*. Antimicrob Agents Chemother 41, 2057–2058 (1997)
- 11. Sreevatsan S, Pan X, Zhang Y, Deretic V, Musser JM: Analysis of the oxyR-ahpC region in isoniazid-resistant and -susceptible *Mycobacterium tuberculosis* complex organisms recovered from diseased humans and animals in di-

verse localities. Antimicrob Agents Chemother 41, 600–606 (1997)

- 12. Hillemann D, Rusch-Gerdes S, Richter E: Feasibility of the GenoType MTBDRsl assay for fluoroquinolone, amikacin– capreomycin, and ethambutol resistance testing of *Mycobacterium tuberculosis* strains and clinical specimens. J Clin Microbiol 47, 1767–1772 (2009)
- 13. Ani A, Bruvik T, Okoh Y, Agaba P, Agbaji O, Idoko J, Dahle UR: Genetic diversity of *Mycobacterium tuberculosis* complex in Jos, Nigeria. BMC Infect Dis 10, 189 (2010)
- 14. Barrett A, Magee JG, Freeman R: An evaluation of the BD ProbeTec ET system for the direct detection of *Mycobacterium tuberculosis* in respiratory samples. J Med Microbiol 51, 895–898 (2002)
- 15. Ling DI, Flores LL, Riley LW, Pai M: Commercial nucleicacid amplification tests for diagnosis of pulmonary tuberculosis in respiratory specimens: meta-analysis and meta-regression. PLoS ONE 3, e1536 (2008)