



Published in final edited form as:

Curr Protoc Immunol. ; 106: 5.10.1–5.10.12. doi:10.1002/0471142735.im0510s106.

Measurement of peptide binding to MHC class II molecules by fluorescence polarization

Liusong Yin¹ and Lawrence J. Stern^{1,2,3}

¹Department of Pathology, University of Massachusetts Medical School, Worcester, MA01605

²Department of Biochemistry and Molecular Pharmacology, University of Massachusetts Medical School, Worcester, MA01605

Abstract

Peptide binding to major histocompatibility complex class II (MHCII) molecules is a key process in antigen presentation and CD4⁺ T cell epitope selection. This unit describes a fairly simple but powerful fluorescence polarization-based binding competition assay to measure peptide binding to soluble recombinant MHCII molecules. The binding of a peptide of interest to MHCII molecules is assessed based on its ability to inhibit the binding of a fluorescence-labeled probe peptide, with the strength of binding characterized as IC₅₀ (concentration required for 50% inhibition of probe peptide binding). Data analysis related to this method is discussed. In addition, this unit includes a support protocol for fluorescence labeling peptide using an amine-reactive probe. The advantage of this protocol is that it allows simple, fast, and high throughput measurements of binding for a large set of peptides to MHCII molecules.

Keywords

MHCII; Peptide binding; Fluorescence polarization; binding competition; IC₅₀

INTRODUCTION

In this unit, we describe a simple and high-throughput fluorescence polarization-based assay to measure peptide binding to MHCII molecules (Basic Protocol). A probe peptide which has a high affinity to MHCII is labeled with a fluorochrome. Unlabeled target peptide is co-incubated with the probe peptide competing for the binding to MHCII. Incubation is typically 72 hours at 37 °C. The binding of labeled probe peptide to MHCII is measured by fluorescence polarization (FP). A Support Protocol is provided to enable labeling of the probe peptide with Alexa488 fluorochrome. IC₅₀, the concentration of unlabeled target peptide required to inhibit the binding of labeled probe peptide by 50% can be obtained by plotting relative binding versus unlabeled target peptide concentration. IC₅₀ represents a reasonable measurement of peptide binding affinity to MHCII molecules.

³Correspondence: Dr. Lawrence J. Stern (Lawrence.Stern@umassmed.edu).

BASIC PROTOCOL

MEASUREMENT OF PEPTIDE BINDING TO MHCII MOLECULES AS IC₅₀ USING A FLUORESCENCE POLARIZATION-BASED BINDING COMPETITION ASSAY

—This protocol describes an assay to measure the binding of peptides to MHCII molecules. Fluorescence-labeled probe peptide is incubated with MHCII molecules, together with serial dilutions of unlabeled target peptide at 37 °C for 72 hours. Binding affinity of target peptide is assessed by IC₅₀, which is the concentration of unlabeled target peptide required to inhibit the binding of labeled probe peptide by 50%. IC₅₀ can be determined by plotting the binding competition curve. A flow chart schematizing the assay is presented in Figure 1.

Materials: Alexa488 labeled influenza hemagglutinin-derived probe peptide (Alexa488-FRR-HA₃₀₆₋₃₁₈; see support protocol for the labeling)

Unlabeled target peptide Double-distilled water (ddH₂O; Millipore)

Dimethylsulfoxide (DMSO)

Purified MHCII molecules (HLA-DR1 in this protocol) (see Unit 18.3) [*Copy ED: Query as to whether HLA-DR1 is available commercially – RC]

Phosphate-buffered saline (PBS), pH7.4 (see recipe)

Binding buffer, pH5.5 (see recipe)

Protease inhibitor cocktail (prepare fresh; see recipe)

1.5 ml black microtube

1.5 ml clear microtube

96-well Polypropylene flat-bottom non-binding black Microplate (Greiner)

Aluminum Sealing Foil for 96-well Microplate (USA Scientific)

37 °C incubator

VICTOR X5 Multilabel Plate Reader (PerkinElmer), POLARstar OPTIMA Plate Reader (BMG LABTECH) or other fluorescence microplate reader equipped for fluorescence polarization

GraphPad Prism 6 graphing and data analyzing software

Prepare peptides: 1. Solubilize lyophilized Alexa488-FRR-HA₃₀₆₋₃₁₈ in a 1.5 ml black microtube with ddH₂O to the stock concentration of 200 μM.

See support protocol for the sequence information and labeling procedure of Alexa488-FRR-HA₃₀₆₋₃₁₈. It is extremely important to keep the fluorescence-labeled peptide in dark and minimize the exposure to light.

2. Solubilize each lyophilized unlabeled target peptide in a 1.5 ml clear microtube with 100% DMSO to the concentration of 20 mg/ml by adding 50 μ l 100% DMSO to 1mg peptide. Then, dilute the peptide into 50% DMSO to the stock concentration of 400 μ M.

The unlabeled target peptides were synthesized at 21st Century Biochemicals with acetylated N-termini and shipped to us as 1mg lyophilized powder in 1.5ml clear microtube. If test peptides are soluble in ddH₂O, one could use ddH₂O to solubilize each peptide, instead of DMSO here. We have verified that the DMSO contained in peptide stock solution up to 2.5% does not interfere with following peptide binding or FP reading.

Prepare binding competition reaction mix: 3. In the 96-well plate, prepare 200 μ l binding competition reaction mixture in each well, containing 100 nM MHCII, 25 nM Alexa488-FRR-HA₃₀₆₋₃₁₈, and series dilution of target peptide, starting at 20 μ M, 1 to 5 dilution, leaving the last well without competing target peptide for determination of FP value for no inhibition. Control wells with 25 nM Alexa488-FRR-HA₃₀₆₋₃₁₈ in buffer but without MHCII are also included for determination of FP value for free peptide.

The reaction is performed in pH5.5 binding buffer, with freshly added protease inhibitor cocktail and incubated at 37 °C for 72 hours. Recombinant MHCII is purified as previously described (Frayser et al., 1999; Stern and Wiley, 1992) and stored in PBS. An example of plate set up can be found in Table 1.

The concentrations for MHCII (100 nM) and labeled probe peptide (25 nM) have been optimized for this assay. The labeled probe peptide is set at lowest concentration for which reliable fluorescence polarization values can be obtained. The MHCII concentration is set by titrating MHCII against fixed labeled peptide concentration, with the selected MHCII concentration which allows for ~50-75% of the maximum binding as judged by polarization assay. It is crucial to optimize the concentrations of MHCII and labeled probe peptide for appropriate FP reading. The starting concentration and dilution factor for serial dilution of target peptides can be changed depending on specific assay configuration.

Detect amount of MHCII-Alexa488-FRR-HA₃₀₆₋₃₁₈ complex by FP: 4. After 72 hr incubation, read FP of the plate using appropriate microplate reader at excitation 485 nm and emission 520 nm.

The key point of this assay is to use FP to measure the fraction of labeled probe peptide bound in complex with MHCII. The FP value for free Alexa488-FRR-HA₃₀₆₋₃₁₈ is ~35 mP, and for purified MHCII- Alexa488-FRR-HA₃₀₆₋₃₁₈ complex is ~350. It is important that the plate reader has been calibrated and tested for polarization measurements. Usually, a standard stock solution of fluorescein is used to calibrate the instrument following the manufacturer's recommendations, to give a FP value of 27 mP in aqueous alkali buffer at room temperature. To simulate the FP of the labeled probe peptide in the bound form with MHCII, we solubilize the labeled probe peptide in different percentage of glycerol. In general, 0% glycerol should give the FP value for free form, ~35 mP; 20% glycerol should give the FP value for 50% bound, ~200 mP; and 50% glycerol should give the FP value for

100% bound, ~350 mP. Alternately, purified MHCII- Alexa488-FRR-HA₃₀₆₋₃₁₈ complex can be used.

Data analysis and calculation of IC₅₀: 5. Use the average of the FP values of wells with control free Alexa488-FRR-HA₃₀₆₋₃₁₈ (no MHCII) as the FP_{free}. Use the average of the FP values of wells without competing target peptide as FP_{no_comp}.

6. Calculate Relative binding for the other wells as $(FP_{\text{sample}} - FP_{\text{free}}) / (FP_{\text{no_comp}} - FP_{\text{free}})$

7. Plot Relative binding versus concentration of unlabeled target peptide, and fit the curve into equation $y = 1 / (1 + [pep] / IC_{50})$, where [pep] is the concentration of unlabeled target peptide; and IC₅₀ is the 50% inhibition concentration.

An example of the binding competition curve and calculation of IC₅₀ for peptides with different binding affinities is shown in Figure 1C.

SUPPORT PROTOCOL

LABELING OF THE PROBE PEPTIDE WITH A FLUOROCHROME—This protocol describes the labeling and purification of a probe peptide FRR-HA₃₀₆₋₃₁₈ with Alexa488 fluorochrome through amine-reactive probe. FRR-HA₃₀₆₋₃₁₈ with single lysine at position 310 (K310) is incubated with Alexa488-TFP at room temperature for 1 hour in 150 mM pH9.0 sodium bicarbonate buffer. The labeled peptide is purified using Jupiter 5u C18 column and lyophilized to store as peptide powder. The procedure can be adapted to label other probe peptides with single lysine in the sequence. The labeled probe peptide can be used to study peptide association and dissociation to MHCII molecules, and to measure IC₅₀ for other unlabeled peptides using as probe peptide in binding competition assay. A flow chart schematizing the protocol is presented in Figure 2.

As an alternate approach, peptides containing a single cysteine residue can be labeled efficiently using thiol-specific reagents, such as Alexa488-maleimide. In this case the binding reaction is done at pH7.5 (for example conveniently in phosphate-buffered saline). The presence of free amino groups (N-terminus and lysine) in the peptide or amine-containing compounds in the reaction buffer is immaterial; however it is crucial that the peptide thiol group should be reduced before reaction. In most cases purified synthetic peptides stored after lyophilization from acidic buffer will be fine, but oxidized peptides will need to be reduced (for example using 50 mM dithiothreitol) and repurified before reaction.

Materials: Influenza hemagglutinin (HA) analog peptide FRR-HA₃₀₆₋₃₁₈ (Ac-PRFVKQNTLRLAT) Alexa Fluor 488 Carboxylic Acid, 2,3,5,6-Tetrafluorophenyl Ester (Alexa488-TFP)

1.5 ml black microtube

150 mM sodium bicarbonate buffer, pH9.0 (or other non-amine buffer) (made fresh)

Double-distilled water (ddH₂O)

Trifluoroacetic acid (TFA)

Acetonitrile

Buffer A (5% TFA in ddH₂O)

Buffer B (5% TFA in acetonitrile)

Jupiter 5u C18 300A 250x4.6mm column (C18 column; Phenomenex)

Agilent 1200 series Liquid Chromatography system (LC; Agilent Technologies) Lyophilizer

Prepare labeling reaction mix: 1. Dissolve 1 mg FRR-HA₃₀₆₋₃₁₈ in 500 µl 150 mM pH9.0 sodium bicarbonate buffer, and then transfer to a 1.5 ml black microtube containing 1 mg Alexa488-TFP. Reserve 5 µg of FRR-HA₃₀₆₋₃₁₈ before transfer to Alexa488-TFP for analysis. Also reserve some Alexa488-TFP by dipping a pipette tip into the fluorochrome powder.

It is extremely important that the labeling reaction mix be free of amine-containing substances such as Tris, glycine, ammonium ions. In order to label each peptide with only one fluorochrome, it is crucial that the peptide has only one active amine group, which is the K310 on FRR-HA₃₀₆₋₃₁₈. If the peptide being labeled or the fluorochrome being used is not soluble in sodium bicarbonate buffer, one can first dissolve it in DMSO or DMF at 10 mg/ml stock solution, and then dilute it into sodium bicarbonate buffer. We observe efficient labeling with molar ratio of fluorochrome to peptide of 2:1 to 5:1. In the protocol described, the ratio is ~2:1. Other fluorescent probes such as fluorescein or tetramethylrhodamine can be used instead of Alexa488, but Alexa488 fluorescence is relatively pH-independent, environment-independent, and relatively stable to photobleaching (Molecular Probes Handbook) that bound and free labeled peptides have clearly distinct FP values. As such it is nearly ideal for this purpose. Other amine-specific modifying groups can be used, but TFP (trifluorophenyl ester) is more stable to hydrolysis in aqueous solution than other commonly used reactive groups such as N-hydroxysuccinimides or imidoesters.

2. Incubate the labeling reaction mix at room temperature for 1 hour with mixing every 10 minutes.

Analysis of labeling by reverse phase LC: 3. Run unlabeled probe peptide FRR-HA₃₀₆₋₃₁₈, free Alexa488-TFP, and a small sample (0.2%) of the labeling reaction mixture over C18 column measuring absorbance at wavelengths to detect both the peptide and the fluorochrome. We use 214 nm, 280 nm and 450 nm. Each sample should be diluted into 100 µl Buffer A before injection onto an equilibrated column.

For unlabeled FRR-HA₃₀₆₋₃₁₈, dilute the reserved 5 µg into 100 µl Buffer A. For free Alexa488-TFP, transfer the reserved powder into 100 µl Buffer A.

To check the progress of the labeling reaction and to identify the labeled peptide peak from unlabeled peptide and free fluorochrome peaks, we aliquot 1 µl of the 500 µl labeling reaction mix and dilute it into 100 µl Buffer A and inject samples of unreacted peptide,

unreacted fluorochrome, and reaction mix. The unlabeled FRR-HA₃₀₆₋₃₁₈ should have strong absorbance at 214 nm, weak at 280 nm and not at 450 nm. The free Alexa488-TFP should have strong absorbance at 214 nm, 280 nm and 450 nm, with elution time distinct from that of unlabeled FRR-HA₃₀₆₋₃₁₈. The labeling reaction mix should show a new peak, having absorbance at 214 nm, 280 nm, and 450 nm, which corresponds to the Alexa488-labeled FRR-HA₃₀₆₋₃₁₈. If so, we can proceed to the purification of Alexa488-FRR-HA₃₀₆₋₃₁₈. If a substantial peak for unlabeled FRR-HA₃₀₆₋₃₁₈ still shows up in the labeling reaction mix, we should increase the ratio of fluorochrome to peptide and/or incubate the labeling reaction mix for a longer time. If the peak for Alexa488-FRR-HA₃₀₆₋₃₁₈ is so close to that of unlabeled FRR-HA₃₀₆₋₃₁₈, we should adjust the elution gradient until they are separate enough for purification. The finalized Buffer B gradient for elution is listed in Table II, and examples of elution traces for unlabeled FRR-HA₃₀₆₋₃₁₈, free Alexa488-TFP, reaction mix and Alexa488-labeled FRR-HA₃₀₆₋₃₁₈ are shown in Figure 2B.

Purification of Alexa488-labeled FRR-HA₃₀₆₋₃₁₈ by reverse phase LC: 4. Inject the remaining ~500 µl of the labeling reaction mix and collect fractions for Alexa488-FRR-HA₃₀₆₋₃₁₈.

We collect 50 fractions during a 10 minutes period, and each vial has 200 µl. One can adjust the number of fractions collected and the amount in each fraction based on the elution profiles.

5. Reinject 1 µl of the pooled fractions of Alexa488-FRR-HA₃₀₆₋₃₁₈ into C18 column to check the purity.

It is recommended to check the labeled probe peptide by mass spectrometry and functional analysis, such as peptide binding in our case.

6. Measure the concentration of Alexa488-FRR-HA₃₀₆₋₃₁₈ by measuring the absorbance of the fluorochrome at its excitation maximum and store as lyophilized powder in -20 °C in dark.

We measure the absorbance of Alexa488-FRR-HA₃₀₆₋₃₁₈ from 200 nm to 600 nm and find that the maximal absorbance is at 495 nm. [Alexa488-FRR-HA₃₀₆₋₃₁₈]=Absorbance₄₉₅/71000 (M), where 71000 is the extinction coefficient at maximal wavelength for fluorochrome Alexa488.

REAGENTS AND SOLUTIONS

Binding buffer, pH5.5: 100 mM sodium citrate (mix of 0.46% w/v anhydrous citric acid and 2.23% w/v sodium citrate dihydrate)

50 mM sodium chloride

5 mM EDTA

0.1% (w/v) octylglucoside

0.1% (w/v) iodoacetamide

0.1% (w/v) sodium azide

Adjust the pH if necessary to 5.5 using sodium hydroxide and hydrochloric acid Store up to 3 months at room temperature

Phosphate-buffered saline: 137 mM sodium chloride

2.7 mM potassium chloride

10 mM disodium hydrogen phosphate

1.8 mM potassium dihydrogen phosphate

Adjust the pH to 7.4 using hydrochloric acid

Store up to 6 months at room temperature

Proteinase inhibitor cocktail (100X): AEBSF (4-(2-Aminoethyl)benzenesulfonyl fluoride hydrochloride)

Aprotinin

Bestatin hydrochloride

E-64 (N-(trans-Epoxy succinyl)-L-leucine 4-guanidinobutylamide)

EDTA

Leupeptin hemisulfate salt

This protease inhibitor cocktail is bought from Sigma-Aldrich as a lyophilized powder (catalog number P2714). We dissolve it into 10 ml ddH₂O and use it as 100X. Store up to 1 year in small aliquots in -20 °C. Use a fresh vial each time.

COMMENTARY—Background Information [*Copy Ed: Please inform the authors that the editorial board suggested that it would be valuable to have at least one or sentences somewhere in the commentary about the potential applicability (or not) to an MHC I peptide binding assay. Since the off rates of MHC I are longer, it would not be directly comparable. – RC]

Antigen presentation to CD4⁺ T cells by MHCII molecules plays a central role in the activation and regulation of adaptive immunity. Characterization of peptide binding to MHCII molecules is crucial in antigen presentation and CD4⁺ T cell epitope selection (Germain, 1994; Roche and Cresswell, 1990b). Measurement of peptide binding is also important in developing, training and evaluation of epitope prediction algorithms based on binding affinity (Greenbaum et al., 2011; Hammer et al., 1992; Peters et al., 2005). Therefore, developing a method capable of measuring binding affinity for a large set of

peptides has been one of the key efforts in understanding the molecular mechanism of MHCII antigen presentation, in identifying CD4+ T cell epitopes, and in improving the efficiency of CD4+ T cell epitope prediction algorithms.

The protocol described in this unit allows for a reliable, fast, simple and high throughput measure of peptide binding to MHCII molecules. This FP-based method uses a well-characterized and widely-used influenza hemagglutinin-derived probe peptide (Ferrante et al., 2008; Ferrante and Gorski, 2010; Joshi et al., 2000; Narayan et al., 2007; Roche and Cresswell, 1990a; Stern et al., 1994; Zhou et al., 2009), to study the binding competition by unlabeled target peptides. First introduced in 2006 (De Wall et al., 2006; Nicholson et al., 2006), the FP assay has now become the standard for measuring MHCII-peptide interactions, in studies of T cell epitope selection (Nastke et al., 2012; Yin et al., 2012) and HLA-DM-mediated peptide exchange (Guce et al., 2013; Painter et al., 2011), and in routine use in our laboratory as well as others in the field (Anders et al., 2011; Ferrante and Gorski, 2010; Ferrante and Gorski, 2012; Pos et al., 2012; Zhou et al., 2009).

Other methods have been used previously in MHCII binding experiments to quantify the amount of bound complex, including fluorescence resonance energy transfer (Joshi et al., 2000; Zarutskie et al., 2001), direct fluorescence measurement (Ferrante and Gorski, 2012; Narayan et al., 2007; Rothbard and Busch, 2001), gel filtration using radioactive peptide (Kim et al., 2013; Sidney et al., 2013; Stern and Wiley, 1992), SDS-PAGE gel shift (Ferrante and Gorski, 2012; Kim et al., 2013; Stratikos et al., 2004), ELISA with biotinylated peptide (Tompkins et al., 1993) and surface plasmon resonance (Kim et al., 2013; Narayan et al., 2009). All these methods allow for a direct measurement of binding affinity and estimation of K_d from the association curves. However, the fact that each target peptide has to be labeled for detection in the association reaction and that multiple time points have to be collected to plot the association curves limits the application for measuring binding affinity for a large set of peptides. Compared with these previously used methods, the major advantage of the protocol described in this Unit is that a single labeled probe peptide can be used for binding determination of many different unlabeled peptides and that only single final read is needed to calculate IC_{50} . The protocol we described here uses 96-well microtiter assay plates, and could easily be applied to 384-well or 1536-well format because the basic FP measurement is concentration- and volume-independent.

One potential disadvantage of this protocol is that it measures only relative affinity by comparison with probe peptide, instead of the absolute equilibrium dissociation constant K_d that describes the binding interaction in energetic terms. Under some circumstances, principally simple equilibrium binding reactions where the competing ligands are in large excess over the binding partner, IC_{50} values can be related to K_d values using the Cheng-Prusoff equation $K_d = IC_{50} / (1 + [probe\ peptide] / K_{d,probe\ peptide})$ (Cheng and Prusoff, 1973). In our case, the concentration of Alexa488-FRR-HA₃₀₆₋₃₁₈ is 25 nM, and the K_d is expected to be ~10 nM based on measurements of the unmodified peptide (Roche and Cresswell, 1990a), but the MHCII concentration is not negligible, even considering the fact that some of the protein does not participate in the binding reaction (Joshi et al., 2000). In addition, there is considerable kinetic complexity in the interaction of MHCII proteins with peptides (de Kroon and McConnell, 1994; Joshi et al., 2000; Narayan et al., 2007; Pos et al., 2012),

including steps not likely to reach equilibrium even during the extended 72 hr incubation time. Thus, relating IC_{50} and K_d values is not as straightforward in this system. However, in most cases K_d is expected to scale with the IC_{50} , so that the relative K_d and certainly the K_d rank are reflected in the IC_{50} values. If necessary, the FP approach described here can be adapted to more elaborate assays of the binding reaction, including kinetic studies and HLA-DM dependence (Anders et al., 2011; Yin et al., 2012).

In summary, in this Unit we describe a protocol for fast, simple and high throughput measurement of peptide binding affinity to MHCII molecules, which will facilitate our understanding of the antigen presentation process and improve our ability to screen and predict CD4+ T cell epitopes.

Critical Parameters: Acquisition of the equipment (mainly the microplate reader and LC system) and reagents (mainly purified MHCII, synthesized peptides and fluorochrome) should be considered before the implementation of this procedure. Many microplate readers can be used to measure absorbance, fluorescence, FP and sometimes luminescence and time-resolved fluorescence by changing the detection modes. HPLC systems are widely used in purifying and analyzing biological chemicals by changing the columns. This increased use of microplate reader and LC system in biochemical and biophysical applications has made them a ready tool for immunologists. For measuring FP, it is important to calibrate the instrument before use as described in this protocol and instrument manual. The success of this protocol depends on appropriate measurement of FP. Therefore, it is crucial to optimize the concentrations of MHCII and labeled probe peptide in the competition assay by titration. Because of the usage of fluorescence-labeled material, the entire procedure should be performed with minimum exposure to light.

The MHCII preparation is a critical aspect of this experiment. The reaction depends on labeled probe peptide binding to MHCII protein, either by direct binding to empty molecules or by displacement of weakly associated peptides. Recombinant MHCII proteins produced in insect cells usually are isolated as a mixture of empty molecules and complexes with weakly associated endogenous peptides, depending on the MHCII allele and culture medium (Sloan et al., 1995; Stern and Wiley, 1992; Vollers and Stern, 2008; Zarutskie et al., 1999). By comparison, MHCII proteins isolated from B cells or other native sources where they usually are occupied with a mixture of tight-binding peptides, and less suitable for peptide binding/exchange assays, although in many cases they can be used but at higher concentrations. We have found recombinant soluble DR1 (HLADRB1*01:01), DR2b (HLA-DRB1*15:01), DR3 (HLA-DRB1*03:01), DR4 (HLADRB1*04:01), DR52 (HLA-DRB3*01:01) and DR2a (HLA-DRB5*01:01) suitable for peptide binding experiments (Anders et al., 2011; Hall et al., 2002; Narayan et al., 2007; Yin et al., 2012), and other MHCII proteins expressed in insect cells such as HLA-DQ (Hou et al., 2011; Sidney et al., 2010), also will be suitable. MHCII proteins produced by expression in *E. coli* inclusion bodies and folded in vitro (Frayser et al., 1999) can be made completely devoid of peptide and work very well for *in vitro* binding assays, but so far this technique appears to be restricted to HLA-DR1 (DRB1:01). Due to the residual proteolytic activities in many preparations of MHCII, it is extremely important to use proteinase inhibitor cocktail in the binding buffer to prevent degradation of peptides during incubation.

In the support protocol to label probe peptide, the key step is to optimize the elution gradient to separate the labeled from unlabeled peptides. Also important is to make sure the probe peptide being labeled has only one active amine group and the buffer is free of amine-containing substances.

Troubleshooting: It is essential to perform pilot experiments to test whether FP can be properly measured for both free unbound probe peptide and probe peptide in complex with MHCII. Fluorescence-labeled probe peptide in different percentages of glycol or purified MHCII-probe peptide can help evaluate the FP measurement system and determine the FP values for bound labeled probe peptide.

An unexpected problem that we encountered during development of this assay was the dependence of the bound peptide FP value on the position in the peptide used for labeling. In initial studies we used a modified CLIP peptide labeled at the P5 position (De Wall et al., 2006) or the FRR-HA₃₀₆₋₃₁₈ peptide labeled at the single amino group present lysine at the P3 position (Painter et al., 2011; Yoon et al., 2012). Peptide positions are numbered according to the pocket that they occupy in the bound structure, with P1 corresponding to a large hydrophobic group near the N-terminus (Stern et al., 1994). For most peptides, positions P-1 to P10 make substantial contacts with MHCII residues when bound in the peptide binding groove. Another group used a modified MPB peptide labeled at an introduced cysteine residue at the P5 position (Call et al., 2009; Nicholson et al., 2006) for binding to DR2b (HLA-DRB1*15:01). In extending this approach to other peptides, we labeled a different CLIP variant at its amino terminus (corresponding to the P-5) position, but this peptide showed only a small polarization change when bound to DR1 relative to the value for the free peptide. Moving the Alexa probe to the P-1 position in the same sequence provided the expected polarization shift. FP measures the rotational mobility of the fluorochrome on the labeled probe peptide. If the fluorochrome is free to move, it will give a low FP value and if it is fixed, it will give a high FP value. Apparently, fluorochrome linked to peptide positions outside the binding groove are not sufficiently immobilized relative to unbound peptide to provide a useful FP shift. Thus, in adapting this method to new peptide sequences, it is important to know the binding register of the probe peptide and label it accordingly to give maximum difference in FP values between bound and free states of the peptide, or to screen several positions if the register is not known.

The binding reaction is sensitive to pH and care should be taken to the pH of binding buffer to pH5.5. It should be noted that some preparations of peptide contain residual TFA which can cause a significant and peptide-dose dependent pH shift which in some cases can artificially reduce peptide binding, particularly at very high peptide concentration. The high citrate concentration in the binding buffer is intended to guard against such pH shifts, but the pH of peptide dilutions can be spot checked with pH paper or a pH microelectrode if problems are encountered, and if necessary the peptide dissolved in water and lyophilized to reduce residual TFA.

Maintenance of the plate reader and LC systems is a prerequisite for accuracy and duplication of this protocol. The plate reader must be calibrated regularly for polarization

measurements. The C18 column needs to be washed each time after purification and stored in solution recommended by the manufacturer.

Anticipated Results: In the basic protocol, after 72 hours incubation at 37°C, the wells with free Alexa488-FRR-HA₃₀₆₋₃₁₈ should have a FP value ~35 and the wells with no competing target peptide should have a FP value ~350. The unlabeled target peptide should show a dose-dependent inhibition on the binding of Alexa488-FRR-HA₃₀₆₋₃₁₈ to MHCII. IC₅₀ for unlabeled FRR-HA₃₀₆₋₃₁₈ (or HA₃₀₆₋₃₁₈) should be approximately 25 nM.

In the support protocol, we start with 1 mg of unlabeled FRR-HA₃₀₆₋₃₁₈ and probably obtain 0.23 mg of Alexa488-FRR-HA₃₀₆₋₃₁₈.

The reproducibility of these results in the basic protocol depends on the quality of purified MHCII, labeled probe peptide and unlabeled target peptides. Using the same batches of those materials, uncertainties between replicate experiments performed on different days usually do not exceed 10%-15% of the calculated IC₅₀.

Time Considerations: The preparation of the binding competition assay in the basic protocol takes about 30 minutes and the incubation is 72 hours, followed by 2 minutes of FP reading of the plate. The fluorescence-labeling of the probe peptide in the support protocol takes about 1 hour for labeling, and 3 hours for purification.

Acknowledgments

This work was supported by NIH grants R37AI38996, R01AI48833, and U19AI57319.


Literature Cited

- Anders AK, Call MJ, Schulze MS, Fowler KD, Schubert DA, Seth NP, Sundberg EJ, Wucherpfennig KW. HLA-DM captures partially empty HLA-DR molecules for catalyzed removal of peptide. *Nat Immunol.* 2011; 12:54–61. [PubMed: 21131964]
- Call MJ, Xing X, Cuny GD, Seth NP, Altmann DM, Fugger L, Krogsgaard M, Stein RL, Wucherpfennig KW. In vivo enhancement of peptide display by MHC class II molecules with small molecule catalysts of peptide exchange. *J Immunol.* 2009; 182:6342–6352. [PubMed: 19414787]
- Cheng Y, Prusoff WH. Relationship between the inhibition constant (K₁) and the concentration of inhibitor which causes 50 per cent inhibition (I₅₀) of an enzymatic reaction. *Biochemical pharmacology.* 1973; 22:3099–3108. [PubMed: 4202581]
- de Kroon AI, McConnell HM. Kinetics and specificity of peptide-MHC class II complex displacement reactions. *J Immunol.* 1994; 152:609–619. [PubMed: 8283041]
- De Wall SL, Painter C, Stone JD, Bandaranayake R, Wiley DC, Mitchison TJ, Stern LJ, DeDecker BS. Noble metals strip peptides from class II MHC proteins. *Nature chemical biology.* 2006; 2:197–201.
- Ferrante A, Anderson MW, Klug CS, Gorski J. HLA-DM mediates epitope selection by a “compare-exchange” mechanism when a potential peptide pool is available. *PLoS One.* 2008; 3:e3722. [PubMed: 19005572]
- Ferrante A, Gorski J. Cutting edge: HLA-DM-mediated peptide exchange functions normally on MHC class II-peptide complexes that have been weakened by elimination of a conserved hydrogen bond. *J Immunol.* 2010; 184:1153–1158. [PubMed: 20038641]
- Ferrante A, Gorski J. A Peptide/MHCII conformer generated in the presence of exchange peptide is substrate for HLA-DM editing. *Sci Rep.* 2012; 2:386. [PubMed: 22545194]

- Frayser M, Sato AK, Xu L, Stern LJ. Empty and peptide-loaded class II major histocompatibility complex proteins produced by expression in *Escherichia coli* and folding in vitro. *Protein Expr Purif.* 1999; 15:105–114. [PubMed: 10024477]
- Germain RN. MHC-dependent antigen processing and peptide presentation: providing ligands for T lymphocyte activation. *Cell.* 1994; 76:287–299. [PubMed: 8293464]
- Greenbaum J, Sidney J, Chung J, Brander C, Peters B, Sette A. Functional classification of class II human leukocyte antigen (HLA) molecules reveals seven different supertypes and a surprising degree of repertoire sharing across supertypes. *Immunogenetics.* 2011; 63:325–335. [PubMed: 21305276]
- Guce AI, Mortimer SE, Yoon T, Painter CA, Jiang W, Mellins ED, Stern LJ. HLA-DO acts as a substrate mimic to inhibit HLA-DM by a competitive mechanism. *Nat Struct Mol Biol.* 2013; 20:90–98. [PubMed: 23222639]
- Hall FC, Rabinowitz JD, Busch R, Visconti KC, Belmares M, Patil NS, Cope AP, Patel S, McConnell HM, Mellins ED, Sonderstrup G. Relationship between kinetic stability and immunogenicity of HLA-DR4/peptide complexes. *Eur J Immunol.* 2002; 32:662–670. [PubMed: 11857340]
- Hammer J, Takacs B, Sinigaglia F. Identification of a motif for HLA-DR1 binding peptides using M13 display libraries. *J Exp Med.* 1992; 176:1007–1013. [PubMed: 1402647]
- Hou T, Macmillan H, Chen Z, Keech CL, Jin X, Sidney J, Strohman M, Yoon T, Mellins ED. An insertion mutant in DQA1*0501 restores susceptibility to HLA-DM: implications for disease associations. *J Immunol.* 2011; 187:2442–2452. [PubMed: 21775680]
- Joshi RV, Zarutskie JA, Stern LJ. A three-step kinetic mechanism for peptide binding to MHC class II proteins. *Biochemistry.* 2000; 39:3751–3762. [PubMed: 10736175]
- Kim A, Ishizuka I, Hartman I, Poluektov Y, Narayan K, Sadegh-Nasseri S. Studying MHC class II peptide loading and editing in vitro. *Methods in molecular biology.* 2013; 960:447–459. [PubMed: 23329506]
- Narayan K, Chou CL, Kim A, Hartman IZ, Dalai S, Khoruzhenko S, Sadegh-Nasseri S. HLA-DM targets the hydrogen bond between the histidine at position beta81 and peptide to dissociate HLA-DR-peptide complexes. *Nat Immunol.* 2007; 8:92–100. [PubMed: 17143275]
- Narayan K, Su KW, Chou CL, Khoruzhenko S, Sadegh-Nasseri S. HLA-DM mediates peptide exchange by interacting transiently and repeatedly with HLA-DR1. *Mol Immunol.* 2009; 46:3157–3162. [PubMed: 19647320]
- Nastke MD, Becerra A, Yin L, Dominguez-Amorocho O, Gibson L, Stern LJ, Calvo-Calle JM. Human CD4+ T cell response to human herpesvirus 6. *J Virol.* 2012; 86:4776–4792. [PubMed: 22357271]
- Nicholson MJ, Moradi B, Seth NP, Xing X, Cuny GD, Stein RL, Wucherpfennig KW. Small molecules that enhance the catalytic efficiency of HLA-DM. *J Immunol.* 2006; 176:4208–4220. [PubMed: 16547258]
- Painter CA, Negroni MP, Kellersberger KA, Zavala-Ruiz Z, Evans JE, Stern LJ. Conformational lability in the class II MHC 310 helix and adjacent extended strand dictate HLA-DM susceptibility and peptide exchange. *Proc Natl Acad Sci U S A.* 2011; 108:19329–19334. [PubMed: 22084083]
- Peters B, Sidney J, Bourne P, Bui HH, Buus S, Doh G, Fleri W, Kronenberg M, Kubo R, Lund O, Nemazee D, Ponomarenko JV, Sathiamurthy M, Schoenberger S, Stewart S, Surko P, Way S, Wilson S, Sette A. The immune epitope database and analysis resource: from vision to blueprint. *PLoS Biol.* 2005; 3:e91. [PubMed: 15760272]
- Pos W, Sethi DK, Call MJ, Schulze MS, Anders AK, Pyrdol J, Wucherpfennig KW. Crystal structure of the HLA-DM-HLA-DR1 complex defines mechanisms for rapid peptide selection. *Cell.* 2012; 151:1557–1568. [PubMed: 23260142]
- Roche PA, Cresswell P. High-affinity binding of an influenza hemagglutinin-derived peptide to purified HLA-DR. *J Immunol.* 1990a; 144:1849–1856. [PubMed: 2307844]
- Roche PA, Cresswell P. Invariant chain association with HLA-DR molecules inhibits immunogenic peptide binding. *Nature.* 1990b; 345:615–618. [PubMed: 2190094]
- Rothbard JB, Busch R. Binding of biotinylated peptides to MHC class II proteins on cell surfaces. *Current protocols in immunology / edited by John E. Coligan ... [et al.].* 2001 Chapter 18:Unit 18 11.

- Sidney J, Southwood S, Moore C, Oseroff C, Pinilla C, Grey HM, Sette A. Measurement of MHC/peptide interactions by gel filtration or monoclonal antibody capture. *Current protocols in immunology* / edited by John E. Coligan ... [et al.]. 2013 Chapter 18:Unit 18 13.
- Sidney J, Steen A, Moore C, Ngo S, Chung J, Peters B, Sette A. Divergent motifs but overlapping binding repertoires of six HLA-DQ molecules frequently expressed in the worldwide human population. *J Immunol*. 2010; 185:4189–4198. [PubMed: 20810981]
- Sloan VS, Cameron P, Porter G, Gammon M, Amaya M, Mellins E, Zaller DM. Mediation by HLA-DM of dissociation of peptides from HLA-DR. *Nature*. 1995; 375:802–806. [PubMed: 7596415]
- Stern LJ, Brown JH, Jardetzky TS, Gorga JC, Urban RG, Strominger JL, Wiley DC. Crystal structure of the human class II MHC protein HLA-DR1 complexed with an influenza virus peptide. *Nature*. 1994; 368:215–221. [PubMed: 8145819]
- Stern LJ, Wiley DC. The human class II MHC protein HLA-DR1 assembles as empty alpha beta heterodimers in the absence of antigenic peptide. *Cell*. 1992; 68:465–477. [PubMed: 1371238]
- Stratikos E, Wiley DC, Stern LJ. Enhanced catalytic action of HLA-DM on the exchange of peptides lacking backbone hydrogen bonds between their N-terminal region and the MHC class II alpha-chain. *J Immunol*. 2004; 172:1109–1117. [PubMed: 14707085]
- Tompkins SM, Rota PA, Moore JC, Jensen PE. A europium fluoroimmunoassay for measuring binding of antigen to class II MHC glycoproteins. *J Immunol Methods*. 1993; 163:209–216. [PubMed: 8354890]
- Vollers SS, Stern LJ. Class II major histocompatibility complex tetramer staining: progress, problems, and prospects. *Immunology*. 2008; 123:305–313. [PubMed: 18251991]
- Yin L, Calvo-Calle JM, Dominguez-Amoroch O, Stern LJ. HLA-DM constrains epitope selection in the human CD4 T cell response to vaccinia virus by favoring the presentation of peptides with longer HLA-DM-mediated half-lives. *J Immunol*. 2012; 189:3983–3994. [PubMed: 22966084]
- Yoon T, Macmillan H, Mortimer SE, Jiang W, Rinderknecht CH, Stern LJ, Mellins ED. Mapping the HLA-DO/HLA-DM complex by FRET and mutagenesis. *Proc Natl Acad Sci U S A*. 2012; 109:11276–11281. [PubMed: 22733780]
- Zarutskie JA, Busch R, Zavala-Ruiz Z, Rushe M, Mellins ED, Stern LJ. The kinetic basis of peptide exchange catalysis by HLA-DM. *Proc Natl Acad Sci U S A*. 2001; 98:12450–12455. [PubMed: 11606721]
- Zarutskie JA, Sato AK, Rushe MM, Chan IC, Lomakin A, Benedek GB, Stern LJ. A conformational change in the human major histocompatibility complex protein HLA-DR1 induced by peptide binding. *Biochemistry*. 1999; 38:5878–5887. [PubMed: 10231540]
- Zhou Z, Callaway KA, Weber DA, Jensen PE. Cutting edge: HLA-DM functions through a mechanism that does not require specific conserved hydrogen bonds in class II MHC-peptide complexes. *J Immunol*. 2009; 183:4187–4191. [PubMed: 19767569]

A Set up binding competition assay



B Read fluorescence polarization (FP) and calculate Relative binding

$$\text{Relative binding} = \frac{\text{FP}_{\text{sample}} - \text{FP}_{\text{free}}}{\text{FP}_{\text{no_comp}} - \text{FP}_{\text{free}}}$$

where FP_{sample} is the FP value for the sample well,
 FP_{no_comp} is the FP value for Alexa488-FRR-HA₃₀₆₋₃₁₈ bound with MHCII without competitor peptide,
 and FP_{free} is the FP value for free Alexa488-FRR-HA₃₀₆₋₃₁₈.

C Data analysis and calculation of IC₅₀




Figure 1. A schematic overview of the steps involved in measuring peptide binding to MHCII molecules by fluorescence polarization

(A) Set up the binding competition assay as shown in Table I and incubate at 37 °C for 72 hours. (B) Read FP on a 96-well microplate reader, which has been calibrated, and convert the FP values into Relative binding using the equation indicated. (C) Fit the competition curve to obtain IC₅₀. Peptide 1 represents a high affinity peptide, as IC₅₀ is 139nM; peptide 2 has no binding affinity as no competition is observed; peptide 3 has an intermediate affinity with an IC₅₀ of 1709nM.




Figure 2. A schematic overview of the steps involved in labeling the probe peptide with a fluorochrome

(A) Set up the labeling reaction mix and incubate at room temperature for 1 hour. (B) Representative elution profile of fluorochrome Alexa488-TFP (first panel, peak (a)), unlabeled FRR-HA₃₀₆₋₃₁₈ (second panel, peak (b)), labeling reaction mix (third panel), and Alexa488-labeled FRR-HA₃₀₆₋₃₁₈ (fourth panel, peak (c)). Peak (a) in the second panel comes from the residual fluorochrome Alexa488-TFP in the last run. The fifth panel shows the Buffer B gradient.

Table 1

Typical set up of the binding competition assay in a 96-well Microplate for measuring peptide binding to MHCII molecules

Target peptide (nM) ^a	Target peptide 1 + MHCII + Alexa488-FRR-HA ₃₀₆₋₃₁₈			Target peptide 2 + MHCII + Alexa488-FRR-HA ₃₀₆₋₃₁₈			Target peptide 3 + MHCII + Alexa488-FRR-HA ₃₀₆₋₃₁₈			Target peptide 4 + MHCII + Alexa488-FRR-HA ₃₀₆₋₃₁₈		
	1	2	3	4	5	6	7	8	9	10	11	12
20000												
4000												
800												
160												
32												
6.4												
1.28												
0												
	FP_no_comp ^b						FP_free ^c					

^aTarget peptide concentrations starts at 20000 nM with a dilution factor of 5. Assay for each peptide is performed with three replicates. The last row H has no competing target peptide.

^bRow H, columns 1-6 contain only MHCII and labeled probe peptide, no competing target peptide, allowing for measurement of FP for labeled probe peptide bound with MHCII without competitor peptide (FP_no_comp).

^cRow H, columns 7-12 contain only labeled probe peptide, no MHCII and competing target peptide, allowing for measurement of FP for free unbound labeled probe peptide (FP_free).

Table 2Elution gradient for purification of Alexa488-FRR-HA₃₀₆₋₃₁₈

Time (minute)	% Buffer B ^a	Flow Rate (ml/minute)
0	0.0	1.0
5 ^b	2.0	1.0
25 ^c	98.0	1.0
28 ^d	98.0	1.0
30 ^e	2.0	1.0
35 ^f	2.0	1.0

^a Buffer B is 5% TFA in acetonitrile, and the counterpart Buffer A is 5% TFA in ddH₂O.

^b Use 5 minutes to equilibrium the C18 column in 2% Buffer B.

^c Run the gradient of Buffer B from 2% to 98% for 20 minutes, which will elute Alexa488-FRR-HA₃₀₆₋₃₁₈ from the C18 column.

^d Keep the gradient of Buffer B at 98% for 3 minutes to wash off any remaining substances from the column.

^e Equilibrium the C18 column back into 2% Buffer B.

^f Stop the current run at 35 minutes and keep the C18 column in 2% Buffer B for the next run.