



Published in final edited form as:

Proc Int Soc Magn Reson Med Sci Meet Exhib Int Soc Magn Reson Med Sci Meet Exhib. 2013 ; 21: 0473.

Accelerated, motion-corrected high-resolution intravascular MRI at 3T

Shashank Sathyanarayana Hegde¹, Yi Zhang², and Paul A Bottomley¹

¹Radiology, Johns Hopkins University, Baltimore, Maryland, United States

²Electrical and Computer Engineering, Johns Hopkins University, Baltimore, Maryland, United States

Audience

Interventionalists, and those interested in atherosclerosis and intravascular MR imaging.

Purpose

Current speeds for intravascular (IV) MRI and MRI endoscopy¹ are limited to ~ 2 frames/s at 3T, rendering high-resolution ($\sim 100\mu\text{m}$) images susceptible to degradation by physiological motion on the order of mm/ms. Here, using projection reconstruction we: (A) reduce sensitivity to motion from the time-scale of individual images, to the time-frame of each projection (TR) by frame-shifting each projection to the antenna, prior to reconstruction. In addition: (B) we apply compressed sensing to provide acceleration factors of up to four-fold. We present data acquired in phantoms (fruit), human vessel specimens and/or apply the methods to retro-actively acquired data as we move toward prospective acquisitions *in vivo*.

Methods

IV MRI with and without mechanical motion, is performed on a Philips 3T scanner using a 2mm diameter 3T loopless antenna receiver, and radial k-space traversal. For motion correction (A), we note that in each projection, there is intense signal surrounding the probe, but the probe itself produces no signal. Further, there is a phase reversal that occurs at the probe (Fig. 1 a, d). These amplitude and phase singularities at the probe's location are detected using a signal derivative algorithm, and used to align all the projections (Fig. 1f). Images reconstructed from these, always have the probe at the center of the field-of-view. Compressed sensing (B), is performed on projection images using uniform under-sampling², while variable-density random under-sampling is used on previously-acquired *in vivo* Cartesian data¹. Images are reconstructed using " ℓ_1 -norm" minimization and wavelet transform^{2,3}.

Results

Motion correction significantly reduces motion artefact compared to conventional reconstruction (Fig. 1b vs. 1c). Radial and Cartesian compressed sensing produced virtually indistinguishable images with only $1/4^{\text{th}}$ to $1/3^{\text{rd}}$ of the original data (Fig. 2, 3). Since the

motion correction algorithm acts on each projection, it was also applied to a radially under-sampled data set (not shown).

Conclusions

3T IV MRI detectors are ideally suited to compressed sensing and motion correction strategies based on their intrinsically radial and sparsely-localized sensitivity profiles and high signal-to-noise ratios. The benefits are much faster IV MRI—approaching real-time (~ 10 fr/s) and reduced motion sensitivity, while retaining the high-resolution (80-300 μ m) image information.

Acknowledgments

Support: NIH R01 HL090728.

References

1. Sathyanarayana S, et al. JACC Cardiovasc Imaging. 2010; 3(11):1158–1165. [PubMed: 21071004]
2. Block, et al. Magn Reson Med. 2007; 57:1086–1098. [PubMed: 17534903]
3. Lustig, et al. Magn Reson Med. 2007; 58(6):1182–1195. [PubMed: 17969013]




Figure 1.

(a) Transverse field of a loopless antenna detector p shows decreasing B_1 with r and azimuthal variation in phase. (b) MRI of an orange shaken ± 3 mm (2D radial GRE; 200 spokes spanning 180° ; $250\mu\text{m}$ in-plane resolution; $\text{TR}/\text{TE}=15/6$ ms) shows debilitating motion artifacts. (c) Projection shifting all but removes streaking, revealing the fruit's underlying structure. A $1/r$ intensity filter has been applied to aid visualization. (d-f) The motion correction algorithm consists of re-aligning every azimuthal projection on p .




Figure 2. Fruit morphology using the complete data set (**a**), is retained in a four-fold under-sampled radial-compressed sense reconstruction (**b**).




Figure 3. (a) Regular Cartesian MRI endoscopy of a rabbit aorta *in vivo* (3D GRE; TR/TE=250/12 ms; in-plane resolution 80 μ m; 3.1 min/5 contiguous slices). (b) Three-fold under-sampling yields a virtually indistinguishable image (cropped for visualization) after compressed-sense reconstruction.