



Published in final edited form as:

*Proc Int Soc Magn Reson Med Sci Meet Exhib Int Soc Magn Reson Med Sci Meet Exhib.* 2013 ; 21: 0530.

## **Spectroscopy with linear algebraic modeling (SLAM): speed and quantification in brain tumor studies**

Yi Zhang<sup>1,2</sup>, Refaat E Gabr<sup>1</sup>, Jinyuan Zhou<sup>1,3</sup>, Robert G Weiss<sup>1,4</sup>, and Paul A Bottomley<sup>1,2</sup>

<sup>1</sup>Division of MR Research, Department of Radiology, Johns Hopkins University, Baltimore, Maryland, United States

<sup>2</sup>Electrical and Computer Engineering, Johns Hopkins University, Baltimore, Maryland, United States

<sup>3</sup>F. M. Kirby Research Center for Functional Brain Imaging, Kennedy Krieger Institute, Baltimore, Maryland, United States

<sup>4</sup>Division of Cardiology, Department of Medicine, Johns Hopkins University, Baltimore, Maryland, United States

---

### **Target audience**

Scientists and clinicians studying large lesions or pathologies with magnetic resonance spectroscopy.

### **Purpose**

Potentially useful clinical MRS chemical shift imaging (CSI) acquisitions are often precluded by long study times. The newly-proposed Spectroscopy with Linear Algebraic Modeling (SLAM) method<sup>1</sup> can provide compartmental-average spectra many-fold faster than CSI and at little cost to signal-to-noise ratio, by incorporating information from scout MRI<sup>1,3</sup>. Prior work demonstrated quantitative agreement in <sup>31</sup>P 1D CSI heart studies<sup>1</sup>. Here, SLAM with a 6 to 14-fold speed-up is extended to 2D <sup>1</sup>H MRS studies of patients with brain tumors. Quantitative results are compared with CSI.

### **Methods**

Proactive SLAM is implemented as follows: (i) Acquire MRI; (ii) Apply a fraction of the central k-space CSI phase encodes based on the number of compartments and speed-up factor; (iii) Co-register the CSI grid with the MRI; (iv) Segment the CSI voxels into the desired MRI-based compartments; (v) Reconstruct compartmental spectra using SLAM. SLAM is tested retroactively on CSI data, by using only the central k-space phase-encoded signals in Step (ii).

SLAM reconstruction was implemented on <sup>1</sup>H CSI data from 15 brain tumor patients undergoing clinical/research MRI and MRS in a 3T Philips scanner, using 1/6 of the acquired (central k-space) CSI data. Five compartments—tumor, contralateral, brain (excluding tumor and contralateral), scalp and background—were segmented for SLAM

reconstruction. Cho/Cr ratios were quantified in compartment spectra and compared with the compartmental average CSI spectra.

The proactive SLAM protocol was applied to volunteers after CSI protocol, using a SLAM speedup factor of ~14. Four compartments—a user-defined “lesion”, brain, scalp and background—were assigned and results compared to CSI.

## Results

Fig. 1 shows segmentation in a brain tumor study (a), and the corresponding spectra (b-d) reconstructed from CSI (blue) and SLAM (red). Fig. 2 shows that SLAM Cho/Cr in tumor (a) and contralateral (b) compartments agree with CSI (vertical) in 15 patients. Fig. 3 shows proactive implementation of SLAM on a normal volunteer.

## Conclusion


SLAM was successfully applied retroactively to patients with brain tumors. It is also implemented proactively for brain studies. SLAM yielded spectra and quantitative results almost identical to CSI, but much faster ( 14 times). SLAM could enable whole-brain MRS in clinical studies that were otherwise precluded by long MRS acquisition times.

## Acknowledgments


Grant support: EB007829, HL61912, EB009731

## References


- [1]. Zhang Y, et al. JMR. 2012; 218:66–76. [PubMed: 22578557]
- [2]. Brown TR, et al. PNAS. 1982; 79:3523–3526. [PubMed: 6954498]
- [3]. Zhang Y, et al. Proc. ISMRM. 2012:0709.



**Fig. 1.** Retroactive SLAM using central  $12 \times 13$  phase encodes out of original  $27 \times 33$  CSI phase encodes. (a) Brain segmentation into 5 compartments. SLAM (red) spectra overlaid on CSI (blue) spectra for the tumor compartment #1 (b), the contralateral compartment #2 (c) and brain compartment #3 (d).



**Fig. 2.** SLAM vs. CSI choline/creatine (Cho/Cr) ratios in 15 brain tumor patients. SLAM used 1/6 of CSI data.



**Fig. 3.** (a) Brain segmentation from normal volunteer with 4 compartments. (b) Proactive 0.8 min SLAM spectrum (red) from compartment #1 compared to 12 min CSI acquisition (blue). (c) 0.8 min SLAM and 12 min CSI spectra from brain compartment #2 (excluding compartment #1).