ORIGINAL ARTICLE



Development of menaquinone-7 enriched nutraceutical: inside into medium engineering and process modeling

Rishipal Singh · Alka Puri · Bibhu Prasad Panda

Revised: 3 August 2014 / Accepted: 2 October 2014 / Published online: 19 October 2014 © Association of Food Scientists & Technologists (India) 2014

Abstract Menaguinone 7 (MK-7) is nutritionally important metabolite found by fermentation mainly using *B. subtilis* species. In this study, soybean medium was modified to improve the MK-7 production using Bacillus subtilis NCIM 2708 under solid state fermentation. The objective of this study was to produce large amount of MK-7 within a short period of time. Nine nutritional components viz. glycerol, mannitol, dextrose, sucrose, yeast extract, malt extract, K₂HPO₄, MgSO₄.7H₂O and CaCl₂ were investigated to obtain the maximum MK-7 concentration. The highest MK-7 concentration 39.039 µg/g was obtained after 24 h of fermentation in the following optimised medium components: soybean 20 g, glycerol 40 ml/kg, mannitol 60 g/kg, yeast extract 4 g/kg, malt extract 8 g/kg and calcium chloride 4 g/kg. The maximum production of MK-7 56.757 µg/g was predicted by point prediction tool of Design Expert 7.1 software (Statease Inc. USA). This data shows 68.78 % validity of the predicted model.

Keywords Soybean · Solid state fermentation · Menaquinone-7 · Medium engineering

Introduction

Vitamin K is fat soluble and occurs naturally in two forms vitamin K1 and vitamin K2. Vitamin K1 is major dietary source of vitamin K, exist naturally in green plants, algae, and photosynthetic bacteria while vitamin K2 includes chicken egg yolk,

Rishipal Singh, Alka Puri and Bibhu Prasad Panda have contributed equally to this paper

R. Singh · A. Puri · B. P. Panda (🖂)

Microbial and Pharmaceutical Biotechnology Laboratory, Centre for Advanced Research in Pharmaceutical Science, Faculty of Pharmacy, Jamia Hamdard, New Delhi 110062, India e-mail: bibhu panda31@rediffmail.com butter, cow liver, certain cheeses and in fermented soybeans (Vermeer and Schurgers 2000; Fujita et al. 2011). Vitamin K3 (menadione) is a synthetic compound that can be converted to K1 in the intestinal tract. It is only used in animal nutrition. Vitamin K2 family has a methylated naphthoquinones ring structure with a variable side chain length of 4 to 13 isoprene units. The compounds in this series are referred to as menaquinone-n, where n denotes the number of isoprene units (Shearer 1990). The level of menaquinone-7 (MK-7) in different food products is very low except for fermented soybean foods such as natto (9.39 µg/g), Hikiwari natto (chopped natto, 8.27 μ g/g) and black soybean *natto* (796 μ g/g) (Kamao et al. 2007). Recently, research has demonstrated that MK-7 consumption significantly reduces the risk of bone fractures (Schurgers et al. 2007; Truong and Booth 2011) and cardiovascular disorders (Gast et al. 2009; Theuwissen et al. 2012). European experts suggested that Recommended Daily Intake of vitamin K, preferably in the form of vitamin K₂ is 200-500 µg/day, which is required for optimal carboxylation of extrahepatic GLA proteins. Natto, a Japanese fermented food, has traditionally been obtained from cooked soybean fermented with Bacillus species under SSF (Sakano et al. 1988). In Japan, natto has received approval as a nutritional food containing high amount of MK-7. Natto is required in a large quantity to meet the therapeutic dose of Vitamin K2 for treatment of osteoporosis and cardiovascular disease. Hence, it is advantageous to develop vitamin K rich *natto* as a dietary supplement.

For fermentation medium engineering, the screening of important parameters influencing fermentation productivity is initially carried out by barrowing method or by Plackett-Burman design and then, RSM is employed to develop, improve and optimise the process using the screened parameters (Kennedy and Krouse 1999). Berenjian et al. (2011a, b) reported that nitrogen source (yeast extract, soy peptone mixture) and carbon source (glycerol) were the most effective nutrients for enhancing MK-7 production. The presence of mono, di and polysaccharides and glycerol in the fermented medium can be used in glycolysis process by *B. subtilis* strains. Menaquinone production is the very complex systematic process. The pathway for isoprene side chain and quinone skeleton (1, 4-naphthoquinone) production is dependent on the presence of carbon sources such as glucose, fructose and glycerol in the fermentation media (Sonenshein et al. 2002). Nutritional requirements of *Bacillus* species play an important role for synthesizing their metabolic products; therefore fermentation medium designing is an important process for production of any secondary metabolite. MK-7 is produced mainly by submerged fermentation using *Bacillus subtilis* species (Sato et al. 2001b; Beulens et al. 2009).

Researchers have reported that genetic manipulation of microbial strains can improve MK-7 productivity during fermentation process (Sato et al. 2001a, b; Tsukamoto et al. 2001). Several researchers have studied the mechanism of MK biosynthesis in Bacillus subtilis (Rowland et al. 1995; Rowland and Taber 1996; Seto et al. 2008). Other studies focused on increasing the MK-7 productivity by modifying the fermentation medium and process parameters. Nevertheless, there are only few examples in literature that involve response surface methodology (RSM), a strategic experimental tool, for enhancing the production of menaquinone-7 (Mahanama et al. 2011; Wu and Ahn 2011; Mahanama et al. 2012a; Berenjian et al. 2011b). Mahanama et al. (2012b) have reported the effect of enzymatic pretreatment on MK-7 production in SSF. Fermentation media and operating conditions have a significant impact on increasing MK-7 production. No research till now conducted for optimization of fermentation medium and followed by process optimization for production of MK-7 rich soya nutraceutical.

Thus, the research was carried out to investigate the effect of different nutritional components and to optimise the significant medium parameters and further optimization of the solid state fermentation condition (process parameters) for maximum MK-7 production in soybeans.

Materials and methods

Materials

Chemicals and microbial strain

All the chemicals used were purchased from Hi-media Laboratories Limited (Mumbai, India) and were of analytical grade. Soybean variety (PUSA-9814) was used as substrate for MK-7 production under solid state fermentation. The soybean variety PUSA-9814 was collected from the Pulse Lab of Indian Agricultural Research Institute (IARI), PUSA, New Delhi. Culture of *Bacillus subtilis* NCIM 2708 was obtained from National Collection of Industrial Microorganisms, National Chemical Laboratory, Pune, Maharashtra, India. The culture was maintained on nutrient agar slants at 4 °C and sub cultured at every 30 days.

Methods

Preparation of seed culture

Spore suspension of *Bacillus subtilis* NCIM 2708 was prepared from actively growing slants using sterile water. Spore suspension was inoculated to sterilized medium: soybean powder 6 % (soybean variety, SL- 525), sodium chloride 0.5 % and distilled water adjusted to pH 7.0 with 0.1 N HCl or NaOH (Hu et al. 2000). These cultures were incubated at 37 °C for 16 h in a shaker incubator at 180 rpm.

Table 1Plackett-Burman De-
sign of twelve experimental trials
generated for nine different me-
dium and two dummy compo-
nents by software Design Expert7.1.6 and Menaquinone-7 (MK-
7) production

		X_3	X_4	X_5	X_6	X_7	X_8	X9	X_{10}	X ₁₁	MK-7 (µg/g)
+	+	_	+	+	+	_	_	_	+	-	26.4
-	+	+	-	+	+	+	-	-	-	+	7.08
+	-	+	+	-	+	+	+	-	-	-	13.14
-	+	-	+	+	-	+	+	+	-	-	2.88
_	-	+	-	+	+	-	+	+	+	-	2.91
-	-	_	+	-	+	+	-	+	+	+	2.82
+	-	_	-	+	-	+	+	-	+	+	2.85
+	+	-	-	-	+	-	+	+	-	+	2.94
+	+	+	-	-	-	+	-	+	+	-	2.85
-	+	+	+	-	-	-	+	-	+	+	2.97
+	-	+	+	+	-	-	-	+	-	+	2.85
-	-	_	-	-	-	-	-	-	-	-	2.88
	 + + + +	- + + - - + + - + + + + + + + +	$\begin{array}{cccccccccccccccccccccccccccccccccccc$								

Designation	Variable	ΣH (+1)	ΣL (-1)	Mean Square	Main Effect	F-value	% Contribution
X ₁	Glycerol	51.03	21.54	72.47	4.91	10.73	13.09
X_2	Mannitol	45.12	27.45	25.49	2.94	3.77	4.59
X ₃	Sucrose	31.8	40.77	71.88	-1.49	10.64	12.94
X_4	Dextrose	51.06	21.51	120.78	4.92	17.89	21.73
X5	Yeast Extract	44.97	27.6	73.66	2.56	10.91	13.25
X ₆	Malt Extract	55.29	17.28	120.02	6.33	17.78	21.6
X ₇	Dummy 1	31.62	40.95	6.44	-1.55	0.95	1.16
X ₈	CaCl ₂	27.69	44.88	7.07	2.86	1.04	1.27
X ₉	Dummy 2	17.25	55.32	7.16	6.34	1.06	1.29
X ₁₀	K ₂ HPO ₄	40.8	31.77	25.14	1.5	3.72	4.52
X ₁₁	MgSO ₄ .7H ₂ O	21.51	51.06	25.32	4.9	3.75	4.56

Table 2 Contribution of different carbon, nitrogen and ions on MK-7 production in experimental trials carried out by Plackett-Burman Design

Solid -state fermentation of soybean

Soybean seeds (20 g) were washed, soaked in water for 12 h, boiled and dehulled. This soybean-based medium was used directly for production of MK-7 after adjusting the pH to 7. Sterilized medium was inoculated with seed culture (1 ml/g of substrate) and incubated for 24 h at 37 °C. After 24 h, fermented soybeans were removed from incubator and kept at 4 °C. The samples were processed for extraction and analysed by HPLC (Shimadzu, Japan).

Optimization of medium

The optimization of medium components for maximization of MK-7 production by *B. subtilis* under solid state fermentation was evaluated in two stages. The first stage, screening of medium component by Plackett-Burman experimental design was done. The second stage was optimization of key nutrients concentration influencing MK-7 production under solid state fermentation by central composite design of response surface methodology (Kennedy and Krouse 1999). Percentage contribution was calculated with the help of design expert 7.1 software (Statease Inc. USA).

Screening of medium components influencing menaquinone-7 production

The Plackett-Burman experimental design (PBD) for eleven variables with nine nutritional components (X_1 to X_9) (independent variables) and two dummy variables (D_1 and D_2) were used to evaluate the relative importance of the nutrients to enhance the production of MK-7 under SSF condition. These variable represented at two-levels, high concentrations (+) and low concentrations (-) and two dummy variables in eleven trials as shown in Table 1. The incubation period, temperature, medium pH, Inoculum volume was maintained

constant throughout the fermentation process (Kennedy and Krouse 1999).

Optimization of medium components

The optimum levels of screened medium components (Table 2) affecting MK-7 production were determined by the Central Composite Design (CCD) of RSM. The five independent variables contributing maximum to MK-7 production such as glycerol, mannitol, malt extract, yeast extract and CaCl₂ were studied at five different levels $(-1, -\alpha, 0, +1, +\alpha)$ with the 50 treatment combinations. The Experimental range and levels of each independent variable is presented in Table 3. Data collected for MK-7 concentrations in each run were analysed using the design expert 7.1.6 software.

Optimization of process parameters

Box-Behnken Design of RSM was used to determine the optimum level of key process variables to obtain the maximum MK-7 concentration. In this study, three independent variables effecting MK-7 production such as fermentation medium pH, fermentation temperature (°C)

 Table 3
 Experimental range and levels of the five significant independent variables used in RSM for optimization of fermentation medium parameters

Variables (g/kg)	Levels	Levels							
	-α	-1	0	+1	+α				
Glycerol	40	30	50	60	70				
Mannitol	40	30	50	60	70				
Malt Extract	4	2	6	8	10				
Yeast Extract	4	2	6	8	10				
CaCl ₂	4	2	6	8	10				

 Table 4
 The predicted and experimental MK-7 concentrations of fermentation trials carried out by CCD of response surface methodology for optimization of fermentation medium parameters

140404448.37 6.78 260404443.03 5.47 34060444 2.94 7.85 54040844 2.94 7.85 54040844 2.94 7.85 54040844 2.94 7.85 66040844 3.48 6.24 86060844 3.48 6.24 86060844 9.97 6.96 74060844 9.87 9.13 114060484 23.55 18.92 126060484 3.66 5.11 146040884 3.66 5.11 146040884 7.77 9.82 154060884 7.77 9.82 1640448 8.01 7.64 194060448 8.01 7.64 194060848 9.78 6.29 214040848 8.62 2.92 214040848 8.62 226040 <th>Run</th> <th>Glycerol (ml/kg)</th> <th>Mannitol (g/kg)</th> <th>Yeast Extract (g/kg)</th> <th>Malt Extract (g/kg)</th> <th>CaCl₂ (g/kg)</th> <th>Actual MK-7 (µg/g)</th> <th>Predicted MK-7 (µg/g)</th>	Run	Glycerol (ml/kg)	Mannitol (g/kg)	Yeast Extract (g/kg)	Malt Extract (g/kg)	CaCl ₂ (g/kg)	Actual MK-7 (µg/g)	Predicted MK-7 (µg/g)
2 60 40 4 4 4 3.03 5.47 3 40 60 4 4 4 2.94 7.85 5 40 40 8 4 4 2.94 7.85 5 40 40 8 4 4 2.94 7.85 7 40 60 8 4 4 9.9 6.96 7 40 60 8 4 4 3.48 6.24 8 60 60 8 4 4 9.87 9.13 10 60 40 4 8 4 9.87 9.13 11 40 60 4 8 4 23.55 18.92 12 60 60 4 8 4 3.66 5.11 14 60 40 8 8 4 3.84 6.23 15 40 60 8 8 4 7.77 9.82 17 40 40 4 8 8.01 7.64 19 40 60 4 4 8 9.03 8.63 20 60 60 4 4 8 9.29 1.25 21 40 40 8 4 8.79 10.25 22 60 40 8 4 8.64 10.33 25 40 40 8 8 9.60 13.58 28 60 60 8 8 $9.$	1	40	40	4	4	4	8.37	6.78
3 40 60 4 4 4 2.94 7.85 5 40 40 8 4 4 2.95 4.19 6 60 40 8 4 4 2.85 4.19 6 60 40 8 4 4 9.90 6.96 7 40 60 8 4 4 10.8 7.56 9 40 40 4 8 4 9.36 12.10 10 60 40 4 8 4 9.35 18.92 12 60 60 4 8 4 2.3.55 18.92 13 40 40 8 8 4 3.66 5.11 14 60 60 8 8 4 9.00 10.16 16 60 60 8 8 4 9.00 8.63 20 60 60 4 4 8 9.03 8.63 21 40 60				4	4			
460604442.947.85540408442.854.19660408449.96.967406084410.87.56940404849.3612.101060404849.879.1311406048423.5518.921260604843.665.111340408843.665.111460408843.645.331540608849.0010.161660608849.786.291740404488.498.531860404489.398.582060604489.398.582140408489.398.582260408489.398.582340608489.4110.33254040889.2410.4826604889.099.58294040889.099.5829								
5 40 40 8 4 4 2.85 4.19 6 60 8 4 4 9.9 6.96 7 40 60 8 4 4 3.48 6.24 8 60 60 8 4 4 10.8 7.56 9 40 40 4 8 4 9.36 12.10 10 60 40 4 8 4 9.35 18.92 12 60 60 4 8 4 23.55 18.92 13 40 40 8 8 4 3.66 5.11 14 60 40 8 8 4 3.66 5.11 14 60 40 8 8 4 7.77 9.82 17 40 40 4 4 8 8.49 8.53 18 60 40 4 4 8 8.90 8.63 20 60 60 4 4 8 9.78 6.29 21 40 40 8 4 8.79 10.25 22 60 40 8 4 $8.9.79$ 10.25 22 60 40 8 4 $8.9.46$ 7.94 23 40 60 8 4 8 9.60 23 40 60 8 4 8 9.60 24 60 60 8 4 8 9.60 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
6 60 40 8 4 4 9.9 6.96 7 40 60 8 4 4 3.48 6.24 8 60 60 8 4 4 10.8 7.56 9 40 40 4 8 4 9.36 12.10 10 60 40 4 8 4 9.35 18.92 12 60 60 4 8 4 23.55 18.92 13 40 40 8 8 4 3.66 5.11 14 60 40 8 8 4 3.84 6.23 15 40 60 8 8 4 7.77 9.82 17 40 40 4 4 8 8.49 8.53 18 60 40 4 4 8 8.63 20 60 60 4 4 8 9.78 6.29 21 40 40 8 4 8.79 10.25 22 60 40 8 4 8 9.79 8.58 24 60 60 8 4 8 9.60 13.58 25 40 40 4 8 9.60 13.58 26 60 40 8 8 9.09 9.58 28 60 60 4 8 $8.9.02$ 9.13 35 50 6 6 6 9.27 <								
74060844 3.48 6.24 86060844 10.8 7.56 94040484 9.36 12.10 106040484 9.35 18.92 126060484 23.55 18.92 126060484 3.66 5.11 146040884 3.66 5.11 146040884 9.00 10.16 166060884 7.77 9.82 174040448 8.49 8.53 186040448 8.01 7.64 194060448 9.03 8.63 206060448 9.78 6.29 214040848 9.78 6.29 214040848 9.78 6.29 226040848 9.78 6.29 234060848 9.79 8.58 246060848 9.24 10.48 254040488 9.60 13.58 286060488 9.02 9.13 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
8606084410.87.56940404849.3612.101060404849.879.1311406048423.5518.9212606048415.3614.501340408843.665.111460408843.846.231540608847.779.821740404488.498.531860404488.017.641940604489.038.632060604489.786.29214040848.9786.29214040848.9398.58226040848.9398.582460608489.992540404889.602660404889.602740604889.092850506669.213330506669.213330506669.21<								
940404849.3612.101060404849.879.1311406048423.5518.9212606048415.3614.501340408843.665.111460408843.846.231540608847.779.821740404488.498.531860404488.017.641940604489.038.632060604489.786.292140408489.786.292140408488.7910.252260408489.398.582460608488.6410.332540404889.2410.482660404889.099.582940408889.029.13306066889.029.133140608889.029.13326060888								
10 60 40 4 8 4 9.87 9.13 11 40 60 4 8 4 23.55 18.92 12 60 60 4 8 4 15.36 14.50 13 40 40 8 8 4 3.66 5.11 14 60 40 8 8 4 3.84 6.23 15 40 60 8 8 4 9.00 10.16 16 60 60 8 8 4 7.77 9.82 17 40 40 4 4 8 8.49 8.53 18 60 40 4 4 8 8.91 7.64 19 40 60 4 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 22 60 40 8 4 8 9.78 6.29 23 40 60 8 4 8 9.78 6.29 24 60 60 8 4 8 9.79 8.58 24 60 60 8 4 8 9.60 13.58 28 60 60 4 8 8 9.02 7.81 30								
114060484 23.55 18.92 126060484 15.36 14.50 134040884 3.66 5.11 146040884 3.84 6.23 154060884 9.00 10.16 166060884 7.77 9.82 174040448 8.49 8.53 186040448 8.01 7.64 194060448 9.03 8.63 206060448 9.78 6.29 214040848 9.78 6.29 214040848 9.79 10.25 226040848 9.99 8.58 246060848 9.24 10.48 254040488 9.24 10.48 266040488 9.09 9.58 286060488 9.02 9.13 3060488 9.02 9.13 314060888 9.27 9.36 314060888 9.27 9.36								
12 60 60 4 8 4 15.36 14.50 13 40 40 8 8 4 3.66 5.11 14 60 40 8 8 4 3.84 6.23 15 40 60 8 8 4 9.00 10.16 16 60 60 8 8 4 7.77 9.82 17 40 40 4 4 8 8.49 8.53 18 60 40 4 4 8 9.03 8.63 20 60 60 4 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.18 13.46 23 40 60 8 4 8 9.39 8.58 24 60 60 8 4 8 9.24 10.48 25 40 40 4 8 8 9.27 9.36 31 40 60 4 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.22 9.13 32 60 60 8 8 8 9.27 9.36 31<								
134040884 3.66 5.11 146040884 3.84 6.23 154060884 9.00 10.16 166060884 7.77 9.82 174040448 8.49 8.53 186040448 9.03 8.63 206060448 9.03 8.63 206060448 9.78 6.29 214040848 9.78 6.29 214040848 9.18 13.46 234060848 9.39 8.58 246060848 9.24 10.48 254040488 9.24 10.48 266040488 9.27 9.36 314060488 9.27 9.36 314060888 9.27 9.36 314060888 9.22 9.13 326060888 9.27 9.36 314060888 9.22 9.13 326060888 9.22 9.13 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
146040884 3.84 6.23 1540608849.0010.161660608847.779.821740404488.498.531860404489.038.632060604489.786.292140408489.786.292140408489.1813.462340608489.398.582460608489.2410.482540404889.2410.482660404889.6013.582860604889.099.58294040889.029.133060408889.029.133260608889.029.133260608889.039.223330506669.6610.623550306669.267.293650706669.367.66395050669.367.6								
1540608849.0010.161660608847.779.821740404488.498.531860404489.038.632060604489.786.292140408489.786.292140408489.786.292140408489.398.582260408489.398.582460608489.398.582460608489.2410.482540404889.2410.482660404889.099.58294040889.279.363140608889.279.363140608889.229.133260608889.299.133260608889.229.133330506669.6910.99375050269.6610.6238505010669.367.66 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
16 60 60 8 8 4 7.77 9.82 17 40 40 4 4 8 8.49 8.53 18 60 40 4 4 8 9.03 8.63 20 60 60 4 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.78 6.29 21 40 40 8 4 8 9.79 8.58 22 60 40 8 4 8 9.18 13.46 23 40 60 8 4 8 9.39 8.58 24 60 60 8 4 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 26 60 4 8 8 9.09 9.58 29 40 40 8 8 9.09 9.58 29 40 40 8 8 9.02 9.13 32 60 60 8 8 8 9.22 33 30 50 6 6 9.21 8.53 35 50 30 6 6 9.66 10.62 38 50 50 6 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>								
174040448 8.49 8.53 186040448 9.03 8.63 194060448 9.03 8.63 206060448 9.78 6.29 214040848 9.78 6.29 226040848 9.18 13.46 234060848 9.39 8.58 246060848 9.39 8.58 246060848 9.44 10.48 266040488 9.60 13.58 286060488 9.09 9.58 29404088 8 9.22 7.81 30604088 9.02 9.13 32606088 8 9.22 33305066 9.21 8.53 35503066 6 9.66 385050106 9.36 7.66 39505066 9.36 7.66 39505066 9.24 9.14 44505066 6 9.24 9.14 44505066 6 <								
1860404488.01 7.64 1940604489.038.632060604489.78 6.29 2140408489.78 6.29 2260408489.18 13.46 2340608489.398.582460608489.24 10.48 2540404889.24 10.48 2660404889.60 13.58 2860604889.099.582940408889.279.363140608889.229.133260608889.299.133260608889.039.223330506669.218.533550306669.6610.6238505010669.367.0340505066109.429.7243505066109.429.724450506669.819.1444505066 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
1940604489.03 8.63 206060604489.78 6.29 2140408489.78 6.29 2260408489.18 13.46 2340608489.39 8.58 2460608489.39 8.58 2460608489.24 10.48 2540404889.24 10.48 2660404889.60 13.58 2860604889.09 9.58 294040888 9.27 9.36 314060888 9.27 9.36 314060888 9.22 9.13 326060888 9.03 9.22 33305066 6 9.66 10.62 3850501066 9.66 10.62 385050669.03 7.03 4050506610 9.42 9.72 415050669.81 9.14 445050666 9.24 9.14								
20 60 60 4 4 8 9.78 6.29 21 40 40 8 4 8 8.79 10.25 22 60 40 8 4 8 9.18 13.46 23 40 60 8 4 8 9.39 8.58 24 60 60 8 4 8 9.39 8.58 24 60 60 8 4 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 27 40 60 4 8 8 9.29 7.81 30 60 4 8 8 9.09 9.58 29 40 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 9.66 10.62 38 50 50 6 6 9.66 10.62 38 50 50 6 6 6 9.81 9.14 41 <								
21 40 40 8 4 8 8.79 10.25 22 60 40 8 4 8 9.18 13.46 23 40 60 8 4 8 9.39 8.58 24 60 60 8 4 8 9.39 8.58 24 60 60 8 4 8 8.64 10.33 25 40 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 26 60 40 8 8 9.24 10.48 27 40 60 4 8 8 9.09 9.58 29 40 40 8 8 9.09 9.58 29 40 40 8 8 9.09 9.58 29 40 40 8 8 9.02 7.81 30 60 4 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 32 60 60 6 6 6 9.66 10.99 35 50 70								
22 60 40 8 4 8 9.18 13.46 23 40 60 8 4 8 9.39 8.58 24 60 60 8 4 8 8.64 10.33 25 40 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.60 13.58 28 60 60 4 8 8 9.09 9.58 29 40 40 8 8 8 9.29 7.81 30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 35 50 70 6 6 6 9.66 10.62 </td <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>								
234060848 9.39 8.58 246060848 8.64 10.33 254040488 9.24 10.48 266040488 9.24 10.48 274060488 9.60 13.58 286060488 9.09 9.58 294040888 10.29 7.81 306040888 9.27 9.36 314060888 9.92 9.13 326060888 9.03 9.22 333050666 10.32 9.75 347050666 9.21 8.53 355030666 9.66 10.62 3850501066 9.36 7.66 3950506106 13.05 11.25 415050666 9.81 9.14 445050666 9.24 9.14								
24 60 60 8 4 8 8.64 10.33 25 40 40 4 8 8 9.24 10.48 26 60 40 4 8 8 9.24 10.48 27 40 60 4 8 8 9.60 13.58 28 60 60 4 8 8 9.09 9.58 29 40 40 8 8 8 9.09 9.58 29 40 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.22 9.13 32 60 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.66 10.62 38 50 50 10 6 6 9.36 7.66 39 50 50 6 6 10 9.42 9.72 41 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.24 9.14								
25 40 40 4 8 8 9.24 10.48 26 60 40 4 8 8 8.46 7.94 27 40 60 4 8 8 9.60 13.58 28 60 60 4 8 8 9.09 9.58 29 40 40 8 8 9.09 9.58 29 40 40 8 8 8 9.02 7.81 30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.22 9.13 32 60 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.66 10.62 38 50 50 2 6 6 9.66 10.62 38 50 50 6 6 2 9.90 8.56 42 50 50 6 6 10 9.42 9.72 43 50 50 6 6 9.24 9.14 44 50 50 6 6 9.24 9.14 45 50 50 6 6 9.72 9.14								
26 60 40 4 8 8 8.46 7.94 27 40 60 4 8 8 9.60 13.58 28 60 60 4 8 8 9.09 9.58 29 40 40 8 8 8 9.09 9.58 29 40 40 8 8 8 9.02 7.81 30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 33 30 50 6 6 6 9.21 8.53 35 50 70 6 6 6 9.66 10.62 38 50 50 6 10 6 13.05 11.25 41 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.24								
27 40 60 4 8 8 9.60 13.58 28 60 60 4 8 8 9.09 9.58 29 40 40 8 8 8 10.29 7.81 30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.69 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 6 10 6 9.36 7.66 39 50 50 6 10 6 13.05 11.25 41 50 50 6 6 9.81 9.14 44 50 50 6 6 9.24 9.14 44 50 50 6 6 6 9.24 9.14								
28 60 60 4 8 8 9.09 9.58 29 40 40 8 8 8 10.29 7.81 30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 8.85 7.29 36 50 70 6 6 6 9.66 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 6 2 6 9.03 7.03 40 50 50 6 10 6 13.05 11.25 41 50 50 6 6 9.81 9.14 44 50 50 6 6 6 9.24 9.14 44 50 50 6 6 6 9.72 9.14								
29 40 40 8 8 8 10.29 7.81 30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.21 8.53 35 50 30 6 6 6 9.66 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 6 2 6 9.03 7.03 40 50 50 6 10 6 13.05 11.25 41 50 50 6 6 10 9.42 9.72 43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.72 9.14								
30 60 40 8 8 8 9.27 9.36 31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.92 9.13 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.69 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 10 6 9.36 7.66 39 50 50 6 10 6 13.05 11.25 41 50 50 6 6 9.81 9.14 44 50 50 6 6 9.24 9.14 44 50 50 6 6 9.72 9.14								
31 40 60 8 8 8 9.92 9.13 32 60 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.21 8.53 35 50 30 6 6 6 9.69 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 10 6 6 9.36 7.66 39 50 50 6 10 6 13.05 11.25 41 50 50 6 6 2 9.90 8.56 42 50 50 6 6 9.81 9.14 44 50 50 6 6 9.24 9.14 45 50 50 6 6 9.72 9.14								
32 60 60 8 8 8 9.03 9.22 33 30 50 6 6 6 10.32 9.75 34 70 50 6 6 6 9.21 8.53 35 50 30 6 6 6 9.21 8.53 35 50 30 6 6 6 9.21 8.53 36 50 70 6 6 6 9.69 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 10 6 6 9.36 7.66 39 50 50 6 2 6 9.03 7.03 40 50 50 6 10 6 13.05 11.25 41 50 50 6 6 10 9.42 9.72 43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.72 9.14								
$\begin{array}{cccccccccccccccccccccccccccccccccccc$								
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	32	60	60	8	8	8	9.03	9.22
$\begin{array}{cccccccccccccccccccccccccccccccccccc$		30	50	6	6	6	10.32	9.75
36 50 70 6 6 6 9.69 10.99 37 50 50 2 6 6 9.66 10.62 38 50 50 10 6 6 9.36 7.66 39 50 50 6 2 6 9.03 7.03 40 50 50 6 10 6 13.05 11.25 41 50 50 6 6 2 9.90 8.56 42 50 50 6 6 10 9.42 9.72 43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.72 9.14	34	70	50	6	6	6	9.21	8.53
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	35	50	30	6	6	6	8.85	7.29
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	36	50	70	6	6	6	9.69	10.99
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	37	50	50	2	6	6	9.66	10.62
40 50 50 6 10 6 13.05 11.25 41 50 50 6 6 2 9.90 8.56 42 50 50 6 6 10 9.42 9.72 43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.24 9.14 45 50 50 6 6 6 9.72 9.14	38	50	50	10	6	6	9.36	7.66
41 50 50 6 6 2 9.90 8.56 42 50 50 6 6 10 9.42 9.72 43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.24 9.14 45 50 50 6 6 6 9.72 9.14	39	50	50	6	2	6	9.03	7.03
42 50 50 6 6 10 9.42 9.72 43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.24 9.14 45 50 50 6 6 6 9.72 9.14	40	50	50	6	10	6	13.05	11.25
43 50 50 6 6 6 9.81 9.14 44 50 50 6 6 6 9.24 9.14 45 50 50 6 6 6 9.72 9.14	41	50	50	6	6	2	9.90	8.56
44 50 50 6 6 6 9.24 9.14 45 50 50 6 6 6 9.72 9.14	42	50	50	6	6	10	9.42	9.72
44 50 50 6 6 6 9.24 9.14 45 50 50 6 6 6 9.72 9.14	43	50	50	6	6	6	9.81	9.14
45 50 50 6 6 6 9.72 9.14	44		50	6	6			
	45	50		6	6	6		
	46	50	50	6	6	6	10.32	9.14

 Table 4 (continued)

Run	2	Mannitol (g/kg)	Extract		-		MK-7
47	50	50	6	6	6	10.08	9.14
48	50	50	6	6	6	10.26	9.14
49	50	50	6	6	6	9.87	9.14
50	50	50	6	6	6	9.24	9.14

and Inoculum Volume (ml/100 g) were studied at different levels (-1, 0, +1) with 17 treatment combinations as presented in Table 6. Data collected for MK-7 concentrations in each run was analysed using the design expert 7.1.6 software.

Extraction of menaquinone-7

Menaquinone-7 was extracted from the fermented soybean samples using solvents propan-2-ol: n-hexane (1:2 v/v) as described by Berenjian et al. (2011b). Crushed samples (5 g) were dispersed in 15 mL of each solvent and vigorously mixed at room temperature for 10 min. This dispersion was

 Table 5
 Analysis of variance values for model of MK-7 production by

 Bacillus subtilis NCIM 2708

Source	Sum of Squares	DF	Mean Square	F value	p-value	
Model	301.53	15	20.1	3.61	0.0009	Significant
A-Glycerol	3.74	1	3.74	0.67	0.4184	
B-Mannitol	34.24	1	34.24	6.14	0.0183	
C-Yeast Extract	21.8	1	21.8	3.91	0.0561	
D-Malt Extract	44.33	1	44.33	7.96	0.0079	
E-CaCl ₂	3.31	1	3.31	0.59	0.4461	
AB	4.23	1	4.23	0.76	0.3899	
AC	33.44	1	33.44	6	0.0196	
AD	5.49	1	5.49	0.98	0.3281	
AE	0.37	1	0.37	0.066	0.7988	
BC	6.33	1	6.33	1.14	0.2941	
BD	17.9	1	17.9	3.21	0.082	
BE	27.81	1	27.81	4.99	0.0322	
CD	38.65	1	38.65	6.94	0.0126	
CE	37.3	1	37.3	6.69	0.0141	
DE	22.6	1	22.6	4.05	0.052	
Residual	189.47	34	5.57			
Adequate precision	11.035					
Lack of Fit				40.79		Significant

R²=0.61, Adj R²=0.4439

Variables	Levels				
	-1	0	+1		
Medium pH	4	7	10		
Temperature (°C)	27	37	47		
Inoculum Volume (ml/100 g)	5	20	35		

tion process parameters used in Box-Behnken's Design

centrifuged at 3,000 rpm for 5 min. The organic layers were separated and concentrated up to 1 mL and filtered through 0.45 μ m membrane filter unit prior to HPLC analysis.

Analysis of menaquinone-7

Menaquinone-7 analysis was performed by HPLC (SHIMADZU, Japan) using a 250 mm x 4.6 mm ID Lichrosper[®] 100 C18 column containing 5 μ m sized particles and a 20 ml loop injector. The optimum separation condition of





Fig. 1 Percentage contribution of different medium parameters and their effects on MK-7 productivity during solid state fermentation process of soybean obtained by Plackett-Burman Design

Fig. 2 The HPLC chromatogram of standard MK-7 (a) and extracted MK-7 (b) from fermented soybean




HPLC was carried out with a multi-step gradient elution at a flow-rate of 1.2 ml/min with the following system: B: acetonitrile, C: Mixture of water and methanol (1:1) acidified to pH 3.0 by ortho-phosphoric acid. After sample injection (20 μ L), the system was maintained at 80 % B for 3.5 min and then increased to 100 % in 4.5 min, held for 2 min and at the end of 10 min, the system was recycled back to 80 % B. The detection was carried out by UV detector at wavelength 254 nm.

Results and discussion

Effect of medium components

The effect of carbon source (glycerol, dextrose, mannitol, sucrose), nitrogen source (malt extract, yeast extract) was investigated in this experiment. Among the entire carbon sources, it was found that dextrose has high contribution, but dextrose was not selected for the RSM study because recommended dietary adequate intake (RDAI) of sugar may exceed above the normal level in patients suffering from diabetic disease. So the glycerol and mannitol was found to have contributing factor 13.09 %, 4.59 %, respectively. The malt extract is better nitrogen source than yeast extract and showed 21.6 % contribution in MK-7 production. Similarly out of all the micro and macro nutrients (K₂HPO₄, MgSO₄, CaCl₂) CaCl₂ was found to be best nutrient for MK-7 production with 1.27 % contribution. All the results of Plackett-Burman experimental design with % contribution towards MK-7 production are shown in Table 2. Sucrose has negative effect on menaquinone-7 production while rest of nutrients including glycerol, mannitol, dextrose, malt extract, yeast extract, K₂HPO₄, CaCl₂ and MgSO₄ has significant positive effect on menaquinone-7 (MK-7) production from soybean by



Bacillus subtilis NCIM 2708 (Fig. 1). The HPL chromatograms of standard MK-7 and extracted MK-7 is shown in Fig. 2.

The PBD is the preliminary technique for rapid illustration of the effects of various medium constituents. It tests each variable at two levels only; hence it cannot give exact idea regarding the optimum level of constituent required in the

 Table 7
 The experimental MK-7 concentrations of fermentation trials carried out by Box-Behnken's Design of response surface methodology for optimization of fermentation process parameters

Run	Mediun pH	Fermentation Temperature (°C)	Inoculum Volume (% v/w)	Menaquinone- 7 (µg/g)
1	4	27	20	7.35
2	10	27	20	0
3	4	47	20	5.76
4	10	47	20	6.27
5	4	37	5	5.28
6	10	37	5	0
7	4	37	35	5.79
8	10	37	35	5.67
9	7	27	5	4.65
10	7	47	5	0
11	7	27	35	4.44
12	7	47	35	3.78
13	7	37	20	3.69
14	7	37	20	5.4
15	7	37	20	5.49
16	7	37	20	3.75
17	7	37	20	5.07

medium. Therefore, further optimization of selected nutrients for menaquinone-7 production is necessary.

Optimization of selected medium component

After screening the medium components viz. mannitol, glycerols, yeast extract, malt extract, $CaCl_2$ were chosen as best C, N and ion sources for further optimization using RSM and CCD. The predicted and experimental MK-7 concentrations were obtained from CCD response surface design in each run are shown in Table 4. Results were analysed using the software Design Expert 7.1.6 and fitted into a multiple nonlinear regression model, resulting in the following equation for MK-7 production.

MK-7(µg/g) =9.14–0.31 A+0.93 B-0.74 C+1.05 D+0.29 E-0.36 AB+1.02 AC-.41 AD+0.11AE - 0.44BC+0.75BD-0.93BE-1.10CD+1.08CE-0.84DE

Where A, B, C, D and E represents glycerol, mannitol, yeast extract, malt extract and calcium chloride respectively in g/kg. The effects of all nutrient parameters on MK-7 production can be compared with the help of response surface plots. Nonlinear regression models resulted in ten response graphs for MK-7 production, shown in Fig. 3. The result adequacies were analysed using the analysis of variance (Table 5) and determination of coefficient R^2 . In this study, the value R^2 is 0.6141 that indicates 61 % variability in the response could be explained by the model. According to ANOVA of nonlinear regression model, F-value" of 3.61 implies the model is significant and a very low probability value (P < 0.0500). There is only a 0.09 % chance that a "model F-value" this large could occur due to noise. "Adequate precision" measures the signal to noise ratio and ratio greater than 4 is desirable. The Adequate precision value 11.035 indicates an adequate signal

 Table 8
 Analysis of variance values for fermentation process parameters of MK-7 production by *Bacillus subtilis* NCIM 2708

Source	Sum of Squares	DF	Mean Square	F Value	<i>p</i> -value	
Model	58.74	6	9.46	4.05	0.0255	Significant
Medium pH	18.73	1	18.73	8.01	0.0178	
Temperature	0.05	1	0.05	0.021	0.8871	
Inc. Volume	11.88	1	11.88	5.08	0.0478	
AB	15.44	1	15.44	5.08	0.0279	
AC	6.66	1	6.66	6.61	0.1224	
BC	3.98	1	3.98	2.85	0.2212	
Residual	23.38	10	2.34	1.7		
Lack of Fit				4.25		Significant
Adequate precision						7.759

R² =0.708, Adj R² =0.5331

which implies that this model can be used to navigate the design space. A negative "Pred R-Squared" implies that the overall mean is a better predictor of your response than the current model ". The "Lack of Fit F-value" of 40.79 implies the Lack of Fit is significant. There is only a 0.01 % chance that a "Lack of Fit F-value" this large could occur due to noise.

Bacillus subtilis NCIM 2708 produce higher amount of MK-7 in glycerol containing media which is in agreement with Berenjian et al. 2011a, b. Yeast extract and malt extract act as undefined nitrogen source results higher MK-7 production since complex nitrogen source contain amino acids, vitamins and provide better genetic requirement for protein & enzyme synthesis.

Optimization of fermentation process


In the study of process parameters, three independent variables (e.g. Temperature, pH, Inoculum volume) were tested. The results of experimental run with predicted and actual MK-7 level in fermentation medium is presented in Table 7, calculation for the analysis of variance is presented in Table 8. The Model F-value of 4.05 and very low probability value (P < 0.0500) implies the model is significant. There is only a 2.55 % chance that a "Model F-Value" this large could occur due to noise. The "Lack of Fit F-value" of 4.25 implies there is a 9.15 % chance that a "Lack of Fit F-value" this large could occur due to noise. "Adequate precision" measures the signal to noise ratio. Adequate precision grater than 4 is desirable and it is 7.759 indicates an adequate signal. Results were analysed using the software Design Expert 7.1.6 and fitted into a multiple nonlinear regression model, resulting in the following equation for MK-7 production.

MK-7 (µg/g) =4.26–1.53 A - 0.079 B +1.22 C +1.97 AB + 1.29 AC +1.00 BC

The effects of all process parameters on MK-7 production can be compared with the help of response plots. Nonlinear regression models resulted in three response graphs for MK-7 production, shown in Fig. 4.

It was found that low pH of medium has positive effect on MK-7 production. The results agree with Hill et al. 1990. Fermentation temperature is inversely proportional to MK-7 production. Optimum temperature for *Bacillus subtilis* was 37 °C, but MK-7 production was found to be decreased with increase in temperature. No measurable effects of Inoculum volume was found on MK-7 production. There was only

Fig. 4 Response surface plots showing the effect of process parameters (pH, temperature and inoculum volume) on MK-7 production during solid state fermentation of soybean



slightly increase in MK-7 production with increase in Inoculum volume. Several research reports maximum MK-7 production as high as 62.32 mg/L under submerged fermentation at 6th day (Berenjian et al. 2011b). However for The present research reports a new optimized fermentation medium and process condition which yield higher MK-7 production (39.039 μ g/g) under solid state fermentation condition only after 24 h.

Conclusions

The *Bacillus subtilis* NCIM 2708 was used for production of MK-7 by solid state fermentation. Statistical experimental design proved to be valuable tools in optimizing a solid state fermentation medium for MK-7 production by *Bacillus subtilis* NCIM 2708. Five medium component viz. glycerol 40 ml/kg, mannitol 60 g/kg, yeast extract 4 g/kg, malt extract 8 g/kg, and calcium chloride 4 g/kg were identified and optimised by RSM. The predicted MK-7 was 56.75 μ g/g by using above optimised medium combination. Practically MK-7 concentration 39.039 μ g/g was achieved in the optimized fermented media. This shows 68.78 % validity of the predicted model.

Acknowledgments The authors gratefully Professor S. K. Lal, IARI, PUSA, New Delhi for providing soybeans variety. All the authors acknowledge the Department of Science and Technology (DST), Govt. of India for providing fellowship for scholar under INSPIRE Programme.

References

- Berenjian A, Mahanama R, Talbot A, Biffin R, Regtop H, Kavanagh J, Dehghani F (2011a) The effect of amino-acids and glycerol addition on MK-7 production. Proc World Congr Eng Comput Sci 11:19–21
- Berenjian A, Mahanama R, Talbot A, Biffin R, Regtop H, Valtchev P, Kavanagh J, Dehghani F (2011b) Efficient media for high menaquinone-7 production: response surface methodology approach. New Biotechnol 28(6):665–672
- Beulens JWJ, Bots ML, Atsma F, Bartelink MLEL, Prokop M, Geleijnse JM, Witteman JCM, Grobbee DE, Van Der Schouw YT (2009) High dietary menaquinone intake is associated with reduced coronary calcification. Atherosclerosis 203(2):489–493
- Fujita Y, Iki M, Tamaki J, Kouda K, Youra A, Kadowaki E, Sato Y, Moon JS, Tomiaka N, Okamoto N, Kurumatani N (2011) Association between vitamin K intake from fermented soybeans, *natto*, and bone mineral density in elderly Japanese men: the Fujiwara-kyo Osteoporosis Risk in Men (FORMEN) study. Osteoporosis Int 23(2):705–714
- Gast GCM, De Roos NM, Sluijs I, Bots ML, Beulens JWJ, Geleijnse JM, Witteman JC, Grobbee DE, Peeters PHM, Van Der Schouw YT (2009) A high menaquinone intake reduces the incidence of coronary heart disease. Nutr Metab Cardiovas 19(7):504–510

- Hill KF, Mueller JP, Taber HW (1990) The *Bacillus subtilis men* CD promoter is responsive to extracellular pH. Arch Microbiol 153: 355–359
- Hu H, Yao S, Mei L, Zhu Z, Hur B (2000) Partial purification of nattokinase from *Bacillus subtilis* by expanded bed adsorption. Biotechnol Lett 22(17):1383–1387
- Kamao M, Suhara Y, Tsugawa N, Uwano M, Yamaguchi N, Uenishi K, Ishida H, Saski S, Okano T (2007) Vitamin K content of foods and dietary vitamin K intake in japanese young women. J Nutr Sci Vitaminol 53(6):464–470
- Kennedy M, Krouse D (1999) Strategies for improving fermentation medium performance: a review. J Ind Microbiol Biot 23(6):456–475
- Mahanama R, Berenjian A, Valtchev P, Talbot A, Biffin R, Regtop H, Dehghani F, Kavanagh JM (2011) Enhanced production of menaquinone 7 via solid substrate fermentation from *Bacillus subtilis*. Int J Food Eng 7(5): DOI:10.2202/1556-3758.2314
- Mahanama R, Boyce D, Berenjian A, Dehghani F, Kavanagh J (2012a) The effect of enzymatic pre- treatment on solid state fermentation of vitamin K2. J Am Sci 8(1):648–654
- Mahanama R, Berenjian A, Dehghani F, Kavanagh J (2012b) Modelling the effect of bed height and particle size for vitamin K2 production in a static bed fermenter. Eng Lett 20(1):16–17
- Rowland B, Taber H (1996) Duplicate isochorismate synthase genes of *Bacillus subtilis*: regulation and involvement in the biosyntheses of menaquinone and 2,3-dihydroxybenzoate. J Bacteriol 178(3):854– 861
- Rowland B, Hill K, Miller P, Driscoll J, Taber H (1995) Structural organization of a Bacillus subtilis operon encoding menaquinone biosynthetic enzymes. Gene 167(1–2):105–109
- Sakano T, Notsumoto N, Nagaoka N (1988) Measurement of K vitamins in food by high performance liquid chromatography with fluorometric detection. Vitam (Jpn) 62(8):393–398
- Sato T, Yamada Y, Ohtani Y, Mitsui N, Murasawa H, Araki S (2001a) Efficient production of menaquinone (Vitamin K2) by minadioneresistant mutant of *Bacillus subtilis*. J Ind Microbiol Biot 26(3):115– 120
- Sato T, Yamadaa Y, Ohtania Y, Mitsuib N, Murasawab H, Araki S (2001b) Production of Menaquinone (vitamin K2)-7 by *Bacillus* subtilis. J Biosci Bioeng 91(1):16–20
- Schurgers LJ, Teunissen KJF, Hamulyák K, Knapen MHJ, Vik H, Vermeer C (2007) Vitamin K-containing dietary supplements: comparison of synthetic vitamin K1 and natto-derived menaquinone-7. Blood 109(8):3279–3283
- Seto H, Jinnai Y, Hiratsuka T, Fukawa M, Furhata K, Itoh N, Dairi T (2008) Studies on a new biosynthetic pathway for menaquinone. J Am Chem Soc 130(17):5614–5615
- Shearer MJ (1990) Vitamin K and vitamin K-dependent proteins. Brit J Haematol 75(2):156–162
- Sonenshein AL, Hoch JA, Losick R (eds) (2002) *Bacillus subtilis* and its closest relatives: from genes to cells. ASM Press, Washington DC
- Theuwissen E, Smit E, Vermeer C (2012) The role of vitamin K in softtissue calcification. Adv Nutr 3:166–173
- Truong JT, Booth SL (2011) Emerging issues in vitamin K research. J Evid Based Complement Altern Med 16(2):73–79
- Tsukamoto Y, Kasai M, Kakuda H (2001) Construction of a *Bacillus subtilis* (natto) with high productivity of vitamin K2 (menaquinone-7) by analog resistance. Biosci Biotechnol Biochem 65(9):2007–2015
- Vermeer C, Schurgers LJ (2000) A comprehensive review of vitamin K and vitamin K antagonists. Hematol Oncol Clin North Am 14(2): 339–353
- Wu WJ, Ahn BY (2011) Improved menaquinone (vitamin K2) production in cheonggukjang by optimization of the fermentation conditions. Food Sci Biotechnol 20(6):1585–1591