

Fuji apple storage time rapid determination method using Vis/NIR spectroscopy

Fuqi Liu and Xuxiang Tang*

Office of Laboratory and Assets Management; Zhejiang Gongshang University; Hangzhou, China

Keywords: Fuji apple, principal component analysis, storage time, signal-to-noise ratio, Vis/NIR spectroscopy

Fuji apple storage time rapid determination method using visible/near-infrared (Vis/NIR) spectroscopy was studied in this paper. Vis/NIR diffuse reflection spectroscopy responses to samples were measured for 6 days. Spectroscopy data were processed by stochastic resonance (SR). Principal component analysis (PCA) was utilized to analyze original spectroscopy data and SNR eigen value. Results demonstrated that PCA could not totally discriminate Fuji apples using original spectroscopy data. Signal-to-noise ratio (SNR) spectrum clearly classified all apple samples. PCA using SNR spectrum successfully discriminated apple samples. Therefore, Vis/NIR spectroscopy was effective for Fuji apple storage time rapid discrimination. The proposed method is also promising in condition safety control and management for food and environmental laboratories.

Introduction

Apples possess various nutrients including vitamins, dietary fibers, amino acids, etc.^{1–3} In recent years, apple has been one of the most popular fruits among customers for its anti-oxidation, anti-senescence, immunity strengthens, and other beneficial functions.^{3–5} So apple has been a valuable nutrition resource in human's daily life. However, apple easily deteriorates after post-harvest due to its high water content and abundant nutrients. Living cells still conduct normal respiration metabolism by consuming large amount of nutrients. In addition, fungi invasion also plays an important role in apple quality decrease.⁶ Therefore, there is a great demand for a novel rapid and nondestructive apple quality discrimination technique.

Vis/NIR spectroscopy, a fast-developed nondestructive detecting technique, has been widely applied in medicine,^{7,8} chemistry^{9,10} and food research, such as meat,^{11,12} drinks,^{13,14} grains,^{15,16} fruits,^{17–19} etc. Jamshidi et al managed to conduct taste characterization of Valencia oranges using reflectance Vis/NIR spectroscopy.²⁰ Magwaza et al applied Vis/NIRS and chemometric analysis to predict fruit defects and postharvest behavior of 'mules Clementine' mandarin fruit.²¹ Zhang et al realized subtle bruises identification on fresh jujube based on NIR spectroscopy.²² So, Vis/NIR spectroscopy technique is frequently used in fruit quality rapid and nondestructive assessment.

In this paper, Fuji apple quality discrimination using nondestructive Vis/NIR spectroscopy were performed. Spectroscopy data were recorded and processed by nonlinear SR method. PCA method was used to analyze original spectroscopy data and SNR eigen values. Vis/NIR spectroscopy was effective for Fuji apple storage time rapid discrimination. The proposed method is promising in fruit freshness rapid detection.

Results and Discussion

Eigen peaks extraction

The original Vis/NIR diffuse reflection spectroscopy to apple samples is displayed in **Figure 1**. The main absorption band ranges from 400 to 1000 nm. The highest absorbance peak appears at about 600 nm and its absorbance reaches near 4000 counts. Eigen peaks corresponding to the wavelengths of 550 nm, 607 nm, 670 nm, 738 nm, 802 nm, and 789 nm are taken for subsequent principal component analysis of apple samples during storage.

SNR spectrum analysis result

SNR spectrum calculated by SR as a function of noise intensity is displayed in **Figure 2**. Each curve first rises with increase of stimulating noise intensity and then declines to a stable value. The initial SNR value in each curve ranges between -75 dB and -65 dB. Different valleys and peaks appear before eigen peak formation. And each SNR curve has its eigen peak at 190 stimulating noise intensity. SNR spectrum clearly discriminates Fuji apples with different storage days. And the maximal SNR value in each curve is taken for subsequent analysis.

PCA results

PCA results based on original spectroscopy data and SNR eigen value are displayed in **Figure 3**. As it shown in **Figure 3(a)**, total contribution of the first 2 components is 84.99%. The first principal component (PC1) under different storage days presents a decreasing trend; the second principal component (PC2) gradually declines. There are some significant overlaps, especially from day0 to day2. So, PCA method could not

*Correspondence to: Xuxiang Tang; Email: txx@zjgsu.edu.cn

Submitted: 02/24/2015; Revised: 03/25/2015; Accepted: 03/27/2015

<http://dx.doi.org/10.1080/21655979.2015.1038001>




Figure 1. Original Vis/NIR diffuse reflection spectroscopy.




Figure 2. SNR spectrum.

totally discriminate Fuji apples with different storage time just by original Vis/NIR spectroscopy data. However, the first 2 principal components include a much higher contribution of 95.77% (see Fig. 3(b)). Components from different storage days have no overlaps between each other. All apple samples in different storage days could be successfully discriminated. Results indicate that PCA on SR processed Vis/NIR spectroscopy data was more effective in Fuji apple storage time rapid discrimination.

In this paper, Fuji apple quality discrimination using non-destructive Vis/NIR spectroscopy were performed. Spectroscopy data were recorded and processed by non-linear SR method. PCA was used to analyze original spectroscopy data and SNR eigen values. Vis/NIR spectroscopy was effective for Fuji apple storage time rapid discrimination.

Materials and Methods

Samples

Twenty fresh Fuji apples in the same bath with near size, ripeness and color (without any pretreatments or mechanical injury) were selected and purchased from a local supermarket in Hangzhou, China. After transporting to our lab within an hour, all samples were washed with distilled water. After draining, all samples were stored at room temperature in dark place. All experiments were performed under room temperature.

Vis/NIR spectroscopy system and measurement

Diagram structure of Vis/NIR spectroscopy system is displayed in Figure 4. The system consists of 3 main parts including sampling device (light controller, sampling pedestal, optical fiber and special halogen lamp), detecting device (visible/NIR spectrometer, USB 2000+, America Ocean Photology Company), and desktop.




Figure 3. PCA results using (a) original spectroscopy data and (b) SNR eigen value.




Figure 4. Diagram structure of Visible/NIR spectroscopy system.

In spectroscopy measurement, sample was fixed in the sampling pedestal by the optical fiber with a dip angle of 75 degrees. One end of the optical fiber was connected with the spectroscopy, while the other end was connected with the special halogen lamp. Sampling pool was a sealed and light-tight cube. By replacing the sampling terminal in the sampling pool, 5 spectroscopy data from different directions could be collected in each sample, which eliminated errors from sampling direction differences. Spectroscopy analysis software (SpectraSuite, American Ocean Photology Company, America) was applied to analyze the data. Samples were analyzed every 2 d. Five samples were randomly taken in each test. All measurements were performed under room temperature.

SR

SR is a typical nonlinear model and proposed by Benzi for the explanation for Earth climate periodic changes.²³⁻²⁵ The schematic diagram of SR is displayed in Figure 5. SR phenomenon has 3 elements: a bistable system, a coherent input, and a noise source,²³ which can be described as

$$\frac{dx}{dt} = -\frac{dV(x)}{dx} + A\sin(2\pi ft + \phi) + I_g(MI(t)) + D\zeta(t) \quad (1)$$

Where x is the position of the Brownian particle, t is the time, A is periodical signal intensity, f is signal frequency. D is external noise intensity. M and D are adjustable parameters, $I(t) = EN(t) + N(t)$ denotes an EN sensor signal $EN(t)$ and




Figure 5. Schematic of SR analysis procedure.

intrinsic noise $N(t)$, $\xi(t)$ is the external noise, and $V(x)$ is the simplest double-well potential with the constants a and b characterizing the system. Noise intensity is a parameter of SR model. SR model is used as a data processing method in this research.

$$V(x) = -\frac{1}{2}ax^2 + \frac{1}{4}bx^4 \quad (2)$$

Eq. (1) can be written as

$$\frac{dx}{dt} = ax - bx^3 + A\sin(2\pi ft + \phi) + I_g(MI(t)) + D\zeta(t) \quad (3)$$

The minima of $V(x)$ are located at $\pm x_m$, where $x_m = (a/b)^{1/2}$. A potential barrier separates the minima with the height given by $\Delta U = a^2/4b$. The barrier top is located at $x_b = 0$. When three elements of SR interact coherently, the potential barrier can be reduced and the Brownian particle may surmount the energy barrier and enter another potential well.²³ The intensity of signals will increase, which makes it possible that the weak signal can be detected from noise background. SNR is the common quantifier for SR and it can be approximately described as²³

$$\text{SNR} = \sqrt{2}\Delta U \left(\frac{A}{D}\right)^2 e^{-\Delta U/D}. \quad (4)$$

Conclusions

In this paper, Fuji apple storage time rapid determination method using Vis/NIR spectroscopy was proposed. Vis/NIR spectroscopy responses to Fuji apples were examined during 6 days' storage. Spectroscopy data were processed by SR technique. And PCA method was used to analyze original spectroscopy data and SNR eigen values. The results demonstrated that PCA method could not totally discriminate Fuji apples with different storage time just via original Vis/NIR spectroscopy. However, SNR spectrum calculated by SR showed different maximal values under the same noise intensity. PCA successfully discriminated all samples using SNR maximal values. This method was effective for Fuji apple storage time rapid discrimination. It presents some advantages including non-destructive, rapid response, high accuracy, etc., and is also promising in condition safety control and management for food and environmental laboratories.

Disclosure of Potential Conflicts of Interest

No potential conflicts of interest were disclosed.

References

- Ma BQ, Chen J, Zheng HY, Fang T, Ogutu C, Li SH, Han YP, Wu BH. Comparative assessment of sugar and malic acid composition in cultivated and wild apples. *Food Chem* 2015; 172:86–91; PMID:25442527; <http://dx.doi.org/10.1016/j.foodchem.2014.09.032>
- Condezo-Hoyos L, Mohanty IP, Noratto GD. Assessing non-digestible compounds in apple cultivars and their potential as modulators of obese faecal microbiota *in vitro*. *Food Chem* 2014; 161:208–15; PMID:24837942; <http://dx.doi.org/10.1016/j.foodchem.2014.03.122>
- Šavikin K, Živković J, Zdunić G, Gođevac D, Đorđević B, Dojčinović B, Đorđević N. Phenolic and mineral profiles of four Balkan indigenous apple cultivars monitored at two different maturity stages. *J Food Compos Anal* 2014; 35:101–11; <http://dx.doi.org/10.1016/j.jfca.2014.05.004>
- Briones-Labarca V, Venegas-Cubillos G, Ortiz-Portilla S, Chacana-Ojeda M, Maureira H. Effects of high hydrostatic pressure (HHP) on bioaccessibility, as well as antioxidant activity, mineral and starch contents in Granny Smith apple. *Food Chem* 2011; 128:520–9; PMID:25212164; <http://dx.doi.org/10.1016/j.foodchem.2011.03.074>
- Hui GH, Wu YL, Ye DD. Fuji apple storage time predictive method using electronic nose. *Food Anal Methods* 2013; 6:82–8; <http://dx.doi.org/10.1007/s12161-012-9414-6>
- Lu HP, Lu LF, Zeng LZ, Fu D, Xiang HL, Yu T, Zheng XD. Effect of chitin on the antagonistic activity of *Rhodosporidium paludigenum* against *Penicillium expansum* in apple fruit. *Postharvest Biol Technol* 2014; 92:9–15; <http://dx.doi.org/10.1016/j.postharvbio.2014.01.009>
- Chiang W-L, Lin T-T, Sureshbabu R, Chia W-T, Hsiao H-C, Liu H-Y, Yang C-M, Sung H-W. A rapid drug release system with a NIR light-activated molecular switch for dual-modality photothermal/antibiotic treatments of subcutaneous abscesses. *J Control Release* 2015; 199:53–62; PMID:25499553; <http://dx.doi.org/10.1016/j.jconrel.2014.12.011>
- Porfire A, Rus L, Loredana AL, Tomuta I. High-throughput NIR-chemometric methods for determination of drug content and pharmaceutical properties of indapamide power blends for tableting. *J Pharmaceut Biomed Anal* 2012; 70:301–9; PMID:22884785; <http://dx.doi.org/10.1016/j.jpba.2012.07.026>
- Zhang Z, Feng WX, Su PY, Liu L, Lü XQ, Song JR, Fan DD, Wong W-K, Jones RA, Su CY. Near-infrared (NIR) luminescent Zn(II)-Ln(III)-containing (Ln=Nd, Yb or Er) Wolf Type II metallopolymer hybrid materials. *Syn Metals* 2015; 199:128–38; <http://dx.doi.org/10.1016/j.synthmet.2014.11.019>
- Alexandrino GL, Poppi RJ. NIR imaging spectroscopy for quantification of constituents in polymers thin films loaded with paracetamol. *Analytica Chimica Acta* 2013; 765:37–44; PMID:23410624; <http://dx.doi.org/10.1016/j.aca.2012.12.017>
- Alamprese C, Casale M, Sinelli N, Lanteri S, Casiraghi E. Detection of minced beef adulteration with turkey meat by UV-vis, NIR and MIR spectroscopy. *LWT-Food Sci Technol* 2013; 53:225–32; <http://dx.doi.org/10.1016/j.lwt.2013.01.027>
- Ma J, Pu HB, Sun DW, Gao WH, Qu JH, Ma KY. Application of Vis-NIR hyperspectral imaging in classification between fresh and frozen-thawed pork *Longissimus Dorsi* muscles. *Int J Refrig* 2015; 50:10–8; <http://dx.doi.org/10.1016/j.ijrefrig.2014.10.024>
- Giovenzana V, Beghi R, Guidetti R. Rapid evaluation of craft beer during fermentation process by vis/NIR spectroscopy. *J Food Eng* 2014; 142:80–6; <http://dx.doi.org/10.1016/j.jfoodeng.2014.06.017>
- Yu HY, Niu XY, Lin HJ, Ying YB, Li BB, Pan XX. A feasibility study on on-line determination of rice wine composition by Vis-NIR spectroscopy and least-squares support vector machines. *Food Chem* 2009; 113:291–6; <http://dx.doi.org/10.1016/j.foodchem.2008.06.083>
- Verdú S, Ivorra E, Sánchez AJ, Barat JM, Grau R. Study of high strength wheat flours considering their physicochemical and rheological characterization as well as fermentation capacity using SW-NIR imaging. *J Cereal Sci* 2015; 62:31–7; <http://dx.doi.org/10.1016/j.jcs.2014.11.002>
- Chen H-Z, Song Q-Q, Tang G-Q, Xu L-L. An optimization strategy for waveband selection in FT-NIR quantitative analysis of corn protein. *J Cereal Sci* 2014; 60:595–601; <http://dx.doi.org/10.1016/j.jcs.2014.07.009>
- Zhang LR, Xu HR, Gu MY. Use of signal to noise and area change rate of spectra to evaluate the Visible/NIR spectral system for fruit internal quality detection. *J Food Eng* 2014; 139:19–23; <http://dx.doi.org/10.1016/j.jfoodeng.2014.04.009>
- He H-J, Sun D-W, Wu D. Rapid and real-time prediction of lactic acid bacteria (LAB) in farmed salmon flesh using near-infrared (NIR) hyperspectral imaging combined with chemometric analysis. *Food Res Int* 2014; 62:476–83; <http://dx.doi.org/10.1016/j.foodres.2014.03.064>
- Giovenzana V, Beghi R, Buratti S, Civelli R, Guidetti R. Monitoring of fresh-cut *Valerianella locusta* Laterr. shelf life by electronic nose and VIS-NIR spectroscopy. *Talanta* 2014; 120:368–75; PMID:24468384; <http://dx.doi.org/10.1016/j.talanta.2013.12.014>
- Jamshidi B, Minaei S, Mohajerani E, Ghasseman H. Reflectance Vis/NIR spectroscopy for nondestructive taste characterization of Valencia oranges. *Comp Electron Agr* 2012; 85:64–9; <http://dx.doi.org/10.1016/j.compag.2012.03.008>
- Magwaza LS, Landahl S, Cronje PJR, Nieuwoudt HH, Mouazen AM, Nicolai BM, Terry LA, Opara UL. The use of Vis/NIRS and chemometric analysis to predict fruit defects and postharvest behaviour of 'Nules Clementine' mandarin fruit. *Food Chem* 2014; 163:267–74; PMID:24912725; <http://dx.doi.org/10.1016/j.foodchem.2014.04.085>
- Zhang SJ, Zhang HH, Zhao YR, Guo W, Zhao HM. A sample identification model for subtle bruises on the fresh jujube based on NIR spectroscopy. *Math Comp Model* 2013; 58:545–50; <http://dx.doi.org/10.1016/j.mcm.2011.10.067>
- Benzi R, Suter A, Vulpiana A. The mechanism of stochastic resonance. *J Phys A* 1981; 14:L453–6; <http://dx.doi.org/10.1088/0305-4470/14/11/006>
- Gammaitoni L, Hanggi P, Jung P, Marchesoni F. Stochastic resonance. *Rev Modern Phys* 1998; 70:223–87; <http://dx.doi.org/10.1103/RevModPhys.70.223>
- Dutta R, Das A, Stocks NG, Morgan D. Stochastic resonance-based electronic nose: a novel way to classify bacteria. *Sensors Actuat B* 2006; 115:17–27; <http://dx.doi.org/10.1016/j.snb.2005.08.033>