Published in final edited form as: *Methods*. 2014 March 1; 66(1): 106–119. doi:10.1016/j.ymeth.2013.10.008.

# Advances in allergen-microarray technology for diagnosis and monitoring of allergy: The MeDALL allergen-chip

Christian Lupinek<sup>a</sup>, Eva Wollmann<sup>a</sup>, Alexandra Baar<sup>a</sup>, Srinita Banerjee<sup>a</sup>, Heimo Breiteneder<sup>b</sup>, Barbara M. Broecker<sup>c</sup>, Merima Bublin<sup>b</sup>, Mirela Curin<sup>a</sup>, Sabine Flicker<sup>a</sup>, Tetiana Garmatiuk<sup>a</sup>, Heidrun Hochwallner<sup>a</sup>, Irene Mittermann<sup>a</sup>, Sandra Pahr<sup>a</sup>, Yvonne Resch<sup>a</sup>, Kenneth H. Roux<sup>d</sup>, Bharani Srinivasan<sup>a</sup>, Sebastian Stentzel<sup>c</sup>, Susanne Vrtala<sup>a</sup>, LeAnna N. Willison<sup>d</sup>, Magnus Wickman<sup>e,f</sup>, Karin C. Lødrup-Carlsen<sup>g</sup>, Josep Maria Antó<sup>h</sup>, Jean Bousquet<sup>i</sup>, Claus Bachert<sup>j</sup>, Daniel Ebner<sup>k</sup>, Thomas Schlederer<sup>k</sup>, Christian Harwanegg<sup>k</sup>, and Rudolf Valenta<sup>a,\*</sup>

<sup>a</sup>Division of Immunopathology, Department of Pathophysiology and Allergy Research, Center for Pathophysiology, Infectiology and Immunology, Medical University of Vienna, Austria

<sup>b</sup>Division of Medical Biotechnology, Department of Pathophysiology and Allergy Research, Center for Pathophysiology, Infectiology and Immunology, Medical University of Vienna, Austria

<sup>c</sup>Department of Immunology, University of Greifswald, Greifswald, Germany

<sup>d</sup>Department of Biological Science, Florida State University, Tallahassee, FL, USA

eSachs' Children's Hospital, Södersjukhuset, Stockholm, Sweden

<sup>f</sup>Institute of Environmental Medicine, Karolinska Institutet, Stockholm, Sweden

<sup>g</sup>Department of Pediatrics, Oslo University Hospital, The University of Oslo, Norway

<sup>h</sup>Centre for Research in Environmental Epidemiology (CREAL), IMIM (Hospital del Mar Research Institute), Universitat Pompeu Fabra, Departament de Ciències Experimentals i de la Salut, CIBER Epidemiología y Salud Pública (CIBERESP), Barcelona, Spain

<sup>i</sup>University Hospital of Montpellier, Hôpital Arnaud de Villeneuve, Montpellier, INSERM 1018, Villejuif, France

<sup>j</sup>Upper Airways Research Laboratory, University Hospital Ghent, Belgium

<sup>k</sup>Phadia Multiplexing, Thermo Fisher Scientific, Vienna, Austria

# Abstract

Allergy diagnosis based on purified allergen molecules provides detailed information regarding the individual sensitization profile of allergic patients, allows monitoring of the development of allergic disease and of the effect of therapies on the immune response to individual allergen molecules. Allergen microarrays contain a large variety of allergen molecules and thus allow the simultaneous detection of allergic patients' antibody reactivity profiles towards each of the

<sup>\*</sup>Corresponding author. Address: Division of Immunopathology, Dept. of Pathophysiology and Allergy Research, Center for Pathophysiology, Infectiology and Immunology, Medical University of Vienna, Waehringer Guertel 18-20, 3Q, 1090 Vienna, Austria. Rudolf.valenta@meduniwien.ac.at (R. Valenta).

allergen molecules with only minute amounts of serum. In this article we summarize recent progress in the field of allergen microarray technology and introduce the MeDALL allergen-chip which has been developed for the specific and sensitive monitoring of IgE and IgG reactivity profiles towards more than 170 allergen molecules in sera collected in European birth cohorts. MeDALL is a European research program in which allergen microarray technology is used for the monitoring of the development of allergic disease in childhood, to draw a geographic map of the recognition of clinically relevant allergens in different populations and to establish reactivity profiles which are associated with and predict certain disease manifestations. We describe technical advances of the MeDALL allergen-chip regarding specificity, sensitivity and its ability to deliver test results which are close to *in vivo* reactivity. In addition, the usefulness and numerous advantages of allergen microarrays for allergy research, refined allergy diagnosis, monitoring of disease, of the effects of therapies, for improving the prescription of specific immunotherapy and for prevention are discussed.

#### Keywords

Allergen-microarray; Sensitivity; Blocking antibodies; Allergy-diagnosis; Recombinant allergen; Immunoglobulin E; Immunoglobulin G

### 1. Introduction

Due to the continuing rise of the worldwide prevalence of allergic diseases the development of strategies for early diagnosis, prevention and treatment of allergic sensitizations becomes an increasingly important issue [1]. The development of IgE responses to allergen molecules, a process termed "allergic sensitization" seems to occur early in life whereas IgE-reactivity profiles of adult allergic patients do not change in a relevant manner [2,3]. This implies that there exists a period early in life during which allergic sensitization takes place and IgE-reactivity profiles are established. This period thus represents a window of opportunity for preventing allergic sensitization and early intervention strategies but the exact duration of this period of plasticity of the allergic immune response has not yet been studied in detail. Likewise, it remains to be investigated whether the development of certain clinical manifestations of allergy (e.g., asthma) and their severity [4] are linked to distinct allergen recognition profiles, levels of pathogenic (i.e., IgE) and protective (i.e., IgG) antibody responses and epitope recognition profiles which would allow early prediction of the clinical outcome. In order to address these questions the European Union-funded project MeDALL (Mechanisms for the development of allergies) (http://medall-fp7.eu/) has been established in which several birth cohorts in allergy and asthma are pooled and serum samples from mothers and children, from birth up to adolescence, are analyzed for allergenspecific IgE and IgG reactivity profiles [5].

Since the amounts of serum samples available from these birth cohorts are very limited a highly sensitive assay allowing to measure IgE and IgG reactivities to a large number of different allergen molecules with minute serum samples is needed. Allergen microarrays containing large numbers of different allergen molecules have been developed already in 2002 [6] with the goal to allow simultaneous detection of IgE reactivities to a large number

of allergen molecules with small serum volumes. These allergen microarrays are comprised of chips containing spotted allergens and were suggested as tools for monitoring the development of allergic sensitization, for the establishment of comprehensive reactivity profiles and to measure eventual changes during the course of disease and various forms of treatment. Today, an allergen microarray termed ImmunoCAP ISAC (Phadia AB, Uppsala, Sweden) is already available for chip-based allergy diagnosis and can be used for determining IgE and IgG responses towards 112 allergen molecules.

The ImmunoCAP ISAC has been used in a large number of clinical studies for a number of years [7], in epidemiologic studies [8], longitudinal follow up of children [9] and to find *de novo* sensitizations [3]. It has many advantages and results are close but not identical to classical immunoassays [10] (WAO Position Paper, in press). However, for some allergens [11-13] and for low IgE levels [14] it may not be sufficiently sensitive. Moreover, important allergens are not included on the ImmunoCAP ISAC chip. The MeDALL allergen-chip was designed to increase the number of allergens on the chip and to improve sensitivity of the test.

We discuss specific features and advantages of allergen microarrays and describe the MeDALL allergen-chip which has been developed for the analysis of sera from birth cohorts pooled in the MeDALL project. The MeDALL allergen-chip represents a comprehensive collection of more than 170 allergen molecules for the reliable detection of allergen-specific antibody signatures. It should facilitate to study the early evolution of the allergic immune response and the generation of a European map of molecular sensitization profiles as a prerequisite for intervention and treatment strategies.

A second important aspect of this study is to demonstrate and clarify major differences between serological test results obtained by microarray and by commonly used systems for measuring quantitative IgE responses. In fact, chip-based IgE measurements provide fundamentally different information compared to quantitative IgE measurements using allergen excess and therefore cannot be directly compared. Both assays are equally sensitive and specific but measurements performed under conditions of low allergen amount allow for the estimation of the effect of allergen-specific blocking IgG on IgE binding to the allergen and thus will more closely resemble biological responses under conditions of natural exposure. These differences imply that chip-based IgE measurements should provide information of higher clinical relevance.

#### 2. Material and methods

#### 2.1. Design of the MeDALL allergen-chip

The development of the allergen microarray platform used for preparing the MeDALL allergen-chip has been described earlier [6]. In order to improve the diagnostic spectrum of allergens on the array and to avoid redundant markers, the allergen panel of the current version of ImmunoCAP ISAC (Phadia AB) has been modified by adding more than 70 allergen molecules and by removing cross-reactive allergens which were present in different versions on the chip (Table 1). Particular care was taken to increase the repertoire of food and respiratory allergens. One version of the MeDALL allergen-chip contained an extensive

collection of bacterial superantigens whereas in another version an extended panel of food allergens was added to study hypersensitivity to wheat. For each protein/allergen, 50–200 fg, corresponding to 1–5 attomol, was spotted in triplicates onto the chip. Customized spotting of the MeDALL-array was done by Phadia Austria GmbH (Vienna, Austria). The newly added allergens are described in the following section.

#### 2.2. Characterization of newly added allergens

**2.2.1. Peanut allergens**—Natural Ara h 1, nAra h 2 and nAra h 6 were purified as described [15,16].

**2.2.2. Almond, cashew and pistachio allergens**—nPru du 6 (prunin) was purified from aqueous almond extract by column chromatography [17] and natural Ana o 2 was purified from aqueous cashew nut protein extract [18] as previously described. Recombinant allergens were generated as previously described for Ana o 1 [19], Ana o 2 [20], Ana o 3 [21], Pru du 4 [22], Pru du 6 [23] and Pis v 3 [24].

Recombinant Pru du 3 was amplified directly from the almond cDNA library (library generation described in detail [23]). The sequence for the lipid transfer protein from almond (*Prunus dulcis*) was available in the NCBI database, accession number CAA65477. The LTP-encoding cDNA was ligated into the pMAL-c4X expression vector, transformed into Rosetta gamiB(DE3)pLysS cells (Novagen/EMD Millipore, Billerica, MA, USA), expressed, and purified as described in detail for Ana o 1 [19].

**2.2.3. Cow's milk allergens**—nBos d 4, nBos d 5, nBos d 8, recombinant  $\alpha$ S1-casein, r $\alpha$ S2-casein, r $\kappa$ -casein, r $\beta$ -casein, transferrin and natural BSA were purchased from Sigma Aldrich (St Louis, MO, USA) or expressed in *Escherichia coli* as described [25-27].

**2.2.4. Wheat allergens**—A wheat cDNA library was screened with sera from wheat food allergic patients as described [28]. Recombinant Tri a 36 and Tri a 36 fragment (amino acids 191-369) [29], clones 37 (thioredoxin h), 38 (glutathione transferase), 112 (1-Cysperoxiredoxin), 123 (profilin) and 126 (dehydrin) as well as clone 10 (serine proteinase inhibitor) were cloned, expressed and characterized as described [28,30]. cDNA-sequences of IgE-reactive clones 79, 85 and 110 that showed partial sequence identities to gamma gliadins, respectively, cDNAs coding for Tri a 37, LTP, beta amylases and avenin-like protein were codon optimized for E. coli expression and cloned with an additional sequence coding for a hexa-histidine tag into the expression vector pET17b (Novagen). Proteins were expressed in E. coli BL21 (DE3) (Stratagene, La Jolla, CA, USA), purified by nickel affinity chromatography (Qiagen, Hilden, Germany) and refolded by dialysis against 10 mM NaH<sub>2</sub>PO<sub>4</sub>, pH 4 (all proteins except LTP, pH 6, and beta amylases, pH 7). Using the same methodology, cDNAs of clones m43 and m82 that showed partial sequence identities to the high molecular weight glutenin x-type subunit Bx7 precursor were cloned into expression vector pMAL-c4X. Proteins were purified and dialyzed against 10 mM Na-Acetate, pH 4 (m43 and m82) as described above.

All purified proteins were tested for IgE-reactivity by non-denaturing, RAST-based IgE dotblotting.

**2.2.5. Antigens important for celiac disease**—The cDNA sequence of gamma gliadin 1 (rGG1) was cloned into pET27b vector (Novagen). The rGG1-protein was expressed with a hexa-histidine tag in *E. coli* BL21 (DE3) (Stratagene) and purified by nickel affinity chromatography (Qiagen). Purified proteins were dialyzed against 50 mM acetic acid and stored at -20 °C (manuscript in preparation). rGG1 derived peptides 4 and 18 were synthesized and purified as described [31]. Tissue transglutaminase protein (expressed in Baculovirus-infected Sf9 insect cells) was purchased from DIARECT AG (Freiburg, Germany).

**2.2.6. Olive pollen allergens**—cDNAs coding for Ole e 5, Ole e 6, Ole e 8 and Ole e 10 were codon-optimized for expression in *E. coli* and ligated into the pET17b vector (Novagen) with a C-terminal hexa-histidine tag. Proteins were expressed in BL21 (DE3) cells (Stratagene), purified by nickel affinity chromatography under native conditions and stored in 10 mM NaH<sub>2</sub>PO<sub>4</sub> (pH 8). IgE-reactivity was assessed by ELISA using sera from patients allergic to olive pollen.

**2.2.7. Mite allergens**—Recombinant Der p 5 (Genbank accession number, S76337), rDer p 7 (U37044), rDer p 21 (DQ354124), rDer p 23 (EU414751) and the clone 16-derived allergen (Vrtala, unpublished) were expressed in *E. coli* as soluble, non-fusion proteins using the pET17 expression system. rDer p 5, rDer p 21 and rDer p 23 were purified as previously described [32-34]. Recombinant Der p 7 was purified by hydrophobic interaction chromatography and hydroxyapatite chromatography while the clone 16-derived protein was purified by anion and cation exchange chromatography. rDer p 4 (AF144060), rDer p 11 (AY189697), rDer p 14 (AF373221), rDer p 15 (DQ078740) and rDer p 18 (DQ078739) were expressed as hexa-histidine-tagged proteins in the inclusion body fraction of *E. coli* extracts and purified by nickel affinity chromatography.

IgE-reactivity was tested for all recombinant HDM allergens with sera from HDM-allergic patients in dot-blot experiments and the allergenic activity was evaluated with basophil activation (CD203c up-regulation) experiments.

**2.2.8. Dog allergens**—Genes coding for the mature forms of Can f 4 (ACY38525.1), Can f 5 (P09582) and Can f 6 (E2QYS2) were inserted into expression vector pET-27b (Novagen). All three genes contained sequences coding for a hexa-histidine tag at the C terminus of the protein and the gene sequences were optimized for *E. coli* expression [35]. Proteins were expressed in *E. coli* BL21 (DE3) (Stratagene), purified by nickel affinity chromatography (Qiagen) and refolded by dialysis against 10 mM Tris–HCl (pH 9), followed by dialysis against 50 mM bicarbonate buffer (pH 9.6). Protein fold was determined by circular dichroism on a JASCO (Tokyo, Japan) J-810 spectropolarimeter.

**2.2.9. Insect venom allergens**—Recombinant Api m 1, rApi m 2 and rVes v 5 were cloned as described [36] and expressed in *E. coli*. Melittin (Api m 4) was purchased from Sigma Aldrich (St Louis, MO, USA). The cDNA coding for Ves v 1 with a 3' sequence coding for a hexa-histidine tag was codon-optimized for the expression in *E. coli* and inserted into plasmid pET17b (Novagen). Recombinant Ves v 1 was expressed in *E. coli* BL21 (DE3) cells (Stratagene) and purified from inclusion bodies by nickel affinity

chromatography (Qiagen). After dialysis, IgE-reactivity of rVes v 1 was assessed by dot-blot assays using sera from wasp-venom allergic patients.

**2.2.10. Staphylococcus aureus toxins**—Sequences encoding *S. aureus* toxins were ligated into expression vector pPR-IBA1 (IBA, Göttingen, Germany) and expressed with a C-terminal Strep-Tag in *E. coli* strain BL21. Cells were lysed by ultrasound and proteins were purified by affinity chromatography using Strep-Tactin columns (IBA). Purified proteins were stored in  $0.5 \times$  PBS.

**2.2.11. Maltose binding protein (MBP)**—The sequence encoding MBP was ligated into pMAL-c4x vector (New England Biolabs, Ipswich, MA, USA), transfected and expressed in *E. coli*.

#### 2.3. Calibration of the MeDALL allergen-chip for different isotypes

Serum pools were analyzed for reactivity to several different allergens on the microarray which were also available as ImmunoCAP tests in order to cover a broad spectrum of antibody specificities and fluorescence intensities. The respective allergen-specific antibodies of the isotype for which the calibration curve was generated were quantified by ImmunoCAP and calibration-curves – fluorescence intensities on the MeDALL allergenchip versus antibody levels as detected by ImmunoCAP – were established.

#### 2.4. Analysis of serum samples on the MeDALL allergen-chip

The microarrays were washed in washing buffer (Phadia AB) for 1 min by stirring. After drying the chips by centrifugation at 1000g for 1 min, aliquots of 35 µl per serum sample that have been centrifuged for one minute were applied and incubated for 2 h at gentle rocking. For detection of IgE, sera were applied undiluted, for detection of IgG and IgG<sub>4</sub>, samples were diluted 1:50 in sample diluent (Phadia AB), the dilution factor for IgA-detection was 1:10. After rinsing the arrays using a spray bottle, another washing step was done as described above, followed by drying by centrifugation. Thirty-five µl of detection antibody were then applied and incubated for 30 min at gentle rocking. After rinsing, washing and drying the chips as described above, the arrays were scanned using a confocal laser scanner (LuxScan-10 K microarray scanner, Capital-Bio, Beijing, People's Republic of China) and evaluated by Microarray Image Analyzer v3.1.2 software (Phadia AB). All incubation steps were performed at room temperature. For calibration and detection of background signals a calibrator serum and sample diluent (Phadia AB), respectively, were included in each run of serum analysis.

# 2.5. Sensitivity of the microarray as determined with monoclonal Bet v 1-specific human IgE and effects of blocking IgG on IgE-binding

Bip1, a Bet v 1-specific monoclonal mouse  $IgG_1$  antibody [37], and chimeric Bip1, a chimeric monoclonal antibody that consists of the variable domains of Bip1 and a human Fce-portion [38] have been described elsewhere.

First, a series of dilutions of chimeric Bip1 was prepared in sample diluent, comprising concentrations of 640, 320, 32, 3.2, 1.6, 0.8, 0.4, 0.2, 0.1 and 0.05 ng/ml of antibody. Each

dilution was measured on the MeDALL allergen-chip in triplicates. To study the effects of blocking IgG on IgE-signals, 3 different concentrations of chimeric Bet v 1-specific IgE, i.e., 312.5, 6.25 and 0.8 ng/ml, were prepared and mixed only with buffer, a 1-, 10-, 100-, 1000-, 10,000-fold (for 6.25 and 0.8 ng/ml IgE) and a 100,000-fold (for 0.8 ng/ml IgE) amount of Bip1. In fact, Bip1 binds to the very same epitope on Bet v 1 with the same affinity as the chimeric IgE. Again, IgE binding of each mixture to Bet v 1 was measured in triplicates using the MeDALL allergen-chip.

Control experiments were performed by mixing Bet v 1-specific IgE at a concentration of 312.5 ng/ml with an excess of two monoclonal IgG<sub>1</sub> antibodies specific for two distinct Bet v 1-peptide epitopes (peptide 2 and peptide 6, respectively) [39], or with an excess of an antibody specific for an unrelated allergen (i.e., anti-Phl p 2 IgG<sub>1</sub>).

The same series of experiments was done with the ImmunoCAP system (Phadia AB) for detection of Bet v 1-specific IgE by testing a dilution series comprising 312.5, 39.1, 19.5, 7.8, 3.9, 2, 1, 0.5 and 0.2 ng/ml of chimeric IgE. Likewise, measurement of Bet v 1-specific chimeric IgE was performed at two different concentrations (312.5 and 20 ng/ml) in the presence of only buffer, 1-fold up to 1000-fold (both for 312.5 and 20 ng/ml of IgE) or 10,000-fold (only for 20 ng/ml IgE) amount of Bip1 using the Bet v 1-ImmunoCAP.

# 2.6. Sensitivity and specificity of IgG measurements exemplified with a monoclonal PhI p 2-specific human $IgG_1$

The human monoclonal Phl p 2-specific  $IgG_1$  was described[40,41]. Such as for IgE, a dilution series of 1 µg/ml, 500 ng/ml, 100, 10, 1, 0.4, 0.3, 0.2 and 0.1 ng/ml was prepared of the monoclonal Phl p 2-specific  $IgG_1$  and measured in triplicates using the MeDALL allergen-chip. An IgG-specific detection antibody conjugate was used that binds to each of the human IgG subclasses.

#### 2.7. Importance of using protein-containing diluents for highly diluted antibodies

In order to investigate the possible binding of highly diluted antibodies to plastic tubes and/or the surface of the chip affecting sensitivity, chimeric Bip1 was diluted to a concentration of 6.25 ng/ml by either using PBS without any proteins or blocking agents or sample diluent containing additives which block non-specific interactions. IgE reactivities of these samples were measured in triplicates on the MeDALL allergen-chip.

# 2.8. Examples of detecting early allergen-specific IgE responses in samples from children of two different birth cohorts

The usefulness of the MeDALL allergen-chip for the evaluation of the development of allergic sensitizations in children is exemplified for sera from two population based birth cohorts, BAMSE from Sweden [42] and ECA from Norway [43]. In the BAMSE cohort, serum samples were analyzed from the age 4, 8 and 16 years, in the ECA cohort at 10 and 16 years.

### 3. Results

# 3.1. The MeDALL allergen-chip contains a large number of clinically relevant allergen molecules

Table 1 provides an overview of the molecules contained on the MeDALL allergen-chip and information regarding allergen-sources, if proteins were purified from extracts (natural allergens) or produced in expression systems (recombinant), if they are glycosylated and to which protein family they belong. Two versions of the chip have been made, one which contains a collection of enterotoxins from *Staphylococcus aureus* involved in respiratory allergies [44] and another comprising an extended panel of antigens relevant for hypersensitivity to wheat [29,30,45-47]. Considering the poor quality of natural house dust mite allergen extracts, we aimed to implement a panel of purified house dust mite allergens that is as complete as possible [48-50]. Furthermore, several cockroach allergens were included [51,52]. The panel of pollen allergens contained allergens useful for the identification of the genuine sensitizing pollen species as well as cross-reactive allergens [53,54]. Furthermore, attempts were made to include relevant food allergens from, e.g., nuts, peanuts, milk and several other allergen sources [55-59]. Column "CCD" (cross-reactive carbohydrate determinants) in Table 1 indicates glycosylated allergen molecules with confirmed reactivity to CCD-specific IgE.

#### 3.2. Technical features of the MeDALL allergen-chip

**3.2.1. Calibration for IgE, IgG, IgG<sub>4</sub> and IgA**—Calibration for IgE, IgG, IgG<sub>4</sub> and IgA has been established for the MeDALL allergen-chip by relating antibody levels to several allergens detected by microarray with the respective results from ImmunoCAP measurements. Screenshots showing the detection of allergen-specific IgE and IgG in the same serum sample illustrate the specificity of the detection systems as demonstrated by distinct reactivity patterns (Fig. 1).

**3.2.2.** Background, cut-off and inter-assay variation of the MeDALL allergenchip for IgE in comparison to ImmunoCAP, and for IgG-detection—Triplicates of a serial dilution of a monoclonal human IgE antibody specific to Bet v 1 were measured on the MeDALL allergen-chip, starting from 640 ng/ml down to 0.05 ng/ml (Fig. 2A–C). An IgE-concentration of 2.42 ng/ml corresponds approximately to 1UA/ml in the ImmunoCAP [60]. Up to a concentration of 32 ng/ml of the mAb (23.3ISU), a linear correlation between Ab-concentration in ng/ml and the signal displayed in ISU was found (Fig. 2B and C). At concentrations higher than 30 ng/ml the curve begins to approach saturation (Fig. 2A). However, high IgE-levels were measured with high reproducibility (mean signal at 640 ng/ml = 132.02ISU, SD = 2.26ISU, CV = 1.7%, Table 2). The mean signal detected at the lowest IgE concentration (0.05 ng/ml) was 0.14ISU (SD = 0.01ISU) with a signal-to-noise ratio of 14. The mean background signal detected with sample diluent alone that was used for dilution of the monoclonal IgE was 0.1ISU (SD = 0.01ISU). Table 2 summarizes results for each of the applied concentrations. To determine background signals for all allergens on the chip, sample diluent was applied as negative control, yielding relevant background signals only for the bacterial toxin Plc (Table 1), i.e., 0.22ISU (SD = 0.28ISU) calculated after 32 measurements of sample diluent.

For comparison, a dilution series of the same monoclonal human IgE specific to Bet v 1 was measured by ImmunoCAP (Fig. 2D-F, Table 2). IgE-dilutions from 0.2–312.5 ng/ml were tested in triplicates. A linear correlation between IgE concentration in ng/ml and UA/ml was demonstrated up to the highest IgE concentration (312.5 ng/ml) applied (Fig. 2D). At levels higher than 100UA/ml the calibration-curve begins to approach saturation-level (calibration-curve not shown). The mean background signal assessed using sample diluent was 0.09UA/ml (SD = 0.01UA/ml). Signal-to-noise ratios and %CV for all dilutions applied are shown in Table 2.

A dilution series of a monoclonal human  $IgG_1$  antibody specific to the timothy grass pollen allergen Phl p 2 with a concentration range from 0.1 ng/ml to 1 µg/ml was also measured in triplicates on the MeDALL allergen-chip (Fig. 3A–C, Table 2). From 0 to 10 ng/ml of the mAb a linear correlation between Ab-concentration in ng/ml and in ISU-G was found (Fig. 3B and C). At higher concentrations, the curve approached saturation-level. At the lowest concentration measured (i.e., 0.1 ng/ml) the signal-to-noise ratio was 4, indicating that even very low IgG-concentrations can be reliably measured by the MeDALL allergen-chip. The mean background signal was 0.02ISU-G (SD = 0.01ISU-G). Results for all dilutions tested are shown in Table 2.

Next, mean variation coefficients and signal-to-noise ratios were calculated for results from triplicate-measurements of IgE- and (only for the microarray) IgG-dilutions up to 1 ng/ml or greater than 1 ng/ml (Table 2). Microarrays and ImmunoCAPs from the same lot were used for testing antibody dilutions. The results for IgE-detection using the MeDALL allergenchip show higher CVs for low concentrations (mean CV for concentrations up to 1 ng/ml = 14.1%, SD = 6.8%) than for higher concentrations (mean CV = 7.7%, SD = 5.5%), but the mean signal-to-noise ratio for low concentrations was sufficient (mean SNR = 8.5; SD = 3.9) to obtain reliable results even for very low IgE concentrations such as 0.05 ng/ml. For IgG detection by the MeDALL allergen-chip, both for high and for low antibody-concentrations, mean %CVs of 14.9% or 16.3% (SD = 7.9% or 10.4%), respectively, were calculated. Signal-to-noise ratios for the two intervals were 8.1 or 21, hence, also for IgG reliable results were obtained both for high and for low antibody concentrations.

Detection of the monoclonal Bet v 1-specific IgE-antibody was also performed by ImmunoCAP for concentrations from 0.2 to 312.5 ng/ml. At the lowest concentration tested (0.2 ng/ml), signals (0.16 UA/ml) were already close to the background level (0.09 UA/ml) whereas with the MeDALL allergen-chip IgE concentrations down to 0.05 ng/ml could be reliably detected. However, in ImmunoCAP, variation coefficients of 2.7% and 2.3% (SD = 1.4% and 2%) were calculated for high (i.e., greater than 1 ng/ml) and for low (i.e., up to 1 ng/ml) IgE-concentrations. Likewise, mean signal-to-noise ratios were 50.1 (SD = 33.6) for concentrations higher than 1 ng/ml and 29 (SD = 2.7) for concentrations up to 1 ng/ml (Table 2, right part).

3.2.3. Simultaneous presence of blocking IgG antibodies inhibits IgE binding to allergens on the chip but not in the ImmunoCAP—In order to test the influence of IgG antibodies on IgE binding to allergens in chip-based and ImmunoCAP assays, we used a monoclonal IgE antibody specific for the major birch pollen allergen Bet v 1 and an IgG antibody recognizing the very same epitope with similar affinity. This experimental set up allowed for a comparison of the effects of blocking IgG on allergen-specific IgE binding for the two technologies. For the MeDALL allergen-chip, three different concentrations of monoclonal Bet v 1-specific human IgE (312.5 ng/ml, 6.25 ng/ml, 0.8 ng/ml) were mixed with the monoclonal  $IgG_1$  (Bip1) that contains the same variable domains as the IgEantibody and hence binds to the same epitope of Bet v 1. The respective IgE-dilutions were incubated with increasing concentrations of Bip1, ranging from a 1:1 ratio of IgE to  $IgG_1$  up to a 1000-fold (for 312.5 ng/ml IgE), 10,000-fold (for 6.25 ng/ml IgE) or a 100,000-fold (for 0.8 ng/ml IgE) excess of Bip1 as compared to the IgE antibody (Fig. 4A–C). At the highest IgE-concentration (312.5 ng/ml) a strong inhibition of the signal was found already at a 1:1 ratio of IgE and IgG (mean inhibition 43.8%, SD = 6.6%), increasing to an almost complete blocking of IgE-binding at a 100-fold and 1000-fold excess of Bip1 (Fig. 4A, Table 3). For lower IgE-concentrations we observed up to a 10-fold excess of Bip1 only a moderate (18.2% for 6.25 ng/ml IgE) or no (for 0.8 ng/ml) decrease of IgE binding (Fig. 4B and C). From a 100-fold excess of competitive IgG on, strong inhibition of the IgE-signal was observed (-71% for 6.25 ng/ml IgE, -28.3% for 0.8 ng/ml IgE) that increased to almost 100% inhibition at 1000- or 10,000-fold excess of Bip1, respectively. However, for the complete range of IgE concentrations applied, Bet v 1-specific IgE could be measured at more than 100-fold excess of competitive IgG (Table 3).

When control experiments were performed by incubating Bet v 1-specific IgE with an excess of an antibody specific for an unrelated allergen (i.e., anti-Phl p 2 IgG<sub>1</sub>) no inhibition of IgE binding was found, whereas mixing with an excess of two monoclonal IgG<sub>1</sub> antibodies specific for two distinct Bet v 1-peptide epitopes (peptide 2 and peptide 6, respectively) gave different levels of inhibition (Fig. 5). Monoclonal IgG<sub>1</sub> specific to Bet v 1-peptide 2 caused 75.5% inhibition, whereas inhibition by the antibody specific for peptide 6 was 47.1% (Fig. 5).

The same inhibitions were performed using the ImmunoCAP system by mixing two different IgE-dilutions (312.5 ng/ml, 20 ng/ml) with up to a 1000-fold (for 312.5 ng/ml IgE) or 10,000-fold (for 20 ng/ml IgE) excess of Bip1. Levels of Bet v 1-specific IgE were then measured using the Bet v 1 ImmunoCAP (Fig. 4 D-E). At high IgE concentrations (312.5 ng/ml), up to a 100-fold excess of Bip1 did not result in detectable inhibition and at low concentrations (20 ng/ml), even 1000-fold excess did not inhibit. At the highest IgG to IgE-ratios applied, an inhibition of 76% (for 312.5 ng/ml IgE) or 57% (for 20 ng/ml IgE) of IgE binding was observed (Table 3). An overlay plot of relative inhibitions of IgE binding as detected by the MeDALL allergen-chip and by ImmunoCAP is shown in Fig. 4F. The levels of inhibition are displayed in Table 3 and show that the MeDALL allergen-chip is much more sensitive (>100-fold) in detecting blocking IgG antibodies than the ImmunoCAP.

Fig. 6 illustrates an important detail when it comes to measuring low concentrations of highly diluted antibodies. Purified monoclonal human Bet v 1-specific IgE was diluted in a

sample diluent which contains additives that block non-specific interaction of the antibody with the surfaces of the test tube and of the microarray, or in PBS, respectively, to a concentration of 6.25 ng/ml. Then, Bet v 1-specific IgE levels were measured by the MeDALL allergen-chip. Using sample diluent containing blocking reagents, a mean IgE-level of 5.66ISU (SD = 0.62ISU) was measured, whereas using PBS without blocking reagents, 2.48ISU (SD = 0.19ISU), i.e., a 56.1% lower level was measured. This result shows that non-specific binding to plastic or other surfaces decreases the availability of antibodies for binding to their specific target. Such adsorption effects may be misinterpreted as instability of the purified reactant and can influence results in all types of assays where highly diluted purified reactants are used (e.g., skin prick testing, cellular assays).

#### 3.2.4. Examples of different progression patterns of allergic sensitization-

Table 4 shows examples for the development of IgE reactivities during early childhood as detected in serum samples from children of the BAMSE cohort from Sweden [42] at the age of 4, 8 and 16 years. Displayed are patterns of IgE responses towards respiratory and food allergens which are detectable early in life, i.e., before the age of 4 years (BAMSE 1), IgE responses which become detectable between 4 and 8 years (BAMSE 2), later, i.e., between 8 and 16 years (BAMSE 3 and 4) and mixed patterns of IgE responses appearing early and late in the same child (BAMSE 5).

Screenshots of IgE detection by microarray in serum samples obtained from one study participant of the Norwegian birth cohort ECA [43] at the age of 10 and 16 years are shown in Fig. 7. At 16 years, the serum sample exhibits reactivities to some additional allergen components, mainly from the same allergen source and IgE levels were higher as illustrated by "warmer" colours (yellow to red).

# 4. Discussion and conclusions

In the present study we demonstrated recent advances in allergen microarray technology using the MeDALL allergen-chip as a model. The MeDALL allergen-chip is a more developed version of the ImmunoCAP ISAC platform in that it contains an expanded repertoire of allergens that can be used for the detection of allergen-specific IgE as well as other immunoglobulin isotypes. The MeDALL allergen-chip has been developed as a research tool for the sensitive detection of IgE-reactivities in samples from children of different birth cohorts of various European countries to study the evolution of the allergenspecific antibody responses in childhood. The MeDALL allergen-chip allows the sensitive detection of allergen-specific IgE towards a large number of clinically relevant allergen molecules as has been demonstrated by IgE titration experiments and comparison with one of the most sensitive tools for measuring allergen-specific IgE, the ImmunoCAP assay. To establish a large panel of clinically relevant allergens on the chip we assembled purified recombinant and natural allergen molecules from different respiratory and food allergen sources, such as house dust mite, peanut and wheat, to complement the allergen panel of ImmunoCAP ISAC (Table 1). Certain molecules were spotted as recombinant, others as natural proteins purified from allergen extracts. In addition to IgE, calibration for IgG, IgG<sub>4</sub> and IgA has been established for the MeDALL allergen-chip and specific detection of other isotypes is being developed. This allows us not only to study IgE reactivities in children's

sera but should also assess the eventual protective effects of maternal allergen-specific IgG transferred to the child or of spontaneously developing IgG responses in study subjects.

A series of experiments was performed using well defined monoclonal antibodies with known specificities to different allergens. In some of these experiments measurements were done with the MeDALL allergen-chip and by ImmunoCAP in parallel in order to study characteristics and special features of these two systems. The ImmunoCAP system was chosen as a reference because it has been widely used by clinicians for measuring allergenspecific IgE responses in the birth cohorts. Both, ImmunoCAP and the MeDALL allergenchip were found to have almost equal background signals for IgE detection, which are in the range of 0.1 (SD = 0.01) ISU or UA/ml, respectively. At the lowest IgE-concentrations applied, i.e., 0.05 ng/ml for the MeDALL allergen-chip and 0.2 ng/ml for ImmunoCAP, signal-to-noise ratios were 14 and 27.1, relative variation coefficients were 7.1% and 3.7%, respectively. These results imply that with both technologies, low IgE-concentrations can be reliably measured, and, with the MeDALL allergen-chip, even at higher sensitivity (Table 2). When detecting IgG using the microarray, almost no relevant background signal (0.02ISU-G, SD = 0.01ISU-G) and a sufficiently high SNR (4 for 0.1 ng/ml IgG) were observed, allowing detection of IgG-concentrations even lower than the smallest concentration tested in the present series. Likewise, when looking at higher concentrations of IgE or IgG, low inter-assay variations were observed (Table 2). However, Fig. 2A demonstrates that at IgE concentrations higher than approximately 30 ng/ml, the linear relation between antibody concentration and ISU is no longer observed. Therefore, IgE levels higher than 30ISU measured with the chip correspond to increasingly higher concentrations of IgE in ng/ml (Fig. 2A). In the ImmunoCAP there is still a linear relation between antibody concentration in ng/ml and UA/ml up to 300 ng/ml of IgE (Fig. 2D). This difference between the two systems is due to the fact that on the microarray, only 50-200 fg of protein are immobilized per spot whereas in ImmunoCAP,  $1-2 \mu g$  of protein, i.e., 10,000,000 times more, are coupled. Consequently, the quality of the protein preparations used and the assessment of antibody-binding capacity to the allergen spots are critically important issues for microarray technology since signal intensities decrease in relation to the proportion of spotted allergen molecules that are denatured.

Apart from the detection of IgE reactivities to defined allergen molecules for diagnosis of allergic sensitization, the measurement of blocking IgG, or antibodies of any other isotype than IgE, is an important aspect, e.g., for the assessment of the immunological efficacy of specific immunotherapy or of potentially protective effects of maternal IgG that are transferred to the child via the placenta or through breast milk or of spontaneously developing allergen-specific IgG responses. Quantification of IgG that binds to a given allergen does not provide information as to whether it inhibits IgE binding to the same allergen. In order to study how different systems for measuring allergen-specific IgE (i.e., microarray versus ImmunoCAP) can provide information regarding the presence of blocking IgG, IgE responses were detected before and after addition of competing IgG. Using ImmunoCAP, at 312.5 ng/ml of Bet v 1-specific IgE, no effect was observed on detected allergen-specific IgE-levels with up to a 100-fold excess of blocking IgG (Fig. 4D). This is due to the fact that a relative excess of allergen is used and therefore, all antibody isotypes specific for Bet v 1 can bind independently from each other. Therefore, the potential

therapeutic effect of allergen-specific IgG cannot be judged by using ImmunoCAP, neither by measuring IgE before and after induction of IgG, nor by direct assessment of IgG. Performing the same experiment using the MeDALL allergen-chip, a 1:1 ratio of IgE (312.5 ng/ml) and competing IgG showed a decrease of the IgE-signal of 44% (Fig. 4A, Table 3). This means that the allergen microarray integrates all antibodies binding to the same epitope into a net IgE value (or value of any other isotype detected). Since blocking IgE from binding to allergen reduces the allergic reaction, it may be assumed that IgE levels measured by microarray should better reflect the potency of IgE to mediate an allergic reaction upon allergen contact. However, the allergen-specific IgE signal is completely abolished only at a 1000-fold or higher excess of blocking antibodies, depending on the absolute allergenspecific IgE level. Relative inhibition of IgE signals detected for all IgE-concentrations investigated and for both CAP and MeDALL allergen-chip are plotted in Fig. 4F showing that at an IgE-concentration of 0.8 ng/ml (microarray) or 312.5 ng/ml (ImmunoCAP), similar findings are obtained for both systems. Moreover, the figure illustrates that at higher IgE concentrations, IgE levels measured are more strongly decreased by the presence of blocking IgG (regarding the relative excess of IgG) than at lower concentrations. This observation becomes clear when considering that, assuming a given number of epitopes available for antibody-binding, at high antibody-concentrations, the saturation-level is approached, especially when the number of epitopes is restricted as in the microarray.

The results obtained in this study illustrate that there is a fundamental difference between allergen microarrays and the ImmunoCAP system for detecting and measuring allergen-specific IgE. In the traditionally used system for detection of allergen-specific IgE, the ImmunoCAP system, there is a relative excess of allergen molecules  $(1-2 \mu g)$  of the respective allergen molecule) as compared to the average concentrations of allergen-specific IgE antibodies in serum. Since there is a large amount of allergen coupled to a large surface, antibodies of all isotypes can bind simultaneously. By contrast, allergen microarrays contain only approximately 5 attomol of protein per spot, corresponding to, e.g., 100 fg in case of a protein of 20 kDa. This means that, compared to ImmunoCAP, one spot of the microarray contains 10,000,000 times less protein.

This difference in immobilized allergen amounts between the two test systems has an important impact on the test results and must be considered for their interpretation. As long as antibody levels are relatively low, both systems yield comparable results but as soon as allergen-specific antibody levels rise, the relatively limited number of epitopes present per spot on the microarray prevents the binding of some of the antibodies that are specific to the respective allergen. As a consequence, in the microarray saturation levels are already reached at lower antibody concentrations which may get close to the *in vivo* situation of allergic patients where only a part of the circulating IgE is bound to the effector cells of the allergic reaction (e.g., mast cells, basophils, antigen presenting cells) via their specific receptors. This assumption needs to be confirmed in clinical studies comparing chip results with *in vivo* sensitivity in patients (e.g., titrated skin tests). Importantly, in case allergen-specific antibodies other than IgE are simultaneously present in serum samples or other biological fluids (e.g., tears, breast milk), competition between IgE and other isotypes will become detectable as a reduction of IgE binding when these isotypes represent true blocking

antibodies recognizing the same epitope or affect IgE binding by other mechanisms such as steric hindrance. On the microarray, like in the *in vivo* situation where only small amounts of allergen enter the body, several parameters will determine which isotype will bind better to the allergen, among them concentration and affinity of the respective antibodies. Regarding allergy diagnosis, it is thus reasonable to assume that the microarray will more closely reflect the biological situation in the patient after incorporation of allergen. Likewise, for successful specific immunotherapy (SIT), higher titers of allergen-specific IgG, rather than IgE antibodies will be induced, and therefore IgG will occupy more epitopes on the allergen molecule than the corresponding IgE. If allergen-specific antibodies in such patients are measured by microarray, one will observe increases of allergen-specific IgG levels and, in those cases where these allergen-specific IgG antibodies can block IgE binding to the allergen, a reduction of IgE binding will be detectable (Lupinek & Valenta, unpublished). Therefore, allergen-microarrays may be useful for the measurement of the successful induction of therapeutically active IgG by SIT and should eventually become a new surrogate marker for therapeutic success of SIT, provided that clinical studies confirm the theoretical assumption based on the *in vitro* results. Test systems such as the ImmunoCAP, which, due to the presence of large amounts of allergen on the solid phase, can accommodate the simultaneous binding of IgE and other isotypes and therefore do not mimic the biological situation (i.e., clinical sensitivity of the patient upon allergen encounter). Yet, they offer the possibility to quantify allergen-specific antibodies which allows, e.g., IgE quantification in order to assess the effect of allergen avoidance or therapeutic IgE depletion on allergen-specific IgE levels.

#### Acknowledgments

#### Funding sources

This work was supported by the European Commission's Seventh Framework Programme under grant agreement No. 261357 and in part by project F4605 of the Austrian Science Fund (FWF) and the Christian Doppler Research Association, Austria.

### Abbreviations

CV	variation coefficient
HDM	house dust mite
IgE	immunoglobulin E
IgG	immunoglobulin G
ISU	ISAC standardized units (for IgE-detection)
ISU-G	ISAC standardized units (for IgG-detection)
MeDALL	mechanisms for the development of allergies
SNR	signal-to-noise ratio
UA/ml	units antigen per milliliter

## References

- Wuthrich B, Schmid-Grendelmeier P, Schindler C, Imboden M, Bircher A, Zemp E, Probst-Hensch N. Int. Arch. Allergy Immunol. 2013; 162:143–148. [PubMed: 23921456]
- [2]. Kulig M, Bergmann R, Klettke U, Wahn V, Tacke U, Wahn U. J. Allergy Clin. Immunol. 1999; 103:1173–1179. [PubMed: 10359902]
- [3]. Lupinek C, Marth K, Niederberger V, Valenta R. J. Allergy Clin. Immunol. 2012; 130:1418–1420. e1414. [PubMed: 22867692]
- [4]. W.H.O.C.C.f. Asthma, Rhinitis. Bousquet J, Anto JM, Demoly P, Schunemann HJ, Togias A, Akdis M, Auffray C, Bachert C, Bieber T, Bousquet PJ, Carlsen KH, Casale TB, Cruz AA, Keil T, Lødrup Carlsen KC, Maurer M, Ohta K, Papadopoulos NG, Roman Rodriguez M, Samolinski B, Agache I, Andrianarisoa A, Ang CS, Annesi-Maesano I, Ballester F, Baena-Cagnani CE, Basagana X, Bateman ED, Bel EH, Bedbrook A, Beghe B, Beji M, Ben Kheder A, Benet M, Bennoor KS, Bergmann KC, Berrissoul F, Bindslev Jensen C, Bleecker ER, Bonini S, Boner AL, Boulet LP, Brightling CE, Brozek JL, Bush A, Busse WW, Camargos PA, Canonica GW, Carr W, Cesario A, Chen YZ, Chiriac AM, Costa DJ, Cox L, Custovic A, Dahl R, Darsow U, Didi T, Dolen WK, Douagui H, Dubakiene R, El-Meziane A, Fonseca JA, Fokkens WJ, Fthenou E, Gamkrelidze A, Garcia-Aymerich J, Gerth van Wijk R, Gimeno-Santos E, Guerra S, Haahtela T, Haddad H, Hellings PW, Hellquist-Dahl B, Hohmann C, Howarth P, Hourihane JO, Humbert M, Jacquemin B, Just J, Kalayci O, Kaliner MA, Kauffmann F, Kerkhof M, Khayat G, Koffi N'Goran B, Kogevinas M, Koppelman GH, Kowalski ML, Kull I, Kuna P, Larenas D, Lavi I, Le LT, Lieberman P, Lipworth B, Mahboub B, Makela MJ, Martin F, Martinez FD, Marshall GD, Mazon A, Melen E, Meltzer EO, Mihaltan F, Mohammad Y, Mohammadi A, Momas I, Morais-Almeida M, Mullol J, Muraro A, Naclerio R, Nafti S, Namazova-Baranova L, Nawijn MC, Nyembue TD, Oddie S, O'Hehir RE, Okamoto Y, Orru MP, Ozdemir C, Ouedraogo GS, Palkonen S, Panzner P, Passalacqua G, Pawankar R, Pigearias B, Pin I, Pinart M, Pison C, Popov TA, Porta D, Postma DS, Price D, Rabe KF, Ratomaharo J, Reitamo S, Rezagui D, Ring J, Roberts R, Roca J, Rogala B, Romano A, Rosado-Pinto J, Ryan D, Sanchez-Borges M, Scadding GK, Sheikh A, Simons FE, Siroux V, Schmid-Grendelmeier PD, Smit HA, Sooronbaev T, Stein RT, Sterk PJ, Sunyer J, Terreehorst I, Toskala E, Tremblay Y, Valenta R, Valeyre D, Vandenplas O, van Weel C, Vassilaki M, Varraso R, Viegi G, Wang DY, Wickman M, Williams D, Wohrl S, Wright J, Yorgancioglu A, Yusuf OM, Zar HJ, Zernotti ME, Zidarn M, Zhong N, Zuberbier T. Severe chronic allergic (and related) diseases: a uniform approach-a MeDALL-GA2LEN-ARIA position paper. Int. Arch. Allergy Immunol. 2012; 158:216–231. [PubMed: 22382913]
- [5]. Bousquet J, Anto J, Sunyer J, Nieuwenhuijsen M, Vrijheid M, Keil T, Me DSG, C.S. Group; E.S. Group; G.S. Group. Int. Arch. Allergy Immunol. 2013; 161:1–10. [PubMed: 23258290]
- [6]. Hiller R, Laffer S, Harwanegg C, Huber M, Schmidt WM, Twardosz A, Barletta B, Becker WM, Blaser K, Breiteneder H, Chapman M, Crameri R, Duchene M, Ferreira F, Fiebig H, Hoffmann-Sommergruber K, King TP, Kleber-Janke T, Kurup VP, Lehrer SB, Lidholm J, Muller U, Pini C, Reese G, Scheiner O, Scheynius A, Shen HD, Spitzauer S, Suck R, Swoboda I, Thomas W, Tinghino R, Van Hage-Hamsten M, Virtanen T, Kraft D, Muller MW, Valenta R. FASEB J. 2002; 16:414–416. [PubMed: 11790727]
- [7]. Patelis A, Gunnbjornsdottir M, Malinovschi A, Matsson P, Onell A, Hogman M, Alving K, Janson C. J. Allergy Clin. Immunol. 2012; 130:397–402. e392. [PubMed: 22633327]
- [8]. Onell A, Hjalle L, Borres MP. Clin. Transl. Allergy. 2012; 2:24. [PubMed: 23254184]
- [9]. Cabrera-Freitag P, Goikoetxea MJ, Beorlegui C, Gamboa P, Gastaminza G, Fernandez-Benitez M, Ferrer M, Blanca M, Sanz ML. Clin. Exp. Allergy. 2011; 41:1440–1446. [PubMed: 21749500]
- [10]. Canonica GW, Ansotegui IJ, Pawankar R, Schmid-Grendelmeier P, van Hage M, Baena-Cagnani CE, Melioli G, Nunes C, Passalacqua G, Rosenwasser L, Sampson H, Sastre J, Bousquet J, Zuberbier T, WAO-ARIA-GA2LEN Task Force. Allen K, Asero R, Bohle B, Cox L, de Blay F, Ebisawa M, Maximiliano-Gomez R, Gonzalez-Diaz S, Haahtela T, Holgate S, Jakob T, Larche M, Matricardi PM, Oppenheimer J, Poulsen LK, Renz HE, Rosario Ne. Rothenberg M, Sanchez-Borges M, Scala E, Valenta R. World Allergy Organ J. 2013; 6:17. [PubMed: 24090398]
- [11]. Javaloyes G, Goikoetxea MJ, Garcia Nunez I, Sanz ML, Blanca M, Scheurer S, Vieths S, Ferrer M. J. Investig. Allergol. Clin. Immunol. 2012; 22:508–513.

- [12]. Gadisseur R, Chapelle JP, Cavalier E. Clin. Chem. Lab. Med. 2011; 49:277–280. [PubMed: 21143018]
- [13]. Asarnoj A, Moverare R, Ostblom E, Poorafshar M, Lilja G, Hedlin G, van Hage M, Ahlstedt S, Wickman M. Allergy. 2010; 65:1189–1195. [PubMed: 20146729]
- [14]. Melioli G, Bonifazi F, Bonini S, Maggi E, Mussap M, Passalacqua G, Rossi ER, Vacca A, Canonica GW, I. Italian Board. Clin. Biochem. 2011; 44:1005–1011. [PubMed: 21627961]
- [15]. Marsh J, Rigby N, Wellner K, Reese G, Knulst A, Akkerdaas J, van Ree R, Radauer C, Lovegrove A, Sancho A, Mills C, Vieths S, Hoffmann-Sommergruber K, Shewry PR. Mol. Nut. Food Res. 2008; 52(Suppl. 2):S272–285.
- [16]. Koppelman SJ, de Jong GA, Laaper-Ertmann M, Peeters KA, Knulst AC, Hefle SL, Knol EF. Clin. Exp. Allergy. 2005; 35:490–497. [PubMed: 15836759]
- [17]. Sathe SK, Wolf WJ, Roux KH, Teuber SS, Venkatachalam M, Sze-Tao KW. J. Agric. Food Chem. 2002; 50:4333–4341. [PubMed: 12105967]
- [18]. Sathe SK, Sze-Tao KW, Wolf WJ, Hamaker BR. J. Agric. Food Chem. 1997; 45:2854–2860.
- [19]. Wang F, Robotham JM, Teuber SS, Tawde P, Sathe SK, Roux KH. J. Allergy Clin. Immunol. 2002; 110:160–166. [PubMed: 12110836]
- [20]. Wang F, Robotham JM, Teuber SS, Sathe SK, Roux KH. Int. Arch. Allergy Immunol. 2003; 132:27–39. [PubMed: 14555856]
- [21]. Robotham JM, Wang F, Seamon V, Teuber SS, Sathe SK, Sampson HA, Beyer K, Seavy M, Roux KH. J. Allergy Clin. Immunol. 2005; 115:1284–1290. [PubMed: 15940148]
- [22]. Tawde P, Venkatesh YP, Wang F, Teuber SS, Sathe SK, Roux KH. J. Allergy Clin. Immunol. 2006; 118:915–922. [PubMed: 17030246]
- [23]. Willison LN, Tripathi P, Sharma G, Teuber SS, Sathe SK, Roux KH. Int. Arch. Allergy Immunol. 2011; 156:267–281. [PubMed: 21720172]
- [24]. Willison LN, Tawde P, Robotham JM. Clin. Exp. Allergy. 2008; 38:1229–1238. [PubMed: 18479490]
- [25]. Hochwallner H, Schulmeister U, Swoboda I, Balic N, Geller B, Nystrand M, Harlin A, Thalhamer J, Scheiblhofer S, Niggemann B, Quirce S, Ebner C, Mari A, Pauli G, Herz U, van Tol EA, Valenta R, Spitzauer S. Clin. Exp. Allergy. 2010; 40:1809–1818. [PubMed: 20860558]
- [26]. Hochwallner H, Schulmeister U, Swoboda I, Twaroch TE, Vogelsang H, Kazemi-Shirazi L, Kundi M, Balic N, Quirce S, Rumpold H, Froschl R, Horak F, Tichatschek B, Stefanescu CL, Szepfalusi Z, Papadopoulos NG, Mari A, Ebner C, Pauli G, Valenta R, Spitzauer S. Allergy. 2011; 66:1201–1207. [PubMed: 21575008]
- [27]. Schulmeister U, Hochwallner H, Swoboda I, Focke-Tejkl M, Geller B, Nystrand M, Harlin A, Thalhamer J, Scheiblhofer S, Keller W, Niggemann B, Quirce S, Ebner C, Mari A, Pauli G, Herz U, Valenta R, Spitzauer S. J. Immunol. 2009; 182:7019–7029. [PubMed: 19454699]
- [28]. Constantin C, Quirce S, Grote M, Touraev A, Swoboda I, Stoecklinger A, Mari A, Thalhamer J, Heberle-Bors E, Valenta R. J. Immunol. 2008; 180:7451–7460. [PubMed: 18490745]
- [29]. Baar A, Pahr S, Constantin C, Scheiblhofer S, Thalhamer J, Giavi S, Papadopoulos NG, Ebner C, Mari A, Vrtala S, Valenta R. J. Immunol. 2012; 189:3018–3025. [PubMed: 22904302]
- [30]. Pahr S, Constantin C, Mari A, Scheiblhofer S, Thalhamer J, Ebner C, Vrtala S, Mittermann I, Valenta R. Clin. Exp. Allergy. 2012; 42:597–609. [PubMed: 22417217]
- [31]. Focke M, Linhart B, Hartl A, Wiedermann U, Sperr WR, Valent P, Thalhamer J, Kraft D, Valenta R. Clin. Exp. Allergy. 2004; 34:1525–1533. [PubMed: 15479266]
- [32]. Weghofer M, Dall'Antonia Y, Grote M, Stocklinger A, Kneidinger M, Balic N, Krauth MT, Fernandez-Caldas E, Thomas WR, van Hage M, Vieths S, Spitzauer S, Horak F, Svergun DI, Konarev PV, Valent P, Thalhamer J, Keller W, Valenta R, Vrtala S. Allergy. 2008; 63:758–767. [PubMed: 18445190]
- [33]. Weghofer M, Grote M, Resch Y, Casset A, Kneidinger M, Kopec J, Thomas WR, Fernandez-Caldas E, Kabesch M, Ferrara R, Mari A, Purohit A, Pauli G, Horak F, Keller W, Valent P, Valenta R, Vrtala S. J. Immunol. 2013; 190:3059–3067. [PubMed: 23460742]

- [34]. Weghofer M, Grote M, Dall'Antonia Y, Fernandez-Caldas E, Krauth MT, van Hage M, Horak F, Thomas WR, Valent P, Keller W, Valenta R, Vrtala S. Int. Arch. Allergy Immunol. 2008; 147:101–109. [PubMed: 18520154]
- [35]. Curin M, Swoboda I, Wollmann E, Lupinek C, Spitzauer S, van Hage M, Valenta R. Microarrayed dog, cat and horse allergens show weak correlation between allergen-specific IgE and IgG responses. J. Allergy Clin. Immunol. 2013
- [36]. Mittermann I, Zidarn M, Silar M, Markovic-Housley Z, Aberer W, Korosec P, Kosnik M, Valenta R. J. Allergy Clin. Immunol. 2010; 125:1300–1307. e1303. [PubMed: 20466415]
- [37]. Laffer S, Vangelista L, Steinberger P, Kraft D, Pastore A, Valenta R. J. Immunol. 1996; 157:4953–4962. [PubMed: 8943401]
- [38]. Laffer S, Hogbom E, Roux KH, Sperr WR, Valent P, Bankl HC, Vangelista L, Kricek F, Kraft D, Gronlund H, Valenta R. J. Allergy Clin. Immunol. 2001; 108:409–416. [PubMed: 11544461]
- [39]. Gieras A, Cejka P, Blatt K, Focke-Tejkl M, Linhart B, Flicker S, Stoecklinger A, Marth K, Drescher A, Thalhamer J, Valent P, Majdic O, Valenta R. J. Immunol. 2011; 186:5333–5344. [PubMed: 21451110]
- [40]. Flicker S, Steinberger P, Norderhaug L, Sperr WR, Majlesi Y, Valent P, Kraft D, Valenta R. Eur. J. Immunol. 2002; 32:2156–2162. [PubMed: 12209627]
- [41]. Padavattan S, Flicker S, Schirmer T, Madritsch C, Randow S, Reese G, Vieths S, Lupinek C, Ebner C, Valenta R, Markovic-Housley Z. J. Immunol. 2009; 182:2141–2151. [PubMed: 19201867]
- [42]. Wickman M, Kull I, Pershagen G, Nordvall SL. Pediatr. Allergy Immunol. 2002; 13(Suppl. 15): 11–13. [PubMed: 12688617]
- [43]. Hovland V, Riiser A, Mowinckel P, Carlsen KH, Lødrup KC. Eur. Respir. J. 2013; 41:838–845. [PubMed: 22903966]
- [44]. Huvenne W, Hellings PW, Bachert C. Int. Arch. Allergy Immunol. 2013; 161:304–314.[PubMed: 23689556]
- [45]. Pahr S, Constantin C, Papadopoulos NG, Giavi S, Makela M, Pelkonen A, Ebner C, Mari A, Scheiblhofer S, Thalhamer J, Kundi M, Vrtala S, Mittermann I, Valenta R. J. Allergy Clin. Immunol. 2013
- [46]. Srinivasan B, Focke-Tejkl M, Swoboda I, Constantin C, Mittermann I, Pahr S, Vogelsang H, Huber WD, Valenta R. Amino Acids. 2013
- [47]. Ebisawa M, Shibata R, Sato S, Borres MP, Ito K. Int. Arch. Allergy Immunol. 2012; 158:71–76. [PubMed: 22212744]
- [48]. Casset A, Mari A, Purohit A, Weghofer M, Ferrara R, Thomas WR, Alessandri C, Chen K-W, De Blay F, Valenta R, Vrtala S. Int. Arch. Allergy Immunol. 2012
- [49]. Vrtala S, Huber H, Thomas WR. Methods. 2013
- [50]. Hales BJ, Elliot CE, Chai LY, Pearce LJ, Tipayanon T, Hazell L, Stone S, Piboonpocanun S, Thomas WR, Smith WA. Int. Arch. Allergy Immunol. 2013; 160:233–240. [PubMed: 23075813]
- [51]. Barbosa MC, Santos AB, Ferriani VP, Pomes A, Chapman MD, Arruda LK. Int. Arch. Allergy Immunol. 2013; 161:213–219. [PubMed: 23549028]
- [52]. Pomes A, Arruda LK. Methods. 2013
- [53]. Douladiris N, Savvatianos S, Roumpedaki I, Skevaki C, Mitsias D, Papadopoulos NG. Int. Arch. Allergy Immunol. 2013; 162:163–172. [PubMed: 23921568]
- [54]. Villalba M, Rodriguez R, Batanero E. Methods. 2013
- [55]. Willison LN, Sathe SK, Roux KH. Methods. 2013
- [56]. Asarnoj A, Glaumann S, Elfstrom L, Lilja G, Lidholm J, Nilsson C, Wickman M. Int. Arch. Allergy Immunol. 2012; 159:209–212. [PubMed: 22677622]
- [57]. Berneder M, Bublin M, Hoffmann-Sommergruber K, Hawranek T, Lang R. Int. Arch. Allergy Immunol. 2013; 161:229–233. [PubMed: 23548307]
- [58]. Hochwallner H, Schulmeister U, Swoboda I, Spitzauer S, Valenta R. Methods. 2013
- [59]. Caballero ML, Umpierrez A, Perez-Pinar T, Moneo I, de Burgos C, Asturias JA, Rodriguez-Perez R. Int. Arch. Allergy Immunol. 2012; 158:232–240. [PubMed: 22398334]
- [60]. Bazaral M, Hamburger RN. J. Allergy Clin. Immunol. 1972; 49:189–191. [PubMed: 4622102]



#### Fig. 1.

Detection of allergen-specific IgE and IgG by the MeDALL allergen-chip. A serum sample from an allergic individual was tested for the presence of allergen-specific IgE (A) and IgG (B) using the MeDALL allergen-chip. Each of the allergen molecules on the chip was spotted in triplicates (groups of three dots aligned vertically), signal strength is visualized in false colours with increasing fluorescence intensities from blue to red/white.

Lupinek et al.



#### Fig. 2.

Measurement of monoclonal human Bet v 1-specific IgE antibody levels using the MeDALL allergen-chip and by ImmunoCAP. Different concentrations of a monoclonal human IgE antibody specific to Bet v 1 (*x*-axes) were tested. Bet v 1-specific IgE-levels measured by the MeDALL allergen-chip (A–C) are shown on the y-axes in ISAC standardized units (ISU). Bet v 1-specific IgE-levels measured by Bet v 1-ImmunoCAP (D–F) are shown in units antigen per milliliter (UA/ml). Mean levels +/– standard deviations of triplicate measurements are shown.





#### Fig. 3.

Measurement of different concentrations of a monoclonal human Phl p 2-specific IgG<sub>1</sub> antibody. Different concentrations of a monoclonal human IgG<sub>1</sub> antibody specific to Phl p 2 (*x*-axes) were measured and are shown on the *y*-axes in ISAC standardized units for IgG (ISU-G). Graphs for IgG-concentrations of 0–1 µg/ml (A), 0–10 ng/ml (B) and 0–0.3 ng/ml (C) are shown. Error bars indicate mean levels +/– standard deviations of triplicate measurements.

Lupinek et al.



#### Fig. 4.

Simultaneous presence of blocking IgG reduces detected Bet v 1-specific IgE-levels. Relative excess of Bip1 compared to the respective IgE concentration and (in parenthesis) absolute concentration of Bip1 in ng/ml or  $\mu$ g/ml is shown on the *x*-axes. Bet v 1-specific IgE-levels detected by the MeDALL allergen-chip (A–C) or by ImmunoCAP (D–E) are displayed on the y-axes in ISU or in UA/ml after subtraction of background signals. IgE-concentrations applied are mentioned in the respective diagrams. Relative inhibition (F) of IgE-signals (*y*-axis) is plotted against levels of relative Bip1-excess (*x*-axis). The respective

IgE-concentrations applied for detection by MeDALL allergen-chip or ImmunoCAP are displayed in the legend. Mean results of triplicates +/- standard deviations are shown.



### Fig. 5.

Bet v 1-specific IgE (312.5 ng/ml) was co-incubated with an excess of two different monoclonal antibodies binding to two distinct Bet v 1-peptides, or of a monoclonal antibody specific for the grass pollen allergen Phl p 2 (*x*-axis). Degree of inhibition of IgE-binding to chip-coupled Bet v 1 is shown on the *y*-axis. Bars show mean values of triplicates +/- standard deviation.



#### Fig. 6.

Effect of different diluents on detected IgE-levels. Monoclonal Bet v 1-specific IgE was diluted to a concentration of 6.25 ng/ml in a sample diluent containing blocking agents or in PBS (*x*-axis) and measured by the MeDALL allergen-chip. IgE-levels are displayed in ISU (*y*-axis). Mean results of triplicate determinations +/- standard deviations are shown.



#### Fig. 7.

Screenshots illustrating the evolution of an IgE-reactivity profile over time. Serum samples from a study participant of the Norwegian birth cohort ECA, obtained at (A) 10 years and (B) 16 years, were analyzed for IgE-reactivities using the MeDALL allergen-chip. Every allergen-molecule was spotted in triplicates (groups of three dots aligned vertically), signal strength is visualized in false colours with increasing fluorescence intensities from blue to red/white.

List of allergens contained on the MeDALL allergen-chip. Allergen-sources, allergen-names according to IUIS-standard, purification from extracts (natural allergens) or production in expression systems (recombinant allergens), glycosylation and allocation to different allergen-families are shown. Proteins that were provided by collaborating groups are highlighted by coloured boxes (column "Allergen"). For wheat allergens, # indicates the number of the respective clones. Some allergen families represented on the chip by several allergen molecules are highlighted by coloured boxes (column "Function of the protein/Allergen family"). In the separate part of the list, allergens that replaced proteins from the first variant of the microarray (Proteins No. 164–176) are listed.

	Version 1 of the MeDALL-Chip									
No.	Allergen-source	Allergen	Rec. or natural	CCD	Route of exposure	Function of the protein/Allergen family				
1		Bet v 1	R		inh	PR-10				
2	Birch	Bet v 2	R		inh	Profilin				
3		Bet v 4	R		inh	Polcalcin				
4	Alder	Aln g 1	R		inh	PR-10				
5		Ole e 1	R		inh	Ole e 1-related protein				
6	Olive	Ole e 7	Ν		inh	nsLTP, type 1				
7		Ole e 9	R		inh	Glucanase				
8	Japanese cedar	Cry j 1	Ν	CCD	inh	Pectate lyase				
9	Cypress	Cup a 1	Ν	CCD	inh	Pectate lyase				
10		Pla a 1	R		inh	Invertase Inhibitor				
11	Plane tree	Pla a 2	Ν	CCD	inh	Polygalacturonases				
12		Pla a 3	R		inh	nsLTP, type 1				
13		Phl p 1	R		inh	Grass group 1 (Beta-Expansin)				
14		Phl p 2	R		inh	Grass group 2/3				
15		Phl p 4	Ν	CCD	inh	Grass group 4 (Berberine bridge enzyme)				
16	Timothy grass	Phl p 5b	R		inh	Grass group 5				
17	Timoiny grass	Phl p 6	R		inh	Grass group 5/6				
18		Phl p 7	R		inh	Calcium-binding proteins (Polcalcin)				
19		Phl p 11	R		inh	Ole e 1-related protein				
20		Phl p 12	R		inh	Profilin				
21	Bermuda grass	Cyn d 1	Ν	CCD	inh	Grass group 1 (Beta-Expansin)				
22	Ragweed	Amb a 1	Ν		inh	Pectate lyase				
23	Muguort	Art v 1	Ν		inh	Defensin-like protein				
24	Mugwon	Art v 3	Ν		inh	nsLTP, type 1				
25	Goosefoot	Che a 1	R		inh	Ole e 1-related protein				
26	Annual mercury	Mer a 1	R		inh	Profilin				
27	Plantain	Pla 1 1	R		inh	Ole e 1-related protein				
28	Wall pellitory	Par j 2	R		inh	nsLTP, type 1				
29	Saltwort	Sal k 1	Ν		inh	Pectin methylesterase				

No.       Allergen-source       Allergen       Rec. or human       CCD       Route of exposure       Function of the protein/Allergen family         30       Instantian       R       inh.skin       Rubber elongation factor         31       Hev b 1       R       inh.skin       Rubber elongation factor         33       Hev b 5       R       inh.skin       Rubber elongation factor         34       Hev b 5       R       inh.skin       Mathematic         35       Alternaria       Alt a 1       R       inh.skin       Portin         36       Alternaria       Alt a 6       R       inh       Matoglycoprotein         38       Aspergillus       Asp f 1       R       inh       Mitoglin family (Ribonucless)         38       Aspergillus       Asp f 6       R       inh       Manantol dehydrogenase         40       Cludosporium       Cla 8       R       inh       Group 1 mite allergens (Cysteine protease)         41       Der p 1       N       inh       Group 1 mite allergens (NPC2 family)         44       Der p 2       R       inh       Group 2 mite allergens (NPC2 family) <td< th=""><th></th><th colspan="8">Version 1 of the MeDALL-Chip</th></td<>		Version 1 of the MeDALL-Chip							
30   Image: Hev b 1   R   inh. skin   Rubber elongation factor     31   Hev b 3   R   inh. skin   Small rubber particle protein     32   Hev b 6.01   R   inh. skin   Acidic protein     33   Hev b 6.01   R   inh. skin   Acidic protein     34   Hev b 6.01   R   inh. skin   Hevin Prohevein)     35   Alternaria   Alt a 1   R   inh. Skin   Portin     36   Alternaria   Ats 1   R   inh   Acidic glycoprotein     37   Aspergillus   Asp f 3   R   inh   Mitogillin family (Ribonuclease)     38   Aspergillus   Asp f 3   R   inh   Manitol dehydrogenase     40   Cladosporum   Cla 8   R   inh   Manitol dehydrogenase     41   Der p 1   N   inh   Group 1 mite allergens (Cysteine protease)     42   Der p 1   N   inh   Group 2 mite allergens (Cysteine protease)     44   Der p 1   N   inh   Group 2 mite allergens (Cysteine protease)     45   Der p 1   R   inh   Group 1 mite allergens (Cysteine protease)     46   Der p 1   R   inh   Group 2 mite allergens     47	No.	Allergen-source	Allergen	Rec. or natural	CCD	Route of exposure	Function of the protein/Allergen family		
31Hev b 3Rinh. skinSmall rubber particle protein32LatexHev b 5Rinh. skinAcidic protein33Hev b 6.01Rinh. skinHevein (Protevcin)34Hev b 8Rinh. skinBofilin35AlternariaAlt a 1RinhAcidic glycoprotein36AlternariaAlt a 6RinhMitogillin family (Ribonuclease)38AspergillusAsp f 1RinhMitogillin family (Ribonuclease)39Asp f 6RinhMinoid dehydrogenase40CladosporiumCla h 8RinhManioid dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (NPC2 family)45Der p 2RinhGroup 2 mite allergens (NPC2 family)46Der p 10RinhGroup 1 mite allergens (NPC2 family)50Der p 10RinhGroup 2 mite allergens (NPC2 family)51Der p 11RinhGroup 2 mite allergens (Parimyosin)52Der p 13RinhChitn-binding domain53Der p 14RinhChitn-binding domain54Der p 15RinhChitn-binding domain55Der p 14RinhChitn-binding domain56Storage miteL	30		Hev b 1	R		inh, skin	Rubber elongation factor		
32LatexHev b 5Rinh, skinAcidic protein33Hev b 6.01Rinh, skinHevein (Prohevein)34Hev b 8Rinh, skinProfilin35Alt a 1RinhAcidic glycoprotein36AlternariaAlt a 6RinhEnclase37AspergillusAsp f 1RinhMitogillin family (Ribonuclease)38AspergillusAsp f 6RinhMinogulin family (Ribonuclease)40CladoporiumCla h 8RinhMinoid dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der p 1NinhGroup 1 mite allergens (NPC2 family)43Der f 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (Alha-amylase)45Der p 7RinhGroup 1 mite allergens (Alha-amylase)46Der p 10RinhGroup 1 mite allergens (Alha-amylase)50Der p 11RinhChritte-binding domain51Der p 15RinhChritte-binding domain52Der p 18RinhChritte-binding domain53Der p 18RinhChritte-binding domain54Der p 18RinhChritte-binding domain55Der p 18RinhChritte-binding domain56Storage miteLep d 2Rinh<	31		Hev b 3	R		inh, skin	Small rubber particle protein		
33Hev b 6.01Rinh, skinHevein (Prohevein)34Hev b 8Rinh, skinProfilin35AlternariaAlt a 1RinhAcidic glycoprotein36AlternariaAlt a 6RinhEnclase37AspergillusAsp f 1RinhMitogillin family (Ribonuclease)38AspergillusAsp f 3RinhMutogillin family (Ribonuclease)39OcladosporiumClah 8RinhMusperoxide dismutase40CladosporiumClah 8RinhManitol dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 2 mite allergens (NPC2 family)43Der f 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 1 mite allergens (NPC2 family)45Der p 7RinhGroup 1 mite allergens (NPC2 family)46Der p 7RinhGroup 1 mite allergens (NPC2 family)50Der p 10RinhGroup 1 mite allergens (NPC2 family)51Der p 11RinhGroup 1 mite allergens (NPC2 family)52Der p 11RinhGroup 1 mite allergens (NPC2 family)53Der p 14RinhGroup 1 mite allergens (NPC2 family)54Der p 15RinhChrin-binding domain55Storage miteLep d 2Rinh <t< td=""><td>32</td><td>Latex</td><td>Hev b 5</td><td>R</td><td></td><td>inh, skin</td><td>Acidic protein</td></t<>	32	Latex	Hev b 5	R		inh, skin	Acidic protein		
34Hev b 8Rinh, skinProfilin35AlternariaAlt a 1RinhAcidic glycoprotein36AlternariaAlt a 6RinhEnolase37Asp f 1RinhMitogilli family (Ribonuclease)38AspergillusAsp f 3RinhPeroxysomal protein39Asp f 6RinhMa superoxide dismutase40CladosporiumCla h 8RinhManitol dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cystein protease)42Der f 1NinhGroup 2 mite allergens (Cystein protease)43Der f 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (Alpha-amylase)45Der f 7RinhGroup 1 mite allergens (Alpha-amylase)46Der p 5RinhGroup 1 mite allergens (Alpha-amylase)47Der p 7RinhGroup 1 mite allergens (Paramyosin)48Der p 11RinhGroup 1 mite allergens (Paramyosin)50Der p 11RinhChritin-binding domain51Der p 11RinhChritin-binding domain52Der p 13RinhGroup 2 mite allergens (Paramyosin)53Der p 14RinhChritin-binding domain54Der p 13RinhChritin-binding domain55Der p 14Rinh <td>33</td> <td></td> <td>Hev b 6.01</td> <td>R</td> <td></td> <td>inh, skin</td> <td>Hevein (Prohevein)</td>	33		Hev b 6.01	R		inh, skin	Hevein (Prohevein)		
35 36AlternariaAlt a 1RinhAcidic glycoprotein36AlternariaAlt a 6RinhEnclase37 37AspergillusAsp f 1RinhMitogillin family (Ribonuclease)38 39AspergillusAsp f 3RinhPeroxysomal protein39CladosporiumCla h 8RinhManuitol dehydrogenase40CladosporiumCla h 8RinhGroup 1 mite allergens (Cystein protease)41Der p 1NinhGroup 2 mite allergens (Cystein protease)42Der p 1NinhGroup 2 mite allergens (NPC2 family)43Der p 2RinhGroup 2 mite allergens (NPC2 family)44Der p 3RinhGroup 1 mite allergens (NPC2 family)45Der p 4RinhGroup 1 mite allergens (NPC2 family)46Der p 5RinhGroup 2 mite allergens (NPC2 family)47Der p 10RinhGroup 1 mite allergens (NPC2 family)48House dust miteDer p 10Rinh49Der p 11RinhGroup 1 mite allergens (NPC2 family)50Der p 13RinhGroup 1 mite allergens (NPC2 family)51Der p 14RinhChtin-binding domain52Der p 15RinhChtin-binding domain53Der p 18RinhChtin-binding domain54Der p 19RinhChtin-binding	34		Hev b 8	R		inh, skin	Profilin		
36AlternationAlt a 6RinhEnclase37AspergillusAsp f 1RinhMitogillin family (Ribonuclease)38AspergillusAsp f 3RinhPeroxysomal protein39CladosporiumCla h 8RinhMannitol dehydrogenase40CladosporiumCla h 8RinhGroup 1 mite allergens (Cysteine protease)41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 2 mite allergens (NPC2 family)43Der f 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (NPC2 family)44Der p 3RinhGroup 1 mite allergens (NPC2 family)45Der p 4RinhGroup 1 mite allergens (NPC2 family)46Der p 7RinhGroup 1 mite allergens (NPC2 family)47Der p 10RinhGroup 1 mite allergens (NPC2 family)48House dust miteDer p 10Rinh49Der p 11RinhGroup 1 mite allergens (NPC2 family)50Der p 12RinhGroup 1 mite allergens (NPC2 family)51Der p 13RinhGroup 2 mite allergens52Der p 14RinhChitin-binding domain53Der p 18RinhChitin-binding domain54Der p 19RinhChitin-binding domain <td>35</td> <td>Altownania</td> <td>Alt a 1</td> <td>R</td> <td></td> <td>inh</td> <td>Acidic glycoprotein</td>	35	Altownania	Alt a 1	R		inh	Acidic glycoprotein		
37AspergillusAsp f 1RinhMitogillin family (Ribonuclease)38AspergillusAsp f 3RinhPeroxysomal protein39Asp f 6RinhMn superoxide dismutase40CladosportumCla h 8RinhMannitol dehydrogenase41Per p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 2 mite allergens (Cysteine protease)43Der f 1NinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (NPC2 family)45Der f 2RinhGroup 7 mite allergens (NPC2 family)46Der p 3RinhGroup 7 mite allergens (Alpha-amylase)47Der p 4RinhGroup 7 mite allergens (Toponyosin)48Der p 7RinhGroup 7 mite allergens (Toponyosin)49Der p 10RinhGroup 7 mite allergens (Toponyosin)50Der p 11RinhChroup 7 mite allergens (Toponyosin)51Der p 18RinhChrolipophorins)53Der p 18RinhChroup 7 mite allergens (Paramyosin)54Der p 18RinhChroup 1 mite allergens (NPC2 family)55Der p 18RinhChritin-binding domain56Storage miteLep d 2Rinh57Blonia tropicalisBlo 15Rinh58Group	36	Allemaria	Alt a 6	R		inh	Enolase		
38AspergillusAsp f 3RinhPeroxysomal protein39Asp f 6RinhMn superoxide dismutase40CladosporiumCla h 8RinhMannitol dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 2 mite allergens (Cysteine protease)43Der f 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (NPC2 family)45Der f 2RinhGroup 4 mite allergens (Alpha-amylase)46Der f 3RinhGroup 7 mite allergens47House dust miteDer f 1Rinh48House dust miteDer f 1Rinh49Der f 1RinhGroup 7 mite allergens51Der p 10RinhGroup 11 mite allergens (Paramyosin)51Der p 14RinhChtitn-binding domain52Der p 15RinhChtitn-binding domain53Der p 21RinhChtitn-binding domain54CockroachBlo 1 5Rinh57Blomia tropicalisBlo 1 5Rinh58CockroachBla g 1Rinh59CockroachBla g 1Rinh60AnisakisAni 1Rinh61Si a 1Rinh62AnisakisAni 1<	37		Asp f 1	R		inh	Mitogillin family (Ribonuclease)		
39Asp f 6RinhMn superoxide dismutase40CladosporiumCla h 8RinhMannitol dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 1 mite allergens (Cysteine protease)43Der p 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (NPC2 family)45Der p 7RinhGroup 5 mite allergens (Alpha-amylase)46Der p 7RinhGroup 7 mite allergens (Alpha-amylase)47House dust miteDer p 10RinhGroup 7 mite allergens (Tropomyosin)50Der p 11RinhGroup 1 mite allergens (Paramyosin)51Der p 11RinhGroup 1 mite allergens (Paramyosin)53Der p 14RinhChilin-binding domain54Der p 12RinhChilin-binding domain55Der p 23RinhChilin-binding domain54Der p 23RinhChilin-binding domain55Storage miteLep d 2Rinh58GockroachBla g 1Rinh59CockroachBla g 1Rinh60CatorageBla g 7Ninh61AnisakisAni s 1Rf62CatFel d 1Rinh64Fel d 1RinhUbroglo	38	Aspergillus	Asp f 3	R		inh	Peroxysomal protein		
40CladosporiumCla h 8RinhMannitol dehydrogenase41Der p 1NinhGroup 1 mite allergens (Cysteine protease)42Der f 1NinhGroup 2 mite allergens (Cysteine protease)43Der p 2RinhGroup 2 mite allergens (NPC2 family)44Der f 2RinhGroup 2 mite allergens (NPC2 family)45Der p 2RinhGroup 2 mite allergens (NPC2 family)46Der p 7RinhGroup 5 mite allergens (Alpha-amylase)47Der p 7RinhGroup 7 mite allergens48House dust miteDer p 10RinhGroup 1 mite allergens (Alpha-amylase)50Der p 10RinhGroup 1 mite allergens (Paramyosin)51Der p 11RinhGroup 11 mite allergens (Paramyosin)52Der p 14RinhChitin-binding domain53Der p 15RinhChitin-binding domain54Der p 23RinhChitin-binding domain55Der p 23RinhChitin-binding domain56Storage miteLep d 2RinhGroup 2 mite allergens58SCockroachBla g 1RinhGroup 5 mite allergens60CockroachBla g 7NinhGroup 5 mite allergens61GaAni s 1RinhGroup 5 mite allergens62Ani sakisAni s 1RfTropom	39		Asp f 6	R		inh	Mn superoxide dismutase		
41Der p 1NinhGroup 1 mite allergens (Cysteine protease)4243444444444546464746474849495050515253545556575854565758596061626162636464656677777788999999999999999999999999999999999999999999999999999 </td <td>40</td> <td>Cladosporium</td> <td>Cla h 8</td> <td>R</td> <td></td> <td>inh</td> <td>Mannitol dehydrogenase</td>	40	Cladosporium	Cla h 8	R		inh	Mannitol dehydrogenase		
42 $42$ $43$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$ $24$	41		Der p 1	N		inh	Group 1 mite allergens (Cysteine protease)		
43 $44$ $44$ $45$ $2mr p 2$ $R$ $mh$ $3roup 2 mite allergens (NPC2 family)444546464616mh3roup 2 mite allergens (NPC2 family)46474616mh3roup 4 mite allergens (Alpha-amylase)4748478mh1nhh3roup 5 mite allergens (Alpha-amylase)484716mh3roup 7 mite allergens (Tropomyosin)499mr p 10R1nhh3roup 10 mite allergens (Tropomyosin)509mr p 14R1nhh3roup 11 mite allergens (Paramyosin)519mr p 14R1nhh3roup 11 mite allergens (Paramyosin)529mr p 14R1nhh3roup 11 mite allergens (Paramyosin)539mr p 14R1nhh3roup 21 mite allergens549mr p 15R1nhh3roup 21 mite allergens559mr p 14R1nhh3roup 21 mite allergens563rorage mite1ep d 2R1nhh578lomia tropicalisBlo t 5R1nhh583roup 28la g 1R1nhh59CockroachBla g 1R1nhh616mr a 15R1nhh3roup 2 mite allergens616mr a 15R1nhh3roup 2 mite allergens62AnisakisAnis 1Rf7ropomyosin636nr a 15$	42		Der f 1	N		inh	Group 1 mite allergens (Cysteine protease)		
44 $45$ $46$ $2mic$ $R$ $mh$ $Group 2$ mite allergens (NPC2 family) $46$ $47$ $48$ $6roup$ $6roup$ $6roup$ $7mite$ $8roup$ <td< td=""><td>43</td><td></td><td>Der p 2</td><td>R</td><td></td><td>inh</td><td>Group 2 mite allergens (NPC2 family)</td></td<>	43		Der p 2	R		inh	Group 2 mite allergens (NPC2 family)		
45 $46$ $47$ $Ber p 4$ $R$ $inh$ $Group 4$ mite allergens (Alpha-amylase)47 $46$ $47$ $Ber p 5$ $R$ $inh$ $Group 5$ mite allergens48 $Ber p 7$ $R$ $inh$ $Group 10$ mite allergens (Tropomyosin)49 $Der p 10$ $R$ $inh$ $Group 10$ mite allergens (Paramyosin)50 $Der p 11$ $R$ $inh$ $Group 11$ mite allergens (Paramyosin)51 $Der p 14$ $R$ $inh$ $Chitin-binding domain52Der p 15RinhChitin-binding domain53Der p 21RinhChitin-binding domain54Der p 23RinhChitin-binding domain55Der p 14RinhChitin-binding domain56Storage miteLep d 2RinhChitin-binding domain56Storage miteLeg d 2RinhGroup 2 mite allergens (NPC2 family)57Blomia tropicalisBlo t 5RinhGroup 5 mite allergens58SeroachBla g 1RinhAspartic protease60CackroachBla g 7NinhSerine protease inhibitor61AnisakisAni s 1RfSerine protease inhibitor63CatFel d 1RinhUteroglobin64Fel d 4RinhSerum Albumin$	44		Der f 2	R		inh	Group 2 mite allergens (NPC2 family)		
46 47 $46$ 47 $47$ 48 $47$ House dust mite $Der p 5$ $R$ inh $Group 5$ mite allergens $49$ $49$ $Der p 10$ $R$ inh $Group 10$ mite allergens (Tropomyosin) $50$ $Der p 11$ $R$ inh $Group 10$ mite allergens (Paramyosin) $51$ $Der p 14$ $R$ inh $Chin-binding domain52Der p 15RinhChiin-binding domain53Der p 18RinhChiin-binding domain54Der p 23RinhChiin-binding domain55Der p 23RinhChiin-binding domain56Storage miteLep 42RinhGroup 2 mite allergens57Blomia tropicalisBlo t 5RinhGroup 5 mite allergens58CockroachBla g 1RinhGoup 5 mite allergens60CockroachBla g 5RinhGlutathione S-transferase61AnisakisAnis 1RfSerine protease inhibitor63Anis 3RfTropomyosin64Fel d 1RinhUteroglobin65CatFel d 2NinhSerum Albumin66Fel d 4RinhLipocalin$	45		Der p 4	R		inh	Group 4 mite allergens (Alpha-amylase)		
47House dust miteDer p 7RinhGroup 7 mite allergens48House dust miteDer p 10RinhGroup 10 mite allergens (Tropomyosin)50Der p 11RinhGroup 11 mite allergens (Paramyosin)50Der p 11RinhVitellogenin (Apolipophorins)51Der p 14RinhChitin-binding domain52Der p 15RinhChitin-binding domain53Der p 18RinhChitin-binding domain54Der p 21RinhGroup 21 mite allergens54Der p 23RinhChitin-binding domain55Storage miteLep d 2RinhGroup 2 mite allergens (NPC2 family)57Blomia tropicalisBlot 5RinhGroup 5 mite allergens58CockroachBla g 1RinhCockroach group 160Bla g 7NinhGlutathione S-transferase61AnisakisAni s 1RfSerine protease inhibitor63AnisakisFel d 1RinhUteroglobin64Fel d 1RinhUteroglobin65CatFel d 2NinhSerum Albumin	46		Der p 5	R		inh	Group 5 mite allergens		
48House dust mireDer p 10RinhGroup 10 mite allergens (Tropomyosin)49950505151525152525354545455656Storage mite57Blomia tropicalis58Bla g 160Cockroach61Bla g 262Anisakis63Anisakis64Fel d 165Cat66606162646465666666	47		Der p 7	R		inh	Group 7 mite allergens		
49Der p 11RinhGroup 11 mite allergens (Paramyosin)5050Der p 14RinhVitellogenin (Apolipophorins)5151Der p 15RinhChitin-binding domain52Der p 18RinhChitin-binding domain53Der p 18RinhChitin-binding domain54Der p 21RinhGroup 21 mite allergens54Der p 23RinhChitin-binding domain55clone 16RinhChitin-binding domain56Storage miteLep d 2RinhGroup 2 mite allergens (NPC2 family)57Blomia tropicalisBlo t 5RinhGroup 5 mite allergens58 $\mathcal{O}$ $\mathcal{O}$ Bla g 1RinhCockroach group 159 $\mathcal{O}$ Bla g 2RinhGlutathione S-transferase61 $\mathcal{O}$ $\mathcal{A}$ nis 3 $\mathcal{R}$ fSerine protease inhibitor62 $\mathcal{A}$ nisakis $\mathcal{A}$ ni s 3 $\mathcal{R}$ fTropomyosin64 $\mathcal{C}$ $\mathcal{F}$ $\mathcal{I}$ $\mathcal{R}$ inhUteroglobin65 $\mathcal{C}$ $\mathcal{F}$ $\mathcal{I}$ $\mathcal{R}$ inhSerine protease inhibitor65 $\mathcal{L}$ $\mathcal{R}$ inh $\mathcal{L}$ $\mathcal{R}$ inh66 $\mathcal{F}$ $\mathcal{I}$ $\mathcal{R}$ inh $\mathcal{L}$ 67 $\mathcal{I}$ $\mathcal{I}$ $\mathcal{R}$ inh $\mathcal{L}$ 68 $\mathcal{I}$ $\mathcal{I}$ $\mathcal{R}$ $$	48	House dust mite	Der p 10	R		inh	Group 10 mite allergens (Tropomyosin)		
50Der p 14RinhVitellogenin (Apolipophorins) $51$ $51$ $Der p 15$ $R$ inhChitin-binding domain $52$ $Der p 18$ $R$ inhChitin-binding domain $53$ $Der p 21$ $R$ inhChitin-binding domain $54$ $Der p 21$ $R$ inhChitin-binding domain $55$ $Der p 23$ $R$ inhChitin-binding domain $56$ $Storage mite$ Lep d 2 $R$ inhGroup 2 mite allergens (NPC2 family) $57$ $Blomia tropicalis$ Blo t 5 $R$ inhGroup 5 mite allergens $58$ $Bla g 1$ $R$ inhCockroach group 1 $59$ $Cockroach$ Bla g 2 $R$ inhGlutathione S-transferase $61$ $Anisakis$ $Ani s 1$ $R$ fSerine protease inhibitor $63$ $Anisakis$ $Ani s 3$ $R$ fTropomyosin $64$ $Fel d 1$ $R$ inhUteroglobin $65$ $Cat$ $Fel d 1$ $R$ inhSerum Albumin	49		Der p 11	R		inh	Group 11 mite allergens (Paramyosin)		
51Der p 15RinhChitin-binding domain $52$ $54$ $Der p 18$ RinhChitin-binding domain $53$ $Der p 21$ RinhGroup 21 mite allergens $54$ $Der p 23$ RinhChitin-binding domain $55$ $Der p 23$ RinhChitin-binding domain $56$ $Storage mite$ Lep d 2RinhChitin-binding domain $56$ $Storage mite$ Lep d 2RinhGroup 2 mite allergens (NPC2 family) $57$ $Blomia tropicalis$ Blo t 5RinhGroup 5 mite allergens $58$ $Pockroach$ Bla g 1RinhCockroach group 1 $59$ $Cockroach$ Bla g 2RinhAspartic protease $61$ $Poly 23$ RinhGlutathione S-transferase $61$ $Anis align 3$ RfSerine protease inhibitor $63$ $Anis 3$ RfTropomyosin $64$ $Cat$ Fel d 1RinhUteroglobin $65$ $Cat$ Fel d 2NinhSerum Albumin	50		Der p 14	R		inh	Vitellogenin (Apolipophorins)		
52Der p 18RinhChitin-binding domain53 $Der p 21$ RinhGroup 21 mite allergens54 $Der p 23$ RinhChitin-binding domain55 $Der p 23$ RinhChitin-binding domain56 $Storage mite$ $Lep d 2$ RinhGroup 2 mite allergens (NPC2 family)57 $Blomia tropicalis$ $Blo t 5$ RinhGroup 5 mite allergens58 $R$ $Inh$ $Cockroach$ $Bla g 1$ Rinh59 $Cockroach$ $Bla g 2$ RinhAspartic protease60 $Bla g 7$ NinhGlutathione S-transferase61 $Bla g 7$ NinhTropomyosin62 $Anisakis$ $Ani s 1$ RfSerine protease inhibitor63 $Cat$ $Fel d 1$ RinhUteroglobin64 $Fel d 1$ RinhUteroglobin65 $Cat$ $Fel d 4$ RinhLipocalin	51		Der p 15	R		inh	Chitin-binding domain		
53Der p 21RinhGroup 21 mite allergens $54$ Der p 23RinhChitin-binding domain $55$ clone 16RinhChitin-binding domain $56$ Storage miteLep d 2RinhGroup 2 mite allergens (NPC2 family) $57$ Blomia tropicalisBlo t 5RinhGroup 5 mite allergens $58$ $58$ Bla g 1RinhCockroach group 1 $59$ CockroachBla g 2RinhAspartic protease $60$ CockroachBla g 5RinhGlutathione S-transferase $61$ Anis 31RfSerine protease inhibitor $63$ Anis 33RfTropomyosin $64$ Fel d 1RinhUteroglobin $65$ CatFel d 1RinhLipocalin	52		Der p 18	R		inh	Chitin-binding domain		
54Der p 23RinhChitin-binding domain $55$ $clone 16$ RinhChitin-binding domain $56$ $Storage mite$ Lep d 2RinhGroup 2 mite allergens (NPC2 family) $57$ $Blomia tropicalis$ Blo t 5RinhGroup 5 mite allergens $58$ $R$ Bla g 1RinhCockroach group 1 $59$ $Cockroach$ Bla g 2RinhAspartic protease $60$ $Cockroach$ Bla g 5RinhGlutathione S-transferase $61$ $R$ inhTropomyosin $62$ $Anisakis$ Ani s 1Rf $63$ $Cat$ Fel d 1RinhUteroglobin $64$ $Fel d 2$ NinhSerum Albumin $65$ $Cat$ Fel d 4RinhLipocalin	53		Der p 21	R		inh	Group 21 mite allergens		
55clone 16RinhChitin-binding domain56Storage miteLep d 2RinhGroup 2 mite allergens (NPC2 family)57Blomia tropicalisBlo t 5RinhGroup 5 mite allergens58 $\mathbb{P}$ Bla g 1RinhCockroach group 159 $\mathbb{P}$ $\mathbb{P}$ Bla g 2RinhAspartic protease60 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 61 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 62 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 63 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 64 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 65 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 66 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ 66 $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$ $\mathbb{P}$	54		Der p 23	R		inh	Chitin-binding domain		
56Storage miteLep d 2RinhGroup 2 mite allergens (NPC2 family) $57$ Blomia tropicalisBlo t 5RinhGroup 5 mite allergens $58$ $80 + 5$ RinhCockroach group 1 $59$ $Cockroach$ Bla g 2RinhAspartic protease $60$ $Cockroach$ Bla g 5RinhGlutathione S-transferase $61$ Bla g 7NinhTropomyosin $62$ $Anisakis$ Ani s 1RfSerine protease inhibitor $63$ $Cat$ Fel d 1RinhUteroglobin $64$ $Fel d 1$ RinhSerum Albumin $65$ $Cat$ Fel d 4RinhLipocalin	55		clone 16	R		inh	Chitin-binding domain		
57Blomia tropicalisBlo t 5RinhGroup 5 mite allergens58 $Bla g 1$ RinhCockroach group 159 $Cockroach$ $Bla g 1$ RinhAspartic protease60 $Bla g 5$ RinhGlutathione S-transferase61 $Bla g 7$ NinhTropomyosin62 $Anisakis$ Ani s 1RfSerine protease inhibitor63 $Ani s 3$ RfTropomyosin64 $El d 1$ RinhUteroglobin65 $Cat$ $Fel d 2$ NinhSerum Albumin66 $Fel d 4$ RinhLipocalin	56	Storage mite	Lep d 2	R		inh	Group 2 mite allergens (NPC2 family)		
$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	57	Blomia tropicalis	Blot 5	R		inh	Group 5 mite allergens		
$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	58		Bla g 1	R		inh	Cockroach group 1		
60CockroachBla g 5RinhGlutathione S-transferase61Bla g 7NinhTropomyosin62AnisakisAni s 1RfSerine protease inhibitor63Anis 3RfTropomyosin64Fel d 1RinhUteroglobin65CatFel d 2NinhSerum Albumin66Fel d 4RinhLipocalin	59		Bla g 2	R		inh	Aspartic protease		
	60	Cockroach	Bla g 5	R		inh	Glutathione S-transferase		
62 63Ani s 1RfSerine protease inhibitor63Ani s 3RfTropomyosin64 65Fel d 1RinhUteroglobin66Fel d 2NinhSerum Albumin66Fel d 4RinhLipocalin	61		Bla g 7	N		inh	Tropomyosin		
63AnisakisAnis 3RfTropomyosin64Fel d 1RinhUteroglobin65CatFel d 2NinhSerum Albumin66Fel d 4RinhLipocalin	62		Ani s 1	R		f	Serine protease inhibitor		
64   Fel d 1   R   inh   Uteroglobin     65   Cat   Fel d 2   N   inh   Serum Albumin     66   Fel d 4   R   inh   Lipocalin	63	Anisakis	Ani s 3	R	1	f	Tropomyosin		
65   Cat   Fel d 2   N   inh   Serum Albumin     66   Fel d 4   R   inh   Lipocalin	64		Fel d 1	R		inh	Uteroglobin		
66 Fel d 4 R inh Lipocalin	65	Cat	Fel d 2	N	1	inh	Serum Albumin		
	66		Fel d 4	R		inh	Lipocalin		
67 Dog Can f 1 R inh Lipocalin	67	Dog	Can f 1	R		inh	Lipocalin		

No.Allergen-sourceAllergenRc. or nuturalCCDRoute of exposureFunction of the protein/Allergen family66Can f 2RinhLipocalin70Can f 3NinhSerum Albounin70Can f 4RinhLipocalin (Odorant-binding protein)71Can f 5RinhArginine Esterase (Trysin-like serine protease)72Can f 5RinhLipocalin73HorseEqu c 1RinhLipocalin74HorseEqu c 3NinhLipocalin75MouseMus m 1NinhLipocalin76MouseMus m 1NinhLipocalin77AnuseMus m 1NinhLipocalin78MouseMus m 1NinhLipocalin79Aru h 3NIStorage protein, 75 globulin (Vicilins)78Aru h 3RIStorage protein, 75 globulin (Vicilins)79Aru h 3RIStorage protein, 75 globulin (Ligumins)80Aru h 3RIStorage protein, 75 globulin (Ligumins)81Aru h 3RIStorage protein, 75 globulin (Ligumins)82Aru h 3RIStorage protein, 15 globulin (Ligumins)83Aru h 3RIStorage protein, 15 globulin (Ligumins)84Aru h 9RII85Cora aRII86<	Version 1 of the MeDALL-Chip								
68Can f 2RinhLipocalin70Can f 3NinhSerun Albumin71Can f 3RinhLipocalin (Odorant-binding protein)72Can f 5RinhArginine Esterase (Trypsin-like serine protease)73Can f 5RinhLipocalin74HorseEqu c 1RinhLipocalin75HorseEqu c 1RinhLipocalin76MouseMus m 1NinhLipocalin77Amb 1NinhLipocalin78MouseMus m 1NinhLipocalin79Ara h 1RfStorage protein, 7S globulin (Vicilins)78Ara h 1RfStorage protein, 2S globulin (Vicilins)79Ara h 3NfStorage protein, 1S globulin (Vicilins)78Ara h 3NfStorage protein, 2S albumin79Ara h 6NfStorage protein, 1S globulin (Legumins)80Ara h 6NfStorage protein, 2S albumin81PearutCor a 1.0401RfStorage protein, 1S globulin (Legumins)83Ara h 6NfStorage protein, 1S globulin (Legumins)84Ara h 9RfMath 9R85Cor a 1.0401RfStorage protein, 1S globulin (Legumins)86Ara h 9RfStorage protein, 1S globulin (Legumins)87Hazehut<	No.	Allergen-source	Allergen	Rec. or natural	CCD	Route of exposure	Function of the protein/Allergen family		
69 70Can f 3NinhSerun Albumin71Can f 3RinhLipocalin (Odorant-binding protein)71Can f 5RinhArginine Esterase (Trypsin-like serine protease)72Can 6RinhLipocalin73HorseEqu c 1RinhLipocalin74HorseEqu c 1RinhLipocalin75MouseMus m 1NinhLipocalin76MouseMus m 1NinhLipocalin77Ara h 1NfStorage protein, 75 globulin (Vicilins)78Ara h 1RfStorage protein, 25 globulin (Vicilins)79Ara h 1RfStorage protein, 11 globulin (Legumins)80Ara h 2RfStorage protein, 11 globulin (Legumins)81PeanutAra h 3RfStorage protein, 25 albumin82Ara h 6NfStorage protein, 25 albumin83Ara h 6NfStorage protein, 11 globulin (Legumins)84Ara h 6NfStorage protein, 25 albumin85Cor a 1.0401Rfstorage protein, 11 globulin (Legumins)86Cor a 1.0401Rfstorage protein, 11 globulin (Legumins)87HazelnutCor a 8Rfaltr.type 188Cor a 1.0401Rfstorage protein, 11 globulin (Legumins)91Jug r 1Rfstorage prote	68		Can f 2	R		inh	Lipocalin		
7071727374781mhLipocalin (Odorant-binding protein)71Can f 5RinhArginine Esterase (Trypsin-like serine protease)73Can f 6RinhLipocalin74HorseEqu c 1RinhLipocalin75HorseEqu c 1RinhLipocalin76MouseMus m 1NinhLipocalin77Arah 1NinhLipocalin78Arah 1NfStorage protein, 78 globulin (Vicilins)79Arah 1RfStorage protein, 28 globulin (Vicilins)79Arah 1RfStorage protein, 28 globulin (Vicilins)79Arah 1RfStorage protein, 28 globulin (Vicilins)80Arah 3NfStorage protein, 28 globulin (Legunins)81PeanutArah 3RfStorage protein, 28 albumin82Arah 6NfStorage protein, 28 albumin83Arah 6NfStorage protein, 28 albumin84Arah 8RfGarage protein, 28 albumin85Arah 8RfStorage protein, 28 albumin86Cor a 1.0401RfStorage protein, 118 globulin (Legunins)87HazelnutCor a 8RfGarage protein, 18 globulin (Legunins)88Jug r 1RfStorage protein, 28 albumin99ManuetJug r 3Rf <td>69</td> <td></td> <td>Can f 3</td> <td>Ν</td> <td></td> <td>inh</td> <td>Serum Albumin</td>	69		Can f 3	Ν		inh	Serum Albumin		
71 $72$ $73$ $74$ $75$ $R$ $1$ inhArginine Estense (Trypsin-like serine protease) $73$ $74$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ $1$ <	70		Can f 4	R		inh	Lipocalin (Odorant-binding protein)		
72Can f SRinhArginine Esterase (Trypsin-like serine protease)73Can f GRinhLipocalin74HorseEqu c 1RinhLipocalin75HorseMus n 1NinhSerum Alkumin76MouseMus n 1NinhLipocalin77Ara h 1NinhLipocalin78Ara h 1RfStorage protein, 75 globulin (Vicilins)79Ara h 1RfStorage protein, 75 globulin (Vicilins)80Ara h 2RfStorage protein, 115 globulin (Legumins)81PeanutAra h 3RfStorage protein, 25 albumin83Ara h 3RfStorage protein, 25 albumin84Ara h 6NfStorage protein, 25 albumin85Ara h 6NfStorage protein, 25 albumin84Ara h 7RfStorage protein, 25 albumin85Ara h 9RfBall Ptypei86Ara h 9RfStorage protein, 25 albumin87HazelnutCor a 8RfStorage protein, 115 globulin (Legumins)88Jug r 1RfStorage protein, 115 globulin (Legumins)89Jug r 1RfStorage protein, 115 globulin (Legumins)91Jug r 1RfStorage protein, 115 globulin (Legumins)92Jug r 1RfStorage protein, 115 globulin (Legumins) <td>71</td> <td></td> <td>Can f 5</td> <td>R</td> <td></td> <td>inh</td> <td>Arginine Esterase (Trypsin-like serine protease)</td>	71		Can f 5	R		inh	Arginine Esterase (Trypsin-like serine protease)		
73Can f 6RinhLipocalin74HorseEqu c 1RinhLipocalin75HorseEqu c 3NinhSerum Albumin76MouseMus m 1NinhLipocalin77MouseMus m 1NinhLipocalin78Ara h 1RfStorage protein, 75 globulin (Vicilins)78Ara h 1RfStorage protein, 75 globulin (Vicilins)79Ara h 1RfStorage protein, 75 globulin (Legumins)79Ara h 3RfStorage protein, 25 albumin80PeanutAra h 3RfStorage protein, 25 albumin81PeanutAra h 3RfStorage protein, 25 albumin82Ara h 6NfStorage protein, 25 albumin83Ara h 6NfStorage protein, 25 albumin84Ara h 9Rfmt1P.type185Ara h 9Rfmt1P.type186Ara h 9Rfstorage protein, 11S globulin (Legumins)87HazelnutCor a 8Rfmt1P.type188Jug r 1RfStorage protein, 25 albumin90WahnutJug r 2NCCDfStorage protein, 75 globulin (Legumins)91Jug r 3Rfstorage protein, 75 globulin (Legumins)92Pro du 4RfStorage protein, 11S globulin (Legumins)93	72		Can f 5	R		inh	Arginine Esterase (Trypsin-like serine protease)		
74 HorseEqu c 1RinhLipocalin75HorseEqu c 3NinhSerum Albumin76MouseMus m 1NinhLipocalin77MouseMus m 1NinhLipocalin78Ara h 1NfStorage protein, 75 globulin (Vicilins)78Ara h 1RfStorage protein, 75 globulin (Vicilins)79Ara h 2RfStorage protein, 75 globulin (Vicilins)70Ara h 3RfStorage protein, 118 globulin (Legumins)80Ara h 3RfStorage protein, 118 globulin (Legumins)81PeanutAra h 6NfStorage protein, 22 albumin82Ara h 6NfStorage protein, 23 albumin83Ara h 6NfStorage protein, 23 albumin84Ara h 6NfBLTP, type 185Ara h 8RfmLTP, type 186Cor a 1.0401RfStorage protein, 118 globulin (Legumins)87HazelnutCor a 8RfmLTP, type 188Qara hfstorage protein, 118 globulin (Legumins)91Jug r 1RfStorage protein, 118 globulin (Legumins)92AnnoutJug r 2NCCDfStorage protein, 118 globulin (Legumins)93Jug r 3RfndLTP, type 194AlmondFStorage protein, 118 globulin (Legumins) <tr< td=""><td>73</td><td></td><td>Can f 6</td><td>R</td><td></td><td>inh</td><td>Lipocalin</td></tr<>	73		Can f 6	R		inh	Lipocalin		
75 $Iab 3e$ Equ c 3NinhSerum Albumin76MouseMus m 1NinhLipocalin77MuseMus m 1NinhLipocalin78Ara h 1RfStorage protein, 75 globulin (Vicilins)78Ara h 1RfStorage protein, 75 globulin (Vicilins)79Ara h 2RfStorage protein, 75 globulin (Uicilins)79Ara h 3RfStorage protein, 118 globulin (Legumins)80Ara h 3RfStorage protein, 118 globulin (Legumins)81PeanutAra h 6NfStorage protein, 25 albumin82Ara h 6NfStorage protein, 25 albumin83Ara h 6NfStorage protein, 25 albumin84Ara h 8RfPR-1085Ara h 9Rfal.TP, type 186Ara h 9Rfstorage protein, 118 globulin (Legumins)87HazelnutCor a 1.0401Rfstorage protein, 118 globulin (Legumins)88MunutIug r 1RfStorage protein, 118 globulin (Legumins)90WalnutIug r 1Rfstorage protein, 78 globulin (Vicilins)91Jug r 3RfndLTP, type 192Pu du 4RfPordini93AnnodRfStorage protein, 118 globulin (Legumins)94AlmondPu du 601RfStora	74	Horse	Equ c 1	R		inh	Lipocalin		
76MouseMus m 1NinhLipocalin77Ara h 1NfStorage protein, 7S globulin (Vicilins)78Ara h 1RfStorage protein, 7S globulin (Vicilins)79Ara h 2RfStorage protein, 7S globulin (Vicilins)80Ara h 2RfStorage protein, 1S globulin (Legumins)81PeanutAra h 3RfStorage protein, 2S albumin82Ara h 6NfStorage protein, 2S albumin83Ara h 6NfStorage protein, 2S albumin84Ara h 6NfStorage protein, 2S albumin85Ara h 6NfStorage protein, 2S albumin86Ara h 9RfBLTP, type i87HazelnutCor a 1.0401Rf, inh88Cor a 1.0401RfStorage protein, 1IS globulin (Legumins)89HazelnutGor a 9NfStorage protein, 2S albumin90WalnutJug r 1RfStorage protein, 1IS globulin (Legumins)91Jug r 2NCCDfStorage protein, 7S globulin (Vicilins)92Pur du 3Rfnd.TP, type i93AlmondRfStorage protein, 1IS globulin (Legumins)94AlmondRfStorage protein, 1IS globulin (Legumins)95Pur du 6.01RfStorage protein, 1IS globulin (Legumins)96Pur du 6.01R	75	norse	Equ c 3	Ν		inh	Serum Albumin		
77Ara h 1NfStorage protein, 7S globulin (Vicilins)78Ara h 1RfStorage protein, 7S globulin (Vicilins)79Ara h 1RfStorage protein, 7S globulin (Vicilins)80Ara h 2RfStorage protein, 7S globulin (Vicilins)81PeanutAra h 3RfStorage protein, 1IS globulin (Legumins)82Ara h 3RfStorage protein, 1IS globulin (Legumins)83Ara h 6NfStorage protein, 2S albumin84Ara h 6NfStorage protein, 2S albumin85Ara h 8RfPR-1086Ara h 9Rfnd. IP, type I87HazehutCor a 1.0401Rf, inh88Cor a 1.0401RfStorage protein, 1S globulin (Legumins)89HazehutCor a 8Rfnd. IP, type I90WahutJug r 1RfStorage protein, 7S globulin (Vicilins)91Jug r 3Rfnd. IP, type I92Pru du 4RfStorage protein, 7S globulin (Vicilins)93Pru du 60NfStorage protein, 1IS globulin (Legumins)94AlmondFStorage protein, 1IS globulin (Legumins)95Pru du 6.01RfStorage protein, 1IS globulin (Legumins)96Pru du 6.02RfStorage protein, 1IS globulin (Vicilins)97Ana o 1Rf <td< td=""><td>76</td><td>Mouse</td><td>Mus m 1</td><td>Ν</td><td></td><td>inh</td><td>Lipocalin</td></td<>	76	Mouse	Mus m 1	Ν		inh	Lipocalin		
78 79Ara h 1RfStorage protein, 7S globulin (Vicilins) $79$ Ara h 2RfStorage protein, 2S albumin $80$ Ara h 2RfStorage protein, 1S globulin (Ucgumins) $81$ PeanutAra h 3RfStorage protein, 1S globulin (Legumins) $81$ Ara h 3RfStorage protein, 1S globulin (Legumins) $82$ Ara h 6NfStorage protein, 2S albumin $83$ Ara h 6NfStorage protein, 2S albumin $84$ Ara h 6NfStorage protein, 2S albumin $84$ Ara h 8RfPR-10 $85$ Ara h 9Rfind.TP, type I $86$ Ara h 9Rfstorage protein, 2S albumin $87$ HazehutCor a 1.0401Rfind.TP, type I $88$ Cor a 1.0401Rfstorage protein, 1S globulin (Legumins) $89$ HazehutGor a 9NfStorage protein, 1S globulin (Vicilins) $90$ WalnutJug r 2NCCDfStorage protein, 1S globulin (Vicilins) $91$ Jug r 3Rfnsl.TP, type I $92$ Jug r 3Rfstorage protein, 1S globulin (Vicilins) $93$ AlmondRfStorage protein, 1IS globulin (Legumins) $94$ AlmondRfStorage protein, 1IS globulin (Legumins) $95$ Pru du 6.01RfStorage protein, 1IS globulin (Vicilins)	77		Arah 1	Ν		f	Storage protein, 7S globulin (Vicilins)		
79Ara h 2RfStorage protein, 2S albumin80818182838483848485868788888990899191929293949595969798989991919192939495959697989991919192939495959697989991919293949595969798989999919192939494959596969798989999999999999999999991919293949495<	78		Ara h 1	R		f	Storage protein, 7S globulin (Vicilins)		
80 81 82PeanutAra h 3NfStorage protein, 11S globulin (Legumins)81 82Ara h 3RfStorage protein, 2S albumin83Ara h 6NfStorage protein, 2S albumin84Ara h 6NfStorage protein, 2S albumin84Ara h 8RfPR-1085Ara h 9RfnsLTP, type 186Ara h 9RfnsLTP, type 187HazelnutCor a 1.0401Rfstorage protein, 11S globulin (Legumins)88Cor a 1.0401RfnsLTP, type 189MahutJug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 1S globulin (Legumins)91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 6NfStorage protein, 1S globulin (Vicilins)94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 2RfStorage protein, 11S globulin (Legumins)98Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins) <td>79</td> <td></td> <td>Ara h 2</td> <td>R</td> <td></td> <td>f</td> <td>Storage protein, 2S albumin</td>	79		Ara h 2	R		f	Storage protein, 2S albumin		
81 82PeanutAra h 3RfStorage protein, 11S globulin (Legumins)83Ara h 6NfStorage protein, 2S albumin84Ara h 6NfStorage protein, 2S albumin84Ara h 8RfPR-1085Ara h 9RfnLTP, type 186Ara h 9RfnLTP, type 187HazelnutCor a 1.0401Rfstorage protein, 11S globulin (Legumins)88Cor a 9NfStorage protein, 11S globulin (Legumins)89PatherJug r 1RfStorage protein, 2S albumin90WahutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192Pru du 4Rfstorage protein, 11S globulin (Legumins)91Pru du 6.01RfStorage protein, 11S globulin (Legumins)92Pru du 6.01RfStorage protein, 11S globulin (Legumins)93Pru du 6.02RfStorage protein, 11S globulin (Legumins)94AlmondFStorage protein, 11S globulin (Legumins)95Pru du 6.02RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 11S globulin (Legumins)98PuAna o 2RfStorage protein, 11S globulin (Legumins	80		Arah 3	Ν		f	Storage protein, 11S globulin (Legumins)		
82Ara h 6NfStorage protein, 2S albumin83Ara h 6NfStorage protein, 2S albumin84Ara h 8RfPR-1085Ara h 9RfnsLTP, type 186Ara h 9RfnsLTP, type 187HazelnutCor a 1.0401RfnsLTP, type 188Cor a 9NfStorage protein, 11S globulin (Legumins)89Jug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 11S globulin (Legumins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)910FFF<	81	Peanut	Arah 3	R		f	Storage protein, 11S globulin (Legumins)		
83Ara h 6NfStorage protein, 2S albumin84Ara h 8RfPR-1085Ara h 9RfnsLTP, type 186Ara h 9RfnsLTP, type 186Cor a 1.0401Rf, inhPR-1087HazelnutCor a 8RfnsLTP, type 188Cor a 9NfStorage protein, 11S globulin (Legumins)89Jug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 11S globulin (Legumins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100PistachioPis v 3RfStorage protein, 75 globulin (Vicilins)101PistachioPis v 3Rf	82		Ara h 6	Ν		f	Storage protein, 2S albumin		
84Ara h 8RfPR-1085Ara h 9Rfisl.TP. type 186Ara h 9Rf, inhPR-1087HazelnutCor a 1.0401Rf, inhPR-1088Cor a 9NfStorage protein, 11S globulin (Legumins)89Jug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDf91Jug r 3Rfisl.TP. type 192Pru du 3Rfisl.TP. type 193Pru du 4RfProfilin94AlmondPru du 6Nf95Pru du 6NfStorage protein, 11S globulin (Legumins)96Pru du 6.01RfStorage protein, 11S globulin (Legumins)97Aano 1RfStorage protein, 11S globulin (Legumins)98Gashew nutAna o 1RfStorage protein, 11S globulin (Legumins)97Ana o 2RfStorage protein, 11S globulin (Legumins)98OAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100PistachioPis v 3RfStorage protein, 11S globulin (Legumins)101PistachioPis v 3RfStorage protein, 2S albumin102PistachioPis v 3RfStorage protein, 2S albumin <td>83</td> <td></td> <td>Arah 6</td> <td>Ν</td> <td></td> <td>f</td> <td>Storage protein, 2S albumin</td>	83		Arah 6	Ν		f	Storage protein, 2S albumin		
85Ara h 9RfnsLTP, type 186 $R$ $R$ $f$ , inhPR-1087HazelnutCor a 1.0401R $f$ , inhPR-1088 $Cor a 8$ R $f$ nsLTP, type 188 $Cor a 9$ N $f$ Storage protein, 11S globulin (Legumins)89 $Pur ur 1$ R $f$ Storage protein, 2S albumin90 $Walnut$ $Jug r 1$ R $f$ Storage protein, 7S globulin (Vicilins)91 $Jug r 3$ R $f$ nsLTP, type 192 $Pur ur 3$ R $f$ nsLTP, type 193 $Pur ur 4$ R $f$ Profilin94AlmondPru $ur 6$ N $f$ Storage protein, 11S globulin (Legumins)95 $Pur ur 6$ N $f$ Storage protein, 11S globulin (Legumins)96Pru $du 6.01$ R $f$ Storage protein, 11S globulin (Legumins)97 $Ana o 1$ R $f$ Storage protein, 11S globulin (Legumins)98 $Ana o 2$ R $f$ Storage protein, 11S globulin (Legumins)99 $Cashew nut$ $Ana o 2$ R $f$ Storage protein, 11S globulin (Legumins)100 $Ana o 2$ R $f$ Storage protein, 11S globulin (Legumins)101 $Ana o 3$ R $f$ Storage protein, 11S globulin (Legumins)102PistachioPis v 3R $f$ Storage protein, 7S globulin (Vicilins)103Brazil nutBer e 1R $f$ Storage pro	84		Ara h 8	R		f	PR-10		
86 87HazelnutCor a 1.0401Rf, inhPR-1087HazelnutCor a 8RfnsLTP, type 188Cor a 9NfStorage protein, 11S globulin (Legumins)89 90ManutJug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 60.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 11S globulin (Legumins)98Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)101PistachioPis v 3RfStorage protein, 11S globulin (Legumins)102PistachioPis v 3RfStorage protein, 7S globulin (Vicilins)103Brazil nutBer e 1RfStorage protein, 7S globulin (Vicilins)	85		Ara h 9	R		f	nsLTP, type 1		
87HazelnutCor a 8RfnsLTP, type 188Cor a 9NfStorage protein, 11S globulin (Legumins)89Jug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.01RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 11S globulin (Legumins)98Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)101PistachioPis v 3RfStorage protein, 11S globulin (Legumins)102PistachioPis v 3RfStorage protein, 75 globulin (Vicilins)103Brazil nutBer e 1RfStorage protein, 75 globulin (Vicilins)	86		Cor a 1.0401	R		f, inh	PR-10		
88Cor a 9NfStorage protein, 11S globulin (Legumins)89Jug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 7S globulin (Vicilins)98Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100PistachioPis v 3RfStorage protein, 11S globulin (Legumins)101PistachioPis v 3RfStorage protein, 11S globulin (Legumins)103Brazil nutBer e 1RfStorage protein, 2S albumin	87	Hazelnut	Cor a 8	R		f	nsLTP, type 1		
89 90WalnutJug r 1RfStorage protein, 2S albumin90WalnutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 19293Pru du 3RfnsLTP, type 19394Pru du 4RfProfilin94Pru du 60NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 11S globulin (Legumins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100PistachioPis v 3RfStorage protein, 11S globulin (Legumins)101PistachioPis v 3RfStorage protein, 11S globulin (Legumins)103Brazil nutBer e 1RfStorage protein, 2S albumin	88		Cor a 9	Ν		f	Storage protein, 11S globulin (Legumins)		
90WahutJug r 2NCCDfStorage protein, 7S globulin (Vicilins)91Jug r 3RfnsLTP, type 192 $Pru du 3$ RfnsLTP, type 193 $Pru du 4$ RfnsLTP, type 194 $Almond$ $Pru du 6$ NfStorage protein, 11S globulin (Legumins)95 $Pru du 6.01$ RfStorage protein, 11S globulin (Legumins)96 $Pru du 6.02$ RfStorage protein, 11S globulin (Legumins)97 $Pru du 6.02$ RfStorage protein, 11S globulin (Legumins)98 $Ana o 1$ RfStorage protein, 11S globulin (Legumins)99 $Cashew nut$ $Ana o 2$ RfStorage protein, 11S globulin (Legumins)100 $Ana o 3$ RfStorage protein, 11S globulin (Legumins)101 $Pis 43$ RfStorage protein, 11S globulin (Legumins)103 $Brazil nut$ Ber e 1RfStorage protein, 2S albumin	89		Jug r 1	R		f	Storage protein, 2S albumin		
91Jug r 3RfnsLTP, type 192Pru du 3RfnsLTP, type 193Pru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Pru du 6.02RfStorage protein, 7S globulin (Vicilins)98Ana o 1RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100PistachioPis v 3RfStorage protein, 11S globulin (Legumins)103Brazil nutBer e 1RfStorage protein, 2S albumin	90	Walnut	Jug r 2	Ν	CCD	f	Storage protein, 7S globulin (Vicilins)		
92Pru du 3RfnsLTP, type I9394AlmondPru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)9595Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)9697Ana o 1RfStorage protein, 7S globulin (Vicilins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100Ana o 2RfStorage protein, 11S globulin (Legumins)101PistachioPis v 3RfStorage protein, 2S albumin103Brazil nutBer e 1RfStorage protein, 2S albumin	91		Jug r 3	R		f	nsLTP, type 1		
93 94AlmondPru du 4RfProfilin94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Pru du 6.02RfStorage protein, 7S globulin (Vicilins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2NfStorage protein, 11S globulin (Legumins)100Ana o 2RfStorage protein, 11S globulin (Legumins)101Ana o 3RfStorage protein, 11S globulin (Legumins)102PistachioPis v 3RfStorage protein, 7S globulin (Vicilins)103Brazil nutBer e 1RfStorage protein, 2S albumin	92		Pru du 3	R		f	nsLTP, type 1		
94AlmondPru du 6NfStorage protein, 11S globulin (Legumins)9596Pru du 6.01RfStorage protein, 11S globulin (Legumins)9696Pru du 6.02RfStorage protein, 11S globulin (Legumins)9797Ana o 1RfStorage protein, 7S globulin (Vicilins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100Ana o 2RfStorage protein, 11S globulin (Legumins)101Ana o 3RfStorage protein, 11S globulin (Legumins)102PistachioPis v 3RfStorage protein, 2S albumin103Brazil nutBer e 1RfStorage protein, 2S albumin	93	1	Pru du 4	R		f	Profilin		
95Pru du 6.01RfStorage protein, 11S globulin (Legumins)96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 7S globulin (Vicilins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100Ana o 2RfStorage protein, 11S globulin (Legumins)101Ana o 3RfStorage protein, 11S globulin (Legumins)102PistachioPis v 3RfStorage protein, 7S globulin (Vicilins)103Brazil nutBer e 1RfStorage protein, 2S albumin	94	Almond	Pru du 6	Ν		f	Storage protein, 11S globulin (Legumins)		
96Pru du 6.02RfStorage protein, 11S globulin (Legumins)97Ana o 1RfStorage protein, 7S globulin (Vicilins)98Ana o 2RfStorage protein, 11S globulin (Legumins)99Cashew nutAna o 2RfStorage protein, 11S globulin (Legumins)100Ana o 2NfStorage protein, 11S globulin (Legumins)101Ana o 2RfStorage protein, 11S globulin (Legumins)101Ana o 3RfStorage protein, 11S globulin (Legumins)102PistachioPis v 3RfStorage protein, 7S globulin (Vicilins)103Brazil nutBer e 1RfStorage protein, 2S albumin	95	1	Pru du 6.01	R		f	Storage protein, 11S globulin (Legumins)		
97     Ana o 1     R     f     Storage protein, 7S globulin (Vicilins)       98     Ana o 2     R     f     Storage protein, 11S globulin (Uicilins)       99     Cashew nut     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       100     Ana o 2     N     f     Storage protein, 11S globulin (Legumins)       101     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       101     Ana o 3     R     f     Storage protein, 2S albumin       102     Pistachio     Pis v 3     R     f     Storage protein, 7S globulin (Vicilins)       103     Brazil nut     Ber e 1     R     f     Storage protein, 2S albumin	96		Pru du 6.02	R		f	Storage protein, 11S globulin (Legumins)		
98     98     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       99     Cashew nut     Ana o 2     N     f     Storage protein, 11S globulin (Legumins)       100     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       100     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       101     Ana o 3     R     f     Storage protein, 2S albumin       102     Pistachio     Pis v 3     R     f     Storage protein, 7S globulin (Vicilins)       103     Brazil nut     Ber e 1     R     f     Storage protein, 2S albumin	97		Ana o 1	R		f	Storage protein, 7S globulin (Vicilins)		
99     Cashew nut     Ana o 2     N     f     Storage protein, 11S globulin (Legumins)       100     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       101     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       101     Ana o 3     R     f     Storage protein, 2S albumin       102     Pistachio     Pis v 3     R     f     Storage protein, 7S globulin (Vicilins)       103     Brazil nut     Ber e 1     R     f     Storage protein, 2S albumin	98		Ana o 2	R		f	Storage protein, 11S globulin (Legumins)		
100     Ana o 2     R     f     Storage protein, 11S globulin (Legumins)       101     Ana o 3     R     f     Storage protein, 2S albumin       102     Pistachio     Pis v 3     R     f     Storage protein, 7S globulin (Vicilins)       103     Brazil nut     Ber e 1     R     f     Storage protein, 2S albumin	99	Cashew nut	Ana o 2	N		f	Storage protein, 11S globulin (Legumins)		
101   Ana o 3   R   f   Storage protein, 2S albumin     102   Pistachio   Pis v 3   R   f   Storage protein, 7S globulin (Vicilins)     103   Brazil nut   Ber e 1   R   f   Storage protein, 2S albumin	100		Ana o 2	R		f	Storage protein, 11S globulin (Legumins)		
102   Pistachio   Pis v 3   R   f   Storage protein, 7S globulin (Vicilins)     103   Brazil nut   Ber e 1   R   f   Storage protein, 2S albumin	101		Ana o 3	R		f	Storage protein, 2S albumin		
103 Brazil nut Ber e 1 R f Storage protein, 2S albumin	102	Pistachio	Pis v 3	R		f	Storage protein, 7S globulin (Vicilins)		
	103	Brazil nut	Ber e 1	R		f	Storage protein, 2S albumin		
104 Sesame Ses i 1 N f Storage protein, 2S albumin	104	Sesame	Ses i 1	Ν		f	Storage protein, 2S albumin		
105 Say Gly m 4 R f PR-10	105	Sov	Gly m 4	R		f	PR-10		

			Version	1 of the	MeDALL-C	hip
No.	Allergen-source	Allergen	Rec. or natural	CCD	Route of exposure	Function of the protein/Allergen family
106		Gly m 5	Ν		f	Storage protein, 7S globulin (Vicilins)
107		Gly m 6	N		f	Storage protein, 11S globulin (Legumins)
108		Tri a 14	R		f, inh	nsLTP, type 1
109		Tri a 14	R		f, inh	nsLTP, type 1
110		Tri a 19.0101	R		f	Storage protein, Omega 5 gliadin
111		Tri a aA TI	Ν		f, inh	Alpha-Amylase / Trypsin Inhibitors
112		m43	R		f	Fragment of HMW Glutenin Bx7
113		m82	R		f	Fragment of HMW Glutenin Bx7
114		Tri a 36 (191-369)	R		f, inh	Fragment of LMW Glutenin
115		Tri a 36	R		f, inh	LMW Glutenin
116	Wheat	Tri a 37	R		f	Thionins (Alpha-purothionin)
117		Tri a 39 (#10)	R		inh	Serine protease inhibitors
118		Tri a Trx (#37)	R		inh	Oxidoreductases (Thioredoxin)
119		Tri a GST (#38)	R		inh	Glutathione S-transferases
120		Tri a 32 (#112)	R		inh	Peroxiredoxins (1-Cys-peroxiredoxin)
121		Tri a 12 (#123)	R		f, inh	Profilin
122		Tri a DH (#126)	R		inh	Dehydrins
123	Buckwheat	Fag e 2	N		f	Storage protein, 2S albumin
124		Act d 1	N		f	Cysteine protease
125	Kiwi	Act d 2	N		f	Thaumatin-like protein
126	Kimi	Act d 5	Ν		f	Kiwellin
127		Act d 8	R		f	PR-10
128	Apple	Mal d 1	R		f	PR-10
129	Peach	Pru p 1	R		f	PR-10
130	1 Euch	Pru p 3	R		f, skin	nsLTP, type 1
131	Celery	Api g 1	R		f	PR-10
132	Pineapple	Ana c 2.0101	Ν	CCD	f, inh	Bromelain (MUXF3, CCD-Marker)
133		Bos d 4	N		f	Alpha-lactalbumin
134		Bos d 5	N		f	Beta-lactoglobulin
135		Bos d 6	N		f, inh	Serum Albumin
136	Cow's Milk	Bos d 8	d 8 N f Caseins		Caseins	
137		Bos d Lactofer	rin N		f	Transferrin
138		Bos d 4	Ν		f	Alpha-lactalbumin
139		Bos d 5	Ν		f	Lipocalin (Beta-lactoglobulin)

	Version 1 of the MeDALL-Chip								
No.	Allergen-source	Allergen	Rec. or natural	ССД	Route of exposure	Function of the protein/Allergen family			
140		Bos d 6	N		f, inh	Serum Albumin			
141		Bos d 8	N		f	Caseins (α-, β-, κ-casein)			
142		Bos d 9	R		f	Caseins (aS1-casein)			
143		Bos d 10	R		f	Caseins (aS2-casein)			
144		Bos d 11	R		f	Caseins (β -casein)			
145		Bos d 12	R		f	Caseins (ĸ-casein)			
146		Bos d Lactofer	rin N		f	Transferrin			
147		Gal d 1	N		f	Ovomucoid			
148	Egg white	Gal d 2	N		f	Ovalbumin			
149		Gal d 3	N		f	Ovotransferrin (Conalbumin)			
150	Egg yolk/chicken	Gal d 5	N		f, inh	Serum Albumin			
151	Codfish	Gad c 1	R		f	Parvalbumin			
152		Pen m 1	N		f	Tropomyosin			
153	Shrimp	Pen m 2	N		f	Arginine Kinase-like protein			
154		Pen m 4	N		f	Calcium-binding proteins (Sarcoplasmic)			
155		Api m 1	R		syst	Insect venom, Phospholipase A2			
156		Api m 1	R		syst	Insect venom, Phospholipase A2			
157	Bee	Api m 2	R		syst	Insect venom, Hyaluronidase			
158		Api m 4	N		syst	Melittin			
159		Api m 4	N		syst	Melittin			
160		Ves v 1	R		syst	Insect venom, Phospholipase A1			
161	Wasp	Ves v 5	R		syst	Insect venom, Antigen 5			
162		Ves v 5	R		syst	Insect venom, Antigen 5			
163	Paper wasp	Pol d 5	R		syst	Insect venom, Antigen 5			
164		Sem	R		skin, inh				
165		SspB	R		skin, inh				
166		Sec	R		skin, inh				
167		HlgC	R		skin, inh				
168		Tsst	R		skin, inh				
169		Hlb	R		skin, inh				
170	Bacterial Toxins	Plc	R		skin, inh				
171		GlpQ	R		skin, inh				
172		Hla	R		skin, inh				
173		Ssl11	R		skin, inh				
174		Sei	R		skin, inh				
175		Seb	R		skin, inh				
176	Control-protein	MBP	R		-	Maltose binding protein			

Lupinek et al.

		Version 2 of the Me	eDALL-Chip (	exc	hange	d allergens)
164		Ole e 5	R		inh	Superoxide dismutases
165	Olive	Ole e 6	R		inh	Oleaceae group 6
166	Olive	Ole e 8	R		inh	Calcium-binding proteins (Polcalcin)
167		Ole e 10	R		inh	Glycosyl hydrolases
168		Tri a gamma gliadin (GG1)	R		f	Gamma-gliadins
169		Peptide 4	R (peptide)		f	Peptides from gamma-gliadin 1
170		Peptide 18	R (peptide)		f	Peptides from gamma-gliadin 1
171	Wheat	Clone 79	R		f	Gamma-gliadins
172	wneai	Clone 85	R		f	Gamma-gliadins
173		Clone 110	R		f	Gamma-gliadins (Fragment)
174		Tri a bA	R		f	beta-Amylases
175		Avenin-like protein	R		f	Avenin
176	Human	htTG2	R		-	Human tissue transglutaminase

Results of the detection of dilution series of monoclonal IgE and IgG<sub>1</sub>. Dilution series of a monoclonal human IgE specific to Bet v 1 were tested by the MeDALL allergen-chip (left top) or Bet v 1-ImmunoCAP (right). A dilution series of a monoclonal human IgG<sub>1</sub> specific to Phl p 2 was tested by the MeDALL allergen-chip (left bottom). Mean results and standard deviations of triplicate measurements are shown in ISU and ISU-G, respectively, for the microarray and in UA/ml for ImmunoCAP. For each dilution, variation coefficients (%CV) and signal-to-noise ratios (SNR) are shown. Orange boxes at the bottom of every table display results obtained for the sample diluent used for diluting the monoclonal antibodies.

	MeDALL-chip											
	Bet v 1-specific monoclonal IgE											
IgE-conc.	IS	U	%CV	mean	SD	SND	mean	SD				
(ng/ml)	Mean	SD	70C V	%CV 9	%CV	SINK	SNR	SNR				
640	132.02	2.26	1.7			58.5						
320	102.88	16.36	15.9			6.3						
32	23.30	1.83	7.9	7.7	5.5	12.7	22.7	21.2				
3.2	3.12	0.29	9.3			10.8						
1.6	1.67	0.07	4.0			25.1						
0.8	0.97	0.24	24.1			4.1						
0.4	0.51	0.05	9.3			10.7						
0.2	0.28	0.04	12.4	14.1	6.8	8.1	8.5	3.9				
0.1	0.18	0.03	17.5			5.7						
0.05	0.14	0.01	7.1			14.0						
0	0.10	0.01	10.0			10.0						

	Phl p 2-specific monoclonal IgG <sub>1</sub>											
IgG-conc.	ISU-G		%CV	mean	SD	SND	mean	SD				
(ng/ml)	Mean	SD	70C V	%CV	%CV	SINK	SNR	SNR				
1000	181.22	15.12	8.3			12.0						
500	161.01	16.26	10.1	14.9	7.0	9.9	0 1	26				
100	117.53	17.70	15.1		1.9	6.6	8.1	5.0				
10	31.61	8.21	26.0			3.8						
1	4.82	1.25	25.8			3.9	21.0	34.4				
0.4	2.33	0.03	1.2			82.4						
0.3	1.64	0.31	18.6	16.3	10.4	5.4						
0.2	1.11	0.12	10.8			9.3						
0.1	0.59	0.15	25.0			4.0						
0	0.02	0.01	50.0			2.0						

Lupinek et al.

			Bet v 1	Immuno(	CAP						
Bet v 1-specific monoclonal IgE											
IgE-conc.	UA/	ml	9/ CV	mean	SD	SND	mean	SD			
(ng/ml)	Mean	SD	70C V	%CV	%CV	SINK	SNR	SNR			
312.5	91.20	3.35	3.7			27.3					
39.1	10.57	0.29	2.7		1.4	36.6		33.6			
19.5	5.49	0.26	4.7	27		21.4	50.1				
7.8	2.16	0.05	2.1	2.7	1.4	47.1					
3.9	1.14	0.01	0.9			114.0					
2.0	0.63	0.01	1.8			54.3					
1.0	0.36	0.01	3.2			30.9					
0.5	0.22	0.00	0.0	2.3	2.0	-	29.0	2.7			
0.2	0.16	0.01	3.7			27.1					
0	0.09	0.01	6.2			16.2					

Inhibition of IgE binding as measured by microarray and by ImmunoCAP after addition of blocking IgG. Different concentrations of monoclonal human Bet v 1-specific IgE (chim.IgE; concentrations indicated in coloured boxes) were incubated with increasing concentrations of blocking IgG (relative excess compared to IgE and, in parentheses, absolute concentrations of IgG indicated in the left columns) and measured by microarray (left tables) or ImmunoCAP (right tables). Absolute responses after subtraction of background signals are indicated in ISU for the MeDALL allergen-chip and in UA/ml for ImmunoCAP. Relative inhibition was calculated for every IgE-concentration applied (baseline responses indicated in red numbers).

MeDALL-chip									
chim. IgE $\rightarrow$		312.5 ng	/ml						
	Absolute r	esponse	% Inhi	bition					
Blocking IgG↓	Mean	SD	Mean	SD					
0× (0 ng/ml)	92.86	2.72	0.00	2.92					
1× (312.5 ng/ml)	52.15	6.13	43.84	6.60					
10× (3.1 μg/ml)	20.05	1.39	78.41	1.50					
100× (31.2 µg/ml)	2.60	0.42	97.20	0.46					
1000× (312.5 μg/ml)	0.32	0.02	99.66	0.02					

chim. IgE $\rightarrow$	6.25 ng/ml						
	Absolute re	esponse	% Inhibition				
Blocking IgG↓	Mean	SD	Mean	SD			
0× (0 ng/ml)	5.58	0.62	0.00	11.17			
1× (6.25 ng/ml)	5.04	0.46	9.56	8.28			
10× (62.5 ng/ml)	4.56	0.74	18.23	13.27			
100× (625 ng/ml)	1.62	0.20	70.95	3.52			
1000× (6.25 μg/ml)	0.18	0.03	96.83	0.62			
10,000× (62.5 μg/ml)	0.02	0.02	99.58	0.27			

chim. IgE $\rightarrow$	0.8 ng/ml			
	Absolute r	esponse	% Inhibition	
Blocking IgG↓	Mean SD		Mean	SD
0× (0 ng/ml)	0.65	0.10	0.00	16.14
1× (0.8 ng/ml)	0.65	0.06	0.00	8.61
10× (8 ng/ml)	0.68	0.13	-4.64	20.10
100× (80 ng/ml)	0.46	0.14	28.35	22.37
1000× (800 ng/ml)	0.09	0.04	85.57	6.25

chim. IgE $\rightarrow$	0.8 ng/ml			
	Absolute re	esponse	% Inhi	bition
Blocking IgG $\downarrow$	Mean	SD	Mean	SD
10,000× (8 µg/ml)	0.03	0.00	95.88	0.00
100,000× (80 µg/ml)	0.00	0.01	100.00	0.89

Bet v 1 ImmunoCAP					
chim. IgE $\rightarrow$	312.5 ng/ml				
	Absolute response		% Inhibition		
Blocking IgG↓	Mean	SD	Mean	SD	
0× (0 ng/ml)	91.11	3.35	0.00	3.67	
1× (312.5 ng/ml)	83.14	2.30	8.74	2.53	
10× (3.1 µg/ml)	90.91	3.95	0.22	4.34	
100× (31.2 μg/ml)	91.74	5.41	-0.70	5.94	
1000× (312.5 μg/ml)	21.54	1.16	76.36	1.27	

chim. IgE $\rightarrow$	20 ng/ml				
	Absolute re	esponse	% Inhibition		
Blocking IgG $\downarrow$	Mean	SD	Mean	SD	
0× (0 ng/ml)	5.39	0.26	0.00	4.76	
1× (20 ng/ml)	5.59	0.10	-3.71	1.90	
10× (200 ng/ml)	5.84	0.06	-8.34	1.02	
100× (2 µg/ml)	6.05	0.20	-12.11	3.75	
1000× (20 µg/ml)	6.17	0.08	-14.34	1.52	
10,000× (200 μg/ml)	2.29	0.26	57.48	4.83	

Detection of early and late sensitization patterns in birth cohorts. IgE-reactivity profiles of five children from the BAMSE-cohort at the age of 4, 8 and 16years, assessed by the MeDALL allergen-chip are displayed. White boxes show IgE-levels below 0.3ISU, green boxes display IgE-levels 0.3 and <1ISU, 1 and 10ISU or higher than 10ISU (from bright to dark green). The first column shows codes of study participants.

	Allergen	IgE detected at			
Individual		4 years	8 years	16 years	
	Bet v 1	8.35	34.7	59.44	
	Act d 8	0	0.73	0.46	
	Aln g 1	4.28	0	1.63	
	Ara h 8	0.22	2.55	1.82	
	Cor a 1.0401	15.9	11.12	4.04	
	Cor a1.01	1.23	6.61	1.56	
	Mal d 1	23.42	15.3	11.83	
BAMSE 1	Pru p 1	0.97	3.12	4.2	
	Cyn d 1	37.14	97.66	101.27	
	Cup a 1	0	2.38	2.15	
	Phl p 1	81.63	96.09	106.65	
	Phl p 2	0	4.19	15.91	
	Phl p 4	4.8	30.31	4.03	
	Phl p 5	0.19	0.31	16.43	
	Can f 2	3.22	2	1.17	
	Bet v 1	0	38.61	161.55	
	Act d 8	0	0.27	2.75	
	Aln g 1	0	0	0.85	
	Api g 1	0	0	9.53	
	Ara h 8	0	0.72	23.78	
	Cor a 1.0401	0	6.96	77.46	
	Cor a1.01	0		20.11	
	Mal d 1	0	11.73	144.86	
	Pru p 1	0	0	8.43	
BAMSE 2	Bet v 2	0	3.39	3.62	
	Phl p 12	0	1.37	1.24	
	Hev b 8	0	7.05	16.08	
	Mer a 1	0	6.29	2.69	
	Pru du 4	0	2.3	8.41	
	Tri a 12	0	1.45	4.29	
	Can f 5	0	1.53	61.3	
	Can f 5	0	0.36	3.6	
	Can f 6	0	0	44.24	

Lupinek et al.

T., 19-4 J., 1	Allongon	IgE detected at			
Individual	Allergen	4 years	8 years	16 years	
	Cyn d 1	0	25.94	51.17	
	Phl p 1	0	60.17	151.32	
	Phl p 2	0	3.54	8.2	
	Phl p 4	0	0.36	3.79	
	Cup a 1	0	0.57	0.65	
	Equ c 1	0	0	17.53	
	Fel d 1	0	77.94	26.9	
	Fel d 2	0	0.98	1.5	
	Fel d 4	0	0	22.19	
	Mus m 1	0	0	2.99	
	Bet v 1	0	0.15	7.08	
	Ara h 8	0	0	1.02	
	Cor a 1.0401	0	0.11	1.01	
	Mal d 1	0	0.1	1.08	
	Gly m 4	0	0	0.49	
	Cyn d 1	0	0.11	3.9	
BAMSE 3	Phl p 1	0	0.99	27.78	
	Equ c 1	0	0	3.09	
	Can f 5	0	0	0.97	
	Can f 6	0	0.15	5.03	
	Fel d 1	0	0	14.51	
	Fel d 4	0	0.13	29.79	
	Mus m 1	0	0	8.79	
	Cyn d 1	0	0.17	47.92	
	Phl p 1	0	0.19	129.4	
BAMSE 4	Phl p 4	0	0.11	25.78	
	Phl p 5	0	0	74.52	
	Phl p 6	0	0	16.13	
	Bet v 1	0	0	10.85	
	Cor a 1.0401	0	0	5.31	
	Cyn d 1	2.71	51.41	14.89	
BAMSE 5	Phl p 1	16.61	52.08	28.42	
	Phl p 2	0	1.58	3.66	
	Phl p 5	0	21.06	11.85	
	Fel d 1	0	0.14	4.81	