

# An Automated Miniaturized Method to Perform and Analyze Antimicrobial Drug Synergy Assays

Peter Chase,\* Imarhia Enogieru,\* Franck Madoux, Eric Bishop, Jacob Beer, Louis Scampavia, and Timothy Spicer

Lead Identification, Department of Molecular Therapeutics, The Scripps Research Institute, Scripps Florida, Jupiter, Florida.

\*These authors contributed equally to this work.

## ABSTRACT

*In the light of emerging antibiotic resistance mechanisms found in bacteria throughout the world, discovery of drugs that potentiate the effect of currently available antibiotics remains an important aspect of pharmaceutical research in the 21st century. Well-established clinical tests exist to determine synergy in vitro, but these are only optimal for low-throughput experimentation while leaving analysis of results and interpretation of high-throughput microscale assays poorly standardized. Here, we describe a miniaturized broth microdilution checkerboard assay and data analysis method in 384-well plate format that conforms to the Clinical Laboratory and Standards Institute (CLSI) methods. This method has been automated and developed to rapidly determine the synergism of current antibiotics with various beta-lactamase inhibitors emerging from our antimicrobial research efforts. This technique increases test throughput and integrity of results, and saves test compound and labor. We facilitated the interpretation of results with an automated analysis tool allowing us to rapidly qualify inter- and intraplate robustness, determine efficacy of multiple antibiotics at the same time, and standardize the results of synergy interpretation. This procedure should enhance high-throughput antimicrobial drug discovery and supersede former techniques.*

## INTRODUCTION

The global increase in antibiotic resistance is occurring at an alarming rate such that many of the currently available FDA-approved drugs are no longer effective.<sup>1-3</sup> In 2013, the Centers for Disease Control reports that ~2 million illnesses and ~23,000 deaths are caused by antibiotic-resistant bacteria annually in the United States.<sup>4</sup> The majority of these resistant organisms are well characterized, both phenotypically and genotypically, and whole cell antimicrobial testing is routine in clinics around

the world.<sup>3</sup> However, it is clear that even when armed with such information physicians are still challenged with new and emerging resistant strains entering their clinical settings. The European Centre for Disease Prevention and Control reports that between 2009 and 2012 resistance to third-generation cephalosporins in *Klebsiella pneumonia* and *Escherichia coli* has increased significantly within the EU and European Economic Area.<sup>5</sup> This is further complicated by the fact that bacteria acquire resistance through mechanisms such as intra- and interstrain transmission by transposable elements.<sup>3,6</sup> Much of the resistance is conferred to the bacteria by beta-lactamases, enzymes that are effective against the cell wall synthesis inhibitors, commonly referred to as beta-lactam molecules.<sup>3,7</sup> Current combination therapies do exist such as with amoxicillin and clavulanic acid, which inhibit extended-spectrum beta-lactamases, and, when used properly, render these antibiotics effective again.<sup>8</sup> However, this combination (Augmentin; GlaxoSmithKline, Philadelphia, PA) and other novel antibiotic combination therapies, such as NXL104 (AstraZeneca, London, United Kingdom) and ceftazidime (GlaxoSmithKline), are less effective than expected against other strains, in particular, those exhibiting class B metallo-beta-lactamases and class D beta-lactamases, for which there are no clinically relevant inhibitors.<sup>9,10</sup>

Complicating this issue is the fact that at the present time, there is a lack of consistent, rapid, cost-effective, and non-controversial testing methods for predicting synergy of current antibiotic combinations with new antimicrobial agents against novel or unknown clinical isolates. Although strip-based Etests (bioMérieux USA, Durham, NC) and various other methods have been established by Clinical Laboratory and Standards Institute (CLSI) for lower throughput testing on agar plates, in tubes, and in 96-well plates, there are currently no agreed-upon testing and analysis designs for use with novel compounds in checkerboard microplate testing.<sup>11</sup> Comparison studies have been performed and results may vary widely based on the method of interpretation used.<sup>12</sup> Lacking a better method for interpreting results lends itself to controversy and demonstrates the need for a more standardized approach to understanding antibiotic efficacy in combination with current beta-lactamase inhibitors (BLIs).

Currently, many laboratories perform drug discovery-based screens in high-density formats. This allows for rapid production of reliable results. Here, we have miniaturized a standard antimicrobial synergy assay to 384-well plate format (wpf). The automation and instrumentation necessary for this task are fairly common as described in the methods. Implementing synergy assays into 384 wpf results in a major increase in throughput and data acquisition compared to the typical 96 wpf. This results in concomitant difficulty in interpreting results, for which we generated a novel, cus-

tomizable, in-house analysis tool known as the Synergy RunTool.

### MATERIALS AND METHODS

#### Minimum Inhibitory Concentration and Synergy Protocols

We use a broth microdilution technique with further miniaturization. All additions are automated using a Biomek FX (Beckman Coulter, Inc., Brea, CA) and a Multidrop Dispenser with stackers (Titertek Instruments, Inc., Huntsville, AL). All tests were performed using cation-adjusted Mueller-Hinton

**Table 1. Antimicrobial Synergy Assay Protocol**

Step	Plate	Parameter	Value	Description
1	Antibiotic	Dispense CAMHB	35 $\mu$ L	Into each well of rows B–P
2	Antibiotic	Antibiotic addition; twofold, horizontal serial dilution	70 $\mu$ L of highest (antibiotic)	Into each well of row A
3	Assay	Dispense CAMHB	15 $\mu$ L	Into each well of all columns except 3 and 13
4	Assay	Test BLI compound addition; twofold, vertical serial dilution	30 $\mu$ L of highest (test BLI)	Into each well of columns 3 and 13
5	Antibiotic into assay	Transfer antibiotic dilutions into assay plates	15 $\mu$ L into each well of the plate	Each well in the assay plate now contains a unique combination of Test BLI and antibiotic
6		Prepare bacteria to 1E+06 cfu/mL	10 mL of bacterial solution per (1/2) plate	Bacteria prepared externally in a polypropylene tube
7	Assay	Dispense bacteria into plate	30 $\mu$ L in each well of columns 1–23	If testing two different bacteria on the same plate, place bacteria 1 in columns 1–12 and bacteria 2 in columns 13–23
8	Assay	Dispense highest (antibiotic) into background control column 24	30 $\mu$ L into each well	Use the same antibiotic solution used in step 2
9	Assay	Incubate	37°C for 18–20 h	Incubate in ample humidification
10	Assay	Read	OD 600	Determine OD 600 using Tecan Safire II
11	Assay	Data analysis	Synergy RunTool	Determination of synergistic effect of antibiotic and test compound


#### Step Notes

- Using the 384-well liquid handle (aka Multidrop) add media to Corning 384-well translucent plates for the purpose of antibiotic dilution.
- Using the Beckman FX serial dilution technique perform serial twofold dilutions down the rows of the same plate used in step 1.
- Using the 384-well liquid handle (aka Multidrop) add media to Corning 384-well translucent plates for the purpose of test inhibitor dilution.
- Using the Beckman FX serial dilution technique perform serial twofold dilutions across the columns of the same plate used in step 3.
- Transfer samples from the antibiotic source plate prepared in steps 1 and 2 into the assay plate prepared in steps 3 and 4.
- Fresh bacteria is diluted in CAMHB.
- Using the Multidrop transfer bacteria into appropriate wells; column 24 being the exception.
- Prepare column 24 as "No Bacteria" control.
- Plates are covered and placed at 37°C overnight in a humidified chamber.
- Read the plates using ABS 590 modality.
- Use the automated Synergy RunTool to import data, analyze, and archive.

BLI, beta-lactamase inhibitor; CAMHB, cation adjusted Mueller-Hinton broth.

broth (CAMHB) (Part No. 297963; Becton, Dickinson and Co., Franklin Lakes, NJ). *Enterobacter cloacae* BAA-1143, a control strain of bacteria known to produce a high level of AmpC beta-lactamase, was obtained from the ATCC. YMC07/8/B3323 *Acinetobacter* species transformed with VIM-2 was a generous gift from Kyungwon Lee and Yunsop Chong at the Yonsei University College of Medicine, Seoul, South Korea. Test inhibitor NXL104, a known AmpC inhibitor, was synthesized in-house. Ceftazidime is commercially available (Part No. C3809; Sigma-Aldrich, St. Louis, MO). Imipenem is commercially available (Part No. 1337809; USP, Rockville, MD). We used 384-well plates (Part No. 3701; Corning, Inc., Corning, NY) with a final assay volume of 60  $\mu$ L, derived from the combination of diluted antibiotic (15  $\mu$ L), test inhibitor (15  $\mu$ L), and bacteria (30  $\mu$ L). As outlined in *Table 1* and depicted in *Figure 1*, CAMHB is dispensed into two separate plates, named the antibiotic plate and the assay plate. The highest concentration of antibiotic is added to the first row of the antibiotic plate using the Beckman Biomek FX, and antibiotic is titrated down the columns using twofold serial dilutions to yield a 13-point curve. Separately, the highest concentration of test inhibitor is added to columns 3 and 13 of

the assay plate and diluted across the rows of the plate to yield a 10-point curve. This is followed by a 15  $\mu$ L transfer of the now diluted antibiotic from the antibiotic plate into the assay plate, making the volume in each well 30  $\mu$ L. Each well thus contains a unique combination of antibiotic and test inhibitor, which is all tracked and registered through barcode within the Scripps database repository (Assay Explorer; Symyx, Santa Clara, CA). Bacteria, either lab adapted strains obtained from various vendors or clinical isolates obtained from investigators around the world, are diluted from fresh stocks to a working stock concentration of  $1 \times 10^6$  cfu/mL, and 30  $\mu$ L of the bacteria is added to the assay plate using a Multidrop Dispenser (Thermo Scientific, Waltham, MA) contained in a BSL-2 enclosure. Each individual well has a final bacterial concentration of  $5 \times 10^5$  cfu/mL, determined as per CLSI methods.<sup>13</sup> The plates are lidded and placed in a humidified chamber within a 37°C incubator for 18–20 h under aerobic conditions. Turbidity results are obtained on the same plate by measuring at OD 590 nm on a Tecan SpectraFluor Plus Reader (Tecan Systems, Inc., San Jose, CA). Plate data are then analyzed using an in-house derived software application called the Synergy RunTool.



**Fig. 1.** Flow diagram of a 384-wpf synergy assay. Separate plates are used to dilute the test inhibitor/s (assay plate) and antibiotics (antibiotic plate), which are then combined using a 384-well transfer device into the assay plate. Test inhibitors (BLIs) are not added to the first two columns of assay plate, just antibiotic and bacteria, allowing for the MIC of antibiotic alone to be measured while also leaving several wells available with no drug serving as positive controls. Bacteria (30  $\mu$ L) are added to columns 1–23. In column 24, 30  $\mu$ L of the highest concentration of antibiotic but, no bacteria, is added serving as the negative growth control. A color diagram is included along with reference legend to help identify wells of interest. ABX, antibiotic; BLIs, beta-lactamase inhibitors; MIC, minimum inhibitory concentration; wpf, well plate format.

### Synergy RunTool

Minimum inhibitory concentrations (MICs) and synergy results were determined as per CLSI methods by incorporating the Synergy RunTool. The Synergy RunTool facilitates the interpretation of large amounts of data quickly and consistently. A proprietary Excel macro™ spreadsheet was created, in-house, to implement these requirements. A file-cycling system functionality was added to quickly load data files one at a time for viewing at roughly 1-s intervals. For each file loaded, all vital criteria are immediately calculated and a final output is created that easily allows the user to identify compounds of interest. The first step consisted of normalizing the plate to the user-selected high and low controls. In this case, these are bacteria in the presence of dimethyl sulfoxide (DMSO) only and media plus DMSO, that is, no test inhibitor and no bacteria, respectively. The MIC determined by the Synergy RunTool is the lowest concentration of the antimicrobial agent that prevents visible growth of bacteria, which is exactly as defined by CLSI.<sup>13</sup> In addition, the Synergy RunTool mathematically determines a unique value to be applied as the MIC cutoff for each compound, which is designated as 3 × standard deviation (SD) + average (Avg) of the “No Bacteria” control. Any well with an OD 590 value greater than the cutoff value has observable growth, whereas those below the cutoff value are translucent and thus determine a unique value for the MIC. The spreadsheet is then able to parse the sample grid to find all relevant MIC values, including the MIC of antibiotic (ABX) alone as well as the MIC of the test inhibitor alone as defined using the following equations:

MIC of Antibiotic Alone = ABX MIC =  $[ABX]_{\min}$  without test inhibitor (TI), where % Growth Bacteria is  $< Avg + 3 \times SD$  of No Bacteria wells.

MIC of Test Inhibitor = TI MIC =  $[TI]_{\min}$  without ABX, where % Growth Bacteria is  $< Avg + 3 \times SD$  of No Bacteria wells.

MIC of Antibiotic in Combination = ABX MIC + TI =  $[ABX]_{\min}$  in  $[TI]_n$ , where Growth Bacteria is  $< Avg + 3 \times SD$  of No Bacteria wells.

MIC of Test Inhibitor in Combination = TI MIC + ABX =  $[TI]_{\min}$  in  $[ABX]_n$ , where Growth Bacteria is  $< Avg + 3 \times SD$  of No Bacteria wells.

A primary metric for determining a successful experiment is whether the measured ABX MIC is within the expected range for the test bacterial strain. An acceptable  $Z' > 0.5$  is also essential for determining a successful experiment.<sup>14</sup> If these

two criteria were met, data were processed to determine whether the inhibitors synergize with the antibiotic. Importantly, the Synergy RunTool can also identify the MIC of the combination of test BLI compound and ABX, which is the basis of determining synergy. With these values determined, the Synergy RunTool measures the Total Fractional Inhibitory Concentration ( $\Sigma$ FIC), which is a standard guide for predicting if a compound has a synergistic effect.<sup>15</sup> This is calculated with the following equation:

$$\Sigma FIC = FIC \text{ Antibiotic} + FIC \text{ Test Inhibitor}$$

where:

$$FIC \text{ Antibiotic} = \frac{MIC \text{ of Antibiotic in Combination}}{MIC \text{ of Antibiotic Alone}}$$

$$FIC \text{ Test Inhibitor} = \frac{MIC \text{ of Test Inhibitor in Combination}}{MIC \text{ of Test Inhibitor Alone}}$$

A  $\Sigma$ FIC value is calculated for every well in the plate because every well has a different combination of antibiotic and test inhibitor concentrations. The test inhibitor is considered to have a synergistic effect when  $\Sigma$ FIC  $\leq 0.5$  or the value that ultimately equates to the minimum test inhibitor concentration (Min [Test Inhibitor] Restore) necessary to cause a fourfold decrease in ABX MIC of the tested strain. The compound is considered indifferent with the antibiotic when the  $\Sigma$ FIC is  $> 0.5$  to  $< 2$  and antagonistic when the  $\Sigma$ FIC is  $\geq 2$ .<sup>15</sup> Figure 2 shows a populated RunTool chart that shows the calculated  $\Sigma$ FICs for each well used in a test case.

Describing a test inhibitors' ability to synergize or potentiate the effect of an antibiotic in this manner can be difficult to translate verbally as well as numerically. To facilitate communication of results to investigators and in particular medicinal chemists, we added a feature to help analyze synergy, which we term potentiation. That's the ability of the test inhibitor, which has less effect on its own, to potentiate or allow an antibiotic to be more effective than when not in combination. In this sense, the effect of simultaneous application of antibiotics and test inhibitor is quantified according to the following equation:

$$\text{Potentiation} = \frac{MIC \text{ antibiotic in resistant bacteria}}{MIC \text{ of antibiotic in the combination with test inhibitor in resistant bacteria}}$$

By definition, an MIC can vary by twofold day to day but not more. With that and similar to the  $\Sigma$ FICs method, a result of a fourfold increase is the minimum threshold to indicate a





synergistic effect. A result of full restoration of antibiotic susceptibility is dependent on the change of the MIC found in resistant versus nonresistant bacteria with the MIC of antibiotic in combination with test inhibitor in resistant bacteria reverting to an MIC found for the antibiotic alone in the nonresistant bacteria.

## RESULTS AND DISCUSSION

By implementing the 384 wpf, we enable the determination of both MIC and synergy on two bacterial strains and at the same time for two antibiotics and two test inhibitors all on the same microtiter plate. This provides 10 twofold serial dilutions of test inhibitor and 13 for the antibiotic for 130 combinations in one-half of a plate or 260 combinations per plate. This is in addition to determining the MICs and  $Z'$  from controls on the same plate. All are done while conforming to the CLSI-recommended procedures.<sup>13</sup> Bacterial strains tested may be laboratory-adapted strains or clinical isolates. The only requirement is that one uses isolates that grow according to CLSI procedures. In addition, the number of serial dilution points for test inhibitor and antibiotic is predicated by the fact that one needs to test and observe at least a fourfold change in MIC in order for synergy to be declared for a test inhibitor and antibiotic combination. In doing so, one can determine synergy with as few as three dilution points, but the number of dilutions to determine restoration or full potentiation may vary from strain to strain and antibiotic tested. As shown in *Figure 2*, we test and identify the concentration required to restore the antibiotic to full effect against the bacteria, a result that would be missed otherwise. Obviously, this affects the throughput and amount of information able to be obtained on one plate; a situation that is flexible but must be assessed and designed based on user input and their knowledge of the bacteria resistance.

Synergy testing requires the ability to discern the effects of, in this case, two drugs in combination versus alone. Previous methods of MIC measurement (single drug testing) and synergy testing included low-throughput testing procedures utilizing test tubes and 96 wpf.<sup>16–18</sup> These methods are more time, labor, and resource intensive and lack uniformity in analysis. The miniaturization to 384 wpf facilitates faster data collection while still conforming to CLSI standards. Traditional tube synergy testing requires a minimum of 36 tubes to test 25 combinations in a 1 mL volume to achieve reliable interpretation of synergy.<sup>13,19</sup> A 96-wpf synergy assay achieves 77 combinations per plate in 0.2 mL volume/well.<sup>13</sup> In contrast to both methods, the 384-wpf synergy test is able to obtain 260 test inhibitor combina-

tions per plate using 0.06 mL/well. Both tube and 96-wpf formats dictate lower throughput and less quantitative interpretation, whereas our method allows for the least amount of materials, the smallest assay volume, and a rapid quantitative assessment. The Synergy RunTool generates all the following data in a matter of seconds (*Fig. 2*): (1) raw data import; (2) SD; (3) median and average value determination of all sample and control regions; (4) normalization of all wells to user selected high and low controls; (5)  $Z'$ ; (6) signal-to-noise ratio; (7) %CV of the “No Bacteria” control (Percent Coefficient of Variation derived by  $[SD/Avg] \times 100$ ); (8) MIC cutoff based on the “No Bacteria” control ( $Avg + 3 \times SD$ ); (9) MIC for every concentration of test inhibitor ( $n = 10$ ) and ABX ( $n = 13$ ); (10)  $\Sigma$ FIC of every possible combination of ABX and test inhibitor; (11) resistant versus nonresistant fold difference for ABX and test inhibitor MICs; (12) minimum test inhibitor and ABX concentrations necessary for 2 $\times$ , 4 $\times$ , 8 $\times$ , and 16 $\times$  fold increase in MIC versus control; and (13) minimum BLI and ABX concentrations necessary to restore resistant bacteria MICs to nonresistant bacteria MICs. From these statistics, we found the following to be the most important:  $Z'$  determination using “No Drug” and “No Bacteria” (assay quality,  $Z' \geq 0.5$ ); ABX MIC-fold difference between nonresistant and resistant bacterial strains (assay quality,  $Fold \geq 8$ );  $Avg + 3 \times SD$  of “No Bacteria” region (determines MIC cutoff through mathematical interpretation confirming no-visible growth);  $\Sigma$ FIC at  $[BLI]_{MIN}$  needed to restore the ABX MIC to its full potential as determined in nonresistant bacteria. These values are essential as they create a standard, noncontroversial method of interpreting quality of synergy assays and results of assays.


To compare our methods with current research, the compound NXL104, a novel but well-known B-lactamase inhibitor, was tested in synergy with ceftazidime, a third generation cephalosporin as shown in *Figure 2*. In this figure, the test inhibitor (BLI) is NXL104 for both the left and right side of the plate, the antibiotic tested in columns 1–12 is ceftazidime, and the antibiotic tested in columns 13–23 is imipenem. ABX MIC is measured in columns 2 and 23. The Synergy RunTool estimated an MIC of 192  $\mu$ g/mL (*i.e.*, the concentration of antibiotic in row A is 192  $\mu$ g/mL) for ceftazidime against ATCC BAA1143; *E. cloacae* positive for AmpC resistance.<sup>20</sup> In column 3, where the concentration of test compound is highest, the MIC of ceftazidime drops to 0.75  $\mu$ g/mL (twofold titrations down the plate). Livermore *et al.* synergy tested *Enterobacter* spp. with AmpC resistance and found that 80% of the isolates had an MIC of 256  $\mu$ g/mL.<sup>17</sup> CLSI predicts the actual MIC for dilution test methods to be

within a twofold range of that measurement (128–256  $\mu\text{g}/\text{mL}$ ). Therefore, the Synergy RunTool's calculation of 192  $\mu\text{g}/\text{mL}$  is within the expected CLSI range.<sup>13</sup> Livermore *et al.* also measured the MIC for the combination of NXL104/ceftazidime when NXL104 remained constant at 4  $\mu\text{g}/\text{mL}$ ; it ranged from 0.25 to 1  $\mu\text{g}/\text{mL}$ .<sup>17</sup> When tested using our 384-well methods and analyzed with the Synergy RunTool, the highest concentration tested of NXL104 was 50  $\mu\text{M}$  or 1.44  $\mu\text{g}/\text{mL}$ . At that concentration, the MIC of ceftazidime was 0.75  $\mu\text{g}/\text{mL}$ , also within the range measured and suggested by Livermore *et al.* We also show that imipenem, as tested on the right half of the plate in *Figure 2*, is an antibiotic already known to be potent against AmpC-producing strains such as BAA-1143, and as such demonstrated the expected result in terms of potency and synergy with NXL104.<sup>21</sup> In this case, it also serves as a reference standard against BAA-1143 bacteria and hence a baseline result for comparison of NXL104's ability to fully restore (potentiate) ceftazidime's efficacy. As shown with data from the experiment above, the Synergy RunTool generates a visual "heat" map (*Fig. 2B*), demonstrating that as test inhibitor concentration decreased along with a concomitant decrease in the ABX MIC, one creates a pattern of steps. The statistical information generated, as seen in *Figure 2E*, includes the MIC of antibiotic as well as the concentration of test inhibitor required to decrease, potentiate, and the MIC of antibiotic by 2-, 4-, 8-, and 16-fold. Taken together, we show NXL104 in combination with ceftazidime to demonstrate the ability of our method to determine the fold reduction in MIC, which can also be taken as

the amount of potentiation a test inhibitor confers to an antibiotic to restore its function against the bacteria being tested.

Compared to existing literature, our miniaturized broth microdilution checkerboard assay produces accepted results, whereas, importantly, the Synergy RunTool accurately interprets those results.

In addition to validating our method using the known synergistic pair described above, we provide further substantiation of our methods ability to determine synergy with different strains of bacteria and various combinations of test inhibitor and imipenem. This work was done in support of medicinal chemistry to determine the potential of test inhibitors to restore antibiotic efficacy against metallo beta-lactamases strains, which ultimately contributed to the identification of two molecular probes, ML121 and ML302.<sup>22,23</sup> We tested various lab-adapted and clinical isolates of bacteria from around the world versus numerous test inhibitors (138) in combination with imipenem. In particular, we tested IMP1 and VIM2 metallo-beta-lactamase-expressing strains, including *E. coli*, *Acinetobacter baumannii*, and *Pseudomonas aeruginosa*. The fold potentiation of test inhibitor was different for each test strain and was reported accordingly. The fold restoration can be predicated using nonmetallo-enzyme transformed or wild-type bacteria of the same species, which are known to be sensitive to imipenem. This was done, in particular, for the lab-adapted and transformed strains; e.g., BL21 *E. coli* (wild type) versus BL21 VIM2 *E. coli*. *Figure 3* represents an example of the outcome found when ML302



**Fig. 3.** (A) The statistics generated using the Synergy RunTool corresponding to the data generated as part of the discovery of probe ML302 (CID53362017). The bacteria tested was clinical isolate YMC07/8/B3323, an *Acinetobacter* spp. expressing VIM2. Fold reduction of the MIC of imipenem by ML302 is shown against this strain as 2 $\times$ , 4 $\times$ , 8 $\times$ , and 16 $\times$  fold. (B) These data were then archived in the Scripps database and also published to the PubChem web site (<https://pubchem.ncbi.nlm.nih.gov/bioassay/624096#section=BioAssay-Target>), which is available to the public. Note the maximum potentiation is listed for this and all other compounds published to the link above.



(CID 53362017) was tested in combination with imipenem to determine its effect on the clinical isolate YMC07/8/B3323; an *Acinetobacter* species transformed with VIM-2. A 16-fold potentiation of imipenem's effect was found at 6.25  $\mu$ M test concentration of ML302. Ultimately, ML302 was able to achieve a 32-fold potentiation of the effect of imipenem in this strain; correlating to a 32-fold reduction in its MIC (Fig. 3B). All synergy data and further information regarding technique and results can be found on the publically available PubChem BioAssay web site ([www.ncbi.nlm.nih.gov/pcassay](http://www.ncbi.nlm.nih.gov/pcassay)) using the following record identifiers; AID 463099, 504620, 624080, 624082, 624095, and 624096.

Notably, the methods described here utilize more than one plate to process the synergy assay, which does require cost and time. The automated equipment we used may not be widely available. However, the gain in efficiency when investigating large numbers of test inhibitors over time, such as that required for medicinal chemistry efforts should provide the impetus to move in this direction. In addition, acoustic transfer devices should handle this type of compound combination procedure and may limit the number of plates and manipulations with tips and other plastic ware. However, this too requires more expensive equipment and skilled engineers and in some cases very specialized source plates.

This methodology has been successfully applied to over 138 separate experiments in our lab, many times with multiple bacteria and test inhibitor/antibiotic combinations, equating to 465 synergy test results. We have successfully moved away from lower throughput methods, thus allowing more compounds to be tested simultaneously using reduced resources while generating a larger matrix of data. The creation of the Synergy RunTool easily handles this increased volume of data, which allows for standardization of results' dissemination, aiding communication between chemists and microbiologists in a manner that ensures results are properly archived, accurate, and reproducible. However, future testing using acoustic transfer devices will be done to compare and contrast the outcomes not only for antibiotic combinations but also for drug-drug combination and repurposing studies.

## ACKNOWLEDGMENTS

The National Institutes of Health (U54MH074404) supported this work. We thank Peter Hodder for his guidance during this effort. Pierre Baillargeon and Lina DeLuca (Lead Identification Division, Translational Research Institute, Scripps Florida) are thanked for their assistance with compound management. We thank Yang-Bo Feng for the synthesis of NXL104. Kyungwon Lee and Yunsop Chong are thanked for their donation of YMC07/8/B3323.

## DISCLOSURE STATEMENT

No competing financial interests exist.

## REFERENCES

- Bush K: Antibacterial drug discovery in the 21st century. *Clin Microbiol Infect*. 2004;10 Suppl 4:10-17.
- Blizzard TA, Chen H, Kim S, et al.: Side chain SAR of bicyclic beta-lactamase inhibitors (BLIs). 1. Discovery of a class C BLI for combination with imipenem. *Bioorg Med Chem Lett* 2010;20:918-921.
- Walsh TR, Toleman MA, Poirel L, Nordmann P: Metallo-beta-lactamases: the quiet before the storm? *Clin Microbiol Rev* 2005;18:306-325.
- CDC: *Antibiotic resistance threats in the United States, 2013*. In: Services USDoHaH (ed.), pp. 1-114. CDC, Atlanta, Georgia, 2013.
- European Centre for Disease Prevention and Control: *Summary of the Latest Data on Antibiotic Resistance in the European Union*, pp. 1-10, Stockholm, 2013.
- Normark BH, Normark S: Evolution and spread of antibiotic resistance. *J Intern Med* 2002;252:91-106.
- Bebrone C: Metallo-beta-lactamases (classification, activity, genetic organization, structure, zinc coordination) and their superfamily. *Biochem Pharmacol* 2007;74:1686-1701.
- Paterson DL, Bonomo RA: Extended-spectrum beta-lactamases: a clinical update. *Clin Microbiol Rev* 2005;18:657-686.
- Livermore DM, Mushtaq S, Warner M, Miossec C, Woodford N: NXL104 combinations versus Enterobacteriaceae with CTX-M extended-spectrum beta-lactamases and carbapenemases. *J Antimicrob Chemother* 2008;62:1053-1056.
- Miossec C: NXL104  $\beta$ -lactamase inhibitor. Paper presented at the Challenge of Antibacterial Drug Development, San Diego, CA, 2007.
- Jorgensen JH, Ferraro MJ: Antimicrobial susceptibility testing: a review of general principles and contemporary practices. *Clin Infect Dis* 2009;49:1749-1755.
- Bonapace CR, Bosso JA, Friedrich LV, White RL: Comparison of methods of interpretation of checkerboard synergy testing. *Diagn Microbiol Infect Dis* 2002;44:363-366.
- CLSI: *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically*. 7th ed., Clinical Laboratory Standards Institute, Wayne, PA, 2006.
- Zhang JH, Chung TD, Oldenburg KR: A simple statistical parameter for use in evaluation and validation of high throughput screening assays. *J Biomol Screen* 1999;4:67-73.
- Orhan G, Bayram A, Zer Y, Balci I: Synergy tests by E test and checkerboard methods of antimicrobial combinations against *Brucella melitensis*. *J Clin Microbiol* 2005;43:140-143.
- Moody J: *Synergism Testing: Broth Microdilution Checkerboard and Broth Macrodilution Methods*. 3rd ed. ASM Press, Washington, DC, 2010.
- Livermore DM, Mushtaq S, Warner M, et al.: Activities of NXL104 combinations with ceftazidime and aztreonam against carbapenemase-producing Enterobacteriaceae. *Antimicrob Agents Chemother* 2011;55:390-394.
- Bonapace CR, White RL, Friedrich LV, Bosso JA: Evaluation of antibiotic synergy against *Acinetobacter baumannii*: a comparison with Etest, time-kill, and checkerboard methods. *Diagn Microbiol Infect Dis* 2000;38:43-50.
- Minond D, Saldanha SA, Subramaniam P, et al.: Inhibitors of VIM-2 by screening pharmacologically active and click-chemistry compound libraries. *Bioorg Med Chem* 2009;17:5027-5037.
- Gootz TD, Sanders CC, Goering RV: Resistance to cefamandole: derepression of beta-lactamases by cefoxitin and mutation in *Enterobacter cloacae*. *J Infect Dis* 1982;146:34-42.
- Jacoby GA: AmpC beta-lactamases. *Clin Microbiol Rev* 2009;22:161-182, Table of Contents.
- Minond D, Saldanha SA, Spicer T, et al.: *HTS Assay for Discovery of Novel Metallo-Beta-lactamase (MBL) Inhibitors*. Probe Reports from the NIH Molecular Libraries Program, Bethesda, MD, 2010.



23. Spicer T, Minond D, Enogieru I, et al.: *ML302, a Novel Beta-Lactamase (BLA) Inhibitor*. Probe Reports from the NIH Molecular Libraries Program, Bethesda, MD, 2010.

Address correspondence to:  
*Timothy Spicer*  
*Lead Identification*  
*Department of Molecular Therapeutics*  
*The Scripps Research Institute*  
*Scripps Florida*  
*130 Scripps Way*  
*Jupiter, FL 33458*  
*E-mail: spicert@scripps.edu*

#### Abbreviations Used

%CV = percent coefficient of variation  
[BLI]<sub>MIN</sub> = minimum BLI compound concentration  
ABX = Antibiotic  
Avg = Average  
BLI(s) = beta-lactamase inhibitor(s)  
CLSI = Clinical and Laboratory Standards Institute  
DMSO = dimethyl sulfoxide  
MIC = minimum inhibitory concentration  
SD = standard deviation  
TI = test inhibitor  
wpl = well plate format  
ΣFIC = Total Fractional Inhibitory Concentration