



HHS Public Access

Author manuscript

J Biolom Screen. Author manuscript; available in PMC 2016 August 06.

Published in final edited form as:

J Biolom Screen. 2015 March ; 20(3): 341–349. doi:10.1177/1087057114557233.

ERK and β -arrestin interaction: a converging-point of signaling pathways for multiple types of cell-surface receptors

Haifeng Eishingdrelo^{1,*}, Wei Sun², Hua Li¹, Li Wang¹, Alex Eishingdrelo¹, Sheng Dai², John C. McKew², and Wei Zheng^{2,*}

¹BioInvenu Corp. East Hanover, NJ 07936, United States

²National Center for Advancing Translational Sciences, National Institutes of Health, Bethesda, MD 20892, United States

Abstract

β -arrestin, a signal adaptor protein, mediates intracellular signal transductions through protein-protein interactions by bringing two or more proteins in proximity. Extracellular signal-regulated kinase (ERK), a protein kinase in the family of mitogen-activated protein kinases (MAPKs), is involved in various receptor signal pathways. Interaction of ERK with β -arrestin or formation of ERK/ β -arrestin signal complex occurs in response to activation of a variety of cell-surface receptors. The ERK/ β -arrestin signal complex may be a common transducer to converge a variety of extracellular stimuli to similar downstream intracellular signaling pathways. By using a cell based protein-protein interaction LinkLight assay technology, we demonstrate a direct interaction between ERK and β -arrestin in response to extracellular stimuli, which can be sensitively and quantitatively monitored. Activations of G-protein coupled receptors (GPCRs), receptor tyrosine kinases (RTKs) and cytokine receptors promote formation of the ERK/ β -arrestin signal complex. Our data indicate that the ERK/ β -arrestin signal complex is a common transducer participated in a variety of receptor signaling pathways. Furthermore, we demonstrate that receptor antagonists or kinase inhibitors can block the agonist induced ERK and β -arrestin interaction. Thus, the ERK/ β -arrestin interaction assay is useful for screening of new receptor modulators.

Keywords

β -arrestin; extracellular signal-regulated kinases; ERK; protein-protein interactions

Introduction

β -arrestin is a signal adaptor protein that brings its interacting partners together to form signal complexes.¹ β -arrestin had been initially considered as G-protein signaling terminator as it promotes GPCR desensitization and internalization. Recently, much broader roles of β -

*To whom correspondence should be addressed: Haifeng Eishingdrelo, Ph.D., BioInvenu Corp., 50 Williams Parkway, East Hanover, NJ 07936, USA. Tel. 973-585-6777; haifeng.eishingdrelo@bioinvenu.com. Wei Zheng, Ph.D., NCATS/NIH, 9800 Medical Center Drive, MSC: 3375, Bethesda, MD 20892, USA. Tel. 301-217-5251; wzheng@mail.nih.gov.

Declaration of Conflicting Interests

H.E. has the commercial interests of the patented LinkLight technology. The other authors have declared that no competing interests exist.

arrestin have been recognized. In addition to interacting with GPCRs, β -arrestin has been indicated to participate in signal transductions of various receptors including receptor tyrosine kinases (RTKs), ion channels and cytokine receptors.² Agonist binding to RTKs, such as epidermal growth factor receptor (EGFR) and fibroblast growth factor receptor (FGFR), leads to activation of β -arrestin-dependent MAPK signaling pathways.^{3, 4} Cytokines, such as TNF α and IL1 β , also activate signaling pathways involving β -arrestins.^{5, 6}

ERK, a MAPK kinase family member, is also one of β -arrestin mediated signaling components.⁷ GPCRs can activate ERK via G-protein-dependent pathways and/or β -arrestin-dependent pathway. G-protein-dependent ERK activation results in the translocation of active ERK proteins into the nucleus, whereas the arrestin-dependent ERK activation retains active ERK proteins in the cytosol. The onset ERK activation and duration of active ERK are also different. G-protein-dependent ERK activation is transient and occurs within a few minutes, while β -arrestin-dependent ERK activation is slower in onset and lasts approximately 60 minutes.⁸⁻¹⁰ Activation of GPCRs is traditionally measured by G-protein-dependent second messengers including cAMP, IP₃, and intracellular calcium mobilization. G-protein-independent signaling pathways for GPCRs have been exploited for drug discovery using β -arrestin based assays.^{11, 12} Several biased GPCR ligands, which selectively activate one pathway over the others, have been identified using β -arrestin based assays.^{13, 14} The biased ligands with increased target specificity and reduced off-target effects may have great potential to become next generation GPCR drugs.

Crosstalk between different types of receptors has been recognized.¹⁵ Shared signal complexes may be the basis of cross-communication between different types of receptors that enable cells to integrate a multitude of signals from the environment. Since both β -arrestins and MAPKs commonly participate in the downstream signaling pathways of different types of receptors, the β -arrestin/ERK signal complex could be one of the signal integration nodes. However, little is known about the regulation of β -arrestin-dependent ERK signaling or the effect of ERK/ β -arrestin interaction on receptor signaling pathways.

Although association of β -arrestins with various kinases has been described in literature, a robust cell-based assay to assess the ERK/ β -arrestin interaction directly is not available. Previously, we have developed a cell-based protein-protein interaction LinkLight assay for assessing GPCR and β -arrestin interaction.¹² Here we have applied the LinkLight technology to determine ERK and β -arrestin interaction. We have found that the technology is capable of capturing transient, phosphorylation-regulated protein-protein interactions and generating stable signals without involving reporter gene transcription and translation. Ability to assess transient protein-protein interactions in cells is particularly useful since many protein interactions are regulated by phosphorylation and dephosphorylation in response to environmental stimuli. The interaction of ERK/ β -arrestin is mediated by activation of GPCRs, RTKs and cytokine receptors. Receptor antagonists and kinase inhibitors can block the interaction. Our data demonstrate that the ERK/ β -arrestin signal complex is a common mediator in signaling pathways of diverse cell-surface receptors.

Materials and Methods

Materials

The compounds and chemicals were purchased from Sigma-Aldrich and Tocris Biosciences.

Cell line and cell culture

U2OS cells were purchased from ATCC and were cultured with McCoy's 5A medium (Gibco catalog # 16600-082) supplemented with 10% FBS (Gibco catalog # 26140-079) and 1X Pen/Strep (Gibco catalog # 15070). Cells were cultured in 37°C incubator with 5% CO₂. Cell culture medium was replaced every three to four days and cells were passed at 90% confluence.

Plasmid construction and stable cell line generation

Full-length cDNA of human ERK2 without stop codon was subcloned in frame with the TEV protease vector described in our previous publication.¹² The β -arrestin-2-permuted luciferase (β -arr-2-pLuc) expression plasmid was also described in the previous publication.¹² Transfection of U2OS cells was performed with Lipofectamine 2000 according to the manufacturer's instructions (Life Technology, catalog # 11668-027). Monoclonal cell lines stably expressing both ERK2-TEV and β -arr-2-pLuc fusion protein were selected with 400 μ g/ml G418 (Invitrogen/Life Technology, catalog # 10131-027) and 300 μ g/ml Hygromycin B (Life Technology, catalog # 10687-010). Multiple stable clones were selected for evaluation.

Luciferase assay

The LinkLight U2OS cells stably expressing both ERK2-TEV and β -arr-2-pLuc were seeded into a 384-well white plate (Becton Dickinson) at 20,000 cells per well with 40 μ l regular medium. Cells were cultured overnight. The next day, the culture medium was replaced with 20 μ l serum-free DMEM medium followed by addition of 5 μ l agonists or stimulants. Agonists were added as serial dilutions with a serum-free DMEM medium from 10 mM DMSO stocks. For antagonist inhibition assays, the medium was replaced with 15 μ l serum-free DMEM medium. Antagonists or inhibitors was added at 5 μ l/well and incubated with the cells for 15 min. Then, the stimulant or agonist was added at 5 μ l/well and incubated for 90 min. An equal volume of Bright-Glo (Promega) or One-Glo (Promega) or BriteLite (Perkin Elmer) or NeoLite (Perkin Elmer) luciferase detection reagent was added for a 10 minute incubation (following manufacturer's instructions). The luminescence signals were determined using a luminiscence plate reader. Alternatively, after 90 min compound incubation, the culture medium was replaced with 10 μ l of luciferase detection reagent and luminescence signals were recorded.

Measurement of FBS response

The LinkLight U2OS cells were seeded in a white 384-well plate and cultured overnight with medium containing 10% FBS. The cell culture medium was replaced with fresh media containing FBS concentrations in a serial dilution (final concentration of FBS from 20 to

0.01%). After 90 min incubation, the luciferase detection reagent was added and luminescence signals were determined.

Measurement of GPCR agonist responses

The LinkLight U2OS cells were seeded in white 384-well plates. After overnight culture, the culture medium was replaced with a serum-free DMEM and a serial dilution of GPCR ligands were added to the cells. After 90 min incubation, a luciferase detection reagent was added to the cells and luminescent signals were recorded. The antagonist effects were measured after cells were pre-incubated with the compounds for 15 min followed by stimulation with 1 μ M agonists for 90 min before proceeding to the luciferase assay.

Measurement of RTK responses

The LinkLight U2OS cells were cultured in a white 384-well plate for overnight. The culture medium was then replaced with a serum-free DMEM. Solutions of EGF or FGF-basic in a serial dilution were added to cells. After 90 min incubation with the agonists, the luciferase assay was performed. For antagonist and kinase inhibitor assays, cells were pre-incubated with antagonists or kinase inhibitors for 15 min followed by agonist stimulation for 90 min then the luciferase assay.

Determination of cytokine receptor response

The LinkLight U2OS cells were cultured in a white 384-well plate for overnight. The culture medium was then replaced with 20 μ l serum-free DMEM. A serial dilution of 5 μ l TNF α or IL1 β were added to the cells for 90 min incubation, then followed by the luciferase assay. For the measurement of response to kinase inhibitor, the cells were treated with staurosporine for 15 min prior to incubation with 25 ng/ml of TNF α or 50 ng/ml IL1 β for 90 min. The luciferase assay was then performed.

Data analysis

Concentration-response curves were generated by using the Prism software (Graphpad Software, Inc. San Diego, CA). All values were expressed as the mean \pm SD (n=3).

Results

Cell-based assay to measure the ERK/ β -arrestin interaction

Upon activation, ERK is phosphorylated and released from the RAF/MEK/ERK signal complex.¹⁶ The phosphorylated ERK binds to other cellular proteins and translocates to new cellular destinations. Binding of ERK to β -arrestin results in cytosolic retention of the ERK/ β -arrestin complex, while its association with nuclear importins results in translocation of ERK into nucleus. To assess ERK and β -arrestin interaction, we constructed the ERK-TEV expression plasmid by fusing ERK2 with a TEV protease at its C-terminus. U2OS cells were co-transfected with the ERK-TEV plasmid and the β -arrestin-2 linked permuted luciferase (β -arrestin-pLuc) expression plasmid. Upon the ERK/ β -arrestin interaction, the proximity of TEV and permuted luciferase leads to cleavage of the permuted luciferase that contains a TEV cleavage site. The cleaved luciferase fragment spontaneously refolds to form an active

luciferase, driven by the self-complementation force. Activity of the reconstituted luciferase can be detected by luciferase detection reagents (Fig. 1A).

To determine if the ERK/ β -arrestin interaction is inducible in the U2OS cell line, we first looked the response to fetal bovine serum (FBS), which contains a number of growth factors, hormones, and ligands. Many of these factors and hormones bind and activate cell-surface receptors. The activated receptors transduce extracellular stimuli through phosphorylation of non-receptor kinases, such as ERK. We found that treatment of cells co-expressing ERK-TEV and β -arrestin-pLuc with FBS dose-dependently increased luciferase signals (Fig. 1B). We also looked at the time course of FBS stimulation. ERK/ β -arrestin interaction signals were reached plateau around 120 min. The result indicated the interaction between ERK-TEV and β -arrestin-pLuc or formation of the ERK/ β -arrestin signal complex. The data also suggested that in the regular cell culture condition with 10% FBS, ERK could be constitutively phosphorylated, or in a dynamic equilibrium between phosphorylation and dephosphorylation states. The U2OS cells transiently expressing only one plasmid (β -arrestin-pLuc) did not show FBS stimulated signals (Fig. S1).

Since cellular signaling pathways involve highly coordinated protein-protein interactions that are regulated by protein phosphorylation and dephosphorylation, we looked whether the ERK and β -arrestin interaction could be blocked by kinase inhibitors. We found that staurosporine (pan-kinase inhibitor) dose-dependently inhibited ERK/ β -arrestin interaction signals (Fig. 1C), indicating that the ERK/ β -arrestin interaction is mediated by protein phosphorylation. To further investigate which kinase may be involved, we examined the inhibitory effect of some kinase inhibitors. Src kinase up-regulates the ERK cascade through direct phosphorylation of RAF and initiates RAF/MEK/ERK signaling cascade.¹⁷ We found that Src kinase inhibitor dasatinib blocked the FBS-induced ERK/ β -arrestin interaction (Fig. 1D), suggesting an involvement of Src kinase in ERK activation. However, c-RAF inhibitor GW5074 or MEK kinase inhibitor CI-1040 did not block the ERK/ β -arrestin interaction signals (Fig. 1E), although RAF and MEK are considered to be immediately upstream of ERK in the RAF/MEK/ERK signaling pathway. Most kinase inhibitors were discovered by cell-free enzymatic assays. However, paradoxical behaviors of these kinase inhibitors in cells were observed.^{18–20} To test whether the crosstalk or feedback mechanism is involved in the RAF/MEK/ERK signaling pathway, we treated the cells with a combination of both kinase inhibitors. We found that the combination of GW5074 and CI1040 inhibited the ERK/ β -arrestin interaction signals in a dose-dependent manner (Fig. 1E). The data indicated that inhibition of either MEK1 or RAF is necessary but not sufficient alone to disrupt the ERK and β -arrestin interaction. The results also suggest that cell-based functional assays may be more physiological relevant for identification and validation of kinase inhibitors.

ERK/ β -arrestin interaction in response to activation of G-protein coupled receptors

Upon activation of GPCRs, GPCRs signal through one or more G protein-mediated signaling pathways. G-protein receptor kinases phosphorylate ligand-occupied GPCRs, the phosphorylation facilitates recruitment of β -arrestins to receptors, and promotes receptor internalization. β -arrestin has also been implicated to link GPCR-mediated MAPK signaling pathways.¹ To examine if activation of GPCRs promotes ERK/ β -arrestin interaction, we

utilized some of endogenously expressed GPCRs in U2OS cells. We found that GPCR agonists, including isoproterenol, dopamine, S1P and adrenaline, dose-dependently stimulated the ERK/ β -arrestin interaction signals with the EC₅₀ values of 71.9, 140, 337, and 918 nM, respectively (Fig. 2A). We also performed HTRF cAMP assays with the ERK/ β -arrestin LinkLight cells. The results demonstrated that exogenous expression of ERK2-TEV and β -arrestin-pLuc did not affect endogenous GPCR signaling (Fig. S2). In addition, we determined the response of exogenously expressed ADRB2 to its agonist in the ERK/ β -arrestin LinkLight cells. We observed a robust agonist-response with an EC₅₀ of 23 nM (Fig. S3), similar to that previously reported.⁵

We then determined whether GPCR antagonists could block the ERK/ β -arrestin interaction induced by GPCR agonists. ICI118551 (β -adrenergic receptor-2 antagonist), and clozapine (dopamine receptor antagonist) dose-dependently inhibited the agonist-induced ERK/ β -arrestin interaction signals (Fig. 2B). The IC₅₀ values for ICI118551 and clozapine were 37.6 and 211 nM, respectively, although they were 3 to 8 folds less potent in comparison to those determined in the classic G-protein signaling based assays.^{19, 21, 22}

ERK/ β -arrestin interaction in response to activation of receptor tyrosine kinases

β -arrestin has recently been implicated to be part of the downstream MAPK signaling cascades elicited by receptor tyrosine kinases, such as insulin-like growth factor and epidermal growth factor (EGF) receptors.² Binding of ligands to RTKs results in RTK auto-phosphorylation and recruits adaptor proteins to propagate signal transductions including the MAPK signaling cascade.^{3, 4} To probe the ERK/ β -arrestin interaction in response to activations of RTKs, we determined the ERK/ β -arrestin interaction signals stimulated by endogenous RTKs expressed in U2OS cells in response to two growth factors EGF and FGF-b.^{18, 23} We found that both EGF and FGF-b dose-dependently stimulated luciferase signals. The EC₅₀ values for EGF and FGF-b were 3.75 and 3.87 ng/ml, respectively (Fig. 3A and B). These values were similar to those of 0.5 to 5.0 ng/ml previously reported.²⁴

We then investigated whether specific RTK inhibitors could inhibit the agonist-induced ERK2 and β -arrestin-2 interaction signals. We found that afatinib (EGFR inhibitor) and PD166285 (FGFR inhibitor) blocked EGF and FGF-basic induced ERK/ β -arrestin interaction signals (Fig. 3C and D). The IC₅₀ values of afatinib and PD166285 were 585 and 338 nM, respectively. In addition, the pan-kinase inhibitor staurosporine inhibited the interaction signal (Fig. S4). Taken together, the results indicate that binding of agonists to RTKs promotes formation of the ERK/ β -arrestin signal complex. This RTK agonist-induced ERK/ β -arrestin interaction signals can be blocked by corresponding RTK inhibitors.

ERK/ β -arrestin interaction in response to activation of cytokine receptors

Tumor necrosis factor alpha (TNF α) is a master cytokine that mediates inflammatory responses and innate immunity.²⁰ TNF α has been implicated in the disease pathogenesis including cancer, sepsis, rheumatoid arthritis, diabetes and inflammatory bowel disease. Major pathways activated by TNF α include caspases, NF κ B, and MAP kinases. The early phase of TNF α signaling induces expression of inflammatory cytokines including TNF α that results in a secondary cytokine expression and secretion. Secondary cytokines contribute to

most of biological functions of TNF α . Functional Involvement of MAP kinases have been indicated in both upstream (TNF α -producing cells) and downstream (TNF α -acting cells) of TNF α signalling.²⁵

To examine the role of the ERK/ β -arrestin signal complex in cytokine receptor signaling, we applied TNF α to the U2OS cells that express endogenous TNF α receptor 1 (TNFR1). TNF α induced robust luciferase signals (Fig. 4A) with an EC₅₀ of 5.65 ng/ml, indicating the participation of the ERK/ β -arrestin signal complex in TNFR1 signaling. We also found that another cytokine IL1 β stimulated ERK and β -arrestin interaction (Fig. 4B) with an EC₅₀ of 2.8 ng/ml that is compatible to the reported value.²⁶ U2OS cells are known to respond to IL1 stimulation.

Since the small molecule TNFR1 and IL1R inhibitors are currently unavailable, we used staurosporine (pan-kinase inhibitor) to block TNF α and IL1 β induced ERK/ β -arrestin interaction signals. We found that staurosporine dose-dependently inhibited the ERK/ β -arrestin interaction signals in response to TNF α or IL1 β stimulation (Fig. 4C and D). The results suggest that the ERK/ β -arrestin interaction has been involved in cytokine receptor signalling pathways, which can be blocked by kinase inhibitors.

Discussion

We established a cell based assay that can quantitatively assess the interaction or formation of ERK/ β -arrestin signal complex. This assay has been employed to investigate the relationship of the signal complex in response to activations of diverse cell surface receptors. We demonstrate that binding of a GPCR, cytokine or RTK agonist to its corresponding cell surface receptor induces the interaction of ERK and β -arrestin or formation of the ERK/ β -arrestin signal complex. The interaction is agonist concentration-dependent and can be blocked by a specific receptor antagonist or kinase inhibitors. Thus, the ERK/ β -arrestin signal complex functions as a convergence point for diverse cell-surface receptor signal transduction.

In the past decade, β -arrestins have been found to play much broader roles in cell biology than previously thought. The β -arrestins appear to be important mediators of a broad range of cellular activities such as cell proliferation, differentiation, migration, and apoptosis. Growing evidence indicates that β -arrestins interact with non-receptor kinases to form signaling complexes for signal propagation and transduction. Activation of diverse cell-surface receptors can lead to intracellular kinase activation. One of the major cellular signaling mediators activated by diverse types of cell-surface receptors is ERK. The ERK protein mediates multiple signaling pathways, which transduces extracellular signals to intracellular signaling pathways by phosphorylating a large number of protein substrates located in both the cytosol and nucleus.

Many cellular signal transducers are shared by diverse signaling pathways.²⁷ For example, the signal adaptors are shared by diverse biological inputs. The shared pathways enable cells to integrate signals of diverse external stimuli to the similar cellular responses. These signal adaptors are regulated by highly coordinated phosphorylation-mediated protein-protein

interactions. Protein phosphorylation induces protein conformational changes, enables phosphorylated proteins to interact with other proteins, influences protein stability, translocates protein to different subcellular locations, and defines downstream signaling cascades. Accordingly, MAPKs including ERK participate in cellular signal transductions elicited by diverse cell-surface receptors.²⁸ The ERK/ β -arrestin signal complex could be one of the central cross-communication integrating nodes.

The mechanisms of how β -arrestins participating in the formation of individual MAPK module/ β -arrestin complexes may vary among different types of receptors. Thus, GPCRs, RTKs and cytokine receptors are intimately involved in multilayer forms of cross-talk that influences a variety of cellular signaling processes. We hypothesized that the formation of ERK/ β -arrestin signal complex could be a common signal integration node in signaling pathways of diverse cell-surface receptors. By utilizing the ERK/ β -arrestin interaction as a functional signal readout, we demonstrated that the ERK/ β -arrestin signal complex is formed after activation of endogenous GPCRs, RTKs, and cytokine receptors in U2OS cells. We also demonstrated that GPCR antagonists, specific RTK inhibitors and kinase inhibitors blocked the ERK/ β -arrestin interaction signals induced by the agonists. Thus, our data indicate that different types of receptors, including GPCRs, RTKs, and cytokine receptors, transduce signals at least partially through the ERK/ β -arrestin signal complex.

By showing that the ERK/ β -arrestin signal complex is a common transducer in multiple receptor signaling pathways, a question of signaling specificity has arisen. β -arrestin binds to two or more cellular proteins in a signaling cascade and brings them in close proximity. This not only facilitates signal propagation and amplification, but also insulates a specific signaling pathway from other signals. Multiple non-receptor kinases, including PI3K, AKT, JNK, p38, and SRC, have also been implicated to interact with β -arrestins.²⁹ In addition, extracellular stimuli may have unequal impacts on formation of various kinase/ β -arrestin signal complexes and their roles in signal transduction need to be investigated further. Signaling specificity can also be achieved by involvement of discrete protein-protein interactions which may have different responsive time courses. Specific subcellular locations of interacting proteins may also define signaling specificity.

Defining mechanisms for receptor signaling pathways is essential for understanding many aspects of cell biology, cancer biology, and immunology, as well as for effectively targeting signaling pathways for drug discovery. Kinase inhibitors are traditionally discovered by biochemical enzymatic assays, binding assays, and antibody-based protein phosphorylation assays. However, functions of kinases in cells are carried out not only by their intrinsic enzymatic activity but also by cellular protein-protein interaction networks. In cells, kinase and phosphatase activities are dynamically regulated in response to environmental stimuli which may be very different from those in cell-free enzymatic activity assays. The cellular feedback and crosstalk mechanisms could compromise effects of kinase inhibitors discovered from cell-free biochemical assays.

Our data indicated that the RAF and MEK kinase inhibitor used individually did not block the ERK/ β -arrestin interaction. However, the interaction was inhibited by a combination of RAF and MEK inhibitors. This could be explained by the crosstalk and feedback

mechanisms in the RAF/MEK/ERK signaling cascade. The result also suggests that the cell-based kinase signaling assay could be a better alternative for identification of kinase inhibitors. However, the precise mechanism of this phenomena is still to be identified.

Receptors including GPCRs, cytokine receptors and RTKs account for over a half of confirmed drug targets. Assessing receptor activity by measurements of its functions, such as signaling pathways, is commonly used in the drug discovery process. The traditional RTK kinase assays used antibody based or reporter-gene based assays to monitor the signaling pathway indirectly. Instead, our assay directly measures the ERK/ β -arrestin interaction which provides more precise spatio-temporal information. We have demonstrated that the ERK/ β -arrestin interaction can be induced by activation of diverse GPCRs, cytokine receptors and RTKs. Selective small molecule antagonists modulating the ERK/ β -arrestin interaction could be discovered using this cell based assay through blockade of the agonist response of either an endogenous receptor or an exogenously expressing receptor. Targeting GPCR's ERK/ β -arrestin signaling pathway for drug discovery may identify new generation of drug candidates that may have better selectivity with fewer adverse effects. RTKs function as key regulators of cell growth, proliferation, and survival. The uncontrolled RTK signaling that leads to cancer has provided the rationale for anti-RTK drug development. Currently, small molecule RTK antagonists are rarely available. Screening of small molecule compound library via RTK-stimulated ERK/ β -arrestin interaction assay may result in discovery of new functional RTK inhibitors. To identify specificity of hits from screening using this assay, a secondary assay is needed to further demonstrate specificity of hits within the interested signaling pathway.

Additionally, the assay time for assessment of the ERK/ β -arrestin interaction is relative short and the assay signals are independent of reporter gene transcription and translation. The feature is particularly useful for cytokine receptor signaling pathways in which the signal transduction is primarily dependent on the transient protein-protein interactions. Unlike RTKs, cytokine receptors such as TNF α receptors have no intrinsic enzymatic activity. Currently, TNF α receptor activity is measured by secreted secondary cytokines or reporter-gene assay driven by the TNF α -response element.³⁰ However, the biphasic nature of TNF α signaling complicates results obtained from those assays. Our assay directly measures the upstream ERK/ β -arrestin interaction that has obvious advantages over those antibody-based and reporter-gene based assays for the measurements of cytokine receptor activity. The cells treated with cytokines produced much robust signals than the cells treated with GPCR and RTK agonists. Potentially, other signaling mechanisms may be participated. Additionally, the S/B ratios and Z' factors of ERK/ β -arrestin LinkLight assays were determined for four agonists. The S/B ratio is 3.66 fold and Z' factor is 0.69 for 10 μ M dopamine; the S/B ratio is 3.78 fold and Z' factor is 0.82 for 10 μ M adrenaline; the S/B ratio is 8.99 fold and Z' factor is 0.78 for 0.1 μ g/ml TNF α ; the S/B ratio is 14.5 fold and Z' factor is 0.81 for 0.1 μ g/ml IL1 β (Fig. S5). Thus this assay is robust and amenable to high throughput screening for drug candidate discovery.

Supplementary Material

Refer to Web version on PubMed Central for supplementary material.

Acknowledgments

The authors thank Dr. Paul S. Wright from sanofi for suggestions and carefully reading the manuscript and Mr. Paul Shinn from the National Center for Advancing Translational Sciences for preparing compound plates. The authors also thank Nicolas Pierre at cisbio for providing HTRF cAMP kit.

Funding

This work was supported by the NIH/NIMH SBIR grant (R43MH101903) awarded to H.E. and was supported by the Intramural Research Programs of the National Center for Advancing Translational Sciences, NIH.

References

1. Lefkowitz RJ, Whalen EJ. beta-arrestins: traffic cops of cell signaling. *Curr Opin Cell Biol.* 2004; 16:162–8. [PubMed: 15196559]
2. Hupfeld CJ, Olefsky JM. Regulation of receptor tyrosine kinase signaling by GRKs and beta-arrestins. *Annu Rev Physiol.* 2007; 69:561–77. [PubMed: 17002595]
3. Akhtar S, Yousif MH, Chandrasekhar B, et al. Activation of EGFR/ERBB2 via pathways involving ERK1/2, P38 MAPK, AKT and FOXO enhances recovery of diabetic hearts from ischemia-reperfusion injury. *PLoS One.* 2012; 7:e39066. [PubMed: 22720029]
4. Lanner F, Rossant J. The role of FGF/Erk signaling in pluripotent cells. *Development.* 2010; 137:3351–60. [PubMed: 20876656]
5. Hwang YS, Jeong M, Park JS, et al. Interleukin-1 beta stimulates IL-8 expression through MAP kinase and ROS signaling in human gastric carcinoma cells. *Oncogene.* 2004; 23:6603–6611. [PubMed: 15208668]
6. Schrofelbauer B, Hoffmann A. How do pleiotropic kinase hubs mediate specific signaling by TNFR superfamily members? *Immunol Rev.* 2011; 244:29–43. [PubMed: 22017429]
7. Luttrell LM. ‘Location, location’: activation and targeting of MAP kinases by G protein-coupled receptors. *J Mol Endocrinol.* 2003; 30:117–126. [PubMed: 12683936]
8. Yao Z, Seger R. The ERK signaling cascade-Views from different subcellular compartments. *Biofactors.* 2009; 35:407–416. [PubMed: 19565474]
9. Tohgo A, Pierce KL, Choy EW, et al. beta-Arrestin scaffolding of the ERK cascade enhances cytosolic ERK activity but inhibits ERK-mediated transcription following angiotensin AT1a receptor stimulation. *J Biol Chem.* 2002; 277:9429–36. [PubMed: 11777902]
10. Eishingdrelo H, Kongsamut S. Minireview: Targeting GPCR Activated ERK Pathways for Drug Discovery. *Curr Chem Genomics Transl Med.* 2013; 7:9–15. [PubMed: 24396730]
11. Hammer MM, Wehrman TS, Blau HM. A novel enzyme complementation-based assay for monitoring G-protein-coupled receptor internalization. *Faseb J.* 2007; 21:3827–34. [PubMed: 17942829]
12. Eishingdrelo H, Cai J, Weissensee P, et al. A cell-based protein-protein interaction method using a permuted luciferase reporter. *Curr Chem Genomics.* 2011; 5:122–8. [PubMed: 22207892]
13. Violin JD, DeWire SM, Yamashita D, et al. Selectively Engaging beta-Arrestins at the Angiotensin II Type 1 Receptor Reduces Blood Pressure and Increases Cardiac Performance. *J Pharmacol Exp Ther.* 2010; 335:572–579. [PubMed: 20801892]
14. Boerrigter G, Whalen EJ, Lark M, et al. Cardiorenal Actions of TRV120027, a Novel, beta-Arrestin-Biased Ligand at the Angiotensin II Type I Receptor, in Healthy and Heart Failure Canines: A Novel Therapeutic Strategy for Acute Heart Failure. *J Card Fail.* 2010; 16:S73–S73.
15. Smith NJ, Luttrell LM. Signal switching, crosstalk, and arrestin scaffolds: novel G protein-coupled receptor signaling in cardiovascular disease. *Hypertension.* 2006; 48:173–9. [PubMed: 16818801]
16. Kolch W. Meaningful relationships: the regulation of the Ras/Raf/MEK/ERK pathway by protein interactions. *Biochem J.* 2000; 351:289–305. [PubMed: 11023813]
17. Lin A, DeFea KA. beta-Arrestin-kinase scaffolds: turn them on or turn them off? *Wires Syst Biol Med.* 2013; 5:231–241.

18. Peters MF, Vaillancourt F, Heroux M, et al. Comparing label-free biosensors for pharmacological screening with cell-based functional assays. *Assay Drug Dev Technol.* 2010; 8:219–27. [PubMed: 20085460]
19. Bilski AJ, Halliday SE, Fitzgerald JD, et al. The pharmacology of a beta 2-selective adrenoceptor antagonist (ICI 118,551). *J Cardiovasc Pharmacol.* 1983; 5:430–7. [PubMed: 6191142]
20. Schall TJ, Lewis M, Koller KJ, et al. Molecular cloning and expression of a receptor for human tumor necrosis factor. *Cell.* 1990; 61:361–70. [PubMed: 2158863]
21. Canton H, Verrielle L, Colpaert FC. Binding of typical and atypical antipsychotics to 5-HT1C and 5-HT2 sites: clozapine potently interacts with 5-HT1C sites. *Eur J Pharmacol.* 1990; 191:93–6. [PubMed: 1982659]
22. Matsumoto RR, Bowen WD, Tom MA, et al. Characterization of two novel sigma receptor ligands: antidystonic effects in rats suggest sigma receptor antagonism. *Eur J Pharmacol.* 1995; 280:301–10. [PubMed: 8566098]
23. McGuinness RP, Proctor JM, Gallant DL, et al. Enhanced Selectivity Screening of GPCR Ligands Using a Label-Free Cell Based Assay Technology. *Comb Chem High T Scr.* 2009; 12:812–823.
24. Miyatake M, Rubinstein TJ, McLennan GP, et al. Inhibition of EGF-induced ERK/MAP kinase-mediated astrocyte proliferation by mu opioids: integration of G protein and beta-arrestin 2-dependent pathways. *J Neurochem.* 2009; 110:662–74. [PubMed: 19457093]
25. Aggarwal BB. Signalling pathways of the TNF superfamily: A double-edged sword. *Nat Rev Immunol.* 2003; 3:745–756. [PubMed: 12949498]
26. Han J, Xie Y, Lan F, et al. Additive effects of EGF and IL-1beta regulate tumor cell migration and invasion in gastric adenocarcinoma via activation of ERK1/2. *Int J Oncol.* 2014
27. Song XF, Coffa S, Fu HA, et al. How Does Arrestin Assemble MAPKs into a Signaling Complex? *Journal of Biological Chemistry.* 2009; 284:685–695. [PubMed: 19001375]
28. Delcourt N, Bockaert J, Marin P. GPCR-jacking: from a new route in RTK signalling to a new concept in GPCR activation. *Trends Pharmacol Sci.* 2007; 28:602–607. [PubMed: 18001849]
29. Luttrell LM, Roudabush FL, Choy EW, et al. Activation and targeting of extracellular signal-regulated kinases by beta-arrestin scaffolds. *Proc Natl Acad Sci U S A.* 2001; 98:2449–54. [PubMed: 11226259]
30. Leister KP, Huang R, Goodwin BL, et al. Two High Throughput Screen Assays for Measurement of TNF-alpha in THP-1 Cells. *Curr Chem Genomics.* 2011; 5:21–9. [PubMed: 21643507]




Figure 1.

Design of the ERK/ β -arrestin interaction LinkLight assay. (A) A schematic of a canonical RAF/MEK/ERK signaling LinkLight assay. Cell surface receptors including GPCRs, RTKs, and cytokine receptors in response to extracellular stimuli lead to activation of the RAF/MEK/ERK signaling pathway. Activated ERK releases from the signal complex tiered by scaffold and anchor proteins and binds to β -arrestin to form ERK/ β -arrestin signal complex. The close proximity between TEV-tagged ERK and pLuc-tagged β -arrestin leads to a proteolytic cleavage of permutated luciferase (pLuc). The cleaved luciferase fragments spontaneously refold to reconstitute an active luciferase, driving by the high affinity of fragment self-complementation. The luminescent signals are then detected after addition of luciferase detection reagent. (B) Dose-response curves generated using the ERK/ β -arrestin interaction assay in response to different concentrations of FBS. After overnight culture,

regular culture medium was replaced with serum-free medium, and cells (20,000 cells/well in a 384-well plate) were stimulated with different concentrations of FBS for 30, 60, 90, 120, and 150 min at 37°C. (C) Dose response of an pan-kinase inhibitor staurosporine. Cells were treated with varying doses of staurosporine in serum-free DMEM for 15 min and then treated with 5% FBS for 90 min. (D) Dose response of the anti-cancer drug Dasatinib (SRC inhibitor). Cells were treated with varying doses of dasatinib in serum-free DMEM for 15 min and then treated with 5% FBS for 90 min. (E) Dose response of individual and combination of GW5074 (c-RAF inhibitor) and CI1040 (MEK inhibitor). Cells were treated with varying doses of kinase inhibitor alone or in combination in serum-free DMEM for 15 min and then stimulated with 5% FBS for 90 min. All errors bars represent SD.

Author Manuscript


Author Manuscript

Author Manuscript

Author Manuscript

2a

2b

**Figure 2.**

ERK/ β -arrestin interaction signals in response to activation of endogenous GPCRs in U2OS cells. (A) Dose-response curves generated using the ERK/ β -arrestin interaction assay in response to varying doses of GPCR agonists adrenaline, isoproterenol, dopamine, and S1P. Cells were stimulated with varying doses of GPCR agonists in serum-free DMEM for 90 min at 37°C. (B) Antagonists inhibit GPCR agonist-induced ERK/ β -arrestin interaction signals. Cells were treated with varying doses of GPCR antagonists BD-1047, ICI118551, and Clozapine in serum-free DMEM for 15 min and then treated with 1 μ M adrenaline, isoproterenol, and dopamine respectively for 90 min.




Figure 3.

ERK/β-arrestin interaction signals in response to activation of endogenous RTKs in U2OS cells. (A and B) Dose-response curves generated using the ERK/β-arrestin interaction LinkLight assay in response to varying doses of RTK ligands EGF and FGF-basic. Cells (20,000/well) were stimulated with varying doses of RTK ligands in serum-free DMEM for 90 min at 37°C. (C and D) Inhibition of growth factor-induced ERK/β-arrestin interaction signals. Cells were treated with varying doses of EGFR inhibitor Afatinib, and FGFR inhibitor PD166285 in serum-free DMEM for 15 min and then stimulated with 25 ng/ml EGF and FGF-basic respectively for 90 min.




Figure 4.

ERK/β-arrestin interaction signals in response to activation of endogenous cytokine receptors in U2OS cells. (A and B) Dose-response curves generated using the ERK/β-arrestin interaction assay in response to varying doses of recombinant human TNFα and IL1β. Cells were stimulated with varying doses of TNFα and IL1β in serum-free DMEM for 90 min at 37°C. (C and D) Inhibition of TNFα and IL1β-induced signals by staurosporine. Cells were treated with varying doses of staurosporine in serum-free DMEM for 15 min and then stimulated with TNFα (25 ng/ml) and IL1β (50 ng/ml) respectively for 90 min.