

VIROLOGY



Multicenter Evaluation of the Cepheid Xpert Hepatitis C Virus Viral Load Assay

M. P. McHugh,^a A. H. B. Wu,^b S. Chevaliez,^c J. M. Pawlotsky,^c M. Hallin,^d K. E. Templeton^a

Specialist Virology Centre, Royal Infirmary of Edinburgh, Edinburgh, United Kingdom^a; Department of Laboratory Medicine, University of California, San Francisco, San Francisco, California, USA^b; National Reference Center for Viral Hepatitis B, C and delta, Department of Virology, Henri Mondor Hospital, Créteil, France^c; Center for Molecular Diagnostics, Laboratoire Hospitalier Universitaire de Bruxelles, Brussels, Belgium^d

ABSTRACT Viral load monitoring for hepatitis C virus (HCV) is necessary to diagnose infection and monitor response to therapy, but the tests involved are currently confined to specialist institutions. There is a need for a fast, accurate assay with limited operator input to enhance the access to viral load monitoring. We evaluated the quantification of HCV RNA in serum and plasma by the Cepheid Xpert HCV Viral Load assay in comparison to the Abbott RealTime HCV assay. Serum and plasma samples were gathered from HCV-infected individuals at four international sites. These were tested with the Xpert HCV Viral Load assay, and results were compared to quantification by the Abbott RealTime HCV assay. An external quality assessment panel of eight samples was also tested. In total, 614 samples were analyzed in the study, and the qualitative results agreed on the two platforms for 588 (95.8%) samples. Further analysis of 396 samples quantified by both tests showed strong correlation (correlation coefficient r = 0.99) across the quantifiable range, with Bland-Altman plot data showing a mean difference (\pm 1.96 standard deviation) of 0.03 \pm 0.44 log₁₀ IU/ml. In the external quality assessment panel, the Xpert HCV Viral Load assay results (quantified in log₁₀ IU per milliliter) were within 1 standard deviation of the target value for all but one sample, which was also similarly misquantified by the Abbott Real-Time HCV assay. The Xpert HCV Viral Load assay performs well compared to a market-leading HCV viral load test and should be considered for instances where rapid near-to-patient testing is required.

KEYWORDS HCV RNA, hepatitis C virus, viral load

epatitis C virus (HCV) is a leading burden on human health worldwide. Recent estimates suggest that there are over 185 million individuals living with chronic HCV infection and that HCV causes over 350,000 deaths every year (1, 2). HCV is a blood-borne virus that infects the liver; chronic infection causes liver disease and extrahepatic manifestations, cirrhosis, hepatocellular carcinoma, and death.

Detection of HCV RNA is necessary to diagnose infection. In addition, determination of HCV viral load is recommended prior to beginning antiviral therapy to establish a baseline, and subsequent monitoring is recommended to ascertain the response to treatment, particularly to assess the sustained virological response (3, 4). Currently, HCV viral load monitoring must be done in specialist centers, with a turnaround time of up to 1 week due to the need to batch samples (5, 6). The Xpert HCV Viral Load assay (Cepheid, Sunnyvale, CA) is a recently CE-marked *in vitro* diagnostic (IVD) test for HCV viral load determination in plasma and serum. Assay components are contained within a single-use cartridge which performs RNA extraction, reverse transcription, and real-time PCR targeting the 5' untranslated region (UTR) of the HCV genome. The test cartridge contains internal controls to ensure accurate test performance and to quantify HCV viral load; results are interpreted by proprietary software. We describe a multi-

Received 16 December 2016 Returned for modification 3 January 2017 Accepted 3 March 2017

Accepted manuscript posted online 8 March 2017

Citation McHugh MP, Wu AHB, Chevaliez S, Pawlotsky JM, Hallin M, Templeton KE. 2017. Multicenter evaluation of the Cepheid Xpert Hepatitis C Virus Viral Load assay. J Clin Microbiol 55:1550–1556. https://doi.org/ 10.1128/JCM.02460-16.

Editor Yi-Wei Tang, Memorial Sloan-Kettering Cancer Center

Copyright © 2017 American Society for Microbiology. All Rights Reserved. Address correspondence to K. E. Templeton, Kate.Templeton@nhslothian.scot.nhs.uk.

	No. of samples w RealTi <i>m</i> e HCV re			
Assay and result	Virus detected, ≥12 IU/ml	Virus detected, <12 IU/ml	Virus not detected	Total no. of samples
Xpert HCV Viral Load				
Virus detected, \geq 10 IU/ml	396	4	0	400
Virus detected, <10 IU/ml	8	24	13	45
Virus not detected	0	1	168	169
Total	404	29	181	614

TABLE 1 Qualitative results of Xpert HCV Viral Load and Abbott RealTime HCV assays

center study assessing the performance of the Cepheid Xpert HCV Viral Load assay compared to that of the reference Abbott RealTime HCV assay.

(These results were presented as a poster at the European Society for Clinical Virology Meeting on 9 to 12 September 2015 in Edinburgh, United Kingdom [abstract 1530] [7]).

RESULTS

Testing of clinical samples. In total, 636 samples from HCV-infected individuals were tested in this study. Thirty-three (5.2%) samples generated indeterminate results with the Xpert HCV Viral Load assay; no particular association was seen with genotype (6 genotype 1, 1 genotype 2, 1 genotype four, and 25 of unknown genotype). Assays of 11 samples were repeated and produced valid results; the remaining 22 samples had insufficient volume for repeat assays and so were excluded from analysis. Of the 614 samples analyzed, 588 (95.8%) had qualitative results that agreed on the two platforms (disagreement was seen in 8 genotype 1, 4 genotype 3, 3 genotype 4, and 11 unknown-genotype samples) (Table 1). The eight samples quantified by the Abbott RealTime HCV assay but not by the Xpert HCV Viral Load assay had mean \pm standard deviation (SD) viral load values of 16.0 \pm 4.6 IU/ml, and the four quantified by the Xpert HCV Viral Load assay but not by the Abbott RealTime HCV assay had mean \pm SD viral load values of 16.0 \pm 3.5 IU/ml. There were also 14 samples with detectable but unquantifiable HCV on one platform and with HCV not detected on the other; these discrepancies most likely represent Poisson distributions of low-level target nucleic acid.

Further analysis of 396 quantified samples (including 169 with known genotype) by Deming regression (Fig. 1A) showed that the results from the two platforms correlated strongly across the quantifiable range (correlation coefficient [r], 0.99; Deming regression equation, Y = 1.02X + 0.09). The correlation between the two assays was strong in analyzing for genotype 1 and 3 samples separately (for genotype 1, r = 0.99, Deming regression equation Y = 1.03X + 0.05; for genotype 3, r = 0.99, Deming regression equation Y = 1.02X + 0.04) (Fig. 1C and E). The level of HCV quantification was higher with the Xpert HCV Viral Load assay in 232 (59.6%) samples, and the overall means \pm 1.96 SD difference of results between the assays was 0.03 \pm 0.44 log₁₀ IU/ml, with a range of -2.64 to $+ 1.02 \log_{10}$ IU/ml (Fig. 1B). A discrepancy in quantification between the two assays was detected for 13 (2.1%) samples that fell outside the estimated agreement interval, defined as \pm 1.96 SD from the overall mean difference (Fig. 1B and Table 2). There was no remaining sample available for these to be investigated further. Further details of the discordant samples are shown in Table 2. The agreement between the results of the assays for genotypes 1 and 3 alone was also acceptable (Fig. 1D and F). Additionally, 100 HCV-negative blood donor samples all produced undetectable results when tested with the Xpert HCV Viral Load assay, indicating 100% specificity (95% confidence interval, 97.0% to 100.0%).

EQA panel testing. The external quality assessment (EQA) panel contained one HCV-negative sample, two HCV genotype 1b samples, and five HCV genotype 3a samples (Table 3). These were run twice using the Xpert HCV Viral Load assay and once




FIG 1 Deming regression and Bland-Altman plot analysis of HCV quantification by the Abbott RealTime HCV and Xpert HCV Viral Load assays. (A) Deming regression of 396 samples quantified by both assays. (B) Bland-Altman plot of 396 samples quantified by both assays. (C) Deming regression of 92 HCV genotype 1 samples quantified by both assays. (D) Bland-Altman plot of 92 HCV genotype 1 samples quantified by both assays. (E) Deming regression of 51 HCV genotype 3 samples quantified by both assays. (F) Bland-Altman plot of 51 HCV genotype 3 samples quantified by both assays. (F) Bland-Altman plot of 51 HCV genotype 3 samples quantified by both assays. (F) Bland-Altman plot of 51 HCV genotype 3 samples quantified by both assays. (F) dualt here are a samples quantified by both assays. (F) Bland-Altman plot of 51 HCV genotype 3 samples quantified by both assays. The Deming regression plot shows the identity line (dashed line) and the Deming fit with a 95% confidence interval (black line and two dotted lines, respectively). The Bland-Altman plot shows the difference between the two assays in HCV quantification results plotted as a function of the mean of the two; the mean difference is highlighted by the black line, and dotted lines indicate \pm 1.96 standard deviation.

using the Abbott RealTime HCV assay. A single sample tested using one Xpert HCV Viral Load assay run and the Abbott RealTime HCV assay run was quantified outside the target range, due to underquantification of genotype 1b HCV. All other samples were called within the target range by both platforms (Table 3).

TABLE 2 Samples with discrepancy between assays

			Xpert HCV Viral Load result	Abbott RealTime HCV result	Difference
Study no.	Storage	HCV Subtype(s)	(log ₁₀ IU/ml)	(log ₁₀ IU/ml)	(log ₁₀ IU/ml)
HCC229065	Frozen	1b	5.96	4.94	1.02
HCC207168	Frozen	1a	5.36	4.61	0.75
HCC111006	Frozen	Unknown	5.81	5.06	0.75
HCC229207	Fresh	1b	4.48	3.80	0.68
HCV167007	Frozen	5a	6.74	6.18	0.56
W7897	Frozen	3a	6.07	5.53	0.54
HCC229071	Frozen	1b	5.36	4.88	0.48
W7930	Frozen	1b	1.57	2.02	-0.45
HCC207027	Frozen	Unknown	1.36	1.82	-0.46
W7820	Frozen	2a/2c	4.49	5.02	-0.53
HCC248052	Frozen	Unknown	6.31	6.89	-0.58
HCC111176	Frozen	Unknown	1.58	2.42	-0.84
HCC248043	Frozen	Unknown	4.13	6.77	-2.64

DISCUSSION

This study showed that the Xpert HCV Viral Load assay accurately quantifies HCV viral load compared to a leading commercial assay used worldwide. Also, the Xpert HCV Viral Load assay performed strongly in a well-respected EQA panel. These performance characteristics, aligned with limited hands-on time (5 min), short run time (105 min), random access testing, and the uncomplicated operator input of the Xpert HCV Viral Load assay, suggest that it can have an important role in diagnosis of HCV infection and HCV viral load monitoring. Our findings confirm those of Gupta et al. (8), who tested 118 plasma samples seropositive for HCV with the Xpert HCV Viral Load and Abbott RealTime HCV assays, Bland-Altman analysis showed a mean difference of 0.04 \log_{10} IU/ml (1.96 SD, -0.42 to 0.49 \log_{10} IU/ml), a result very similar to those reported in our multicenter study. We believe that the Xpert HCV Viral Load assay could be used for routine viral load monitoring as well as in situations such as organ donor testing or needle stick injuries, where a rapid evaluation of HCV status is required, and could even be used in a near-patient capacity in HCV clinics or prisons (9).

Although we found the Xpert HCV Viral Load assay result to be slightly higher overall than the Abbott RealTime HCV result ($\leq 2 \log_{10} IU/mI$), the Xpert HCV Viral Load assay result appears to have been lower (Fig. 1B), with 3 of 29 samples falling outside a 1.96 SD from the mean difference. Also, the Abbott RealTime HCV assay quantified HCV in 8 samples that were detected by the Xpert HCV Viral Load assay at levels below the limit of quantification, compared to 4 samples quantified by the Xpert HCV Viral Load assay (Table 1). However, no samples were quantified by the Abbott RealTime HCV assay but not detected by the Xpert HCV Viral Load assay. We believe that this finding will have limited impact in clinical use, as individuals generally have viral loads of $>2 \log_{10} IU/mI$

TABLE 3 Performance of Xpert HCV Viral Load and Abbott RealTime HCV assays with EQA samples

				Abbott RealTime HCV		Xpert HCV Viral Load result			
		Target result		result		Run 1		Run 2	
Sample no.	Sample content or result	Value ^a	SD ^a	Value ^a	Within target range	Value ^a	Within target range	Value ^a	Within target range
HCVRNA14B-05	HCV type 1b	4.25	0.34	4.16	Yes	4.24	Yes	4.46	Yes
HCVRNA14B-02	HCV type 1b	3.24	0.24	2.97	No	2.93	No	3.20	Yes
HCVRNA14B-01	HCV type 3a	3.96	0.20	3.84	Yes	3.96	Yes	3.96	Yes
HCVRNA14B-08	HCV type 3a	2.93	0.25	2.82	Yes	2.81	Yes	2.99	Yes
HCVRNA14B-06	HCV type 3a	2.10	0.28	1.92	Yes	2.09	Yes	2.00	Yes
HCVRNA14B-07	HCV type 3a	3.82	0.23	3.79	Yes	3.89	Yes	3.84	Yes
HCVRNA14B-03	HCV type 3a	2.86	0.26	2.81	Yes	2.92	Yes	2.89	Yes
HCVRNA14B-04	Negative	ND ^b		ND		ND		ND	

^aLog₁₀ IU/ml. ^bND, not detected.

TABLE 4 Known HCV genotypes tested

Genotype and subgenotype	п	%
1	110	51.4
1a	47	22.0
1b	33	15.4
1d	1	0.5
Subtype unknown	29	13.6
2	15	7.0
2a/c	6	2.8
2b	3	1.4
21	1	0.5
Subtype unknown	5	2.3
3	61	28.5
3a	47	22.0
Subtype unknown	14	6.5
4	21	9.8
4f	1	0.5
4g	1	0.5
4h	1	0.5
4r	1	0.5
Subtype unknown	17	7.9
5	2	0.9
5a	2	0.9
6	4	1.9
Subtype unknown	4	1.9
1 or 6	1	0.5
Total	214	100.0

at diagnosis (8, 10), and the Xpert HCV Viral Load assay did not report any samples falsely negative in our sample collection. Overall, the Xpert HCV Viral Load assay detected HCV in 445 samples compared to 433 samples with HCV detected by the Abbott RealTime HCV assay (Table 1).

Results from one study site showed HCV viral loads significantly lower than those from the other sites by analysis of variance (ANOVA). That site was in France, where the level of access to antiviral treatment is high; the lower viral loads seen at that site were likely due to higher treatment rates as opposed to a methodological difference (11, 12). One drawback of this work was that most samples were tested retrospectively after being frozen, although samples were tested on each platform from the same freeze-thaw cycle. Another drawback was that the genotypes were not available for all samples, although we included 214 samples of 7 different genotypes covering those seen throughout the world (Table 4) (13). Eight of these samples, including one of the two genotype 5a samples, fell outside the agreement interval by Bland-Altman analysis (Fig. 1B and Table 2). Prior to implementation, laboratories should validate that the Xpert HCV Viral Load assay reliably quantifies the HCV types circulating locally. Further studies investigating the performance of the Xpert HCV Viral Load assay in samples with known genotypes, particularly genotype 5a, are warranted.

Disadvantages of the Xpert HCV Viral Load assay include the high cost per test cartridge (\in 35), the lack of external positive and negative controls in the reagent kit, and the relatively high sample input volume (1 ml). Also, one test cartridge in our study leaked after sample addition. Although a rare occurrence, this is a potentially serious issue. If the Xpert HCV Viral Load assay is to be used outside the laboratory setting, it will still be necessary to ensure that operators understand and adhere to appropriate risk assessments for handling infectious samples.

Overall, the Xpert HCV Viral Load assay had performance comparable to that of a

market-leading assay in testing samples from four countries and an international EQA panel. The Xpert HCV Viral Load assay has the potential to change how viral load monitoring is performed for HCV in the future.

MATERIALS AND METHODS

Clinical samples. The study was a method comparison of the Xpert HCV Viral Load assay (Cepheid, Sunnyvale, CA) to the Abbott RealTime HCV assay (Abbott Molecular, Des Plaines, IL) in fresh and frozen human plasma (EDTA) and serum specimens. Remnant samples were collected at four health care sites in Europe and the United States after standard-of-care HCV viral load testing. Frozen (\leq 70°C) samples were tested within 6 months of freezing and were tested on both platforms from the same freeze-thaw cycle; after thawing, samples were held at 4°C for a maximum of 6 h prior to testing. Fresh samples were held at 4°C and tested within 72 h of collection. Samples were from known HCV-positive adults (18 to 89 years old) and were collected per the assay manufacturer's instructions. Different technicians performed sample randomization and testing, and the individuals running the Xpert HCV Viral Load assay were blind to the Abbott RealTime HCV assay result. The study methods were in accordance with local ethical guidelines, and written informed consent was not deemed necessary at each study site. Samples were processed with the Abbott RealTime HCV assay on an Abbott m2000 instrument according to the manufacturer's instructions; all sites are accredited to perform this test. For the Xpert HCV Viral Load assays, 1 ml of specimen was added to a test cartridge and the cartridge was loaded into a GeneXpert instrument. Three sites used the GeneXpert Dx XVI system, and one site used a GeneXpert Infinity 80 system. Acrometrix (Benicia, CA) external controls (HCV high-level, mid-level, and negative controls) were tested each day. The linear range of the Xpert HCV Viral Load assay is 10 IU/ml to 10⁸ IU/ml; the linear range of the Abbott RealTime HCV assay is 12 IU/ml to 10⁸ IU/ml. Results were reported as follows: HCV detected (with the associated quantitation reported in IU per milliliter); HCV not detected; HCV detected at <12 IU/ml (Abbott RealTime HCV) or <10 IU/ml (Xpert HCV Viral Load); and HCV quantified at \geq 12 IU/mI (Abbott RealTime HCV) or ≥10 IU/mI (Xpert HCV Viral Load). Samples with error, invalid, or no-result outputs on the Xpert HCV Viral Load were classed as indeterminate.

A total of 589 samples were collected between August 2013 and October 2014. Of these, four were not collected in line with protocol, one leaked after loading of the Xpert HCV Viral Load assay cartridge, and one frozen sample was thawed for >24 h before Xpert testing; all these were excluded, leaving 583 samples (558 frozen, 25 fresh) for testing. The HCV genotype classification was available for samples from two participating centers; a further 53 frozen plasma samples of various genotypes were obtained from Diagnostic Laboratory Services (DLS; Los Osos, CA) and from one participating center (Table 4). Genotype were determined in independent accredited centers by different methods (Versant HCV Genotype Assay [Innogenetics, Ghent, Belgium], Trugene HCV Genotyping kit [Siemens, Erlangen, Germany], RealTime HCV Genotype [Abbott], Hepatitis C Viral RNA Genotype LiPA [Quest Diagnostics, London, United Kingdom], AmpliSens HCV Genotype-FRT [InterLabService, Moscow, Russia], and in-house reverse transcription-PCR [RT-PCR] and sequencing). The 53 additional samples were collected between 2003 and 2016 and were stored as frozen samples as described above. Also, specimens were collected from 100 blood donors who were confirmed to be HCV seronegative by U.S. Food and Drug Administration-licensed antibody and nucleic acid tests.

External quality assessment panel testing. Quality Control for Molecular Diagnostics (Glasgow, United Kingdom) is an external quality assessment (EQA) provider for routine clinical molecular diagnostic laboratories. The HCV panel contains eight simulated plasma samples containing different genotypes and target RNA concentrations (Table 3). The target range was defined as the mean \pm one standard deviation (SD) of results from all respondents (n = 100) corresponding to the distribution (reference number HCVRNA14B).

Statistical analysis. Descriptive statistics are presented as means \pm SD. We used Deming regression to compare the results of quantification by the two assays, using the Method Comparison Regression (MCR) R package (Roche, Basel, Switzerland). Deming regression takes account of measurement errors for both methods and so is suited to comparing laboratory assays. We also use Bland-Altman plot analysis to highlight the differences between the assays, where the differences between the two assays are plotted against the mean of the results from the two assays. Bland-Altman analysis was performed with the BlandAltmanLeh R package (University of Greifswald, Greifswald, Germany).

ACKNOWLEDGMENT

We thank Cepheid for financially supporting this study.

REFERENCES

- Perz JF, Armstrong GL, Farrington LA, Hutin YJ, Bell BP. 2006. The contributions of hepatitis B virus and hepatitis C virus infections to cirrhosis and primary liver cancer worldwide. J Hepatol 45:529–538. https://doi.org/10.1016/j.jhep.2006.05.013.
- Mohd Hanafiah K, Groeger J, Flaxman AD, Wiersma ST. 2013. Global epidemiology of hepatitis C virus infection: new estimates of agespecific antibody to HCV seroprevalence. Hepatology 57:1333–1342. https://doi.org/10.1002/hep.26141.
- EASL. 2017. European Association for the Study of the Liver recommendations on treatment of hepatitis C. 2016. J Hepatol 66:153–194. https:// doi.org/10.1016/j.jhep.2016.09.001.
- AASLD-IDSA. 27 September 2016. Recommendations for testing, managing, and treating hepatitis C. www.hcvguidelines.org. Accessed 25 November 2016.
- 5. Michelin BDA, Muller Z, Stelzi E, Marth E, Kessler HH. 2007. Evaluation of the Abbott RealTime HCV assay for quantitative detection of hepatitis C

virus RNA. J Clin Virol 38:96–100. https://doi.org/10.1016/j.jcv.2006.11 .007.

- Vermehren J, Yu ML, Monto A, Yao JD, Anderson C, Bertuzis R, Scheider G, Sarrazin C. 2011. Multi-center evaluation of the Abbott RealTime HCV assay for monitoring patients undergoing antiviral therapy for chronic hepatitis C. J Clin Virol 52:133–137. https://doi.org/10.1016/j.jcv.2011.07 .007.
- McHugh M, Wu A, Pawlotsky JM, Chevaliez S, Hallin M, Templeton K. 2015. Evaluation of the Cepheid Xpert[®] HCV Viral Load assay, abstr 1530, poster 2. Abstr Eur Soc Clin Virol Meet, Edinburgh, United Kingdom, 9 to 12 September 2015.
- Gupta E, Agarwala P, Kumar G, Maiwall R, Sarin SK. 2017. Point-of-care testing (POCT) in molecular diagnostics: performance evaluation of GeneXpert HCV RNA test in diagnosing and monitoring of HCV infection. J Clin Virol 88:46–51. https://doi.org/10.1016/j.jcv.2017.01.006.
- Rich JD, Beckwith CG, Macmadu A, Marshall BDL, Brinkley-Rubinstein L, Amon JJ, Milloy MJ, King MRF, Sanchez J, Atwoli L, Altice FL. 2016. Clinical care of incarcerated people with HIV, viral hepatitis, or tuberculosis. Lancet 388:1103–1114. https://doi.org/10.1016/S0140-6736(16)30379-8.
- Alonso R, Pérez-García F, Ampuero D, Reigada E, Bouza E. 2017. New direct-acting antivirals for patients with chronic HCV infection: can we monitor treatment using an HCV core antigen assay? Diagn Microbiol

Infect Dis 87:243–246. https://doi.org/10.1016/j.diagmicrobio.2016.11 .010.

- Geri G, Maynard M, Rosenthal E, Fontaine H, Lacombe K, Slama L, Goujard C, Loustaud-Ratti V, Bergmann JF, Morlat P, Vittecoq D, Alric L, Cacoub P, for the GERMIVIC Group. 2014. Care of hepatitis C virus infection in France: modifications in three consecutive surveys between 1995 and 2010. Liver Int 34:1349–1357. https://doi.org/10.1111/liv .12388.
- 12. Razavi H, Waked I, Sarrazin C, Myers RP, Idilman R, Calinas F, Vogel W, Mendes Correa MC, Hézode C, Lázaro P, Akarca U, Aleman S, Balık I, Berg T, Bihl F, Bilodeau M, Blasco AJ, Brandão Mello CE, Bruggmann P, Buti M, Calleja JL, Cheinquer H, Christensen PB, Clausen M, Coelho HS, Cramp ME, Dore GJ, Doss W, Duberg AS, El-Sayed MH, Ergör G, Esmat G, Falconer K, Félix J, Ferraz ML, Ferreira PR, Frankova S, García-Samaniego J, Gerstoft J, Giria JA, Gonçales FL, Jr, Gower E, Gschwantler M, Guimarães Pessôa M, Hindman SJ, Hofer H, Husa P, Kåberg M, et al. 2014. The present and future disease burden of hepatitis C virus (HCV) infection with today's treatment paradigm. J Viral Hepat 21:34–59. https:// doi.org/10.1111/jvh.12248.
- Messina JP, Humphreys I, Flaxman A, Brown A, Cooke GS, Pybus OG, Barnes E. 2015. Global distribution and prevalence of hepatitis C virus genotypes. Hepatology 61:77–87. https://doi.org/10.1002/hep.27259.