



Published in final edited form as:

Curr Protoc Nucleic Acid Chem. ; 69: 10.17.1–10.17.30. doi:10.1002/cpnc.31.

A High-Throughput Process for the Solid-Phase Purification of Synthetic DNA Sequences

Andrzej Grajkowski¹, Jacek Cie lak¹, and Serge L. Beaucage¹

¹Laboratory of Biological Chemistry, Food and Drug Administration, Silver Spring, Maryland

Abstract

An efficient process for the purification of synthetic phosphorothioate and native DNA sequences is presented. The process is based on the use of an aminopropylated silica gel support functionalized with aminoxyalkyl functions to enable capture of DNA sequences through an oximation reaction with the keto function of a linker conjugated to the 5'-terminus of DNA sequences. Deoxyribonucleoside phosphoramidites carrying this linker, as a 5'-hydroxyl protecting group, have been synthesized for incorporation into DNA sequences during the last coupling step of a standard solid-phase synthesis protocol executed on a controlled pore glass (CPG) support. Solid-phase capture of the nucleobase- and phosphate-deprotected DNA sequences released from the CPG support is demonstrated to proceed near quantitatively. Shorter than full-length DNA sequences are first washed away from the capture support; the solid-phase purified DNA sequences are then released from this support upon reaction with tetra-*n*-butylammonium fluoride in dry dimethylsulfoxide (DMSO) and precipitated in tetrahydrofuran (THF). The purity of solid-phase-purified DNA sequences exceeds 98%. The simulated high-throughput and scalability features of the solid-phase purification process are demonstrated without sacrificing purity of the DNA sequences.

Keywords

Synthetic DNA sequences; Solid-phase purification; High-throughput capability; Large scale purification; Cost-effective process

INTRODUCTION

The use of synthetic DNA sequences and their analogues for recognition and binding to mRNAs encoding disease-causing proteins has led to the development of nucleic acid-based drugs capable of inhibiting the expression of those proteins (Crooke, 2004). This approach to control the expression of disease-causing proteins requires the production of synthetic nucleic acid sequences in sufficient quantities and of high purity to support preclinical and clinical studies. In contrast, total gene synthesis for synthetic biology applications demands minute amounts of numerous highly pure synthetic DNA sequences. Although the phosphoramidite approach to the chemical synthesis of nucleic acid sequences (Beaucage and Caruthers, 1981; Sinha et al., 1984; UNIT 3.3) is efficient and amenable to pharmaceutical scale production, the purification of those sequences is a daunting challenge. Despite the fact that the coupling efficiency of phosphoramidite monomers is near

quantitative, when performed on a controlled pore glass support, the full-length nucleic acid sequences are mixed with impurities consisting of shorter sequences, resulting from incomplete phosphoramidite coupling at each cycle of the nucleic acid sequence assembly. Other process-related impurities include deletion sequences caused by failure to quantitatively prevent the growth of shorter than full-length sequences or to completely remove the 5'-hydroxyl protecting group at each step of the nucleic acid sequence synthesis protocol. It should also be noted that longer than full-length nucleic acid sequences can occur when the activation of phosphoramidite monomers by a weak acid causes, to a limited extent, the premature cleavage of the acid-labile group protecting the 5'-hydroxyl function of the newly extended nucleic acid sequence (Eleuteri et al., 2000; Wei, 2013). Although longer than full-length nucleic acid sequences are produced in minute amounts, the physicochemical similarity of these impurities to the desired nucleic acid sequence makes their removal very challenging. In the context of large-scale production of nucleic acid-based drugs, HPLC-based methods (UNIT 10.5) including reversed-phase and/or anion-exchange HPLC are currently the preferred techniques for purifying nucleic acid sequences. These techniques require high capacity instruments and accessories (e.g., preparative columns) in addition to large volumes of consumables (e.g., buffered aqueous and organic elution solvents). The HPLC-based purification of nucleic acid sequences is also time-consuming given that, depending on the nature of individual nucleic acid sequence, more than one purification run may be required to achieve the level of purity needed for pharmaceutical applications. Furthermore, a major technical drawback associated with any HPLC purification process is the tedious removal of relatively large volumes of aqueous solvents that are generated during the course of the purification procedure. Thus, HPLC-based methods for the purification of nucleic acid sequences are neither cost-effective nor amenable to high-throughput purification processes; only a single nucleic acid sequence can be purified per run unless numerous instruments are dedicated for this purpose.

This unit describes an innovative solid phase method for the purification of synthetic unmodified and modified (i.e., phosphorothioate) DNA sequences. Specifically, the preparation of a capture support that is critically needed for the proposed purification process is detailed in Basic Protocol 1. The chemical synthesis of linkers for solid-phase capture of DNA sequences and for determining the surface density of aminoxy functions covalently bound to the capture support is delineated in Basic Protocol 2. The preparation of 5'-functionalized deoxyribonucleosides and deoxyribonucleoside phosphoramidites for solid-phase capture of DNA sequences is described in Basic Protocol 3. The solid-phase synthesis of 5'-functionalized DNA sequences for solid-phase purification is reported in Support Protocol 1. Capture of crude 5'-functionalized DNA sequences by the capture solid support and release of the purified sequences from the support is outlined in Basic Protocol 4. A scale up process for the solid-phase purification of phosphorothioate DNA sequences is detailed in Alternate Protocol 1.

BASIC PROTOCOL 1. GENERAL PROCEDURE FOR THE PREPARATION OF A CAPTURE SUPPORT

As illustrated in Fig. 1, this protocol provides a general method for the preparation of capture support **3** from commercial 3-aminopropylated silica gel, 7-oxooctanoic acid (**1**) and *O,O'*-1,3-propanediylbishydroxylamine dihydrochloride.

Materials

7-Oxooctanoic acid (**1**, Aldrich)
O,O'-1,3-Propanediylbishydroxylamine dihydrochloride (Aldrich)
1,1'-Carbonyldiimidazole (CDI, Aldrich)
3-Aminopropyl silica gel (~1 mmol NH₂, Aldrich)
Tetrahydrofuran, dry (THF, Acros)
Acetonitrile (MeCN, Aldrich)
Triethylamine (Et₃N, Aldrich)
Cap Mix A (Glen research)
Cap Mix B (Glen Research)
10-mL screw cap glass vials with caps
25- and 50-mL round-bottom flasks
100-mL Erlenmeyer flasks
50-mL glass cylinder
Rubber septa for 14/20-glass joints
Magnetic stir bars
Magnetic stirrer
8-mesh Drierite with indicator (Aldrich)
Drierite guard tube
Reflux condenser (LabGlass)
Vacuum line (ChemGlass)
High-vacuum oil pump (Edwards)
Temperature controlled orbital shaker (Thomas Scientific)
15-mL, 30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity (ChemGlass)
5- and 10-mL Luer-tipped glass syringes
21-G stainless steel syringe needles

Prepare the solid support 2

- 1 Place 7-oxooctanoic acid (**1**, 0.58 g, 3.0 mmol) and a magnetic stir bar into a flame-dried 25-mL round bottom flask. Add 1,1'-carbonyldiimidazole (CDI, 0.49 g, 3.0 mmol). Seal the flask with a reflux condenser and a Drierite drying tube.
- 2 Remove the drying tube. Add 4-mL dry THF through the reflux condenser using a 5-mL glass syringe. Connect the drying tube to the reflux condenser. Stir the solution for 2 hr at ~25 °C using a magnetic stirrer.
- 3 Place commercial 3-aminopropyl silica gel (1.0 g, ~1 mmol NH₂) into a 30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Add 20-mL of 1:5 (v/v) Et₃N/MeCN to the funnel using a 50-mL glass cylinder. Allow the filtrate to drip into the Erlenmeyer. Connect the Buchner funnel to an in-house vacuum line. Apply vacuum to air-dry the support.
- 4 Remove the reflux condenser/drying tube set from the flask of step 2. Add the dry support of step 3 to the solution. Securely seal the flask with a rubber septum. Agitate the suspension at 65 °C for 24 hr using an orbital shaker.
- 5 Transfer the suspension to a 30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Connect the Buchner funnel to an in-house vacuum line. Apply vacuum and wash the support successively with 20-mL THF and 20-mL MeCN using a 50-mL glass cylinder. Air-dry the support.
- 6 Place the support in a flame dried 50-mL round bottom flask. Seal the flask with a rubber septum. Add 20-mL of a commercial solution of acetic anhydride, 1-methylimidazole and pyridine in THF [Cap Mix A:Cap Mix B (1:1 v/v)] through the septum using a 10-mL glass syringe and a 21-G syringe needle as a vent. Remove the needle. Mechanically agitate the solution for 30 min at ~25 °C using an orbital shaker.

This step is required to inactivate amino functions left unreacted on solid support **2**.
- 7 Repeat step 5 with the following exception: Apply vacuum and wash the support twice with 20-mL MeCN using a 50-mL glass cylinder. Air dry the support. Transfer the support to a 25-mL round bottom flask connected to a high vacuum line/oil vacuum pump system. Apply vacuum for 1 hr at ~25 °C. Store the solid support **2** (981 mg) at -20°C until needed.

Loss of solid support **2** occurred when transferring the support from the fritted glass vacuum filtration Buchner funnel to a round-bottom flask (i.e., steps 3 to 4, 5 to 6 and step 7).

Prepare the capture support 3

- 8 Place the solid support **2** (1.0 g) into a 10-mL screw cap glass vial. Add a solution of *O,O'*-1,3-propanediylbishydroxylamine dihydrochloride (537 mg, 3.00 mmol) in 4-mL H₂O using a 5-mL glass syringe. Seal the glass vial and mechanically agitate the suspension for 6 hr at 65 °C using an orbital shaker.
- 9 Repeat step 5 with the following exception: Apply vacuum and wash the support successively with 20-mL DMF and 20-mL MeCN using a 50-mL glass cylinder. Transfer the support to a 25-mL round bottom flask connected to a high vacuum line/oil vacuum pump system. Apply vacuum over a period of 1 hr at ~25 °C. Store the solid support **3** (995 mg) at -20°C until needed.

BASIC PROTOCOL 2. PREPARATION OF CHEMICAL LINKERS FOR THE SOLID-PHASE CAPTURE OF DNA SEQUENCES AND FOR DETERMINING THE SURFACE DENSITY OF AMINOXY FUNCTIONS COVALENTLY BOUND TO SUPPORT 3

Functionalization of the terminal 5'-hydroxyl of DNA sequences with carbonyl-derived linkers is necessary for capture by the solid support **3** upon reaction with its highly reactive aminoxy groups. Fig. 2 depicts the reaction of *N,N'*-dimethylethylenediamine with γ -butyrolactone (**4a**) or 5,5-dimethyldihydrofuran-2-one (**4b**) to provide the aminated amido alcohol **5a** or **5b**. The reaction of CDI-activated 7-oxooctanoic acid **1** with **5a** or **5b** gives the diamidated keto alcohol **6a** or **6b**; monomethoxytritylation of **6a** affords the linker **6c** for measuring the surface density of aminoxy functions covalently attached to the capture support **3**.

Materials

N,N'-Dimethylethylenediamine (Aldrich)
 γ -Butyrolactone (Aldrich)
5,5-Dimethyldihydrofuran-2-one (ACS Scientific)
7-Oxooctanoic acid (**1**, Aldrich)
1,1'-Carbonyldiimidazole (CDI, Aldrich)
4-Methoxytrityl chloride (Aldrich)
Chloroform (CHCl₃, Fisher Scientific)
Dichloromethane (CH₂Cl₂, Fisher Scientific)
Methanol (MeOH, FisherScientific)
Acetonitrile (MeCN, Aldrich)
Tetrahydrofuran, dry (THF; Acros)
Pyridine, dry (Acros)

Toluene (Fisher Scientific)

Triethylamine (Et₃N, Aldrich)

Dimethyl sulfoxide (DMSO, Aldrich)

Deuterated dimethylsulfoxide (DMSO-*d*₆, Aldrich)

Phosphomolybdic acid (Aldrich)

3% Trichloroacetic acid in dichloromethane (Deblock solution, Glen Research)

Anhydrous sodium sulfate (Aldrich)

8-mesh Drierite with indicator (Aldrich)

Drierite guard tube

4- and 25-mL screw cap glass vials with caps

25- 50- and 100-mL round-bottom flasks (Kontes)

Rubber septa for 14/20 glass joints (Aldrich)

100-mL Separatory funnel

50- and 100-mL Erlenmeyer flasks

Reflux condenser (LabGlass)

50-mL glass cylinder

Pipettor (Corning)

1-mL pipette tips (Fisher Scientific)

5- and 10-mL Luer-tipped glass syringes

21-G stainless steel syringe needles

10-mL volumetric flask (Thomas Scientific)

13×100-mm disposable glass tubes

5 × 20-cm glass chromatography columns (LabGlass)

Silica gel (60-Å, 230 to 400 mesh; EMD)

2.5 × 7.5-cm TLC plates precoated with a 250-μm layer of silica gel 60 F254 (EMD)

Pasteur pipettes (VWR)

Vacuum line (LabGlass)

15-mL, 30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity (ChemGlass)

Magnetic stir bars

Magnetic stirrer

Magnetic hot-plate stirrer (Corning)

Dry heat block (VWR)
High-vacuum oil pump (Edwards)
Temperature controlled orbital shaker (Thomas Scientific)
Rotary evaporator connected to a vacuum pump (Büchi)
Hand-held UV₂₅₄ lamp (UVP)
UV/vis spectrophotometer (Agilent Technologies)
NMR spectrometer (Bruker)
Additional reagents and equipment for column chromatography (*APPENDIX 3E*)
and TLC (*APPENDIX 3D*)

Synthesize 4-hydroxy-N-methyl-N-(2-(methylamino)ethyl)butanamide (**5a**)

- 1 Place *N,N'*-Dimethylethylenediamine (3.25 g, 40.0 mmol) and commercial γ -butyrolactone (**4a**, 1.72 g, 20.0 mmol) in a 25-mL screw-cap glass vial. Seal the vial and heat it at 90 °C for 16 hr using a dry heat block.
- 2 Pour the liquid into a 25-mL round-bottom flask. Rinse the vial with 5-mL MeCN and pour the liquid into the round bottom flask. Concentrate to an oil using a rotary evaporator connected to a high vacuum oil pump.

Purify and isolate **5a**

- 3 Use a Pasteur pipette to spread the oily material on top of a 5 × 20-cm chromatography glass column packed with silica gel (~40 g) that has been equilibrated in 9:1 (v/v) CHCl₃/MeOH.
- 4 Elute the column using a gradient of methanol (10% → 20%) in chloroform and collect 6-mL fractions in 13 × 100 mm disposable glass tubes.
- 5 Analyze fractions by TLC (*APPENDIX 3D*) on a 2.5 × 7.5-cm EMD silica gel 60 F₂₅₄ TLC plate using 9:1 (v/v) CHCl₃/MeOH as the eluent. Pool appropriate fractions in an appropriately sized round-bottom flask and rotoevaporate under reduced pressure (~20 mmHg) to afford the linker **5a** (2.32 g, 14.4 mmol) as an oil in a yield of 72%.

The product is visualized on TLC [*R*_f (9:1 (v/v) CHCl₃/MeOH)= 0.2] by soaking the silica gel glass plate in a commercial phosphomolybdic acid solution. The product appears as a dark blue spot upon heating the glass plate on a magnetic hot plate stirrer.

4-Hydroxy-N-methyl-N-(2-(methylamino)ethyl)butanamide (**5a**). ¹H NMR (300 MHz, DMSO-*d*₆): δ 3.39 (dt, *J* = 6.6, 2.0 Hz, 2H), 3.32 (dt, *J* = 6.6, 2.0 Hz, 2H), 2.95 (s, 1.5H), 2.79 (s, 1.5H), 2.59 (t, *J* = 6.6 Hz, 1H), 2.53 (t, *J* = 6.6 Hz, 1H), 2.34 (t, *J* = 7.4 Hz, 1H), 2.29 (t, *J* = 7.4 Hz, 1H), 2.28 (s, 1.5H), 2.25 (s, 1.5H), 1.63 (m, 2H). ¹³C NMR (75 MHz, DMSO-*d*₆): δ 171.92, 171.91, 60.3, 49.7, 49.0, 48.9, 46.6,

36.2, 36.0, 35.4, 33.1, 29.3, 28.6, 28.3, 28.0. +ESI-TOF-HRMS calcd for $C_8H_{18}N_2O_2$ $[M + H]^+$ 175.1400, found 175.1405.

Synthesize 4-hydroxy-N,4-dimethyl-N-(2-(methylamino)ethyl)-pentanamide (5b)

- 6 Repeat step 1 with the exception of replacing **4a** with commercial 5,5-dimethyldihydrofuran-2-one (**4b**, 2.28 g, 20.0 mmol).
- 7 Repeat step 2.

Purify and isolate 5b

- 8 Repeat steps 3 through 5.

The linker **5b** is isolated as an oil (2.99 g, 14.8 mmol) in a yield of 74%. The product is visualized on TLC [R_f (9:1 (v/v) $CHCl_3/MeOH$)=0.2] as reported in the annotation of step 5.

4-Hydroxy-N,4-dimethyl-N-(2-(methylamino)ethyl)-pentanamide (**5b**). 1H NMR (300 MHz, $DMSO-d_6$): δ 4.23 (br s, 1H), 3.33 (q, J = 6.6 Hz, 2H), 2.97 (s, 1.5H), 2.79 (s, 1.5H), 2.61 (t, J = 6.6 Hz, 1H), 2.53 (t, J = 6.6 Hz, 1.5H), 2.37-2.27 (m, 3H), 2.28 (s, 1.5H), 2.25 (s, 1.5H), 1.57 (m, 2H), 1.08 (s, 3H), 1.07 (s, 3H). ^{13}C NMR (75 MHz, $DMSO-d_6$): δ 172.5, 172.4, 79.1, 49.7, 49.0, 48.9, 46.6, 38.7, 38.3, 36.2, 35.9, 33.1, 28.9, 28.1, 27.4. +ESI-TOF-HRMS calcd for $C_{10}H_{22}N_2O_2$ $[M + H]^+$ 203.1754, found 203.1763.

Synthesize N-(2-(4-hydroxy-N-methylbutanamido)ethyl)-N-methyl-7-oxooctanamide (6a)

- 9 Add 1,1'-carbonyldiimidazole (2.52 g, 15.5 mmol) and a magnetic stir bar to a solution of 7-oxooctanoic acid (**1**, 2.45 g, 15.5 mmol) in 15-mL dry THF in a 50-mL round-bottom flask. Seal the flask with a reflux condenser and a Drierite drying tube. Stir the solution using a magnetic stirrer for 2 hr at ~ 25 °C.
- 10 Remove the reflux condenser/drying tube set. Add **5a** (2.32 g, 14.4 mmol) to the solution using a pipettor. Seal the flask with a rubber septum. Place the sealed flask on a magnetic hot plate stirrer. Stir the reaction at 65 °C over a period of 24 hr.
- 11 Remove the flask from the magnetic hot plate stirrer. Remove the rubber septum. Rotoevaporate the liquid to a thick oil under reduced pressure (~ 20 mmHg) using a vacuum pump.
- 12 Add 40-mL $CHCl_3$ using a 50-mL glass cylinder to dissolve the oily material. Pour the solution into a 100-mL separatory funnel.
- 13 Add 20-mL water using a 50-mL glass cylinder. Seal the separatory funnel and shake vigorously. Allow phase separation. Collect the organic phase into a 100-

mL round-bottom flask and rotoevaporate to an oil under reduced pressure (~20 mmHg) using a vacuum pump.

Purify and isolate **6a**

- 14** Use a Pasteur pipette to add 4-mL CHCl₃ and spread the solution on top of a 5 × 20-cm chromatography glass column packed with silica gel (~40 g) that has been equilibrated in CHCl₃.
- 15** Elute the product off the column using a gradient of methanol (0→4%) in chloroform and collect 6-mL fractions in 13 × 100-mm disposable glass tubes.
- 16** Analyze fractions by TLC (APPENDIX 3D) on a 2.5 × 7.5-cm EMD silica gel 60 F₂₅₄ TLC plate using 9:1 (v/v) CHCl₃/MeOH as the eluent. Pool appropriate fractions in an appropriately sized round-bottom flask and rotoevaporate under reduced pressure (~20 mmHg) to afford the linker **6a** (3.57g, 11.3 mmol) as an oil in a yield of 78%.

The linker **6a** is visualized on TLC [*R_f* (9:1 (v/v) CHCl₃/MeOH)= 0.6; **5a**: *R_f* (9:1 (v/v) CHCl₃/MeOH)= 0.2] as reported in the annotation of step 5.

N-(2-(4-Hydroxy-N-methylbutanamido)ethyl)-N-methyl-7-oxooctanamide (**6a**). ¹H NMR (300 MHz, DMSO-d₆): δ 4.22 (m, 1H), 3.47 (s, 0.4H), 3.39 (m, 3H), 3.33 (s, 0.4H), 2.98 (s, 0.6H), 2.96 (s, 0.6H), 2.94 (s, 1H), 2.91 (s, 1H), 2.82 (d, J = 2.1 Hz, 0.6H), 2.80 (s, 1H), 2.40 (t, J = 7.3 Hz, 2H), 2.33-2.21 (m, 3H), 2.17 (t, J = 7.3 Hz, 0.8H), 2.06 (s, 3H), 1.56-1.52 (m, 1.8H), 1.50-1.38 (m, 4.2H), 1.22 (m, 2H), 1.08 (s, 1.4H), 1.06 (s, 3.8H); ¹³C NMR (75 MHz, DMSO-d₆): δ 208.4, 172.9, 172.7, 172.5, 172.4, 172.2, 172.0, 171.8, 171.7, 79.1, 68.25, 68.21, 47.3, 47.2, 46.6, 46.4, 45.8, 45.6, 44.1, 44.0, 42.6, 38.7, 38.2, 36.0, 35.8, 35.1, 35.0, 33.40, 33.37, 33.05, 33.01, 32.4, 32.3, 31.4, 29.6, 29.2, 28.4, 28.3, 28.2, 28.14, 28.06, 27.2, 24.7, 24.6, 24.3, 24.2, 23.1. +ESI-TOF-HRMS calcd for C₁₈H₃₄N₂O₄ [M + H]⁺ 343.2591, found 343.2623.

Synthesize N-(2-(4-hydroxy-N,4-dimethylpentanamido)ethyl)-N-methyl-7-oxooctanamide (**6b**)

- 17** Repeat step 9.
- 18** Repeat step 10 with the exception of replacing **5a** with **5b** (2.91 g, 14.4 mmol).
- 19** Repeat step 11 through 16.

The linker **6b** is isolated as an oil (3.97 g, 11.6 mmol) in a yield of 80%. The product is visualized on TLC [*R_f* (9:1 (v/v) CHCl₃/MeOH)= 0.6; **5b**: *R_f* (9:1 (v/v) CHCl₃/MeOH)= 0.2] as reported in the annotation of step 5.

N-(2-(4-Hydroxy-N,4-dimethylpentanamido)ethyl)-N-methyl-7-oxooctanamide (**6b**). ¹H NMR (300 MHz, DMSO-d₆): δ 4.22 (m, 1H), 3.47 (s, 0.4H), 3.39 (m, 3H), 3.33 (s, 0.4H), 2.98 (s, 0.6H), 2.96 (s, 0.6H), 2.94 (s, 1H), 2.91 (s, 1H), 2.82 (d, J = 2.1 Hz, 0.6H), 2.80 (s, 1H), 2.40 (t, J = 7.3 Hz, 2H), 2.33-2.21 (m, 3H), 2.17 (t, J = 7.3 Hz, 0.8H), 2.06 (s, 3H), 1.56-1.52 (m, 1.8H), 1.50-1.38 (m, 4.2H), 1.22 (m, 2H), 1.08 (s, 1.4H), 1.06 (s, 3.8H). ¹³C NMR (75 MHz, DMSO-d₆): δ 208.4, 172.9, 172.7, 172.5, 172.4, 172.2, 172.0, 171.8, 171.7, 79.1, 68.25, 68.21, 47.3, 47.2, 46.6, 46.4, 45.8, 45.6, 44.1, 44.0, 42.6, 38.7, 38.2, 36.0, 35.8, 35.1, 35.0, 33.40, 33.37, 33.05, 33.01, 32.4, 32.3, 31.4, 29.6, 29.2, 28.4, 28.3, 28.2, 28.14, 28.06, 27.2, 24.7, 24.6, 24.3, 24.2, 23.1. +ESI-TOF-HRMS calcd for C₁₈H₃₄N₂O₄ [M + H]⁺ 343.2591, found 343.2623.

Synthesize N-(2-(4-((4-methoxyphenyl)diphenylmethoxy)-N-methylbutanamido)ethyl)-N-methyl-7-oxooctanamide (**6c**)

- 20 Place linker **6a** (628 mg, 2.00 mmol) and a magnetic stir bar in a 100-mL round bottom flask. Seal the flask with a rubber septum.
- 21 Add sequentially 10-mL dry pyridine and a solution of 4-methoxytrityl chloride (927 mg, 3.00 mmol) in 10-mL dry pyridine through the rubber septum using a 10-mL glass syringe and a 21-G syringe needle as a vent. Remove the syringe needle vent and stir the reaction mixture at ~25 °C over a period of 4 h using a magnetic stirrer.
- 22 Remove the rubber septum. Quench the reaction by adding 20-mL water using a 50-mL glass cylinder.

Extract the crude product

- 23 Add 30-mL CHCl₃ using a 50-mL glass cylinder. Pour the liquid into a 100-mL separatory funnel. Seal the funnel and shake vigorously. Allow phase separation.
- 24 Collect the organic phase and mix with 5 g anhydrous sodium sulfate. Filter the suspension through a 30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL round bottom flask. Connect the Buchner funnel to an in-house vacuum line and apply vacuum. Remove the funnel and rotoevaporate the solvent to an oil under reduced pressure (~20 mmHg) using a vacuum pump.

Purify and isolate the crude linker **6c**

- 25 Use a Pasteur pipette to add 4-mL CHCl₃ and spread the solution on top of a 5 × 20-cm chromatography glass column packed with silica gel (~40 g) that has been equilibrated in a 99.5:0.5 (v/v) CHCl₃/pyridine.
- 26 Elute the product from the column using a gradient of 0 → 2% methanol in a 99.5:0.5 (v/v) CHCl₃/pyridine. Collect 6-mL fractions in 13 × 100-mm disposable glass tubes.

- 27 Analyze fractions by TLC (APPENDIX 3D) on a 2.5 × 7.5-cm EMD silica gel 60 F₂₅₄ TLC plate using 9:1 (v/v) CHCl₃/MeOH as the eluent and UV₂₅₄ for detection. Pool appropriate fractions [**6c**: *R_f* = 0.8; **6a**: *R_f* (9:1 (v/v) CHCl₃/MeOH) = 0.6] in an appropriately sized round-bottom flask and rotoevaporate under reduced pressure (~20 mmHg) using a vacuum pump.

The linker **6c** is needed to determine the concentration or surface density of aminoxy functions covalently bound to the capture support **3**.

- 28 Dissolve the purified material in 5-mL toluene. Rotoevaporate the solvent to dryness using an oil vacuum pump. Repeat this step twice.

Pure linker **6c** (1.01 g, 1.7 mmol) is isolated as an oil in a yield of 85%.

N-(2-(4-((4-Methoxyphenyl)diphenylmethoxy)-N-methylbutanamido)ethyl)-N-methyl-7-oxooctanamide (**6c**). ¹H NMR (300 MHz, DMSO-d₆): δ 7.39-2.29 (m, 8H), 7.25-7.20 (m, 4H), 6.89 (d, *J* = 8.8 Hz, 2H), 3.74 (s, 3H), 3.42 (s, 0.5H), 3.36 (s, 2H), 2.98 (m, 2H), 2.93 (s, 0.7H), 2.90 (s, 0.8H), 2.88 (s, 1H), 2.87 (s, 1H), 2.79 (m, 2H), 2.35 (m, 3H), 2.26-2.12 (m, 3H), 2.04 (s, 1H), 2.03 (s, 1H), 2.02 (s, 1H), 1.77 (m, 2H), 1.42 (m, 4H), 1.19 (m, 2H). ¹³C NMR (75 MHz, DMSO-d₆): δ 208.3, 172.1, 172.0, 171.9, 171.7, 171.6, 171.5, 171.4, 158.1, 144.59, 144.57, 135.4, 129.8, 128.8, 128.1, 127.9, 127.8, 126.7, 125.2, 113.1, 85.5, 62.6, 62.5, 62.4, 55.0, 47.3, 47.2, 46.4, 46.3, 45.8, 45.6, 44.13, 44.10, 42.6, 35.9, 35.8, 35.1, 35.0, 33.4, 33.3, 33.0, 32.34, 32.28, 31.4, 29.58, 29.56, 29.3, 29.2, 28.38, 28.35, 28.3, 28.2, 25.4, 25.1, 24.7, 24.5, 24.3, 24.2, 23.1, 21.0. +ESI-TOF-HRMS calcd for C₃₆H₄₆N₂O₅Cs [M + Cs]⁺ 719.2456, found 719.2458.

Determine the surface density of aminoxy functions covalently bound to capture support **3**

- 29 Place 150 mg of solid support **3** into a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 50-mL Erlenmeyer. Pour 5-mL of 3:7 (v/v) Et₃N/MeCN into the funnel. Allow the filtrate to drip into the Erlenmeyer.
- 30 Pour 10-mL MeCN into the funnel. Connect the Erlenmeyer to an in-house vacuum line. Apply vacuum and air-dry the support. Transfer 20 mg of the support to a 4-mL screw cap glass vial.

Functionalize the capture support **3** with the linker **6c**

- 31 Add a solution of linker **6c** (60 mg, 100 μmol) in 300-μL 5:1 (v/v) DMSO/H₂O. Seal the vial and mechanically agitate it at ~25 °C using an orbital shaker over a period of 24 hr

- 32** Pour the suspension into a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Connect the funnel to an in-house vacuum line and apply vacuum.
- 33** Pipet 10 mL MeCN into the glass funnel. Apply vacuum. Repeat this step twice and air dry the support. Remove the funnel from the Erlenmeyer.

Measure the surface density of the covalently attached aminoxy functions

- 34** Pipet 1-mL 3% trichloroacetic acid in dichloromethane in the glass funnel. Allow the yellow-colored liquid to drip into a 10-mL volumetric flask. Repeat this step nine times. Add dichloromethane as needed to bring the volume to the mark of the 10-mL volumetric flask.
- 35** Measure the absorbance of the solution at 478 nm ($\epsilon = 56 \text{ mL}\cdot\text{cm}^{-1}\cdot\mu\text{mol}^{-1}$) against a blank solution of 3% trichloroacetic acid in dichloromethane using a UV/vis spectrophotometer.

The surface density of aminoxy functions is $146 \pm 7 \mu\text{mol}$ per gram of capture support **3**. This measurement serves as a reference standard to ensure consistency/reproducibility in the preparation of support **3**. The importance of this measurement relates to the determination of optimal amount of support **3** needed to ensure complete capture of the DNA sequences to be purified.

BASIC PROTOCOL 3. PREPARATION OF 5'-FUNCTIONALIZED DEOXYRIBONUCLEOSIDES AND DEOXYRIBONUCLEOSIDE PHOSPHORAMIDITES FOR SOLID-PHASE CAPTURE OF DNA SEQUENCES

The conjugation of the linker **6b** to the 5'-hydroxy function of 2'-deoxythymidine and *N*-protected 2'-deoxyribonucleosides (**7a-d**) upon reaction with dichlorodiisopropylsilane in the presence of imidazole is detailed in this protocol. As schematically shown in Figure 3, conversion of the 5'-functionalized deoxyribonucleosides **8a-d** to their corresponding deoxyribonucleoside phosphoramidite derivatives **9a-d** is also described.

Materials

N-(2-(4-hydroxy-*N*,4-dimethylpentanamido)ethyl)-*N*-methyl-7-oxooctanamide (**6b**, *Basic Protocol 2*)

1*H*-Imidazole (Aldrich)

N,N-Dimethylformamide, dry (DMF, Acros)

N,N-Diisopropylethylamine (Aldrich)

Dichlorodiisopropylsilane (Aldrich)

*N*⁶-Benzoyl-2'-deoxyadenosine (**7a**, ChemGenes)

*N*⁴-Benzoyl-2'-deoxycytidine (**7b**, ChemGenes)

*N*²-Isobutyryl-2'-deoxyguanosine (**7c**, ChemGenes)
2'-Deoxythymidine (**7d**, ChemGenes)
Acetone (Fisher Scientific)
Benzene (C₆H₆, Aldrich)
Deuterated benzene (C₆D₆, Aldrich)
Deuterated dimethyl sulfoxide (DMSO-*d*₆, Aldrich)
Triethylamine (Et₃N, Aldrich)
Ethyl acetate (EtOAc, VWR)
Dichloromethane, dry (CH₂Cl₂, Acros)
Chloroform (CHCl₃, Fisher Scientific)
Methanol (MeOH, Fisher Scientific)
Isopropyl ether (Aldrich)
2-Cyanoethyl *N,N*-diisopropylchlorophosphoramidite (Aldrich)
Sodium bicarbonate (NaHCO₃, Fisher Scientific)
Anhydrous sodium sulfate (Fisher Scientific)
Argon or nitrogen sources
30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity (ChemGlass)
50- and 100-mL round-bottom flasks (Kontes)
100-mL Erlenmeyer flask (Kontes)
Rubber septa for 14/20 glass joints (Aldrich)
50-mL glass cylinder
100-mL separatory funnel
Chromatography glass column (LabGlass)
Silica gel (60-Å, 230 to 400 mesh; EMD)
13×100-mm disposable glass tubes
Pasteur pipettes (VWR)
1- and 10-mL glass syringes
21-G stainless steel syringe needles
Magnetic stir bars
Magnetic stirrer
Pipettor (Corning)
1-mL pipet tips (Fisher Scientific)

Dry ice-isopropyl ether bath
Dry ice-acetone bath
Wet ice bath
Rotary evaporator connected to a vacuum pump (Büchi)
Hand-held UV₂₅₄ lamp (UVP)
UV/vis spectrophotometer (Agilent Technologies)
NMR spectrometer (Bruker)
Lyophilizer

Synthesize the 5'-functionalized deoxyribonucleosides 8a-d

- 1 Place the linker **6b** (780 mg, 2.28 mmol), imidazole (184 mg, 2.70 mmol) and a magnetic stir bar in a flame-dried 50-mL round bottom flask.
- 2 Add, under an argon atmosphere, 5-mL dry DMF using a 10-mL syringe and pipet 2.35 mL *N,N*-diisopropylethylamine (13.5 mmol) using a pipettor. Seal the flask with a rubber septum and cool the solution to ~0 °C using a crushed wet ice bath.
- 3 Add 730- μ L dichlorodiisopropylsilane (4.50 mmol) through the rubber septum using a 1-mL glass syringe. Stir the cold solution for 1 hr at ~0 °C using a magnetic stirrer.
- 4 Allow the reaction mixture to warm up to room temperature while stirring over a period of 4 hr. Cool the reaction mixture to ~ -60 °C in a dry ice-isopropyl ether bath.
- 5 Add dropwise, over a period of 1 min, a 5-mL solution of 2'-deoxythymidine (**7a**, 1.92 g, 5.40 mmol) and imidazole (368 mg, 5.40 mmol) in dry DMF, through the rubber septum using a 10-mL glass syringe. Stir the solution for 1 hr at ~ -60 °C. Allow the reaction mixture to warm up to ~0 °C in a crushed wet ice bath while keeping stirring it at this temperature for 3 h.

Quench the reaction

- 6 Remove the rubber septum and pour 35-mL of a cold (~0 °C) 5% aqueous solution NaHCO₃.

Extract the product

- 7 Pour the aqueous solution into a 100-mL separatory funnel. Add 40-mL EtOAc into the separatory funnel using a glass cylinder. Seal the funnel and shake it vigorously. Allow complete phase separation.
- 8 Collect the organic layer into a 100-mL round bottom flask and rotoevaporate to an oil under reduced pressure (~ 20 mmHg) using a vacuum pump.

Purify and isolate the 5'-functionalized deoxyribonucleosides 8a

- 9 Use a Pasteur pipette to add 4-mL CHCl₃ and spread the solution on top of a 5 × 20-cm chromatography glass column packed with silica gel (~40 g) that has been equilibrated in CHCl₃.
- 10 Elute the product from the column using a gradient of 0 → 6% methanol in chloroform. Collect 6-mL fractions in 13 × 100-mm disposable glass tubes.
- 11 Analyze the fractions by TLC (APPENDIX 3D) on a 2.5 × 7.5-cm EMD silica gel 60 F₂₅₄ TLC plate using 9:1 (v/v) CHCl₃/MeOH as the eluent and UV₂₅₄ for detection. Pool appropriate fractions [*R_F*(9:1 (v/v) CHCl₃/MeOH) = 0.7] in an appropriately sized round-bottom flask and rotoevaporate the volatiles under reduced pressure (~20 mmHg) using a vacuum pump.

The purified 5'-functionalized deoxyribonucleoside **8a** (1.32 g, 1.59 mmol) is isolated as a solid in a yield of 70%.

8a. ¹H NMR (300 MHz, DMSO-d₆): δ 11.23 (s, 1H), 8.74 (s, 1H), 8.62 (s, 0.7H), 8.59 (s, 0.3H), 8.32 (s, 1H), 8.07 (d, J = 7.8 Hz, 2H), 7.65 (dd, J = 7.6, 7.4 Hz, 1H), 7.56 (dd, J = 7.6, 7.4 Hz, 2H), 6.51 (t, J = 6.8 Hz, 1H), 5.45 (m, 1H), 4.54 (br s, 1H), 3.98 (m, 2H), 3.85 (m, 1H), 3.38 (m, 4H), 2.92 (m, 4H), 2.79 (m, 2H), 2.38 (t, J = 7.3 Hz, 2H), 2.31-2.13 (m, 4H), 2.05 (s, 3H), 1.66 (m, 2H), 1.43 (m, 4H), 1.21 (m, 7H), 0.97 (m, 16H). ¹³C NMR (75 MHz, DMSO-d₆): δ 208.3, 172.4, 172.2, 172.1, 172.0, 171.9, 171.8, 171.71, 171.66, 165.6, 151.9, 151.4, 150.3, 143.0, 142.8, 133.4, 132.4, 128.5, 128.4, 125.9, 125.8, 87.2, 83.74, 83.66, 79.2, 73.2, 73.1, 70.1, 70.4, 68.3, 68.2, 62.9, 62.7, 47.2, 46.5, 46.4, 45.8, 45.7, 44.1, 44.0, 42.6, 38.5, 35.9, 35.8, 35.0, 33.5, 33.3, 33.2, 33.0, 32.4, 32.3, 31.5, 29.6, 29.3, 29.2, 28.4, 28.3, 28.2, 28.1, 28.0, 27.3, 27.1, 24.7, 24.6, 24.3, 24.2, 23.1, 17.60, 17.52, 17.47, 17.43, 17.3, 17.21, 17.18, 17.11, 17.0, 12.9, 12.7, 12.6, 12.32, 12.29. +ESI-TOF-HRMS calcd for C₄₁H₆₃N₇O₈SiNa [M + Na]⁺ 832.4399, found 832.4411.

- 12 Repeat steps 1 through 11 with the following exception:

Replace the deoxyribonucleoside **7a** in step 5 with **7b** (1.79 g, 5.40 mmol), **7c** (1.82 g, 5.40 mmol) or **7d** (1.31 g, 5.40 mmol). The purified 5'-functionalized deoxyribonucleoside **8b** (1.16 g, 1.48 mmol), **8c** (890 mg, 1.13 mmol) or **8d** (1.08 g, 1.54 mmol) is isolated as a solid in a yield of 65%, 50%, and 68% respectively.

8b. ¹H NMR (300 MHz, DMSO-d₆): δ 11.25 (s, 1H), 8.27 (d, J = 7.3 Hz, 1H), 8.01 (d, J = 7.3 Hz, 2H), 7.63 (m, 1H), 7.52 (m, 2H), 7.37 (d, J = 7.3 Hz, 1H), 6.16 (t, J = 6.2 Hz, 1H), 5.37 (d, J = 4.4 Hz, 1H), 4.28 (m, 1H), 4.01-3.89 (m, 3H), 3.47 (s, 0.4H), 3.38 (m, 3H), 2.98 (s, 0.7H), 2.95 (s, 0.7H), 2.94 (s, 1.1H), 2.90 (s, 1.1H), 2.81 (m, 2H), 2.38 (m, 5H), 2.23 (t, J = 7.3 Hz, 2H), 2.15 (m, 2H), 2.05 (s, 3H), 1.70 (m, 2H), 1.43 (m, 4H), 1.27 (s, 6H), 1.20 (m, 2H), 1.02 (m, 14H). ¹³C

NMR (75 MHz, DMSO- d_6): δ 208.29, 208.27, 172.4, 172.2, 172.0, 171.9, 171.8, 171.70, 171.66, 167.3, 163.0, 154.3, 144.3, 133.1, 132.7, 128.41, 128.38, 95.8, 87.1, 86.0, 73.34, 73.30, 69.3, 62.3, 47.3, 47.2, 46.5, 46.4, 45.8, 45.7, 44.1, 44.0, 42.6, 40.8, 35.9, 35.8, 35.1, 35.0, 33.5, 33.3, 33.2, 33.0, 32.4, 32.3, 31.5, 29.6, 29.39, 29.37, 28.4, 28.3, 28.2, 28.13, 28.07, 27.3, 27.1, 24.7, 24.6, 24.31, 24.25, 23.1, 17.61, 17.59, 17.51, 17.48, 12.8, 12.6. +ESI-TOF-HRMS calcd for $C_{40}H_{63}N_5O_9SiNa$ [M + Na]⁺ 808.4287, found 808.4298.

8c. ¹H NMR (300 MHz, DMSO- d_6): δ 12.07 (s, 1H), 11.68 (s, 1H), 8.16 (s, 1H), 6.22 (t, J = 6.6 Hz, 1H), 5.37 (d, J = 4.2 Hz, 1H), 4.42 (m, 1H), 3.88 (d, J = 7.8 Hz, 1H), 3.83 (m, 2H), 3.44-3.34 (m, 4H), 2.95 (s, 0.8H), 2.94 (s, 0.8H), 2.91 (s, 1.2H), 2.89 (s, 1.2H), 2.77 (m, 3H), 2.63 (m, 2H), 2.38 (t, J = 7.3 Hz, 2H), 2.33-2.13 (m, 5H), 2.05 (s, 3H), 1.66 (m, 2H), 1.43 (m, 4H), 1.21 (m, 7H), 1.12 (d, J = 6.8 Hz, 6H), 0.95 (m, 14H). ¹³C NMR (75 MHz, DMSO- d_6): δ 208.3, 180.1, 172.4, 172.2, 172.1, 172.0, 171.9, 171.8, 171.69, 171.66, 154.8, 148.3, 148.0, 137.2, 128.3, 120.3, 87.1, 82.8, 73.2, 73.1, 70.0, 62.9, 47.3, 47.2, 46.5, 46.3, 45.7, 45.6, 44.1, 44.0, 42.6, 35.9, 35.8, 35.05, 35.02, 34.7, 33.5, 33.3, 33.2, 33.0, 32.4, 32.3, 31.4, 29.6, 29.3, 28.4, 28.3, 28.2, 28.1, 28.0, 27.2, 27.1, 24.7, 24.5, 24.3, 24.2, 23.1, 18.83, 18.78, 17.55, 17.52, 17.48, 17.43, 17.1, 12.7, 12.6. +ESI-TOF-HRMS calcd for $C_{38}H_{65}N_7O_9SiNa$ [M + Na]⁺ 814.4505, found 814.4518.

8d. ¹H NMR (300 MHz, DMSO- d_6): δ 11.33 (s, 1H), 7.42 (s, 1H), 6.17 (t, J = 6.8 Hz, 1H), 5.31 (d, J = 4.5 Hz, 1H), 4.25 (m, 1H), 3.85 (m, 3H), 3.46 (s, 0.4H), 3.39 (s, 2 H), 3.36 (m, 1.1H), 2.97 (s, 0.8H), 2.95 (s, 0.7H), 2.92 (s, 1H), 2.90 (s, 1H), 2.80 (m, 2H), 2.39 (t, J = 7.3 Hz, 2H), 2.13 (m, 1H), 2.22 (m, 2H), 2.12 (m, 3H), 2.06 (s, 3H), 1.76 (s, 3H), 1.67 (m, 2H), 1.44 (m, 4H), 1.21 (m, 8H), 1.00 (m, 13H), 0.95 (s, 1H), 0.93 (s, 1H). ¹³C NMR (75 MHz, DMSO- d_6): δ 208.3, 172.4, 172.2, 172.0, 171.9, 171.8, 171.71, 171.67, 163.1, 150.3, 135.6, 109.4, 86.5, 83.6, 73.3, 73.2, 70.0, 62.8, 47.3, 47.2, 46.5, 46.4, 45.8, 45.7, 44.1, 44.0, 42.6, 35.9, 35.8, 35.1, 35.0, 33.5, 33.3, 33.1, 33.0, 32.4, 32.3, 31.5, 29.6, 29.4, 28.4, 28.3, 28.2, 28.1, 28.0, 27.3, 27.1, 24.7, 24.6, 24.3, 24.2, 23.1, 17.60, 17.56, 17.50, 17.48, 17.3, 12.8, 12.63, 12.59, 12.7. +ESI-TOF-HRMS calcd for $C_{34}H_{60}N_4O_9SiNa$ [M + Na]⁺ 719.4022, found 719.4031.

Synthesize the 5'-functionalized deoxyribonucleoside phosphoramidites 9a-d

- 13** Place a magnetic stir bar into a flame-dried 100-mL round-bottom flask. Seal the vessel with a rubber septum.
- 14** Add a solution of the 5'-functionalized 2'-deoxythymidine (**8a**, 810 mg, 1.00 mmol) in 10-mL anhydrous CH_2Cl_2 through the rubber septum using a 10-mL glass syringe and a 21-G syringe needle as a vent.

- 15 Add 335- μ L 2-cyanoethyl *N,N*-diisopropylchlorophosphoramidite (1.50 mmol) using a 1-mL glass syringe. Stir the reaction mixture at $\sim 25^{\circ}\text{C}$ for 3 hr using a magnetic stirrer.

Quench the reaction

- 16 Add 10-mL water through the rubber septum using a 10-mL glass syringe and a 21-G syringe needle as a vent. Remove the rubber septum and add 25-mL CH_2Cl_2 in the flask using a 50-mL glass cylinder.

Extract the product

- 17 Pour the liquid into a 100-mL separatory funnel. Seal the funnel and shake vigorously. Allow phase separation.
- 18 Collect the organic layer in a 100-mL Erlenmeyer and mix with 5 g anhydrous sodium sulfate. Filter the suspension through a 30-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL round bottom flask. Connect the funnel to an in-house vacuum line. Apply vacuum.
- 19 Remove the funnel and rotoevaporate the solvent to an oil under reduced pressure (~ 20 mmHg) using a vacuum pump.

Purify and isolate the 5'-functionalized deoxyribonucleoside phosphoramidites **9a**

- 20 Use a Pasteur pipette to add 4-mL 9:1 (v/v) $\text{C}_6\text{H}_6/\text{Et}_3\text{N}$ and disperse the solution on top of a 5×20 -cm chromatography glass column packed with ~ 25 g silica gel that has been equilibrated in 9:1 (v/v) $\text{C}_6\text{H}_6/\text{Et}_3\text{N}$.
- 21 Elute the product from the column using 9:1 (v/v) $\text{C}_6\text{H}_6/\text{Et}_3\text{N}$ as the eluent. Collect 6-mL fractions in 13×100 -mm disposable glass tubes.
- 22 Analyze the fractions by TLC (APPENDIX 3D) on a 2.5×7.5 -cm EMD silica gel 60 F₂₅₄ TLC plate using 9:1 (v/v) $\text{C}_6\text{H}_6/\text{Et}_3\text{N}$ as the eluent and UV₂₅₄ for detection. Pool appropriate fractions [$R_f(9:1 \text{ (v/v) } \text{C}_6\text{H}_6/\text{Et}_3\text{N}) = 0.7$] into a 100-mL round bottom flask and rotoevaporate the solvents to a colorless oil under reduced pressure (~ 20 mmHg) using a vacuum pump.

Lyophilize the 5'-functionalized deoxyribonucleoside phosphoramidite **9a**

- 23 Dissolve the purified phosphoramidite **9a** in 10-mL dry benzene. Freeze the solution by manually swirling it in a dry-ice/acetone bath. Lyophilize the frozen solution under high vacuum using an oil vacuum pump connected to a glass trap cooled in liquid nitrogen-filled Dewar flask.

Triethylamine-free **9a** (837 mg, 830 μ mol) is isolated as a colorless oil in a yield of 83%.

9a. ^{31}P NMR (121 MHz, C_6D_6): δ 149.0, 148.8, 148.76, 148.72. +ESI-TOF-HRMS calcd for $\text{C}_{50}\text{H}_{80}\text{N}_9\text{O}_9\text{PSi}$ [$\text{M} + \text{H}$]⁺ 1010.5659, found 1010.5675.

24 Repeat steps 13 through 23 with the following exception:

Replace the 5'-functionalized deoxyribonucleoside **8a** in step 14 with **8b** (785 mg, 1.00 mmol), **8c** (791 mg, 1.00 mmol) or **8d** (696 mg, 1.00 mmol).

The purified triethylamine-free 5'-functionalized deoxyribonucleoside phosphoramite **9b** (788 mg, 800 μ mol), **9c** (743 mg, 795 μ mol) or **9d** (717 mg, 800 μ mol) is isolated as a colorless oil in a yield of 80%, 75%, and 80% respectively.

9b. ^{31}P NMR (121 MHz, C_6D_6): δ 148.3, 148.25, 148.22, 148.1, 147.91, 147.89. +ESI-TOF-HRMS calcd for $\text{C}_{49}\text{H}_{80}\text{N}_7\text{O}_{10}\text{PSi}$ [$\text{M} + \text{H}$] $^+$ 986.5546, found 986.5548.

9c. ^{31}P NMR (121 MHz, C_6D_6): δ 148.77, 148.73, 148.0, 147.92, 147.88. +ESI-TOF-HRMS calcd for $\text{C}_{47}\text{H}_{82}\text{N}_9\text{O}_{10}\text{PSi}$ [$\text{M} + \text{H}$] $^+$ 992.5764, found 992.5773.

9d. ^{31}P NMR (121 MHz, C_6D_6) δ 148.7, 148.6, 148.49, 148.46, 148.37. +ESI-TOF-HRMS calcd for $\text{C}_{43}\text{H}_{77}\text{N}_6\text{O}_{10}\text{PSiCs}$ [$\text{M} + \text{Cs}$] $^+$ 1029.4257, found 1029.4266.

SUPPORT PROTOCOL 1. SOLID-PHASE SYNTHESIS OF 5'-FUNCTIONALIZED PHOSPHOROTHIOATE AND NATIVE DNA SEQUENCES

This protocol delineates the synthesis and deprotection of 5'-functionalized phosphorothioate and native DNA sequences (Fig. 4).

Additional Materials

Long-chain alkylamine controlled-pore glass (500 Å LCAA-CPG or 2000 Å LCAA-CPG) support functionalized with 1- μ mol 5'-*O*-(4,4'-dimethoxytrityl)-2'-deoxythymidine as the leader nucleoside (Glen Research).

5'-*O*-(4,4'-dimethoxytrityl)-dA^{Bz}, -dG^{iBu}, -dC^{Bz} and -dT phosphoramidite monomers (ChemGenes)

5'-Functionalized deoxyribonucleoside phosphoramidites **9a-d** (*Basic Protocol 3*)

Anhydrous acetonitrile diluent (MeCN, Glen Research)

Reagents for oligonucleotide synthesis (all available from Glen Research):

0.45 M 1*H*-tetrazole or 0.25 M 5-ethylthio-1*H*-tetrazole in MeCN

Cap A solution: acetic anhydride in tetrahydrofuran (THF)/pyridine

Cap B solution: 1-methylimidazole in THF

Oxidation solution: 0.05 M 3*H*-1,2-benzodithiol-3-one 1,1-dioxide in MeCN or 0.02 M iodine in THF/pyridine/water

Deblocking solution: 3% trichloroacetic acid in CH_2Cl_2

Concentrated aqueous ammonia (NH₃, Aldrich)
4-mL screw cap glass vial with caps (Fisher Scientific)
1-mL glass syringe
1-mL plastic syringe
21-G stainless steel syringe needles
Heat block (VWR)
DNA/RNA synthesizer (394 DNA/RNA synthesizer, Applied Biosystems)

Synthesize the DNA Sequences

- 1 Perform a 1- μ mol scale solid-phase synthesis of the phosphorothioate DNA sequence **10a–d** (Fig. 4) using an Applied Biosystems 394 DNA/RNA synthesizer in trityl-on mode according to the manufacturer's instructions (also see *APPENDIX 3C*).

*Each synthesis is performed using a long-chain alkylamine controlled-pore glass (500 Å LCAA-CPG or 2000 Å LCAA-CPG for the native DNA sequence **10f**) support functionalized with 5'-O-(4,4'-dimethoxytrityl)-2'-deoxythymidine as the leader nucleoside. The syntheses are carried out using a DNA/RNA synthesizer and commercial appropriately protected deoxyribonucleoside phosphoramidites, which are each dissolved in dry MeCN to give a 0.1 M solution. Each 5'-functionalized deoxyribonucleoside phosphoramidite **9a–d** is dissolved in dry MeCN to provide a 0.15 M solution. Each solution is placed in a distinct vial connected to the DNA/RNA synthesizer through an additional delivery port. Commercial 1H-tetrazole solution was used for phosphoramidite activation in the solid-phase synthesis of DNA sequences **10a–f**. The reaction times for the coupling, capping, oxidation and 5'-deblocking steps in the synthesis of native and phosphorothioate DNA sequences are 120, 60, 60 s and 60 s, respectively. However, the capping step in the synthesis of phosphorothioate DNA sequences (**10a–d**) is performed after the oxidative sulfuration step, which is effected using 0.05 M 3H-1,2-benzodithiol-3-one 1,1-dioxide in MeCN (Iyer et al., 1990); the standard 0.02 M iodine solution in THF/pyridine/water is employed in the oxidation step of native DNA sequences. The last coupling reaction of each synthesis is performed using any of the activated 5'-functionalized deoxyribonucleoside phosphoramidites **9a–d** over a period of 180 s.*

Deprotect the 5'-functionalized DNA sequences

- 2 Transfer the LCAA-CPG-linked 5'-functionalized DNA sequence **10a** from the synthesis column to a 4-mL screw cap glass vial. Add 1-mL of concentrated

aqueous ammonia using a 1-mL glass syringe. Tightly seal the vial and place it in a heat block kept at ~65 °C over a period of 16 hr.

- 3 Allow the vial to cool to ~25 °C. Carefully open the vial and transfer the ammoniacal solution to an empty 4-mL screw cap glass vial using a 1-mL plastic syringe. Evaporate the solution to approximately half its original volume using a stream of air.

The deprotected DNA sequence **10a** can be used immediately for solid-phase purification or kept frozen at -20 °C in the sealed glass vial for subsequent solid-phase purification.

- 4 Repeat steps 2 and 3 for deprotection of each crude DNA sequence **10b–f**.

BASIC PROTOCOL 4. GENERAL PROCEDURE FOR THE SOLID-PHASE CAPTURE AND RELEASE OF PHOSPHOROTHIOATE AND NATIVE DNA SEQUENCES

The details of the capture of crude 5'-functionalized DNA sequences **10a–f** by the solid support **3** are provided in this protocol. As shown in Fig 5, capture of the DNA sequences occurs from the reaction of the aminoxy groups of support **3** with the keto function of each 5'-functionalized DNA sequences **10a–f** and formation of stable oxime conjugates on the surface of solid support **11a–f**. Fig. 5 also shows the release of the captured DNA sequences from support **11a–f** upon treatment with tetra-*n*-butylammonium fluoride to provide the highly pure DNA sequences **12a–f** (Fig. 6).

Materials

Capture solid support **3** (*Basic Protocol 1*)

5'-Functionalized DNA sequences **10a–f** (*Basic Protocol 3*)

Acetonitrile (MeCN, Aldrich)

Dimethyl sulfoxide (DMSO, Aldrich)

Peroxide-free tetrahydrofuran (THF, Aldrich)

Triethylamine (Aldrich)

Methoxytrimethylsilane

2 M Triethylammonium acetate solution (Applied Biosystems)

Tetra-*n*-butylammonium fluoride hydrate (Aldrich)

Tetra-*n*-butyl ammonium chloride (Aldrich)

Concentrate aqueous ammonia (NH₃, Aldrich)

4-mL screw cap glass vial with caps (Fisher Scientific)

1.5-mL microcentrifuge tubes (Fisher Scientific)

15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity (ChemGlass)
100-mL Erlenmeyer flask (Fisher Scientific)
25-mL round-bottom flask
1- and 10-mL plastic syringes (Fisher Scientific)
1- and 10-mL glass syringes (Fisher Scientific)
21-G stainless steel syringe needles
Pasteur pipettes (VWR)
Pipettor (Corning)
1-mL pipet tips (Fisher Scientific)
Spatula
Temperature-controlled orbital shaker (Thomas Scientific)
Heath block (VWR)
RP-HPLC (Agilent Technologies)
5 μ m Supelcosil 25 cm \times 4.6 mm LC-18S column (Supelco)
Eppendorf microcentrifuge
Rotary evaporator connected to a vacuum pump (Büchi)
High vacuum oil pump

Capture the crude 5'-functionalized phosphorothioate or native DNA sequences

- 1** Capture support **3** (150 mg) is added to a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Add a 10-mL 1:5 (v/v) Et₃N/MeCN using a 10-mL plastic syringe. Allow the filtrate to drip into the Erlenmeyer. Connect the funnel to an in-house vacuum line and apply vacuum to air-dry the capture support.
- 2** Transfer the dry support to a 4-mL screw cap glass vial. Add successively solid tetra-*n*-butylammonium chloride (14 mg, 50 μ mol) using a spatula and the solution (~ 500 μ L in water) of crude 5'-functionalized DNA sequence **10a** (170 OD₂₆₀) obtained from step 3 of Support Protocol 1 using a 1-mL plastic syringe. Seal the flask and agitate the suspension at 65 °C over a period of 3 hr using an orbital shaker.
- 3** Stop the shaker. Allow the suspension to settle. Verify the near completeness of the capture reaction by RP-HPLC analysis of a 1- μ L aliquot of the supernatant.

The completeness of the capture reaction is made based on a side-by-side comparison of the RP-HPLC chromatographic profile of crude 5'-functionalized DNA sequence **10a**, prior to capture, with the RP-HPLC

profile of the 1- μ L supernatant aliquot recorded after capture of **10a**. RP-HPLC analyses have been performed using a 5 μ m Supelcosil LC-18S column (25 cm \times 4.6 mm) under the following conditions: starting from 0.1 M (pH 7.0), a linear gradient of 2.5% MeCN/min is pumped at a flow rate of 1 mL/min for 40 min.

- 4 Repeat step 1 through 3 for solid-phase capture of each 5'-functionalized DNA sequence **10b–f**.
- 5 Transfer the suspended solid support **11a** (Fig. 5) to a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Add 10-mL of 1:1 (v/v) CH₃CN/H₂O. Connect the funnel to an in-house vacuum line. Apply vacuum to air-dry the support.
- 6 Transfer the support to a 4-mL screw cap glass vial. Add a 1-mL of 1:1 (v/v) concentrated aqueous NH₃/CH₃CN using a 1-mL glass syringe. Seal the vial and place it in a heat block for 30 min at 65 °C. Allow the flask to cool to ~25 °C.
- 7 Transfer the suspended solid support **11a** to a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Connect the funnel to an in-house vacuum line and apply vacuum to air-dry the support.
- 8 Repeat steps 6 and 7 once.
- 9 Add 1-mL DMSO to the Buchner funnel using a 1-mL glass syringe. Apply vacuum. Repeat step 9 four more times and air-dry the support **11a**.

The washing steps 5 through 9 are necessary to wash away the shorter than full length DNA sequences that might have adsorbed to the support **11a**.

Release of the DNA sequence **12a** from the Support **11a**

- 10 Transfer the air-dried solid support to a 4-mL screw cap glass vial. Add 0.5-mL of a 1.0 M tetra-*n*-butylammonium fluoride in dry DMSO using a 1-mL glass syringe. Seal the flask and place it in a heat block over a period of 3 hr at 65 °C. Allow the glass vial to cool to ~25 °C.

1.0 M tetra-*n*-butylammonium fluoride in dry DMSO is prepared by mixing solid tetra-*n*-butylammonium fluoride hydrate (262 mg) in 1.0-mL dry DMSO.
- 11 Open the glass vial. Add 200 μ L methoxytrimethylsilane and 200 μ L MeCN to the suspension. Seal the vial and keep it at ~25 °C for 30 min.

Methoxytrimethylsilane is used to consume excess fluoride ions.
- 12 Pour the suspension into a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 25-mL round-bottom flask. Add a 10-mL of 1:5 (v/v) Et₃N/MeCN using a 10-mL glass syringe. Allow the filtrate to drip into

the flask. Connect the funnel to an in-house vacuum line. Apply vacuum to air-dry the support.

- 13 Repeat step 6.
- 14 Transfer the suspended solid support **11a** to a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to the 25-mL round bottom flask used in step 12. Allow the filtrate to drip into the flask. Apply vacuum to air-dry the support.
- 15 Repeat successively step 6 and 14. Remove the Buchner funnel. Concentrate the filtrates collected in the 25-mL round bottom flask under reduced pressure to a volume of ~100 μ L using a high vacuum oil pump.

The washing steps 12 through 15 are necessary to ensure that the full length DNA sequence **12a** is completely washed away from the support **11a**.

Precipitate the DNA sequence **12a**

- 16 Add 1-mL peroxide-free THF to the concentrated material. Swirl the flask to yield a uniform precipitate of the solid-phase purified DNA sequence **12a**.
- 17 Pipet the precipitate into a 1.5-mL microcentrifuge tube using a Pasteur pipette. Centrifuge at 16,000*g* for 15 min at 25 °C.
- 18 Carefully remove the supernatant using a pipettor. Resuspend the DNA pellet in 1-mL peroxide-free THF. Centrifuge at 16,000*g* for 15 min at 25 °C.
- 19 Repeat step 18 two more times.
- 20 Carefully remove the supernatant using a pipettor. Dry the DNA pellet under reduced pressure using a SpeedVac system. Store the dried DNA pellet at -20 °C until further use.
- 21 Repeat steps 5 through 20 with each of the solid support **11b–f** to provide the solid-phase purified DNA sequences **12b–f**.

The solid-phase purification of DNA sequence **10a** demonstrates that the sequences **10b–f** can each be, technically, solid-phase purified in a parallel or high-throughput manner. The yields of solid-phase purified DNA sequences depend of several factors including: (i) composition (purine or pyrimidine rich; modified or unmodified) and length of the DNA sequence; (ii) solubility of the sequence in the capture solution, (iii) quality of the phosphoramidite **9a–d** and (iv) performance of the DNA synthesizer. Under optimal conditions, the yield and purity of solid-phase purified DNA sequences are comparable or better than that obtained by RP-HPLC, as discussed in the background section of this report.

ALTERNATE PROTOCOL 1. SCALABLE SOLID-PHASE PURIFICATION OF A RELEVANT DNA SEQUENCE

This protocol provides evidence that the solid-phase purification of synthetic DNA sequences can be qualified as a scalable process. This is demonstrated through the solid-phase synthesis of ten identical 5'-functionalized phosphorothioate DNA sequences, each on a 1 μ mol scale. These individual sequences will be purified together using a 10-fold increase of the amount of capture support **3** typically needed for a 1 μ mol scale solid-phase purification process.

Materials

5'-Functionalized phosphorothioate DNA sequences (**10a**, *Support Protocol 1*)
Capture support **3** (*Basic Protocol 1*)
Acetonitrile (MeCN, Aldrich)
Dimethyl sulfoxide (DMSO, Aldrich)
Dimethyl sulfoxide, dry (Acros)
Tetra-*n*-butylammonium chloride (Aldrich)
Tetra-*n*-butylammonium fluoride hydrate (Aldrich)
Methoxytrimethylsilane (Aldrich)
Peroxide-free tetrahydrofuran (THF, Aldrich)
Triethylamine (Et₃N, Aldrich)
25- and 100-mL round-bottom flasks
100-mL Erlenmeyer flask
20-mL screw cap glass vial with caps (Kimble)
10-mL glass cylinder
15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity (ChemGlass)
5-mL plastic syringe
5- and 10-mL glass syringes
21-G stainless steel syringe needles (Fisher Scientific)
Pipettor (Corning)
1-mL pipette tips (Fisher Scientific)
Heat block (VWR)
DNA/RNA synthesizer (394 DNA/RNA synthesizer, Applied Biosystems)
Rotary evaporator connected to a vacuum pump (Büchi)
Temperature-controlled orbital shaker (Thomas Scientific)

Vortexer (VWR)

Synthesize the DNA Sequence 10a

- 1 Repeat step 1 of Support Protocol 1.

The solid-phase synthesis of each of the ten phosphorothioate DNA sequences (**10a**) has been achieved on a 1- μ mole scale using a DNA/RNA synthesizer capable of assembling four DNA sequences in a single run. Thus, two synthesis runs of four identical sequences and a single run of two identical sequences have been performed to provide **10a** from ten individual sequences.

Deprotect the 5'-functionalized DNA sequence 10a

- 2 Repeat step 2 of Support Protocol 1 for each of the ten individual sequences.
- 3 Allow each of the ten vials to cool to ~ 25 °C. Carefully open each vial and pool each ammoniacal solution together into a 25-mL round-bottom flask using a 5-mL plastic syringe. Rotoevaporate to half (~ 5 mL) of the original volume (~ 10 mL) under reduced pressure (~ 20 mmHg).
- 4 Add 5.0-mL of a 0.1 M solution of tetra-*n*-butylammonium chloride in 1:1 (v/v) DMSO/H₂O to the crude DNA solution of step 3 using a 5-mL glass syringe.

Capture the crude 5'-functionalized phosphorothioate DNA sequence 10a

- 5 Repeat step 1 of Basic Protocol 4 using 1.50 g of capture support **3**.

The amount of capture support **3** is a 10-fold increase of the amount of support used (150 mg) for the solid-phase purification of **10a** that was synthesized on a 1- μ mol scale (see step 1 of Support Protocol 1)
- 6 Place the capture support **3** of step 5 into a 20-mL screw cap glass vial. Pour the crude DNA solution of step 4 into the vial. Tightly seal the vial and agitate the suspension over 3 hr at 65 °C using an orbital shaker.
- 7 Repeat step 3 of Basic Protocol 4 using 0.5 μ L aliquot of the supernatant.
- 8 Transfer the suspended solid support **11a** to a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL Erlenmeyer. Connect the funnel to an in-house vacuum line. Apply vacuum to air-dry the support.
- 9 Transfer the support to a 20-mL screw cap glass vial. Add a 10-mL of 10% Et₃N in 1:1 (v/v) CH₃CN/H₂O using a 10-mL glass syringe. Seal the vial and place it in a heat block for 30 min at 65 °C. Allow the flask to cool at ~ 25 °C.
- 10 Repeat successively steps 8 and 9 two more times.
- 11 Repeat step 8.
- 12 Add 5-mL DMSO to the Buchner funnel using a 5-mL glass syringe. Apply vacuum. Repeat step 12 four more times and air-dry the support **11a**.

Release of the DNA sequence 12a from the Support 11a

- 13 Transfer the air-dried solid support to a 20-mL screw cap glass vial. Add 5-mL of a 1.0 M tetra-*n*-butylammonium fluoride in dry DMSO using a 5-mL glass syringe. Seal the flask and place it in a heat block over a period of 3 hr at 65 °C. Allow the glass vial to cool to ~25 °C.
- 14 Open the glass vial. Add 2 mL methoxytrimethylsilane and 2 mL MeCN to the suspension using a 5-mL glass syringe. Seal the vial and keep it at ~25 °C for 30 min.
- 15 Pour the suspension into a 15-mL fritted glass vacuum filtration Buchner funnel of coarse porosity connected to a 100-mL round bottom flask. Allow the filtrate to drip into the flask. Connect the funnel to an in-house vacuum line. Apply vacuum to air-dry the support.
- 16 Transfer the air-dried solid support to a 20-mL screw cap glass vial. Add 10-mL of 10% Et₃N in 1:1 (v/v) MeCN/H₂O. Seal the flask and place it in a heat block over a period of 30 min at 65 °C. Allow the glass vial to cool to ~25 °C.
- 17 Repeat step 15.
- 18 Repeat successively steps 16 and 17. Remove the funnel from the 100-mL round-bottom flask containing all filtrates collected in step 15, 17 and 18. Concentrate the filtrates to ~ 1-mL by rotoevaporation under reduced pressure (~20 mmHg) using a vacuum pump.

Precipitate the DNA sequence 12a

- 19 Add 19-mL peroxide-free THF to the concentrated material. Seal the vial. Manually swirl the vial to yield a uniform precipitate of the solid-phase-purified DNA sequence **12a**. Allow the product to settle.
- 20 Open the vial and carefully pipet most of the clear supernatant using a pipettor. Add 10-mL peroxide-free THF to the vial using a 10-mL glass cylinder. Seal the vial and vortex it for 30 sec, using a vortexer. Allow the product to settle.
- 21 Repeat step 20 three times.
- 22 Open the vial. Pipet approximately half of the clear supernatant using a pipettor. Swirl the vial and pour the suspension into a 25-mL-round bottom flask. Rinse the vial by adding 5-mL peroxide-free THF using a 10-mL glass cylinder. Swirl and pour the wash into the 25-mL round-bottom flask.
- 23 Rotoevaporate the suspension to dryness under reduced pressure (~20 mmHg). Store indefinitely the purified nucleic acid sequence **12a** at -20 °C.

COMMENTARY

Background Information

The ability to design and synthesize DNA and RNA sequences has had a massive impact on the rapidly expanding fields of synthetic biology and nucleic acid-based drug development.

Typically, DNA/RNA sequences and their analogues have been used for recognition and binding to mRNAs encoding disease-causing proteins and have led to the production of nucleic acid-based drugs capable of inhibiting the expression of these proteins through either an antisense (Crook, 2004) or an RNA interference (Dorsett and Tuschl, 2004) pathway. Both approaches to silencing gene expression require the manufacture of synthetic nucleic acid sequences in large quantities (e.g., millimoles) and high purity for preclinical and clinical investigations. Conversely, total gene synthesis for synthetic biology applications requires small amounts (e.g., nanomoles) of a large number of highly pure synthetic DNA sequences. Although the automated chemical synthesis of nucleic acid sequences is efficient and can be scaled up for pharmaceutical production, the purification of these sequences is challenging; the full-length nucleic acid sequences are mixed with shorter sequences, due to incomplete coupling of activated deoxyribonucleoside phosphoramidites at each chain extension step of the nucleic acid sequence assembly. Even though these shorter than full-length sequences are produced in small amounts, the physicochemical similarity of these sequences to the desired nucleic acid sequence makes them very difficult to chromatographically separate them. HPLC-based methods including reversed-phase HPLC and anion-exchange HPLC are currently the preferred techniques for either small or large-scale purification of nucleic acid sequences. However, these methods require expensive instruments and accessories; considerable volumes of buffered aqueous and organic solvents are also needed for large-scale purification of DNA sequences, thereby contributing to high operational cost. In addition to not being cost-effective, HPLC-based purification methods are not amenable to parallel purification processes; only a single nucleic acid sequence can be purified per run and depending on the nature of each nucleic acid sequence, more than one purification run may be required to achieve the level of sequence purity required for pharmaceutical or synthetic biology applications. One important limitation of any large- or small-scale HPLC purification of DNA sequences is the burdensome removal of relatively large volumes of aqueous solvents generated during the purification process; this operation requires pricey equipment as well. HPLC-based purification methods are time-consuming and often do not completely resolve shorter than full-length sequences from the desired DNA sequence. Alternatively, polyacrylamide gel electrophoresis (PAGE, see UNIT 10.4) can efficiently separate shorter nucleic acid sequences from full-length DNA sequences; however the high-throughput or large scale purification of DNA sequences by PAGE, has to the best of our knowledge not been reported in the scientific literature. Several methods orthogonal to reversed-phase HPLC or electrophoretic techniques have, however, been proposed for the purification of nucleic acid sequences. These methods are based on ion-pair (Swiderski et al., 1994), hydrophobic (Sproat et al., 1999) or affinity (Fang and Bergstrom, 2003a, 2003b, 2003c, 2004; Pearson et al., 2005; Beller and Bannwarth, 2005; Dandapani, 2006; Mishra et al., 2006) chromatography. All of these purification techniques are either not amenable to either highly parallel or large-scale purification of nucleic acid sequences. A conceptually different approach to the purification of synthetic DNA sequences has been proposed (Fang et al., 2011; Pokharel et al., 2014). The method consists of catching shorter than full-length sequences through a polymerization process; the full-length sequence is then extracted from the polymer and isolated by precipitation in *n*-butanol. This process has been claimed to be potentially useful for large-scale and high-throughput purification of DNA sequences (Pokharel et al., 2014). In this context, an innovative solid-phase-based

The solid-phase synthesis of 5'-functionalized phosphorothioate and native DNA sequences (**10a-d** and **10e-f**, respectively) has been conducted using a commercial long-chain alkylamine controlled-pore glass support (LCAA-CPG) and standard deoxyribonucleoside phosphoramidites according to standard protocols. The 5'-functionalized phosphoramidites **9a-d** are used in the last coupling cycle of each DNA sequence assembly. It is critically important that phosphoramidite coupling efficiency and the capping of unreacted 5'-hydroxy functions be optimal for solid-phase purification of DNA sequences; less than optimal phosphoramidite coupling and capping reactions will result in poorer recovery of solid-phase-purified DNA sequences. In this regard, the coupling time of 5'-functionalized phosphoramidites **9a-d** has been extended to 180 s to ensure the highest coupling efficiency of these 5'-sterically-demanding phosphoramidite monomers. Post-synthesis deprotection and release of the DNA sequences **10a-f** from LCAA-CPG have been performed under basic conditions according to standard protocols (Ellington and Pollard, 2000). Each aqueous ammonia solution containing the crude DNA sequence **10a-f** (Figs. 4 and 5) and shorter than full-length sequences is evaporated to one-half of its original volume. Solid tetra-*n*-butylammonium chloride (TBACl) and the capture support **3** are then sequentially added to the aqueous solution of DNA sequences; the suspension is kept at 65 °C over a period of 3 hr. TBACl is used to exchange the ammonium counterions of thiophosphate or phosphate diester functions with tetrabutylammonium ions to enhance the solubility of phosphorothioate DNA sequences in aqueous solutions. RP-HPLC analysis of the pre- and post-capture solutions of each DNA sequence synthesized using the 5'-functionalized deoxyribonucleoside phosphoramidites **9a-d** indicates that the oximation reaction resulting in the capture of each 5'-functionalized DNA sequence **10a-f** is in all cases near complete (Grajkowski et al. 2016). The solid supports **11a-f** (Fig. 5) are then treated twice with a warm (55 °C) solution of aqueous ammonia in acetonitrile to wash off unbound shorter than full-length DNA sequences by filtration; a final wash with anhydrous DMSO is necessary to ensure efficient release of the DNA sequences **12a-f** from their respective supports (**11a-f**) upon treatment with 1.0 M tetra-*n*-butylammonium fluoride (TBAF) in dry DMSO for 3 hr at 65 °C. The solid-phase purified DNA sequences **12a-f** (Fig. 6) are isolated by precipitation in dry THF and characterized by ESI-TOF-MS. The purity of these sequences has been evaluated by RP-HPLC and by PAGE under denaturing conditions (Fig. 7, Fig. 8, Fig. 9, Fig. 10 and Fig. 11).

Although the RP-HPLC chromatograms show the purity of solid-phase-purified DNA sequences is essentially 100%, PAGE analyses of the same sequences show minute amounts of shorter than full-length DNA sequences, indicating that the coupling efficiency of commercial phosphoramidites and that of **9a-d** is less than 100%. The presence of tiny amounts of shorter than full-length sequences also underscores the affinity of those sequences for silica gel, which makes it difficult to quantitatively wash them off the supports **11a-f** prior to release.

With the objective of demonstrating that the solid-phase purification process can be qualified as highly parallel and scalable, the solid-phase synthesis of ten (10) identical phosphorothioate DNA sequences (**10a**) is carried out, each on a scale of 1 μmol using the conditions described in *Support Protocol 1* and *Alternate Protocol 1*. Upon completion of the

syntheses, deprotection, and release of each DNA sequence from the CPG support, the ammoniacal solution of each sequence are pooled together and concentrated under reduced pressure. The amounts of the capture support **3**, reagents, and solvent required for the capture of **10a** are increased by 10-fold while keeping the final concentration of the reagents the same as that reported for individual syntheses. The capture reaction is performed under conditions identical to those described in *Basic Protocol 4* for individual syntheses in terms of reaction time and temperature. As anticipated, the capture reaction resulted in the near complete (>98%) disappearance of the DNA sequence **10a** (Figs. 12A–B). Release of the DNA sequence **12a** from the support **11a** (Fig. 12C) has been done using a 10-fold increase of 1.0 M TBAF/DMSO solution that had been required for a 1 μ mol scale reaction, while keeping reaction time and temperature conditions the same. A 10-fold increase in the volume of THF is necessary to precipitate **12a** (Fig. 12D). The yield of the solid-phase-purified DNA sequence has been measured by UV spectroscopy at 260 nm and found to be 995 ODs; this yield is nearly proportional to that of solid-phase-purified **12a** (105 ODs) isolated from a 1 μ mol process scale. These results conclusively demonstrate that the solid-phase purification of nucleic acid sequences can be achieved in a highly parallel and scalable process. Thus, the proposed solid-phase purification process provides a cost-effective high-throughput alternative to conventional chromatographic methods for the purification and production of phosphorothioate or native DNA sequences with a comparable, if not better, level of purity.

Critical Parameters and Troubleshooting

When preparing the capture support **3**, it is critically important to adequately execute step 6 of *Basic Protocol 1* because any free amine functions, left unreacted on the solid-support **2** and/or **3**, has the potential of forming relatively unstable imino-conjugates with 5'-functionalized DNA sequences **10a–f** during the capture process. This can lead to premature hydrolysis of the imino-conjugates and ultimately result in poorer recovery of the solid-phase-purified DNA sequences **12a–f**.

The 5'-functionalized deoxyribonucleoside phosphoramidites **9a–d** (see Basic Protocol 3) must be prepared under strictly anhydrous conditions and must not be exposed to any acidic reagents or conditions, which would make these compounds prone to hydrolysis. With the purpose of preventing this unfavorable event, readers are referred to the Critical Parameters and Troubleshooting sections of *UNITS 2.7 & 3.17*, which address those critical issues associated with the preparation and use of nucleoside phosphoramidites. Precipitation of the 5'-functionalized deoxyribonucleoside phosphoramidites **9a–d** from cold hexane after silica-gel purification is helpful in rendering those phosphoramidites free of *H*-phosphonate impurities and residual triethylamine; if not completely removed, triethylamine can neutralize 1*H*-tetrazole or other acidic activators and reduce the coupling efficiency of any phosphoramidites during automated solid-phase DNA synthesis. The presence of moisture can also compromise the coupling efficiency of phosphoramidites and should be taken into consideration; this concern has also been addressed in *UNIT 2.7*.

The capture support **3** is acidic. It is therefore crucial to neutralize this support by exposing it to a solution of triethylamine in acetonitrile (see step 1 of *Basic Protocol 4*) prior to

capture of the 5'-functionalized DNA sequences **10a–f** to avoid depurination of these sequences.

The purity of solid-phase-purified DNA sequences **12a–f** is highly dependent on the efficiency of eluting off shorter than full-length DNA sequences from the silica-based support **11a–f**. The affinity of shorter than full-length DNA sequences for silica-based supports increases with the size and polarity of these sequences. Steps 5 through 9 of *Basic Protocol 4* have been optimized to wash off most of those shorter than full length DNA sequences. The same rationale applies when releasing the full-length DNA sequences (**12a–f**) from the supports **11a–f**; steps 12 through 15 of Basic Protocol 4 have also been optimized to isolate the DNA sequences with the highest recovery yields and purity (ca. 98+ %).

Anticipated Results

The two-step functionalization of commercial 3-aminopropylated silica gel with CDI-activated 7-oxooctanoic acid and *O,O'*-1,3-propanediylbishydroxylamine dihydrochloride should provide the capture support **3** with a surface density of ~ 150 μmol of aminoxy functions per gram of support. The reaction of *N,N'*-dimethylethylenediamine with 4,4-dimethyl- γ -butyrolactone affords the aminated amido alcohol **5b** in a yield of 74%; amidation of **5b** with CDI-activated 7-oxooctanoic acid gives the ketodiamidated alcohol **6b** in a yield of 80%. Conjugation of the linker **6b** to the 5'-hydroxyl of each deoxyribonucleoside **7a–d**, to provide 5'-functionalized deoxyribonucleosides (**8a–d**) with yields in the range of 50–70%. Phosphitylation of **8a–d** under standard conditions leads to the corresponding phosphoramidite derivatives **9a–d** in yields of 75–83%. The coupling efficiency of the 5'-functionalized phosphoramidites **9a–d**, used in the last solid-phase synthesis coupling cycle of each DNA sequences listed in Figure 4, should be ~99%. Capture of the 5'-functionalized DNA sequences **10a–f** by the solid support **3** is near complete, as shown in Figs. 12A–B. The high-throughput synthesis and 10-fold scale-up solid-phase purification process of the 5'-functionalized DNA sequence **10a** is performed under conditions highly similar to those used for the 1- μmole scale process; release of the captured **10a** from support **11a** yields 995 OD₂₆₀ of solid-phase purified DNA sequence **12a** after precipitation in dry THF (Figs. 12C–D). This yield should be nearly proportional to that measured from the solid-phase-purified **12a** (105 OD₂₆₀) isolated from a 1 μmol -scale purification process.

Time Considerations

Preparation of the intermediate solid support **2** will take ~26 hr to complete. The preparation of capture solid support **3** from support **2** can be achieved within 17 hr. Preparation and purification of each of the amidated alcohols **5a** and **5b** is accomplished within 24 hr. Synthesis and purification of the ketodiamidated alcohols **6a** and **6b** are each achieved within 35 hr. Synthesis and purification of the ketodiamidated alcohol **6c** needed for measuring the surface density of aminoxy functions conjugated to the capture support **3** will take 10 hr to complete. Measurement of the surface density of aminoxy functions per gram of support **3** is completed ~26 hr. Preparation and purification of each 5'-functionalized deoxyribonucleoside **8a–d** is achieved within 16 hr. Phosphitylation of **8a–d**

to the corresponding phosphoramidite derivatives **9a–d** and purification of the monomeric phosphoramidites including lyophilization can each be done within 25 hr. Solid phase synthesis and deprotection of each DNA sequence **10a–f** are completed within 24 hr. Solid-phase capture of each DNA sequence including washing off shorter than full-length DNA sequences from the support **11a–f** is done within 5 hr. Release of solid-phase purified DNA sequences **12a–f** from the support **11a–f** including precipitation, centrifugation and washes are each done with 6 hr. A 10-fold scale-up solid-phase purification of any DNA sequences can be achieved essentially within the time required for a small scale purification process. Many of the above tasks can be performed in a parallel manner.

Acknowledgments

This work was supported through FDA intramural funds.

Literature Cited

- Beaucage SL, Caruthers MH. Deoxynucleoside phosphoramidites -- A new class of key intermediates for deoxypolynucleotide synthesis. *Tetrahedron Lett.* 1981; 22:1859–1862.
- Beaucage, SL. Oligodeoxyribonucleotides synthesis -- Phosphoramidite Approach. In: Agrawal, S., editor. *Methods in Molecular Biology, Vol.20: Protocols for Oligonucleotides and Analogs.* Humana Press, Inc; Totowa: 1993. p. 33-61.
- Beller C, Bannwarth W. Noncovalent attachment of nucleotides by fluorine-fluorine interactions: Application to a simple purification principle for synthetic DNA fragments. *Helv Chim Acta.* 2005; 88:171–179.
- Cielak J, Grajkowski A, Ausin C, Gapeev A, Beaucage SL. Permanent or reversible conjugation of 2'-O- or 5'-O-aminooxymethylated nucleosides with functional groups as a convenient and efficient approach to the modification of RNA and DNA sequences. *Nucleic Acids Res.* 2012; 40:2312–2329. [PubMed: 22067450]
- Cielak J, Ausin C, Grajkowski A, Beaucage SL. The 2-cyano-2,2-dimethylethanamine-*N*-oxymethyl group for the 2'-hydroxyl protection of ribonucleosides in the solid-phase synthesis of RNA sequences. *Chem Eur J.* 2013; 19:4623–4632. [PubMed: 23417977]
- Cooke ST. Antisense strategies. *Curr Mol Med.* 2004; 4:465–487. [PubMed: 15267220]
- Dandapani S. Recent applications of fluorine separation methods in organic and bioorganic chemistry. *QSAR Comb Sci.* 2006; 25:681–688.
- Defrancq E, Lhomme J. Use of an aminoxy linker for the functionalization of oligodeoxyribonucleotides. *Bioorg Med Chem Lett.* 2001; 11:931–933. [PubMed: 11294394]
- Dorsett Y, Tuschl T. siRNAs: applications in functional genomics and potential as therapeutics. *Nat Rev Drug Discovery.* 2004; 3:318–329. [PubMed: 15060527]
- Eleuteri A, Capaldi DC, Krotz AH, Cole DL, Ravikumar VT. Pyridinium trifluoroacetate/*N*-methylimidazole as an efficient activator for oligonucleotide synthesis via the phosphoramidite method. *Org Process Res Develop.* 2000; 4:182–189.
- Ellington, A., Pollard, JD, Jr. Introduction to the synthesis and purification of oligonucleotides. In: Beaucage, SL, Bergstrom, DE, Glick, GD., Jones, RA., editors. *Current Protocols in Nucleic Acid Chemistry.* Vol. I. John Wiley & Sons, Inc; Hoboken: 2000. p. A.3C.1-A.3C.22.
- Fang S, Bergstrom DE. Reversible biotinylation phosphoramidite for 5'-end-labeling, phosphorylation, and affinity purification of synthetic oligonucleotides. *Bioconjug Chem.* 2003a; 14:80–85. [PubMed: 12526696]
- Fang S, Bergstrom DE. Fluoride cleavable biotinylation phosphoramidite for 5'-end-labeling and affinity purification of synthetic oligonucleotides. *Nucleic Acids Res.* 2003b; 31:708–715. [PubMed: 12527780]
- Fang, S., Bergstrom, DE. Reversible Biotinylation of the 5'-Terminus of Oligodeoxyribonucleotides and its Application in Affinity Purification. In: Beaucage, SL, Bergstrom, DE, Glick, GD., Jones,

- RA., editors. *Current Protocols in Nucleic Acid Chemistry*. Vol. I. John Wiley & Sons, Inc; Hoboken: 2003c. p. 4.20.1-4.20.17.
- Fang S, Bergstrom DE. Reversible 5'-end biotinylation and affinity purification of synthetic RNA. *Tetrahedron Lett.* 2004; 45:7987-7990.
- Fang S, Fueangfung S. Scalable synthetic oligodeoxynucleotide purification with use of a catching by polymerization, washing, and releasing approach. *Org Lett.* 2010; 12:3720-3723. [PubMed: 20704419]
- Fang S, Fueangfung S, Lin X, Zhang X, Mai W, Bi L, Green SA. Synthetic oligodeoxynucleotide purification by polymerization of failure sequences. *Chem Commun.* 2011; 47:1345-1347.
- Fina NJ, Edwards JO. The alpha effect. A review. *Int J Chem Kinet.* 1973; 5:1-26.
- Forget D, Boturyn D, Defrancq E, Lhomme J, Dumy P. Highly efficient synthesis of peptide-oligonucleotide conjugates: Chemoselective oxime and thiazolidine formation. *Chem Eur J.* 2001a; 7:3976-3984. [PubMed: 11596939]
- Forget D, Renaudet O, Defrancq E, Dumy P. Efficient preparation of carbohydrate-oligonucleotide conjugates (COCs) using oxime bond formation. *Tetrahedron Lett.* 2001b; 42:7829-7832.
- Forget D, Renaudet O, Boturyn D, Defrancq E, Dumy P. 3'-Oligonucleotides conjugation via chemoselective oxime bond formation. *Tetrahedron Lett.* 2001c; 42:9171-9174.
- Grajkowski A, Cieslak J, Beaucage SL. Solid-phase purification of synthetic DNA sequences. *J Org Chem.* 2016; 8:6165-6175.
- Iyer RP, Phillips LR, Egan W, Regan JB, Beaucage SL. The automated synthesis of sulfur-containing oligodeoxyribonucleotides using 3*H*-1,2-benzodithiol-3-one 1,1-dioxide as a sulfur-transfer reagent. *J Org Chem.* 1990; 55:4693-4698.
- Katajisto J, Virta P, Lönnberg H. Solid-phase synthesis of multiantennary oligonucleotide glycoconjugates utilizing on-support oximation. *Bioconjug Chem.* 2004; 15:890-896. [PubMed: 15264878]
- Kawasaki AM, Casper MD, Prakash TP, Manalili S, Sasmor H, Manoharan M, Cook PD. Synthesis, hybridization, and nuclease resistance properties of 2'-*O*-aminoxyethyl (2'-*O*-AOE) modified oligonucleotides. *Tetrahedron Lett.* 1999; 40:661-664.
- McBride LJ, McCollum C, Davidson S, Efcavitch JW, Andrus A, Lombardi SJ. A new, reliable cartridge for the rapid purification of synthetic DNA. *BioTechniques.* 1988; 6:362-367. [PubMed: 3273859]
- Mishra R, Mishra S, Misra K. Synthesis and application of fluorouracil-tagged oligonucleotides. *Chem Lett.* 2006; (10):1184-1185.
- Morvan F, Sanghvi YS, Perbost M, Vasseur JJ, Bellon L. Oligonucleotide mimics for antisense therapeutics: Solution phase and automated solid-support synthesis of MMI linked oligomers. *J Am Chem Soc.* 1996; 118:255-256.
- Pearson WH, Berry DA, Stoy P, Jung K-Y, Sercel AD. Fluorouracil affinity purification of oligonucleotides. *J Org Chem.* 2005; 70:7114-7122. [PubMed: 16122230]
- Pokharel D, Yuan Y, Fueangfung S, Fang S. Synthetic oligodeoxynucleotide purification by capping failure sequences with a methacrylamide phosphoramidite followed by polymerization. *RSC Adv.* 2014; 4:8746-8757.
- Salo H, Virta P, Hakala H, Prakash TP, Kawasaki AM, Manoharan M, Lönnberg H. Aminoxy functionalized oligonucleotides: Preparation, on-support derivatization, and postsynthetic attachment to polymer support. *Bioconjug Chem.* 1999; 10:815-823. [PubMed: 10502348]
- Sinha ND, Biernat J, McManus J, Köster H. Polymer support oligonucleotide synthesis XVIII: Use of β -cyanoethyl-*N,N*-dialkylamino-*N*-morpholino phosphoramidite of deoxynucleosides for the synthesis of DNA fragments simplifying deprotection and isolation of the final product. *Nucl Acids Res.* 1984; 12:4539-4557. [PubMed: 6547529]
- Sproat BS, Rupp T, Menhardt N, Keane D, Beijer B. Fast and simple purification of chemically modified hammerhead ribozymes using a lipophilic capture tag. *Nucleic Acids Res.* 1999; 27:1950-1955. [PubMed: 10101206]
- Swiderski PM, Bertrand EL, Kaplan BE. Polystyrene reverse-phase ion-pair chromatography of chimeric ribosymes. *Anal Biochem.* 1994; 216:83-88. [PubMed: 8135369]

- Trader DJ, Carlson EE. Siloxyl ether functionalized resins for chemoselective enrichment of carboxylic acids. *Org Lett.* 2011; 13:5652–5655. [PubMed: 21957886]
- Trevisiol E, Renard A, Defrancq E, Lhomme J. The oxyamino-aldehyde coupling reaction: An efficient method for the derivatization of oligonucleotides. *Tetrahedron Lett.* 1997; 38:8687–8690.
- Urdea MS, Horn T. Solid-supported synthesis, deprotection and enzymatic purification of oligodeoxyribonucleotides. *Tetrahedron Lett.* 1986; 27:2933–2936.
- Wei X. Coupling activators for the oligonucleotide synthesis via phosphoramidite approach. *Tetrahedron.* 2013; 69:3615–3637.
- Wilk, A., Grajkowski, A., Chmielewski, MK., Phillips, LR., Beaucage, SL. Deoxyribonucleoside Phosphoramidites. In: Beaucage, SL, Bergstrom, DE, Glick, GD., Jones, RA., editors. *Current Protocols in Nucleic Acid Chemistry*. Vol. I. John Wiley & Sons, Inc; Hoboken: 2001. p. 2.7.1-2.7.12.





Figure 1.
Preparation of the capture support **3**. CDI, 1,1'-carbonyldiimidazole; Si, 3-aminopropyl silica gel.

**Figure 2.**

Preparation of the diamidated keto alcohol linkers **6a–c**. Conditions: (i) MMTr-Cl, dry pyridine, 25 °C, 4 hr. CDI, 1,1'-carbonyldiimidazole; MMTr, 4-monomethoxytrityl.

Oligomer	DNA Sequence ^a
10a	5'-CLF-d(A _{PS} C _{PS} A _{PS} C _{PS} T _{PS} G _{PS} T _{PS} G _{PS} A _{PS} A _{PS} T _{PS} C _{PS} G _{PS} A _{PS} T _{PS} G _{PS} C _{PS} C _{PS} A _{PS} T)
10b	5'-CLF-d(C _{PS} T _{PS} C _{PS} C _{PS} G _{PS} T _{PS} A _{PS} C _{PS} C _{PS} T _{PS} T _{PS} A _{PS} C _{PS} G _{PS} T _{PS} C _{PS} T _{PS} T _{PS} G _{PS} T)
10c	5'-CLF-d(G _{PS} T _{PS} G _{PS} A _{PS} G _{PS} T _{PS} A _{PS} G _{PS} C _{PS} G _{PS} A _{PS} A _{PS} C _{PS} G _{PS} T _{PS} G _{PS} A _{PS} A _{PS} G _{PS} T)
10d	5'-CLF-d(T _{PS} A _{PS} T _{PS} C _{PS} C _{PS} G _{PS} T _{PS} A _{PS} G _{PS} C _{PS} T _{PS} A _{PS} A _{PS} C _{PS} G _{PS} T _{PS} C _{PS} A _{PS} G _{PS} T)
10e	5'-CLF-d(A _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P C _P A _P T)
10f	5'-CLF- d(T _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P A _P A _P T _P T _P G _P C _P A _P C _P T _P G _P T _P G _P - A _P A _P T _P C _P G _P A _P T _P G _P C _P C _P A _P T _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P C _P A _P T)

^a CLF, capture linker function; PS, phosphorothioate diester; P, phosphate diester

Figure 4.

5'-Functionalized DNA sequences for solid-phase capture and purification.

Oligomer	DNA Sequence ^a
12a	5'-d(A _{PS} C _{PS} A _{PS} C _{PS} T _{PS} G _{PS} T _{PS} G _{PS} A _{PS} A _{PS} T _{PS} C _{PS} G _{PS} A _{PS} T _{PS} G _{PS} C _{PS} C _{PS} A _{PS} T)
12b	5'-d(C _{PS} T _{PS} C _{PS} C _{PS} G _{PS} T _{PS} A _{PS} C _{PS} C _{PS} T _{PS} T _{PS} A _{PS} C _{PS} G _{PS} T _{PS} C _{PS} T _{PS} T _{PS} G _{PS} T)
12c	5'-d(G _{PS} T _{PS} G _{PS} A _{PS} G _{PS} T _{PS} A _{PS} G _{PS} C _{PS} G _{PS} A _{PS} A _{PS} C _{PS} G _{PS} T _{PS} G _{PS} A _{PS} A _{PS} G _{PS} T)
12d	5'-d(T _{PS} A _{PS} T _{PS} C _{PS} C _{PS} G _{PS} T _{PS} A _{PS} G _{PS} C _{PS} T _{PS} A _{PS} A _{PS} C _{PS} G _{PS} T _{PS} C _{PS} A _{PS} G _{PS} T)
12e	5'-d(A _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P A _P T)
12f	5'-d(T _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P A _P A _P T _P G _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P A _P T _P G _P T _P G _P - A _P A _P T _P C _P G _P A _P T _P G _P C _P A _P T _P C _P A _P C _P T _P G _P T _P G _P A _P A _P T _P C _P G _P A _P T _P G _P C _P A _P T)

^a PS, phosphorothioate diester; P, phosphate diester; A, 2'-deoxyadenosine; C, 2'-deoxycytidine; G, 2'-deoxyguanosine; T, 2'-deoxythymidine

Figure 6.
Solid-phase purified DNA sequences released from the supports **11a–f**.




Figure 7. Solid-phase purification of the phosphorothioate DNA sequence **10a**. RP-HPLC analysis of solid-phase-purified **12a** that was released from the support **11a**. Inset: Purity analysis of the solid-phase-purified **12a** by PAGE. Chromatographic conditions are reported in the annotation of step 3 of *Basic Protocol 4*; electrophoretic conditions are reported in Grajkowski et al. 2016. Abbreviation: BP, bromophenol blue dye.




Figure 8. Solid-phase purification of the phosphorothioate DNA sequence **10b**. RP-HPLC analysis of solid-phase-purified **12b** that was released from the support **11b**. Inset: Purity analysis of the solid-phase-purified **12b** by PAGE. Abbreviation: BP, bromophenol blue dye.




Figure 9. Solid-phase purification of the phosphorothioate DNA sequence **10c**. RP-HPLC analysis of solid-phase-purified **12c** that was released from the support **11c**. Inset: Purity analysis of the solid-phase-purified **12c** by PAGE. Abbreviation: BP, bromophenol blue dye.




Figure 10. Solid-phase purification of the phosphorothioate DNA sequence **10d**. RP-HPLC analysis of solid-phase-purified **12d** that was released from the support **11d**. Inset: Purity analysis of the solid-phase-purified **12d** by PAGE. Abbreviation: BP, bromophenol blue dye.







Figure 11. Solid-phase purification of the native DNA sequence **10f**. RP-HPLC analysis of solid-phase-purified **12f** that was released from the support **11f**. Inset: Purity analysis of the solid-phase-purified **12f** by PAGE. Abbreviation: BP, bromophenol blue dye.





D



Figure 12.

10-Fold scale up solid-phase purification of the phosphorothioate DNA sequence **12a**. (A) RP-HPLC profile of unpurified **10a**. (B) RP-HPLC profile of unpurified **10a** after capture by the support **3**. (C) RP-HPLC analysis of solid-phase purified **12a** that has been released from the support **11a**. (D) Photograph of precipitated **12a** accounting for 995 OD₂₆₀ of the DNA sequence. Inset: Purity analysis of the solid-phase purified **12a** by PAGE. Abbreviation: BP, bromophenol blue dye.

Author Manuscript

Author Manuscript

Author Manuscript

Author Manuscript