

RESEARCH ARTICLE

A validated LC-MS/MS method for cellular thyroid hormone metabolism: Uptake and turnover of mono-iodinated thyroid hormone metabolites by PCCL3 thyrocytes

Keith H. Richards¹, Nancy Schanze¹, Ray Monk¹, Eddy Rijntjes¹, Daniel Rathmann², Josef Köhrle^{1*}

1 Charité-Universitätsmedizin Berlin, Institut für Experimentelle Endokrinologie, Augustenburger Platz 1, Berlin, Germany, **2** Chromicent GmbH, Johann-Hittorf-Strasse 8, Berlin, Germany

* josef.koehrle@charite.de



OPEN ACCESS

Citation: Richards KH, Schanze N, Monk R, Rijntjes E, Rathmann D, Köhrle J (2017) A validated LC-MS/MS method for cellular thyroid hormone metabolism: Uptake and turnover of mono-iodinated thyroid hormone metabolites by PCCL3 thyrocytes. PLoS ONE 12(8): e0183482. <https://doi.org/10.1371/journal.pone.0183482>

Editor: Michelina Plateroti, University Claude Bernard Lyon 1, FRANCE

Received: January 27, 2017

Accepted: August 5, 2017

Published: August 24, 2017

Copyright: © 2017 Richards et al. This is an open access article distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

Data Availability Statement: All relevant data are within the paper and its Supporting Information files.

Funding: This work was supported by the Deutsche Forschungsgemeinschaft DFG (KO922-16/1-2 to J.K. and KO922-17/1-2 to J.K.) in the framework of the DFG funded Priority Programme SPP1629 THYROIDTRANSACT and a Charité thesis scholarship (to N.S.). The company Chromicent GmbH played no part in any capacity in the project

Abstract

Tyrosine and phenolic ring de-iodination of thyroid hormones (TH) is crucial for regulating their physiological activity. Furthermore, reactions such as de-carboxylation to thyronamines (TAM) and de-amination to thyroacetic acids (TAc) produce TH metabolites (THM) with distinct biological properties. This needs to be considered when studying effects of TH and THM. The accurate and precise quantitative analysis of TH and THM in cell culture supernatants and cell lysates are key procedures required for studying the *in vitro* metabolism of TH. We report here the development of a liquid-liquid extraction/isotope dilution-liquid chromatography-electrospray tandem mass spectrometry (LC-MS/MS) method for the quantification of 9 thyronines (TN) and 6 TAM in human hepatocellular carcinoma Hep G2 cell lysate extracts. In addition, we adapted the method to quantify TH, TAM and TAc, in cell lysates of FBS-depleted rat thyroid epithelium PCCL3 cells. The methods for both cell lines were validated by rigorous assessment of linearity, limits of quantification and detection (LLOQ and LLOD respectively), intra- and inter-day accuracy, precision, process efficiency (PE), matrix effect (ME) and relative recovery (RE). Calibration curves covering 11 concentrations (based on 400 µl of lysate) were linear in the range 0.016–50 nM and 0.010–50 nM for Hep G2 and PCCL3 cells respectively. The lower limits of quantification were in the range 0.031 to 1 nM. We applied the PCCL3 version of the LC-MS/MS method to the analysis of lysed cell extracts from PCCL3 cells that had been incubated with 3-iodo-L-thyronine (T₁), 3-iodothyronamine (3-T₁AM) and 3-iodothyroacetic acid (3-T₁Ac). Over the course of 30 minutes incubation 3-T₁AM was de-iodinated to 4-[4-(2-aminoethylphenoxy)]phenol (thyronamine, T₀AM) and de-aminated to 3-T₁Ac respectively, whilst T₁ underwent de-iodination to T₀. This data indicates avid metabolism of these mono-iodinated compounds and the utility of LC-MS/MS to quantify such cellular metabolism.


funding, intellectual input or experimentation, nor in the decision to publish/preparation of the manuscript. Daniel Rathmann's role (who now works for Chromicent GmbH) was carried out whilst still an employee at Charité Universitätsmedizin, Berlin. The DFG solely provided financial support in the form of authors' salaries and/or research materials, they had no other involvement.

Competing interests: We have no commercial affiliation to Chromicent GmbH including consultancy, patents, products in development, or marketed products, etc. Dr Rathmann's position as an employee of Chromicent GmbH does not alter our adherence to PLOS ONE policies on sharing data and materials. (as detailed online in the guide for authors <http://journals.plos.org/plosone/s/competing-interests>).

Introduction

The thyroid hormone (TH) 3,5,3'-triiodo-L-thyronine (T_3) regulates a variety of processes that ensure proper development, growth and metabolism. Most of the circulating T_3 is generated by de-iodination of the phenolic ring of the less active TH 3,5,3',5'-tetraiodo-L-thyronine (T_4)—a reaction catalysed by deiodinases 1 and 2 *in vivo* [1, 2]. Inactivation of T_4 is also accomplished by de-iodination, and leads to the formation of 3,3',5'-triiodothyronine (reverse T_3 , rT_3); similarly, de-iodination of T_3 generates either the active 3,5-diodothyronine (3,5- T_2) or the inactive 3,3'-diodothyronine 3,3'- T_2), [3]. Besides de-iodination reactions, several other pathways of TH metabolism are possible. TH metabolites (THM, see Fig 1.) include thyronamines (TAM), resulting from TH de-carboxylation, and thyroacetic acids (TAc) resulting from the deamination of TAM. Some of these THM are endogenous and possess biological activity [4–7]. For example, 3,5-diiiodothyronine (3, 5- T_2) exerts thyromimetic action in rodents [8, 9] and treatment with 3-iodothyronamine (3- T_1 AM) or 4-[4-(2-aminoethylphenoxy)]phenol (T_0 AM) produces partially TH antagonistic effects such as hypothermia in mice and Djungarian hamsters [10, 11]. The mechanisms of action of TH and THM in cell culture systems *in vitro* are of high scientific interest; however, uptake, release and intracellular metabolism affect their bioavailability or may lead to the formation of products with their own distinct biological activity in the experimental system. To elucidate how TH and THM are utilised by cell types derived from different tissues can help clarify their mode(s) of action. Hence, the development of a validated and convenient analytical method for TH, TAM and TAc in cell extracts *in vitro* is of major importance.

We recently published a validated analytical method based on liquid-liquid extraction and isotope dilution high performance liquid chromatography/tandem mass spectrometry (LC-MS/MS) for the determination of 15 TH/THM (T_0 to T_4 thyronines (TN) and TAM, see Fig 1 for a complete list of compounds) in cell culture media extracts [12]. The method demonstrated the accurate, reproducible quantification of THM within a single 10 min analysis, with



L-Thyroxine (T_4): $R^1, R^2, R^3, R^4 = I$	3,5,3',5'-Tetraiodothyronamine (T_4 AM): $R^1, R^2, R^3, R^4 = I$	3-Iodothyroacetic acid (T_1 Ac)
3,5,3'-Triiodo-L-Thyronine (T_3): $R^1 = H, R^2, R^3, R^4 = I$	3,5,3'-Triiodothyronamine (T_3 AM): $R^1 = H, R^2, R^3, R^4 = I$	Thyroacetic acid (T_0 Ac): $R^1 =$
3,3',5'-L-Thyronine (rT_3): $R^2 = H, R^1, R^3, R^4 = I$	3,3',5'-Tetraiodothyronamine (rT_3 AM): $R^2 = H, R^1, R^3, R^4 = I$	
3,5-Diiodo-L-thyronine (3,5- T_2): $R^2, R^3 = I, R^1, R^4 = H$	3,5-Diiodothyronamine (3,5- T_2 AM): $R^2, R^3 = I, R^1, R^4 = H$	
3,3'-Diiodo-L-thyronine (3,3'- T_2): $R^1, R^2 = H, R^3, R^4 = I$	3-iodothyronamine (3- T_1 AM): $R^1, R^2, R^4 = H, R^3 = I$	
3',5'-Diiodo-L-thyronine (3',5'- T_2): $R^2, R^3 = H, R^1, R^4 = I$	Thyronamine (T_0 AM): $R^1, R^2, R^3, R^4 = H$	
3-Iodo-L-thyronine (3- T_1): $R^1, R^2, R^4 = H, R^3 = I$		
3'-Iodo-L-thyronine (3'- T_1): $R^1, R^2, R^3 = H, R^4 = I$		
L-Thyronine (T_0): $R^1, R^2, R^3, R^4 = H$		

Fig 1. Molecular structures of TH, TAM and TAc.

<https://doi.org/10.1371/journal.pone.0183482.g001>

lower limits of quantification (LLOQ) in the range 0.078–0.234 nM. We applied the method to quantify the de-iodination metabolites 3,3'-T₂, 3-T₁ and T₀ that were detected in DMEM medium when T₃ was incubated with primary hepatocytes [12]. We recently reported a preliminary adaption of the above method to analyse a limited number of TN (T₄, T₃ and rT₃) in Madin-Darby canine kidney 1 cell lysate extracts as part of a study on molecular characterization of monocarboxylate transporters involved in cellular TH uptake and efflux [13].

We now describe the extension of this method to enable the analysis and quantification of TH, THM and TAc in extracted, lysed cells. The method has been validated for the 15 TH/THM, using eight inter day replicates, for the human hepatocellular carcinoma cell line Hep G2, applying US Federal Drug Administration guidelines [14]. In addition, the method was validated (with 3-inter day repeats) for rat thyroid epithelium PCCL3 cells [15, 16] with the inclusion of 3-iodothyroacetic acid (3-T₁Ac) and thyroacetic acid (T₀Ac). We applied the latter method to the elucidation of potential metabolism pathways in cultured rat thyroid PCCL3 cells by incubating with 3-T₁, 3-T₁AM and 3-T₁Ac.

Materials and methods

Chemicals, reagents, buffers, calibrators and quality control samples

The internal standards (IS) ¹³C₆-T₄, ¹³C₆-T₃, ¹³C₆-rT₃, and ¹³C₆-3,3'-T₂ were obtained from Isosciences LLC (King of Prussia, PA, USA). Further IS, ²H₄-3-T₁AM, ²H₄-3-T₁Ac and ²H₄-T₀Ac, were a generous gift from Dr. T.S. Scanlan (Portland, OR, USA).

T₄AM, T₃AM and rT₃AM, were provided by Dr. R. Thoma (Formula GmbH, Berlin, D), who also purified 3-T₁AM by preparative HPLC. Dr. R. Smits (ABX Advanced Biochemical Compounds, Radeberg, Germany) performed synthesis of high-purity 3-T₁AM and T₀AM. 3,5-T₂AM, T₁Ac and T₀Ac were also kindly provided by T.S. Scanlan; T₄, T₃, 3,3'-T₂, 3,5-T₂, 3',5'-T₂, 3-T₁, 3'-T₁ and T₀ were sourced from Henning Arzneimittel GmbH (Berlin, Germany). ¹³C- and ¹H-labelled internal standards contained no measurable concentrations of unlabelled analogues and no cross-contamination from other ¹³C- or ¹H-labelled TH. TN, TAM and TAc were of minimum 99.6% purity. No corrections were made for potential cross-contaminating compounds.

Certified reference standards of T₃, 3,3'-T₂ and T₄ at 100 µg/mL in 0.1 N ammonium hydroxide/methanol were purchased from Qmx Laboratories (Thaxted, UK).

CHROMASOLV® Plus 2-Propanol, *t*-butyl methyl ether, methanol and dimethyl sulfoxide (DMSO) were from Sigma-Aldrich Chemie GmbH (Munich, Germany) as was formic acid, hydrochloric acid (37%), glacial acetic acid, sodium hydroxide and EDTA. Sucrose was purchased from Carl Roth (Karlsruhe, Germany) and Hepes from Merck (Darmstadt, Germany). Water was purified by a Milli-Q water purification system from Millipore (Billerica, MA, USA) and was of 18.2 mΩ quality.

Coon's F12 medium and fetal bovine serum (FBS) were purchased from Biochrom (Berlin, Germany); bovine TSH, insulin, transferrin, hydrocortisone, glycyl-histidyl-lysine somatostatin and DMEM (Dulbecco's Modified Eagle's medium) were obtained from Sigma Aldrich (Munich, Germany).

Individual stock solutions of TN, TAM and TAc were prepared at 10 mM in DMSO. A 20 µM mixed TN/TAM (15 compounds) working solution of T₄, T₃, rT₃, 3,3'-T₂, 3,5-T₂, 3',5'-T₂, 3-T₁, 3'-T₁, T₀, T₄AM, T₃AM, rT₃AM, 3,5-T₂AM, 3-T₁AM and T₀AM was prepared in 50% v/v methanol/water containing 0.1% formic acid (re-suspension buffer, RS). For PCCL3 analyses an additional mixed TN/TAM/TAc (17 compounds) working solution, also at 20 µM, was prepared containing 3-T₁Ac and T₀Ac in addition to the above THM. A working solution of certified reference standards containing a mixture of T₄, T₃ and 3,3'-T₂ at 20 µM

was prepared in RS. Further working solutions at 1000, 100 and 10 nM were prepared in RS and used to spike cell lysates to obtain calibrator and quality control samples. For Hep G2 cell lysates in 24-well cell culture plates, calibrators were spiked at 0.016, 0.031, 0.063, 0.125, 0.25, 0.5, 1.25, 2.5, 5, 20 and 50 nM. In-house and external spiked quality-control samples (QQAL) at low (LQAL, 0.10 nM), medium (MQAL, 12.5 nM) and high (HQAL, 40 nM) concentrations were prepared by spiking cell lysates with either the appropriate in-house or certified reference working standard solutions respectively. For PCCL3 cell lysates in 6-well cell culture plates, calibrators were spiked at 0.01, 0.05, 0.1, 0.5, 1.0, 2.5, 12.5, 25, 50, 125, 250 and 500 nM. In-house QQALs at low (LQAL, 0.10 nM), medium (MQAL, 1.0 nM) and high (HQAL, 2.5 nM) concentrations were prepared by spiking cell lysates with in-house working standard solutions.

IS stock solutions for Hep G2 analyses were prepared in DMSO and contained $^{13}\text{C}_6\text{-T}_4$, $^{13}\text{C}_6\text{-T}_3$, $^{13}\text{C}_6\text{-rT}_3$, $^{13}\text{C}_6\text{-3,3'-T}_2$ and $^2\text{H}_4\text{-3-T}_1\text{AM}$ (all at 2 μM). For PCCL3 analyses an additional IS stock solution was prepared containing the same 5 stable-isotopically labelled standards plus $^2\text{H}_4\text{-T}_0\text{Ac}$ and $^2\text{H}_4\text{-T}_1\text{Ac}$ (also at 2 μM in DMSO). All stock solutions were stored at -20°C in brown glass vials with Teflon-lined screw tops.

Cell culture

All cell culture work was performed under standard culturing conditions (5% CO_2 and 37°C) with standard aseptic techniques used throughout. Hep G2 medium was used without phenol red indicator.

Hep G2 cells were cultured in DMEM culture medium supplemented with 10% FBS. Cells were seeded in 24-well plates and grown confluent. At 24 h post plating, the medium was changed. After a further 24 h the medium was again changed, but this time for culture media without FBS. This state of cellular starvation was maintained for another 24 h, after which plates destined for standard curve production were then taken to ice, the media aspirated and the wells washed twice with ice cold PBS. The plates were stored dry at -20°C .

Rat thyroid PCCL3 cells were cultured in Coon's F12 medium supplemented with 5% FBS, 1 mU / ml bovine TSH and 5 hormones (10 μg / ml insulin, 5 μg / ml transferrin, 10 nM hydrocortisone, 20 ng / ml glycyl-histidyl-lysine and 10 ng / ml somatostatin). After seeding in 6-well plates at a density of 400,000 cells per well, the cells were grown for 24 h in medium as described above. 24 h before experiments the cells were washed twice with PBS and grown in medium without FBS to deplete them from TH and their metabolites. 24 h fasted PCCL3 cells that were to be used to investigate TH metabolism were subjected to incubations with TH (or just vehicle) in Krebs Ringer Buffer for various time points before being washed and stored in the same manner as those destined for use in making standard curves.

Cell lysis and sample extraction

Homogenisation Buffer pH 7.4 (250 mM D-(+)Sucrose, 20 mM Hepes, 1 mM EDTA; HB) was made and stored in aliquots at -20°C for single use. HB was mixed 1:1 with 0.1 N NaOH to produce Lysis Buffer 1 (LB1). Lysis Buffer 2 (LB2) consisted of 30% v/v glacial acetic acid in HB. Both LB1 and LB2 were made freshly for each cell lysis.

The protocol stated was applied to both 6- (PCCL3) and 24- (Hep G2) well cell culture plates. Cells were lysed at room temperature by adding 200 μl LB1 per well with orbital shaking (600 rpm) for 3 min followed by the addition of 100 μl LB2. Calibrator and quality control designated wells were spiked with 95 μl of appropriate TN/TAM or TN/TAM/TAc working solution; those wells designated as unspiked or matrix controls received 95 μl RS. All wells except those designated as matrix controls were then spiked with 5 μl of the appropriate IS mix.

The resulting 400 µl cell lysates were transferred to 2.0 ml Eppendorf tubes and incubated in the dark at 37°C for 1 h. 1 ml of 30% v/v 2-propanol /*t*-butyl methyl ether was added to the lysate and the mixture was vortexed for 5 min. The phases were separated via centrifugation at 2000 rcf for 5 min, at room temperature. The upper (non-aqueous) phase was transferred to a new 2.0 ml Eppendorf tube, and the extraction process was repeated with a further 1 ml of 2-propanol /*t*-butyl methyl ether. The combined organic phases were dried to a pellet in an Eppendorf 5301 Concentrator.

Samples, calibrators and quality controls were re-suspended in 100 µl of RS. Samples designated as matrix controls were re-suspended in 95 µl of TN/TAM or TN/TAM/TAc working solution at the appropriate concentration, and then spiked with 5 µl of the appropriate IS Mix. All samples were stored at -20°C for LC-MS/MS analysis.

Incubation of PCCL3 cells with mono-iodinated TH, TAM and TAc

After 24 h under FBS depletion, PCCL3 cells were incubated +/- 500 nM 3-T₁AM, 3-T₁, 3-T₁Ac or DMSO in Krebs Ringer buffer (119 mM NaCl, 4.74 mM KCl, 1.19 mM KH₂PO₄, 25 mM NaHCO₃, 2.54 mM CaCl₂, 1.19 mM MgCl₂, 10mM Hepes, 11 mM glucose, pH 7.4, 1.5mL) for 5, 10 or 30 min. After removal of the incubation buffer, cells were washed 3 times with PBS and frozen at -20°C to await cell lysis and extraction as described in detail in the previous section. Supernatants were also stored at -20°C; 400 µL aliquots were extracted as previously described [12] and analysed with the LC-MS/MS method described in the next section.

HPLC–MS/MS analysis

Sample analyses were performed by auto-injection of samples with a PAL HTC-xt auto-sampler (CTC Analytics AG, Zwingen, Switzerland) onto an HSS PFP 2.5 µm 3.0x100 mm column (Waters, Milford, MA, USA) maintained at 40°C and served by a 1260 binary HPLC system (Agilent technologies GmbH, Waldbronn, Germany). The HPLC column outlet was linked to a Sciex API 6500 QTRAP mass spectrometer (SCIEX Germany GmbH, Darmstadt, Germany) fitted with a Turbo Spray IonDrive electrospray source. The positive ion electrospray LC-MS/MS method used for the analysis of TN, TAM and accompanying IS has been described in detail [12] and was applied without modification to the analysis of extracted Hep G2 cell lysates.

For the analysis of extracted PCCL3 cell lysates the electrospray method was extended to facilitate positive/negative ion switching. The source temperature was reduced to 450°C and the ionspray voltage was -4500 V. The additional MRM transitions listed in Table 1 were incorporated into the negative ion experiment. Dwell times of 10 ms were used throughout. The complete duty cycle (positive plus negative ion) amounted to 0.73 sec.

Table 1. Negative ion electrospray tandem mass spectrometry and retention time (R_t) parameters for thyroacetic acids and internal standards.

Compound	R _t (min)	PI (m/z)	Pro 1 (m/z)	Pro 2 (m/z)	DP (V)	CE (V)	CXP (V)
² H ₄ -T ₀ Ac	4.71	247.0	203.1	105.9	-15	-14 (-26)	-9 (-5)
² H ₄ -T ₁ Ac	5.15	372.9	329.0	126.8	-5	-8 (-16)	-15 (-13)
T ₀ Ac	4.72	243.0	198.9	106	-10	-12 (-24)	-9 (-5)
3-T ₁ Ac	5.16	369.0	324.9	126.7	-15	-8 (-26)	-13 (-19)

PI: precursor ion; Pro 1: quantifier product ion; Pro 2: confirmation product ion; DP: declustering potential; CE: collision energy; m/z: mass to charge ratio; CXP: collision exit potential; V: volts. Values in parentheses refer to Pro 2.

<https://doi.org/10.1371/journal.pone.0183482.t001>

For both Hep G2 and PCCL3 experiments the relative recoveries (RE), process efficiencies (PE) and matrix effects (ME) were determined. The ME represents the impact of interfering compounds (from the cell lysates and buffers) on analyte ionisation [17] and was calculated as the percentage of [post extraction spiked MRM peak area]/[spiked solvent MRM peak area] [18]. The PE, representing the amount of substance recovered after undergoing the extraction process, is defined as the percentage of the [pre extraction spiked MRM peak area]/[post extraction spiked MRM peak area]. A poor PE thus indicated the analyte(s) of interest were lost during the extraction process [19]. The RE describes the percentage of the [pre extraction spiked MRM peak area]/[spiked solvent MRM peak area]. Lower limits of quantification (LLOQ) and detection (LLOD) were determined as described in reference [12] for calibrator sample chromatographic peaks passing precision and accuracy tolerances $\leq 20\%$ and showing signal:noise $> 10:1$ and $> 6:1$ respectively.

Results

Chromatographic peak measurement and identification

Representative chromatograms for TAM and TN in Hep G2 and for TAM, TN and TAc in PCCL3 cell lysate extracts are displayed in Figs 2 and 3 respectively. For both Hep G2 and




Fig 2. T₁AM, T₀AM and T₁Ac chromatograms from a PCCL3 cell lysate extract after 5 min incubation with T₁AM. 2 MRM Traces per analyte, overlaid. Upper trace: T₁AM, middle trace: T₀AM (4.56 min), lower trace: T₁Ac.

<https://doi.org/10.1371/journal.pone.0183482.g002>

PCCL3 methods the lowest chromatographic peak width was 12 sec. (obtained for 3,3'-T₂); the MS duty cycle allowed for a minimum of 15 data points to describe each peak to ensure reproducible peak integration. For each analyte and internal standard two MRM transitions were monitored; in most cases (except for 3-T₁AM where background peaks interfered at the LLOQ) the most intense peak was used as quantification ion and the less intense peak as confirmation ion. To confirm the identity of metabolites arising from the incubation of PCCL3 cells with TH, the observed quantification ion/confirmation ion MRM peak area ratios (MRM_{qc}) were compared with expected values (calculated from a standard solution of TN/TAM/TAc in RS) applying a ± 20% tolerance.




Fig 3. Thyronamine and thyronine chromatograms from Hep G2 cell lysate extracts spiked at 50 nM. 2 MRM Traces per analyte, overlaid. Upper trace: T₀ (4.12 min); 3-T₁ (4.54 min); 3-T₁' (4.80 min); 3,5-T₂ (4.60 min); 3,3'-T₂ (5.16 min); 3',5'-T₂ (5.35 min); T₃ (5.01 min); rT₃ (5.64 min); T₄ (5.36 min). **Lower trace:** T₀AM (4.65 min); T₁AM (5.30 min); 3,5-T₂AM (5.34 min); T₃AM (5.97 min); rT₃AM (7.10 min); T₄AM (inset, 6.56 min).

<https://doi.org/10.1371/journal.pone.0183482.g003>

Table 2. Linearity and lower limits of detection and quantification for TH and TAM spiked lysed Hep G2 and PCCL3 cell extracts in 24- and 6-well plates respectively.

Compound	Hep G2			PCCL3		
	R ² Range	LLOQ [nM]	LLOD [nM]	R ² Range	LLOQ [nM]	LLOD [nM]
T ₀	0.995–0.999	*	*	0.996–1.000	0.050	0.025
3-T ₁	0.996–1.000	0.063	0.031	0.997–1.000	0.050	0.025
3'-T ₁	0.995–1.000	0.031	0.016	0.997–1.000	0.050	0.025
3,5-T ₂	0.996–1.000	0.125	0.063	0.996–1.000	0.050	0.025
3,3'-T ₂	0.996–1.000	0.063	0.031	0.998–0.999	0.100	0.025
3',5'-T ₂	0.995–0.999	0.031	0.016	0.997–1.000	0.050	0.025
T ₃	0.995–1.000	0.063	0.016	0.997–0.999	0.040	**
rT ₃	0.997–1.000	0.063	0.031	0.998–0.999	0.100	0.010
T ₄	0.996–1.000	0.063	0.016	0.997–0.999	0.040	**
T ₀ AM	0.995–1.000	0.125	0.063	0.996–0.999	0.500	0.025
3-T ₁ AM	0.996–1.000	0.031	0.016	0.997–0.999	0.100	0.020
3,5-T ₂ AM	0.996–1.000	0.063	0.031	0.998–1.000	0.050	0.025
T ₃ AM	0.994–1.000	0.063	0.031	0.997–0.999	0.050	0.010
rT ₃ AM	0.995–0.999	0.063	0.031	0.997–0.999	0.050	0.010
T ₄ AM	0.996–0.998	0.250	0.125	0.995–0.999	0.500	0.100
T ₀ Ac	-	-	-	0.995–0.999	1.000	0.330
3-T ₁ Ac	-	-	-	0.996–0.999	1.000	0.600

R²: Pearson coefficient, 1/x weighted linear regression. LLOQ: lower limit of quantification LLOD: lower limit of detection.

*LLOD/LLOQ not quantifiable (lysed cell extracts contained ca. 0.10 nM T₀)

** LLOD not quantifiable (lysed cell extracts contained ca. 0.010nM T₄ and T₃)

<https://doi.org/10.1371/journal.pone.0183482.t002>

Excellent linearity for most THM in both cell lines and acceptable LLOD/LLOQ

1/x Weighted linear regression calibration curves for each analyte were obtained by plotting MRM peak area ratios (analyte/IS) against concentration over the working range 0.016–50 nM for Hep G2 cell lysate extracts (except for T₀: 0.50–50 nM) and 0.01–125 nM for PCCL3 cell lysate extracts. 1/x weighted linear correlation coefficients (r²) for all compounds are displayed in Table 2, and excellent linearity across the dynamic range of both methods for all analytes is evidenced by a minimum value for r² of 0.994.

The LLOD were in the range of 0.016–0.125 nM for Hep G2 cell lysate extracts in 24-well format and 0.01–0.6 nM for PCCL3 cell lysate extracts in 6-well format. The LLOQ ranged from 0.031–0.250 nM for Hep G2 and 0.04–1.0 nM for PCCL3. LLOD and LLOQ could not be determined for T₀ in Hep G2 cell lysate extracts, due to a background concentration of ca. 0.06 nM. Similarly, LLOD could not be calculated for T₃ and T₄ in PCCL3 cell lysate extracts because of ca. 0.01 nM background concentration for both TN.

Reasonable accuracy and precision for spiked quality-control samples (QQAL) in both cell lines

For Hep G2 inter- and intra-batch precision and accuracy were determined at low (LQAL), medium (MQAL) and high (HQAL) concentrations of TN and TAM. Certified reference solutions of three TN (3,3'-T₂, T₃ and T₄) were purchased from a commercial supplier, diluted accordingly and used alongside in-house prepared solutions for independent comparison. 8-Day inter-batch CVs for in-house QQALs (Table 3) were all < 13.5% across the entire

concentration range, with the exception of T₀ (not quantifiable for LQAL, otherwise < 5.7%), 3,5-T₂, T₀AM and T₄AM (LLOQ < LQAL, otherwise < 12.5%) and < 8.7% for certified reference QQALs (with the exception of LQAL for 3,3'-T₂ and T₃ which were < 17.8%). Intra-batch CVs were < 15.0% (in-house solutions only, Table 4). 8-Day inter-batch accuracy for in house QQALs were between 92.5–109.4% and for certified reference QQALs between 94.0–114.5% (Table 3). Intra-batch accuracies were determined for in-house solutions only and exhibited values in the range 92.0–120.0% (Table 4).

For PCCL3 inter-batch precision and accuracy were determined for QQALs using in-house prepared spike-solutions at 3 concentrations spanning the appropriate range for the occurrence of THM investigated in the accompanying application (0.10, 1.0 and 2.5 nM). 3-Day inter-batch CVs were all < 10.5% with the exception of the T₂, T₃ and T₄ TAM (< 23.1% for LQAL, < 14.2% for MQAL and < 16.4% for HQAL) (Table 5). Intra-batch CVs (Table 6), measured at 0.10 (LQAL), 12.5 (MQAL) and 40 (HQAL) nM were < 6.1% for all THM for HQAL, < 4.4% at MQAL (with the exception of rT₃AM, T₃AM and T₄AM which were below 19.4%) and < 9.1% at LQAL (with the exception of 3-T₁, 3'-T₁ and T₀AM, which were < 14%). 3-Day inter-batch accuracy for all QQALs (Table 5) were between 92.6–118.7% for all TN and TAM, with the exception of T₀AM, T₄AM, T₀Ac and T₁Ac for LQAL which was <LLOQ for these compounds. Intra-batch accuracies for all QQAL were in the range 88.7–108.8% (Table 6). The necessity for using weighted linear regression for calibration curves that span a large concentration range was evidenced by comparison of the mean accuracy value obtained at the LQAL (0.1 nM) concentration for PCCL3 cells with or without 1/x weighting across a range of TN/TAM (14 compounds) for a single day. The mean % error using 1/x weighted calibration curves was 12%, whereas for non-weighted curves the mean % error was 233%.

Table 3. Inter-day (n = 8) accuracy and precision for TH and TAM spiked lysed Hep G2 cell extracts in 24-well plates.

Compound	LQAL 0.10 nM*			MQAL 12.5 nM*			HQAL 40.0 nM*		
	Mean conc. [nM]	CV [%]	Accuracy [%]	Mean conc. [nM]	CV [%]	Accuracy [%]	Mean conc. [nM]	CV [%]	Accuracy [%]
T ₀	-	-	-	12.4	4.2	99	40.6	5.7	102
3-T ₁	0.10	9.9	104	12.5	8.3	100	41.1	5.7	103
3'-T ₁	0.11	8.7	107	12.6	5.9	100	41.1	7.1	103
3,5-T ₂	<LLOQ			12.6	7.2	101	41.9	6.8	105
3,3'-T ₂	0.11 [0.11]	8.7 [14.2]	109 [114]	12.6 [13.2]	6.6 [7.2]	100 [106]	41.5 [41.5]	6.9 [6.6]	104 [104]
3',5'-T ₂	0.11	10.8	106	12.6	5.6	101	40.8	5.5	102
T ₃	0.11 [0.09]	9.3 [17.8]	105 [94]	12.4 [11.8]	8.3 [8.7]	99 [95]	42.8 [39.4]	7. [6.0]	107 [99]
rT ₃	0.10	7.5	104	12.3	8.2	98	41.0	6.9	102
T ₄	0.09 [0.10]	13.5 [8.1]	108 [103]	9.8 [12.1]	8.3 [6.1]	98 [96]	35.0 [39.4]	7.1 [5.0]	109 [98]
T ₀ AM	<LLOQ			12.8	5.7	102	41.6	6.5	104
3-T ₁ AM	0.11	12.0	107	12.5	6.4	100	41.3	5.7	103
3,5-T ₂ AM	0.10	6.1	103	12.0	5.8	96	41.3	5.8	103
T ₃ AM	0.10	10.1	100	12.3	8.6	98	41.0	8.4	102
rT ₃ AM	0.10	11.4	104	12.1	11.1	97	41.2	8.8	103
T ₄ AM	<LLOQ			11.6	12.5	92	40.1	11.0	100

CV [%]: Percentage coefficient of variation. LQAL: low concentration spiked quality-control sample. MQAL: medium concentration spiked quality-control sample. HQAL: high concentration spiked quality-control sample. Values in parentheses were obtained by spiking samples with externally purchased TH standard solutions (see experimental section); all other values were from internally prepared TH standard solutions. Lysed cell extracts contained ca. 0.060 nM T₀ making quantification at 0.1 nM unreliable.

*All spiked concentrations for internally prepared TH standard solutions as labelled except: T₄ LQAL (0.08 nM), MQAL (10 nM), HQAL (32 nM).

<https://doi.org/10.1371/journal.pone.0183482.t003>

Table 4. Intra-day (n = 8) accuracy and precision for TH and TAM spiked lysed Hep G2 cell extracts in 24-well plates.

Compound	LQAL 0.10 nM			MQAL 12.5 nM			HQAL 40.0 nM		
	Mean conc. [nM]	CV [%]	Accuracy [%]	Mean conc. [nM]	CV [%]	Accuracy [%]	Mean conc. [nM]	CV [%]	Accuracy [%]
T ₀	-	-	-	12.1	6.9	97	38.9	6.4	97
3-T ₁	0.10	12.6	104	12.6	5.7	101	40.4	5.5	101
3'-T ₁	0.09	15.0	92	12.1	6.3	97	40.0	6.6	100
3,5-T ₂	<LLOQ			12.5	5.8	100	40.0	5.4	100
3,3'-T ₂	0.10	8.5	99	12.3	5.0	99	39.5	4.8	99
3',5'-T ₂	0.10	11.3	95	12.3	5.3	98	40.2	5.9	100
T ₃	0.10	10.2	104	12.2	5.1	98	38.7	7.6	97
rT ₃	0.12	11.9	116	11.5	3.1	92	38.2	8.0	96
T ₄	0.11	12.8	108	12.0	4.5	96	39.9	5.1	100
T ₀ AM	<LLOQ			12.2	5.1	98	38.7	7.6	97
3-T ₁ AM	0.11	9.6	108	12.5	5.4	100	40.1	6.1	100
3,5-T ₂ AM	0.12	6.3	120	12.2	6.1	97	40.7	5.1	102
T ₃ AM	0.12	11.2	116	11.9	6.5	95	39.6	7.3	99
rT ₃ AM	0.11	9.7	110	11.9	5.7	95	38.9	7.8	97
T ₄ AM	<LLOQ			12.9	4.3	103	38.0	7.6	95

Lysed cell extracts contained ca. 0.10 nM T₀ making quantification at 0.10 nM unreliable.

<https://doi.org/10.1371/journal.pone.0183482.t004>

Recoveries, matrix effects and process efficiencies as analysed in two cell lines

Recoveries, process efficiencies and matrix effects determined for extracted Hep G2 cell lysates in 24-well plate format are displayed in [Table 7](#). The recoveries for all TN, TAM and IS across

Table 5. Inter-day (n = 3) accuracy and precision for TH and TAM spiked lysed PCCL3 cell extracts in 6 well plates.

Compound	LQAL 0.10 nM			MQAL 1.0 nM			HQAL 2.5 nM		
	Mean Conc. [nM]	CV [%]	Accuracy [%]	Mean Conc. [nM]	CV [%]	Accuracy [%]	Mean Conc. [nM]	CV [%]	Accuracy [%]
T ₀	0.11	14.8	111	1.1	5.0	107	2.7	3.6	107
3-T ₁	0.11	9.6	112	1.1	1.9	105	2.7	1.5	107
3'-T ₁	0.11	7.0	106	1.0	1.4	104	2.6	1.6	105
3,5-T ₂	0.10	5.3	103	1.1	1.9	111	2.8	1.5	113
3,3'-T ₂	0.11	10.4	106	1.0	2.3	104	2.6	1.6	106
3',5'-T ₂	0.11	10.5	106	1.0	4.9	104	2.7	5.0	108
T ₃	0.10	1.7	96	1.0	2.6	100	2.6	2.6	106
rT ₃	0.11	10.1	107	1.0	1.8	100	2.6	1.6	104
T ₄	0.11	10.4	108	1.0	3.2	100	2.6	2.0	104
T ₀ AM	<LLOQ			1.0	7.9	101	2.6	7.2	106
3-T ₁ AM	0.11	8.1	111	1.0	1.9	105	2.7	1.6	109
3,5-T ₂ AM	0.11	13.1	110	1.0	3.1	102	2.7	4.0	108
T ₃ AM	0.12	21.8	119	1.0	9.1	95	2.6	9.1	104
rT ₃ AM	0.11	23.1	107	0.9	13.2	93	2.6	13.8	102
T ₄ AM	<LLOQ			0.9	14.2	87	2.4	16.4	95
T ₀ Ac	<LLOQ			0.9	5.9	90	2.7	2.5	107
3-T ₁ Ac	<LLOQ			1.0	1.8	101	2.5	1.6	100

<https://doi.org/10.1371/journal.pone.0183482.t005>

Table 6. Intra-day (n = 8) accuracy and precision for TH, TAM and TAc in spiked lysed PCCL3 cell extracts in 6-well plates.

Compound	0.10			12.5			40.0		
	Mean Calculated Concentration	% Relative Standard Deviation	% Accuracy	Mean Calculated Concentration	% Relative Standard Deviation	% Accuracy	Mean Calculated Concentration	% Relative Standard Deviation	% Accuracy
T0	0.10	8.6	103	12.7	3.8	102	39.8	1.4	99
3T1	0.10	12.6	99	12.5	3.2	100	39.3	2.2	98
3'T1	0.09	12.7	94	12.8	2.8	103	40.1	1.1	100
3,5T2	0.10	5.3	100	12.9	2.4	103	39.3	2.2	98
3,3'T2	0.10	9.1	98	12.8	1.2	102	40.5	1.3	101
3',5'-T2	0.11	6.6	108	12.6	2.4	101	39.9	1.9	100
T3	0.11	8.6	106	12.7	1.9	102	40.4	1.8	101
rT3	0.11	7.0	106	12.2	1.5	97	38.7	3.6	97
T4	0.10	4.5	103	12.6	1.6	101	39.9	1.6	100
T0AM	0.09	14.0	89	12.5	1.2	100	39.0	1.7	98
3T1AM	0.10	5.6	95	12.4	0.8	99	39.9	1.0	100
3,5T2AM	0.09	7.9	94	12.1	3.2	97	39.5	1.5	99
T3AM	0.10	8.8	104	11.1	17.0	89	39.4	3.2	98
rT3AM	0.10	6.3	101	11.1	18.8	89	38.4	4.2	96
T4AM	<LLOQ			11.4	19.4	91	36.9	6.1	92
T0Ac	<LLOQ			12.6	2.5	101	39.3	3.5	98
T1Ac	<LLOQ			13.6	4.4	109	40.2	3.7	101

<https://doi.org/10.1371/journal.pone.0183482.t006>

the range of QQALs (0.1 to 40 nM) were between 73–93%, the process efficiencies between 51–95%, and the matrix effects between 70–106%.

The same QA parameters were determined at a single concentration (1 nM) for extracted PCCL3 cell lysates in 6-well plate format. Recoveries for all TN, TAM and TAc, including IS varied from 75–101%, with % CV that varied between 2 and 17%. Process efficiencies were in the range 58–89% with % CV from 1–17%; matrix effects ranged from 77–98% with % CV between 1 and 14%. This data is displayed in [Table 8](#).

Application

Identification of TN, TAM and TAc metabolites generated by PCCL3 cells. Identification of THM in extracts of lysed PCCL3 cells and supernatants was confirmed or rejected by both R_t and MRM_{qc} as follows: T₀: confirmed (R_t expected: 4.13, found: 4.10 min; MRM_{qc} within -4/+7% and -5/+14% of expected values in cell lysates and supernatants respectively); T₀AM: confirmed (R_t expected: 4.51, found: 4.51 min; MRM_{qc} within -9/+4% and -14/+12% of expected values in cell lysates and supernatants respectively). 3-T₁Ac: confirmed in cell lysates (R_t expected: 5.14, found: 5.13 min; MRM_{qc} within -14/+16% of expected values); confirmed only in cell supernatants with > 10 min incubation time (MRM_{qc} -11/+15% of expected values, for incubation times < 10 min the MRM_{qc} exceeded +/- 20%). T₀Ac: unconfirmed in cell lysates (R_t expected: 4.71, found: 4.73 min; MRM_{qc} irreproducible and in many cases > ±20% of expected value); confirmed in supernatants with > 10 min incubation time (MRM_{qc} -18/+3% of expected values)).

Method application to cell extracts and supernatants of rat thyrocytes after incubation with mono-iodinated thyroid hormone metabolites. PCCL3 thyrocytes were incubated with 500 nM of mono-iodinated representatives of the TAM, the TN and the TAc family for 5, 10 and 30 min. After incubation with 3-T₁AM the un-metabolised parent substance was

Table 7. Inter-day (n = 8) process efficiency, matrix effect and recovery for TH and TAM spiked lysed Hep G2 cell extracts in 24-well plates.

Compound	0.10 nM						12.5 nM						40.0 nM														
	PE	±	% RSD	ME	±	% RSD	RE	±	% RSD	PE	±	% RSD	ME	±	% RSD	RE	±	% RSD	PE	±	% RSD	ME	±	% RSD	RE	±	% RSD
² H ₄ -3-T ₁ AM	69	±	12	84	±	11	82	±	10	69	±	10	82	±	9	81	±	13	71	±	10	86	±	6	83	±	13
¹³ C ₆ -3,3'T ₂	73	±	8	87	±	8	85	±	8	74	±	8	83	±	5	85	±	11	75	±	9	87	±	6	86	±	11
¹³ C ₆ -T ₃	78	±	8	91	±	7	86	±	8	80	±	8	88	±	5	89	±	11	78	±	10	88	±	8	88	±	12
¹³ C ₆ -rT ₃	67	±	10	84	±	11	80	±	7	69	±	8	79	±	8	85	±	10	70	±	10	84	±	8	84	±	12
¹³ C ₆ -T ₄	67	±	13	84	±	13	80	±	10	69	±	11	80	±	10	84	±	11	68	±	10	81	±	10	85	±	12
T ₀	148	±	17	170	±	15	87	±	9	71	±	8	91	±	7	76	±	12	72	±	11	92	±	6	78	±	13
3-T ₁	78	±	10	85	±	8	91	±	8	74	±	10	88	±	8	83	±	13	75	±	10	90	±	7	84	±	11
3'-T ₁	74	±	9	84	±	12	88	±	10	70	±	9	84	±	7	81	±	10	73	±	10	87	±	7	84	±	13
3,5-T ₂	<LLOQ			<LLOQ			<LLOQ			81	±	10	96	±	8	83	±	12	80	±	11	94	±	6	86	±	10
3,3'-T ₂	75	±	9	82	±	10	92	±	10	73	±	9	87	±	6	84	±	10	76	±	13	88	±	6	86	±	11
3',5'-T ₂	70	±	13	80	±	13	88	±	11	67	±	9	83	±	9	81	±	11	71	±	10	86	±	7	82	±	11
T ₃	81	±	13	89	±	10	91	±	10	72	±	11	87	±	8	81	±	12	75	±	11	90	±	7	84	±	11
rT ₃	75	±	10	86	±	11	88	±	10	68	±	8	85	±	7	83	±	10	71	±	10	87	±	6	82	±	11
T ₄	72	±	20	83	±	18	86	±	12	64	±	15	83	±	13	79	±	10	68	±	15	86	±	9	79	±	12
T ₀ AM	<LLOQ			<LLOQ			<LLOQ			74	±	13	87	±	9	83	±	13	77	±	13	90	±	7	84	±	12
3-T ₁ AM	75	±	13	87	±	15	87	±	14	69	±	11	87	±	8	78	±	13	72	±	12	90	±	5	80	±	13
3,5-T ₂ AM	71	±	13	84	±	12	85	±	14	64	±	14	85	±	11	74	±	12	67	±	14	86	±	8	78	±	15
T ₃ AM	64	±	26	77	±	33	85	±	19	61	±	15	79	±	16	80	±	11	66	±	13	84	±	9	78	±	16
rT ₃ AM	63	±	14	80	±	14	79	±	19	59	±	15	79	±	11	79	±	11	62	±	16	81	±	8	76	±	18
T ₄ AM	<LLOQ			<LLOQ			<LLOQ			51	±	21	70	±	11	79	±	10	60	±	21	74	±	15	73	±	24

Values were obtained by spiking samples with in-house prepared TH standard solutions. Lysed cell extracts contained ca. 0.10 nM T₀, thus results at 0.10 nM are skewed.

<https://doi.org/10.1371/journal.pone.0183482.t007>

detectable in cell lysates with a clear increase in concentration over time (Fig 4A) indicating its cellular uptake and intracellular accumulation. The extracellular concentration of 3-T₁AM remained almost constant over the course of 30 min (Fig 4C). The de-iodination product T₀AM and the de-amination product 3-T₁Ac were detected in all cell extract samples (Fig 4A), and were below LLOD in DMSO treated negative controls and in the parent stock solution. These metabolites accumulated in the supernatants as evidenced (Fig 4C) by their increase in concentration over time. 3-T₁ was detectable in cell lysates after incubation, indicating cellular uptake or binding of the parent substance (Fig 4B). Compared to 3-T₁AM incubation, binding or uptake of 3-T₁ occurred to a substantially lesser extent. Furthermore, in contrast to 3-T₁AM, no clear increase in 3-T₁ concentration was observed with incubation time in cell lysates. The T₁ concentration in supernatants remained almost constant (Fig 4D).

The de-iodination product T₀ was detected at low concentrations in all lysates after incubation with 3-T₁ and in low but moderately increasing concentration in cell supernatants (Fig 4D), but not in DMSO controls or the 3-T₁ stock solution. After incubation with 3-T₁Ac, only small amounts of the parent substance were detectable in cell lysates for each time point (highest concentration approximately 5 nM, data not shown), the bulk of the 3-T₁Ac remained unchanged in the cell supernatant (data not shown). As mentioned above, the putative de-iodination product T₀Ac was detected in cell lysates, but (even after 30 min incubation time) at concentrations close to the limit of quantification making the confirmation of its identity irreproducible. However, T₀Ac was detected in supernatants after 10 and 30 min incubation at mean concentrations of 5.6 and 10 nM respectively (data not shown).

Table 8. Inter-day (n = 3) process efficiency, matrix effect and recovery for TH and TAM spiked lysed PCCL3 cell extracts in 6 well plates.

Compound	QQAL 1 nM (n = 3)								
	PE	±	SD	ME	±	SD	RE	±	SD
² H ₄ -T ₀ Ac	89	±	12	89	±	8	100	±	17
² H ₄ -T ₁ Ac	80	±	7	98	±	7	82	±	6
² H ₄ -3-T ₁ AM	76	±	6	90	±	7	85	±	7
¹³ C ₆ -3,3'T ₂	77	±	5	94	±	7	82	±	4
¹³ C ₆ -T ₃	80	±	6	93	±	3	85	±	7
¹³ C ₆ -rT ₃	66	±	7	88	±	6	75	±	6
¹³ C ₆ -T ₄	72	±	10	92	±	4	79	±	11
T ₀	72	±	3	87	±	2	82	±	3
3-T ₁	81	±	5	90	±	3	90	±	3
3'-T ₁	78	±	6	93	±	3	84	±	4
3,5-T ₂	88	±	9	96	±	4	91	±	6
3,3-T ₂	79	±	4	93	±	3	85	±	2
3',5'-T ₂	75	±	6	94	±	1	80	±	7
T ₃	74	±	2	91	±	2	81	±	4
rT ₃	70	±	5	89	±	4	79	±	5
T ₄	68	±	10	94	±	3	78	±	3
T ₀ AM	84	±	8	97	±	2	87	±	8
3-T ₁ AM	79	±	4	90	±	1	88	±	5
3,5-T ₂ AM	73	±	5	88	±	3	83	±	3
T ₃ AM	69	±	1	82	±	11	84	±	9
rT ₃ AM	64	±	4	81	±	10	80	±	10
T ₄ AM	58	±	10	79	±	13	80	±	14
T ₀ Ac	77	±	17	77	±	14	101	±	12
3-T ₁ Ac	82	±	8	96	±	13	86	±	5

<https://doi.org/10.1371/journal.pone.0183482.t008>

Table 9 shows mass balance data for the above incubation experiments. For all three incubation substances over all time points, mass balance totals were in the range of 66–84%. For 3-T₁ incubation, 4–5% of the parent substance was accounted for in the extracts of cell lysates, 62–63% was evident in the supernatant. The de-iodinated metabolite T₀ accounted for only 0.2% of the total mass balance in cell lysates after 30 min and in marginally higher amounts (0.4%) in supernatants. The percentage balance of 3-T₁AM remained constant at ca. 40% in supernatants over the entire course of 3-T₁AM incubation, but increased from 29% after 5 min to 41% after 30 min in cell lysates. The percentage balance of the de-iodination metabolite T₀AM, which was present at low concentrations, increased in supernatants from 0.1 (10 min) to 0.8% (30min), but stayed constant at 0.2–0.4% in cell lysates. The de-amination product 3-T₁Ac constituted only 0.1–0.2% of the total in cell lysates for all time points and was detected only after 30 min in supernatants at 1.2%. Finally, T₁Ac was not substantially accumulated by PCCL3 cells and remained between 77–81% of the total accountable pmol amount in supernatants. The % pmol of the de-iodinated product T₀Ac increased from 1.1 to 1.9% of total pmol from 10 to 30 min.

Discussion

LC-MS/MS method validation

We applied our previously published LC-MS/MS method [12], with minor modifications, to the analysis of Hep G2 cell lysate extracts. In brief, the method features a 10 min analysis time




Fig 4. Application of the LC-MS/MS method to extracts of PCCL3 rat thyrocytes. THM in extracts of PCCL3 (A, B) cell and supernatants (C, D) after incubation with 500 nM 3-T₁AM (A, C) or 3-T₁ (B, D) in Krebs-Ringer buffer as measured by LC-MS/MS. Bars represent means of 3 independent experiments carried out on different days +/- SEM.

<https://doi.org/10.1371/journal.pone.0183482.g004>

Table 9. Mass balance in cells/supernatants for THM after incubation of mono-iodinated THM with PCCL3 rat thyrocytes.

Incubation compound (pmol)	Incubation time	5min						10min						30min					
		Cells		Supernatants		Total pmol %	Cells		Supernatants		Total pmol %	Cells		Supernatants		Total pmol %			
		pmol (%)	(%)	pmol (%)	(%)		pmol (%)	(%)	pmol (%)	(%)		pmol (%)	(%)	pmol (%)	(%)				
T ₁ (775.9)	T ₁	32 (4%)	(4%)	498 (62%)	(62%)	66%	34 (4%)	(4%)	502 (63%)	(63%)	67%	39 (5%)	(5%)	507 (63%)	(63%)	69%			
	T ₀	1 (0.1%)	(0.1%)	0.3 (0.03%)	(0.03%)		1 (0.2%)	(0.2%)	1 (0.1%)	(0.1%)		2 (0.2%)	(0.2%)	3 (0.4%)	(0.4%)				
3-T ₁ AM (812.7)	3-T ₁ AM	238 (29%)	(29%)	323 (40%)	(40%)	69%	276 (34%)	(34%)	333 (41%)	(41%)	76%	334 (41%)	(41%)	329 (40%)	(40%)	84%			
	T ₀ AM	1 (0.2%)	(0.2%)	1 (0.1%)	(0.1%)		3 (0.4%)	(0.4%)	2 (0.2%)	(0.2%)		3 (0.3%)	(0.3%)	6 (0.8%)	(0.8%)				
	3-T ₁ Ac	1 (0.1%)	(0.1%)	nd (0%)	(0%)		1 (0.2%)	(0.2%)	nd (0%)	(0%)		1 (0.1%)	(0.1%)	9 (1.2%)	(1.2%)				
3-T ₁ Ac (777.6)	3-T ₁ Ac	4 (0.6%)	(0.6%)	596 (77%)	(77%)	77%	4 (0.5%)	(0.5%)	599 (77%)	(77%)	79%	3 (0.3%)	(0.3%)	632 (81%)	(81%)	84%			
	T ₀ Ac	1 (0.1%)	(0.1%)	nd (0%)	(0%)		1 (0.2%)	(0.2%)	8 (1.1%)	(1.1%)		1 (0.1%)	(0.1%)	14 (1.9%)	(1.9%)				

pmol: absolute amount in cell lysate (0.4 mL) or supernatant (1.5 mL), (%); relative % pmol for individual compound with respect to pmol of incubation compound, **Total pmol %:** summation of % pmol for all compounds detected in cell lysates or supernatants with respect to pmol of incubation compound. % pmol values were not adjusted to account for relative recoveries.

<https://doi.org/10.1371/journal.pone.0183482.t009>

for the 9 TN and 6 TAM in positive ion mode and complete baseline chromatographic separation for the three isobaric T_2 isomers (due to identical precursor and product ion selections, these isomers must be separated chromatographically for individual quantification). Wherever possible, structurally identical, stable isotopically labelled IS for quantification were used for each analyte to enhance the accuracy of the method [20]. In cases where no identical IS was commercially available, an IS with the closest structural similarity was applied.

We also adapted this method for the analysis of PCCL3 cell lysate extracts by using positive/negative electrospray ionisation switching to further incorporate negative ion electrospray MRM transitions for T_0Ac and T_1Ac and their respective ($^2H_4-T_0Ac$ and $^2H_4-T_1Ac$) IS.

Across the range of QQAL concentrations and cell lines PE were ca. 17–24% lower than for those measured in cell culture media [12]. For the most part RE were not substantially worse for cell lysates extracts (4–13% less across both cell lines and QQAL concentrations) compared to culture media extracts. ME for Hep G2 in contrast to cell culture media tended to decrease substantially with increasing QQAL concentration and were 22% higher for the LQAL (0.10 nM) than for the culture media LQAL (0.75 nM), 15% higher at MQAL (10 vs 12.5 nM) but comparable at HQAL (40 vs 80 nM) concentrations. These results would seem to be consistent with greater LC-MS/MS ion suppression effects caused by the dirtier cell lysate matrix. The marginally lower RE could arise from hormone adhesion to cell culture plates.

Application of the LC-MS/MS method to extracts of PCCL3 rat thyrocytes

We successfully applied the described LC-MS/MS method to investigate THM in PCCL3 thyrocytes after incubation with 3- T_1AM , 3- T_1 and 3- T_1Ac . The data obtained demonstrate that 3- T_1AM is de-iodinated to T_0AM and de-aminated to 3- T_1Ac , while 3- T_1 is de-iodinated to T_0 by PCCL3 cells. This indicates 3- T_1AM and 3- T_1 uptake and intracellular metabolism by the cell line and confirms tyrosyl ring de-iodination activity as well as activity of an amino oxidase in PCCL3 cells as previously suggested [21]. Thus, the interpretation of 3- T_1AM effects on the hypothalamus-pituitary-thyroid (HPT) axis in mice treated with this biogenic amine [21] and actions observed in thyrocyte cell culture, should consider the possibility of direct biological action of its metabolites. Our findings of *in vitro* 3- T_1AM de-iodination to T_0AM and de-amination to 3- T_1Ac are in line with previous data obtained in different cell lines by others ([22], [23]). However, in the related thyrocyte model FRTL-5 the catabolite T_0Ac was detected additionally. In this study, virtually all of the produced T_0Ac was detected in the incubation medium [22].

PCCL3 cells accumulate T_1 and T_1AM and T_0AM , but not T_1Ac or T_0Ac . We observed a time-dependent increase in recovery, probably due to cellular accumulation of TH/THM, which may indicate some adhesion of TH/THM to cell culture plastic in the absence of TH/THM binding proteins in incubation medium. Cells de-iodinate T_1 to T_0 , which mainly accumulates in the supernatant, suggesting uptake of T_1 and export of T_0 . Cells also accumulate T_1AM and de-iodinate it to T_0AM , the majority of which is again exported; T_1Ac also appears in supernatants, but only at later time points. In contrast, T_1Ac is not accumulated by PCCL3 cells, but rather exported over time and also de-iodinated to T_0Ac , which appears in the supernatant.

Generally, only limited metabolism of these TH/THM occurs in PCCL3 over the first 30 min. Less than 1% of T_1 is de-iodinated to T_0 , around 5% of 3 T_1AM is metabolized to T_0AM and 3 T_1Ac and less than 2.5% of 3 T_1Ac is de-iodinated to T_0Ac when offered as the sole substrate to PCCL3 cells. Recoveries of TH and THM carrying only one or no iodine atom ranges between 66 and 84% (values were not corrected for any losses in the sample extraction process,

the relative recoveries of the method would account for the fact that mass balances do not total 100%). These recoveries are comparable to those of the classical TH, T₄ and T₃, in this and other cell lines during uptake and metabolism experiments ([12], [13]). Interestingly, 3T₁Ac exhibits the highest recovery and lowest accumulation in cell extracts, which might indicate that both its precursors 3T₁ and 3T₁AM are more avidly bound by cells, and metabolized intracellularly. Some (irreversible) binding of both substrates to cellular targets cannot be excluded. However, the time dependent increase in recovery may also indicate that initially these TH/THM compounds may adsorb to cell culture plastic surfaces and thus escape complete extraction by our currently employed procedure.

The observation of significant amounts of iodine-free T₀ in extracts of HepG2 cells cannot yet be explained satisfactorily. As blank solutions, reagents and standards used during the pre-analytical workup and LC-MS/MS analysis do not contain measurable amounts of T₀, we suppose that it might be carried over from T₀-containing FCS, which is required for the initial steps of cell culture. Alternatively, HepG2 cells may transiently accumulate T₀ from serum in the culture medium or eventually generate it from accumulated TH and THM. Up to now, fate and metabolism of T₀ have rarely been studied in cell culture experiments, whilst formation and renal secretion into urine has been reported using various methods ([24], [25], [26]). In humans, T₀ has been described and quantified by immunoassays and mass spectrometry at low nanomolar concentrations in urine, indicating that not all iodothyronines are fully de-iodinated to recover their content of the essential trace element iodine for its recycling by the thyroid gland.

In conclusion, we have shown that TH and THM can be quantified precisely and accurately to a LLOQ of 0.031–0.125 nM for TN (except for T₀ in HepG2 cells, which contained baseline concentrations) and 0.031–0.05 nM for TAM in Hep G2 and PCCL3 cell lysate extracts. ME for all TN, TAM and TAc were acceptable for both cell lines (>70%) and RE were >73% in all cases. For both cell lines, PE, an indication of the overall efficacy of pre-analytic extraction and LC-MS/MS measurement for TH and THM in the sample matrix, were substantially less than for culture media [12]. This is due to increased ME from the cell sample matrix. PE were however reproducible for all TN, TAM and TAc in PCCL3 cells (10% CV or better) with the exception of T₀Ac (17% CV) and better than 20% for all TN and TAM in Hep G2 cell lysate extracts with the exception of T₃AM and T₄AM (26 and 21% CV respectively). We therefore expect the reported methods to find wide applicability for the quantitative analysis of TH and THM in cell lysates.

By applying our positive/negative electrospray switched tandem LC-MS/MS method to the analysis of cell lysate and supernatant extracts recovered from the 30 min incubation of mono-iodinated TN, TAM and TAc with PCCL3 rat thyrocytes, we were able to follow the progress of limited de-iodination of T₁, T₁AM and T₁Ac and the de-amination of T₁AM to T₁Ac.

Acknowledgments

This work was supported by the Deutsche Forschungsgemeinschaft DFG (KO922-16/1-2 to J.K. and KO922-17/1-2 to J.K.) in the framework of the DFG funded Priority Programme SPP1629 THYROIDTRANSACT and a Charité thesis scholarship (to N.S.).

Author Contributions

Conceptualization: Keith H. Richards, Nancy Schanze, Ray Monk, Eddy Rijntjes, Daniel Rathmann, Josef Köhrle.

Data curation: Keith H. Richards, Nancy Schanze.

Formal analysis: Keith H. Richards, Nancy Schanze.

Funding acquisition: Josef Köhrle.

Investigation: Keith H. Richards, Nancy Schanze, Ray Monk.

Methodology: Keith H. Richards, Nancy Schanze, Ray Monk, Daniel Rathmann.

Project administration: Josef Köhrle.

Resources: Josef Köhrle.

Supervision: Josef Köhrle.

Validation: Keith H. Richards.

Writing – original draft: Keith H. Richards, Nancy Schanze, Ray Monk, Eddy Rijntjes, Daniel Rathmann, Josef Köhrle.

Writing – review & editing: Keith H. Richards, Nancy Schanze, Ray Monk, Eddy Rijntjes, Daniel Rathmann, Josef Köhrle.

References

1. Pilo A, Iervasi G, Vitek F, Ferdeghini M, Cazzuola F, Bianchi R. Thyroidal and peripheral production of 3,5,3'-triiodothyronine in humans by multicompartmental analysis. *Am J Physiol*. 1990; 258(4 Pt 1): E715–26. PMID: [2333963](#).
2. Köhrle J. Local activation and inactivation of thyroid hormones: the deiodinase family. *Mol Cell Endocrinol*. 1999; 151(1–2):103–19. PMID: [10411325](#).
3. Dentice M, Salvatore D. Deiodinases: the balance of thyroid hormone: local impact of thyroid hormone inactivation. *J Endocrinol*. 2011; 209(3):273–82. <https://doi.org/10.1530/JOE-11-0002> PMID: [21398344](#).
4. Hoefig CS, Köhrle J, Brabant G, Dixit K, Yap B, Strasburger CJ, et al. Evidence for extrathyroidal formation of 3-iodothyronamine in humans as provided by a novel monoclonal antibody-based chemiluminescent serum immunoassay. *J Clin Endocrinol Metab*. 2011; 96(6):1864–72. Epub 2011/04/15. <https://doi.org/10.1210/jc.2010-2680> PMID: [21490071](#).
5. Lehmpfuhl I, Brabant G, Wallaschofski H, Ruchala M, Strasburger CJ, Köhrle J, et al. Detection of 3,5-diiiodothyronine in sera of patients with altered thyroid status using a new monoclonal antibody-based chemiluminescence immunoassay. *Thyroid*. 2014; 24(9):1350–60. Epub 2014/06/27. <https://doi.org/10.1089/thy.2013.0688> PMID: [24967815](#).
6. Visser TJ. Pathways of thyroid hormone metabolism. *Acta Med Austriaca*. 1996; 23(1–2):10–6. PMID: [8767510](#).
7. Wu SY, Green WL, Huang WS, Hays MT, Chopra IJ. Alternate pathways of thyroid hormone metabolism. *Thyroid*. 2005; 15(8):943–58. <https://doi.org/10.1089/thy.2005.15.943> PMID: [16131336](#).
8. Jonas W, Lietzow J, Wohlgemuth F, Hoefig CS, Wiedmer P, Schweizer U, et al. 3,5-Diiodo-L-thyronine (3,5-t2) exerts thymimetic effects on hypothalamus-pituitary-thyroid axis, body composition, and energy metabolism in male diet-induced obese mice. *Endocrinology*. 2015; 156(1):389–99. <https://doi.org/10.1210/en.2014-1604> PMID: [25322465](#); PubMed Central PMCID: [PMC4272399](#).
9. Padron AS, Neto RA, Pantaleao TU, de Souza dos Santos MC, Araujo RL, de Andrade BM, et al. Administration of 3,5-diiiodothyronine (3,5-T2) causes central hypothyroidism and stimulates thyroid-sensitive tissues. *J Endocrinol*. 2014; 221(3):415–27. <https://doi.org/10.1530/JOE-13-0502> PMID: [24692290](#); PubMed Central PMCID: [PMC4045230](#).
10. Scanlan TS, Suchland KL, Hart ME, Chiellini G, Huang Y, Kruzich PJ, et al. 3-Iodothyronamine is an endogenous and rapid-acting derivative of thyroid hormone. *Nat Med*. 2004; 10(6):638–42. Epub 2004/05/18. <https://doi.org/10.1038/nm1051> PMID: [15146179](#).
11. Bräulke LJ, Klingenspor M, DeBarber A, Tobias SC, Grandy DK, Scanlan TS, et al. 3-Iodothyronamine: a novel hormone controlling the balance between glucose and lipid utilisation. *J Comp Physiol B*. 2008; 178(2):167–77. Epub 2007/10/04. <https://doi.org/10.1007/s00360-007-0208-x> PMID: [17912534](#).
12. Rathmann D, Rijntjes E, Lietzow J, Köhrle J. Quantitative Analysis of Thyroid Hormone Metabolites in Cell Culture Samples Using LC-MS/MS. *Eur Thyroid J*. 2015; 4(Suppl 1):51–8. <https://doi.org/10.1159/000430840> PMID: [26601073](#); PubMed Central PMCID: [PMC4640291](#).

13. Johannes J, Braun D, Kinne A, Rathmann D, Köhrle J, Schweizer U. Few Amino Acid Exchanges Expand the Substrate Spectrum of Monocarboxylate Transporter 10. *Mol Endocrinol*. 2016; 30(7):796–808. <https://doi.org/10.1210/me.2016-1037> PMID: 27244477.
14. Committee BC. Guidance for industry: Bioanalytical Method Validation. In: CVM C, editor. <http://www.fda.gov/downloads/Drugs/Guidances/ucm070107.pdf>; FDA; 2001. p. 1–25.
15. Fusco A, Berlingieri MT, Di Fiore PP, Portella G, Grieco M, Vecchio G. One- and two-step transformations of rat thyroid epithelial cells by retroviral oncogenes. *Mol Cell Biol*. 1987; 7(9):3365–70. PMID: 3670314; PubMed Central PMCID: PMC367981.
16. Fusco A, Berlingieri MT, di Fiore PP, Grieco M, Portella G, M S, et al. Epithelial Transformation by Newly Isolated Murine Retroviruses. *From Oncogenes to Tumor Antigens*. Elsevier Science Publishers B, Amsterdam; 1985.
17. Vogeser M, Seger C. Pitfalls Associated with the Use of Liquid Chromatography-Tandem Mass Spectrometry in the Clinical Laboratory. *Clinical Chemistry*. 2010; 56(8):1234–44. <https://doi.org/10.1373/clinchem.2009.138602> PMID: 20511452
18. Matuszewski BK, Constanzer ML, Chavez-Eng CM. Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC-MS/MS. *Analytical chemistry*. 2003; 75(13):3019–30. PMID: 12964746.
19. Buhrman DL, Price PI, Rudewiczcor PJ. Quantitation of SR 27417 in human plasma using electrospray liquid chromatography-tandem mass spectrometry: A study of ion suppression. *J Am Soc Mass Spectrom*. 1996; 7(11):1099–105. [https://doi.org/10.1016/S1044-0305\(96\)00072-4](https://doi.org/10.1016/S1044-0305(96)00072-4) PMID: 24203071.
20. Richards K, Rijntjes E, Rathmann D, Köhrle J. Avoiding the pitfalls when quantifying thyroid hormones and their metabolites using mass spectrometric methods: the role of quality assurance. *Molecular and Cellular Endocrinology*. 2017;in press.
21. Schanze N, Jacobi SF, Rijntjes E, Mergler S, Del Olmo M, Hoefig CS, et al. 3-Iodothyronamine Decreases Expression of Genes Involved in Iodide Metabolism in Mouse Thyroids and Inhibits Iodide Uptake in PCCL3 Thyrocytes. *Thyroid*. 2017; 27(1):11–22. Epub 2016 Dec 21. <https://doi.org/10.1089/thy.2016.0182> PMID: 27788620; PubMed Central PMCID: PMC27788620.
22. Agretti P, De Marco G, Russo L, Saba A, Raffaelli A, Marchini M, et al. 3-Iodothyronamine metabolism and functional effects in FRTL5 thyroid cells. *J Mol Endocrinol*. 2011; 47(1):23–32. <https://doi.org/10.1530/JME-10-0168> PMID: 21511808.
23. Ghelardoni S, Chiellini G, Frascarelli S, Saba A, Zucchi R. Uptake and metabolic effects of 3-iodothyronamine in hepatocytes. *J Endocrinol*. 2014; 221(1):101–10. <https://doi.org/10.1530/JOE-13-0311> PMID: 24627446.
24. Pittman CS, Buck MW, Chambers JB Jr., Urinary metabolites of 14 C-labeled thyroxine in man. *J Clin Invest*. 1972; 51(7):1759–66. <https://doi.org/10.1172/JCI106977> PMID: 5032524; PubMed Central PMCID: PMC367981.
25. Willetts P, Crossley DN, Ramsden DB, Hoffenberg R. The role of thyronine in thyroid hormone metabolism. *J Clin Endocrinol Metab*. 1979; 49(4):658–60. <https://doi.org/10.1210/jcem-49-4-658> PMID: 479355.
26. Ramsden DB, Farmer M, Mohammed MN, Willetts P. Mass spectral properties of volatile derivatives of thyronine (T0) and use of these in the study of thyronine excretion in eu-, hyper- and hypothyroidism. *Biomed Mass Spectrom*. 1984; 11(4):193–8. <https://doi.org/10.1002/bms.1200110411> PMID: 6733255.