

Prospects and Frontiers of Stem Cell Toxicology

Shuyu Liu,^{1,2} Nuoya Yin,^{1,2} and Francesco Faiola^{1,2}

The development of stem cell biology has revolutionized regenerative medicine and its clinical applications. Another aspect through which stem cells would benefit human health is their use in toxicology. In fact, owing to their ability to differentiate into all the lineages of the human body, including germ cells, stem cells, and, in particular, pluripotent stem cells, can be utilized for the assessment, *in vitro*, of embryonic, developmental, reproductive, organ, and functional toxicities, relevant to human physiology, without employing live animal tests and with the possibility of high throughput applications. Thus, stem cell toxicology would tremendously assist in the toxicological evaluation of the increasing number of synthetic chemicals that we are exposed to, of which toxicity information is limited. In this review, we introduce stem cell toxicology, as an emerging branch of *in vitro* toxicology, which offers quick and efficient alternatives to traditional toxicology assessments. We first discuss the development of stem cell toxicology, and we then emphasize its advantages and highlight the achievements of human pluripotent stem cell-based toxicity research.

Keywords: stem cell toxicology, developmental toxicity, pluripotent stem cells, *in vitro* toxicology, embryonic stem cells

Introduction

TOXICOLOGY, THE STUDY of the adverse effects of chemicals on living organisms, serves human society in many ways, to not only protect humans and the environment from the deleterious effects of toxicants but also facilitate the development of safer chemicals to be employed as clinical drugs, pesticides, food additives, and so on. Toxicity can be acute or chronic, comes from different routes of exposure, varies from one organ to another, as well as yields to different outcomes according to age, genetic background, gender, diet, physiological conditions, and the health status of the organism.

Many *in vivo* and *in vitro* testing methods are available in toxicology. Although animal tests have played a critical role in toxicology, a trend for an alternative *in vitro* toxicology, mostly consistent in cell-based assays, started in the 1950s with the publication of the 3R principle (Replacement, Reduction, and Refinement) [1]. In fact, *in vitro* toxicity tests take advantage of well-developed cell culture protocols and pre-validated cell-based toxicity tests with multiple biological endpoints. Besides, cell-based toxicity tests usually require shorter time and still allow for more replicates, than *in vivo* assays do [2].

Currently, the list of potentially hazardous chemicals we can get exposed to is continuously and rapidly updated, due

to the exponential discovery and production of artificial materials. Therefore, it is extremely urgent to develop high-throughput experimental systems that will let us screen for the potential toxicity of all these substances.

Hence, in this review, we introduce stem cell toxicology, an emerging branch of *in vitro* toxicology, which offers effective and efficient alternatives to traditional toxicology assessments. We first discuss the development of stem cell toxicology, and then emphasize its advantages and highlight the achievements of human pluripotent stem cell (hPSC)-based toxicity research. We hope this review will be inspiring and innovative, and nurture more efforts into stem cell-based toxicology research.

Development of Stem Cell Toxicology

Stem Cell Toxicology, originally inspired by the mouse embryonic stem cell test (EST), includes toxicology research on stem cells and stem cell-derived differentiated cells. Stemness allows for *in vitro* cell maintenance for many generations without loss of differentiation potential and without genetic manipulations. Moreover, contrary to somatic or cancer cells, stem cell-based toxicity models let us perform developmental toxicity assays, in particular with mouse pluripotent stem cell- and hPSC-based models. Therefore, those unique advantages of

¹State Key Laboratory of Environmental Chemistry and Ecotoxicology, Research Center for Eco-Environmental Sciences, Chinese Academy of Sciences, Beijing, P.R. China.

²College of Resources and Environment, University of Chinese Academy of Sciences, Beijing, P.R. China.

© Shuyu Liu et al. 2017; Published by Mary Ann Liebert, Inc. This is an Open Access article distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

stem cell-based toxicity models call for further development and utilization of this new crossover field.

Mouse stem cell-based toxicology

After the successful derivation of mouse embryonic stem cells (mESCs) [3] and induction of mouse-induced PSCs (iPSCs) [4], the rapid development of stem cell technology also allowed for mESC-based toxicology research.

The first toxicity test based on mESCs was designed to test the cytotoxicity of potential teratogens [5]. It not only showed overall consistent results with *in vivo* data but also suggested that mESCs were more sensitive to tested chemicals, compared with mouse fibroblasts [5]. After this encouraging attempt, the mouse EST was proposed as the first complete and accepted *in vitro* toxicity test and has been continually improved and refined ever since [6]. The original EST consisted of, on 10-day toxicant exposure, MTT cytotoxicity assays on D3 mESCs and 3T3 mouse fibroblasts for the detection of IC50s (half maximal inhibitory concentrations), as well as a cardiogenic differentiation test with D3 mESCs, to determine the ID50 (the equivalent of IC50 for differentiating cells) [6].

Several efforts have been addressed to the optimization of the EST, particularly regarding its efficiency in experimental operation, accuracy, and objectiveness of endpoint determination [7], as well as toward high-throughput applications and extensive utilization (the main refinements of the mESC-based EST are summarized in Table 1). One early improvement was the introduction of molecular endpoints for gene and protein expression. For instance, the expression of the cardiac marker myosin heavy chain, quantified by flow cytometry, helped judge to which extent cardiogenesis was affected by tested chemicals [26–28] in a more objective and accurate way, compared with the original scoring of beating clusters of cardiomyocytes by microscopic observation [6]. In another example, groups of genes representative of the whole genome and that were involved in toxicity responses were selected for more predictive, reliable, and yet less laborious assays [10,20,29]. More recently, a shortened EST based on the detection of the cardiac and neural crest *Hand1* gene expression, through a luciferase reporter assay, showed good reproducibility [22]. In addition, studies on the effects of flusilazole and monophthalates on mESC cardiogenesis convinced that low chemical concentrations, although not cytotoxic, compromised mESC cardiogenesis by down-regulating the expression of related genes in a dose-dependent manner [20,30]. These refinements allowed for the investigation of the underlying molecular events triggered by chemical exposure, especially for changes in molecular levels that could be relevant for development, instead of merely cell viability.

Another major strategy to improve the original EST has taken advantage of differentiation procedures besides cardiogenesis, such as neural differentiation, endothelial differentiation, and osteogenesis. Such strategies were motivated by the fact that in some cases the traditional EST underestimated the possibility that different developmental lineages might show different responses on exposure to the same chemicals [8,14,31]. The use of cardiogenic induction as the only differentiation procedure would not provide a comprehensive toxicity assessment. Therefore, to address this concern, the neural EST was introduced and successfully validated with six classical neurotoxicants with diverse mechanisms of toxicities [16]. Similar refinements were obtained by other groups

[13,29,32]. There were also EST-inspired tests based on the differentiation of mESCs into endothelial cells [33] and osteoblasts [14,31], which showed to be effective and responsive to novel embryotoxic compounds [33], strongly validating the choice of using additional differentiation protocols in the EST.

Currently, mESC-based toxicity assays represent a powerful tool for the high-throughput toxicity screening of emerging organic chemicals [25], for chemical use and exposure guideline purposes. For instance, the US Environmental Protection Agency (EPA) and others developed a toxicity assay based on the monolayer differentiation of mESCs into cardiomyocytes, called mESC adherent cell differentiation and cytotoxicity (ACDC) assay [34], enabling the high-throughput screening of environmental pollutants, which also gives information on the underlying molecular mechanisms [35].

Human stem cell-based toxicology

Both mESCs and human ESCs (hESCs) are necessary for stem cell toxicology research, but with different roles. In the perspective of stem cell biology, hESCs differ from mESCs in not only their morphology and expression of surface antigens but also gene markers and cytokines [3,36,37]. In addition, when derived from the inner cell mass of the blastocyst, mESCs are generally in the naive state (the highest pluripotency level) whereas hESCs are in the primed state (a lower pluripotency level) [38]. Thus, a rational corollary of these differences is that in toxicological terms, hESCs may be more appropriate than mESCs in terms of the representation of toxicity effects toward humans. In addition, hESC-based toxicity tests may be more reliable and indicative, since they avoid interspecies differences. Studies have found out hESCs were more sensitive to the adverse effects of a broad range of chemicals [39–42], and could be impaired by certain chemicals that were not necessarily toxic to mice, such as thalidomide [43], or chemicals, including ethanol and caffeine, which have raised health concerns [44,45].

The first hESC-based toxicity tests took advantage of the neuronal induction from hESCs. Exposure of hESC-derived neurons to specific dopaminergic neurotoxins resulted in apoptotic cell death, production of reactive oxygen species, and loss of neuronal functions [46]. Two years after this study, the European Center for the Validation of Alternative Methods (ECVAM) proposed the establishment of the human EST [47]. In this case, H1 hESCs and hMRC-5 embryonic lung fibroblasts were exposed to two well-known developmental toxicants, all-trans retinoic acid (RA) and 5-fluorouracil (5-FU), for either 4 or 10 days. MTT cytotoxicity tests and gene expression analyses for cardiogenic markers demonstrated that the toxic responses caused by RA and 5-FU were compatible with previous *in vivo* data. These two studies demonstrated already that hESCs could be employed in cytotoxicity, developmental toxicity, and functional toxicity assays.

To assess the vulnerability of the human nervous and cardiovascular systems, neonatally and postnatally, hESC- or human iPSC (hiPSC)-based toxicity research has been embraced with a variety of artificial chemicals, such as bisphenol A, perfluorooctane sulfonate, trichloroethylene, trimethyltin chloride, ibuprofen, and paraquat [48–50], as well as more complex materials including silver nanoparticles [51,52] and gold nanoparticles [53], for their effects on cardiogenesis and neural development.

TABLE 1. PRIMARY REFINEMENTS OF THE EMBRYONIC STEM CELL TEST

<i>Chemicals</i>	<i>Cell line</i>	<i>Exposure duration</i>	<i>Differences from Spielmann's original EST</i>	<i>Findings</i>	<i>Reference</i>
15 pharmaceuticals	DBA/11acJ murine derived ESCs, 3T3 A31	Cytotoxicity assay: 7 days Differentiation assay: 10 days	Different ESC line	Other murine ESC lines can replace D3 in the EST; false positives and negatives can still be significant	7
Penicillin G, 5-FU, RA, diphenylhydantoin, valproic acid, thalidomide	Balb/c 3T3, D3	Cell viability: 10 days Differentiation assays (contraction assay and collection of RNA samples): 10 days	Involvement of multiple marker genes of other lineages (neurogenesis, osteogenesis, and chondrogenesis markers)	First time to include multiple gene expressions in the EST; results are promising and open to further improvement and optimization	8
6-Aminonicotinamide, 5-FU, methylmercury chloride, hydroxyl urea, valproic acid, boric acid, methoxyacetic acid, lithium chloride, penicillin G, saccharin, diphenhydramine, acrylamide	D3	Cell viability assay: 6 days Differentiation assay: 10 days	Instead of being seeded in 24-well plates, mESCs are seeded and cultured in low attachment 96-well plate to facilitate high-throughput tests. Introduce relative embryotoxicity potency values that are derived from the ranking of the embryotoxic potential of test compounds relative to positive controls	No need for 3T3 to provide reference toxicity information. Toxicity ranking seems more reliable than the absolute classification of toxicity	9
Monobutyl phthalate, methoxyacetic acid, valproic acid, RA, 5-FU, penicillin G	D3	Cell viability assay: 3 days Differentiation assay: 10 days	Use resazurin proliferation assay; involve whole-genome gene expression profiling	Several genes are very sensitive and significantly dysregulated after compound exposure, which are useful as additional endpoints in the EST	10
Triazoles, flusilazole, hexaconazole, cyproconazole, myclobutaniol, triadimefon, triticonazole, 5-FU	D3	Cell viability assay: 5 days Differentiation assay: 10 days	Compare the EST with in vivo tests, the rat postimplantation WEC, and zebrafish embryotoxicity test	The EST is able to represent developmental toxicity, and partially mimic the processes in utero	11
9-cis-RA, 13-cis-RA, acitretin, ATRA, TTNPB, etretinate, retinol	D3	Cell viability: 10 days Cell differentiation: 10 days	Compare the EST with in vivo tests, the rat limb bud micromass test, and the postimplantation rat WEC	The EST, although it could not completely imitate toxicological kinetics in vivo, was overall a promising in vitro alternative in toxicology	12
Methylmercury, valproic acid, AsV, AsIII, saccharin, isoniazid, ascorbic acid	D3, 3T3	Cytotoxicity: 10 days Neural differentiation assay: 12 days	Neural differentiation instead of cardiac induction with exposure duration of 12 days. Use flow cytometry and detection of neural gene expression as endpoints	Neural differentiation could serve as an alternative process in the EST. The molecular endpoints introduced in the neural differentiation-based EST were <i>Tuj1</i> mRNA and protein levels	13

(continued)

TABLE 1. (CONTINUED)

<i>Chemicals</i>	<i>Cell line</i>	<i>Exposure duration</i>	<i>Differences from Spielmann's original EST</i>	<i>Findings</i>	<i>Reference</i>
Mono-ethylhexyl phthalate, valproic acid, methotrexate, 6-aminonicotinamide, methoxyacetic acid, penicillin G	D3	Cell viability: 5 days Osteoblast differentiation assay: 21 days. Cardiac differentiation assay: 10 days	Involvement of osteoblast differentiation and molecular endpoints to evaluate it. Comparison between osteoblast and cardiomyocyte differentiations on exposure to same chemicals	Osteoblast differentiation can be an alternative to cardiogenesis in the EST, and may give different results	14
Phenol, p-fluorophenol, p-heptyloxyphenol, p-mercaptophenol, p-methylketophenol	D3	Cell differentiation assay: 10 days	Compare the EST with <i>in vivo</i> tests and the WEC assay	The EST gives toxicity rankings of tested phenols that are different from the rankings given by <i>in vivo</i> tests and the WEC assay; exposure doses in the EST have to consider the kinetics of <i>in vivo</i> absorption, metabolism, elimination, and excretion	15
Acetaldehyde, carbamazepine, flusilazole, monoethylhexylphthalate, penicillin G sodium salt, phenytoin	D3	Cell viability test: 48 h Morphological scoring: 72 h Whole-genome expression profiling: 24 h	Neural differentiation Different exposure durations Resazurin cell viability assay Include genome profiling	The neural differentiation-modified EST is valid; transcriptomics provides mechanistic information	16
MeHgCl, monosodium L-glutamate, penicillin G, poly-L-ornithine, sodium arsenite, sodium valproate, chlorpyrifos-ethyl, parathion-ethyl	D3	Cell viability: 4 or 5 days Differentiation: 2 or 3 days	Differentiation to neural cells Cell proliferation tests are based on ELISA. Cell viability tests are based on CellTiter-Blue Cell Viability Assay. Involvement of β III-Tubulin enzyme-linked immunosorbent assay	This method is suitable for high-throughput screening but does not necessarily represent relevant concentrations <i>in vivo</i> and is not applicable for acute and chronic toxicities	17
Bisphenol A, genistein, as well as combined with bisphenol A and 5-FU	D3, 3T3	Cell viability test: 10 days Differentiation assay: 10 days	Cell Titer 96 Aqueous One Solution Cell Proliferation Assay for cell viability test; cells are exposed to two chemicals	Bisphenol A and genistein, to which we are exposed daily unintentionally, have combined embryotoxic effects that become synergistic at low concentrations	18
38 teratogens 39 nonteratogens	D3	Cell viability test: 72 h Cell differentiation assay: 96 h	Shorter exposure times; include gene expression analysis for 12 potential molecular endpoints	The Molecular Embryonic Stem Cell Developmental Toxicity Assay facilitates high-throughput screenings of potential teratogens with good predictivity and concordance with <i>in vivo</i> data	19

(continued)

TABLE 1. (CONTINUED)

Chemicals	Cell line	Exposure duration	Differences from Spielmann's original EST	Findings	Reference
Monobutyl phthalate, monobenzyl phthalate, mono-(2-ethylhexyl) phthalate, monomethyl phthalate	D3	Cell viability: 5 days Differentiation assay: 10 days	Incorporate RNA microarray analyses as additional endpoints	A total of 668 commonly expressed genes are altered after exposure, proving the validity of transcriptomics in the EST	20
5-FU, hydroxyurea, saccharin; silver nanomaterial, coated and uncoated zinc oxide, titanium and silica nanomaterials	D3, 3T3	Cell viability: 10 days Cell differentiation: 10 days	Skip the step of EB formation in petri dishes and transfer EBs directly to 24-well plates. Add nanomaterial once to avoid continuous accumulation in cells	This simplified protocol shows to be more suitable to facilitate nanotoxicity research for medical or therapeutic nanomaterial uses	21
6-aminonicotinamide, all-trans RA, 5-bromo-2'-deoxyuridine, dexamethasone, methoxyacetic acid, salicylic acid sodium salt, ascorbic acid, acrylamide, D-(+)-camphor, 5-FU	Linearized Hand1-promoter-Luc plasmid transfected C57BL/6 mice derived ESCs	Cell viability: 5 days Differentiation assay: 5 days	Monitor <i>Hand1</i> expression via Luciferase reporter assay, which at the same time indicates both proliferation and differentiation	The expression of <i>Hand1</i> by Luciferase reporter gene assay is reproducible and relatively accurate	22
Simvastatin	D3, 3T3	Cytotoxicity: 10 days. Differentiation assay (both hanging drop method and monolayer differentiation): 10 days	Include both EB hanging drop method and monolayer differentiation. Molecular endpoints are maker genes for each germ layer	Genes of the mesodermal lineage are most sensitive to the two drugs; the hanging drop method and monolayer differentiation give rise to consistent results	23
Chinese herbal extracts from <i>Atractylodes macrocephala</i> , <i>Coptis chinensis</i> , <i>Radix isatidis</i> , <i>Flos genkwa</i>	OG2 mESCs, BALB/c 3T3	Cell viability: 10 days Differentiation assay: 10 days	Cell viability assay: CCK8 assay Differentiation assay based on myosin heavy chain gene expression	<i>Atractylodes macrocephala</i> and <i>Radix isatidis</i> are non-embryotoxic, <i>Coptis chinensis</i> is weakly embryotoxic whereas <i>Flos genkwa</i> is strongly embryotoxic	24
5-FU, RA, valproic acid, diphenhydramine, LiCl, saccharin, penicillin G	D3, 3T3	Cell viability: 5 days Differentiation assay: 5 days	Based on only monolayer culture with 5-day exposure. Check 16 genes for the three germ layers as endpoints for differentiation	Monolayer culture is applicable in the EST with gene expression detection for three germ layers	25

ATRA, all trans retinoic acid; EB, embryoid body; ESC, embryonic stem cell; EST, embryonic stem cell test; 5-FU, 5-fluorouracil; mESCs, mouse embryonic stem cells; RA, retinoic acid; TTNPB, (E)-4[2-(5,6,7,8-tetrahydro-5,5,8,8-tetramethyl-2-naphthalenyl)-1-propenyl]benzoic acid; WEC, whole embryo culture.

In summary, stem cell toxicology offers great potential and is open to improvements and refinements in parallel with the development of stem cell biology. It would allow for toxicity tests to be designed for many different kinds of cells, with the ultimate goal of obtaining trustable toxic information relative to human health.

Advantages of Stem Cell-Based Toxicity Tests

The unique properties of stem cells qualify them as one of the best materials for *in vitro* tests. Unlike primary cells, stem cells, such as ESCs and iPSCs, are able to proliferate extensively *in vitro* under defined conditions. Theoretically, ESCs and iPSCs can be passaged for indefinite generations without transformation, as established cell lines. Conversely, cancer or immortalized cell lines, although they can be also maintained *in vitro* indefinitely, are prone to queries as to whether or not they represent healthy or normal cells and are sensitive enough in toxicity assays. Stem cell-based toxicity tests avoid these issues and, consequently, stand out as promising alternatives to animal tests owing to their recapitulating of many key *in vivo* features.

Stem cells also hold the capacity to differentiate into other cell types. In other words, with proper differentiation conditions, it is possible to generate cells of all the different tissues of the body, including germ cells, in particular when starting from ESCs. These cells can then be used in toxicity studies without *ex vivo* derivation or ethical concerns. For example, hESCs are able to differentiate into cardiomyocytes, allowing for toxicity assays either during the process of cardiac differentiation or with the terminally differentiated cells [54], therefore permitting developmental and functional toxicity tests, respectively. Thus, stem cells offer the distinct possibility to study cellular, developmental, functional, and reproductive toxicities with various types of stem cells and stem cell-derived offspring, particularly when primary cell types are impossible to be derived directly.

Consistency between stem cell toxicity tests and live animal assays is necessary before concluding that by using stem cells in toxicity assays we can get rid of animal experiments altogether [55]. Summarized in Table 1 are the consistency assessments between EST-based stem cell toxicity tests and *in vivo* or other *in vitro* tests, such as the WEC (whole embryo culture) [11,12,15]. They revealed a general consistency between *in vivo* and *in vitro* results, advocating for the use of stem cell toxicity tests that would also allow for tests with multiple chemicals, doses, and durations at one time, with considerably less consumption of time, money, and labor.

Therefore, the advantages of stem cells mentioned earlier would help face the tremendous challenge of the numerous environmental pollutants, drugs, and industrial chemicals with vague toxicology information [11,39], by employing reliable and high-throughput experimental systems that will eventually contribute to the rise of stem cell toxicology [56].

Major Achievements and Future Prospects of Human Stem Cell Toxicology

In this section, we focus only on the growing number of research studies recently reported for hESCs and hiPSCs related teratogen and drug screenings, with an emphasis on their potentiality for more comprehensive toxicity assays.

hESC systems as powerful tools for developmental toxicity research

hESCs are pluripotent, as they have the capacity to differentiate into cells of all three germ layers, as well as germ cells. Thus, toxicity assays with hESCs can provide toxicity information at very early stages after fertilization, since differentiating ESCs mimic the dynamic process of embryonic development. Besides, hESC-derived cells can also serve as a useful alternative to study toxicity effects on a variety of somatic cells.

Cytotoxicity assays on hESCs involve the use of the chemical(s) of interest, to study not only cell death but also whether pluripotency of hESCs would be affected by the chemicals. For example, the embryotoxic compound 5-FU downregulated the expression of *Oct4* and *Nanog* in hESCs, as well as significantly dysregulated *HDAC9* (involved in neuronal, skeletal muscle and adipocyte differentiation), *DLK1* (involved in embryogenesis and differentiation of MSCs to chondrocytes), and *NFE2L3* (involved in differentiation, inflammation, and carcinogenesis) [57], indicating that early exposure to 5-FU may result in multiple malformations. A similar study focused on how teratogens affected the hESC pluripotent state by observing changes in gene networks, and it proposed the explanation that adverse effects occurring during embryonic development might be originated from exposure during very early stages, even before ESCs in the blastocyst start differentiating [58]. Therefore, such a study helps understand how chemicals can influence the process of lineage fate decision in pluripotent cells.

Developmental toxicity tests with the hESC system benefit from the process of embryoid body (EB) formation, which recapitulates the key features of gastrulation [59]. Therefore, toxicity tests during EB formation yield information about which germ layer(s) will be affected and, therefore, which lineage differentiation will be impaired. For example, several natural and synthetic estrogens that readily cross the placenta, such as estradiol, estriol, and bisphenol A, negatively affected endodermal, mesodermal, and ectodermal differentiation during EB formation [49]. Moreover, a common component in toothpaste, fluoride, which had been shown to affect murine and human early embryogenesis, suppressed hESC proliferation and induced apoptosis, but also stimulated the generation of ectoderm and mesoderm at the expenses of endoderm, when administered at high doses [60]. Studies like these have drawn attention to the lack of toxicity information for emerging chemicals that have been widely used.

Toxicity assays for a variety of pharmaceuticals, materials, and environmental pollutants with hESC-derived somatic cells have been performed and validated, implying a great potentiality for these systems in functional toxicology. Cardiomyocytes are one of the earliest cell types being successfully derived from hESCs and have been already used in toxicity tests, for a few common drugs [54,61–63], and in a high-throughput screening for teratogens [8,44]. For instance, the chemotherapeutic agent doxorubicin caused acute cytotoxicity in hESC-derived cardiomyocytes, due to membrane damage [61,62]. Further, trichloroethylene, a ubiquitously detected industrial chemical, disrupted the transition from cardiac progenitor cells to cardiomyocytes [48]. Neural differentiation of hESCs, as another well-established *in vitro* process, has been introduced in stem cell toxicology [64,65]. Several studies

showed that a variety of chemicals or materials might impair the development of different types of neural cells [51,53,66,67], and they also deciphered the underlying mechanisms. For example, gold nanoparticles were toxic to neural hESC derivatives, due to the alteration of DNA methylation and hydroxymethylation levels during neural differentiation [53]. In addition, an immediate and significant increase in oxidative stress caused by silver nanoparticle on hESC-derived neural progenitor cells was due to dysregulation at both transcriptional and post-transcriptional levels [52].

There are also studies on more complex neural cell populations, such as neural rosette [68,69] and artificial neural constructs [50,70]. For example, highly consistent responses toward RA between in vitro neural rosette and in vivo systems confirmed the potentiality of neural rosette generation for developmental neurotoxicity tests [69]. One promising study focused on the generation of neural constructs, comprising different kinds of cells, in 3D structures mimicking the developing brain. That study obtained highly accurate toxicity results with known neural toxicants and controls [70], and together with another toxicity study based on 3D neural spheres [50], demonstrated the potential of 3D cell constructs or even organoids in toxicology research.

With the establishment of additional hESC-based differentiation protocols, other hESC-based models in toxicology have been validated, and several different kinds of toxicities from a variety of chemicals, materials, and environmental pollutants have been measured. For example, hESC-derived germ cell-like cells have been utilized in reproductive toxicity tests that facilitated and accelerated the screening of the remarkable number of potential endocrine disruptors to which we are continuously exposed [71]. A report on the perturbation of low-dose bisphenol A on the direct differentiation of hESCs into prostate organoids suggested that prenatal exposure could perturb the morphogenesis of the prostate [72].

Besides, systems biology methodologies, such as transcriptomics, proteomics, and metabolomics, have been also incorporated into stem cell toxicology research. For example, one group exposed hESCs to several teratogenic drugs and detected small molecules whose levels were significantly altered in response to those drugs, proving the convenience of metabolomics in toxicity assays [73]. A follow-up toxicity study based on metabolomics achieved an 83% predictive accuracy [74]. Another investigation using ornithine and cysteine, as indicators of developmental toxicity, identified novel potential teratogens with 77% accuracy [75], providing evidence for metabolomics effectiveness as an alternative method in toxicology. Optimally, integrating metabolomics, transcriptomics, and proteomics assays in stem cell toxicology would guarantee a more comprehensive analysis of the toxicity mechanisms.

In conclusion, from a toxicological point of view, hESCs have great potential. Novel and versatile toxicity screening platforms based on hESCs are advocated to advance the field of in vitro toxicology.

hiPSCs as promising tools for personalized toxicology analyses

iPSCs have drawn attention and raised expectations in regenerative medicines, since their discovery. They are

equivalent to ESCs in potency, although they are not derived from embryos. In fact, they are generated by reprogramming somatic cells, generally obtained from an adult individual, into the pluripotent state. Thus, this method could provide a large quantity of patient-specific pluripotent cells with no immune response issues for potential cell base transplantations and other clinical applications. Moreover, hiPSCs come with almost no ethical concerns compared with the derivation and use of hESCs. Besides regenerative medicine, iPSCs can potentially be employed for drug development and toxicology studies, as previously reviewed [76–79]. A study comparing the neurotoxicity effects of KG-501, an inhibitor of the CREB pathway, on hESC- and hiPSC-derived neurons and glia, verified the potentialities of hiPSCs in neural developmental toxicity assays [80]. Cardiomyocyte and hepatogenic differentiation from hiPSC has been also employed for drug screenings. Table 2 summarizes the findings in several drug screenings based on iPSC-derived cardiomyocytes and hepatocytes, proving the robustness of these models.

Reprogramming enables generation of iPSCs from patients as well as healthy people. In this respect, toxicology research can target different populations. For example, iPSC-derived cardiomyocytes from patients whose genetic backgrounds made them more prone to cardiac diseases resulted in more sensitivity to perturbations by cardio-toxic drugs, whereas cardiomyocytes from healthy people were more tolerant [85]. One investigation compared doxorubicin-induced toxicity in iPSC-derived cardiomyocytes from breast cancer patients who either suffered from doxorubicin-induced cardiotoxicity or did not. Results showed that cardiomyocytes derived from patients who experienced doxorubicin-induced cardiotoxicity were more sensitive to doxorubicin [89].

Further, iPSC-derived neural cells provide great research material, especially for patient-specific neurotoxicity information. Efforts on modeling neurotoxicity with iPSC-derived neural cell types are emphasized in one study about chemotherapy-induced toxicity [90], as well as in two others on embracing hiPSC neural differentiation in high-throughput screenings [91,92]. Those examples also clearly demonstrate that hiPSC-based assays are open to improvement and adaptation to personalized toxicology analysis for diagnosis and therapies.

Taking everything together, hiPSC applications in medical research and drug screenings support the idea of employing them as well as iPSC-derived cells in stem cell toxicology. In fact, those methods offer the possibility of a more personalized toxicology in that they allow to investigate whether patients or specific groups of healthy people with different genetic backgrounds respond differently to toxicants. However, there have not yet been a lot of persuasive reports on toxicity assays of environmental chemicals based on hiPSCs [93].

Conclusive Remarks

Currently, toxicology still heavily relies on live animal tests, which are prone to errors due to interspecies inconsistencies [94] when we applied them to predict human health risks. Unlike drug screenings, which always include clinical trials with human volunteers, toxicology tests cannot rely on human subjects directly. In addition, we are continuously exposed to a variety of different chemicals,

TABLE 2. DRUG SCREENING AND TOXICITY ASSAYS BASED ON HUMAN-INDUCED PLURIPOTENT STEM CELL-DERIVED CARDIOMYOCYTES OR HEPATOCYTES

<i>Differentiation lineage</i>	<i>Drug</i>	<i>Findings</i>	<i>Reference</i>
Hepatocyte	Amiodarone Alatoxin B1 Troglitazone Ximelagatran	Amiodarone, Alatoxin B1 and Troglitazone cause dose-dependent toxic effects on hiPSC-derived hepatocytes after 14 days of exposure. There is specific induction of phospholipidosis and steatosis after exposure. hiPSC-derived hepatocytes show time-dependent toxicity effects from exposure to the drugs, suggesting that the model is applicable in long-term toxicity test and that certain drugs may have accumulative toxic effects on the liver	62
Hepatocyte	24 drugs	3D spheroid culture of hepatocytes not only produces high efficiencies but also allows cells to grow with better oxygen conditions, which gives rise to more functional hepatocytes with enhanced performance in drug screenings	81
Hepatocyte	238 compounds from the Screen-Well™ Hepatotoxicity Library (Enzo Life Sciences)	High-content automated screening assays based on hiPSC-derived hepatocytes provide information on cell viability, nuclear shape and intensity, cytoskeleton integrity, mitochondrial potential, autophagy, and lipid accumulation by different staining methods	82
Hepatocyte	Staurosporine Acetaminophen	hiPSC-derived hepatocytes are highly similar to human primary hepatocytes in regulating and executing apoptosis after drug exposure, compared with two other human liver cell lines, suggesting that hiPSC-derived hepatocytes are a good alternative to primary hepatocytes	83
Cardiomyocyte	24 drugs	hiPSC-derived cardiomyocytes show impedance changes after drug exposure.	84
Cardiomyocyte	Cisapride Nicorandil Alfuzosin Verapamil	Although iPSC-derived cardiomyocytes are not as mature as cardiomyocytes derived from the human body, based on their gene expression, iPSC-derived cardiomyocytes can be utilized in drug screenings, especially suitable for toxicity tests for high-risk populations	85
Cardiomyocyte	131 drugs consisting of both cardiotoxic and cardio-safe ones	384-well plate-based toxicity assays with good accuracy in classification of cardiotoxicity of drugs	86
Cardiomyocyte	Ponatinib	Ponatinib induces cell death, troponin secretion, and reactive oxygen species and lipid formation, inhibits ABL activation and survival pathways, and disrupts actin cytoskeleton structures and beating of hiPSC-derived cardiomyocytes, at medical significant doses	87
Cardiomyocyte	23 drugs	Toxicity endpoints for hiPSC-derived cardiomyocytes cell viability and function provide robust evaluation of drug toxicity	88
Cardiomyocyte	Doxorubicin	iPSC-derived cardiomyocytes from patients experiencing doxorubicin-induced cardiotoxicity are more sensitive than the ones from patients who do not experience cardiotoxicity	89

hiPSC, human-induced pluripotent stem cell.

including emerging environmental pollutants, which present potentially high risk for our health. Thus, reliable toxicity models based on the human physiology are needed to assess the adverse effects of all these chemicals on humans [95].

More than 20 years of exploration of ESC-based in vitro toxicity assays confirmed the importance of stem cells in toxicology research. Human PSCs enable toxicologists to study developmental toxicity, as well as functional toxicities on mature terminally differentiated cell types that are difficult to be derived directly from humans. The development of human stem cell toxicology benefited from all the efforts in the design and improvement of the mouse EST, which was the first stem cell-based developmental toxicity test to eliminate the use of pregnant animals. Nowadays, hPSC-derived

cardiomyocytes, neural cells, and hepatocytes are major lineages being employed. However, more stem cell-derived lineages, such as intestinal, gastric, and pulmonary cells, should be implemented in developmental and functional toxicity tests. In particular, differentiation of hESCs/iPSCs into primordial germ cells and gametes would broaden the scope of the in vitro reproductive toxicity investigations. In addition, advances in the synthesis of biomaterials, and conditions for 3D culture and differentiation of hESCs, have allowed the generation of more complex cell constructs, which provide a dynamic cell system resembling the in vivo situation more comprehensively than 2D systems. This would also make toxicology more predictive. The breakthrough of the generation of hiPSCs has also permitted drug screenings and toxicity

assays with cells derived from patients, building the foundations for a more personalized toxicology.

In conclusion, stem cell toxicology, although still developing, has already been proved to be a very comprehensive toxicity system, which allows for the assessment of embryonic, developmental, organ, reproductive, and functional toxicities. When fully evolved, it will definitely represent the gold standard for in vitro toxicology.

Acknowledgments

The authors were funded by grants from the Chinese Academy of Sciences Strategic Priority Research Program (XDB14040301), the National Natural Science Foundation of China (21577166), the Chinese Academy of Sciences Hundred Talent Program (29[2015]30), and the Key Research Program of Frontier Sciences, CAS (QYZDJ-SSW-DQC017).

Author Disclosure Statement

No competing financial interests exist.

References

- Russel WMS and RL Burch, ed. (1959). *The Principles of Humane Experimental Technique*. Methuen & Co., London. Special edition published by Universities Federation for Animal Welfare (UFAW), 1992.
- Schwetz BA, RE Morrissey, F Welsch and RA Kavlock. (1991). In vitro teratology. *Environ Health Perspect* 94: 265–268.
- Evans MJ and MH Kaufman. (1981). Establishment in culture of pluripotential cells from mouse embryos. *Nature* 292:154–156.
- Takahashi K and S Yamanaka. (2006). Induction of pluripotent stem cells from mouse embryonic and adult fibroblast cultures by defined factors. *Cell* 126:663–676.
- Laschinski G, R Vogel and H Spielmann. (1991). Cytotoxicity test using blastocyst-derived euploid embryonic stem cells—a new approach to in vitro teratogenesis screening. *Reprod Toxicol* 5:57–64.
- Spielmann H, I Pohl, B Döring, M Liebsch and F Moldenhauer. (1997). The embryonic stem cell test, an in vitro embryotoxicity test using two permanent mouse cell lines: 3T3 fibroblasts and embryonic stem cells. *In Vitro Toxicol* 10:119–127.
- Paquette JA, SW Kumpf, RD Streck, JJ Thomson, RE Chapin and DB Stedman. (2008). Assessment of the embryonic stem cell test and application and use in the pharmaceutical industry. *Birth Defects Res B Dev Reprod Toxicol* 83:104–111.
- zur Nieden NI, G Lempka and HJ Ahr. (2004). Molecular multiple endpoint embryonic stem cell test—a possible approach to test for the teratogenic potential of compounds. *Toxicol Appl Pharmacol* 194:257–269.
- Peters AK, M Steemans, E Hansen, N Mesens, GR Verheyen and P Vanparys. (2008). Evaluation of the embryotoxic potency of compounds in a newly revised high throughput embryonic stem cell test. *Toxicol Sci* 105:342–350.
- van Dartel DA, JL Pennings, LJ de la Fonteyne, MH van Herwijnen, JH van Delft, FJ van Schooten and AH Piersma. (2010). Monitoring developmental toxicity in the embryonic stem cell test using differential gene expression of differentiation-related genes. *Toxicol Sci* 116:130–139.
- de Jong E, M Barenys, SA Hermsen, A Verhoef, BC Ossendorp, JG Bessems and AH Piersma. (2011). Comparison of the mouse embryonic stem cell test, the rat whole embryo culture and the zebrafish embryotoxicity test as alternative methods for developmental toxicity testing for six 1,2,4-triazoles. *Toxicol Appl Pharmacol* 253:103–111.
- Louisse J, S Gönen, IMCM Rietjens and M Verwei. (2011). Relative developmental toxicity potencies of retinoids in the embryonic stem cell test compared with their relative potencies in in vivo and two other in vitro assays for developmental toxicity. *Toxicol Lett* 203:1–8.
- Baek DH, TG Kim, HK Lim, JW Kang, SK Seong, SE Choi, SY Lim, SH Park, BH Nam, et al. (2011). Embryotoxicity assessment of developmental neurotoxicants using a neuronal endpoint in the embryonic stem cell test. *J Appl Toxicol* 32:617–626.
- de Jong E, L van Beek and AH Piersma. (2012). Osteoblast differentiation of murine embryonic stem cells as a model to study the embryotoxic effect of compounds. *Toxicol In Vitro* 26:970–978.
- Strikwold M, RA Woutersen, B Spenkeliink, A Punt and IMC Rietjens. (2012). Relative embryotoxic potency of p-substituted phenols in the embryonic stem cell test (EST) and comparison to their toxic potency in vivo and in the whole embryo culture (WEC) assay. *Toxicol Lett* 213:235–242.
- Theunissen PT, JF Robinson, JL Pennings, MH van Herwijnen, JC Kleinjans and AH Piersma. (2012). Compound-specific effects of diverse neurodevelopmental toxicants on global gene expression in the neural embryonic stem cell test (ESTn). *Toxicol Appl Pharmacol* 262:330–340.
- Hayess K, C Riebeling, R Pirow, M Steinfath, D Sittner, B Slawik, A Luch and AE Seiler. (2013). The DNT-EST: a predictive embryonic stem cell-based assay for developmental neurotoxicity testing in vitro. *Toxicology* 314: 135–147.
- Kong D, L Xing, R Liu, J Jiang, W Wang, L Shang, X Wei and W Hao. (2013). Individual and combined developmental toxicity assessment of bisphenol A and genistein using the embryonic stem cell test in vitro. *Food Chem Toxicol* 60:497–505.
- Panzica-Kelly JM, KC Brannen, Y Ma, CX Zhang, OP Flint, LD Lehman-McKeeman and KA Augustine-Rauch. (2013). Establishment of a molecular embryonic stem cell developmental toxicity assay. *Toxicol Sci* 131:447–457.
- Schulpen SHW, JF Robinson, JLA Pennings, DAM van Dartel and AH Piersma. (2013). Dose response analysis of monophthalates in the murine embryonic stem cell test assessed by cardiomyocyte differentiation and gene expression. *Reprod Toxicol* 35:81–88.
- Corradi S, E Dakou, A Yadav, LC Thomassen, M Kirsch-Volders and L Leyns. (2015). Morphological observation of embryoid bodies completes the in vitro evaluation of nanomaterial embryotoxicity in the embryonic stem cell test (EST). *Toxicol In Vitro* 29:1587–1596.
- le Coz F, N Suzuki, H Nagahori, T Omori and K Saito. (2015). Hand1-Luc embryonic stem cell test (Hand1-Luc EST): a novel rapid and highly reproducible in vitro test for embryotoxicity by measuring cytotoxicity and differentiation toxicity using engineered mouse ES cells. *J Toxicol Sci* 40:251–261.
- Jyoti S and S Tandon. (2015). Genetic basis for developmental toxicity due to statin intake using embryonic stem cell differentiation model. *Hum Exp Toxicol* 34:965–984.

24. Li LY, FF Cao, ZJ Su, QH Zhang, XY Dai, X Xiao, YD Huang, Q Zheng and H Xu. (2015). Assessment of the embryotoxicity for four Chinese herbal extracts using embryonic stem cell test. *Mol Med Rep* 12:2348–2354.
25. Romero AC, E del Rio, E Vilanova and MA Sogorb. (2015). RNA transcripts for the quantification of differentiation allow marked improvements in the performance of embryonic stem cell test (EST). *Toxicol Lett* 238:60–69.
26. Seiler AEM and H Spielmann. (2011). The validated embryonic stem cell test to predict embryotoxicity in vitro. *Nat Protoc* 6:961–978.
27. Seiler A, A Visan, R Buesen, E Genschow and H Spielmann. (2004). Improvement of an in vitro stem cell assay for developmental toxicity: the use of molecular endpoints in the embryonic stem cell test. *Reprod Toxicol* 18:231–240.
28. Buesen R, E Genschow, B Slawik, A Visan, H Spielmann, A Luch and A Seiler. (2009). Embryonic stem cell test remastered: comparison between the validated EST and the new molecular FACS-EST for assessing developmental toxicity in vitro. *Toxicol Sci* 108:389–400.
29. Pennings JLA, DAM van Dartel, JF Robinson, TE Pronk and AH Piersma. (2012). Gene set assembly for quantitative prediction of developmental toxicity in the embryonic stem cell test. *Toxicology* 284:63–71.
30. van Dartel DA, JL Pennings, LJ de la Fonteyne, KJ Brauers, S Claessen, JH van Delft, JC Kleinjans and AH Piersma. (2011). Concentration-dependent gene expression responses to flusilazole in embryonic stem cell differentiation cultures. *Toxicol Appl Pharmacol* 251:110–118.
31. Chen X, DK Hansen, G Merry, C DeJarnette, G Nolen, D Sloper, JE Fisher, W Harrouk, MS Tassinari and AL Inselman. (2015). Developing osteoblasts as an endpoint for the mouse embryonic stem cell test. *Reprod Toxicol* 53: 131–140.
32. van Dartel DA, SH Schulpen, PT Theunissen, A Bunschoten, AH Piersma and J Keijer. (2014). Dynamic changes in energy metabolism upon embryonic stem cell differentiation support developmental toxicant identification. *Toxicology* 324:76–87.
33. Festag M, B Viertel, P Steinberg and C Sehner. (2007). An in vitro embryotoxicity assay based on the disturbance of the differentiation of murine embryonic stem cells into endothelial cells. II. Testing of compounds. *Toxicol In Vitro* 21:1631–1640.
34. Barrier M, S Jeffay, HP Nichols, KJ Chandler, MR Hoopes, K Slentz-Kesler and ES 3rd Hunter. (2011). Mouse embryonic stem cell adherent cell differentiation and cytotoxicity (ACDC) assay. *Reprod Toxicol* 31:383–391.
35. Chandler KJ, M Barrier, S Jeffay, HP Nichols, NC Kleinstreuer, AV Singh, DM Reif, NS Sipes, RS Judson, et al. (2011). Evaluation of 309 environmental chemicals using a mouse embryonic stem cell adherent cell differentiation and cytotoxicity assay. *PLoS One* 6:e18540.
36. Thomson JA, J Itskovitz-Eldor, SS Shapiro, MA Waknitz, JJ Swiergiel, VS Marshall and JM Jones. (1998). Embryonic stem cell lines derived from human blastocysts. *Science* 282:1145–1147.
37. Ginis I, Y Luo, T Miura, S Thies, R Brandenberger, S Gerecht-Nir, M Amit, A Hoke, MK Carpenter, J Itskovitz-Eldor and MS Rao. (2004). Differences between human and mouse embryonic stem cells. *Dev Biol* 269:360–380.
38. Hanna J, AW Cheng, K Saha, J Kim, CJ Lengner, F Soldner, JP Cassidy, J Muffat, BW Carey and R Jaenisch. (2010). Human embryonic stem cells with biological and epigenetic characteristics similar to those of mouse ESCs. *Proc Natl Acad Sci (USA)* 107:9222–9227.
39. Flora SJS and A Mehta. (2009). Monoisoamyl dimercaptosuccinic acid abrogates arsenic-induced developmental toxicity in human embryonic stem cell-derived embryoid bodies: comparison with in vivo studies. *Biochem Pharmacol* 78:1340–1349.
40. He X, S Imanishi, H Sone, R Nagano, XY Qin, J Yoshinaga, H Akanuma, J Yamane, W Fujibuchi and S Ohsako. (2012). Effects of methylmercury exposure on neuronal differentiation of mouse and human embryonic stem cells. *Toxicol Lett* 212:1–10.
41. Lin S, V Tran and P Talbot. (2009). Comparison of toxicity of smoke from traditional and harm-reduction cigarettes using mouse embryonic stem cells as a novel model for preimplantation development. *Hum Reprod* 24:386–397.
42. Lin S, S Fonteno, J Weng and P Talbot. (2010). Comparison of the toxicity of smoke from conventional and harm reduction cigarettes using human embryonic stem cells. *Toxicol Sci* 118:202–212.
43. Meganathan K, S Jagtap, V Wagn, J Winkler, JA Gaspar, D Hildebrand, M Trusch, K Lehamann, J Hescheler, H Schlüter and A Sachinidis. (2012). Identification of thalidomide-specific transcriptomics and proteomics signatures during differentiation of human embryonic stem cells. *PLoS One* 7:e44228.
44. Mayshar Y, O Yanuka and N Benvenisty. (2011). Teratogen screening using transcriptome profiling of differentiating human embryonic stem cells. *J Cell Mol Med* 15: 1393–1401.
45. Nash R, M Krishnamoorthy, A Jenkins and M Csete. (2012). Human embryonic stem cell model of ethanol-mediated developmental toxicity. *Exp Neurol* 234:127–135.
46. Zeng X, J Chen, X Deng, Y Liu, MS Rao, JL Cadet and WJ Freed. (2006). An in vitro model of human dopaminergic neurons derived from embryonic stem cells: MPP+ toxicity and GDNF neuroprotection. *Neuropsychopharmacology* 31:2708–2715.
47. Adler S, C Pellizzer, L Hareng, T Hartung and S Bremer. (2008). First steps in establishing a developmental toxicity test method based on human embryonic stem cells. *Toxicol In Vitro* 22:200–211.
48. Jiang Y, D Wang, G Zhang, J Tong and T Chen. (2015). Disruption of cardiogenesis in human embryonic stem cells exposed to trichloroethylene. *Environ Toxicol* 31:1372–1380.
49. Kim H, YY Kim, SY Ku, SH Kim, YM Choi and SY Moon. (2013). The effect of estrogen compounds on human embryoid bodies. *Reprod Toxicol* 20:661–669.
50. Sandström J, E Eggermann, I Charvet, A Roux, N Toni, C Greggio, A Broyer, F Monnet-Tschudi and L Stoppini. (2017). Development and characterization of a human embryonic stem cell-derived 3D neural tissue model for neurotoxicity testing. *Toxicol In Vitro* 38:124–135.
51. Begum AN, JS Aguilar, L Elias, and Y Hong. (2016). Silver nanoparticles exhibit coating and dose-dependent neurotoxicity in glutamatergic neurons derived from human embryonic stem cells. *Neurotoxicology* 57:45–53.
52. Oh JH, MY Son, MS Choi, S Kim, AY Choi, HA Lee, KS Kim, J Kim, CW Song and S Yoon. (2016). Integrative analysis of genes and miRNA alterations in human embryonic stem cell-derived neural cells after exposure to silver nanoparticles. *Toxicol Appl Pharmacol* 299:8–23.
53. Senut M, Y Zhang, F Liu, A Sen, DM Ruden and G Mao. (2016). Size-dependent toxicity of gold nanoparticles on

- human embryonic stem cells and their neural derivatives. *Small* 12:631–646.
54. Braam SR, L Tertoolen, A van de Stolpe, T Meyer, R Passier and CL Mummery. (2010). Prediction of drug-induced cardiotoxicity using human embryonic stem cell-derived cardiomyocytes. *Stem Cell Res* 4:107–116.
 55. Trosko JE and C Chang. (2010). Factors to consider in the use of stem cells for pharmaceutical drug development and for chemical safety assessment. *Toxicology* 270:18–34.
 56. Faiola F, N Yin, X Yao and G Jiang. (2015). The rise of stem cell toxicology. *Environ Sci Technol* 49:5847–5848.
 57. Jung EM, YU Choi, HS Kang, H Yang, EJ Hong, BS An, JY Yang, KH Choi and EB Jeung. (2015). Evaluation of developmental toxicity using undifferentiated human embryonic stem cells. *J Appl Toxicol* 35:205–218.
 58. Yamane J, S Aburatani, S Imanishi, H Akanuma, R Nagano, T Kato, H Sone, S Ohsako and W Fujibuchi. (2016). Prediction of developmental chemical toxicity based on gene networks of human embryonic stem cells. *Nucleic Acids Res* 44:5515–5528.
 59. Mehta A, VBR Konala, A Khanna and AS Majumdar. (2008). Assessment of drug induced developmental toxicity using human embryonic stem cells. *Cell Biol Int* 32:1412–1424.
 60. Fu X, FN Xie, P Dong, QC Li, GY Yu and R Xiao. (2016). High-dose fluoride impairs the properties of human embryonic stem cells via JNK signaling. *PLoS One* 11:e0148819.
 61. Andersson H, D Steel, J Asp, K Dahlenborg, M Jonsson, A Jeppsson, A Lindahl, B Kagedal, P Sartipy and C Mandenius. (2010). Assaying cardiac biomarkers for toxicity testing using biosensing and cardiomyocytes derived from human embryonic stem cells. *J Biotechnol* 150:175–181.
 62. Holmgren G, AK Sjögren, I Barragan, A Sabirsh, P Sartipy, J Synnergren, P Björquist, M Ingelman-Sundberg, TB Andersson and J Edsbacke. (2014). Long-term toxicity testing using human pluripotent stem cell-derived hepatocytes. *Drug Metab Dispos* 42:1401–1406.
 63. Zhao Q, X Wang, S Wang, Z Song, J Wang and J Ma. (2017). Cardiotoxicity evaluation using human embryonic stem cells and induced pluripotent stem cell-derived cardiomyocytes. *Stem Cell Res Ther* 8:54.
 64. Schulpen SHW, E de Jong, LJJ de la Fonteyne, A de Klerk and AH Piersma. (2015). Distinct gene expression responses of two anticonvulsant drugs in a novel human embryonic stem cell based neural differentiation assay protocol. *Toxicol In Vitro* 29:449–457.
 65. Schulpen SHW, JLA Pennings and AH Piersma. (2015). Gene expression regulation and pathway analysis after valproic acid and carbamazepine exposure in a human embryonic stem cell-based neurodevelopmental toxicity assay. *Toxicol Sci* 146:311–320.
 66. Ehashi T, N Suzuki, S Ando, K Sumida and K Saito. (2014). Effects of valproic acid on gene expression during human embryonic stem cell differentiation into neurons. *J Toxicol Sci* 39:383–390.
 67. Stummann TC, L Hareng and S Bremer. (2009). Hazard assessment of methylmercury toxicity to neuronal induction in embryogenesis using human embryonic stem cells. *Toxicology* 257:117–126.
 68. Colleoni S, C Galli, SG Giannelli, MT Armentero, F Blandini, V Broccoli and G Lazzari. (2010). Long-term culture and differentiation of CNS precursors derived from anterior human neural rosettes following exposure to ventralizing factors. *Exp Cell Res* 316:1148–1158.
 69. Colleoni S, C Galli, JA Gaspar, K Meganathan, S Jagtap, J Hescheler, A Sachinidis and G Lazzari. (2011). Development of a neural teratogenicity test based on human embryonic stem cells: response to retinoic acid exposure. *Toxicol Sci* 124:370–377.
 70. Schwartz MP, Z Hou, NE Propson and J Thomson. (2015). Human pluripotent stem cell-derived neural constructs for predicting neural toxicity. *Proc Natl Acad Sci U S A* 112:12516–12521.
 71. West FD, WM Henderson, P Yu, JY Yang, SL Stice and MA Smith. (2012). Metabolomic response of human embryonic stem cell-derived germ-like cells after exposure to steroid hormones. *Toxicol Sci* 129:9–20.
 72. Calderon-Gierszal EL and GS Prins. (2015). Directed differentiation of human embryonic stem cells into prostate organoids in vitro and its perturbation by low-dose bisphenol A exposure. *PLoS One* 10:e0133238.
 73. West PR, Weir AM, Smith AM, Donley ELR and GG Cezar. (2010). Predicting human developmental toxicity of pharmaceuticals using human embryonic stem cells and metabolomics. *Toxicol Appl Pharmacol* 247:18–27.
 74. Kleinstreuer NC, AM Smith and PR West. (2011). Identifying developmental toxicity pathways for a subset of ToxCast chemicals using human embryonic stem cells and metabolomics. *Toxicol Appl Pharmacol* 257:111–121.
 75. Palmer JA, AM Smith, LA Egnash, KR Conard, PR West, RE Burrier, EL Donley and FR Kirchner. (2013). Establishment of assessment of a new human embryonic stem cell-based biomarker assay for developmental toxicity screening. *Birth Defects Res B Dev Reprod Toxicol* 98:343–363.
 76. Baxter MA, C Rowe, J Alder, S Harrison, KP Hanley, BK Park, NR Kitteringham, CE Goldring and NA Hanley. (2010). Generating hepatic cell lineages from pluripotent stem cells for drug toxicity screening. *Stem Cell Res* 5:4–22.
 77. Sartipy P and P Björquist. (2011). Concise review: Human pluripotent stem cell-based models for cardiac and hepatic toxicity assessment. *Stem Cells* 29:744–748.
 78. Hou Z, J Zhang, MP Schwartz, R Stewart, CD Page, WL Murphy and JA Thomson. (2013). A human pluripotent stem cell platform for assessing developmental neural toxicity screening. *Stem Cell Res Ther* 4 (Suppl1):S12.
 79. Sinnecker D, K Laugwitz and A Moretti. (2014). Induced pluripotent stem cell-derived cardiomyocytes for drug development and toxicity testing. *Pharmacol Ther* 143: 246–252.
 80. Pistollato F, J Louise, B Scelfo, M Menecozzi, B Accordi, G Basso, JA Gaspar, D Zagoura, M Barilari, et al. (2014). Development of a pluripotent stem cell derived neuronal model to identify chemically induced pathway perturbations in relation to neurotoxicity: Effects of CREB pathway inhibition. *Toxicol Appl Pharmacol* 280:378–388.
 81. Takayama K, K Kawabata, Y Nagamoto, K Kishimoto, K Tashiro, F Sakurai, M Tachibana, K Kanda, T Hayakawa, MK Furue and H Mizuguchi. (2013). 3D spheroid culture of hESC/hiPSC-derived hepatocyte-like cells for drug toxicity testing. *Biomaterials* 34:1781–1789.
 82. Sirenko O, J Hesley, I Rusyn and EF Cromwell. (2014). High-content assays for hepatotoxicity using induced pluripotent stem cell-derived cells. *Assay Drug Dev Technol* 12:43–54.
 83. Sjögren AK, M Liljevald, B Glinghammar, J Sagemark, XQ Li, A Jonebring, I Cotgreave, G Brolén and TB Andersson. (2014). Critical differences in toxicity mechanisms in induced pluripotent stem cell-derived hepatocytes,

- hepatic cell lines and primary hepatocytes. *Arch Toxicol* 88:1427–1437.
84. Guo L, RM Abrams, JE Babiarz, JD Cohen, S Kameoka, MJ Sanders, E Chiao and KL Colaja. (2011). Estimating the risk of drug-induced proarrhythmia using human induced pluripotent stem cell-derived cardiomyocytes. *Toxicol Sci* 123:281–289.
 85. Liang P, F Lan, AS Lee, T Gong, V Sanchez-Freire, Y Wang, S Diecke, K Sallam, JW Knowles, et al. (2013). Drug screening using a library of human induced pluripotent stem cell-derived cardiomyocytes reveals disease-specific patterns of cardiotoxicity. *Circulation* 127:1677–1691.
 86. Sirenko O, EF Cromwell, C Crittenden, JA Wignall, FA Wright and I Rusyn. (2013). Assessment of beating parameters in human induced pluripotent stem cells enables quantitative in vitro screening for cardiotoxicity. *Toxicol Appl Pharmacol* 273:500–507.
 87. Talbert DR, KR Doherty, PB Trusk, DM Moran, SA Shell and S Bacus. (2015). A multi-parameter in vitro screen in human stem cell-derived cardiomyocytes identifies ponatinib-induced structural and functional cardiac toxicity. *Toxicol Sci* 143:147–155.
 88. Doherty KR, DR Talbert, PB Trusk, DM Moran, SA Shell and S Bacus. (2015). Structural and functional screening in human induced-pluripotent stem cell-derived cardiomyocytes accurately identifies cardiotoxicity of multiple drug types. *Toxicol Appl Pharmacol* 285:51–60.
 89. Burridge PW, YF Li, E Matsa, H Wu, SG Ong, A Sharma, A Holmström, AC Chang, MJ Coronado, et al. (2016). Human induced pluripotent stem cell-derived cardiomyocytes recapitulate the predilection of breast cancer patients to doxorubicin-induced cardiotoxicity. *Nat Med* 22:547–556.
 90. Wheeler HE, C Wing, SM Delaney, MM Komatsu and E Dolan. (2015). Modeling chemotherapeutic neurotoxicity with human induced pluripotent stem cell-derived neuronal cells. *PLoS One* 10:e0118020.
 91. Pei Y, J Peng, M Behl, NS Sipes, KR Shockley, MS Rao, RR Tice and X Zeng. (2016). Comparative neurotoxicity screening in human iPSC-derived neural stem cells, neurons and astrocytes. *Brain Res* 1638:57–73.
 92. Ryan KR, O Sirenko, F Parham, JH Hsieh, EF Cromwell, RR Tice and M Behl. (2016). Neurite outgrowth in human induced pluripotent stem cell-derived neurons as a high-throughput screen for developmental neurotoxicity or neurotoxicity. *Neurotoxicology* 53:271–281.
 93. Yamada S, M Asanagi, N Hirata, H Itagaki, Y Sekino and Y Kanda. (2016). Tributyltin induces mitochondrial fission through Mfn1 degradation in human induced pluripotent stem cells. *Toxicol In Vitro* 34:257–263.
 94. Jennings P. (2014). The future of in vitro toxicology. *Toxicol In Vitro* 29:1217–1221.
 95. La Farré M, S Pérez, L Kantiani and D Barceló. (2008). Fate and toxicity of emerging pollutants, their metabolites and transformation products in the aquatic environment. *TrAC Trends Anal Chem* 27:991–1007.

Address correspondence to:

Francesco Faiola
State Key Laboratory of Environmental Chemistry
and Ecotoxicology
Research Center for Eco-Environmental Sciences
Chinese Academy of Sciences
Beijing 100085
P.R. China

E-mail: faiola@rcees.ac.cn

Received for publication July 20, 2017

Accepted after revision September 5, 2017

Prepublished on Liebert Instant Online September 6, 2017