



Assays of D-Amino Acid Oxidase Activity

Elena Rosini^{1,2}, Laura Caldinelli^{1,2} and Luciano Piubelli^{1,2*}

¹ Department of Biotechnology and Life Sciences, University of Insubria, Varese, Italy, ² The Protein Factory Research Center, Politecnico di Milan and University of Insubria, Milan, Italy

OPEN ACCESS

Edited by:

Jumpei Sasabe,
Keio University, Japan

Reviewed by:

Mariarita Bertoldi,
University of Verona, Italy
Argyrides Argyrou,
AstraZeneca (United Kingdom),
United Kingdom

*Correspondence:

Luciano Piubelli
luciano.piubelli@uninsubria.it

Specialty section:

This article was submitted to
Structural Biology,
a section of the journal
Frontiers in Molecular Biosciences

Received: 17 October 2017

Accepted: 28 December 2017

Published: 18 January 2018

Citation:

Rosini E, Caldinelli L and Piubelli L
(2018) Assays of D-Amino Acid
Oxidase Activity.
Front. Mol. Biosci. 4:102.
doi: 10.3389/fmolb.2017.00102

D-amino acid oxidase (DAAO) is a well-known flavoenzyme that catalyzes the oxidative FAD-dependent deamination of D-amino acids. As a result of the absolute stereoselectivity and broad substrate specificity, microbial DAAOs have been employed as industrial biocatalysts in the production of semi-synthetic cephalosporins and enantiomerically pure amino acids. Moreover, in mammals, DAAO is present in specific brain areas and degrades D-serine, an endogenous coagonist of the N-methyl-D-aspartate receptors (NMDARs). Dysregulation of D-serine metabolism due to an altered DAAO functionality is related to pathological NMDARs dysfunctions such as in amyotrophic lateral sclerosis and schizophrenia. In this protocol paper, we describe a variety of direct assays based on the determination of molecular oxygen consumption, reduction of alternative electron acceptors, or α -keto acid production, of coupled assays to detect the hydrogen peroxide or the ammonium production, and an indirect assay of the α -keto acid production based on a chemical derivatization. These analytical assays allow the determination of DAAO activity both on recombinant enzyme preparations, in cells, and in tissue samples.

Keywords: D-amino acids, D-amino acid oxidase, flavoproteins, enzymatic activity, analytical detection


INTRODUCTION

D-Amino acid oxidase (DAAO, EC 1.4.3.3) is a flavoenzyme containing a non-covalently bound FAD molecule per 40 kDa monomer that belongs to the dehydrogenase/oxidase class of flavoproteins (Pollegioni et al., 2007a). DAAO catalyzes the oxidative deamination of the D-isomer of α -amino acids to the corresponding α -keto acids, according to the scheme reported in **Figure 1**. It shows a strict selectivity toward D-isomers (i.e., it does not oxidize L-amino acids) and possesses a broad substrate specificity: it oxidizes aliphatic, aromatic and polar α -D-amino acids, whereas D-aspartate and D-glutamate are not substrates for DAAO (Pollegioni et al., 2008).

The reaction catalyzed by DAAO can be divided into (Pollegioni et al., 2008):

- 1) a reductive half-reaction, in which the enzyme catalyzes the dehydrogenation of α -D-amino acids to their imino acid counterparts with concomitant reduction of FAD to FADH₂ [reaction (1) in **Figure 1**];

Abbreviations: 4-AAP, 4-aminoantipyrine; CNS, central nervous system; DAAO, D-amino acid oxidase; DCPIP, 2,6-Dichlorophenol-indophenol; DNaseI, desossiribonuclease I; DNP, 2,4-dinitrophenylhydrazine; EDTA, ethylenediaminetetraacetic acid; FAD/FADH₂, flavin adenine dinucleotide oxidized/reduced form; GDH, glutamate dehydrogenase; hDAAO, human D-amino acid oxidase; HRP, horseradish peroxidase; MG, methylene green; NAD⁺/NADH, nicotinamide adenine dinucleotide oxidized/reduced form; NMDARs, N-methyl-D-aspartate receptors; o-DNS, o-dianisidine; RgDAAO, D-amino acid oxidase from *Rhodotorula gracilis*; Thi, thionine.



- 2) an oxidative half-reaction, in which FADH₂ is reoxidized by molecular oxygen (the physiological electron acceptor) to produce hydrogen peroxide [reaction (2a) in **Figure 1**]; FADH₂ can also be reoxidized *in vitro* by artificial electron acceptors, such as some redox dyes [reaction (2b) in **Figure 1**].
- 3) The spontaneous hydrolysis of the imino acid to give the corresponding α -keto acid and ammonium [reaction (3) in **Figure 1**].

DAAO is present in all eukaryotic taxa with the exception of plants (Pollegioni et al., 2007a). All known DAAOs share the same catalytic mechanism, but show important differences in biochemical and structural properties, such as catalytic efficiency, substrate specificity, oligomeric state, stability, kinetic mechanism, and FAD binding. The diverse properties of DAAOs from different sources reflect their physiological role (for an exhaustive review, see Pollegioni et al., 2007a). DAAO is an efficient catabolic enzyme in yeasts where it oxidizes α -D-amino acids, making it possible to use them as carbon, nitrogen, or energy sources. The main role of DAAO in mammals (including humans) is related to its presence in selected brain areas where it is devoted to the catabolism of D-serine. This α -D-amino acid is a neuromodulator acting as a coagonist of the N-methyl-D-aspartate receptors (NMDARs). Dysfunction of this receptor has been correlated to different neurological or psychiatric diseases: an overstimulation of these receptors is involved in stroke, epilepsy, in neurodegenerative pathologies such as Parkinson's and Alzheimer's diseases and amyotrophic lateral sclerosis. An hypostimulation of NMDARs, in contrast, is involved in psychiatric diseases such as schizophrenia, attention-deficit hyperactivity disorder, and chronic depression (Ross et al., 2006; Wu et al., 2007; Mitchell et al., 2010; Pollegioni and Sacchi, 2010; Nagasawa et al., 2012; Sacchi et al., 2012; Zhou and Sheng, 2013). Thus, by modulating D-serine levels, DAAO plays a key role in regulating NMDARs activation state in mammals. Owing

to of its role in human brain, molecules that modulate or inhibit DAAO can act as innovative drugs for the treatment of the many diseases linked to NMDARs dysfunction (Sacchi et al., 2013).

Microbial DAAOs possess properties that render them suitable for industrial biotechnological applications: for example, they are stable enzymes and show broad substrate specificity, high turnover number (**Table 1**), and a tight binding with the FAD cofactor. The production of DAAOs in large amounts as recombinant proteins together with the availability of the 3D-structure (e.g., Mattevi et al., 1996; Mizutani et al., 1996; Umhau et al., 2000; Kawazoe et al., 2006), allowed the design and production by protein engineering techniques of enzyme variants with new and evolved properties (Pollegioni et al., 2007b; Rosini et al., 2008, 2009, 2010; Wang et al., 2008; Wong et al., 2010; Pollegioni and Molla, 2011; Golubev et al., 2014).

The main biotechnological applications of DAAO have been reviewed in detail in Pollegioni et al. (2007b, 2008) and Pollegioni and Molla (2011) and comprise: (i) production of 7-amino cephalosporanic acid from cephalosporin C, the starting molecule for the production of semi-synthetic cephalosporins; (ii) resolution of racemic mixtures of amino acids to produce enantiomerically pure natural or synthetic amino acids as fine chemicals; (iii) detection and quantification of D-amino acids in biological samples and foodstuff—D-amino acids are components of the bacterial cell wall and thus their presence in food could indicate bacterial contamination and can be used as an indicator of foodstuff aging; and (iv) use as selective marker in plants.

DAAO activity can be determined using different assays based on: (i) molecular oxygen consumption, directly measured by a Clark-type oxygen electrode (Sacchi et al., 2004; Molla et al., 2006); (ii) reduction of redox dyes used as alternative electron acceptors (instead of molecular oxygen), followed by a colorimetric assay (Brugger et al., 2014); (iii) hydrogen peroxide production, indirectly detected by enzyme-coupled assays in

TABLE 1 | Comparison of kinetic properties of DAAO from different organisms.

	Yeast DAAOs		Mammalian DAAOs		
	RgDAAO ^a	TvDAAO ^a	pkDAAO ^a	rDAAO ^b	hDAAO ^c
Kinetic parameters on D-alanine^d					
$k_{cat, app}$ (s ⁻¹)	81 ± 5	46 ± 3	7.3 ± 0.6	27 ± 1	5.2 ± 0.1
K_m, app (mM)	1.0 ± 0.2	7.0 ± 0.9	1.7 ± 0.3	140 ± 20	1.3 ± 0.2

RgDAAO, DAAO from *Rhodotorula gracilis*; TvDAAO, DAAO from *Trigonopsis variabilis*; pkDAAO, DAAO from pig kidney; rDAAO, DAAO from rat; hDAAO, DAAO from *Homo sapiens*.

^a(Pollegioni et al., 2004), ^b(Frattini et al., 2011), ^c(Molla et al., 2006). ^dThe activity was determined using the oxygen-consumption assay, at 25°C, pH 8.5, and at air saturation ($[O_2] = 0.253$ mM).

which a chromogenic (*o*-dianisidine or 4-aminoantipyrine; Job et al., 2002; Sacchi et al., 2004; Rosini et al., 2008) or a fluorogenic (Amplex[®] UltraRed, Invitrogen; Sacchi et al., 2011; Hopkins et al., 2013) substrate is used in combination with horseradish peroxidase, or in which a chromogenic reagent (Purpald[®]: 4-amino-3-hydrazino-5-mercapto-1,2,4-triazole, Sigma-Aldrich) is employed in combination with catalase (Dickinson and Jacobsen, 1970; Watanabe et al., 1978; Sasabe et al., 2014); (iv) α -keto acid formation, directly detected by measuring the increase in absorbance in the ultraviolet range (Tedeschi et al., 2012) or indirectly detected following the reaction of the α -keto acid produced with 2,4-dinitrophenylhydrazine, to give a chromogenic 2,4-dinitrophenylhydrazone derivative (Nagata et al., 1988); and (v) ammonium production, indirectly detected using an enzyme-coupled assay in which α -ketoglutarate is converted to L-glutamate by the NAD⁺-dependent L-glutamate dehydrogenase (Holme and Goldberg, 1975).

Here, we describe in detail the analytical assays listed above; full lists of materials and equipment required and step-by-step procedures are given, together with examples of application and anticipated results.

MATERIALS AND EQUIPMENT

Unless otherwise stated, all chemicals and reagents were of analytical grade and were purchased from Sigma-Aldrich (Milan, Italy). Solutions were prepared at room temperature using ultrapure water (e.g., Milli-Q grade water) and assays were carried out at 25°C.

One DAAO unit is defined as the amount of enzyme that converts 1 μ mol of D-amino acid per minute at 25°C (Molla et al., 1998). O_2 concentration at 25°C and at air saturation is 0.253 mM. Unless otherwise stated, in order to obtain reliable specific activity values on different D-amino acids, a substrate concentration 10-fold higher than the K_m value is used.

Enzymes

Recombinant RgDAAO and hDAAO were expressed and purified from *Escherichia coli* cells as previously described (Fantinato et al., 2001; Molla et al., 2006). The concentration of purified enzyme preparations was determined by using the extinction coefficient at ~ 450 nm (12.6 and 12.2 mM⁻¹cm⁻¹ for RgDAAO and hDAAO, respectively).

Determination of Molecular Oxygen Consumption: Polarographic Assay

DAAO activity can be assayed polarographically using an oxygen electrode (Sacchi et al., 2004; Molla et al., 2006).

- Oxygraph system (Hansatech Instr. Ltd, Pentney, Norfolk, UK): a highly sensitive S1 Clark-type polarographic oxygen electrode disc is mounted in a DW1/AD electrode chamber and connected to the Oxygraph electrode control unit. The electrode disc comprises a central platinum cathode and a concentric circular silver anode. By using the electrode chamber, dissolved oxygen can be measured in liquid-phase samples in the 0.2–2.5 mL volume range. Precise temperature control of the sample and electrode disc (the oxygen solubility decreases as temperature increases) is achieved by connecting the water jacket of the electrode chamber to a thermoregulated circulating water bath. The sample is continuously stirred in order to ensure that the dissolved oxygen is kept evenly distributed throughout the reaction vessel.
- Electrolyte: 50% saturated solution of potassium chloride.
- Cigarette Rizla+ Blue regular rolling paper (Rizla).
- Thin layer of polytetrafluoroethylene (PTFE) membrane, selectively permeable to molecular oxygen (Hansatech Instruments Ltd).
- Sodium dithionite crystals.
- Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 75 mM (final concentration) disodium pyrophosphate buffer, pH 8.5.

Determination of the Redox Dye Reduction

Different redox dyes are known for their ability to directly react with the reduced form of the cofactor of flavin-dependent oxidases and dehydrogenases and thus can be used as artificial electron acceptors for the assay of these enzymes. 2,6-Dichlorophenol-indophenol (DCPIP), methylene green (MG), and thionine (Thi) show absorption changes in the visible range between their oxidized and reduced forms and were thus selected to directly detect different enzymatic activities, including DAAO (Brugger et al., 2014).

- Microtiter plate reader (Infinite 200, Tecan, Cernusco sul Naviglio, Milan, Italy).
- 96-well microplates.
- 50 mM phosphate buffer, pH 8.0.
- Redox dye (electron acceptor) solutions: 3 mM DCPIP or 0.5 mM MG or 1 mM Thi, dissolved in water.
- Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 50 mM phosphate buffer, pH 8.0.

Determination of Hydrogen Peroxide Production

Hydrogen peroxide can be detected by a second enzyme-coupled reaction in which the chromogenic *o*-dianisidine (*o*-DNS) (Sacchi et al., 2004; Rosini et al., 2008) or 4-aminoantipyrine (4-AAP) (Job et al., 2002) or the fluorogenic Amplex[®] UltraRed

Reagent (Invitrogen, Thermo-Fisher Sci. Co., Carlsbad, CA, USA; Sacchi et al., 2011; Hopkins et al., 2013) reagent substrates are used together with horseradish peroxidase (HRP, EC 1.11.1.7; Roche, Mannheim, Germany; specific activity: ~ 225 U/mg lyophilized. One unit of HRP is defined as the amount of enzyme that converts 1 μ mol of hydrogen peroxide per minute at 25°C in the presence of guaiacol; **Figure 2**), or the chromogenic Purpald[®] reagent is used together with catalase (EC 1.11.1.6; Genencor, Palo Alto, CA, USA; one unit of catalase is defined as the amount of enzyme which produces 1 μ mol of hydrogen peroxide per minute at 30°C, pH 6.8; **Figure 3**; Dickinson and Jacobsen, 1970; Watanabe et al., 1978; Sasabe et al., 2014).

(A) Horseradish Peroxidase and *o*-Dianisidine Coupled Assay

1. UV/Vis spectrophotometer (V-560, Jasco, Cremella, Lecco, Italy) connected to a water bath for temperature control or a microtiter plate reader (Infinite 200, Tecan).
2. 1.5 mL plastic disposable microcuvettes or 96-well microplates.
3. Automated liquid handling system (epMotion[®] 5075, Eppendorf, Milan, Italy).
4. 100 mM disodium pyrophosphate buffer, pH 8.5.
5. 10 mM *o*-DNS dihydrochloride in water. The solution must be freshly prepared every day and kept on ice.
6. 100 U/mL HRP in 100 mM disodium pyrophosphate buffer, pH 8.5. A concentrated stock solution (e.g., 2,500 U/mL) can be stored at 4°C for several months; the diluted working solution must be freshly prepared every day and kept on ice.
7. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 100 mM disodium pyrophosphate buffer, pH 8.5.
8. Lysis buffer: 50 mM disodium pyrophosphate buffer, pH 8.5, 1 mM EDTA, 100 mM NaCl, 40 μ g/mL lysozyme (from chicken egg white), and 1 μ g/mL DNaseI (from bovine pancreas; Roche), all final concentrations. Lysozyme and DNaseI must be freshly added from appropriate stock solutions.

(B) Horseradish Peroxidase and 4-Aminoantipyrine Coupled Assay

1. UV/Vis spectrophotometer (V-560, Jasco) connected to a water bath for temperature control or a microtiter plate reader (Infinite 200, Tecan).
2. 1.5 mL plastic disposable microcuvettes or 96-well microplates.
3. Automated liquid handling system (epMotion[®] 5075, Eppendorf).
4. 100 mM disodium pyrophosphate buffer, pH 8.5.
5. 15 mM 4-AAP solution in water. The solution must be freshly prepared every day and kept on ice.
6. 200 mM phenol.
7. 250 U/mL HRP in 100 mM disodium pyrophosphate buffer, pH 8.5. A concentrated stock solution (e.g., 2,500 U/mL) can be stored at 4°C for several months; the diluted working solution must be freshly prepared every day and kept on ice.


8. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 100 mM disodium pyrophosphate buffer, pH 8.5.
9. Lysis buffer: 50 mM disodium pyrophosphate buffer, pH 8.5, 1 mM EDTA, 100 mM NaCl, 40 μ g/mL lysozyme, and 1 μ g/mL DNaseI (all final concentrations). Lysozyme and DNaseI must be freshly added from appropriate stock solutions.

(C) Horseradish Peroxidase and Amplex[®] UltraRed Coupled Assay

1. Microtiter plate reader (Infinite 200, Tecan).
2. Black 96-well microplates.
3. 50 mM sodium phosphate buffer, pH 7.4.
4. 10 U/mL HRP in 50 mM sodium phosphate buffer, pH 7.4. A concentrated stock solution (e.g., 2,500 U/mL) can be stored at 4°C for several months; the diluted working solution must be freshly prepared every day and kept on ice.
5. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 50 mM sodium phosphate buffer, pH 7.4.
6. Inhibitor stock solutions: different compounds dissolved in dimethyl sulfoxide at 10 mM final concentration and stored at 4°C.
7. Lysis buffer: 50 mM sodium phosphate buffer, pH 7.4, 1 μ M pepstatin, 2 μ M leupeptin, 500 μ M phenylmethylsulfonyl fluoride (PMSF), 10 μ M FAD (all final concentrations). Pepstatin, leupeptin, PMSF, and FAD must be freshly added from appropriate stock solutions.
8. 0.5 M sodium azide in 50 mM sodium phosphate buffer, pH 7.4. This solution can be stored at room temperature for several months.
9. 500 mM sodium benzoate in 50 mM sodium phosphate buffer, pH 7.4.
10. 1 mM and 10 mM FAD solutions in 50 mM sodium phosphate buffer, pH 7.4.
11. Amplex[®] UltraRed Reagent solution: allow this reagent to equilibrate at room temperature before use. Prepare a 10 mM stock solution adding 340 μ L of dimethyl sulfoxide to the content of one vial (each vial contains 1 mg of the lyophilized powder) and vortex. Prepare 25 μ L aliquots and store in the dark at -20°C .
12. Amplex[®] UltraRed Stop Reagent solution: dissolve the content of the Stop Reagent vial in 1.45 mL of ethanol. Vortex briefly and mix 1 mL of this solution with 1 mL of water. The solution must be freshly prepared every day and stored in the dark at 4°C. The solution is colorless; the appearance of dark coloration is indicative of decomposition.
13. Working solutions (WS): a different working solution is required for each application:

(a) Detection of DAAO activity of purified protein preparations on different D-amino acids

WS_A: mix 50 μ L of 10 mM Amplex[®] UltraRed Reagent solution and 100 μ L of 10 U/mL HRP in 5 mL (final volume) of 50 mM sodium phosphate buffer, pH 7.4;



(b) Detection of DAAO activity in U87 cells. For this application, three different WS are required:

WS_{B1}: add 150 μ L of 400 mM D-serine, 75 μ L of 1 mM FAD and 30 μ L of 500 mM sodium azide to 1.2 mL of 50 mM sodium phosphate buffer, pH 7.4. Just before aliquoting the working solution in the wells, add 15 μ L of 10 mM Amplex® UltraRed Reagent solution and 30 μ L of 10 U/mL HRP. The final volume is 1.5 mL;

WS_{B2}: to be used for the negative control without substrate. This solution has the same composition as WS_{B1}, except that D-serine is replaced by 150 μ L of buffer (50 mM sodium phosphate buffer, pH 7.4);

WS_{B3}: to be used for the negative control in the presence of benzoate (a well-known DAAO inhibitor): 9 μ L of


500 mM sodium benzoate are added to 1.5 mL of WS_{B1} (final sodium benzoate concentration: 3 mM).

(c) Identification of hDAAO Inhibitors

WS_C: mix 220 μ L of 1 M D-serine, 4 μ L of 10 mM FAD, 100 μ L of 10 mM Amplex® UltraRed Reagent solution, and 200 μ L of 10 U/mL HRP in 10 mL (final volume) of 50 mM sodium phosphate buffer, pH 7.4.

- Hydrogen peroxide stock solution 30% (w/w) in water (9.8 M); the diluted solutions (10–100 μ M, to be used for the calibration curves) must be freshly prepared every day.

(D) Catalase and Purpald® Coupled Assay



1. UV/Vis spectrophotometer (V-560, Jasco) connected to a water bath for temperature control.
2. 1.5 mL plastic disposable microcuvettes.
3. 5 mL plastic tubes.
4. 133 mM sodium pyrophosphate buffer, pH 8.3.
5. 34 mM Purpald® in 0.5 M HCl.
6. 0.1 mM FAD in 133 mM sodium pyrophosphate buffer, pH 8.3.
7. 0.75% (w/v) KIO₄ in 0.2 M KOH.
8. 70% (v/v) methanol.
9. 10% (w/v) trichloroacetic acid.
10. 5 M KOH.
11. 700 U/mL catalase in 133 mM sodium pyrophosphate buffer, pH 8.3. The solution must be freshly prepared every day and kept on ice.
12. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 133 mM sodium pyrophosphate buffer, pH 8.3.
13. Human post mortem tissue samples from different regions of the central nervous system (CNS).

Determination of α -Keto Acid Production

α -Keto acids concentration can be determined by measuring their absorbance in the ultraviolet range (Tedeschi et al., 2012) or following their chemical derivatization with 2,4-dinitrophenylhydrazine (DNP) to produce the corresponding 2,4-dinitrophenylhydrazone derivatives, which show an absorption maximum at 445 nm (Figure 4; Nagata et al., 1988).

A) Direct Spectrophotometric Assay

1. UV/Vis spectrophotometer (V-560, Jasco) connected to a water bath for temperature control.
2. 1.5 mL quartz microcuvettes.
3. 75 mM disodium pyrophosphate buffer, pH 8.5.
4. α -Keto acid solution in 75 mM disodium pyrophosphate buffer, pH 8.5 (for the calibration curve).
5. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 75 mM disodium pyrophosphate buffer, pH 8.5.



B) Indirect Assay with 2,4-Dinitrophenylhydrazine

1. UV/Vis spectrophotometer (V-560, Jasco) connected to a water bath for temperature control.
2. 1.5 mL glass or plastic disposable microcuvettes.
3. 1 mM DNP solution in 1 M HCl.
4. 0.6 M NaOH.
5. 75 mM disodium pyrophosphate buffer, pH 8.5.
6. α -Keto acid solution in 75 mM disodium pyrophosphate buffer, pH 8.5 (for the calibration curve).
7. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 75 mM disodium pyrophosphate buffer, pH 8.5.

Determination of Ammonium Production: Glutamate Dehydrogenase and α -Ketoglutarate Coupled Assay

The ammonium produced by DAAO reaction on a D-amino acid can be detected by a second enzyme-coupled reaction in which α -ketoglutarate is converted to L-glutamate by L-glutamate dehydrogenase from bovine liver (GDH, EC 1.4.1.3; specific activity ≥ 35 U/mg protein. One GDH unit is defined as the amount of enzyme that reduces 1 μ mol of α -ketoglutarate to L-glutamate per minute at pH 7.3 and 25°C, in the presence of ammonium ions) and NADH is oxidized to NAD⁺ (Figure 5; Holme and Goldberg, 1975). The decrease of absorbance at 340 nm, due to the conversion of NADH to NAD⁺, is followed.

1. UV/Vis spectrophotometer (V-560, Jasco) connected to a water bath for temperature control.
2. 1.5 mL quartz microcuvettes.
3. 1,000 U/mL GDH in 75 mM disodium pyrophosphate buffer, pH 8.5.
4. 75 mM disodium pyrophosphate buffer, pH 8.5.
5. 5 mM α -ketoglutarate in 75 mM disodium pyrophosphate buffer, pH 8.5.
6. 0.25 mM NADH.



7. Substrates: D-amino acid solutions (at different concentrations, depending on the corresponding K_m values) in 75 mM disodium pyrophosphate buffer, pH 8.5.

STEPWISE PROCEDURES

Determination of Molecular Oxygen Consumption: Polarographic Assay

1. Preparation of the electrode disc: before use, a small drop of electrolyte is placed on top of the dome of the central platinum cathode of the electrode disc and a 1.5 cm² cigarette Rizla+ Blue rolling paper spacer is placed over the electrolyte. The paper spacer is covered with a similar sized piece of thin, oxygen-permeable polytetrafluoroethylene membrane. Once prepared, the disc is mounted in the base of the electrode chamber and 1 mL of deionized water is added in the chamber.
2. Calibration: the disc is connected to the electrode control unit, which applies a polarizing voltage (0.7 V) between the platinum and silver electrodes and generates a small current proportional to oxygen concentration in the sample. The amperometric value measured with 1 mL of stirred air-saturated water is considered as 100% of oxygen content, and the amperometric value measured after adding a few crystals of sodium dithionite corresponds to 0% of oxygen content.
3. Measurement of DAAO activity: 1 mL of the substrate solution is placed in the electrode chamber, thermostated at 25°C, and the registration of amperometric signal is started. Once the recorded signal is stable, 10 μ L of DAAO solution (\sim 0.05 U) is added. The DAAO activity is calculated from the initial (1–2 min) oxygen consumption rate (μ mol/min) using the

following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{\mu\text{mol}_{\text{O}_2}/\text{min}}{\text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where mL_{TOT} is the final volume of the reaction mixture (1.01 mL). To accurately determine the DAAO activity, the initial oxygen consumption rate should be in the 0.02–0.05 μ mol/min range. The detection range is reported in **Table 3**.

Application: Determination of the Steady-State Kinetic Parameters and Substrate Specificity

The apparent kinetic parameters are investigated by this assay, employing different concentrations of each D-amino acid as substrate (Sacchi et al., 2004; Molla et al., 2006). The initial reaction rates (μ mol/min) at different substrate concentrations (in the 0.1–10-fold K_m range) are fitted using the Michaelis-Menten equation:

$$v = \frac{V_{\text{max}} \times [S]}{K_m + [S]}$$

where V_{max} and K_m are the apparent (i.e., at fixed O_2 concentration, 0.253 mM) maximal velocity and Michaelis-Menten constant for the substrate tested.

Determination of the Redox Dye Reduction

This assay is based on the reduction of one of the redox dyes DCPIP, MG, and Thi by DAAO. These dyes show absorption changes between their oxidized (colored) and reduced (colorless) forms in the visible range, with absorption maxima at 520, 655, and 600 nm, respectively (Brugger et al., 2014).

1. The reaction mixture is prepared in a 96-well microplate, mixing the electron acceptor solution (450 μM DCPIP, 75 μM MG, or 100 μM Thi, all final concentrations) and the D-amino acid solution in 50 mM phosphate buffer, pH 8.0, in a total volume of 100 μL .
2. 10–20 μL of DAAO are added (~ 0.2 U).
3. The electron acceptor reduction is followed at 520 nm for DCPIP, 655 nm for MG, and 600 nm for Thi ($\text{Abs}_{x \text{ nm}}$).
4. The DAAO activity is calculated from the initial absorbance variation at the corresponding wavelength per minute ($\Delta\text{Abs}_{x \text{ nm}}/\text{min}$) using the following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{\Delta\text{Abs}_{x \text{ nm}}/\text{min}}{\varepsilon_{x \text{ nm}} (\text{redox dye}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where the molar extinction coefficients ($\varepsilon_{x \text{ nm}}$) of DCPIP, MG, and Thi are 6.6, 46.6, and 55.4 $\text{mM}^{-1}\text{cm}^{-1}$, at 520, 655, and 600 nm, respectively. mL_{TOT} is the final volume of the reaction mixture (0.1 mL). For an accurate determination of the DAAO activity, the $\Delta\text{Abs}_{x \text{ nm}}/\text{min}$ should be in the 0.1–0.3 range. The detection range is reported in **Table 3**.

Determination of Hydrogen Peroxide Production

A) Horseradish Peroxidase and *o*-Dianisidine Coupled Assay

The hydrogen peroxide produced by the DAAO reaction is reduced by HRP to water and *o*-DNS is simultaneously oxidized to give a colored compound showing an absorption maximum at 440 nm with an extinction coefficient of 13 $\text{mM}^{-1}\text{cm}^{-1}$ (see **Figures 1, 2A**).

1. The reaction mixture is prepared in a plastic disposable microcuvette or in a 96-well microplate, mixing the D-amino acid solution, 1 mM *o*-DNS, and 1 U/mL of HRP, in 75 mM disodium pyrophosphate buffer, pH 8.5 (all final concentrations), in a total volume of 1 mL or 300 μL , for microcuvettes or 96-well microplates, respectively.
2. 10–20 μL of DAAO are added (~ 0.01 U).
3. The absorbance at 440 nm ($\text{Abs}_{440 \text{ nm}}$) is monitored.
4. The DAAO activity is calculated from the initial absorbance variation at 440 nm per minute ($\Delta\text{Abs}_{440 \text{ nm}}/\text{min}$) using the following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{\Delta\text{Abs}_{440 \text{ nm}}/\text{min}}{\varepsilon_{440 \text{ nm}} (\textit{o}-\text{DNS}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where the molar extinction coefficient at 440 nm of *o*-DNS is 13 $\text{mM}^{-1}\text{cm}^{-1}$ and mL_{TOT} is the final volume of the reaction mixture (1 or 0.3 mL). To accurately determine the DAAO activity, the $\Delta\text{Abs}_{440 \text{ nm}}/\text{min}$ should be in the 0.1–0.3 range. The detection range is reported in **Table 3**.

Application: Screening for RgDAAO Variants with Broader Substrate Specificity

The reaction catalyzed by RgDAAO has been exploited in the analytical determination of the D-amino acid content in biological samples (Rosini et al., 2008). Since the enzyme is not active on acidic D-amino acids, it cannot be used to detect

the total amount of D-amino acids. In an attempt to obtain DAAO variants to detect and quantify the total amount of D-amino acids, a directed evolution approach was performed; the wild-type DAAO cDNA was randomly mutated by error-prone PCR and a colorimetric screening procedure was set up to determine any alteration in the substrate specificity of RgDAAO (Sacchi et al., 2004). This screening procedure is performed in an automated liquid handling system (see **Figure 6**):

1. To the 1 mL saturated *E. coli* culture 1 mM isopropyl β -D-thiogalactopyranoside (IPTG, final concentration) is added and then incubated at 30°C for 2 h.
2. 50 μL of the cells culture are transferred to different wells of a 96-well microplate (one microplate for each substrate tested).
3. 150 μL of lysis buffer are added in each well and the plates incubated at 37°C for 30 min.
4. DAAO activity is assayed on the crude extracts by adding 100 μL of developing solution (100 mM D-amino acid, 1 mM *o*-DNS and 1 U/mL HRP, in 100 mM disodium pyrophosphate buffer, pH 8.5, all final concentrations).
5. The time course of the absorbance change at 440 nm is followed at room temperature by using a microtiter plate reader and comparing it with that of cells expressing wild-type RgDAAO as control.

(B) Horseradish Peroxidase and 4-Aminoantipyrine Coupled Assay

This assay determines the amount of hydrogen peroxide produced by the reaction catalyzed by DAAO coupled to the oxidation of 4-AAP by HRP in the presence of phenol. The resulting quinoneimine product shows an absorption maximum at 505 nm with an extinction coefficient of 6.58 $\text{mM}^{-1}\text{cm}^{-1}$ (see **Figure 2B**).

1. The reaction mixture is prepared in a plastic disposable microcuvette or in a 96-well microplate, mixing the D-amino acid solution, 1.5 mM 4-AAP, 2 mM phenol, and 2.5 U/mL HRP, in 75 mM disodium pyrophosphate buffer, pH 8.5 (all final concentrations), in a total volume of 1 mL or 300 μL , for microcuvettes or 96-well microplates, respectively.
2. 10–20 μL of DAAO are added (~ 0.05 U).
3. The absorbance at 505 nm ($\text{Abs}_{505 \text{ nm}}$) is monitored.
4. The DAAO activity is calculated from the initial absorbance variation at 505 nm per minute ($\Delta\text{Abs}_{505 \text{ nm}}/\text{min}$) using the following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{2 \times \Delta\text{Abs}_{505 \text{ nm}}/\text{min}}{\varepsilon_{505 \text{ nm}} (\text{quinoneimine}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where the molar extinction coefficient at 505 nm of quinoneimine is 6.58 $\text{mM}^{-1}\text{cm}^{-1}$ and mL_{TOT} is the final volume of the reaction mixture (1 or 0.3 mL). To accurately determine the DAAO activity, the $\Delta\text{Abs}_{505 \text{ nm}}/\text{min}$ should be in the 0.1–0.3 range. The detection range is reported in **Table 3**. Owing to the stoichiometry of the reaction, in which two molecules of hydrogen peroxide react with 4-AAP to produce one molecule of quinoneimine (see **Figure 2B**), a factor 2 is added at the numerator.




FIGURE 6 | Schematic representation of the screening procedure on crude extract of *E. coli* cells expressing RgDAAO variants.

C) Horseradish Peroxidase and Amplex[®] UltraRed Coupled Assay

This assay is based on the reaction catalyzed by HRP on the fluorogenic substrate Amplex[®] UltraRed with hydrogen peroxide in a 1:1 stoichiometric ratio to give the brightly fluorescent product Amplex[®] UltroRed (see Figure 2C).

General procedure:

1. The reaction mixture, as described for each application below, is prepared and set up in a 96-well microplate.
2. The microplate is incubated for 30–60 min in the dark at 25°C under agitation.
3. The Stop Reagent solution is added to each well and the microplate is incubated for 5 min in the dark at 25°C under agitation. This step allows the fluorescence signal-generating reaction to be terminated at a user-determined time point, ensuring that the timing of the standard and unknown samples measurements is the same. After reagent addition, the fluorescence signal is stable for at least 3 h.
4. The fluorescence intensity is recorded using a fluorimeter microplate reader: $\lambda_{exc} = 535$ nm (slit = 25 nm); $\lambda_{em} = 595$ nm (slit = 20 nm).
5. The enzymatic activity is calculated on the basis of a calibration curve obtained using known amounts of purified recombinant DAAO (0.02–0.4 mU range). The detection range is reported in Table 3. Due to the response variability

of the fluorescence signal, the calibration curve must be carried out on each 96-well microplate, and triplicates for each standard and sample are recommended.

Applications

(a) Detection of DAAO Activity on Different D-Amino Acids of Purified Protein Preparations

1. Preparation of the reaction mixture: 50 μ L of WS_A (0.1 mM Amplex[®] UltraRed Reagent, 0.5 U/mL HRP, in 50 mM sodium phosphate buffer, pH 7.4), 45 μ L of D-amino acid solution, and 5 μ L of DAAO (~1 mU) are mixed in each well;
2. Reaction is stopped by adding 20 μ L of Stop Reagent solution after 30 min.

(b) Detection of DAAO Activity in U87 Cells

The enzymatic assay is performed on U87 glioblastoma control (not transfected) cells and on the same cells transfected with the plasmid expressing hDAAO (Sacchi et al., 2011).

1. The U87 cells are resuspended in lysis buffer (4×10^5 cells per mL) and disrupted by sonication.
2. The soluble fraction is collected by centrifugation at $16,000 \times g$ for 30 min at 4°C.
3. Preparation of the reaction mixture: 50 μ L of cell extract (the supernatant obtained in the previous step) and 50 μ L of WS_{B1} (40 mM D-serine, 50 μ M FAD, 10 mM sodium azide,

0.1 mM Amplex[®] UltraRed Reagent, 0.2 U/mL HRP, in 50 mM sodium phosphate buffer, pH 7.4) are mixed in each well.

On each microplate, three negative controls must be set up:

- *control I* (without cell extract): 50 μ L of lysis buffer and 50 μ L of WS_{B1} are mixed in each well;
 - *control II* (without substrate): 50 μ L of cell extract and 50 μ L of WS_{B2} are mixed in each well;
 - *control III* (with a DAAO inhibitor): 50 μ L of cell extract and 50 μ L of WS_{B3} are mixed in each well.
4. Reaction is stopped by adding 20 μ L of Stop Reagent solution after 60 min.
 5. The fluorescence intensity values should be similar for all three controls (background value) and should be 2–3% of the values of the samples. The fluorescence values of the samples (with both cell extract and substrate) are subtracted from those recorded for *control II* (with cell extract and without substrate).

(c) Identification of hDAAO Inhibitors

The effect of small molecules on the enzymatic activity of hDAAO is evaluated by measuring the activity of purified recombinant hDAAO at different concentrations of inhibitors, in the presence of D-serine as substrate (Hopkins et al., 2013; Terry-Lorenzo et al., 2014, 2015).

1. The preincubation mixtures are prepared in Eppendorf-type test tubes mixing 25 μ L of 0.01 mg/mL hDAAO with 625 μ L of different inhibitors (in the 0–1.25 mM final concentration range), in 50 mM sodium phosphate buffer, pH 7.4, 4 μ M FAD.
2. The tubes are incubated for 30 min at room temperature.
3. Preparation of the reaction mixture: 130 μ L of preincubation mixture are transferred to each well of a microplate (four wells for each inhibitor's concentration tested) and 70 μ L of WS_C (22 mM D-serine, 4 μ M FAD, 0.1 mM Amplex[®] UltraRed Reagent, 0.2 U/mL HRP, in 50 mM sodium phosphate buffer, pH 7.4) are added.
4. Reaction is stopped by adding 40 μ L of Stop Reagent solution after 30 min.
5. The fluorescence value determined for the control in the absence of the compound is used as reference (activity = 100%).
6. Data are fit to a standard, four-parameter equation to determine curve top, bottom, inhibitor concentration giving 50% inhibition (IC₅₀), and Hill Slope:

$$y = \text{Bottom} + \frac{(\text{Top} - \text{Bottom})}{(1 + 10^{((\text{LogIC}_{50} - x) \times \text{Hill slope}))}}$$

where Bottom is the fluorescence value at the bottom plateau and Top is the fluorescence value at the top plateau (Terry-Lorenzo et al., 2015).

(D) Catalase and Purpald[®] Coupled Assay

The hydrogen peroxide produced by the DAAO reaction is reduced to water by catalase in the presence of methanol

as hydrogen donor, giving formaldehyde that, in turn, reacts with the chromogenic compound Purpald[®] to give the colored product 6-mercapto-*s*-triazolo-[4,3-*b*]-*s*-tetrazine showing an absorption maximum at 550 nm with an extinction coefficient of 7.74 mM⁻¹cm⁻¹ (Dickinson and Jacobsen, 1970; Watanabe et al., 1978; Sasabe et al., 2014; see **Figure 3**).

1. The reaction mixture is prepared in a plastic tube mixing 300 μ L of the D-amino acid solution, 200 μ L of 0.1 mM FAD, 300 μ L of 700 U/mL catalase, in 133 mM sodium pyrophosphate buffer, pH 8.3, and 100 μ L of 70% methanol (final volume: 900 μ L).
2. 100 μ L of DAAO solution are added (\sim 0.02 U).
3. The reaction mixture is incubated at 37°C for 15 min under agitation.
4. The reaction is stopped by adding 1 mL of 10% (w/v) trichloroacetic acid (final volume: 2 mL) followed by centrifugation at 700 \times g for 20 min.
5. A blank solution is prepared, adding 1 mL of 10% (w/v) trichloroacetic acid before the DAAO solution is added.
6. To 1 mL of the supernatant from step 4, 1 mL of 5 M KOH and 1 mL of 34 mM Purpald[®] in 0.5 M HCl are added.
7. The solution from step 6 is incubated at room temperature for 15 min, followed by the addition of 1 mL of 0.75% (w/v) KIO₄ in 0.2 M KOH with vigorous shaking.
8. 1 mL of the solution from step 7 is transferred in a plastic disposable microcuvette and the absorbance at 550 nm (Abs_{550 nm}) is measured.
9. The absorbance at 550 nm is corrected by subtracting the absorbance value of a blank solution and the DAAO activity is calculated from the Δ Abs_{550 nm} using the following equation (Watanabe et al., 1978):

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{2 \times \Delta \text{Abs}_{550 \text{ nm}} / 15 \text{ min}}{\epsilon_{550 \text{ nm}} (\text{chromogen}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where the molar extinction coefficient at 550 nm of 6-mercapto-*s*-triazolo-[4,3-*b*]-*s*-tetrazine is 7.74 mM⁻¹cm⁻¹ and mL_{TOT} is the final volume of the solution from step 7 (4 mL). The detection range is reported in **Table 3**. In this assay an average velocity corresponding to the first 15 min of the reaction is determined. A factor 2 is added at the numerator because only half of the initial sample (from step 4) is reacted with Purpald[®].

Application: Detection of DAAO Activity in Tissues

The enzymatic assay is performed on post mortem human brain tissues dissected after perfusion from the thoracic aorta with ice-cold phosphate buffer saline (PBS), pH 7.4 (Sasabe et al., 2014).

1. The tissue is homogenized in 7 mM sodium pyrophosphate buffer, pH 8.3, at 3,500 rpm for 2 min and centrifuged at 5,500 \times g for 10 min at 4°C. The supernatant (tissue lysate) is recovered and the pellet discarded.
2. 50 μ L tissue lysate is added to the reaction mixture containing 150 μ L of 100 mM D-alanine, 100 μ L 0.1 mM FAD, 150 μ L of

- 700 U/mL catalase in 133 mM sodium pyrophosphate buffer, pH 8.3, and 50 μ L of 70% (v/v) methanol.
- The reaction mixture is incubated at 37°C for 60 min under stirring.
 - Reaction is stopped by adding 500 μ L of 10% (w/v) trichloroacetic acid and the reaction mixture is centrifuged at 700 \times g for 20 min. The supernatant is recovered and the pellet discarded.
 - To 250 μ L of the supernatant solution 250 μ L of 5 M KOH and 250 μ L of 34 mM Purpald[®] in 0.5 M HCl are added.
 - The reaction mixture is incubated at room temperatures for 15 min to which 250 μ L 0.75% (w/v) KIO₄ in 0.2 M KOH is added under vigorous shaking.
 - The DAAO activity is measured on the supernatant and expressed as the amount of D-alanine oxidized per min per milligram of total proteins. Total protein concentration is determined by the method of (Lowry et al., 1951).

Determination of α -Keto Acid Production

A) Direct Spectrophotometric Assay

- Preparation of reaction mixture: 990 μ L of D-amino acid solution in 75 mM disodium pyrophosphate buffer, pH 8.5, are transferred into a quartz microcuvette.
- The initial absorbance value at the wavelength corresponding to the maximum absorption peak of the α -keto acid produced from the starting amino acid is recorded.
- Measurement of DAAO activity: 10 μ L of DAAO (\sim 0.01 U) are added and the increase of absorbance at the selected wavelength, due to α -keto acid production, is recorded. The DAAO activity is calculated from the Δ Abs/min using the following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{\Delta \text{Abs}_{x \text{ nm}} / \text{min}}{\epsilon_{x \text{ nm}} (\alpha\text{-keto acid}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

- where $\epsilon_{x \text{ nm}}$ is the extinction coefficient at the selected wavelength of the α -keto acid produced and mL_{TOT} is the final volume of the reaction mixture (1 mL). The detection range is reported in **Table 3**.
- When the extinction coefficient of the α -keto acid is unknown, a calibration curve is set up by recording the absorbance values at known concentrations of the α -keto acid.

(B) Indirect Assay with 2,4-Dinitrophenylhydrazine

- Preparation of DAAO reaction mixture: 10 μ L of DAAO (\sim 0.05 U) is added to the D-amino acid solution in 75 mM disodium pyrophosphate buffer, pH 8.5 (final volume = 0.3 mL), and the reaction mixture is incubated at 25°C for 10 min.
- Derivatization of α -keto acid produced by DAAO reaction with DNP (**Figure 4**): 0.15 mL of 1 mM DNP, solubilized in 1 M HCl, are added to the reaction mixture.
- The reaction is stopped by adding 1.05 mL of 0.6 M NaOH, followed by incubation at 25°C for 5 min to allow color development.
- A blank solution is prepared using the same procedure, but without substrate (D-amino acid).

- The absorbance at 445 nm of the reaction mixture is corrected by subtracting the absorbance value of blank solution and the DAAO activity is calculated from the Δ Abs_{445 nm} using the following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{\Delta \text{Abs}_{445 \text{ nm}} / 10 \text{ min}}{\epsilon_{445 \text{ nm}} (\alpha\text{-keto acid derivative}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where $\epsilon_{445 \text{ nm}}$ is the molar extinction coefficient at 445 nm of the α -keto acid derivative and mL_{TOT} is the final volume of the reaction mixture. In this assay an average velocity corresponding to the first 10 min of the reaction is determined. The detection range is reported in **Table 3**. If the extinction coefficient of the 2,4-dinitrophenylhydrazone derivative of the α -keto acid is unknown, a calibration curve is set up by recording the absorbance values obtained upon derivatization of different known concentrations of the α -keto acid with DNP under the same experimental conditions.

Determination of Ammonium Production: Glutamate Dehydrogenase and α -Ketoglutarate Coupled Assay

- Preparation of the reaction mixture: D-amino acid solution, 5 mM α -ketoglutarate, 0.25 mM NADH, and 10 U/mL of GDH, in 75 mM disodium pyrophosphate buffer, pH 8.5 (all final concentrations).
- 1 mL of this solution is transferred to a quartz microcuvette and the absorbance at 340 nm of the reaction mixture is recorded.
- Measurement of DAAO activity: 10 μ L of DAAO (\sim 1 U) are added and the decrease in absorbance at 340 nm, due to the conversion of NADH to NAD⁺ by GDH reaction, is proportional to the amount of ammonium produced by DAAO reaction. The DAAO activity is calculated from the Δ Abs_{340 nm}/min using the following equation:

$$\frac{U}{\text{mL}_{\text{DAAO}}} = \frac{\Delta \text{Abs}_{340 \text{ nm}} / \text{min}}{\epsilon_{340 \text{ nm}} (\text{NADH}) \times \text{mL}_{\text{DAAO}}} \times \text{mL}_{\text{TOT}}$$

where the molar extinction coefficient at 340 nm of NADH is 6.3 mM⁻¹cm⁻¹ and mL_{TOT} is the final volume of the reaction mixture. The detection range is reported in **Table 3**.

ANTICIPATED RESULTS

Determination of the Steady-State Kinetic Parameters and Substrate Specificity

The apparent (i.e., at fixed O₂ concentration) kinetic parameters ($k_{\text{cat,app}}$ and $K_{\text{m,app}}$ values) at air saturation ([O₂] = 0.253 mM at 25°C) are determined for the wild-type RgDAAO enzyme on various D-amino acids using the oxygen consumption assay (Sacchi et al., 2004). As shown in **Figure 7**, the chemical nature of the substrates profoundly influences the kinetic parameters; a 6- and 60-fold lower activity is apparent on basic and acidic

substrates (D-lysine and D-glutamate, respectively) in comparison to the value obtained on the neutral D-alanine (see Table 2).

Screening for RgDAAO Variants with Broader Substrate Specificity

In an attempt to evolve the substrate specificity of RgDAAO, five rounds of random mutagenesis (under different PCR reaction conditions) are performed (Sacchi et al., 2004): a rapid and sensitive assay is required to select the clone expressing the variant with the activity of interest. By using the colorimetric assay based on the reaction of RgDAAO on the D-amino acids of interest coupled to HRP and *o*-DNS, RgDAAO variants showing an altered substrate specificity could be identified (Sacchi et al., 2004). The crude extract of *E. coli* culture expressing clones 1-7 and 2-41, deriving from the first and the second set of PCR conditions, respectively, showed a similar activity on D-alanine, but an altered activity on D-arginine with respect to the wild-type enzyme (Figure 8). Furthermore, clones 4-903 and 5-249, deriving from the fourth and the fifth set of PCR conditions, respectively, expressed RgDAAO variants showing a higher activity on at least two of the three substrates tested (D-alanine, D-aspartate, and D-arginine) compared to the wild-type enzyme (see Figure 8). These four RgDAAO variants possess a broader substrate specificity and thus are useful for analytically determining the total content of D-amino acids.




FIGURE 7 | Michaelis-Menten plot of the kinetics of wild-type RgDAAO on D-alanine (black), D-lysine (gray), and D-glutamate (light gray and inset) as substrate, at pH 8.5 and 25°C. All measurements were performed at air saturation ($[O_2] = 0.253$ mM) using the oxygen consumption assay.

TABLE 2 | Comparison of the apparent kinetic parameters of RgDAAO on different substrates.

	V_{max} (U/mg)	K_m (mM)	V_{max}/K_m
D-Alanine	100 ± 3	1.3 ± 0.2	76.9
D-Lysine	16.4 ± 0.5	14.6 ± 1.4	1.1
D-Glutamate	1.7 ± 0.1	76.4 ± 6.4	0.02

Identification of hDAAO Inhibitors

The Amplex[®] UltraRed-based fluorescence enzyme assay is a useful and adequate screening technique to identify novel inhibitors of hDAAO. As stated in the Introduction section, this flavoenzyme controls the levels of D-serine in the brain, the coagonist of the NMDARs, a central regulator of synaptic plasticity, learning, and memory (Panatier et al., 2006; Zhou and Sheng, 2013); the inhibitors of hDAAO have the potential to treat nervous system disorders such as schizophrenia (Sacchi et al., 2013). Computational tools are used to identify libraries of compounds, then screened for hDAAO inhibition. The affinity of the different molecules is based on IC_{50} values (i.e., the compound concentration at which hDAAO activity is halved). As shown in Figure 9, all compounds tested inhibited the enzymatic




FIGURE 8 | Comparison of DAAO activity determined by the coupled reaction with HRP and *o*-DNS and expressed as absorbance change at 440 nm of wild-type and variants RgDAAO identified following five rounds of random mutagenesis and determined on D-alanine (black), D-aspartate (gray), and D-arginine (light gray) as substrates (30 mM final concentration). Figure obtained on the basis of the results published in Sacchi et al. (2004).





FIGURE 9 | Enzyme inhibition assays performed by the Amplex[®] UltraRed method using recombinant hDAAO, D-serine as substrate and three different compounds (1, black; 2, gray; and 3, light gray) as potential inhibitors. The plot displays the mean values ± SD, $n = 4$.

activity, although compound 1 was a more potent inhibitor compared to compounds 2 and 3 ($IC_{50} = 1.7 \pm 0.1$, 10.6 ± 0.9 and $43.5 \pm 2.3 \mu\text{M}$, respectively).

Detection of DAAO Activity in Tissues

Owing to the involvement of hDAAO in neuropsychiatric disorders (see Introduction), the investigation of its distribution based on its activity in the human CNS is of utmost importance. Thus, DAAO activity in post mortem samples from different regions of the CNS of mice and humans was determined using a coupled assay with catalase and Purpald[®] (Sasabe et al., 2014). DAAO activity was detected in all the human specimens examined; significant activity was found along the corticospinal tract, rubrospinal tract, nigrostriatal system, ponto-/olivo-cerebellar fibers, and in the anterolateral system. DAAO activity was also detected in human primary somatosensory cortex (PSC), in primary motor cortex (PMC), and in the internal capsule, whereas in mice there was no such activity in the cerebral cortex and in the striatum (Figure 10).



CONCLUSIONS

It is useful to determine DAAO activity for different biotechnological and medical applications (Pollegioni et al., 2007b, 2008; Pollegioni and Molla, 2011): (i) investigation of kinetic parameters on different substrates in order to identify the optimized enzyme for a specific industrial process (i.e., synthesis of semi-synthetic antibiotics, resolution of racemic mixtures, etc.); (ii) identification of novel inhibitors of the human DAAO to be used as potential drugs for the treatment of neuropsychiatric disorders; (iii) evolution of DAAO variants more active at low oxygen concentration for the application in anticancer enzyme therapy; and (iv) development of enzyme-based biosensors for the analytical determination of D-amino acids content in biological fluids.

DAAO activity can be determined using direct or coupled assays which differ in terms of sensitivity, required instrumentation and cost, as described in this paper and summarized in Table 3.

NOTES

- Assay buffers must be pre-equilibrated at the working temperature.
- Instruments with similar specifications of the ones listed in “Materials and Equipment” section can be used.
- The reaction rate can be increased by incubating the reaction mixtures at higher temperature (e.g., 37°C instead of 25°C), taking in account the thermal stability of the enzymes in the experimental conditions.
- In absorbance and fluorescence assays, a sample containing buffer instead of the substrate and another one containing buffer instead of the enzyme must be performed as negative controls to exclude the presence of artifacts. The presence of contaminants in the sample might affect the absorbance or fluorescence signals (e.g., the presence of nucleic acids and/or proteins can increase the absorbance intensity in the UV region).
- In coupled assays with HRP or GDH the production of the expected α -keto acid is not directly monitored but it can be verified by different analytical methods (e.g. HPLC, NMR, and MS).

TABLE 3 | Comparison of different assays of DAAO activity.

Assay	Type of assay	Type of detection	Sensitivity range (mU _{DAAO} /mL)	Cost per assay
Molecular oxygen consumption	Direct, continuous	Amperometry	10–100	\$
Redox dye reduction	Direct, continuous	Abs at different wavelengths	500–5000	\$
o-DNS	Coupled (HRP), continuous	Abs at 440 nm	5–25	\$\$
4-AAP	Coupled (HRP), continuous	Abs at 505 nm	25–100	\$\$
Amplex [®] UltraRed	Coupled (HRP), discontinuous	Fluorescence (λ_{exc} at 535 nm, λ_{em} at 590 nm)	0.25–2.5	\$\$\$
Purpald [®]	Coupled (catalase), discontinuous	Abs at 550 nm	100–500	\$\$
α -Keto acid	Direct, continuous	Abs at λ_{max}	5–25	\$
DNP	Chemical derivatization, discontinuous	Abs at 445 nm	100–500	\$
Ammonium production	Coupled (GDH), continuous	Abs at 340 nm	500–2,500	\$\$

- As a general rule, continuous assays should be preferred (if possible), since the initial rate estimation is easier, due to the fact that the shape of the progress curve is visible all along the assay time and any discrepancy from the linearity is thus promptly detected. On the contrary, in discontinuous assays, to assess the linearity of the initial rate, withdrawal of different samples at fixed incubation times is necessary. Furthermore, the termination of the reaction and withdrawal of samples at fixed-time can introduce time and volume inaccuracy.

REFERENCES

- Brugger, D., Krondorfer, I., Zahma, K., Stoisser, T., Bolivar, J. M., Nidetzky, B., et al. (2014). Convenient microtiter plate-based, oxygen-independent activity assays for flavin-dependent oxidoreductases based on different redox dyes. *Biotechnol. J.* 9, 474–482. doi: 10.1002/biot.201300336
- Dickinson, R. G., and Jacobsen, N. W. (1970). A new sensitive and specific test for the detection of aldehydes: formation of 6-mercapto-3-substituted-s-triazolo-[4,3-b]-s-tetrazines. *J. Chem. Soc. D. Chem. Commun.* 24, 1719–1720. doi: 10.1039/c29700001719
- Fantinato, S., Pollegioni, L., and Pilone, M. S. (2001). Engineering, expression and purification of a His-tagged chimeric D-amino acid oxidase from *Rhodotorula gracilis*. *Enzyme Microb. Technol.* 29, 407–412. doi: 10.1016/S0141-0229(01)00400-8
- Fratini, L. F., Piubelli, L., Sacchi, S., Molla, G., and Pollegioni, L. (2011). Is rat an appropriate animal model to study the involvement of D-serine catabolism in schizophrenia? Insights from characterization of D-amino acid oxidase. *FEBS J.* 278, 4362–4373. doi: 10.1111/j.1742-4658.2011.08354.x
- Golubev, I. V., Komarova, N. V., Ryzhenkova, K. V., Chubar, T. A., Savin, S. S., and Tishkov, V. I. (2014). Study of the structure-function-stability relationships in yeast D-amino acid oxidase: hydrophobization of alpha-helices. *Acta Nat.* 6, 76–88.
- Holme, D. J., and Goldberg, D. M. (1975). Coupled optical rate determination of amino acid oxidase activity. *Biochim. Biophys. Acta* 377, 61–70. doi: 10.1016/0005-2744(75)90286-7
- Hopkins, S. C., Heffernan, M. L., Saraswat, L. D., Bowen, C. A., Melnick, L., Hardy, L. W., et al. (2013). Structural, kinetic, and pharmacodynamic mechanisms of D-amino acid oxidase inhibition by small molecules. *J. Med. Chem.* 56, 3710–3724. doi: 10.1021/jm4002583
- Job, V., Marcone, G. L., Pilone, M. S., and Pollegioni, L. (2002). Glycine oxidase from *Bacillus subtilis*. Characterization of a new flavoprotein. *J. Biol. Chem.* 277, 6985–6993. doi: 10.1074/jbc.M111095200
- Kawazoe, T., Tsuge, H., Pilone, M. S., and Fukui, K. (2006). Crystal structure of human D-amino acid oxidase: context-dependent variability of the backbone conformation of the VAAGL hydrophobic stretch located at the si-face of the flavin ring. *Protein Sci.* 15, 2708–2717. doi: 10.1110/ps.062421606
- Lowry, O. H., Rosebrough, N. J., Farr, A. L., and Randall, R. J. (1951). Protein measurement with the Folin phenol reagent. *J. Biol. Chem.* 193, 265–275.
- Mattevi, A., Vanoni, M. A., Todone, F., Rizzi, M., Teplyakov, A., Coda, A., et al. (1996). Crystal structure of D-amino acid oxidase: a case of active site mirror-image convergent evolution with flavocytochrome b2. *Proc. Natl. Acad. Sci. U.S.A.* 93, 7496–7501. doi: 10.1073/pnas.93.15.7496
- Mitchell, J., Paul, P., Chen, H. J., Morris, A., Payling, M., Falchi, M., et al. (2010). Familial amyotrophic lateral sclerosis is associated with a mutation in D-amino acid oxidase. *Proc. Natl. Acad. Sci. U.S.A.* 107, 7556–7561. doi: 10.1073/pnas.0914128107
- Mizutani, H., Miyahara, I., Hirotsu, K., Nishina, Y., Shiga, K., Setoyama, C., et al. (1996). Three-dimensional structure of porcine kidney D-amino acid oxidase at 3.0 Å resolution. *J. Biochem.* 120, 14–17. doi: 10.1093/oxfordjournals.jbchem.a021376
- Molla, G., Sacchi, S., Bernasconi, M., Pilone, M. S., Fukui, K., and Pollegioni, L. (2006). Characterization of human D-amino acid oxidase. *FEBS Lett.* 580, 2358–2364. doi: 10.1016/j.febslet.2006.03.045

AUTHOR CONTRIBUTIONS

All authors listed have made a substantial, direct and intellectual contribution to the work, and approved it for publication.

ACKNOWLEDGMENTS

The work of LP was supported by Fondo di Ateneo per la Ricerca.

- Molla, G., Vegezzi, C., Pilone, M. S., and Pollegioni, L. (1998). Overexpression in *Escherichia coli* of a recombinant chimeric *Rhodotorula gracilis* D-amino acid oxidase. *Protein Expr. Purif.* 14, 289–294. doi: 10.1006/prep.1998.0956
- Nagasawa, M., Ogino, Y., Kurata, K., Otsuka, T., Yoshida, J., Tomonaga, S., et al. (2012). Hypothesis with abnormal amino acid metabolism in depression and stress vulnerability in Wistar Kyoto rats. *Amino Acids* 43, 2101–2111. doi: 10.1007/s00726-012-1294-y
- Nagata, Y., Shimojo, T., and Akino, T. (1988). Two spectrophotometric assays for D-amino acid oxidase: for the study of distribution patterns. *Int. J. Biochem.* 20, 1235–1238. doi: 10.1016/0020-711X(88)90225-X
- Panatier, A., Theodosis, D. T., Mothet, J. P., Touquet, B., Pollegioni, L., Poulain, D. A., et al. (2006). Glia-derived-serine controls NMDA receptor activity and synaptic memory. *Cell* 125, 775–784. doi: 10.1016/j.cell.2006.02.051
- Pollegioni, L., and Molla, G. (2011). New biotech applications from evolved D-amino acid oxidases. *Trends Biotechnol.* 29, 276–283. doi: 10.1016/j.tibtech.2011.01.010
- Pollegioni, L., and Sacchi, S. (2010). Metabolism of the neuromodulator D-serine. *Cell. Mol. Life Sci.* 67, 2387–2404. doi: 10.1007/s00018-010-0307-9
- Pollegioni, L., Caldinelli, L., Molla, G., Sacchi, S., and Pilone, M. S. (2004). Catalytic properties of D-amino acid oxidase in cephalosporin C bioconversion: a comparison between proteins from different sources. *Biotechnol. Prog.* 20, 467–473. doi: 10.1021/bp034206q
- Pollegioni, L., Molla, G., Sacchi, S., Rosini, E., Verga, R., and Pilone, M. S. (2008). Properties and applications of microbial D-amino acid oxidases: current state and perspectives. *Appl. Microbiol. Biotechnol.* 78, 1–16. doi: 10.1007/s00253-007-1282-4
- Pollegioni, L., Piubelli, L., Sacchi, S., Pilone, M. S., and Molla, G. (2007a). Physiological functions of D-amino acid oxidases: from yeast to humans. *Cell. Mol. Life Sci.* 64, 1373–1394. doi: 10.1007/s00018-007-6558-4
- Pollegioni, L., Sacchi, S., Caldinelli, L., Boselli, A., Pilone, M. S., Piubelli, L., et al. (2007b). Engineering the properties of D-amino acid oxidases by a rational and a directed evolution approach. *Curr. Protein Pept. Sci.* 8, 600–618. doi: 10.2174/138920307783018677
- Rosini, E., Molla, G., Ghisla, S., and Pollegioni, L. (2010). On the reaction of D-amino acid oxidase with dioxygen: O₂ diffusion pathways and enhancement of reactivity. *FEBS J.* 278, 482–492. doi: 10.1111/j.1742-4658.2010.07969.x
- Rosini, E., Molla, G., Rossetti, C., Pilone, M. S., Pollegioni, L., and Sacchi, S. (2008). A biosensor for all D-amino acids using evolved D-amino acid oxidase. *J. Biotechnol.* 135, 377–384. doi: 10.1016/j.jbiotec.2008.06.001
- Rosini, E., Pollegioni, L., Ghisla, S., Orrù, R., and Molla, G. (2009). Optimization of D-amino acid oxidase for low substrate concentrations – towards a cancer enzyme therapy. *FEBS J.* 276, 4921–4932. doi: 10.1111/j.1742-4658.2009.07191.x
- Ross, C. A., Margolis, R. L., Reading, S. A., Pletnikov, M., and Coyle, J. T. (2006). Neurobiology of schizophrenia. *Neuron* 52, 139–153. doi: 10.1016/j.neuron.2006.09.015
- Sacchi, S., Caldinelli, L., Cappelletti, P., Pollegioni, L., and Molla, G. (2012). Structure–function relationships in human D-amino acid oxidase. *Amino acids* 43, 1833–1850. doi: 10.1007/s00726-012-1345-4
- Sacchi, S., Cappelletti, P., Giovannardi, S., and Pollegioni, L. (2011). Evidence for the interaction of D-amino acid oxidase with pLG72 in a glial cell line. *Mol. Cell. Neurosci.* 48, 20–28. doi: 10.1016/j.mcn.2011.06.001

- Sacchi, S., Rosini, E., Molla, G., Pilone, M. S., and Pollegioni, L. (2004). Modulating D-amino acid oxidase substrate specificity: production of an enzyme for analytical determination of all D-amino acids by directed evolution. *Prot. Eng. Des. Sel.* 17, 517–525. doi: 10.1093/protein/ghz064
- Sacchi, S., Rosini, E., Pollegioni, L., and Molla, G. (2013). D-Amino acid oxidase inhibitors as a novel class of drugs for schizophrenia therapy. *Curr. Pharm. Des.* 3119, 2499–2511. doi: 10.2174/1381612811319140002
- Sasabe, J., Suzuki, M., Imanishi, N., and Aiso, S. (2014). Activity of D-amino acid oxidase is widespread in the human central nervous system. *Front. Synaptic Neurosci.* 6:14. doi: 10.3389/fnsyn.2014.00014
- Tedeschi, G., Pollegioni, L., and Negri, A. (2012). “Assays of D-amino acid oxidases” in *Unnatural Amino Acids*, eds L. Pollegioni and S. Servi (New York, NY: Humana Press), 381–395.
- Terry-Lorenzo, R. T., Chun, L. E., Brown, S. P., Heffernan, M. L., Fang, Q. K., Orsini, M. A., et al. (2014). Novel human D-amino acid oxidase inhibitors stabilize an active-site lid-open conformation. *Biosci. Rep.* 34, 487–499. doi: 10.1042/BSR20140071
- Terry-Lorenzo, R. T., Masuda, K., Sugao, K., Fang, Q. K., Orsini, M. A., Sacchi, S., et al. (2015). High-throughput screening strategy identifies allosteric, covalent human D-amino acid oxidase inhibitor. *J. Biomol. Screen.* 20, 1218–1231. doi: 10.1177/1087057115600413
- Umhau, S., Pollegioni, L., Molla, G., Diederichs, K., Welte, W., Pilone, M. S., et al. (2000). The X-ray structure of D-amino acid oxidase at very high resolution identifies the chemical mechanism of flavin-dependent substrate dehydrogenation. *Proc. Natl. Acad. Sci. U.S.A.* 97, 12463–12468. doi: 10.1073/pnas.97.23.12463
- Wang, S. J., Yu, C. Y., Lee, C. K., Chern, M. K., and Kuan, I. C. (2008). Subunit fusion of two yeast D-amino acid oxidases enhances their thermostability and resistance to H₂O₂. *Biotechnol. Lett.* 30, 1415–1422. doi: 10.1007/s10529-008-9694-5
- Watanabe, T., Motomura, Y., and Suga, T. (1978). A new colorimetric determination of D-amino acid oxidase and urate oxidase activity. *Anal. Biochem.* 86, 310–315. doi: 10.1016/0003-2697(78)90347-0
- Wong, K. S., Fong, W. P., and Tsang, P. W. K. (2010). A single Phe54Tyr substitution improves the catalytic activity and thermostability of *Trigonopsis variabilis* D-amino acid oxidase. *New Biotechnol.* 27, 78–84. doi: 10.1016/j.nbt.2009.11.002
- Wu, S., Basile, A. S., and Barger, S. W. (2007). Induction of serine racemase expression and D-serine release from microglia by secreted amyloid precursor protein (sAPP). *Curr. Alzheimer Res.* 4, 243–251. doi: 10.2174/156720507781077241
- Zhou, Q., and Sheng, M. (2013). NMDA receptors in nervous system diseases. *Neuropharmacology* 74, 69–75. doi: 10.1016/j.neuropharm.2013.03.030

Conflict of Interest Statement: The authors declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

Copyright © 2018 Rosini, Caldinelli and Piubelli. This is an open-access article distributed under the terms of the Creative Commons Attribution License (CC BY). The use, distribution or reproduction in other forums is permitted, provided the original author(s) or licensor are credited and that the original publication in this journal is cited, in accordance with accepted academic practice. No use, distribution or reproduction is permitted which does not comply with these terms.