


RESEARCH

Open Access



An innovative diagnostic technology for the codon mutation C580Y in *kelch13* of *Plasmodium falciparum* with MinION nanopore sequencer

Kazuo Imai^{1,2}, Norihito Tarumoto^{1,2}, Lucky Ronald Runtuwene³, Jun Sakai^{1,2}, Kyoko Hayashida⁴, Yuki Eshita^{4,5,6,7}, Ryuichiro Maeda⁸, Josef Tuda⁹, Hideaki Ohno¹⁰, Takashi Murakami^{2,11}, Shigefumi Maesaki^{1,2}, Yutaka Suzuki³, Junya Yamagishi^{4,12} and Takuya Maeda^{2,11*} 

Abstract

Background: The recent spread of artemisinin (ART)-resistant *Plasmodium falciparum* represents an emerging global threat to public health. In Southeast Asia, the C580Y mutation of *kelch13* (*k13*) is the dominant mutation of ART-resistant *P. falciparum*. Therefore, a simple method for the detection of C580Y mutation is urgently needed to enable wide-spread routine surveillance in the field. The aim of this study is to develop a new diagnostic procedure for the C580Y mutation using loop-mediated isothermal amplification (LAMP) combined with the MinION nanopore sequencer.

Results: A LAMP assay for the *k13* gene of *P. falciparum* to detect the C580Y mutation was successfully developed. The detection limit of this procedure was 10 copies of the reference plasmid harboring the *k13* gene within 60 min. Thereafter, amplicon sequencing of the LAMP products using the MinION nanopore sequencer was performed to clarify the nucleotide sequences of the gene. The C580Y mutation was identified based on the sequence data collected from MinION reads 30 min after the start of sequencing. Further, clinical evaluation of the LAMP assay in 34 human blood samples collected from patients with *P. falciparum* malaria in Indonesia revealed a positive detection rate of 100%. All LAMP amplicons of up to 12 specimens were simultaneously sequenced using MinION. The results of sequencing were consistent with those of the conventional PCR and Sanger sequencing protocol. All procedures from DNA extraction to variant calling were completed within 3 h. The C580Y mutation was not found among these 34 *P. falciparum* isolates in Indonesia.

Conclusions: An innovative method combining LAMP and MinION will enable simple, rapid, and high-sensitivity detection of the C580Y mutation of *P. falciparum*, even in resource-limited situations in developing countries.

Keywords: Malaria, *Plasmodium falciparum*, Artemisinin resistance, LAMP, Nanopore sequencer, MinION™, *kelch 13*

*Correspondence: t_maeda@saitama-med.ac.jp

¹¹ Department of Microbiology, Saitama Medical University, 38 Morohongo, Moroyama-machi, Iruma-gun, Saitama 350-0495, Japan
Full list of author information is available at the end of the article



Background

Malaria is a mosquito-borne infectious disease and one of the most severe and prevalent public health problems in tropical and subtropical areas worldwide. In recent years, the spread of artemisinin (ART)-resistant *Plasmodium falciparum* malaria, which is defined as delayed parasite clearance from the blood following appropriate treatment with ART monotherapy or ART-based combination therapy (ACT) [1], has gained global attention, particularly in Southeast Asia and southern China [2–4]. The World Health Organization (WHO) currently recommends ACT instead of artesunate monotherapy as the first-line treatment for suspected *P. falciparum* malaria infections, and has urgently proposed a global plan for the monitoring and surveillance of ART-resistant *P. falciparum* to assess the threat of emerging resistant strains [5].

A molecular marker of ART resistance in *P. falciparum* has not been completely identified or validated; however, mutations in the *kelch13* (*k13*)-propeller gene are thought to be associated with ART resistance based on in vitro and in vivo investigations [1]. Although amino acid substitutions N458Y, Y493H, R539T, I543T, R561H, and C580Y could all be associated with ART-resistance and reduce the cure rate following ACT treatment, C580Y is thought to be the dominant mutation and a potential drug resistance marker, especially in Southeast Asia [6–9]. Recently, ART-resistant strains of *P. falciparum* with C580Y mutation have spread across Cambodia, north-eastern Thailand, and southern Laos. Simultaneously, the mutation could have developed resistance to piperaquine—one of the ACT partner drugs—because this drug is related to high rates of failure in the treatment of malaria [8]. Thus, there is an urgent need to strengthen the surveillance and elimination of ART-resistant *P. falciparum* C580Y mutation in Southeast Asia.

Point mutations and their prevalence can be confirmed by polymerase chain reaction (PCR) amplification and Sanger sequencing, which are widely used and well established. However, these techniques are unsuitable in malaria-endemic areas due to their high cost and need for sophisticated laboratories, and they cannot keep up with the current demand for field-based assays. Clinical samples need to be stored by appropriate protocols and transferred to laboratories from remote locations, which creates several disadvantages, including the cost of transportation, transfer time, and potential loss of clinical samples. Therefore, a simple and easy-to-perform molecular genotyping method is urgently required to accurately assess the distribution of ART-resistant *P. falciparum* in malaria-endemic areas.

Loop-mediated isothermal amplification (LAMP) is an isothermal nucleic acid amplification technique that

has some operational advantages over conventional PCR procedures because it is fast, easy to use, and cost-effective [10]. Amplification is mediated by a loop structure and relies on an auto-cycling procedure performed with strand-displacement of *Bst* DNA polymerase under isothermal conditions. Furthermore, the amplification products in LAMP reactions can be directly judged by turbidity or changes in fluorescence with the naked eye. Given these properties, the LAMP method is expected to become an important and widely-used clinical diagnostic technique in point-of-care testing for infectious diseases and will require only limited equipment and labor. The advantages of LAMP are particularly relevant in areas with limited resources [11, 12]. Recently, Loopamp™ MALARIA Pan/Pf Detection Kit (Eiken Chemical Co., Ltd., Tokyo, Japan) has become commercially available for the detection of *Plasmodium* parasites.

The MinION nanopore sequencer is a pocket-sized and USB-connected portable real-time sequencer developed by Oxford Nanopore Technologies (ONT; Oxford, UK). The advantages of MinION include the simple and rapid preparation of samples, portability, real-time sequencing, relatively low cost, and the need for minimal equipment and personnel. Recently, MinION has been shown to be an alternative method for whole genome sequencing that can be used for genomic surveillance and infection control in hospitals by enabling the rapid identification of pathogens, including Ebola virus, in areas with limited resources [13, 14]. In addition, an innovative diagnostic technology using the LAMP assay combined with MinION, which is rapid, simple, highly sensitive, and cost-effective, has been introduced for the serotyping of dengue virus and genotyping of *Plasmodium* parasites [15, 16].

In this study, a single nucleotide polymorphism genotyping method for C580Y in *k13* of *P. falciparum* with clinical specimens using the LAMP assay combined with MinION nanopore sequencing was developed. This is a quick, sensitive, and simple method that requires minimal equipment and personnel.

Methods

Design of LAMP primers

The oligonucleotide LAMP primers for the detection of the codon mutation of C580Y were designed based on the sequences of *k13* to encompass C580Y within the F1–B1 primer pairs via online LAMP primer design software (PrimerExplorer 5, <http://primerexplorer.jp/index.html>; Eiken Chemical). Based on in silico analysis of the *k13* sequence in six human *Plasmodium* parasites (*P. falciparum*, *Plasmodium vivax*, *Plasmodium ovale curtisi*, *Plasmodium ovale wallikeri*, *Plasmodium knowlesi*, and *Plasmodium malariae*), these primers were located on

diverse sequences for species-specific LAMP amplification (Additional file 1: Fig. S1). The sequences of each selected primer are given in Table 1 and their positions are shown in Fig. 1.

Plasmid construction

The templates used for the analytical LAMP reactions were the pEX-A2J1 plasmids harbouring partial sequences within the *k13* of *P. falciparum*, which were constructed by Eurofines Genomics Co., Ltd. (Tokyo, Japan) based on the reference sequences (GenBank accession numbers: *P. falciparum*, KT956001.1 and KR537460.1). The sequences of the constructed plasmids are summarized in Fig. 1. The plasmids were serially diluted tenfold and adjusted from 1.0×10^1 to 1.0×10^5 copies/ μ L to determine the detection threshold and specificity of the LAMP reactions. The sequences of the constructed plasmids were confirmed by the Sanger method using F3 and B3 primers (Table 1).

LAMP reactions

The reactions were performed using a Loopamp DNA Amplification Reagent Kit (Eiken Chemical). Detection of the LAMP amplicons was performed by real-time measurements of turbidity using a Loopamp real-time turbidimeter LA-200 (Eiken Chemical) and visual observations of the colour changes with the naked eye under natural light. The final reaction volume was 25 μ L, comprising 40 pmol of the FIP primer and 40 pmol of the BIP primer, 20 pmol of each loop primer, 5 pmol of F3 and B3 primers, 1 μ L of the Bst DNA polymerase, reaction buffer (20 mM Tris-HCl, 10 mM KCl, 8 mM MgSO₄, 10 mM (NH₄)₂SO₄, 0.1% Tween 20, 0.8 M betaine, and 1.4 mM each dNTP), and 1 μ L of plasmid DNA or 8 μ L of genomic DNA samples. Reactions were carried out at 62,

63, 64, and 65 °C in duplicate to find the optimum temperature for the LAMP amplification.

The detection limit and specificity of the LAMP assay

Ten-fold serial dilutions of each plasmid DNA ranging from 1.0×10^1 to 1.0×10^5 copies/ μ L were produced to determine the detection limit of the LAMP primers, and stored at -20 °C in duplicate until use. Additionally, to evaluate the specificity of the test, seven different protozoa (*Toxoplasma gondii*, *Cryptosporidium parvum*, *Giardia intestinalis*, *Entamoeba histolytica*, *Leishmania donovani*, *Trypanosoma brucei rhodesiense*, and *Trypanosoma cruzi*) isolated by the National BioResource Project (<http://www.nbrp.jp/>) were also examined.

Clinical samples and ethics

Clinical samples were obtained from Sam Ratulangi University in Manado and Bitung, North Sulawesi, Indonesia from August to December 2010. All patients were diagnosed based on symptoms, blood smears stained with Giemsa stain, and nested PCR according to previously reported protocols [17–19]. Peripheral blood was collected by FTA Elute cards (GE Healthcare Life Sciences, Little Chalfont, UK) from each patient and stored at room temperature.

Both the design and protocol of this study conformed to the Helsinki Declaration and were approved by the Institutional Ethics Committee. All samples were collected after written informed consent has been obtained. The methods of collection and analysis of the human samples were approved and cleared by the Institutional Ethical Review Board of Sam Ratulangi University and the University of Tokyo.

DNA extraction from clinical samples

Genomic DNA was extracted following the manufacturer’s instructions for FTA Elute cards. Briefly, each disk with a diameter of 3.0 mm was cut from the blood spot areas and washed three times with 200 μ L of distilled water. DNA was eluted in 30 μ L of TE buffer (10 mM Tris-HCl, 0.1 mM EDTA, pH 8.0) with a heat block at 95 °C for 30 min. Finally, 8.0 μ L of supernatant was used as template DNA for both LAMP assays.

Sequencing of LAMP products by MinION

To confirm the presence of codon mutation C580Y in *k13*, amplicon sequencing of LAMP reactions by MinION sequencer was performed. LAMP amplicons derived from plasmid DNA harboring partial sequences of *k13* with or without C580Y mutation and clinical samples were purified with Agencourt AMPure XP beads (Beckman Coulter Inc., Brea, CA). For the multiplex and real-time MinION sequencing, a Rapid Barcoding

Table 1 Nucleotide sequences of the LAMP primers constructed for *kelch13* of *Plasmodium falciparum*

Primer	Sequence (5’-3’)	Length
LAMP primers		
F3	TGGGGGATATGATGGCTC	18
B3	ATTATCAATACCTCCAACACAT	23
FIP	AGCTGATGATCTAGGGGTATTCAA-TTC TATTATACCGAATGTAGAAGCA	49
BIP	AATGGGAACAATTTCCATATGCCT-GAT TAAGGTAATTAAGCTGCTC	48
LF	CCCATGCTTTCATACGATGATCATA	25

F3 and B3, outer primers; FIP and BIP, inner primers; LF, loop primer
The FIP primer consisted of F2 and the complementary strand (F1c)
The BIP primer consisted of B2 and the complementary strand (B1c)




Fig. 1 Alignment of the partial sequences of the *kelch13* of *Plasmodium falciparum*, which was constructed within the pEX-A2J1 plasmid. The constructed sets of LAMP primers are shown as lines and boxes. Asterisks show the specific sequences located at the codon position of C580Y in artemisinin-resistant *P. falciparum*

Sequencing Kit (SQK-RBK001, ONT) was used to sequence up to 12 samples on a single flow cell. LAMP amplicons (200 ng each) were measured by Qubit 3.0 Fluorometer (Thermo Fisher Scientific) and simultaneously processed for the barcoding and library preparation using the Rapid Barcoding Sequencing Kit, according to the manufacturer’s instructions. After the MinION Platform QC run, the DNA library was loaded into

MinION Flow Cell (FLO-MIN107 R9.5 Version) and the “NC_48Hr_sequencing_FLO-MIN107_SQK-RBK001_plus_basecaller” protocol was initiated using MinKNOW software (ONT, v1.4.2).

Analysis of MinION data

Figure 2 shows the workflow of the data analysis in this study. Local basecalling was performed using




Fig. 2 Schematic of the data analysis workflow in this study

MinKNOW in real time and automatically. FAST5 reads were collected at 30 min and 48 h after the start of sequencing. FASTQ files were extracted from each collected FAST5 read via Poretools (v 0.6.0). Demultiplexing and adapter trimming was performed using Porechop (v 0.2.2), and each FASTQ per ONT-barcode data was mapped by BWA-MEM (v 0.7.15) using the sequence of *k13* without the C580Y codon mutation between a pair of LAMP primers, F1c (5'-TTG AAT ACC CCT AGA TCA TCA GC-3') and B1c (5'-AAT GGG AAC AAT TTC CAT ATG CCT-3'). The mapped data with a "Map-Quality" < 60 were discarded via Samtools (v 1.5.0). The single nucleotide variant calling from mapped data was obtained via Samtools and Bcftools (v 1.5.0), and cutoff values of 50 for single nucleotide variant quality and read depth were determined. Mapped data were visualized by IGV software (v 2.3.8). FAST5 reads and mapped data were analysed via Poretools and Samtools, respectively.

k13 amplification and Sanger sequencing

The Sanger method was used to perform PCR amplicon sequencing to confirm and validate the sequence results of LAMP amplicons via the MinION sequencer for the clinical specimens. For the specific amplification of *k13*, nested PCR according to previously reported protocols were used [7]. Briefly, the first- and second-round PCR amplicons were generated by PrimeSTAR GXL DNA Polymerase (Takara Bio Inc., Kusatsu, Japan) with the following primers: first round forward primers, 5'-CGG AGT GAC CAA ATC TGG GA-3' and reverse primers,




Fig. 3 The detection limit of LAMP for *kelch13* of *Plasmodium falciparum* with tenfold serial dilutions of plasmid DNA

5'-GGG AAT CTG GTG GTA ACA GC-3'; second round forward primers, 5'-GCC AAG CTG CCA TTC ATT TG-3' and reverse primers, 5'-GCC TTG TTG AAA GAA GCA GA-3'. Thermal cycling was carried out under the following conditions: 94 °C for 2 min, followed by 40 cycles at 98 °C for 10 s, 55 °C for 30 s, 68 °C for 1 min, with a final extension at 68 °C for 7 min in both the first- and second-round amplifications. The 1000-fold dilution products of the first round of PCR were used as templates for the second round. All second-round PCR products were analyzed using 1.2% (w/v) agarose gel, stained with ethidium bromide, and purified by Agencourt AMPure XP beads (Beckman Coulter). Sanger sequencing of these purified PCR amplicons was performed by Eurofines Genomics Co., Ltd. (Tokyo, Japan).

Results

Sensitivity and specificity of LAMP

A tenfold serial dilution of each plasmid DNA was amplified to determine the lower detection limit of the constructed LAMP assay and the optimal temperature was observed at 62 °C. Figure 3 shows the results of detection based on real-time turbidity; amplification of the target DNA is indicated by the rising curve. The minimum amounts of plasmid DNA on real-time turbidities with LA-200 were 1.0×10^1 /reaction within 60 min (Fig. 2). These results were similar to those obtained by observing visual changes within 60 min. Additionally, 1.0 ng of genomic DNA derived from seven other protozoa (*Toxoplasma gondii*, *Cryptosporidium parvum*, *Giardia intestinalis*, *Entamoeba histolytica*, *Leishmania donovani*, *Trypanosoma brucei rhodesiense*, and *Trypanosoma cruzi*) were not amplified by the LAMP procedure.

Clinical evaluation of LAMP assays

Parasitemia was confirmed in 40 blood samples by nested-PCR (33 *P. falciparum*, 6 *P. vivax*, and 1

Table 2 Summary of the results of nested PCR and LAMP amplifications, and the identification of *kelch13* mutations of *P. falciparum*

		LAMP reaction		Identification of mutation			
		Positive	Negative	PCR and Sanger		LAMP and MinION	
				C580Y	Wild type	C580Y	Wild type
Nested PCR							
<i>P. f.</i>	33	33	0	0	33	0	33
<i>P. f./P. v.</i>	1	1	0	0	1	0	1
<i>P. v.</i>	6	0	6	–	–	–	–

P. f., *P. falciparum*; *P. v.*, *P. vivax*; *P. f./P. v.*, *P. falciparum* and *P. vivax* co-infection

co-infection with *P. falciparum* and *P. vivax*). All 33 *P. falciparum* and the 1 co-infection cases were positive ($n = 34$; 100%). None of the 6 *P. vivax* blood samples were amplified with the LAMP primers (Table 2).

Sequencing of LAMP products by MinION

Figure 4 and Additional file 2: Table S1 show the results of FAST5 reads analysis generated 30 min after the start of MinION sequencing. The mean collected read numbers and read bases per run were 10,052 and 2215, respectively. Of the total collected reads, 1.57–7.75% were classified into ONT-barcodes after barcode demultiplexing; however, reads from 35.8 to 70.6% were unclassified. Based on the analysis of the mapped data, a depth of coverage above 50 \times was achieved at each base on the reference sequence. Because LAMP amplicons contain random repeating structures, a read was mapped as a reference sequence more than once by BWA-MEM.

Table 3 showed the results of mapped data obtained from the analysis workflow using LAMP amplicons with plasmids harbouring the partial sequences of *k13* with or without the C580Y mutation. Variant callings from collected reads were consistent with the plasmid sequences. The average accuracy of the mapped reads of LAMP amplicons generated from plasmids with the codon mutation C580Y was 89.0%, and at the base of the mutation C580Y (position 6/from “G” to “A”), 86.4% reads were correctly matched with “A”, but 13.6% reads were mismatched (Additional file 3: Fig. S2, Additional file 2: Table S2).

None of the 34 LAMP amplicons from clinical specimens showed the C580Y mutation based on variant calling from collected reads at 30 min. These results were consistent with the results of PCR amplicon sequencing by conventional Sanger methods (Table 2 and Additional

file 2: Table S1). In addition, the same analysis of reads collected 48 h after the start of MinION sequencing were also performed (Additional file 4: Fig. S3 and Additional file 2: Table S3). The results of variant calling were consistent with the results at 30 min.

Discussion

Although several sophisticated LAMP methods have been developed for the detection of *Plasmodium* parasites from clinical specimens [11, 12, 20–33], none have been designed to detect single nucleotide polymorphisms (SNPs) associated with ART resistance in *P. falciparum*. The LAMP method combined with MinION sequencing method was capable of simultaneously detecting both *P. falciparum* parasitaemia and the codon mutation in *k13* in up to 12 samples at a high level of sensitivity without expensive or large experimental devices. The practical advantage of this method is that it enables real-time monitoring of ART resistance in patients diagnosed with *P. falciparum* malaria, even in areas with limited resources, because the samples can be prepared and sequenced without storage or transfer to a fully-equipped laboratory. Another advantage of this protocol is that sequence data can be stored in the cloud and shared at anytime, anywhere in the world.

The gold standard diagnostic method for malaria to date is the microscopic examination of Giemsa-stained thin or thick blood smears. Several rapid diagnostic tests (RDT), such as the rapid immunochromatographic test, are already manufactured and used throughout the world. However, misdiagnosis is common in microscopic examinations or rapid immunochromatographic tests of mixed infections when a patient has low-density parasitemia or is asymptomatic [34–36]. Consequently, low-density parasitaemia or asymptomatic carriers of *P.*

(See figure on next page.)

Fig. 4 The results of FAST5 reads analysis collected 30 min from the start of MinION sequencing. **a** Histogram of FAST5 read sizes from each MinION sequencing run. **b** Collector’s curve reflecting the total base pairs of sequencing yield over time for each MinION sequencing run. **c** Depth of coverage for each ONT-barcode number and MinION sequencing run




Table 3 Summary of the MinION sequencing of LAMP amplicons from plasmids

Reference plasmid	Total reads	Mapped reads > MQ60	Coverage	% bases depth > 50	Accuracy (%)	Sequencing result
C580Y	389	793	748.10	100	89.54	C580Y
Wild type	476	822	761.83	100	89.9	Wild type

falciparum could be reservoirs for malaria infections [37, 38]. Therefore, highly sensitive diagnostic tools that can detect even low-density parasitaemia are essential for the elimination and accurate epidemiological surveillance of ART-resistant *P. falciparum*. Molecular diagnostic tests such as PCR or LAMP methods are highly sensitive and can detect plasmodium infections even at very low levels of parasitaemia [20, 39–42]. In particular, the detection limit of standard nested-PCR [18], real-time PCR [43], and commercial LAMP assay for malaria (Loopamp™ MALARIA Pan/Pf Detection Kit, Eiken Chemical, Japan) are respectively 6.0, 0.7, and 25.0 parasites/mL, and exceed the detection limit of both microscopy (50–500 parasites/mL) [44, 45] and RDT for *P. falciparum*-specific histidine-rich protein II, *Plasmodium* lactate dehydrogenase, or aldolase (100 parasites/mL) [44]. In this study, the detection limit of the LAMP method for *k13* gene was 1.0×10 copies/reaction (corresponding to 12.5 parasite/mL), and the levels of sensitivity and specificity were 100% in the clinical samples. The sensitivity of the LAMP assay is equivalent to previously described LAMP assays (sensitivity ~100% and specificity ~100%) [11, 12, 20–33].

Using the high specificity of the LAMP method, techniques for genotyping gene polymorphisms without sequence analysis have been developed, including LAMP with FIP and/or BIP primers, which are designed to contain a single-nucleotide polymorphism at each 5' end (SNPs-LAMP) (<http://loopamp.eiken.co.jp/>) or insert artificial mutations using FIP and/or BIP primers (ARMS-LAMP) [46, 47]. However, the design of LAMP primers for the above-mentioned protocol is not only quite difficult and limited, but also less sensitive than conventional LAMP assays [48]. In contrast, the LAMP method combined with MinION sequencing can achieve unlimited LAMP primer designs for genotyping polymorphisms.

A diagnostic method for dengue virus via the real-time LAMP method combined with MinION sequencer has been previously developed [16]. However, the multiplex sequence methods for that protocol use barcoded-LAMP FIP primers, which insert a unique 24-base ONT-barcode tag within the F1c sites. In this study, the detection limit of “Barcode-LAMP FIP primers” was one-hundredth that of the “Non-barcoded LAMP FIP primers” (1.0×10^3

copies/reaction vs 1.0×10 copies/reaction), and the “Barcode-LAMP FIP primers” can detect only 85.7% of *P. falciparum* parasitaemia cases (Additional file 5: Fig. S4, Additional file 2: Tables S4, S5). In the study by Yamagishi et al., the sensitivity of real-time LAMP in the clinical cases was only 80%, and the LAMP reaction took 90 min. The sensitivity of the above-mentioned LAMP assay was inferior to previously reported LAMP assays for dengue virus (sensitivity ~100%) [49–51]. In addition, the protocol took over 2 h to prepare the library using the Ligation Sequencing Kit (SQK-LSK002). In this study, all procedures from DNA extraction to variant calling were completed within 3 h.

It is well documented that MinION reads have a much higher error rate, especially via 1D-chemistry, than conventional Sanger sequencing or other next generation sequencing platforms, including 454, Illumina, and IonTorrent technology [52–55]. The rate of accuracy of 1D-read chemistry with R9.4 flow cell is about 90% (<https://nanoporetech.com/>). To generate accurate sequences, genomic regions must be read multiple times, with errors eliminated through consensus averaging. Because of the high error rate of MinION, a 50-fold read coverage of the genome positions was sufficient to accurately determine the genotypes. LAMP amplicons are constructed of complicated repeating structures; therefore, a single read generated from MinION contains several multiple base sequences of the target region. This LAMP method combined with MinION sequencing could collect sufficient sequence data within 30 min of starting MinION sequencing. The results of variant calling by this method were completely consistent with the results of Sanger sequencing. Therefore, this procedure is a highly accurate method of variant calling. However, it did not achieve sufficient accuracy for deep sequencing to detect low-frequency mutations within a sample because the average rate of accuracy was 89.4%.

Therefore, if the proportion of the C580Y mutant strain is smaller than that of the wild type strain in the same clinical specimen, the C580Y mutation may not be accurately detected. It is thought that this disparity occurred for the following reasons: the high error rate of MinION reads, the fidelity of the LAMP enzyme, and the mapping algorithm of BWA-MEM. ONT has introduced 1D2-read chemistry with higher precision than 1D-read chemistry

(<https://nanoporetech.com/>), which will improve the precision of the combined LAMP and MinION sequencer method in the future.

This study was limited due to the small number of malaria samples; hence, the performance of the sensitivity and specificity in clinical samples have not been demonstrated. The specificity in clinical samples that contained only *P. vivax* were evaluated. In addition, there were no C580Y mutations in *k13* found among 34 *P. falciparum* isolates collected in 2010 in Indonesia. An accurate evaluation of the sensitivity and specificity of this method requires further investigation in a clinical setting.

The LAMP method combined with MinION sequencing detected SNPs only in the 95 bp region incorporating the C580Y codon mutation in *k13* amplified by the LAMP primers. The other codon mutations in *k13* are reportedly associated with ART resistance. This procedure will be used to detect other mutations for the surveillance of ART resistance in the next study. Due to the flexibility afforded in the design of LAMP primers, this procedure can easily be used to detect other SNPs associated with drug-resistance genes not only in *P. falciparum*, but also in bacteria or fungi.

Conclusions

In this study, an innovative diagnostic technology to detect the codon mutation C580Y in *k13* of *P. falciparum* was demonstrated. The method, which uses the LAMP assay combined with MinION sequencer, is rapid, simple, and highly sensitive. This procedure could contribute to the epidemiological surveillance of ART-resistant *P. falciparum* and be applied to the analyses of sequence polymorphisms or genotyping. The next step for refining this procedure will be to conduct a clinical evaluation to verify that the sensitivity and specificity are sufficient and consistent in resource-limited endemic regions.

Additional files

Additional file 1: Fig. S1. Sequence alignment of the kelch-propeller domain of human *Plasmodium* parasites and the primer locations for LAMP. Pf, *P. falciparum*; Pv, *P. vivax*; Poc, *P. ovale curtisi*; Pow, *P. ovale wallikeri*; Pm, *P. malariae*; Pk, *P. knowlesi*.

Additional file 2: Table S1. The analysis of FAST5 reads collected 30 min after the start of MinION sequencing. **Table S2.** The results of the visualization analysis and the variant calling using the reads collected 30 min and 48 h after the start of MinION sequencing. **Table S3.** The analysis of FAST5 reads collected 48 h after the start of MinION sequencing. **Table S4.** Nucleotide sequences of the Barcoded-LAMP primers constructed for *kelch 13* of *Plasmodium falciparum*. **Table S5.** Summary of the results of nested PCR and Barcoded-LAMP amplifications in clinical samples.

Additional file 3: Fig. S2. A comparison of visualizations of mapped MinION reads with Sanger sequencing trace data. a: The upper image shows mapped reads from LAMP amplicons generated from plasmid DNA with

C580Y as the reference sequence (*kelch13* of wild type; KT956001.1) visualized by igvtools. The lower image is a Sanger sequencing trace of the C580Y allele using the same sample. b: The upper image shows mapped reads from plasmid DNA with the wild type as the reference sequence. The lower image is the Sanger sequencing trace of the wild type allele. Asterisks show the specific sequences located at the codon position of C580Y in *kelch13*.

Additional file 4: Fig. S3. The results of FAST5 reads analysis collected 48 h from the start of MinION sequencing. a: Histogram of FAST5 read sizes from each MinION sequencing run. b: Collector's curve reflecting the total base pairs of the sequencing yield over time for each MinION sequencing run. c: Depth of coverage for each ONT-barcode number and MinION sequencing run.

Additional file 5: Fig. S4. The detection limit of Barcode-LAMP assay for *kelch13* of *Plasmodium falciparum* with tenfold serial dilutions of plasmid DNA.

Authors' contributions

JT, YS, JY, and TMa conceived and designed the study; KI, NT, LRR, JS, and TMa gathered and analysed the data, and drafted and edited the manuscript; KH, YE, RM, HO, TMu, SM, YS, and JY supervised the study and revised the manuscript; YS, JY, and TMa performed project administration. All authors read and approved the final manuscript.

Author details

¹ Department of Infectious Disease and Infection Control, Saitama Medical University, 38 Morohongo, Moroyama-machi, Iruma-gun, Saitama 350-0495, Japan. ² Center for Clinical Infectious Diseases and Research, Saitama Medical University, 38 Morohongo, Moroyama-machi, Iruma-gun, Saitama 350-0495, Japan. ³ Department of Computational Biology and Medical Science, Graduate School of Frontier Sciences, The University of Tokyo, 5-1-5 Kashiwanoha, Kashiwa, Chiba 277-8562, Japan. ⁴ Research Center for Zoonosis Control, Hokkaido University, North 20, West 10 Kita-ku, Sapporo, Hokkaido 001-0020, Japan. ⁵ Faculty of Medicine, Oita University, 1-1 Hasama-machi, Yufu, Oita 879-5593, Japan. ⁶ Department of Medical Entomology, Faculty of Tropical Medicine, Mahidol University, 420/6 Ratchawithi Road, Thung Phaya, Ratchathewi, Bangkok 10400, Thailand. ⁷ Research Institute for Microbial Diseases, Osaka University, 3-1 Yamadaoka, Suita, Osaka 565-0871, Japan. ⁸ Division of Biomedical Sciences, Department of Basic Veterinary Medicine, Obihiro University of Agriculture and Veterinary Medicine, 2-11 Inada-cho, Obihiro, Hokkaido 080-8555, Japan. ⁹ Department of Parasitology, Faculty of Medicine, Sam Ratulangi University, Kampus Unsrat, Bahu Manado 95115, Indonesia. ¹⁰ Department of Infectious Diseases and Infection Control, Saitama Medical Center, Saitama Medical University, 1981 Kamoda, Kawagoe, Saitama 350-8550, Japan. ¹¹ Department of Microbiology, Saitama Medical University, 38 Morohongo, Moroyama-machi, Iruma-gun, Saitama 350-0495, Japan. ¹² Global Station for Zoonosis Control, GI-CoRE, Hokkaido University, North 20, West 10 Kita-ku, Sapporo, Hokkaido 001-0020, Japan.

Acknowledgements

This work was supported in part by the National BioResource Project in Japan. (<http://www.nbrp.jp/>).

Competing interests

The authors declare that they have no competing interests.

Availability of data and materials

The datasets supporting the conclusions of this article are included within the article.

Consent for publication

Not applicable.

Ethics approval and consent to participate

Ethical clearances were obtained from Institutional Ethical Review Board of Sam Ratulangi University and The University of Tokyo (Approval Number

10–49). The purpose of the study was explained in detail and written, informed consent was obtained from all participants. Confidentiality was maintained by avoiding the use of names and other identifiers.

Funding

This research did not receive any specific grant from funding agencies in the public, commercial, or not-for-profit sectors.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Received: 7 October 2017 Accepted: 22 May 2018

Published online: 29 May 2018

References

- WHO. Global Malaria Programme Artemisinin and artemisinin-based combination therapy resistance. Status report. Geneva: World Health Organization; 2016.
- Dondorp AM, Nosten F, Yi P, Das D, Phyo AP, Tarning J, et al. Artemisinin resistance in *Plasmodium falciparum* malaria. *N Engl J Med*. 2009;361:455–67.
- Hawkes M, Conroy AL, Kain KC. Spread of artemisinin resistance in malaria. *N Engl J Med*. 2014;371:1944–5.
- Noedl H, Se Y, Schaefer K, Smith BL, Socheat D, Fukuda MM. Evidence of artemisinin-resistant malaria in western Cambodia. *N Engl J Med*. 2008;359:2619–20.
- WHO. Guidelines for the treatment of malaria. 3rd ed. Geneva: World Health Organization; 2015.
- Miotto O, Amato R, Ashley EA, MaLinnis B, Almagro-Garcia J, Amaratunga C, et al. Genetic architecture of artemisinin-resistant *Plasmodium falciparum*. *Nat Genet*. 2015;47:226–34.
- Ariey F, Witkowski B, Amaratunga C, Beghain J, Langlois AC, Khim N, et al. A molecular marker of artemisinin-resistant *Plasmodium falciparum* malaria. *Nature*. 2014;505:50–5.
- Imwong M, Suwannasin K, Kunasol C, Sutawong K, Mayxay M, Rekol H, et al. The spread of artemisinin-resistant *Plasmodium falciparum* in the Greater Mekong subregion: a molecular epidemiology observational study. *Lancet Infect Dis*. 2017;17:491–7.
- Imwong M, Jindakhad T, Kunasol C, Sutawong K, Vejakkama P, Dondorp AM. An outbreak of artemisinin resistant falciparum malaria in Eastern Thailand. *Sci Rep*. 2015;5:17412.
- Mori Y, Notomi T. Loop-mediated isothermal amplification (LAMP): a rapid, accurate, and cost-effective diagnostic method for infectious diseases. *J Infect Chemother*. 2009;15:62–9.
- Morris U, Khamis M, Aydin-Schmidt B, Abass AK, Msellem MI, Nassor MH, et al. Field deployment of loop-mediated isothermal amplification for centralized mass-screening of asymptomatic malaria in Zanzibar: a pre-elimination setting. *Malar J*. 2015;14:205.
- Aydin-Schmidt B, Xu W, Gonzalez IJ, Polley SD, Bell D, Shakely D, et al. Loop mediated isothermal amplification (LAMP) accurately detects malaria DNA from filter paper blood samples of low density parasitaemias. *PLoS One*. 2014;9:e103905.
- Quick J, Ashton P, Calus S, Chatt C, Gossain S, Hawker J, et al. Rapid draft sequencing and real-time nanopore sequencing in a hospital outbreak of Salmonella. *Genome Biol*. 2015;16:114.
- Kugelman JR, Wiley MR, Mate S, Ladner JT, Beitzel B, Fakoli L, et al. Monitoring of Ebola Virus Makona evolution through establishment of advanced genomic capability in Liberia. *Emerg Infect Dis*. 2015;21:1135–43.
- Elia MY, Isra W, Mochammad H. The evaluation on molecular techniques of reverse transcription loop-mediated isothermal amplification (RT-LAMP), reverse transcription polymerase chain reaction (RT-PCR), and their diagnostic results on MinION™ Nanopore sequencer for the detection of dengue virus serotypes. *Am J Microbiol Res*. 2015;3:118–24.
- Yamagishi J, Runtuwene LR, Hayashida K, Mongan AE, Thi LAN, Thuy LN, et al. Serotyping dengue virus with isothermal amplification and a portable sequencer. *Sci Rep*. 2017;7:3510.
- Iwasaki M, Toshihiro Y, Otsuka K, Suzuki W, Nagamine K, Hase T, et al. Validation of the loop-mediated isothermal amplification method for single nucleotide polymorphism genotyping with whole blood. *Genome Lett*. 2003;2:119–26.
- Singh B, Bobogare A, Cox-Singh J, Snounou G, Abdullah MS, Rahman HA. A genus- and species-specific nested polymerase chain reaction malaria detection assay for epidemiologic studies. *Am J Trop Med Hyg*. 1999;60:687–92.
- Singh B, Kim Sung L, Matusop A, Radhakrishnan A, Shamsul SS, Cox-Singh J, et al. A large focus of naturally acquired *Plasmodium knowlesi* infections in human beings. *Lancet*. 2004;363:1017–24.
- Cook J, Aydin-Schmidt B, Gonzalez IJ, Bell D, Edlund E, Nassor MH, et al. Loop-mediated isothermal amplification (LAMP) for point-of-care detection of asymptomatic low-density malaria parasite carriers in Zanzibar. *Malar J*. 2015;14:43.
- Cuadros J, Perez-Tanoira R, Prieto-Perez L, Martin-Martin I, Berzosa P, Gonzalez V, et al. Field evaluation of malaria microscopy, rapid malaria tests and loop-mediated isothermal amplification in a rural hospital in South Western Ethiopia. *PLoS One*. 2015;10:e0142842.
- Hayashida K, Simukoko H, Simmunza M, Ndebe J, Chota A, Namangala B, et al. Direct detection of falciparum and non-falciparum malaria DNA from a drop of blood with high sensitivity by the dried-LAMP system. *Parasit Vectors*. 2017;10:26.
- Hopkins H, Gonzalez IJ, Polley SD, Angutoko P, Ategeka J, Asiimwe C, et al. Highly sensitive detection of malaria parasitemia in a malaria-endemic setting: performance of a new loop-mediated isothermal amplification kit in a remote clinic in Uganda. *J Infect Dis*. 2013;208:645–52.
- Lau YL, Fong MY, Mahmud R, Chang PY, Palaeya V, Cheong FW, et al. Specific, sensitive and rapid detection of human *Plasmodium knowlesi* infection by loop-mediated isothermal amplification (LAMP) in blood samples. *Malar J*. 2011;10:197.
- Lau YL, Lai MY, Fong MY, Jelip J, Mahmud R. Loop-mediated isothermal amplification assay for identification of five human plasmodium species in Malaysia. *Am J Trop Med Hyg*. 2016;94:336–9.
- Lee PW, Ji DD, Liu CT, Rampao HS, do Rosario VE, Lin IF, et al. Application of loop-mediated isothermal amplification for malaria diagnosis during a follow-up study in Sao Tome. *Malar J*. 2012;11:408.
- Lucchi NW, Demas A, Narayanan J, Sumari D, Kabanywany A, Kachur SP, et al. Real-time fluorescence loop mediated isothermal amplification for the diagnosis of malaria. *PLoS One*. 2010;5:e13733.
- Lucchi NW, Ljolje D, Silva-Flannery L, Udhayakumar V. Use of malachite green-loop mediated isothermal amplification for detection of *Plasmodium* spp. parasites. *PLoS One*. 2016;11:e0151437.
- Modak SS, Barber CA, Geva E, Abrams WR, Malamud D, Ongagna YS. Rapid point-of-care isothermal amplification assay for the detection of malaria without nucleic acid purification. *Infect Dis (Auckl)*. 2016;9:1–9.
- Mohon AN, Elahi R, Khan WA, Haque R, Sullivan DJ Jr, Alam MS. A new visually improved and sensitive loop mediated isothermal amplification (LAMP) for diagnosis of symptomatic falciparum malaria. *Acta Trop*. 2014;134:52–7.
- Oriero EC, Okebe J, Jacobs J, Van Geertruyden JP, Nwakanma D, D'Alessandro U. Diagnostic performance of a novel loop-mediated isothermal amplification (LAMP) assay targeting the apicoplast genome for malaria diagnosis in a field setting in sub-Saharan Africa. *Malar J*. 2015;14:396.
- Patel JC, Lucchi NW, Srivastava P, Lin JT, Sug-Aram R, Aruncharus S, et al. Field evaluation of a real-time fluorescence loop-mediated isothermal amplification assay, RealAmp, for the diagnosis of malaria in Thailand and India. *J Infect Dis*. 2014;210:1180–7.
- Sirichaisinthop J, Buates S, Watanabe R, Han ET, Suktawonjaroenpon W, Krasaesub S, et al. Evaluation of loop-mediated isothermal amplification (LAMP) for malaria diagnosis in a field setting. *Am J Trop Med Hyg*. 2011;85:594–6.
- Cheng Q, Cunningham J, Gatton ML. Systematic review of sub-microscopic *P. vivax* infections: prevalence and determining factors. *PLoS Negl Trop Dis*. 2015;9:e3413.
- Okell LC, Ghani AC, Lyons E, Drakeley CJ. Submicroscopic infection in *Plasmodium falciparum*-endemic populations: a systematic review and meta-analysis. *J Infect Dis*. 2009;200:1509–17.
- Falade CO, Ajayi IO, Nsungwa-Sabiiti J, Siribie M, Diarra A, Serme L, et al. Malaria rapid diagnostic tests and malaria microscopy for guiding malaria

- treatment of uncomplicated fevers in Nigeria and prereferred cases in 3 African countries. *Clin Infect Dis*. 2016;63:S290–7.
37. Lindblade KA, Steinhart L, Samuels A, Kachur SP, Slutsker L. The silent threat: asymptomatic parasitemia and malaria transmission. *Expert Rev Anti Infect Ther*. 2013;11:623–39.
 38. Bousema T, Drakeley C. Epidemiology and infectivity of *Plasmodium falciparum* and *Plasmodium vivax* gametocytes in relation to malaria control and elimination. *Clin Microbiol Rev*. 2011;24:377–410.
 39. Baltzell KA, Shakely D, Hsiang M, Kemere J, Ali AS, Bjorkman A, et al. Prevalence of PCR detectable malaria infection among febrile patients with a negative *Plasmodium falciparum* specific rapid diagnostic test in Zanzibar. *Am J Trop Med Hyg*. 2013;88:289–91.
 40. Kattenberg JH, Tahita CM, Versteeg IA, Tinto H, Traore Coulibaly M, D'Alessandro U, et al. Evaluation of antigen detection tests, microscopy, and polymerase chain reaction for diagnosis of malaria in peripheral blood in asymptomatic pregnant women in Nanoro, Burkina Faso. *Am J Trop Med Hyg*. 2012;87:251–6.
 41. Kyabayinze DJ, Zongo I, Cunningham J, Gatton M, Angutoko P, Ategeka J, et al. HRP2 and pLDH-based rapid diagnostic tests, expert microscopy, and PCR for detection of malaria infection during pregnancy and at delivery in areas of varied transmission: a prospective cohort study in Burkina Faso and Uganda. *PLoS One*. 2016;11:e0156954.
 42. Wu L, van den Hoogen LL, Slater H, Walker PG, Ghani AC, Drakeley CJ, et al. Comparison of diagnostics for the detection of asymptomatic *Plasmodium falciparum* infections to inform control and elimination strategies. *Nature*. 2015;528:S86–93.
 43. Perandin F, Manca N, Calderaro A, Piccolo G, Galati L, Ricci L, et al. Development of a real-time PCR assay for detection of *Plasmodium falciparum*, *Plasmodium vivax*, and *Plasmodium ovale* for routine clinical diagnosis. *J Clin Microbiol*. 2004;42:1214–9.
 44. Wongsrichanalai C, Barcus MJ, Muth S, Sutamihardja A, Wernsdorfer WH. A review of malaria diagnostic tools: microscopy and rapid diagnostic test (RDT). *Am J Trop Med Hyg*. 2007;77:119–27.
 45. Milne LM, Kyi MS, Chiodini PL, Warhurst DC. Accuracy of routine laboratory diagnosis of malaria in the United Kingdom. *J Clin Pathol*. 1994;47:740–2.
 46. Tamura S, Maeda T, Misawa K, Osa M, Hamamoto T, Yuki A, et al. Development of a highly resolved loop-mediated isothermal amplification method to detect the N526K ftsI mutation of beta-lactamase-negative ampicillin-resistant *Haemophilus influenzae*. *J Microbiol Methods*. 2017;141:108–14.
 47. Lee D, Kim EJ, Kilgore PE, Takahashi H, Ohnishi M, Tomono J, et al. A novel loop-mediated isothermal amplification assay for serogroup identification of *Neisseria meningitidis* in cerebrospinal fluid. *Front Microbiol*. 2015;6:1548.
 48. Torigoe H, Seki M, Yamashita Y, Sugaya A, Maeno M. Detection of *Haemophilus influenzae* by loop-mediated isothermal amplification (LAMP) of the outer membrane protein P6 gene. *Jpn J Infect Dis*. 2007;60:55–8.
 49. Parida M, Horioko K, Ishida H, Dash PK, Saxena P, Jana AM, et al. Rapid detection and differentiation of dengue virus serotypes by a real-time reverse transcription-loop-mediated isothermal amplification assay. *J Clin Microbiol*. 2005;43:2895–903.
 50. Hu SF, Li M, Zhong LL, Lu SM, Liu ZX, Pu JY, et al. Development of reverse-transcription loop-mediated isothermal amplification assay for rapid detection and differentiation of dengue virus serotypes 1–4. *BMC Microbiol*. 2015;15:265.
 51. Teoh BT, Sam SS, Tan KK, Johari J, Danlami MB, Hooi PS, et al. Detection of dengue viruses using reverse transcription-loop-mediated isothermal amplification. *BMC Infect Dis*. 2013;13:387.
 52. Schirmer M, Ijaz UZ, D'Amore R, Hall N, Sloan WT, Quince C. Insight into biases and sequencing errors for amplicon sequencing with the Illumina MiSeq platform. *Nucleic Acids Res*. 2015;43:e37.
 53. Shao W, Boltz VF, Spindler JE, Kearney MF, Maldarelli F, Mellors JW, et al. Analysis of 454 sequencing error rate, error sources, and artifact recombination for detection of low-frequency drug resistance mutations in HIV-1 DNA. *Retrovirology*. 2013;10:18.
 54. Liu L, Li Y, Li S, Hu N, He Y, Pong R, et al. Comparison of next-generation sequencing systems. *J Biomed Biotechnol*. 2012;2012:251364.
 55. Jain M, Tyson JR, Loose M, Ip CLC, Eccles DA, O'Grady J, et al. MinION analysis and reference consortium: phase 2 data release and analysis of R9.0 chemistry. *F1000Res*. 2017;6:760.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

