

# SCIENTIFIC REPORTS



OPEN

## Engineering of Primary Human B cells with CRISPR/Cas9 Targeted Nuclease

Matthew J. Johnson<sup>1,2,3</sup>, Kanut Laoharawee<sup>1,2,3</sup>, Walker S. Lahr<sup>1,2,3</sup>, Beau R. Webber<sup>1,2,3</sup> & Branden S. Moriarity<sup>1,2,3</sup>

B cells offer unique opportunities for gene therapy because of their ability to secrete large amounts of protein in the form of antibody and persist for the life of the organism as plasma cells. Here, we report optimized CRISPR/Cas9 based genome engineering of primary human B cells. Our procedure involves enrichment of CD19<sup>+</sup> B cells from PBMCs followed by activation, expansion, and electroporation of CRISPR/Cas9 reagents. We are able to expand total B cells in culture 10-fold and outgrow the IgD<sup>+</sup> IgM<sup>+</sup> CD27<sup>-</sup> naïve subset from 35% to over 80% of the culture. B cells are receptive to nucleic acid delivery via electroporation 3 days after stimulation, peaking at Day 7 post stimulation. We tested chemically modified sgRNAs and Alt-R gRNAs targeting *CD19* with Cas9 mRNA or Cas9 protein. Using this system, we achieved genetic and protein knockout of *CD19* at rates over 70%. Finally, we tested sgRNAs targeting the *AAVS1* safe harbor site using Cas9 protein in combination with *AAV6* to deliver donor template encoding a splice acceptor-*EGFP* cassette, which yielded site-specific integration frequencies up to 25%. The development of methods for genetically engineered B cells opens the door to a myriad of applications in basic research, antibody production, and cellular therapeutics.

B cells and their downstream effectors, plasma blasts and plasma cells, are central to the humoral arm of the adaptive immune system and are the only cell lineage that secretes antibodies<sup>1,2</sup>. These antigen-specific antibodies help to protect the host from infection via neutralization and opsonization of pathogens and toxins<sup>3,4</sup>. B cells are also considered professional antigen presenting cells (APCs) with the ability to present exogenous antigen to naïve T cells on MHC class II<sup>5</sup>. Additionally, B cells play a critical role in the development and maintenance of immunological memory through the generation of memory B cells capable of rapidly reinitiating an antigen-specific immune response upon reencountering their cognate antigen, as well as long-lived plasma cells, which passively maintain low levels of antigen-specific antibodies in the plasma and mucosal surfaces<sup>6</sup>.

Several unique features make B cells and plasma cells an attractive target for genome engineering. B cells are easily isolated in large number from the peripheral blood and can be activated, grown, and expanded in culture<sup>7-9</sup>. Plasma cells specifically upregulate multiple pathways to increase their ability to produce massive amounts of protein and receive pro survival signals from bone marrow stromal cells to extend their longevity<sup>10</sup>. Thus, a system where B cells are isolated from the peripheral blood, engineered to express a specific gene, matured to plasma cells, and then reintroduced back to the host would be an excellent gene therapy platform for long term, high titer expression of protein in the serum.

Previous work on genome engineering B cells has focused on editing primary human peripheral B cells, or editing hematopoietic stem cells and then maturing them to B cells *in vitro*, to express a specific monoclonal antibody (mAb) using lentiviral vectors or transposons<sup>11,12</sup>. However, the non-site specific integration of lentiviral vectors and transposons and their predilection for integrating into transcriptionally active regions carries the inherent risk of harmful mutations which greatly limits any future clinical applications of such approaches<sup>13</sup>.

The site-specific nature of the CRISPR/Cas9 nuclease system makes it ideal for precision engineering of B cells and recent advancements in CRISPR/Cas9 technologies have drastically improved editing efficiency in primary human cells<sup>14</sup>. The use of high quality Cas9 mRNA or protein and chemically modified guide RNAs (gRNAs) were demonstrated to increase the editing efficiency to over 75% in primary cell types<sup>14</sup>. Additional work has

<sup>1</sup>Department of Pediatrics, University of Minnesota, Minneapolis, MN, 55455, USA. <sup>2</sup>Center for Genome Engineering, University of Minnesota, Minneapolis, MN, 55455, USA. <sup>3</sup>Masonic Cancer Center, University of Minnesota, Minneapolis, MN, 55455, USA. Correspondence and requests for materials should be addressed to B.S.M. (email: [mori0164@umn.edu](mailto:mori0164@umn.edu))

demonstrated that Adeno Associated Virus Serotype 6 (AAV6) functions as a highly efficient DNA template donor for homologous recombination (HR) in both hematopoietic stem cells and T cells<sup>15,16</sup>. Therefore, we have developed methods for using the CRISPR/Cas9 system with chemically modified sgRNA to knockout endogenous genes and combined this with AAV6 to engineer primary human B cells.

Here, we report the site-specific engineering of primary human B cells using the CRISPR/Cas9 system. The use of chemically modified sgRNA is compared with the two-component AltR system<sup>17,18</sup> and the use of purified Cas9 protein is compared with chemical modified Cas9 mRNA. Finally, we demonstrate the use of Cas9 protein, AAVS1-targeting sgRNA, and AAV6 as a DNA donor for HR to introduce a DNA sequence encoding splice acceptor-EGFP into the AAVS1 locus of B cells.

## Martials and Methods

**Peripheral Blood Mononuclear Cells (PBMCs) isolation.** Human PBMCs from de-identified, normal, healthy donors were obtained by automated leukapheresis (Memorial Blood Centers) and further depleted of red blood cells by lysis with ACK buffer (Thermo Fisher Scientific) for 3 minutes at room temperature. PBMCs were then cultured in RPMI 1640 (Thermo Fisher Scientific) supplemented with 2mM L-glutamine (Invitrogen), 1% streptomycin and penicillin (Invitrogen), and 10% Fetal Bovine Serum (Gibco) at a density of  $1 \times 10^6$  cells/ml. Signed informed consent was obtained from all donors and the study was approved by the University of Minnesota Institutional Review Board (IRB study number 1602E84302). All methods were performed in accordance with relevant the guidelines and regulations.

**Isolation and expansion of B cells.** B cells were isolated from PBMCs by immunomagnetic negative selection using EasySep Human B Cell Isolation Kit (Stemcell Technologies) in accordance with the manufacturer's instructions. B cells were cultured in StemMACS HCS Expansion Media XF (Miltenyi Biotec) supplemented with 1% streptomycin and penicillin (Invitrogen), 5% Human AB Serum (Valley Biomedical), and 125 IU/mL of IL4 (Miltenyi Biotec) at a density of  $5 \times 10^5$  cells/ml. B cells were expanded by crosslinking CD40 using Human CD40-Ligand Multimer Kit (Miltenyi Biotec) at a concentration of 8 U/ml in accordance with the manufacturer's instructions. Media, IL4, and multimeric CD40L were refreshed every 3–4 days throughout the duration of all experiments.


**Electroporations.** For electroporations,  $3 \times 10^5$  B cells were added to a combination of 1  $\mu$ g of sgRNA, 1.5  $\mu$ g of Cas9 mRNA (TriLink), 1  $\mu$ g of EGFP mRNA (TriLink), or 1  $\mu$ g of Cas9 protein (Integrated DNA Technologies) and brought up in 10  $\mu$ l of T Buffer (Thermo Fisher Scientific) for electroporation in the Neon Transfection System (Thermo Fisher Scientific). Cells were then loaded in 10  $\mu$ l tips and electroporated in accordance with the manufacturer's instructions using settings of 1400 volts, 10 ms width, and 3 pulses unless indicated otherwise.

**Editing with CRISPR/Cas9.** gRNAs targeting *CD19* and *AAVS1* genes were designed using the Crispr MIT website (<http://crispr.mit.edu/>) along with Cas-OFFinder which together determined the most efficient and specific guides. crRNAs were ordered from Integrated DNA Technologies for testing of the Alt-R system. Tracr and crRNA were combined according to the manufacturer's protocol and these crRNA: tracrRNA duplexes were then mixed with Cas9 protein before being electroporated into primary human B cells<sup>17,18</sup>. For experiments using chemically modified sgRNAs, guides were ordered from Synthego Biosciences with 2'-O-methyl and 3' phosphorothioate internucleotide modified linkages at the 3' and 5' ends of sgRNAs. sgRNAs were complexed with Cas9 protein for 20 minutes similar to the two component Alt-R system. Guide RNA sequences for *CD19* and *AAVS1* were CUGUGCUGCAGUGCCUCAA and GUCACCAAUCCUGUCCUAG, respectively. Analysis of Cas9 activity was performed by amplifying a region of genomic DNA around the site of interest. The primer pair AAATTCAGGAAAGGTTGGAAG and GCGGACCTCTTCTGTCCATG was used to amplify a region flanking the *CD19* cut site and primers GAGATGGCTCCAGGAAATGG and ACCTCTCACTCCTTTCATTGG were used to amplify the *AAVS1* region. These PCR amplicons were then submitted for Sanger sequencing and subsequent chromatographs were uploaded to the TIDE website (<https://tide-calculator.nki.nl>) with their respective non-edited control chromatograms<sup>19</sup>.

For experiments inducing a t(8;14) translocation, sgRNAs targeting known *MYC* and *IGH* Burkitt Lymphoma translocation regions were designed<sup>20</sup>. The sgRNA sequence AUAAAGCAGGAAUGUCCGAC was used for *MYC* and AUAUUCCACCCAGGUAGUGG was used for *IGH*. The primer pair AACCAGGTAAGCACCGAATCC and CTGGCTCACACAGGCGATATGC was used to amplify a region flanking the *MYC* cut site and the primer pair AGCATCTCAATACCCTCCTCTTGG and CAGCCCCAGTTCAGCCTTGTTTAG was used to amplify the *IGH* region and the TIDE website was used to quantify indel formation. The CTGGCTCACACAGGCGATATGC *MYC* reverse primer and AGCATCTCAATACCCTCCTCTTGG *IGH* forward primer which flank the translocation site were used to confirm the t(8;14) translocation had occurred.

The GGACGAGCTGTACAAGTAACG EGFP primer and GAGACAGTGACCAACCATCC AAVS1 reverse primer was used to amplify a region flanking the 3' end of EGFP and the AAVS1 cutsite to confirm site specific integration of the donor template had occurred.

**B cell subset analysis.** PBMCs or sorted B cells were washed with PBS and incubated with Fixable Viability Dye eFluor780 (eBioscience) for 10 mins at room temperature. Cells were then washed with a buffer containing 0.5% BSA and stained with fluorescently labeled antibodies to CD10 (Biolegend), CD19 (BD Horizon), CD20 (Biolegend), CD21 (BD OptiBuild), CD27 (Biolegend), CD38 (BD Horizon), CD138 (Biolegend), CXCR5 (BD Horizon), IgD (BD Pharmingen), IgG (Biolegend), and IgM (Biolegend) (Supplementary Table 1) for 20 mins at room temperature. Samples were then washed again with a buffer containing 0.5% BSA, fixed with 1% PFA, and analyzed on a LSRFortessa flow cytometer (BD Biosciences). Data analysis was performed using FlowJo version 9.9.3 (FlowJo LLC).



**Figure 1.** B cells are efficiently expanded with CD40L and IL4. Total cell numbers (A) and viability (B) of isolated B cells in culture and during expansion (n = 3). (C) Percentage of Naïve B cells (IgD+ IgM+ CD27−, blue line), Pro or Pre B cells (IgD+ IgM− CD27−, red line), Natural Effector B cells (IgD+ IgM+ CD27+, green line), Memory B cells (IgD− CD27+ CXCR5+, orange line), and Activated Memory B cells (IgD− CD27+ CXCR5−, purple line) of total live CD19+ B cells in culture over time. (n = 3).

**rAAV6 construction and transduction.** Splice acceptor-EGFP rAAV virus with 500 base pair homology arms targeting the AAVS1 site was produced using a custom ordered gBlock Gene Fragment (Integrated DNA Technologies), the AAV backbone, and Gibson assembly. After the vector was constructed it was sent to Vigene Biosciences for AAV6 packaging and concentration. Primary human B cell cultures were transfected with AAVS1 targeting sgRNA pre-complexed to Cas9 protein immediately before being transduced with AAV6 at the indicated MOIs.

**Statistical Analysis.** Data analyses were performed using FlowJo version 9.9.3. Statistical analyses were performed using one-way ANOVA tests with GraphPad Prism version 7.0c software (GraphPad Software). Error bars represent standard error in all cases. Measurements were considered statistically different at  $p < 0.05$ .


**Data Availability.** All data generated and analyzed in this study are included in this published article and supplementary files.

## Results

**Naïve B cells are activated and expanded in culture with CD40L and IL4.** To determine the optimal timing and culture conditions for transfections, sorted CD19+ B cells were placed into culture and activated by CD40L crosslinking and IL4 for 18 days<sup>21</sup>. B cell numbers remain roughly constant for 7 days before beginning to rapidly expand (Fig. 1A). Similarly, the viability of B cells reduced to below 60% initially on days 1 ( $57 \pm 8\%$  (mean  $\pm$  SD)) and 3 ( $56 \pm 6\%$ ), before rebounding on day 7 ( $90 \pm 2\%$ ) and remaining over 85% for the duration of the experiment (Fig. 1B).


Flow cytometry was used to determine the frequency of commonly described subsets of B cells present in the peripheral blood, including IgD+ IgM+ CD27− naïve B cells, IgD+ IgM− CD27− pro and pre B cells, IgD+ IgM+ CD27+ natural effector B cells, and IgD− CD27+ CXCR5+ memory and IgD− CD27+ CXCR5− activated memory B cells (Supplementary Figure 1)<sup>22</sup>. Plasmablasts and Plasma cells were undetectable at all time points (data not shown). The naïve B cell population represented  $35 \pm 6\%$  of the total purified CD19+ B cell culture, but expanded to  $81 \pm 2\%$  by day 7 and remained over 80% throughout the duration of the experiment (Fig. 1C).

**Multiple electroporation conditions allow for efficient transfection of activated B cells.** To determine the ideal electroporation condition for CD19+ B cells in culture, we measured the efficiency of EGFP mRNA transfection across 13 conditions of varying voltage, pulse-width, and pulse number (Fig. 2A). While all conditions tested achieved over 50% transfection efficiency as measured by GFP expression and over 80% post electroporation cell viability, Condition 12 resulted in the optimal combination of EGFP+ frequency, cell viability, and post-electroporation cell expansion (Fig. 2B). Therefore, Condition 12, a combination of 1400 volts, 10 ms pulse-width, and 3 pulses, was selected for all subsequent experiments.




**Figure 2.** Electroporation conditions affect B cell transfection efficiency, viability, and post transduction expansion. (A) Table of electroporation conditions tested. (B) The percentage of B cells transfected with chemically modified *EGFP* coding mRNA (*upper panel*) as well as post transfection cell viability (*middle panel*) and expansion (*lower panel*) ( $n = 2$ ). Dashed lines indicate pre-electroporation viability and cell count.


**B cells become permissive to transfection with mRNA by day 3.** To determine the optimal time point for transfection following activation, CD19+ B cells were transfected with mRNA encoding EGFP at multiple time points between day 0 and day 10 post-activation, and the frequency of EGFP+ cells was assayed 48hrs after transfection by flow cytometry (Fig. 3). No EGFP expression was detected in any B cell subset following transfection on day 0 and day 1 post-activation, indicating a lack of permissiveness to mRNA transfection at these time points. However, mRNA transfection on days 3, 7, and 10 post-activation resulted in EGFP expression in all



**Figure 3.** B cells become permissive to electroporation in culture. Frequency of EGFP<sup>+</sup> cells 2 days after electroporation with mRNA encoding EGFP (left panel) at the indicated time points (n = 6).



**Figure 4.** B cells are efficiently edited by CRISPR/Cas9. (A) B cells were enriched and isolated, electroporated with either chemically modified sgRNA or Alt-R gRNA in combination with either Cas9 protein or chemically modified mRNA encoding Cas9 protein, and then assessed for indel formation (n = 6). (B) Frequency of indel formation as measured by TIDE analysis (left panel) and loss of CD19 expression as measured by flow cytometry (right panel), viability (left panel) and total cell numbers (right panel) (\*indicate p < 0.05) (n = 6).



**Figure 5.** Simultaneous application of RNPs targeting *MYC* and *IGH* creates a t(8;14) Burkitt Lymphoma translocation. (A) Diagram of relative positions of *MYC* on Chromosome 8 (left diagram), the *IGH* locus on Chromosome 14 (center diagram) and the t(8;14) translocation (right diagram) associated with Burkitt Lymphoma. Arrows indicate primer direction and triangles indicate gRNA target sites. (B) Indel pattern and efficiency of RNPs targeting *MYC* (top panel) and *IGH* (bottom panel). (C) PCR amplification using *MYC* (left panel), *IGH* locus (middle panel), or a combination of *MYC/IGH* locus primers (right panel) of samples engineered with RNPs targeting *MYC*, the *IGH* locus, or a combination of *MYC* and *IGH* locus, as well as no gRNA (represented by '-') and H<sub>2</sub>O controls.

subsets. The naïve B cell subset, which at these time points represents the majority of the culture, was transfected at  $41 \pm 8\%$  on day 3, and over 50% on days 7 ( $58 \pm 9\%$ ) and 10 ( $52 \pm 7\%$ ). These frequencies are below those observed in Fig. 2B, possibly due to reduced cell activation due to the disruption caused by frequent resampling throughout the time-course experiment. All subsequent experiments in this study were performed on Day 7.


**Assessment of optimal gRNA and Cas9 format for engineering B cells.** Both chemically modified sgRNAs or two component trRNA/crRNA CRISPR/Cas9 systems have been successfully used in primary cells<sup>14,18,23</sup>. Therefore, to examine efficacy between two component and sgRNA systems in B cells we targeted *CD19*, which is expressed on all peripheral blood B cells, with the two component Alt-R system or chemically modified sgRNA in combination with either Cas9 protein or chemically modified Cas9 mRNA (Fig. 4A). Chemically modified sgRNA in combination with either Cas9 protein or Cas9 mRNA was able to induce over 70% indel formation (Cas9 protein  $74 \pm 19\%$  and Cas9 mRNA  $72 \pm 4\%$ ), resulting in  $>70\%$  loss of cell surface CD19 protein (Cas9 protein  $75 \pm 11\%$  and Cas9 mRNA  $72 \pm 13\%$ ) in expanded B cells (Fig. 4B, upper panels). Comparably high editing efficiency was observed with Alt-R gRNA when pre-complexed with Cas9 protein (indel formation  $65 \pm 9\%$  and CD19 protein knockout  $72 \pm 11\%$ ), however no editing was observed following co-delivery of Alt-R gRNA and Cas9 mRNA.

While a high editing frequency was demonstrated in both conditions using chemically modified sgRNA, there were significantly different impacts on the post editing expansion of the cells. Cells treated with chemical modified sgRNA and Cas9 mRNA had significantly lower cell numbers ( $130,000 \pm 111,000$  cells and  $59 \pm 18\%$  viability) compared to cultures treated with chemically modified sgRNA and Cas9 protein ( $529,000 \pm 150,000$  cells and  $75 \pm 4\%$  viability) (Fig. 4B lower panels). In addition, we tested two sgRNAs targeting the *BCL2* gene to further validate our findings beyond *CD19*. When used as part of an RNP these sgRNAs achieved  $91 \pm 13\%$  and  $69 \pm 19\%$  indel formation (Supplementary Figure 2). Therefore, the combination of chemically modified sgRNA and Cas9 protein (henceforth referred to as RNP) was chosen for use in all subsequent experiments.

**CRISPR/Cas9 RNPs can be used to induce chromosomal translocations.** One potential application of this technology is the induction of chromosomal translocations in a site-specific manner for the study of B cell malignancies. As a proof of concept we modeled the t(8;14) translocation characteristic of most instances of Burkitt Lymphoma (Fig. 5A)<sup>24,25</sup>. In this disease, the gene encoding the transcription factor c-Myc (*MYC*) is translocated to the immunoglobulin heavy chain (*IGH*) locus where enhancers and other elements are thought to strongly up regulate its expression<sup>26–28</sup>. RNPs targeting *MYC* and *IGH* were used individually and in combination to create DSBs in B cells. Both *MYC* and *IGH* RNPs created high frequency DSB individual at 78.0% and 56.6% respectively (Fig. 5B). As expected, PCR amplification products using *MYC* and *IGH* locus primer pairs were detectable in all conditions (Fig. 5C, left and center panels). However, the hallmark t(8;14) translocation, amplified using *IGH* locus forward primer and the *MYC* reverse primer, was only detected in samples engineered with a combination of RNPs targeting *MYC* and the *IGH* locus (Fig. 5C, right panel). The observed band matches the size of the predicted 295 bp amplification product.

**AAV6 effectively delivers DNA template for HR in B cells.** Adeno Associated Virus serotype 6 (AAV6) has been successfully used as a homologous recombination DNA template donor in both CD34+ hematopoietic stem cells and T cells<sup>15,29</sup>. Given their developmental relationship to B cells as members of the lymphohematopoietic hierarchy, we speculated that AAV6 might also function as a DNA template donor in B cells. To test the efficacy of AAV6 as a DNA template donor, a splice acceptor-based system was designed and constructed that would express EGFP only if integration into the *AAVS1* 'safe harbor' site was achieved<sup>30,31</sup>, but not if the vector remained





**Figure 6.** B cells can be engineered using homologous recombination with rAAV6 donor and CRISPR/Cas9. (A) A schematic design of the DNA donor template for a splice acceptor-EGFP system targeting the *AAVS1* site by HR used in this study. (B) Percentage of EGFP+ cells represent integration frequency following transfection with rAAV6 containing a splice acceptor-EGFP system at the indicated MOIs ( $n = 3$ ). (C) PCR amplification using a combination of EGFP and *AAVS1* primers of samples engineered with RNPs targeting *AAVS1* and exposed to rAAV6 containing a splice acceptor-EGFP system at the indicated MOIs, as well as no gRNA (represented by ‘-’) and H<sub>2</sub>O controls.

episomal (Fig. 6A). Activated primary human CD19+ B cells were electroporated with an RNP targeting the *AAVS1* locus immediately prior to exposure to AAV6 vector at the indicated MOIs (Fig. 6B). Site specific donor DNA integration as evidenced by EGFP expression was observed at all MOIs tested, peaking at  $18 \pm 4\%$  with an MOI of 250,000. Control samples, transfected with Cas9 protein, but not sgRNA, showed no EGFP expression above background indicating no episomal expression of the EGFP transcript. Integration was further confirmed using junction PCR amplification of the region spanning the 3' end of EGFP insert and *AAVS1* site (Fig. 6C). The observed bands match the size of the predicted 1.5 kb amplification product. Absence of signal in the both no AAV6 and the no gRNA controls confirm the specificity of this primer set.

## Discussion

B cells are an attractive candidate for use in gene therapy due to their ease of isolation from peripheral blood, rapid growth, and their ability to be converted to long lived plasma cells which produced high levels of protein. The ability to genetically engineer B cells opens a new avenue for using B cells to treat various diseases and to study their basic biology in a human context. Here, we report that RNP, chemically modified sgRNA co-transfected with Cas9 mRNA, and the two-component Alt-R system all efficiently engineer B cells. However, delivery of RNP results in significantly higher post-editing expansion. Furthermore, we show AAV6-based donor templates can be successfully used in conjunction with the CRISPR/Cas9 system for targeted integration of a gene of interest at the AAVS1 'safe harbor' site with integration efficiency up to 25%.

While this manuscript was in preparation a group from the University of Washington published a related manuscript on the precision engineering of B cells with CRISPR/Cas9 in *Molecular Therapy*<sup>32</sup>. Many of our overall findings—that B cells can be efficiently edited with CRISPR/Cas9 ribonucleoprotein and that AAV6 is an efficient DNA template donor—are in agreement with their data. In addition, we present data on the subsets of B cells in expanded/activated CD19+ B cell populations, clearly demonstrating that the naïve subset outgrows other subsets in culture and that all subsets are efficiently electroporated after Day 3. We also closely examined alternative methods of CRISPR/Cas9 engineering, demonstrating that the use of Cas9 mRNA, instead of Cas9 protein, has a negative impact on post-engineering B cell health and expansion. This observation may be due to activation of cytosolic RNA sensors such as TLR7, TLR8 or RIG1, which are present in immune cells, including B cells, and have been shown to impact activation, growth regulation, and apoptosis<sup>33,34</sup>. In addition, we demonstrate that CRISPR/Cas9 reagents can be used to model the translocation characteristic of Burkitt Lymphoma. In contrast, their manuscript further describes the *in vitro* maturation of B cells into plasma cells and demonstrates that B-cell activation factor (BAFF) helps to successfully engraft plasma cells into immunodeficient mice. Other, smaller differences in the exact methodologies used, such as differing culture media and B cell isolation kits, appear to have limited impact. Together, these papers collectively reinforce the viability of using CRISPR/Cas9 engineered B cells as a tool for basic B cell research and as a novel platform for gene therapy.

This technology could be adapted for use in research into the basic biology of B cells. Researchers can knock-out or overexpress genes for use in both *in vitro* and *in vivo* humanized mouse models. For example, it has been reported that expression of BCL2, BCMA, and MCL1 are essential for the long-term survival of plasma cells *in vivo* in mice<sup>35,36</sup>. However, their role in the long-term survival of human plasma cells *in vivo* remains difficult to study. CRISPR/Cas9 editing of primary human B cells will make it possible to knockout or over express these gene *in vitro* before maturing cells to plasma cells and adoptively transferring them into a humanized mouse model to assess the gene of interest's role in *in vivo* plasma cell longevity. Furthermore, we demonstrate the potential for creating chromosomal translocations present in Burkitt Lymphoma. However, an initial characterization of these cell cultures with this translocation showed no obvious transformed properties, i.e. increased proliferation or cytokine independent growth. Moving forward, we will study the process of transformation in greater depth, including determining the impact of additional secondary mutations known to frequently occur in Burkitt Lymphoma<sup>28</sup>.

This technology could also be adapted for the treatment of enzymopathies, a type of genetic disorder that results in missing or defective enzymes which are usually treated by costly periodic enzyme replacement therapies. Autologous B cells could be isolated, engineered to express the deficient enzyme, and matured into plasma cells before being adoptively transferred back to the original donor where they can express the enzyme. The high protein production potential and longevity of plasma cells could potentially allow for large amounts of the enzyme to collect in the serum for decades after the adoptive transfer. Additional transplants could be applied if serum enzyme levels begin to wane over time.

Treatment using certain monoclonal antibody (mAb) therapies, particularly those used to treat autoimmunity, work by antagonistically binding cytokine or cytokine receptors thereby inhibiting inflammation<sup>37–39</sup>. The biological half-life of these mAbs generally ranges from one to three weeks so biweekly or monthly injections of the mAb are required to keep titers at an efficacious level. Autologous B cells could be engineered to express therapeutic mAbs, before being matured into plasma cells and adoptively transferred back to the host where they would secrete large amounts of mAb into the serum. This could potentially eliminate the need for the repeated injections currently needed in mAb therapies. This could also be further adapted as an alternative method to vaccination for diseases where broadly neutralizing antibodies are protective, but very difficult to elicit using traditional methods of vaccination. Further work is required to fully realize the potential of the engineered B cell/plasma cell lineage for therapeutic applications.

## References

1. Radbruch, A. *et al.* Competence and competition: the challenge of becoming a long-lived plasma cell. *Nat. Rev. Immunol.* **6**, 741–750 (2006).
2. Slifka, M. K., Antia, R., Whitmire, J. K. & Ahmed, R. Humoral immunity due to long-lived plasma cells. *Immunity* **8**, 363–372 (1998).
3. Hessel, A. J. *et al.* Fc receptor but not complement binding is important in antibody protection against HIV. *Nature* **449**, 101–104 (2007).
4. Forthall, D. N. Functions of Antibodies. *Microbiol. Spectr.* **2**, 1–17 (2014).
5. Rodríguez-Pinto, D. B cells as antigen presenting cells. *Cell. Immunol.* **238**, 67–75 (2005).
6. Amanna, I. J. & Slifka, M. K. Mechanisms that determine plasma cell lifespan and the duration of humoral immunity. *Immunol. Rev.* **236**, 125–138 (2010).
7. Spriggs, M. K. *et al.* Recombinant human CD40 ligand stimulates B cell proliferation and immunoglobulin E secretion. *J. Exp. Med.* **176**, 1543–1550 (1992).
8. Liebig, T. M., Fiedler, A., Zoghi, S., Shimabukuro-Vornhagen, A. & von Bergwelt-Baildon, M. S. Generation of Human CD40-activated B cells. *J. Vis. Exp. JoVE*, <https://doi.org/10.3791/1373> (2009).
9. Naito, M. *et al.* CD40L-Tri, a novel formulation of recombinant human CD40L that effectively activates B cells. *Cancer Immunol. Immunother.* **62**, 347–357 (2013).
10. Cassese, G. *et al.* Plasma cell survival is mediated by synergistic effects of cytokines and adhesion-dependent signals. *J. Immunol. Baltim. Md 1950* **171**, 1684–1690 (2003).



11. Fusil, F. *et al.* A Lentiviral Vector Allowing Physiologically Regulated Membrane-anchored and Secreted Antibody Expression Depending on B-cell Maturation Status. *Mol. Ther.* **23**, 1734–1747 (2015).
12. Luo, X. M. *et al.* Engineering human hematopoietic stem/progenitor cells to produce a broadly neutralizing anti-HIV antibody after *in vitro* maturation to human B lymphocytes. *Blood* **113**, 1422–1431 (2009).
13. Staunstrup, N. H. *et al.* Hybrid Lentivirus-transposon Vectors With a Random Integration Profile in Human Cells. *Mol. Ther. J. Am. Soc. Gene Ther.* **17**, 1205–1214 (2009).
14. Hendel, A. *et al.* Chemically modified guide RNAs enhance CRISPR-Cas genome editing in human primary cells. *Nat. Biotechnol.* **33**, 985–989 (2015).
15. Wang, J. *et al.* Highly efficient homology-driven genome editing in human T cells by combining zinc-finger nuclease mRNA and AAV6 donor delivery. *Nucleic Acids Res.* **44**, e30–e30 (2016).
16. Wang, J. *et al.* Homology-driven genome editing in hematopoietic stem and progenitor cells using zinc finger nuclease mRNA and AAV6 donors. *Nat. Biotechnol.* **33**, 1256–1263 (2015).
17. Al Abdallah, Q., Ge, W. & Fortwendel, J. R. A Simple and Universal System for Gene Manipulation in *Aspergillus fumigatus*: *In Vitro*-Assembled Cas9-Guide RNA Ribonucleoproteins Coupled with Microhomology Repair Templates. *mSphere* **2** (2017).
18. Mikheikin, A. *et al.* DNA nanomapping using CRISPR-Cas9 as a programmable nanoparticle. *Nat. Commun.* **8**, 1665 (2017).
19. Brinkman, E. K., Chen, T., Amendola, M. & van Steensel, B. Easy quantitative assessment of genome editing by sequence trace decomposition. *Nucleic Acids Res.* **42**, e168 (2014).
20. Blum, K. A., Lozanski, G. & Byrd, J. C. Adult Burkitt leukemia and lymphoma. *Blood* **104**, 3009–3020 (2004).
21. Rush, J. S. & Hodgkin, P. D. B cells activated via CD40 and IL-4 undergo a division burst but require continued stimulation to maintain division, survival and differentiation. *Eur. J. Immunol.* **31**, 1150–1159 (2001).
22. Allman, D. & Pillai, S. Peripheral B cell subsets. *Curr. Opin. Immunol.* **20**, 149–157 (2008).
23. Eyquem, J. *et al.* Targeting a CAR to the TRAC locus with CRISPR/Cas9 enhances tumour rejection. *Nature* **543**, 113 (2017).
24. Berger, R. & Bernheim, A. Cytogenetics of Burkitt's lymphoma-leukaemia: a review. *IARC Sci. Publ.* 65–80 (1985).
25. Boerma, E. G., Siebert, R., Kluin, P. M. & Baudis, M. Translocations involving 8q24 in Burkitt lymphoma and other malignant lymphomas: a historical review of cytogenetics in the light of today's knowledge. *Leukemia* **23**, 225 (2009).
26. Adams, J. M. *et al.* The c-myc oncogene driven by immunoglobulin enhancers induces lymphoid malignancy in transgenic mice. *Nature* **318**, 533–538 (1985).
27. Sander, S. *et al.* Synergy between PI3K signaling and MYC in Burkitt lymphomagenesis. *Cancer Cell* **22**, 167–179 (2012).
28. Love, C. *et al.* The genetic landscape of mutations in Burkitt lymphoma. *Nat. Genet.* **44**, 1321 (2012).
29. Ling, C. *et al.* High-Efficiency Transduction of Primary Human Hematopoietic Stem/Progenitor Cells by AAV6 Vectors: Strategies for Overcoming Donor-Variation and Implications in Genome Editing. *Sci. Rep.* **6**, 35495 (2016).
30. DeKelver, R. C. *et al.* Functional genomics, proteomics, and regulatory DNA analysis in isogenic settings using zinc finger nuclease-driven transgenesis into a safe harbor locus in the human genome. *Genome Res.* **20**, 1133–1142 (2010).
31. Hockemeyer, D. *et al.* Highly efficient gene targeting of expressed and silent genes in human ESCs and iPSCs using zinc finger nucleases. *Nat. Biotechnol.* **27**, 851–857 (2009).
32. Hung, K. L. *et al.* Engineering Protein-Secreting Plasma Cells by Homology-Directed Repair in Primary Human B Cells. *Mol. Ther.* **26**, 456–467 (2017).
33. Masumi, A. *et al.* Enhanced RIG-I expression is mediated by interferon regulatory factor-2 in peripheral blood B cells from hepatitis C virus-infected patients. *Biochem. Biophys. Res. Commun.* **391**, 1623–1628 (2010).
34. Jensen, S. & Thomsen, A. R. Sensing of RNA Viruses: a Review of Innate Immune Receptors Involved in Recognizing RNA Virus Invasion. *J. Virol.* **86**, 2900–2910 (2012).
35. Peperzak, V. *et al.* Mcl-1 is essential for the survival of plasma cells. *Nat. Immunol.* **14**, 290–297 (2013).
36. O'Connor, B. P. *et al.* BCMA Is Essential for the Survival of Long-lived Bone Marrow Plasma Cells. *J. Exp. Med.* **199**, 91–98 (2004).
37. Rau, R. Adalimumab (a fully human anti-tumour necrosis factor  $\alpha$  monoclonal antibody) in the treatment of active rheumatoid arthritis: the initial results of five trials. *Ann. Rheum. Dis.* **61**, ii70–ii73 (2002).
38. Brekke, O. H. & Sandlie, I. Therapeutic antibodies for human diseases at the dawn of the twenty-first century. *Nat. Rev. Drug Discov.* **2**, 52 (2003).
39. Welch, B. Adalimumab (Humira) for the Treatment of Rheumatoid Arthritis. *Am. Fam. Physician* **78**, 1406 (2008).

## Acknowledgements

The authors are grateful to Mitchell Kluesner for his careful review of the manuscript. This study was supported by NIH/NIAID 5 R21 AI128087-02.

## Author Contributions

M.J.J., K.L., W.S.L., B.R.W., and B.S.M. designed research; M.J.J., K.L., and W.S.L. performed research; M.J.J., K.L., W.S.L., and B.S.M. analyzed data; and M.J.J. and B.S.M. wrote the paper.

## Additional Information

**Supplementary information** accompanies this paper at <https://doi.org/10.1038/s41598-018-30358-0>.

**Competing Interests:** The authors declare they have submitted a patent application for the methods and use of engineered B cells based on the work published in this manuscript.

**Publisher's note:** Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.



**Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit <http://creativecommons.org/licenses/by/4.0/>.

© The Author(s) 2018