PET of Adoptively Transferred Chimeric Antigen Receptor T Cells with ⁸⁹Zr-Oxine

Michael R. Weist^{1,2}, Renate Starr³, Brenda Aguilar³, Junie Chea¹, Joshua K. Miles¹, Erasmus Poku¹, Ethan Gerdts³, Xin Yang³, Saul J. Priceman³, Stephen J. Forman³, David Colcher¹, Christine E. Brown³, and John E. Shively¹

¹Department of Molecular Immunology, Beckman Research Institute, City of Hope Medical Center, Duarte, California; ²Irell and Manella Graduate School of Biological Sciences, City of Hope Medical Center, Duarte, California; and ³Department of Hematology, City of Hope Medical Center, Duarte, California

Chimeric antigen receptor (CAR) T cell therapy is a promising clinical approach for reducing tumor progression and prolonging patient survival. However, improvements in both the safety and the potency of CAR T cell therapy demand quantitative imaging techniques to determine the distribution of cells after adoptive transfer. The purpose of this study was to optimize 89Zr-oxine labeling of CAR T cells and evaluate PET as a platform for imaging adoptively transferred CAR T cells. Methods: CAR T cells were labeled with 0-1.4 MBq of 89Zroxine per 10⁶ cells and assessed for radioactivity retention, viability, and functionality. In vivo trafficking of 89Zr-oxine-labeled CAR T cells was evaluated in 2 murine xenograft tumor models: glioblastoma brain tumors with intracranially delivered IL13Rq2-targeted CAR T cells, and subcutaneous prostate tumors with intravenously delivered prostate stem cell antigen (PSCA)-targeted CAR T cells. Results: CAR T cells were efficiently labeled (75%) and retained more than 60% of the 89Zr over 6 d. In vitro cytokine production, migration, and tumor cytotoxicity, as well as in vivo antitumor activity, were not significantly reduced when labeled with 70 kBg/10⁶ cells. IL13Rα2-CAR T cells delivered intraventricularly were detectable by PET for at least 6 d throughout the central nervous system and within intracranial tumors. When intravenously administered, PSCA-CAR T cells also showed tumor tropism, with a 9-fold greater tumor-to-muscle ratio than for CARnegative T cells. Conclusion: 89Zr-oxine can be used for labeling and imaging CAR T cells while maintaining cell viability and function. On the basis of these studies, we conclude that 89Zr-oxine is a clinically translatable platform for real-time assessment of cell therapies.

Key Words: ⁸⁹Zr-oxine; positron emission tomography; chimeric antigen receptor T cells; prostate cancer; glioblastoma

J Nucl Med 2018; 59:1531–1537 DOI: 10.2967/jnumed.117.206714

The genetic modification of T cells to express chimeric antigen receptors (CARs), which enable these cells to specifically target and kill cancer cells, has become an important immunotherapeutic approach in the treatment of cancer (*1*–*3*). CAR T cells combine the target specificity of antibodies with the intracellular signaling components of the T cell receptor and costimulatory molecules to redirect T cell–mediated immune responses (*4*). Significant progress

Received Dec. 5, 2017; revision accepted Apr. 17, 2018. For correspondence or reprints contact: John E. Shively, City of Hope Medical Center, 1500 E. Duarte Rd., Duarte, CA 91010.

E-mail: jshively@coh.org Published online May 4, 2018.

COPYRIGHT © 2018 by the Society of Nuclear Medicine and Molecular Imaging.

in CAR T cell therapy has been made over the past several decades, as led by the clinical success of CD19-CAR T cells against refractory B-cell malignancies, in which complete responses are achieved in most patients (5–8). Despite ongoing progress, challenges remain in improving CAR T cell therapy, including poor tumor targeting (9), normal-tissue toxicity (10), and lack of CAR T cell persistence (11). Real-time imaging to determine the dynamics of CAR T cell trafficking in both patients and animals can help address these issues. Because CAR T cells are manufactured ex vivo, they can readily be labeled for imaging before infusion. PET offers significant advantages over other imaging modalities because it is quantitative, highly sensitive, and subject to minimal tissue attenuation. Optimal nuclear imaging is achieved when the radionuclide half-life approximates the biologic action of the imaging target. Of the commercially available PET radionuclides, 89Zr has a relatively long half-life (78.4 h) suitable for imaging CAR T cells over the days required for trafficking and antitumor responses.

⁸⁹Zr-oxine has previously been used for labeling and imaging murine lymphocytes (12–15) but not human CAR T cells. ⁸⁹Zr-oxine has preparation and functional characteristics similar to those of its SPECT counterpart, ¹¹¹In-oxine, which has been used for many types of cells (16). ¹¹¹In-oxine has been used to visualize the trafficking and tumor accumulation of γ-δ-T cells and HER2-primed T cells (17,18). In comparison to SPECT, PET has higher resolution (19), making ⁸⁹Zr-oxine an attractive choice for imaging cells. Another PET-based approach for labeling cells is ¹⁸F-FDG, which is inadequate for tracking CAR T cells because of its short half-life (109.7 min) and relatively poor radioactivity retention (20). On the basis of the existing clinical applications of ⁸⁹Zr (21) and its oxine complex, we selected ⁸⁹Zr-oxine for labeling CAR T cells and imaging with PET.

We present an alternative method for labeling cells with ⁸⁹Zr-oxine in which ⁸⁹Zr-oxine is prepared from ⁸⁹Zr-oxalate and retained within CAR T cells. Before cells were imaged in vivo, we identified an optimal labeling activity that maintained cell viability and function while maximizing the radioactive signal for imaging. To further assess the robustness of this PET imaging platform, we tracked ⁸⁹Zr-oxine–labeled CAR T cells administered via 3 different delivery routes: intravenous, intratumoral, and intraventricular. These data demonstrate the tumor tropism of adoptively transferred CAR T cells and their distribution based on administration routes.

MATERIALS AND METHODS

Labeling Cells with 89Zr-Oxine

⁸⁹Zr-oxalate from 3D Imaging was aliquoted into borosilicate vials. On the basis of the specific activity of aqueous ¹¹¹In-oxine preparations (22), a 50-μL quantity of 6.9 mM 8-hydroxyquinoline (oxine) in

0.1 M 4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid (HEPES) was added to each 37 MBq of ⁸⁹Zr-oxalate, pH was adjusted to 7.4 with 1 M HEPES, and the solution was incubated at room temperature for 15 min. After 2 washings with Hanks buffered salt solution (HBSS) (without calcium, magnesium, or phenol red), CAR T cells (20 × 10⁶/mL) were incubated with 0–1,400 kBq of unpurified ⁸⁹Zr-oxine solution per 10⁶ cells for 30 min at 37°C. Past ⁸⁹Zr- and ¹¹¹Inoxine reports indicate that purification of the oxine complex is not necessary (*14*,*15*,*22*) and that oxine is rapidly lost from labeled cells after washing (*23*). The radiolabeled cells were washed twice with T cell medium (IL15 [0.5 ng/mL] and IL2 [50 U/mL]) before in vitro experiments or with HBSS plus 1% human serum albumin before in vivo experiments.

Labeling Efficiency of 89Zr-Oxine

The yield of the ^{89}Zr -oxine synthesis was determined by instant thin-layer chromatography (92% \pm 0.1%, n=35) and reversed-phase high-performance liquid chromatography (82%) (Supplemental Fig. 1; supplemental materials are available at http://jnm.snmjournals.org). ^{89}Zr -oxine remained at the origin, whereas unincorporated ^{89}Zr migrated to the solvent front developed with 0.9% NaCl. A gradient of 0%–100% acetonitrile with 0.1% trifluoroacetic acid on a SunFire C18 column (100 Å, 3.5 μm , 4.6 \times 150 mm; Waters) separated ^{89}Zr -oxine complexes from unincorporated ^{89}Zr . Given this high yield, the ^{89}Zr -oxine solution was immediately used for cell labeling.

Cell Lines

PBT030-2 is a low-passage, patient-derived IL13R α 2-positive tumor line derived from glioblastoma tissue that was obtained in accordance with City of Hope Institutional Review Board-approved protocols and cultured as described previously (24). To generate the PBT030-2 ffLuc-positive line for visualization with luminescence, the cells were modified by lentiviral transduction to express enhanced green fluorescent protein and firefly luciferase (ffLuc) reporter genes (24). The adherent glioblastoma cell line, U251T, was cultured as described previously (25). Prostate cancer PC3 cells (CRL-1435; ATCC) were genetically modified to express prostate stem cell antigen (PC3-PSCA) by lentiviral transduction and enriched for antigen expression by flow cytometry.

IL13Rα2-CAR T cells were generated from T cells enriched for central memory as described by Brown et al. (1). This second-generation IL13Rα2-CAR consists of a membrane-bound interleukin 13 (IL13) mutein with a mutation at the E13 site, an extracellular IgG4-Fc spacer containing L235E and N297Q point mutations (26), a CD4 transmembrane domain, and intracellular signaling domains of 4-1BB and CD3ζ. The lentiviral cassette also includes a truncated CD19 separated from the CAR via a T2A ribosomal skip sequence (27). The PSCA-CAR is identical to the IL13Rα2-CAR except that the extracellular tumor-targeting domain comprises a single-chain variable fragment derived from the humanized anti-PSCA antibody clone A11 (kindly provided by Drs. Anna Wu and Robert Reiter, UCLA) (28). For CAR T cell engineering, the cells were stimulated with CD3/ CD28 Dynabeads (Invitrogen) at a ratio of 3:1 (T cell to bead) and genetically modified via lentiviral transduction. T cells were maintained in culture at $0.6-1 \times 10^6$ cells/mL in X-VIVO 15 medium (Bio Whittaker) plus 10% fetal calf serum (Hyclone) and supplemented with T cell medium 3 times per week. Both T cell populations were enriched for CAR expression by selecting for the truncated CD19 component coexpressed with CAR using an anti-CD19 immunomagnetic isolation kit (Stem Cell Technologies).

Cell-Labeling Efficiency and Radioactivity Retention

Cell-labeling efficiency was determined as the percentage radioactivity in the final cell pellet over the total radioactivity added. For retention experiments, CAR T cells were cultured in T cell medium

and sampled at the indicated time points. Radioactivity within the cell pellet and supernatant was assayed in a γ -counter (Packard Cobra Auto-Gamma Counter) to determine retention (% retention = cell pellet radioactivity/total radioactivity \times 100%). For subcellular fractionation experiments, 5×10^6 cells were processed per the manufacturer's instructions using a subcellular fractionation kit (catalog no. 78840; ThermoFisher Scientific). Each fraction was subsequently counted for radioactivity.

Assays for CAR T Cell Function

After labeling and incubation under standard culture conditions, the cytotoxicity, cytokine production, and migration potential of labeled IL13R α 2-CAR T cells were evaluated. For all experiments, cell viability was assessed with 0.1% propidium iodide on a MoxiFlow cassette-based flow cytometer (Orflo Technologies) before labeling, after labeling, and at the indicated time points in culture. For the nonradioactive ^{nat}Zr-oxine preparation, 10 μ L of a 17 μ M ^{nat}ZrCl₄ solution in 1 M oxalic acid was compounded with 260 μ L of 6.8 mM oxine in 0.1 M HEPES and adjusted to pH 7.4 with 1 M HEPES. Equivalent volumes of ^{nat}Zr-oxine were added at specific activities toxic to CAR T cells (1.5 and 12 MBq/10⁶ cells).

The cytotoxicities of $^{89}\mathrm{Zr}\text{-}\mathrm{oxine}\text{--}\mathrm{labeled}$ CAR T cells were also monitored using biophotonic luciferase assays as described by Brown et al. (25). A 5×10^4 PBT030-2 quantity of ffLuc-positive cells was coincubated with $5\text{--}50\times10^4$ $^{89}\mathrm{Zr}\text{-}\mathrm{oxine}\text{--}\mathrm{labeled}$ IL13R α 2-CAR T cells in 0.2 mL at the specified time points for 7 hours at 5% CO₂ and $37^{\circ}\mathrm{C}$. The plates were then supplemented with luciferin to a final concentration of 0.14 mg/mL, incubated for 10 min in the dark, and measured for luminescence (1-s exposure) with the Victor3 plate reader (PerkinElmer). Wells with the addition of 1% sodium dodecyl sulphate were used as a reference for maximum killing, and wells without CAR T cells were used as a reference for minimum killing. For all in vitro comparisons, 2-way ANOVA was performed to determine the significance of differences between labeling conditions.

 $^{89}\text{Zr-oxine-labeled IL13R}\alpha\text{2-CAR T}$ cells were stimulated for 2 h by plate-bound IL13R\$\alpha\text{2-Fc}\$ (614-IR-100; R&D Systems [0–1,000 ng/well]), and the secreted interferon (IFN)-\$\gamma\$ was quantified by an enzyme-linked immunosorbent assay (catalog no. 430108; BioLegend) per the manufacturer's instructions. Labeled and unlabeled cells were also stimulated with antigen-positive tumor cells in an in vitro stimulation assay. A 2.5 × 10^5 quantity of T cells was incubated overnight with 5 × 10^4 PBT030-2 ffLuc-positive cells in 0.5 mL of T cell medium in 24-well plates. After a 16-h incubation, the supernatants were removed, and a cytokine 10-plex human panel (catalog no. LHC0001M; ThermoFisher) was used to quantify human IFN-\$\gamma\$ within the supernatant, per the manufacturer's instructions.

The migration potential of 89 Zr-oxine–labeled IL13R α 2-CAR T cells was assessed in 96-well ChemoTx plates with a 5- μ m pore diameter (catalog no. 106-5; NeuroProbe). In the lower chambers, 30 μ L of unconditioned Dulbecco modified Eagle medium with 2% fetal calf serum or U251T-conditioned Dulbecco modified Eagle medium with 2% fetal calf serum were added (prepared as described previously (29)). The upper chambers were loaded with 7.5×10^4 CAR T cells/30 μ L (2.5 \times 106 cells/mL). After 4 h at 37°C, residual cells were scraped off the polycarbonate filter, and the plate was centrifuged for 2 min at 400g. The filter was removed, and cells in the lower chamber were counted with the MoxiFlow. Percentage migration was calculated as the number of cells in the lower chamber divided by the total number of cells plated per well.

Tumor Targeting of CAR T Cells

All mouse studies were approved by the City of Hope Institutional Animal Care and Use Committee. Under anesthesia from an intraperitoneal injection of ketamine (132 mg/kg) and xylazine (8.8 mg/kg), 6- to 8-wk-old male NOD. Cg-Prkdc^{scid}Il2rg^{tm1Wjl}/SzJ (NSG)

mice were stereotactically injected with PBT030-2 ffLuc-positive glioma cells into the left forebrain (30). A 5-mm skin incision was made along the sagittal suture, and a burr hole was drilled into the skull 2 mm lateral and 0.5 mm anterior to the bregma. A 1×10^5 quantity of cells in 2 µL of Dulbecco modified Eagle medium/F-12 (50:50) plus B-27 and heparin supplement was injected as follows: 1 μL each at depths of 2.5 and 2.25 mm from the dura. Eight days after tumor implantation, the mice were treated intratumorally at the tumor coordinates (i.e., intracranially) or in the right lateral ventricle (intraventricularly) with 89Zr-oxine-labeled (70 kBq/106 cells) or unlabeled IL13Rα2-CAR T cells. Mice for therapy and PET studies were treated with 0.2×10^6 and 2.0×10^6 cells, respectively. Control mice with tumors were treated intratumorally with unlabeled mock-transduced T cells. For intratumoral injections, T cells in 3 µL of phosphatebuffered saline were injected as follows: 0.5 µL each at a depth of 3.0, 2.75, 2.50, 2.35, 2.25, and 2.15 mm from the dura. For intraventricular injections, a 5-mm skin incision was made along the sagittal suture, a burr hole drilled into the skull 0.9 mm lateral and 0.3 mm caudal to the bregma, and T cells in 3 µL of phosphate-buffered saline were injected at a depth of 2.5 mm from the dura. Tumor growth was monitored and quantified by luminescence imaging on an AmiX optical imaging system and AMIView software (Spectral Instruments), and survival was assessed according to the guidelines of the American Veterinary Medical Association (30). Tukey multiple-comparison tests and log-rank tests were used to compare differences in flux




FIGURE 1. ⁸⁹Zr-oxine does not affect cell viability or function at optimized cell labeling activity. After ⁸⁹Zr-oxine labeling at specified specific activities, IL13Rα2-CAR T cells were cultured under standard conditions for up to 144 h. (A) Viability of untransduced (mock) and IL13Rα2-CAR T cells was monitored over 144 h of cell culture. (B) IL13Rα2-CAR T cell killing of PBT030-2 ffLuc-positive cells was assessed in 7-h luciferase cytotoxicity assay at 10:1 T cell-to-tumor ratio (n = 6) at specified time points after ⁸⁹Zr-oxine labeling and subsequent cell culturing. (C) Enzyme-linked immunosorbent assay for IFN-γ cytokine production of ⁸⁹Zr-oxine-labeled IL13Rα2-CAR T cells after incubation with plate-bound IL13Rα2 for 4 h (n = 2). (D) IFN-γ production determined by cytokine bead analysis of ⁸⁹Zr-oxine-labeled IL13Rα2-CAR T cells after overnight incubation with PBT030-2 ffLuc-positive cells. (E) Migration of IL13Rα2-CAR T cells to U251T conditioned medium or control medium in 4-h Boyden chamber chemotaxis assay (n = 4). All data are expressed as mean ± SE.

and survival, respectively. In the subcutaneous prostate PC3-PSCA tumor model, the mice were injected with 2.5×10^6 PC3-PSCA cells in HBSS into the right flank as described previously (28). After 28 d, the mice were injected via the tail vein with 7.5×10^6 T cells (70 kBq/ 10^6 cells) or 0.37 MBq of cell-free 89 Zr-oxine.

PET Imaging and 89Zr-Biodistribution Studies

PET scans were acquired with an Inveon small-animal PET/CT system (Siemens Medical Solutions) at an energy window of 350–605 keV with a time resolution of 6 ns. In the first acquisition session for the glioblastoma model, the mice remained under anesthesia from the ketamine and xylazine received before CAR T cell implantation, further supplemented by 2%–4% isoflurane in oxygen. At subsequent time points, 2%–4% isoflurane in oxygen alone was used for anesthesia. All PET scans were 10–20 min long. The data were sorted into 2-dimensional sinograms using the Fourier rebinning method and corrected for intrascan radiodecay, detector nonuniformity, and random coincidence noise. Images were reconstructed using iterative 3-dimensional ordered-subsets expectation maximization (4 iterations). CT scans were also acquired at the last imaging time point.

For biodistribution studies with ⁸⁹Zr-oxine, major organs were removed (heart, lung, liver, spleen, stomach, large intestine, small intestine, kidneys, right quadriceps, right femur, left stifle joint, and lumbar vertebrae) after euthanasia by cardiocentesis. Brains were washed in phosphate-buffered saline and sectioned with a razor blade to separate

the right and left cortices, olfactory bulbs, and cerebellum. Tissues were weighed and measured for radioactivity using an automated γ-counter (Wallac Wizard 3; PerkinElmer) along with 3 dose standards. After correcting for background and radioactive decay, the percentage injected dose (%ID) and %ID/g of each tissue were calculated. Region-of-interest tracing of subcutaneous tumors was performed in the IRW software suite (Siemens). Vivo-Quant (Invicro) was used to analyze PET/CT images of the intracranial model, and the 3D Brain Atlas tool of VivoQuant was used to quantify radioactivity in the right and left cortices of the brain.

Histology

After the terminal PET imaging time point, brains or subcutaneous tumors were collected and fixed in 4% paraformaldehyde. Immunohistochemical staining for human CD3 T cells was performed on 10-µm-thick sections from paraffin-embedded tissue, which were deparaffinized in xylene followed by graded ethanol. The samples were then quenched in 3% hydrogen peroxide and incubated at high pH. The slides were washed in buffer, incubated for 30 min with anti-human CD3 (Dako M7254 and clone F7.2.38, 1:200) in dilution buffer, rewashed in buffer, and incubated with secondary antibody. After further washes in buffer, the slides were incubated with the chromogen diaminobenzidine tetrahydrochloride, rinsed in distilled H2O, counterstained with hematoxylin, rinsed in distilled H2O, rehydrated through grades of ethanol and xylene, and mounted. Adjacent slides were also stained by a standard hematoxylin and eosin

RESULTS

Optimization of 89Zr-Oxine Labeling for CAR T Cells

To ensure that in vivo signal was representative of cell localization, we assessed radionuclide retention and subcellular localization of 89 Zr after labeling. The efficiency of 89 Zr-oxine labeling of T cells was $75\% \pm 17\%$ (n=33). By 144 h (6 d), over 60% of the radioactivity was cell-associated (Supplemental Fig. 2A). Subcellular fractionation of the cells indicated that radioactivity was mostly cytoplasmic, which decreased over time (Supplemental Fig. 2B). These results demonstrate that 89 Zr-oxine is retained within CAR T cells for an extended period.

We next determined the maximum ^{89}Zr -oxine labeling specific activity that maintains CAR T cell viability and function. Using IL13R α 2-CAR T cells as a model platform, we labeled cells at specific activities of 0–1,400 kBq/10⁶ cells and monitored for viability and cell function over 144 h (6 d) (Fig. 1). Cells labeled at 70 kBq/10⁶ cells had the smallest decrease in viability among the activities tested (6.6% reduction) (Fig. 1A). Because labeled CAR T cells must also kill antigen-expressing tumor cells, the cytotoxicities of ^{89}Zr -oxine–labeled T cells were assessed with in vitro killing assays up to 145 h after labeling. When T cells were cocultured with PBT030-2 ffLuc-positive cells, CAR T cells labeled with 70 kBq/10⁶ cells were not significantly affected by

the 89 Zr-oxine labeling (P = 0.24) (Fig. 1B and Supplemental Fig. 3A). The chemical burden of an equivalent preparation of natZroxine to 1,500 kBq/106 cells did not significantly affect killing (Supplemental Fig. 3B). To test the effects of 89Zr-oxine labeling on antigen-dependent cytokine production, IL13Rα2-CAR T cells were activated with plate-bound IL13Rα2 and evaluated for IFN-γ production. Immediately after labeling, IFN-y production of cells labeled with 280 kBq/106 cells was not significantly reduced in comparison to unlabeled cells (P = 0.06) (Fig. 1C). Because cells were imaged over multiple days after radiolabeling, we also evaluated cytokine production over time in an in vitro stimulation assay. After 89Zr-oxine labeling and culturing, samples of IL13Rα2-CAR T cells were incubated with PBT030-2 ffLuc-positive cells overnight, and the supernatants were analyzed for cytokine production (Fig. 1D). Production of IFN-γ by ⁸⁹Zr-oxine–labeled CAR T cells was unaffected at 70 kBq/10⁶ cells. Finally, we evaluated the chemotaxis of ⁸⁹Zr-oxine-labeled IL13Rα2-CAR T cells to U251T-conditioned medium. 89Zr-oxine labeling at 70 kBq/106 cells did not significantly reduce the migration of CAR T cells (Fig. 1E). We observed both time- and dose-dependent effects of ⁸⁹Zr-oxine on cell viability and function. Taken together, we identified a maximum labeling activity of 70 kBq/10⁶ cells for maintaining CAR T cell function over 6 d in vitro, and this activity was used for subsequent in vivo studies. 89Zr-Oxine-Labeled Cells Retain




FIGURE 2. ⁸⁹Zr-oxine–labeled IL13Rα2-CAR T cells are therapeutic in orthotopic glioma model. (A) Schema of bioluminescence and survival analysis of NSG mice engrafted with PBT030-2 glioblastoma multiforme tumors treated (Tx) on day 8 with ⁸⁹Zr-oxine–labeled (70 kBq/10⁶ cells) or unlabeled IL13Rα2-CAR T cells, or unlabeled mock T cells injected intracranially (I.C.T.) or intraventricularly (I.C.V.). (B) Mean luminescence and SE of quantified tumor growth for each treatment group over 64 d (P = 0.006 for CAR T cells [⁸⁹Zr-labeled and unlabeled] vs. tumor only, n = 6). (C) Kaplan–Meier survival analysis of treatment groups over 120 d (P = 0.0007, n = 6). (D) Representative bioluminescent images of each treatment group.

⁸⁹Zr-Oxine–Labeled Cells Retain Potent Antitumor Activity

For 89Zr-oxine-labeled CART cells to be relevant for clinical use, the antitumor activity of the labeled cells should be comparable to that of unlabeled CAR T cells. Thus, we compared the therapeutic potential of 89Zr-oxine-labeled to unlabeled IL13Rα2-CAR T cells in an orthotopic glioblastoma tumor model with IL13Rα2-positive PBT030-2 ffLuc-positive cells in NSG mice (Fig. 2A). We observed an equivalent reduction in tumor size for mice treated intratumorally with labeled or unlabeled IL13Rα2-CAR T cells, as compared with mock-transduced T cells or tumor only (P = 0.006 each) over 64 d (Figs. 2B and)2D). A similar survival benefit was also observed (Fig. 2C); mice treated with either labeled or unlabeled CAR T cells survived to the end of the study (day 120), whereas none of the mock or untreated mice survived past 83 d (P = 0.0007 each).

PET Imaging of ⁸⁹Zr-Oxine–Labeled IL13Rα2-CAR T Cells After Intracranial Adoptive Transfer

Our team is currently evaluating the locoregional delivery of IL13Rα2-CAR T cells in early-phase clinical trials for the treatment of glioblastoma (1,31); therefore, we evaluated the utility of ⁸⁹Zr-oxine for tracking intracranially injected CAR T cells. Cohorts of mice were imaged with PET/CT and used for ⁸⁹Zr-biodistribution studies after CAR T cell administration at the tumor site or into the cerebrospinal fluid via the

ventricular system. Independent of the presence of tumor, intraventricularly injected CAR T cells were initially localized to the cranium, including the olfactory bulbs at 0.5 h, and then became distributed throughout the spinal column over 6 d (Supplemental Fig. 4A, top and middle rows). PET imaging also revealed that intraventricularly delivered CAR T cells trafficked out of the central nervous system to the liver and spleen by day 6. By contrast, intratumorally injected CAR T cells remained at their injection site throughout the 6 d of PET imaging (Supplemental Fig. 4A, bottom row). Biodistribution studies also confirmed that radioactivity went to the spine, spleen, blood, and liver of mice injected intraventricularly, whereas most of the radioactivity remained at the injection site for intratumorally treated mice (Supplemental Figs. 4B and 4C; Supplemental Table 1). Tumor-bearing mice treated intraventricularly with ⁸⁹Zr-oxine–labeled IL13Rα2-CAR T cells showed a leftfrontal-cortex focus corresponding to the tumor site (Fig. 3, top row) that was not observed in mice without tumors (Fig. 3, middle row). The same anatomic area was seen in mice treated intratumorally (Fig. 3, bottom row). Because radioactivity was observed in the olfactory bulbs of mice treated intraventricularly, we affirmed that the focus corresponded to the tumor in the cortex by fitting the




FIGURE 3. Tumor trafficking of 89 Zr-oxine–labeled IL13Rα2-CAR T cells in brain. (A) Representative PET/CT scans of NSG mice with or without PBT030-2 ffLuc-positive tumors in left cortex at 140 h after treatment with 89 Zr-oxine–labeled IL13Rα2-CAR T cells (70 kBq/10 6 cells). In mice receiving intraventricular injections (I.C.V.), labeled CAR T cells were injected into cerebrospinal fluid via right ventricle. Intratumoral injections were administered intracranially (I.C.T.) into left cortex. Arrows represent tumor location. (B) Immunohistochemistry for human CD3 (brown) showing trafficking of I.C.V.-delivered IL13Rα2-CAR T cells to contralateral PBT030-2 ffLuc-positive tumors. Scale bar is 250 μm.

VivoQuant 3D Brain Atlas to coregistered PET/CT scans (Supplemental Fig. 5). Tumor infiltration of human CD3 T cells, as assessed by immunohistochemistry, also confirmed the tumor trafficking of intraventricularly delivered unlabeled and $^{89}\text{Zr-oxine-labeled CAR}$ T cells (Fig. 3B). By injecting the hemisphere contralateral to the tumor and measuring the radioactivity in the tumor-bearing hemisphere, we could monitor CAR T cell migration. We observed a 2-fold increase (1.01 vs. 0.51) in the ratio of left-hemisphere–to–right-hemisphere radioactivity in tumor-bearing mice, compared with tumor-free mice (Supplemental Fig. 4C; Supplemental Table 1). Image quantification for the left and right cortices, as determined by the 3D Brain Atlas, yielded similar ratios. These data indicate that intraventricularly injected IL13R α 2-CAR T cells traffic to contralateral tumors and become distributed throughout the central nervous system and major organ systems.

PET Imaging of ⁸⁹Zr-Oxine–Labeled PSCA-CAR T Cells After Intravenous Injection

To further assess 89Zr-oxine as an imaging platform for CAR T cells, we evaluated the intravenous delivery of PSCA-CAR T cells in a subcutaneous prostate tumor model of PC3-PSCA cells (Fig. 4). ⁸⁹Zr-oxine-labeled PSCA-CAR T cells, ⁸⁹Zr-oxine-labeled untransduced T cells (mock), or cell-free 89Zr-oxine was injected into mice bearing subcutaneous PC3-PSCA xenografts and imaged over 7 d. Radioactivity within the first 4 h was observed in the lungs of mice treated with T cells, whereas the radioactivity of the cell-free 89Zroxine was rapidly distributed throughout the blood (Fig. 4A). Over 6 d, the signal in mice treated with ⁸⁹Zr-oxine-labeled T cells shifted from the lungs to the liver and spleen whereas the signal in cell-free control mice rapidly shifted to the liver and bone (Fig. 4A; Supplemental Table 2). Radioactivity in subcutaneous tumors was detected in mice treated with CAR but not with mock T cells (Fig. 4C; Supplemental Table 2). PET scans revealed a steady accumulation of signal in subcutaneous tumors over 7 d (Fig. 4D). Tumor radioactivity was much higher in mice treated with CAR T cells (2.44 %ID/g) than in mice treated with mock cells (0.31 %ID/g), indicating enhanced CAR T cell trafficking to the tumor. The tumorto-blood ratios of mice treated with 89Zr-oxine-labeled CAR T cells was, on average, 20.0, in comparison to 3.54 with 89Zr-oxinelabeled mock cells. Radioactivity was also present in tumors from the cell-free 89Zr-oxine control but was due to residual blood activity (Fig. 4A, bottom row); the average tumor-to-muscle ratio was greater for the CAR T cells (14.7) than for the cell-free control mice (2.59). CD3-positive cells were observed in the tumors of mice treated with 89Zr-oxine-labeled and unlabeled PSCA-CAR T cells but not in the mock cell-treated tumors (Fig. 4B). These data show that ⁸⁹Zr-oxine can visualize the tumor targeting and whole-body distributions of intravenously injected CAR T cells.

DISCUSSION

One of the primary goals of this study was to provide proof-of-principle data for a clinical imaging platform using $^{89}\mathrm{Zr}\text{-}\mathrm{oxine}$ to track CAR T cells after adoptive transfer in patients. Our team has a longstanding interest in the locoregional delivery of CAR T cells for the treatment of brain tumors and other solid tumors (1,24,27,28,31). In our ongoing clinical trial (NCT02208362), we are evaluating the safety and antitumor effects of up to 200×10^6 IL13R α 2-CAR T cells for the treatment of recurrent glioblastoma via intratumoral and intraventricular routes of administration. We recently reported a case study involving a patient with recurrent




FIGURE 4. Visualizing in vivo tumor targeting of ⁸⁹Zr-oxine–labeled PSCA-CAR T cells. (A) Representative PET (left) and PET/CT (right) maximum-intensity projections at 0.5 and 162 h, respectively. NSG mice with PC3-PSCA tumors in right flank (arrow) were intravenously treated with ⁸⁹Zr-oxine PSCA-CAR T cells (top row), ⁸⁹Zr-oxine–labeled mock-transfected T cells (middle row), or ⁸⁹Zr-oxine alone (bottom row). (B) Antihuman CD3 (brown) immunostaining for T cells in tumors of labeled or unlabeled mock or PSCA-CAR T cells. (C) Axial section showing PET signal of ⁸⁹Zr-oxine PSCA-CAR T cells at subcutaneous tumor (arrow) of representative mouse at 162 h. (D) Quantification of PET images in subcutaneous tumors of mice treated with ⁸⁹Zr-oxine PSCA-CAR T cells over 7 d.

multifocal glioblastoma multiforme who received multiple infusions of IL13Rα2-CAR T cells delivered first intratumorally after tumor resection and then intraventricularly to the cerebral spinal fluid. Although intratumoral therapy appeared to control local tumor recurrence, complete elimination of multifocal tumors, including lesions in the spine, was observed only after intraventricular delivery. The presented data further suggest that intraventricular but not intratumoral delivery yields cell distributions capable of surveying the central nervous system. When labeling the clinical product at 70 kBq/106 cells for imaging, we expect doses of up to 14 MBq per intracranial infusion. Although 40% less than intravenous 89Zrantibodies (32), the smaller volume of distribution within the cerebral spinal fluid should yield comparable activities per unit volume. This experience highlights the importance of methods, such as ⁸⁹Zroxine cell imaging, that enable the tracking of adoptively transferred CAR T cells for optimizing routes of delivery and dosing schedules.

We demonstrated that human CAR T cells can be labeled and imaged with ⁸⁹Zr-oxine. The successful labeling of CAR T cells was in part achieved by generating ⁸⁹Zr-oxine from ⁸⁹Zr-oxalate at a reproducibly high efficiency without further purification. Previously reported methods for generating ⁸⁹Zr-oxine have required extraction with chloroform (*12,13*) or conversion to ⁸⁹ZrCl₄ and addition of polysorbate-80 (*14,15*), which was excluded to avoid any unanticipated cell toxicity and only marginally increases the radiolabeling efficiency (92%–97%). Cell labeling with ⁸⁹Zr-oxine was more efficient than with the previous methods (75% vs. 44% (*14*) or 61% (*13*)), possibly because of the cell type or the higher cell concentration used in our labeling. Importantly, we showed that ⁸⁹Zr-oxine

labeling at 70 kBq/10⁶ cells did not affect CAR T cell in vitro function or antitumor efficacy. The toxicity of higher radioactive doses on cell viability and function was not due to chemical burden and did not present until several days after labeling. Previously published studies also used ⁸⁹Zr-oxine at specific activities of 32–334 kBq/10⁶ murine bone marrow and CD8-positive T cells (14,15), an activity range that approximates the optimal 70 kBq/10⁶ human CAR T cells determined in our study. These studies provide evidence that the tracking of ⁸⁹Zr-oxine–labeled cells by PET imaging is representative of their unlabeled counterparts.

We also visualized the trafficking of intravenous 89Zr-oxine-labeled PSCA-CAR T cells to subcutaneous tumors with PET imaging, demonstrating the utility of this platform for tracking systemically delivered cells. As a control, we found that intravenously injected 89Zr-oxine rapidly bound serum proteins and was mostly bone- and liverassociated after 6 d. This bone-seeking phenomenon is similar to that reported previously for weakly chelated ⁸⁹Zr (33). It follows that 89Zr lost from cells because of death or efflux can be estimated from the activity of the bone. In a previous study, Sato et al. (14) showed tumor trafficking and therapy of intravenous 89Zr-oxine-labeled CD8-positive T cells to subcutaneous tumors; however, T cells in the tumor were

not confirmed, and ⁸⁹Zr-oxine–labeled cell therapy was not compared with an unlabeled control. A small percentage of our PSCA-CAR T cell dose accumulated in tumors over 7 d (0.48% and 0.70%), which is equivalent to 37,000 and 54,000 cells of the initial cell dose. This finding could be due to a limited number of cells reaching arterial circulation because of pulmonary entrapment (*34*), inadequate levels of chemotactic stimuli and adhesion molecules (*35*), or downregulation of the tumor antigen (*36*). By visualizing and quantifying cell uptake, future studies will investigate the properties of CAR T cells and tumors that affect optimal tumor trafficking and infiltration.

An inherent limitation of 89Zr-oxine or other ex vivo labeling methods is that the measured signal reflects the number of injected cells and does not account for cell division or death. For instance, a minority of injected cells may proliferate and initiate the majority of the antitumor response (37), and 89Zr-oxine would detect only the preliminary migration. Modifying CAR T cells with reporter genes would account for proliferation after injection; in fact, our IL13Rα2-CAR T cells have been modified to express herpes simplex virus type 1 thymidine kinase and imaged with PET in patients (38). This approach enables longitudinal imaging if transfected cells are present but is limited by reduced specificity, an additional ex vivo transfection, and potential immunogenicity. If both imaging strategies are used, improvements in CAR T cell therapy can be evaluated, such as enhanced ex vivo expansion methods (39), cell engineering (40), or tumor invasion (41). Thus, both ⁸⁹Zr-oxine cell labeling and reporter genes are complementary methods for evaluating the percentage of injected dose and the total number of effector cells at a given time within tumors, respectively.

CONCLUSION

89Zr-oxine was generated using a rapid, biocompatible method that efficiently labeled CAR T cells with minimal loss of viability or function at the optimal specific activity of 70 kBq/10⁶ cells. Preclinically, this method can be used to optimize the tumor targeting of different CAR T cell therapies, such as those with newly engineered CARs, tumor models, or routes of injection. This method can also be used to evaluate the effect of immunomodulatory compounds or radiation therapy on the migration of ⁸⁹Zr-oxine-labeled cells. Because both ⁸⁹Zr and oxine formulations are already used for clinical imaging, ⁸⁹Zr-oxine labeling should have immediate clinical utility for evaluating the dynamics of CAR T cell trafficking. This information would provide real-time evaluation of cell therapy, offering clinicians the opportunity to intervene or continue therapy if tumor targeting is observed.

DISCLOSURE

Research reported in this publication was supported by the National Cancer Institute of the National Institutes of Health under awards P30CA033572 and R01FD005129, the Prostate Cancer Foundation Challenge Award, Gateway for Cancer Research (G-14-600), and the California Institute for Regenerative Medicine (CIRM; TR3-05641). No other potential conflict of interest relevant to this article was reported. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

REFERENCES

- Brown CE, Alizadeh D, Starr R, et al. Regression of glioblastoma after chimeric antigen receptor T-cell therapy. N Engl J Med. 2016;375:2561–2569.
- Priceman SJ, Forman SJ, Brown CE. Smart CARs engineered for cancer immunotherapy. Curr Opin Oncol. 2015;27:466–474.
- Debinski W, Gibo DM, Hulet SW, Connor JR, Gillespie GY. Receptor for interleukin 13 is a marker and therapeutic target for human high-grade gliomas. Clin Cancer Res. 1999;5:985–990.
- Thomas C, Tampe R. Proofreading of peptide-MHC complexes through dynamic multivalent interactions. Front Immunol. 2017;8:65.
- Maus MV, Grupp SA, Porter DL, June CH. Antibody-modified T cells: CARs take the front seat for hematologic malignancies. Blood. 2014;123:2625–2635.
- 6. Ramos CA, Savoldo B, Dotti G. CD19-CAR trials. Cancer J. 2014;20:112-118.
- Batlevi CL, Matsuki E, Brentjens RJ, Younes A. Novel immunotherapies in lymphoid malignancies. Nat Rev Clin Oncol. 2016;13:25–40.
- Gross G, Waks T, Eshhar Z. Expression of immunoglobulin-T-cell receptor chimeric molecules as functional receptors with antibody-type specificity. *Proc Natl Acad Sci USA*. 1989;86:10024–10028.
- Kakarla S, Gottschalk S. CAR T cells for solid tumors: armed and ready to go? Cancer J. 2014;20:151–155.
- Bonifant CL, Jackson HJ, Brentjens RJ, Curran KJ. Toxicity and management in CAR T-cell therapy. Mol Ther Oncolytics. 2016;3:16011.
- Sadelain M, Brentjens R, Riviere I. The promise and potential pitfalls of chimeric antigen receptors. Curr Opin Immunol. 2009;21:215–223.
- Ferris TJ, Charoenphun P, Meszaros LK, Mullen GE, Blower PJ, Went MJ. Synthesis and characterisation of zirconium complexes for cell tracking with Zr-89 by positron emission tomography. *Dalton Trans.* 2014;43:14851–14857.
- Charoenphun P, Meszaros LK, Chuamsaamarkkee K, et al. [89Zr]oxinate4 for long-term in vivo cell tracking by positron emission tomography. Eur J Nucl Med Mol Imaging. 2015;42:278–287.
- Sato N, Wu H, Asiedu KO, Szajek LP, Griffiths GL, Choyke PL. 89Zr-oxine complex PET cell imaging in monitoring cell-based therapies. *Radiology*. 2015; 275:490–500.
- Asiedu KO, Koyasu S, Szajek LP, Choyke PL, Sato N. Bone marrow cell trafficking analyzed by ⁸⁹Zr-oxine positron emission tomography in a murine transplantation model. *Clin Cancer Res.* 2017;23:2759–2768.

- Roca M, de Vries EF, Jamar F, Israel O, Signore A. Guidelines for the labelling of leucocytes with ¹¹¹In-oxine: Inflammation/Infection Taskgroup of the European Association of Nuclear Medicine. *Eur J Nucl Med Mol Imaging*. 2010;37:835– 841
- Stanton SE, Disis ML. Clinical significance of tumor-infiltrating lymphocytes in breast cancer. J Immunother Cancer. 2016;4:59.
- Nicol AJ, Tazbirkova A, Nieda M. Comparison of clinical and immunological effects of intravenous and intradermal administration of alpha-galactosylceramide (KRN7000)-pulsed dendritic cells. Clin Cancer Res. 2011;17:5140–5151.
- Livieratos L. Technical pitfalls and limitations of SPECT/CT. Semin Nucl Med. 2015;45:530–540.
- Botti C, Negri DR, Seregni E, et al. Comparison of three different methods for radiolabelling human activated T lymphocytes. Eur J Nucl Med. 1997;24:497–504.
- Jauw YW, Menke-van der Houven van Oordt CW, Hoekstra OS, et al. Immunopositron emission tomography with zirconium-89-labeled monoclonal antibodies in oncology: what can we learn from initial clinical trials? Front Pharmacol. 2016;7:131.
- Kelly JD, inventor; GE Healthcare Ltd., assignee. Indium-111 oxine complex composition. U.S. patent US4335095. June 15, 1982.
- Thakur ML, Segal AW, Louis L, Welch MJ, Hopkins J, Peters TJ. Indium-111labeled cellular blood components: mechanism of labeling and intracellular location in human neutrophils. *J Nucl Med.* 1977;18:1022–1026.
- Brown CE, Starr R, Martinez C, et al. Recognition and killing of brain tumor stem-like initiating cells by CD8+ cytolytic T cells. Cancer Res. 2009;69:8886– 8893
- Brown CE, Wright CL, Naranjo A, et al. Biophotonic cytotoxicity assay for highthroughput screening of cytolytic killing. J Immunol Methods. 2005;297:39–52.
- Jonnalagadda M, Mardiros A, Urak R, et al. Chimeric antigen receptors with mutated IgG4 Fc spacer avoid fc receptor binding and improve T cell persistence and antitumor efficacy. Mol Ther. 2015;23:757–768.
- Hong H, Brown CE, Ostberg JR, et al. L1 cell adhesion molecule-specific chimeric antigen receptor-redirected human T cells exhibit specific and efficient antitumor activity against human ovarian cancer in mice. *PLoS One*. 2016;11: e0146885.
- Priceman SJ, Gerdts EA, Tilakawardane D, et al. Co-stimulatory signaling determines tumor antigen sensitivity and persistence of CAR T cells targeting PSCA+ metastatic prostate cancer. Oncoimmunology. 2017;7:e1380764.
- Brown CE, Vishwanath RP, Aguilar B, et al. Tumor-derived chemokine MCP-1/CCL2 is sufficient for mediating tumor tropism of adoptively transferred T cells. J Immunol. 2007;179:3332–3341.
- Kahlon KS, Brown C, Cooper LJ, Raubitschek A, Forman SJ, Jensen MC. Specific recognition and killing of glioblastoma multiforme by interleukin 13-zetakine redirected cytolytic T cells. *Cancer Res.* 2004;64:9160–9166.
- Brown CE, Badie B, Barish ME, et al. Bioactivity and safety of IL13Ralpha2redirected chimeric antigen receptor CD8+ T cells in patients with recurrent glioblastoma. Clin Cancer Res. 2015;21:4062–4072.
- Deri MA, Zeglis BM, Francesconi LC, Lewis JS. PET imaging with ⁸⁹Zr: from radiochemistry to the clinic. *Nucl Med Biol.* 2013;40:3–14.
- Abou DS, Ku T, Smith-Jones PM. In vivo biodistribution and accumulation of 89Zr in mice. Nucl Med Biol. 2011;38:675–681.
- Fischer UM, Harting MT, Jimenez F, et al. Pulmonary passage is a major obstacle for intravenous stem cell delivery: the pulmonary first-pass effect. Stem Cells Dev. 2009;18:683–692.
- Slaney CY, Kershaw MH, Darcy PK. Trafficking of T cells into tumors. Cancer Res. 2014;74:7168–7174.
- Gardner R, Wu D, Cherian S, et al. Acquisition of a CD19-negative myeloid phenotype allows immune escape of MLL-rearranged B-ALL from CD19 CAR-T-cell therapy. *Blood*. 2016;127:2406–2410.
- Sathaliyawala T, Kubota M, Yudanin N, et al. Distribution and compartmentalization of human circulating and tissue-resident memory T cell subsets. *Immunity*. 2013;38:187–197.
- Keu KV, Witney TH, Yaghoubi S, et al. Reporter gene imaging of targeted T cell immunotherapy in recurrent glioma. Sci Transl Med. 2017;9:eaag2196.
- Song GY, Hong L, Zhang H, Zhang J, Li Y, Feng H. Clinical outcomes of combined lateral extra-articular tenodesis and intra-articular anterior cruciate ligament reconstruction in addressing high-grade pivot-shift phenomenon. *Arthroscopy.* 2016;32:898–905.
- Jackson HJ, Rafiq S, Brentjens RJ. Driving CAR T-cells forward. Nat Rev Clin Oncol. 2016;13:370–383.
- Scarfò I, Maus MV. Current approaches to increase CAR T cell potency in solid tumors: targeting the tumor microenvironment. *J Immunother Cancer*. 2017:5:28.