

Received: 21 January 2019 Accepted: 27 March 2019 Published online: 10 April 2019

OPEN Loop-mediated isothermal amplification (LAMP) approach for detection of heat-resistant Talaromyces flavus species

Jacek Panek 🗅 & Magdalena Frąc

Talaromyces flavus is a soilborne fungus that can contaminate fruits. It constitutes serious influence on heat-processed food spoilage, as T. flavus belongs to the heat-resistant fungi group, which are able to survive the pasteurization process. Moreover T. flavus has been reported to be capable of mycotoxigenicity, therefore they have a serious threat to human health. To maintain the safety of food production, sensitive method for T. flavus detection was developed. The loop mediated amplification, abbreviated LAMP, reactions were designed as specific for detection of DNA replication licensing factor gene of T. flavus. The specificity of assay was confirmed by use of 5 T. flavus strains and 35 other fungal isolates. The achieved limit of detection was 1fq of T. flavus genomic DNA and 64 ascospores in 1 q of strawberry fruits or soil samples.

Fungal spoilage of heat-processed food emerges as a major problem in food industry. Estimated losses caused to fruit juice production by heat-resistant fungi (HRF) cost food industry millions of USD solely in United States of America^{1,2}. HRF were reported to inhabit soil and contaminate agricultural products in countries around the world, making them a worldwide threat^{1,3-8}. Aforementioned spoilage origins from HRF that contaminate raw agricultural products by either anemochory or physical contact with soil that is inhabited by HRF such as Talaromyces flavus^{3,5,9}.

T. flavus is one of the most common fungi among HRF group. The HRF are able to withstand high temperature treatment, such as pasteurization process 10,11. Heat resistance of T. flavus labelled as decimal reduction time (D-value) was determined to achieve D_{90} °C = 6 minutes and D_{95} °C = 1 minute¹². Furthermore, *T. flavus* was reported to possess capability of producing plenty of mycotoxins and to grow under microaerobic and aerobic conditions¹³⁻¹⁶. What is also significant, other HRF Neosartorya fischeri was reported to be resistant to such growth control factors as fungicides or antibiotics¹⁷. These findings and these traits further establish *T. flavus* and HRF as a potent risk factor to the food industry and a risk for consumer health. Moreover, currently the best remedy to provide HRF-free heat-processed products is detection of the HRF in raw products or in soil under plantations. To achieve this, successful methods based on PCR and qPCR techniques to detect HRF of genus Neosartorya, Byssochlamys and species Talaromyces flavus were developed 18-20. Yet, emergence of isothermal amplification techniques may provide the possibility to develop HRF detection methods that can be regarded as easier and cheaper assay to deploy within food industry and on the field. Among the isothermal amplification methods, the ones that are used the most commonly, due to necessity to use only one enzyme, are Rolling circle amplification (RCA), Cross-priming amplification (CPA), Polymerase spiral reaction (PSR) and Loop-mediated isothermal amplification (LAMP)^{21–24}.

The aim of this study was to develop LAMP method for detection of Talaromyces flavus that provides specificity and rapidity of achieving results. LAMP is a technique that due to use of 6 primers that are complementary to 8 regions of DNA sequence assures very high specificity^{24,25}. Moreover, the isothermal character of reaction and presence of many priming sites for DNA synthesis provides possibility to observe results earlier than in competing techniques such as PCR or qPCR. Apart from that, LAMP reaction provides possibility to detect occurrence of amplification product by either measurement fluorescence of minor-groove dsDNA binding dyes, by measurement of increasing turbidity of reaction or by colorimetric detection using pH sensitive dyes. It is important to mention, that described detection methods can render the results visible to the naked eye. Moreover, laboratory

Institute of Agrophysics, Polish Academy of Sciences, Doświadczalna 4, 20-290, Lublin, Poland. Correspondence and requests for materials should be addressed to J.P. (email: j.panek@ipan.lublin.pl)

Species	Targeted sequence	Primer name	Primer sequences 5'→3'	Final concentration
		Tf_F3	CTGCTGTTATGCGAGACC	0.25 μΜ
	DNA	Tf_B3	CATACAGAGGGTTGGCTG	0.25 μΜ
Talaromyces flavus	replication	Tf_FIP	CATCGATRCAGCAGATACCGTGATGAGATGGTRCTGGAAG	0.5 μΜ
Tataromyces jiavus	licensing factor gene	Tf_BIP	CGACAAGATGGACGACTCRGCCGGCCTTGGAGATAGAA	0.25 μΜ
	lactor gene	Tf_LoopF	TATCAGCGAGCACAAGASC	0.5 μΜ
		Tf_LoopB	TCACGAAGTTATGGAACAGCA	0.25 μΜ

Table 1. Primers for the LAMP assay.

KX587558.1 Talaromyces flavus strain DSM 63536 DNA replication licensing factor gene, partial cds




Figure 1. Positioning and orientation of regions used to design LAMP primers within the nucleotide sequence of the DNA replication licensing factor gene of *T. flavus* DSM 63536 (GenBank accession no. KX587558).

equipment such as qPCR thermocyclers, fluorimeters or turbidimeters may be used to visualise the results^{26–31}. The target of developed reaction was *Talaromyces flavus DNA replication licensing factor* (*rlf*) gene. The DNA licensing factor is a protein complex that allows replication to begin at the replication origin site and is a part of the replication licensing system, which regulates the replication and ensures that no section of genome is replicated more than once in a single cell cycle³². The *DNA replication licensing factor* has already been reported as a possible marker and target gene for detection of *T. flavus* by qPCR reaction²⁰.

Material and Methods

Development of the LAMP method. In order to design LAMP primers, 10 nucleotide sequences of *T. flavus rlf gene* (accession number GU396448.1 - GU396452.1, KX587558-KX587562) were obtained from GenBank database and aligned with MEGA6 software³³. After that, possible regions for primer design were found and checked by BLAST performed within GenBank database to ensure the uniquity and specificity of them³⁴. Then, the primers were designed with help of LAMP Designer software (Optigene) and once again confirmed by BLAST to reassure their specificity (Table 1 and Fig. 1). During the assessment of specificity of the designed LAMP reaction, DNA extracted from 40 fungi isolates (Table 2) was used. Employed strains were identified prior by sequencing of internal transcribed spacer 1 (ITS1) and/or D2 region of large subunit of rRNA gene (D2 LSU rRNA), with comparison to MicroSeq, GenBank and Mycobank databases. The studied pool of microorganisms consisted of 5 *T. flavus* strains commercially available from Leibniz Institute DSMZ-German Collection of Microorganisms and Cell Cultures (Braunschweig, Germany) and from NITE Biological Resorce Center (Tokyo, Japan), and collection of 35 different isolates containing other HRF, cosmopolite fungi and fungi isolated from soil and strawberries belonging to genus *Talaromyces, Penicillium, Aspergillus, Neosartorya* and *Byssochlamys*.

The tested fungal species were grown on Potato Dextrose Agar (PDA) at 27 °C for 7 days. The DNA isolation and quality assessment was performed as in previous work²⁰. The mycelium was homogenized in sterile conditions

Fungi	Reference	Isolation source	Accession number of sequences in GenBank (*RLF, bD2 LSU, cTTS1)	LAMP results
Talaromyces flavus	DSM 63536	Nepal; soil	^a KX587558, ^b KX587552, ^c MG878427	+
Talaromyces flavus	NBRC 31879	Unknown	^a KX587559, ^b KX587553, ^c MG878428	+
Talaromyces flavus	NBRC 30574	Nepal; soil	^a KX587560, ^b KX587554, ^c MG878429	+
Talaromyces flavus	NBRC 7233	Japan; unknown	^a KX587561, ^b KX587555, ^c MG878430	+
Talaromyces flavus	NBRC 102293	Japan; grassland	^a KX587562, ^b KX587556, ^c MG878431	+
Talaromyces macrosporus	NBRC 30070	Unknown; photo album	^b KX587557, ^c MG878432	_
Byssochlamys nivea	LMEM G49/11	Poland; fruit-based silage	^b MG877737	_
Neosartorya fischeri	LMEM G48/12	Poland; strawberry jam	^b KC179765	-
Byssochlamys fulva	DSM 62097	Germany; bottled apple juice	^b MG877729	-
Byssochlamys fulva	DSM 1808	Unknown; grape juice	^b MG877730	-
Neosartorya glabra	LMEM G5/14	Poland; soil	^b MG877735	-
Talaromyces ruber	LMEM G141/14	Poland; soil	^b KM822632	-
Talaromyces islandicus	LMEM G142/14	Poland; soil	^b KM822633	-
Penicillium camembertii	LMEM G77/11	Poland; grass silage	^b MG877736	_
Aspergillus niger	LMEM G99/13	Poland; soil	^b MG877738	-
Aspergillus niger	LMEM G100/13	Poland; soil	^b MG877739	-
Aspergillus niger	LMEM G101/13	Poland; soil	^b MG877740	-
Aspergillus niger	LMEM G102/13	Poland; soil	^b MG877741	-
Aspergillus niger	LMEM G105/13	Poland; soil	^b MG877742	-
Penicillium miczynskii	LMEM G4/18	Poland; soil	^b MG877731	-
Talaromyces angelicus	LMEM G5/18	Poland; soil	^b MG877732	-
Aspergillus fumigatus	LMEM G6/18	Poland; soil	^b MG877733	-
Talaromyces angelicus	LMEM G7/18	Poland; strawberries	^b MG877734	-
Talaromyces angelicus	LMEM G8/18	Poland; strawberries	^b MG881801, ^c MG878420	-
Talaromyces angelicus	LMEM G9/18	Poland; strawberries	^b MG881802, ^c MG878421	-
Talaromyces kabodanensis	LMEM G10/18	Poland; strawberries	^b MG881803, ^c MG878422	-
Talaromyces kabodanensis	LMEM G11/18	Poland; strawberries	^b MG881804, ^c MG878423	-
Talaromyces kabodanensis	LMEM G12/18	Poland; strawberries	^b MG881805, ^c MG878424	-
Talaromyces kabodanensis	LMEM G13/18	Poland; strawberries	^b MG881806, ^c MG878425	-
Talaromyces aculeatus	LMEM G14/18	Poland; strawberries	^b MG881807, ^c MG878426	-
Talaromyces kabodanensis	LMEM G15/18	Poland; soil	^b MG881808, ^c MG878410	-
Talaromyces kabodanensis	LMEM G16/18	Poland; soil	^b MG881809, ^c MG878411	-
Talaromyces kabodanensis	LMEM G17/18	Poland; soil	^b MG881810, ^c MG878412	-
Talaromyces kabodanensis	LMEM G18/18	Poland; soil	^b MG881811, ^c MG878413	-
Talaromyces kabodanensis	LMEM G19/18	Poland; soil	^b MG881812, ^c MG878414	-
Talaromyces kabodanensis	LMEM G20/18	Poland; soil	^b MG881813, ^c MG878415	-
Talaromyces kabodanensis	LMEM G21/18	Poland; soil	^b MG881814, ^c MG878416	-
Talaromyces kabodanensis	LMEM G22/18	Poland; soil	^b MG881815, ^c MG878417	-
Talaromyces kabodanensis	LMEM G23/18	Poland; soil	^b MG881816, ^c MG878418	-
Talaromyces kabodanensis	LMEM G24/18	Poland; soil	^b MG881817, ^c MG878419	_

Table 2. List of fungal strains used and LAMP specificity results. DSM – Leibniz Institute DSMZ-German Collection of Microorganisms and Cell Cultures, Braunschweig, Germany; NBRC – National Institute of Technology and Evaluation, Tokyo, Japan; LMEM – Laboratory of Molecular and Environmental Microbiology, Institute of Agrophysics, Polish Academy of Sciences, Lublin, Poland.

with using of FastPrep-24 instrument (MP Biomedicals) under the following settings: $6\,\text{m/s}$, $40\,\text{seconds}$, and with the matrix containing glass beads with diameter $3.15\,\text{mm}$ ($500\,\text{mg}$) and $1.45\,\text{mm}$ ($250\,\text{mg}$), and reagents: $400\,\mu\text{l}$ of Lyse F buffer (EURx), $10\,\mu\text{l}$ of Proteinase K (EURx) and $3\,\mu\text{l}$ of RNase A (EURx). Next, the mycelium was sonicated in $240\,\text{W}$ ultrasonic cleaner (Polsonic SONIC-6) for $10\,\text{minutes}$, After disintegration the fungal genomic DNA was extracted according to manufacturer protocol with Plant & Fungi DNA Purification Kit (EURx) and suspended in $200\,\mu\text{l}$ of Tris-HCl buffer ($10\,\text{mM}$ Tris-HCl, pH 8.5). The quality and quantity of the DNA was determined spectrophotometrically by NanoDrop instrument (Thermo Fisher Scientific) and fluorometrically with Quantus instrument (Promega) and deemed as good for further analysis ($A_{260}/A_{280} = 2.05\,\text{SD} = 0.08$ while $A_{260}/A_{230} = 1.94\,\text{SD} = 0.06$).

In order to determine if developed method can be used with quick DNA extraction methods, PrepMan[™] Ultra Sample Preparation Reagent (Thermo Fisher Scientific) was used. The mycelium of tested fungal species, and also soil samples or strawberry fruits were transferred with using 10 µl inoculation loop into 1.5 ml Eppendorf

tube containing 100 μ l of PrepMan reagent, thoroughly vortexed and heated to 100 °C for 10 minutes. After that, samples were cooled to the room temperature for 2 minutes and centrifuged (120 s, 20 000 xg). 50 μ l of supernatant was preserved (-20 °C) for reactions.

The LAMP reactions were performed with 7500 Fast thermocycler (Applied Biosystems). The reaction contained 6 µl of Isothermal MasterMix (ISO-001, OptiGene), 3 µl of PrimerMix (Tables 1) and 1 µl of the DNA sample. Total reaction volume was 10 µl. Reactions were conducted in constant temperature of 65 °C through 90 min with fluorescence data collection every 1 minute. Isothermal MasterMix contained GspSSD polymerase with strand displacement activity and fluorescent dye intercalating to dsDNA minor-groove compatible with FAM channel. To determine the detection limit of developed reactions serial 10-fold dilutions of *T. flavus* DSM 63536 genomic DNA solution obtained with EURx method (100 ng/µl) and PrepMan method (0.1 ng/µl) were prepared. Concentrations of 10 ng/µl, 1 ng/µl, 100 pg/µl, 10 pg/µl, 10 pg/µl, 100 fg/µl, 10 fg/µl and 1 fg/µl (EURx method) and 100 pg/μl, 10 pg/μl, 10 pg/μl, 100 fg/μl, 10 fg/μl and 1 fg/μl (PrepMan method) were obtained and 1 μl of each was used as a template for LAMP reaction. Sterilized DirectQ water was used as non-template control (NTC). All reactions were performed in biological triplicates. The results were visualized as amplification plots (ΔR) expressed as raw fluorescence units (RFU) in time. To obtain ΔR , on the basis of the raw fluorescence values, a baseline was subtracted from the raw data. Presence or absence of amplification plot determined the positivity or negativity of tested reactions. Moreover, cycle thresholds were calculated with assumption that each data collection step corresponds to 1 cycle. Each reaction was followed with melting curve analysis (65 °C to 95 °C, $\Delta 0.016$ °C/s).

Application of LAMP method using environmental samples. To determine applicability of developed LAMP method for detection of *T. flavus* in food and soil samples the LAMP assays were performed with DNA isolated from artificially contaminated environmental samples - strawberry fruits and soil, which is a source of fungal contamination of agricultural products and food.

In this study, soil samples, which were obtained from strawberry plantations in Poland, Lubelskie voivodeship and strawberries fruits of Fragaria × ananassa 'Senga Sengana' cultivar were contaminated with ascospores of T. flavus. To confirm absence of natural occurrence of T. flavus in examined environmental samples the method of Beuchat and Pitt⁴ was used for fungal strains isolation. Achieved isolates were additionally analysed and identified by sequencing of D2LSU rRNA gene. The results confirmed, that studied samples were free of natural occurrence of T. flavus. Environmental samples were not sterilised prior to study to ensure, that designed method is able to detect T. flavus DNA in presence of other, autochthonic DNA. T. flavus DSM 63536 ascospores were suspended in sterile water and vortexed for 30 seconds. The concentration of suspension was measured microscopically with Thoma's counting chamber and resulted in 64000 ascospores per 1 ml of sterile water. Then, serial 10-fold dilutions (64000, 6400, 640 and 64 ascospores per 1 ml of water) of suspension were prepared. After that, 1 ml of each prepared inoculum was added to 1 g of strawberry fruits and soil samples and then, mixtures were vortexed for 15 seconds. Next, 200 ml of contaminated sample was transferred to sterile tubes. NTC samples were prepared by addition of sterile water to soil and strawberry fruits samples instead of inoculum. To isolate DNA from contaminated strawberry fruits, already described above procedure involving EURx Plant & Fungi DNA Purification kit and also PrepMan Ultra Sample Preparation Reagent were used. Extraction from contaminated soil samples was conducted with Soil DNA Purification Kit (EURx) and PrepMan Ultra Sample Preparation Reagent according to manufacturer's protocols.

Validation of developed method. In order to further determine the potential applicability of developed approach, validation tests were performed on strawberry fruits and soil samples obtained from strawberry plantations located in Lublin Voivodeship, Poland. Samples were prior examined with use of standard method for detection of HRF⁴ with addition of sequencing of D2LSU rRNA gene. DNA was extracted from 10 soil and 10 strawberry fruits samples with above mentioned EURx and PrepMan methods.

Results

Specificity and sensitivity of LAMP method. Developed LAMP reaction for detection of *T. flavus* provided specific results. Positive amplification results were observed with all five *T. flavus* isolates: *T. flavus* DSM 63536, *T. flavus* NBRC 31879, *T. flavus* NBRC 30574, *T. flavus* NBRC 7233 and *T. flavus* NBRC 102293, and with *T. flavus* DSM 63536 DNA mixed with *Byssochlamys nivea* LMEM G49/11 DNA in ratio of 100fg:10 ng of genomic DNA, while no amplification was observed with DNA templates isolated from other tested 35 fungal isolates (Fig. 2a). Calculated Ct of positive reactions was 19.47, SD ± 2.72 . Figure 2b presents melting curves of obtained amplification products. The fact that proposed reaction did not amplify DNA isolated from fungi commonly found in environmental samples belonging to genus *Penicillium* and *Aspergillus*, nor from other common HRF like *Byssochlamys*, *Neosartorya* and other species of *Talaromyces* sp., isolated from soil or strawberry fruits samples, suggest that developed reaction is very specific. What is more, the same results were obtained for fungal DNA isolated with PrepMan extraction method. Positive reaction observed in samples containing mixed DNA of *T. flavus* and *B. nivea* in ratio of 1:100000 suggest, that presence of external DNA does not interfere with amplification.

Sensitivity of LAMP method for detection of *T. flavus*. The amplification plots visible at the Fig. 3 represent the results of LAMP sensitivity assay. In order to assess the sensitivity of developed LAMP approach, reactions were prepared with serially diluted DNA matrices. With use of EURx DNA extraction method, positive results were observed for all tested DNA matrix concentrations ranging from 10 ng/µl to 1 fg/µl, yet no amplification was observed in non-template controls (Fig. 3a). Obtained results show that the limit of detection for proposed reaction is 1 fg/µl. Results for DNA matrix concentrations of 10 ng/µl to 100 fg/µl were observed in about




Figure 2. The specificity of *Talaromyces flavus* LAMP detection method. (a) The amplification plots (ΔR) of *T. flavus* and 35 other studied isolates (Table 2) in time. Legend describes tested *T. flavus* isolates, while rest of tested strains are described in Table 2. NTC stands for non-template control. (b) Melting curve analysis of obtained products. DNA isolated with EURx Plant&Fungi DNA Purification Kit.




Figure 3. Sensitivity study of *Talaromyces flavus* LAMP detection method. (a) The amplification plots (ΔR) of pure genomic *T. flavus* DSM 63536 DNA added into reaction in different concentrations in time. (b) The amplification plot (ΔR) of *T. flavus* detected in strawberries contaminated by different number of ascospores added to 1 g of sample (a/g) prior the DNA isolation reaction. DNA isolated with EURx Plant&Fungi DNA Purification Kit. (c) The amplification plot (ΔR) of *T. flavus* detected in soil contaminated by different number of ascospores added to 1 g of sample prior the DNA isolation reaction. NTC stands for non-template control. DNA isolated with EURx Soil DNA Purification Kit.

30 minutes (Ct of 23.29, SD \pm 2.97), while positive results of reactions with DNA matrices of 10 fg/µl and 1 fg/µl were observed after about 45 (Ct of 37.42, SD \pm 1.29) and 55 minutes (Ct of 44.68, SD \pm 4.62), respectively. With deployment of PrepMan DNA extraction method, the positive results were observed in range of 100 pg/µl to 10 fg/µl

in all replications (Ct 23.03, SD \pm 1.35; Ct 35.99, SD \pm 1.81; Ct 38.6, SD \pm 2.46; Ct 55.27, SD \pm 3.48; Ct 71.2, SD \pm 3.52). As for matrix containing 1 fg/µl, weak positive results were observed in only one repetition (Ct 88.72).

The sensitivity assessments of *T. flavus* detection from contaminated strawberries and soil are presented at Fig. 3b,c, respectively. LAMP reactions resulted in positive amplifications plots. The assay was able to detect occurrence of DNA isolated with EURx methods from samples inoculated with 64 ascospores per 1 g of each tested environmental sample types: soil and strawberry fruits. Time needed to observe positive results was similar for both kinds of samples, and resulted in about 90 minutes, about 75 minutes, about 40 minutes and about 30 minutes for DNA isolated from 64, 640, 6400 and 64000 ascospores per 1 gram of sample, respectively. Simultaneously no amplification was observed in NTC. These results indicate that the developed LAMP method can accurately detect *T. flavus* ascospores in strawberry fruits and soil samples. Moreover, developed assay was able to positively amplify DNA isolated with PrepMan method from soil samples contaminated with 640 or more ascospores per 1 g. However, in case of DNA extracted with the same method from contaminated strawberry fruits, the detection was possible only in the highest concentration of ascospores (64 000 per 1 g).

Validation of developed method. Out of 10 tested soil samples 4 were prior identified as spoiled with *T. flavus*. Both tested method of DNA extraction allowed to detect *T. flavus* in those samples, however with EURx method, one more sample was positive. In case of strawberry fruits samples, out of 5 samples identified prior as contaminated with *T. flavus*, positive reaction results were observed in all of these 5 samples, when DNA was extracted with EURx method and only in one sample extracted with PrepMan method.

Discussion

Studies towards identification and preservation of fungal contamination in food products are important to reassure high quality and health safety of food products. Heat-resistant fungi play important role in endangering food production as they are able to survive the pasteurization process. Moreover, as they in general are not considered plant pathogens, symptoms of their occurrence on fruit may not be observed. These attributes give rise to claim that one of the best methods to provide HRF-free food production and safety involves early detection of HRF contamination in raw products9. It is important, that food processing industry often requires quality control methods that are rapid, sensitive and specific. The common, traditional agar plate based microbiological methods do not provide, especially in case of quickness and specificity³⁵. What is more, the other methods, such as liquid chromatography-mass spectrometry methods developed in recent years can determine mycotoxins to evaluate postharvest quality of fruits³⁶ by they don't detect the presence of fungi at the very early stage before contamination of samples by toxins. It is important to say, that presence of HRF contamination, even without mycotoxin secretion, presents threat to the food production, as it may spoil the final product during its shelf life. To fulfil these requirements, molecular biology techniques were vastly employed, especially to detect bacteria and fungi causing human diseases, yet the HRF occurrence detection seem to be neglected^{37,38}. In recent years, PCR and qPCR based methods for detection of HRF genus Neosartorya¹⁹, genus Thermoascus³⁸, genus Byssochalmys and Hamigera¹⁸ and species Talaromyces flavus²⁰. Due to emergence of isothermal techniques such as LAMP, researchers developed methods to utilize these techniques to detect microorganisms and viruses important for the animal and plant food processing and consumer health 30,39-44 and in general public health sensu largo 45-51. However, LAMP methods for detection of fungi are scarce, while there are no isothermal methods concentrated on detecting any of HRF, yet alone Talaromyces flavus.

In this study a loop-mediated isothermal amplification method for detection of *T. flavus* was developed. The developed assay was evaluated for detection of *T. flavus* in environmental samples, such as strawberry fruits and soil samples. Designed method and primers were confirmed to be specific and to discriminate *T. flavus* from other *Talaromyces* species, HRF and common cosmopolite fungi. Designed method targeted *Talaromyces flavus DNA replication licensing factor gene*, which was previously reported as a possible and promising marker gene for detection of *Talaromyces flavus*²⁰. The replication licensing system acts to ensure that the no section of the genome is replicated more than once in a single cell cycle, while the licensing factor is protein complex that allows an origin of replication to begin DNA replication at that site³².

Proposed LAMP assay for *T. flavus* detection was specific and sensitive. During studies amplification was observed with DNA isolated from 5 *Talaromyces flavus* strains, while no amplification was seen with DNA isolated from 35 different fungal isolates: *Talaromyces, Byssochlamys, Neosartorya, Penicillium* and *Aspergillus* genera. Detection limits of developed method were: 1 fg of genomic *T. flavus* DNA, DNA isolated from soil sample contaminated with 64 *T. flavus* ascospores per one gram of sample and DNA isolated from strawberries sample contaminated with 64 *T. flavus* ascospores per one gram of sample. Moreover, proposed method remained specific in occurrence of *Byssochlamys nivea* DNA and DNA of microorganisms naturally inhabiting soil and strawberries samples. It is important to state, that proposed LAMP method turned out to be 10 times more sensitive to detection of *T. flavus* in environmental samples then the qPCR method when the same gene was used as a target for amplification²⁰. Achieved limits of detection were significantly better than already published methods for detection of other fungi^{42,43}.

Conclusion

Summarizing, a sensitive and specific LAMP assay for heat resistant fungi of *T. flavus* species was developed. Proposed assay is suitable to detect referenced HRF effectively in both soil and strawberry samples. This method asserts higher sensitivity than qPCR method for detection of the same fungus, while maintain almost the same reagent cost per reaction (approximately 0.5\$/per reaction)²⁰, allowing to detect 1 fg of genomic *T. flavus* DNA and DNA isolated from soil or strawberry samples contaminated with 64 *T. flavus* ascospores per one gram of sample. This method contains the same drawback as many other methods based on molecular biology techniques – requires prior extraction of DNA, however very quick and simple extraction protocols, such as PrepMan

Ultra Sample Preparation Reagent DNA extraction method may be employed with lower, yet satisfying detection limits of 10 fg of genomic *T. flavus* DNA and 640 ascospores per gram of soil. Moreover, thanks to potential use of different, non-instrumental detection methods such as turbidimetric, intercalating dyes or pH sensitive dyes, this method may be deployed to detect presence of *T. flavus* under field conditions. The results of this study and proposed method may be applicable and adopted for quality monitoring of *T. flavus* contamination in agriculture and food-processing industry.

Data Availability

The datasets generated and analysed during the current study are available from the corresponding author on reasonable request.

References

- Dijksterhuis, J. Heat-Resistant Ascospores. In Food Mycology: A Multifaceted Approach to Fungi and Food (eds Dijksterhuis, J. & Samson, R. A.) 101–117 (CRC Press, 2007).
- 2. Tournas, V. Heat-Resistant Fungi of Importance to the Food and Beverage Industry. Crit. Rev. Microbiol. 20, 243-263 (1994).
- 3. Amaeze, N. J., Ugwuanyi, J. O. & Obeta, J. A. N. Studies Of Heat Resistant Fungi In The Soil: Talaromyces flavus Isolated In Nigerian Soils. New York Sci. J. 3, 8–14 (2010).
- 4. Beuchat, L. R. & Pitt, J. L. Detection and enumeration of heat-resistant molds. In Compendium of the methods for the microbiological examination of food (eds Downes, F. P. & Ito, K.) 217–222 (APHA, 2001).
- 5. Dethoup, T., Manoch, L., Visarathanonth, N., Chamswarng, C. & Kijjoa, A. Morphology and Distribution of *Talaromyces flavus* from Soil and Potential Use as a Biological Control Agent against Plant Pathogenic Fungi. *Thai J. Agric. Sci.* 40, 37–50 (2007).
- 6. Manoch, L., Dethoup, T., Yilmaz, N., Houbraken, J. & Samson, R. A. Two new *Talaromyces* species from soil in Thailand. *Mycoscience* 54, 335–342 (2013).
- 7. Naraghi, L., Heydari, A., Rezaee, S., Razavi, M. & Afshari-Azad, H. Biological control of Verticillium wilt of greenhouse cucumber by *Talaromyces flavus*. *Phytopathol. Mediterr.* **49**, 321–329 (2010).
- 8. Fravel, D. Ř. & Adams, P. B. Estimation of United States and World Distribution of Talaromyces flavus. Mycologia 78, 684 (1986).
- 9. Frac, M., Jezierska-Tys, S. & Yaguchi, T. Occurrence, Detection, and Molecular and Metabolic Characterization of Heat-Resistant Fungi in Soils and Plants and Their Risk to Human Health. *Adv. Agron.* 132, 161–204 (2015).
- 10. Jesenská, Z., Piecková, E. & Bernát, D. Heat resistance of fungi from soil. Int. J. Food Microbiol. 19, 187-192 (1993).
- 11. Samson, R. A. et al. Phylogeny and nomenclature of the genus *Talaromyces* and taxa accommodated in *Penicillium* subgenus. *Biverticillium*. Stud. Mycol. **70**, 159–183 (2011).
- 12. King, A. D. Heat resistance of *Talaromyces flavus* ascospores as determined by a two phase slug flow heat exchanger. *Int. J. Food Microbiol.* 35, 147–151 (1997).
- 13. Proksa, B. Talaromyces flavus and its metabolites. Chem. Pap. 64, 696-714 (2010).
- 14. Ayer, W. A. & Racok, J. S. The metabolites of *Talaromyces flavus*: Part 1. Metabolites of the organic extracts. *Can. J. Chem.* 11, 2085–2094 (1990).
- 15. Frisvad, J. C., Filtenborg, O., Samson, R. A. & Stolk, A. C. Chemotaxonomy of the genus *Talaromyces*. *Antonie Van Leeuwenhoek* 57, 179–189 (1990).
- 16. Zhai, M.-M. et al. The Bioactive Secondary Metabolites from Talaromyces species. Nat. Products Bioprospect. 6, 1-24 (2016).
- 17. Panek, J., Frac, M. & Bilińska-Wielgus, N. Comparison of Chemical Sensitivity of Fresh and Long-Stored Heat Resistant Neosartorya fischeri Environmental Isolates Using BIOLOG Phenotype MicroArray System. *PLoS One* 11, e0147605 (2016).
- 18. Nakayama, M. et al. A Rapid Method for Identifying Byssochlamys and Hamigera. J. Food Prot. 73, 1486-1492 (2010).
- 19. Yaguchi, T. et al. Method for Identifying Heat-Resistant Fungi of the Genus. Neosartorya. J. Food Prot. 75, 1806–1813 (2012).
- 20. Panek, J. & Frac, M. Development of a qPCR assay for the detection of heat-resistant *Talaromyces flavus*. Int. J. Food Microbiol. 270, 44–51 (2018).
- 21. Liu, W. et al. Polymerase Spiral Reaction (PSR): A novel isothermal nucleic acid amplification method. Sci. Rep. 5, 12723 (2015).
- Fang, R. et al. Cross-Priming Amplification for Rapid Detection of Mycobacterium tuberculosis in Sputum Specimens. J. Clin. Microbiol. 47, 845–847 (2009).
- 23. Banér, J., Nilsson, M., Mendel-Hartvig, M. & Landegren, U. Signal amplification of padlock probes by rolling circle replication. *Nucleic Acids Res.* 26, 5073–8 (1998).
- 24. Notomi, T. et al. Loop-mediated isothermal amplification of DNA. Nucleic Acids Res. 28, e63 (2000).
- 25. Nagamine, K., Hase, T. & Notomi, T. Accelerated reaction by loop-mediated isothermal amplification using loop primers. *Mol. Cell. Probes* 16, 223–229 (2002).
- 26. Njiru, Z. K. et al. African trypanosomiasis: Sensitive and rapid detection of the sub-genus *Trypanozoon* by loop-mediated isothermal amplification (LAMP) of parasite DNA. *Int. J. Parasitol.* **38**, 589–599 (2008).
- Mori, Y., Nagamine, K., Tomita, N. & Notomi, T. Detection of Loop-Mediated Isothermal Amplification Reaction by Turbidity Derived from Magnesium Pyrophosphate Formation. Biochem. Biophys. Res. Commun. 289, 150–154 (2001).
- Mori, Y., Kitao, M., Tomita, N. & Notomi, T. Real-time turbidimetry of LAMP reaction for quantifying template DNA. J. Biochem. Biophys. Methods 59, 145–157 (2004).
- 29. Tomita, N., Mori, Y., Kanda, H. & Notomi, T. Loop-mediated isothermal amplification (LAMP) of gene sequences and simple visual detection of products. *Nat. Protoc.* 3, 877–882 (2008).
- Woźniakowski, G., Kozdruń, W. & Samorek-Salamonowicz, E. Loop-mediated isothermal amplification for the detection of goose circovirus. Virol. J. 9, 110 (2012).
- 31. Tanner, N. A., Zhang, Y. & Evans, T. C. Visual detection of isothermal nucleic acid amplification using pH-sensitive dyes. *Biotechniques* 58, 59–68 (2015).
- 32. Tada, S. & Blow, J. J. The Replication Licensing System. *Biol. Chem.* **379**, 941–949 (1998).
- 33. Tamura, K., Stecher, G., Peterson, D., Filipski, A. & Kumar, S. MEGA6: Molecular Evolutionary Genetics Analysis Version 6.0. *Mol. Biol. Evol.* 30, 2725–2729 (2013).
- 34. Altschul, S. F., Gish, W., Miller, W., Myers, E. W. & Lipman, D. J. Basic local alignment search tool. J. Mol. Biol. 215, 403–10 (1990).
- 35. Arbefeville, S., Harris, A. & Ferrieri, P. Comparison of sequencing the D2 region of the large subunit ribosomal RNA gene (MicroSEQ®) versus the internal transcribed spacer (ITS) regions using two public databases for identification of common and uncommon clinically relevant fungal species. *J. Microbiol. Methods* 140, 40–46 (2017).
- Fornal, E., Parfieniuk, E., Czeczko, R., Bilinska-Wielgus, N. & Frac, M. Fast and easy liquid chromatography–mass spectrometry method for evaluation of postharvest fruit safety by determination of mycotoxins: Fumitremorgin C and verruculogen. *Postharvest Biol. Technol.* 131, 46–54 (2017).
- 37. Kordalewska, M., Brillowska-Dąbrowska, A., Jagielski, T. & Dworecka-Kaszak, B. PCR and real-time PCR assays to detect fungi of *Alternaria alternata* species. *Acta Biochim. Pol.* 62, 707–12 (2015).
- 38. Hosoya, K. et al. Risk analysis and rapid detection of the genus Thermoascus, food spoilage fungi. Food Control 41, 7-12 (2014).

- 39. Woźniakowski, G., Wencel, P. & Samorek-Salamonowicz, E. Detection of CoHV-1 by loop-mediated amplification (LAMP). Application of LAMP for CoHV-1 incidence monitoring in domestic pigeons. *Lett. Appl. Microbiol.* **59**, 610–4 (2014).
- 40. Woźniakowski, G., Samorek-Salamonowicz, E. & Kozdruń, W. Comparison of loop-mediated isothermal amplification and PCR for the detection and differentiation of Marek's disease virus serotypes 1, 2, and 3. Avian Dis. 57, 539–43 (2013).
- 41. Kursa, O., Woźniakowski, G., Tomczyk, G., Sawicka, A. & Minta, Z. Rapid detection of *Mycoplasma synoviae* by loop-mediated isothermal amplification. *Arch. Microbiol.* **197**, 319–25 (2015).
- 42. Almasi, M. A. et al. Development and Application of Loop-Mediated Isothermal Amplification Assay for Rapid Detection of Fusarium oxysporum f. Sp. lycopersici. J. Plant Pathol. Microbiol. 04 (2015).
- 43. Niessen, L. & Vogel, R. F. Detection of Fusarium graminearum DNA using a loop-mediated isothermal amplification (LAMP) assay. *Int. J. Food Microbiol.* **140**, 183–91 (2010).
- 44. Yan, M. et al. Direct detection of various pathogens by loop-mediated isothermal amplification assays on bacterial culture and bacterial colony. *Microb. Pathog.* **102**, 1–7 (2017).
- 45. Sattabongkot, J., Tsuboi, T., Han, E.-T., Bantuchai, S. & Buates, S. Loop-mediated isothermal amplification assay for rapid diagnosis of malaria infections in an area of endemicity in Thailand. J. Clin. Microbiol. 52, 1471–7 (2014).
- Geojith, G., Dhanasekaran, S., Chandran, S. P. & Kenneth, J. Efficacy of loop mediated isothermal amplification (LAMP) assay for the laboratory identification of Mycobacterium tuberculosis isolates in a resource limited setting. J. Microbiol. Methods 84, 71–73 (2011).
- 47. Curtis, K. A., Rudolph, D. L. & Owen, S. M. Rapid detection of HIV-1 by reverse-transcription, loop-mediated isothermal amplification (RT-LAMP). J. Virol. Methods 151, 264–270 (2008).
- Njiru, Z. K. et al. Loop-mediated isothermal amplification (LAMP) method for rapid detection of Trypanosoma brucei rhodesiense. PLoS Negl. Trop. Dis. 2, e147 (2008).
- Poon, L. L. M. et al. Sensitive and Inexpensive Molecular Test for Falciparum Malaria: Detecting Plasmodium falciparum DNA Directly from Heat-Treated Blood by Loop-Mediated Isothermal Amplification. Clin. Chem. 52, 303 LP–306 (2006).
- 50. Lee, P.-W. et al. Application of loop-mediated isothermal amplification for malaria diagnosis during a follow-up study in São Tomé. Malar. J. 11, 408 (2012).
- 51. Cook, J. et al. Loop-mediated isothermal amplification (LAMP) for point-of-care detection of asymptomatic low-density malaria parasite carriers in Zanzibar. Malar. J. 14, 43 (2015).

Acknowledgements

The research is funded by the Polish Ministry of Science and Higher Education under the Diamond Grant program, project No.: DI2012 024042.

Author Contributions

J.P. and M.F. contributed to the design, planning and carrying of the research, to the analysis of the results and to the writing of the manuscript.

Additional Information

Competing Interests: The authors declare no competing interests.

Publisher's note: Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Open Access This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit https://creativecommons.org/licenses/by/4.0/.

© The Author(s) 2019