

Evaluation of CEDIA and DRI Drugs of Abuse Immunoassays for Urine Screening on a Thermo Indiko Plus Analyzer

Katharina M. Köhler, Raija Hammer, Kathrin Riedy, Volker Auwärter, and Merja A. Neukamm*

Institute for Forensic Medicine, Medical Center – University of Freiburg, Faculty of Medicine, University of Freiburg, Germany

Background: For analysis of urine samples during abstinence control for driving ability assessment (medical and psychological assessment, MPA), a reliable screening method for ethyl glucuronide and drugs of abuse (cannabinoids, opiates, cocaine, amphetamines, methadone, and benzodiazepines) is needed. **Methods:** In this study CEDIA and DRI immunoassays were applied on a Thermo Fisher Scientific Indiko Plus analyzer. Precision and accuracy as well as sensitivity

and specificity at the required cut-offs for the MPA were evaluated. **Results:** The specificity was satisfactory and ranged from 91% for methamphetamine to 100% for opiates, cocaine metabolite, amphetamine, EDDP, and benzodiazepines. Moreover, sensitivity was 100% for all assays except for cannabinoids (91%). **Conclusion:** The presented method can therefore be recommended for abstinence control. *J. Clin. Lab. Anal.* 31:e22021, 2017. © 2016 Wiley Periodicals, Inc.

Key words: abstinence control; diagnostic tests; ethyl glucuronide; forensic toxicology

INTRODUCTION

Abstinence control for drugs of abuse and ethanol is often demanded for driving license regranting during medical and psychological assessment (MPA). The MPA guidelines require sensitive polytoxicological analyses of urine samples including the parameters and the quantitative MPA cut-offs given in Table 1 (1).

At the MPA cut-offs, the sensitivity of the applied screening method should be close to 100%, and preferably, positive samples with concentrations slightly below the MPA cut-off should be detected as well. According to the MPA guidelines these samples have to be reported as positives because the consumption of the drug is proven. In addition, false-positive samples should be avoided for economic reasons, as samples with a positive screening result have to be reanalyzed by chromatographic-mass spectrometric methods. Typically, immunoassays are used for screening purposes. Thus, at the defined immunoassay cut-offs, a reasonable compromise has to be achieved regarding sensitivity and specificity. Immunoassays applied for forensic analysis in the context of abstinence control have to be validated by measuring authentic samples to show that the method is capable of detecting the target analytes in at least 9 of 10 samples with concentrations in the range of the MPA cut-off (2).

Most of the urine samples collected during abstinence control for MPA are negative, so an efficient analysis procedure usually consists of a fast and easy screening method followed by confirmation of positive screening results using a quantitative method. Musshoff et al. (3) evaluated CEDIA and DRI drugs of abuse immunoassays for abstinence control in the context of MPA. The tests were carried out on an Olympus AU 400. ELISA manual microtiter tests were also evaluated for this purpose (4, 5).

The aim of this study was to evaluate the performance of CEDIA and DRI drugs of abuse and ethyl glucuronide tests on the Thermo Scientific Indiko Plus Analyzer, which is a fully automated bench-top analyzer for clinical chemistry analytics for the purpose of abstinence control according to the MPA and German society of Toxicological and Forensic Chemistry (GTFCh) guidelines. The main

*Correspondence to: Merja A. Neukamm, Institute for Forensic Medicine, Forensic Toxicology, Medical Center – University of Freiburg, Albertstr. 9, 79104 Freiburg, Germany. E-mail: merja.neukamm@uniklinik-freiburg.de

Received 2 February 2016; Accepted 7 June 2016

DOI 10.1002/jcla.22021

Published online in Wiley Online Library (wileyonlinelibrary.com).

TABLE 1. Parameters and MPA Cut-Offs in Urine for Abstinence Control (1)

Substance class Target analyte	Urine [ng/ml]
Cannabinoids THC-COOH	10 (after hydrolysis)
Opiates Morphine (codeine, dihydrocodeine)	25 (after hydrolysis)
Cocaine (metabolite) Benzoylecgonine	30
Amphetamines Amphetamine, methamphetamine, MDMA, MDEA	50
Methadone EDDP	50
Benzodiazepines Nordiazepam, oxazepam, hydroxy-alprazolam, hydroxy-bromazepam, 7-amino-flunitrazepam, lorazepam	50
Ethyl glucuronide	100

requirement is sufficient sensitivity of the assays to reliably detect the low concentrations required in the context of the MPA. In addition, the precision of the assays was assessed.

MATERIALS AND METHODS

Samples

Drug and ethyl glucuronide (EtG)-positive urine samples obtained in our laboratory for routine drug-testing were selected. Only positive samples with concentrations up to twice the MPA cut-off were used (for opiates up to triple the MPA cut-off). If necessary, the samples were diluted with blank urine samples. For some analytes spiked urine samples were used additionally to increase the number of positive samples (analytes and respective numbers of spiked samples: methamphetamine: 15; MDMA: 15; benzoylecgonine: 9; morphine: 10). For benzodiazepines, positive samples of each of the target analytes according to the MPA parameters were included. As expected, all of the oxazepam-positive samples included also nordiazepam and temazepam. Two of the lorazepam-positive samples included also less than 2 ng/ml of oxazepam. One of the two alprazolam-positive samples included also 17 ng/ml nordiazepam and 28 ng/ml oxazepam. One bromazepam positive sample was used. For opiates, all authentic samples contained mainly morphine next to low concentrations of normorphine and codeine. The blank urine samples were tested negative with confirmatory methods. Sample numbers and concentration ranges are given in Table 2.

TABLE 2. Overview of Authentic Samples Used for Evaluation of Cut-Offs, Sensitivity, and Specificity

Immunoassay	Number of samples		Concentration range [ng/ml]
	Negative	Positive	
Cannabinoids	28	28	2.0–19
Opiates (morphine)	22	26	1.0–75
Cocaine (metabolite)	27	42	2.0–40
Amphetamine	67	32	2.0–92
Methamphetamine	22	25	2.0–55
Ecstasy	23	33	2.0–61
EDDP	23	16	1.0–85
Benzodiazepines	31	48	5.0–100
Nordiazepam	(31)	17	6.5–26
Oxazepam	(31)	33	1.7–100
Hydroxy-alprazolam	(31)	2	18–35
7-Aminoflunitrazepam	(31)	4	30–55
Lorazepam	(31)	10	7.9–66
Ethyl glucuronide	34	27	80–190

Positive samples show concentrations above the limit of detection of the confirmative method. The concentration ranges of respective target analytes in positive samples are given.

Chemicals and Reagents

Analytical standards were purchased from LGC Standards GmbH (Wesel, Germany). All solvents and chemicals were at least of analytical or HPLC grade. Reagents, calibration, and control materials for the immunoassay testing were obtained from Thermo Fisher Scientific (Passau, Germany). The following CEDIA tests were used: THC (cannabinoids), Opiate, Cocaine, Amphetamine OFT, Methamphetamine OFT, EDDP (methadone metabolite), and Benzodiazepine. The Benzodiazepine CEDIA included a beta-glucuronidase protocol. The following DRI tests were used: XTC (ecstasy) and EtG (ethyl glucuronide). Beta-glucuronidase/arylsulfatase (*E. coli*) was obtained from Roche (Mannheim, Germany).

Immunoassay Testing

Four calibration points were applied for the CEDIA tests and three calibration points were used for the DRI tests. For all tests, a calibration point of zero (CEDIA-Negative Cal or DRI-Negative Cal, respectively) was applied. The calibration samples were prepared by dilution of the following standard calibration materials with CEDIA- or DRI-Negative Cal (matrix: urine), if necessary: cannabinoids: THC 25 ng/ml Cal, THC 50 ng/ml Cal; opiates, cocaine (metabolite), methamphetamine, EDDP, and benzodiazepines: Multidrug Optional Cut-Off Cal; amphetamine: Amphetamine 200 ng/ml Serum Cal (matrix: serum); ecstasy: Ecstasy 250 ng/ml Cal; and ethyl glucuronide: DRI EtG CE calibrator. The following concentrations (ng/ml) were used for calibration:

cannabinoids: 12.5, 25, and 50; opiates: 37.5, 100, and 300; cocaine: 37.5, 100, and 300; amphetamine: 25, 67, and 200; methamphetamine: 37.5, 100, and 300; EDDP: 25, 50, and 100; benzodiazepines: 50, 100, and 200; ecstasy: 50 and 250; and ethyl glucuronide: 100 and 500. The calibration matrix was urine for all parameters except for amphetamine, where a mixture of serum and urine was applied. All calibration samples were analyzed in duplicate. A point-to-point calibration model was applied. Calibration was performed on every working day.

The analyses were performed using standard settings as recommended by the manufacturer. Shortly, the volume of reagents 1 and 2 was 55 µl, the first incubation time was 300 sec, and the second incubation time was 240 sec for CEDIA tests and 90 sec for DRI tests. The main detection wavelength was 575 nm for CEDIA tests and 340 nm for DRI tests. The following sample volumes were used: THC 9 µl, Opiate 6 µl, Cocaine 9 µl, Amphetamine OFT 6 µl, Methamphetamine OFT 8 µl, EDDP 5 µl, Benzodiazepine 6 µl, XTC 22 µl, and EtG 23 µl.

Quality controls (QCs) in the low, the middle, and the high concentration range were used. The QC samples were prepared by dilution of standard control or calibration material with CEDIA- or DRI-Negative Cal (matrix: urine), if necessary. Different lots were used for dilution of QC and calibration samples. Concentrations and dilutions are given in Table 3.

The matrix of the standard control material was urine for all parameters except for amphetamine, where a serum standard was applied and diluted with Negative Cal (matrix: urine), to reach the desired concentrations of 20, 25, and 50 ng/ml. One of the QC concentrations per parameter was in the range of the MPA cut-off. A negative urine control (material: DAU Neg Ctrl) was applied after each calibration. The intra-day precision and accuracy were calculated for all three QC levels with 20 samples. Inter-day precision and accuracy were calculated for the low and the high QCs with 7–13 samples. For calculating the accuracy, the manufacturer's stated concentrations and respective dilutions were applied as target concentrations. For amphetamine and methamphetamine, OFT assays, which were originally designed for the analysis of oral fluid samples, were applied, as their sensitivity is higher than the sensitivity of the respective standard clinical urine assays.

Confirmation Analysis

All urine samples were confirmed with standard chromatographic-mass spectrometric methods, which were fully validated according to the guidelines of the

TABLE 3. Quality Control Target Concentrations and Respective Standard Material Dilution

Immunoassay and QC	Target concentration	Material and dilution	
<i>Cannabinoids</i>			
Low	9.4	THC 25 Ctrl Low	1 + 1
Med	12.5	THC 25 Cal	1 + 1
High	16	THC 25 Ctrl High	1 + 1
<i>Opiates</i>			
Low	28	MD Optional Ctrl Low	1 + 7
Med	37.5	MD Optional CutOff Cal	1 + 7
High	47	MD Optional Ctrl High	1 + 7
<i>Cocaine (metabolite)</i>			
Low	28	MD Optional Ctrl Low	1 + 7
Med	37.5	MD Optional CutOff Cal	1 + 7
High	47	MD Optional Ctrl High	1 + 7
<i>Amphetamine</i>			
Low	20	Amphetamine 200 Serum Cal	1 + 9
Med	25	Amphetamine 200 Serum Cal	1 + 7
High	50	Amphetamine 200 Serum Cal	1 + 3
<i>Methamphetamine</i>			
Low	28	MD Optional Ctrl Low	1 + 7
Med	37.5	MD Optional CutOff Cal	1 + 7
High	47	MD Optional Ctrl High	1 + 7
<i>Ecstasy</i>			
Low	37.5	MD Select Ctrl Low	1 + 9
Med	50	Ecstasy 250 Cal	1 + 4
High	75	MD Select Ctrl High	1 + 7
<i>EDDP</i>			
Low	37.5	MD Optional Ctrl Low	1 + 1
Med	50	MD Optional CutOff Cal	1 + 3
High	75	MD Optional Ctrl High	1 + 1
<i>Benzodiazepines</i>			
Low	38	MD Optional Ctrl Low	1 + 3
Med	50	MD Optional CutOff Cal	1 + 3
High	63	MD Optional Ctrl High	1 + 3
<i>Ethyl glucuronide</i>			
Low	75	EtG 375 Ctrl	1 + 4
Med	100	EtG 100 Cal	Undiluted
High	375	EtG 375 Ctrl	Undiluted

German society of Toxicological and Forensic Chemistry (GTFCh; 2) and accredited according to DIN EN ISO 17025. For THC-carboxylic acid (THC-COOH), 1 ml urine per sample was spiked with deuterated THC-COOH and hydrolyzed with 10 M sodium hydroxide (NaOH) for 15 min at 60°C. The sample was extracted via automated solid phase extraction (SPE) on Chromabond C18 cartridges, derivatized with N-methyl-N-trimethylsilyltrifluoroacetamide (MSTFA), and analyzed with gas chromatography-mass spectrometry (GC-MS) (GC: 6890N GC Agilent, Waldbronn, Germany, GC-column: J&W HP5-MS, Agilent, Waldbronn, Germany, MS: 5973N MSD, Agilent, Waldbronn, Germany). For ethyl glucuronide, 100 µl urine per sample was spiked with deuterated

ethyl glucuronide, diluted with methanol, and analyzed with liquid chromatography-tandem mass spectrometry (LC-MS/MS) (LC: LC20AD, Shimadzu, Duisburg, Germany, LC-column: Hypercarb, Thermo Scientific, Dreieich, Germany, MS: QTrap 4000, ABSciex, Toronto, Kanada). For all other drugs, sulfates or glucuronides in the samples were cleaved with beta-glucuronidase/arylsulfatase for 2 hr at 45°C. For benzodiazepines, 1 ml urine per sample was spiked with a mixture of deuterated analytes, extracted with 1-chlorobutane (liquid-liquid extraction) at pH 9, and analyzed with LC-MS/MS (LC: G1312, Agilent, Waldbronn, Germany, LC-column: Polar RP, Phenomenex, Aschaffenburg, Germany, MS: QTrap2000, ABSciex, Toronto, Canada).

For the remaining analytes (basic drugs), 1 ml urine per sample was spiked with a mixture of deuterated analytes, extracted via automated SPE on Chromabond Drug cartridges, and analyzed with LC-MS/MS (LC: G1312, Agilent, Waldbronn, Germany, LC-column: Luna PFP, Phenomenex, Aschaffenburg, Germany, MS: QTrap2000, ABSciex, Toronto, Canada). The limits of detection (LOD) and limits of quantification (LOQ) are given in Table 4.

Data Analysis

Receiver operating characteristic (ROC) plots (6) were generated by SPSS software (version 22, IBM, Armonk, NY). Sensitivity and specificity were calculated from the ROCs with Microsoft Excel 2010 (Redmont, WA).

Sensitivity and Specificity

The sensitivity of an assay was defined as the ability of the assay to correctly identify samples with analyte concentrations above the MPA cut-off. The sensitivity was calculated as the number of immunoassay-positive samples above the respective assay cut-off minus the

number of false-positive samples then divided by the number of samples with analyte concentrations above the MPA cut-off. A high sensitivity therefore indicates a low rate of false-negative samples. For calculating the sensitivity, the result of the mass spectrometric analysis was regarded as positive, if the measured concentration reached or exceeded the required MPA cut-off, thus samples with a negative immunoassay result and a concentration between the MPA cut-off and the LOD are not counted as false negatives. According to the MPA guidelines, drug-positive samples with concentrations below the MPA cut-off have to be reported as well because the consumption of the respective drug is proven. This has to be taken into account when evaluating the sensitivity of an assay. Samples with concentrations below the MPA cut-off, which are identified as positive by the respective immunoassay, are therefore regarded as true-positive samples.

The specificity of an assay was defined as the ability of the assay to correctly identify samples with analyte concentrations below the LOD of the confirmation method. The specificity was calculated as the number of immunoassay-negative samples below the respective assay cut-off minus the number of false-negative samples then divided by the number of samples with analyte concentrations below the LOD of the confirmation method. A high specificity therefore indicates a low rate of false-positive samples. For calculating the specificity, the result of the mass spectrometric analysis was regarded as positive, if the measured concentration exceeded the LOD of the method. Thus, samples with a positive immunoassay result and a concentration between the MPA cut-off and the LOD are not counted as false positives but true positives and, according to the MPA guidelines, they have to be reported as positive.

RESULTS AND DISCUSSION

Precision and Accuracy

Precision, accuracy, and target concentrations of the control samples are given in Table 4. The relative standard deviation (RSD %) was 11% for the low control of EDDP and below 10% for all other concentrations and parameters. Thus, the precision is high in the low and the high concentration range as well (Table 5).

Methamphetamine showed the highest intra-day and inter-day bias of all parameters (25–32 % and 17–23%, respectively), the bias for cocaine, ecstasy, and benzodiazepines was in the $\pm 25\%$ range. The bias for cannabinoids, opiates, amphetamine, EDDP, and ethyl glucuronide was in the $\pm 15\%$ range.

TABLE 4. Limits of Quantification and Limits of Detection in Urine of the Confirmatory Methods

Target analyte	Limit of detection [ng/ml]	Limit of quantification [ng/ml]
THC-COOH	2.0	5.0
Morphine (opiates)	1.0	2.5
Benzoylcegonine	2.0	5.0
Amphetamine	2.0	2.0
Methamphetamine	2.0	2.0
MDMA, MDA (ecstasy)	2.0	2.0
EDDP	1.0	2.5
Benzodiazepines	0.5–5.0	0.5–5.0
Ethyl glucuronide	75	75

TABLE 5. Precision and Accuracy of Quality Controls

Immunoassay	MPA Cut-off	Target concentration [ng/ml]			Precision (RSD %)					Accuracy (bias, %)				
					Intra-day			Inter-day		Intra-day			Inter-day	
		Low	Middle	High	Low	Middle	High	Low	High	Low	Middle	High	Low	High
Cannabinoids	10	9.4	12.5	16	2.4	1.0	0.9	5.2	5.1	-1	13	1	-12	0
Opiates	25	28	37.5	47	2.0	1.4	1.4	6.0	3.5	2	4	7	1	7
Cocaine (metabolite)	30	28	37.5	47	1.1	2.9	1.5	3.1	3.2	18	17	25	15	21
Amphetamine	50	20	25	50	3.8	6.5	2.5	9.9	6.6	13	10	11	10	8
Methamphetamine	50	28	37.5	47	4.4	6.4	4.3	5.3	4.5	25	26	32	17	23
Ecstasy	50	37.5	50	75	7.6	1.5	0.7	4.4	3.3	-19	3	2	-21	-2
EDDP	50	37.5	50	75	6.1	4.1	2.8	11	3.6	3	-8	7	10	11
Benzodiazepines	50	38	50	63	2.4	2.3	3.5	4.7	4.5*	3	0	1	16	13*
Ethyl glucuronide	100	75	100	375	4.8	4.5	0.9	9.7	2.0	-7	3	-6	-1	-6

*N = 4.

Sensitivity and Specificity

For evaluating the sensitivity and specificity, only samples with concentrations up to twice (opiates: triple) the MPA cut-off were applied (see Table 2) because these are the crucial concentrations for the evaluation of the performance of an immunoassay used for abstinence screening.

Sensitivity and specificity diagrams are given in Figure 1. Possible assay cut-offs and respective sensitivities and specificities are given in Table 6. The assays for opiates, cocaine (metabolite), EDDP, and benzodiazepines showed 100% specificity and sensitivity, simultaneously. At these assays the cut-off could also be varied in a certain range without generating




Fig. 1. Specificity and sensitivity diagrams calculated from ROC data.

TABLE 6. Assay Cut-Offs and Respective Sensitivity and Specificity

Immunoassay	Assay cut-off		Specificity [%]
	[IA-units]	Sensitivity [%]	
Cannabinoids	4.0	91	96
Opiates	8	100	100
Cocaine (metabolite)	10	100	100
Amphetamine	8	100	97
Methamphetamine	36	100	91
Ecstasy	33	100	96
EDDP	9	100	100
Benzodiazepines	13	100	100
Ethyl glucuronide	90	100	97

false negatives, e.g., for opiates from 8 to 16, for cocaine (metabolite) from 10 to 32, for EDDP from 9 to 46, and for benzodiazepines from 13 to 30.

For cannabinoids, amphetamine, methamphetamine, ecstasy, and ethyl glucuronide, the specificity declines with increasing assay cut-off before the sensitivity reaches 100%. For those assays several false-positive samples have to be accepted at the low necessary cut-off required to fulfill the MPA guidelines. The lack of specificity of the methamphetamine and ecstasy assays at low cut-offs results from cross-reactivity with endogenous phenethylamines which can be present in urine samples. Only the cannabinoids assay showed a sensitivity below 100% as there was one sample with a confirmed concentration of 13 ng/ml THC-COOH and a cannabinoid immunoassay value of 0.4 (assay cut-off: 4.0). This sample can be regarded as an outlier, resulting from a potential adulteration of the sample with, e.g., an oxidant during collection (leading to a reduction in enzyme activity). The resulting sensitivity for cannabinoids of 91% is acceptable as the GTFCh guidelines (2) demand a rate of detection of at least 9 of 10 positive samples (sensitivity $\geq 90\%$) for the validation of immunoassays for MPA cut-offs.

CONCLUSION

The evaluation of CEDIA and DRI drugs of abuse and ethyl glucuronide immunoassays for abstinence

control on the Thermo Indiko Plus Analyzer showed good results. The method showed sufficient precision and accuracy for urine screening. The sensitivity and specificity are mostly 100% at the assay cut-offs required to meet the MPA criteria. For the majority of analytes false-positive immunoassay results are rare, thus only few unnecessary confirmation analyses have to be carried out. More importantly, the only false-negative result occurred at the cannabinoids assay. So the presence of a drug or metabolite in urine in the concentration range of the MPA cut-off would be revealed in almost every case. The presented method can therefore be recommended for urine screening for abstinence control on drugs of abuse or ethanol.

ACKNOWLEDGEMENTS

The authors would like to thank Thermo Fisher Scientific for providing reagents and equipment as well as technical and scientific support. The Institute for Clinical Chemistry and Laboratory Medicine, Medical Center – University of Freiburg is acknowledged for their collaboration.

REFERENCES

1. Schubert W, Dittmann V, Brenner-Hartmann J. 2013. Beurteilungskriterien: Urteilsbildung in der Fahreignungsbeurteilung, 3rd edn. Bonn: Kirschbaum Verlag.
2. Peters FT, Hartung M, Herbold M, Schmitt G, Daldrup T, Musshoff F. Anhang B zu den Richtlinien der GTFCh zur Qualitätssicherung bei forensisch-toxikologischen Untersuchungen; Anforderungen an die Validierung von Analysemethoden. *Toxicchem Krimtech* 2009;76:185–199.
3. Musshoff F, Wolters T, Lott S, Ippisch J, Gradl S, Madea B. Optimization and validation of CEDIA drugs of abuse immunoassay tests in serum and urine on an Olympus AU 400. *Drug Test Anal* 2013;5:366–371.
4. Kirschbaum KM, Musshoff F, Wilbert A, Röhrich J, Madea B. Direct ELISA kits as a sensitive and selective screening method for abstinence control in urine. *Forensic Sci Int* 2011;207:66–69.
5. Agius R, Nadulski T, Moore C. Validation of LUCIO((R))-Direct-ELISA kits for the detection of drugs of abuse in urine: Application to the new German driving licence re-granting guidelines. *Forensic Sci Int* 2012;215:38–45.
6. Zweig MH, Campbell G. Receiver-operating characteristic (ROC) plots: A fundamental evaluation tool in clinical medicine. *Clin Chem* 1993;39:561–577.