#### **METHODOLOGY ARTICLE**

**Open Access** 

# A T-DNA mutant screen that combines high-throughput phenotyping with the efficient identification of mutated genes by targeted genome sequencing



Ulrike Frank<sup>1</sup>, Susanne Kublik<sup>2</sup>, Dörte Mayer<sup>1</sup>, Marion Engel<sup>3</sup>, Michael Schloter<sup>2</sup>, Jörg Durner<sup>1,4</sup> and Frank Gaupels<sup>1\*</sup>

#### **Abstract**


**Background:** Nitrogen dioxide (NO<sub>2</sub>) triggers hypersensitive response (HR)-like cell death in *Arabidopsis thaliana*. A high-throughput mutant screen was established to identify genes involved in this type of programmed cell death.

**Results:** Altogether 14,282 lines of SALK T-DNA insertion mutants were screened. Growing 1000 pooled mutant lines per tray and simultaneous NO<sub>2</sub> fumigation of 4 trays in parallel facilitated high-throughput screening. Candidate mutants were selected based on visible symptoms. Sensitive mutants showed lesions already after fumigation for 1 h with 10 ppm (ppm) NO<sub>2</sub> whereas tolerant mutants were hardly damaged even after treatment with 30 ppm NO<sub>2</sub>. Identification of T-DNA insertion sites by adapter ligation-mediated PCR turned out to be successful but rather time consuming. Therefore, next generation sequencing after T-DNA-specific target enrichment was tested as an alternative screening method. The targeted genome sequencing was highly efficient due to (1.) combination of the pooled DNA from 124 candidate mutants in only two libraries, (2.) successful target enrichment using T-DNA border-specific 70mer probes, and (3.) stringent filtering of the sequencing reads. Seventy mutated genes were identified by at least 3 sequencing reads. Ten corresponding mutants were re-screened of which 8 mutants exhibited NO<sub>2</sub>-sensitivity or -tolerance confirming that the screen yielded reliable results. Identified candidate genes had published functions in HR, pathogen resistance, and stomata regulation.

**Conclusions:** The presented  $NO_2$  dead-or-alive screen combined with next-generation sequencing after T-DNA-specific target enrichment was highly efficient. Two researchers finished the screen within 3 months. Moreover, the target enrichment approach was cost-saving because of the limited number of DNA libraries and sequencing runs required. The experimental design can be easily adapted to other screening approaches e.g. involving high-throughput treatments with abiotic stressors or phytohormones.

**Keywords:** *Arabidopsis thaliana*, T-DNA insertion, Mutant screen, High-throughput phenotyping, Target enrichment, Next generation sequencing, Adapter ligation-mediated PCR

<sup>&</sup>lt;sup>1</sup>Institute of Biochemical Plant Pathology, Helmholtz Zentrum München, Neuherberg, Germany





<sup>\*</sup> Correspondence: frank.gaupels@gmx.de

Frank et al. BMC Plant Biology (2019) 19:539 Page 2 of 12

#### **Background**

Forward genetics by mutant screens is a widely used approach to experimentally assign biological functions to genes or characterize physiological processes through the identification of involved genes [1]. Thousands of mutant lines were generated in Arabidopsis thaliana, Oryza sativa, Solanum lycopersicum, Medicago truncatula, and other plant species by insertional mutagenesis of genes using retrotransposons or Agrobacterium tumefaciens-derived T-DNA [2-5]. For Arabidopsis, the SALK collection represents the largest source of T-DNA insertion mutants covering almost the complete genome [1]. The fact that the T-DNA sequence is known facilitates convenient mapping of the insertion site. For instance, several PCR-based methods such as the adapter ligation-mediated PCR rely on insert-specific primers to selectively amplify the T-DNA/genomic DNA junction before identification of the mutated gene by sequencing [6]. Alternatively, insert specific biotinylated 70mer probes can be employed for the isolation of T-DNAcontaining DNA fragments by streptavidin beads [7]. Combining this target enrichment with next generation sequencing (NGS) allowed the simultaneous identification of multiple insertion sites in a complex pool of DNA from as much as 64 different mutants [7]. Both, adapter ligation-mediated PCR and target enrichment followed by NGS (i.e. targeted genome sequencing) are suitable for high-throughput mappings of mutation sites [6, 7].

Nitrogen dioxide (NO<sub>2</sub>) is a notorious toxic air pollutant but also an upcoming endogenous signal in plant cells where it arises from the oxidation of NO and nitrite or decomposition of peroxynitrite [8, 9]. Fumigation experiments revealed that NO2 has beneficial as well as detrimental effects on plant performance dependent on the applied gas concentration. Short-term treatments of Arabidopsis for 1 h with 10 ppm NO<sub>2</sub> did not cause visible leaf damage but triggered transient salicylic acid (SA) accumulation and basal pathogen resistance against Botrytis cinerea and Pseudomonas syringae [10]. Treatments of Arabidopsis seedlings for 8 h with 10 ppm or 1 h with 30 ppm NO<sub>2</sub> caused apparent lesion formation [11, 12]. NO<sub>2</sub>-induced cell death was reminiscent of the hypersensitive response (HR) of resistant plants towards avirulent pathogens, which ultimately cumulates in programmed cell death (PCD) [13-15]. NO<sub>2</sub>-induced cell death - like HR-PCD - was found to be dependent on simultaneous signaling by H<sub>2</sub>O<sub>2</sub> and NO, and was enhanced in SA deficient mutants [11, 12, 16, 17]. Moreover, both cell death events were accompanied by the accumulation of fluorescent phenolic compounds (Additional file 1: Figure S1) and oxylipins including jasmonates [11]. Collectively, these findings suggest that NO<sub>2</sub> exposure leads to the dose-dependent induction of basal pathogen resistance and HR-like cell death. How  $\mathrm{NO}_2$  exerts these different effects in plants remains ambiguous.

Therefore, a high-throughput genetic screen was developed to identify new genes involved in NO2 sensitivity and -tolerance. Altogether 14,282 individual lines of SALK T-DNA insertion mutants were screened for their NO<sub>2</sub>-induced leaf symptoms. Sensitive mutants displayed lesions already after fumigation with 10 ppm NO<sub>2</sub> whereas the other tested plants were unaffected. Tolerant mutants survived even exposure to 30 ppm NO<sub>2</sub> without severe damage. Subsequently, T-DNA-mutated genes were identified by adapter ligation-mediated PCR or genome sequencing after T-DNA-specific target enrichment. The latter technique proved to be particularly efficient facilitating the identification of 162 genes with putative functions in NO<sub>2</sub> sensitivity or tolerance. An initial re-screen of corresponding mutants confirmed that 8 of 10 tested candidate mutants were altered in NO<sub>2</sub>-induced leaf damage as compared to wild-type (WT) plants.

Overall, the current paper describes an experimental set-up that facilitates the efficient genome-wide screening, identification and investigation of T-DNA insertion mutants. The experimental design can be adapted to other high-throughput treatments e.g. with abiotic stressors or phytohormones. Although the screen was optimized for Arabidopsis, it would also be possible to screen insertion mutants in other plant species.

#### Challenges/tasks

Whole genome screen → many mutant lines must be handled

NGS → large datasets must be handled → advanced strategy for data handling/filtering required

Screen should be cost-, time, and labor-saving

Advantages and disadvantages of published screens and mutant collections

#### **Results**

# High-throughput screening for NO<sub>2</sub>-sensitive and -tolerant T-DNA insertion mutants

Altogether 14,282 SALK T-DNA insertion mutants ("confirmed homozygous" mutant sets) were screened for their sensitivity or tolerance towards  $\mathrm{NO}_2$ . In order to limit the investment of hands-on time and climate chamber space the different mutant lines were not grown separately but were pooled (Fig. 1). Each  $\sim 1000$  mutant lines (2–4 seeds per line) were combined resulting in only 14 trays of plants (Fig. 1, Fig. 2a). Fourteenday-old mutant seedlings were short-term fumigated for 1 h with 10 ppm  $\mathrm{NO}_2$  and (2 days later) 30 ppm  $\mathrm{NO}_2$ . The exposure chamber allowed  $\mathrm{NO}_2$  treatment of 4 trays in parallel, i.e. fumigation of all 14 trays with one

Frank et al. BMC Plant Biology (2019) 19:539 Page 3 of 12

#### 14,282 T-DNA insertion mutants



14 seed stocks of ~1,000 pooled mutant lines each, 2-4 seeds/line → ~3,000 seeds/stock



14 trays with mutant plants



Sequential NO<sub>2</sub> fumigations with 10 and 30 ppm NO<sub>2</sub>



Sampling of NO<sub>2</sub>-sensitive and -tolerant mutants



Identification of mutants by adapter ligation-mediated PCR or genome sequencing after target enrichment

**Fig. 1** Strategy of the NO<sub>2</sub> screen of Arabidopsis T-DNA insertion mutants

concentration of  $\mathrm{NO}_2$  took only 4 h. Sensitive plants developed lesions already after exposure to 10 ppm  $\mathrm{NO}_2$  whereas the other plants were not affected (Fig. 2b). Subsequent treatment with 30 ppm  $\mathrm{NO}_2$  caused severe leaf damage in almost all plants (Fig. 2c) whereas tolerant mutants displayed only weak symptoms (Fig. 2d). Altogether, 124 SALK T-DNA insertion lines showing distinct  $\mathrm{NO}_2$  phenotypes were sampled and DNA was extracted.

#### Identification of mutated genes by adapter ligationmediated PCR

Adapter ligation-mediated PCR is the standard method employed by the SALK institute to define T-DNA insertion sites [6]. It is based on digestion of the DNA by the restriction enzyme *Ase*1, Ase adapter ligation, and selective amplification of T-DNA-containing fragments using the T-DNA left border-specific primer LBb1 and the adapter-specific primer AP2 before sequencing. Hence, adapter ligation-mediated PCR promotes the targeted amplification and sequencing of T-DNA/genomic DNA

junctions. Here, this approach was tested with three NO<sub>2</sub>-tolerant candidate mutants detected in the screen.


Selective amplification of the T-DNA-containing fragments and separation in an agarose gel resulted in a clear band for mutant 1 but weaker bands for mutants 2 and 3 (Fig. 3a). The low DNA concentrations in the excised bands prevented successful sequencing. Therefore, multiple PCR reactions were separated in an agarose gel, and the major band for mutant 1 and mutant 2 was cut (Fig. 3b). Furthermore, the two strongest bands for mutant 3 were analyzed to check whether they represent two independent T-DNA insertion sites (Fig. 3b). Sequencing of the eluted DNA revealed that in all cases the investigated DNA fragments contained the T-DNA left border (Fig. 3c). The DNA from mutant 1 additionally included a 726 bases-long sequence with 100% homology to the Arabidopsis gene AT5G55620. For mutant 2 the 592 bases adjacent to the T-DNA border shared significant homology with the gene AT1G13860/QUASIMODO-LIKE1/QUL1. Both analyzed mutant 3 bands contained a short DNA stretch of only 33 bases next to the T-DNA border that was homologous to the gene AT2G16630/FUSED OUTER CUTICU-LAR LEDGE1/FOCL1. The bases adjacent to this sequence did not result in a BLAST (Basic Local Alignment Search Tool) hit. Hence, both analyzed bands for mutant 3 represented the same T-DNA insertion. In sum, adapter ligation-mediated PCR successfully identified T-DNA insertion sites in all three investigated mutants.

# Identification of mutated genes by genome sequencing after T-DNA target enrichment

In a second approach (Fig. 4a) DNA aliquots of the 124 sampled mutant seedlings were combined in two pools each consisting of 1 µg DNA from 59 or 65 mutant seedlings. The preparation of the two libraries involved digestion of the DNA, adapter ligation, and amplification. At this stage T-DNA-containing DNA fragments were hidden within the complex total DNA pools. Therefore, 6 T-DNA border-specific biotinylated 70mer probes served to semi-purify the T-DNA-containing DNA fragments (Fig. 4a and b) [7]. Successful purification was confirmed by quantitative polymerase chain reaction (qPCR) with T-DNA-specific primers before identification of the T-DNA insertion sites by next generation sequencing. It is important to note that the target enrichment approach required processing of only two DNA pools and only two runs of next generation sequencing.

High proportions of long reads (400–700 bases) and Phred scores around 30 (i.e. 99.9% probability of a correct base call) indicated a good sequencing quality for both libraries. Further analysis of the sequencing results

Frank et al. BMC Plant Biology (2019) 19:539 Page 4 of 12



**Fig. 2** Detection of NO<sub>2</sub>-sensitive and -tolerant mutants. **a.** Approximately 3000 mutant plants per tray were grown for 2 weeks. Plants were fumigated for 1 h with 10 ppm NO<sub>2</sub> and 2 days later with 30 ppm NO<sub>2</sub>. **b** Example of a sensitive mutant showing lesions after fumigation with 10 ppm NO<sub>2</sub> while neighbouring plants were unaffected. Pictures were taken at 48 h after fumigation. **c** Most plants were severely damaged after exposure to 30 ppm NO<sub>2</sub>. **d** However, tolerant mutants exhibited only weak symptoms. Pictures were taken at 48 h – 96 h after fumigation. Note that the highlighted (red arrows) mutants exhibited altered NO<sub>2</sub>-induced symptoms compared to all neighboring plants

revealed that in both libraries approx. 75% of the reads contained T-DNA proving again the successful T-DNA enrichment (Table 1). By comparison, in another study based on target enrichment 40% of the reads mapped to the T-DNA [7]. More stringent filtering criteria were applied to extract high-quality reads as a pre-requisite for the correct determination of T-DNA insertion sites. The raw datasets were filtered by BLASTN searches against the T-DNA containing vector pBIN-pROK2 with highscoring segment pair (HSP) length set to be between 30 and 343 bp and an Expect (E) value threshold of below 5.72\*E<sup>-11</sup>. This procedure removed > 97% of the total sequencing reads (Table 1). Additional trimming of the data for read length (>40 base pairs (bp)) and quality (low quality limit 0.05, ambiguous bases ≤2) as well as disposal of adapter-only- and vector-only sequences reduced the number of remaining reads to 264 for library I and 374 for library II corresponding to 0.1 and 0.2% of the original sequencing reads, respectively. 255 (97%) and 367 (98%) of these reads mapped to the Arabidopsis genome suggesting that the bioinformatic processing of the datasets yielded meaningful results (Table 1).


The altogether 622 reads mapped to 162 Arabidopsis genes. All 3 T-DNA insertion sites defined by adapter ligation-mediated PCR were also detected by the target enrichment approach which can be interpreted as internal cross-validation of both methods. Although *AT2G16630/FOCL1* had a T-DNA insertion as shown

by both approaches, the corresponding mutants were not present in the SALK mutant collection used for the screen. A second site mutation of AT2G16630 seems unlikely because adapter ligation-mediated PCR detected only a single T-DNA insertion. More likely, the used seed stocks included an AT2G16630 mutant that was incorrectly annotated or labelled. For instance, the SALK\_042357 line has a T-DNA mutation in the position determined by the targeted genome sequencing (Additional file 3: Figure S2) but was not included in the screened mutant sets. Another remarkable finding of the target enrichment approach was that the characterization of 124 NO2-sensitive and -tolerant mutant seedlings led to the identification of 162 T-DNA-containing genes. This could be related to unnoticed secondary T-DNA insertions in the tested mutants or experimental artefacts such as sequencing errors. Seventy mutated genes were identified by ≥3 reads, 23 genes by 2 reads, and 69 genes by 1 read. The most solid candidate genes were selected for further characterization based on a read number ≥ 3 and presence of the corresponding mutant line in the screened mutant collection.

# Confirmation of the NO<sub>2</sub>-induced mutant phenotypes by a re-screen

A re-screen was initiated to confirm the altered NO<sub>2</sub>-induced cell death phenotypes of the identified candidate

Frank et al. BMC Plant Biology (2019) 19:539 Page 5 of 12



**Fig. 3** Identification of candidate genes by adapter ligation-mediated PCR. **a** DNA of the mutants 1, 2, and 3 (M1, M2, M3) was extracted and T-DNA containing fragments were selectively amplified by nested PCR. DNA of a wild-type plant (WT, Col-0) served as T-DNA negative control. For each mutant line 3 and 9 ng DNA was loaded on an agarose gel. DNA standards are indicated. **b** Multiple PCR reactions (9 ng DNA per lane) were loaded on an agarose gel to obtain sufficient amounts of DNA for sequencing. Major bands (blue arrow heads) were cut, and the eluted DNA was sequenced. **c** Results of sequencing, homology search and candidate gene alignment are shown. For mutant 3 both cut bands contained DNA with the displayed consensus sequence. Homology of mutant DNA to the pBIN-pROK2 vector (T-DNA) or an Arabidopsis gene is highlighted in yellow and red, respectively. Asterisks indicate matching bases. Only 100 bases around the T-DNA-gene junctions are shown. BLAST hit for Mutant 1: NM\_124944.4, AT5G55620, E value 0.0, 100% identities (726/726). BLAST hit for Mutant 2: CP002684.1, AT1G13860, QUL1, E value 0.0, 89% identities (527/592). BLAST hit for Mutant 3: NM\_127215.5, AT2G16630, FOCL1, E value 2e-06, 100% identities (33 of 33). The T-DNA-flanking genomic sequences used for the BLAST searches are given in Additional file 2: Table S1

mutants. Visible symptoms are not a suitable estimate of cell death because they can be influenced e.g. by altered leaf toughness or enhanced chlorosis in certain mutants. Therefore, cell death was quantified by ion leakage measurements after fumigation for 40 min with 30 ppm  $NO_2$ . In addition, the basal stomatal conductance, which is a determining factor of  $NO_2$  uptake via the stomatal pores, was assessed in untreated leaves using a porometer [11]. The 10 tested mutants were identified by  $\geq 3$  reads in the target enrichment approach (Table 2).

Additionally, the mutants *qul1* and *at5g55620* were also identified by adapter ligation-mediated PCR.

All tested mutants had a basal ion leakage similar to WT confirming that they did not develop spontaneous lesions (Additional file 4: Figure S3a). Eight of 10 investigated mutants showed different extents of HR-like cell death after NO<sub>2</sub> fumigation compared to WT plants (Fig. 5; Additional file 4: Figure S3b and c). The mutants far1, erdj3b, qul1, and cyp76c1 were sensitive but bub3.2, at1g04930, at5g55620, and cesa7 rather tolerant

Frank et al. BMC Plant Biology (2019) 19:539 Page 6 of 12




Fig. 4 The T-DNA target enrichment approach. a Overview of the experimental procedure. b Sequences of the T-DNA border-specific 70mer probes used for target enrichment. Partial sequences of the pBIN-pROC2 vector containing the T-DNA left and right borders are shown. Binding sites of the 70mer probes specific for the T-DNA left border (LB1 (yellow background color), LB2 (red frame), LB3 (blue background)) and right border (RB1 (green background), RB2 (blue frame), RB3 (orange background)) are indicated

to NO<sub>2</sub>. *qul1* was found to be tolerant in the original screen but sensitive in the re-screen. The reason for this inconsistency is unknown. Two mutants had an enhanced stomatal conductance phenotype (Fig. 5). For *cyp76c1* the increased stomatal uptake of NO<sub>2</sub> correlated with stronger symptom development whereas *ddi1* showed similar symptoms like WT. By contrast, the strong NO<sub>2</sub> tolerance of *at5g55620* and *cesa7* was linked to reductions of the stomatal conductance by 42 and 65%, respectively. Three mutants exhibited elevated and 2 mutants decreased NO<sub>2</sub>-induced ion leakage but no stomata phenotype. Ion leakage was 18–37% higher in *far1*, *erdj3b*, and *qul1* but 16 and 47% lower in *bub3.2* and *at1g04930* compared to WT (Fig. 5).

The  $NO_2$ -sensitive or -tolerant mutants described above are candidates for further characterization. Although SALK mutants with confirmed T-DNA insertion

sites were purchased for the study, homozygosity of the mutants must be confirmed by genotyping (Fig. 6). After the re-screen, genotyping will be necessary only for the most promising mutants. Mutant backcrossing or the investigation of two independent mutant lines per candidate gene prevents artefacts e.g. due to unnoticed secondary T-DNA insertions. Ultimately, it will be investigated whether knockout- and over-expressor lines of selected genes show specific phenotypes related to NO<sub>2</sub>-induced cell death, HR-PCD and pathogen resistance.

#### Discussion

The project aimed at establishing an efficient genome-wide screen of Arabidopsis mutants to identify genes involved in  $\mathrm{NO}_2$  sensitivity and tolerance. The mutant collection investigated so far comprised 14,282 SALK T-DNA insertion lines, where altogether 10,849 individual genes were mutated. Hence, the collection covered more

**Table 1** Filtering of sequencing reads improves mapping to the Arabidopsis genome

|   | Library I    |       | Library II   |       |
|---|--------------|-------|--------------|-------|
|   | No. of reads | %     | No. of reads | %     |
| Total sequencing reads  | 207,785      | 100.0 | 154,847      | 100.0 |
| 1. Filter: T-DNA containing reads   | 163,166      | 78.5  | 116,160      | 75.0  |
| 2. Filter <sup>a</sup> : HSP length, e-value  | 1536         | 0.7   | 3932         | 2.5   |
| 3. Filter <sup>b</sup> : Read length and quality, removal of adapter-only and vector-only reads | 264          | 0.1   | 374          | 0.2   |
| Mapped to the Arabidopsis genome  | 255          | 0.1   | 367          | 0.2   |

a BLASTN search against the pBIN-PROK2 sequence with the settings: High-scoring segment pair (HSP) length  $\geq$  30 bp to  $\leq$ 343 bp, E value  $\leq$ 5.72\*E $^{-11}$ ; b read length > 40 bp, low quality limit 0.05, ambiguous bases  $\leq$ 2, removal of the 5' terminal nucleotide, removal of adapter-only and vector-only sequences

Frank et al. BMC Plant Biology (2019) 19:539 Page 7 of 12


**Table 2** Description of re-screened candidate mutants

| Mutated gene | Gene name | Reads | SALK line    |
|--------------|-----------|-------|--------------|
| AT4G15090    | FAR1      | 4     | SALK_031652C |
| AT3G62600    | ERDJ3B    | 10    | SALK_055599C |
| AT1G13860    | QUL1      | 15    | SALK_048823C |
| AT2G45560    | CYP76C1   | 10    | SALK_010566C |
| AT1G49910    | BUB3.2    | 4     | SALK_151687C |
| AT1G04930    | AT1G04930 | 14    | SALK_015636C |
| AT5G55620    | AT5G55620 | 5     | SALK_018370C |
| AT5G17420    | CESA7     | 3     | SALK_029940C |
| AT1G52810    | AT1G52810 | 5     | SALK_050959C |
| AT3G13235    | DDI1      | 6     | SALK_066713C |

than one third of the complete Arabidopsis genome. High-throughput screening was facilitated by pooling of each  $\sim 1000$  mutant lines ( $\sim 3000$  seeds) resulting in only 14 trays of mutant seedlings to be handled (Fig. 1). Growing the plants on soil within a climate chamber prevented exposure of the plants to artificial conditions prevailing e.g. in petri dishes or 96-well plates. Additionally, this setup allowed the efficient treatment of many mutant lines in parallel. Other main features of the screen included short-term NO $_2$  fumigations for only 1 h, NO $_2$  treatments of 4 trays in parallel, and simple selection of candidate mutants

for their visible  $\mathrm{NO}_2$  phenotypes (Fig. 2). This way, 124 candidate mutants displaying distinct  $\mathrm{NO}_2$  phenotypes were collected. A single person could carry out the  $\mathrm{NO}_2$  screen within one-month period.

Adapter ligation-mediated PCR identified T-DNA insertion sites in all three investigated mutants (Fig. 3). However, at least in our hands and without robotics the method was suitable only for low numbers of mutants because it was rather time consuming. Moreover, costly enzymes, chemicals and sequencing were required for every single mutant analysis. By contrast, targeted genome sequencing required the processing of just 2 libraries containing the pooled DNA of all 124 candidate mutants (Fig. 4). The identification of mutated genes further involved target enrichment using T-DNA border-specific 70mer probes (Fig. 4), next generation sequencing, and stringent filtering of the sequencing reads by bioinformatic tools (Table 1). This strategy led to the identification of 70 T-DNA insertion sites by at least 3 sequencing reads. It is important to note that all 3 gene mutations detected by adapter ligation-mediated PCR were also found by targeted genome sequencing indicating that both approaches yielded similar results. Overall, the one-by-one analysis of mutants by adapter ligation-mediated PCR reduced the risk that a mutant remained unidentified without notice whereas NGS after



**Fig. 5** Re-screen of candidate mutants. Ten mutants defective in genes identified by the target enrichment approach were re-screened. The genes AT1G13860 (QUL1) and AT5G55620 were also identified by adapter ligation-mediated PCR. The tested mutants were not genotyped. Ion leakage was determined upon fumigation with 30 ppm NO<sub>2</sub> whereas stomatal conductance was measured in untreated plants. Asterisks indicate statistically significant differences ( $p \le 0.05$ ) between WT and tested mutants according to one-way ANOVA and Holm-Sidak posthoc test (ion leakage: n = 15-28; stomatal conductance: n = 12-23)

Frank et al. BMC Plant Biology (2019) 19:539 Page 8 of 12

### NO<sub>2</sub> dead-or-alive-screen



# Mutant identification by targeted genome sequencing



Re-screen, ion leakage and stomatal cond. measurements



# Genotyping and functional characterization of candidate mutants

**Fig. 6** Summary of the T-DNA mutant screen. Genotyping and functional characterization of candidate mutants represents future work not described in the current report

target enrichment had the advantage of being time and cost saving.

Ten mutants defective in identified candidate genes were re-screened. Eight of these mutants showed altered NO2-induced cell death compared to WT plants confirming that results of the screen were reproducible. The NO<sub>2</sub> sensitivity of cyp76c1 mutant correlated with an enhanced leaf uptake of NO2 due to an increased stomatal conductance whereas the reduced stomatal conductance of cesa7 and at5g55620 probably accounted for the strong NO<sub>2</sub> tolerance of these mutants (Fig. 5). The CELLULOSE SYNTHASE 7-deficient mutant cesa7 is disturbed in xylem-mediated water transport, which triggers a constitutive down-regulation of the stomatal aperture to reduce water loss due to transpiration [18]. Moreover, inhibition of the cell wall biosynthesis in cesa7 leads to enhanced abscisic acid-dependent resistance against Ralstonia solanacearum [19]. AT5G55620 is an ethylene-responsive gene that codes for an uncharacterized protein [20]. Functions of CYTOCHROME P450 76C1 (CYP76C1) are also unknown. However, the data presented here argue for a role of AT5G55620 and CYP76C1 in stomatal regulation. It will be interesting to learn whether pathogen-induced HR-PCD is altered in cyp76c1, at5g55620, and cesa7.

Other mutants showed a distinct NO<sub>2</sub>- but no stomata phenotype (Fig. 5). Hence, the respective genes code for

proteins that probably function in cell death induction or protection rather than stomata regulation. far1, erdj3b, and qul1 were sensitive but bub3.2 and at1g04930 tolerant towards NO2. The transcription factor FAR-RED IMPAIRED RESPONSE 1 (FAR1) prevents light-induced oxidative stress and SA-dependent cell death [21], and thus might also protect cells from NO<sub>2</sub>induced cell death and HR-PCD. ERDJ3B is a chaperone involved in the assembly of the immune receptor kinase EFR [22]. Its role in basal pathogen resistance could be connected to cell death protection by unknown mechanisms. The ubiquitin ligase BUB3.2 acts in cell division. BUB3.2 gene expression was down-regulated in Arabidopsis pistils after infection with the fungal pathogen Fusarium graminearum [23]. The fact that ion leakage is decreased in the knockout mutant argues for BUB3.2 being a positive regulator of NO2-induced HR-like cell death. qul1 showed 37% increased but at1g04930 47% decreased ion leakage after NO<sub>2</sub> fumigation compared to the WT. Functions of the mutated genes were not yet investigated. Therefore, the current results represent initial evidence that the genes could be involved in cell death regulation.

 $\mathrm{NO}_2$  triggers basal pathogen resistance or HR-like cell death in a dose-dependent manner [10, 11]. Accordingly, the results described here suggest that the  $\mathrm{NO}_2$  dead-oralive screen is a useful tool to identify genes related to HR (–like)-PCD, pathogen resistance, and stomata regulation.

#### **Conclusions**

The described experimental system allows for the efficient screening and identification of Arabidopsis T-DNA insertion mutants. Main features of the NO2 dead-oralive screen included (1.) the growth and treatment of 1000 pooled mutant lines per tray and (2.) candidate mutant selection based on obvious damage symptoms. Main features of the mutant identification by the target enrichment approach included (1.) pooling of DNA from ~ 60 candidate mutants per library, (2.) T-DNA-specific target enrichment before next-generation sequencing, and (3.) stringent filtering of the sequencing reads. Mutant identification by adapter ligation-mediated PCR might be advisable only if robotics is available. A re-screen revealed that 8 of ten tested mutants showed NO<sub>2</sub> phenotypes thereby confirming that the results of the screen were reliable. The described experimental system can be applied for any screening of Arabidopsis T-DNA insertion mutants based on high-throughput treatments inducing obvious phenotypes. Moreover, the screening procedure could even be adapted to insertion mutants in other plant species such as rice and tomato.

Frank et al. BMC Plant Biology (2019) 19:539 Page 9 of 12

#### **Methods**

#### Plants and growth conditions

SALK T-DNA insertion mutants were used in this study. Generation, formal identification, and submission of the mutant lines to public seed stocks was done by the SALK Institute Genomic Analysis Laboratory (SIGnAL) [24]. The four sets of confirmed homozygous lines with the NASC (Nottingham Arabidopsis Stock Centre) order numbers N27941, N27951, N27942, and N27952 comprised altogether 14,282 T-DNA mutants. All lines are available from NASC (http://arabidopsis.info). Fourteen seed pools of approx. 1000 mixed mutant lines each were prepared. To this end, 2-4 seeds per line were flipped out of every mutant stock tube by a bended needle forming a little hook. Thus, one seed pool contained approx. 3000 seeds, and in sum > 42,000 seeds were used in the screen. Fourteen plant trays (size:  $0.6 \text{ m} \times 0.4 \text{ m} \times$ 0.06 m) with inlets (tray-sized, holes in the bottom to prevent water logging) were filled with soil (4 parts Floradur propagation substrate (Floragard) mixed with 1 part quartz sand). Seeds of every pool were evenly spread in each tray and vernalized for 2 days at 4°C in the dark.

#### NO<sub>2</sub> fumigations

The mutant screen was done in walk-in size chambers housing 4 air tight fumigation chambers (www.helmholtz -muenchen.de/eus/facilities/phytotron), which allowed simultaneous fumigation of 4 plant trays in parallel. The growth conditions were set to 250 µmol m<sup>-2</sup> s<sup>-1</sup> light intensity, 14 h light-10 h dark cycle, 23 °C/18 °C (day/night), and 70%/90% (day/night) relative humidity. The plants were grown for 2 weeks in the chambers, and mutants exhibiting chlorosis or lesions were removed before the NO<sub>2</sub> treatments. NO<sub>2</sub> was generated by the reaction of 15% NO with 100%  $O_2$  in mixing vessels containing Raschig glass rings. The NO<sub>2</sub> concentrations were adjusted by regulating the NO flux rate. The concentrations of NO2 were monitored with an AC32M chemiluminescent oxides of nitrogen analyzer (Environnement S.A.). Generally, the fumigations started in the morning at approx. 2 h after onset of the light period. Plants were fumigated with 10 ppm NO<sub>2</sub> for 1 h and 2 days later with 30 ppm NO<sub>2</sub>. 10 ppm NO<sub>2</sub> had no visible effect on Col-0 (mutant background) and most mutants. However, exposure to 30 ppm NO<sub>2</sub> caused dead leaf areas or complete leaf collapse in most of the tested plants. Sensitive plants displayed symptoms already after 10 ppm NO<sub>2</sub> whereas tolerant plants were hardly affected even by 30 ppm NO<sub>2</sub>.

#### **Imaging**

UV-induced fluorescence was detected using a handheld UV lamp (Blak-Ray B-100AP; UVP). Plants were photographed with a Nikon DC300 digital camera.

#### Sampling of candidate mutants, and DNA extraction

Mutants showing NO<sub>2</sub> phenotypes were sampled at 48 h after fumigation when symptoms became clearly visible due to bleaching of the dead leaf areas. Mutants that showed obvious phenotypes compared to the overall appearance of the fumigated seedlings were selected based on the assumption that the vast majority of mutants does not have an NO<sub>2</sub> phenotype different from the Col-0 background line. NO2 sensitive and tolerant mutant seedlings without roots were sampled into polypropylene tubes containing 10 to 12 glass beads (1.7-2.0 mm, Roth) and were immediately frozen in liquid N2. The sampled seedlings were homogenized twice for 10 s in a Silamat S6 bead mill (Ivoclar Vivadent). DNA was extracted using the DNeasy Plant Mini kit (Qiagen) and quantified by Quant-iT Picogreen dsDNA (Life Technologies) according to the manufacturer's instructions.

#### Adapter ligation-mediated PCR

Genes containing T-DNA insertions were identified by adapter ligation-mediated PCR following a published protocol [6]. In brief, genomic DNA from mutant plants was digested with the *Ase*1 restriction enzyme, and Ase adapters (long strand adapter 2 plus short strand of adapter Ase) were ligated to the cut sites. Fragments containing T-DNA were selectively amplified by LBa1/AP1 primer pairs binding T-DNA as well as adapter sequences. Subsequently, a nested PCR with LBb1/AP2 primers was run to selectively amplify T-DNA/gDNA junctions before sequencing. According to O'Malley et al. (2007) this step can be omitted when using Ase adapters but we found that the additional PCR enhances sequencing quality and success rate.

# T-DNA-specific target enrichment and next generation sequencing

Genomic DNA from 59 (Pool I) or 65 (Pool II) plants was pooled to give two samples each containing 1 µg DNA. The DNA was precipitated by mixing with 1/10 sample volume of 3 M sodium acetate, pH 5.2, and 2 volumes of 100% cold ethanol. Samples were then incubated at -20 °C for at least 20 min. After centrifugation at maximum speed the resulting pellet was washed once with 1 mL 70% ethanol, briefly air-dried, and finally resuspended in 100 µl of Tris-EDTA (TE) buffer. Two libraries were prepared from the DNA pools I and II using the GS-FLX+ Rapid Library (Roche) according to the manufacturer's instructions, which involved shearing of the DNA, fragment end repair, and adapter ligation. Small fragments were removed by AMPure XP beads (Beckman), and the appropriate length distribution of DNA fragments between 500 and 1500 bp was confirmed by measurements with the Bioanalyzer 2100 (Agilent). Libraries I and II were amplified by ligation Frank et al. BMC Plant Biology (2019) 19:539 Page 10 of 12

mediated PCR using adapter-specific Rapid-A and -B primers as described in the NimbleGen SeqCap EZ Library LR User's Guide v2.0 [25], and DNA was quantified with the Quant-iT Picogreen dsDNA assay (Life Technologies).

DNA fragments containing T-DNA were isolated from the amplified libraries using T-DNA border-specific 70mer probes, the SegCap EZ Hybridization and Wash Kit (Roche Nimblegen), and SeqCap EZ Developer Reagent (Roche Nimblegen). The following biotinylated 70mer probes designed to bind T-DNA left- and right border (LB/RB) sequences of the Agrobacterium tumefaciens vector pROK2 were employed for target enrichment: LB1: GAG CTG TTG GCT GGC TGG TGG CAG GAT ATA TTG TGG TGT AAA CAA ATT GAC GCT TAG ACA ACT TAA TAA C; LB2: TAT ATT GTG GTG TAA ACA AAT TGA CGC TTA GAC AAC TTA ATA ACA CAT TGC GGA CGT TTT TAA TGT ACT G; LB3: CAA CAG CTG ATT GCC CTT CAC CGC CTG GCC CTG AGA GAG TTG CAG CAA GCG GTC CAC GCT GGT TTG CCC C; RB1: CTC CGC TCA TGA TCA GAT TGT CGT TTC CCG CCT TCA GTT TAA ACT ATC AGT GTT TGA CAG GAT ATA TTG G; RB2: CAG TGT TTG ACA GGA TAT ATT GGC GGG TAA ACC TAA GAG AAA AGA GCG TTT ATT AGA ATA ATC GGA TAT T; RB3: GGG TAA ACC TAA GAG AAA AGA GCG TTT ATT AGA ATA ATC GGA TAT TTA AAA GGG CGT GAA AAG GTT TAT C.

Following the NimbleGen SeqCap EZ Library LR User's Guide v2.0 a hybridization mixture was prepared including among others 1 µg amplified library DNA and the six different 70mer probes. The latter were adjusted to  $3.75 \times 10^6$  molecules of each oligonucleotide in a total volume of 4.5 µl water, replacing the SeqCap EZ Library in the User's Guide [7]. The hybridization mixture was incubated at 95 °C for 10 min and then at 47.5 °C for 40 h [7]. DNA fragments bound to the biotinylated probes were captured with Dynabeads M-270 streptavidin (Invitrogen) and amplified using 454 Rapid-A and -B primers. Target enrichment was verified by qPCR with the T-DNA right border-specific primers T-DNA\_R\_Rev CTG TGG TTG GCA TGC ACA TAC and T-DNA\_R\_For AGA TTG TCG TTT CCC GCC TT. Small fragments, primers, and primer dimers were removed with AMPure XP beads, the fragment length distribution was checked by Bioanalyzer 2100 (Agilent), and DNA concentration was determined with the Quant-iT Picogreen dsDNA assay (Life Technologies). Emulsion PCR, emulsion breaking, and sequencing using a 454 Genome Sequencer FLX instrument were performed as described in the NimbleGen SeqCap EZ Library LR User's Guide v2.0 [26].

As Roche stopped its service for GS FLX+ sequencing in 2017, future studies will be performed using the Illumina MiSeq platform (Illumina). Here, the TruSeq DNA sample prep kit (Illumina) will be applied for sample preparation as described by the manufacturer. The read length of Illumina sequencing is in the range of 250 bases per read, if the paired end modus is used. Illumina sequencing has previously been reported to be applicable for targeted genome sequencing [7].

# Bioinformatic analysis and T-DNA insertion line identification

Bioinformatic analysis was performed using the CLC Genomics Workbench 11 software (Qiagen). Two datasets of 207,785 reads (dataset 1) and 154,847 reads (dataset 2) were used for the identification of T-DNA insertion sites. To identify T-DNA-containing sequences all reads were blasted against the T-DNA insertion vector pBIN-pROK2 using the BLASTN function under default settings. By filtering the data for high-scoring segment pair (HSP) length  $\geq 30$  and  $\leq 343$  bp and an evalue ≤5,72\*E<sup>-11</sup> unspecific hits and parts of "only T-DNA" hits were filtered out. Then the reads were further trimmed with the following settings in the CLC Genomic Workbench: read length > 40 bp, low quality limit of 0.05, max. 2 ambiguous bases. Finally, one nucleotide on the 5' terminus and adapter sequences were cleaved off and vector sequences from the vector pBIN-pROK2 were labeled as trimmed and therefore ignored for the mapping. Afterwards the remaining reads were mapped against the Arabidopsis genome sequence. Sequence alignments were performed with the ApE-Plasmid Editor (http://jorgensen.biology.utah.edu/wayned/ape/).

# Re-screen using ion leakage- and stomatal conductance measurements

Selected candidate mutants and WT plants were grown in pots at 65–85 µmol m<sup>-2</sup> s<sup>-1</sup> light intensity, 14 h light– 10 h dark cycle, 20 °C/18 °C (day/night), and 65-68% relative humidity. Three-week-old plants were used for the ion leakage measurements whereas basal stomatal conductance was determined in 4-week-old plants. A fumigation chamber housing individual plants up to one tray of plants was used for re-screening of selected candidate mutants as described recently [10, 27]. Plants were fumigated for 40 min with 30 ppm NO<sub>2</sub>, which caused visible symptoms in 40-60% of the leaf area in WT plants. Immediately after fumigation 2 seedlings were collected into 30 mL of deionized water, and the background water conductivity (μS cm<sup>-1</sup>) was determined using a conductivity meter (GLM 020A, Greisinger Electronic). After 24 h the sample conductivities were measured and the background values were subtracted. The resulting corrected sample conductivities

Frank et al. BMC Plant Biology (2019) 19:539 Page 11 of 12

were normalized to their respective conductivities measured after freezing and reheating to RT (i.e. 100% conductivity). Results are given as relative ion leakage. Basal stomatal conductance was measured in untreated plants using the Leaf Porometer Model SC-1 (Decagon Devices). The measurements took place in the growth chamber at  $2\,h-4\,h$  after start of the light period.

#### **Supplementary information**

**Supplementary information** accompanies this paper at https://doi.org/10. 1186/s12870-019-2162-7.

**Additional file 1: Figure S1.**  $NO_2$ -induced cell death is associated with UV-induced green-blue fluorescence. **a.** Fumigation with 30 ppm  $NO_2$  for 1 h caused partial leaf collapse as visualized under white light (lower panel) and UV-induced emission of green-blue fluorescence in dying leaf areas (upper panel). Pictures were taken immediately after the fumigation. **b.** Untreated plants show blue fluorescence under UV illumination.

**Additional file 2: Table S1.** Results of sequencing after adapter ligation-mediated PCR.

**Additional file 3: Figure S2.** The identified T-DNA insertion site in the gene *AT2G16630* corresponds to the mutant line SALK\_042357. **a.** Sequence reads from dataset 1 (red color) map to a T-DNA insertion in the 5' UTR region of *AT2G16630* (yellow color represents the coding sequence (CDS); blue color represents the gene sequence). One sequence read mapped to the left and 2 reads mapped to the right border of the T-DNA insertion (bold red lines). The T-DNA left border/genomic DNA junction of the mutant line SALK\_042357 (see http://signal.salk.edu/cgi-bin/tdnaexpress? JOB = TEXT&TYPE = DATA&QUERY=SALKseq\_042357.1) is indicated in green color. **b.** The SALK\_042357 sequence and the mapped sequencing read from our study are identical, except for one base. Screenshots from CLC Genomics workbench.

**Additional file 4: Figure S3.** Visible  $NO_2$ -induced symptoms of the rescreened mutants. **a.** Basal ion leakage is similar between WT and the rescreened candidate mutants (n = 3). **b.** Plants before and **c.** at 72 h after treatment with 30 ppm  $NO_2$  for 40 min.

#### Acknowledgements

We thank Dr. Normand Brisson for helpful suggestions on the design of the 70mer probes. We thank Dr. Rafał Maksym and Dr. Anton Schäffner for managing the SALK mutant collection. The support of Andreas Albert and Hans Lang with the set-up of the fumigation chamber is acknowledged.

#### Authors' contributions

FG and JD planned the research. FG carried out the  $NO_2$  screen. DM performed the adapter ligation-mediated PCR. SK, ME, and MS did the next generation sequencing. UF analyzed the sequencing results. FG and UF wrote the manuscript. All authors have read and approved the manuscript.

#### **Funding**

Not applicable.

#### Availability of data and materials

The datasets used and/or analysed during the current study are available from the corresponding author on reasonable request.

#### Ethics approval and consent to participate

Not applicable.

#### Consent for publication

Not applicable.

#### Competing interests

The authors declare that they have no competing interests.

#### **Author details**

<sup>1</sup>Institute of Biochemical Plant Pathology, Helmholtz Zentrum München, Neuherberg, Germany. <sup>2</sup>Research Unit for Comparative Microbiome Analysis, Helmholtz Zentrum München, Neuherberg, Germany. <sup>3</sup>Scientific Computing Research Unit, Helmholtz Zentrum München, Neuherberg, Germany. <sup>4</sup>Chair of Biochemical Plant Pathology, Technische Universität München, Freising, Germany.

Received: 11 September 2019 Accepted: 26 November 2019 Published online: 04 December 2019

#### References

- Alonso JM, Ecker JR. Moving forward in reverse: genetic technologies to enable genome-wide phenomic screens in Arabidopsis. Nat Rev Genet. 2006;7(7):524–36.
- Tadege M, Wen J, He J, Tu H, Kwak Y, Eschstruth A, et al. Large-scale insertional mutagenesis using the Tnt1 retrotransposon in the model legume Medicago truncatula. Plant J. 2008;54(2):335–47.
- Wang N, Long T, Yao W, Xiong L, Zhang Q, Wu C. Mutant resources for the functional analysis of the rice genome. Mol Plant. 2013;6(3):596–604.
  Available from. https://doi.org/10.1093/mp/sss142.
- Lo SF, Fan MJ, Hsing YI, Chen LJ, Chen S, Wen IC, et al. Genetic resources offer efficient tools for rice functional genomics research. Plant Cell Environ. 2016;39(5):998–1013.
- Pérez-Martín F, Yuste-Lisbona FJ, Pineda B, Angarita-Díaz MP, García-Sogo B, Antón T, et al. A collection of enhancer trap insertional mutants for functional genomics in tomato. Plant Biotechnol J. 2017;15(11):1439–52.
- O'Malley RC, Alonso JM, Kim CJ, Leisse TJ, Ecker JR. An adapter ligationmediated pcr method for high-throughput mapping of t-dna inserts in the arabidopsis genome. Nat Protoc. 2007;2(11):2910–7.
- Lepage É, Zampini É, Boyle B, Brisson N. Time- and Cost-Efficient Identification of T-DNA Insertion Sites through Targeted Genomic Sequencing. PLoS One. 2013;8(8):e70912.
- Pryor WA. Free radical biology and medicine: it's a gas, man! AJP Regul Integr Comp Physiol. 2006;291(3):R491–511 Available from: http://ajpregu. physiology.org/cgi/doi/10.1152/ajpregu.00614.2005.
- Groß F, Durner J, Gaupels F. Nitric oxide, antioxidants and prooxidants in plant defence responses. Front Plant Sci. 2013;4(Oct):419.
- Mayer D, Mithöfer A, Glawischnig E, Georgii E, Ghirardo A, Kanawati B, et al. Short-term exposure to nitrogen dioxide provides basal pathogen resistance. Plant Physiol. 2018;178(1):468–87.
- Kasten D, Mithöfer A, Georgii E, Lang H, Durner J, Gaupels F. Nitrite is the driver, phytohormones are modulators while NO and H2O2 act as promoters of NO2-induced cell death. J Exp Bot. 2016;67(22):6337–49 Available from: http://jxb.oxfordjournals.org/lookup/doi/10.1093/jxb/erw401.
- Qu Y, Wang Y, Yin Q, Huang L, Jiang Y, Li G, et al. Multiple biological processes involved in the regulation of salicylic acid in Arabidopsis response to NO2 exposure. Environ Exp Bot. 2018;149:9–16.
- Desaki Y, Miya A, Venkatesh B, Tsuyumu S, Yamane H, Kaku H, et al. Bacterial lipopolysaccharides induce defense responses associated with programmed cell death in rice cells. Plant Cell Physiol. 2006;47(11):1530–40.
- Naito K, Taguchi F, Suzuki T, Inagaki Y, Toyoda K, Shiraishi T, et al. Amino acid sequence of bacterial microbe-associated molecular pattern flg22 is required for virulence. Mol Plant-Microbe Interact [Internet]. 2008;21(9):1165–74 Available from: http://apsjournals.apsnet.org/doi/10.1094/MPMI-21-9-1165.
- Thomma BPHJ, Nürnberger T, Joosten MHAJ. Of PAMPs and effectors: the blurred PTI-ETI dichotomy. Plant Cell. 2011;23(1):4–15 Available from: http:// www.plantcell.org/lookup/doi/10.1105/tpc.110.082602.
- Delledonne M, Zeier J, Marocco A, Lamb C. Signal interactions between nitric oxide and reactive oxygen intermediates in the plant hypersensitive disease resistance response. Proc Natl Acad Sci U S A. 2001;98(23):13454–9 Available from: http://www.pnas.org/content/98/23/13454.abstract.
- Gaupels F, Kuruthukulangarakoola GT, Durner J. Upstream and downstream signals of nitric oxide in pathogen defence. Curr Opin Plant Biol. 2011;14(6):707–14.
- Liang YK, Xie X, Lindsay SE, Wang YB, Masle J, Williamson L, et al. Cell wall composition contributes to the control of transpiration efficiency in Arabidopsis thaliana. Plant J. 2010;64(4):679–86.
- Hernandez-Blanco C, Feng DX, Hu J, Sanchez-Vallet A, Deslandes L, Llorente F, et al. Impairment of cellulose synthases required for Arabidopsis secondary cell

Frank et al. BMC Plant Biology (2019) 19:539 Page 12 of 12

- wall formation enhances disease resistance. Plant Cell Online. 2007;19(3):890–903 Available from: http://www.plantcell.org/cgi/doi/10.1105/tpc.106.048058.
- García MJ, Lucena C, Romera FJ, Alcántara E, Pérez-Vicente R. Ethylene and nitric oxide involvement in the up-regulation of key genes related to iron acquisition and homeostasis in Arabidopsis. J Exp Bot. 2010;61(14):3885–99.
- 21. Ma L, Tian T, Lin R, Deng XW, Wang H, Li G. Arabidopsis FHY3 and FAR1 regulate light-induced myo-inositol biosynthesis and oxidative stress responses by transcriptional activation of MIPS1. Mol Plant. 2016;9(4):541–57. Available from. https://doi.org/10.1016/j.molp.2015.12.013.
- Nekrasov V, Li J, Batoux M, Roux M, Chu ZH, Lacombe S, et al. Control of the pattern-recognition receptor EFR by an ER protein complex in plant immunity. EMBO J. 2009;28(21):3428–38. Available from. https://doi.org/10. 1038/emboj.2009.262.
- Mondragón-Palomino M, John-Arputharaj A, Pallmann M, Dresselhaus T. Similarities between reproductive and immune pistil Transcriptomes of Arabidopsis species. Plant Physiol. 2017;174(3):1559–75.
- Alonso JM, Stepanova AN, Leisse TJ, Kim CJ, Chen H, Shinn P, et al. Genome-wide insertional mutagenesis of *Arabidopsis thaliana*. Science (80-). 2003;301(5633):653–7.
- NimbleGen R. NimbleGen SeqCap EZ Library SR User's Guide v1.0; 2011. p. 1–52. Available from: papers3://publication/uuid/23462899-D688-49F9-A6BC-C40ES28ECE6A
- Pilloni G, Granitsiotis MS, Engel M, Lueders T. Testing the limits of 454 pyrotag sequencing: Reproducibility, quantitative assessment and comparison to T-RFLP fingerprinting of aquifer microbes. PLoS One. 2012; 7(7):e40467.
- 27. Kasten D, Durner J, Gaupels F. Gas alert: The NO<inf>2</inf> pitfall during NO fumigation of plants. Front Plant Sci. 2017;8(January):85.

#### **Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

#### Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

#### At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

