



## Isolation and culture of oculomotor, trochlear, and spinal motor neurons from prenatal *Is<sup>lMN</sup>:GFP* transgenic mice

Ryosuke Fujiki<sup>1,2,3,4,9</sup>, Joun Y. Lee<sup>1,2,10</sup>, Julie A. Jurgens<sup>1,2,3,7</sup>, Mary C. Whitman<sup>2,5,6</sup>, Elizabeth C. Engle<sup>1,2,3,4,5,6,7,8</sup>

<sup>1</sup>Department of Neurology, Boston Children's Hospital, Massachusetts 02115, USA.

<sup>2</sup>FM Kirby Neurobiology Center, Boston Children's Hospital, Massachusetts 02115, USA.

<sup>3</sup>Department of Neurology, Harvard Medical School, Boston, Massachusetts 02115, USA.

<sup>4</sup>Medical Genetics Training Program, Harvard Medical School, Boston, Massachusetts 02115, USA.

<sup>5</sup>Department of Ophthalmology, Boston Children's Hospital, Boston, Massachusetts 02115, USA.

<sup>6</sup>Department of Ophthalmology, Harvard Medical School, Boston, Massachusetts 02115, USA.

<sup>7</sup>Broad Institute of M.I.T. and Harvard, Cambridge, Massachusetts 02142, USA.

<sup>8</sup>Howard Hughes Medical Institute, Chevy Chase, Maryland 20815, USA.

<sup>9</sup>Present address: Department of Neurology, Kokura Memorial Hospital, Kitakyushu, Fukuoka 802-0001, Japan

<sup>10</sup>Present address: Department of Genetics, Albert Einstein College of Medicine, Bronx, New York 10461, USA

### Abstract

Oculomotor neurons (CN3s) and trochlear neurons (CN4s) exhibit remarkable resistance to degenerative motor neuron diseases such as amyotrophic lateral sclerosis (ALS) when compared to spinal motor neurons (SMNs). The ability to isolate and culture primary mouse CN3s, CN4s and SMNs would provide an approach to study mechanisms underlying this selective vulnerability. To date, most protocols have been for heterogeneous cell cultures, which can confound interpretation of experimental outcomes. To minimize the problems associated with mixed-cell populations, pure cultures are indispensable. Here, we describe in detail the first protocol to efficiently purify and cultivate CN3s/CN4s alongside SMNs counterparts from the same embryos using embryonic day 11.5 (E11.5) *Is<sup>lMN</sup>:GFP* transgenic mouse embryos. We provide details on tissue dissection and dissociation, FACS-based cell isolation, and *in vitro* cultivation of cells from CN3/CN4 and SMN nuclei. This protocol adds a novel *in vitro* CN3/CN4 culture system to existing protocols and simultaneously provides a pure species- and age-matched SMN culture for comparison. Analyses focusing on the morphological, cellular, molecular, and electrophysiological

---

Corresponding Author: Elizabeth C. Engle, [elizabeth.Engle@childrens.harvard.edu](mailto:elizabeth.Engle@childrens.harvard.edu), Tel: (617) 919-4030, Fax: (617) 730-4834.

#### DISCLOSURES

The authors declare no conflict of interest.

characteristics of motor neurons are feasible in this culture system. We hope this protocol will enable research into the mechanisms that define motor neuron development, selective vulnerability, and disease.

## SUMMARY

Here, we present a protocol to yield homogeneous cell cultures of primary oculomotor, trochlear, and spinal motor neurons. These cultures can be used for comparative analyses of the morphological, cellular, molecular, and electrophysiological characteristics of ocular and spinal motor neurons.

## Keywords

Motor neuron; oculomotor neuron; trochlear neuron; primary culture; mouse embryonic motor neuron culture; FACS; *Is<sup>MN</sup>:GFP* transgenic mouse; cell purification; cell isolation

## INTRODUCTION

The culture of primary motor neurons is a powerful tool which enables the study of neuronal development, function, and susceptibility to exogenous stressors. Motor neuron cultures are particularly useful for the study of neurodegenerative diseases such as amyotrophic lateral sclerosis (ALS)<sup>1,2</sup>, whose disease mechanisms are incompletely understood. Interestingly, despite significant cell death of spinal motor neurons (SMNs) in both ALS patients and ALS model mice, oculomotor neurons (CN3s) and trochlear neurons (CN4s) are relatively spared<sup>1,3-9</sup>. Therefore, comparative analyses of pure cultures of CN3s/CN4s and SMNs could give important clues about mechanisms underlying relative vulnerability. Unfortunately, a major barrier to such analyses has been the inability to grow purified cultures of these motor neurons.

Many protocols have been described for the purification of SMNs from animal models. Most of these protocols use density gradient centrifugation<sup>10-12</sup> and/or p75<sup>NTR</sup>-antibody-based cell-sorting panning techniques<sup>13-16</sup>. Density gradient centrifugation exploits the larger size of SMNs relative to other spinal cells, whereas p75<sup>NTR</sup> is an extracellular protein expressed exclusively by SMNs in the spinal cord. Nearly 100% pure SMN cultures have been generated by one or both of these protocols<sup>11,12,14</sup>. These protocols have not been successful for generating CN3/CN4 cultures because CN3s/CN4s do not express p75<sup>NTR</sup>, are smaller than SMNs and therefore more difficult to isolate based on size, and other specific CN3/CN4 markers have not been identified. Instead, *in vitro* studies of CN3s or CN4s have relied on dissociated<sup>17-21</sup>, explant<sup>17,22-26</sup>, and slice<sup>27,28</sup> cultures, which are composed of heterogeneous cell types, and no protocols have existed for isolation and culture of primary CN3s or CN4s.

Here, we describe a protocol for the visualization, isolation, purification, and cultivation of CN3, CN4, and spinal motor neurons from the same embryonic day 11.5 (E11.5) *Is<sup>MN</sup>:GFP* transgenic mice<sup>29</sup> (Figures 1 and 2A). *Is<sup>MN</sup>:GFP* specifically labels motor neurons with a farnesylated GFP that localizes to the cell membrane. This protocol enables species- and

age-matched comparison of multiple types of motor neurons in order to elucidate pathological mechanisms in motor neuron disease.

## PROTOCOL

All experiments utilizing laboratory animals were performed in accordance with NIH guidelines for the care and use of laboratory animals, and with approval of the Animal Care and Use Committee of Boston Children's Hospital.

### 1. Setting up Timed Matings Prior to Dissection

- 1.1. To generate prenatal embryonic mice for motor neuron harvest, set up timed matings between adult *Isl<sup>MN</sup>:GFP* transgenic mice 11.5 days prior to the day of neuron isolation. Before setting up timed matings, weigh each female mouse. For the purpose of developing this protocol, 129S1/C57BL/6J *Isl<sup>MN</sup>:GFP* mice aged 2–9 months were used and timed matings were set up in the evening.
- 1.2. Examine female mice for vaginal plugs the following morning. The date on which the plug is identified is embryonic day 0.5.
- 1.3. Weigh female mice and examine for pups using ultrasound (see Table of Materials) between E8.5–11. In addition to recognition of vaginal plugs, successful mating can be confirmed by detection of weight gain in female mice (usually >1.5g on E9.5, if there are more than 5–6 embryos) and visual confirmation of embryos under ultrasound. Embryos are easily detectable by ultrasound after E9.5. Females that gain weight as described above are more often pregnant than those that do not gain weight, so ultrasounds are conducted only on females that have gained weight. However, female mice can gain weight for reasons other than pregnancy, so weight gain alone is not a reliable indicator of pregnancy. Ultrasound confirmation prevents unnecessary sacrifice of females that are not pregnant, but is not crucial if unavailable.

### 2. Dissection Conditions and Preparation of Instruments

- 2.1. Perform all coating (except acid-cleaning of coverslips), media preparation, tissue dissociation (except centrifugation and incubation), and culture work in a laminar flow hood to ensure the sterility of media and embryonic motor neurons.
- 2.2. With careful attention to sterile technique, tissue dissection can be conducted outside of a laminar flow hood with minimal risk of contamination.
- 2.3. Sterilize one dissection plate, one pair of micro-dissecting scissors, one pair of thumb dressing forceps, two pairs of Dumont #5 tweezers, one micro dissection knife, and one Moria mini perforated spoon by immersing in 70% ethanol prior to use.

### 3. PDL/Laminin Coating of Dishes/Coverslips

Note 1. Culture dissociated primary motor neurons in 96-well or 24-well plates, depending on the number of cells required for the application. Cells can be imaged directly in the tissue

culture plate without the use of coverslips if the wells are optically transparent and of thicknesses compatible with imaging.

- 3.1.** For applications requiring coverslips, prepare acid-cleaned, sterilized, and air-dried coverslips at least two days prior to neuron isolation as previously described<sup>30</sup>. Batches of coverslips can be prepared in this manner well in advance of experiments and can be stored up to 6 months without impact on experimental quality.
- 3.2.** Prepare a working solution of 20 µg/ml poly D-lysine (PDL) in phosphate buffered saline (PBS) two days prior to neuron isolation.
  - 3.2.1.** Aliquot PDL (1mg/ml) in advance and store at –20 °C as a stock solution.
- 3.3.** Cover the surface of each coverslip or well of the tissue culture plate with sufficient PDL solution working solution (e.g. 100 µl per well on 96-well plates or 500 µl per well on 24-well plates with or without coverslips). Incubate overnight at 37 °C.
- 3.4.** The following day, wash ×3 with sterilized water.
  - 3.4.1.** Parafilm-sealed and washed PDL-coated plates can be stored at 4 °C for up to 1 month if they are dried completely after the final wash.
- 3.5.** Prepare a working solution with 10 µl laminin (1.1–1.2mg/ml) in 1.2 ml PBS.
  - 3.5.1.** Aliquot laminin stocks (1.1–1.2mg/ml) in advance and store at –80 °C.
- 3.6.** Cover the surface of each coverslip or well with sufficient laminin solution to evenly coat the surface. Incubate for at least 2 h at 37 °C prior to use. PDL/laminin-coated plates and coverslips can be stored up to 1 week at 37 °C, but freshly coated laminin is preferred. Remove laminin directly before plating.<sup>30</sup> Laminin should not be allowed to dry out. If the plates are to be stored for more than several hours, they should be wrapped in parafilm.

#### **4. Preparation of Dissection and Motor Neuron Culture Medium and Dissociation Solutions**

- 4.1.** Prepare dissection medium.
  - 4.1.1.** Thaw heat-inactivated horse serum overnight at 4 °C, prepare 1 ml aliquots, and store at –20 °C. Thaw aliquots at RT directly before use.
  - 4.1.2.** Store B27-supplement (50X) in 1 ml aliquots at –20 °C. Thaw aliquots at RT immediately before use. Avoid freeze/thaw cycles.
  - 4.1.3.** Store GlutaMAX supplement (100X) at 4 °C or –20 °C. Divide into 0.5 ml aliquots if storing at –20 °C.

- 4.1.4.** Store penicillin-streptomycin (10,000 U/ml) in 1 ml aliquots at  $-20^{\circ}\text{C}$ . Thaw aliquots at RT immediately before use.
- 4.1.5.** To make dissection medium, mix 9.4 ml Hibernate E with 200  $\mu\text{l}$  horse serum (2%), 200  $\mu\text{l}$  of 50X B27 supplement (1X), 100  $\mu\text{l}$  of 100X GlutaMAX supplement (1X), and 100  $\mu\text{l}$  of 10,000 U/ml penicillin-streptomycin (100 U/ml). All parentheses indicate final concentrations of each reagent. Use this medium for collecting dissected tissues and making the final suspension of dissociated cells.
- 4.1.6.** Add 500  $\mu\text{l}$  dissection medium to individual 1.7 mL microcentrifuge tubes for tissue collection the day before neuron isolation. Prepare one tube for each tissue type that will be collected (e.g. positive control, negative control, CN3/CN4, SMN) and store at  $4^{\circ}\text{C}$  prior to use.
- 4.1.7.** Combine 49 ml Hibernate E low fluorescence media with 1 ml B27 supplement (1X final concentration) and fill a 24-well plate with this medium (2 ml/well) the day before neuron isolation. This dish will be used for collecting mouse embryos. Store at  $4^{\circ}\text{C}$  prior to use.
- 4.2.** Prepare motor neuron culture medium.
- 4.2.1.** Prepare 250  $\mu\text{l}$  aliquots of 25 mM 2-mercaptoethanol in Leibovitz's L15 medium. Store at  $-20^{\circ}\text{C}$ .
- 4.2.2.** Generate 5  $\mu\text{l}$  aliquots of 100  $\mu\text{g}/\text{ml}$  solutions of BDNF, CNTF, and GDNF diluted in sterilized water. Store at  $-80^{\circ}\text{C}$  and thaw aliquots at RT immediately before use.
- 4.2.3.** Prepare 10 mM forskolin solution by adding DMSO (64  $\mu\text{l}$ ) to 5 mg (1.0670  $\mu\text{mol}$ ) of forskolin and vortex well to dissolve completely. Then add sterilized water (1.003 ml) to the DMSO solution and vortex well. Store 12  $\mu\text{l}$  aliquots of 10 mM forskolin at  $-20^{\circ}\text{C}$  and thaw aliquots at RT immediately before use.
- 4.2.4.** Prepare 12  $\mu\text{l}$  aliquots of 100 mM isobutylmethylxanthine (IBMX) diluted in DMSO. Store at  $-20^{\circ}\text{C}$  and thaw aliquots at RT immediately before use.
- 4.2.5.** To make motor neuron culture medium, mix 9.4 ml neurobasal medium with 200  $\mu\text{l}$  horse serum (2%), 200  $\mu\text{l}$  of 50X B27 supplement (1X), 100  $\mu\text{l}$  of 100X GlutaMAX supplement (1X), 100  $\mu\text{l}$  of 10,000 U/ml penicillin-streptomycin (100 U/ml), and 20  $\mu\text{l}$  of 25 mM 2-mercaptoethanol (50  $\mu\text{M}$ ), preferably directly before use. This step can be performed up to 1 day before neuron isolation. All parentheses indicate final concentrations of each reagent.

- 4.2.6.** Just before use, add 1  $\mu$ l each of 100  $\mu$ g/ml BDNF (10 ng/ml), CNTF (10 ng/ml), and GDNF (10 ng/ml); and 10  $\mu$ l of 10 mM forskolin (10  $\mu$ M); and 10  $\mu$ l of 100 mM IBMX (100  $\mu$ M) to motor neuron culture medium. Pre-warm medium to 37 °C. All parentheses indicate final concentrations of each reagent.
- 4.3.** Prepare dissociation solutions.
- 4.3.1.** Prepare a papain solution (20 units/ml papain and 0.005% DNase) and an albumin-ovomucoid inhibitor solution (1mg/ml ovomucoid inhibitor, 1mg/ml albumin, and 0.005% DNase) following manufacturer's instruction.
- 4.3.2.** Prepare 500  $\mu$ l aliquots of each of ovomucoid inhibitor and papain solutions and store at -80 °C. Thaw aliquots at 37 °C immediately before use.

## 5. Ventral Midbrain and Spinal Cord dissection

Note: Perform all of the following except for steps 5.1.1–5.1.3 and 5.1.5–5.1.6 under a fluorescence dissection stereomicroscope. Total dissection time per experiment is typically 3–5 h, depending on the proficiency of dissection technique and the number of motor neurons required for each experiment.

- 5.1.** Ventral midbrain dissection
- 5.1.1.** Euthanize a pregnant mouse approximately 11.5 days post-fertilization by carbon dioxide gas and cervical dislocation.
- 5.1.2.** Spray the abdomen thoroughly with ethanol and remove the uterus using sterile micro-dissecting scissors and thumb dressing forceps.
- 5.1.3.** Wash the uterus briefly in sterile PBS, then transfer to the dissection plate filled with pre-chilled sterile PBS.
- 5.1.4.** Remove the *Isl<sup>MN</sup>:GFP*-positive embryos carefully from the uterus using sterile micro-dissecting scissors, thumb dressing forceps, and Dumont #5 tweezers in ice-cold sterile PBS under the bright light of the microscope.
- 5.1.5.** Using a sterile Moria mini-perforated spoon, transfer each embryo to a separate well of a 24-well plate filled with pre-chilled Hibernate-E low fluorescence medium supplemented with 1 $\times$  B27. Keep the 24-well plate on ice.
- 5.1.6.** Transfer one embryo to a sterile dissection plate and cover completely with ice-cold sterile Hank's balanced salt solution (HBSS).
- 5.1.7.** All of the following dissection steps are performed under fluorescein isothiocyanate (FITC) illumination of the microscope.

Using tweezers, remove the tail and face of the embryo without damaging midbrain (Figure 2Ba).

- 5.1.8.** Place the embryo prone with limbs straddled underneath and tail pointing toward the front of the microscope (toward the dissector which is indicated by \*) (Figures 2Bb).
- 5.1.9.** Using tweezers, slit open the roof of the fourth ventricle in order to generate a small opening. Use this opening to hook tweezers into the space created between the fourth ventricle and its roof, and dissect along the dorsal surface of the embryo rostral to the cortex and lateral to the floor plate and motor column (Figure 2Ca and b).
- 5.1.10.** Open the dissected tissue in an open-book manner to reveal the GFP-positive CN3 and CN4 nuclei. A small piece of tissue from the ventral midbrain containing mesenchyme, CN3, and CN4 will now be exposed.
- 5.1.11.** Carefully separate the ventral midbrain from the embryo and remove meningeal tissue using tweezers and a micro-dissecting knife.
- 5.1.12.** Dissect the bilateral GFP-positive CN3 and CN4 nuclei away from the floor plate and other GFP-negative surrounding tissue using tweezers and a micro-dissecting knife (Figure 2D). Maximize GFP-positive motor neurons in the excised tissue but take care not to damage them (avoid touching).
- 5.1.13.** If collection of separate CN3 and CN4 nuclei is desired, cut along the midline of these two nuclei (yellow dotted line in Figure 2D).
- 5.1.14.** Using a P1000 pipette, collect the dissected ventral midbrain tissue with minimal HBSS and place in the CN3/CN4-labelled 1.7 mL microcentrifuge tube filled with Dissection Media (4.1.5.). Store on ice until dissociation.
- 5.1.15.** Continue pooling ventral midbrains from additional embryos in the same tube until the total number meets the experimental requirement (refer to “8. Culture of Purified Primary Motor Neurons” for ideal cell numbers). Since tissues are subject to stress during dissociation and sorting, pooled collection of at least 10 ventral midbrains yielding approximately  $1 \times 10^4$  CN3/CN4 motor neurons is recommended.

## **5.2.** Ventral spinal cord dissection

- 5.2.1.** Keep the embryo prone with the head facing the front of the microscope (toward the dissector). Hold the embryo with one pair of tweezers and insert the tip of the other pair of tweezers into the unopened caudal part of the fourth ventricle.

- Author Manuscript
- Author Manuscript
- Author Manuscript
- Author Manuscript
- Author Manuscript
- 5.2.2. Open the rest of the hindbrain and spinal cord dorsally over the whole rostro-caudal extent of the embryo. Open by cutting dorsal tissue, starting from the fourth ventricle and working toward the central canal of the caudal spinal cord using forceps as scissors (Figure 2Ca and b). Take care not to damage the ventral spinal cord (avoid contact) during this procedure.
  - 5.2.3. Hold the embryo with one pair of tweezers and pinch off the flap of dorsal tissue on each side with the other pair of tweezers (Figures 2Ea and b).

Note: Excised dorsal tissues contain dorsal skin, mesenchyme, dorsal root ganglia (DRGs), dorsal hindbrain, and spinal cord. Remove as much of these tissues as possible without damaging SMN nuclei, since they are adhesive and can trap SMNs during filtering or cause clogging during FACS sorting.
  - 5.2.4. Remove the ventral spinal cord by using the micro-dissection knife to pierce directly below the GFP-positive SMN. Lift up the ventral spinal cord with saw-like movements on both sides (Figures 2Fa and b).
  - 5.2.5. Cut the floating ventral spinal cord transversely directly above C1 where the first GFP-positive anterior horn projects (Figure 2G). Also cut transversely at the upper boundary of the lower limb (Figure 2G). Remove the cervical (C1) - lumbar (L2-L3) portion of the ventral spinal cord after this procedure.
  - 5.2.6. Place the ventral spinal cord dorsal side up and hold by pressing the GFP-negative tissue between the GFP-positive SMN columns with one pair of tweezers. Remove the remaining attached mesenchyme, DRGs, and dorsal spinal cord by trimming both sides of the GFP-positive SMN column with the micro-dissection knife (Figure 2H). Take care to maximize GFP-positive motor neurons without damaging them.
  - 5.2.7. Using a P1000 pipette, collect the dissected ventral spinal cord tissue with minimal HBSS and place in the SMN-labelled 1.7 mL microcentrifuge tube filled with Dissection Media. Store on ice until dissociation.
  - 5.2.8. Continue pooling ventral spinal cords from additional embryos in the same tube until the total number meets the experimental requirements. Since tissues are subject to stress during dissociation and sorting, pooled collection of at least 3 ventral spinal cords yielding approximately  $2.1 \times 10^4$  SMN is recommended.
  - 5.2.9. Collect facial motor neurons and extremities of *Isl<sup>MN</sup>:GFP* mouse embryos as GFP-positive and -negative controls for fluorescence-activated cell sorting (FACS), respectively. Extremities are GFP-negative because the GFP-positive axons of the spinal motor neurons have not yet extended into the extremities at this embryonic age.

## 6. Tissue Dissociation

Note: Total dissociation time is typically 1.5 h per experiment.



- 6.1.** Warm papain and albumin-ovomuroid inhibitor solution aliquots to 37 °C 30 min prior to dissociation.
- 6.2.** Briefly spin down microdissected tissues at a low speed.
- 6.3.** Using a P100 pipette, carefully remove as much Hibernate E as possible without aspirating tissues.  
  
Note: Be sure to remove all residual Hibernate E after this step to avoid reducing the efficacy of papain dissociation in the next step.
- 6.4.** Add the appropriate volume of papain solution (Table 1) to each of the 1.7 mL microcentrifuge tubes containing microdissected tissue samples. The appropriate volume of papain for dissociation was determined in order to maximize effective dissociation while minimizing stress on the cells.
- 6.5.** Gently triturate  $\times 8$  with a P200 pipette. Perform all trituration steps gently to preserve motor neuron viability.
- 6.6.** Incubate the tubes containing the tissues for 30 min at 37 °C, agitating by finger flicking  $\times 10$  every 10 min.
- 6.7.** Gently triturate each suspension  $\times 8$  with a P200 pipette after incubation.
- 6.8.** Spin down the cells at 300  $\times g$  for 5 min.
- 6.9.** To ensure the efficacy of ovomucoid inhibition in the next step, use a P1000 pipette to remove and discard as much supernatant as possible without aspirating the tissues.
- 6.10.** Resuspend pellets in the appropriate volume of albumin-ovomuroid inhibitor solution (Table 1) by gently triturating  $\times 8$  with a P200 pipette.
- 6.11.** Wait 2 min to allow any remaining pieces of undissociated tissue to settle to the bottom of the tube.
- 6.12.** Collect as much supernatant as possible without aspirating undissociated tissues using a P200 pipette. Transfer the supernatant to fresh 1.7 mL microcentrifuge tubes.
- 6.13.** Typically, some chunks of tissue will remain undissociated after step 6.12. To maximize the final yield of dissociated cells while minimizing stress on the cells dissociated previously (contained in the supernatant of step 6.12), repeat steps 6.10–6.12 for the undissociated tissues that remain in the original 1.7 mL microcentrifuge tubes.
- 6.14.** Spin down the cells at 300  $\times g$  for 5 min.
- 6.15.** Carefully remove and discard the supernatant using a P1000 pipette.
- 6.16.** Resuspend the pellet in the appropriate volume of Dissection Media (Table 1) by pipetting  $\times 8$  using a P1000 pipette. The appropriate volume of final suspension was determined so that cell density does not exceed  $10^7$  cells/ml, which can

block the stream of the flow cytometry machine, but also so that cells are not excessively dilute, which results in a slowed sorting speed.

- 6.17. Filter the suspensions through 70  $\mu\text{m}$  cell strainers to eliminate any large clumps or undigested tissue. Transfer the suspensions into 5 ml round bottom polystyrene test tubes and store on ice until required.

## 7. Fluorescence-activated Cell Sorting (FACS)

Note 1: This protocol was optimized using a FACS sorter equipped with a 15 mw 405 nm violet laser, a 100 mw 488 nm blue laser, a 75 mw 594 nm orange laser, and a 40 mw 640 nm red laser. Cells were sorted as sheath fluid in sterile PBS under aseptic conditions through a 100  $\mu\text{m}$  nozzle. In order to minimize cell stress, the flow rate was set to a sample pressure of 1–3, such that a maximum of 1000–4000 events per second were acquired.

Note 2: Total FACS time is typically 1–2 h per experiment.

- 7.1. Set up voltages for forward and side scatter so that the cell population can be visualized properly. Setting up appropriate voltages for cell sorting is complex and requires an experienced FACS operator.
- 7.2. First plot cells based on size [Forward Scatter Area (FSC-A)] versus internal complexity [Side Scatter Area (SSC-A)] to distinguish different cell populations. Draw a gate around live cells as indicated in Figure 3Aa and Ba to exclude debris and dead cells. Group the cells within the gated region as population 1 (P1).
- 7.3. To exclude cell clumps and doublets, plot P1 cells next based on Side Scatter Width (SSC-W) versus SSC-A. Gate the population of single cells as population 2 (P2) (Figure 3Ab and Bb).
- 7.4. Plot P2 cells based on Forward Scatter Width (FSC-W) versus FSC-A and gate the population of single cells as population 3 (P3) (Figure 3Ac and Bc).

Note: Using two consecutive gates in 7.3 and 7.4 excludes cell clumps and doublets (high FSC-W and high FSC-A).

- 7.5. Gate P3 cells based on GFP versus allophycocyanin (APC). The APC channel detects auto-fluorescence and gating on this channel avoids capturing auto-fluorescent cells. Use GFP-negative cells to adjust the voltage for FITC/GFP fluorescent channels. Ideally, position gates for these cell populations around  $10^2$ . Select gate thresholds for GFP-positive population 4 (P4) individually for each type of motor neuron (Figure 3Ad and Bd).

Note: Set the GFP gate much higher for SMNs than for CN3s/CN4s in order to obtain a pure culture (Figure 3Ad and Bd). A lower GFP gate for SMN cultures leads to contamination of the cultures by glia and non-motor neurons. This is likely because there is low-level GFP expression in some glia and non-motor neurons due to a leaky promoter. The percentage of GFP-positive cells as compared to total cells is typically 0.5–1.5% for CN3s/CN4s and 1.5–2.5% for SMNs. If dissection was successful, these numbers can be used as a benchmark

to determine the appropriate position for the GFP-positive gate (Figures 3Ae and Be).

- 7.6. Perform FACS according to manufacturer's protocol. Collect P4 cells into fresh 1.7 mL microcentrifuge tubes filled with 500  $\mu$ l motor neuron culture medium. Store on ice until plating.

Note: Although cells can be sorted directly into wells, this results in an uneven number of cells per well. Sort the cells into 1.75 mL microcentrifuge tubes and then plate manually in order to achieve more even plating distribution.

## 8. Culture of Purified Primary Motor Neurons

- 8.1. Dilute FACS-isolated CN3/CN4 and SMN suspensions with prewarmed (37 °C) motor neuron culture medium to densities of  $5 \times 10^3$  and  $1 \times 10^4$  cells/ml, respectively.

Note: One E11.5 embryo yields approximately  $1 \times 10^3$  CN3/CN4 and  $7 \times 10^3$  SMN. However, these yields rely heavily on the purity of dissected tissues, the thoroughness of cell dissociation, and the appropriate thresholding of GFP gates during FACS.

- 8.2. Transfer 96-well plates precoated with PDL and laminin from the 37 °C tissue culture incubator to the laminar flow hood and aspirate laminin from each well. Use plates and coverslips immediately without washing.

- 8.3. Add 200  $\mu$ l of diluted CN3/CN4 and SMN suspensions into the each well of PDL/laminin-coated 96-well plates. Final cell densities should be  $1 \times 10^3$  and  $2 \times 10^3$  cells/well for CN3/CN4 and SMN, respectively.

Note: Initial plating density of SMNs in 96-well plates ( $2 \times 10^3$  cells/well) is double that of CN3s/CN4s ( $1 \times 10^3$  cells/well) in order to obtain similar final motor neuron numbers and densities at 2 and 9 DIV ( $4-6 \times 10^2$  and  $2-4 \times 10^2$  cells per well, respectively).

- 8.4. Culture neurons in a 37 °C, 5% CO<sub>2</sub> incubator.

- 8.5. Feed neurons every 5 days by removing half of the old media (100  $\mu$ l) and replacing with the same volume of fresh motor neuron culture medium. Neuronal processes should become visible on 1 DIV and become thicker and longer by 14 DIV (Figure 4). Neuronal cell bodies become enlarged and tend to aggregate in long-term cultures, particularly for SMNs (Figure 4).

Note: Perform all medium and solution changes by leaving half of the original medium volume in order to avoid detaching cultured cells. This includes fixation and immunocytochemistry (ICC) steps. If all media is removed, regardless of how gently, most of the cells will detach and be washed away.

## REPRESENTATIVE RESULTS

The aim of this protocol was to highly purify and culture both primary CN3s/CN4s and SMNs long-term to enable comparative analyses of the mechanisms underlying motor neuron disorders (see Figures 1 and 2 for overview).

Once neurons have been successfully isolated and grown in culture, nearly pure primary CN3/CN4 and SMN cultures can be obtained (Figure 5A and B) and maintained for at least 14 DIV (Figures 4 and 6). Purities of CN3/CN4 and SMN cultures at 2 DIV were  $93.5 \pm 2.2\%$  and  $86.7 \pm 4.7\%$ , respectively, when assessed by ICC using the motor neuron marker Islet1 and neuronal marker TUJ1 (Figure 5B). However, these high purities rely heavily on the age of the embryos and on setting appropriate thresholds for GFP gates during FACS (Figure 3). Dissection of embryos at E10.5 is more difficult than dissection at E11.5 due to increased softness and adhesiveness of tissues, resulting in decreased motor neuron yields. However, the purities of E10.5 CN3s/CN4s and SMNs are comparable to those for E11.5 embryos (92.8% and 82.2% at 2 DIV, respectively; data obtained from a single experiment). However, the purities of CN3s/CN4s and SMNs dramatically decrease when E13.5 embryos are used, even if only the highest GFP-positive population is collected (20.7% and 7.4% at 2 DIV, respectively; data obtained from a single experiment), probably due to the expression of GFP in non-motor neurons (Figure 7). This same tendency also holds for E12.5 cultures, although it is much less dramatic. Therefore, embryos at E12.5 or older are inappropriate for use in the purification of motor neurons using this protocol.

Pure motor neuron cultures are valuable for understanding isolated growth patterns, behaviors, and vulnerabilities of motor neurons. This example demonstrates how these cultures can be used to test motor neuron responses to chemical treatment. To determine if primary CN3s/CN4s and SMNs show differential responses to endoplasmic reticulum (ER) stressors, primary monocultures of CN3s/CN4s and SMNs were obtained using this protocol and treated with varying concentrations of an ER stressor, cyclopiazonic acid (CPA). Neurons were treated with cyclopiazonic acid (CPA; 5, 10, 15, 20, 25, or 30  $\mu\text{M}$ ) or vehicle control [dimethyl sulfoxide (DMSO)] at 2 DIV and fixed 3 days later for ICC to evaluate survival ratios (Figure 8A). The number of viable neurons in each sample was counted and survival ratios were calculated as the number of viable cells in drug-treated wells divided by the number of viable cells in the wells treated with DMSO. CN3/CN4 monocultures were significantly more resistant to CPA treatment (10–25  $\mu\text{M}$ ) as compared to SMN monocultures (Figure 9 and Figure 8B)<sup>31</sup>.

In conclusion, this protocol allows for the generation of highly purified primary mouse embryonic CN3/CN4 and SMN cultures that provide a powerful and reliable system for the investigation of neuronal behavior.

## DISCUSSION

Historically, *in vitro* studies of CN3 and/or CN4 motor neurons have relied on heterogeneous cultures such as dissociated<sup>17–21</sup>, explant<sup>17,22–26</sup>, and slice<sup>27,28</sup> cultures because these cells cannot be distinguished from surrounding cells based on size, and

specific markers for these cells have not been reported. The present protocol is a comprehensive method for the isolation and culture of primary E11.5 murine CN3s/CN4s and SMN from the same embryos and confirm the high purity of the cultures. By generating pure SMN and CN3/CN4 cultures from the same mouse embryos, the protocol enables well controlled comparisons of the *in vitro* behaviors of CN3s/CN4s versus SMNs isolated from wildtype as well as mutant embryos.

The pure cultures of CN3s/CN4s and SMNs generated by this protocol allow comparative studies of morphological, cellular, molecular, and electrophysiological characteristics of these motor neurons. In theory, because other cranial motor neuron populations can be visualized and dissected from this *Is<sup>SMN</sup>:GFP* transgenic mouse line (including abducens, motor trigeminal, facial, and hypoglossal), this protocol could be expanded for their isolation and culture as well, provided that FACS GFP gates are adjusted appropriately. Finally, the FACS-sorted motor neurons derived from this protocol can be subjected to genomic [eg. Assay for Transposase-Accessible Chromatin using sequencing (ATAC-seq)] and/or transcriptomic (eg. RNA sequence<sup>40</sup>) analyses to study normal development and the selective vulnerability of specific motor neuron subtypes in neurodegenerative disorders.<sup>31</sup>

There are multiple steps in this protocol that are critical to maximize the number of pure, healthy motor neurons in the isolated culture system. (1) During the dissection, the GFP-negative tissues (eg. mesenchyme and DRGs) should be maximally removed without damaging the motor neurons, since these tissues are adhesive and can trap SMNs during filtering or cause clogging during FACS sorting. (2) During tissue dissociation, the minimal essential volume of papain should be used and the cells must be treated gently with minimal but sufficient trituration. Papain was used for the tissue dissociation step in this protocol, since preliminary data indicated that it is less destructive than trypsin to both CN3/CN4 and SMN. Survival ratios based on plated numbers of CN3s/CN4s and SMNs at 2 days *in vitro* (DIV) increased from 39.5% to 52.7% and from 52.3% to 58.4%, respectively, when papain was used instead of trypsin (0.25%, 4 min incubation). Although these numbers are derived from a single experiment performed before full optimization, additional reports also suggest that trypsin is suboptimal for cell extraction from nervous system tissues<sup>12,34-36</sup>. (3) During FACS, fluorescent vital dyes (propidium iodide and calcein blue) and small sorting nozzles (eg. 70  $\mu$ m) should not be used since they are deleterious to motor neuron survival. Use of large sorting nozzles (100  $\mu$ m or larger) is highly recommended since SMN cell death increases significantly when a 70  $\mu$ m nozzle is used. Setting the appropriate gating thresholds for GFP-positive cells in FACS is a critical step in order to obtain pure cultures. (4) Motor neuron cultures are supplemented with forskolin, IBMX, and growth factors (BDNF, CNTF, and GDNF). Forskolin and IBMX have been reported to additively promote SMN survival<sup>37,38</sup>. Preliminary data from the present studies suggest that forskolin and IBMX also additively increase CN3/CN4 survival. Survival ratio based on plated numbers of CN3s/CN4s at 2 DIV increased from 17.5% to 26.9%, 31.9%, and 37.0% when IBMX, forskolin, and IBMX+forskolin were added, respectively (numbers are based on a single experiment performed prior to full optimization of cell culture conditions). Perform all medium changes and washes of cultured cells by leaving half of the original volume to avoid detaching cultured cells. (5) Finally, reduce the time spent between dissection and plating of

motor neurons (e.g. by shortening dissection time using multiple dissectors) to improve the viability of the cultures.

There are four major potential problems that may arise when following this protocol. (1) Low yield of motor neurons after FACS. Potential causes include using young embryos (e.g. E10.5), which have fewer motor neurons; insufficient removal of adhesive GFP-negative tissues during dissection (e.g. mesenchyme and DRGs), which can trap motor neurons and lead to their removal during filtration; insufficient papainization/trituration during dissociation; and/or setting the GFP-positive gate too high during FACS. (2) Low purity of the motor neuron cultures, which most likely arises from the use of older embryos (eg. E12.5) and/or from setting the GFP-positive gate too low during FACS. (3) Low number of attached motor neurons in culture. Potential causes for this include inappropriate FACS sorting and/or inadequate PDL/laminin coating of plates/coverlips. (4) Low viability of cultured motor neurons. Potential causes include rough and/or prolonged dissection; excessive papainization/trituration during dissociation; inappropriate handling of cells throughout the protocol (e.g. rough pipetting of cells, failure to place cells on ice, failure to pre-chill PBS and HBSS); and/or excessive time between euthanization of pregnant mice and final plating of the cells. Use of reagents which are not fresh and/or inappropriate concentrations can also impair experimental outcomes.

There are three major limitations of this protocol. (1) *Isl<sup>MN</sup>:GFP* transgenic mice and FACS sorting are both fairly expensive. They are, however, crucial for this protocol as there is currently no alternative method capable of generating highly purified CN3s/CN4s in a more economical fashion. (2) The small E10.5-E12.5 age window for the embryonic mice, and the difficulty of confirming that appropriately aged embryos are present, especially if an ultrasound machine is not available. If only pure SMNs are required, they can be derived from E12.5–15.0 mouse embryos using methods such as gradient centrifugation<sup>10–12</sup> and/or p75<sup>NTR</sup>-antibody-based cell-sorting panning techniques<sup>13–16</sup>. (3) Protein-based assays that require a large amount of starting material (e.g. western blot analysis) thus are not feasible from these cultures due to the small yield of motor neurons (especially CN3s/CN4s). Stem cell-derived motor neurons<sup>31,39</sup>, which can be generated limitlessly, could in theory be substituted for this purpose.

## ACKNOWLEDGEMENTS


We thank Brigitte Pettmann (Biogen, Cambridge, MA, USA) for instruction in SMN dissection techniques; the Dana Farber Cancer Institute Flow Cytometry Facility, the Immunology Division Flow Cytometry Facility of Harvard Medical School, The Joslin Diabetes Center Flow Cytometry Core, Brigham and Women's Hospital Flow Cytometry Core, and Boston Children's Hospital Flow Cytometry Research Facility for FACS isolation of primary motor neurons; A.A. Nugent, A.P. Tenney, A.S. Lee, E.H. Nguyen, M.F. Rose, additional Engle laboratory members, and Project ALS consortium members for technical assistance and thoughtful discussion. This study was supported by Project ALS. In addition, R.F. was funded by the Japan Heart Foundation / Bayer Yakuin Research Grant Abroad and NIH Training grant in Genetics T32 GM007748; J.J. was supported by the NIH/NEI training program in the Molecular Bases of Eye Diseases (5T32EY007145-16) through Schepens Eye Research Institute and by the Developmental Neurology Training Program Postdoctoral Fellowship (5T32NS007473-19) through Boston Children's Hospital; M.C.W was supported by NEI (5K08EY027850) and Children's Hospital Ophthalmology Foundation (Faculty Discovery Award); and E.C.E. is a Howard Hughes Medical Institute Investigator.

## REFERENCE

1. Kiernan MC et al. Amyotrophic lateral sclerosis. *Lancet*. 377 (9769), 942–55 (2011). [PubMed: 21296405]
2. Wood-Allum C, Shaw PJ Motor neurone disease: a practical update on diagnosis and management. *Clinical medicine (London, England)*. 10 (3), 252–8 (2010).
3. Nijssen J, Comley LH, Hedlund E Motor neuron vulnerability and resistance in amyotrophic lateral sclerosis. *Acta neuropathologica*. 133 (6), 863–85 (2017). [PubMed: 28409282]
4. Gizzi M, DiRocco A, Sivak M, Cohen B Ocular motor function in motor neuron disease. *Neurology*. 42 (5), 1037–46 (1992). [PubMed: 1579227]
5. Kanning KC, Kaplan A, Henderson CE Motor neuron diversity in development and disease. *Annual review of neuroscience*. 33, 409–40 (2010).
6. Nimchinsky EA et al. Differential vulnerability of oculomotor, facial, and hypoglossal nuclei in G86R superoxide dismutase transgenic mice. *The Journal of comparative neurology*. 416 (1), 112–25 (2000). [PubMed: 10578106]
7. Angenstein F et al. Age-dependent changes in MRI of motor brain stem nuclei in a mouse model of ALS. *Neuroreport*. 15 (14), 2271–4 (2004). [PubMed: 15371748]
8. Niessen HG et al. In vivo quantification of spinal and bulbar motor neuron degeneration in the G93A-SOD1 transgenic mouse model of ALS by T2 relaxation time and apparent diffusion coefficient. *Experimental neurology*. 201 (2), 293–300 (2006). [PubMed: 16740261]
9. Spiller KJ et al. Selective Motor Neuron Resistance and Recovery in a New Inducible Mouse Model of TDP-43 Proteinopathy. *The Journal of neuroscience*. 36 (29) 7707–17 (2016). [PubMed: 27445147]
10. Graham JM Isolation of a mouse motoneuron-enriched fraction from mouse spinal cord on a density barrier. *ScientificWorldJournal*. 2, 1544–6 (2002). [PubMed: 12806137]
11. Gingras M, Gagnon V, Minotti S, Durham HD, Berthod F Optimized protocols for isolation of primary motor neurons, astrocytes and microglia from embryonic mouse spinal cord. *Journal of neuroscience methods*. 163 (1), 111–8 (2007). [PubMed: 17445905]
12. Beaudet MJ et al. High yield extraction of pure spinal motor neurons, astrocytes and microglia from single embryo and adult mouse spinal cord. *Scientific reports*. 5, 16763 (2015). [PubMed: 26577180]
13. Camu W, Henderson CE. Purification of embryonic rat motoneurons by panning on a monoclonal antibody to the low-affinity NGF receptor. *Journal of neuroscience methods*. 44 (1), 59–70 (1992). [PubMed: 1434751]
14. Arce V et al. Cardiotrophin-1 requires LIFRbeta to promote survival of mouse motoneurons purified by a novel technique. *Journal of neuroscience research*. 55 (1), 119–26 (1999). [PubMed: 9890440]
15. Wiese S et al. Isolation and enrichment of embryonic mouse motoneurons from the lumbar spinal cord of individual mouse embryos. *Nature protocols*. 5 (1), 31–8 (2010). [PubMed: 20057379]
16. Conrad R et al. Lectin-based isolation and culture of mouse embryonic motoneurons. *Journal of visualized experiments*. (55), doi: 10.3791/3200 (2011).
17. Lerner O et al. Stromal cell-derived factor-1 and hepatocyte growth factor guide axon projections to the extraocular muscles. *Developmental neurobiology*. 70 (8), 549–64 (2010). [PubMed: 20506246]
18. Ferrario JE et al. Axon guidance in the developing ocular motor system and Duane retraction syndrome depends on Semaphorin signaling via alpha2-chimaerin. *Proceedings of the National Academy of Sciences of the United States of America*. 109 (36), 14669–74 (2012).
19. Clark C, Austen O, Poparic I, Guthrie S  $\alpha$ 2-Chimaerin regulates a key axon guidance transition during development of the oculomotor projection. *The Journal of neuroscience*. 33 (42), 16540–51 (2013). [PubMed: 24133258]
20. Theofilopoulos S et al. Cholestenic acids regulate motor neuron survival via liver X receptors. *The Journal of clinical investigation*. 124 (11), 4829–42 (2014). [PubMed: 25271621]


21. Montague K, Guthrie S, Poparic I In Vivo and In Vitro Knockdown Approaches in the Avian Embryo as a Means to Study Semaphorin Signaling. *Methods in molecular biology*. 1493, 403–16 (2017). [PubMed: 27787867]
22. Porter JD, Hauser KF Survival of extraocular muscle in long-term organotypic culture: differential influence of appropriate and inappropriate motoneurons. *Developmental biology*. 160 (1), 39–50 (1993). [PubMed: 8224547]
23. Serafini T et al. Netrin-1 is required for commissural axon guidance in the developing vertebrate nervous system. *Cell*. 87 (6), 1001–14 (1996). [PubMed: 8978605]
24. Varela-Echavarría A Tucker, A. Püschel, A.W. Guthrie, S. Motor axon subpopulations respond differentially to the chemorepellents netrin-1 and semaphorin D. *Neuron*. 18 (2), 193–207 (1997). [PubMed: 9052791]
25. Irving C, Malhas A, Guthrie S, Mason I Establishing the trochlear motor axon trajectory: role of the isthmic organiser and Fgf8. *Development*. 129 (23), 5389–98 (2002). [PubMed: 12403710]
26. Chen J, Butowt R, Rind HB, von Bartheld CS GDNF increases the survival of developing oculomotor neurons through a target-derived mechanism. *Molecular and cellular neurosciences*. 24 (1), 41–56 (2003). [PubMed: 14550767]
27. Whitman MC, et al. Loss of CXCR4/CXCL12 Signaling Causes Oculomotor Nerve Misrouting and Development of Motor Trigeminal to Oculomotor Synkinesis. *Investigative ophthalmology & visual science*. 59, (12), 5201–5209 (2018). [PubMed: 30372748]
28. Whitman MC, Bell JL, Nguyen EH, Engle EC Ex Vivo Oculomotor Slice Culture from Embryonic GFP-Expressing Mice for Time-Lapse Imaging of Oculomotor Nerve Outgrowth. *Journal of visualized experiments*. (Pending Publication), e59911, In-press (2019).
29. Lewcock JW, Genoud N, Lettieri K, Pfaff SL The ubiquitin ligase Phr1 regulates axon outgrowth through modulation of microtubule dynamics. *Neuron*. 56: 604–20 (2007). [PubMed: 18031680]
30. Bibel M, Richter J, Lacroix E, Barde YA, Generation of a defined and uniform population of CNS progenitors and neurons from mouse embryonic stem cells. *Nature protocols*. 2 (5), 1034–43 (2007). [PubMed: 17546008]
31. An D, Fujiki R, et al. Stem cell-derived cranial and spinal motor neurons reveal proteostatic differences between ALS resistant and sensitive motor neurons. *eLife*. 8:e44423 (2019)
32. Huber AB, et al. Distinct roles for secreted semaphorin signaling in spinal motor axon guidance. *Neuron* 48 (6), 949–64 (2005). [PubMed: 16364899]
33. Plachta N, et al. Identification of a lectin causing the degeneration of neuronal processes using engineered embryonic stem cells. *Nature neuroscience*. 10 (6), 712–9 (2007). [PubMed: 17486104]
34. Eide L, McMurray CT Culture of adult mouse neurons. *BioTechniques*. 38 (1), 99–104 (2005). [PubMed: 15679091]
35. Brewer GJ, Torricelli JR Isolation and culture of adult neurons and neurospheres. *Nature protocols*. 2 (6), 1490–8 (2007). [PubMed: 17545985]
36. Seibenhener ML, Wotten MW Isolation and culture of hippocampal neurons from prenatal mice. *Journal of visualized experiments*. (65), doi: 10.3791/3634 (2012).
37. Hanson MG Jr., Shen S, Wiemelt AP, McMorris FA, Barres BA Cyclic AMP elevation is sufficient to promote the survival of spinal motor neurons in vitro. *The Journal of neuroscience*. 18 (18), 7361–71 (1998). [PubMed: 9736656]
38. Lamas NJ, et al. Neurotrophic requirements of human motor neurons defined using amplified and purified stem cell-derived cultures. *Plos one*. 9 (10), doi: 10.1371 (2014).
39. Mazzoni EO, et al. Synergistic binding of transcription factors to cell-specific enhancers programs motor neuron identity. *Nature neuroscience* 16 (9), 1219–27 (2013). [PubMed: 23872598]





**Figure 1: Scheme for preparation of mouse embryonic motor neurons.**


The schematic illustrates the steps involved in the isolation and culture of mouse embryonic motor neurons and the approximate time in hours or days for each step. The order of the dissection procedure for CN3/CN4 and for SMN are each labelled sequentially 1 through 4. Abbreviations: CN3/CN4: oculomotor neuron/trochlear neuron; SMN: spinal motor neuron; FACS: fluorescence-activated cell sorting; h: hour; d: day.



**Figure 2: Dissection of the ventral midbrain and the cervical (C1) - lumbar (L2-L3) portion of the ventral spinal cord.**


A) Lateral (a) and dorsal (b) views of GFP-positive motor neurons in an E11.5 *Isl<sup>MN</sup>:GFP* transgenic mouse embryo under fluorescein isothiocyanate (FITC) illumination. A whole mount E11.5 embryo was prepared as previously described<sup>32</sup> in order to make the embryo transparent. Subsequently, the embryo was analyzed by immunofluorescence labeling with anti-GFP staining (green). Images were captured under a confocal microscope. Scale bars: 200  $\mu$ m (lateral view) and 400  $\mu$ m (dorsal view). Abbreviations: S: superior; I: inferior; V:

ventral; D: dorsal. **B-H**) Dissection steps highlighted on images of E11.5 ventral midbrain and ventral spinal cord tissues taken with an equipped camera under bright light (Bb) or FITC illumination using a fluorescence dissection stereomicroscope. Scale bars: 200  $\mu\text{m}$  (D) and 1 mm (A-C, E-H). **B**) **a**) Remove the face and tail of the embryo by cutting along the red lines. **b**) Embryo positioned for dissection. Positioning of the front of the microscope is indicated by \*. **C**) Cut along the solid red line in order to slit open the roof of the fourth ventricle [**a**] lateral view and **b**] dorsal view]. Use this opening to cut along the surface of the embryo dorsal to the brain (trajectory indicated by dashed red arrow). This will expose the tissue containing mesenchyme, CN3, and CN4, which can be lifted out of the cranium. For SMN dissection, insert forceps into the same opening between the fourth ventricle and its roof, then cut toward the caudal side of the embryo (trajectory indicated by dashed yellow arrow). **D**) Final view of the ventral midbrain containing bilateral GFP-positive CN3 and CN4 nuclei. The edges of the tissue are highlighted by a red rectangle. Cut along yellow dotted line to collect CN3 and CN4 nuclei separately, if desired. **E**) After opening the rest of the hindbrain and spinal cord, pinch off flapping dorsal tissues above the red lines on both sides with tweezers [**a**] before and **b**] after]. **F**) Bilaterally remove excess tissue ventral to the spinal cord along the red line [**a**] before and **b**] after]. **G**) Cut the ventral spinal cord at the two locations indicated by the red lines. On the rostral side, cut the floating ventral spinal cord transversely above C1 where the first GFP-positive anterior horn projects. Cut the caudal end of the spinal cord transversely at the upper boundary of the lower limb. Once these cuts are made, the cervical (C1) through lumbar (L2-L3) portion of the ventral spinal cord can be dissected away. **H**) Final view of the ventral spinal cord containing GFP-positive SMN columns.




**Figure 3: Representative sort plots of ventral midbrains (A) and ventral spinal cords (B).** (Aa and Ba). Forward Scatter Area (FSC-A) versus Side Scatter Area (SSC-A) sort plot before exclusion of debris and dead cells. Ab, Bb, Ac, Bc) Sort plots for exclusion of cell clumps (b) and doublets (c) based on Width (SSC-W) versus SSC-A and Forward Scatter Width (FSC-W) versus FSC-A, respectively. Ad and Bd) Sort plots to isolate *Is<sup>SMN</sup>:GFP*-positive motor neurons. Set the GFP gate higher for SMNs (Bd) than for CN3s/CN4s (Ad) in order to obtain a pure culture. Ae and Be) Percentages of cells gated for collection by FACS sorting. %Parent represents the percentage of cells in the current gated population

relative to the number of cells in the previous gated cell population, whereas %Total represents the percentage of gated cells relative to total cells. Expected percentages of GFP-positive cells as compared to total cells (boxed in red) are 0.5–1.5 % for CN3/CN4 and 1.5–2.5% for SMN. If dissection was performed successfully, these percentages can be used as a benchmark to set up the GFP-positive gate in (**Ad and Bd**).



**Figure 4: Phase-contrast images of primary CN3/CN4 and SMN monocultures at 2, 7, and 14 DIV.**

Representative differential interference contrast images of primary CN3/CN4 and SMN cultures were captured at 2, 7, and 14 DIV with inverted fluorescence microscope using corresponding image acquisition and processing software and 40× objectives. Neuronal processes become thicker and longer by 14 DIV. Neuronal cell body sizes become enlarged and tend to aggregate in long-term cultures, especially for SMNs. Both cultures can be maintained at least 14 DIV. Scale bar: 50  $\mu$ m.




**Figure 5: Characterizations of isolated E11.5 mouse CN3/CN4 and SMN cultures.**

**A)** Representative immunocytochemistry images of E11.5 mouse CN3s/CN4s (top) and SMNs (bottom) cultured for 2 DIV. Neurons were analyzed by immunofluorescence labeling with the neuronal marker TUJ1 (green) and the motor neuron marker Islet1 (red), and nuclei were counterstained with DAPI (blue). Almost all of the cultured cells were motor neurons (TUJ1<sup>+</sup>, Islet1<sup>+</sup>). Images were captured with inverted fluorescence microscope using corresponding image acquisition and processing software and 20× objectives. Samples were imaged and processed to achieve maximum signal intensity without saturated pixels. All of


the microscopic work and image processing in the following figures were performed in these conditions unless otherwise specified. Scale bar: 100  $\mu\text{m}$ . **B)** The purities of E11.5 mouse CN3/CN4 and SMN cultures at 2 DIV. The purities of CN3/CN4 and SMN cultures were  $93.5 \pm 2.2\%$  and  $86.7 \pm 4.7\%$ , respectively. Dead neuronal cell bodies were assessed by screening for pyknotic nuclear morphology and membrane swelling. Neuronal processes that showed signs of beading and swelling were classified as degenerating processes. Cells with neither cell body death nor degenerating processes were counted as viable non-motor neurons (TUJ1<sup>+</sup>, Islet1<sup>-</sup>) or viable motor neurons (TUJ1<sup>+</sup>, Islet1<sup>+</sup>).<sup>33</sup> The purities of motor neuron cultures were calculated as the number of viable motor neurons divided by the total number of viable non-motor neurons plus viable motor neurons. Values represent the mean  $\pm$  SEM of 3 separate experiments. Not significant ( $P > 0.05$ ) by Student's t test. Cell counting was performed manually under 20 $\times$  magnification.






**Figure 6: Representative immunocytochemistry of primary CN3/CN4 and SMN monocultures at 2, 7, and 14 DIV.**

Primary CN3/CN4 and SMN cultures were analyzed at 2, 7, 14 DIV by immunofluorescence labeling with TUJ1 (green) and nuclei were counterstained with DAPI (blue). Neuronal processes become thicker and longer by 14 DIV. Neuronal cell body sizes become enlarged and tend to aggregate in long-term cultures, particularly for SMNs. Both CN3/CN4 and SMN cultures can be maintained at least 14 DIV. Images were captured under 10 $\times$  magnification. Scale bar: 200  $\mu$ m.



**Figure 7: Characterization of E13.5 mouse CN3/CN4 and SMN isolated cultures.**


E13.5 CN3/CN4 and SMN were isolated and cultured using this protocol and analyzed at 2 DIV by immunofluorescence labeling with TUJ1 (green) and Islet1 (red), and nuclei were counterstained with DAPI (blue). Many non-motor-neuronal cells (TUJ1<sup>+</sup>, Islet1<sup>-</sup>), as indicated by arrows, were present, resulting in a drastic decrease in both CN3/CN4 and SMN purity, with the decrease more pronounced in SMN cultures. Images were captured under 20× magnification. Scale bar: 100 μm.



**Figure 8: Representative application of primary motor neuron culture demonstrating that CN3s/CN4s are selectively resistant to ER stress induced by CPA.**

**A)** Experimental outline: primary CN3/CN4 and SMN monocultures were treated with CPA or vehicle control (DMSO) at 2 DIV and cell viabilities were evaluated through immunocytochemistry analysis after 3 days of treatment. This outline has been modified from published work<sup>31</sup>. **B)** Quantification of survival ratios of CN3s/CN4s and SMNs treated with 5–30  $\mu\text{M}$  CPA for 3 days from 2 DIV. Neurons were analyzed by immunofluorescent labeling of cells with TUJ1, and nuclei were counterstained with DAPI.

Survival ratios were calculated as the number of viable cells (refer to Figure 5B legend) in drug-treated wells divided by the number of viable cells in wells containing vehicle alone (DMSO). Cell counting was performed manually under 20× magnification. Values represent the mean ± SEM of 4 separate experiments. \* $P < 0.05$ ; \*\*\* $P < 0.005$  by Student's  $t$  test. This figure has been modified from published work<sup>31</sup>.



**Figure 9. Representative immunocytochemistry of primary CN3/CN4 and SMN monocultures after 3-day exposure to increasing concentrations of CPA beginning at 2 DIV.**

Neurons were analyzed by immunofluorescent labeling of cells with TUJ1 (green) and nuclei were counterstained with DAPI (blue). Primary CN3s/CN4s were more resistant to CPA treatment than primary SMNs. Images were captured under 10 $\times$  magnification. Scale bar: 200  $\mu\text{m}$ .

**Table 1:**  
**Appropriate volumes of papain, albumin-ovomuroid, and final suspension (Hibernate E)**  
**used in dissociation steps.**

The appropriate volumes of papain and albumin-ovomuroid to be used with various numbers of ventral midbrain and ventral spinal cord tissues were modified from the manufacturer's instruction after several rounds of optimization. Since tissues are subject to stress during dissociation and sorting, pooled collection of more than 10 ventral midbrains and more than 3 ventral spinal cords is recommended. The volume of papain was determined by considering the balance between effective dissociation and the stress of this procedure. The volume of albumin-ovomuroid inhibitor solution is half of that of papain. The appropriate volume of Hibernate E final suspension was determined such that cell density does not exceed  $10^7$  cells/ml, but also such that the cells do not become excessively diluted.

Number of midbrains (X)	Papain	Albumin-ovomuroid	Hibernate E
$10 \leq X \leq 20$	200 $\mu$ l	100 $\mu$ l	600 $\mu$ l
$20 < X \leq 30$	300 $\mu$ l	150 $\mu$ l	700 $\mu$ l
$30 < X \leq 40$	400 $\mu$ l	200 $\mu$ l	800 $\mu$ l
Number of spinal cords (Y)	Papain	Albumin-ovomuroid	Hibernate E
$3 \leq Y \leq 5$	200 $\mu$ l	100 $\mu$ l	500 $\mu$ l
$5 < Y \leq 10$	400 $\mu$ l	200 $\mu$ l	800 $\mu$ l
$10 < Y \leq 15$	600 $\mu$ l	300 $\mu$ l	1200 $\mu$ l

Table of Materials

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Alexa Fluor 2488-conjugated goat anti-mouse IgG (H+L)	Thermo Fisher Scientific	A-11001	1:400
Alexa Fluor 594-conjugated F(ab') <sub>2</sub> goat anti-rabbit IgG (H+L)	Thermo Fisher Scientific	A-11072	1:400
B27 Supplement (50X), serum free	Thermo Fisher Scientific	17504-044	This has 4 laser system equipped with 405, 488, 594, and 640 nm lasers. For filtering the dissociating cells before FACS.
BD FACSAria III SORP Flow Cytometer	BD Bioscience	-	
BD Falcon 70µm Nylon Cell Strainers	CORNING	352350	
BD Falcon Round Bottom Test Tubes With Snap Cap	CORNING	352054	
BDNF Human	ProSpec-Tany TechnoGene, Ltd.	CYT-207	
Cell Culture microplate, 96 well, PS, F-bottom (Chimney Well)	Greiner Bio-One International	655090	We tried multiple 96-well dishes and this was the best one for culture and analyses after ICC
Circular Cover Glasses for microscopy	Karl Hecht & Assistent	1001/14	We used this coverslip since the area was large (diamater: 14 mm).
CNTF Human	ProSpec-Tany TechnoGene, Ltd.	CYT-272	
Cycloiazonic acid from <i>Penicillium cyclospium</i>	Sigma-Aldrich	C1530	CPA. One of ER stressors.
4',6-diamidino-2-phenylindole (DAPI)	Thermo Fisher Scientific	D1306	
Dimethyl sulfoxide	Sigma-Aldrich	D2650	DMSO
Dumont #5 Forceps Inox Tip Size .05 × .01 mm Biologie Tips	Roboz Surgical Instrument	RS-5015	
Forskolin	Thermo Fisher Scientific	BP25205	
GDNF Human	ProSpec-Tany TechnoGene, Ltd.	CYT-305	
GlutaMAX supplement	Thermo Fisher Scientific	35050-061	
Hanks' Balanced Salt Solution (HBSS)	Thermo Fisher Scientific	14175-095	
HibernatE	BrainBits	HE	
HibernatE low fluorescence	BrainBits	HELF	Fluorescence which hinders observation of embryo's GFP expressions should be low.
Horse serum, heat inactivated, New Zealand origin	Thermo Fisher Scientific	26050-070	
IBMX	Tocris Cookson	2845	Isobutylmethylxanthine
Laminin	Thermo Fisher Scientific	23017-015	
Leibovitz's L15 medium	Thermo Fisher Scientific	11415064	
2-Mercaptoethanol	Sigma-Aldrich	M6250	
Micro Dissecting Scissors	Roboz Surgical Instrument	RS-5913	
Micro Knife 4.75" 1.7 × 27 mm blade	Roboz Surgical Instrument	RS-6272	

Name of Material/ Equipment	Company	Catalog Number	Comments/Description
Moria Mini Perforated Spoon	Fine Science Tools	10370-19	
mouse monoclonal antibody to neuronal class III $\beta$ -tubulin (TUBB3)	BioLegend	801202	1:500, TUJ1
Nikon Perfect Focus Eclipse Ti live cell fluorescence microscope and Elements software	Nikon	-	Differential interference contrast images and immunocytochemistry images of the cell cultures were captured with these equipments
Nitric Acid 90%, Fuming (Certified ACS)	Fisher Scientific	A202-212	For rinsing coverslips
Olympus 1.7ml Microtubes, Clear	Genesee Scientific	22-281	These are the tubes that we described "1.7 mL microcentrifuge tubes" in the context.
Papain Dissociation System	Worthington Biochemical Corp	LK003150	Papain solution and albumin-ovomucoid inhibitor solution are prepared from this kit.
Penicillin-streptomycin (10,000 U/ml)	Thermo Fisher Scientific	15140-122	
Phosphate buffered saline (PBS)	Thermo Fisher Scientific	10010-023	
Poly D-lysine (PDL)	MilliporeSigma	A-003-E	
rabbit monoclonal antibody to Islet1	Abcam	ab109517	1:200
SMZ18 and SMZ1500 zoom stereomicroscopes with DS-R1 camera	Nikon	-	Dissection was performed and images of dissected embryos and tissues are captured under these fluorescence microscopes.
Sylgard 170 Black Silicone Encapsulant - A+B 0.9 Kg kit	Dow Corning	1696157	We make dissection dishes using this kit.
TC treated Dishes, 100 $\times$ 20 mm	Genesee Scientific	25-202	We make dissection dishes using this dish.
Thum Dressing Forceps 4.5" Serrated 2.2 mm Tip Width	Roboz Surgical Instrument	RS-8100	
Transducer for LOGOQ e VET	GE Healthcare	L8-18i-RS	For ultrasound on female mice
Veterinary ultrasound machine	GE Healthcare	LOGOQ e VET	For ultrasound on female mice
Zeiss LSM 700 series laser scanning confocal microscope and Zen Software	Carl Zeiss	-	Confocal image of the embryo was captured with these equipments