

Amplex Red Assay for Measuring Hydrogen Peroxide Production from *Caenorhabditis elegans*

Ozgur Karakuzu¹, Melissa R. Cruz¹, Yi Liu^{1, 2} and Danielle A. Garsin^{1, 2, *}

¹Department of Microbiology and Molecular Genetics, The University of Texas Health Science Center at Houston, Houston TX, USA; ²MD Anderson Cancer Center UTHealth Graduate School of Biomedical Sciences, Houston TX, USA

*For correspondence: Danielle.A.Garsin@uth.tmc.edu

[Abstract] Reagents such as Amplex[®] Red have been developed for detecting hydrogen peroxide (H₂O₂) and are used to measure the release of H₂O₂ from biological samples such as mammalian leukocytes undergoing the oxidative burst. *Caenorhabditis elegans* is commonly used as a model host in the study of interactions with microbial pathogens and releases reactive oxygen species (ROS) as a component of its defense response. We adapted the Amplex[®] Red Hydrogen Peroxide/Peroxidase Assay Kit to measure H₂O₂ output from live *Caenorhabditis elegans* exposed to microbial pathogens. The assay differs from other forms of ROS detection in the worm, like dihydrofluorescein dyes and genetically encoded probes such as HyPer, in that it generally detects released, extracellular ROS rather than intracellular ROS, though the distinction between the two is blurred by the fact that certain species of ROS, including H₂O₂, can cross membranes. The protocol involves feeding *C. elegans* on a lawn of the pathogen of interest for a period of time. The animals are then rinsed off the plates in buffer and washed to remove any microbes on their cuticle. Finally, the animals in buffer are distributed into 96-well plates and Amplex[®] Red and horseradish peroxidase (HRP) are added. Any H₂O₂ released into the buffer by the worms will react with the Amplex[®] Red reagent in a 1:1 ratio in the presence of HRP to produce the red fluorescent excitation product resorufin that can be measured fluorometrically or spectrophotometrically, and the amount of H₂O₂ released can be calculated by comparison to a standard curve. The assay is most appropriate for studies focused on released ROS, and its advantages include ease of use, the ability to use small numbers of animals in a plate reader assay in which measurements can be taken either fluorometrically or spectrophotometrically.

Keywords: Amplex, Hydrogen peroxide, *C. elegans*, ROS, NADPH oxidase, Oxidative burst, Pathogen

[Background] Interest in developing assays to study the amount and type of reactive oxygen species produced by *C. elegans* was primarily driven by the use of this model organism in studies focused on aging and infection. Oxidative damage caused by ROS has long been hypothesized to contribute to aging, and therefore measuring the intracellular oxidative environment and how it changes over time and in various aging mutants is relevant. Likewise, ROS are generated during the *C. elegans* immune response where they are capable of acting as signaling or effector molecules. Understanding the sources of ROS and being able to measure their magnitude is therefore of key importance (reviewed by [Labuschagne and Brenkman, 2013; McCallum and Garsin, 2016; Miranda-Vizuete and Veal, 2017]). Commonly used methods to measure ROS in *C. elegans* include dihydrofluorescein dyes, the

genetically encoded HyPer probe, and the Amplex[®] Red assay. Additional methods with more specific applications include MioSox, a modified dihydroethidium (DHE) that is specific for measuring ROS in the mitochondria, and the Grx1-roGFP2 genetically-encoded probe that measures changes in redox potential (Labuschagne and Brenkman, 2013). All current methodologies that measure ROS have caveats that must be carefully considered when deciding what method to use and how to interpret the results.

The dihydrofluorescein dyes, including 2,7-dichlorodihydrofluorescein diacetate (H₂DCFDA) commonly used in *C. elegans*, are colorless and membrane permeable in the reduced form, but become fluorescent and non-membrane permeable when oxidized. Therefore, these dyes can penetrate the cytoplasmic membrane and provide a decent measure of intracellular ROS (Labuschagne and Brenkman, 2013). However, to measure the fluorescence from *C. elegans* using a plate reader, a large number of animals (1,000) are required (Schulz *et al.*, 2007; Ewald *et al.*, 2017). Alternatively, animals can be visualized using fluorescence microscopy (Sem and Rhen, 2012). Two disadvantages to H₂DCFDA include the fact that it is not specific to H₂O₂, and it may exert free radical properties, contributing to the generation of ROS in addition to its measurement (Labuschagne and Brenkman, 2013; Dikalov and Harrison, 2014).

HyPer, in contrast, is considered very specific to H₂O₂. It consists of a cpYFP (circularly permuted yellow fluorescent protein) inserted into the regulatory domain of OxyR-RD and placed under the control of a *C. elegans* promoter. Both tissue-specific promoters and non-specific promoters have been utilized in the worm depending on the application (Back *et al.*, 2012; Knoefler *et al.*, 2012). The sensor works based on the formation of an intramolecular disulfide bond that occurs upon interaction with H₂O₂. The strength of the two fluorescent emission peaks released by HyPer changes upon formation of the disulfide bond; one decreases and the other increases in a proportionate manner (Belousov *et al.*, 2006). The advantage of such a ratiometric approach is that quantification can be done independently of the expression levels of the protein and therefore no normalization is required (Labuschagne and Brenkman, 2013). An additional advantage is the ability to monitor hydrogen peroxide levels over time and observe both increases and decreases. Disadvantages include the fact that HyPer is sensitive to pH changes and overexpression of these probes can significantly affect H₂O₂ redox signaling of the cell potentially confounding results. Furthermore, measuring the emission peaks requires sophisticated microscopy of individual animals or large numbers of animals (1,000/sample) to obtain accurate results from a plate reader (Labuschagne and Brenkman, 2013; Dikalov and Harrison, 2014; Ewald *et al.*, 2017).

The Amplex[®] Red assay is based on the oxidation of 10-acetyl-3,7-dihydroxyphenoxazine, which is catalyzed by horseradish peroxidase (HRP) in the presence of H₂O₂ to produce a red fluorescent oxidation product, resorufin, at a 1:1 ratio (Figure 1). Resorufin can be measured both fluorometrically and spectrophotometrically, though the fluorescent measurements are more sensitive and can detect lower levels of H₂O₂ (Mohanty *et al.*, 1997). The Amplex[®] Red Hydrogen Peroxide/Peroxidase Assay Kit was first adapted to *C. elegans* by our group to measure the output of ROS following exposure to microbial pathogens (Chavez *et al.*, 2007; Chavez *et al.*, 2009; Tiller and Garsin, 2014; van der Hoeven, *et al.*, 2015; Liu *et al.*, 2019). The assay measures the amount of H₂O₂ excreted by the worm into the

surrounding medium. The major source of ROS during this response appears to be an NADPH oxidase called BLI-3, which is found in the hypodermis (skin), intestine and pharynx (Edens *et al.*, 2001; Chavez *et al.*, 2007 and 2019; van der Hoeven *et al.*, 2015). BLI-3 is a member of the dual oxidase (DUOX) family of proteins, which all feature a peroxidase domain in addition to an NADPH oxidase domain that generates H₂O₂ rather than superoxide (Aguirre and Lambeth, 2010). Bioinformatic predictions and experimental evidence indicate that BLI-3 is located in the cytoplasmic membranes of the tissue in which it is produced and oriented in a manner such that H₂O₂ will be excreted into extracellular environments (Edens *et al.*, 2001; Chavez *et al.*, 2009; Aguirre and Lambeth, 2010; Moribe *et al.*, 2012; van der Hoeven *et al.*, 2015). Based on the localization of BLI-3, H₂O₂ is released from the hypodermal surface and into the lumen of the GI tract whereupon it is excreted due to the defecation cycle that the animals continuously undergo. When measuring H₂O₂ from animals exposed to pathogen, the levels can be significantly, but not completely, suppressed by the addition of an NADPH oxidase inhibitor (Chavez *et al.*, 2007 and 2009). These results suggest that there are additional sources of ROS during infection that are contributing the pool of H₂O₂ that is being measured. A likely source is the mitochondria, which are known to become dysfunctional during pathogen assault (Pellegrino *et al.*, 2014), and ROS generated from this source in the form of H₂O₂ would be able to cross membranes and contribute to the extracellular pools. Advantages to using the Amplex[®] Red assay for measuring ROS in the worm include the fact that it is specific to H₂O₂, it can be done on samples containing as few as 30 worms in a plate reader, and it does not require sophisticated genetics or microscopy (Labuschagne and Brenkman, 2013; Tiller and Garsin, 2014). Disadvantages include the fact that there is no tissue specificity and one must be careful to maintain samples in the dark, as resorufin is light sensitive (Labuschagne and Brenkman, 2013). Herein, we present a detailed protocol of our adaption of the Amplex[®] Red Hydrogen Peroxide/Peroxidase Assay Kit to measuring H₂O₂ output from the worm following pathogen exposure.




Figure 1. Conversion of Amplex Red to Resorufin in the presence of H₂O₂ and horseradish peroxidase (HRP)

Materials and Reagents

1. 96-well microplate, black, clear bottom with lid (Corning costar, catalog number: 3603)
2. 60 mm Petri dish (Fisher Scientific, catalog number: FB0875713A)
3. 1.7 ml microcentrifuge tube
4. 15 ml conical tube
5. Micropipette tips
6. PluriStrainer 10 μ m (pluriSelect, catalog number: 43-50010-01)
7. Serological pipettes
8. *Caenorhabditis elegans* wild-type N2 strain (*Caenorhabditis* Genetics Center, University of Minnesota)
9. *Escherichia coli* OP50 strain (*Caenorhabditis* Genetics Center, University of Minnesota)
10. *Enterococcus faecalis* OG1RF strain (ATCC 47077)
11. Autoclaved distilled water
12. Cholesterol (Sigma-Aldrich, catalog number: C-8503)
13. Absolute ethanol (Fisher, catalog number: BP2818-4)
14. Sodium chloride (NaCl) (Fisher, catalog number: S271-500)
15. Peptone (Becton, Dickinson and Company, catalog number: 211677)
16. Agar (Becton, Dickinson and Company, catalog number: 214530)
17. Calcium chloride (CaCl₂) (EMD, catalog number: 3000)
18. Magnesium sulfate (MgSO₄) (EMD, catalog number: MX0070-1)
19. Potassium phosphate, dibasic (K₂HPO₄) (Sigma-Aldrich, catalog number: P3786)
20. Potassium phosphate, monobasic (KH₂PO₄) (Sigma-Aldrich, catalog number: 795488)
21. Sodium phosphate dibasic (Na₂HPO₄) (Sigma-Aldrich, catalog number: S3264)
22. Nystatin (EMD Millipore, catalog number: 475914-1GM)
23. Gentamicin sulfate (Sigma-Aldrich, catalog number: G3632-5G)
24. Horse radish peroxidase (part of kit: Invitrogen, catalog number: A22188)
25. Dimethyl sulfoxide
26. Amplex[®] Red hydrogen peroxide/peroxidase assay kit (Invitrogen, catalog number: A22188)
27. Hydrogen peroxide (VWR, catalog number: VW3540-2)
28. Brain heart infusion powder (Becton, Dickinson and Company, catalog number: 211059)
29. Luria-Bertani powder (Becton, Dickinson and Company, catalog number: 244620)
30. Brain heart infusion (BHI) agar plates with gentamicin (50 μ g/ml), made from BHI powder with agar added. Autoclaved, cooled, gentamycin added before pouring
31. Brain heart infusion (BHI) broth, made from BHI powder and autoclaved.
32. Nematode growth medium agar plates (see Recipes)
33. M9W buffer (see Recipes)
34. 50 mg/ml gentamicin sulfate (see Recipes)
35. 10 mg/ml nystatin (see Recipes)

36. 5 mg/ml cholesterol (see Recipes)
37. Potassium Phosphate Buffer (1 M, pH 6.0) (see Recipes)

Equipment

1. Micropipettes
2. Fluorescence and absorbance microplate reader (Biotek, Cytation 5 multimode reader)
3. Incubators for stable temperatures; 20 °C, 25 °C, 37 °C
4. Tabletop centrifuge
5. Stereo microscope (Zeiss, Stemi 2000)
6. Vortex mixer
7. Autoclave

Software

1. Gen5 (Biotek, <https://www.biotek.com/products/software-robotics-software/gen5-microplate-reader-and-imager-software/>)
2. Microsoft Excel (Microsoft, <https://products.office.com/en-us/excel>)
3. GraphPad Prism (GraphPad, <https://www.graphpad.com/scientific-software/prism/>)

Procedure

A. Worm synchronization and growth

1. Transfer 10-15 N2 wild-type adult worms from a crowded non-starved NGM plate to a fresh 100 mm NGM plate.
2. Let the worms lay eggs overnight at 20 °C.
3. Next day, collect all the worms in the plate in 5 ml of M9W buffer.
4. Filter the worm suspension through a 10 µm pluriStrainer into a 50 ml conical tube. Only L1-stage worms can pass through the filter.
5. Count the number of L1-stage worms that pass through the filter.
6. Put approximately 500 of L1-staged worms on a fresh 100 mm NGM agar plate with an OP50 bacterial lawn.
7. Allow the worms to develop to the late L3 to early L4 stage at 20 °C, about 36 h, before starting the next step at which you will need worms at the young-adult stage.

B. Preparation of *E. faecalis* plates and control *E. coli* plates for worm exposure

1. The night before the worms reach the proper stage (young-adult), start an overnight culture of *E. faecalis* OG1RF in 2 ml of BHI broth and *E. coli* OP50 in 2 ml of LB broth at 37 °C with shaking.

2. Add 400 μ l of the *E. faecalis* culture to 100 mm BHI agar plates containing 50 μ g/ml of gentamycin and spread to create a lawn. Gentamycin prevents the growth of any *E. coli* transferred over with the worms, as *E. faecalis* is inherently resistant. Grow at 37 °C for 4-6 h. Additionally, add 400 μ l of the *E. coli* culture to 100 mm NGM agar plates and grow at 37 °C overnight. Make one plate each for every condition you wish to test.
3. Remove the *E. faecalis* plates from the incubator and allow them to cool to room temperature.
4. Wash the worms (young-adult stage) off from the NGM agar plate by adding 5 ml of M9W buffer. Transfer the suspended worms into 15 ml conical tubes. Wash the remaining worms off the plates using extra M9W buffer, if necessary.
5. Either let the worms settle to the bottom of the tube or spin them down for one minute at 200 x *g*. Aspirate the liquid, being careful not to disturb the collected worms. Wash the worms twice with M9W buffer. Transfer worms to *E. faecalis* (OG1RF) or control *E. coli* (OP50) plates.
6. Let the plates incubate overnight at 25 °C (about 16 h). This is the standard temperature used when exposing *C. elegans* to pathogenic bacteria.

C. Preparation of reagents for Amplex Red Assay (enough for 20 reactions)

1. Dilute the 5x reaction buffer that comes with the Amplex[®] Red assay kit to make a total of 25 ml of working buffer.
2. Thaw a tube of the horse radish peroxidase (HRP) and aliquot 100 μ l into 1.7 ml microcentrifuge tubes. One tube should then be aliquoted further into 20 μ l volumes, also in microcentrifuge tubes. Keep one 20 μ l volume out for immediate use and freeze the rest at -20 °C. Avoid repeated freezing and thawing.
3. Using the dimethyl sulfoxide (DMSO) that comes with the kit, add 60 μ l to the Amplex[®] Red reagent provided in the kit (154 μ g) and dissolve. Aliquot the dissolved reagent into 11 μ l volumes in 1.7 ml microcentrifuge tubes. Leave one tube out for immediate use and freeze the rest at -20 °C.
4. Right before you are ready to begin the assay, add the reagents together with the following volumes: 970 μ l 1x reaction buffer, 20 μ l HRP, 10 μ l Amplex[®] Red. Mix well and keep in the dark until use.

D. Amplex[®] Red assay for *C. elegans*

1. Wash the worms off the *E. faecalis* and control plates and use 1x reaction buffer to transfer.
2. Collect the worms in 15 ml tubes suspended in reaction buffer.
3. Let the worms settle to the bottom of the tube or spin them down for one minute at 200 x *g*. Aspirate the liquid, being careful not to disturb the collected worms.
4. Wash the worms 3-4 times in 1 ml of reaction buffer to get rid of the bacteria.
5. For each sample, transfer the 1 ml of worms and reaction buffer to a 60 mm empty Petri plate. Cut the end of a pipette tip to make the opening larger and, looking under the stereomicroscope, suck up about 15 worms in 25 μ l of reaction buffer and transfer to a well of a 96-well plate. Do

this twice to obtain 30 +/- 2 worms per 50 µl of reaction buffer/well. Count animals and adjust if necessary. Repeat twice more to generate technical triplicates of the sample.

6. Prepare dilutions in 1x reaction buffer of H₂O₂ for the standard curve as described in the Amplex[®] Red kit literature. We use 0 µM, 1 µM, 2 µM, 3 µM, 4 µM and 5 µM as final concentrations. Dispense 50 µl of 2x concentrated H₂O₂ standards to 96-well plate in duplicate. Desired final concentrations of H₂O₂ standards will be achieved after adding the same volume of Amplex[®] red reagent mix. It is recommended that the H₂O₂-standard wells should be as far as possible from the sample wells in order to prevent interference of volatile H₂O₂ from high-concentration standard wells into the sample wells.
7. Mix the Amplex[®] Red reagent, HRP and buffer (Step C4) and add 50 µl to each well of the 96-well plate with the worm samples, controls and standards. Record the time assay begins.

E. Reading/recording assay

Using a microplate reader, read the fluorescence using 540 +/-20 nm excitation and 590 +/-20 nm emission, bottom optics position and “extended” gain. We use a Biotek Cytation 5 microplate reader with Gen5 software to obtain the data. Measurements are expressed as **Relative Fluorescence Units (RFU)**. See Figure 2 for plate organization and fluorescence readings of a sample experiment designed to test the difference of H₂O₂ production by worms in the presence of non-pathogenic (*E. coli* OP50 strain) and pathogenic (*E. faecalis* OG1RF strain) bacteria as food. Take time points at 30, 60, 120 and 180 min. Generally, the 60-min reading demonstrates the best dynamic range.

	1	2	3	4	5	6	7	8	9	10	11	12
A	7633	7603	14391	14868	35386	35925	101726	101705	74698	78419	102865	100803
B	0uM	0uM	1uM	1uM	2uM	2uM	3uM	3uM	4uM	4uM	5uM	5uM
H ₂ O ₂ standards												
C												
D	<i>E. coli</i> OP50			<i>E. faecalis</i> OG1RF								
	Rep-1	Rep -2	Rep-3	Rep-1	Rep -2	Rep-3						
E	11321	12312	11851		16036	15981	16367	← Samples with worms				
F												
G												
H												

Figure 2. Sample plate organization and fluorescence readings. Row “A” has duplicates of increasing concentrations (0-5 µM) of H₂O₂ for creating the standard curve. Row “E” has technical triplicates of the samples with the worms exposed to *E. coli* OP50 and *E. faecalis* OG1RF bacteria. Note that the samples are placed in the wells relatively far from the H₂O₂ standards. The reading was taken 60 min after mixing the samples with the Amplex Red reagent.

Data Analysis

A. Calculation of the amount of H₂O₂ produced

1. Calculate the **refined RFU (rRFU)** of the samples by subtracting the measurement of the 0 μM H₂O₂ standard (the average of the technical duplicates) from the readings. The results represent the RFU generated by the worms in the samples during the measured time period.
2. Calculate the **rRFU** of the standards by subtracting the measurement of 0 μM H₂O₂ standard (background fluorescence of the reaction buffer) from the measurements of all of the standards.
3. Plot the rRFU data of the standards using GraphPad Prism or similar software to create a linear standard curve. See Figure 3 for a sample standard curve created using GraphPad Prism.
4. Calculate the **slope (m)** and the **y-intercept (b)** of the standard curve equation (rRFU = m[H₂O₂] + b). GraphPad Prism software automatically calculates the values.
5. Use the formula below to calculate the amount of H₂O₂ produced per minute for each sample.
H₂O₂ production:

$$\left(\frac{pmoles}{minutes}\right) = \frac{(\text{refined RFU} - \text{y intercept}) \times 100}{\text{Slope} \times \text{Time (min)}}$$

Note: Below is the detailed derivation of the formula above from H₂O₂ concentration.

H₂O₂ concentration in each well:

$$[\text{H}_2\text{O}_2] \text{ micro molar (}\mu\text{M)} = \frac{(\text{rRFU} - \text{y intercept})}{\text{Slope}}$$

Total H₂O₂ amount produced in each well with 100 μl volume:

$$\text{micro moles} = \text{volume (liters)} \times [\text{H}_2\text{O}_2]$$

$$\text{micro moles} = 100 \times 10^{-6} \times [\text{H}_2\text{O}_2]$$

$$\text{pico moles} = 100 \times 10^{-6} \times [\text{H}_2\text{O}_2] \times 10^6$$

$$\text{pico moles} = 100 \times [\text{H}_2\text{O}_2]$$

H₂O₂ production rate

$$\left(\frac{pmoles}{minutes}\right) = \frac{100 \times [\text{H}_2\text{O}_2]}{\text{time (minutes)}}$$

H₂O₂ production rate

$$\left(\frac{pmoles}{minutes}\right) = \frac{100 \times (\text{rRFU} - \text{y intercept})}{\text{Slope} \times \text{Time (min)}}$$


A

X	Group A	
H ₂ O ₂ Concentration(μM)	Refined RFUs	
X	A:Y1	A:Y2
0	0	0
1	6758	7265
2	27753	28322
3	44093	49403
4	67065	70816
5	95232	93200

B

Interpolation		A
		Refined RFUs
		Y
1	Line	
2	Best-fit values	
3	YIntercept	-7430
4	Slope	19302
5	95% CI (asymptotic)	
6	YIntercept	-13494 to -1366
7	Slope	17299 to 21305
8	Goodness of Fit	
9	Degrees of Freedom	10

C



$$y = mx + b$$

$$rRFU = 19302 [H_2O_2] - 7430$$

Figure 3. Sample standard curve and equation. A. Refined relative fluorescence units (rRFU) were plotted against concentration. B and C. Best-fit linear curve (B) and equation (C) were generated by Graphpad Prism.

B. Statistical analysis of the data

- Each sample has three technical replicates. GraphPad Prism or similar software is used to analyze the data. Create a bar graph by averaging the H₂O₂ production ($\frac{pmoles}{minutes}$) for each sample. Error bars represent the standard error of mean (SEM) for the three replicates. See Figure 4 for a sample bar graph of H₂O₂ production in the presence of non-pathogenic and pathogenic bacteria; *E. coli* (OP50 strain) and *E. faecalis* (OG1RF), respectively.
- The unpaired two-tailed *t*-test is used to test the statistical significance of the data.
- To verify results with independent biological replicates, the experiment should be replicated a total of at least three times on different days with independent batches of worms. In our experience, the exact amount of ROS produced by a given batch of worms on a given day differs depending on slight differences in developmental timing and other variables. Therefore, it is not feasible to average experimental values obtained on different days. However, the PATTERNS are consistent. For example, in this experiment the animals exposed to *E. faecalis* ALWAYS produce significantly more ROS than animals exposed to non-pathogenic *E. coli*. For publication purposes, a representative experiment can be shown as primary data with biological replicates added to a supplement, if necessary.




Figure 4. Sample bar graph showing the rate of H₂O₂ production. By each sample as picomoles per minute. Error bars represent the standard error of mean (SEM) from technical triplicates of each sample ($P = 0.0002$).

Recipes

1. 50 mg/ml gentamicin sulfate
50 mg gentamicin sulfate is dissolved in 1 ml of distilled water and stored at -20 °C
2. 10 mg/ml nystatin
10 mg nystatin is suspended in 1 ml 100% ethanol and stored at 4 °C
3. 5 mg/ml cholesterol
5 mg cholesterol is dissolved in 1 ml of 100% ethanol and stored at room temperature
4. Nematode growth medium (NGM) agar plates (per liter)
 - a. Autoclave and cool to 50 °C
 - 3 g NaCl
 - 2.5 g peptone
 - 18 g agar
 - 975 ml distilled water
 - b. Add following while stirring with the same order:
 - 25 ml of potassium phosphate buffer (1 M, pH 6.0)
 - 1 ml of MgSO₄ (1 M)
 - 1 ml of CaCl₂ (1 M)
 - 1 ml of nystatin (10 mg/ml)
 - 1 ml of cholesterol (5 mg/ml)
 - c. Pour 10 ml or 25 ml of NGM into 60 mm or 100 mm Petri dish, respectively. Let plates solidify overnight

5. M9W buffer
 - 3 g of KH_2PO_4
 - 6 g of Na_2HPO_4
 - 5 g of NaCl
 - 1 L distilled water
 - Autoclave
 - Add 1 ml of sterile MgSO_4 (1 M)
6. Potassium Phosphate Buffer (1 M, pH 6.0) per liter
 - 30.2 g of K_2HPO_4
 - 118.2 g of KH_2PO_4
 - Complete volume to 1 L with distilled water
 - Autoclave and dispense in 50 ml conical tubes for future use

Acknowledgments

This work was supported by the National Institute of Allergy and Infectious Diseases of the National Institutes of Health under award number R01AI076406 to D. A. G. The content is solely the responsibility of the authors and does not represent the official views of the National Institutes of Health. Also acknowledged, are the previous publications from on which this protocol was based (Chavez *et al.*, 2007 and 2009; Liu *et al.*, 2019; Tiller and Garsin, 2014; van der Hoeven *et al.*, 2015).

Competing Interests

The authors declare no competing interests.

References

1. Aguirre, J. and Lambeth, J. D. (2010). [Nox enzymes from fungus to fly to fish and what they tell us about Nox function in mammals](#). *Free Radic Biol Med* 49(9): 1342-1353.
2. Back, P., De Vos, W. H., Depuydt, G. G., Matthijssens, F., Vanfleteren, J. R. and Braeckman, B. P. (2012). [Exploring real-time *in vivo* redox biology of developing and aging *Caenorhabditis elegans*](#). *Free Radic Biol Med* 52(5): 850-859.
3. Belousov, V. V., Fradkov, A. F., Lukyanov, K. A., Staroverov, D. B., Shakhbazov, K. S., Terskikh, A. V. and Lukyanov, S. (2006). [Genetically encoded fluorescent indicator for intracellular hydrogen peroxide](#). *Nat Methods* 3(4): 281-286.
4. Chavez, V., Mohri-Shiomi, A. and Garsin, D. A. (2009). [Ce-Duox1/BLI-3 generates reactive oxygen species as a protective innate immune mechanism in *Caenorhabditis elegans*](#). *Infect Immun* 77(11): 4983-4989.

5. Chavez, V., Mohri-Shiomi, A., Maadani, A., Vega, L. A. and Garsin, D. A. (2007). [Oxidative stress enzymes are required for DAF-16-mediated immunity due to generation of reactive oxygen species by *Caenorhabditis elegans*](#). *Genetics* 176(3): 1567-1577.
6. Dikalov, S. I. and Harrison, D. G. (2014). [Methods for detection of mitochondrial and cellular reactive oxygen species](#). *Antioxid Redox Signal* 20(2): 372-382.
7. Edens, W. A., Sharling, L., Cheng, G., Shapira, R., Kinkade, J. M., Lee, T., Edens, H. A., Tang, X., Sullards, C., Flaherty, D. B., Benian, G. M. and Lambeth, J. D. (2001). [Tyrosine cross-linking of extracellular matrix is catalyzed by Duox, a multidomain oxidase/peroxidase with homology to the phagocyte oxidase subunit gp91phox](#). *J Cell Biol* 154(4): 879-891.
8. Ewald, C. Y., Hourihan, J. M., Bland, M. S., Obieglo, C., Katic, I., Moronetti Mazzeo, L. E., Alcedo, J., Blackwell, T. K. and Hynes, N. E. (2017). [NADPH oxidase-mediated redox signaling promotes oxidative stress resistance and longevity through memo-1 in *C. elegans*](#). *Elife* 6: e19493.
9. Knoefler, D., Thamsen, M., Konieczek, M., Niemuth, N. J., Diederich, A. K. and Jakob, U. (2012). [Quantitative *in vivo* redox sensors uncover oxidative stress as an early event in life](#). *Mol Cell* 47(5): 767-776.
10. Labuschagne, C. F. and Brenkman, A. B. (2013). [Current methods in quantifying ROS and oxidative damage in *Caenorhabditis elegans* and other model organism of aging](#). *Ageing Res Rev* 12(4): 918-930.
11. Liu, Y., Kaval, K. G., van Hoof, A. and Garsin, D. A. (2019). [Heme peroxidase HPX-2 protects *Caenorhabditis elegans* from pathogens](#). *PLoS Genet* 15(1): e1007944.
12. McCallum, K. C. and Garsin, D. A. (2016). [The role of reactive oxygen species in modulating the *Caenorhabditis elegans* immune response](#). *PLoS Pathog* 12(11): e1005923.
13. Miranda-Vizuete, A. and Veal, E. A. (2017). [Caenorhabditis elegans as a model for understanding ROS function in physiology and disease](#). *Redox Biol* 11: 708-714.
14. Mohanty, J. G., Jaffe, J. S., Schulman, E. S. and Raible, D. G. (1997). [A highly sensitive fluorescent micro-assay of H₂O₂ release from activated human leukocytes using a dihydroxyphenoxazine derivative](#). *J Immunol Methods* 202(2): 133-141.
15. Moribe, H., Konakawa, R., Koga, D., Ushiki, T., Nakamura, K. and Mekada, E. (2012). [Tetraspanin is required for generation of reactive oxygen species by the dual oxidase system in *Caenorhabditis elegans*](#). *PLoS Genet* 8(9): e1002957.
16. Pellegrino, M. W., Nargund, A. M., Kirienko, N. V., Gillis, R., Fiorese, C. J. and Haynes, C. M. (2014). [Mitochondrial UPR-regulated innate immunity provides resistance to pathogen infection](#). *Nature* 516(7531): 414-417.
17. Schulz, T. J., Zarse, K., Voigt, A., Urban, N., Birringer, M. and Ristow, M. (2007). [Glucose restriction extends *Caenorhabditis elegans* life span by inducing mitochondrial respiration and increasing oxidative stress](#). *Cell Metab* 6(4): 280-293.
18. Sem, X. and Rhen, M. (2012). [Pathogenicity of *Salmonella enterica* in *Caenorhabditis elegans* relies on disseminated oxidative stress in the infected host](#). *PLoS One* 7(9): e45417.

19. Tiller, G. R. and Garsin, D. A. (2014). [The SKPO-1 peroxidase functions in the hypodermis to protect *Caenorhabditis elegans* from bacterial infection](#). *Genetics* 197(2): 515-526.
20. van der Hoeven, R., Cruz, M. R., Chavez, V. and Garsin, D. A. (2015). [Localization of the dual oxidase BLI-3 and characterization of its NADPH oxidase domain during infection of *Caenorhabditis elegans*](#). *PLoS One* 10(4): e0124091.