

Patients with Common Cold Coronaviruses Tested Negative for IgG Antibody to SARS-CoV-2

Stephen M. Brecher,^{a,b} Jonathan Dryjowicz-Burek,^a Hongbo Yu,^{a,c} Sheldon Campbell,^{d,e} Nora Ratcliffe,^{f,g} Kalpana Gupta^{h,i}

^aDepartment of Pathology and Laboratory Medicine, VA Boston Healthcare System, West Roxbury, Massachusetts, USA

^bDepartment of Pathology and Laboratory Medicine, Boston University School of Medicine, Boston, Massachusetts, USA

^cDepartment of Pathology, Harvard Medical School, Boston, Massachusetts, USA

Journal of

MICROBIOLOGY Clinical Microbiology®

AMERICAN SOCIETY FOR

^dDepartment of Pathology and Laboratory Medicine, West Haven VA Healthcare System, West Haven, Connecticut, USA

^eDepartment of Laboratory Medicine, Yale School of Medicine, New Haven, Connecticut, USA

Department of Pathology and Laboratory Medicine, White River Junction VA Healthcare System, White River Junction, Vermont, USA

9Department of Pathology and Laboratory Medicine, Dartmouth School of Medicine, Hanover, New Hampshire, USA

^hDepartment of Medicine, VA Boston Healthcare System, West Roxbury, Massachusetts, USA

ⁱDepartment of Medicine, Boston University School of Medicine, Boston, Massachusetts, USA

KEYWORDS Abbott IgG, cross-reacting antibody, SARS-CoV-2, coronavirus

he need for accurate antibody testing in patients following symptomatic or asymptomatic infections with severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is well documented. How we best utilize the data obtained from these antibody studies is still a haunting question (1). There are now over 160 serologic kits on the market for the detection of antibodies to this virus (2). Most are not Emergency Use Authorization (EUA)/FDA approved, and due to the poor performance of some of these assays, the FDA has requested more stringent data (2). In any case, a validation study is necessary prior to reporting patient results. One early issue in the validation/ evaluation of antibody tests for evidence of SARS-CoV-2 infection is the possibility of cross-reacting antibodies from the plasma of patients who had been infected with one or more of the common cold coronaviruses (coronavirus 229E, HKU1, NL63, and OC43). Antibody testing for SARS-CoV-2 in these patients could result in reduced specificity of the SARS-CoV-2 antibody assays due to false-positive results from cross-reacting antibodies (3). This is a problem for numerous reasons, especially in a low-prevalence population where there could be more false positives than true positives. It also ties in with issues outlined in the Infectious Diseases Society of America (IDSA) guidelines (4) and in a commentary in Lancet (5) on how to best utilize antibody test data, especially when there could be false-positive results, including cross-reacting antibodies to the four common cold coronaviruses.

In the course of performing a validation study, plasma from three distinct groups of patients was selected for immunoglobulin G (lgG) antibody testing. The validation samples were from patients with previous exposure to SARS-CoV-2, as determined by a positive PCR test (Xpert Xpress SARS-CoV-2; Cepheid, Sunnyvale, CA; or cobas SARS-CoV-2 assay; Roche Molecular Systems, Branchburg, NJ); those with negative SARS-CoV-2 PCR tests; and patients with previous non-SARS-CoV-2 respiratory infections. This last group included plasma from 20 patients who had positive viral respiratory panel PCRs (FilmArray respiratory panel 2; BioFire Diagnostics, LLC, Salt Lake City, UT) for one of the four common cold coronaviruses. The plasma from these 20 patients was collected more than 4 weeks after the positive PCR. This would allow enough time for the synthesis of IgG antibodies in these patients (Table 1).

Testing was performed on an Abbott Architect i2000SR (Abbott Park, IL) automated analyzer using the SARS-CoV-2 immunoglobulin G (IgG) assay. The assay is designed to

Citation Brecher SM, Dryjowicz-Burek J, Yu H, Campbell S, Ratcliffe N, Gupta K. 2020. Patients with common cold coronaviruses tested negative for IgG antibody to SARS-CoV-2. J Clin Microbiol 58:e01029-20. https://doi.org/10 .1128/JCM.01029-20.

Editor Alexander J. McAdam, Boston Children's Hospital

Copyright © 2020 American Society for Microbiology. All Rights Reserved.

Address correspondence to Stephen M. Brecher, stephen.brecher@va.gov.

Accepted manuscript posted online 27 May 2020

Published 23 July 2020

| Patient ^a | RP2 ^b PCR result | RP2 ^b date | Serum date ^e | SARS-CoV-2 lgG result ^c | SARS-CoV-2 PCR result ^d |
|----------------------|-----------------------------|------------------------|-------------------------|------------------------------------|------------------------------------|
| 1 | CV OC43 | 1/28/2019 ^e | 4/21/2020 | Negative | Negative |
| 2 | CV NL63 | 12/29/2019 | 4/22/2020 | Negative | Negative |
| 3 | CV HKU1 | 1/2/2020 | 4/4/2020 | Negative | Not tested |
| 4 | CV HKU1 | 1/11/2020 | 4/20/2020 | Negative | Negative |
| 5 | CV HKU1 | 1/12/2020 | 4/22/2020 | Negative | Not tested |
| 6 | CV NL63 | 2/7/2020 | 4/9/2020 | Negative | Negative |
| 7 | CV HKU1 | 2/11/2020 | 4/21/2020 | Negative | Not tested |
| 8 | CV 229E | 2/18/2020 | 4/20/2020 | Negative | Negative |
| 9 | CV NL63 | 3/2/2020 | 4/22/2020 | Negative | Not tested |
| 10 | CV NL63 | 3/4/2020 | 4/22/2020 | Negative | Not tested |
| 11 | CV NL63 | 3/4/2020 | 4/22/2020 | Negative | Not tested |
| 12 | CV NL63 | 3/4/2020 | 5/1/2020 | Negative | Not tested |
| 13 | CV NL63 | 3/5/2020 | 4/22/2020 | Negative | Not tested |
| 14 | CV HKU1 | 3/6/2020 | 5/4/2020 | Negative | Negative |
| 15 | CV NL63 | 3/6/2020 | 4/22/2020 | Negative | Not tested |
| 16 | CV NL63 | 3/9/2020 | 4/21/2020 | Negative | Negative |
| 17 | CV HKU1 | 3/9/2020 | 4/18/2020 | Negative | Negative |
| 18 | CV HKU1 | 3/9/2020 | 4/30/2020 | Negative | Not tested |
| 19 | CV OC43 | 3/11/2020 | 4/29/2020 | Negative | Not tested |
| 20 | CV NL63 | 3/23/2020 | 4/22/2020 | Negative | Negative |

aMales (17) and females (3); age range, 28-88.

^bRespiratory panel 2 film array; BioFire Diagnostics, LLC, Salt Lake City, UT.

^cAbbott Architect SARS-CoV-2 IgG antibody test.

^dEither Cepheid Xpert Express SARS-CoV-2 PCR or Roche cobas PCR assay (see text).

^eAll dates are mo-day-year.

detect IgG antibodies to the nucleocapsid protein of SARS-CoV-2. The antibody binds to SARS-CoV-2 antigen-coated microparticles and undergoes a chemiluminescent reaction, producing a direct relationship between the amount of IgG and relative light units. The presence of antibody is determined by comparing the relative light units in the reaction to the relative light units in the calibrator. The presence of antibody above the quantitative cutoff of 1.4 (index sample calibrator) is interpreted as positive.

All 20 patients tested negative for IgG antibody to SARS-CoV-2 (Table 1). Although the sample size was minimal, these data are reassuring that at least for the Abbott Architect SARS-CoV-2 antibody test, plasma from patients with documented positive PCRs for these four common cold coronaviruses did not test positive for the SARS-CoV-2 IgG antibody. This small study does not rule out that possibility. It does provide data that in our small study population, cross-reacting antibodies were not detected. Conclusions are limited by the small sample size of a predominantly elderly male population, consistent with the veteran population we studied. However, this multisite study, including data from 3 regional Veterans Affairs (VA) institutions (MA, CT, and VT) suggests that cross-reacting antibodies are not detected when testing for SARS-CoV-2 IgG antibody.

REFERENCES

- Theel ES, Slev P, Wheeler S, Coururier MR, Wong SJ, Kadkhora K. 2020. The role of antibody testing for SARS-CoV-2: is there one? J Clin Microbiol https://doi.org/10.1128/JCM.00797-20.
- Burton TM. 2020. FDA sets standards for coronavirus antibody tests in crackdown of fraud. https://www.wsj.com/articles/fda-sets-standards-for -coronavirus-antibody-tests-in-crackdown-on-fraud-11588605373.
- 3. Gorse GJ, Patel GB, Vitale JN, O'Connor TZ. 2010. Prevalence of antibodies

to human coronaviruses is lower in nasal secretions than in serum. Clin Vaccine Immunol 17:1875–1880. https://doi.org/10.1128/CVI.00278-10.

- IDSA. 2020. IDSA COVID-19 antibody primer. https://www.idsociety.org/ news--publications-new/articles/2020/emphasizing-need-for-more -information-idsa-releases-antibody-testing-primer2/.
- 5. Petherick A. 2020. Developing antibody tests for SARS-CoV-2. Lancet 395:1101–1102. https://doi.org/10.1016/S0140-6736(20)30788-1.