REVIEW ARTICLE

Update on Current Technologies for Deep Brain Stimulation in Parkinson's Disease

Michelle Paff, Aaron Loh, Can Sarica, Andres M. Lozano, Alfonso Fasano^{2,3,4}

ABSTRACT

Deep brain stimulation (DBS) is becoming increasingly central in the treatment of patients with Parkinson's disease and other movement disorders. Recent developments in DBS lead and implantable pulse generator design provide increased flexibility for programming, potentially improving the therapeutic benefit of stimulation. Directional DBS leads may increase the therapeutic window of stimulation by providing a means of avoiding current spread to structures that might give rise to stimulation-related side effects. Similarly, control of current to individual contacts on a DBS lead allows for shaping of the electric field produced between multiple active contacts. The following review aims to describe the recent developments in DBS system technology and the features of each commercially available DBS system. The advantages of each system are reviewed, and general considerations for choosing the most appropriate system are discussed.

Key Words Deep brain stimulation; Parkinson disease; Surgery; Technology; Therapy.

The development of deep brain stimulation (DBS) over 30 years ago has provided a means of relieving the conditions of patients with a variety of movement disorders and has since become increasingly important in the treatment of Parkinson's disease (PD).^{1,2} The two most common targets for DBS for PD are the globus pallidus internus (GPi) and the subthalamic nucleus (STN). While both targets are effective for reducing levodopa-induced dyskinesias and improving mobility, STN-DBS has the added benefit of reducing the daily levodopa requirement; on the other hand, the GPi is considered a safer target from a psychiatric and cognitive perspective.³⁻⁵

DBS surgery has been shown to be effective for improving the motor symptoms and quality of life of PD patients; however, there are some symptoms that may be resistant or refractory to stimulation. For example, gait, speech and postural dysfunction tend to deteriorate independently and at a faster rate than appendicular symptoms.6 STN-DBS in particular is known to exacerbate stuttering and hypophonia in PD patients. Additionally, stimulation at contacts residing within specific regions of the STN may have differential effects on axial and appendicular symptoms. For example, stimulation of the dorsal half of the lateral STN been shown to improve step velocity, length and balance compared to stimulation of ventral regions.8 The potential of STN-DBS to adversely affect axial symptoms is exemplified by a case report in which de novo freezing of gait (FOG) was produced with stimulation of the anteromedial region of the STN.9 By contrast, GPi DBS is associated with fewer axial adverse events although this observation mainly comes from uncontrolled open studies.⁶ These findings demonstrate potential challenges that may arise when optimizing stimulation settings for individual patients.

The therapeutic potential of DBS relies heavily on ideal placement of the DBS electrode within the target nucleus. While mul-

Received: April 26, 2020 Revised: June 22, 2020 Accepted: July 5, 2020

Corresponding author: Alfonso Fasano MD, PhD Edmond J. Safra Program in Parkinson's Disease, Morton and Gloria Shulman Movement Disorders Centre, Toronto Western Hospital, 399 Bathurst St, 7McL410, Toronto, ON, M5T 2S8 Canada / Tel: +1-416-603-5800 / Fax: +1-416-603-5004 / E-mail: alfonso.fasano@uhn.ca

© This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (https://creativecommons.org/licenses/by-nc/4.0) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

¹Division of Neurosurgery, University Health Network, University of Toronto, Toronto, Canada

²Edmond J. Safra Program in Parkinson's Disease, Morton and Gloria Shulman Movement Disorders Centre, Toronto Western Hospital, UHN, Division of Neurology, University of Toronto, Toronto, Canada

³Krembil Brain Institute, Toronto, Canada

⁴Center for Advancing Neurotechnological Innovation to Application (CRANIA), Toronto, Canada



tiple localization techniques, such as image guidance, microelectrode recordings, and clinical evaluation, may be implemented during the surgical procedure, small deviations from the planned trajectory can result in a reduced therapeutic window for stimulation. The STN is a small, lens-shaped nucleus measuring approximately 8 mm in length and 6 mm in width that is divided into dorsolateral sensorimotor, associative and medioventral limbic regions. Studies have shown that optimal implantation of the electrode into the dorsolateral region of the STN is necessary to produce improvements in motor symptoms while avoiding stimulation-related side effects. DBS lead revisions are uncommon due to the risk of surgical complications; thus, clinicians may have to rely on programming strategies to optimize therapeutic benefit in the case of suboptimally placed DBS electrodes.

Obtaining an optimal clinical response while minimizing stimulation-related side effects can be achieved in the majority of cases by manipulating stimulation parameters, such as electrode polarity, frequency, pulse width, and voltage as well as by shaping the volume of tissue activation (VTA). For example, bipolar stimulation can narrow the shape of the VTA to avoid the fibers of the internal capsule or other eloquent structures. Furthermore, some studies have reported improvements in gait with switching to low-frequency (60 Hz) high-voltage stimulation 13,14 as opposed to standard protocols implementing high frequency stimulation (HFS). Recently, technological advances in DBS systems have given rise to additional programming strategies that allow for further flexibility in manipulating stimulation parameters. Such advancements may aid clinicians in optimizing program settings for patients with complex Parkinsonian symptoms or suboptimally positioned DBS leads. The present review will discuss new techniques for stimulation implemented by the latest generation of DBS systems with a focus on the unique features of each system.

UPDATE ON STIMULATION

The goals of effective DBS programming in order of importance are to maximize clinical benefit, avoid stimulation-related side effects, and minimize current consumption. In general, stimulation through an optimally positioned DBS electrode should produce a significant clinical benefit with relatively low energy consumption. Monopolar stimulation, in which one (or two in case of 'double monopolar') of the electrode contacts serves as a cathode against the implantable pulse generator (IPG) case, produces a radial diffusion of current around the active contact. The radius of the VTA increases with increasing amplitude (voltage or current). Conversely, in bipolar stimulation, one contact serves as the cathode while another serves as the anode. This produces a narrower VTA with the greatest field strength centered around the cathode. Bipolar stimulation may be helpful in cases where

the VTA produced by monopolar stimulation violates surrounding eloquent structures leading to side effects. Interestingly, bipolar stimulation has been found to be associated with longer battery life than double monopolar and even monopolar stimulation using a voltage-constant IPG. ¹⁶ The therapeutic benefit from STN-DBS in PD patients has been observed with stimulation frequencies greater than 50 Hz, and clinical improvements in rigidity and bradykinesia are maximized at approximately 130 Hz. ¹⁷ Conversely, stimulation at 5 Hz has been observed to worsen symptoms of rigidity and bradykinesia. ¹⁷

Constant current

Until recently, all DBS IPGs utilized single-source, constant voltage-driven stimulation. The consequence of this design is that the impedance of the DBS electrode at the contact-tissue interface determines the magnitude of current flow to the neural elements, and this impedance may change over time (e.g., after IPG replacements). Current Medtronic IPGs can provide either voltage-constant or current-constant stimulation, but the former is usually preferred as more programming options are available. Currently, the majority of manufacturers only produce current-constant IPGs.

Pulse width

Pulse width is usually kept at the lowest duration of $60~\mu s$ unless there is lack of clinical benefit with increasing amplitude. While the amount of current required to excite neural tissue decreases with increased pulse width, the amplitude can be more easily manipulated with most IPGs. Shortening the pulse width $<60~\mu s$ has shown some utility in increasing the therapeutic window of stimulation, namely, the range of amplitudes that produce clinical benefit without side effects. Pulse widths $<60~\mu s$ are only possible with use of newer IPGs, some of which are compatible with electrodes produced by different manufacturers. Shortening the pulse width below $60~\mu s$ was shown to reduce stimulation-related side effects and increase the therapeutic window of stimulation in two patients with essential tremor (ET) when they were switched from a Medtronic IPG to a Boston Scientific Vercise IPG. 20

Frequency

Varying stimulation frequency has recently been explored for treating gait and axial symptoms in PD that responded poorly to the traditional HFS of 130 Hz. DBS targeting the pedunclopontine nucleus is an experimental procedure previously investigated to address FOG in PD patients, and its optimal stimulation settings utilized lower frequencies than other targets, most commonly frequencies < 60 Hz.²¹ Low-frequency stimulation (LFS) was hypothesized to provide better control of axial symp-

toms, including FOG, even when delivered to the STN.13 In a small study of 13 PD patients with FOG and prominent axial symptoms, LFS was found to improve gait and reduce freezing episodes while maintaining a reduction in the Unified PD Rating Scale (UPDRS) motor score although the majority of patients had to increase their levodopa dose to mitigate worsening Parkinsonian symptoms.¹³ Similarly, LFS was found to improve gait and axial symptoms in a cohort of 14 patients, and the optimal contacts for LFS were located ventral to the contacts for optimal HFS of 130 Hz.²² Variable frequency stimulation (VFS) involves alternating HFS with LFS. Jia et al.²³ proposed that VFS may serve as an effective strategy for treating FOG and axial symptoms in PD patients while maintaining optimal treatment of appendicular symptoms. This was exemplified in a case of one PD patient with FOG who improved with VFS. Subsequently, a small case study of 4 PD patients with FOG symptoms demonstrated that VFS improved gait speed and reduced episodes of FOG.²⁴ Additionally, an acute open-label study found that high stimulation frequency of 10 kHz was reported to yield good mobility with fewer stimulation-induced paresthesias and speech disturbances than the traditional stimulation frequency of 130 Hz.²⁵

Temporal fractionation

Temporal fractionation involves using two separate programs in an alternating fashion through a single DBS electrode. The two programs may have different amplitudes, polarities, and pulse widths; however, the frequency must be the same, with a maximum allowed frequency of 125 Hz (i.e., the maximum IPG output is 250 Hz). This method of stimulation allows for different areas of the target nucleus to be stimulated simultaneously. This technology is generally known as 'interleaving stimulation', as introduced by Medtronic; a similar technology is called 'multi-stim set' in Abbott devices (see below).

Case reports have described interleaving stimulation to be effective for treating refractory dystonia²⁶ and for treating tremor predominant symptoms in a patient with comorbid PD and ET.²⁷ This technology may be particularly helpful in cases in which a DBS electrode is suboptimally placed. Miocinovic et al.²⁸ described 3 patients, each of whom had one DBS lead that was suboptimally positioned within the STN. Stimulation using various monopolar and double monopolar programs improved some symptoms but worsened others and caused stimulation-related side effects. Use of interleaving programs was effective for providing stimulation of the desired areas while avoiding surrounding structures. Similarly, Zhang et al.²⁹ reported a series of 12 patients in whom multiple trials of conventional stimulation were ineffective in controlling various symptoms. Prior to implementing interleaving stimulation, the patients' complaints included dysarthria, stimulation-related dyskinesias, gait disturbances, and incomplete control of Parkinsonian symptoms. MRI evaluation revealed suboptimal placement of the affected lead in most patients in this series; however, adequate control of symptoms was able to be achieved in all with interleaving stimulation. In addition to avoiding stimulation-related side effects, interleaving stimulation has been shown to improve dyskinesias in a subset of PD patients. Interestingly, this effect was achieved with activation of a dorsal contact residing in the zona incerta, just outside the STN.³⁰ These reports exemplify how interleaving provides greater flexibility in programming, enabling the clinician to better balance Parkinsonian symptoms and stimulation-related side effects.

Current fractionation

The concept of current fractionation involves multiple independent current sources that apply constant current through each electrode, enabling the clinician to directly control the amount of current delivered to the surrounding tissues. Control of current flow through each contact individually is termed multiple independent current control (MICC) in Boston Scientific IPGs. Utilizing this paradigm, current can be distributed in a controlled fashion between multiple contacts, allowing for the VTA to be adjusted to a desired shape to fit the target region.³¹ Current fractionation is different than temporal fractionation, which uses technologies such as interleaving (see above).

MICC is not possible on Abbott devices, which use a different—less versatile—method called coactivation. Coactivation allows for multiple contacts to be stimulated as if they were a single electrode (see below). In a modeling study, coactivation was shown to be associated with lower power consumption than MICC due to lowered impedance at the electrode-tissue interface.³²

Computational models have demonstrated that direct current steering may have the potential to stimulate desired neuronal populations even in the case of suboptimally placed DBS electrodes. In fact, based on these models, selective therapeutic stimulation may be achieved across a range of electrode locations. The utility of this method of stimulation has been demonstrated in some case series and case reports. For example, current steering to shift the stimulation toward a more proximal contact, changing the shape of the VTA to a tear-drop-shaped distribution has been used in a patient who developed dyskinesias with increasing amplitude at his most used therapeutic contact. This effectively alleviated the patient's Parkinsonian symptoms without eliciting dyskinesias. Parkinsonian symptoms without eliciting dyskinesias.

However, targeting with such accuracy would require detailed information on the location of the electrode within the target nucleus of the individual and consequently may complicate the programming process. The authors suggest that visualization software could help guide clinicians by means of patient-specific compu-



tational models during the programming process.³³ For example, the GUIDE system by Boston Scientific (Valencia, CA, USA) is a three-dimensional visual software that provides information about the location of the DBS lead based on preoperative MRI and postoperative CT images. The VTA can be modeled when the stimulation settings are selected, allowing clinicians to visualize the resulting VTA in relation to surrounding structures. Use of this programming software was shown to shorten the duration of the initial programming session by a mean of 75%.³⁵

Anodic stimulation

Traditionally, cathodic stimulation has been favored, since activation of myelinated fibers requires 3-8 times more stimulation strength when the active contact acts as an anode.36 For this reason, first-generation DBS systems were designed such that the IPG case served as the anode while the active contact served as the cathode. Theoretically, the active contact on a DBS electrode can act as either a cathode or an anode. Cathodic and anodic stimulation are believed to differentially affect neural tissue depending on the orientation of the fibers, such that passing fibers are activated by cathodal stimulation and orthogonally oriented fibers are activated by anodic stimulation.³⁷ With the recent development of second-generation DBS systems that allow for increased flexibility in stimulation parameters, the effects of anodic stimulation in comparison to cathodal stimulation have been studied.38,39 Kirsch et al.38 studied anodic stimulation in 10 PD patients and reported higher side effect thresholds and lower UP-DRS III motor scores with anodic stimulation than with cathodic stimulation. Another study involving 10 PD patients found an increased therapeutic window and higher side effect threshold with anodic stimulation.³⁹ However, anodic stimulation was not able to control tremors in two patients, and anodic stimulation is also known to result in higher battery consumption.³⁹

Semibipolar stimulation

In semibipolar stimulation, the anode is divided equally between the IPG case and another contact on the DBS electrode. Based on clinical experience, this configuration is believed to be useful for avoiding stimulation-related side effects. Soh et al.³⁹ compared the effects of semibipolar stimulation to bipolar, monopolar cathodal and monopolar anodic stimulation in a study group of 10 PD patients. Semibipolar stimulation was found to have a significantly higher side effect threshold than cathodal stimulation; furthermore, semibipolar stimulation was found to have a lower battery consumption than bipolar stimulation.³⁹

Directional leads

DBS lead design has exhibited significant growth in recent years with the development of commercially available directional leads. The traditional DBS lead structure consists of 4 contacts of 1.5 mm in length separated from one another by interspaces of 1.5 mm in the case of the Medtronic 3387 model and 0.5 mm in the case of the Medtronic 3389 model lead. As described previously, monopolar stimulation produces a spherical electric field that diffuses radially from the active contact. Directional leads, on the other hand, contain multiple contacts in a radial distribution around the shaft of the DBS lead. Each contact may be stimulated individually or in combination by means of temporal fractionation (Abbott devices) or current fractionization (MICC, Boston Scientific). Passing current through individual contacts enables current steering along a vector perpendicular to the lead. Commercially available directional leads have a circumferential ring configuration made up of 3 individual contacts. Stimulation through all 3 contacts simultaneously produces an electric field equivalent to that of monopolar stimulation through a traditional DBS lead. Importantly, the smaller size of each individual contact leads to a greater electric field density at the contact-tissue interface; therefore, less current is required to produce a therapeutic effect. This feature may potentially contribute to preserved battery life.

The advantage of the directional lead design was first proposed by computational models. 40,41 One of these models suggested that current steering through segmented leads could shift the center of the VTA as much as 1.0–1.3 mm. 42 While this distance may seem insufficient to compensate for a poorly positioned DBS lead, studies in human subjects have demonstrated that directional stimulation through a most favorable contact could increase the therapeutic window for stimulation and increase the current threshold for producing side effects. 43-45

Two initial studies performed in PD patients were intraoperative alone. 44,46 Pollo et al. 44 calculated that the average therapeutic window was 41.3% wider than that of omnidirectional stimulation, and the average therapeutic current was 43% less for the best directional stimulation than for omnidirectional stimulation. Contarino et al. 46 found similar results in a double-blind intraoperative study in eight PD patients by means of a noncommercialized 32-contact electrode. Dembek et al. 43 investigated the clinical benefit of directional stimulation in a chronic double-blind crossover study in which patients received a trial of best directional stimulation compared to omnidirectional stimulation. This group found no difference in motor outcomes between the two conditions other than a significant positive difference in hand rotation in the directional condition.

The utility of directional leads for STN-DBS was demonstrated in the PROGRESS study in which 234 PD patients were implanted with Abbott directional leads targeting the STN.⁴⁷ For the first three months, patients received traditional omnidirectional stimulation, and during the following 3-month period, pa-

tients were switched to directional stimulation. The results of this prospective, single-arm crossover study yielded a significant mean 41% increase in the therapeutic window for stimulation with directional stimulation. Additionally, the amount of current required to produce therapeutic effects was reduced by 39%. 47 The capability and limitation of directional leads to steer current in the case of a misplaced DBS lead is depicted in Figure 1.

One potential disadvantage of directional leads is the increased complexity of programming due to the addition of multiple contacts directed in different orientations. Fortunately, monopolar review can be performed in a similar fashion as with traditional lead types, and directional stimulation may be implemented only if needed to address side effects or complex symptoms. Algorithms implementing anatomical models based on MR and diffusion tensor imaging to automate programming of directional lead have been proposed. 48 However, access to such technology is likely to be a limiting factor for many community neurologists.

UPDATE ON DEVICES

DBS system manufacturers

Medtronic

Medtronic moved its headquarters from Minneapolis, MN, USA to Dublin, Ireland, after acquiring Covidien, a rival medical device company, in 2015. The Medtronic 3382 DBS lead and the ITREL IPG was the first DBS system to gain U.S. Food and Drug Administration (FDA) approval in the United States for treatment of ET and tremor-dominant PD in 1997. In 2002, this system became FDA-approved for treatment of bradykinesia and rigidity in PD. The early IPGs were single channel devices, such that bilateral DBS required implantation of 2 IPGs (e.g., Soletra was approved in South Korea in 2001). The Kinetra IPG by Medtronic, which became available in 1998, was the first dual-channel IPG (never launched in South Korea). In 2008, the Kinetra IPG was replaced by the Activa Primary Cell (PC), which was approved in South Korea in 2010. An IPG serving a single channel (Activa SC) was also subsequently produced and approved in many countries, including in South Korea in 2011. In 2015, the FDA approved Medtronic DBS systems for full-body MRI when certain safety conditions were met. In addition to the commercially available Activa IPGs, Medtronic launched a novel system named Activa PC+S for research use in 2013. This system enables the sensing and recording of select brain activity while simultaneously providing targeted DBS therapy. Following the success of this sensing-enabled system, a commercially available form named Percept PC was released in January 2020. This new IPG is using the BrainSense technology that enables physicians to track patient-specific brain signals and correlate these with patient-recorded symptoms and treatment side effects. Additionally, the Percept PC system allows for 3T full-body MRI scan has a lower pulse width range lower limit and an improved battery longevity than the Activa PC⁴⁹).

Abbott

St. Jude Medical (St. Paul, MN, USA) was the first medical device company to develop constant-current DBS systems, the Libra and Libra XP. The Libra DBS system was tested for efficacy in PD, with 136 patients implanted between October 2005 and April 2009 in a multicenter trial involving 15 centers.⁵⁰ Additionally, the Libra system was used to treat ETs in 127 patients at 12 investigational sites.⁵¹ These trials were designed to establish the safety and efficacy of this DBS system for the purposes of gaining FDA-approval and did not compare constant-current to con-




Figure 1. A: Stimulation of the dorsolateral STN through an ideally positioned DBS lead. B: Current steering through a directional DBS lead is unable to compensate for a misplaced lead that resides outside the target nucleus. C: Use of current steering through a directional lead is the compensate for a misplaced lead that resides outside the target nucleus. C: Use of current steering through a directional lead that is positioned slightly posterior to the STN is capable of steering current toward the dorsolateral STN. STN: subthalamic nucleus, DBS: deep brain stimulation.



stant-voltage devices. The Libra system was FDA-approved as an investigational device in 2012 to investigate the efficacy of subcallosal cingulate gyrus stimulation for the treatment of refractory depression in a multicenter trial.⁵² Subsequently, the rechargeable Brio DBS system was granted FDA approval to treat PD and ETs in 2015.

Despite gaining FDA approval, the Libra, Libra XP, and Brio systems were never commercialized and are not available for clinical use in the United States. They are instead used in many other countries (e.g., in South Korea, Libra was launched in 2013 and Brio in 2016), St. Jude Medical was later acquired by Abbott Laboratories (Abbott Park, IL, USA). The Infinity DBS system was the first St. Jude/Abbott DBS system to be commercially available in the United States. The Infinity system was also the first FDA-approved system to utilize a directional lead, which was granted in 2016 (same year in South Korea). This system can be programmed using an iPad Mini (Apple, Cupertino, CA, USA), and patients can connect the implant with a paired iPod Touch (Apple). Communication between the IPG device and these programmers occurs via Bluetooth, and no physical contact is necessary. More recently, Abbott received FDA approval for a software upgrade allowing patients implanted with the Infinity DBS system to receive full-body MRI.

Boston Scientific

Boston Scientific launched the Vercise DBS system in Europe in 2012. This system was evaluated in the VANTAGE trial, a multicenter, open-label, nonrandomized study conducted at 6 centers in 6 European countries.⁵³ Forty patients received bilateral STN-DBS using an eight-contact DBS lead and a multisource, constant-current, rechargeable DBS IPG. Measures of clinical improvement and the incidence of stimulation-related side effects were similar to those in previous trials of DBS for PD.53 The Vercise DBS system was granted FDA approval in December 2017 based on the results of the INTREPID study.⁵⁴ In the multicenter, double-blind, randomized, controlled trial, 160 subjects received bilateral STN-DBS with the Vercise system and then were randomized to receive 12 weeks of active or sham stimulation. The study achieved its primary endpoint of change in the number of hours spent with good mobility based on patient diaries, and the results at one year of follow up reported an improvement in the UPDRS III motor score of 49.2% compared to baseline.⁵⁵ A few years later in Europe in January of 2019, the Vercise Gevia IPG with the Cartesia directional DBS lead became FDA approved, and later that year, the FDA granted labeling of the ImageReady MRI Vercise Gevia DBS system, enabling patients implanted with the device to qualify for full-body MRI. In South Korea, the Vercise PC was approved in 2017, and Gevia was approved one year later.

PINS Medical

PINS Medical Co., Ltd. (Beijing, China) was established in 2008, and the PINS DBS system was developed as less-costly alternative to Medtronic devices. As of October 2016, 5,000 patients have been implanted with PINS DBS systems throughout China. The PINS IPG models resemble those of Medtronic, with single-channel, dual-channel, and dual-channel rechargeable options. The PINS DBS electrodes also resemble those of Medtronic with four 1.5 mm-long contacts spaced 1.5 mm or 0.5 mm apart. The PINS DBS system, utilizing a single-channel IPG, was evaluated in 40 patients in two centers in China in an open-label trial⁵⁶ in which 35 patients received bilateral and 5 patients received unilateral STN-DBS between November 2009 and December 2011. The PINS IPG was the first device to allow for VFS, in which stimulation at two different frequencies may be interleaved. VFS is proposed to be helpful in treating axial symptoms in PD patients. In a small case series of 4 PD patients with FOG, VFS was observed to improve the UPDRS III motor score by a mean of 14% more than HFS, and gait speed was increased by a mean of 45% in 3 of 4 patients.57

SceneRay

SceneRay Co., Ltd. (Suzhou, China) was initiated in 2015, and in 2019, the SceneRay DBS system obtained registration with the Chinese Food and Drug Administration. The SceneRay model 1180 IPG is a dual-channel device that is capable of remote wireless programming. The currently available electrode models contain four contacts of 1.5 mm in length and are separated from one another by 0.5, 1.0, or 1.5 mm intervals. SceneRay has also designed an electrode to specifically to target the nucleus accumbens and anterior limb of the internal capsule.⁵⁸ The lead has four contacts of 3.0 mm in length; the two ventral-most contacts are separated by 2 mm, and the more dorsal contacts are separated by 4 mm. Both SceneRay and PINS promote web-based, remote, wireless DBS programming systems in which patients may have their DBS settings adjusted at home by a clinician remotely located in a hospital or clinic. While the advantages of remote programming capabilities include increased access to neurologists and reduced travel costs for patients, the programming sessions tend to be lengthier and susceptible to network instability and communication errors.59

Commercially available DBS pulse generators

Medtronic Activa series

The Activa series is currently the most commonly implanted IPG. The Activa PC is smaller in size (volume of 37 cm³) and was proposed to have a longer battery life; however, studies have found a significant shorter life span for the Activa PC than for

the Kinetra. 60,61 Subsequent additions to the Activa series included the rechargeable dual-channel Activa RC (22 cm³, approved in South Korea in 2010) and the single-channel nonrechargeable Activa SC (27 cm³). These IPGs can be used as either constantvoltage or constant-current devices. These devices are capable of interleaving stimulation, and patients who are implanted with these devices may receive full-body MRI when certain conditions are met ('conditionally' safe). The Activa PC and RC have a frequency range of 2-250 Hz, while the Activa SC has a frequency range of 3-250 Hz. Each device is capable of pulse widths ranging from 60-450 µs. These Medtronic devices are FDA approved for bilateral STN and GPi stimulation for PD, unilateral thalamic stimulation for ETs, and unilateral or bilateral stimulation of the GPi or STN for treatment of chronic, drug-refractory segmental or generalized dystonia. The same indications are approved in South Korea.

Abbott Infinity

The Abbott Infinity 5 and Infinity 7 IPGs are currently commercially available in the United States. Both models are constant-current, dual-channel, nonrechargeable devices. The infinity 5 IPG has a volume of 30.4 cm³ and has been designed for centers that utilize a single lead, similar to the Medtronic Activa SC. The Infinity 7 is a larger IPG designed for bilateral stimulation, with a volume of 38.6 cm³ and a predicted life span of 5 years. These devices have a frequency range of 2-240 Hz and pulse width range of 20-500 μs. Temporal and current fractionation is possible by means of multi-stim set and coactivation, respectively. Multi-stim set involves alternating stimulation of two separate contacts, and coactivation involves current fractionation to two separate contacts. Patients implanted with Infinity DBS systems can receive a full-body MRI when certain conditions are met. Currently, the Infinity system is FDA approved for bilateral STN and GPi stimulation for PD and for bilateral thalamic stimulation for ETs. In South Korea, these devices are also approved for dystonia.

Boston Scientific Vercise

Current commercially available IPGs include the Vercise PC, Vercise RC and Gevia. The Vercise PC is a constant-current, dual-channel, nonrechargeable device with a volume of 33 cm³. The Vercise RC is a rechargeable IPG with a volume of 22.7 cm³. The Vercise PC and RC are not MRI safe. The Gevia IPG is a dual-channel, rechargeable device with a volume of 19.8 cm³ that is MRI conditional. These devices have a frequency range of 2–255 Hz and a pulse width range of 20–450 µs. Temporal fractionation is accomplished using different 'areas' (2 per channel), and current fractionation may be accomplished via MICC independently per area. These systems are currently FDA approved for bilateral

STN stimulation for PD. In South Korea, these devices are approved for PD, dystonia and medication-refractory tremors.

PINS Medical

The available IPGs manufactured by PINS Medical include the G102, G102R, and G102A. The G102 is a constant-voltage, nonrechargeable, dual-channel IPG similar to the Medtronic Activa PC with a volume of $51.2~\rm cm^3$. The G102R is a dual-channel, rechargeable device with a volume of $31.5~\rm cm^3$, and the G102A is a single-channel, nonrechargeable IPG with a volume of $35.9~\rm cm^3$. The frequency range of these devices is $2-250~\rm Hz$, and the range of pulse widths is $30-450~\rm \mu s$. While these devices are not capable of current steering, similar to the Medtronic devices, they have been used for VFS and remote wireless programming. These devices are not FDA approved in the United States, but they currently are widely implanted throughout China.

SceneRay

The SceneRay model 1180 IPG is a constant-voltage, dual-channel, nonrechargeable IPG with a volume of 51 cm 3 . The frequency range of this device is 1–1,600 Hz, and the range of possible pulse widths is 60–960 μ s. Similar to PINS Medical, use of wireless remote programming has been reported with this device.

The features of each of these devices are summarized in Table 1.

Commercially available DBS leads

Medtronic

Current commercially available DBS leads from Medtronic include the models 3387, 3389, and 3391. Each lead is a traditional quadripolar DBS lead of 1.27 mm in diameter. The contacts of the 3387 and 3389 leads are 1.5 mm in length. The 3387 lead has a contact spacing of 1.5 mm, and the 3389 lead has a contact spacing of 0.5 mm. The 3391 DBS lead is specifically designed for targeting the anterior limb of the internal capsule for treatment of refractory obsessive-compulsive disorder. The contacts of the 3391 lead are 3 mm in length and separated from one another by 4 mm. The Medtronic DBS leads are MRI conditional.

Abbott Infinity

The Abbott DBS leads include both 4-contact quadripolar electrodes and 8-contact directional electrodes. Additionally, Abbott manufactures adaptors for Medtronic DBS electrodes that allow for compatibility with Abbott Infinity IPGs. The quadripolar leads for the Infinity system have a diameter of 1.29 mm and a total length of either 30 cm (models 6166 and 6167) or 40 cm (models 6168 and 6169). All have 4 contacts of 1.5 mm in length with spacing of either 0.5 mm (models 6166 and 6168) or 1.5 mm (models 6167 and 6169). These DBS leads are not MRI compat-



Table 1. Features of currently available implantable pulse generators

Neuromodulation	No. of	Volume (cm³)	RC	Frequency range (Hz)	VFS	Pulse widths (µs)	Fractionation		Directional	MRI	Remote
system cha	channels						Temporal	Current	lead	safety	programming
Abbott (St. Jude Me	dical) Infinit	у									
5	2	30.4	No	2–240	No	20-500	MSS	CA	Yes	С	No
7	2	38.6	No	2–240	No	20-500	MSS	CA	Yes	С	No
Boston Scientific Vercise											
PC	2	33.0	No	2–255	No	10–450	Α	MICC	Yes	U	No
RC	2	22.7	Yes	2–255	No	10-450	Α	MICC	No	U	No
Gevia	2	19.8	Yes	2–255	No	20-450	Α	MICC	Yes	С	No
Medtronic Activa											
PC	2	37.0	No	2–250	No	60–450	IL	No	No	С	No
RC	2	22.0	Yes	2–250	No	60–450	IL	No	No	С	No
SC	1	27.0	No	3–250	No	60–450	IL	No	No	С	No
PINS Medical											
G102	2	51.2	No	2, 5–130, 132–250	Yes	30–450	No	No	No	UK	Yes
G102R	2	31.5	Yes	2, 5–130, 132–250	Yes	30–450	No	No	No	UK	Yes
G101A	1	35.9	No	2, 3–130, 132–250	Yes	30-450	No	No	No	UK	Yes
SceneRay Co., Ltd.											
1180	2	51	No	1–1,600	No	60–960	No	No	No	UK	Yes

A: areas, C: conditional, CA: coactivation, IL: interleaving, MICC: multiple independent current control, MSS: multi-stim set, PC: primary cell, RC: rechargeable, SC: single cell, U: unsafe, UK: unknown, VFS: variable frequency stimulation.

ible. Similar to the quadripolar leads, Abbott directional leads have a diameter of 1.29 mm and come in lengths of 30 cm (models 6170 and 6171) or 40 cm (models 6172 and 6173). The center contacts, 2 and 3, consist of 3 radial contacts (named 'a,' 'b,' and 'c,' e.g., 2a) that allow for directional current steering, and contacts 1 and 4 are single contacts similar to those of the typical quadripolar leads. The directional leads include a radio-opaque marker proximal to contact 4 that can be used to determine the orientation of the directional contacts (direction of segment 'a'). The contact spacing of models 6170 and 6172 is 0.5 mm, and it is 1.5 mm for models 6171 and 6173, making these leads ideal for directional stimulation in targets such as the GPi or the Vim. The Abbott directional DBS leads are MRI conditional.

Boston Scientific

Boston Scientific manufactures a linear, 8-contact lead as well as the Vercise Cartesia 8-contact directional lead, both of which have a diameter of 1.3 mm. Each lead type may have a length of 30 cm (models DB-2201-30DC and DB-2202-30) or 45 cm (DB-2201-45DC and DB-2202-45). Models DB-2201-30DC and DB-2201-45DC contain 8 contacts of 1.5 mm in length spaced by 0.5 mm. The Vercise Cartesia directional leads (models DB-2202-30 and DB-2202-45) contain two levels of 3 radially arranged directional contacts as well as a distal active tip and a proximal omnidirectional contact. The nonsegmented contacts are present at the bottom and top level of the lead (contact 1 and 8 for one channel and 9 and 16 for the other); segmented contacts are 2,

3, and 4 for the second level (10, 11, and 12 for the second channel), and 5, 6, and 7 are for the third level (13, 14, and 15 for the second channel). Similar to the Abbott directional lead, there is a proximal radio-opaque marker that enables determination of the orientation of the directional contacts (indicating the direction of the anterior contacts 2 and 5 or 10 and 13 for the other channel). All leads have a contact length of 1.5 mm and intercontact spacing of 0.5 mm. With a contact span of 15.5 mm, the linear 8-contact lead may be advantageous if a longer length of stimulation is required, such as when targeting the anterior limb of the internal capsule.

Similarly to Abbott, Boston Scientific manufactures adaptors for Medtronic DBS electrodes that allow for compatibility with their IPGs (23).²⁰

PINS Medical

The DBS leads manufactured by PINS Medical closely resemble those of Medtronic. The models L302 and L302S are quadripolar leads with contacts of 1.5 mm in length separated from one another by 1.5 mm, similar to the 3397 Medtronic lead. Resembling the Medtronic 3389 lead, the model L301 and L301S leads include 1.5 mm contacts separated by 0.5 mm intervals. The PINS DBS leads have a diameter of 1.3. There is currently no available data to suggest that these leads are MRI conditional.

SceneRay

SceneRay currently offers three types of DBS electrodes, which

differ in the spacing between contacts. All contacts are 1.5 mm in length, and each electrode is 1.27 mm in diameter. The model 1200 lead has 0.5 mm spacing between contacts, the model 1210 lead has 1.5 mm spacing, and the model 1211 has 1.0 mm spacing. Each type of lead is available in lengths of either 300 mm or 400 mm. One unique feature of the SceneRay DBS system is the TouchLoc lead anchoring device, which is an internationally

patented lead fixture device that has been reported to result in less DBS lead drift and require significantly fewer adjustments than traditional lead anchors.⁶³

The features of each of the DBS leads described above are summarized in Table 2, and each electrode is visually presented in Figure 2.

Table 2. Features of currently available leads (see Figure 2)

Manufacturer and lead model	Lead diameter (mm)	Contact length (mm)	Contact spacing (mm)	No. of contacts	Directional	Active tip	MR safety
Abbott (St. Jude Medical) Infinity							
6166 & 6168	1.29	1.5	0.5	4	No	No	U
6167 & 6169	1.29	1.5	1.5	4	No	No	U
6170 & 6172	1.29	1.5	0.5	8	Yes	No	С
6171 & 6173	1.29	1.5	1.5	8	Yes	No	С
Boston Scientific							
Linear 8-contact (DB-2201-30DC/DB-2201045DC)	1.3	1.5	0.5	8	No	No	С
Cartesia directional (DB-2202-30/DB-2202-45)	1.3	1.5	0.5	8	Yes	Yes	С
Medtronic							
3387	1.27	1.5	1.5	4	No	No	С
3389	1.27	1.5	0.5	4	No	No	С
3391	1.27	3.0	4.0	4	No	No	С
PINS Medical							
L301 & L301S	1.3	1.5	0.5	4	No	No	UK
L302 & L302S	1.3	1.5	1.5	4	No	No	UK
SceneRay Co., Ltd.							
1200	1.27	1.5	0.5	4	No	No	UK
1210	1.27	1.5	1.5	4	No	No	UK

C: conditional, U: unsafe, UK: unknown.




Figure 2. Scaled visual representation of currently commercially available deep brain stimulation electrodes. The features and measurements of each lead type are summarized in Table 2.



CHOOSING THE BEST SYSTEM

Many factors may contribute to the choice of system, including the anatomical target, patient factors, and surgeon and neurologist preferences. The choice to use traditional quadripolar leads versus directional leads is most applicable when considering the size of the target nucleus and the potential for stimulation-related side effects. Conversely, patient factors are important when considering using a rechargeable versus a nonrechargeable IPG. While rechargeable devices are proposed to reduce the number of battery changes, sparing the patient repeated surgical procedures, patients must be reliable and capable of recharging their battery on a consistent basis. While the Medtronic Activa RC and Boston Scientific Gevia require charging only once weekly, patients who are disabled or have lack of social support may have difficulty assuming this responsibility. Additionally, the choice of DBS system relies heavily on the preferences and comfort of the treating surgeons and neurologists. Many practitioners have extensive experience using the Medtronic system since Medtronic has dominated the market for the past several decades. Patients living in remote regions or who are likely to relocate may find that Medtronic devices are more widely available.

Anatomical considerations

The STN is a small nucleus, measuring only 6-8 mm in its largest dimension. A typical trajectory of a DBS lead through the STN begins at the dorsolateral border and projects medially and ventrally while ideally avoiding the associative and limbic regions. Surrounding the nucleus laterally and, to some extent, anteriorly is the internal capsule. A lead positioned too lateral within the STN may produce contralateral face and limb contractions via current spread to internal capsule fibers. A lead that is positioned too medial may stimulate fibers of the adjacent oculomotor tract, resulting in eye deviation and diplopia. Additionally, there have been reports of stimulation of ventrally located contacts causing acute hypomanic episodes⁶⁴ as well as stimulation of the underlying substantia nigra resulting in hypomania and depression. 65,66 Anteriorly positioned leads may project into the lateral hypothalamic region, leading to stimulation-related autonomic effects, such sweating and nausea.

The GPi, in contrast to the STN, is a larger target associated with a lower incidence of side effects. The internal capsule borders the GPi anteriorly and medially, and the optic tract is located ventrally. Also located ventral to the GPi is the ansa lenticularis, which separates the GPi from the nucleus basalis of Meynert and the amygdala. A typical DBS lead trajectory extends along the posterolateral aspect of the GPi and just dorsal to the optic tract. While side effects are generally rare, medial or posterior placement of the DBS lead can result in capsular stimulation and

contralateral face and limb contractions, and ventral deviation can produce phosphene and visual phenomena.

The potential for stimulation-related side effects with STN-DBS may warrant consideration of directional DBS leads. One technical point to consider when implanting directional leads is that the second and third levels of contacts contain the directional contacts such that the lead may need to be implanted deeper than usual in order to place the directional contacts within the target region.

Another potential strategy for increasing the therapeutic window for stimulation while avoiding stimulation-induced side effects is the use of shorter pulse widths < 60 µs. Experimental studies in animals have demonstrated that axons are more excitable than cell bodies, 18 and the threshold for activation of neural elements varies with stimulus strength and pulse width. Longer pulse widths may facilitate stimulation of larger diameter, myelinated axons, such as those arising from pyramidal tracts, and shorter pulse widths are believed to target current to cell bodies and smaller axons. Reich et al. 19 observed an inverse relationship between pulse width duration and the threshold for side effects in a study group of four patients with STN-DBS. The authors implemented the Boston Scientific Vercise neurostimulation system for delivery of pulse widths < 60 µs and observed a twofold increase in the therapeutic window of stimulation with a pulse width of 30 µs compared with the standard pulse width of 60 µs. Currently, both the Boston Scientific Vercise and Abbott Infinity systems allow for pulse widths shorter than 60 µs.

Rechargeable vs. nonrechargeable IPGs

The most common nonrechargeable IPG, the Medtronic Activa PC, has an average battery lifespan ranging from 2.6 to 4.5 years, 60,61,67,68 which is approximately 2.5 years shorter than its predecessor, the Kinetra. This is considerably shorter than the marketed lifespan of the PC, which is advertised as approximately 5 years. The lifespan of an IPG is largely dependent on the current demands of the patient's stimulation settings, with high-voltage stimulation necessitating more frequent battery changes. Rechargeable IPGs were first introduced in 2008, and the predicted lifespan of rechargeable IPGs according the manufacturers ranges from 10 for Brio (Abbott/St. Jude) to 15 (Medtronic), up to 25 (Boston Scientific) years, which has yet to be confirmed.

The goal of rechargeable IPGs is to reduce the number of battery changes, which may reduce costs associated with multiple surgical procedures and minimize the chances of hardware infections. Some studies have reported an increased incidence of infections with repeated IPG replacements, ^{69,70} but one study did not find such an association.⁷¹ Nevertheless, the use of a rechargeable IPG over a 9-year period (the initially approved lifespan of Activa RC) has been estimated to lead to savings of \$60,900 based

Table 3. Advantages and disadvantages of rechargeable IPGs

Advantages	Disadvantages					
Smaller size • Better suited for women and children	Inconvenience • Need for recharging on a regular basis					
Patient autonomy Less frequent battery changes Less frequent follow-up appointments Lower infection risk Fewer surgical procedures	Possibility for technical difficulties • Poor connection between the IPG and charger • Risk of sudden therapy withdraw if charging fails Not favorable for all patients • Patients with cognitive issues					
	Patients who are not technology savvyPatients without caregivers					
Cost-effectiveness Fewer surgical procedures Reduced use of hospital resources	Cost ineffectiveness Patients with a reduced lifespan Patients requiring IPG removal (e.g., infection, lack of efficacy)					

IPG: implantable pulse generator.

on a single-center study involving 206 patients. 72 Since rechargeable IPGs became available, many patients have been switched to rechargeable devices, and some have even received rechargeable IPGs during initial DBS implantation.^{72,73} Patient surveys have demonstrated that the vast majority of patients find the recharging process to be easy, and the incidence of adverse events and temporary interruptions in therapy is low. 67,72,73

The decision to use a rechargeable versus a nonrechargeable IPG is generally left to the discretion of the treating physician and may be influenced by the current requirements of the stimulation settings as well as the reliability of the patient to recharge their device (Table 3). The Multi-Recharge Trial, conducted at four neurosurgical centers across Germany, reported on the survey results of 56 patients switched to rechargeable IPGs and 139 patients who were implanted with rechargeable IPGs during DBS implantation.⁶⁷ According to survey results, patients charged their IPGs once every 10 ± 8 days and spent a mean of 122 ± 175 minutes per week recharging their device. Dystonia patients found the recharging process to be significantly more convenient than PD patients. Interestingly, the Activa RC was found to require significantly more charging time per week than the Boston Scientific Vercise and Abbott (St. Jude) Brio. While charging failures were uncommon, they were more likely to occur with patients who were switched from a nonrechargeable IPG. While the decision to use a rechargeable IPG should take patient factors into consideration, the results of these survey studies suggest that rechargeable IPGs are effective and easy to manage even in older patients.67,73

CONCLUSIONS

Technological advances in DBS hardware and software presently comprise the greatest area of growth within the field of functional neurosurgery. While studies evaluating the features of new DBS systems have found gains in measures of effective stimulation, further research is needed to characterize the impact of these technologies on patient outcomes and quality of life. In a smallsample, cross-sectional, questionnaire-based study, it was found that patients treated with an established DBS system (Medtronic) were less inclined to think that additional programming would lead to further improvement, thus suggesting that the additional features of newer devices may complicate the DBS programming of patients, increasing the number of visits necessary to convince patients that stimulation is optimized and thus possibly increasing physicians' and patients' frustrations.74

Many of the features of these new systems currently focus on minimizing stimulation-related side effects. Further understanding of the pathophysiology of these disorders—and PD in particular—is required to determine if DBS therapy is capable of addressing symptoms refractory to standard medications or DBS, such as postural instability, ataxia and FOG. Exploration of other surgical targets, such as the substantial nigra pars reticulata^{75,76} and pedunculopontine nucleus,77,78 have sought to address levodopa-resistant FOG in PD. With continued research, deeper insights into the neural networks that contribute to these disorders in synergy with advancements in stimulation technology will hopefully yield new therapies that can improve the quality of life of these patients.

Conflicts of Interest

AML has received honoraria from Abbott, Boston Scientific, and Medtronic. AF has received honoraria from Abbott, Boston Scientific, Brainlab and Medtronic.

Acknowledgments

Authors are grateful to Desiree Burdyshaw, Binith Cheeran (Abbott) and Benjamin Nagy (Boston Scientific) for having provided some technical information about the devices produced by their companies.

Author Contributions

Conceptualization: Alfonso Fasano. Data curation: Michelle Paff. Methodology: Michelle Paff. Supervision: Andres M. Lozano, Alfonso Fasano. Visualization: Aaron Loh. Writing-original draft: Michelle Paff, Aaron Loh, Can Sarcia. Writing-review & editing: Andres M. Lozano, Alfonso Fasano. Approval of final manuscript: all authors.



ORCID iDs

 Michelle Paff
 https://orcid.org/0000-0002-6780-2855

 Aaron Loh
 https://orcid.org/0000-0001-5297-6775

 Can Sarica
 https://orcid.org/0000-0001-8419-7426

 Andres M. Lozano
 https://orcid.org/0000-0001-8257-3694

 Alfonso Fasano
 https://orcid.org/0000-0001-5346-0180

REFERENCES

- Siegfried J, Lippitz B. Bilateral chronic electrostimulation of ventroposterolateral pallidum: a new therapeutic approach for alleviating all parkinsonian symptoms. Neurosurgery 1994;35:1126-1129. discussion 1129-1130.
- Limousin P, Krack P, Pollak P, Benazzouz A, Ardouin C, Hoffmann D, et al. Electrical stimulation of the subthalamic nucleus in advanced Parkinson's disease. N Engl J Med 1998;339:1105-1111.
- Williams A, Gill S, Varma T, Jenkinson C, Quinn N, Mitchell R, et al. Deep brain stimulation plus best medical therapy versus best medical therapy alone for advanced Parkinson's disease (PD SURG trial): a randomised, open-label trial. Lancet Neurol 2010;9:581-591.
- 4. Odekerken VJJ, van Laar T, Staal MJ, Mosch A, Hoffmann CF, Nijssen PC, et al. Subthalamic nucleus versus globus pallidus bilateral deep brain stimulation for advanced Parkinson's disease (NSTAPS study): a randomised controlled trial. Lancet Neurol 2013;12:37-44.
- Deep-Brain Stimulation for Parkinson's Disease Study Group; Obeso JA, Olanow CW, Rodriguez-Oroz MC, Krack P, Kumar R, Lang AE. Deepbrain stimulation of the subthalamic nucleus or the pars interna of the globus pallidus in Parkinson's disease. N Engl J Med 2001;345:956-963.
- Fasano A, Aquino CC, Krauss JK, Honey CR, Bloem BR. Axial disability and deep brain stimulation in patients with Parkinson disease. Nat Rev Neurol 2015;11:98-110.
- Tsuboi T, Watanabe H, Tanaka Y, Ohdake R, Yoneyama N, Hara K, et al. Distinct phenotypes of speech and voice disorders in Parkinson's disease after subthalamic nucleus deep brain stimulation. J Neurol Neurosurg Psychiatry 2015;86:856-864.
- Johnsen EL, Sunde N, Mogensen PH, Ostergaard K. MRI verified STN stimulation site--gait improvement and clinical outcome. Eur J Neurol 2010;17:746-753.
- Tommasi G, Lopiano L, Zibetti M, Cinquepalmi A, Fronda C, Bergamasco B, et al. Freezing and hypokinesia of gait induced by stimulation of the subthalamic region. J Neurol Sci 2007;258:99-103.
- Parent A, Hazrati LN. Functional anatomy of the basal ganglia. II. The place of subthalamic nucleus and external pallidum in basal ganglia circuitry. Brain Res Brain Res Rev 1995;20:128-154.
- Herzog J, Fietzek U, Hamel W, Morsnowski A, Steigerwald F, Schrader B, et al. Most effective stimulation site in subthalamic deep brain stimulation for Parkinson's disease. Mov Disord 2004;19:1050-1054.
- Godinho F, Thobois S, Magnin M, Guenot M, Polo G, Benatru I, et al. Subthalamic nucleus stimulation in Parkinson's disease: anatomical and electrophysiological localization of active contacts. J Neurol 2006;253:1347-1355
- Moreau C, Defebvre L, Destée A, Bleuse S, Clement F, Blatt JL, et al. STN-DBS frequency effects on freezing of gait in advanced Parkinson disease. Neurology 2008;71:80-84.
- Xie T, Kang UJ, Warnke P. Effect of stimulation frequency on immediate freezing of gait in newly activated STN DBS in Parkinson's disease. J Neurol Neurosurg Psychiatry 2012;83:1015-1017.
- Volkmann J, Moro E, Pahwa R. Basic algorithms for the programming of deep brain stimulation in Parkinson's disease. Mov Disord 2006;21 Suppl 14:S284-S289.
- Almeida L, Rawal PV, Ditty B, Smelser BL, Huang H, Okun MS, et al. Deep brain stimulation battery longevity: comparison of monopolar versus bipolar stimulation modes. Mov Disord Clin Pract 2016;3:359-366.
- 17. Moro E, Esselink RJ, Xie J, Hommel M, Benabid AL, Pollak P. The impact on Parkinson's disease of electrical parameter settings in STN stimula-

- tion. Neurology 2002;59:706-713.
- Ranck JB Jr. Which elements are excited in electrical stimulation of mammalian central nervous system: a review. Brain Res 1975;98:417-440.
- Reich MM, Steigerwald F, Sawalhe AD, Reese R, Gunalan K, Johannes S, et al. Short pulse width widens the therapeutic window of subthalamic neurostimulation. Ann Clin Transl Neurol 2015;2:427-432.
- Soh D, Lozano AM, Fasano A. Hybrid deep brain stimulation system to manage stimulation-induced side effects in essential tremor patients. Parkinsonism Relat Disord 2019;58:85-86.
- Sidiropoulos C, Walsh R, Meaney C, Poon YY, Fallis M, Moro E. Low-frequency subthalamic nucleus deep brain stimulation for axial symptoms in advanced Parkinson's disease. J Neurol 2013;260:2306-2311.
- Khoo HM, Kishima H, Hosomi K, Maruo T, Tani N, Oshino S, et al. Lowfrequency subthalamic nucleus stimulation in Parkinson's disease: a randomized clinical trial. Mov Disord 2014;29:270-274.
- Jia F, Guo Y, Wan S, Chen H, Hao H, Zhang J, et al. Variable frequency stimulation of subthalamic nucleus for freezing of gait in Parkinson's disease. Parkinsonism Relat Disord 2015;21:1471-1472.
- Jia F, Hu W, Zhang J, Wagle Shukla A, Almeida L, Meng FG, et al. Variable frequency stimulation of subthalamic nucleus in Parkinson's disease: rationale and hypothesis. Parkinsonism Relat Disord 2017;39:27-30.
- Harmsen IE, Lee DJ, Dallapiazza RF, De Vloo P, Chen R, Fasano A, et al. P.089 ultra-high frequency deep brain stimulation at 10,000 Hz improves motor function (abstracts: 54th Annual Congress of the Canadian Neurological Sciences Federation). Can J Neurol Sci 2019;46(Suppl 1):S37.
- Kovács N, Janszky J, Nagy F, Balás I. Changing to interleaving stimulation might improve dystonia in cases not responding to pallidal stimulation. Mov Disord 2012;27:163-165.
- Baumann CR, Imbach LL, Baumann-Vogel H, Uhl M, Sarnthein J, Sürücü
 O. Interleaving deep brain stimulation for a patient with both Parkinson's
 disease and essential tremor. Mov Disord 2012;27:1700-1701.
- Miocinovic S, Khemani P, Whiddon R, Zeilman P, Martinez-Ramirez D, Okun MS, et al. Outcomes, management, and potential mechanisms of interleaving deep brain stimulation settings. Parkinsonism Relat Disord 2014;20:1434-1437.
- Zhang S, Zhou P, Jiang S, Wang W, Li P. Interleaving subthalamic nucleus deep brain stimulation to avoid side effects while achieving satisfactory motor benefits in Parkinson disease: a report of 12 cases. Medicine (Baltimore) 2016;95:e5575.
- Kern DS, Picillo M, Thompson JA, Sammartino F, di Biase L, Munhoz RP, et al. Interleaving stimulation in Parkinson's disease, tremor, and dystonia. Stereotact Funct Neurosurg 2018;96:379-391.
- Butson CR, McIntyre CC. Current steering to control the volume of tissue activated during deep brain stimulation. Brain Stimul 2008;1:7-15.
- Zhang S, Silburn P, Pouratian N, Cheeran B, Venkatesan L, Kent A, et al. Comparing current steering technologies for directional deep brain stimulation using a computational model that incorporates heterogeneous tissue properties. Neuromodulation 2020;23:469-477.
- Chaturvedi A, Foutz TJ, McIntyre CC. Current steering to activate targeted neural pathways during deep brain stimulation of the subthalamic region. Brain Stimul 2012;5:369-377.
- Barbe MT, Maarouf M, Alesch F, Timmermann L. Multiple source current steering--a novel deep brain stimulation concept for customized programming in a Parkinson's disease patient. Parkinsonism Relat Disord 2014;20:4 71-473.
- Pavese N, Tai YF, Yousif N, Nandi D, Bain PG. Traditional trial and error versus neuroanatomic 3-dimensional image software-assisted deep brain stimulation programming in patients with Parkinson disease. World Neurosurg 2020;134:e98-e102.
- BeMent SL, Ranck JB Jr. A quantitative study of electrical stimulation of central myelinated fibers. Exp Neurol 1969;24:147-170.
- Anderson DN, Duffley G, Vorwerk J, Dorval AD, Butson CR. Anodic stimulation misunderstood: preferential activation of fiber orientations with anodic waveforms in deep brain stimulation. J Neural Eng 2019;16: e016026.

- Kirsch AD, Hassin-Baer S, Matthies C, Volkmann J, Steigerwald F. Anodic versus cathodic neurostimulation of the subthalamic nucleus: a randomized-controlled study of acute clinical effects. Parkinsonism Relat Disord 2018;55:61-67.
- Soh D, Ten Brinke TR, Lozano AM, Fasano A. Therapeutic window of deep brain stimulation using cathodic monopolar, bipolar, semi-bipolar, and anodic stimulation. Neuromodulation 2019;22:451-455.
- Keane M, Deyo S, Abosch A, Bajwa JA, Johnson MD. Improved spatial targeting with directionally segmented deep brain stimulation leads for treating essential tremor. J Neural Eng 2012;9:e046005.
- Martens HCF, Toader E, Decré MMJ, Anderson DJ, Vetter R, Kipke DR, et al. Spatial steering of deep brain stimulation volumes using a novel lead design. Clin Neurophysiol 2011;122:558-566.
- Teplitzky BA, Zitella LM, Xiao Y, Johnson MD. Model-based comparison of deep brain stimulation array functionality with varying number of radial electrodes and machine learning feature sets. Front Comput Neurosci 2016;10:e58
- Dembek TA, Reker P, Visser-Vandewalle V, Wirths J, Treuer H, Klehr M, et al. Directional DBS increases side-effect thresholds—A prospective, double-blind trial. Mov Disord 2017;32:1380-1388.
- Pollo C, Kaelin-Lang A, Oertel MF, Stieglitz L, Taub E, Fuhr P, et al. Directional deep brain stimulation: an intraoperative double-blind pilot study. Brain 2014;137:2015-2026.
- Steigerwald F, Müller L, Johannes S, Matthies C, Volkmann J. Directional deep brain stimulation of the subthalamic nucleus: a pilot study using a novel neurostimulation device. Mov Disord 2016;31:1240-1243.
- Contarino MF, Bour LJ, Verhagen R, Lourens MA, de Bie RM, van den Munckhof P, et al. Directional steering: a novel approach to deep brain stimulation. Neurology 2014;83:1163-1169.
- 47. Schnitzler AS, Mir PM, Brodsky MB, Verhagen LV, Groppa SG, Cheeran BC, et al. Directional versus omnidirectional Deep Brain Stimulation for Parkinson's disease: Results of a prospective, blinded, multi center, single-arm crossover study [abstract no. 201]. Mov Disord 2019;34(Suppl 2).
- Anderson DN, Osting B, Vorwerk J, Dorval AD, Butson CR. Optimized programming algorithm for cylindrical and directional deep brain stimulation electrodes. J Neural Eng 2018;15:e026005.
- Young DT. Medtronic receives CE mark approval for the Percept[™] PC neurostimulator DBS system with BrainSense[™] technology [Internet]. Dublin: Medtronic. c2020. [cited 2020 May 24]. Available from: http://newsroom.medtronic.com/node/30866/pdf.
- Okun MS, Gallo BV, Mandybur G, Jagid J, Foote KD, Revilla FJ, et al. Subthalamic deep brain stimulation with a constant-current device in Parkinson's disease: an open-label randomised controlled trial. Lancet Neurol 2012;11:140-149.
- Wharen RE Jr, Okun MS, Guthrie BL, Uitti RJ, Larson P, Foote K, et al. Thalamic DBS with a constant-current device in essential tremor: a controlled clinical trial. Parkinsonism Relat Disord 2017;40:18-26.
- Holtzheimer PE, Husain MM, Lisanby SH, Taylor SF, Whitworth LA, Mc-Clintock S, et al. Subcallosal cingulate deep brain stimulation for treatment-resistant depression: a multisite, randomised, sham-controlled trial. Lancet Psychiatry 2017;4:839-849.
- 53. Timmermann L, Jain R, Chen L, Maarouf M, Barbe MT, Allert N, et al. Multiple-source current steering in subthalamic nucleus deep brain stimulation for Parkinson's disease (the VANTAGE study): a non-randomised, prospective, multicentre, open-label study. Lancet Neurol 2015;14:693-701.
- 54. Vitek J, Jain R, Starr P. INTREPID trial: a prospective, double blinded, multi-center randomized controlled trial evaluating deep brain stimulation with a new multiple-source, constant-current rechargeable system in Parkinson's disease (P5.016). Neurology 2017;88(Suppl 16):P5.016.
- 55. Vitek J, Jain R, Starr P; INTREPID Clinical Study Group. Two year outcomes: a prospective, double blinded, multicenter randomized controlled trial evaluating deep brain stimulation with a new multiple source, constant current rechargeable system in Parkinson's disease (INTREPID) (P1.8-026). Neurology 2019;92(Suppl 15):P1.8-026.
- 56. Liu HG, Ma Y, Zhang K, Ge M, Meng FG, Feng T, et al. Subthalamic deep

- brain stimulation with a new device in Parkinson's disease: an open-label trial. Neuromodulation 2013;16:212-218. discussion 218.
- 57. Jia F, Wagle Shukla A, Hu W, Almeida L, Holanda V, Zhang J, et al. Deep brain stimulation at variable frequency to improve motor outcomes in Parkinson's disease. Mov Disord Clin Pract 2018;5:538-541.
- Zhang C, Li D, Jin H, Zeljic K, Sun B. Target-specific deep brain stimulation of the ventral capsule/ventral striatum for the treatment of neuropsychiatric disease. Ann Transl Med 2017;5:e402.
- Zhang C, Zhang Y, Zhan S, Li D, Jin H, Denys D, et al. Telemedical deep brain stimulation: merits and limitations. Stereotact Funct Neurosurg 2018;96:272-273.
- 60. Fisher B, Kausar J, Garratt H, Hodson J, White A, Ughratdar I, et al. Battery longevity comparison of two commonly available dual channel implantable pulse generators used for subthalamic nucleus stimulation in Parkinson's disease. Stereotact Funct Neurosurg 2018;96:151-156.
- Sette AL, Seigneuret E, Reymond F, Chabardes S, Castrioto A, Boussat B, et al. Battery longevity of neurostimulators in Parkinson disease: a historic cohort study. Brain Stimul 2019;12:851-857.
- 62. Li D, Zhang C, Gault J, Wang W, Liu J, Shao M, et al. Remotely programmed deep brain stimulation of the bilateral subthalamic nucleus for the treatment of primary Parkinson disease: a randomized controlled trial investigating the safety and efficacy of a novel deep brain stimulation system. Stereotact Funct Neurosurg 2017;95:174-182.
- Wang T, Pan Y, Zhang C, Zhan S, Sun B, Li D. Lead fixation in deep brain stimulation: comparison of three lead anchoring devices in China. BMC Surg 2019;19:e92.
- 64. Mallet L, Schüpbach M, N'Diaye K, Remy P, Bardinet E, Czernecki V, et al. Stimulation of subterritories of the subthalamic nucleus reveals its role in the integration of the emotional and motor aspects of behavior. Proc Natl Acad Sci U S A 2007;104:10661-10666.
- 65. Ulla M, Thobois S, Llorca PM, Derost P, Lemaire JJ, Chereau-Boudet I, et al. Contact dependent reproducible hypomania induced by deep brain stimulation in Parkinson's disease: clinical, anatomical and functional imaging study. J Neurol Neurosurg Psychiatry 2011;82:607-614.
- Blomstedt P, Hariz MI, Lees A, Silberstein P, Limousin P, Yelnik J, et al. Acute severe depression induced by intraoperative stimulation of the substantia nigra: a case report. Parkinsonism Relat Disord 2008;14:253-256.
- 67. Jakobs M, Helmers AK, Synowitz M, Slotty PJ, Anthofer JM, Schlaier JR, et al. A multicenter, open-label, controlled trial on acceptance, convenience, and complications of rechargeable internal pulse generators for deep brain stimulation: the Multi Recharge Trial. J Neurosurg 2019 Aug 16 [epub]. https://doi.org/10.3171/2019.5.JNS19360.
- Helmers AK, Lübbing I, Deuschl G, Witt K, Synowitz M, Mehdorn HM, et al. Comparison of the battery life of nonrechargeable generators for deep brain stimulation. Neuromodulation 2018;21:593-596.
- 69. Fytagoridis A, Heard T, Samuelsson J, Zsigmond P, Jiltsova E, Skyrman S, et al. Surgical replacement of implantable pulse generators in deep brain stimulation: adverse events and risk factors in a multicenter cohort. Stereotact Funct Neurosurg 2016;94:235-239.
- Thrane JF, Sunde NA, Bergholt B, Rosendal F. Increasing infection rate in multiple implanted pulse generator changes in movement disorder patients treated with deep brain stimulation. Stereotact Funct Neurosurg 2014;92: 360-364.
- 71. Frizon LA, Hogue O, Wathen C, Yamamoto E, Sabharwal NC, Jones J, et al. Subsequent pulse generator replacement surgery does not increase the infection rate in patients with deep brain stimulator systems: a review of 1537 unique implants at a single center. Neuromodulation 2017;20:444-449
- Hitti FL, Vaughan KA, Ramayya AG, McShane BJ, Baltuch GH. Reduced long-term cost and increased patient satisfaction with rechargeable implantable pulse generators for deep brain stimulation. J Neurosurg 2018; 131:799-806.
- Jakobs M, Kloß M, Unterberg A, Kiening K. Rechargeable internal pulse generators as initial neurostimulators for deep brain stimulation in patients with movement disorders. Neuromodulation 2018;21:604-610.



- Soh D, Ten Brinke TR, Lozano AM, Fasano A. Novel deep brain stimulation technologies for Parkinson's disease: more expectations, more frustrations? Mov Disord Clin Pract 2019;7:113-114.
- 75. Weiss D, Wächter T, Meisner C, Fritz M, Gharabaghi A, Plewnia C, et al. Combined STN/SNr-DBS for the treatment of refractory gait disturbances in Parkinson's disease: study protocol for a randomized controlled trial. Trials 2011;12:e222.
- 76. Valldeoriola F, Muñoz E, Rumià J, Roldán P, Cámara A, Compta Y, et al. Simultaneous low-frequency deep brain stimulation of the substantia nig-
- ra pars reticulata and high-frequency stimulation of the subthalamic nucleus to treat levodopa unresponsive freezing of gait in Parkinson's disease: a pilot study. Parkinsonism Relat Disord 2019;60:153-157.
- Moro E, Hamani C, Poon YY, Al-Khairallah T, Dostrovsky JO, Hutchison WD, et al. Unilateral pedunculopontine stimulation improves falls in Parkinson's disease. Brain 2010;133:215-224.
- 78. Thevathasan W, Coyne TJ, Hyam JA, Kerr G, Jenkinson N, Aziz TZ, et al. Pedunculopontine nucleus stimulation improves gait freezing in Parkinson disease. Neurosurgery 2011;69:1248-1253. discussion 1254.