

EDITORIAL

Open Access



# Clinical application of esophageal manometry: how I do it

Elias Baedorf Kassis<sup>1\*</sup>  and Daniel Talmor<sup>2</sup>

Our group uses esophageal manometry routinely to personalize mechanical ventilation in patients with acute respiratory distress syndrome (ARDS) [1, 2]. Esophageal pressures ( $P_{es}$ ) allow for differentiation of chest wall, lung and respiratory system mechanics, and we use this for PEEP titration [1, 2], monitoring of parenchymal lung stress, limiting peak end-inspiratory transpulmonary pressures and monitoring for ventilator synchrony [3, 4].

We find that esophageal manometry is straightforward in the majority of patients although proper training and application are important. The initial step is to assure correct placement with insertion of stand-alone catheters or feeding tubes with integrated esophageal balloons which are similar to routine gastric tubes. Typical depth of insertion ranges from 33 to 40 cm, depending on body size and we assure proper placement through functional bedside assessment. First, we look for the presence of cardiac oscillations to assure correct position posterior to the heart. If absent, this suggests the balloon is too deep or shallow and we incrementally adjust while monitoring for these oscillations. Next we perform expiratory breath holds, with changes in  $P_{es}$ , airway ( $P_{ao}$ ) and transpulmonary pressure ( $P_L = P_{ao} - P_{es}$ ) monitored during gentle chest pushes. Proper position is confirmed when  $P_{es}$  and  $P_{ao}$  increase in equal measure, with no change in the calculated  $P_L$ . If  $P_{ao}$  increases more  $P_{es}$ , this suggests that position is too deep and the balloon is adjusted incrementally with repeat chest pushes. This may be confirmed with gentle abdominal pushes (with  $P_{es}$  increasing more than  $P_{ao}$ ). (Table 1).

Using a balloon with a consistent working range of inflation volume is helpful for obtaining consistent and accurate measurements. While optimal inflation volume can be confirmed based upon the pressure–volume characteristics of the balloon itself [5], this is time-consuming and not required in practice when using a balloon with a known acceptable range. Overinflation results in inaccurately high measured pressures secondary to the compliance of the balloon, while underinflation causes dampening of waveform variation. Visualization and interpretation of data are facilitated by integrated pressure sensors within the ventilator or can be recorded using stand-alone devices as we used in the EPVent and EPVent2 studies [1, 2].

One of our primary applications of esophageal manometry is for titration of positive end-expiratory pressure (PEEP). Critically ill patients frequently exhibit increased chest wall weight and elevated basal end-expiratory pleural pressures secondary to edema, effusions, abdominal hypertension and other causes that may lead to derecruitment, increased lung elastance and hypoxemia. We measure  $P_{es}$  as a surrogate for pleural pressure [6] and if the pleural pressure is larger than the measured airway/alveolar pressure ( $PL = P_{ao} - P_{es}$ ), these collapsing pressures can be countered with the application of PEEP. Our EPVent [2] and EPvent2 [1] studies investigated the use of esophageal manometry to titrate PEEP and while the latter study did not show clear benefit compared with empiric high-PEEP, further analysis suggested a benefit when end-expiratory  $P_L$  were maintained in a tight physiological range of  $-2$  to  $+2$  cmH<sub>2</sub>O with PEEP adjustment (publication under review) which is how we practice clinically. We aim for an end-expiratory  $P_L$  of zero *regardless of the FiO<sub>2</sub>* which is distinct from the original sliding-scale protocols [1, 2] and is in part secondary to the slight

\*Correspondence: enbaedor@bidmc.harvard.edu

<sup>1</sup> Division of Pulmonary and Critical Care and Harvard Medical School, Beth Israel Deaconess Medical Center, Boston, MA 02215, USA  
Full list of author information is available at the end of the article



**Table 1 Tricks and troubleshooting**

	How we do it	Troubleshooting
Proper placement	<ol style="list-style-type: none"> <li>(1) Placement depth: Usual depth is 33–40 cm (a good starting point is 37 cm)</li> <li>(2) Balloon inflation: Use a balloon with a consistent working volume. Optimization of volume otherwise will need to be done by measuring the pressure–volume characteristics of the balloon itself which is not always feasible</li> <li>(3) Cardiac oscillations: Cardiac oscillations should be present to confirm placement posterior to the heart above the diaphragm</li> <li>(4) End-Expiratory Hold Chest Pushes: Pes and Pao should increase in equal measures with chest push resulting in no change in <math>P_L</math></li> </ol>	<ol style="list-style-type: none"> <li>(1) Depth is incrementally adjusted while looking for oscillations and doing chest pushes</li> <li>(2) Depending on body habitus and the unique patient, the amplitude of oscillations may be widely variable</li> <li>(3) If patient is ACTIVELY breathing, expiratory breath holds can still confirm placement. Pao and Pes DECREASES in this case without change in <math>P_L</math></li> </ol>
PEEP titration	<ol style="list-style-type: none"> <li>(1) Measurement of Pes and PL during expiratory holds</li> <li>(2) Adjust PEEP until PL at end –expiration = zero</li> <li>(3) We use in most patients with moderate-severe ARDS</li> <li>(4) Especially useful with obesity or abdominal hypertension</li> </ol>	<ol style="list-style-type: none"> <li>(1) We no longer use the sliding scale <math>FiO_2</math>—<math>P_L</math> used in the EPVent studies, targeting <math>P_L</math> = zero with acceptable range from -2 to +2 cmH<sub>2</sub>O</li> </ol>
Monitoring cyclic and total lung stress	<ol style="list-style-type: none"> <li>(1) Measure end-inspiratory PL: This measurement is obtained when the plateau pressure is measured during an inspiratory breath hold. We keep the end-inspiratory PL &lt; 20cmH<sub>2</sub>O and ideally aim for &lt; 15cmH<sub>2</sub>O to provide additional safety</li> <li>(2) Measure <math>\Delta P_L</math>: Calculated as the end-inspiratory PL minus the end expiratory PL. This provides a more targeted driving pressure measurement than the respiratory system values, and we aim for less than 10–12cmH<sub>2</sub>O</li> <li>(3) Targeted titration of tidal volume if above target values if allowable with ventilation requirements</li> </ol>	<ol style="list-style-type: none"> <li>(1) With large cardiac oscillations, use the diastole phase for measurements to be consistent</li> <li>(2) We recommend PEEP titration/optimization to maximize compliance prior to targeted tidal volume reduction</li> </ol>
Dyssynchrony and neuromuscular blockade	<ol style="list-style-type: none"> <li>(1) This is easiest with systems that integrate xy plots</li> <li>(2) This is more advanced level application and beyond routine use as above</li> </ol>	<ol style="list-style-type: none"> <li>(1) Not recommended for routine use as requires more specialty equipment and training</li> </ol>

benefit gained from mediastinal artifact in the “actual” vs. measured  $P_L$  [7].


We find esophageal manometry particularly useful with morbid obesity and ARDS [8, 9], allowing for measurement of elevated pleural pressures and safe application of high PEEP levels (~20–30cmH<sub>2</sub>O) to offload the weight of the chest wall. Conversely, esophageal manometry is also useful in determining when applied PEEP is *too high* allowing for targeted titration to lower PEEP which may prevent the harmful effects of overdistension (Fig. 1).

As a more specific measure of lung stress, we routinely monitor the cyclic distending pressures across the lungs (transpulmonary driving pressure [ $\Delta P_L$ ]) [10]. While the respiratory system driving pressure correlated with mortality in patients with ARDS [11], we believe it is inadequate due to the inherent variability and heterogeneity of the chest wall which we can directly measure using Pes. We target a  $\Delta P_L$  of less than 10–12cmH<sub>2</sub>O due to lung inhomogeneity and local stress raisers [12], which could prevent lung injury [13] and is in agreement with our

retrospective mortality data [10].  $\Delta P_L$  is easily measured as the end-inspiratory PL (plateau pressure equivalent) minus the end-expiratory PL (total PEEP equivalent).

In addition to the cyclic lung stress, we use  $P_{es}$  to measure the total lung stress (end-inspiratory  $P_L$ ) in addition to plateau pressure. A plateau pressure < 30cmH<sub>2</sub>O represents a widely varying level of lung stress depending on the chest wall mechanics. While safe levels have not been clearly defined, we have extrapolated thresholds from our understanding of the relationship between stress, strain and specific elastance, and data suggesting high  $\Delta P_L$  and end-inspiratory  $P_L$  can bring the lung to total lung capacity and lead to lethal ventilator induced lung injury [14]. As such, our practice is to limit end-inspiratory  $P_L$  to less than 20cmH<sub>2</sub>O (and ideally even lower, to < 15cmH<sub>2</sub>O), to decrease overdistension and improve the margin of safety.

If we identify a patient with elevated  $\Delta P_L$  or total end-inspiratory  $P_L$ , this data is used to facilitate targeted tidal volume reduction, to bring these values within safer



limits. We recognize that further prospective investigation of these limits is warranted to better clarify targets, but we synthesize these data with other clinical data to help inform our bedside care. Importantly, with widely variable chest wall pressures and elastance, *we cannot predict if we are reaching these thresholds* of cyclic and total stress without the use of an esophageal balloon.

We also use esophageal manometry in the bedside assessment of patient-ventilator synchrony using the chest wall pressure volume loops for identification of passive ventilator delivered breaths, spontaneous breaths, dyssynchrony [3, 4], for titration of neuromuscular blockade [6] and for direct measurement of work of breathing, inspiratory muscle efforts and lung-directed mechanical power to assess when levels of effort may be harmful [15]. In conclusion, esophageal balloon catheters are easily placed and interpreted. Measured esophageal pressures and calculation of transpulmonary pressures have broad applications for personalized care of mechanically ventilated patients with PEEP titration, measurement of lung stress and assessment for ventilator synchrony.

#### Authors' contributions

Article was written and edited by both authors.

#### Funding

Dr Baedorf Kassis' work was conducted with the support of a KL2 award (an appointed KL2 award) from Harvard Catalyst | The Harvard Clinical and Translational Science Center (National Center for Advancing Translational Sciences, National Institutes of Health Award KL2 TR002542). The content is solely the responsibility of the authors and does not necessarily represent the official views of Harvard Catalyst, Harvard University and its affiliated academic healthcare centers, or the National Institutes of Health.

#### Availability of data and material

Not applicable. If your manuscript does not contain any data, please state not applicable in this section.

#### Consent for publication

Not applicable. No consent needed for this editorial.

#### Ethics approval and consent to participate

Not applicable. There was no ethics approval needed for this editorial.

#### Competing interests

The authors declare that there are no competing interests regarding publication of this editorial. Of note both Dr Baedorf Kassis and Dr Talmor receive educational speaking fees from Hamilton Medical for unrelated work.

**Author details**

<sup>1</sup> Division of Pulmonary and Critical Care and Harvard Medical School, Beth Israel Deaconess Medical Center, Boston, MA 02215, USA. <sup>2</sup> Department of Anesthesia, Critical Care and Pain Medicine and Harvard Medical School, Beth Israel Deaconess Medical Center, Boston, MA 02215, USA.

Received: 8 December 2020 Accepted: 23 December 2020

Published online: 05 January 2021

**References**

1. Beitler JR, Sarge T, Banner-Goodspeed V, Gong MN, Cook DJ, Novack V, et al. Lung mechanics to guide positive end-expiratory pressure in acute respiratory distress syndrome: The EPVent-2 Randomized Clinical Trial. 2019.
2. Talmor D, Sarge T, Malhotra A, O'Donnell CR, Ritz R, Lisbon A, et al. Mechanical ventilation guided by esophageal pressure in acute lung injury. *N Engl J Med*. 2008;359:2095–104.
3. Baedorf Kassis E, Su HK, Graham AR, Novack V, Loring SH, Talmor DS. Reverse trigger phenotypes in acute respiratory distress syndrome. *Am J Respir Crit Care Med*. 2020;
4. Baedorf Kassis E, Loring SH, Talmor D. Lung volumes and transpulmonary pressure are decreased with expiratory effort and restored with passive breathing in ARDS: a reapplication of the traditional Campbell diagram. *Intensive Care Med*. 2018;44.
5. Mojoli F, Iotti GA, Torriglia F, Pozzi M, Volta CA, Bianzina S, et al. In vivo calibration of esophageal pressure in the mechanically ventilated patient makes measurements reliable. *Crit Care*. 2016.
6. Baedorf Kassis E, Train S, MacNeil B, Loring SH, Talmor D. Monitoring of neuromuscular blockade: a comparison of train-of-four and the Campbell diagram. *Intensive Care Med*. 2018;44:2305–6.
7. Baedorf Kassis E, Loring SH, Talmor D, Terragni P, Mascia L, Ranieri VM. A fixed correction of absolute transpulmonary pressure may not be ideal for clinical use: discussion on “Accuracy of esophageal pressure to assess transpulmonary pressure during mechanical ventilation.” *Intensive Care Med*. 2017;43.
8. De Santis Santiago R, Teggie Droghi M, Fumagalli J, Marrazzo F, Florio G, Grassi LG, et al. High pleural pressure prevents alveolar overdistension and hemodynamic collapse in ARDS with Class III obesity. *Am J Respir Crit Care Med*. 2020.
9. Fumagalli J, Santiago RRS, Teggie Droghi M, Zhang C, Fintelmann FJ, Trotschel FM, et al. Lung recruitment in obese patients with acute respiratory distress syndrome. *anesthesiology*. 2019.
10. Baedorf Kassis E, Loring SH, Talmor D. Mortality and pulmonary mechanics in relation to respiratory system and transpulmonary driving pressures in ARDS. *Intensive Care Med*. 2016.
11. Amato MBP, Meade MO, Slutsky AS, Brochard L, Costa ELV, Schoenfeld DA, et al. Driving pressure and survival in the acute respiratory distress syndrome. *N Engl J Med*. 2015.
12. Cressoni M, Amint M, Cadringer P, Chiurazzi C, Febres D, Gallazzi E, et al. Quantification of stress raisers in ARDS. *Crit Care. BioMed Central*; 2013;17:P105.
13. Mauri T, Yoshida T, Bellani G, Goligher EC, Carteaux G, Rittayamai N, et al. Esophageal and transpulmonary pressure in the clinical setting: meaning, usefulness and perspectives. *Intensive Care Med*. 2016.
14. Protti A, Cressoni M, Santini A, Langer T, Miletto C, Febres D, et al. Lung stress and strain during mechanical ventilation: Any safe threshold? *Am J Respir Crit Care Med*. 2011.
15. Schaefer M, Loring SH, Talmor D, Baedorf Kassis E. Comparison of mechanical power estimations in mechanically ventilated patients with the acute respiratory distress syndrome: a secondary data analysis from the EPVent study. *Intensive Care Med*. 2020;Under Revi.

**Publisher's Note**

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more [biomedcentral.com/submissions](https://biomedcentral.com/submissions)

