



# Improved method for isolation of high-quality total RNA from *Agave tequilana* Weber roots

Luis F. Maceda-López<sup>1</sup> · José L. Villalpando-Aguilar<sup>2</sup> · Eleazar García-Hernández<sup>1</sup> · Emmanuel Ávila de Dios<sup>3</sup> · Silvia B. Andrade-Canto<sup>4</sup> · Dalia C. Morán-Velázquez<sup>1</sup> · Lorena Rodríguez-López<sup>1</sup> · Demetrio Hernández-Díaz<sup>1</sup> · Manuel A. Chablé-Vega<sup>1</sup> · Laura Trejo<sup>5</sup> · Elsa Góngora-Castillo<sup>6</sup> · Itzel López-Rosas<sup>7</sup> · June Simpson<sup>8</sup> · Fulgencio Alatorre-Cobos<sup>7</sup>

Received: 1 February 2020 / Accepted: 26 December 2020 / Published online: 15 January 2021  
© King Abdulaziz City for Science and Technology 2021

## Abstract

Together with their undeniable role in the ecology of arid and semiarid ecosystems, *Agave* species are emerging as a model to dissect the relationships between crassulacean acid metabolism and high efficiency of light and water use, and as an energy crop for bioethanol production. Transcriptome resources from economically valuable *Agaves* species, such as *Agave tequilana* and *A. salmiana*, as well as hybrids for fibers, are now available, and multiple gene expression landscape analyses have been reported. Key components in molecular mechanisms underlying drought tolerance could be uncovered by analyzing gene expression patterns of roots. This study describes an efficient protocol for high-quality total RNA isolation from phenolic compounds-rich *Agave* roots. Our methodology involves suitable root handling and collecting in the field and using saving-time commercial kits available. RNA isolated from roots free of lignified out-layers and clean cortex showed high values of quality and integrity according to electrophoresis and microfluidics-based platform. Synthesis of long full-length cDNAs and PCR amplification tested the suitability for downstream applications of extracted RNA. The protocol was applied successfully to *A. tequilana* roots but can be used for other *Agave* species that also develop lignified epidermis/exodermis in roots.

**Keywords** *Agave* · Full-length cDNA · Lignin · RNA extraction · Roots

## Introduction

*Agave* species (Agavaceae) are xeric monocots considered keystone species in arid and semiarid regions and endemic to America. The genus *Agave* groups about 166 species, most of them are native to Mexico where they grow in different habits, from sea level to an elevation of 3400 m (Good-Avila et al. 2006). Traditionally, *Agave* plants have been used for beverages, food, fiber, medicines, shelter, and as ornamentals (Niechayev et al. 2019). Because of their remarkable tolerance to long drought periods, high light efficiency, minimal water needs, high biomass productivity, and high contents of soluble sugar and lignocellulose, *Agave* species are considered an emerging biofuel crop (Davis et al. 2011; Escamilla-Treviño 2012; Gross et al. 2013; Corbin et al. 2015).

Additionally, analysis of the crassulacean acid metabolism (CAM) in *Agave* plants is opening up a new avenue of research to create more efficient crops based on CAM traits (Abraham et al. 2016). This growing interest, which began in the last decade, is contrasted with the relatively scarce transcriptomic resources available for some *Agave* species (Table 1), and with the fact that there are no sequenced reference genomes thus far. Due to a well-developed tequila industry and a potential for bioethanol production, transcriptomics data have mostly been reported for *A. tequilana* leaves, and recently also for the *Agave* hybrid No. 11648 [(*A. amaniensis* × *A. angustifolia*) × *A. amaniensis*], used for fiber production (Table 1). Transcriptomics landscapes for other organs like roots have been less explored; however, *Agave* tolerance to drought could be directly controlled by physiology and morpho-anatomical traits of roots, modulated by gene expression. The *Agave* root system is fairly shallow (Nobel 1976) (Fig. 1a–c). This root system is originated from a short-lived primary root. At maturity, the system is composed by multiple shallow

✉ Fulgencio Alatorre-Cobos  
fulgencio@colpos.mx

Extended author information available on the last page of the article

**Table 1** Protocols used for total RNA isolation in several Agave species


Agave specie/tissue/analyses	Protocol used	References
<i>A. tequilana</i> Weber var. azul/pines and plants/gene cloning	Aurum Total RNA Mini kit (Bio-Rad, 7326820)	Ávila-Fernández et al. (2007)
<i>A. tequilana</i> /leaf, pine, and floral tissues/transcriptome sequencing	TRIzol (Thermo Fisher Scientific, 15596026) Concert (Thermo Fisher Scientific)	Martínez et al. (2010)
<i>A. tequilana</i> and <i>A. deserti</i> /roots, pine and leaf/transcriptome sequencing	TRIzol (Thermo Fisher Scientific, 15596026)	Gross et al. (2013)
Agave hybrid No. 11648/leaf/RT-PCR analysis	E.Z.N.A. <sup>®</sup> Plant RNA (Omega, R6827-01)	Gao et al. (2014)
<i>A. tequilana</i> , <i>A. victoriae-reginae</i> , and <i>A. striata</i> /roots, pine, leaf and floral tissues/transcriptome sequencing and RT-qPCR analysis	TRIzol (Thermo Fisher Scientific, 15596026)	Avila-de Dios et al. (2015)
<i>A. tequilana</i> /bulbs/RT-qPCR analysis	TRIzol (Thermo Fisher Scientific, 15596026)	Abraham et al. (2015)
<i>Agave americana</i> /meristem, rhizome, root, and pine/transcriptome sequencing	Spectrum™ Plant Total RNA (Sigma-aldrich, STRN10)	Abraham et al. (2016)
<i>Agave salmiana</i> /callus, leaves, and root/transcriptome sequencing and RT-qPCR analysis	TRIzol (Thermo Fisher Scientific, 15596026)	Cervantes-Pérez et al. (2018)
Agave hybrid No. 11648, <i>A. americana</i> , <i>A. deserti</i> , and <i>A. tequilana</i> /leaf/transcriptome sequencing	RNAprep Pure Plant Kit (Tiangen, DP432)	Huang et al. (2019a; b)
<i>Agave sisalana</i> /leaf/transcriptome sequencing	Phenol-based RNA isolation protocol	Batcho et al. (2019)
<i>Agave sisalana</i> /leaf/transcriptome sequencing	TRIzol (Thermo Fisher Scientific, 15596026)	Sarwar et al. (2019)
<i>A. tequilana</i> /shoot apical meristems and leaves/transcriptome sequencing and RT-qPCR analysis	TRIzol (Thermo Fisher Scientific, 15596026)	Avila-de Dios et al. (2019)
<i>A. tequilana</i> /leaf, roots, and pine/gene cloning	Trizol (Sigma-aldrich, T9424)	Sierra et al. (2019)
Agave hybrid No. 11648/pine and leaf/RT-qPCR analysis	RNAprep Pure Plant Kit (Tiangen, DP432)	Deng et al. (2019)
<i>Agave angustifolia</i> /leaf/RT-qPCR analysis	Trizol (Sigma-aldrich, T9424)	Sara et al. (2020)

roots that grow vertically and give rise to usually horizontal lateral roots, which are extended in a radial pattern (Dubrovsky 1997; Cruz-Ramírez et al. 2009). Connected to this network, it is usual to find rhizomes, which are not roots but underground stems. In *A. deserti*, alternated periods of wet conditions and drought induce suberization of the epidermis, hypodermis, and endodermis of young lateral roots (Huang and Nobel 1992), and lipid accumulation and lignification of the epidermis/exodermis in older roots (North and Nobel 1991, 1995). Cardiac and saponin glycosides have also been reported in *A. americana* roots (Kadam et al. 2012). Although RNA isolation is considered a routine protocol, for Agave roots, this step could be a bit tricky. Our laboratory experience with total RNA isolation from mature lateral roots of *A. tequilana* shows that lignin accumulates on out-layers and cortex (Fig. 1a–l) and can co-precipitate and yield low-quality RNA if they are not removed. Here, we report a detailed procedure for collecting and handling Agave root samples from soil and a protocol for total RNA extraction at micro-scale. This protocol is based on the standard guanidine thiocyanate-based method widely used for total RNA isolation in Agave species (Table 1), and yields total RNA with suitable quality for downstream applications such as gene expression analysis, cloning or sequencing.

## Materials and methods

### Plant material and lignin deposition analysis

Samples of lateral roots of four 6-year-old plants of *A. tequilana* Weber var. Azul were collected at field. Mature root samples for lignin deposition analysis or morphology by electronic scanning microscopy (SEM) was collected and kept in water or fixed in 4% paraformaldehyde solution, respectively. Lignin deposition was revealed by basic fuchsin staining (Sigma-Aldrich, 215597) as reported (Kapp et al. 2015) and imaged under a stereomicroscope (Leica M165 C, Leica, Germany) coupled with a Leica DFC290 camera, and confocal microscopy (Zeiss AXIO Imager.Z2, Carl Zeiss, Germany). For SEM analysis, root samples were gradually dehydrated with ethanol and frozen using liquid nitrogen. Micrographs were obtained using a JEOL JSM-6360LV Scanning Electron Microscope (Japan), equipped with a cold plate and an integrated digital image acquisition system. For RNA extraction, root samples collected in the field were immediately frozen in liquid nitrogen and stored at  $-80^{\circ}\text{C}$  until processed. All



**Fig. 1** Root system of *A. tequilana*. Root system of 2-year-old plant (a) and 6-year-old plants (b, c). In mature roots (red roots, d), high accumulation of lignin in cortex (e–g) and epidermis/exodermis (e, h, i) was revealed by staining with basic Fuchsin, and imaged under light microscopy (f) and confocal microscopy (g–i). This gradual

lignin deposition from young (j) to mature roots (k) alters viability out-layers cells. A slight lignin deposition (g) is also observed in vascular tissues (g, l). Diagram of a cross section of agave root was redrawn of North and Nobel 1995

common solutions and standard lab ware were prepared as recommended for routine RNA isolation protocols, and reagent setup as suggested by manufacture’s instructions.

### Recommendations for collecting and handling Agave root samples from soil

1. Agave plants are very susceptible to fungal pathogens if soil is not well drained; be careful to sample only healthy plants.

2. Clean harvested roots using dry paper towels and then cut them quickly into 5-cm long fragments using pruning shears.
3. When collecting mature roots, eliminate lignified epidermis/exodermis using tweezers as possible, and remove the red layer deposited on cortex using a paper towel soaked with absolute ethanol (clean roots). This red layer is highly soluble in absolute ethanol (Figs. 1, 2).
4. Place clean root fragments into a Falcon 50 ml conical tube and freeze using liquid nitrogen. Store at  $-80^{\circ}\text{C}$  until use.
8. Add 500  $\mu\text{l}$  isopropanol per 1 ml QIAzol reagent used in step 1. Mix thoroughly by inverting the tube 4–6 times. Do not mix by pipetting or vortexing.
9. Let the tube sit on the benchtop at room temperature for 10 min for RNA precipitation. Optional: RNA precipitation can occur at  $-20^{\circ}\text{C}$  overnight.
10. Centrifuge at  $12,000\times g$  for 10 min at  $4^{\circ}\text{C}$ .
11. Decant the supernatant and wash the pellet adding at least 1 ml of 75% ethanol, and vortex briefly (not more than 5 seg).
12. Centrifuge at  $7500\times g$  for 5 min at  $4^{\circ}\text{C}$ .
13. Decant the supernatant, spin the tube and remove the ethanol completely by pipetting. If necessary, repeat the spin and pipetting. Caution: be careful not to touch the pellet.
14. Dry the pellet by air for 15–20 min.
15. Re-suspend the pellet with DEPC-water.

### Total RNA extraction protocol (a workflow diagram is shown in Fig. 2)

This protocol enables up to 8–10 mini-extractions to be performed simultaneously. We set up the protocol using QIAzol lysis reagent (Qiagen, 79306), but other kits such as TRIzol reagent (Thermo Fisher Scientific) also work:

#### Total RNA isolation

1. Pour liquid nitrogen into a sterile porcelain mortar and pestle to pre-cool them. Repeat this step as necessary. Mortar and pestle will be ready for use when a small volume of liquid nitrogen remains at the bottom of the mortar.
2. Transfer the root sample collected to the pre-cooled mortar and grind to get a homogenous and fine powder. Tip: do not grind large volume of samples. Small samples will yield better results.
3. Transfer the powder to 1.5 ml RNase-free microcentrifuge tubes. If possible, weigh 100 mg aliquot per tube. Tips: the powder can be transferred using a sterile small spatula or 1 ml RNA-free tips; in both cases, spatula or tips have to be previously pre-cooled. If using tips, use one tip per tube to reduce contamination risks. Add powder to fill  $\frac{3}{4}$  volume of the conical space of the tube. This volume will be enough to get approximately 100 mg.
4. Add 1 ml of QIAzol lysis reagent per 100 mg of powder. Cap the microcentrifuge tube and immediately vortex until getting a defrosted homogenate. Keep at room temperature for 5 min. Stop: if necessary, keep the tubes at  $-80$  for maximum 1 week and then go to step 5.
5. Add 200  $\mu\text{l}$  chloroform. Cap the tube and vortex briefly. Place the tube on the benchtop at room temperature for 2–3 min.
6. Centrifuge at  $12,000\times g$  for 15 min at  $4^{\circ}\text{C}$ .
7. Transfer the upper aqueous phase into a fresh tube. Be careful to avoid the interphase to avoid contamination.

#### Quantification and analysis of quality of total RNA


Concentration and purity parameters (A260/280 and A260/230) of isolated RNA can be assessed using NanoDrop microvolume spectrometer and before standard electrophoresis analysis.

#### Total RNA purification by silica-based spin column

To remove DNA and other contaminants that could be co-precipitated during RNA isolation, a purification procedure is a must for downstream applications. Here, we purified Agave total RNA samples using the RNeasy<sup>®</sup> MinElute<sup>®</sup> Cleanup Kit (Qiagen, 74204) that is based on RNA-adsorption in silica spin columns. Below, the standard protocol is described and some suggestions are included to improve the performance of RNA purification and quality.

1. According to the manufacturer's suggestion, up to 45  $\mu\text{g}$  RNA can be purified per column. Tip: mix the previously isolated RNA stocks and up to 100  $\mu\text{l}$  with DEPC-water. Use preferentially DEPC-water from double-autoclaved aliquots.
2. Add 100  $\mu\text{l}$  cold 70% ethanol and mix well by pipetting 4–6 times.
3. Transfer immediately the mix to the RNeasy mini spin column (pink column) placed in a 2 ml collection tube (supplied by the manufacturer).
4. Centrifuge for 15 s at  $8000\times g$ . Discard the flow-through.
5. Add 700  $\mu\text{l}$  buffer RW1 to the spin column. Close the lid and centrifuge for 15 at  $8000\times g$ . Discard the flow-through.

**Fig. 2** Workflow diagram of the principal steps for total RNA isolation. Contamination by phenolic compounds through the protocol is illustrated for control roots (red roots)



6. Add 500  $\mu$ l buffer RPE to the spin column. Close the lid and centrifuge for 15 at  $8000 \times g$ . Discard the flow-through.

7. Add 500  $\mu$ l buffer RPE to the spin column. Close the lid and centrifuge for 1 min at  $8000 \times g$ . Discard the flow-through. Place the spin column in a fresh 2 ml

collection tube. CAUTION: some mini spin columns can have an internal circle shoulder just above the silica filter. Remove any ethanol or wash buffer trapped in this shoulder by pipetting.

8. Dry the spin column by centrifuging for 2 min at full speed. Keep the tube open during centrifuging. Optional: cut the cap if needed.
9. Discard the collection tube and place the spin column in new fresh 1.5 ml tube. Add 30–50  $\mu$ l RNase-free water directly to the spin column membrane. Let rest for 1 min on the benchtop. Tip: the manufacturer supplies RNase-free water or you could use previously prepared DEPC-water.
10. To elute the RNA, centrifuge for 1 min at full speed.
11. Hereafter, keep the RNA samples in ice. Quantify RNA concentration using a microvolume spectrometer and check RNA integrity by electrophoresis. Store purified total RNA samples at  $-20^{\circ}\text{C}$  until use.

### Evaluation of amenability for downstream applications

After spin column purification, 1  $\mu$ g total RNA from clean roots was used to synthesize cDNA using reverse transcriptase and oligo (dT) (10  $\mu$ M), following the manufacturer's protocol (SuperScript III, Invitrogen, 18080093). Specific primers for *CELLULOSE SYNTHASE7*-like gene (CESA7) forward: 5'-ATGTCCTTGATGGAGTCTGGTCTG-3' and reverse: 5'-TTAGCACTCCACTCCACACTGTTAAG-3') and for *SECONDARY WALL-ASSOCIATED NAC DOMAIN PROTEIN1*-like gene (SND1) (forward: 5'-ATGAGCATCACTGTGAATGGCCA-3' and reverse: 5'-TTATATTGAGGCGTGGGGCACG-3') were used to amplify full-length cDNAs using a proofreading DNA polymerase (Platinum SuperFi, Invitrogene, 12351010). The PCR program for CESA7-like gene was used as follows: 98  $^{\circ}\text{C}$  for 3 min; 35 cycles of 98  $^{\circ}\text{C}$  for 30 s, 62  $^{\circ}\text{C}$  for 30 s, and 72  $^{\circ}\text{C}$  for 1.5 min, and final extension at 72  $^{\circ}\text{C}$  for 5 min. The PCR program for SND1-like gene was used as follows: 98  $^{\circ}\text{C}$  for 3 min; 35 cycles of 98  $^{\circ}\text{C}$  for 30 s, 62  $^{\circ}\text{C}$  for 30 s, and 72  $^{\circ}\text{C}$  for 30 s, and final extension at 72  $^{\circ}\text{C}$  for 5 min. PCR products (5  $\mu$ l) were analyzed in agarose gels and visualized as described above.

### Statistical analyses and images processing


All data were analyzed as completely randomized factorial design. One-way independent ANOVA analyses were carried out test using SAS software (SAS Institute, USA) and Tukey's honestly significant difference (HSD) tests were applied to find out which treatment's means (compared with each other) are different ( $P \leq 0.05$ ). All images were prepared using Remove.bg (Kaleido AI GmbH, Austria),

Microsoft PowerPoint (Microsoft, USA) and GIMP 2.10.12 free software (The GIMP Development Team, USA).

## Results and discussion

High-quality RNA is a key prerequisite for successful gene cloning and expression analysis during molecular studies of plant development or responses to environmental clues. Although reliable protocols and a range of commercial kits are available, isolation of high-quality RNA is still challenging when plants are enriched with polysaccharides, polyphenols and others secondary metabolites. Most of these contaminants can degrade or co-precipitate with RNA due to physical interactions or similar chemical properties (Alemzadeh et al. 2005; Guan et al. 2019). This has led to optimize protocols previously reported according to plant species or even tissue-specific (Rubio-Piña and Zapata-Pérez 2011; Guan et al. 2019). In the case of Agave species, roots are organs with secondary metabolites, most likely polyphenols, and even lignified (Fig. 1). Our protocol for total RNA extraction of such roots simplifies the procedure of eliminating double precipitation by high LiCl concentrations or using polyvinylpyrrolidone (PVP) or polyethylene glycol (PEG), by doing just a suitable handling of Agave roots previous RNA isolation (Fig. 2). Removing lignified epidermis/exodermis and metabolites accumulated in the root cortex (clean roots) substantially reduced phenolic compounds that remained in the supernatant and precipitated with RNA isolated, as is shown by extractions of clean or control (red) roots processed with the protocol proposed here (Fig. 2). A comparison of RNA precipitation during 10 min at root temperature (RT) or at  $-20^{\circ}\text{C}$  overnight showed that RNA precipitation from polyphenol-rich samples provides better RNA yield when is carried out for short time (10 min) (Fig. 3a, b). Moreover, quality of isolated RNA can be increased if polyphenolic compounds are previously eliminated of the processed root samples (Fig. 3a, b). This statement is supported by our RNA quality data from red roots (control) or clean roots, which were determined by spectrophotometry (Fig. 3b, left), and confirmed by high-resolution automated electrophoresis (Agilent 2100 Bioanalyzer system) (Fig. 3c). Although purification by silica-based column can reduce the level of polyphenols/carbohydrates, in our case, RNAs isolated from lignified root (red) samples kept low levels of purity and quality even after purification as confirmed by RNA integrity number (RIN) values obtained by Bioanalyzer analysis (Fig. 3c). Taken together, these data indicate that removing lignified epidermis and metabolites accumulated in epidermis/exodermis and cortex is an essential step to get high-quality RNA from Agave roots, as reported in other tropical crops such as coffee (Huded et al. 2018). In our protocol, electrophoretic


**Fig. 3** Removing lignin and phenolic compounds in Agave roots is an essential step to get high-quality RNA. **a** Electrophoretic resolution patterns from ethidium bromide-stained gels and **b** yield parameters (quantity and quality) of total RNA samples isolated from control roots (red roots) and clean roots. Two RNA precipitation conditions (10 min at room temperature and  $-20^{\circ}\text{C}$  overnight) were also contrasted. Box and whisker plots represent mean value  $\pm$  SD. Asterisks indicate significant statistical differences between treatments determined by Tukey's HSD test ( $P \leq 0.05$ ). **c** Agilent bio-analyzer gel-like images of total RNA purified by spin column



patterns of RNAs, visualized by standard electrophoresis and microfluidics-based platform, suggest no RNA degradation but likely an irreversible binding of polyphenol-derived

quinones to RNA molecules as previously reported (Graham et al. 1993; Huded et al. 2018).

In *Arabidopsis thaliana*, *CELLULOSE SYNTHASE7* or *AtCESA7* (At5g17420) is involved in cellulose biosynthesis,



**Fig. 4** PCR amplification of full-length cDNAs of two agave genes related to cell wall using cDNA from total RNA of clean roots. Unique and abundant PCR amplicons were observed for the (a) CESA7-like gene and (b) SDN1-like gene, indicating no DNA con-

tamination in isolated RNA samples. c 1 kb plus (Thermo Fisher Scientific, 10787018) is shown to help to determine the PCR product length

while *SECONDARY WALL-ASSOCIATED NAC DOMAIN PROTEIN1* or *AtSDN1* (At1g32770) is a master modulator of lignin formation (reviewed in Handakumbura and Hazen 2012; McCahill and Hazen 2019). Here, we evaluated the suitability of RNA isolated from clean agave roots by amplifying the full-length cDNAs of CESA7-like gene and SDN1-like genes, found after mining Agave transcriptomes previously reported (Gross et al. 2013; Avila de Dios et al. 2019). Electrophoresis analysis showed a specific amplified PCR product of 3108 pb for CESA7, and a 1227 pb-length product for the like-SDN1 gene, as expected for amplification from full-length cDNAs (Fig. 4). No extra re-amplification PCR rounds were needed for the subsequent cDNA cloning.

## Conclusion

Our improved protocol for handling and RNA extraction allowed the isolation of high-quality total RNA from lignified and high-polyphenol accumulating roots from *A. tequilana*. This effective and simple protocol yielded suitable

RNA for synthesis of long cDNAs. Removing outer dead epidermis of roots did not compromise the root transcriptomic landscape isolated, as we were able to amplify full-length cDNAs related to secondary cell wall formation. The protocol developed was used successfully in *A. tequilana* roots but can be used for other Agave species that also develop lignified epidermis/exodermis in roots.

**Acknowledgements** This work was funded by CONACyT (Grant 1049) and Colpos. Luis F. Maceda-Lopez is grateful for a M.Sc. fellowship from CONACyT. Eleazar Garcia-Hernández and Dalia C. Moran Velazquez were supported by an undergraduate scholarship from CONACyT (Grant 1049). Authors especially grateful to Victor H. Fernandez-Carrillo, manager of Santa Genoveva farm, for supporting during Agave samples collecting.

**Author contributions** LFML, JLVA, DCMV, and FAC designed and performed the RNA extraction, quantification, and quality RNA analysis. JLVA, EGH, LRL, DHD, and MACV performed RT-PCR analysis. Lignin deposition and SEM analysis were carried out by LFML, SBAC, LRL, and FAC. Data mining was performed by LFML, EAD, JS, LT, ILR, EGC, and FAC. LFML, LT, EGC, ILR, JS, and FAC discussed results. LFML and FAC conceived and designed the research and wrote the paper. All authors read and approved the final manuscript.



## Compliance with ethical standards

**Conflict of interest** The authors declare that they have no conflict of interest in the publication.







**Ethical approval** Each of the authors confirms that this manuscript is original, has not been previously published, and is not currently under consideration by any other journal. Additionally, all of the authors have approved the contents of this paper and have agreed to the 3 Biotech's submission policies.

## References

- Abraham PE, Yin H, Borland AM et al (2016) Transcript, protein and metabolite temporal dynamics in the CAM plant *Agave*. *Nat Plants* 2:16178. <https://doi.org/10.1038/nplants.2016.178>
- Abraham-Juárez MJ, Hernández-Cárdenas R, Santoyo-Villa JN et al (2015) Functionally different PIN proteins control auxin flux during bulbil development in *Agave tequilana*. *J Exp Bot* 66:3893–3905. <https://doi.org/10.1093/jxb/erv191>
- Alemzadeh A, Fujie M, Usami S, Yamada T (2005) Isolation of high-quality RNA from high-phenolic tissues of eelgrass (*Zostera marina* L.) by keeping temperature low. *Plant Mol Biol Rep* 23:421. <https://doi.org/10.1007/BF02788892>
- Avila-de Dios E, Gomez-Vargas AD, Damian-Santos ML, Simpson J (2015) New insights into plant glycoside hydrolase family 32 in *Agave* species. *Front Plant Sci* 6:594. <https://doi.org/10.3389/fpls.2015.00594>
- Avila-de Dios E, Delaye L, Simpson J (2019) Transcriptome analysis of bolting in *A. tequilana* reveals roles for florigen, MADS, fructans and gibberellins. *BMC Genomics* 20:473. <https://doi.org/10.1186/s12864-019-5808-9>
- Ávila-Fernández Á, Olvera-Carranza C, Rudiño-Piñera E et al (2007) Molecular characterization of sucrose: sucrose 1-fructosyltransferase (1-SST) from *Agave tequilana* Weber var. azul. *Plant Sci* 173:478–486. <https://doi.org/10.1016/j.plantsci.2007.07.009>
- Batcho-Agossa A, Sarwar-Muhammad B, Tariq L et al (2019) Identification and characterisation of heat shock protein gene (HSP70) family and its expression in *Agave sisalana* under heat stress. *J Hort Sci Biotechnol* 1–13:470–482. <https://doi.org/10.1080/14620316.2019.1685412>
- Cervantes-Pérez SA, Espinal-Centeno A, Oropeza-Aburto A et al (2018) Transcriptional profiling of the CAM plant *Agave salmiana* reveals conservation of a genetic program for regeneration. *Dev Biol* 442:28–39. <https://doi.org/10.1016/j.ydbio.2018.04.018>
- Corbin-Kendall R, Byrt-Caitlin S, Bauer S et al (2015) Prospecting for energy-rich renewable raw materials: agave leaf case study. *PLoS ONE* 10:e0135382. <https://doi.org/10.1371/journal.pone.0135382>
- Cruz-Ramírez A, Calderón-Vázquez C, Herrera-Estrella L (2009) Effect of nutrient availability on root system development. *Annual Plant Reviews Volume 37: Root Development* 288–324. <https://doi.org/10.1002/9781444310023.ch11>
- Davis-Sarah C, Dohleman-Frank G, Long-Stephen P (2011) The global potential for *Agave* as a biofuel feedstock. *Gcb Bioenergy*. <https://doi.org/10.1111/j.1757-1707.2010.01077.x>
- Deng G, Huang X, Xie L et al (2019) Identification and expression of SAUR Genes in the CAM Plant *Agave*. *Genes* 10:555. <https://doi.org/10.3390/genes10070555>
- Dubrovsky JG (1997) Determinate primary-root growth in seedlings of Sonoran desert Cactaceae; its organization, cellular basis, and ecological significance. *Planta* 203:85–92. <https://doi.org/10.1007/s00050168>
- Escamilla-Treviño LL (2012) Potential of plants from the genus *Agave* as bioenergy crops. *BioEnergy Res* 5:1–9. <https://doi.org/10.1007/s12155-011-9159-x>
- Gao J, Yang F, Zhang S et al (2014) Expression of a hevein-like gene in transgenic *Agave* hybrid No. 11648 enhances tolerance against zebra stripe disease. *Plant Cell Tissue Organ Cult* 119:579–585. <https://doi.org/10.1007/s11240-014-0557-6>
- Good-Avila SV, Souza V, Gaut BS, Eguiarte LE (2006) Timing and rate of speciation in *Agave* (Agavaceae). *Proc Natl Acad Sci USA* 103:9124–9129. <https://doi.org/10.1073/pnas.0603312103>
- Graham GC (1993) A method for extraction of total RNA from *Pinus radiata* and other conifers. *Plant Mol Biol Rep* 11:32–37. <https://doi.org/10.1007/BF02670557>
- Gross SM, Martin JA, Simpson J et al (2013) *De novo* transcriptome assembly of drought tolerant CAM plants, *Agave deserti* and *Agave tequilana*. *BMC Genomics* 14:563. <https://doi.org/10.1186/1471-2164-14-563>
- Guan L, Ma X, Zhou X et al (2019) An optimized method to obtain high-quality RNA from cassava storage root. *3 Biotech* 9:118. <https://doi.org/10.1007/s13205-019-1608-9>
- Handakumbura PP, Hazen SP (2012) Transcriptional regulation of grass secondary cell wall biosynthesis: playing catch-up with *Arabidopsis thaliana*. *Front Plant Sci* 3:74. <https://doi.org/10.3389/fpls.2012.00074>
- Hernández-Castellano S, Garruña-Hernández R, Us-Camas R et al (2020) *Agave angustifolia* albino plantlets lose stomatal physiology function by changing the development of the stomatal complex due to a molecular disruption. *Mol Genet Genomics*. <https://doi.org/10.1007/s00438-019-01643-y>
- Huang B, Nobel PS (1992) Hydraulic conductivity and anatomy for lateral roots of *Agave deserti* during root growth and drought-induced abscission. *J Exp Bot* 43:1441–1449. <https://doi.org/10.1093/jxb/43.11.1441>
- Huang X, Xiao M, Xi J et al (2019a) *De novo* transcriptome assembly of *Agave* H11648 by Illumina sequencing and identification of cellulose synthase genes in *Agave* species. *Genes* 10:103. <https://doi.org/10.3390/genes10020103>
- Huang X, Xie L, Gbokie T et al (2019b) Transcriptome dataset of leaf tissue in *Agave* H11648. *Brown Univ Dig Addict Theory Appl* 4:62. <https://doi.org/10.3390/data4020062>
- Huded AKC, Jingade P, Mishra MK (2018) A rapid and efficient SDS-based RNA isolation protocol from different tissues of coffee. *3 Biotech* 8:183. <https://doi.org/10.1007/s13205-018-1209-z>
- Kadam PV, Yadav KN, Deoda RS et al (2012) Pharmacognostic and phytochemical studies on roots of *Agave Americana* (Agavaceae). *Intern J Pharmacogn Phytochem Res* 4:92–96
- Kapp N, Barnes WJ, Richard TL, Anderson CT (2015) Imaging with the fluorogenic dye Basic Fuchsin reveals subcellular patterning and ecotype variation of lignification in *Brachypodium distachyon*. *J Exp Bot* 66:4295–4304. <https://doi.org/10.1093/jxb/erv158>
- Martínez-Hernández A, Mena-Espino M, Herrera-Estrella AH et al (2010) Construcción de bibliotecas de ADNc y análisis de expresión génica por RT-PCR en agaves. *Revista latinoamericana de química* 38:21–44
- McCahill IW, Hazen SP (2019) Regulation of cell wall thickening by a medley of mechanisms. *Trends Plant Sci*. <https://doi.org/10.1016/j.tplants.2019.05.012>
- Niechayev NA, Jones AM, Rosenthal DM, Davis SC (2019) A model of environmental limitations on production of *Agave americana* L. grown as a biofuel crop in semi-arid regions. *J Exp Bot* 70:6549–6559. <https://doi.org/10.1093/jxb/ery383>
- Nobel PS (1976) Water relations and photosynthesis of a desert CAM plant, *Agave deserti*. *Plant Physiol* 58:576–582. <https://doi.org/10.1104/pp.58.4.576>
- North GB, Nobel PS (1991) Changes in hydraulic conductivity and anatomy caused by drying and rewetting roots of

- Agave deserti* (Agavaceae). *Am J Bot* 78:906–915. <https://doi.org/10.1002/j.1537-2197.1991.tb14494.x>
- North GB, Nobel PS (1995) Hydraulic conductivity of concentric root tissues of *Agave deserti* Engelm. under wet and drying conditions. *New Phytol* 130:47–57. <https://doi.org/10.1111/j.1469-8137.1995.tb01813.x>
- Rubio-Piña JA, Zapata-Pérez O (2011) Isolation of total RNA from tissues rich in polyphenols and polysaccharides of mangrove plants. *Electron J Biotechnol*
- Sara HC, René GH, Rosa UC et al (2020) *Agave angustifolia* albino plantlets lose stomatal physiology function by changing the development of the stomatal complex due to a molecular disruption. *Mol Genet Genomics* 295:787–805. <https://doi.org/10.1007/s00438-019-01643-y>
- Sarwar MB, Ahmad Z, Rashid B et al (2019) *De novo* assembly of *Agave sisalana* transcriptome in response to drought stress provides insight into the tolerance mechanisms. *Sci Rep* 9:396. <https://doi.org/10.1038/s41598-018-35891-6>
- Sierra-Gómez Y, Rodríguez-Hernández A, Cano-Sánchez P et al (2019) A biophysical and structural study of two chitinases from *Agave tequilana* and their potential role as defense proteins. *FEBS J* 286:4778–4796. <https://doi.org/10.1111/febs.14993>

## Authors and Affiliations

Luis F. Maceda-López<sup>1</sup>  · José L. Villalpando-Aguilar<sup>2</sup> · Eleazar García-Hernández<sup>1</sup> · Emmanuel Ávila de Dios<sup>3</sup> · Silvia B. Andrade-Canto<sup>4</sup> · Dalia C. Morán-Velázquez<sup>1</sup> · Lorena Rodríguez-López<sup>1</sup> · Demetrio Hernández-Díaz<sup>1</sup> · Manuel A. Chablé-Vega<sup>1</sup> · Laura Trejo<sup>5</sup>  · Elsa Góngora-Castillo<sup>6</sup>  · Itzel López-Rosas<sup>7</sup>  · June Simpson<sup>8</sup>  · Fulgencio Alatorre-Cobos<sup>7</sup> 

Luis F. Maceda-López  
maceda.luis@outlook.com

José L. Villalpando-Aguilar  
nachfornez@gmail.com

Eleazar García-Hernández  
garcia.eleazar@colpos.mx

Emmanuel Ávila de Dios  
aviladedios@gmail.com

Silvia B. Andrade-Canto  
sbac@cicy.mx

Dalia C. Morán-Velázquez  
moran.dalia@colpos.mx

Lorena Rodríguez-López  
lodoyarena96@gmail.com

Demetrio Hernández-Díaz  
demet2422@gmail.com

Manuel A. Chablé-Vega  
chable.manuel@colpos.mx

Laura Trejo  
laura.trejo@st.ib.unam.mx

Elsa Góngora-Castillo  
elsa.gongora@cicy.mx

Itzel López-Rosas  
itzel.rosas@colpos.mx

June Simpson  
june.simpson@cinvestav.mx

- <sup>1</sup> Colegio de Postgraduados Campus Campeche, Carretera Haltunchen-Edzna km 17.5, Sihochac, 24450 Campeche, Mexico
- <sup>2</sup> Tecnológico Nacional de México, Instituto Tecnológico de Chiná, Calle 11 entre 22 y 28, 24050 Chiná, Mexico
- <sup>3</sup> Departamento de Ingeniería Genética, Centro de Investigación y Estudios Avanzados, 3681 Irapuato, Mexico
- <sup>4</sup> Centro de Investigación Científica de Yucatán, A.C., Calle 43 No.130 x 32 y 34, Chuburna de Hidalgo, 97205 Mérida, Yucatán, Mexico
- <sup>5</sup> CONACYT Research Fellow, Laboratorio de Biodiversidad y Cultivo de Tejidos Vegetales, Instituto de Biología, UNAM, 90640 Santa Cruz Tlaxcala, Tlaxcala, Mexico
- <sup>6</sup> CONACYT Research Fellow, Centro de Investigación Científica de Yucatán, A.C., Calle 43 No.130 x 32 y 34, Chuburna de Hidalgo, 97205 Mérida, Yucatán, Mexico
- <sup>7</sup> CONACYT Research Fellow, Colegio de Postgraduados Campus Campeche, Carretera Haltunchen-Edzna km 17.5, Sihochac, 24450 Campeche, Mexico
- <sup>8</sup> Departamento de Ingeniería Genética, Centro de Investigación y Estudios Avanzados, 3681 Irapuato, Mexico