



# Complete Genome Sequence of Neonatal Clinical Group B Streptococcal Isolate CJB111

Brady L. Spencer,<sup>a</sup> Anushila Chatterjee,<sup>a</sup> Breck A. Duerkop,<sup>a</sup> Carol J. Baker,<sup>b,c</sup> Kelly S. Doran<sup>a</sup>

<sup>a</sup>Department of Immunology and Microbiology, School of Medicine, University of Colorado-Anschutz Medical Center, Aurora, Colorado, USA

<sup>b</sup>Department of Pediatrics, McGovern Medical School, University of Texas Health Science Center, Houston, Texas, USA

<sup>c</sup>Children's Memorial Hermann Hospital, Houston, Texas, USA

**ABSTRACT** Group B *Streptococcus* (GBS) is an asymptomatic colonizer of the female reproductive tract but can cause maternal and neonatal infections and adverse pregnancy outcomes. Here, we closed the genome sequence of strain CJB111, a neonatal GBS clinical isolate from a case of late-onset bacteraemia without focus (Houston, TX; 1990).

The Gram-positive beta-hemolytic bacterium *Streptococcus agalactiae* (group B *Streptococcus* [GBS]) asymptotically colonizes the gastrointestinal and female genital tracts of healthy adults but can cause neonatal infections (pneumonia, bacteraemia, meningitis [1]) and adverse pregnancy outcomes (2). Serotype V GBS isolates are emerging among adults and infants (3–6), and serotype V isolate CJB111 exhibits hypervirulence and vaginal persistence in murine models of GBS infection and colonization (7, 8).

CJB111 (ATCC BAA23) was isolated by Carol J. Baker from the blood of a female infant with late-onset sepsis on 18 July 1990, grown in Todd Hewitt broth (THB), and stored in glycerol at –90°C. The patient received intravenous ampicillin and gentamicin 1 day post onset of illness. Upon GBS isolation, therapy was switched on day 3 to intravenous penicillin (10-day treatment total). Following clean cerebrospinal fluid tests, she was diagnosed with bacteraemia without focus at age 55 days and discharged without apparent sequelae. While CJB111's sequence is currently available in 155 contigs (GenBank accession no. AAJQ01000000), a closed genome sequence may ease future genomic analyses.

CJB111 was grown statically overnight at 37°C in THB, genomic DNA was purified (Genta PureGene Yeast/Bact kit), and concentration and quality were confirmed by NanoDrop spectroscopy. The Microbial Genome Sequencing Center (MiGS; Pittsburgh, PA) performed short- and long-read sequencing (Illumina and Oxford Nanopore technologies [ONT], respectively) and *de novo* assembly. Default parameters were used except where otherwise noted. Short reads were obtained using the Illumina Nextera kit and NextSeq 550 platform (9). For ONT sequencing, libraries were prepared using kit SQK-LSK109 to the manufacturer's specifications (no DNA size selection/shearing), sequencing was performed on a MinION R9 flow cell, and base calling was performed using Guppy v4.2.2 (GPU mode) (10). Illumina paired-end reads (2 × 150 bp) and ONT long reads were provided as fastq files (Illumina: 7,410,044 reads, 989,364,400 bases, 472× coverage; ONT: 175,394 reads, 650,701,562 bases, 310× coverage,  $N_{50}$  value of 4,577 bp). bcl2fastq v2.20.0.422 was used for demultiplexing, quality control, and trimming of the Illumina reads (11) and Porechop v0.2.4 for quality trimming and removing adapters for ONT sequencing (12). Hybrid assembly via Unicycler v0.4.8 with a verbosity value of 2 (13) yielded six contigs, which were further assembled into three nonoverlapping contigs upon mapping to CJB111 contigs (AAJQ01000000) in Geneious v11.1.5 (14). The genome

**Citation** Spencer BL, Chatterjee A, Duerkop BA, Baker CJ, Doran KS. 2021. Complete genome sequence of neonatal clinical group B streptococcal isolate CJB111. *Microbiol Resour Announce* 10:e01268-20. <https://doi.org/10.1128/MRA.01268-20>.

**Editor** Julie C. Dunning Hotopp, University of Maryland School of Medicine

**Copyright** © 2021 Spencer et al. This is an open-access article distributed under the terms of the [Creative Commons Attribution 4.0 International license](#).

Address correspondence to Kelly S. Doran, kelly.doran@cuanschutz.edu.

**Received** 3 November 2020

**Accepted** 21 December 2020

**Published** 14 January 2021

**TABLE 1** Primers used in this study

PCRs flanking nonoverlapping contigs	Primer name <sup>a</sup>	Sequence (5'-3')	PCR $T_m$ (°C)	PCR product length (bp)	PCR extension time (min:s)
PCR 1	5121	TATCAATAACGATAGTATGCCAGCG	65	4,712	2:21
	3121	TTCCAATAGGTCTGATAGTGAGGTG			
PCR 2	5122	GTTTGTGCAGTCGTCTTATCTC	63	3,392	1:42
	3122	CGTCGGAATTAAATCTTGGAAATACC			
PCR 3	5123	GGCATCAGGAATGATCTGATTACAC	65	2,064	1:02
	3123	TGCCTCCATTGGATTACTGTATAC			
PCR 4	5124	GACTCGATAGGGTATATGGTGC	65	5,126	2:34
	3124	GGTTCGATTGCGTTACTGCG			
Sequencing primers					
PCR 1	5131	GTGACATAGATTGGAATAGGGTTAGC			
	5132	TATTCTCAGTGTCTGTACTACTGC			
	5133	AAATCTTGGCAGACAGTGGTTATC			
	5134	CAACAGGAGGAACCTGAGAAGTTC			
	5135	TACAATCCATCTCTGAAATTAC			
	3131	GTGACATAGATTGGAATAGGGTTAGC			
	3132	ATAATAAGGTGTCAGACAAACTCGC			
	3133	GGTTCGTCATTTATGAATGGTGTAC			
	3134	TTGACTATGGTTATGCTTCAGG			
	3135	TTCTCAACCTTGATTCTCTTTGG			
	3136	GTGCCGTTCAAAGGTGCGT			
	3137	CCGGGCTCGCTCATAGATAAG			
	5136	TATGCTCTCATAGGTAAACACCACC			
	5137	AACGATCACCTAACATTAGTACCTGC			
PCR 2	5138	TCTATCTGTTCTGTTCCTTG			
	5139	TTTAGGTTAGAAAGGAGACTGCC			
	5140	TACTTCAAATGGTATGCAAGCTATGG			
	3138	GCTGAACAAGCTGCTTATTGC			
	3139	TTTAGTTGAGGATGCTTATCGAG			
	3140	AGTTATCTGTCTATAAGGAATGTCG			
	3141	AAGCTATGGTTGAAGCTGTTG			
	5141	TTAAATTAACTCCTGAAGTACTCCG			
	3142	AGGTAATTCCATTCTCACCTGAAG			
	3143	TTTCGGGACAATTCTTGAACGTGAG			
PCR 4	5125	GTAACTAGTTATCTAGCCTGTAGC			
	5126	CACGAAAGCACTTAATCCGTG			
	5127	CCCTTGACTACATAAGTACTAACCC			
	5128	CTGTTATAATCAGCTCCATGAGC			
	5129	TTTCCCCTGCATTCCCATAGACC			
	5130	GCCTATCCAATTATTGCTTGGAG			
	3125	TTTACCTCTGTCATCCACAATC			
	3126	GCAAAGCAATTGATTCCGTCTT			
	3127	AAAGTGTGTTACCAACTCTGAAG			
	3128	AAATTATGAATCAGGCATGCTCCTGG			
	3129	AATAAGCCTGAAACCAGTCAGAG			
	3130	CATCACTTGGCCTCTATTATT			

<sup>a</sup> Primer names beginning with 5 indicate forward primers. Primer names beginning with 3 indicate reverse primers.

sequence was closed via PCR using primers flanking the nonoverlapping contigs (Table 1) and Phusion high-fidelity polymerase/buffer (New England Biolabs) under the following cycling conditions on a Bio-Rad T100 thermal cycler: 98°C, 2-min hot start; 34 cycles (98°C, 10<sup>seconds</sup>,  $T_m$  °C, 20<sup>seconds</sup>, 72°C, 30- second extension/kb); and 72°C, 10-min extension. Purified PCR amplicons (Qiagen) were Sanger sequenced using Applied Biosystems 3730/3500xl genetic analyzers, yielding 2× sequencing in both directions. Contigs were assembled *de novo*, overlapping ends were trimmed, and the genome sequence was rotated manually in Geneious v11.1.5 to start with *dnaA*.

The CJB111 sequence was deposited at GenBank as one circular contig (2,093,987 bp;

GC content, 35.52%). BUSCO\_v1 and CheckM v1.0.18 confirmed the genome completeness (15, 16). GenBank annotated the CJB111 genome sequence using PGAP v4.13 (17).

**Data availability.** The CJB111 sequence is available in GenBank under accession no. [CP063198](#). The raw sequence reads are accessible under Sequence Read Archive accession no. [SRX9273111](#), [SRX9273112](#), and [SRX9273113](#); BioProject accession no. [PRJNA663970](#); and BioSample accession no. [SAMN16191206](#).

## ACKNOWLEDGMENTS

This work was funded by NIH/NIAID F32 AI143203 to B.L.S., NIH/NINDS R01 NS116716 to K.S.D., and NIH/NIAID R01 AI141479 to B.A.D.

We acknowledge Dan Snyder and the Microbial Genome Sequencing Center (Pittsburgh, PA) for Illumina and Nanopore sequencing and Mark Farrell and the Barbara Davis Center Molecular Biology Core at University of Colorado-Anschutz Medical Campus (Aurora, CO) for Sanger sequencing.

## REFERENCES

- Verani JR, McGee L, Schrag SJ. 2010. Prevention of perinatal group B streptococcal disease: revised guidelines from CDC, 2010. *MMWR Recomm Rep* 59:1–32.
- Seale AC, Bianchi-Jassir F, Russell NJ, Kohli-Lynch M, Tann CJ, Hall J, Madrid L, Blencowe H, Cousens S, Baker CJ, Bartlett L, Cutland C, Gravett MG, Heath PT, Ip M, Le Doare K, Madhi SA, Rubens CE, Saha SK, Schrag SJ, Sobanjo-ter Meulen A, Vekemans J, Lawn JE. 2017. Estimates of the burden of group B streptococcal disease worldwide for pregnant women, stillbirths, and children. *Clin Infect Dis* 65:S200–S219. <https://doi.org/10.1093/cid/cix664>.
- A'Hearn-Thomas B, Khatami A, Randis TM, Vurayai M, Mokomane M, Arscott-Mills T, Banda FM, Mazhani T, Lepere T, Gaolebale P, Nchingane S, Chamby A, Gegick M, Suzman E, Steenhoff AP, Ratner AJ. 2019. High rate of serotype V *Streptococcus agalactiae* carriage in pregnant women in Botswana. *Am J Trop Med Hyg* 100:1115–1117. <https://doi.org/10.4269/ajtmh.18-0847>.
- Burcham LR, Spencer BL, Keeler LR, Runft DL, Patras KA, Neely MN, Doran KS. 2019. Determinants of group B streptococcal virulence potential amongst vaginal clinical isolates from pregnant women. *PLoS One* 14: e0226699. <https://doi.org/10.1371/journal.pone.0226699>.
- Yoon IA, Jo DS, Cho EY, Choi EH, Lee HJ, Lee H. 2015. Clinical significance of serotype V among infants with invasive group B streptococcal infections in South Korea. *Int J Infect Dis* 38:136–140. <https://doi.org/10.1016/j.ijid.2015.05.017>.
- Flores AR, Galloway-Pena J, Sahasrabhojane P, Saldana M, Yao H, Su X, Ajami NJ, Holder ME, Petrosino JF, Thompson E, Margarit YR, Rosini R, Grandi G, Horstmann N, Teatero S, McGeer A, Fittipaldi N, Rappuoli R, Baker CJ, Shelburne SA. 2015. Sequence type 1 group B *Streptococcus*, an emerging cause of invasive disease in adults, evolves by small genetic changes. *Proc Natl Acad Sci U S A* 112:6431–6436. <https://doi.org/10.1073/pnas.1504725112>.
- Patras KA, Doran KS. 2016. A murine model of group B *Streptococcus* vaginal colonization. *J Vis Exp* e54708. <https://doi.org/10.3791/54708>.
- Faralla C, Metruccio MM, De Chiara M, Mu R, Patras KA, Muzzi A, Grandi G, Margarit I, Doran KS, Janulczyk R. 2014. Analysis of two-component systems in group B *Streptococcus* shows that RgfAC and the novel FspSR modulate virulence and bacterial fitness. *mBio* 5:e00870-14. <https://doi.org/10.1128/mBio.00870-14>.
- Baym M, Kryazhimskiy S, Lieberman TD, Chung H, Desai MM, Kishony R. 2015. Inexpensive multiplexed library preparation for megabase-sized genomes. *PLoS One* 10:e0128036. <https://doi.org/10.1371/journal.pone.0128036>.
- Wick RR, Judd LM, Holt KE. 2019. Performance of neural network basecalling tools for Oxford Nanopore sequencing. *Genome Biol* 20:129. <https://doi.org/10.1186/s13059-019-1727-y>.
- Illumina. bcl2fastq: a proprietary Illumina software for the conversion of bcl files to basecalls. [https://support.illumina.com/sequencing/sequencing\\_software/bcl2fastq-conversion-software.html](https://support.illumina.com/sequencing/sequencing_software/bcl2fastq-conversion-software.html).
- Wick RR. 2017. Porechop. <http://github.com/rwick/Porechop>.
- Wick RR, Judd LM, Gorrie CL, Holt KE. 2017. Unicycler: resolving bacterial genome assemblies from short and long sequencing reads. *PLoS Comput Biol* 13:e1005595. <https://doi.org/10.1371/journal.pcbi.1005595>.
- Kearse M, Moir R, Wilson A, Stones-Havas S, Cheung M, Sturrock S, Buxton S, Cooper A, Markowitz S, Duran C, Thaler T, Ashton B, Meintjes P, Drummond A. 2012. Geneious Basic: an integrated and extendable desktop software platform for the organization and analysis of sequence data. *Bioinformatics* 28:1647–1649. <https://doi.org/10.1093/bioinformatics/bts199>.
- Simao FA, Waterhouse RM, Ioannidis P, Kriventseva EV, Zdobnov EM. 2015. BUSCO: assessing genome assembly and annotation completeness with single-copy orthologs. *Bioinformatics* 31:3210–3212. <https://doi.org/10.1093/bioinformatics/btv351>.
- Parks DH, Imelfort M, Skennerton CT, Hugenholtz P, Tyson GW. 2015. CheckM: assessing the quality of microbial genomes recovered from isolates, single cells, and metagenomes. *Genome Res* 25:1043–1055. <https://doi.org/10.1101/gr.186072.114>.
- Tatusova T, DiCuccio M, Badretdin A, Chetvernin V, Nawrocki EP, Zaslavsky L, Lomsadze A, Pruitt KD, Borodovsky M, Ostell J. 2016. NCBI Prokaryotic Genome Annotation Pipeline. *Nucleic Acids Res* 44:6614–6624. <https://doi.org/10.1093/nar/gkw569>.