


ORIGINAL ARTICLE

Comparison of the analytical performance between the Oncomine Dx Target Test and a conventional single gene test for epidermal growth factor receptor mutation in non-small cell lung cancer

Tadashi Sakaguchi¹ , Akemi Iketani², Kazuki Furuhashi¹, Yuki Nakamura¹, Yuta Suzuki¹, Kentaro Ito¹, Kentaro Fujiwara¹, Yoichi Nishii¹, Koji Katsuta², Osamu Taguchi¹ & Osamu Hataji¹

¹ Respiratory Center, Matsusaka Municipal Hospital, Matsusaka, Japan

² Pathology Department, Matsusaka Municipal Hospital, Matsusaka, Japan

Keywords

EGFR mutations; next-generation sequencing; non-small cell lung cancer; Oncomine Dx Target Test; PNA-LNA PCR clamp.

Correspondence

Tadashi Sakaguchi, Respiratory Center, Matsusaka Municipal Hospital, 1550, Tonomachi, Matsusaka, Mie 515-0073, Japan.
Tel: +81 598 23 1515
Fax: +81 598 21 8790
Email: mchsakaguchi@city-hosp.matsusaka.mie.jp

Received: 16 October 2020;

Accepted: 17 November 2020.

doi: 10.1111/1759-7714.13767

Thoracic Cancer **12** (2021) 462–467

Abstract

Background: Next-generation sequencing (NGS) has been implemented in clinical oncology to analyze multiple genes and to guide targeted therapy; however, little is known about the performance of the Oncomine Dx Target Test compared with conventional single gene tests for detecting *EGFR* mutations. The objective of this study was to evaluate the performance of the Oncomine Dx Target Test compared with a PNA-LNA PCR clamp test to detect *EGFR* mutations.

Methods: We retrospectively reviewed consecutive patients with non-small cell lung cancer (NSCLC) from whom FFPE samples were simultaneously submitted for the Oncomine Dx Target Test, and a PNA-LNA PCR clamp test using the same specimen. We subsequently compared the analysis success rates and detection rates between the two tests.

Results: A total of 116 samples were identified. The success rates and detection rates of *EGFR* mutations in the total number of samples were 90% and 28%, respectively for the Oncomine Dx Target Test, and 100% and 35% for the PNA-LNA PCR clamp test. The Oncomine Dx Target Test was unable to analyze three samples (2%) due to the samples not passing the nucleic acid concentration threshold, and nine (8%) samples had invalid results. The exon 19 deletion was not detected by the Oncomine Dx Target Test in four cases (4%).

Conclusions: The analytical performance of the Oncomine Dx Target Test analysis for *EGFR* mutations may not be comparable with conventional single gene tests due to both invalid and false-negative results.

Key points

Significant findings of the study

- The success rate of the Oncomine Dx Target Test was significantly lower than the PNA-LNA PCR clamp test.
- Among the samples successfully analyzed, four exon 19 deletions were not detected by the Oncomine Dx Target Test.

What this study adds

- The analytical performance of the Oncomine Dx Target Test may not be comparable with conventional single gene tests.
- We should revise the sampling procedures, and review the sample quality assessment methods, to improve the analytical performance.

Introduction

Targeted therapies for non-small cell lung cancer (NSCLC) patients harboring driver oncogene alterations have been proven to have promising antitumor activities, and are generally recommended in the clinical guidelines.^{1–3} Epidermal growth factor receptor-tyrosine kinase inhibitors (EGFR-TKIs) play an important role in effective therapy, especially for NSCLC patients in East Asia.⁴

Conventional single-gene tests, such as the cobas *EGFR* assay, the therascreen *EGFR* assay as an in vitro diagnostic (IVD) test, and the peptide nucleic acid-locked nucleic acid polymerase chain reaction (PNA-LNA PCR) clamp assay as a laboratory-developed test (LDT) for *EGFR* mutations, have been conducted to select suitable patients responsive to EGFR-TKI.^{5–7} Recently, as more gene targets have been identified and recognized in clinical settings for NSCLC, more biomarker tests have been required to make appropriate treatment decisions. As the number of single-gene tests has increased, tissue consumption for analysis of multiple single-gene tests has also increased, and the completion rates of the ordered tests has decreased.⁸

Next-generation sequencing (NGS) can detect multiple gene variants simultaneously enabling comprehensive genetic testing. The Oncomine Dx Target Test (Ion Torrent PGM Dx Sequencer; Thermo Fisher Scientific) is an NGS panel for NSCLC, approved by the US Food and Drug Administration in June 2017.⁹ It is a qualitative, in vitro diagnostic test that uses high-throughput parallel sequencing technology to detect sequence variations in 46 genes on DNA and RNA isolated from formalin-fixed and paraffin-embedded (FFPE) specimens using the Ion PGM Dx System. In February 2019, the test was approved for use in Japan as a companion diagnostic for targeted therapies on four driver mutations; *EGFR*, anaplastic lymphoma kinase (*ALK*), ROS proto-oncogene 1, receptor tyrosine kinase (*ROS1*), and B-Raf proto-oncogene, serine/threonine kinase (*BRAF*) (p.V600E).

The Oncomine Dx Target Test is expected to analyze multiple gene targets reliably and simultaneously; however, the performance of the Oncomine Dx Target Test compared with conventional single gene tests for detecting *EGFR* mutations has not previously been fully evaluated. Therefore, in this study, we retrospectively evaluated the performance of the Oncomine Dx Target Test compared with the PNA-LNA PCR clamp test, in detecting *EGFR* mutations in patients with NSCLC.

Methods

Patient selection

This retrospective study was conducted at Matsusaka Municipal Hospital, Japan. We reviewed the electronic data

from consecutive NSCLC patients whose FFPE samples were simultaneously submitted for an Oncomine Dx Target Test, and a PNA-LNA PCR clamp test, using the same specimen, from August 2019 to February 2020. Samples collected in other hospitals, and archived samples were excluded. Clinical data assessments included patient characteristics, sampling methods, staging, histology, pathological findings, and the results of genetic tests. This study was approved by the institutional review board of Matsusaka Municipal Hospital (Approval date: 20 April 2020; IRB number J-76-200 410-5-2). Informed consent was obtained by the opt-out method.

Sample processing and genetic tests

Small tissue samples collected by endobronchial biopsy/transbronchial biopsy (EBB/TBB), computed tomography-guided percutaneous needle biopsy (CTNB) and fine needle aspiration (FNA) were immediately placed in 10% neutral buffered formalin (NBF) and fixed over 12 to 18 hours at room temperature. In cases of surgical lung resection, the samples from limited resections, including lung segmentectomy and wedge resections, were immediately placed in 10% NBF after sampling for intraoperative rapid diagnosis (IRD) for 24 to 48 hours at room temperature. The specimens from lung lobectomies performed up to January 2020 were stored in a refrigerator at 4°C for less than three hours, after sampling for IRD, then placed in 10% NBF for 24 to 48 hours at room temperature. Meanwhile, using samples from lung lobectomies performed after January 2020, a tissue sample of the correct size (10 mm × 10 mm), and enriched in tumor cells was prepared for the Oncomine Dx Target Test concurrently when sampling for IRD, and placed in 10% NBF immediately after sampling. Formalin-fixed tissues underwent serial processing and were then embedded in paraffin to create FFPE blocks. The amount of tumor cells, and tumor content of the sample stained with hematoxylin and eosin, were evaluated by skilled cytopathologists. Trimming, which eliminates samples containing little or no tumor cells in small biopsy samples, was performed, but macro- and microdissection were not performed in our institution. If the tumor content was <20% after trimming of small biopsy samples, or the amount of tumor cells was insufficient, the sample was not submitted for the Oncomine Dx Target Test. For the Oncomine Dx Target Test, 10 to 15 slide-mounted 5 µm sections of small biopsy samples and 5 to 10 slide-mounted 5 µm sections of surgical resection samples were submitted to LSI Medience Laboratories (Tokyo, Japan). For the PNA-LNA PCR clamp test, five slides of 5 µm sections from small biopsy samples and surgical resection samples were submitted. LSI Medience Laboratories performed the Oncomine Dx Target Tests

based on Thermo Fisher's Ion AmpliSeq technology, and the PNA-LNA PCR clamp tests were performed using the PNA-LNA PCR clamp assay.^{7,9} The Oncomine Dx Target Test panel included 46 genes (Table S1).

Specific EGFR mutations detectable by each test

Specific EGFR mutations detectable by the Oncomine Dx Target Test were identified in exon 18 (p.G709X and p.G719A/C/S/D), exon 19 (deletion), exon 20 (p.S768I and p.T790M), and exon 21 (p.L858M, p.L858R and p.L861Q/R) and other mutations in exon 18 to 21. Exon 20 insertion has been included in Oncomine Dx Target Test reports since May 2020 in Japan. Specific EGFR mutations detectable by PNA-LNA PCR clamp test were identified in exon 18 (p.G709X and p.G719A/C/S), exon 19 (deletion and insertion), exon 20 (p.S768I, p.T790M and insertion), and exon 21 (p.L833X, p.L858R and p.L861Q) and other mutations in exon 18 to 21 (p.V769M, p.V834L, p.K860I, etc.). p.S768I, exon 20 insertion and p.L833X have been included in PNA-LNA PCR clamp test reports since October 2019.

Outcomes

We evaluated the success rate and detection rate of EGFR mutations for the Oncomine Dx Target Test compared with those of the PNA-LNA PCR clamp test. The analysis results were regarded as successful if all the results in the following EGFR mutations reported for each time period were completely available; exon 18 (p.G719A/C/S), exon 19 (deletion), exon 20 (p.S768I, p.T790M, insertion), and exon 21 (p.L858R and p.L861Q), which are the mutations considered required for detection due to clinical

implications by the Japanese Lung Cancer Society. Meanwhile, the analysis results were regarded as unsuccessful if the sample did not pass the nucleic acid concentration threshold, or if one or more EGFR sequence results mentioned above were invalid due to a failure to meet the DNA sample quality control (QC) metrics, or no call. The detection rates of EGFR mutation were calculated as the rates of samples with detected EGFR mutations in all samples, or samples diagnosed with adenocarcinoma. We also evaluated concordance rate, and the clinical settings of patients with the discordant results between the two tests. The concordance rate was calculated as follows: concordance rate = N_c/N_t , where N_c was the number of samples showing the concordant results, and N_t was the total number of successfully analyzed samples.

Statistical analysis

Statistical analyses were performed using Pearson's Chi-squared test for comparison of analysis success rates and detection rates. *P*-values less than 0.05 indicated statistical significance. Statistical analyses were performed using SPSS software, version 26.0 (SPSS Inc., Chicago, USA).

Results

Sample characteristics

A total of 116 samples were identified for comparison analysis. The sample characteristics for the analysis are shown in Table 1. The main sampling methods were EBB/TBB (60%), surgical resection (21%), and CTNB (16%). More than half of the tumor histology was adenocarcinoma (64%), followed by squamous cell carcinoma (29%).

Success and detection rates for EGFR mutations for each test

The success rates for the Oncomine Dx Target Test and PNA-LNA PCR clamp test are shown in Table 2. Although the success rate of the Oncomine Dx Target Test was 90%, the success rate of the PNA-LNA PCR clamp test was 100% ($P < 0.01$). EGFR mutation results from the Oncomine Dx Target Test could not be assessed in three samples (2%) due to a failure to pass the nucleic acid concentration threshold. Nine samples returned invalid results; seven samples (6%) generated invalid results due to a failure to meet the DNA sample QC metrics, or no call for all EGFR mutations specified above, and two samples (2%) generated invalid results for the subset of the mutations due to no call.

The detection rates for EGFR mutations of each test for all samples, and adenocarcinoma, are shown in Fig 1.

Table 1 Sample characteristics

	Total samples	
	<i>N</i> = 116	(%)
Sampling method		
EBB/TBB	69	60%
Surgical resection	24	21%
CTNB	19	16%
Others	4	3%
Histology		
ADC	74	64%
Sq	34	29%
Non-Sq Non-ADC	7	6%
NSSC NOS	1	1%

ADC, adenocarcinoma; CTNB, computed tomography-guided needle aspiration; EBB, endobronchial biopsy; NSSC NOS, non-small cell carcinoma, not otherwise specified; Sq, squamous cell carcinoma; TBB, transbronchial biopsy.

Table 2 Analysis success rates of Oncomine Dx Target Test and PNA-LNA PCR clamp test

	Total samples <i>N</i> = 116			
	Oncomine Dx Target Test	(%)	PNA-LNA PCR clamp	(%)
Success of analysis	104	90%	116	100%
Not passing the nucleic acid concentration threshold	3	3%	0	0%
Invalid results for all <i>EGFR</i> mutations	7	6%	0	0%
Invalid results for subset of <i>EGFR</i> mutations	2	2%	0	0%

EGFR, epidermal growth factor receptor.

Although the detection rate of the Oncomine Dx Target Test was 28%, the rate of the PNA-LNA PCR clamp test was 35% for all samples ($P = 0.32$). For adenocarcinoma, the detection rates were 43% for the Oncomine Dx Target Test, and 53% for the PNA-LNA PCR clamp test ($P = 0.32$). Among all samples, 36 common mutations and six uncommon mutations, including two compound mutations were detected by the PNA-LNA PCR clamp test. Three samples with exon 21 p.L858R could not be detected by the Oncomine Dx Target Test due to a failure to pass

the nucleic acid concentration threshold and invalid results. Among the samples successfully analyzed by each test, five discordant results were reported, and the concordance rate was 95%. The Oncomine Dx Target Test failed to detect four exon 19 deletion mutations that were detected with the PNA-LNA PCR clamp test. Although the Oncomine Dx Target Test reported a negative for *EGFR* mutation in one sample where exon 20 insertion was detected by the PNA-LNA PCR clamp test, the Oncomine Dx Target Test did not include exon 20 insertion as a




Figure 1 The detection rates for *EGFR* mutations. (a) Detection rate of the Oncomine Dx Target Test by all histology and adenocarcinoma. (b) Detection rate of the PNA-LNA PCR clamp test by all histology and adenocarcinoma.

target mutation to be reported when the sample was analyzed.

Clinical information of discordant cases

The four discordant cases of exon 19 deletion between the Oncomine Dx Target Test and the PNA-LNA PCR clamp test are detailed in Table 3. Although all samples were small biopsy samples collected by EBB/TBB, these had a sufficient amount of tumor cells and content to be viable. Tumor histology of three cases were adenocarcinoma, and one was large cell carcinoma (LCC). One patient had disease progression early after the administration of EGFR-TKI, and another patient who was administered EGFR-TKI could not have their response assessed due to an adverse event of pneumonitis. Two early stage patients underwent surgery and were not administered EGFR-TKI.

Discussion

This is the first report to indicate that the Oncomine Dx Target Test may not be comparable to conventional single gene testing for the detection of *EGFR* mutations. Among patients harboring *EGFR* mutations detected by the PNA-LNA PCR clamp test, the Oncomine Dx Target Test was unable to detect seven *EGFR* mutations (6%). All were common mutations, specifically; one harboring exon 21 p.L858R due to not passing the nucleic acid concentration threshold, two harboring exon 21 p.L858R due to invalid results, and four harboring exon 19 deletion due to false negative results. In the current Japanese medical insurance system, the samples which fail analysis by the Oncomine Dx Target Test due to not passing the nucleic acid concentration threshold can be submitted for conventional single gene tests; however, the samples which fail due to invalid results cannot be submitted for conventional single gene tests using medical insurance in clinical practice. Missing

suitable patients responsive to EGFR-TKIs should be avoided, especially in a high prevalence area for *EGFR* mutations.⁴

With regard to the cases with discordant results, excluding the sample where exon 20 insertion was detected by the PNA-LNA PCR clamp test, all samples were collected by EBB/TBB and exon 19 deletion could not be detected by the Oncomine Dx Target Test in our study. Two factors were considered to cause the difference in results between the tests; one is the lower sensitivity of the Oncomine Dx Target Test compared with the PNA-LNA PCR clamp test, and the other is the difference in the range of variants reported by each test. The estimated limit of detections (LODs) of the Oncomine Dx Target Test has been previously reported to be 6% allele frequency (AF) for *EGFR* exon 19 deletion.¹⁰ A clinical bridging study to establish the assurance of the Oncomine Dx Target Test compared with the theascreen *EGFR* assay, whose LODs were 1%–2% for *EGFR* common mutations, was performed by Thermo Fisher Scientific and showed good concordance.¹⁰ Meanwhile, the LODs for the PNA-LNA PCR clamp test have been estimated to be 0.1%–1%⁷; therefore, it is expected that the PNA-LNA PCR clamp test is able to detect a lower concentration of *EGFR* mutation. In cases 1 and 2, the Catalogue of Somatic Mutations in Cancer (COSMIC) IDs of the exon 19 deletions identified by the PNA-LNA PCR clamp test, were variants that are detectable by the Oncomine Dx Target Test. Therefore, the cause of the discordant results for these two cases would likely be the lower sensitivity of the Oncomine Dx Target Test. The COSMIC IDs in cases 3 and 4 could not be identified; therefore, the cause of the discordant results was unknown. When the results of the Oncomine Dx Target Test for samples with high tumor content were false negative due to low sensitivity, the clinical meaning of detecting these genetic abnormalities was unclear because the tumors harboring the mutation may not have been dominant in the

Table 3 Clinical information of discordant cases

No	Age	Sex	Sampling method	Histology	Result of <i>EGFR</i> mutation		Stage (UICC-eighth)	Treatment	Response to EGFR-TKI
					Oncomine Dx Target Test	PNA-LNA PCR clamp			
1	47	M	TBB	ADC	Negative	Detected exon 19 deletion	IVB	Osimertinib	PD
2	70	M	EBB	LCC	Negative	Detected exon 19 deletion	Postoperative recurrence	Osimertinib	NA due to pneumonitis
3	68	M	TBB	ADC	Negative	Detected exon 19 deletion	IA3	Operation	Not administered
4	68	F	TBB	ADC	Negative	Detected exon 19 deletion	IB	Operation	Not administered

ADC, adenocarcinoma; EBB, endobronchial biopsy; EGFR, epidermal growth factor receptor; LCC, large cell carcinoma; NA, not assessed; PD, progression disease; TBB, transbronchial biopsy; TKI, tyrosine kinase inhibitor; UICC, Union for International Cancer Control.

whole tumor.¹¹ Meanwhile, when the false negative results due to low sensitivity come from low tumor content, trimming or macro- and microdissection may be useful to eliminate nontumor cells, and EGFR-TKI would be effective on the tumors. Although the sample in Case 1 had enough tumor content, Case 1 did not respond to EGFR-TKI. We speculate that this was due to the sample being taken from a heterogeneous tumor, in which the areas that harbored the mutation did not form the dominant portion of the tumor. This resulted in the tumor sample failing to meet the sensitivity limit of the Oncomine Dx Target Test, the discordant results between the Oncomine Dx Target Test and the PNA-LNA PCR clamp test, and the lack of response to EGFR-TKI.

There were several limitations to this study. First, this was a relatively small retrospective study. Although the results of this study suggested that the Oncomine Dx Target Test had a risk of false negatives for *EGFR* mutations, especially exon 19 deletion variants, more samples would be needed to evaluate the frequency of false negatives for each variant and the clinical meaning of detecting the mutations. Second, this study was conducted in a single institute, and therefore the results of this study may not be applicable to other institutions because the methods of tissue sampling, sample preparation process, and judgment of whether or not to submit a sample for Oncomine Dx Target Test, vary in each institution. Finally, we did not perform macro- and microdissection; however, macro- and microdissection would further improve the tumor content and decrease the false negative rate.

In conclusion, the success rate of Oncomine Dx Target Test analysis is generally well tolerated; however, we should consider the risk of missing *EGFR* mutations compared with conventional single gene testing due to invalid and false-negative results.

Acknowledgments

This research did not receive any specific grant from funding agencies in the public, commercial, or not-for-profit sectors. The authors thank the staff of Matsusaka Municipal Hospital for their advice and assistance. The authors also thank Mr Shane LeGros for reviewing the language of this article.

Disclosure

Matsusaka Municipal Hospital received research grant funding from Novartis, GlaxoSmithKline, AstraZeneca, Daiichi Sankyo, Bayer, and Boehringer Ingelheim. K. Ito has received speaker fees as honoraria from Eli Lilly Japan, Chugai, AstraZeneca, MSD, Boehringer Ingelheim Japan, Ono, and Pfizer Japan. O. Taguchi received speaker fees as

honoraria from AstraZeneca. O. Hataji received speaker fees as honoraria from Novartis Pharma, AstraZeneca, and Boehringer Ingelheim Japan. The remaining authors declare no conflict of interest.

References

- 1 Ettinger DS, Wood DE, Aggarwal C *et al.* NCCN clinical practice guidelines in oncology: Non-small cell lung cancer. Version 7.2019. 2019. To view the most recent version, visit NCCN.org Accessed 30 August 2019.
- 2 Planchard D, Popat S, Kerr K *et al.* Metastatic non-small cell lung cancer: ESMO clinical practice guidelines for diagnosis, treatment and follow-up. *Ann Oncol* 2018; **29**: iv192–237.
- 3 Akamatsu H, Ninomiya K, Kenmotsu H *et al.* The Japanese lung cancer society guideline for non-small cell lung cancer, stage IV. *Int J Clin Oncol* 2019; **24**: 731–70.
- 4 Zhang YL, Yuan JQ, Wang KF *et al.* The prevalence of *EGFR* mutation in patients with non-small cell lung cancer: A systematic review and meta-analysis. *Oncotarget* 2016; **7**: 78985–93.
- 5 Kimura H, Ohira T, Uchida O *et al.* Analytical performance of the cobas *EGFR* mutation assay for Japanese non-small-cell lung cancer. *Lung Cancer* 2014; **83**: 329–33.
- 6 Therascreen *EGFR* RGQ PCR Kit Handbook Version 1. <https://www.qiagen.com/ch/resources/download.aspx?id=db794cae-999b-4362-aba3-455ebfd807a5&lang=en>.
- 7 Nagai Y, Miyazawa H, Tanaka T *et al.* Genetic heterogeneity of the epidermal growth factor receptor in non-small cell lung cancer cell lines revealed by a rapid and sensitive detection system, the peptide nucleic acid-locked nucleic acid PCR clamp. *Cancer Res* 2005; **65**: 7276–82.
- 8 Yu TM, Morrison C, Gold EJ *et al.* Multiple biomarker testing tissue consumption and completion rates with single-gene tests and investigational use of Oncomine dx target test for advanced non-small-cell lung cancer: A single-center analysis. *Clin Lung Cancer* 2018; **20**: 20–9.
- 9 Food and Drug Administration. Oncomine™ Dx Target Test Part I: Sample Preparation and Quantification User Guide. Revision C.0 2017. 2017 https://www.accessdata.fda.gov/cdrh_docs/pdf16/P160045C.pdf. Accessed 14 June 2020.
- 10 Food and Drug Administration. Summary of Safety and Effectiveness Data 2017. 2017. https://www.accessdata.fda.gov/cdrh_docs/pdf16/P160045B.pdf. Accessed April 8, 2018.
- 11 Jamal-Hanjani M, Wilson GA, McGranahan N *et al.* Tracking the evolution of non-small-cell lung cancer. *N Engl J Med* 2017; **376**: 2109–21.

Supporting Information

Additional Supporting Information may be found in the online version of this article at the publisher's website:

Table S1. List of genes searched for in the Oncomine Dx Target Test.