



Published in final edited form as:

Methods Mol Biol. 2021 ; 2193: 111–120. doi:10.1007/978-1-0716-0845-6_11.

Wound Matrix Stiffness Imposes on Macrophage Activation

Pu Duann¹, Pei-Hui Lin^{2,3}

¹Research and Development, Salem Veteran Affairs Medical Center, Salem, USA

²Davis Heart and Lung Research Institute, The Ohio State University, Columbus, USA

³Department of Surgery, The Ohio State University, Columbus, USA

Abstract

The immune system depends on two major paths—the innate and the adaptive immunity. Macrophage, with its unique features as the first line of immune defense to engulf and digest invaders, serves as the key effector cells integrating those two paths. The dynamic plasticity of macrophage activation during wound repair, inflammation resolution, and tissue remodeling are emerging biomedical and bioengineering hot topics in immune function studies such as the various secretions of cytokines and chemokines and the signaling pathways with ligands and their cognate receptors. Better knowledge on how physical/mechanical and multicellular microenvironment on the modulation of macrophage functions will create innovative therapies to boost host defense mechanism and assist wound healing. In this, we describe an easy method to measure functions (gene expressions) of human and mouse macrophages in response to mechanical microenvironment changes by embedding isolated macrophages in polymerized hyaluronan gel with different wound matrix stiffness.

Keywords

Hyaluronan; Extracellular matrix; Matrix stiffness; Cytokines; Chemokines; Gene expression; PBMC

1 Introduction

Macrophages contribute to diverse biology including tissue development, angiogenesis, inflammation, organ injury, and repair [1, 2]. These multifunctional macrophages phenotypes are due to their dynamic adoption/polarization in responding to signals from their microenvironments [3]. Macrophages play essential roles in all phases of wound healing process, such as boosting host immune defense from microbes via promoting and resolving inflammation through inflammatory cytokines secretions, dead cells removal, supporting epithelial cell proliferation, and tissue remodeling (scar formation) after a wound occurs [4]. Macrophages respond to chemical and physical cues upon wound healing. As such, dysregulated macrophages-mediated inflammation could also lead to tissue destruction. The tissue microenvironments such as matrix property of wound stiffness have recently opened up new research avenues and therapeutic opportunities in wound healings.

Hyaluronan (HA, hyaluronic acid) is a polyanionic, nonsulfated glycosaminoglycan composed of long chain disaccharide units of GlcNAc and D-glucuronic acid polymer that repeat up to 20,000 or more times to form a linear cable with length > 20 μm and high molecular mass (>8 million Da) [5, 6]. In fact, the wide range of molecular sizes of HA dictates its different biological functions and cellular signaling [7]. HA is a major extracellular matrix (ECM) carbohydrate polymer widely distributed throughout human body and enriched in many tissues including connective, epithelial, and neural tissues, joints, and skin. In normal human body, HA constitutes 0.02% body weight (15 g in 70 kg body weight) with rapid turnover as approximately a third of the body's HA turns over daily [5]. HA is involved in many physiological and pathobiological processes such as fertilization, cell proliferation, tissue development, immune response, wound healing, atherosclerosis, angiogenesis, tissue fibrosis, and cancers [8, 9]. These versatile biological functions are related to its unique biosynthesis process. HA is formed directly within the pericellular undersurface of the plasma membrane and extruded into extracellular space by active plasma membrane-localized hyaluronan synthase (HAS) which is transported as inactive enzyme within vesicles through the Golgi [10].

HA greatly affects tissue mechanics during wound healing as its production could escalate to 80 fold of basal level during injury and inflammation, and decline to normal level after wound remodeling [5]. It is widely accepted that cellular actions of HA are mediated by specific HA surface receptors expressed on immune cells such as CD44 [11]. Exogenous low molecular weight HA (defined as <300 kDa) was shown to enhance wound repair and increases collagen III synthesis in aged murine dermal wounds model [12]. Both matrix-bound HA content and hydrogel stiffness affect mechanical environment [13]. Moreover, high and low molecular weight HA differentially influence macrophage activation [14]. In vitro gel systems were thus designed to assess activation of macrophages infiltrating in wound with different stiffness. For example, HA polymer mixed with formulated supplemental elements to polymerize into gels with defined stiffness was used to simulate mechanical environment of stiffed wound [15].

In this protocol, we have outlined a standard protocol applying HA gel to study macrophages activation over wounds mechanical environment, such as stiffness. We demonstrated how physical changes in the stiffness of extracellular matrix (ECM) alter macrophages functions via changes in expression of inflammatory mediators. Culture of macrophages isolated from peripheral blood mononuclear cells (PBMC) on crosslinked hyaluronan hydrogels with different HA contents (stiffness) posed different inflammatory activation. Our studies confirmed that physical tension independently contributes toward macrophage activation with subsequent differential induction of gene expressions between wildtype and a mutant mouse strain. This study thus provides insight on the role of HA on macrophage function during wound healing or tissue remodeling during disease.

2 Materials

1. Peripheral blood mononuclear cells (PBMC) derived from animals (wildtype or mutant mouse strain) or human origin (*see* Note 1).

2. Ficoll Histopaque-1077 (colloidal silica particles with density 1.077 g/mL, Sigma 10771) for PBMC isolation. Allow to come to room temperature before use.
3. Cell counting (hematocytometer and 0.4% trypan blue).
4. HyStem[®]-C Hydrogel Kits (from Advanced BioMatrix) (*see* Note 2).
5. RNA extraction with TRIzol (Life Technologies) and RNeasy Mini Kit (Qiagen).
6. c-DNA synthesis with reverse transcription and gene expression polymerase chain reaction (PCR) arrays (various vendors).
7. Standard tissue culture media RPMI 1640, fetal bovine serum (FBS), and penicillin and streptomycin.
8. Confocal microscope LSM 780 (Zeiss).

3 Methods

3.1 PBMC and Macrophage Isolation

1. Collect whole blood, roughly 1 mL per mouse, into a heparinized tube (*see* Note 3).
2. In a biosafety hood, add three parts (by volume) sterile RPMI 1640 medium to dilute the whole blood in a 15-mL Falcon conical tube.
3. Pour 4 mL Ficoll Histopaque-1077 into a 15-mL conical tube for each blood sample isolation.
4. Carefully layer 4 mL diluted blood on top of each Ficoll layer, taking care not to mix blood and Ficoll (lay Falcon conical tube almost horizontal and slowly layer the diluted blood over Ficoll).
5. Perform density gradient centrifugation with a horizontal rotor. Centrifuge at $400 \times g$ for 30 min at room temperature (15–20 °C) with deceleration and acceleration both set at 1, brake off.
6. Following centrifugation, carefully remove the top plasma layer using a Pasteur pipette and discard.
7. Gently aspirate the middle PBMC layer using a Pasteur pipette and add to a new Falcon tube.

¹The PBMC could be obtained from either mouse or human origins. PBMCs contain enriched leukocytes and monocytes. This method could be applied to isolate macrophages from wildtype or mutant murine strains with immune dysfunctions [16]. Alternatively, heparinized human whole blood (65 mL) could be obtained from normal healthy donors (sourced from blood bank) or specific patient targets [17].

²HyStem[®]-C Hydrogel Kits contains three biocompatible components: thiol-modified hyaluronan (Glycosil) (100 kDa), thiol-modified denatured collagen (Gelin-S[®]), and thiol-reactive crosslinker, PEGDA (Extralink). Each is packaged as sterile lyophilized powder.

³We demonstrated mouse whole blood PBMC isolation here. For human PBMC isolation, scale up reagents accordingly. In general, routinely $1-2 \times 10^6$ PBMC will be obtained per mL blood. Assuming 10% monocyte content, $2.5-3 \times 10^7$ monocytes is expected to be isolated from 250 mL whole blood. Compared with normal subjects, abnormal macrophages behaviors might be detected from patients with remarkable keloid scar condition.

8. Wash PBMCs by adding 10 mL sterile RPMI 1640 and centrifuging at $300 \times g$ for 10 min at RT.
9. Carefully remove supernatant with a pipette to avoid disturbing the loose pellet.
10. Repeat washing step once.
11. Resuspend cell pellet in 10 mL of complete culture media (RPMI 1640–10% FBS).
12. Enrich macrophages from PBMC preparation by adherence method (**steps 13 and 14**) (*see* Note 4).
13. Resuspend the PBMCs (roughly $10\text{--}15 \times 10^6$ cells) and seed cells into T-75 flask.
14. Allow to adhere in a 5% CO₂ incubator at 37 °C for 2 h in 10 mL of complete culture media (RPMI 1640–10% FBS).
15. Remove nonadherent cells and carefully wash the adherent cells twice with RPMI 1640–10% FBS.
16. Verify cell viability and numbers (by 0.4% trypan blue exclusion method) under microscopy with hemocytometer or alternative cell counter.

3.2 Preparation of Hyaluronan Hydrogel Components with Optimal Stiffness

The desired stiffness of hydrogel is achieved by crosslinking various HA contents. In the current study, we focused on two HA stiffness groups, low (equivalent to Young's modulus 70 Pa) and high (110 Pa) [15]. A blank control was included with macrophages cultured in growth medium only (zero stiffness). Table 1 serves as an example of 24-well tissue culture format. Meanwhile, HA hydrogel preparation should be coordinated with “macrophage plating” (Subheading 3.3) to maximize cell viability and wellness in casted gel.

1. Prepare the working solutions for each of the three components from HyStem[®]-C Hydrogel Kits (*see* Note 5).

Formulate HA gels accordingly (*see* Table 1).

Low stiffness group: 1% Glycosil: 1% Gelin-S: 1% Extralink in volume ratio of 1:1:8.

High stiffness group: 2% Glycosil: 2% Gelin-S: 2% Extralink in volume ratio of 1:1:8.

2. Reconstitute and prepare 1% and 2% Thiolated collagen (Gelin-S) solution in sterile ddH₂O following manufacturer's protocol. Reconstitute each vacuumed vial with lyophilized Gelin-S powder with 0.5- or 1- volume/mL ddH₂O which would give rise to either 2% or 1%, respectively.

⁴A typical yield would be achieved for monocyte (assuming 10–25% in content) by seeding 1.0×10^6 PBMC per cm² of adhesive area. An alternative to isolate monocytes is affinity column purification where monocyte/macrophage will be retained in affinity column, while the rest of the cells elute off [16]. Follow manufacturer's protocol accordingly.

⁵We followed gel preparation method according to Nagy et al. ([15]) to achieve three different (blank control, low and high) gel stiffness by mixing sodium hyaluronate (HA) in thiolated collagen (Gelin-S) and crosslinker PEGDA.

3. Reconstitute and prepare 1% and 2% Sodium HA solution (Glycosil) in sterile ddH₂O by reconstituting each vacuumed vial with lyophilized Glycosil powder with 0.5- or 1- volume/mL ddH₂O which would give rise to either 2% or 1%, respectively.
4. Prepare the two “Glycosil: Gelin-S mixtures” (Table 1) and set aside. First, Mix 1% Gelin-S (from Subheading 3.2, **step 2**) with 1% glycosil (from Subheading 3.2, **step 3**) at 1:1 ratio by volume. Second, mix 2% Gelin-S with 2% glycosil at 1:1 ratio by volume.
5. Reconstitute and prepare 1% and 2% polyethylene glycol diacrylate (PEGDA) and set aside. Reconstitute PEGDA powder in vial with ddH₂O to generate 1% or 2% solutions. It will require four parts (by volume) of PEGDA for each one part (volume) of “Glycosil:Gelin-S mixture.”

3.3 Macrophage Plating and HA Hydrogel Polymerization

1. Prepare macrophage at a cell density of $1-5 \times 10^5$ cells per well (24-well culture plate) (*see* Note 6).
2. The hydrogel which will crosslink its constituents must be in a complete fluidity form when macrophages are mixed in.
3. Gently mix macrophages with the “Glycosil: Gelin-S mixture” (from Subheading 3.2, **step 4**) according to Table 1; ensure well mixing. This final mixture, in volume, contains macrophages, Glycosil, and Gelin-S at a ratio of 2:1:1.
4. Initiate polymerization by adding crosslinker PEGDA (from Subheading 3.2, **step 5**) into above mixture (from Subheading 3.3, **step 3**) with gentle pipetting (*see* Note 7).
5. Gel formation with embedded macrophages will proceed in sterile condition in cell culture incubator at 37 °C, 5% carbon dioxide with saturated humidity for 30 min.
6. Replenish wells with 250 μ L fresh culture medium (RPMI 1640, FBS10%) and culture for 72 h (*see* Note 8 and Fig. 1a).
7. Image macrophages in HA hydrogels with confocal microscope (Fig. 1b, c).

3.4 Gene Expression Analysis with RT-PCR Arrays

1. Remove top layer of culture media.

⁶One mL of mouse peripheral blood would yield $1-7 \times 10^6$ PBMC and 25% of them are macrophage/monocytes. A practical preparation is to pool total recovery macrophages from plastic adherence of mouse PBMC, make a master mixing, and reconstitute to 0.5 mL with culture media (RPMI 1640–10% FBS). Divide those concentrated cells to 6–8 aliquots. Each aliquot, containing approximately $1-5 \times 10^5$ monocyte/macrophages, is subjected to each experimental group analysis.

⁷Gently pipette-mixing the cells in all three hydrogel components would ensure even distribution of cells embedded in gel from bottom to top. This would ensure that macrophages respond to only calculated gel stiffness in all three experimental groups. This step is critical and might require several practices to achieve perfection.

⁸Macrophages are cultured in basal medium (RPMI 1640 with 10% FBS, glutamine, and antibiotics penicillin/streptomycin) to ensure no further activation from growth factors such as M-CSF (myeloid-colony stimulating factor (M-CSF)) and only respond to physical environmental cue.

2. Retrieve gel. Cells retrieving and lysis are done by applying 120 μ L of TRIzol to gently crush the HA hydrogel by applying gentle force and pipetting (*see Note 9*).
3. Proceed with total RNA extraction and purification with TRIzol reagent and RNeasy Mini Kit according to manufacturer's protocols.
4. Measure RNA concentration using a NanoDrop spectrophotometer and use 500–1000 ng RNA to synthesize cDNA using high-capacity RNA to cDNA kit according to manufacturer's protocols. The cDNA can be stored at -20°C until PCR reaction (*see Note 10*).
5. Dilute cDNA 10 \times with PCR-grade water.
6. Set up PCR reaction with RT² PCR array analysis according to manufacturer's protocols (*see Note 11*).
7. Alternatively, PCR could be performed with designed target primers (*see Note 12*).
8. Run RT-PCR reaction with target-specific PCR program.
9. Calculate expression of IL-1 α , IL-1 β , IL-5, MMP2 and MMP3 using the C_T method. Use average expression of the three housekeeping genes (GAPDH, actin, and HSP90) to calculate the relative expression (Fig. 2).

Acknowledgments

This work was supported by an Ohio State University intramural Lockwood Research grant (to P.H.L.).

References

1. Hamidzadeh K, Christensen SM, Dalby E et al. (2017) Macrophages and the recovery from acute and chronic inflammation. *Annu Rev Physiol* 79:567–592. 10.1146/annurev-physiol-022516-034348 [PubMed: 27959619]
2. Vannella KM, Wynn TA (2017) Mechanisms of organ injury and repair by macrophages. *Annu Rev Physiol* 79:593–617. 10.1146/annurev-physiol-022516-034356 [PubMed: 27959618]
3. Murray PJ (2018) Immune regulation by monocytes. *Semin Immunol* 35:12–18. 10.1016/j.smim.2017.12.005 [PubMed: 29290545]
4. Lin PH, Sermersheim M, Li H et al. (2017) Zinc in wound healing modulation. *Nutrients* 10(1). 10.3390/nu10010016
5. Fraser JR, Laurent TC, Laurent UB (1997) Hyaluronan: its nature, distribution, functions and turnover. *J Intern Med* 242(1):27–33. 10.1046/j.1365-2796.1997.00170.x [PubMed: 9260563]
6. Hascall VC (2019) The journey of hyaluronan research in the journal of biological chemistry. *J Biol Chem* 294(5):1690–1696. 10.1074/jbc.TM118.005836 [PubMed: 30710015]

⁹. TRIzol contains corrosive phenol and should be handled with caution in a chemical hood or a biosafety hood.

¹⁰. Use 500–1000 ng total RNA to synthesize cDNA in a final 20- μ L reaction mixture.

¹¹. Several target RT²PCR arrays are commercially available, such as inflammation, fibrosis cytokines, and chemokines arrays.

¹². Alternatively, PCR gene expression analysis could be performed with self-designed target primers. Mix 5 μ L (2 \times) Fast Sybr Green Master Mix, 0.5 μ L forward and reverse primers (from 10 μ M stock, the final primer concentration is 500 nM), and 4 μ L diluted cDNA (equivalent to original 50–100 ng original RNA) for each 384-well final RT-PCR reaction. Final volume is 10 μ L per well.

7. Tavianatou AG, Caon I, Franchi M et al. (2019) Hyaluronan: molecular size-dependent signaling and biological functions in inflammation and cancer. *FEBS J* 286(15):2883–2908. 10.1111/febs.14777 [PubMed: 30724463]
8. Jiang D, Liang J, Noble PW (2011) Hyaluronan as an immune regulator in human diseases. *Physiol Rev* 91(1):221–264. 10.1152/physrev.00052.2009 [PubMed: 21248167]
9. Spinelli FM, Vitale DL, Demarchi G et al. (2015) The immunological effect of hyaluronan in tumor angiogenesis. *Clin Transl Immunology* 4(12):e52. 10.1038/cti.2015.35 [PubMed: 26719798]
10. Weigel PH, Hascall VC, Tammi M (1997) Hyaluronan synthases. *J Biol Chem* 272 (22):13997–14000. 10.1074/jbc.272.22.13997 [PubMed: 9206724]
11. Johnson P, Ruffell B (2009) CD44 and its role in inflammation and inflammatory diseases. *Inflamm Allergy Drug Targets* 8(3):208–220 [PubMed: 19601881]
12. Damodarasamy M, Johnson RS, Bentov I et al. (2014) Hyaluronan enhances wound repair and increases collagen III in aged dermal wounds. *Wound Repair Regen* 22 (4):521–526. 10.1111/wrr.12192 [PubMed: 25041621]
13. Chen JE, Pedron S, Harley BAC (2017) The combined influence of hydrogel stiffness and matrix-bound hyaluronic acid content on Glioblastoma invasion. *Macromol Biosci* 17(8). 10.1002/mabi.201700018
14. Rayahin JE, Buhrman JS, Zhang Y et al. (2015) High and low molecular weight hyaluronic acid differentially influence macrophage activation. *ACS Biomater Sci Eng* 1(7):481–493. 10.1021/acsbiomaterials.5b00181 [PubMed: 26280020]
15. Nagy N, de la Zerda A, Kaber G et al. (2018) Hyaluronan content governs tissue stiffness in pancreatic islet inflammation. *J Biol Chem* 293 (2):567–578. 10.1074/jbc.RA117.000148 [PubMed: 29183997]
16. Zhang X, Goncalves R, Mosser DM (2008) The isolation and characterization of murine macrophages. *Curr Protoc Immunol*. Chapter 14:Unit 14.11. 10.1002/0471142735.im1401s83
17. Kelly A, Grabiec AM, Travis MA (2018) Culture of human monocyte-derived macrophages. *Methods Mol Biol* 1784:1–11. 10.1007/978-1-4939-7837-3_1 [PubMed: 29761383]




Fig. 1.

3D z-stack images of plated macrophages embedded in hyaluronan gel. HA gel could be formulated with various stiffness. **(a)** Schematic diagram depicts three various formulations to prepare three HA hydrogels (blank, low, or high stiffness), respectively. Diagrams shown side (top panels) and top (middle panels) views of cultured macrophages with different HA formulations (bottom panels). **(b and c)** Two representative z-stack microscopic images captured from a high-stiffness gel at different z-levels (11.8-mm and 12-mm for **b** and **c**, respectively). Macrophages were shown evenly embedded in HA gels with nearly identical cell density at different z-levels. Scale bars, 100 μm




Fig. 2. Stiffness of HA hydrogel affects gene expressions of isolated macrophages derived from wildtype and mutant mice. Differential gene expressions (RT² PCR array) after 72-h culture under different HA gel stiffness were detected from macrophages of wildtype and mutant strain, e.g., wound-prone (IL10^{-/-}) mice. Changes of gene expression (relative to wildtype) were indicated by folds ($2^{-\Delta C_t}$, in logarithm). Compared to wildtype macrophages, mutant macrophages displayed aberrant cytokines and MMP gene expressions in response to different microenvironmental HA stiffness

Table 1

Formulations of blank, low, and high stiffness HA gel

Macrophage in RPM11640–10% FBS	100 μ L	100 μ L	100 μ L
Glycosil (100-kDa hyaluronate)	0	50 μ L (1%)	50 μ L (2%)
Gelin-S (thiolated collagen)	0	50 μ L (1%)	50 μ L (2%)
Extralink (^a PEGDA-crosslinker)	0	400 μ L (1%)	400 μ L (2%)
RPMI 1640-10 FBS	750 μ L	250 μ L	250 μ L
Stiffness ^b (pa)	Blank (0 pa)	Low (70 pa)	High (110 pa)

^aPEGDA-Polyethylene Glycol Diacrylate Acrylate

^bStiffness is a mechanical property that a solid material preserves its shape resisting to stress applied on it. It is measured according to Young's modulus, derived as the reciprocal of the stress, presented in unit of Pa. In SI base units, Pa = kg m⁻¹ s⁻²